


S100 Ceiling Mount S100 Wall Mount



Planning Manual

M-30-1385-en

Issue 6.0

Printed on 05. 11. 2007



Key to symbols

Different symbols used in this user's manual draw your attention to safety aspects and useful tips. The symbols are explained in the following.



Warning!

The **warning triangle** indicates potential sources of danger which may constitute a risk of injury for the user or a health hazard.



Caution:

The **square** indicates situations which may lead to malfunction, defects, collision or damage of the instrument.



Note:

The **hand** indicates hints on the use of the instrument or other tips for the user.

OPMI[®]

OPMI[®] is a registered trademark of Carl Zeiss.






Contents

S100 product range	5
Ceiling mount	6
Wall mounts	6
Support units	6
General technical data	7
Minimum distances and pivot points	8
Motion ranges of S100 suspension systems	10
Standby positions	17
Installation sites in the OR or the treatment room	18
Technical data / Dentistry	23
Ceiling mount for dentistry	24
Ceiling mount for dentistry, recommended work ranges	26
Wall mounts for dentistry, dimensions	28
Wall mounts for dentistry, recommended work ranges	30
S100 Centro suspension system for dentistry	32
Technical data / ENT	33
Alternative mounting possibilities for OPMI pico in ENT	34
Ceiling mount for ENT	36
Ceiling mount for ENT, recommended work ranges	38
Wall mounts for ENT, dimensions	40
Wall mounts for ENT, recommended work ranges	42
Technical data / Ophthalmology	45
Ceiling mount for ophthalmology	46
Ceiling mount for ophthalmology, recommended work ranges	48
Wall mounts for ophthalmology	50
Wall mounts for ophthalmology, recommended work ranges	52

Technical data of examination units	55
S100 examination unit, with illumination module	56
S100 examination unit including light guide	57
Customer's preparatory responsibilities	59
Customer's responsibilities	60
Planning the electrical installation	64
Anchoring of the S100 ceiling mount	66
Anchoring of the S100 wall mount	70
Confirmation of the structural calculation and execution of installation	75
Index	77

S100 product range

Ceiling mount	6
Wall mounts	6
Support units	6

Ceiling mount	Cat. No.: 1151-515	
		
Wall mounts	with 950mm suspension arm	with 600 mm suspension arm
	Cat. No.: 1080-347	Cat. No.: 1148-883
		
Support units		
for attachment to external examination unit,	including light guide and CZ illumination module	including light guide
	Cat. No.: 1095-091	Cat. No.: 1095-094
		

General technical data


Minimum distances and pivot points	8
Motion ranges of S100 suspension systems	10
Motion range of S100 ceiling mount	10
Motion range of the wall mount with a 950mm suspension arm	12
Motion range of the wall mount with a 600mm suspension arm	14
Standby positions	17
Possible standby positions of the wall mounts	17
Installation sites in the OR or the treatment room	18
Ceiling mount: ideal installation above the patient	18
Ceiling mount: alternative installation to the patient's side	19
Wall mount with a 950mm suspension arm, possible installation sites	20
Wall mount with a 600mm suspension arm, possible installation sites	21

Minimum distances and pivot points

Minimum distances

For the installation of the wall mount, please observe the necessary minimum distances from walls as a function of the selected standby position; also see the drawing below. Please also make sure that the mount cannot collide with doors or other examination configurations.

Any collisions with doors and other examination configurations must also be avoided in the case of the ceiling mount.




Determining the pivot points

The specifications for the installation areas always refer to the pivot points of the mounts. The illustration below shows you the positions of the pivot points.

Key

- 1 Pivot point of wall mount
- 2 Pivot point of ceiling mount



Motion ranges of S100 suspension systems

Motion range of S100 ceiling mount



Notes:

- The motion and work range of the suspension system is determined by the 45° and 122° positions between the suspension arm and the carrier arm. The system can be easily moved between these two positions, allowing reliable positioning of the surgical microscope. Plan the location of the suspension system relative to the patient in such a way that you normally work exactly in the middle between these two positions. It is possible to move the arms a little beyond the positions given.
- The ceiling mount provides a circular motion and work range of 360°. This is achieved by the swivel range of the carrier arm mounted on the column of 163° and the possibility of also folding the suspension arm toward the carrier arm in the opposite direction.

Key

- 1 45° position of the suspension arm relative to the carrier arm, corresponds to a radius of 700 mm (27.56")
- 2 122° position of the suspension arm relative to the carrier arm, corresponds to a radius of 1,300 mm (51.18")
- 3 Carrier arm, length 500 mm (19.69"), provides a swivel range of ±163° about the column, supports the suspension arm
- 4 Motion and work range of the ceiling mount
- 5 Suspension arm, length 950 mm (37.40"), provides a swivel range of ±163° relative to the carrier arm



Caution:


When planning and installing the mount, make sure that the orientation of the swivel range is correct.



- 6 Center of the carrier arm's swivel range relative to the ±163° stops

Caution:

When planning and installing the mount, make sure that the orientation of the swivel range is correct.



Motion range of the wall mount with a 950mm suspension arm




Notes:

- The motion and work range of the suspension system is determined by the 45° and 122° positions between the suspension arm and the carrier arm. The system can be easily moved between these two positions, allowing reliable positioning of the surgical microscope. Plan the location of the suspension system relative to the patient in such a way that you normally work exactly in the middle between these two positions. It is possible to move the arms a little beyond these positions.
- The wall mount with a 950mm suspension arm provides a semi-circular motion and work range of 180°. This is achieved by the swivel range of the carrier arm of 105° and the possibility of also folding the suspension arm toward the carrier arm in the opposite direction.

Key

- 1** 45° position of the suspension arm relative to the carrier arm, corresponds to a radius of 700 mm (27.56")
- 2** 122° position of the suspension arm relative to the carrier arm, corresponds to a radius of 1,300 mm (51.18")
- 3** Carrier arm, length 500 mm (19.69"), provides a swivel range of $\pm 105^\circ$, supports the suspension arm
- 4** Length of the suspension arm, 950 mm (37.40")
- 5** Center of the carrier arm's swivel range relative to the $\pm 105^\circ$ stops
- 6** Motion and work range of the wall mount



Motion range of the wall mount with a 600mm suspension arm



Notes:

- The motion and work range of the suspension system is determined by the 45° and 122° positions between the suspension arm and the carrier arm. The system can be easily moved between these two positions, allowing reliable positioning of the surgical microscope. Plan the location of the suspension system relative to the patient in such a way that you normally work exactly in the middle between these two positions. It is possible to move the arms a little beyond these positions.
- The wall mount with a 600 mm suspension arm provides a semi-circular motion and work range of 180°. This is achieved by the swivel range of the carrier arm of 105° and the possibility of also folding the suspension arm toward the carrier arm in the opposite direction.

Key


- 1 45° position of the suspension arm relative to the carrier arm, corresponds to a radius of 700 mm (27.56")
- 2 122° position of the suspension arm relative to the carrier arm, corresponds to a radius of 1,300 mm (51.18")
- 3 Carrier arm, length 500 mm (19.69"), provides a swivel range of ±105°, supports the suspension arm
- 4 Suspension arm, length 600 mm (23.62"), provides a swivel range of ±163° relative to the carrier arm



Caution:

When planning and installing the mount, make sure that the orientation of the swivel range is correct.

- 5 Center of the carrier arm's swivel range relative to the ±105° stops
- 6 Motion and work range of the wall mount




Standby positions

Possible standby positions of the wall mounts


The following illustrations show the possible standby positions for the wall mounts with 950mm and 600mm suspension arms.

Wall mount with a 950mm suspension arm



1:20


Wall mount with a 600mm suspension arm



1:20


Installation sites in the OR or the treatment room


Ceiling mount: ideal installation above the patient




1:20

Ceiling mount: alternative installation to the patient's side



Wall mount with a 950mm suspension arm, possible installation sites


Wall mount with a 600mm suspension arm, possible installation sites

Technical data / Dentistry

Ceiling mount for dentistry	24
Dimensions of ceiling mount used in dentistry	24
Determining the ceiling height for dentistry	25
Ceiling mount for dentistry, recommended work ranges	26
Work range recommended for ideal installation site	26
Work range recommended for alternative installation site	27
Wall mounts for dentistry, dimensions	28
Wall mount with a 950mm suspension arm for dentistry	28
Wall mount with a 600mm suspension arm for dentistry	29
Wall mounts for dentistry, recommended work ranges	30
Work ranges recommended for wall mount with a 950 mm suspension arm	30
Work ranges recommended for wall mount with a 600mm suspension arm	31
S100 Centro suspension system for dentistry	32

Ceiling mount for dentistry

Dimensions of ceiling mount used in dentistry



Determining the ceiling height for dentistry



Caution:

- Taking the notes into account, please answer the following questions. Only then will Carl Zeiss be able to process your order properly and to install the ceiling mount correctly.
- Enclose a completed copy of this page with each order.

Order No.:

Customer address / Delivery address:

- Please measure exactly and enter the results here:

Ceiling height

(room height from floor to structural ceiling):

..... mm

The column length of the ceiling mount will be adapted to the ceiling height specified in the order.

For **dentistry**, determine the column length of the ceiling mount as follows:

Ceiling height - 1810 mm (71.26") = column length

..... mm - 1810 mm (71.26") = mm



Note:

The maximum possible column length is 1500 mm (59.06")!

The minimum possible column length is 200 mm (7.87")!


This results in the following limit values for the room height in dentistry: maximum 3310 mm (130.31"), minimum 2010 mm (79.13")

Should you obtain a different value in your calculation, please contact Carl Zeiss.


The ceiling flange must be mounted by the customer. See "Anchoring of the S100 ceiling mount", page 66.

Ceiling mount for dentistry, recommended work ranges

Work range recommended for ideal installation site




Work range recommended for alternative installation site




Wall mounts for dentistry, dimensions

Wall mount with a 950mm suspension arm for dentistry

A No stop




Wall mount with a 600mm suspension arm for dentistry




Wall mounts for dentistry, recommended work ranges

Work ranges recommended for wall mount with a 950 mm suspension arm



Work ranges recommended for wall mount with a 600mm suspension arm




S100 Centro suspension system for dentistry

The S100 Centro suspension system is an S100 suspension system without column. This version has been specially designed for installation on the Centro column from KaVo.

An adapter from KaVo (KaVo Cat. No.: 1.002.0345) permits the S100 Centro suspension system to be installed on the Centro column.

Interface (1) between the Centro column and the S100 Centro suspension system (i.e. the axis of adapter 1.002.0345) must be aligned in a vertical position (max. deviation $\pm 0.5^\circ$).

1 Interface between Centro column and S100 Centro suspension system




Technical data / ENT

Alternative mounting possibilities for OPMI pico in ENT	34
Ceiling mount for ENT	36
Dimensions of ceiling mount for ENT	36
Determining the ceiling height for ENT	37
Ceiling mount for ENT, recommended work ranges	38
Work range recommended for ideal installation site	38
Work range recommended for alternative installation site	39
Wall mounts for ENT, dimensions	40
Wall mount with a 950mm suspension arm for ENT	40
Wall mount with a 600mm suspension arm for ENT	41
Wall mounts for ENT, recommended work ranges	42
Work ranges recommended for wall mount with a 950 mm suspension arm	42
Work ranges recommended for wall mount with a 600mm suspension arm	43


Alternative mounting possibilities for OPMI pico in ENT

- 1 Direct vertical mounting of OPMI pico in the downward direction on the S100 suspension arm via the microscope's support arm.
- 2 Mounting of OPMI pico with a 120° coupling on the S100 suspension arm.



Ceiling mount for ENT

Dimensions of ceiling mount for ENT



Determining the ceiling height for ENT



Caution:

- Taking the notes into account, please answer the following questions. Only then will Carl Zeiss be able to process your order properly and to install the ceiling mount correctly.
- Enclose a completed copy of this page with each order.

Order No.:

Customer address / Delivery address

- Please measure exactly and enter the results here:

Ceiling height (room height from floor to structural ceiling): mm
Admissible limit values for ENT	Maximum height: 3,450 mm (138.98") Minimum height: 2,250 mm (88.58")

The column length of the ceiling mount will be adapted to the ceiling height specified in the order.

For **ENT**, determine the column length of the ceiling mount as follows:

Ceiling height - 1,960 mm (77.17") = length of the column

.....mm - 1,960 mm (77.17") = mm



Note:

The maximum possible column length is 1,500 mm (59.06")!


The minimum possible column length is 200 mm (7.87")!

Should you obtain a different value in your calculation, please contact Carl Zeiss.


The ceiling flange must be mounted by the customer. See "Anchoring of the S100 ceiling mount", page 66.

Ceiling mount for ENT, recommended work ranges

Work range recommended for ideal installation site




Work range recommended for alternative installation site




Wall mount with a 950mm suspension arm for ENT

A No stop




Wall mount with a 600mm suspension arm for ENT




Wall mounts for ENT, recommended work ranges

Work ranges recommended for wall mount with a 950 mm suspension arm



Work ranges recommended for wall mount with a 600mm suspension arm




1:20

Technical data / Ophthalmology

Ceiling mount for ophthalmology	46
Dimensions of ceiling mount in ophthalmology	46
Determining the ceiling height for ophthalmology	47
Ceiling mount for ophthalmology, recommended work ranges	48
Work range recommended for ideal installation site	48
Work range recommended for alternative installation site	49
Wall mounts for ophthalmology	50
Wall mount with a 950mm suspension arm for ophthalmology	50
Wall mount with a 600mm suspension arm for ophthalmology	51
Wall mounts for ophthalmology, recommended work ranges	52
Wall mount with a 950mm suspension arm	52
Wall mount with a 600mm suspension arm	53

Ceiling mount for ophthalmology

Dimensions of ceiling mount in ophthalmology



Determining the ceiling height for ophthalmology



Caution:

- Taking the notes into account, please answer the following questions. Only then will Carl Zeiss be able to process your order properly and to install the ceiling mount correctly.
- Enclose a completed copy of this page with each order.

Order No.:

Customer address / Delivery address

- Please measure exactly and enter results here:

Ceiling height

(room height from floor to structural ceiling):

..... mm

Admissible limit values for ophthalmology:

Maximum height: 3,530 mm
(138.98")
Minimum height: 2,250 mm (88.58")

The column length of the ceiling mount will be adapted to the ceiling height specified in the order.

For **ophthalmology**, determine the column length of the ceiling mount as follows:

Ceiling height - 2,040 mm (80.31") = length of the column

..... mm - 2,040 mm = mm
(80.31")



Note:

The maximum possible column length is 1,500 mm (59.06")!


The minimum possible column length is 200 mm (7.87")!

Should you obtain a different value in your calculation, please contact Carl Zeiss.


The ceiling flange must be mounted by the customer. See "Anchoring of the S100 ceiling mount", page 66.

Ceiling mount for ophthalmology, recommended work ranges

Work range recommended for ideal installation site



Work range recommended for alternative installation site



Wall mount with a 950mm suspension arm for ophthalmology

Technical drawing of a mobile crane showing dimensions in millimeters (mm) and inches (").

Horizontal Dimensions:

- 320 mm / 12.60"
- 500 mm / 19.69"
- 950 mm / 37.40"

Vertical Dimensions:


- 2158 mm / 84.96"
- 1734 mm / 68.27"
- 1624 mm / 63.27"
- 954 mm / 37.56"
- 1650 mm / 64.96"
- 400 mm / 15.75"
- 296 mm / 11.65"
- 400 mm / 15.75"

Other Dimensions:

- ± 105° (Rotation angle)
- A (Point of rotation)

Scale: 1:20


Wall mount with a 600mm suspension arm for ophthalmology



Wall mounts for ophthalmology, recommended work ranges


Wall mount with a 950mm suspension arm

Ideal site of installation and work range for ophthalmology



Wall mount with a 600mm suspension arm

Ideal site of installation and work range for ophthalmology




Technical data of examination units

S100 examination unit, with illumination module	56
--	-----------

S100 examination unit including light guide	57
--	-----------

S100 examination unit, with illumination module


S100 examination unit for attachment to an external system carrier
(incl. light guide and CZ illumination module)



1:20

S100 examination unit including light guide

S100 examination unit for attachment to an external system carrier
(including light guide)



1:20

Customer's preparatory responsibilities

Customer's responsibilities	60
Checking the installation conditions	60
Building vibrations	62
Planning the electrical installation	64
Anchoring of the S100 ceiling mount	66
Mounting the ceiling flange	66
Aligning the column in a vertical position	68
Constructional requirements for the ceiling mount	68
Anchoring of the S100 wall mount	70
S100 wall flange (1244-708)	70
Aligning the wall flange in a vertical position	71
Wall plate for S100 wall mount	72
Constructional requirements for wall mounts	74
Confirmation of the structural calculation and execution of installation	75

Customer's responsibilities

Checking the installation conditions



Note:

Zeiss service staff can mount the ceiling or wall mount only if all points of the following checklist applicable to the relevant installation conditions have been fulfilled.

The actual load on the ceiling or wall depends on a large number of different factors, which must be determined in detail by a structural engineer on a case-to-case basis.

- Constructional requirements for ceiling mounts, see page 68,
- Constructional requirements for wall mounts, see page 74.



Caution:

- Make sure that a structural engineer checks the installation conditions during the planning procedure. Structural verification must be performed prior to the installation of the mount.
We recommend filing the structural verification in the ceiling or wall mount documentation.
- Obtain a written confirmation from a structural engineer stating that the applicable national codes and regulations have been complied with.
- Please add a copy of the "Confirmation of structural calculation" to your order (see page 75).
- If any differences exist between the planning documents and the actual on-site situation, please inform your contact at Carl Zeiss or the planning expert prior to the installation of the ceiling or wall flange or of the pre-installation set.
- The system should be installed on a hard, level surface, e.g. concrete. Do not countersink the ceiling flange, pre-installation set or wall flange, but mount it directly on the respective surface.
- On-site conditions also include building vibrations, which the structural engineer responsible must take into account right during the planning phase (see page 62).
Obtain a written confirmation from your structural engineer stating that possible building vibrations have been taken into account (see page 75).

- Obtain the structural verification for an existing substructure from a structural engineer, i.e. make sure that the engineer confirms the effective strength of the existing anchors.
Structural verification must be performed prior to the installation of the mount.
Please add a copy of the "Confirmation of structural calculation" to your order (see page 75).
- If an existing substructure is to be used, make sure that the maximum inclination of the ceiling anchor plate under load does not exceed 0.5°
- The structural engineer must check that no modifications have been made to the substructure since the original installation.
- When ordering a ceiling mount for installation on an existing substructure, you have to take into account the height of the existing substructure in the calculation of the column length.

**Warning!**

- If an existing flange plate or intermediate piece must be exchanged, never re-use the old anchors. New anchor holes must be drilled.
When calculating the effective strength of the new anchors, the structural engineer must take into account the weakening effect of the old holes in the ceiling or wall.

Planning the installation

- Inspect the mounting components supplied for completeness and damage.
- Prior to the installation of a ceiling or wall mount, you must always check that the actual installation conditions - in particular the room height for ceiling mounts - correspond to the specifications in the drawing.
- Make sure that the anchors calculated by your structural engineer are properly installed, that the nuts and washers required to mount the Zeiss flange are readily available at the installation site and that the maximum tightening torque for the nuts is indicated.
- For new installation, the ceiling flange must be mounted before the false ceiling is attached.

Building vibrations

The on-site requirements also include the low-vibration design of the ceiling in the OR. This must be taken into account right during the planning phase for the ceiling mount.

The following information primarily refers to the ceiling mount, but it can also be applied to a wall mount or floor stand.

Two types of excitation factors must be distinguished:

Single events which excite short-term vibration

Induced by inadvertent knocks against the suspension system or strong impact against ceiling, wall or floor. This is the most frequent, but least critical type of excitation. The ZEISS ceiling mount features excellent damping against this type of vibration and displays a short recovery time. In extreme cases, surgery has to be interrupted for a few moments.

Constant excitation causing sustained vibration

The excitation energy of factors such as elevators, air conditioning systems, construction work, traffic does not easily reach the ceiling mount via the building. This type of excitation is extremely rare, but may lead to permanent vibration of the ceiling mount in extreme cases. This becomes particularly visible at high magnifications as used in the surgical microscope.

The following result of a study is intended to help you understand constant vibration excitation occurring in rare cases in ceiling mounts, eliminate it or prevent it from the outset in new installations.

Like any kinematic system, the Zeiss ceiling mount displays eigenfrequencies (=resonance frequencies) that range between 2 and 80 Hz depending on the position of the system and the accessories attached to it. The ZEISS ceiling mount provides very effective damping in particular against higher frequencies above 10 Hz which are typical of buildings with electrical excitation factors. Nevertheless it may happen in rare cases, when the ceiling mount's suspension arm is in a specific position, that even high frequencies which are critical for this particular constellation lead to vibrations (e.g.: 17 Hz, ..., 19 Hz, ...).

Such excitation factors can usually be suppressed by the following measures:


- elimination of the excitation source (e.g. repair of damping of the air conditioning system)
- constructional damping measures in the ceiling installation

Due to the large number of parameters involved and the variety of potential building/ceiling mount constellations, Carl Zeiss Surgical is unable to give an absolute guarantee for the vibration-free suspension of the ceiling mount, even if the building meets the applicable ISO standards.

However, it is highly improbable that constant vibrations are transferred to the ceiling mount if all requirements regarding vibrations in the OR ceiling are met.

- Max. vibration velocity (RMS) at the installation points for the ceiling mount.
- $V_{\max} < 0.1 \text{ mm / s}$ or $V_{\max} < 4 \text{ milli-inches / s}$ or below the curves (diagram) for the specified frequency range.

Sources: Carl Zeiss Surgical (in-house study), ISO 10811, recommendations for ORs.



Planning the electrical installation

The following work is the customer's responsibility:


- A socket with a properly installed protective ground connection must be provided at the installation site near the ceiling flange. Terminals for power connection are located on the ceiling mount.
- Make sure that the power outlet has a properly connected protective ground connector.
- Potential equalization: take the necessary actions in the building to incorporate the system in protective "potential equalization" measures.

Installation of conduits

- In the planning process, make appropriate allowance for conduits required for additional future applications (communications, video, monitors, etc.).

Power requirements at the installation site

- Please note the applicable national codes and regulations concerning the line cross section and fuse strength.
- Please use the following electrical specifications as a basis for on-site electrical protection measures:

Rated voltage	115 VAC, (100 - 120) VAC \pm 10% 230 VAC, (220 - 240) VAC \pm 10%
Rated frequency	50 - 60 Hz
Current consumption	Current consumption of halogen light source: Max. 2.0 A at 115 VAC Max. 1.0 A at 230 VAC Current consumption of xenon light source: Max. 5.0 A at 115 VAC Max. 2.5 A at 230 VAC
Electrical standard	Complying with IEC 60601-1 / EN 60 601-1; UL 60 601-1, CAN / CSA-C22.2 601.1 -M90 Protection class I, degree of protection IPXO, Type B equipment 

The system has been designed for continuous operation.

Anchoring of the S100 ceiling mount

For recommended installation sites of the ceiling mount, please see the illustrations from page 18.



Caution:

The ceiling mount must only be installed on the original ceiling flange supplied by the company Carl Zeiss.

Mounting the ceiling flange




Caution:

When planning and installing the system, make sure that the orientation of the swivel range is correct:

- The swivel range of the carrier arm about the column is limited by two stops at $\pm 163^\circ$. In accordance with the installation drawings, the swivel range of $\pm 163^\circ$ must be marked on the column and ceiling flange with a waterproof pen.
The two stops determine the range beyond which the carrier arm cannot reach.
- Two stops at $\pm 163^\circ$ determine the range which the suspension arm cannot reach.


Key

- A** Ceiling flange, dia. 300 mm (11.81")
- B** Ceiling height
2,250 mm (88.58")
Maximum 3,300 mm (129.92") in dentistry
Maximum 3,450 mm (135.83") in ENT
Maximum 3,530 mm (138.98") in ophthalmology
- C** Cover, dia. 361 mm (14.21")
- D** Column, dia. 60 mm (2.36")
- E** Column length



1 : 10

Aligning the column in a vertical position



- The column must be aligned in a vertical position, with a maximum admissible deviation of $\pm 0.5^\circ$:
 - Use C-type washers 305497-0002-000 to level out any unevenness of the mounting surface.
t = 1 mm / 0.04", thickness of the C-type washer
 - Or, if the mounting surface is very uneven, you can use longer threaded bolts with additional nuts and washers to align the flange component.

Constructional requirements for the ceiling mount

The actual load on the ceiling depends on a large number of different factors. The requirements to be met by the ceiling or substructure result from the addition of perpendicular forces and torques produced by the suspension system and accessories. These are the forces to be transmitted into the structural ceiling via the ceiling anchors,
 The column must be aligned in a vertical position (max. deviation $\pm 0.5^\circ$).

Forces and torques



Warning!

The structural engineer must ensure in each individual case that the structural ceiling has a sufficient load capacity for the forces and torques listed below. He must also take into account any additional loads on the ceiling and add an appropriate safety margin, and must observe the applicable national codes and regulations.



Caution:

- The perpendicular forces and torques specified below include an additional load of 100 kg (220 lb), which is generated when a person hangs on the end of the suspension arm (or the microscope). Further safety margins have **not** been incorporated.
- The perpendicular forces and torques have been calculated on the basis of the maximum permissible load on the suspension arm.
- For an S100 ceiling mount installed on the ceiling flange, the structural ceiling must have the following load capacity:
 - Perpendicular force: min. 1200 N (min. 270 lbf)
 - Torque: min. 1350 Nm (min. 996 lbf.ft)

Anchoring of the S100 wall mount

For recommended installation sites of the wall mount, please see the illustrations from page 20.



Caution:

The wall mount must only be installed on the original wall flange or, if necessary, on the original wall plate from Carl Zeiss.

S100 wall flange (1244-708)



Caution:

When planning and installing the mount, make sure that the orientation of the swivel range is correct:

- Two stops at $\pm 163^\circ$ determine the range which the carrier arm cannot reach.
- In the wall mount with a 600mm suspension arm, two stops at $\pm 163^\circ$ determine the range which the suspension arm cannot reach.

Key


A Cover, W x H = 213 x 426 mm (8.39 x 16.77")

B S100 wall flange (1244-708)

C Distance of top row of bores from floor
1910 mm (75.20") for dentistry
2060 mm (81.10") for ENT
2140 mm (84.25") for ophthalmology

t = 12 mm / 0.47" thickness of the wall flange plate


Aligning the wall flange in a vertical position



- The wall flange must be aligned in a vertical position, with a maximum admissible deviation of $\pm 0.5^\circ$:

Use C-type washers 305497-0002-000 to level out any unevenness of the mounting surface.

$t = 1 \text{ mm} / 0.04''$, thickness of the C-type washer



Wall plate for S100 wall mount

The wall plate is used to install the wall mount on less stable walls. Compared with direct installation using the wall flange only, the wall plate distributes the torques over a larger area, thus reducing the load on individual points.

The wall plate is the main component of the mounting kit "Adapter for S100, complete" (1277-816).

The mounting kit also includes caps for covering any bores in the wall plate which you do not need.



Caution:


- The original wall flange from Carl Zeiss must always be used to install the wall mount on the wall plate. The wall plate must not be inclined by more than 0.2° under the load of the wall mount.

1 Wall plate for S100 wall mount (1277-816)
shown here with screwed-on S100 wall flange (1244-708)

2 S100 wall flange (1244-708)

C Installation height of wall plate and wall flange
Distance of row of bores from floor:
1910 mm (75.20") for dentistry
2060 mm (81.10") for ENT
2140 mm (84.25") for ophthalmology

t = 12 mm / 0.47" thickness of the wall plate



Constructional requirements for wall mounts

The actual load on the wall depends on a large number of different factors. The requirements to be met by the wall or substructure result from the addition of perpendicular forces and torques produced by the suspension system and accessories. These are the forces to be transmitted into the structural wall via the wall mount.

The customer is responsible for providing a wall with sufficient load capacity. Light-construction walls, for example, can be reinforced by additional, high-stability posts.

The wall flange must be aligned in a vertical position, with a max. deviation of $\pm 0.5^\circ$.

Forces and torques



Warning!

The structural engineer must ensure in each individual case that the structural wall has a sufficient load capacity for the forces and torques listed below. He must also take into account any additional loads on the wall and add an appropriate safety margin, and must observe the applicable national codes and regulations.



Caution:

- The perpendicular forces and torques specified below include an additional load of 100 kg (220 lb), which is generated when a person hangs on the end of the suspension arm (or the microscope). Further safety margins have **not** been incorporated.
- The perpendicular forces and torques have been calculated on the basis of the maximum permissible load on the suspension arm, i.e. a surgical microscope with maximum admissible accessories.
- For an S100 wall mount, the structural wall must be able to accept the following forces and torques:

S100 wall mount with a 950mm suspension arm

- Perpendicular force: min. 1,110 N (min. 250 lbf)
- Torque, vertical: min. 1,690 Nm (min. 1,250 lbf.ft)
- Torque, horizontal: min. 1,340 Nm (min. 990 lbf.ft)

S100 wall mount with a 600mm suspension arm

- Perpendicular force: min. 1,090 N (min. 245 lbf)
- Torque, vertical: min. 1,350 Nm (min. 1,000 lbf.ft)
- Torque, horizontal: min. 1,020 Nm (min. 755 lbf.ft)

Confirmation of the structural calculation and execution of installation

Sales order no.:

Customer address /

Delivery address:

By signing below, the following persons confirm that they have performed their work in a proper and orderly way:

The **structural engineer** for

- the selection and layout of the installation site, taking possible building vibrations into account
- the structural calculation, taking into account the applicable national regulations and the planning manual
- the structural checking of an existing substructure
- the structural calculation of a substructure built on site
- the final checking and release of the structural calculations:

Name and address of the structural engineer:

.....

Date

.....

Signature

The **installer** for

the proper mounting of the pre-installation set or the ceiling or wall flange from Carl Zeiss:

Name and address of the executing company and installer's name:

.....

Date

.....

Signature

Index

A

Adapter for S100 wall mount, compl. (1277-816) 72
Additional load of 100 kg 68, 74
Aligning the column in a vertical position 68
Anchoring of S100 ceiling mount 66
Anchoring of S100 wall mount 70

B

Building vibrations 60, 62

C

Carrier units, technical data 55
Ceiling flange 66
Ceiling flange, mounting of 66
Ceiling flange, original from Carl Zeiss 66
Ceiling height 25, 37, 47, 66
Ceiling mount 6
Ceiling mount for dentistry 24
Ceiling mount for dentistry, recommended work ranges 26
Ceiling mount for ENT 36
Ceiling mount for ENT, recommended work ranges 38
Ceiling mount for ophthalmology 46
Ceiling mount for ophthalmology, recommended work ranges 48
Ceiling mount, alternative installation on patient's side 19
Ceiling mount, ideal installation above the patient 18
Centro suspension system 32
Checking the installation conditions 60
Collision with other examination configurations 8
Column diameter 66
Column length 66
Column length for dentistry 25
Column length for ENT 37
Column length for ophthalmology 47
Column, vertical alignment 68
Conduits, installation of 64
Confirmation of execution of installation 75
Confirmation of structural calculation 75
Constructional requirements 60
Constructional requirements for ceiling mount 68
Constructional requirements for wall mounts 74
Cover 66, 70
C-type washers (305497-0002-000) 68, 71
Customer's responsibilities 60

Customer's preparatory responsibilities 59

D

Dentistry, technical data for 23

Determining the ceiling height for dentistry 25

Dimensions of ceiling mount for dentistry 24

E

ENT, technical data for 33

Execution of installation, confirmation of 75

F

Forces and torques 68, 74

I

Installation conditions, checking of 60

Installation height of wall plate and wall flange 72

Installation site 60

Installation site, layout 75

Installation sites in OR or treatment room 18

Installation, proper 75

K

Key to symbols 2

L

Load on the ceiling or wall 60

Load on the ceiling, actual 68

Load on the wall, actual 74

Loads on the wall, additional 74

M

Minimum distances 8

Motion range of S100 ceiling mount 10

Motion range of wall mount with 600mm suspension arm 14

Motion range of wall mount with 950mm suspension arm 12

Mounting height, determining 25, 37, 47

Mounting possibilities for OPMI pico in ENT 34

O

Ophthalmology, technical data 45

Original ceiling flange 66

Original wall flange 70, 72

Original wall plate 70

P

Perpendicular force 68, 74

Pivot point of ceiling mount 9

Pivot point of wall mount 9

Pivot points 8

Planning the electrical installation 64

Potential equalization 64
Power requirements at installation site 64
Product range 5
Protective ground connection socket 64

R

Reinforcing the light-construction walls 74
Review and release 75
Row of bores, distance from floor 70

S

S100 examination unit incl. light guide 57
S100 examination unit, with illumination module 56
S100 product range 5
Signature, confirm by 75
Standby position 8
Standby positions 17
Structural calculation 75
Structural calculation, confirmation of 75
Structural engineer 60, 68, 74, 75
Structural engineer, obligatory written confirmation 60
Structural verification 60
Substructure, produced on-site 75
Substructures, existing 60
Support units 6
Swivel range, orientation of 66, 70

T

Technical data / Dentistry 23
Technical data / ENT 33
Technical data for ophthalmology 45
Technical data of examination units 55
Threaded bolts, longer 68
Torque 68
Torque, horizontal 74
Torque, vertical 74

V

Vibration, short-term 62
Vibration, sustained 62

W

Wall flange 70
Wall flange, original from Carl Zeiss 72
Wall flange, vertical alignment 71
Wall mount 6
Wall mount for ENT 40
Wall mount on less stable walls 72
Wall mount with 600mm suspension arm, possible installation sites 21

Wall mount with 600mm suspension arm, standby position 17
Wall mount with 950mm suspension arm, possible installation sites 20
Wall mount with 950mm suspension arm, standby position 17
Wall mount, constructional requirements 74
Wall mount, ENT, with 600mm suspension arm 41, 50, 51
Wall mounts for dentistry, dimensions 28
Wall mounts for dentistry, recommended work ranges 30
Wall mounts for ENT, dimensions 40
Wall mounts for ENT, recommended work ranges 42
Wall mounts for ophthalmology 50
Wall mounts for ophthalmology, recommended work ranges 52
Work ranges, recommended, ceiling mount for dentistry 26
Work ranges, recommended, ceiling mount for ENT 38
Work ranges, recommended, ceiling mount for ophthalmology 48
Work ranges, recommended, wall mount for ophthalmology 52
Work ranges, recommended, wall mounts for dentistry 30
Work ranges, recommended, wall mounts for ENT 42

Z

Zeiss service 60



Carl Zeiss Surgical GmbH
A Carl Zeiss Meditec Company

73446 Oberkochen
Germany

Fax: ++49 (0) 7364-20 48 23
Email: surgical@meditec.zeiss.com
Internet: www.meditec.zeiss.com