



iCAM

User Manual

April 2006

Copyright 2006 DTA s.r.l.

All rights reserved. The reproduction of any part of this manual is allowed only with the written authorization by DTA s.r.l..

The contents of this manual may be subject to changes without any warning.

DTA are not responsible for errors that may occur in this manual.

Revision 1.5 dated 05/04/2006

CONTENTS

<i>INTRODUCTION</i>	<i>4</i>
<i>PERSONAL COMPUTER MINIMUM REQUIREMENTS</i>	<i>4</i>
<i>PERSONAL COMPUTER RECOMMENDED REQUIREMENTS</i>	<i>4</i>
<i>SCHEME OF THE INSTALLATION PROCEDURE</i>	<i>5</i>
<i>CONNECTING THE FDL INTERFACE CARD TO THE PC</i>	<i>6</i>
<i>LIBRARY & FDL-PCI INSTALLATION</i>	
<i>WINDOWS XP</i>	<i>7</i>
<i>LINUX</i>	<i>11</i>
<i>WINDRIVER REGISTRATION</i>	<i>13</i>
<i>CONNECTING CABLES AND POWER LEADS</i>	<i>16</i>
<i>9-PIN USER PORT CONNECTOR</i>	<i>17</i>
<i>POWER SUPPLY CONNECTOR</i>	<i>17</i>
<i>44-PIN PARALLEL PORT CONNECTOR</i>	<i>18</i>
<i>OPTICAL WINDOW CLEANING</i>	<i>19</i>
<i>SPECIFICATIONS</i>	<i>20</i>
<i>INDEX</i>	<i>21</i>

INTRODUCTION

The iCam cooled CCD camera series is specifically designed for industrial and scientific applications. This series of cameras aims to making the most of readout speed of Kodak KAF sensors: the goal to reach is the best possible performance by exploiting its high sensitivity in the shortest possible image downloading time.

Every component is integrated in the aluminium head: CCD, electromechanical shutter, Single stage Peltier cooling system, electronics. An external motorized filter wheel is available, as optional.

The camera is connected to the PC by means of a FDL-PCI Interface card (alternatively, a fiber optic connection can be made). Thanks to its cooling method, 12-bit (multirange) A/D Converter and remarkable speed during image acquisition (it can reach up to 12.5 Mpixel/s of readout speed), iCAM is the ideal solution to shoot extremely weaks events, subject to quick changes.

Among the main application fields: Biology, Astronomy, Plasma Physics, Semiconductor Physics



PERSONAL COMPUTER MINIMUM REQUIREMENTS


- CPU Pentium III 500 Mhz .
- 256 Mb RAM.
- Microsoft Windows 98.
- Hard Disk with at least 200 Mb available.

PERSONAL COMPUTER RECOMMENDED REQUIREMENTS

- CPU Pentium IV 2 GHz or higher.
- 512 Mb RAM.
- Colour monitor 19".

SCHEME OF THE INSTALLATION PROCEDURE

When you use the camera for the first time, you must install the library and the FDL-PCI interface. For this purpose, please follow the procedure described in the next pages and schematized in the diagram below:




CONNECTING THE FDL INTERFACE CARD TO THE PC

To install your FDL interface card, do as follow :

- ⇒ Turn your PC off, including any peripheral.
- ⇒ Remove the external case of the PC (please see the relevant instructions in the PC Handbook).
- ⇒ Choose a free PCI expansion slot. Remove the relevant back cover by means of a screwdriver.
- ⇒ Insert the camera interface card into the expansion slot. Please make sure that the card is properly and fast inserted. Fix the card by screwing the relevant screw again.
- ⇒ Reassemble the PC case.

Turn the PC on once again. The interface card installation has been completed.



FDL-PCI interface card

In case of a fiber optic link connection, a dedicated FDL-PCI is provided with the proper connectors.

LIBRARY & FDL-PCI INSTALLATION

The installation of the library is always linked to the installation of the camera, of course. For this purpose you just need to install the PCI interface as well.

Let's analyse the sequence of operations to be carried out. It will be shown the example with WINDOWS XP Operating System. The operation are similar for the other Operating Systems.


Once you turn the system on for the first time with the FDL-PCI interface on, you will be asked to specify where the available drivers are. Follow the steps on the basis of your Operative System.

WINDOWS XP


- 1) Insert the FDL-PCI card into the slot of the motherboard and turn the PC on. As soon as the system gets started, a window (like the one shown below) will appear, noticing you a new hardware was found. You'll be requested to connect for Windows update, let select: "No, not at this time".




- 2) Insert the "ViSTA" CD into the CD-ROM reader, close the "Quick Install" (it may open due to the CD autorun). Select the option "Install from a list or specific location (Advanced)" (recommended choice). Then click on "Next" to continue the installation.



- 3) Check the function selected in the window (recommended choice) below (in this case, E:\ indicates the CD-ROM drive) and click on next to continue. The PC will start to search for the most suitable driver to the peripheral.



- 4) Select the folder that contains drivers for FDL-PCI hardware and click on OK. Then, wait while the wizard installs the software (just few seconds).



- 5) When the installation has been completed, a window like the one below will appear. Click on finish to close the wizard and restart the PC (recommended choice).



LINUX

System Requirements

- Linux 2.0.31 or higher (can be embedded Linux or Linux 2.4x)
- An x86 processor.
- Any 32-bit development environment supporting C (such as GCC).

Preparing the system for installation

In Linux, kernel modules must be compiled with the same header files that the kernel itself was compiled with. Since WinDriver installs the kernel module **windrvr6.o**, it must compile with the header files of the Linux kernel during the installation process.

Therefore, before you install WinDriver for Linux, verify that the Linux source code and the file **versions.h** are installed on your machine:

Install linux kernel source code

- If you have yet to install Linux, please choose **Custom** installation when performing the installation and then choose to install the source code.
- If Linux is already installed on the machine, you must check to see if the Linux source code was installed. You can do this by looking for linux in the **/usr/src** directory. If the source code is not installed, you can either reinstall Linux with the source code, as described above, or you can install the source code by following these steps:

1. Login as super user.
2. Type:
/\$ rpm -i /<source location>/ <Linux distributor>/RPMS/kernel-source-<version number>
(For example: to install the source code from the Linux installation CD-ROM, for RedHat 7.1, type:
/\$ rpm -i /mnt/cdrom/RedHat/RPMS/kernel-source-2.4.2.-2.i386rpm)

TIP!

If you do not have an RPM with the source code, you may download it from:
<http://rpmfind.net/linux/RPM/>.

Install version.h

- The file **version.h** is created when you first compile the Linux kernel source code. Some distributions provide a compiled kernel without the file **version.h**. Look under **/usr/src/linux/include/linux/** to see if you have this file. If you do not, please follow these steps:

1. Type:
/\$ make xconfig
2. Save the configuration by choosing **Save and Exit**.
3. Type:
/\$ make dep.

Before proceeding with the installation, you must also make sure that you have a linux symbolic link. If you do not, please create one by typing:

```
/usr/src$ ln -s <target kernel>/ linux
```

(For example: for Linux 2.4 kernel type:

```
/usr/src$ ln -s linux-2.4/ linux)
```

Installation

1. Insert the **DTA SDK CD** into your Linux machine CD drive or copy the downloaded file to your preferred directory.
2. Change directory to your preferred installation directory (your home directory, for example):
`/$ cd ~`
3. Extract the file **WDxxxLN.tgz** (where xxx is the version number):
`~$ tar xvfz /<file location>/WDxxxLN.tgz`
For example:
 - o From a CD:
`~$ tar xvfz /mnt/cdrom/LINUX/WDxxxLN.tgz`
 - o From a downloaded file:
`~$ tar xvfz /home/username/WDxxxLN.tgz`
4. Change directory to WinDriver (this directory gets created by tar):
`~$ cd WinDriver/`

NOTE:

From version 5.x and above this directory gets created by tar, but in versions preceding 5.x the WinDriver directory does not get created by the extraction. Therefore, when working with versions preceding 5.x (version 4.33, for example) first create a directory (e.g., WinDriver) before proceeding with the installation.

`(/$ mkdir ~/WinDriver)`

5. Install WinDriver:
 - a. `~/WinDriver$ make`
 - b. Become super user:
`~/WinDriver$ su`
 - c. Install the driver:
`~/WinDriver# make install`
6. Create a symbolic link so that you can easily launch the DriverWizard GUI
`~/WinDriver$ ln -s ~/WinDriver/wizard/wdwizard/ usr/bin/wdwizard`
7. Change the read and execute permissions on the file **wdwizard** so that ordinary users can access this program.
8. Change the user and group ids and give read/write permissions to the device file **/dev/windrvr6** depending on how you wish to allow users to access hardware through the device.
9. You can now start using WinDriver to access your hardware and generate your driver code!

Restricting Hardware Access on Linux


CAUTION:

Since **/dev/windrvr6** gives direct hardware access to user programs, it may compromise kernel stability on multi-user Linux systems. Please restrict access to the DriverWizard and the device file **/dev/windrvr6** to trusted users.

For security reasons the WinDriver installation script does not automatically perform the steps of changing the permissions on **/dev/windrvr6** and the DriverWizard executable (**wdwizard**).

WINDRIVER REGISTRATION

- 1) After the installation of the drivers for the FDL-PCI card, we need to record the files to make them work properly. Insert the "ViSTA 3" CD into the CD-ROM reader, wait a few seconds so as to allow the PC to load the "QUICK INSTALL" menu.



- 2) Once the PC has loaded the "QUICK INSTALL" menu, follow the steps below:


- Register the drivers according to your Operating System, clicking on the proper **Drivers Registration** field:
FDL-PCI on XP/2000 (for WINDOWS XP/2000)
FDL-PCI on 98/Me (for WINDOWS 98/Me)
- Install ViSTA 3 just clicking on the proper field

Once you have installed the drivers and the "ViSTA 3" software, restart the PC.

When you restart the PC, it is suggested to check the correct installation of the peripheral, that is if the installation of the necessary drivers for the FDL-PCI card and WINDRIVER has been carried out good.

According to the different Operating Systems, you have to click on "System Properties" by means of the right key of the mouse. The images below refer to the case you use Windows XP as OS.

A menu with different options will open; click on "Hardware" by means of the left key of the mouse and on "Device Manager". Then check if the FDL-PCI peripheral has been installed correctly (see the images below).



The tests have been carried out with the following OS:

WINDOWS 98 SECOND EDITION

WINDOWS 2000 WITHOUT SERVICE PACK

WINDOWS 2000 WITH SERVICE PACK 3/4/5

WINDOWS ME

WINDOWS XP WITHOUT SERVICE PACK


WINDOWS XP WITH SERVICE PACK 1/2

CONNECTING CABLES AND POWER LEADS

Now it is possible to connect the cables between the interface and the camera.

We advise you to carry out this operation with all your **instruments not powered**.

The connectors are fixed by means of screws; these, however, should only be slightly tightened to avoid any damage to them. For the iCAM series, a proper power supply with cable is provided. The camera is turned on at the power supply.




In the figure above, you can find the connectors required to connect the camera to the PC and to the power supply.

Starting from the top left clock-wise, the following connectors are shown:

- USER PORT: 9-pin connector on which you can find synchronization pins and serial communication with the camera
- Power connector
- BNC connector for trigger input (TTL)
- PARALLEL PORT: high density 44-pin connector for parallel communication by means of a given FDL-PCI card.

9-PIN USER PORT CONNECTOR




1	RX	2	TX	3	VSYN	4	HSYN	5	TRG
6	GND	7	SHUT	8	SHTD1	9	SHTD2	-	-

- **RX:** controller serial input
- **TX:** controller serial output
- **VSYN:** vertical synchronism input
- **HSYN:** horizontal synchronism input
- **TRG:** input signal to start the acquisition.
- **GND:** ground.
- **SHUT:** shutter output signal.
- **SHTD1/ SHTD2:** driver signal of the shutter. **WARNING: 65V signal on SHTD1 and SHTD2 pins.**


Note: All signals are TTL compatible, if not differently reported.

POWER SUPPLY CONNECTOR



Frontal View

44-PIN PARALLEL PORT CONNECTOR



1	D0	2	D3	3	D6	4	D9	5	D12
6	D15	7	ACK	8	R/W	9	GND	10	GND
11	D18	12	D21	13	D24	14	D27	15	D30
16	D1	17	D4	18	D7	19	D10	20	D13
21	VCC	22	REQ	23	STB	24	GND	25	D16
26	D19	27	D22	28	D25	29	D28	30	D31
31	D2	32	D5	33	D8	34	D11	35	D14
36	VCC	37	A/D	38	GND	39	GND	40	D17
41	D20	42	D23	43	D26	44	D29	-	-

- **D0 – D31**: data, if the camera being controlled by the HR IV is a Hurricane, all the 32 given bits are significant, otherwise in case of a HiRes only the less significant 16 bits are significant.
- **REQ**: by means of this pin, the platform starts the transfer request to the PC.
- **ACK**: response signal to the REQ signal from the PC.
- **A/D**: by means of this pin, it is indicated the presence of addresses or data on the bus.
- **R/W**: indicates whether the operation to execute is a reading or a writing.
- **STB**: STROBE signal.
- **GND**: ground.
- **VCC**: +5V.

NOTE: all signals are TTL compatible.

By means of this port, it is possible to carry out a parallel connection between CCD Camera and PC.

OPTICAL WINDOW CLEANING

Both the optical window and the CCD cleaning are carried out in the clean room by means of a 30-magnifying power microscope. This procedure removes any dust which can otherwise bring about unmistakable marks on the image you have taken. In particular such spots increase as the focal ratio gets wider. In other words an image may not show any mark at f/5.6 but it can be clearly noted at f/32 because of an obvious geometrical problem of projection. Due to the shutter or elapsed time, the external surface of the optical window may gather dirty particles that may be easily removed.

To perform such task we use a compressed air cylinder specifically designed for optical cleaning.

Be careful !!! There are similar products that, instead of using compressed air, use a liquefiable gas : at ambient pressure it quickly gassifies, thus “triggering” an air-compressed-effect. You do **not** absolutely have to use these products : they may give rise to heavy marks or rings on the windows itself.



A product we can recommend is DUST-OFF provided by EDMUND-OPTICS. Thanks to DUST-OFF (or any other similar product) it is very easy to get rid of any microparticles : keep the shutter open for a few seconds (by the camera control program) and spray some air blast. We kindly advise you against using cloths, optical paper and cleaning liquids because the dirt will be only mixed up or, even worse, increased. The risk is to finally damage the coating of the optical window itself !

SPECIFICATIONS

SHUTTER:

Electromechanical. Exposure time: from 0.005 s to 9999 s

A/D CONVERTER:

12 bit

SETTABLE GAINS:

64

READ OUT SPEED:

up to 12.5 Mpixel/s

INTERFACE:

FDL-PCI 8, 16, 32 bit.

COOLING:

Single stage PELTIER (-40°C below ambient)

CCD TEMPERATURE CONTROL:

± 0.1 °C

OPTICAL WINDOW:

BK7 with double AR coating

FILTER WHEEL:

External (optional)

BACKFOCUS:

17.5 mm

MAX TOTAL NOISE:

25 e⁻

MOUNT:

According to different sensors

BINNING:

From 1 x 1 to 8 x 8 or arbitrary

POWER SUPPLY:

230V

MAX ABSORPTION:

50 W

WEIGHT:

1.5 kg

DIMENSIONS:

95x95x250 mm³

INDEX

9

9-PIN USER PORT CONNECTOR 17

44

44-PIN PARALLEL PORT CONNECTOR 18

C

CONNECTING THE FDL INTERFACE CARD TO THE PC 6

CONNECTING CABLES AND POWER LEADS 16

I

INTRODUCTION 4

L

LIBRARY & FDL-PCI INSTALLATION 7

LINUX 11

O

OPTICAL WINDOW CLEANING 19

P

PERSONAL COMPUTER MINIMUM REQUIREMENTS 4

PERSONAL COMPUTER RECOMMENDED REQUIREMENTS 4

POWER SUPPLY CONNECTOR 17

S

SCHEME OF THE INSTALLATION PROCEDURE 5

SPECIFICATIONS 20

W

WINDOWS XP 7

WINDRIVER REGISTRATION 13



Viale Campania 23
56021 Cascina (Pisa) - ITALY
web site: www.dta.it e-mail: sales@dta.it