

**ADC COMPACT  
AUTOPROCESSING  
SOFTWARE**

---

USER MANUAL



# TABLE OF CONTENTS

## **SWITCHING THE PROCESSING STATION ON AND GETTING THE SYSTEM MAIN SCREEN**

Switching the processing station on	1
Switching the processing station off	3

## **THE SYSTEM MAIN SCREEN**

Selecting and de-selecting images	6
Selecting by criteria	6
Sorting images	7
Printing a selected image	8
Sending images to review stations	11
Archiving images	13
Monitoring the status of the processing station	14
Managing messages	14
Locking Screen	17

## **SYSTEM MONITORING**

Disk management	18
Input queue management	19
Output queue management	21
Setup	23
Queues	27
Devices	28

## **TROUBLESHOOTING**

Automatic anticipation	29
General procedures in case of malfunctions	29

## **PERFORMANCE REQUIREMENTS OF THE PROCESSING STATION 33**

## **PERFORMANCE REQUIREMENTS - MONITORS 35**

# Switching the processing station on and getting the system main screen


Normally, the processing station may be left on day and night. The processing station is equipped with a screen saver to prevent the screen from burning in when it is not used for a long period of time. However, if you need to move the processing station to a different location or if it has to be disconnected from the mains, proceed as follows when switching the processing station on or off.

## Switching the processing station on

Proceed as follows to switch the processing station on.


- 1 Put the main switch in the ON position. The main switch is located on the back of the processing station close to the power supply cable.

The processing station performs its initialization procedure and asks your login name.



- 2 Type your user name and press [ENTER] or click the OK button.

The processing station asks your password.




- 3 Type your password and press [ENTER] or click the OK button.

The processing station automatically loads the application software and displays the system main screen.



- 2 Click 'Stop user interface'.

The processing station displays the Log-out menu.



- 3 Click the Power Off button for a complete shut-down of the processing station.



If you click the Log-out button, the system returns to the Log-in prompt.



- 4 Click the OK button to proceed.





The images displayed on the alphanumeric image list are selected by default in accordance with the selection criteria that have been configured in the Setup module of the system main screen. In the example screen above, the preconfigured selection criteria are 'Patient name', 'Identification date', 'Examination' and the 'Sent' feature.

### Selecting and de-selecting images

To perform an operation on an image such as printing, sending or archiving an image, you first have to select it.

You can select an image by just clicking on it. A selected image appears on a white background.

You can deselect an image by clicking it again, or you can click 'Deselect' to de-select the currently selected image.

### Selecting by criteria

You can also use selection criteria. The following criteria are available:

- by identification date
- by patient ID
- by patient name
- by patient last name
- by radiologist
- by examination type
- by examination subtype
- by RIS ID
- by Digitizer ID
- by User Info (i.e. four user configurable data fields).

You can configure four criteria from this list. For an explanation on how to configure the selection criteria, we refer to the section 'Setup' on page 23. Each of the selected criteria has a matching button in the upper left corner of the screen.

When no selection is entered, all images appear on the screen (on several successive pages).

Example: If you click 'Patient name' then select a given patient from the menu, the screen will display images of this patient only. You can also combine selection criteria: e.g. you may want to select the images of one examination type from a given date. The images that appear on the screen all fulfil those four conditions. If a button is marked with an asterisk '\*' it means that this selection criterion is inactive.

### Sorting images

The default sort algorithm puts the most recently identified images at the head of the list. However, you can change this order by clicking the 'Sort on:' button.

You have 8 possibilities to sort on:

- identification date
- patient name
- radiologist
- examination
- patient ID
- RIS ID
- Digitizer ID
- User Info (i.e. four user configurable data fields).


You can also use the 'Sort on:' button in combination with the selection criteria mentioned in the previous sub-section.

## Printing a selected image

To print an image, proceed as follows:

- 1 Select the image you want to print.
- 2 Click 'Print'.

The processing station returns the print menu:




- 3 Click a button from the top row to select the size of film you want to print onto.

Depending on the configuration of the printer destination and the printer, you have the choice between 3 sizes:

- 8 x 10 inch
- 11 x 14 inch
- 14 x 17 inch


The button of the film size and layout as selected by the ADC Compact ID Software is pressed automatically by default. The printer destination as configured in the ID Station determines which size buttons can be selected. The Scopix laser printers print on all three sizes, Drystar 2000 only on 8 x 10 inch film size and the Drystar 3000 printer on 11 x 14 inch and 14 x 17 inch. As a result only the size selection buttons that are relevant are selectable. In the illustration below, for example, only the 8 x 10 " film size button can be activated, as the printer destination was initially configured to the Drystar 2000 (Di 2000) printer.



- 4 Click a button from the next row to select the basic orientation of the image or image boxes.


You can choose between:

- the portrait format



or


- the landscape format



- 5 The next set of buttons determines which part of the actual film surface is used when a 1 on 1 reproduction needs to be printed.


You can choose either of the following two options:

- Full size, i.e. the full film surface is printed on (while preserving the aspect ratio)



or

- True size, i.e. any cassette size is printed “true size” (=100%) on all film sizes (See standard film division table on pages 22 and following of the ADC Compact System Overview Manual). In a number of cases, for example, when a 35 x 45 cm cassette size is printed on a 8 x 10 inch film size, parts of the image will be clipped at all sides of the image.



- 6 In case you wish to assign the image to be printed the highest priority in the printing queue, click the Emergency button.



7 Select Printer



8 Click 'Print'.

Or click 'Cancel' if you decide not to print the selected image(s).



### Sending images to review stations

The IMPAX product line provides gateways, transmit, archive and review stations (DICOM STORE = optional) and hard copy printing (DICOM PRINT = standard) for digital imaging. The processing station features an interface software module to integrate ADC Compact and IMPAX. You can interactively send 8 or 12-bit images to IMPAX modules using the DICOM protocol. 8-bit images already have a Look-up Table, whereas 12-bit images are MUSICA-processed with separate Look-up Table. To allow for fully-fledged CR functionality also private elements have been defined. Refer to the DICOM Conformance Statement. For more information on the IMPAX product line, consult your local Agfa sales organisation.

Only images that have been given a default destination on the ADC Compact ID Station can be sent to an IMPAX review station.

If you want to send (re-send) a processed image to an IMPAX review station, proceed as follows:


1 Select the image that you want to send to the review station.

The 'Send' button becomes active.

**Note**

You can select an image by just clicking on it. A selected image appears on a white background.

The following picture shows the system main screen with one image selected.



Note You can deselect an image by clicking on it again, or you can click 'Deselect' to deselect the currently selected image.

- 2 Click 'Send'.



## Archiving images

Only images that have been assigned a default destination on the ADC Compact ID Station can be archived to an Impax archive station.


If you want to archive a processed image, proceed as follows:

- 1 Select the image that you want to archive. In case the selected image was previously assigned a default archiving station destination, the 'Archive' button becomes active.

### Note

You can select an image by just clicking on it. A selected image appears on a white background.

The following picture shows the system main screen with one image selected.



### Note

You can deselect an image by clicking on it again, or you can click 'Deselect' to deselect the currently selected image.

- 2 Click 'Archive'.

## Monitoring the status of the processing station

With the dynamic bar graph in the upper right corner of the browser screen you can monitor the status of the processing station.

- D: indicates the degree of occupation of the hard disk
- I: indicates the current input rate
- O: indicates the current output rate



These three meters indicate the activity of the processor and of the incoming and outgoing images. The input rate and output rate indicators show between zero and four jobs. If there are more than four jobs, the meter turns white.

The activity of the background module is indicated as well:

- X: indicates that the background is not running
- \ alternating with /: indicates that the background is running

## Managing messages

If the application wants to inform you on certain events, conditions or problems, it generates a message. A window pops up displaying the message. After you have read the message and clicked 'OK', the message is put in the message queue.

Messages can have different priorities: this is reflected in the status of the message. The status of a message can be one of the following:


- message  
A message with a 'message' status is intended to notify the user of certain events.

- warning  
A 'Warning' status informs the user about abnormal conditions.
- panic message  
These messages inform the operator about failed operations.

To manage the message queue, proceed as follows:

- 1 From any screen, click the button labelled 'message' in the upper right corner.

Clicking this button opens a window with information about the oldest message with the highest priority. An example of this window is shown below:



The following information can be displayed:

- description of the message and the status (message, warning or panic)
- a possible solution
- the application that has generated the message
- date and time the message was generated
- frequency

- message ID-number
- the number of messages as well as the position of the currently selected message in the queue are also displayed (e.g. 2/2 means that this message is the second out of two).

If the problem persists and you contact your Agfa service engineer, be sure that you have this information at hand.

2 You have the following options:

- Click 'Previous' to view the message that comes before the currently selected message in the queue.
- Click 'Next' to view the message that comes after the currently selected message in the queue.
- Click 'View' to pop up a window that lists all messages in the queue. This window looks similar to the following one:



3 Click the 'Cancel' button.

- 4 If you have solved the problem and you want to delete the currently selected message from the queue, click the 'Delete' button.

As long as a panic message is not deleted, it continues to pop up.

- 5 Click the 'Cancel' button.

### Locking Screen

In order to forestall accidental erasure of data whenever the operator temporarily leaves the processing station, the screen will turn blank and be locked **automatically** after a (re-configurable) lapse of time.

The screen can also be locked **interactively**. In order to lock the screen, proceed as follows:

- 1 Click the Controller menu.
- 2 Click Lock Screen.

The processing station returns the following message:



In order to unlock the screen, enter your password.

# System Monitoring

## Disk management

From the general System Monitoring screen, click the button labelled 'Disk' in the upper left corner.



Clicking this button opens a window with information about the disk structure.



The system gives you the number of stored images on your system and the amount (in Megabytes) used. Since, as an option, you can install more than one disk, the window shows the name of the disk.

You can also see the number of images that you can add without overwriting older images.

**Note**

The number of images is an approximation since different images can have different sizes.

Click the 'Cancel' button to return to the general System Monitoring screen.

Input queue management

**Input queue window**

From the general System Monitoring screen, click the button labelled 'Input' in the upper left corner. Clicking this button opens a window with information about the incoming images. An example of this window is shown below:



The screen shows the patient demographics and the examination data of the images that are currently queued for the processing station.

Clicking the 'Cancel' button returns you to the general System Monitoring screen.

### **Changing the job priority**

Normally, every job has a priority of 5. The processing station reads jobs according to the first-in-first-out rule. However, you can change the sequence.

If you want a job to jump the queue, proceed as follows:

- 1 From the general System Monitoring screen, select the button labelled 'Input'.

The Input window pops up.

- 2 Select the job that you want to give a higher priority.

The selected job is highlighted.

- 3 Click the 'Emergency' button

The job moves to the first place in the queue and its priority changes from 5 to 0.

- 4 Click 'Cancel' to return to the general System Monitoring screen.

### **Holding jobs**

The normal status of a job is 'queued'. If the processing station is actually reading the image, the job status changes into 'active'.

If you select a job and click the 'Hold' button, the job status becomes 'hold'. The processing station will not read the image unless you select it again and click the 'Unhold' button.

### **Cancel jobs**

If you want to cancel a job, proceed as follows:

- 1 From the general System Monitoring screen, select the button labelled 'Input'.

The Input window pops up.

- 2 Select the job that you want to remove from the queue.



The background of the job becomes white.

- 3 Click the 'Delete' button.

The status changes in 'Cancelled'.


*Caution* If a job status is 'Cancelled', you can, by no means, retrieve the job.

- 4 Click 'Cancel' to return to the general System Monitoring screen.

Output queue management

### Output queue window

From the general System Monitoring screen, click the button labelled 'Output' in the upper left corner. Clicking this button opens a window with information about the images in the output queue. An example of this window is shown below:



The screen shows the patient demographics of the images that are currently queued for the different output devices such as printers, archive stations and review stations.

Clicking the 'Cancel' button returns you to the general System Monitoring screen.

## **Changing the job priority**

Normally, every job has a priority of 5. Consider, as an example, a Laser Imager. Jobs are printed according to the first-in-first-out rule. However, you can change the sequence.

If you want a job to jump the queue, proceed as follows:

- 1 From the general System Monitoring screen, select the button labelled 'Output'.

The Imager window pops up.

- 2 Select the job that you want to give a higher priority.

The background of the job becomes white.

- 3 Click the 'Emergency' button.

The job moves to the first place in the queue and its priority changes from 5 to 0.

## **Holding jobs**

The normal status of a job is 'queued'. If an output device is actually processing the image, the job status changes into 'active'.

If you select a job and click the 'Hold' button, the job status becomes 'hold'. The output device will not process the image unless you select it again and click the 'Unhold' button.

## **Cancel jobs**

If you want to cancel a job, proceed as follows:

- 1 Select the job that you want to remove from the queue.

The background of the job becomes white.

- 2 Click the 'Delete' button.

The status changes in 'Cancelled'.

*Caution* If a job status is 'Cancelled', you cannot retrieve the job.

## Setup

From the general System Monitoring screen, click the button labelled 'Setup'. Clicking this button opens a window that allows you to configure the system. An example of this window is shown below:



Clicking the 'Save' button stores the currently configured settings whereas clicking 'Cancel' returns to the general System Monitoring screen.

### Selection criteria

The selection criteria that you configure are used by default to select the images in the main screen. Select four different criteria from the following list:

- by identification date
- by patient ID
- by patient name
- by patient last name
- by radiologist
- by examination type
- by examination subtype
- by RIS ID
- by Digitizer ID
- by User Info (i.e. a user configurable data field).
- by “Sent”

In case “Sent” has been configured as a selection criterium in the Set-up menu, the system gives you a choice from among six options:





The “sent” and “archived” options only comprise images that have been sent respectively archived manually, whereas the “printed” option also includes the autorouted images.

### Wait-message Icons


You are offered a choice between 'Simple' and 'Fancy' wait-message icons.

By default the wait-message icon assumes the shape of an egg-timer. By clicking the 'Fancy' button in the Set-up menu, you can turn the icon into a more fanciful shape, like the ones displayed below:



### Panic message level

The panic message levels can be set by clicking the 'More' respectively 'Fewer' buttons in the Panic message level field of the Set-up menu, as shown in the screen below.



A selection can be made between three different levels:

- Level 1  
When the system is set to level 1, any warning or error message irrespective of its importance is displayed and highlighted on the screen.
- Level 2  
When panic message level 2 is selected, the system displays warnings or error messages, whenever errors of a somewhat higher priority than in level 1 occur.
- Level 3  
When level 3 is selected in the Set-up menu, the system only displays and highlights messages of the highest priority.

#### **Max. number images**



The maximum number of storable images is calculated on the basis of the disk space available in the system.

This configurable number of storable images can be raised (never lowered), as the initial system calculations are based on storage capacity required for storing original exposed images and as subsequently processed images account for far less storage space.

## Queues

Using the Queues option in the Set-up menu, you can define periods for transmitting both the Archive and Review/Softcopy queues. In case the queues are not set by the user, the system sends the images to the archive and/or review stations continuously, by default.


For both queues you can pre-set three different off-peak periods during which you want images to be archived or sent to a review station.

In order to do so proceed as follows:

- 1 Click the Queues button in the Set-up menu.
- 2 Click and drag the pointer across the period you want the images to be sent to the Archive respectively Review destinations.

The period can also be indicated by keyboarding the time of day in the appropriate boxes.

The system displays the indicated off-peak period(s) on a black background.



## Devices

Using the Devices option in the Set-up menu, you can consult a table showing the last known status of all devices in the system, as displayed in the following sample screen.



The screenshot shows a window titled "Device status." with a list of devices and their statuses. The list is as follows:

Device	Status
0001	Off
1234-L1	Off
ADC_DR0A	Unknown
ADC_LR11	Unknown
ADC_LR0A	Unknown
ADC_SC0	Unknown
ADC_SC12	Unknown
ADC_T02	Unknown
ADC_A58	Unknown
ADC_R017	Unknown
ADC_DR56	Unknown
ADC_DR512	Unknown
ADC_R01	Unknown
ADC_A58	Unknown
ADC_A017	Unknown
ADC_DR56	Unknown
ADC_DR512	Unknown
ADC_A01	Unknown
ADC_SC0	Unknown
ADC_SC12	Unknown

A "Cancel" button with a close icon is located in the bottom right corner of the window.



# Troubleshooting

## Automatic anticipation

The processing station is programmed to anticipate problems. The actions taken are fully automatic and most of them are transparent to the operator of the processing station. The automatic actions are explained below:

- Every night, the processing station automatically makes a consistency check between the Oracle database of the images and the files of the images to anticipate database errors.
- If an internal software problem occurs during operation, the processing station automatically restarts the MIMOSA software.

## General procedures in case of malfunctions

Although the processing station is programmed to anticipate problems and to correct malfunctions automatically, a user intervention is sometimes required. If you face a problem that is not automatically solved you can take one of the following actions:

- Reset the user interface
- Reset the background
- Reboot the system.

Always try to reset the user interface first. If this does not solve the problem, reset the background. Rebooting the system is the ultimate action you can undertake. If this does not solve the problem, contact a service engineer.

## Reset and stop the user interface

If the software performs slower than normal or faces a deadlock, resetting the user interface is the most appropriate action. To reset the user interface, proceed as follows:

- 1 Make sure the main screen is displayed.

If not, click 'Select' in the upper right corner of the current screen.

- 2 Click 'Controller'.

The processing station returns the controller menu.



- 3 Select 'Reset User interface'

If the previous actions did not solve the problem, then continue with step 4.

- 4 Click 'Controller' again.

- 5 Select 'Stop User Interface'.

The processing station displays the Log-out menu.



- 6 Click either the Log-out button or Power Off button.



If you click the Log-out button and then the OK-button, the system asks you for confirmation and then returns to the Log-in prompt.



If you click the Power Off button and then the OK-button, the system returns the following screen:

OK to Power off Workstation  
[ yes, no, ?, q ]

- 7 Type "yes" and press [ENTER].

The system informs you that the Shut-down procedure has started and then returns the OK prompt. Switch off the computer, or type "boot" to start up the system again.

### **Reset the background**

This action solves typical background problems. To reset the background, proceed as follows:

- 1 Make sure the main screen is displayed.  
  
If not, click 'Select' in the upper right corner of the current screen.
- 2 Click 'Controller' in the right corner of the screen and select 'Reset Background' to reset the software.

### **Check System**

This function clears the system of corrupt data, if any. To check the system, proceed as follows:

- 1 Make sure the main screen is displayed.  
If not, click 'Select' in the upper right corner of the current screen.
- 2 Click 'Controller' in the right corner of the screen and select 'Check System'.

### **Software Version**

When selecting 'Software Version' in the 'Controller' roll-down menu, the system displays the software version and the operating system.

# Performance Requirements of the processing station

## Hardware

- SUN ULTRA 1 model 170
- 128 MB RAM
- 1 internal disk 2.1 GB
- 5 GB DAT tape drive (optional)
- 1 internal CD ROM drive
- 1 floppy drive 3.5 inch
- 1 keyboard (Country specific)
- 1 mouse + mousepad

## Software

- Operating system Unix (Sun Solaris 2.5)
- Graphic User interface: Xwindows
- Relational data base: Oracle 7.2.2
- Motif
- Agfa ADC Compact Autoprocessing software
- Agfa Musica image processing software

### Application software

- Chest
- Abdomen
- Skeleton

### Optional application software

- Uro / Tomo
- Dental (in development)
- Paediatrics (in development)
- Full Leg / Full Spine (in development)

The following applications are not yet supported by the ADC Compact System:

- Angiography
- Mammography

# Performance Requirements - Monitors

## Standard Brightness Monitor

### Synchronization

Horizontal Scan Rate	48 to 108 kHz (automatic adjustment)
Vertical Scan Rate	60 to 80 Hz (automatic adjustment)

### Input Signals

ECL Digital	D-Sub 25 pin connector (1, 2 or 4 bits per pixel)
Analog	3 BNC connector (composite on video or separate sync.)
	13W3 and other connectors optional

### CRT

Screen Size	21" FS
Phosphor	P104 (other options available)
Non-glare AR Panel	Optional
Display Mode	Non-interlaced

### Display Area

Horizontal Width	16" max (user adjustable)
Vertical Height	12" max (user adjustable)

### Video

Video Band Width	200 MHz (up to 300 MHz optional)
------------------	----------------------------------

## Resolution

Horizontal	1024 to 2048
Vertical	768 to 1536
Brightness	65 Footlamberts nominal HRMS™ allows the monitor to adjust automatically to a range of horizontal and vertical scan frequencies, providing compatibility with many hi-res display controllers for PCs, PS/2s, MACs, Sun, Apollo, DEC, Silicon Graphics and other workstations.

## User Controls

Side	Power, Brightness and Contrast Controls
Rear	Height, Vertical Position, Width, Horizontal Position

## Dimensions

Height	16.5" (19" with tilt-swivel base)
Width	19.5"
Depth	19"
Weight	48 lbs. (55 lbs. with tilt-swivel base)

## Power Supply

Power Consumption	100 Watts
Input Voltage	90 to 130 VAC, 220 to 250 VAC
Input Frequency	47 to 63 Hz



## High Brightness Monitor

**Greyscale Monitor** SMM 2183 L

### Power requirements

Input Voltage	90 - 264 V wide range power supply
Power Frequency	47 - 65 Hz
Power Consumption	Max 150 Watts
Power factor control	according to IEC 1000
Power saving	according to VESA DPMS

### CRT specifications

Size	21" Flat Square
Deflection Angle	90°
Light Transmission	30%
Phosphor Type	P45, cadmium free
Surface	Multicoated conductive panel (AR/AS)
Gun System	Dispenser Cathode (long life)
Focusing	static and dynamic

### General performance

Horizontal Frequency	30 to 83 kHz
Vertical Frequency	50 to 120 Hz
Formats	19 Max, self recognizing and autoscaling

### Display performance

Display Area (WxH)	400 mm x 300 mm
Non-Linearity	≤ 2%
Raster stability	0.05 mm Max swim&jitter
Maximum Brightness	≥ 175 ftL @ 1 Vp-p Input Signal
High Voltage Regulation	0,2% Max size change

### Video amplifier, inputs

Bandwidth	140 MHz @ 80 V Modulation
Rise & Fall times	3 ns
Connectors	BNC Type
Impedance	75 ohms
Video Level	0.75 to 1.2 Vp-p
Sync Level	0.1 to 0.6 Vp-p

### Front panel controls

Controls	Power on/off, Contrast/Brightness
Geometry	H/V-phase, H/V-amplitude, pin & barrel, raster rotation
Ambient light sensor	for automatic contrast control

### Operating conditions

Temperature	
- Operating	+ 10 °C...+ 35 °C
- Storage	- 25 °C...+ 70 °C
Humidity	
- Operating	20% - 80% rel. humidity
- Storage	20% - 80% rel. humidity

### Mechanical specifications

Dimensions (WxDxH)	499 x 520 x 476 mm
Weight	approx. 30.5 kg

### Approvals

Safety	UL 1950, CSA, IEC 950, EN 60 950, DHHS
EMC	IEC 601-1-2
Radiation	MPR II
Others	CE mark TÜV-Ergo (optional) ISO 9001 certified plant





Printed in Belgium  
Published by Agfa-Gevaert N.V., B-2640 Mortsel-Belgium  
2202A GB 199801

