




Historage-Elec Controls



Leica @ BRDG

***Caution: Never attempt to hold or grab any component of the machine to stop motion***

# Microtome SOP



Mechanical Controls



American Optical @BRDG

Prepared by: Bob Morrison

STLCC , Instrumentation Specialist

May 2008, updated Dec 2015 for models at BRDG


# Microtome: American Optical, Model 820, SN 70298



*Hotlink to Effective Use and Care of  
Microtome.....(LARGE pdf, 93 pgs)*

*Hotlink to American Optical, Model 820 Reference Manual.....(pdf)*

# Microtome: American Optical, Lube Schedule



Key to oil and grease points on the AO "E20" Microtome

- O<sub>1</sub> 2 drops Pike Oil every month for average usage and 2 drops every two weeks for heavy usage.
- O<sub>2</sub> 2 drops of Pike Oil every month.
- G Grease surface every 6 months for average usage and every 3 months for heavy usage.
- ★ Relubricate knife holder parts and knife holder slideway with thin annealing oil or Oil, or light grease, after each cleaning.

For maintenance instructions see Instruction Manual 820-301 (Effective Use and Proper Care of the Microtome) Chapter III, section D, Page 18.  
NOTE: For best results, do not operate faster than one revolution per second. For difficulties in "Sectioning Difficulties", see manual 820-301, Chapter II, Page 6.

AMERICAN OPTICAL CORPORATION

SCIENTIFIC INSTRUMENT DIVISION

BRIDGEPLAND, N. Y. 14220  
Bridgeland, N. Y. 14220



# Microtome: Leica, Model SM 2000R




[\*Hotlink to Leica Model SM2000R, Reference Manual.....\(pdf\)\*](#)

# Microtome : General Setup and Operating Instructions

1. Prepare specimen in paraffin or desired mounting material
2. Before plugging in the device, mount specimen and the knife and adjust per page 4
3. Plug in the device and set electronic controls per page 5
4. Release the cutting flywheel by moving the locking lever up, page 3
5. Engage the flywheel motion by moving the lever on the flywheel clockwise (cw), page 3
6. Engage knife platform motion by rotating the platform lock lever clockwise (significant force is required), page 3
7. **Cancel or stop using the electronic buttons ONLY. Never attempt to hold or grab any component of the machine to stop motion, page 5.**

# Microtome; Mechanical Controls, Major

Flywheel; Specimen Motion (up, down, lock, unlock) rotate to move specimen above knife area



Free Flywheel lever ; down to lock, up to allow cutting motion

Speciment Adjust (see next slide)

Knife Adjustments (see next slide)

Knife Mount Lock; CW to clamp knife for motion and slicing

Rate of Cut, 0-200 mm/sec  
(speed of up/down motion) CW to increase.

# Microtome; Mechanical Controls, Knife (K) and Specimen (S) (DO NOT PLUG IN DEVICE AT THIS TIME)

S1. Release to remove entire block if needed


S2. Clamp specimen to block, CW

S3: Specimen side control lock

K2. Release Levers; set angle of knife mount


K3. Set height of knife, screw CCW, then slide bar

K1. Release Levers; set angle of knife mount





# Microtome; Historange with Electronic Controls



***Hotlink to Historange model LKB 2218 Manual.....(pdf)***



# Microtome: Other Protocols or Notes

- Future home of other or more details protocols....