

USER'S MANUAL

MODEL: 3473-50

MODEL: 3473-70

150MM ELECTROMAGNET

Date Sold: _____

Serial number: _____

PROPRIETARY

THIS DOCUMENT CONTAINS CONFIDENTIAL INFORMATION PROPRIETARY TO GMW ASSOCIATES. IT MUST NOT BE REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY WAY EXCEPT FOR THE INSTALLATION, OPERATION OR MAINTENANCE OF GMW ASSOCIATES PRODUCTS.

TABLE OF CONTENTS

SPECIFICATIONS	Section 1
Table 1 Model 3473-50 General Specifications	
Table 1 Model 3473-70 General Specifications	
Table 2 Model 3473-50/3473-70 Electrical and Water Connections	
WARNINGS [Refer to this section before operation of Electromagnet]	Section 2
INSTALLATION	Section 3
Unpacking Instructions	
Mounting Position	
Pole Selection and Installation	
Electrical Circuit	
Interlocks	
Cooling	
OPERATION	Section 4
General	
Calibration	
Field Control Operation	
MAINTENANCE	Section 5
STANDARD OPTIONS	Section 6
Motorized Rotating Drive	
Bias Coils	
Pole Spacer [for 100mm Pole Caps]	
Pole Spacer [for 150mm Pole Caps]	
Probe Holder	
CUSTOM OPTIONS	Section 7
EXCITATION CURVES	Section 8
TEST DATA	Section 9
DRAWINGS	Section 10
Elmwood 3450 Thermostats	
Imo Gems Series Flow Switch	

Continued...

DRAWINGS

Drawing 11801281 3473-50 Electromagnet General Assembly
Drawing 11801282 3473-70 Electromagnet General Assembly
Drawing 11900710 3473-50/P62-4050A Electromagnet Electrical Assembly
Drawing 13900210 3473-50/P62-4050A Electromagnet Electrical Wiring
Drawing 11900700 3473-70/P63-60110A Electromagnet Electrical Assembly
Drawing 13900200 3473-70/P63-60110A Electromagnet Electrical Wiring
Drawing 11900220 3473-70/DF858 Electromagnet Electrical Assembly
Drawing 13900100 3473-70/DF858 Electromagnet Electrical Wiring
Drawing 11900230 3473-50/BOP20-40 Electromagnet Electrical Assembly
Drawing 13900080 3473-50/BOP20-40 Electromagnet Electrical Wiring
Drawing 13900090 3473-50/BOP20-20/20-20 Power Supply Electrical Wiring
Drawing 11900110 Electromagnet Assembly Sequence to Rolling/Rotating Base
Drawing 11803230 Electromagnet Assembly to Rolling/Rotating Base (45° Mtg)
Drawing 11803200 Electromagnet Assembly to Rolling Base (45° Mtg)
Drawing 11803210 Electromagnet Assembly to Rotating Base (45° Mtg)
Drawing 11900090 Electromagnet Assembly to Rolling/Rotating Base (Horz Mtg)
Drawing 11900080 Electromagnet Assembly to Rolling Base (Horz Mtg)
Drawing 11900100 Electromagnet Assembly to Rotating Base (Horz Mtg)
Drawing 11803430 Electromagnet Rolling/Rotating Base Assembly
Drawing 11803170 Electromagnet Rolling Base Assembly
Drawing 11802090 Electromagnet Rotating Base Assembly
Drawing 11803250 Electromagnet Assembly to Vertical Mount
Drawing 11900070 Electromagnet Assembly to Horizontal Mount
Drawing 17803180 Electromagnet Vertical Mount Bracket
Drawing 17900170 Electromagnet Horizontal Mount Bracket
Drawing 17800520 Electromagnet 45 Degree Mount Bracket
Drawing 18900020 Electromagnet Tool Kit
Drawing 17801350 Pole Cap (150, 100, 75, 50, 25mm)
Drawing 18800361 Shipping Crate Assembly
Drawing 18800410 Packing Box Pole Cap Pair

**Section 1
SPECIFICATIONS**

Table 1. Model 3473-50 Specifications

Pole Diameter:	150mm (6 inch)
Pole Gap:	0 - 127mm (0 to 5 inch)
Standard Pole Caps:	150mm (6 inch) cylindrical 100mm (4 inch) tapered 75mm (3 inch) tapered 50mm (2 inch) tapered 25mm (1 inch) tapered
Coils (series connection)	
coil resistance (20°C)	0.72 Ohm
max resistance (hot)*	0.87 Ohm
max power (air)	20A/17V (0.5kW)
max power (water)	50A/44V (2.2kW)
Self Inductance	
Water Cooling (18°C)	3 liters/m (0.8 US gpm) 0.8 bar (12 psid)
Overtemperature Interlock	Elmwood 3450G thermostat part number 3450G 611-1 L50C 89/16 mounted on each coil and wired in series. Contact rating 120Vac,0.5A. Closed below 50°C.
Water Flow Interlock	Imo/Gems flow switch part number FS927 Part No.70823 mounted on outlet side of water circuit. Contact rating 0.17A/120Vac (non inductive). Set to open at a flow of less than 2.5 l/min (0.7 USgpm)
Dimensions	Drawing 11801281 686mm W x 405mm D x 570mm H (27.0 inch W x 16.0 inch D x 22.4 inch H)
Weight	600 kg (1320 lb)

***CAUTION - The value of maximum coil resistance given should not be exceeded. At this resistance the coils are at maximum safe temperature for continuous operation.**

**Section 1
SPECIFICATIONS**

Table 1. Model 3473-70 Specifications

Pole Diameter	150mm (6 inch)
Pole Gap	0 - 96mm (0 to 3.8 inch)
Standard Pole Caps	150mm (6 inch) cylindrical 100mm (4 inch) tapered 75mm (3 inch) tapered 50mm (2 inch) tapered 25mm (1 inch) tapered
Coils (series connection)	
coil resistance (20°C)	0.72 Ohm
max resistance (hot)*	0.87 Ohm
max power (air)	20A/17V (0.5kW)
max power (water)	70A/59V (4.1kW)
Self Inductance	
Water Cooling (18°C)	6 liters/m (1.6 US gpm) 2.0 bar (30 psid)
Overtemperature Interlock	Elmwood 3450G thermostat part number 3450G 611-1 L50C 89/16 mounted on each coil and wired in series. Contact rating 120Vac,0.5A. Closed below 50°C.
Water Flow Interlock	Imo/Gems flow switch part number FS927 Part No.70825 mounted on outlet side of water circuit. Contact rating 0.17A/120Vac (non inductive). Set to open at a flow of less than 4 l/min (1.1 USgpm).
Dimensions	Drawing 11801282 686mm W x 405mm D x 570mm H (27.0 inch W x 16.0 inch D x 22.4 inch H)
Weight	610 kg (1340 lb.)

***CAUTION - The value of maximum coil resistance given should not be exceeded. At this resistance the coils are at maximum safe temperature for continuous operation.**

Section 1
SPECIFICATIONS

Table 2. Model 3473-50/3473-70 Electrical and Water Connections

DC Current (as seen from the rear refer to Drawing 11801281/2)

Right hand terminal: Negative
Left hand terminal: Positive

Ground

An M6 screw (Item 40 on drawing 11801281/2) is provided near the Interlock Block connections to enable the magnet frame to be grounded according to local safety regulations. It is normally appropriate to connect the magnet frame to the power supply ground.

Interlocks (refer to Drawing 11801281/2)

1	Water flow
2	Water flow
3	Temperature
4	Temperature
5	No connection
6	No connection
7	Signal ground
8	Spare (No connection)

Water (refer to Drawing 11801281/2)

Outlet: ¼ inch NPT
Inlet: ¼ inch NPT
(mating couplings for ¼ inch hose provided)

CAUTION - Ensure that the high current connections are tight. Loose connections may lead to oxidation and overheating. The field stability may be degraded and the current terminations damaged.

Section 2

WARNINGS

REFER TO WARNINGS BELOW BEFORE OPERATING ELECTROMAGNET

1 Personnel Safety

In operation the magnet fringing field is in excess of 0.5mT (5G). This can cause malfunctioning of heart pacemakers and other medical implants. We recommend that the fringing field should be mapped and warning signs be placed outside the 0.5mT (5G) contour. Entry to this region should be restricted to qualified personnel

2 Ferromagnetic Objects

During operation the magnet exerts strong magnetic attraction towards ferromagnetic objects in the near vicinity of its pole gap or coils. Loose objects can be accelerated to sufficient velocity to cause severe personnel injury or damage to the coils or precision pole faces if struck. Keep ferromagnetic tools clear!

3 Arcing

This magnet stores considerable energy in its field during operation. Do not disconnect any current lead while under load or the magnetic field energy will be discharged across the interruption causing hazardous arcing.

4 Coil Hot Resistance

Do not exceed the maximum coil hot resistance given in the specifications or coil overheating and possible damage may occur.

5 Interlocks

These should *always* be connected if the magnet is operated unattended, to avoid the possibility of coil overheating caused by excessive power dissipation or inadequate cooling.

6 Watches, Credit Cards, and Magnetic Disks

Do not move magnetically sensitive items into the close vicinity of the magnet. Even some anti-magnetic watches can be damaged when placed in close proximity to the pole gaps during operation. Credit cards, and magnetic disks are affected by magnetic fields as low as 0.5mT (5G). Depending on the previous operating field and the pole gap, the remanent field in the gap can be in excess of 50G (5mT) with the magnet power supply off or disconnected.

Section 3

INSTALLATION

Caution: This is a heavy system. All movement, lifting and installation of the 3473 Electromagnet must be under the supervision of an experienced person to prevent the possibility of serious injury or damage to the Electromagnet and associated equipment.

Unpacking Instructions and Damage Inspection

To unpack the electromagnet please use the following procedure (Refer to Drawing 18800361).

1. First remove all of the "Hex Head Screws" located at the lower edge of all the side panels of the "Crate Top Cover".
2. Gently rock the "Crate Top Cover" to work it loose from the shipping crate base.
3. Use one person on each side of the shipping crate grip the side panels of the Crate Top Cover. Lift "Crate Top Cover" high enough to clear top of electromagnet, walk cover sideways to clear area and place on floor.
4. Inspect the magnet to ensure that no damage has occurred to the magnet in shipment. If damage is evident report the damage in detail to the shipper for claim and simultaneously notify GMW in case assessment of the damage must be made. If no damage is found proceed with magnet unpacking and installation.
5. Remove the M16 Hex Bolts that secure the magnet to the steel shipping angle brackets.
6. Remove the hex lag bolts that secure the steel "shipping angle brackets" to shipping crate base, and remove shipping angle brackets.
7. Install M16 lifting eye and washer to top of magnet yoke, screw down firmly.
8. The magnet is now prepared for final installation, follow the appropriate following procedure to install to 45°, vertical or direct mounting.

Direct Mounting

1. With suitable lifting equipment (e.g. 900kg (2000 lb.) minimum safe lifting rating), lift magnet 50mm (2") clear of shipping crate base.
2. Slide shipping crate base clear.
3. Lower magnet to 50mm (2") above floor.
4. Move magnet to final location and secure using the steel shipping angle brackets. The brackets can be modified to suit installation space needs.
- 5.

Rolling or Rolling/Rotating Base Mounting (refer to Drawing 11900110)

Caution do not attempt to move magnet and rolling base or rolling/rotating base until the magnet has been firmly bolted down to the base (refer to figure 6).

1. To mount on rolling base or rolling/rotating base lift magnet from BOTH FRONT EYEBOLTS high enough to clear top of base (refer to figure 5).
2. Slide rolling base or rolling/rotating base underneath, lower magnet to 12mm (0.5") above base top surface (refer to figure 5).
3. Position rolling base or rolling/rotating base so the tapped hole in the base are aligned with the 45° mounting bracket hole (refer to figure 5).
4. Lower magnet onto rolling base or rolling/rotating base assembly (refer to figure 5).
5. Secure magnet and 45° mounting assembly to rolling base or rolling/rotating base with M16 x 25 long Hex Head Bolts (refer to figure 6).
6. Move magnet and rolling base or rolling/rotating base to desired location.

Continued...

Section 3

INSTALLATION

Rolling or Rolling/Rotating Base Mounting (continued)

7. Screw down the four support legs located on each corner of the rolling or rolling/rotating base until the wheels clear the floor by 6mm (.25").
8. Secure the support legs with the locknut.
9. Secure rolling/rotating base to an adequate concrete floor to prevent movement and possible injury to personnel during an earthquake.

Pole Cap Selection and Installation (Refer to drawing 11801291/2)

Using the field uniformity and induction curves determine the most desirable pole cap; cylindrical or tapered. in general:

If a uniform field is required use a cylindrical cap.

If a high field is required use a tapered cap.

Pole cap removal (refer to drawing 11801281/2)

1. Turn off the power supply
2. Draw pole caps about 20mm into the pole sleeves.
3. Loosen the axial draw stud nut (item 35 on drawing 11801281/2).
4. Insert the hex key wrench into the end of the draw stud (item 6 on drawing 11801281/2).
5. Remove draw stud (item 6 on drawing 118801281/2) while supporting the pole cap.

Pole cap fitting.

1. Ensure the pole caps, pole cores, and pole sleeves are clean and free from debris.
2. Reverse the above pole cap removal sequence.

Electrical Circuit

Never connect or remove cables from the magnet with the power supply connected. The stored energy in the magnet can cause arcing resulting in severe injury equipment damage.

The magnet has two coils which are connected in series, (11801281/2). Refer to drawing. The power supply cables should be connected directly to the dc current terminals marked + and -. Recommended current cable for the 3473-50 is stranded copper of 16mm² cross section (4 AWG). For the 3473-50 the cable size should be increased to 25mm² cross section (3 AWG).

Because the magnet stores a significant amount of energy in its magnetic field, special care should be taken to insure that the current terminations are secure and cannot work loose in operation. Local heating at the terminations can cause rapid oxidation leading to a high contact resistance and high power dissipation at the terminals. If left unattended this can cause enough local heating to damage the terminals and the coils.

Section 3

INSTALLATION

The 3473 Interlocks

The Model 3473-50 uses two thermostats, Elmwood 3450G Part Number 3450G611-1 L50C 89/16. They are wired in series and terminated in positions 3 and 4 on the Interlock Terminal block. The thermostats are normally closed, opening when the coil central cooling plate temperature exceeds $50^{\circ}\text{C} +/3^{\circ}\text{C}$. The 3473-70 uses six thermostats. The flow switch is connected to terminals 1 and 2. The contacts are normally open, closing when the water flow exceeds approximate 2.5l/min. for the 3473-50 and 4.0l/min for the 3473-70.

Cooling

The Model 3473 can be operated to an average coil temperature of 70°C . Assuming an ambient environment temperature of 20°C and a temperature coefficient of resistivity for copper of $0.0039/^{\circ}\text{C}$, the hot resistance of the coil should not exceed 20% more than the ambient temperature "cold" resistance. The coil thermostats will open when any coil cooling plate temperature exceeds approximately 50°C . Clean, cool ($16^{\circ}\text{C} - 20^{\circ}\text{C}$) water at 3 l/min at 0.8 bar (12 psid) should be used to cool the 3473-50 magnet, and clean, cool ($16^{\circ}\text{C} - 20^{\circ}\text{C}$) water at 6 l/min at 2.0 bar (30 psid) for the 3473-70.

The cooling copper tubes are electrically isolated from the coils to avoid electrochemical corrosion. A 50 micron filter should be placed before the input to the magnet to trap particulates and avoid unreliable operation of the water flow switch interlock.

For continuous operation of the magnet it may be appropriate to use a recirculating chiller to reduce water and drainage costs. The chiller capacity will depend on whether cooling is required for the magnet alone or magnet and power supply. For the Model 3473-50 Electromagnet alone a suitable chiller is the Bay Voltex Model: RRS-0850 for the Model 3473-70 alone use the Bay Voltex Model: RRS-1650. Use distilled or deionized water with a biocide to prevent bacterial growth and corrosion. Do not use corrosion inhibitors in high quality electrical systems since the water conductivity is increased which can result in increased leakage currents and electrochemical corrosion.

At currents of approximately 20A and below the Model 3473 can be operated safely without water cooling. However the coil temperature will vary with the power dissipation. This results in dimensional and permeability changes of the magnet yoke and air cooling is not suitable when high field stability is required.

Freon, oil, ethylene glycol or other cooling mediums can be used. The flow required will be approximately inversely proportional to their specific heats. An experimental determination of the flow and pressure required will be necessary.

Avoid cooling the magnet below the dew point of the ambient air. Condensation may cause electrical shorts and corrosion.

During operation the resistance can be checked using a voltmeter across each coil. The voltage will rise to a constant value once thermal equilibrium has been reached. If it is desired to save water, the flow can be reduced until the hot resistance is approached. NOTE: This adjustment must be made slowly enough to allow for the thermal inertia of the coils.

Section 4

OPERATION

General

The magnet operates as a conventional electromagnet.

1. Adjust the poles to the desired gap with the poles approximately symmetrical about the center magnet line. To reduce mechanical backlash when the magnetic field is applied, it is best to set the poles by increasing the gap.
2. Adjust the cooling water flow to about 3 liters/min (0.8 USgpm) for the 3472-50. For the 3473-70 set water flow to about 6 liters/min (1.6 US gpm,). For operation at less than maximum power the water flow may be correspondingly reduced. Note that the inlet water temperature will determine the actual flow rate required. The above specified figures were determined with a water inlet temperature of $<18^{\circ}\text{C}$.
3. Turn on the power supply and increase the current until the desired field is reached.

Calibration

The induction curves may be used to estimate the field in the air gap to within four or five percent. More accurate field determination may be obtained by deriving experimentally a calibration curve for the particular pole and air gap combination being used. Magnetic hysteresis in the yoke and poles can cause an error of 30 to 70G (3 to 7mT) with an arbitrary application of such a calibration curve. This effect may be reduced to less than one percent by following a prescribed 'current setting schedule' designed to make the magnet 'forget' its prior magnetic history. The schedule should of course be used both in establishing the calibration curve and in its subsequent use. A possible schedule would be:

From zero current, increase to maximum current and reduce again to zero current. Increase again to maximum current and reduce to the current to give the desired field setting. Approaching the desired field from a higher setting will typically produce better field uniformity. This is because the field changes at the pole edges will normally lag the field change at the center thereby helping to compensate the radial decrease in field.

Greater precision in setting up the calibration curve will be achieved with the use of a digital gaussmeter and by making a numerical table. This table used with an interpolation routine will eliminate the error associated with reading a graph.

In any event, three points need to be remembered:

1. A calibration curve or table is only as good as the precision employed in generating it.
2. The field is defined only at the point it is measured. It will generally be different at a different point in the air gap. For example, the induction curves refer to the field on the pole axis and at the center of the air gap (median plane).

Calibration - continued

3. The field is most directly a function of the current in the magnet coils. Voltage across the coils is not a good measure of field since the electrical resistance of the coils depends on the temperature (about 0.4% per degree celsius).

Section 4

OPERATION

Field Control Operation

The necessity to use calibration curves can be avoided by using a field controller to sense the magnetic field and provide a corresponding power supply control signal through the power supply programming inputs. Contact GMW for suitable instrumentation.

Section 5

MAINTENANCE

Periodically check that the pole adjustment mechanism is clean, properly lubricated and free of grit and dirt, which may cause binding of the mechanism. Otherwise no particular maintenance is required. Be very careful not to damage the relatively soft pole surface since this may degrade the magnetic field uniformity in the gap.

Note that the surface treatments used provide good corrosion protection but in order to maintain the inherent mechanical precision of the magnet, heavy build-up of plating materials is deliberately avoided. As a result, high humidity or otherwise seriously corrosive atmospheres can defeat the protection mechanisms. Check the equipment periodically and use an appropriate corrosion protection when the magnet is stored for an extended period.


Section 6

STANDARD OPTIONS

PROHIBITARY
 THIS DRAWING IS THE PROPERTY OF GMW. IT IS TO BE USED ONLY FOR THE PURPOSES SPECIFIED HEREIN. IT IS NOT TO BE REPRODUCED, COPIED, OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE WRITTEN PERMISSION OF GMW.


TOP VIEW

NOTE: ROTATING BASE SHOWN AT THE 180° POSITION




FRONT VIEW

NOTE: ROTATING BASE SHOWN AT THE 180° POSITION



UNDERNEATH VIEW



STEPPER MOTOR ELECTRICAL CONNECTIONS

SIDE VIEW

MECHANICAL STOPS

HOMING MICROSWITCH CO DEEG POSITION

DISENGAUGE LEVER RELEASES MOTOR DRIVE

RELEASE PINS USED TO REMOVE MOTOR DRIVE ASSY FROM ROLLING/ROTATING BASE ASSY

55 (2.11) SPOOL CLEARANCE

135 (5.31) SERVICE ACCESS

NOTE: ROTATING BASE SHOWN AT THE 180° POSITION

85 (3.00)

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	02/11/01	ENGORZAK


ITEM	QTY	DESCRIPTION	UNIT
13	1	GROMMET 25MM OD X 20MM D	
12	2	SKIRT-ION	
11	1	DN 4.33	
10	5	BN-7.92	
9	6	BN 7.92	
8	1	DN 9.12	
7	6	DN 9.12	
6	1	17.901.230	
5	1	109001.01	
4	2	17.901.020	
3	1	11.900.820	
2	1	11.900.811	
1	1	11.900.800	

DO NOT SCALE
 DIMENSIONS OF THIS DRAWING ARE TO BE USED FOR MANUFACTURE AND ASSEMBLY.
 GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (650)992-8292 Fax: (650)992-8298
 MOTORIZED ROT. DRIVE
 3473/3472/5403
 A1 11900800
 SCALE 1:12 (1:1)

PROVISIONAL
 THIS DRAWING IS THE PROPERTY OF GMW. IT IS TO BE USED ONLY FOR THE PROJECT AND SPECIFICATIONS STATED HEREON. IT IS NOT TO BE REPRODUCED, COPIED, OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE WRITTEN PERMISSION OF GMW.

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	07/07/01	CS0025.AS
B	ADD DIA 17.818, MTD SPRING SHAPE	11/29/01	BL0005.AS


SCALE 4:1



DETAIL "B"



BOTTOM VIEW



WORN GEAR
 BALL PLUNGERS USED FOR OVERLOAD CLUTCH

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
18	1	169004770	LABEL, SPOOL IDENTIFICATION	
17	2	BN 782	WASHER, FLAT M6 X 1.6 S/S	
16	1	DN 1587	NUT, HEX BOWED M8 X 12 S/S	
15	4	BN 782	WASHER, LOCK M8 X 1.2	
14	5	BN 792	WASHER, LOCK M8 X 1.4	
13	1	BN80	BOLT, M6 X 16 HEX N/UN	
12	5	DN 912	SHCS M8 X 16 S/S	
11	4	DN 912	SHCS M8 X 25 S/S	
10	2	DN 913 A2	SHSS M4 X 10 S/S	
9	2	DN 913 A2	SHSS M4 X 12 S/S	
8	4	R 1.0MM	PINNO WIRE 5/5	
7	1	DN 5325	DOWEL PIN, M6 X 80	
6	1	17901090	STOP BAR	
5	1	17901100	STOP BAR GLIDE	
4	6	SRMH-1DN	BALL PLUNGER, VUER 5/5	
3	1	17901130	SPOOL CLAMP RING	
2	1	12900220	WORN GEAR, BRONZE	
1	1	11900830	SPOOL WELDED ASSY	

DO NOT SCALE
 DIMENSIONS
 FINISHES
 MATERIALS
 TOLERANCES
 UNLESS OTHERWISE SPECIFIED

GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (650)902-8292 Fax: (650)902-8298


MOTORIZED ROT. DRIVE
 SPOOL ASSEMBLY

REV A1 11900820

SCALE 1:2 (M L)

SHEET 1 OF 1

PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 PROPRIETARY TO GMW, INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY
 MANNER THAT VIOLATES THE PATENT RIGHTS OR TRADE
 SECRETS OF GMW, INC. OR IN VIOLATION OF ANY
 APPLICABLE LAWS OR REGULATIONS.



BOTTOM VIEW

REAR VIEW

TOP VIEW

REVISIONS			
REV	DESCRIPTION	DRAWN	DATE
A	RELEASE		07/07/97
		G.DOUGLAS	

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
10	6	BN 752	WASHER, LOCK SP/S M2 X 0.5 SP/S	
9	6	BN 792	WASHER, LOCK SP/S M3 X 0.9 SP/S	
8	6	DN 912	BOLT, SHCS M2 X 10 S/S	
7	6	DN 912	BOLT, SHCS M3 X 10 S/S	
6	2	V4NT7	MICROSWTCH, BURGESS	
5	1	V4NT9	SHAFT, ZERO MICROSWTCH	
4	1	17901170	BRACKET, LIMIT MICROSWTCH	
3	2	17901160	BRACKET, LIMIT MICROSWTCH	
2	1	17901150	BRACKET, ZERO MICROSWTCH	
1	1	17901470	STOP BLOCK	

GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (650)902-8292, Fax: (650)902-8298

**MOTORIZED, ROT. DRIVE
 STOP BLOCK ASSY**

DATE: 05/02/97
 DRAWN: G.DOUGLAS
 CHECK: G.DOUGLAS
 ENGINEERING: G.DOUGLAS

DO NOT SCALE FROM DRAWING
 DIMENSIONS & TOLERANCES
 UNLESS OTHERWISE SPECIFIED

LINEAR	METERS	MM
ANGULAR	DEGREES	DEC. MIN.
FINISH	RA	1.6


THIRD ANGLE PROJECTION

SOFTWARE: AUTOCAD 13

SCALE: 1:1 WT kg SHEET 1 OF 1

PROVISIONAL
 THIS DRAWING IS THE PROPERTY OF GMW ELECTRIC INC. IT IS TO BE USED ONLY FOR THE PROJECT AND QUANTITY SPECIFIED THEREIN. IT IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE WRITTEN PERMISSION OF GMW ELECTRIC INC.


TOP VIEW




MOTOR AND LIMIT SWITCH CABLE PART NO 16900400

SERIAL COMMUNICATIONS CABLE PART NO 16900410

AC POWER INPUT



STEPPER MOTOR CONTROLLER



CONTROL COMPUTER

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	02/17/91	EDWARDS

REVISIONS

ITEM	QTY	PART NUMBER	DESCRIPTION	UNIT
1				
2				
3				
4				
5				
6				
7				
8				

DO NOT SCALE	DATE	1/11/91
BY	EDWARDS	
CHECKED	EDWARDS	
APPROVED	EDWARDS	
DATE	1/11/91	
PROJECT	GMW	
DESCRIPTION	MOTORIZED ROT. DRIVE ELECTRICAL ASSY	
REV	A1	11901020
REV	A	
SOFTWARE	AUTOCAD 13	
AUTHOR		
SCALE	1:2	
SHEET	1	OF 1

A

B

C

D




8 7 6 5 4 3 2 1

PROHIBITED: THE INFORMATION CONTAINED HEREIN IS UNCLASSIFIED EXCEPT WHERE SHOWN OTHERWISE. IT IS THE PROPERTY OF GMW AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM.

REV	DESCRIPTION	SHEET	DATE	APPROVED
A	RELEASE	1	08/11/21	BOB/DAK
B	ADD CONTROLLER CONNECTOR LINKING	2	10/19/21	BOB/DAK

MOTOR DRIVE ASSEMBLY

PART NO: 11900800/1



- NOTE:
- 1 USE BLACK WIRES FROM GREEN AND WHITE PAIRS. DO NOT MIX WITH BLACK WIRES USED FOR STEPPER MOTOR.
 2. HOMING SWITCH USES NORMALLY OPEN CONTACT. LIMIT SWITCHES USE NORMALLY CLOSED CONTACTS.


ITEM	QTY	PART NUMBER	DESCRIPTION	UNIT
1	1	11900800	MOTOR DRIVE ASSEMBLY	ASSEMBLY
2	1	16904400	NOTER AND LIMIT SWITCH CABLE	ASSEMBLY
3	1	10000000	CONTROLLER 10 WAY CONNECTOR	CONNECTOR
4	1	90000000	CONTROLLER 9 WAY CONNECTOR	CONNECTOR

GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (650)902-6292 Fax: (650)902-6298


MOTORIZED, ROT DRIVE
 ELECTRICAL WIRING

REV A1 13900350

SCALE 2:1 1st Ed




PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 PROPRIETARY TO GDMV, INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY
 MANNER WITHOUT THE WRITTEN PERMISSION OF GDMV, INC.



NOTE: THIS DRAWING SHOWS REAR INSTALLATION OF PROBE MOUNT
 FOR FRONT INSTALLATION OF PROBE MOUNT SEE DRAWING NO 11901262


SIDE VIEW



HUB ANGLE ADJUSTABLE
 IN 15° INCREMENTS
 FROM -45° to +45°

MPT PROBEHOLDER
 SEPARATE ORDER ITEM

- 1 LOOSEN THUMB NUT 2mm
- 2 PULL HUB FORWARD 2mm
- 3 ROTATE TO ANGLE REQUIRED
- 4 ROTATE SLIGHTLY BACK AND FORTH TO FIND INDEX PIN
- 5 PUSH HUB REARWARDS
- 6 TIGHTEN THUMB NUTS



REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	08/27/98	G.DOUGLAS


ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
21	1	10900320	LABEL, IDENTIFICATION	
20	1	1SBMHB	BALL PLUNGER, M8 S/S VLER	
19	2	VSM 12771B	DOWEL PIN M1 X 5 S/S [Index Pin]	
18	1	BN 1073	SET SCREW, M6 X 5 SLOTTED HD NYLON	
17	4	ISD 738Q	SHCS M4 X 6 BUTTON HD S/S	
16	5	DIN 7991	SHCS, M4 X 6 FLAT HEAD S/S	
15	2	DIN 917	SHSS M4 X 8 CONE POINT S/S	
14	2	Q8M04D070TN	THUMB NUT, NYLON	
13	3	18-830	ITEM PRODUCTS, END CAP, PLASTIC	
12	1	17902010	BASE STUD	
11	1	17902000	HUB STUD	
10	1	17901990	HUB INSERT [for Sentron Hall Probes]	
9	1	17901980	HUB INSERT [for Metrolob NMR probes]	
8	1	17901970	HUB INSERT [for Grp3 MPT Hall Probes]	
7	1	17901960	HUB COVER	
6	1	17901950	HUB BASE	
5	1	17901944	VERTICAL MOUNTING EXTRUSION	
4	1	17901930	BASE NUT	
3	1	17901920	BASE SUPPORT	
2	1	17902050	BASE MOUNTING EXTRUSION	
1	1	17902040	BASE MOUNTING PLATE	

DATE	SCALE	DESCRIPTION	REVISION
08/27/98	1:1	PROBE MOUNT MODEL: 3473	A

GMV
 955 Industrial Rd, San Carlos, CA 94070
 Tel: (650)902-8292. Fax: (650)902-8298.




PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 PROPRIETARY TO GDMV INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED IN ANY MANNER
 WITHOUT THE WRITTEN CONSENT OF GDMV INC.



NOTE: THIS DRAWING SHOWS FRONT INSTALLATION OF PROBE MOUNT
 FOR REAR INSTALLATION OF PROBE MOUNT SEE DRAWING NO 11901261

SIDE VIEW



- HUB ANGLE ADJUSTABLE
 IN 15° INCREMENTS
 FROM -45° to +45°
- 1 LOOSEN THUMB NUT 2mm
 - 2 PULL HUB FORWARD
 - 3 ROTATE TO ANGLE REQUIRED
 - 4 ROTATE SLIGHTLY BACK AND FORTH TO FIND INDEX PIN
 - 5 PUSH HUB REARWARDS
 - 6 TIGHTEN THUMB NUTS

MPT PROBEHOLDER
 SEPARATE ORDER ITEM

MPT HALL PROBE
 SEPARATE ORDER ITEM



REV	DESCRIPTION	DRAFT	DATE	APPROVED
A	RELEASE		08/27/98	G.DOUGLAS

REV	DESCRIPTION	DRAFT	DATE	APPROVED
21	1 10900320			
20	1 SBMHB			
19	2 VSM 12771B			
18	1 BN 1073			
17	4 ISO 7380			
16	5 DIN 7991			
15	2 DIN 917			
14	2 Q8M04D070TN			
13	3 18-830			
12	1 17902010			
11	1 17902000			
10	1 17901990			
9	1 17901980			
8	1 17901970			
7	1 17901960			
6	1 17901950			
5	1 17901944			
4	1 17901930			
3	1 17901920			
2	1 17902030			
1	1 17902020			

GMW
 955 Industrial Rd, San Carlos, CA 94070
 Tel: (650)902-8292. Fax: (650)902-8298.

PROBE MOUNT
 MODEL: 3473

DATE: 08/27/98
 DRAWN BY: G.DOUGLAS
 CHECKED BY: [blank]
 ENGINEERING DATE: [blank]

DO NOT SCALE FROM DRAWING
 DIMENSIONS & TOLERANCES
 UNLESS OTHERWISE SPECIFIED

LINE	UNIT	CONV	CONV
X, X.X, X.XX	IN	MM	25.4
4, 4.0, 4.00	IN	MM	25.4
5, 5.0, 5.00	IN	MM	25.4
6, 6.0, 6.00	IN	MM	25.4
7, 7.0, 7.00	IN	MM	25.4
8, 8.0, 8.00	IN	MM	25.4
9, 9.0, 9.00	IN	MM	25.4
10, 10.0, 10.00	IN	MM	25.4
11, 11.0, 11.00	IN	MM	25.4
12, 12.0, 12.00	IN	MM	25.4
13, 13.0, 13.00	IN	MM	25.4
14, 14.0, 14.00	IN	MM	25.4
15, 15.0, 15.00	IN	MM	25.4
16, 16.0, 16.00	IN	MM	25.4
17, 17.0, 17.00	IN	MM	25.4
18, 18.0, 18.00	IN	MM	25.4
19, 19.0, 19.00	IN	MM	25.4
20, 20.0, 20.00	IN	MM	25.4
21, 21.0, 21.00	IN	MM	25.4

THIRD ANGLE PROJECTION

SCALE 1:1 WT kg SHEET 1 OF 1

Section 7


CUSTOM OPTIONS

Section 8

EXCITATION CURVES

GMW Associates
Electromagnet Excitation Plot
Field Vs Gap

Contract No:	Page: 1 of 1	Date: May 05, 1994
Customer:		Engr: R Yass
Model: 3473-70	Power Supply: D/F 854 100-100	Set Current: 70 Amps
Serial No: 22	Serial No: 9101033	Target Field:
Pole Face: As per table below	Position: X=0, Y=0, Z=0	
Serial No: None	Notes:	
Pole Gap: As per table below		
Pole Spacers: None		




GMW Associates

Electromagnet Excitation Plot

Field Vs Current

Contract No:	Page: 1 of 5	Date: May 05, 94
Customer:		Engr: R Yass
Model: 3473-70	Power Supply: D/F 854 100-100	Set Current:
Serial No: 22	Serial No: 9101033	Target Field:
Pole Face: 150	Position: X=0, Y=0, Z=0	
Serial No: None	Notes:	
Pole Gap: As per table below		
Pole Space: None		



GMW Associates Electromagnet Excitation Plot Field Vs Current

Contract No:

Page: 2 of 5

Date: May 05, 94

Customer:

Engr: R Yass

Model: 3473-70

Power Supply: D/F 854 100-100

Set Current:

Serial No: 22

Serial No: 9101033

Target Field:

Pole Face: 100


Position: X=0, Y=0, Z=0

Serial No: None

Notes:

Pole Gap: As per table below

Pole Spacers: None




GMW Associates

Electromagnet Excitation Plot


Field Vs Current

Contract No:	Page: 3 of 5	Date: May 05, 94
Customer:		Engr: R Yass
Model: 3473-70	Power Supply: D/F 854 100-100	Set Current:
Serial No: 22	Serial No: 9101033	Target Field:
Pole Face: 75	Position: X=0, Y=0, Z=0	
Serial No: None	Notes:	
Pole Gap: As per table below		
Pole Spacers: None		



GMW Associates
Electromagnet Excitation Plot
Field Vs Current

Contract No:	Page: 4 of 5	Date:	May 05, 94
Customer:		Engr:	R Yass
Model: 3473-70	Power Supply: D/F 854 100-100	Set Current:	
Serial No: 22	Serial No: 9101033	Target Field:	
Pole Face: 50	Position: X=0, Y=0, Z=0		
Serial No: None	Notes:		
Pole Gap: As per table below			
Pole Spacers: None			




GMW Associates

Electromagnet Excitation Plot


Field Vs Current

Contract No:	Page: 5 of 5	Date: May 05, 94
Customer:		Engr: R Yass
Model: 3473-70	Power Supply: D/F 854 100-100	Set Current:
Serial No: 22	Serial No: 9101033	Target Field:
Pole Face: 25	Position: X=0, Y=0, Z=0	
Serial No: None	Notes:	
Pole Gap: As per table below		
Pole Spacers: None		



Section 9

TEST DATA



PROPRIETARY

This document contains information proprietary to GMW Inc. It must not be reproduced or disclosed to others or used in any other way, in whole or part, except as authorized in writing by GMW Inc.


MAGNETIC PLOTTING AXIS

80900020

A

SHEET 1 OF 1

1 SHIM SHOWN FITTED TO POLE AT 90 DEG POSITION



LH POLE: CAP REMOVED

RH POLE: CAP REMOVED

LH POLE SHIM DETAILS		
NUMBER	THICKNESS	POSITION
1	_____mm	_____deg
2	_____mm	_____deg
3	_____mm	_____deg
4	_____mm	_____deg

RH POLE SHIM DETAILS		
NUMBER	THICKNESS	POSITION
1	_____mm	_____deg
2	_____mm	_____deg
3	_____mm	_____deg
4	_____mm	_____deg

MAGNET MODEL: _____

DATA LOGGED BY: _____

MAGNET SERIAL NO: _____

DATA LOGGED DATE: _____


PROPRIETARY
This document contains information proprietary to GMW Inc. It must not be reproduced or disclosed to others or used in any other way, in whole or part, except as authorized in writing by GMW Inc.

MAGNET FINAL SHIMMING LOG

3474-0001

C


SHEET 1 OF 1



1. THOROUGHLY CLEAN AND DEGREASE AREA WHERE SHIM IS TO BE FITTED.
2. CUT SHIM STOCK TO DIMENSIONS SHOWN.
3. APPLY DOUBLE SIDED TAPE 0.1mm THICK TO AREA SHOWN.
4. FIT SHIM TO POLE FACE, ENSURE TAPE IS KEPT WITHIN UNDERCUT AREA.
5. REASSEMBLE POLE CAPS ONTO MAGNET.
6. REMAP MAGNET, IF RESULTS WITHIN SPECIFICATION THEN GO TO ITEM 7. IF OUTSIDE SPECIFICATION ADJUST SHIMS, REMAP THEN GO TO ITEM 7.
7. FILL IN SHIMMING DETAILS ON SHEET NO 3474-0001.

PROPRIETARY <small>This document contains information proprietary to GMW Inc. It must not be reproduced or disclosed to others or used in any other way, in whole or part, except as authorized in writing by GMW Inc.</small>	FINAL SHIMMING METHOD		
3474-0002	A	SHEET 1 OF 1	

TAKE MEASUREMENTS OF GAP TAPER AT POINTS MARKED WITH BORE GAUGE. RECORD RESULTS BELOW



MAGNET GAP TAPER LOG			
X + 70		Y + 70	
X - 70		Y - 70	
X DIFF		Y DIFF	

MAGNET MODEL: _____

DATA LOGGED BY: _____

MAGNET SERIAL NO: _____

DATA LOGGED DATE: _____

PROPRIETARY
This document contains information proprietary to GMW Inc. It must not be reproduced or disclosed to others or used in any other way, in whole or part, except as authorized in writing by GMW Inc.

MAGNET GAP TAPER TEST

80900120

A

SHEET 1 OF 1

GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PL



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	16	Pole Gap 19 mm	Date Oct 13, 1992
Coil Set	70A Sn 1469&1470	Pole Shims 0.004 fitted	NMR Signal -650mV

Power Supply 8.2 Amps 082000 ADC 8.2 % Current

Start Time 14:45
 Start Field
 0.5040260

Finish Time 15:10
 Finish Field
 0.5039170

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	0.5040740	0.5039130	0.5037570	0.5036130	0.5035000
10 +	0.5041710	0.5040060	0.5038690	0.5037510	0.5036590
0	0.5043060	0.5041440	0.5040260	0.5039060	0.5038550
10 -	0.5043640	0.5042140	0.5040870	0.5039990	0.5039630
20 -	0.5045000	0.5043600	0.5042610	0.5041920	0.5041540





GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PL

Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	16	Pole Gap 19 mm	Date Oct 13, 1992
Coil Set	70A Sn 1469&1470	Pole Shims 0.004 fitted	NMR Signal -200mV

Power Supply 37.9 Amps 3796000 ADC 37.9 % Current

Start Time 11:35
 Start Field
 1.5047420
 Finish Time 12:15
 Finish Field
 1.5046450

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	n/s	n/s	1.5030600	1.5024200	n/s
10 +	1.5030500	1.5037000	1.5042400	1.5034500	n/s
0	1.5036400	1.5046300	1.5047420	1.5044000	1.5034600
10 -	1.5034100	1.5046500	1.5048250	1.5045230	1.5035300
20 -	n/s	1.5041300	1.5044600	1.5041000	n/s



GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model 3473	Pole Face 150 mm	Engr Greg Douglas
Serial No 16	Pole Gap 19 mm	Date Oct 13, 1992
Coil Set 70A Sn 1469&1470	Pole Shims 0.004 fitted	NMR Signal -650mV

Power Supply 8.2 Amps 082000 ADC 8.2 % Current

Start Time 15:15
 Start Field
 0.5038870

 Finish Time 15:35
 Finish Field
 0.5038540

Plot Z = 0.0					
Y	X (mm)				
	10 -	5 -	0	5 +	10 +
10 +	0.5039080	0.5038400	0.5037730	0.5037140	0.5036580
5 +	0.5039570	0.5038830	0.5038210	0.5037690	0.5037290
0	0.5040180	0.5039480	0.5038870	0.5038500	0.5037930
5 -	0.5040550	0.5039880	0.5039280	0.5038780	0.5038470
10 -	0.5041240	0.5040650	0.5040050	0.5039560	0.5039190



GMW ASSOCIATES LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	16	Pole Gap 19 mm	Date Oct 13, 1992
Coil Set	70A Sn 1469&1470	Pole Shims none	NMR Signal -450mV

Power Supply 8.2 Amps 082000 ADC 8.2 % Current

Start Time 15:50
Start Field
0.5028100

Finish Time 16:05
Finish Field
0.5027600

Plot Z = 0.0					
Y	X (mm)				
	10 -	5 -	0	5 +	10 +
10 +	0.5029950	0.5028690	0.5027380	0.5026100	0.5024780
5 +	0.5030140	0.5028780	0.5027500	0.5026260	0.5025150
0	0.5030570	0.5029260	0.5028100	0.5026690	0.5025420
5 -	0.5030580	0.5029230	0.5027930	0.5026640	0.5025500
10 -	0.5030950	0.5029740	0.5028420	0.5027110	0.5025840



GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model 3473	Pole Face 150 mm	Engr Greg Douglas
Serial No 16	Pole Gap 19 mm	Date Oct 16, 1992
Coil Set 50A Sn 654 & 655	Pole Shims none	NMR Signal -400mV

Power Supply 8.3 Amps 083000 ADC 8.3% Current

Start Time 13:30
 Start Field
 0.5047420

 Finish Time 13:55
 Finish Field
 0.5047240

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	0.5055160	0.5051560	0.5048040	0.5044870	0.5041900
10 +	0.5054800	0.5050930	0.5047410	0.5044160	0.5041430
0	0.5054890	0.5051020	0.5047420	0.5044200	0.5041160
10 -	0.5055200	0.5051300	0.5047730	0.5044360	0.5041290
20 -	0.5055450	0.5051880	0.5048200	0.5044670	0.5041300



GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PL



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	16	Pole Gap 100 mm	Date Oct 15, 1992
Coil Set	50A Sn 654 & 655	Pole Shims none	Mapped with DTM-141

Power Supply 48.1 Amps 481300 ADC 48 % Current

Start Time 15:05
 Start Field
 0.53155

Finish Time 16:20
 Finish Field
 0.53155

Plot Y = 0.0					
Z	X (mm)				
	40 -	20 -	0	20 +	40 +
40 +			0.54814		
20 +			0.53776		
0	0.51588	0.52867	0.53155	0.52671	0.51005
20 -			0.53797		
40 -			0.54859		



GMW ASSOCIATES LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	16	Pole Gap 19 mm	Date Oct 11, 1992
Coil Set	70A Sn 1469 & 1470	Pole Shims none	NMR Signal -120mV

Power Supply 37.9 Amps 3796000 ADC 37.9 % Current

Start Time 17:50
Start Field
1.5037360

Finish Time 18:25
Finish Field
1.5037850

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	N/S	1.5025020	1.5019890	1.5008900	N/S
10 +	1.5031490	1.5037200	1.5035380	1.5024180	1.5002480
0	1.5032360	1.5040800	1.5037360	1.5024040	1.5005940
10 -	1.5038260	1.5040260	1.5036990	1.5026320	1.5010500
20 -	1.5028900	1.5031800	1.5031150	1.5019500	1.5007600



GMW ASSOCIATES

LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model	3473	Pole Face	150 mm	Engr	Greg Douglas
Serial No	16	Pole Gap	19 mm	Date	Oct 11, 1992
Coil Set	70A Sn 1469 & 1470	Pole Shims	none	NMR Signal	-450mV

Power Supply 8.2 Amps 082000 ADC 8.2 % Current

Start Time 17:10
Start Field
0.5033240

Finish Time 17:35
Finish Field
0.5033110

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	0.5038730	0.5035830	0.5032830	0.5029760	0.5026760
10 +	0.5038790	0.5035740	0.5033020	0.5030210	0.5027390
0	0.5039180	0.5036030	0.5033240	0.5030510	0.5028040
10 -	0.5039550	0.5036430	0.5033500	0.5030850	0.5028490
20 -	0.5040320	0.5037320	0.5034460	0.5031790	0.5029250



GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PL



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	18	Pole Gap 130 mm	Date Oct 22, 1992
Coil Set	50A Sn 1989 & 1988	Pole Shims none	Mapped with DTM-141

Power Supply 34.0 Amps 340546 ADC 34 % Current

Start Time 19:45
 Start Field
 0.28940

Finish Time 21:15
 Finish Field
 0.28940

Plot Y = 0.0					
Z	X (mm)				
	40 -	20 -	0	20 +	40 +
40 +	0.30815	0.30590	0.30498	0.30585	0.30864
20 +	0.28558	0.29238	0.29404	0.29223	0.28569
0	0.27730	0.28666	0.28940	0.28662	0.27683
20 -	0.28667	0.29304	0.29436	0.29259	0.28598
40 -	0.30934	0.30633	0.30541	0.30644	0.30918



GMW ASSOCIATES
LABORATORY ELECTROMAGNET FIELD UNIFORMITY PLOT



Model	3473	Pole Face 150 mm	Engr Greg Douglas
Serial No	18	Pole Gap 19 mm	Date Oct 22, 1992
Coil Set	50A Sn 1989 & 1988	Pole Shims none	NMR Signal -580mV

Power Supply 8.3 Amps 083000 ADC 8.3% Current

Start Time 20:50
 Start Field
 0.5035620

Finish Time 21:15
 Finish Field
 0.5035460

Plot Z = 0.0					
Y	X (mm)				
	20 -	10 -	0	10 +	20 +
20 +	0.5041230	0.5038590	0.5036100	0.5034100	0.5032640
10 +	0.5040940	0.5037800	0.5035480	0.5033540	0.5032270
0	0.5041090	0.5038000	0.5035620	0.5033640	0.5032230
10 -	0.5041470	0.5038500	0.5036000	0.5034120	0.5032600
20 -	0.5042370	0.5039750	0.5037280	0.5035320	0.5033600



Section 10

DRAWINGS

SERIES 3450/3450R/3455R/3455RBV 15 AMP THERMOSTATS


Typical Applications:

Power Supplies

Communication
Equipment

Medical Equipment

Computers (Where
High AMP Loads are
Present)




The Series 3450/3455R is a snap-acting, non-adjustable precision thermostat especially suited for industrial and electrical equipment.

The 3450 (.390" or 10mm overall) is ideal for applications that require precision control of high electric loads to 8 Amp resistive.

The 3450R and 3455R have a patented metal insert rivet construction.

The 3455R (.484" or 12.5mm) overall, has higher spacing as required by European approval agencies. Model 3455RBV is an epoxy overmold version of the 3455R, specifically designed for electrical insulation or protection in a high humidity environment. Consult factory for performance qualifications.

To insure that a safe combination of thermostat and application is achieved, the purchaser must determine product suitability for their individual requirements.



*Series 3450/3450R/3455R/3455RBV

MODEL	ELECTRIC LIFE CYCLES	120 VAC	240 VAC	277VAC
3450	100,000	8.0A	-	-
3450R/	100,000	15A	8.3A	7.2A
3455R	100,000	4.4FLA 25.4LRA	2.2FLA 13.2LRA	-
	6,000	5.8FLA 34.8LRA	2.9FLA 17.4LRA	-
3455RBV	100,000	15A	8.3A	-
	6,000	5.8A 34.8LRA	2.9A 17.4LRA	-

A: Amps

FLA: Full Load Amps

LRA: Locked Rotor Amps

Contacts are available for millivolt and milliamp applications.

*Includes UL and CSA ratings.

Consult Elmwood Sensors for additional ratings.

Key Features:

- Electric Rating to 15 Amp 120 VAC Resistive
- Environmental Exposure 0° to 350° F (-18° to 177° C)
- UL recognized and CSA certified and European Approved
- Single-Pole, Single-Throw (SPST)
- Pre-set and Tamperproof
- Variety of Mounting Brackets and Terminals Available

SERIES 3450/3450R/3455R/3455RBV 15 AMP THERMOSTATS

Standard Temperature Characteristics

Operating Temperature Range The tightest specification determines the group	Tolerance Allowable ^a ± at mean temperature set points				Standard Mean Differential Nominal degrees between opening and closing points		Price Group ^a
	Open		Close		°F	°C	
	±°F	±°C	±°F	±°C			
32° to 79°F 0° to 25°C	5	2.8	8	4.4	30-50	16-28	I
	5	2.8	7	3.9	25-29	14-16	II
	5	2.8	6	3.3	20-24	11-13	III
	5	2.8	6	3.3	15-19	8-11	IV
80° to 200°F 25° to 95°C	5	2.8	8	4.4	30-50	16-28	I
	5	2.8	7	3.9	25-29	14-16	II
	5	2.8	6	3.3	20-24	11-14	III
	6	2.2	5	2.8	15-19	8-11	IV
201 to 250°F 96° to 120°C	6	4.4	8	4.4	30-50	16-28	I
	6	3.9	7	3.9	25-29	14-16	II
	6	3.3	6	3.3	20-24	11-14	III
	6	2.8	6	2.8	15-19	8-11	IV
251 to 302°F 121.7° to 148.9°C	7	3.9	8	4.4	30-50	16-28	I
	7	3.9	7	3.9	30-50	16-28	II
	7	3.9	7	3.9	20-29	11-16	III
	6	3.3	7	3.9	15-19	8-11	IV

^aGrouped according to level of accuracy required. Group I with greatest latitude is less expensive than Group II, etc. Please consult factory for temperature ranges, tolerances and differentials not noted. The operating temperature ranges include tolerances.

The ± tolerances shown have been established after careful review of many thermostat applications. Attempts should be made to establish the widest acceptable tolerance possible. For example, the chart may list a tolerance of ±5°F (±2.8°C); however, ±6°F (±3.3°C) may be acceptable for the application at reduced cost.

Note: Temperature checking methods may be slightly different, and allowance for a 1.8°F (1°C) variance should be considered.

See Section B of the Terminal and Bracket Guide for dimensional characteristics.

Operating Parameters

Dielectric Strength	MIL-STD-202 Method 301 -2000 VAC 60 Hz - Terminal to Case
Insulation Resistance	MIL-STD-202 Method 302 Cond. B - 500 Megohms - 500 Volts DC applied
Environmental Exposure	0° to 350°F (-18° to 177°C)
Operating Temp. Range	32° to 302°F (0° to 150°C)
Contact Resistance	MIL-STD-202, Method 307 - 50 Milliohms
Marking	MIL-STD-1285
Weight	6 Grams (Brackets and wire leads not included)
Materials	Base: Phenolic Terminals: Plated Brass or Steel Closure: Aluminum, Stainless Steel, or Brass Brackets: Aluminum, Stainless Steel, or Brass Contacts: Silver

UL and CSA Listings

UL and CSA Listings are for use in equipment where the acceptability of the combination of the thermostat and equipment is determined by Underwriters' Laboratories, Inc. and/or the Canadian Standards Association.

UL File E36103, UL File SA4469 (3455RBV only), UL File MH8267 (3455R only), CSA File 21048.

FS-927 Series – Small Design For Tight Instrumentation Packages

Flow Rate Settings: 0.10 GPM to 1.50 GPM

Port Size: 1/4" NPT

Primary Construction Material: Brass

Setting Type: Fixed

Measuring only 1" x 2-3/4", these compact switches are ideal for use where space is at a premium. Designed for use with water and oil, these switches are suitable for high volume OEM applications. They are ideal for coolant or lubricant flow monitoring in portable equipment and many other applications with space constraints.



U.L. Recognized — File No. E31926

Specifications

Wetted Materials	
Housing and Piston	Brass
Spring	316 Stainless Steel
Other Wetted Parts	Stainless Steel
Operating Pressure, Maximum	1000 PSIG
Operating Temperature	-20°F to +225°F (-29°C to +107°C)
Set Point Accuracy	±15% Maximum
Set Point Differential	20% Maximum
Switch*	SPST, 20 VA
Inlet/Outlet Ports	1/4" NPT
Electrical Termination	No. 18 AWG, 24" L., Polymeric Lead Wires

*See "Electrical Data" on Page 3 for more information.

How to Order – Standard Models

Specify Part Number based on flow setting and switch operation.



Liquids other than water: Special calibration is available from GEMS for media other than water. Please consult factory with your requirements, including flow media, operating pressure, flow set point and liquid viscosity (SSU).

Flow Setting GPM	Part Numbers	
	Normally Open @ No Flow	Normally Closed @ No Flow
0.10	70820	70826
0.25	70821	70827
0.50	70822	70828
0.75	70823	70829
1.00	70824	70830
1.50	70825	70831

Notes:


- Flow settings are calibrated using water @ +70°F on increasing flow, with units in a vertical position (lead wires up).
- Care should be taken by specifiers to ensure fluid compatibility with the above listed wetted materials.
- Use of 50 micron filtration is recommended.

Dimensions





An FS-927 unit is shown silhouetted against the already small FS-925 unit. It illustrates just how little space is required to provide protection to your valuable OEM equipment.

Pressure Drop – Typical



Tests conducted with units in a vertical position (lead wires up), with water at +70°F

REAR VIEW OF MAGNET SHOWN WITH TERMINAL COVER ITEM 10 REMOVED.




MAGNET SPECIFICATIONS

POLE DIAMETER	150 mm [6.1"]
POLE GAP	D-1/2" [1.27"]
FILE GAP	D-1/2" [1.27"]
FILE GAP	100, 70, 50, 30, 20 mm
	[4.1"] [3.7"] [2.1"] [1.1"]
COILS (Series connected)	D, 67 ohm
MAX RESISTANCE	500/4V
MAX CURRENT	3.17 A [10.9 gm] @ D, 9 bar, 112 ps-1d
FIDELITY	3.17 A [10.9 gm] @ D, 9 bar, 112 ps-1d
TEMPERATURE	OPEN CIRCUIT ABOVE 50°C (122°F) [1.7 gm]
FLUID INTERLOCK	OPEN CIRCUIT 1.5 (1/min [0.7 gm])
MECHANICAL	600 HD (1580 US)

NOTE: DO NOT EXCEED THE MAXIMUM SPECIFIED COIL RESISTANCE OR COIL CURRENT AND POSSIBLE DAMAGE MAY OCCUR.

INTERLOCK SCHEMATIC



USER CONNECTION

COOLANT FLOW: DK
 OPEN: BLOW 2.5 l/min
 MAX CURRENT [Non inductive]
 TEMPERATURE INTERLOCK: DK
 OPEN: ABOVE 50°C
 MAX CURRENT 0.5A 120VAC [Non inductive]
 TEMP SWITCH LINK: NO USER CONNECTION
 SIGNAL GROUND

REVISED

REV	DESCRIPTION	SHEET	DATE	APPROVED
1	ISSUED	1		
2	REVISED	1		
3	REVISED	1		
4	REVISED	1		
5	REVISED	1		
6	REVISED	1		
7	REVISED	1		
8	REVISED	1		
9	REVISED	1		
10	REVISED	1		

REV	DESCRIPTION	QTY	UNIT	APPROVED
1	1.00000000	1	PCB	
2	1.00000000	1	PCB	
3	1.00000000	1	PCB	
4	1.00000000	1	PCB	
5	1.00000000	1	PCB	
6	1.00000000	1	PCB	
7	1.00000000	1	PCB	
8	1.00000000	1	PCB	
9	1.00000000	1	PCB	
10	1.00000000	1	PCB	
11	1.00000000	1	PCB	
12	1.00000000	1	PCB	
13	1.00000000	1	PCB	
14	1.00000000	1	PCB	
15	1.00000000	1	PCB	
16	1.00000000	1	PCB	
17	1.00000000	1	PCB	
18	1.00000000	1	PCB	
19	1.00000000	1	PCB	
20	1.00000000	1	PCB	
21	1.00000000	1	PCB	
22	1.00000000	1	PCB	
23	1.00000000	1	PCB	
24	1.00000000	1	PCB	
25	1.00000000	1	PCB	
26	1.00000000	1	PCB	
27	1.00000000	1	PCB	
28	1.00000000	1	PCB	
29	1.00000000	1	PCB	
30	1.00000000	1	PCB	
31	1.00000000	1	PCB	
32	1.00000000	1	PCB	
33	1.00000000	1	PCB	
34	1.00000000	1	PCB	
35	1.00000000	1	PCB	
36	1.00000000	1	PCB	
37	1.00000000	1	PCB	
38	1.00000000	1	PCB	
39	1.00000000	1	PCB	
40	1.00000000	1	PCB	
41	1.00000000	1	PCB	
42	1.00000000	1	PCB	
43	1.00000000	1	PCB	
44	1.00000000	1	PCB	
45	1.00000000	1	PCB	
46	1.00000000	1	PCB	
47	1.00000000	1	PCB	
48	1.00000000	1	PCB	
49	1.00000000	1	PCB	
50	1.00000000	1	PCB	



GMW
 MAGNET ASSEMBLY
 MODEL: 3473-50
 A1 11801281

DO NOT SCALE
 GMW
 465 Inverness Rd, San Jose, CA 94106
 Tel: (650)992-8293 Fax: (650)992-8199

DATE: 04/25/90
 DRAWN: J. J. J.
 CHECKED: J. J. J.
 APPROVED: J. J. J.

REV 1 OF 1

REAR VIEW OF MAGNET SHOWN WITH TERMINAL COVER ITEM 10 REMOVED




MAGNET SPECIFICATIONS

- PILE DIAMETER 150 mm [6"]
- PILE GAP 0-70 mm [3.8"]
- PILE GAPS CALIBRATION 150 mm [6"]
- PILE GAPS 147mm [5.8"] [2"] [1.5"]
- COILS (series connected)
- MAX RESISTANCE 0.87 OHM
- MAX POWER [air] 204/17W
- MAX POWER [water] 704/59W
- COOLING: 5 l./sec/min [1.6 gpm] 2.0 bar [30 psi]
- THERMAL INTERLOCK OPEN CIRCUIT ABOVE 50°C [122°F]
- FLOW INTERLOCK OPEN CIRCUIT BELOW 4 l./sec/min [1.1 gpm]
- WEIGHT: 610 kg [1340 lbs]

NOTE: DO NOT EXCEED THE MAXIMUM SPECIFIED COIL RESISTANCE OR COIL OVERHEATING AND POSSIBLE DAMAGE MAY OCCUR

INTERLOCK SCHEMATIC



REVISED

REV	DESCRIPTION	SHEET	DATE	APPROVED
1	ISSUED	1		
2	REVISED	1		
3	REVISED	1		
4	REVISED	1		
5	REVISED	1		
6	REVISED	1		
7	REVISED	1		
8	REVISED	1		
9	REVISED	1		
10	REVISED	1		

ITEM	QTY	DESCRIPTION	UNIT	DATE	APPROVED
1	1	MAGNET ASSEMBLY			
2	1	TERMINAL COVER			
3	1	WATER PUMP			
4	1	COOLANT FLOW SWITCH			
5	1	TEMPERATURE INTERLOCK			
6	1	FLOW SWITCH			
7	1	TERMINAL BLOCK			
8	1	WATER PUMP			
9	1	COOLANT FLOW SWITCH			
10	1	TEMPERATURE INTERLOCK			
11	1	FLOW SWITCH			
12	1	TERMINAL BLOCK			
13	1	WATER PUMP			
14	1	COOLANT FLOW SWITCH			
15	1	TEMPERATURE INTERLOCK			
16	1	FLOW SWITCH			
17	1	TERMINAL BLOCK			
18	1	WATER PUMP			
19	1	COOLANT FLOW SWITCH			
20	1	TEMPERATURE INTERLOCK			
21	1	FLOW SWITCH			
22	1	TERMINAL BLOCK			
23	1	WATER PUMP			
24	1	COOLANT FLOW SWITCH			
25	1	TEMPERATURE INTERLOCK			
26	1	FLOW SWITCH			
27	1	TERMINAL BLOCK			
28	1	WATER PUMP			
29	1	COOLANT FLOW SWITCH			
30	1	TEMPERATURE INTERLOCK			
31	1	FLOW SWITCH			
32	1	TERMINAL BLOCK			
33	1	WATER PUMP			
34	1	COOLANT FLOW SWITCH			
35	1	TEMPERATURE INTERLOCK			
36	1	FLOW SWITCH			
37	1	TERMINAL BLOCK			
38	1	WATER PUMP			
39	1	COOLANT FLOW SWITCH			
40	1	TEMPERATURE INTERLOCK			
41	1	FLOW SWITCH			
42	1	TERMINAL BLOCK			
43	1	WATER PUMP			
44	1	COOLANT FLOW SWITCH			
45	1	TEMPERATURE INTERLOCK			
46	1	FLOW SWITCH			
47	1	TERMINAL BLOCK			
48	1	WATER PUMP			
49	1	COOLANT FLOW SWITCH			
50	1	TEMPERATURE INTERLOCK			
51	1	FLOW SWITCH			
52	1	TERMINAL BLOCK			


COOLANT FLOW: DK
 OPEN: BELOW 4.5 l./min
 MAX CURRENT [Non Inductive]
 TEMPERATURE INTERLOCK: DK
 OPEN: ABOVE 50°C
 MAX CURRENT 0.5A 120VAC [Non Inductive]
 TEMP SWITCH LINK: NO USER CONNECTION
 SIGNAL GROUND
 SPARE

DO NOT SCALE
 GMW
 965 Invercauld Rd San Diego, CA 94070
 Tel: (650)992-8293 Fax: (650)992-8199
 MAGNET ASSEMBLY
 MODEL: 3473-70
 A1 11801282
 SHEET 1 OF 1

REVISIONS
 1. REVISED BY: [Name]
 DATE: [Date]
 2. REVISED BY: [Name]
 DATE: [Date]

REV	DESCRIPTION	SHEET	DATE	APPROVED
1	REVISED	01/17/01	10:30:24	


MODEL: 3473-50 MAGNET



REAR VIEW TERMINAL COVER REMOVED

*** WARNING ***
 CHECK AC POWER VOLTAGE AND FREQUENCY MATCH POWER SUPPLY
 SPECIFIED REQUIREMENTS BEFORE APPLYING AC INPUT POWER

POWER TEN MODEL: P62B-4050A POWER SUPPLY



POWER SUPPLY REAR VIEW


NOTE

1. POWER SUPPLY SHOWN WITH 2 PHASE 208V AC INPUT
 2. REFER TO TABLE ON DWG 13900210 FOR AC INPUT RATINGS OTHER THAN 2 PHASE 208V
- N/S = NOT SUPPLIED

REV	DATE	DESCRIPTION	BY	CHKD
1	08/23/01	CABLE TIE NYLON 2.5mm WIDE BAR-LOCK		
2	08/23/01	CABLE TIE NYLON ADHESIVE BAR-LOCK		
3	11/15-20/01	POWER CABLE TIE SB 3 CORE 12AWG N/S		
4	11/15-20/01	PLUG 3PHASE/20A NYLON BRANT		
5	11/15-20/01	CABLE CLAMP THOMAS & BETTS		
6	12/20/01/10	AC TERMINAL BOX		
7	11/20/01/80	CURRENT & INTERRUPT CABLE BUL		


DO NOT SCALE
 ELECTRICAL WIRING
 3473-50/P62-4050
 11900710

DATE	08/23/01	BY	AW	CHKD	AW
REV	1	DATE	08/23/01	BY	AW
REV	2	DATE	11/15/01	BY	AW
REV	3	DATE	11/15/01	BY	AW
REV	4	DATE	11/15/01	BY	AW
REV	5	DATE	11/15/01	BY	AW
REV	6	DATE	12/20/01	BY	AW
REV	7	DATE	11/20/01	BY	AW



REVISIONS
 1. REVISED
 DATE 08/17/01
 ENGINEER J.A.

REVISED
 SHEET 1 OF 1
 DATE 08/17/01
 ENGINEER J.A.



- NOTE:
- ① MAGNET WATER FLOW AND TEMP INTERLOCKS LINKED IN SERIES ON MAGNET.
 - ② WHITE AND GREEN WATER FLOW INTERLOCK WIRES TIED TO TERM 7 & TERM 8.
 - ③ REMOVE LINK BETWEEN PIN 3 & PIN 8 ON POWER SUPPLY AT J1.
 - ④ AND CONNECT INTERLOCKS WIRES AS SHOWN.
 - ⑤ INSULATE WHITE AND GREEN WIRES WITH HEAT SHRINK SLEEVING (not used).

ELECTROMAGNET SYSTEM ELECTRICAL REQUIREMENTS

AC INPUT	3 PHASE, 480 VAC	2000V	-	-
AC INPUT	FLAT, 120V CURRENT	15	-	-
RECOMMENDED WIRE SIZE	20	-	-	-
RECOMMENDED WIRE TYPE	1.15-20P	-	-	-
RECOMMENDED AC INPUT SIZE	3 SQ/MM	-	-	-

NOTE: DRAWING SHOWS POWER SUPPLY FOR US 200V AC POWER

DO NOT SCALE

GMW
 P.O. Box 2518, Redwood City, CA 94064
 Tel: (415)802-8200 Fax: (415)802-8296


ELECTRICAL WIRING
 3473-50/P62-4050
 A1 13900210

DATE: 08/17/01
 SHEET 1 OF 1

REVISIONS
 1. REVISED BY: [Name]
 DATE: [Date]
 2. REVISED BY: [Name]
 DATE: [Date]

REV	DESCRIPTION	DATE	BY
1	REVISED	08/20/01	EDWARDS


MODEL: 3473-70 MAGNET



REAR VIEW TERMINAL COVER REMOVED

*** WARNING ***
 CHECK AC POWER VOLTAGE AND FREQUENCY MATCH POWER SUPPLY
 SPECIFIED REQUIREMENTS BEFORE APPLYING AC INPUT POWER

POWER TEN MODEL: P63C-60110A POWER SUPPLY




POWER SUPPLY REAR VIEW

NOTE

1. POWER SUPPLY SHOWN WITH 3 PHASE 208V AC INPUT
 2. REFER TO TABLE ON DWG. 13900200 FOR AC INPUT RATINGS OTHER THAN 3 PHASE 208V
- N/S = NOT SUPPLIED

REV	DATE	DESCRIPTION	BY
1	11/15/01	POWER CORD TYPE SD 4 CORE 12AWG N/S	
2	11/15/01	PLUG 3PHASE/20K NYLON, BR/ANT N/S	
3	1/23/02	CABLE CLAMP HINDAS & BELITS	
4	5/10/04	AC TERMINAL BOX	
5	1/18/00	CABLE TIE NYLON Green WHIT. BAR-LOK	
6	1/18/00	CABLE TIE INTERLOCK CABLE BUN	

REV	DATE	DESCRIPTION	BY
1	04/25/70	3473-70	
2	08/20/01	REVISED	
3	08/20/01	REVISED	




DO NOT SCALE
 ELECTRICAL ASSY
 3473-70/P63-60110
 11900700
 A1


GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (415)802-0202 Fax: (415)802-0296

REVISIONS
 1. REVISED BY: [Name]
 DATE: [Date]
 2. REVISED BY: [Name]
 DATE: [Date]

REAR VIEW TERMINAL COVER REMOVED



MODEL: 3473-70 ELECTROMAGNET



DANFYSIK MODEL: 858 65V/70A POWER SUPPLY

DANFYSIK MODEL: 858 GPIB SWITCH SETTINGS

POWER SUPPLY MOTHER BOARD		GPIB INTERFACE BOARD	
DIP SWITCH SW1	DIP SWITCH SW4	DIP SWITCH SW1	DIP SWITCH SW2
1 OFF	1 OFF	1 OFF	1 OFF
2 ON	2 OFF	2 DN	2 OFF
3 ON	3 OFF	3 DN	3 ON
	4 ON	4 DN	4 OFF
	5 OFF	5 DN	5 ON
	6 OFF	6 OFF	6 ON
	7 DN	7 OFF	7 ON
	8 OFF	8 OFF	8 ON

DIP SW1 AND SW2 LOCATED ON M/B
 NOT USED FOR GPIB INTERFACE
 SW3/SW4 USE NORMAL LOGIC OFF=0
 SW1 SW2 SHOWN SET TO ADDRESS 1 [Default]
 SW1 1-5 USES AVERAGE LOGIC OFF=1
 SW4 USES NORMAL LOGIC OFF=0

NOTE

1. POWER SUPPLY SHOWN WITH US 3 PHASE 208V AC INPUT
2. GPIB INTERFACE IS OPTIONAL EQUIPMENT
3. REFER TO TABLE ON DWG 1390100 FOR AC INPUT RATINGS OTHER THAN 3 PHASE 208V AC
4. GPIB INTERFACE FITTED INTERNALLY ON LATER MODELS OF 858 N/S=NOT SUPPLIED

*** WARNING ***
 CHECK AC POWER VOLTAGE AND FREQUENCY MATCH POWER SUPPLY SPECIFIED REQUIREMENTS BEFORE APPLYING AC INPUT POWER

REV/ISSUE

REV	DESCRIPTION	SHEET	DATE	APPROVED
A	REVISION	08/01	01/02/2014	

1	LEAD Q01-1-3	GPIB INTERFACE CABLE	N/S
2	44-1214	POWER CABLE, TYPE SB 4 CORE 12AWG N/S	N/S
3	1115-20P	PLUG, 3PHASE/20A, NYLON, BRANIT	N/S
4	Q0481	CABLE TIE ADHESIVE MFG. NYL. BAR-10K	N/S
5	10314	CABLE TIE, NYLON 6/6MM WHIT. BAR-10K	N/S
6	1199020	CABLE TIE, NYLON 6/6MM WHIT. BAR-10K	N/S

DO NOT SCALE
 DIMENSIONS
 DIMENSIONS IN PARENTHESIS ARE FOR INFORMATION ONLY
 DIMENSIONS IN SQUARE BRACKETS ARE FOR INFORMATION ONLY
 DIMENSIONS IN BRACKETS ARE FOR INFORMATION ONLY
 DIMENSIONS IN BRACKETS ARE FOR INFORMATION ONLY

GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (650)992-8292 Fax: (650)992-8199

ELECTRICAL ASSEMBLY
 3473-70/DF858
 11990220
 A1

REVISIONS
REV. 1
DATE
BY/APP. ENGINEER

REVISED
SHEET 1 OF 1
DATE
BY/APP. ENGINEER

MODEL: 3473-50 MAGNET

KEPCO MODEL: BOP 20-20M BIPOLAR POWER SUPPLY


REAR VIEW TERMINAL COVER REMOVED

SLAVE SUPPLY REAR VIEW

MASTER SUPPLY REAR VIEW

*** WARNING ***
CHECK AC POWER VOLTAGE AND FREQUENCY MATCH POWER SUPPLY SPECIFIED REQUIREMENTS BEFORE APPLYING AC INPUT POWER


NOTE
1. POWER SUPPLY SHOWN WITH 115V AC INPUT
2. CABLE INTERFACE IS ORIGINAL EQUIPMENT
3. REFER TO TABLE FOR DWG. 13900000 FOR AC INPUT RATINGS OTHER THAN 115V AC INPUT
N/S=NOT SUPPLIED



5	2	KEPCO	AC POWER CORD 115V US TYPE	
4	1	UL-20-15	MASTERS POWER OUTLET STRIP 20A	N/S
3	1	5-2DP	AC POWER TRIG 20A (Part of Item 4)	N/S
2	1	5MD 48B-2	CRB CABLE (2M LMG)	Z
1	1	1B900X0	CABLENT & INTERLOCK CABLE 4M	Z

DO NOT SCALE
ELECTRICAL ASSY
3473-50/BOP20-40
A1 11900230
A

GMW
P.O. Box 2578, Redwood City, CA 94064
Tel: (415)802-8200 Fax: (415)802-8296




REVISIONS

NO.	DESCRIPTION	DATE	BY
1	REVISED		

REVISED

NO.	DESCRIPTION	DATE	BY
1	REVISED		



NOTE:
1. CUT OUT GREEN WIRE. USE WHITE, RED, BLACK ONLY.

WIRING TO MAGNET TEMP INTK

1	1	1	1
2	2	2	2
3	3	3	3
4	4	4	4
5	5	5	5
6	6	6	6
7	7	7	7
8	8	8	8

1	1	1	1
2	2	2	2
3	3	3	3
4	4	4	4
5	5	5	5
6	6	6	6
7	7	7	7
8	8	8	8

DO NOT SCALE
ELECTRICAL WIRING
BOP 2020/2020

GMW
P.O. Box 2578, Redwood City, CA 94064
Tel: (650)992-8290 Fax: (650)992-8299

13
A1 139000990
A

ASSEMBLY SEQUENCE MODEL# 3473 ELECTROMAGNET ON TO 45° MOUNTING AND ROLLING/ROTATING BASE

FIGURE 1




FIGURE 2




FIGURE 3




FIGURE 4




FIGURE 5





FIGURE 6



NOTE: LEFT LAG MOUNTED BY BOTH OTHER FRONT LIFTING SYSTEMS EXCEPT NOT SHOWN

NOTE: LEFT LAG MOUNTED BY BOTH OTHER FRONT LIFTING SYSTEMS EXCEPT NOT SHOWN

NOTE: LEFT LAG MOUNTED AT BOTH OTHER FRONT LIFTING SYSTEMS EXCEPT NOT SHOWN


Min. x 3/16 IN. DIA. 1/4 IN. SPRING WASHER, FLAT WASHER (2 Pkts) MOUNTED TO MOUNT

(1) (2) (3) (4) (5) (6)

Min. x 3/16 IN. DIA. 1/4 IN. SPRING WASHER, FLAT WASHER (2 Pkts) MOUNTED TO MOUNT

PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 PROPRIETARY TO GMM INC. IT MUST NOT BE
 REPRODUCED OR TRANSMITTED IN ANY FORM OR
 BY ANY MEANS, ELECTRONIC OR MECHANICAL,
 INCLUDING PHOTOCOPYING, RECORDING, OR BY
 ANY INFORMATION STORAGE AND RETRIEVAL
 SYSTEMS, WITHOUT THE WRITTEN PERMISSION
 IN WRITING BY GMM INC.

REVISIONS				
REV	DESCRIPTION	DRAFT	DATE	APPROVED
A	RELEASE		11/09/03	A.MARTIN
B	NEW T/B; ROTATE DWG 90°. ADD DIMS		05/11/94	G.DOUGLAS
C	ADD VERT DIM MAG C/L TO ROT BASE TOP		02/18/03	G.DOUGLAS



ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
PARTS LIST				

DRAWN: A.MARTIN DATE: 11/09/03
 CHECK: DATE: (UNLESS OTHERWISE SPECIFIED)
 ENGINEERING DATE:

DO NOT SCALE FROM DRAWING DIMENSIONS & TOLERANCES
 (UNLESS OTHERWISE SPECIFIED)

LINEAR	INCHES	MILLIMETERS	TITLE
XXX	±.001	±0.03	ROLL/ROTATING BASE
XX	±.01	±0.3	MODEL: 3473 45°MTG
X	±.08	±1.3	
DEG.	±.3	±0.5	

FINISH: 63
 THIRD ANGLE PROJECTION


NEXT ASSY SYSTEM: 3473
 SOFTWARE: AUTOCAD 2000

GMM
 955 Industrial Rd, San Carlos, CA 94070
 Tel: (650)802-8292. Fax: (650)802-8298.

SCALE: 1:2.5 WT KG SHEET 1 OF 1

PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 BELONGING TO GMW, INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY
 MANNER THAT VIOLATES THE TERMS OF CONFIDENTIALITY
 AGREEMENTS WITH GMW, INC.

REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	05/11/91	



M16 X 30 HEX. BDL 1, S/S
 SPRING WASHER, FLAT WASHER
 C8 PLSJ. BRACKET TO BASE/MAGNET


ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
1				
2				
3				
4				

DRAWN	DATE	DO NOT SCALE	TITLE
G. DOUGLAS	10/05/94	FROM DRAWING	ROLLING BASE ASSY
CHECK	DATE	REVISIONS & DIMENSIONS	MODEL: 3473/H MTC
ENGINEERING	DATE	WASHER DRAWING REVISIONS	
		LENGTH	
		WIDTH	
		THICKNESS	
		FINISH	
		THIRD ANGLE PROJECTION	
SOFTWARE:	1.3		


GMW	P.O. Box 2578, Redwood City, CA 94064
	Tel: (415)802-8292 Fax: (415)802-8298
SCALE	1:2.5 WT kg
SHEET	1 OF 1

REVISIONS


REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	06/27/04	C DOUGLAS
B	ADD 3 HOLES IN HOLES	09/17/04	C DOUGLAS
C	ADD MOUNTING HOLES TO ROLLING BASE HOLES	07/06/07	C DOUGLAS



TOP VIEW



FRONT VIEW




SIDE VIEW


MOUNTING HOLES	
A-4405	DIRECT MOUNTING
B-3473/3472	DIRECT MOUNTING
C1/C2	54.72 45° MOUNTING
D1/D2	54.73 45° MOUNTING
C1/C3	54.72 HORIZ MOUNTING
D1/D3	54.73 HORIZ MOUNTING
E	MWD MOTOR DRIVE MOUNTING
F	MWD MOTOR DRIVE MOUNTING UNDER ROLLING BASE

REV	DATE	BY	CHKD	DESCRIPTION
1	06/27/04	C DOUGLAS	C DOUGLAS	RELEASE
2	09/17/04	C DOUGLAS	C DOUGLAS	ADD 3 HOLES IN HOLES
3	07/06/07	C DOUGLAS	C DOUGLAS	ADD MOUNTING HOLES TO ROLLING BASE HOLES

REV	DATE	BY	CHKD	DESCRIPTION
1	06/27/04	C DOUGLAS	C DOUGLAS	RELEASE
2	09/17/04	C DOUGLAS	C DOUGLAS	ADD 3 HOLES IN HOLES
3	07/06/07	C DOUGLAS	C DOUGLAS	ADD MOUNTING HOLES TO ROLLING BASE HOLES



DISCLAIMER: THIS DRAWING IS THE PROPERTY OF GMW MAGNET SYSTEMS, INC. AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE WRITTEN PERMISSION OF GMW MAGNET SYSTEMS, INC.



TOP VIEW

SIDE VIEW

FRONT VIEW

MOUNTING HOLES	
A=5403	DIRECT MOUNTING
B=3473/3472	DIRECT MOUNTING
C1/C2=3472	45° MOUNTING
D1/D2=3473	45° MOUNTING
C1/C3=3472	HORZ MOUNTING
D1/D3=3473	HORZ MOUNTING
E=RD14MG	BASE MOUNTING
F=NRD	MOTOR DRIVE MOUNTING

NOTE 1. PARTS FITTED AT GMW

REV	DATE	DESCRIPTION	BY	CHKD
15	1	LABEL, GMW MAGNET SYSTEMS		
14	4	DN 125 A		
13	4	DN 127 B		
12	4	DN 923		
11	10	DN 127 B		
10	10	DN 433		
9	10	DN 433		
8	2	REX CHDS 4RT		
7	2	REX CHDS 4RT		
6	4	17802180		
5	4	17802125		
4	1	17802122		
3	2	17802121		
2	1	17802110		


DO NOT SCALE
 A. MARTINI (R) 06/01
 455 Inlandia Rd San Diego, CA 94070
 (650)802-8280 Fax: (650)802-8289

ROLLING BASE ASSY
 3473/3472/5403
 A1 11803170

REV A1 11803170
 DATE 12/13/10
 BY A1
 CHKD A1

REV	DATE	DESCRIPTION	BY	CHKD
A	1	REVISION		
B	1	REV 7/8 ADD HOLE 2.4 DIA 1" HOLES		
C	1	REV 1/10 ADD HOLES 1.5 DIA 1" HOLES		
D	1	ADD 3473 HOLES		
E	1	ADD 3473 HOLES		

REVISIONS
 1. REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.
 2. REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.
 3. REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.
 4. REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.



REV	DESCRIPTION	DATE	APPROVED
1	REVISION		
2	REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.		
3	REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.		
4	REVISED TO ADD DIMENSIONS FOR THE NEW DESIGN.		

MOUNTING HOLES

A=540S DIRECT MOUNTING
 B=3473/3472 DIRECT MOUNTING
 C1/C2=3472 45° MOUNTING
 D1/D2=3473 45° MOUNTING
 C1/C3=3472 HORIZ MOUNTING
 D1/D3=3473 HORIZ MOUNTING
 E=MOTORIZED ROT BASE SPOOL

- NOTES**
- 1. ADJUST SET SCREW FOR MINIMUM CLEARANCE ALLOWING FOR FULL FREE ROTATION; AND LOCITE
 - 2. FORM DETAIL TO PLATE DIA TO PREVENT ENDS FROM SPRINGING LOOSE
 - 3. GREASE BEARING SURFACES BEFORE ASSEMBLY
 - 4. ITEM 14 AND ITEM 15 ONLY USED IF ROTATING BASE SOLD SEPARATELY. SEE DWG NO 11802430 FOR DETAILS ON MOUNTING ROTATING BASE TO ROLLING BASE

ITEM	DESCRIPTION	QTY	UNIT
15	WASHER M12 SPRING 5/5	4	
14	ROD M12 x 56 HEX HD 5/5	4	
13	SCREW M6 x B SHSS .0VAL P1 5/5	4	
12	LABEL ANGLE GRADUATIONS 0-360°	1	
11	LABEL VERNIER INDEX	1	
10	SPACER BEARING	1	
9	BALL BEARING 114mm I15/32" BA SIF	1	
8	UN 9/12	4	
7	CLAMP RETAINING	4	
6	HANDHELD MIO	2	
5	LOWER THRUST BEARING PLATE	1	
4	UPPER THRUST BEARING PLATE	1	
3	SCREW M12 x 20 SK HD CAP 5/5	4	
2	TRANSITION PLATE	1	
1	TRANSITION PLATE	1	

DO NOT SCALE

GMW
 P.O. Box 2578 Redwood City, CA 94064
 Tel: (650)902-0200 Fax: (650)902-0299


ROTATING BASE
 3473/3472/5403

DATE: 11/18/02
 DRAWN: A1
 CHECKED: A1
 SCALE: 1:2 (1:1)
 SHEET 1 OF 1

PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 BELONGING TO GDMW, INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED TO OTHERS OR USED IN
 ANY MANNER WITHOUT THE EXPRESS WRITTEN
 PERMISSION OF GDMW, INC.

ELECTROMAGNET REAR VIEW

VIEW SHOWS TERMINAL COVER REMOVED



REVISIONS

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	08/12/93	G.DONIGLAS
B	ADD DIMS. AND MOUNTING HOLE DETAILS	05/11/94	G.DONIGLAS
C	SHOW REAR OF MAGNET WITH WATER/ELECT CONNECTIONS	08/25/99	G.DONIGLAS

HOLE PATTERN
 FOR
 VERTICAL MOUNTING

*** WARNING ***
 WHEN THE 3473 ELECTROMAGNET IS VERTICALLY MOUNTED
 IT MUST BE BOLTED TO THE FLOOR SECURELY AS SHOWN.

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
6	4		BOLT, M16 or 5/8" EXPANSION	
5	4	DIN 933	BOLT, M6 x 3D HEX HD S/S	
4	8	DIN 127 B	WASHER, M16 SPRING LOCK, S/S	
3	8	DIN 125 A	WASHER, M16 FLAT, S/S	
2	2	17803180	VERTICAL MOUNTING BRACKET	
1	1	3473	ELECTROMAGNET, 150MM SDA or 70A	


DATE	SCALE	TITLE
08/12/93	1:2	VERTICAL MGT ASSY
05/11/94		MODEL: 3473
08/25/99		

REV	DATE	BY	APP'D
A2	118032250		
C			

SCALE 1:2 WT kg SHEET 1 OF 1

PROPRIETARY
 HAS REMAINING COPY RIGHTS RESERVED. INFORMATION
 IS UNCLASSIFIED BY DATE 08/11/93 BY A.MARTIN
 REASON: UNCLASSIFIED BY DATE 05/08/94 BY G.DONAGLAS
 AUTHORITY: 48 CFR 1.101-11.1
 A WARNING BY CIVIL INC.

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	08/11/93	A.MARTIN
B	NEW T/B INCREASE HOLE SIZE TO 18MM	05/08/94	G.DONAGLAS



Ø18.0 [4 PLS]
 .5 X 45° C/SK
 REMOVE BURR ON INSIDE

- NOTES**
- 1 MATERIAL: 150 X 75 X 6mm M. S CHANNEL
 - 2 PAINT INSTRUMENT TAN TD BSL TP8580010

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
PARTS LIST				
DRAWN	A.MARTIN	DATE	08/11/93	
CHECK		DATE		
ENGINEERING		DATE		
DO NOT SCALE FROM DRAWING				
CHECK DIMENSIONS AGAINST DRAWING DIMENSIONS				
LINKER	ADT/ST	DATE	08/03	
MARK	4/09	DATE	08/03	
DES	5/03	DATE	08/03	
FINISH	1.2	DATE	08/03	
THIRD ANGLE PROJECTION				
SOFTWARE: 13				
AUTOCAD				
SCALE	1:1	WT KG		
SHEET	1	OF	1	




GMW
 P.O. Box 2578, Redwood City, CA 94064
 Tel: (415)802-8292 Fax: (415)802-8298

VERT MTG BRACKET
 MODEL: 3473

PROPRIETARY
 HAS REMAINING COMMON LAW PROTECTION
 RESPECTFULLY TO OWN. ALL RIGHTS RESERVED.
 REPRODUCTION OR DISSEMINATION OF THIS
 DRAWING OR ANY PART THEREOF WITHOUT
 THE WRITTEN PERMISSION OF GMW INC.
 IS PROHIBITED BY LAW.

PART No 17900172 [LH] BRACKET

REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	04/22/94	G.DONIGLAS
B	REVISE BOTTOM MTG HOLES, ADD CHAMFER	05/09/94	G.DONIGLAS



BREAK SHARP EDGES 2mm [2 PLS]
 MACHINE THIS SLOT ON PART No 17900172 ONLY

Ø18.0 [4 PLS]
 0.5 X 45° C/SK
 REMOVE BURR ON INSIDE





- NOTES
- 1 MATERIAL: 150 X 75 X 6mm MS CHANNEL
 - 2 PAINT INSTRUMENT TAN TO BSL TP85B0010
 - 3 BREAK ALL SHARP EDGES 1mm

PART No 17900172 [LH] AS DRAWN
 PART No 17900171 [RH] MIRROR IMAGE

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
PARTS LIST				
DRWN	G.DONIGLAS	DATE	10/22/94	
CHKD		DATE		
ENGR		DATE		
DO NOT SCALE FROM DRAWING				
P.O. Box 2578, Redwood City, CA 94064				
Tel: (415)802-4292 Fax: (415)802-8298				
GMW				
TITLE: HORZ MTG BRACKET				
MODEL: 3473				
DATE	3473	REV	A2	179000170
SCALE	1:2	WT kg	SHEET	1 OF 1
SOFTWARE	AUTOCAD 13			

REVISIONS

REV	DESCRIPTION	DATE	BY	APPROVED
A	ISSUE			
B	ADD BRACKET AND COACH SCREW			
C	ADD BRACKET AND COACH SCREW			
D	ADD NUTS			




- NOTE:
1. THE 3475 SHIPPING CRATE HAS A ONE PIECE COVER
 2. REMOVE ALL SHARP EDGES 0.2MM

- COVER REMOVAL:
1. REMOVE THE COVER SECURING SCREWS
 2. GRIP THE COVER AT THE TOP LH AND RH CORNERS
 3. LIFT THE COVER VERTICALLY HIGH ENOUGH TO CLEAR THE MAGNET
 4. MOVE THE COVER SIDEWAYS AND PLACE ON FLOOR



SHIPPING WEIGHT 660 KG (1,450 lbs)

REV	DATE	DESCRIPTION	BY	APPROVED
A	01/01/01	ISSUE		
B	01/01/01	ADD BRACKET AND COACH SCREW		
C	01/01/01	ADD BRACKET AND COACH SCREW		
D	01/01/01	ADD NUTS		



PROPRIETARY
 THIS DRAWING CONTAINS CONFIDENTIAL INFORMATION
 BELONGING EXCLUSIVELY TO GDMV, INC. IT MUST NOT BE
 REPRODUCED OR DISCLOSED TO OTHERS OR USED IN ANY
 MANNER THAT VIOLATES THE CONFIDENTIALITY RIGHTS
 AS MAINTAINED BY GDMV, INC.


REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
A	RELEASE	08/01/91	A. MARTIN
B	CORRECT TOGGLE LATCH VIEW (INVERT)	08/17/99	G.DONIGLAS



R5 TYPICAL
 ON CORNERS INDICATED

#4 X 12 C/SK POSIDRIVE

TOGGLE LATCH
 DZUS TL803B



12mm CUSTOM BOARD
 GLUE AND PIN NAIL JOINTS

ITEM	QTY	PART NUMBER	DESCRIPTION	NOTE
PARTS LIST				
DO NOT SCALE FROM DRAWING				
CHECK DIMENSIONS AND DIMENSIONS AGAINST DRAWING DIMENSIONS				
DRYAN A. MARTIN	DATE	BB/01/94	DATE	BB/01/94
ENGINEERING	DATE		DATE	
SCALE	1:2	WT	kg	
SOFTWARE	AUTOCAD 13			
TITLE	POLE PACKING BOX			
REV	A2 18800410			
REV	B			

A B C D 4 3 2 1

A B C D 4 3 2 1