
Agilent Technologies

Agilent Stethos Electronic Stethoscope

## Users Guide for Stethos Electronic Stethoscope

Electronic Stethoscope  
Elektronisches Stethoskop  
Stéthoscope électronique  
Estetoscopio electrónico  
Stetoscopio elettronico  
Estetoscópio electrónico  
Elektronische stethoscoop  
Elektroninen stetoskoopi  
Elektroniskt stetoskop  
Elektroniskt stetoskop  
Elektronický fonendoskop  
Elektronikus sztetoszkóp  
Stetoskop elektroniczny  
Elektronik steteskop  
Elektronisk stetoskop  
электронный стетоскоп  
Ηλεκτρονικό Στηθοσκόπιο  
電子聽診器  
电子听诊器

<b>STETHOS Electronic Stethoscope</b>	
<b>Operating Instructions .....</b>	<b>1</b>
<b>STETHOS Elektronisches Stethoskop</b>	
<b>Benutzerhandbuch .....</b>	<b>3</b>
<b>Stéthoscope électronique STETHOS</b>	
<b>Mode d'emploi .....</b>	<b>5</b>
<b>STETHOS Elektronische stethoscoop</b>	
<b>Gebruikershandleiding.....</b>	<b>7</b>
<b>Stetoscopio elettronico STETHOS</b>	
<b>Istruzioni per l'uso .....</b>	<b>9</b>
<b>Estetoscopio electrónico STETHOS</b>	
<b>Instrucciones de funcionamiento .....</b>	<b>11</b>
<b>Estetoscópio Electrónico STETHOS</b>	
<b>Manual do Utilizador.....</b>	<b>13</b>
<b>STETHOS, elektroniskt stetoskop</b>	
<b>Bruksanvisning .....</b>	<b>15</b>
<b>Elektroninen STETHOS stetoskooppi</b>	
<b>Käyttöohjeet.....</b>	<b>17</b>
<b>STETHOS, elektroniskt stetoskop</b>	
<b>Brugervejledning.....</b>	<b>19</b>
<b>Ηλεκτρονικό Στηθοσκόπιο STETHOS, της</b>	
<b>Οδηγίες Λειτουργίας .....</b>	<b>21</b>
<b>STETHOS, elektroniskt stetoskop</b>	
<b>Bruksanvisning .....</b>	<b>23</b>
<b>Электронный стетоскоп STETHOS</b>	
<b>Руководство пользователя .....</b>	<b>25</b>
<b>Stetoskop elektroniczny STETHOS</b>	
<b>Instrukcja obsługi.....</b>	<b>27</b>
<b>STETHOS Elektronický Fonendoskop</b>	
<b>Provozní pokyny .....</b>	<b>29</b>
<b>STETHOS Elektronikus sztetoszkóp</b>	
<b>Felhasználói kézikönyv .....</b>	<b>31</b>
<b>STETHOS Elektronik Stetoskop</b>	
<b>İşletme Talimatları .....</b>	<b>33</b>
<b>STETHOS 電子聽診器</b>	
<b>取扱説明書 .....</b>	<b>35</b>
<b>STETHOS 电子听诊器</b>	
<b>使用说明 .....</b>	<b>37</b>
<b>STETHOS 電子聽診器</b>	
<b>使用說明 .....</b>	<b>39</b>

# STETHOS Electronic Stethoscope

## Operating Instructions

### Introduction

Congratulations! The STETHOS electronic stethoscope you have purchased is the result of years of research and development to provide you the ultimate in medical auscultation.

### Technical Features

#### General Features

With its built-in micro-chip, your new STETHOS is at the cutting edge of technology. Its main features include:

- Fully electronic design eliminates sound loss and resonance effects
- Bell, Diaphragm and Extended Diaphragm modes with 8 levels of sound amplification
- Hand tremor and ambient noise filtering with auto impact noise muting and auto shut-off for extended battery life
- Ergonomically designed to look and feel like your traditional stethoscope


#### Frequency Response

The key to the success of STETHOS is its similarity with high fidelity stereo systems, which include an amplifier and an equalizer that allow the user to obtain sound amplification as well as a controlled frequency response.

This unique feature lets you amplify sounds within the range of the standard Diaphragm and Bell modes, while minimizing undesirable frequencies. Thus, it is possible to reduce interference caused by ambient noise. The Extended Diaphragm mode also enables you to hear high frequency sounds such as those produced by mechanical valvular cardiac prostheses.

The following chart shows the improved frequency response provided by STETHOS over conventional stethoscopes.

The below graph demonstrates how the STETHOS amplifies the preferred low frequency sounds while reducing the undesirable high frequency ambient noise.

### Operating Instructions

#### Headset Fitting

To ensure optimal transmission of sound, STETHOS was developed with the ear morphology in mind. The earphones are thus identified by the letters R (right) and L (left).

#### Tension Adjustment


To ensure maximum comfort for the user's ears, STETHOS is equipped with a roll-on tension adjustment sleeve. Once adjusted, pressure is maintained at a constant level.

#### Power On

To turn on the instrument, press the **•** key located at the center of the chest piece. A red light will flash to indicate that the Diaphragm mode is activated.

#### Mode Selection

To set to the Bell mode, press the **•** key and wait for the green light to flash. The **•** key allows you to alternate between the Diaphragm and Bell modes during auscultation. To access the Extended Diaphragm mode, press and hold the **•** and **+** keys simultaneously. The red light will flash 3 times.

#### Volume Control

To increase the volume level, press the **+** key. To lower it, press the **-** key. STETHOS has 8 volume levels. In the Diaphragm and Extended Diaphragm modes, the red light will flash to indicate a change in volume level. In the Bell mode, the green light will flash to indicate a volume change. If the light (red or green) does not flash, the minimum (1) or maximum (8) level has been reached.

#### Programming Volume Level

To save your preferred volume, simply adjust it to your desired level and press the **•** key and **-** key simultaneously. The light will flash 3 times. Your STETHOS will automatically reach this pre-programmed volume as soon as it is activated. The volume can be reprogrammed by repeating these steps.

#### Power Off

To turn off the instrument, simply press the **•** and **-** keys simultaneously and the red and green lights will flash once. This step can be avoided as the instrument's automatic shut-off feature will turn the device off two minutes after activating the last key.


#### Low Battery Indicator

The simultaneous flashing (3 times) of the green and red lights, when the instrument is turned on, is an alarm indicating that the batteries must be replaced.

#### Replacing Batteries

To replace the batteries in your STETHOS:

1. Unscrew the non-chill ring by pressing it firmly against the palm of your hand.
2. Slide off the battery cover.
3. Remove the batteries. (Please discard in an environmentally safe manner.)

To reassemble the chest piece, repeat these steps in reverse order. However, to

screw the Non-Chill Ring back on, it is preferable to turn the chest piece upside down while pressing on the Battery door.

In view of the low power consumption of your STETHOS, batteries will operate for approximately one year under normal use. Replacement batteries must be type 357.

The batteries are not intended to be recharged. If you are not using your STETHOS for a period, beyond expected life of battery, it is suggested you remove the batteries.

#### Replacing Eartips

To ensure optimal comfort, two types of Eartips are supplied with your STETHOS. To replace them, simply unscrew the tips and replace them with the alternative pair.

#### Maintenance

STETHOS requires no specific maintenance. However, in order to reduce the risk of cross contamination, we recommend wiping the stethoscope clean with alcohol. The eartips should be cleaned periodically to prevent cerumen build-up.

**Avoid exposure to extreme temperatures and any contact with concentrated solvents. Do not immerse your stethoscope in any liquid.**

**There are no user serviceable parts with exception of the batteries inside the STETHOS.**

#### Warranty for your STETHOS

Every STETHOS is warranted against defects in materials or workmanship for a period of one year from purchase date. During the warranty period, Agilent Technologies will, at its option, either repair or replace components or product which prove to be defective.

The foregoing warranty shall not apply to defects resulting from unauthorized modifications or misuse by Buyer. No other warranty is expressed or implied, Agilent Technologies specifically disclaims the implied warranties or merchantability and fitness for a particular purpose.

#### Packaging Symbols

Latex-Free

Consult guide for more information

Complies with CAN\CSA-C22.2 no. 601.1-M90 and UL2601-1 standards

Type BF: Degree of protection against electric shock

#### Parts Ordering and Service

The following replacement parts are available:

- M4504-60050      Replacement Batteries
- M4504-60000      Eartips soft, 12.7mm
- M4504-60010      Eartips firm, 13.3mm
- M4504-60020      Diaphragm
- M4504-60030      Non-chill ring
- M4504-60040      Battery cover
- M4504-91000      User's manual

Use only Agilent Technologies approved parts, using any other parts or materials can degrade performance, including safety.

For parts ordering, technical specifications and service information, contact an Agilent Technologies specialist near you or visit our websites at

[www.agilent.com/healthcare/stethos](http://www.agilent.com/healthcare/stethos)  
[www.agilent.com/healthcare/medsupplies](http://www.agilent.com/healthcare/medsupplies)

See back inside cover for country location and phone number.

#### Technical Specifications

Supply voltage	4.5 VDC
Supply current	1mA (typ)

#### Operating Environment

Temperature	-5°C to 60°C
Humidity	20 to 95 %RH
Altitude	up to 4550m

Degree of protection against ingress of water: IPX0

Degree of safety of application in the presence of a flammable mixture with air or with oxygen or nitrous oxide: Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

#### Storage Environment

Temperature	-40°C to 70°C
Humidity	15 to 95 %RH
Altitude	up to 4550m

#### Radio Frequency Electromagnetic Field

as per EN 60601-1-2 (1993)

#### Intended Use

The Stethos electronic stethoscope is intended for the purpose of electronically amplifying and listening to biological sounds, such as the heart, arteries and veins and other internal organs.

#### Contraindications

Stethos is not intended as a substitute for medical care or for diagnosis and treatment by unlicensed or unsupervised/untrained persons. The Stethos electronic stethoscope is contraindicated for use other than listening to biological sounds.


#### Caution

**Choking Hazard** from small parts included with STETHOS. Keep out of reach of small children.

## Basic Operating Functions

- Power On** red light one time  
**Diaphragm Mode** red light one time  
**Bell Mode** green light one time  
**To increase volume** one time at each level  
**To decrease volume** one time at each level

*Note:* When the instrument is turned on, the red light flashes indicating that the diaphragm mode has been activated.

[www.agilent.com/healthcare/stethos](http://www.agilent.com/healthcare/stethos)

## Advanced Features

- Extended Diaphragm Mode** (red light flashes 3 times)

- Volume Programming** (red or green light flashes 3 times)

*Note:* When turned on, the instrument sets itself automatically to the pre-programmed volume.

- Power Off**

or the instrument will shut off automatically after 2 minutes

## Low Battery Detector

red and green lights flash 3 times only at Power On

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**Canada**  
Agilent Technologies Canada Inc.  
Healthcare Business Centre  
17500 Transcanada Highway  
Kirkland, Québec H9J 2X8

1-800-291-6743 option 1  
Monday to Friday, 8:30 a.m. to 5:00 p.m. EST

**United States**  
Agilent Technologies, Inc.  
Healthcare Solutions Group Headquarters  
3000 Minuteman Road  
Andover, MA 01810

1-800-225-0230  
M-F, 8AM - 7 PM EST

**Latin America, Brazil, Chile and Mexico**  
Agilent Technologies, Inc.  
5200 Blue Lagoon Drive  
9th Floor  
Miami, FL 33126  
(305) 267-4220

**Europe**  
Agilent Technologies  
Deutschland GmbH  
Herrenberger Straße 130  
71034 Boeblingen  
Germany  
fax (+49) 7031 464-2346

**Austria:**  
Phone: (+43) 1 25125 - 600  
Fax: (+43) 1 25125 - 7360  
[office\\_medical@agilent.com](mailto:office_medical@agilent.com)

**Belgium/Luxembourg**  
(+32) 2 778-3085

**Finland**  
(+358) 10 855 2450

**France**  
(+33) 01 69 82 62 16

**Germany**  
0130 860-781

**Italy**  
(+39) 2 926081

**Medical Distribution**  
**Europe, Middle East, Africa**  
(+41) 22 780-6162

**Netherlands**  
(+31) 20 547-6335

**Poland**  
(+48) 22 608 7700

**Russia**  
Phone (+7) 095 797 3979  
Fax: (+7) 095 797 3903


**Spain**  
(+34) 91 631 3100  
902 30 40 50

**Sweden**  
(+46) 8 506 488 12

**Switzerland**  
(+41) 1 735-9225 (German Swiss)  
(+41) 22 780-6111 (Suisse Romande)

**United Kingdom**  
0870 607 76 77  
0870 607 76 88

Copyright © 2000  
Agilent Technologies  
Printed in Canada 3/99  
AW-M4504-91000-2 Rev. D

**Agilent Technologies**

**Agilent Technologies**

\*+H179M4504A3.\*

© Copyright 2000  
Agilent Technologies  
[www.agilent.com/healthcare/stethos](http://www.agilent.com/healthcare/stethos)  
Printed in Canada  
AW-M4504-91000-2 Rev. D  
Printed on recycled paper