# **AMBIVISION**

# digital portable B mode Ultrasound Scanner

AV-3618V

**User Manual for Veterinary use** 

**Version A0** 

AMBISER Technology Corp., Ltd.

# AMBISEN User Manual

# Revision Effective Date Page

Document

No.

AV-3618V A0 2008/02/21

3/68

#### AV-3618V Digital Portable B mode ultrasound scanner

#### **Table of Content**

1 II	NTRODUCTION	
1.1	Copyright	6
1.2	STATEMENT	
1.3	MANUFACTURER'S WARRANTY	
1.4	SYSTEM LABELS DIAGRAM	8
1.5	CONTRAINDICATION OF THE USAGE	10
2 S	UMMARY	1
2.1	Characters	11
2.2	RANGE OF APPLICATION	11
2.3	IMPLEMENT STANDARDS	11
3 S	YSTEM INTRODUCTION	12
3.1	APPEARANCE	12
3.2	Panel	12
3.3	TECHNICAL SPECIFICATIONS	13
3.4	BLOCK SCHEMATIC DIAGRAM	14
3.5	FUNDAMENTAL PRINCIPLE	14
3.6	EQUIPMENT CONFIGURATION	14
3.	6.1 Standard configuration pieces	14
3.	6.2 Optional pieces	15
4 II	NSTALLING THE EQUIPMENT	10
4.1	OPERATING ENVIRONMENTAL REQUIREMENTS	16
4.2	UNPACKING AND ACCEPTANCE INSPECTION	16
4.3	STARTING THE INSTALLATION	16
4.4	PROBE HOOK'S INSTALLATION (SELECTED MATING)	17
4.5	EQUIPOTENTIAL CONNECTION	18
5 K	EYBOARD AND TRACKBALL OPERATION	20
5.1	SCREEN DISPLAY	20
5.2	OPERATING KEYBOARD	20
5.3	TRACK BALL	28
6 O	PERATION PROCEDURES	29
6.1	GETTING STARTED	29
6.2	DIAGNOSIS	29
6.3	MENU OVERVIEW	29
6.4	Measurement	33
	Copyright ©2010 by AMBISEA Technology Corp., Ltd - All Right Reserved	

### AMBISER User Manual

# Document No. AV-3618V Revision A0 Effective Date 2008/02/21 Page 4/68

#### **AV-3618V Digital Portable B mode ultrasound scanner**

	6.4.	1 Distance measurement	33
	6.4.2	2 Circumference and Area Measurement	34
	6.4	3 Volume measurement	36
	6.4.4	4 Cardiac rate measurement (only at B/M modes)	37
	6.5	SAVE AND SWITCH TO STORYBOARD VIEW MODE	38
	6.5.	1 Saving a image	38
	6.5.2	2 View a saved image	38
	6.5	3 Review storyboard mode	39
	6.6	OBSTETRICAL CALCULATION	40
	6.7	DIAGNOSTIC REPORT	46
	6.8	IMAGE PROCESSING	47
	6.9	IMAGE PRINT	48
	6.10	IMAGE UPLOAD TO COMPUTER	49
	6.11	SHUT DOWN THE EQUIPMENT	49
7	CHI	ECK AND MAINTENANCE	50
	7.1	DETECT	50
	7.2	EFFECTIVE SERVICE LIFE	50
	7.3	HOST MAINTENANCE	50
	7.4	PROBE MAINTENANCE	50
	7.5	CLEANSING	52
	7.6	CORRECTLY USE THE PROBE	53
	7.7	EQUIPMENT DETECTION AND CALIBRATION	53
8	TRA	ANSPORTATION AND STORAGE	54
	8.1	ENVIRONMENTAL REQUIREMENTS OF TRANSPORTATION AND STORAGE	54
	8.2	TRANSPORTATION	
	8.3	STORAGE.	
9		IPLE MALFUNCTION EXAMINATION AND TROUBLESHOOTING	
		EXAMINATION	
	9.2	TROUBLESHOOTING	55
<b>A</b> l	PPENI	DIX A LABELS DIAGRAM	56
	1. AV-3	3618V Main unit nameplate diagram	56
	2. AV-3	3618V DEVICE PROBE NAMEPLATE DIAGRAM	56
	3. AV-3	3618V DEVICE OPTIONAL PROBE NAMEPLATE DIAGRAM	56
	4.AV-3	3618V Packing diagram	57
A]	PPENI	DIX B ACOUSTIC OUTPUT PARAMETERS	58
	TEST N	Mode: B-mode Transducer Type: C1-7/60R	58
		MODE: B/M-MODE TRANSDUCER TYPE: C1-7/60R	
	12011	Convigant ©2010 by AMRISEA Technology Corp. Ltd All Right Reserved	

# AMBISEN User Manual

Document No.	AV-3618V
Revision	A0
Effective Date	2008/02/21
Page	5/68

#### AV-3618V Digital Portable B mode ultrasound scanner

TEST MODE: B-MODE	TRANSDUCER TYPE: C1-8/20R	60
TEST MODE: B/M-MODE	TRANSDUCER TYPE: C1-8/20R	61
TEST MODE: B/M-MODE	TRANSDUCER TYPE: EC1-2	62
TEST MODE: B/M-MODE	TRANSDUCER TYPE: EC1-2	63
TEST MODE: B-MODE	TRANSDUCER TYPE: L1-3	64
TEST MODE: B/M-MODE	TRANSDUCER TYPE: L1-3	65
APPENDIX C OBSTETRI	C TABLE	66
GESTATIONAL TABLE 1: E	QUINE	66
GESTATIONAL TABLE 2: B	OVINE	67
GESTATIONAL TABLE 3: SI	TEED	68

AMBISEN User Manual	Document No.	AV-3618V
	Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner	Effective Date	2008/02/21
AV-3010V Digital Portable B filode ditrasoulid Scaliller	Page	6/68

#### 1 Introduction

#### 1.1 Copyright

This publication, including pictures and illustrations, is property of AMBISEA Technology Corp., Ltd. This is under the protection of International Copyright Law.

#### 1.2 Statement

Information in this document is not changed in annotated. The manufactory does not state nor observe any other warranty basing on this point, and definitely give up any implied warranty basing on any special purpose of selling or making benefit.

Without previous written permission from the manufacture, this document must not be photocopied, reproduced or translated into any other languages.

The manufacture preserves the right of revision on this document without notice.

Some pictures in this manual, which are schematic diagrams for indication only, may disaccord with the real object; therefore the real object should be regarded as the final.

#### 1.3 Manufacturer's warranty

AMBISEA Technology Corp., Ltd assumes the responsibility for device security, reliability and performance only under the preconditions that the disassembly, assembly and maintenance of the device are all performed by its assigned professional and the device is used strictly in compliance with the operation manual.

AMBISEA Technology Corp., Ltd ensures a guarantee period within a year and a half since the delivery day and promises there is no problem with the new device within the warranty period, AMBISEA Technology Corp., Ltd will maintain the device and replace the parts of non-man-made damages free. But will not repair or replace the device that the surface is damaged.

This guarantee will only available for failures occurred, which the device is operated in compliance by using the operation manual. And the guaranteed device can only be used in the prescribed range given in the manual.

This guarantee excludes losses or damages caused by external reasons, such as thunder struck, earthquake, theft, unsuitable use or abuse and refitting the device.

AMBISEA Technology Corp., Ltd do not has responsible for the devices' damages, which caused by the arbitrary connection to the different devices.

AMBISEA Technology Corp., Ltd do not takes on the responsible for losses, damages or injuries caused by the delayed service requests.

When there is problem with the device (products), please contact AMBISEA Technology Corp., Ltd in the meantime explains the device model, serial number, date of purchase and the nature of problem.

#### **Announcements**

To ensure the operational safety and the long-term stable performance, please read this operation manual closely and understand the equipment functions, operation and maintenance information before operating the device, especially pay attention to the contents of "Warning", "Caution" and "Note".

Misoperation or inobservance of the instructions, which were given by the manufacturer or its agents, may

### AMBISEN User Manual

Document No.	AV-3618V
Revision	A0
Effective Date	2008/02/21
Page	7/68

#### AV-3618V Digital Portable B mode ultrasound scanner

result in device damage or personal injury.

The following convention works through this manual to declaration the special emphasis of some information.

"Warning": displays the neglect of it will cause the severe personal injury, death or realized property loss.

"Caution": displays the neglect of it will cause slight personal injury or property damage.

"Note": to remind the user of the installation, operation or maintenance information. These information is very significant but without risk. Any warning against dangers shall not be contained in NOTE.

#### Icon description

#### Equipment icon description

Note! Refer to accompanying documents

Turn-on (general supply)

O Disconnect (general supply)

→ Signal output

Equipotent Equipotent

**IPX7** Watertight-proof

Risky voltage

• USB port

RS-232 port

Recycling of the carton

Shipping package icon description

Ţ	Handle with Care	
.c.X	Temperature extremes	
<u>11</u>	Upward	
5	Piling limit	
<b>党</b>	Keep dry	
*	Protect against heat	

AMBISEN User Manual	Document No.	AV-3618V
	Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner	Effective Date	2008/02/21
AV-3010V Digital Portable B filode ditrasoulid Scaliller	Page	8/68

#### 1.4 System labels diagram

Refer to Appendix A for device labels.

#### Common suggestion of the Equipment regulation

- ◆ In operation
- 1. Heat radiation holes are strictly prohibited to be covered.
- 2. After shutdown, do not switch on the device within 2 3 minutes.
- 3. In scanning, if any abnormal cases are found, stop the scanning immediately and shut down the device.
- 4. Patient is not allowed to touch any non-application parts.
- 5. Do not press the keyboard overexertion; otherwise the equipment will be damaged.
- ◆ After operation
- 1. Switch off the device.
- 2. Pull out the plug from the power supply socket instead of pulling the cable.
- 3. Cleaning the coupling gel on the probe with soft medical sterilized tampon.
- 4. Put the probe into its box.

#### **Common safety information**

System approach to the safety and reliability of the device are taken into consideration as a part of designing and producing process, because this will also affect the safety of the operators and patients.

The following safety precaution must be implemented:

- 1. The device shall be operated by qualified operating staff or under their instructions.
- 2. Do not reconfigure the parameters of the device without permission. If necessary, please ask for AMBISEA Technology Corp., Ltd or its authorized agent for helping.
- 3. The device has already been regulated into its optimal performance. Do not adjust any preset control or the switch unless operate access the instructions in the manual.
- 4. If there is fault, please shut down the device at once and contact with AMBISEA Technology Corp., Ltd or their authorized agents.
- 5. If the device is needed to connect with the other company electronic or mechanical devices, please consult AMBISEA Technology Corp., Ltd.



# Document No. AV-3618V Revision A0 Effective Date 2008/02/21 Page 9/68

#### AV-3618V Digital Portable B mode ultrasound scanner

6. Operation, storage and transportation environment of the device

Environmental requirements under the normal operation:

a) Environment temperature range:  $+5^{\circ}$ C $\sim$ +40 $^{\circ}$ C

b) Relative humidity range: ≤80%

c) Atmosphere pressure range: 86KPa~106KPa

Environment requirements of the storage and transportation:

a) Environment temperature range: - 40°C~+55°C

b) Relative humidity range:  $10\% \sim 100\%$ 

c) Atmosphere pressure range: 50KPa~106KPa

- 7. Ultrasound might cause harm on human body, please avoid the long time radiation. Appendix B for information of sound output parameters.
- 8. Please use standard power cable as the input line of the network power supply for the adapter to reduce risk.
- AMBISEA Technology Corp., Ltd shall not take any responsibility for any risk resulted from propelled / unauthorized re-fitment by the users.
- 10. The power supply plug is the standard three pins plug; its protective grounding pin (terminal) should be connected with the protective grounding line from the main power.

#### **EMC** statement

AV-3618V shall not affect the basic performance of radio service and other equipments and can work well in its stated electromagnetic environment.

#### **Warning**

When AV-3618V works in strong electromagnetic environment exceeding its statement, its image might be interfered and the diagnoses could be influenced. At this time, please stop operating until the EMC interference is removed.

#### **Warning**



**AV-3618V Digital Portable B mode ultrasound scanner** 

Document No.	AV-3618V
Revision	A0
Effective Date	2008/02/21
Page	10/68

When AV-3618V works at the state of overlapping or paralleling with other equipments, there might be unexpected EMC problems; If it must work close to other equipments, please observe carefully and check if some equipment is influenced by unexpected EM coupling.

#### Warning

Replacement of non-standard probe parts may cause unexpected EMC problem.

#### NOTE

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (e.g. IEC 60950 for data processing equipment and IEC 60601-1-1 for medical equipment). Furthermore, all configurations shall comply with the valid version of the system standard IEC 60601-1-1. Everybody who connects additional equipment to the signal input or signal output part configures a medical system, and is therefore responsible that the system complies with the requirements of the valid version of the system standard IEC60601-1. IF in doubt, consult the technical service department or your local representative.

#### 1.5 Contraindication of the usage

- ◆ The equipment is not suitable for inspection on organs containing with gas such as lung, etc.
- ♦ It is recommended not to check the parts with wounded or acute inflammation to avoid the cross infection.

Patient in the following situations are not allowed to be checked with vagina and rectum probes: inflammation of the vagina, such as trichomonas virginities, colpomycosis, venereal disease etc: The unmarried, deformity of vagina, menstrual period, postmenopausal vagina atrophy, person with difficulty carry on the vagina ultrasonic examination, colporrhagia, placenta previa patient, etc.

The following patient is prohibited to use the puncture

High blood pressure, coronary heart disease, disturbance of blood coagulation and the bleeding tendency patient

AMBISEN User Manual	Document No.	AV-3618V
/ 11 IDIOCEN User Manual	Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner	Effective Date	2008/02/21
AV-3010V Digital Portable B illoue ultrasound Scailler	Page	11/68

#### 2 Summary

#### 2.1 Characters

This equipment is high resolution linear/convex ultrasound scanning diagnostic equipment. It adopts micro-computer control and digital scan converter (DSC), digital beam-forming (DBF), real time dynamic aperture (RDA), real time dynamic receiving apodization, real time Dynamic receiving focusing (DRF), Digital frequency Scan (DFS), 8 segments TGC, frame correlation technologies to endue its image with clarity, stability and high resolution.

There are five display modes: B, B+B, B+M, M and 4B; And 256 gray scale.

The system can process real time image display, freeze, save, load, zoom, up and down flip, left and right flip, black and white flip, and capacity cine loop; Multi-level scanning depth, angle, dynamic range, acoustic power, frame correlation factor regulation and focus number, focal space, focus position, etc.

Date, clock display; name, sex, age, doctor and hospital annotation; Distance, circumference, area, volume and heart rate measurement; EDD Measurement for equine, bovine, sheep, dog and cat. It offers 12 body marks. And Multiple-choice probes for the clinic diagnosis demands.

PAL-D video output offers connection to external video image printer and large-screen display and other equipments. High speed USB port provides real time image transfer to the PC.

Touch function folding keyboard made by apheliotropic silica gel and the trackball provides immediate, convenient and flexible operation.

The equipment is jet molding enclosure and potable structure. The usage of non-industrial frequency transformer controls power supply, programmability parts and the surface mounting technology (SMT) make the whole unit highly integrated.

The equipment consists with the main unit and electronic convex array probe. Standard configuration is LV2-3/6.5MHz Rectal Linear probe. And C1-7/60R/3.5MHz convex array probe, L1-3/7.5Mhz high frequency linear array probe, C1-8/20R/5.0MHz micro convex probe, EC1-2/13R/6.5MHz Transvaginal and EL1-3/7.5MHz endo-rectal probe for options.

#### 2.2 Range of application

Application in animals' diagnosis, such as horse, moo-cow, sheep, cat and dog, etc.

#### 2.3 Implement Standards

The equipment is designed and manufactured rigidly according to GB10152-1997  $(B \ Ultrasound \ Diagnosis \ Equipment)$ , GB9706.1—1995  $(Medical \ Electrical \ Equipment)$  Part I: Safety General Requirement) and GB9706.9-1997  $(Medical \ Electrical \ Equipment)$  Special Safety Requirement for Medical B-ultrasound Diagnosis and Monitor Equipment). The equipment is C class B-Ultrasound diagnosis equipment and its electricity shock prevention type is Class I B type.

The environment test requirements conform with requirements of climate environment test group II, mechanical environment test group II of GB/T-14710-93  $\langle$  Medical Electrical Equipment Environment Requirements and Test Methods $\rangle$ .

AMBISEN	User Manual	Document No.	AV-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	12/68

#### 3 System Introduction

#### 3.1 Appearance




Figure 3-1 appearance's size sketch

#### 3.2 Panel





Figure 3-2. Front panel sketch



Copyright ©2010 by AMBISEA Technology Corp., Ltd - All Right Reserved

AMBISEN	User Manual	Document No.	AV-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	13/68

Figure 3-3. Back panel sketch

### 3.3 Technical specifications

MC	DDEL			AV-3618V		
		Standard		Option		
Pı	robe	LV2-3/6.5M Hz Rectal Linear array probe	EL1-3/7.5 MHz Endo-Rectal probe L1-3/7.5MHz Linear probe	C1-8/20R/5.0 MHz Micro-convex probe	C1-7/60R/3.5 MHz convex array probe	EC1-2/13 R/6.5MHz Transvagi nal probe
Detect I	Depth(mm)	≥80	≥80	≥80	≥160	≥60
Resoluti on	Lateral	≤1 (Depth≤60)	≤1 (Depth≤60)	≤3 (Depth≤60)	<pre>\$2 (Depth≤80)      ≤3 (80<depth≤ 130)<="" pre=""></depth≤></pre>	≤1 (Depth≤ 40)
(mm)	Axial	≤1 (Depth≤80)	$\leq$ 1 (Depth $\leq$ 80)	≤1 (Depth≤60)	$\begin{array}{c} \leqslant 2 \\ (\text{Depth} \leqslant 80) \\ \leqslant 3 \\ (80 < \text{Depth} \leqslant \\ 130) \end{array}$	≤1 (Depth≤ 40)
Blind z	cone(mm)	€3	€3	≪8	€5	€7
Geometr	Horizontal	€5	€5	≤15	≤15	€10
position precisio n	Vertical	€5	<b>≤</b> 5	≤10	≤10	<b>≤</b> 5
Moni	tor size	10 inch				
Displa	ay mode	B、B+B、B+M、M、4B				
Image §	gray scale			256 level		
Cin	e loop	809 frames (MAX)				
Image	estorage	32 frames				
Scar	n angle	50%~100%				
Scar	n depth	40mm~240mm				
Acoust	tic power			2 steps		
Dynan	nic range		10	00dB~130dB		
Ima	ge flip	Up/down, left/right, black/ white				
Focus	position			Adjustable		

AMBISEN	User Manual	Document No.	AV-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	14/68

Focal space	5 level
Measurement	Distance, circumference, area, volume, heart rate. GA, EDD
Notation	Date, time, name, sex, age, doctor, hospital name. Full screen words edit.
output report	1 type
posture mark	12 types
USB port	High speed USB 2.0 (device)
power consumption (MAX)	100VA

#### 3.4 Block Schematic Diagram




Figure 3-4 Electric principal block diagrams

#### 3.5 Fundamental principle

Digital Ultrasonic Diagnostic System is based on the following procedure: different human organism possess different densities and speeds of transmission of ultrasound, i.e. different acoustic impedance (product of media density and sound speed). When piezo-chip (transductor) gets certainly regulated electric impulse, it will produce ultrasound with certain frequency. When this ultrasound (sound energy) is injected into human body, different organ surfaces will produce reflection echo, the different size reflection is received by the transduction which emitted ultrasound and is changed into electric impulse, when this electric impulse is amplified, demodulated, digital scanned, shifted and some other handling, video standard signal is produced and organ cross-Sectional images are displayed on the monitor.

#### 3.6 Equipment configuration

#### 3.6.1 Standard configuration pieces

- Mainframe
- ◆ LV2-3/6.5MHz Rectal Linear probe

AMBISEN	RISER Manual		AV-3618V for vet
		Revision	A0
AV 2610V Digital D	Effective Date	2008/02/21	
AV-3618V Digital Po	Page	15/68	

- application software suite(CD)
- power cord
- ◆ High speed USB cable
- ◆ 2 pieces of fuse tube F2AL250V
- ◆ Coupling gel 250ml
- User Manual
- ◆ Final examination report
- Packing List

#### 3.6.2 Optional pieces

- ◆ C1-7/60R/3.5MHz convex array probe
- ◆ C1-8/20R/5.0MHz micro convex probe
- ◆ L1-3/7.5MHz Linear probe
- ◆ EL1-3/7.5MHz endo-rectal probe
- ◆ EC1-2/13R/6.5MHz Transvaginal probe
- Video printer
- Trolley
- ◆ Abdomen probe puncture support
- Probe clasp

#### **Warning**

If select spare parts models prescribed above. The manufacturer will not assume the risks such as safety problem, non-expected drop of EMC performance that cause by arbitrary adoption of spare parts out of prescription.

AMBISEN	MRISER User Manual		AV-3618V for vet
		Revision	Α0
AV 2640V Di-it-I D	Effective Date	2008/02/21	
AV-3618V Digital Portable B mode ultrasound scanner		Page	16/68

#### 4 Installing the Equipment

#### 4.1 Operating environmental requirements

- a) Environment temperature range:  $+5^{\circ}$ C $\sim$ +40 $^{\circ}$ C
- b) Relative humidity range: ≤80%
- c) Atmosphere pressure range: 86KPa~106KPa
- d) Power supply: AC110V-10%~130V+10%/220V-10%~240V+10%, 50Hz/60Hz

Using the equipment, avoid strenuous vibration, keep it away from devices with high field, intense magnetic field or high voltage; avoid strong sunlight blazing down on the display; keep the device well-ventilated, moisture proof and dustproof.

#### 4.2 Unpacking and Acceptance Inspection

After unpacking, please check the device according to the "packing List" and there has not been damage or had any fault during transportation. Install it following with the requirements and methods described in "Installation".

#### **Warning**

At acceptance inspection, if the equipment has been damaged, it must be forbidden to use for security reason.

#### Warning

The probe should be protected from felling off or crashing, and the manufacturer assumes no responsibility for this kind of hazard.

#### 4.3 Starting the installation

- 1. Check the power supply to see if it is in the expected range, then connect the equipment and the power socket with the cable (see figure 3-3).
- 2. A 96-core connector is used for electronic convex array probe (see figure 4-1).
- 3. Connection of 96-coreconnector. First, switch the locker on the probe connector to "OPEN", position the connector localization pin with the hole on the socket, then plug in the connector and switch the locker to "LOCK". When the locker is moved to "OPEN", unplug the connector to remove the probe from the main unit. (see figure 4-2).

AMBISEN	User Manual	Document No.	AV-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	17/68




Figure 4-1 96-core probe connection sketch




Figure 4-2. Probe socket sketch

#### **Note**

Avoid by all means unplugging or plugging the probe connector at state of log on in case the probe and main unit be damaged.

Once the probe is connected with the main unit, do not unplug nor plug it at discretion in case poor contact happen.

#### Warning

Must not touch the contact pin in the probe connector.

#### 4.4 Probe hook's installation (selected mating)

Put the probe hook into the hook fixing seat, wind the nut up until the probe hook is fixed in the right position. See figure 4-3, 4-4:

14515.5	44DIC C)		AV-3618V for vet
AMBISER	User Manual	Revision	A0
AV 2C10V Divital Daviah	Effective Date	2008/02/21	
AV-3618V Digital Portabl	Page	18/68	

Fix the probe hook nut



Figure 4-3 Installing the probe hook step 1




Figure 4-4 Installing the probe hook step 2

#### 4.5 Equipotential connection

See the figure:

44516.6	44514.45		AV-3618V for vet
AMBISEN	User Manual	Revision	A0
AV 2610V Digital D	Effective Date	2008/02/21	
AV-3618V Digital Portable B mode ultrasound scanner		Page	19/68



#### Warning

Equipotential: When this equipment is together used with other device, equipotential should be in consideration.

When the equipment is in use, the doctor and patient are under the risk of uncontrollable compensating current influence, which is caused by different electric potential between facilities and tangible current carrying parts. The most safe solution is set up a united equipotential network, and connect the medical equipments with the equipotential network in the treatment room.

AMBISEN	User Manual	Document No.	AV-3618V
ווטכופווויו		Revision	A0
AV-3618V Digital Po	Effective Date	2008/02/21	
AV-3618V Digital Portable B mode ultrasound scanner		Page	20/68

#### 5 Keyboard and Trackball Operation

#### 5.1 Screen display




Figure 5-1 Equipment interface illustration

#### 5.2 Operating keyboard




Figure 4-2. Operating keyboard sketch

AMBISEN	User Manual	Document No.	AV-3618V
רוטכופווויז		Revision	A0
AV-3618V Digital P	Effective Date	2008/02/21	
AV-3618V Digital Portable B mode ultrasound scanner		Page	21/68



Press these keys at annotation mode to put in characters at the cursor position. Among the alphabet keys, there are some keys with dual functions:

SN	Key	Function
1	Report	At freeze, press this key to enter or exit report interface
2	BPower	At real time or freeze, press this key to adjust acoustic power
3	N Now-P	At real time or freeze, press this key to clear notation and report nformation on the image or the characters area
4	G	Total gain control key. There are four steps: 25, 30, 35, and 40.  Press this key to set them circularly



At annotation, press this key to input space after cursor.



At annotation, press this key to delete the input characters.



Numbers are used for time, date settings, age notation and function menu and selection.

There are the following dual functions keys:

SN	Key	Function
1	<b>[-1]</b>	Press this key to move cursor or sampling line left
2	12	Press this key to move cursor up
3	[13]	Press this key to move cursor down
4	-4	Press this key to move cursor or sampling line right
5	0	Control the cursor moving speed by direction keys

AMBISEN	User Manual	Document No.	AV-3618V
רוטכופו ווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
711 00=01 21g.tm. 1		Page	22/68



Press this key to clear off the measuring marks, notation data and measurement results on the image (When some menu is displayed, it can not be cleared until the current task is finished).



Press this key to display the notation menu, and press the above mentioned numbers or letters for notation.



Press this key to display the OB table, and press the number keys to carry out the corresponding function.



Press this key to display image post-process menu, and press the number keys to carry out the corresponding function.




Press this key to restart the system when malfunction and misoperation make the system "DEAD" (unable to start the system by pressing any other key).

# W Volocity Correlation Frame-correlation factors regulation/M mode velocity option

This key has two functions.

At "B/M" or 'M' mode, enter real time state; press this key to adjust the refresh velocity of "M mode image". There are four levels: 3.005 \, 2.505 \, 2.005 \, 1.255. Press this key repeatedly to set the velocity circularly. The current velocity displays on the right of the screen.



At "B", "B/B" and "4B" modes enter real time state; press this key the "frame correlation factors." will display \_ ↑ )or ↓ to adjust the parameters of frame correlation on the right top of the screen, then press the factors. The system provides eight factors: 0.25, 0.35, 0.45, 0.55, 0.65, 0.75, 0.85 and 0.95, . Press the continuous key, the eight factors is setting circularly. The current parameter appears on the right top of the

AMBISEN	1BISEN User Manual	Document No.	AV-3618V
רוטכונויי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	23/68

screen (see figure 4-1).



#### **Probes** conversion key

This equipment is fitted with two probe connectors (A, B). Press this key to convert the probes. The current probe model appears on the right top of the screen (see figure 4-1). The probe model and their codes are:

LV2-3	LV2-3/6.5MHz Rectal Linear probe
C1-8/20R	C1-8/20R/5.0MHz Micro-convex probe
C1-7/60R	C1-7/60R/3.5MHz convex probe
EC1-2/13R	EC1-2/13R/6.5MHz Transvaginal probe
L1-3	L1-3/7.5MHz linear probe or EL1-3/7.5MHz Endo-rectal probe

#### Tips:

- The device can automatically identify probes.
- When it is connected with two probes. The system default working probe is the one connecting to socket A.
- Please shut down the system first before replacing probes. Restart the system, it can realize automatic with default setting.



At real time of freeze, press this key to save the current image: See "5.5.1 Image saving" for details.



Press this key to recall the stored image; See "5.5.2 Image loading" for details.



Press this key to flip image up and down. As given below:




Figure 5-3. Image up and down flip sketch

AMRISEN	User Manual	Document No.	AV-3618V
וויי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	24/68

# Image left-right

Press this key to flip image left and right. As given below:




Figure 5-4. Image left and right flip sketch

Probe scan direction mark (in red circle) is the mark for image left and right flip. The default scan mode is

#### figure A

# Probe frequency conversion key

Press this key to change the working frequency of the probe. The frequency of each probe is given bellow:

En do-rectal probe)

LV2-3 — 5.5MHz、6.5MHz、7.5MHz (LV2-3/6.5MHz Rectal Linear probe)

C1-7/60R — 2.5MHz、3.5MHz、5.0MHz (C1-7/60R/3.5MHz convex probe)

C1-8/20R — 4.5MHz、5.0MHz、5.5MHz (C1-8/20R/5.0MHz micro-convex probe)

EC1-2/13 — 5.5MHz、6.5MHz、7.5MHz (EC1-2/R13/6.5MHz Transvaginal probe)

L1-3 — 6.5MHz、7.5MHz、8.5MHz (L1-3/7.5MHz Linear probe EL1-3/7.5MHz

The current probe frequency displays on the top right of the screen (see figure 4-1).



For volume, fetus weight measurement and histogram (see 5.4 for details).



Preinstall function



At freeze state, press this key to enter cine loop, see "5.5.3 Cine Loop" for details.

#### 



Backlight Keyboard Backlight switch key

Press the Backlight key, switch on or switch off the Backlight of the keyboard



Preinstall function.





For distance, circumference, area, volume measurement operating. See Chapter 5.4.



Operate the trackball to select the option in menu.

Confirm the contents when noting.

At real time or freeze process, press this key to display azimuth mark  $\checkmark$  .

### B model display key

At free position, press this key to enter B mode (default mode is single B at starting up).

# Double B mode display key

Press this key to enter BB mode. And there will be two B mode images displayed on the screen. One of them is a "Frozen" image and the other is a "real-time" image; Repress this key can make the two images switch between those two "Frozen" and "real-time" modes. Press "Freeze" key can freeze these two images.

# BM/M mode display key

Press this key to enter B/M mode. B mode and M mode images will be displayed on the screen ("BM" or "B+M" for shortened form). At the same time, B mode real time image is on the left and M mode real time image is on the right.

Press this key again to enter M mode, then M mode image will be displayed on the screen. (Press this key to convert from BM to M, M to BM.)

# BB 4B mode display key

Under the single B or BB mode, press this key to enter 4B mode and display 4B mode images. One is real time, the other three are all frozen images; Press this key repeatedly to switch each image between "freeze" and

AMRISEN	User Manual	Document No.	AV-3618V
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	26/68

"real time" circularly. Press the freeze key to get four "Frozen" images.








Single B mode

BB mode

BM mode





M mode

4B mode

Figure 5-5. Five display modes sketch





#### Adjusting keys

Adjusting keys are used to adjust focus position, dynamic range, scanning angle, zoom, cine loop See" Chapter 5 Operation Procedures")

#### Fpin )\$(

#### The number of Focal point selection key

Press this key to lighten up the "focus number" on the screen right top, press again to change the focus numbers. There are maximum two focuses. Press this key to switch between one to two focus



#### Focus space selection

When set the focus number as two, their distance can be adjusted. There are 5 levels focus space 2, 3, 4, 5, 6 and press this key repeatedly to switch among them circularly.



#### Scan angle selection key

Press this key to lighten up the "Angle" on the Top right screen, press or to change the scanning angles. The angel range is physical scanning angles among 50%-100%.

Tips: When linear probe or rectal linear probe is used, this function can not be used. Reduce the

AMRISEN	MBISER User Manual	Document No.	AV-3618V
רושכונו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	27/68

#### angle can improve the image's frame rate and increase the cine loop frames

# Dynamic range adjusting key

Press this key to lighten up "DR" on the screen top right corner. Press or key to adjust the range. The equipment's adjust dynamic range among 100dB-131dB. Default is recommended.



Press or keys to adjust scan depth. The current depth displays at the bottom of screen (see figure 5-1). The depth adjusting range of each probe is given below:

		EL1-3/7.5 MHz		
	LV2-3/6.5MHz	Endo-rectal	C1-7/60R/3.5	EC1-2/13R/6.5
Probe Model	Rectal Linear	probe	MHz convex	MHz
	probe	L1-3/7.5MHz	probe	TV probe
		Linear probe		
Depth range(mm)	60~ <b>130</b>	<b>4</b> 0~90	<b>70</b> ~240	<b>50</b> ~120

Table 5-1. Probes and their depth adjusting range table

Zoom key

At real time or freeze state, press this key to display or close sampling windows and zooming window.

Move the trackball can control the "zooming window" and the "sampling window". Press switch the current moving window.




Figure 5-6. Zoom sketch

Tips:

AMBISEN	User Manual	Document No.	AV-3618V
וויי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	28/68

In order to change the sizes of "zooming window" and "sampling window", press or window, press or window, press window, press or window, press window, press or zooming window is a quarter of the image, which is displayed on the area of the screen.



#### **IJImage freeze key**

Press this key to switch between freeze and real time display.

Tips: The screen's right bottom is displayed a "freezen", which means the image, has been freezen.

#### 5.3 Track ball

Track ball is fast, convenient in operation. In this equipment, the trackball functions show as below:

- Move measuring mark during measurement.
- Select these items in the menu during operation.
- Move the cursor at annotation state.
- Move the sampling line at BM mode.
- Display the puncture guide line's and moves it.
- Control movie replay on single frame show at cine loop mode.
- Move the sampling window and zooming window at zooming state.
- •At histogram statistic mode, Move the statistic sampling window.

#### Tips:

- do not press the trackball too heavy when operating the trackball.
- Keep the trackball surface clean.

AMBISEN	User Manual	Document No.	AV-3618V
7111919611		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	29/68

#### 6 Operation Procedures

#### 6.1 Getting started

Switch on the equipment's rear panel power, indicator light shows on the panel and the start interface displays on screen. Press any key (except key) to enter the state of scanning. Adjust the monitor's brightness, contrast knob, TGC and overall gain (GAIN) to acquire a satisfied visual effect.



Figure 6-1 Start interface

#### 6.2 Diagnosis

Apply some coupling gel on the diagnoses part of the patient and press the probe acoustical window close contact with this part. The Section of the tissue image pattern will display on the screen, move the probe slightly to find out the proper depth and optimum position; at the same time, adjust the brightness, contrast, overall gain (GAIN), 8 segment TGC to acquire the best section sonogram of the part.

#### Tips:

Adjust the total gain (GAIN) knob slowly.

#### Note

- 1. Do not press the probe too heavy and avoid breaking the probe or causing discomfort.
- 2. Do use proper probe and frequency to exam the target part.

#### 6.3 Menu overview

Press Annotation key to activate notation menu. Show as below:

AMBISEN	User Manual	Document No.	AV-3618V
7111919611		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	30/68

1. NAME
2. AGE
3. SEX
4. COMMENT
5. TIME
6. HOSP
7. DOCT
8. LANGUAGE
9. ERASE
A. EXIT

Figure 6-2 Notation menu

Press key to select "1. NAME" to enter the patient name, 15 characters maximum can be put in (letter, number or space), press to delete the input mistakes. (Same as the following comment)

The name input box will appear on screen as follows:

PLEASE ENTER NAME:

PLEASE ENTER NAME:

After input, press any key (except letters and number keys) to exit. To cancel input, you can press any key. The name will display in "NAME" on the top of the screen. For example: press the character keys input WANG NING,

then press any non-character key (e.g. ) to finish input.

Press to select "2. AGE" to input patient's age, 3 characters input maximum. The age input box is listed bellow:

PLEASE ENTER AGE: ■

After input, press any key (except letter, number keys) to exit. To cancel input, you can press any key. The age will display after "AGE" on the top of the screen. For example: input 27, then press any non-number key

(e.g. ) to finish input.

Press select "3. SEX", input patient's sex. The input box is listed bellow:

PLEASE ENTER SEX:

1. MALE 2. FEMALE

Copyright ©2010 t, All Right Reserved

AMBISEN	User Manual	Document No.	AV-3618V
רוטכופווויז		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	31/68

Press to select "1. MALE", and press	to select "2	2. FEMALE".	
Press to select "4. COMMENT" to ent Move the trackball cursor, you can annotate ar the screen will display the Chinese character's	nywhere within tl		ne bottom of
cursor will input the same characters.	o shortede keyy pr	cos the corresponding named he	(6 5)) ene
When finishing the comment, press "COMMEN cancel the comment.	NT" key on the b	oard to end this operation or pres	cancel to
	Confirm	on the leastern made then among	Reference
Tips: At real time or freeze state, press change the mark's direction ,use the di		ve the location mark, then press trackball to mark the image ar	
specimen.			
Press to select "5. TIME" to modify the bellow:	e system's time a	nd date. The time and date input bo	x is listed as
	YY-MM-DD		
	HH-MM-SS		

For example, the date and time are: 2006-9-22 9:35:30, then the input should be:

YY-MM-DD 060922 HH-MM-SS 093530


Press to select "6. HOSP" to input hospital name, input letters or numbers maximum 18 characters. The hospital name input box is listed bellow:

PLEASE ENTER HOSP: ■

After input, press any key (except letter and number) to confirm exit. To cancel input, you can press any key. The hospital name will display after "Hospital" at the bottom of the screen. (see figure 4-1 shows)

Press to select "7. DOCT" to put in doctor's name in letter or number maximum 15characters. The doctor name input box is listed bellow:

AMBISEN	User Manual	Document No.	AV-3618V
וואכונו וו ז		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	32/68



Press to select "8. LANGUAGE", choose the system's interface language. The language selection box is given bellow:

PLEASE ENTER LANGUAGE: 1. CHINESE 2. ENGLISH

Press number key to select Chinese and to select English

Press to select "9. ERASE", erase image's storage section. The remind box is given bellow:

ERASE ALL STORAGE?

1. YES 2. NO

Press number key to confirm clearance. When processing, a word "ERASING..." displays on the upper left screen to indicate that system is processing clearing and no other operation can be done at this moment. When the reminder disappears, the image storage is cleared.

Press to give up and exit.

Press to select "A. EXIT" to exit comment menu.

#### Tips:

- ●When image storage erasing is in process (system reminds "Erasing..."), please do not take other operations in case the equipment was damaged.
- ●The trackball can be used for operating the menu: when the menu displays, move the trackball up

and down to select the items, when the selected item is lightened, press "confirm" to enter it.

• key functions:

At the real time or freeze states, press this key to clear all the information about the patient.

N New-P Key functions:

1. at real time or freeze state, press this key to clear off the measure marks, comments, measurement results, doctor, name, age and sex in image area

AMBISEN	User Manual	Document No.	AV-3618V
רושכונו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	33/68

2. At report interface, clear off all the information and measurement results except "Hospital", "doctor" and "comment" information.

#### 6.4 Measurement

Distance, circumference, area and volume can be measured by controlling the direction keys or trackball; Unit of distance and circumference is mm, unit of area is mm<sup>2</sup> and volume is cm<sup>3</sup>.

There are four measuring marks corresponding to measurement display:

- +: D1
- ×: D2
- **%**: D3
- **₩**: D4

In circumference, area measurement, there have two measuring marks corresponding to measurement display as below:

- +: C1, A1
- ×: C2, A2

#### 6.4.1 Distance measurement

#### Operation steps:

- 1. Press "DIST" key, and the first cursor appears on the screen;
- 2. Press "Reference" key, the second cursor appears on the screen;
- 3. Move the cursor to the start point of this measurement by controlling the trackball, press "Reference " key again to confirm the start point;
- 4. Move another cursor to the end of this measurement, press "Cancel" to fulfill the measurement and

exit; (**Tips:** Press "Reference" key repeatedly to switch between the start and end points of the measurement).

5. Continue to measure just repeat step 1-4 to get at most 4 groups of distance measurement. The measuring results display on the right below of the screen as the following figure shows:



Figure 6-3. Distance measurement sketch

AMRISEN	User Manual	Document No.	AV-3618V
רוטכופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	34/68

The four groups of data are D1, D2, D3, D4, among them D1/D2 is the ratio of D1 and D2; D3/D4 is the ratio of D3 and D4.

#### 6.4.2 Circumference and Area Measurement

The circumference and area can be measured in two ways.

Press key to display this menu as below:

PLEASE INPUT: 1. FREEHAND 2.ELLIPSE

Option 1 is FREEHAND measuring method and option 2 is ELLIPSE measuring method.

#### a. Freehand method operation:

- 1. Press key to select freehand method and call out the measuring cursor, move the cursor to the start point of the measurement.
- 2. Press key and move the cursor along the edge of the measurement area to the end point by controlling the trackball;
- 3. Press key again to end the measuring of circumference and area.
- 4. Press key and repeat the above steps 2-3 to process another measurement. Two groups of data can be measured at most. And the measuring results display on the screen right as below:



Figure 6-4 Circumference and area measuring sketch (Freehand method)

AMRISEN	User Manual	Document No.	AV-3618V
ווטכופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	35/68

C1 and A1are the circumference and area of the first group data;

C2 and A2 are the circumference and area of the second group data;

C1/C2 is the ratio of the two circumferences;

A1/A2 is the ratio of the two areas.

#### b. Ellipse measuring method

- 1. Press key to select Ellipse method, an ellipse measuring mark display on the screen, which is named the measured area, can be moved to any place on the image by trackball.
- 2. Press key, then move the trackball to change the size of this measured area.

  Move trackball left and right, the measured area shrinks or enlarges horizontally; Move trackball up and down to shrink or enlarge the measured area vertically;

#### Tips:

- Press key repeatedly to shift between moving the mark and adjust the size by trackball.
- 3. Press key, Move trackball to adjust the angle of the measured area.
- 4. Press key to finish measuring.
- 5. Press key and repeat the above steps 2-4 to process another measurement. Two group's data can be measured at most. And the measuring results display on the screen right as below:



Figure 6-5 Circumference and area measuring sketch (Ellipse method)

C1 and A1 are the circumference and area of the first group data; C2 and A2 are the circumference and area of the second group data; C1/C2 is the ratio of the two circumferences;

AMBISEN	User Manual	Document No.	AV-3618V
וושכונו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	36/68

A1/A2 is the ratio of the two areas.

#### Note:

Ellipse measuring method only achieve under the freeze state.

#### 6.4.3 Volume measurement

#### Volume can be measured in two ways.

1. Volume measurement utilizes 3-axle method to measure 3 groups of distances' dates to calculate from them.

To fulfill volume measurement must measure three times of distances and then press key to get the volume.

If the measured distance data is less than three groups, pressing will not get the volume; If four

distances information are measured, when pressing key, the result is the value of the first three data(D1, D2, and D3) to calculate the results.

Operation procedure: (take kidney as an example)

- 1. Get the kidney cross section image and profile section image respectively and freeze them.
- 2. Measure the long axis and short axis lengths of kidney cross section with distance measurement method,
- 3. Measure the longitudinal profile se of kidney's diameter.
- 4. Press key to fulfill volume measurement. The volume cost displays at the right corner of the screen behind "Vm1" as the following figure shows:



Figure 6-6. volume measurement sketch (3-axil method)

- 2. Utilizes ellipse measuring method to measure 2 groups of areas and then calculate them. Operation procedure: (take kidney as an example)
- 1. Get the kidney cross section image and profile section image respectively and freeze them.
- 2. Measure the kidney cross section circumference and area with circumference and area measuring method.

AMBISEN	User Manual	Document No.	AV-3618V
רושכופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	37/68

3. Press key to fulfill measurement. The volume cost automatically displays at the right corner of the screen behind "Vm1" as the following figure shows:



Figure 6-7 volume measurement sketch (Ellipse method)

Tips: If only one group of circumference and area are measured, Volume will not be displayed.

#### 6.4.4 Cardiac rate measurement (only at B/M modes)

- 1. At B/M modes, move trackball to change the sampling line position and freeze it ,and acquires a satisfying heart beat oscillogram.
- 2. Measure the distance between the three cardiac cycle's wave peaks according to distance measurement method. 4 groups data display at screen lower-right corner, from the left to right, they are: heart rate (HR) (unit: per rhythm /per minute.), rate of slope (EF) (unit: mm/s), time (T) (unit: ms), refresh velocity (ST) (unit: s). as following figure shows:




Figure 6-8. Heart rate measurement

Tips:

AV-3618V Digital Portable B mode ultrasound scanner    Effective Date   2008/02/21   Page   38/68	MBISER User Manual	Document No.	AV-3618V
Page 38/68  *Please identify the functions of			
Please identify the functions of key and key: Press key to delete  measuring marks and results, press key to clear off all annotation notes (including Doctor, name, age and sex) and measuring marks, results as well as report data.  • You can fulfill measurement by operating 1 1 2 1 3 4 direction keys as well as trackball.  6.5 Save and switch to storyboard view mode  The system provides separate memory area for "storyboard view" and "image save".  6.5.1 Saving a image  •When getting a satisfying image. press to save the current image. At the same time, the currestored image serial number automatically displays at left top corner of the image area. Such as "SAVING0! When the image's saving is finished, the serial number will disappear. After the saving, the system enters freeze mode, press key to return to real time state.  •The equipment can store 32 images at most, and these images will be automatically numbered by their savin order. For example: The number 01- 20 images has already been saved. press key, then the currently saved image will be numbered 21; When the memory is full (i.e. Already 32 images in it), the pressing key, the screen will display as below:  STORAGE IS FULL ERASE NO.01?  1. YES 2. NO  Tips replace the stored image number 01; press key 1 to replace it with current image; press key to exit this operation.  Select "2" No" to abandon saving image, then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  Tips:  •When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.	AV-3618V Digital Portable B mode ultrasound scanner		
**Rey and **Rey: Press **Rey to delete measuring marks and results, press **Rey to clear off all annotation notes (including Doctor, name, age and sex) and measuring marks, results as well as report data.  **You can fulfill measurement by operating **-1 12 13 4 direction keys as well as trackball.  6.5 Save and switch to storyboard view mode  The system provides separate memory area for "storyboard view" and "image save".  6.5.1 Saving a image  **When getting a satisfying image. press ** to save the current image. At the same time, the currestored image serial number automatically displays at left top corner of the image area. Such as "SAVING0! When the image's saving is finished, the serial number will disappear. After the saving, the system enters freeze mode, press ** trees* key to return to real time state.  **The equipment can store 32 images at most, and these images will be automatically numbered by their savin order. For example: The number 01- 20 images has already been saved. press ** key, then the currently saved image will be numbered 21; When the memory is full (i.e. Already 32 images in it), the pressing ** key, the screen will display as below:  **STORAGE IS FULL.ERASE NO.01?**  1. YES 2. NO  Tips replace the stored image number 01; press key ** 1 to replace it with current image; press key to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  **Tips:**  **When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.		Page	38/68
Note: The equipment can store 32 images at most, and these stared will be automatically numbered by their saving order. For example: The number 01: 20 images has already been saved.  **STORAGE IS FULL_ERASE NO.01?**  1. YES 2. NO  **When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  **Offer of the saving image in the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  **Offer of the saving image i		🕅	delete
*You can fulfill measurement by operating	-	•	cluding Doctor,
trackball.  6.5 Save and switch to storyboard view mode  The system provides separate memory area for "storyboard view" and "image save ".  6.5.1 Saving a image  •When getting a satisfying image. press by to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING	name, age and sex) and measuring marks, results as well as rep	oort data.	
The system provides separate memory area for "storyboard view" and "image save".  6.5.1 Saving a image  •When getting a satisfying image. press to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING	, , ,	→4 direction keys	as well as
The system provides separate memory area for "storyboard view" and "image save".  6.5.1 Saving a image  •When getting a satisfying image. press to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING	6.5 Save and switch to storyboard view mode		
•When getting a satisfying image, press  when getting a satisfying image, press  to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING	•	l "image save ".	
When getting a satisfying image, press to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING0! When the image's saving is finished, the serial number will disappear. After the saving, the system enters freeze mode, press key to return to real time state.  The equipment can store 32 images at most, and these images will be automatically numbered by their saving order. For example: The number 01- 20 images has already been saved, press key, then the currently saved image will be numbered 21; When the memory is full (i.e. Already 32 images in it), the pressing key, the screen will display as below:  STORAGE IS FULL.ERASE NO.01?  1. YES 2. NO  Tips replace the stored image number 01; press key  to exit this operation.  Select "2 No" to abandon saving image, then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  Tips:  When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image			
•When getting a satisfying image. press  to save the current image. At the same time, the current stored image serial number automatically displays at left top corner of the image area. Such as "SAVING0! When the image's saving is finished, the serial number will disappear. After the saving, the system enters freeze mode, press  key to return to real time state.  •The equipment can store 32 images at most, and these images will be automatically numbered by their savin order. For example: The number 01- 20 images has already been saved, press key, then the currently saved image will be numbered 21; When the memory is full (i.e. Already 32 images in it), the pressing  **STORAGE IS FULL.ERASE NO.01?**  1. YES 2. NO  Tips replace the stored image number 01; press key to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  **Tips:*  •*When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	6.5.1 Saving a image		
Tips replace the stored image number 01; press key to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  Tips:  When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	When the image's saving is finished, the serial number will disappear freeze mode, press key to return to real time state.  •The equipment can store 32 images at most, and these images will be order. For example: The number 01- 20 images has already been currently saved image will be numbered 21; When the memory is	automatically number saved, press	e system enters red by their savin key, then th
Tips replace the stored image number 01; press key to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  Tips:  When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	pressing key, the screen will display as below:		
Tips replace the stored image number 01; press key to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covimage number 02, reason by analogy in turn.  Tips:  •When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	STORAGE IS FULL.ERASE NO.01	?	
to exit this operation.  Select "2 No" to abandon saving image. then "Save" again, the system will remind you if you want to covi image number 02, reason by analogy in turn.  Tips:  •When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	1. YES 2. NO		
•When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	to exit this operation. Select "2 No" to abandon saving image. then "Save" again, the system		
•When the image memory is full, if some saved images are reloaded, then processing the saving operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	Tino		
operation, the system will remind you if you want to delete this loaded image and save the current there from this serial number.  6.5.2 View a saved image	·		
there from this serial number.  6.5.2 View a saved image	<ul> <li>When the image memory is full, if some saved images are r</li> </ul>	eloaded, then proces	ssing the saving
		loaded image and sa	ave the current
At real time or freeze, press "Load" key to display the following dialog box:	6.5.2 View a saved image		
	At real time or freeze, press "Load" key to display the following dialo	g box:	

AMRISEN	User Manual	Document No.	AV-3618V
רוביייי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	39/68

Please put in image number:

Input the saved image number, e.g. 01, then presses any key (except numbers and letters. Press backspace key to clear the mistakes and input again to load the first image). On the screen left bottom, 01/32 means: 01 is the serial number and 32 is the total amount of saved images in the memory. Press key to return to real time or freeze state. Repeat the above to load other images. Press steps to load other images. Tips: • You can not load any image from an empty memory area. • The input number should be from 01 to 32, when the input number is more than 32, the system will not take any operation. 6.5.3 Review storyboard mode After starting up, system has entered the real time state, first collect cine loop images; it will last for 30 seconds. Cine loop Ш key to review images in circulation. 1. Switch the equipment into freeze mode. Press 2. During review, press to enter manual review mode Press to play increase to play decrease by frame. And press image by frame, and press key to return to autoreview mode. 3. When review at "B", "B/B" mode, switch between "B/B" and "B" to play the images in difference windows. 4. To exit **cine loop**, press Tips: •If changing the image scanning mode, probe or display mode, CINE IOOD can not be processed immediately. Normally CiNe IOOD can be done 30 seconds after scanning.

•At CINE 1000 mode, roll trackball can play cine by single frame. When roll trackball up, cine is played from small number to big, roll it down, the reverse. Press "Cine loop" key to go on auto

AMRISEN	AMBISER User Manual	Document No.	AV-3618V
רוטכופווויז		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	40/68

playback.

•Replace probe or change the scan angle, the CİNE IOOD frame may change .The standard results depend on the display data on the left bottom of the screen.

#### 6.6 Obstetrical calculation

The device provides gestation age measurement for horse, cows, sheep, pigs, cats, dogs and other animals through parameters such as GS, BL, HL, SL, USD, HD, BD, CRL, etc. to automatically display the relevant gestation age (gestation weeks).

Operating steps:




Figure 6-9 the obstetrical table sketch

Repeat the above operation to select the required obstetrics table circularly (It can be showed the expected date of childbirth of DOG and CAT).

After that, select the distance parameter according to distance measurement method, the corresponding gestation age data will be automatically showed below "G·A=", the expected date of childbirth will be automatically showed below "EDD=", details are given below:


#### 1. EQUINE: Estimate horse gestation age according to GS

Check routine on horse:

- (1). Clean out the excreta in horse's rectum.
- ②. Palpate by hand and preliminary estimate the gestation state to determine the germinal tissue for B ultrasonic scanning.
- ③. Hold the probe, and insert it into the rectum with hand on the probe closely to ensure that your hand will feel any following situation firstly. Position your hand dorsal part between the rectum wall and the transducer.
- 4. The internal structure of horse is displayed on the screen. Bladder lies in the vertical crossed part, the

AMRISEN	User Manual	Document No.	AV-3618V
ווטכופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	41/68

following parts are cervix tissue and uterine body. Observe horizontally, uterine horn usually looks in round shape. Move the probe inside the rectum until the coupling section of uterine horn and uterine body is displayed, and then diverse the transducer to the direction of uterine horn as the following figure (Figure 6-10) shows:



- 1 Rectum
- 2 Uterine horns
- 3 Uterine bodies
- 4 Ovaries
- 5 Vaginas
- 6 Bladders

Figure 6-10. Probe position in womb and ovary sketch

⑤. The measuring method of GS diameter is given in the following figure (Figure 6-11). It can go in horizontal direction or vertical direction with no difference.




Figure 6-11. Horse's gestation age measurement sketch

- ⑥. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ". Through measurement, a chart can be made to determine growth curve and hereby to judge the gestation age corresponding to the embryo size. Gestation age calculation begins from the copulation day instead of the impregnation day.
- 2. Moo-cow BL: Estimate gestation age according to cow body length.


Check routine on cow:

- 1). Clean out the excreta in cow's rectum.
- ②. Palpate by hand and preliminary estimate the gestation state to determine the germinal tissue for B Copyright ©2010 by AMBISEA Technology Corp., Ltd All Right Reserved

AMRISEN	AMBISER User Manual	Document No.	AV-3618V
רובונויי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	42/68

ultrasonic scanning.

- ③. Hold the probe, and insert it into the rectum with hand on the probe closely to ensure that your hand will feel any following situation firstly. Position your hand dorsal part between the rectum wall and the transducer.
- ④. The internal structure of cow is displayed on the screen. Bladder lies in the vertical crossed part, the following parts are cervix tissue and uterine body. Observe horizontally, uterine horn usually looks in round shape. Move the probe inside the rectum until the coupling Section of uterine horn and uterine body displayed, and then diverse the transducer to the direction of uterine horn as the following figure (Figure 6-12) shows:



- 1 Rectum
- 2 Uterine horns
- 3 Uterine bodies
- 4 Ovaries
- 5 Vaginas
- 6 Bladders

Figure 6-12. Probe position in womb and ovary sketch

⑤. To measure fetal body diameter, a vertical Section should be selected, that is a plane covers from both sides through fetal neck, chest cavity and abdomen. Body diameter could be measured within 60 to 150 days gestation.

Body diameter measuring method is given in the following figure:




Figure 6-13. Cow's body length measurement sketch

⑥. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ".

#### 3. Moo-cow – SL: Estimate gestation age according to cow stomach length.

- (1). Make the cow in a state of stand.
- ②. Put the probe, a little bit left or right of the centre, on the ventral abdominal wall closely. If the is mud on this part, clean with water first incase the abdomen pelvic structure could not be displayed accurately.
- The screen should be big enough to display the maximum stomach longitudinal length. With the growth

AMRISEN	AMBISEN User Manual	Document No.	AV-3618V
רושכונו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	43/68

of gestation age, the fetal stomach macro-axis increases regularly. Measuring method of stomach maximum macro-axis is given in the following figure:




Figure 6-14. Cow's cardiac measurement sketch

④. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ".

#### 4. Moo-cow - HL: Estimate gestation age according to cow heart length.

- ①. Make the cow in a state of stand.
- ②. Put the probe, a little bit left or right of the centre, on the ventral abdominal wall closely. If the is mud on this part, clean with water first incase the abdomen pelvic structure could not be displayed accurately.
- ③. To measure the heart macro-axis, screen should display the maximal longitudinal axis of heart. With the growth of gestation age, the fetal heart macro-axis growths regularly. Measuring method of maximal macro-axis of heart is given in the following figure:




Figure 6-15. Cow's heart macro-axis measurement sketch

④. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ".

#### 5. SWINE-HL: estimate gestation age according to pig heart length

Check routine on pigs:

- ①. Make the cow in a state of stand.
- ②. Put the probe, a little bit left or right of the centre, on the ventral abdominal wall closely along the side of teats and skull to rear leg. If the is mud on this part, clean with water first incase the abdomen pelvic structure could not be displayed accurately.

AMRISEN	AMBISEN User Manual	Document No.	AV-3618V
ווייבופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	44/68




Figure 6-16. Pig's gestation age measurement position sketch

③. To measure the heart macro-axis, screen should display the maximal longitudinal axis of heart. With the growth of gestation age, the fetal heart macro-axis increase regularly. Measuring method is given in the following figure:




Figure 6-17. Swine heart macro-axis measurement sketch

④. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ".

#### 6. SHEEP-USD: Estimate gestation age according to hilum-spine length of sheep

There are two methods to check pregnant ewe with B ultrasound:

Use convex probe or linear probe check the abdomen; and recta probe for rectum. These two methods are both quite useful to diagnose pregnancy. On the base of some publication, we agree that both the methods on pregnancy diagnosis are proved to be effective.

- With 35 days pregnancy, result of rectal way is more precise than that of abdominal way;
- During 35 to 70 days pregnancy, both are equally effective.
- After 70 days pregnancy, abdominal method would be better. This way would be more practical for the Bigger womb. Abdomen check
  - ①. Abdomen check can be done when the ewe is standing, lying or squatting. Put the probe on the middle abdomen without fur and is set to be the proper position for checking in advance.
  - ②. If the is mud on this part, clean with water first incase the abdomen pelvic structure could not be displayed accurately.

AMBISEN	User Manual	Document No.	AV-3618V
ווישכונו		Revision A	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	45/68

- ③. Measure the hilum-spine length.
- ④. Measure selected parameter distance according to distance measurement method, the corresponding gestation age data will automatically shows behind " G·A ".

#### 7. CAT - HD: Estimate gestation age according to cat head diameter

To measure the head diameter is to measure the interior maximum diameter of skull between ventral and dorsal part.

Measuring method of head diameter is given in the following figure:



Figure 6-18. Cat head diameter measurement sketch

#### 8. CAT-BD: Calculate the gestation age according to cat BD

After fetal head formed, binary top diameter measurement becomes a routine in ultrasonic examination. The measuring method is:

According to distance measurement method to measure distance of selected parameters, the corresponding gestation age data will automatically show behind " G·A".

## 9. DOG-GSD: Calculate the gestation age according to dog gestation saccus diameter

The method is the same as that of horse.

A. DOG-CRL: Calculate the gestation age according to dog CRL

The method is the same as the cow's.

#### B.DOG-HD: Calculate the gestation age according to dog HD

The method is the same as the cat's.

# C. DOG-BD: Calculate the gestation age according to BD

The method is the same as the cat's.

After measurement, press "key to clear the screen for the next measurement. Otherwise.

Press "D" to select "D. Exit" to exit OB table.

AMBISEN	User Manual	Document No.	AV-3618V
וועוווווווווווווווווווווווווווווווווווו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	46/68

Tip: the system would not show the relative items when the distance is less than the measurements defined as the table below:

EQUINE - GSD	D1<6mm
BOVINE - BL	D1<8mm
BOVINE - SL	D1<1mm
BOVINE - HL	D1<3mm
SWINE - HL	D1<31mm
SHEEP - USD	D1<15mm
CAT - HD	D1<15mm
CAT - BD	D1<17mm
DOG - GSD	D1<1mm
DOG - CRL	D1<1mm
DOG - HD	D1<14mm
DOG - BD	D1<16mm

## 6.7 Diagnostic Report

The information of patient, diagnostic comments, measuring results as well as hospital, date, time, doctor will be saved in the report. System will generate a routine report page automatically according to different measurement items. Therein, distance, circumference, area, volume, ect. will be saved in Routine Report; measuring results such as gestation age, expected date of delivery, fetus weight, etc. will be saved in Obstetrics Report

●At any of the B, BB, BM, M, 4B mode, when fulfilling patient information comment, diagnoses, and foci

distance, circumference, area, volume measurements, keep the image in freeze state, press key to switch to Routine Report interface as given bellow:




Figure 6-19. Routine report sketch

Tips:

AMRISEN	AMBISER User Manual	Document No.	AV-3618V
ווטכופוווו		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	47/68

• At report interface, press key to add comments in comment frame but the other information can not be amended.

• Press key repeatedly to switch between diagnoses interface and report interface.

## 6.8 Image Processing

Press key to display the image processing menu (as below). It provides image positive-negative flip, GAMA calibration, statistic histogram, edge gain and other image processing functions.



Press to select "1. Posi-Nega" to change the image polarity as given below:





Figure 6-20. Polarity sketch

Press to select "2. GAMA calibration" to lighten up "GAMA", on the top right corner of the screen.

There are 4 levels G1, G2, G3 and OFF, press to select them circularly.

Tips:

• This function is set for video output and external devices (e.g. video printer) adjustment and is recommended to set at "OFF" during operation.

AMRISEN	User Manual	Document No.	AV-3618V
ווייכופוווי		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	48/68

Press to select "3. Histogram" to process histogram statistics as the following procedure:

At freeze, press " key to display the menu, select "3. Histogram" to call out a rectangle window on the image, press "Measurement" key to calculate the gray scale pixel numbers. The result displays on screen right as the following figure shows:




Figure 6-21. Histogram sketch

Horizontal ordinate indicates gray scale degree, vertical ordinate indicates number.

PT indicates total quantity of pixel dots in the rectangular window.

Gm indicates the corresponding gray scale of the curve at the vertical peak.

Pm indicate the quantity of pixel dots at Gm gray scale.

From the above figure, in the rectangular area, the total number of pixel dots is 10000. At dray scale 28, there are 28, the most image pixel dots.

Press to select "4. Edge enhancement" to lighten up the "EDGEEN" on the top right corner of the screen. Then edge enhancement is adjustable. There are 9 levels: OFF, EE1-EE8 which can be adjusted by pressing.

#### Tips:

• This function is a special item and recommended to set at "OFF" during operation.

Press to select "5. Exit" to exit the menu.

## 6.9 Image print

Connect the VIDEO IN port of the video printer to the VIDEO OUT port on rear panel of the equipment, then operate according to the video printer operation manual.

AMBISEN	User Manual	Document No.	AV-3618V
ָר הוביריי			A0
AV-3618V Digital P	Effective Date	2008/02/21	
AV 3010V Digital Fortable B mode altrasound scanner		Page	49/68

## 6.10 Image Upload to computer

Connect the USB communication port of the equipment to the computer's USB port.

The high speed USB2.0 mouthpiece can upload images to the computer at current time.

In the accompanying disc, there are USB device driver, image preview program Amcap and SDK documents for relevant redevelopment.

Install the USB device driver, open Amcap.EXE, and then you can view images according to the instruction files in the disc.

Using SDK documents in the disc, you can process redevelopment as required.

# 6.11 Shut down the equipment

Turn the power off.

Tips: It would be better to unplug the AC supply if do not use the equipment for a long time.

#### Note

Must not plug or unplug the power plug when the equipment is not shut down; If it needs to switch on the equipment immediately after just shutdown, please wait for 2 or 3 minutes to avoid breaking the equipment.

AMBISEN	User Manual	Document No.	Av-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	50/68

#### 7 Check and Maintenance

#### 7.1 Detect

Periodic detect the device power cord and probe cable, if find any damage or breakage, must not use the device and replace it immediately.

#### 7.2 Effective Service life

Bases on the manufacturer's design, production related files, this model's use life is six years. The Product's material will gradually aging, if the product continually use over the designed use life, it may bring the problem of the performance reduced and fault rate raise.

Note

To Discard the device according to the local law.

Do not discard it mixing with other household garbage.

#### Warning

The manufacturer shall not assume the responsibility of risks caused by using the device beyond its service life.

#### 7.3 Host maintenance

Instrumentation environment should accord with "4.1 operation environmental requirement".

If device enclosure needs cleaning, shutdown the device first and then wipe with alcohol sponges. And only at the switch off condition open the device shell to use the blower to get rid off the dust.

When the device does not work for a long time, should pack the device according to the instructions in the package. Store it properly in the warehouse. The storage environment should accord with "8.1 Transportation and storage environmental requirements".

#### 7.4 Probe maintenance

Probe is an expensive and frangible part. Never hit it or drop it on floor. When stop the working, place the probe in the box





key to set the probe at freeze state.

Please use medical ultrasound coupling gel during diagnoses. The water prevention level is IPX7, so do not let water exceed the acoustic window (see figure 7-1, 7-2, 7-3, 7-4, 7-5). Do daily inspection on the probe enclosure to see if it is cracked and avoid liquid leakage to spoil the inner components

AMBISEN	User Manual	Document No.	Av-3618V for vet
ווסוסבוו	Oser manuar	Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	51/68




Figure 7-1.C1-7/60R/3.5MHz convex probe water-proof sketch



Figure 7-2. EC1-2Endocavity probe water-proof sketch




Figure 7-3.L1-3/7.5MHz HF linear probe water-proof sketch




Figure 7-4.C1-8/20R/5.0MHz convex probe water-tight sketch



Figure 7-5 EL1-3/7.5MHz rectal probe water-proof sketch

Once the probe is connected to the main unit, do not remove it at will to avoid poor contact.

AMBISEN	User Manual	Document No.	Av-3618V for vet
		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	52/68

#### Note

Probe might be damaged due to long time covered coupling gel.

#### Note

Clean the probe head after every use.

Do not clean the probe with a surgical brash neither soft brash. Only soft cloth can be used to clean it.

#### Note

Do not press the probe on the patient too long to avoid discomfort.

## 7.5 Cleansing

- 1. When the enclosure needs cleaning, wipe it with soft dry cloths and then wipe gently with sponge dipped with 75% medical alcohol.
- 2. When the inner part of the equipment needs cleaning, power off the equipment first and open the enclosure and vacuum it.

#### Warning

To prevent accidents, please power off the equipment when cleaning it or the probe.

#### Caution

- 1. Please refer to instructions prescribed by the manufacturer closely when using detergents.
- 2. Be careful with cleaning of the display, because it is very easy to scratch and spoil. Please wipe it with dry soft cloth.
- 3. Please do not clean the inner part of the device.
- 4. Please do not place the device in liquid.
- 5. Do not leave any detergent on the device surface.
- 6. Though there will be no chemical reaction between the device enclosure and most of those detergents, We still suggest no detergent in cleaning lest the device surface is spoiled.

AMDIC CN	User Manual	Document No.	Av-3618V for vet
AMBISEN User Manual		Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	53/68

#### Warning

Must not use extender, ethylene oxide or any other organic solvent which end to deface the probe's protective foil.

Must not place the probe in liquid or detergent.

Must keep the equipment or probe from any type of liquid's infiltration.

Must not clean device or probe by airing or heating.

# 7.6 Correctly use the probe

In order to extend probe's service life and obtain optimum performance, follow the instructions as below:

- 1. Periodic inspection on probe cable, socket and acoustic window.
- 2. Shutdown the device first and then connect or disconnect the probe.
- 3. Do not drop probe or flint body, and never hit the probe acoustic window, otherwise probe should be damaged.
- 4. Put the probe in the probe box when it is not in use.
- 5. Never heat the probe.
- 6. Never bend or pull probe cable, otherwise the internal connection should be broken.
- 7. Use coupling gel only on probe header and then clean probe.
- 8. Inspect probe acoustical window, enclosure and cable seriously after probe cleaning. Never use the probe again if any crack or breakage is found.

#### 7.7 Equipment detection and calibration

1. Check the leakage current of the device annually referring to the following data.

Test Items			Standard Requirements
Continuous leakage	Leakage current	Normal	≤100
current under normal	to Shell	Single Malfunction	≤500
temperature (uA)	Leakage current	Normal	≤100
	to Patient	Single Malfunction	≤500
Dielectric endurance temperature (V)	under normal	A- a <sub>2</sub>	4000V/1min No flashover No

2. Test the software of obstetrics, area, and circumference measurement; please refer to Appendix C for detailed data.

AMBISEN	User Manual	Document No.	Av-3618V for vet
		Revision	A0
AV 2610V Digital	Dortoble P mode ultracound scanner	Effective Date	2008/02/21
AV-3618V Digital Portable B mode ultrasound scanner		Page	54/68

# 8 Transportation and Storage

# 8.1 Environmental Requirements of Transportation and Storage

Environmental temperature: - 40°C ~+55°C;

Relative humidity: 10%~100%;

Atmospheric pressure: 50KPa~106KPa;

# 8.2 Transportation

All marks printed on the packing box are in accordance with the requirements of GB191-2000 <Packing, Storage and Transportation Marks>. Simple shockproof installations are made in the packing box, which adapt to voyage, railway, and highway and steamship transportation. Avoid rain, inversion and impact.

# 8.3 Storage

- When the device is stored for more than 6 months, take out the device from the packing box. Electrify the device for 4
  hours, and then put the device into the packing box according to the instructions way on the packing box. Do not
  crisscross/piling devices or keep the device to floor, wall and ceiling.
- Keep the repository ventilative. Avoid direct sunshine and caustic gas.

AMBISEN	User Manual	Document No.	Av-3618V for vet
ן או ופוסבגו	Osei manuai	Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	55/68

# 9 Simple Malfunction Examination and Troubleshooting

## 9.1 Examination

- Examine whether the supply power is in normal state. The power line of the main frame has been properly connected and inserted in the electric socket.
- Examine whether the probe has been connected to the device correctly.


# 9.2 Troubleshooting

## • Replacement of fuse

Position the driver into the slot on the fuse cap and press it, then contra-rotate it to loose the cap. Take off the fuse tube (protective tube) and replace it then put back the cap and take the reverse measure to fix the cap and fasten it. Fuse specification:  $\phi 5 \times 20$ ,F2L250V.

## Troubleshooting (see the following diagram)

No.	Malfunction phenomena	Troubleshooting	
1	Open power switch of the device, with no	1. Examine power supply;	
	signal appeared on the screen and the	2. Examine supply line and plug;	
	indicator light is not on.	3. Examine whether the fuse is burned-out;	
		4. Examine the lightness adjusting knob of the	
		monitor.	
2	Discontinuous stria and snow appear on the	1. Examine power supply: strike fire interfere of the	
	screen.	other devices;	
		2. Environmental examination: Electric and	
		magnetic field interfere;	
		3. Examine power and probe plugs: whether they are	
		well connected.	
3	Unclear image display	1.Adjust the brightness and contrast knobs on	
		the front panel;	
		2.Adjust the 8 level TGC and the total gain knob	
		on the front panel;	
		3. Clean the screen optical filter.	
4	Unclear near field	1. Adjust the total gain knob on the front panel	
		and the 8 level TGC;	
		2. Adjust the focus position to the near field;	
5	Unclear far field	1. Adjust the total gain knob on the front panel	
		and the 8 level TGC;	
		2. Adjust the focus number, space and	
		position to set the focus in the far field.	



# **Appendix A Labels Diagram**

1. AV-3618V Main unit nameplate diagram:

Item	B-Ultrasonic Diagnostic Equipment	Model	AV-3618V
Power supply	AC110-130V/220-240V,50Hz/60Hz	Power Consumption	100VA
Manufactured		Serial No.	
Manufacturer		Safety	Class   Type B
ADD	5A,ShunTian Bldg,62 GuiMiao Road,Nanshan District,Shenzhen City,China		

# AV-3618V

# **B-Ultrasonic Diagnostic Equipment**

2. AV-3618V device probe nameplate diagram:



3. AV-3618V device optional probe nameplate diagram:






L1-3/7.5MHz

EL1-3/7.5MHz



AMBISEN	User Manual	Document No.	Av-3618V for vet
(1) IRIDELI Oser Mandai		Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	57/68

# 4.AV-3618V Packing diagram:



Above contents are subject to change without prior notice.

AMBISEN	User Manual	Document No.	Av-3618V for vet
(1) IBIDELI OSEL Mandai		Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	58/68

# **Appendix B Acoustic Output Parameters**

Test Mode: B-mode Transducer Type: C1-7/60R

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power (mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	72.4	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.72	
3	lob(mW/cm2)	Output beam intensity. IEC 61157 4.2.2c)	18.04	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived Intensity. IEC 61157 4.2.2d)	9.56	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:57% Frame rate:34	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	40	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 1.59 (□): 8.50	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	Not applicable	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	34.13	Test
10	Output beam dimensions (cm2)	Dimensions parallel(  ) and perpendicular(\(\perpendicular\)) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken As the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	4.013	frequen cy: 2.5 MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working Frequency. IEC 61157 4.2.2j)	2.78	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2I)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

AMBISEN	1RISER User Manual		Av-3618V for vet
וון ופוסכרו	Oser Manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	59/68

Test Mode: B/M-mode Transducer Type: C1-7/60R

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power (mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	64.6	
2	P <sub>-</sub> (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.72	
3	I <sub>ob</sub> (mW/cm <sup>2</sup> )	Output beam intensity. IEC 61157 4.2.2c)	16.09	
4	I <sub>spta</sub> (mW/cm²)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	33.89	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:57% Frame ate:38 ST:1.00s	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	40	
7	W <sub>pb6</sub> (mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(∥): 1.59 (⊥): 8.55	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	0.222	Test
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	37.31	frequency: 2.5 MHz
10	Output beam dimensions (cm²)	Dimensions parallel(    ) and perpendicular( \( \) ) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	4.013	
11	f <sub>awf</sub> (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	2.78	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2I)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initializatio n mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

AMBISEN	User Manual		Av-3618V for vet
אסכומו וגי	Oser mailuai	Revision A0	Α0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	60/68

Test Mode: B-mode Transducer Type: C1-8/20R

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power(mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	34.8	
2	P_ (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.82	
3	I <sub>ob</sub> (mW/cm <sup>2</sup> )	Output beam intensity. IEC 61157 4.2.2c)	20.17	
4	I <sub>spta</sub> (mW/cm²)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	25.74	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:49% Frame rate:71	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	33	
7	W <sub>pb6</sub> (mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(    ): 2.55 ( ⊥ ): 3.43	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	Not applicable	Test
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	67.93	frequency: 5.0
10	Output beam dimensions (cm <sup>2</sup> )	Dimensions parallel(  ) and perpendicular(\(\perp \)) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	1.725	MHz
11	f <sub>awf</sub> (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	3.89	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable" (n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

AMBISEN User Manual		Document No.	Av-3618V for vet
וזסכנו ווי)	Oser manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	61/68

Test Mode: B/M-mode Transducer Type: C1-8/20R

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power(mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	31.4	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.82	
3	lob(mW/cm 2)	Output beam intensity. IEC 61157 4.2.2c)	18.20	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	46.85	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:49% Frame rate:66 ST:1.00s	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral.  IEC 61157 4.2.2f)	33	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 2.12 (□): 3.99	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	0.228	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	62.50	
10	Output beam dimensions (cm2)	Dimensions parallel(    ) and perpendicular( \( \_ \)) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	1.725	Test frequency: 5.0 MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	3.89	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a).  IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

Copyright ©2010 by AMBISEA Technology Corp., Ltd - All Right Reserved

AMBISEN	User Manual	Document No.	Av-3618V for vet
ווסוסבוו	Oser manual	Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	62/68

Test Mode: B/M-mode Transducer Type: EC1-2

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power (mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	29.4	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.94	
3	lob(mW/cm2)	Output beam intensity. IEC 61157 4.2.2c)	21.41	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	39.01	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:49% Frame rate:AO	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	26	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 1.47 (□): 2.41	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	Not applicable	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	95.05	Test frequenc
10	Output beam dimensions (cm2)	Dimensions parallel(    ) and perpendicular( \( \) ) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group.  IEC 61157 4.2.2i)	1.373	y: 6.5 MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	4.39	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable" (n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no".  IEC 61157 4.2.20)	Yes	

AMBISEN	User Manual	Document No.	Av-3618V for vet
וזסכנו וני)	Oser manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	63/68

Test Mode: B/M-mode Transducer Type: EC1-2

	·	rest wode: b/w-mode transducer ry	JO. 20. 2	
Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power(mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	25.6	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	1.94	
3	lob(mW/cm 2)	Output beam intensity. IEC 61157 4.2.2c)	18.64	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	59.31	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Angle:49% Frame rate:AO ST:1.00s	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	26	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 1.46 (□): 2.34	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	0.233	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	93.45	Test frequency: 6.5
10	Output beam dimensions (cm2)	Dimensions parallel(  ) and perpendicular(\(\perp \)) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may Be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	1.373	MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	4.40	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

AMBISEN	User Manual		Av-3618V for vet
וזסכנו וני)	USEI Manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	64/68

Test Mode: B-mode Transducer Type: L1-3

Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power (mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	26.4	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	2.63	
3	lob(mW/cm2)	Output beam intensity. IEC 61157 4.2.2c)	17.86	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	11.81	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Frame rate:58	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral.  IEC 61157 4.2.2f)	21	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 0.78 (□): 2.57	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	Not applicable	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	55.24	Test
10	Output beam dimensions (cm2)	Dimensions parallel(  ) and perpendicular(\( \_ \)) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	1.478	frequency: 7.5 MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	5.20	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable" (n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initialization mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no".  IEC 61157 4.2.20)	Yes	

AMBISEN	FR User Manual		Av-3618V for vet
וזשכומו וני	OSEI Manuai	Revision	Α0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	65/68

Test Mode: B/M-mode Transducer Type: L1-3

rest wiode: b/wi-mode transducer type: L1-3				
Nº	Test Item	Requirements and Clauses	Results	Remarks
1	Maximum power (mW)	Maximum temporal-average power output. IEC 61157 4.2.2a)	27.6	
2	P- (MPa)	Peak-negative acoustic pressure. IEC 61157 4.2.2b)	2.63	
3	lob(mW/cm 2)	Output beam intensity. IEC 61157 4.2.2c)	18.67	
4	Ispta (mW/cm2)	Spatial-peak temporal average derived intensity. IEC 61157 4.2.2d)	29.60	
5	System settings	Ultrasound instrument console settings. IEC 61157 4.2.2e)	Focus number:1 Power:80% Frame rate:58 ST:1.00s	
6	Lp (mm)	Distance form the transducer output face to the point of maximum pulse-pressure-squared integral. IEC 61157 4.2.2f)	21	
7	Wpb6(mm)	-6dB pulse-width. IEC 61157 4.2.2g)	(□): 0.98 (□): 4.10	
8	prr(kHz)	Pulse repetition rate. IEC 61157 4.2.2h)	0.232	
9	srr (Hz)	Scan repetition rate. IEC 61157 4.2.2h)	55.55	Test frequenc
10	Output beam dimensions (cm2)	Dimensions parallel( $\parallel$ ) and perpendicular( $\perp$ ) to the reference direction shall be specified. For scanning modes, these shall refer to the central scan line only. In many cases, especially contact systems, these dimensions may be taken as the geometrical dimensions of the ultrasonic transducer or ultrasonic transducer element group. IEC 61157 4.2.2i)	1.478	y: 7.5 MHz
11	fawf (MHz)	Arithmetic-mean acoustic-working frequency. IEC 61157 4.2.2j)	5.20	
12	APF	Acoustic power-up fraction. IEC 61157 4.2.2k)	0%	
13	Power-up mode	In system in which the user defines the power-up mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2l)	B-mode freeze	
14	AIF	Acoustic initialization fraction. IEC 61157 4.2.2m)	100%	
15	Initializatio n mode	If appropriate. in systems in which the user defined the initialization mode, this shall be stated as either "user defined" or "not applicable"(n / a). IEC 61157 4.2.2n)	B-mode	
16	Acoustic output freeze	If the system has acoustic output freeze then this shall be stated as "yes", otherwise it shall be stated as "no". IEC 61157 4.2.20)	Yes	

AMBISEN	User Manual	Document No. Av-3618V for vet	Av-3618V for vet
וזסכנו וני)	Oser manual	Revision	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	66/68

# **Appendix C Obstetric table**

**Gestational Table 1: Equine** 

Measurement(mm)	Week	D
(Gestational Sac Diameter)	week	Day
6	1	4
8	1	4
10	1	5
12	1	6
14	1	6
16	2	0
18	2	0
20	2	1
22	2	2
24	2	3
26	2	5
28	4	1
30	4	2
32	4	3
34	4	4
40	5	0
42	5	2
44	5	3
46	5	4
48	5	5
50	5	6
52	6	1
54	6	2
56	6	3

All measurements +/- 3 days

AMBISER User Manual		Document No.	Av-3618V for vet
וזסכנו וני)	Oser manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	67/68

# **Gestational Table 2: Bovine**

Measurement (mm)	Week	Day
(Body Length)	- week	
8	4	0
10	5	0
12	5	1
14	5	2
16	5	3
18	5	5
20	5	5
22	5	6
24	5	6
26	6	1
28	6	1
30	6	1
32	6	2
34	6	3
36	6	3

All measurements +/- 3 days

AMBISEN	AMRISER User Manual		Av-3618V for vet
וו וסוסבוו	Oser Manual	Revision A0	A0
AV-3618V Digital Portable B mode ultrasound scanner		Effective Date	2008/02/21
		Page	68/68

# **Gestational Table 3: Sheep**

Measurement(mm)	Wools	Dow
(Umbilicus to Spine	Week	Day
15	7	1
18	7	3
21	7	6
24	8	1
27	8	4
30	9	0
33	9	2
36	9	4
39	10	0
42	10	2
45	10	5
48	11	3
51	11	5
54	12	1
57	12	2
60	12	4
63	12	6
66	13	2
69	13	4
72	14	2
75	14	4
78	15	0
81	15	2
84	15	5
87	16	6
90	17	0
93	17	1
96	17	3
99	17	6

All measurements +/- 3 days