

SuperScript[®] III Platinum[®] CellsDirect Two-Step qRT-PCR Kit with SYBR[®] Green

**For two-step real-time quantitative RT-PCR from
cell lysate using SYBR[®] Green I fluorescent dye**

Catalog Nos. 11738-060 and 11738-068

Rev Date: 28 June 2010

Manual part no. 250751

MAN0000472

Table of Contents

Kit Contents and Storage.....v

Introduction.....1

Lysing Cells4

First-Strand cDNA Synthesis8

qPCR — Guidelines and Recommendations.....9

qPCR — Instruments Using PCR Tubes/Plates11

qPCR — Roche LightCycler®13

Troubleshooting.....15

References17

Related Products.....18

Purchaser Notification19

Technical Service22

Kit Contents and Storage

Shipping and Storage Kit components are shipped on dry ice and should be stored at -20°C. ROX Reference Dye should be stored in the dark.

Kit Size and Modules Each kit includes a cDNA Synthesis module and a qPCR module. Each cDNA synthesis reaction provides enough cDNA for multiple qPCR reactions.

<u>Kit Catalog Number</u>	<u>Number of Reactions</u>	
	<u>cDNA Synthesis</u>	<u>qPCR</u>
11737-060	25	100
11737-068	100	500

cDNA Synthesis Module	<u>Component</u>	<u>100-rxn kit</u>	<u>500-rxn kit</u>
	Resuspension Buffer	2 × 250 µl	2 × 1 ml
	DNase I (1 U/µl)	125 µl	500 µl
	10X DNase I Buffer	40 µl	160 µl
	25-mM EDTA	120 µl	400 µl
	RT Enzyme Mix (contains SuperScript® III RT, 100 units/µl; and RNaseOUT™ Recombinant Ribonuclease Inhibitor, 20 units/µl)	50 µl	200 µl
	2X RT Reaction Mix*	500 µl	2 × 1 ml
	<i>E. Coli</i> RNase H (2 U/µl)	30 µl	100 µl
	HeLa Total RNA (10 ng/µl)	10 µl	10 µl
	Lysis Enhancer	125 µl	500 µl
	*Oligo(dT) ₂₀ (2.5 µM), random hexamers (2.5 ng/µl), 10 mM MgCl ₂ , and dNTPs		

qPCR Module	<u>Component</u>	<u>100-rxn kit</u>	<u>500-rxn kit</u>
	Platinum® SYBR® Green qPCR SuperMix-UDG*	2 × 1.25 ml	12.5 ml
	50-mM Magnesium Chloride (MgCl ₂)	1 ml	2 × 1 ml
	20X Bovine Serum Albumin (BSA), UltraPure (1 mg/ml)	300 µl	1.25 ml
	ROX Reference Dye	100 µl	500 µl
	*SYBR® Green I, 60 U/ml Platinum® Taq DNA polymerase, 40-mM Tris-HCl (pH 8.4), 100-mM KCl, 6 mM MgCl ₂ , 400-µM dGTP, 400-µM dATP, 400-µM dCTP, 400-µM dUTP, 40 U/ml UDG, and stabilizers.		

Continued on next page

Kit Contents and Storage, continued

Materials Supplied by the User

The following additional items are required for use with this kit:

- Coulter Counter or hemacytometer
 - Microcentrifuge
 - qPCR instrument
 - Trypsin (for adherent cell cultures only)
 - 1X cold phosphate-buffered saline PBS, without calcium or magnesium
 - 0.2-ml thin-walled PCR tubes or 96-well PCR plates
 - Ice
 - Pipettes
 - Disposable gloves
-

Introduction

System Overview

The SuperScript® III Platinum® CellsDirect Two-Step qRT-PCR Kit with SYBR® Green is an optimized kit for synthesizing first-strand cDNA directly from mammalian cell lysate without first isolating RNA, and then amplifying the cDNA in a real-time quantitative PCR (qPCR) reaction using Platinum® SYBR® Green qPCR SuperMix-UDG.

In traditional qRT-PCR, RNA is first isolated from cells in a time-consuming procedure that can lead to a loss of material. Using the CellsDirect cDNA Synthesis System, the cells are lysed and the cDNA is generated from the lysate in a single tube with minimal handling and no sample loss. DNase I is added to eliminate genomic DNA prior to first-strand synthesis. After synthesis, the first-strand cDNA can be transferred directly to the qPCR reaction without intermediate organic extractions or ethanol precipitations.

This kit has been optimized for small cell samples, ranging from 10,000 cells down to a single cell. The use of SuperScript® III Reverse Transcriptase ensures high specificity and high yields of cDNA from small amounts of starting material—as little as 10 pg total RNA. The use of Platinum® SYBR® Green qPCR SuperMix-UDG ensures optimal qPCR performance using SYBR® Green I dye, with excellent sensitivity and a linear dose response over a wide range of target concentrations.

Advantages of the Kit


This kit has the following advantages:

- Compatible with a wide range of mammalian cell types grown under different treatment conditions
- Cell lysis and first-strand cDNA synthesis in the same tube minimizes reagent loss, sample loss, and handling time
- Total lysate volume is used in the cDNA synthesis reaction, providing greater yields with a limited number of cells and allowing for detection of rare transcripts
- SuperScript® III Reverse Transcriptase, with reduced RNase H activity and higher thermal stability, produces high yields of cDNA in the first-strand synthesis reaction for greater sensitivity and enhanced detection of rare transcripts
- Platinum® SYBR® Green qPCR SuperMix-UDG ensures optimal sensitivity and performance in qPCR using SYBR® Green I fluorescent dye, with built-in carryover contamination protection and a linear dose response over a wide range of target concentrations

Continued on next page

Introduction, Continued

Diagram of cDNA Synthesis from Cell Lysate



SuperScript® III RT

SuperScript® III Reverse Transcriptase is an engineered version of M-MLV RT with reduced RNase H activity and increased thermal stability. The enzyme can be used to synthesize first-strand cDNA at temperatures up to 55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases.

Because SuperScript® III RT is not inhibited significantly by ribosomal and transfer RNA, it can effectively synthesize first-strand cDNA directly from total RNA. The concentration of SuperScript® III RT in this system has been optimized to synthesize first-strand cDNA from total RNA in cell lysate.

Continued on next page

Introduction, Continued

Platinum® SYBR® Green qPCR SuperMix-UDG

Platinum® SYBR® Green qPCR SuperMix-UDG is a ready-to-use reaction cocktail containing all components, except primers, for the amplification and detection of DNA in qPCR. It contains SYBR® Green I fluorescent dye, Platinum® *Taq* DNA polymerase, Mg⁺⁺, uracil-DNA glycosylase (UDG), proprietary stabilizers, and deoxyribonucleotide triphosphates (dNTPs), with dUTP instead of dTTP. The concentration of the SuperMix allows for the addition of primers and template.

SYBR® Green I is a fluorescent dye that binds directly to double-stranded DNA (dsDNA). In qPCR, as dsDNA accumulates, the dye generates a signal that is proportional to the DNA concentration and that can be detected using qPCR instruments. SYBR® Green I in this SuperMix formulation can quantify as few as 10 copies of a target gene in as little as 1 pg of template DNA or RNA. It has a broad dynamic range of six orders of magnitude, and is compatible with melting curve analysis.

Platinum® *Taq* DNA polymerase is precomplexed with specific monoclonal antibodies that inhibit polymerase activity during reaction assembly at room temperature. Full polymerase activity is restored after the denaturation step in PCR cycling, providing an automatic “hot start” in PCR and thereby increasing amplification efficiency, sensitivity, and yield.

UDG and dUTP are included in the mixture to prevent the reamplification of carryover PCR products between reactions. dUTP in the mix ensures that any amplified DNA will contain uracil. UDG, or uracil-N-glycosylase, removes uracil residues from single- or double-stranded DNA, preventing dU-containing DNA from serving as template in future PCRs. Incubation of subsequent PCRs with UDG before cycling destroys any contaminating dU-containing PCR product from previous reactions. After this decontamination step, UDG is inactivated by the high temperatures during normal PCR cycling, thereby allowing the amplification of genuine target sequence(s).

Control RNA

HeLa Total RNA is included in the kit as a control. The concentration of HeLa Total RNA provided (10 ng/μl) is equivalent to 1,000 cells.

Methods

Lysing Cells

Introduction

In this step, you lyse your cells in Resuspension Buffer and Lysis Enhancer and perform a DNase I digestion to remove genomic DNA from the sample.

Cell Types and Density

This kit has been optimized for small cell samples, ranging from 1 to 10,000 cells. This kit is compatible with several different mammalian cell lines, including HeLa, COS-7, 293, Jurkat, CV1, and K562. Cells may be grown under a variety of conditions and treatments. Any type of culture vessel can be used.



Important

- We recommend using a maximum of 10,000 cells per reaction. Higher numbers of cells may inhibit reverse transcription and result in reduced yields and/or truncated cDNA product.
- Make sure that all solutions and equipment that come in contact with the cells are sterile. Always use proper sterile technique and work in a laminar flow hood when handling cells.

Required Materials

The following materials are provided by the user:

- Mammalian cell cultures in growth media
- Coulter Counter or hemacytometer
- Centrifuge (for pelleting cells)
- Incubator, water bath, or thermal cycler preheated to 75°C
- Trypsin (for adherent cell cultures only)
- 1X cold phosphate-buffered saline (PBS), without calcium or magnesium
- 0.2-ml thin-walled PCR tubes or 96-well PCR plates
- Ice
- Pipettes

The following materials are provided in the kit:

- Resuspension Buffer
- Lysis Enhancer
- DNase I, Amplification Grade (1 U/ μ l)
- 10X DNase I Buffer
- EDTA, 25 mM
- Optional: Control HeLa Total RNA



Note

All steps should be performed on ice, and reagents should be chilled and/or thawed immediately prior to use. The incubator should be **preheated** to 75°C.

Continued on next page

Lysing Cells, continued

Lysing Adherent Cells or Cells in Suspension

Use the following lysis procedure for adherent cell cultures in vessels larger than 24-well plates. For cells in suspension, skip Steps 1–4 and proceed to Step 5 below.

1. Add enough trypsin to cover the adherent cells in your tissue culture dish, plate, or flask (e.g., for a 10-cm dish, use ~1 ml; for a T75 flask, use ~3 ml).
2. Incubate for 5 minutes at room temperature or in a 37°C incubator.
3. Check for cell detachment under a microscope. If cells have not detached, gently tap the dish or flask to dislodge the cells, or let the cells incubate longer, checking them every minute under a microscope.
4. When all the cells have detached, add serum-containing media to a final volume of 10 ml (for 6- and 12-well plates, add a 1X–2X volume of media). Note that the media must contain serum to inactivate the trypsin.
5. Pipet the cells gently up and down to mix, and then transfer the cell suspension to a centrifuge tube.
6. Spin the cells at $200 \times g$ for 5 minutes to pellet (or spin at the recommended speed and time for your cell line).
7. Aspirate the media and wash the cell pellet with 5–10 ml of 1X cold PBS.
8. Spin the cells at $200 \times g$ for 5 minutes to pellet.
9. Aspirate the PBS and resuspend the pellet in 500 μ l to 1 ml of 1X cold PBS. Mix the cell solution gently.
10. Collect a small aliquot to verify that the cells are at the desired concentration. Determine cell density electronically using a Coulter Counter or manually using a hemacytometer chamber.
11. Adjust the cell density using cold PBS so that it falls within the range of 1–10,000 cells/ μ l. Count the cells again to verify cell concentration.
12. To a 0.2-ml thin-walled PCR tube or plate well on ice, add 1 μ l of Lysis Enhancer and 10 μ l of Resuspension Buffer. **Note:** A master mix of Lysis Enhancer and Resuspension Buffer may be prepared for multiple reactions.
13. Transfer 1–2 μ l of cells (<10,000 cells) to the PCR tube/well.
Control Reaction: For the control reaction, add 1 μ l of Control HeLa Total RNA to the PCR tube or plate well instead of cell lysate.
14. Transfer the tube/plate to an incubator, water bath, or thermal cycler preheated to 75°C and incubate for 10 minutes.
Control Reaction: For the control reaction, incubate for 3 minutes.
15. After incubation, spin briefly to collect the condensate and proceed to **DNase I Digestion**, page 7.

Continued on next page

Lysing Cells, continued



Note

For adherent cells grown in tissue culture wells, note the following:

- Seed cells in tissue culture wells so that 10 μ l of resuspended cells will yield the desired concentration.
- **Master mix:** Before starting the following procedure, prepare a master mix of Lysis Enhancer and Resuspension buffer for multiple reactions. Add 1 μ l of Lysis Enhancer for every 10 μ l of Resuspension Buffer.



Important

You can order additional CellsDirect Resuspension Buffer and Lysis Enhancer from Invitrogen (Catalog no. 11739-010). Additional buffer and enhancer may be required if you are using 48-well or 24-well plates in your experiments.

Lysing Cells in Tissue Culture Wells

For adherent cells grown in tissue culture wells (i.e., in 24-well, 48-well, or 96-well plates), perform the following lysis procedure.

1. Aspirate the media in each well and wash each well with 1X cold PBS. Aspirate the PBS.
2. Add the Lysis Enhancer/Resuspension Buffer master mix (see **Note** above) to each well. For 96-well plates, add at least 11 μ l of the buffer/enhancer mix to each well. For 24-well plates, add at least 110 μ l of the buffer/enhancer mix to each well. The master mix should cover the cells in the well.
3. Incubate the plates on ice for up to 10 minutes. During that period, tap the plate periodically and check the cells under a microscope every 2–3 minutes to see whether they have detached or burst.
4. After 10 minutes, gently pipet the cells up and down to dislodge the remaining attached cells.
5. Transfer 10 μ l of the cell suspension to a 0.2-ml thin-walled PCR tube or plate well.

Control Reaction: For the control reaction, add 10 μ l of Resuspension Buffer and 1 μ l of Lysis Enhancer to a PCR tube or plate well, and then add 1 μ l of Control HeLa Total RNA.

6. Transfer the tube/plate to an incubator or thermal cycler preheated to 75°C and incubate for 10 minutes.

Control Reaction: For the control reaction, incubate for 3 minutes.

7. After incubation, spin briefly to collect the condensation, and proceed to **DNase I Digestion**, page 7.

Continued on next page

Lysing Cells, continued

DNase I Digestion

In this step, you treat the cell lysate with DNase I to degrade any contaminating DNA.

1. Place each tube/plate from Step 15, page 5, or Step 7, page 6, on ice, and add the following:

<u>Component</u>	<u>Amount</u>
DNase I, Amplification Grade (1 U/ μ l)	5 μ l
10X DNase I Buffer	1.6 μ l

2. Mix by gently pipetting up and down or briefly vortexing, and spin briefly to collect the contents.
3. Incubate the tube/plate at 25°C (or room temperature) for 5 minutes. **Note:** A longer incubation time (up to 10 minutes) may be used for larger samples (>5,000 cells). However, incubation times exceeding 10 minutes can greatly reduce cDNA yield.
4. Spin briefly, and add 4 μ l of 25-mM EDTA to each tube/well on ice. Mix by gently pipetting up and down, and spin briefly to collect the contents.
5. Incubate at 70°C for 10 minutes.
6. Spin briefly and proceed to **First-Strand cDNA Synthesis**, page 8.

First-Strand cDNA Synthesis

Required Materials

The following materials are provided by the user:

- Thermal cycler preheated to 25°C
- Ice
- Pipettes

The following materials are provided in the kit:

- 2X RT Reaction Mix
- RT Enzyme Mix (contains SuperScript® III RT, 100 units/μl; and RNaseOUT™ Recombinant Ribonuclease Inhibitor, 20 units/μl)
- RNase H (2 U/μl)

First-Strand cDNA Synthesis

1. To each tube/plate from **DNase I Digestion**, Step 6, page 7, add the following:

<u>Component</u>	<u>Amount</u>
2X RT Reaction Mix	20 μl
RT Enzyme Mix*	2 μl

*For negative RT controls, use 1 μl of sterile, distilled water and 1 μl of RNaseOUT™ Recombinant Ribonuclease Inhibitor instead of the RT Enzyme Mix.

2. Spin the tube/plate briefly to collect the contents.
3. Transfer the tube/plate to a thermal cycler preheated to 25°C and incubate for 10 minutes.
4. Incubate at 50°C for 20 minutes.
5. Inactivate the reaction at 85°C for 5 minutes.
6. Add 1 μl of RNase H (2 U/μl) to each tube/well and incubate at 37°C for 20 minutes.
7. Chill the reaction on ice, and store at –20°C or proceed directly to qPCR.

qPCR — Guidelines and Recommendations

Introduction

After first-strand cDNA synthesis, you can proceed directly to qPCR without additional purification.

Required Materials

The following materials are provided by the user:

- qPCR instrument
- Appropriate PCR plates/tubes for instrument
- Primers
- Pipettes

The following materials are provided in the kit:

- Components of the qPCR module
-

Instrument Settings

Platinum® SYBR® Green qPCR SuperMix-UDG can be used with a variety of instruments, including the ABI PRISM® 7000/7300/7500/7700/7900 and GeneAmp® 5700, Bio-Rad iCycler™, Stratagene Mx4000® and Mx3000P™, Corbett Research Rotor-Gene™, MJ Research DNA Engine Opticon® and Opticon® 2, Cepheid Smart Cycler®, and Roche LightCycler®. Optimal cycling conditions will vary; refer to your instrument manual for operating instructions.

The protocols on the following pages have been optimized for the ABI PRISM® 7700 and the Roche LightCycler®.

Primers

Primer selection is one of the most important parameters for qPCR when using a SYBR® Green detection system. To design primers, we strongly recommend using a primer design software program such as OligoPerfect™, available on the Web at www.invitrogen.com/oligos. In OligoPerfect™ designer, enter your target sequence and select *PCR: Detection* from the Application pulldown menu. Using primer design software will ensure that primers are specific for the target sequence and free of internal secondary structure, and avoid complementation at 3'-ends within each primer and with each other.

When designing primers, keep in mind that the amplicon length should be approximately 80–250 bp to optimize the efficiency of qPCR. Optimal results may require a titration of primer concentrations between 100 and 500 nM. A final concentration of 200 nM per primer is effective for most reactions.

Amplicon Size

For best results, the amplicon should be limited to 80–250 bp in size.



Note

Since PCR is a powerful technique capable of amplifying trace amounts of DNA, all appropriate precautions should be taken to avoid cross-contamination.

Continued on next page

qPCR — Guidelines and Recommendations, continued

Melting Curve Analysis

Melting curve analysis should always be performed during qPCR to identify the presence of primer dimers and analyze the specificity of the reaction. Melting curve analysis can identify primer dimers by their lower annealing temperature compared to that of the amplicon. The presence of primer dimers in samples containing template decreases PCR efficiency and obscures analysis and determination of cycle thresholds.

The formation of primer dimers most often occurs in no-template controls, where the polymerase enzyme is essentially idle, and in this case the quantitative analysis of the template samples is not affected. Melting curve analysis of no-template controls can discriminate between primer dimers and spurious amplification due to contaminating nucleic acids in reagent components.

Magnesium Concentration

Platinum® SYBR® Green qPCR SuperMix-UDG includes magnesium chloride at a final concentration of 3 mM. Optimal performance for any given target may require adjusting this level of magnesium. If necessary, use the 50-mM magnesium chloride solution included in the kit to increase the magnesium concentration.

ROX Reference Dye

ROX Reference Dye can be used to adjust for non-PCR-related fluctuations in fluorescence between reactions, and provides a stable baseline in multiplex reactions. Its use is optional. It is composed of a glycine conjugate of 5-carboxy-X-rhodamine, succinimidyl ester (25 µM) in 20 mM Tris-HCl (pH 8.4), 0.1 mM EDTA, and 0.01% Tween® 20.

ROX is supplied at 50X concentration. Add 1 µl of ROX for every 50 µl of reaction volume. To prepare a master mix of ROX and Platinum® SYBR® Green qPCR SuperMix-UDG:

1. Add ROX Reference Dye to Platinum® SYBR® Green qPCR SuperMix-UDG, at a ratio of 1 µl of ROX for every 25 µl of SuperMix-UDG.
2. Mix by vortexing for 10 seconds.
3. Store mixture at either -20°C or 4°C in the dark. Use 26 µl of ROX/SuperMix-UDG mixture per 50 µl of reaction volume.

Note: Use of ROX Reference Dye is not supported on the iCycler™, Rotor-Gene™, Opticon®, and LightCycler® platforms.. ROX Reference Dye is not required on the ABI PRISM® 7900.

Bovine Serum Albumin (BSA)

BSA (ultrapure, non-acetylated) is included as a separate tube in each kit for use in LightCycler® reactions.

qPCR — Instruments Using PCR Tubes/Plates

Introduction

This section provides a cycling program, reaction mixture, and protocol for qPCR instruments that use PCR tubes/plates (e.g., ABI PRISM®, Stratagene Mx4000® and Mx3000P™, Corbett Research Rotor-Gene™, MJ Opticon®). For a protocol using the Roche LightCycler®, see page 13.

Note: This cycling program is recommended as a starting point and guideline. Optimal cycling temperatures and times may vary for different target sequences, primer sets, and instruments.

After programming the instrument and preparing the reaction mix, follow the protocol on the following page to perform the reaction.



Note

You can use 1–8 µl of cDNA template in the following protocol, depending on the concentration of the template. Adjust the volume of water in the master mix accordingly for a final reaction volume of 50 µl. We recommend 4 µl of template as a general starting point.

Cycling Program

Program the qPCR instrument as follows:

50°C for 2 minutes hold (UDG incubation)

95°C for 2 minutes hold

50 cycles of:

95°C, 15 seconds

60°C, 30 seconds

Melting Curve Analysis: Refer to instrument documentation.

Master Mix

Use the following table to prepare a master mix of all components except template and water. **Note:** Preparation of a master mix is **crucial** in qPCR to minimize pipetting errors.

<u>Component</u>	<u>1 rxn</u>	<u>50 rxns</u>
Platinum® SYBR® Green qPCR SuperMix-UDG ¹	25 µl	1250 µl
ROX Reference Dye (optional)	1 µl	50 µl
Forward primer, 10 µM	1 µl	50 µl
Reverse primer, 10 µM	1 µl	50 µl
Autoclaved, distilled water ²	to 42–49 µl	to 2100–2550 µl

¹Final concentration: 0.06 U/µl Platinum® Taq DNA polymerase, 20-mM Tris-HCl (pH 8.4), 50-mM KCl, 3-mM MgCl₂, 200-µM dGTP, 200-µM dATP, 200-µM dCTP, 200-µM dUTP, 1 U UDG

²Volume of water used in the master mix will depend on template volume used in the reaction (see steps 3–4 next page).

Continued on next page

qPCR — Instruments Using PCR Tubes/Plates, continued

Protocol

1. Program the qPCR instrument to perform a brief UDG incubation immediately followed by PCR amplification, as shown on the previous page. Optimal cycling temperatures and times may vary for different target sequences, primer sets, and instruments.
 2. Prepare a master mix of all components except template as specified on the previous page.
 3. For each reaction, add 42–49 μl of the master mix (depending on template volume) to a 0.2-ml microcentrifuge tube or each well of a 96-well PCR plate.
 4. Add 1–8 μl of cDNA template from the first-strand synthesis reaction (Step 7, page 8) to each reaction vessel, for a final reaction volume of 50 μl . (Use 4 μl of template as a general starting point.) Cap or seal the tube/plate.
 5. Gently mix and make sure that all components are at the bottom of the reaction vessel. Centrifuge briefly if needed.
 6. Place reactions in a thermal cycler programmed as described above. After cycling, hold the reaction at 4°C until further analysis. Collect and analyze the results.
-

qPCR — Roche LightCycler®

Introduction

This section provides a cycling program, reaction mixture, and protocol for the Roche LightCycler®.

After programming the instrument and preparing the reaction mix, follow the protocol on the following page to perform the reaction.

Cycling Program

Program the LightCycler® as follows:

Program choice: Amplification

Analysis mode: Quantification

50°C for 2 minutes hold (UDG incubation)

92°C for 1 minute hold

50 cycles of:

92°C, 5 seconds

60°C, 30 seconds (single acquire)

Melting Curve Analysis

Program choice: Melting curve

Analysis mode: Melting curves

95°C, 0 seconds (20°C/second transition)

55°C, 15 seconds (20°C/second transition)

92°C, 0 seconds (0.1°C/second transition — continuous acquisition)

40°C, 0 seconds

Master Mix

Use the following table to prepare a master mix of all components except template. **Note:** Preparation of a master mix is **crucial** in qPCR to minimize pipetting errors.

<u>Component</u>	<u>1 rxn</u>	<u>34 rxns</u>
Platinum® SYBR® Green qPCR		
SuperMix-UDG*	10 µl	340 µl
BSA, UltraPure (1 mg/ml)	1 µl	34 µl
Forward primer, 10 µM	1 µl	34 µl
Reverse primer, 10 µM	1 µl	34 µl
Autoclaved, distilled water	to 18 µl	to 612 µl

*Final concentration: 0.06 U/µl Platinum® Taq DNA polymerase, 20-mM Tris-HCl (pH 8.4), 50-mM KCl, 3-mM MgCl₂, 200-µM dGTP, 200-µM dATP, 200-µM dCTP, 200-µM dUTP, 1 U UDG

Continued on next page

Protocol

1. Program the LightCycler® to perform a brief UDG incubation immediately followed by PCR amplification, as shown on the previous page. Optimal cycling temperatures and times may vary for different target sequences and primer sets.
 2. Set the fluorescence on the LightCycler® to the F1 channel.
 3. Prepare a master mix of all components except template as specified on the previous page.
 4. For each reaction, add 18 µl of the master mix to each capillary tube.
 5. Add 2 µl of the cDNA from the first-strand synthesis reaction (Step 7, page 8) to each capillary tube for a final reaction volume of 20 µl, and cap the tube.
 6. Centrifuge tubes at $700 \times g$ for 5 seconds.
 7. Place reaction tubes in the rotor of the LightCycler® and run the program. After cycling, hold the reaction at 4°C until further analysis. Collect and analyze the results.
-

Troubleshooting

Problem	Possible Cause	Suggested Solution
Cells in tissue-culture wells do not detach/burst	Incubation temperature of lysis reaction is too low	Incubate lysis reaction at room temperature instead of on ice.
No amplification curve appears on the qPCR graph	There is no PCR product	Run the PCR product on a gel to determine whether PCR worked. Then proceed to the troubleshooting steps below.
No PCR product is evident, either in the qPCR graph or on a gel	Procedural error	Confirm that all steps were followed. Use the Control RNA to verify the efficiency of the first-strand reaction (see the next page on troubleshooting with the Control RNA).
	RNA is degraded	Add control total HeLa RNA to sample to determine if RNase is present in the first-strand reaction. A longer DNase I digestion can hydrolyze the RNA in the sample. Use a digestion time of <10 minutes. Maintain aseptic conditions to prevent RNase contamination.
	Primer design is suboptimal	Verify your primer selection. We recommend using validated pre-designed primers or designing primers using dedicated software programs or primer databases.
	Target mRNA contains strong transcriptional pauses	Maintain an elevated temperature after the annealing step. Increase the temperature of first-strand reaction (up to 55°C).
PCR product is evident in the gel, but not on the qPCR graph	qPCR instrument settings are incorrect	Confirm that you are using the correct instrument settings (dye selection, reference dye, filters, acquisition points, etc.).
	Problems with your specific qPCR instrument	See your instrument manual for tips and troubleshooting.
Poor sensitivity	Not enough starting template RNA	Increase the number of cells used
Higher than expected signal	Too much first-strand product was used in qPCR	Decrease amount of the first-strand product in qPCR.

Continued on next page

Troubleshooting, continued

Problem	Possible Cause	Suggested Solution
Signals are present in no-template controls, and/or multiple peaks are present in the melting curve graph	Template or reagents are contaminated by nucleic acids (DNA, cDNA)	<p>Use melting curve analysis and/or run the PCR products on a 4% agarose gel after the reaction to identify contaminants.</p> <p>To reduce the risk of contamination, take standard precautions when preparing your PCR reactions. Ideally, amplification reactions should be assembled in a DNA-free environment. We recommend using aerosol-resistant barrier tips.</p>
	Primer dimers or other primer artifacts are present	<p>Use melting curve analysis to identify primer dimers. We recommend using validated pre-designed primer sets or designing primers using dedicated software programs or primer databases.</p> <p>Primer contamination or truncated or degraded primers can also lead to artifacts. Check the purity of your primers by gel electrophoresis. If agarose gels are used, we recommend cooling the gels before visualization with intercalating dyes.</p>
Product detected at higher than expected cycle number	RNA is degraded	<p>Add control total HeLa RNA to sample to determine if RNase is present in the first-strand reaction.</p> <p>A longer DNase I digestion can hydrolyze the RNA in the sample. Use a digestion time of <10 minutes.</p> <p>Maintain aseptic conditions to prevent RNase contamination.</p>
Product detected at lower-than-expected cycle number, and/or positive signal from no-template controls	Template or PCR carry-over contamination	Isolate source of contamination and replace reagent(s). Use separate dedicated pipettors for reaction assembly and post-PCR analysis. Assemble reactions (except for target addition) in a DNA-free area. Use aerosol-resistant pipet tips or positive displacement pipettors.
Unexpected bands after gel analysis	Contamination by genomic DNA	<p>Do not omit the DNase Digestion step on page 7. For larger samples (>1,000 cells), use a longer DNase I incubation time, i.e., up to 10 minutes.</p> <p>Design primers that anneal to sequence in exons on both sides of an intron or exon/exon boundary of the mRNA to allow differentiation between amplification of cDNA and products potential contaminating genomic DNA.</p> <p>To test if products were derived from DNA, prepare a negative RT control.</p>
	Nonspecific annealing of qPCR primers	<p>Vary the annealing conditions.</p> <p>Optimize magnesium concentration for each template and primer combination.</p>

References

- Berger, S.L. and Kimmel, A.R. (1987) *Methods Enzymol* 152, 316.
- Chomczynski, P. and Sacchi, N. (1987) *Anal. Biochem.* 162, 156.
- Chou, Q., Russel, M., Birch, D., Raymond, J., and Bloch, W. (1992) Prevention of pre-PCR mis-priming and primer dimerization improves low-copy-number amplifications. *Nucl. Acids Res.* 20, 1717.
- Compton, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 39, Academic Press, Inc.
- D'Alessio, J. M., Gruber, C. E., Cain, C., and Noon, M. C. (1990) *Focus*[®] 12, 47.
- Frohman, M.A., Dush, M.K, and Martin, G.R. (1988) *Proc. Nat. Acad. Sci USA* 85, 8998.
- Gerard, G.F. (1994) *Focus*[®] 16, 102.
- Gerard, G.F., Schmidt, B.J., Kotewicz, M.L., and Campbell, J.H. (1992) *Focus*[®] 14, 91.
- Higuchi R., Fockler C., Walsh P.S., Griffith R. (1992) Simultaneous amplification and detection of specific DNA sequences. *Bio/Technology* 10, 413–417.
- Ishiguro, T., Saitoh, J., Yawata, H., Yamagishi, H., Iwasaki, S., and Mitoma, Y. (1995) Homogeneous quantitative assay of hepatitis C virus RNA by polymerase chain reaction in the presence of a fluorescent intercalater. *Anal. Biochem.* 229, 207.
- Lee, C.C. and Caskey, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 46, Academic Press, Inc.
- Lindahl, T., Ljungquist, S., Siebert, W., Nyberg, B., and Sperens, B. (1977) DNA N-glycosidases: properties of uracil-DNA glycosidase from *Escherichia coli*. *J. Biol. Chem.* 252, 3286.
- Longo, M., Berninger, M., and Hartley, J. (1990) Use of uracil DNA glycosylase to control carry-over contamination in polymerase chain reactions. *Gene* 93, 125.
- Meateo, R., Sico, E., Wang, L.F., Fung, J. (1994) Uracil DNA Glycosylase: Minimizing Residual Enzyme Activity After Nested PCR. *Focus* 16.4, 104-5.
- Sambrook J., Fritsch, E.F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- Sharkey, D.J., Scalice, E.R., Christy, K.G., Atwood, S.M., and Daiss, J.L. (1994) Antibodies as thermolabile switches: high temperature triggering for the polymerase chain reaction. *BioTechnology* 12, 506-509.
- Simms, D., Cizdziel, P.E., and Chomczynski, P. (1993) *Focus*[®] 15, 99.
- Wittwer C.T., Herrmann M.G., Moss A.A., and Rasmussen R.P. (1997) Continuous fluorescence monitoring of rapid cycle DNA amplification. *BioTechniques* 22, 130-138.
- Takagi, M., Nishioka, M., Kakiyama, H., Kitabayashi, M., Inoue, H., Kawakami, B., Oka, M., and Imanaka, T. (1997) *Appl. Environ. Microbiol.* 63, 4504.
- Schwabe, W., Lee, J.E., Nathan, M., Xu, R.H., Sitaraman, K., Smith, M., Potter, R.J., Rosenthal, K., Rashtchian, A., Gerard, G.F. (1998) *Focus*[®] 20, 30.
-

Related Products

Product	Size	Cat. No.
CellsDirect Resuspension and Lysis Buffer	10 ml Resuspension Buffer and 1 ml Lysis Enhancer	11739-010
E-Gel® Pre-cast Agarose Gels		
0.8% Starter Pak	9 gels and base	G5000-08
1.2% Starter Pak	9 gels and base	G5000-01
2% Starter Pak	9 gels and base	G5000-02
4% Starter Pak	9 gels and base	G5000-04
UltraPure™ Agarose	100 g	15510-019
	500 g	15510-027
UltraPure™ Agarose 1000	100 g	10975-035
100-bp DNA Ladder	50 µg	15628-019
123-bp DNA Ladder	100 µg	15613-011
	250 µg	15613-029
1-Kb Plus DNA Ladder	250 µg	10787-018
	1,000 µg	10787-026

Purchaser Notification

Limited Use Label License No. 1: Thermostable Polymerases

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, and 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008; Phone (760) 603-7200 or e-mail: outlicensing@lifetech.com.

Continued on next page

Purchaser Notification, continued

Limited Use Label License No. 14: Direct Inhibition by Anti-Polymerase Antibodies

Licensed to Life Technologies Corporation, under U.S. Patent Nos. 5,338,671; 5,587,287; and foreign equivalents for use in research only.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity) in a manner consistent with the accompanying product literature. The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents that are owned or controlled by Life Technologies Corporation and/or Molecular Probes, Inc. which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Continued on next page

Purchaser Notification, continued

**Limited Use
Label License
No. 276: Dye-
Intercalation
Detection Assays**

The purchase price of this product includes a limited, non-transferable license under U.S. Patent Nos. 5,994,056 and 6,171,785 and under any corresponding foreign patent claims owned by Roche Molecular Systems, Inc. or F. Hoffmann-La Roche Ltd, to use only this amount of the product for dye-intercalation detection assays and related processes described in said patent solely for the purchaser's own internal research and development activities. No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of this product. Further information on purchasing licenses for these patented processes may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA. No rights are conveyed expressly, by implication or by estoppel under any apparatus or system claim of any patent or patent application.

Trademarks

ABI PRISM® and GeneAmp® are registered trademarks of Applera Corporation.

iCycler™, Mx4000™, Rotor-Gene™, DNA Engine Opticon™, and Smart Cycler® are trademarks of their respective companies.

LightCycler® is a registered trademark of Idaho Technologies, Inc.

Technical Service

Web Resources



Visit the Invitrogen web site at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
 - Complete technical support contact information
 - Access to the Invitrogen Online Catalog
 - Additional product information and special offers
-

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:

5791 Van Allen Way
Carlsbad, CA 92008 USA
Tel: 1 760 603 7200
Tel (Toll Free): 1 800 955 6288
Fax: 1 760 602 6500
E-mail:
tech_support@invitrogen.com

Japanese Headquarters:

LOOP-X Bldg. 6F
3-9-15, Kaigan
Minato-ku, Tokyo 108-0022
Tel: 81 3 5730 6509
Fax: 81 3 5730 6519
E-mail:
jpinfo@invitrogen.com

European Headquarters:

Inchinnan Business Park
3 Fountain Drive
Paisley PA4 9RF, UK
Tel: +44 (0) 141 814 6100
Tech Fax: +44 (0) 141 814 6117
E-mail:
eurotech@invitrogen.com

SDS

SDSs (Safety Data Sheets) are available on our web site at www.invitrogen.com/sds.

Quality Control

The Certificate of Analysis (CofA) provides detailed quality control information for each product. The CofA is available on our website at www.invitrogen.com/cofa, and is searchable by product lot number, which is printed on each box.

Continued on next page

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, report it to our Technical Support Representatives.

Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

©2010 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Notes



Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information visit our web site at www.invitrogen.com