

## PURIFLASH® 800 - User Manual

ITM - 2012/10/31

### Process Chromatography

#### Large Scale Purification

✓ Flash columns up to F1600 format & more



## PURIFLASH® 800

800 ml/min @ up to 10 bar



INTERFINE CHEMICALS



**interchim®**

211 bis, avenue JF Kennedy  
BP 1140 - 03100 Montluçon  
Fax +33 4 70 03 82 60

Hotline +33 4 70 03 73 01 • [interfine@interchim.fr](mailto:interfine@interchim.fr)

# Summary

---

<b>1 &gt; General description</b>	<b>page A. 3</b>
<b>2 &gt; Packing list</b>	<b>page A. 7</b>
<b>3 &gt; Installation</b>	<b>page A. 8</b>
3.1 Assembling the system	
3.2 Install the tubing and prepare the containers	
3.3 Set up the column, the pre-column and the dry-load	
3.4 Prepare the fractions collection	
3.5 Position the stylus	
3.6 Make the connections (electrical and data transmission)	
3.7 Options and accessories	
<b>4 &gt; Operation</b>	<b>page A. 17</b>
4.1 Control the system	
4.2 Make a liquid injection	
4.3 Make a solid injection	
<b>5 &gt; Safety regulations</b>	<b>page A. 20</b>
<b>6 &gt; Maintenance</b>	<b>page A. 21</b>
6.1 Introduction	
6.2 Removing a pump head	
6.3 Replacing a seal	
6.4 Cleaning / replacing a check valve	
6.5 Replacing the detector deuterium lamp	
6.6 Cleaning the detector flow cell	
<b>7 &gt; Preventive measures - Upkeep schedule</b>	<b>page A. 26</b>
<b>8 &gt; Warranty</b>	<b>page A. 27</b>
8.1 Generalities	
8.2 Limited warranty of consumable parts	
<b>9 &gt; List of spare parts + Blow-up views</b>	<b>page A. 30</b>
<b>10 &gt; Specifications</b>	<b>page A. 36</b>
<b>► Appendix 1 : Rack sets available</b>	<b>page A. 37</b>
<b>► Appendix 2 : Service report</b>	<b>page A. 38</b>

## 1 > General description

The **PURIFLASH® 800** system is a modern instrument designed for Flash chromatography which allows to withdraw mobile phases from four containers of solvents, form a gradient, inject liquid or solid samples and provide the solution to a chromatography column.

A detector measures the absorbance and gives a signal called chromatogram which is visualized with the built-in computer and the control software "INTERCHIM Software" installed. Then fractions can be collected with the fraction collector.



The characteristics of a **PURIFLASH® 800** are:

- ✓ **Flow rate:** Up to 800mL/min
- ✓ **Pressure:** Max. 10 bar (145psi)
- ✓ **Gradient:** Binary or Quaternary
- ✓ **Column size:**  
From F4 to F1600 format - according to column trademark
- ✓ **UV Detection:**  
Spectrophotometer dual-wavelength  
(DAD technology, range: 200 - 600nm)
- ✓ **Collection:** Several types of racks available.

### The complete system consists of following modules:


#### The 2-head pump:

As mentioned above, the pump can deliver mobile phase at flow rates up to eight hundred millilitres per minute with low residual pulsation. It is a high performance double-piston pump (reciprocal movement). Thanks to a system "pulleys + belt", a stepper motor drives two special cams on which lean the two pistons. Their motion generates a cycle and then creates the suction.

A static mixing chamber at the outlet ensures a suitable mixing of solvents whatever the working flow rate.

Moreover, an optic sensor can detect leakages in the pump drain pan. The pumping is then interrupted and the user informed via the control software.

The 2 or 4-ways suction valve allows to perform binary or quaternary elution gradients.



## The injection unit and columns holder system:


This device gives the possibility to inject sample manually (liquid or solid via dry-load) directly to the chromatography column installed.

The adjustable height holder of the black POM holder allows to adapt all types of columns. A knurled knob can block this holder on the vertical guidance rail. The bottom union is fixed.



## The UV diode array detector:


This instrument measures the amount of light absorbed by a liquid passing through the cell at a given wavelength. It is equipped with a deuterium lamp. The two wavelengths and the SCAN function available are selected within a large range, from 200 to 600nm or from 190 to 840nm (in option), P/N: JO5140.






**The fraction collector:**

The fraction collector is placed downstream of the detector. After the detection of the signal (chromatogram), the results obtained at the outlet of the column (peaks) can be directed to tubes disposed on racks. A built-in 3-way valve allows to choose between the collection or the rejection of the liquid.

**The computer and its touch screen:**

The system is controlled by the central processing unit via a touch screen, which replace the keyboard and the mouse. The operator uses the control software and can create methods to run the instrument.




This module contains the main ON/OFF switch and centralizes the connection for data and power supply to the other modules.

Four USB connectors are also available, 1 in front of the system under the screen and 3 to the left hand side (cable maximum length = 2m).

### The external connections:

Some others external connectors allow to plug in additional devices like second detector or options. The maximum length of all external cables connected to the system is 2m.





## 2 > Packing list

### Standard system *PURIFLASH*<sup>®</sup> 215

Item:

✓ *PURIFLASH*<sup>®</sup> 800 system composed of:

- Pump module,
- Column holder module,
- Pre-column holder,
- Detector module,
- Fraction collector module,
- Computer & power supply module.

✓ **Tubing with fittings (ETFE):**

- Suction (OD = 3/16" - 1.50m - ID = 3.2mm - x4 or x2 depends to the gradient valve)
- Gradient valve outlet to Tee Pump inlet (OD = 3/16" - 28cm - ID = 3.2mm) (present on the system)
- Pump outlet to Pre-column holder inlet (OD = 1/8" - 60cm - ID = 2.4mm)
- Pre-column holder outlet to column holder inlet (OD = 1/8" - 90cm - ID = 2.4mm)
- Column holder outlet to detector inlet (OD = 1/8" - 70cm - ID = 2.4mm)
- Detector outlet to fraction collector inlet (OD = 1/8" - 30cm - ID = 2.4mm)
- Waste outlet (OD = 1/8" - 1.50m - ID = 2.4mm)
- Pressure test OQ (OD = 1/16" - 50cm - 1/16" - ID = 0.5mm)

✓ **Rack set "DT8250" for 29.5mm diameter tubes (4 racks)**

✓ **Solvent Tray with a drainage system**

✓ **Serial communication cables (2 x Sub-D9 male/female - 1 x Sub-D9 female/female)**

✓ **Distribution power cables (2 x main voltage supply - 1 x 24V)**

✓ **Binding column**

✓ **General power cable (model according to the destination country)**

✓ **USB key (including the file of installation ".exe" of the freeware Interchim<sup>®</sup> Software Demo)**

✓ **Stylus**

✓ **Keyboard**


✓ **Documentation:**

- User Manual for *PURIFLASH*<sup>®</sup> 800 system
- Instruction Manual for "Interchim<sup>®</sup> Software"
- "Adjusting Protocol of Assembled Unit" for "Flash 06 DAD" UV detector
- PDF document: Operational Qualification and test report ("My documents")
- Quick start for "Interchim<sup>®</sup> Software"

## 3 > Installation

### 3.1 Assembling the system

Select a suitable location by ensuring adequate air-space (~15cm around the system) to provide ventilation of the controller, the pump, the detector cooling fans. Its total weight is approximately 60kg.



✓ **Assemble the system in the following order:**

- Fraction collector : ~ 25kg
- Unit control : ~ 15kg
- Detector : ~ 7kg
- Pump : ~ 8kg
- Column Holder: ~ 3kg
- Solvent tray : ~ 1kg

✓ **Put then the stainless steel solvent tray supplied on top to dispose solvents (up to five).** They cover the totality of the system with a drainage system and their design avoids any leakage risk.





### 3.2 Install the tubing and prepare the containers

The tubes situated outside the housing are delivered apart in order to avoid all deterioration during transfer, they are already equipped with the suitable fittings (nut + ferrule) to connect them easily to the system. Regarding the main circuit, two sizes of tube are used:

3/16" - ID = 3.2mm for suction

1/8" - ID = 2.4mm for others.


Follow the instructions below to correctly install the tubing:

#### Suction:


Provide yourself with the two or four "Suction tubes" (l = 1.5m, 3/16" - ID = 3.2mm) supplied, they are more brittle than others so take care not to bend them.

For the tube "A", take the extremity where there is a fitting with a long spring and put the ferrule in the PEEK body of the suction valve. Screw the PEEK nut, tighten it hand-tight and position properly the long spring.

Renew the same operation for "B", "C" and "D" as indicated on the photo.



To ensure the pumping continuity of the solvents and avoid the air entry in the hydraulic circuit, the others extremities can be equipped with ballasts (in option, P/N: FV1290). It consists in fact to put a fitting (nut + set ferrule) and screw it in a stainless steel piece which serve of weight to maintain the tubes at the bottom of the containers.



Position the solvents containers on the solvent trays for instance, and put the ballasted tubes into them.

#### GRADIENT VALVE:

The gradient valve to Tee Pump inlet tube is already connected into the pump (l = 28cm - 3/16" ID = 3.2mm).

### Pump to pre-column holder:

Using the 60cm x 2.4mm ID tube insert one end into the mixing chamber and the other end into the pre-column holder inlet. Position the tubing, seat the ferrules and hand tight each fitting.

### Pre-column holder to column holder:

Using the 90cm x 2.4mm ID tube insert one end into the pre-column holder outlet and the other end into the column holder inlet. Position the tubing, seat the ferrules and hand tight each fitting.

### Column holder to detector:

Using the 70cm x 2.4mm ID tube insert one end into the column holder outlet of the stainless steel bulkhead union 90° and the other end into the detector inlet. Position the tubing, seat the ferrules and hand tight each fitting.

### Detector to fraction collector:

Using the 30cm x 2.4mm ID tube insert one end into the detector outlet and the other end into the fraction collector inlet. Seat the ferrules and hand tighten.

### Waste:

Using the 1.5m x 2.4mm tube attach the end with the fitting into the fraction collector waste port and hand tighten. On the other end, insert into an appropriate waste bottle.

## 3.3 Set up the column, the pre-column and the dry-load

A large range of columns can be used on the column holder, columns packed with normal or reverse phase, with an external diameter up to 120mm.

The "Extension kit for Large columns (P/N: DU9671)" is supplied with this device, this overall dimension can go up to 120mm.

On the pre-column holder, dry-load or pre-column (4g up to 300g) can be used.

### Set up the column:



- 1 -

Unscrew lightly the knurled knob in order to allow its vertical movement.



- 3 -

Move up the large black POM holder on the vertical guidance rail. Put the PTFE outlet adapter (P/N: DZ5200) at the bottom of the large column (800g and 1600g). Place the bottom fitting of the column on the "Adapter ETFE or Stainless steel Luer lock female to 1/4-28 male" located on the stainless steel bulkhead union 90°.



- 2 -

Unlock the black ring of the fitting "Adapter PEEK Luer lock male to 1/4-28 male" located on the column holder.



- 4 -

Move down the large black POM holder. Put the stainless steel inlet adapter (P/N: DZ5200) at the top of the large column (800g and 1600g). Lock again the black ring of the fitting "Adapter PEEK Luer lock male to 1/4-28 male" to the top fitting of the adapter.

**Set up the dry-load and pre-column:**

To inject solid sample, make dry loading injection, or use the pre-column, you have the possibility to use a column of which the maximum dimensions are 27cm in height and 70mm for external diameter. Install it as follows:



- 1 -

Unscrew lightly the knurled knob of the holder in order to allow its vertical movement.



- 3 -

Move up the large black POM holder on the vertical guidance rail and place the bottom fitting of the dry-load or pre-column on the "Adapter ETFE or Stainless steel Luer lock female to 1/4-28 male" located on the stainless steel bulkhead union 90°.



- 2 -

Unlock the black ring of the fitting "Adapter PEEK Luer lock male to 1/4-28 male" of the pre-column holder.



- 4 -

Move down the large black POM holder and lock again the blazzack ring of the fitting "Adapter PEEK Luer lock male to 1/4-28 male" to the top fitting of the dry-load


**3.4 Prepare the fractions collection**

According to the type of collection required (small or big volumes, short or long fractionation, etc...), the equipments to install are different.


The standard rack set (29.5 x 200mm - P/N: PF4440) is supplied by default.



**CAUTION: Before any action, take a diagonal plier and carefully cut the cable tie which holds the arm during the transfer.**



Place first the tubes in the holes and then put the racks on the tray. Thanks to notches machined on racks and positioning guides already mounted on the system, their emplacement is easy and necessarily proper.




You can also equip your collector with an optional Enhancement tube device for 13 x 100m (P/N : DW3360). It is designed for elevating the racks with small height tubes. Its installation is simple too because you have just to emplace it on the collector tray using positioning guides.

**Several models of rack sets are available, find a detailed list on “APPENDIX 1”.**

It is important to maintain appropriate spacing of the dispensing needle above the tubes. A space of 5 - 25mm between the end of the needle to the top of the tube is recommended.

To adjust the needle height turn off the system, remove the racks and gently move the arm into an easy to access position. Do not use the needle to relocate the arm as bending and mis-alignment may occur. Loosen the set screw and adjust the height of the needle and then tighten the set screw. When the system is turned on the arm will return to its “home” position.



### 3.5 Position the stylus

The stylus is supplied to control the system via computer touch screen. A fixed holder allows to easily put it away when it is not used and its extensible cord (linked to the holder) gives the suitable length.



Stick the stylus holder just as you like in order to profit at best of its features (i.e. to the left side of the system as showed on this picture).


### 3.6 Make the connections (electrical and data transmission)

The complete system is supplied with appropriate fuses and is ready to operate at the line voltage of the shipping destination.

The Computer module serves as the primary point for power distribution and data collection. The ports on the modules are labelled to ensure proper connections. Using the cables provided connect each module to the unit control.


**Unit control:**

Back panel of the unit control



**Detector:**

Back panel of the detector






Take first the special cable for 24V supply and screw the round male connector in the suitable socket (label "DET.") of the unit control. Screw the other end (round female connector) in the back panel of the detector.

Link the unit control (label "COM2: DETECTOR") and the detector module using the Sub-D9 female/female cable supplied. Do not forget to set the connectors by fastening the two fixing screws.

## Pump:



Take the power cable for main voltage supply and connect it at the same time in the plugs of the unit control (label "PUMP") and in the sockets of the pump module.

Link the unit control (label "COM1: PUMP") and the pump module using the Sub-D9 male/female cable supplied. Do not forget to set the connectors by fastening the two fixing screws.

## Collector:




Take the power cable for main voltage supply and connect it at the same time in the plugs of the unit control (label "COLL.") and in the sockets of the pump module.

Link the unit control (label "COM3: COLLECTOR") and the collector module using the Sub-D9 male/female cable supplied. Do not forget to set the connectors by fastening the two fixing screws.

Plug the AC power cable supplied in the general switch located to the left side of the system, and then to the mains, knowing that the connector depends on the destination country (comply with country regulations in effect).

For American system: 110-120VAC ; 50-60Hz ; 6.2A (delay action fuses 2 x 6.3A)

For European system: 220-240VAC ; 50-60 Hz ; 3.1A (delay action fuses 2 x 3.15A)



Check that all the ON/OFF switches of secondary modules (pump, detector, collector) are on “1” and switch on the unit with the main ON/OFF switch.

### 3.7 Options and accessories

1. “Solvent level sensors 5-1 (global)” (P/N: FSP720) or “Solvents level sensors 5-5 (Separated)” (P/N: FSP730)
2. “Bar code reader” (P/N: GO2380)
3. “Enhancement tube device for 13 x 100mm tubes” (P/N: DW3360)
4. “Video camera + TLC plates holder” (P/N: GO2390)
5. “Fume enclosure” (P/N: PF4350)
6. “Quad system” (P/N: PF4360)
7. “2<sup>nd</sup> acquisition channel” (P/N: PF4320)
8. “FLASH- ELSD” (P/N: FI0320)
9. “Funnel rack set 1 unit” (P/N : FJ4331) or “Funnel rack set 4 units” (P/N: FJ4330)
10. “Safety solvent caps kit - tubing 3/16” ” (1 unit) (P/N: JO1620)
11. “Safety waste cap (with container 20L + filter)” (P/N: JO4500)
12. “Drainage system of the tray collector” (P/N: DZ7790)
13. “Injection pump + electric switching valve (200mL/mm; 30 bar)” (P/N: FY9960)
14. “CarouXel” (P/N: JO3750)
15. “Autosampler” (P/N: LO8850)
16. “Mass Spectrometer + interface” (P/N: JO1090)
17. “2<sup>nd</sup> fraction collector coupled” (P/N: FK6140)
18. “Ballasting kit for 3/16” and 1/8” tubing” (P/N: FV1290)
19. “Racks set”
20. “Stand Alone column holder for flash” (P/N: PF4530)

## Designation / Items:

### 1. Solvent level sensors 5-1 (global)

- Communication interface 5-1 + cable
- Sensors + cables (x 5)
- Straps (80cm) (x 5)
- Removable cable ties (x 5)
- Power adapter
- Hardware
- Installation procedure

### 1. Solvent level sensors 5-5 (separated)

- Communication interface 5-5 + cable
- Sensors + cables (x 5)
- Straps (80cm) (x 5)
- Removable cable ties (x 5)
- Power adapter
- Hardware
- Installation procedure

### 2. Bar code reader

- Reader with USB connection
- Bracke
- Hardware
- Installation procedure

### 3. Enhancement tube device for 13 x 100mm tubes

### 4. Video camera + TLC plates holder

- Webcam (with USB connection)
- Bracket
- CD for Webcam drivers installation
- TLC plates holder
- Hardware
- Installation procedure

### 5. Fume enclosure

- Side panel
- Transparent door
- 5 screws
- Collar
- Magnetic flap
- Hardware
- Installation procedure

### 6. Quad system

- Module for 4 additional columns installation
- Inlet tubing with stainless steel fitting set (1/8" - ID = 1.6mm - 1.50m)
- Outlet tubing (1/8" - ID = 1.6mm - 1.50m)
- PEEK 1/4-28 long nut for 1/8" tubing (x 3)
- ETFE ferrule for 1/8" tubing (x 3)
- CD for TRP converter drivers installation
- External control cable
- USB cable
- Power cable (model according to the destination country)
- Hardware
- Installation procedure

### 7. 2<sup>nd</sup> Acquisition channel

- Module USB-1408FS + USB cable
- USB key for drivers installation
- Installation procedure

### 8. FLASH-ELSD

- FLASH-ELSD
- 1V Output cable
- Power cable
- Splitter box with micro-needle (option)
- ETFE Tubing (1/16", ID = 0.25mm) ; Splitter Outlet - ELSD Inlet : length = 30cm, fitting 10 - 32 Male (P/N : 780771) (connected to Splitter box)
- ETFE Tubing (1/8", ID = 1.6mm) ; Splitter Outlet - UV Inlet: length = 72cm + Adapter 1/4-28F à 10 - 32 Male (P/N: 737664) +spring 20 cm (P/N : PFS800) (connected to Splitter box)
- ETFE Tubing (1/8", ID = 1.6mm) ; Column Outlet - Splitter Inlet: length = 72cm + Adapter 1/4 - 28F to 10 - 32 Male (P/N: 737664) (connected to Splitter box)
- Manometer ( P/N: FJ6720) (option)
- 2 red or white tubing to connect the gas (6mm OD, 2m + 1 m and fitting)
- 1 PEEK fitting 10 - 32 Male, fitting Inlet ELSD (P/N: 780771)
- Adapter 1/4" - 28 Female - 1/4" - 28 Female Nylon (P/N:187210) (only for the *PURIFLASH® 800*)
- Installation procedure

### 9. Funnel rack set 1 unit or funnel rack set 4 units

### 10. Safety solvent caps kit - tubing 3/16"

### 11. Safety waste cap

- Container 20L
- Funnel
- Filter

### 12. Drainage System of the tray collector

### 13. Injection pump + electric switching valve (200mL/mm; 30 bar)

### 14. CarouXel

- Power cable
- RS 232 cable
- Installation procedure

### 15. Autosampler

- Power cable
- RS 232 cable
- User's manuel
- Installation procedure

### 16. Mass spectrometer + interface

### 17. 2<sup>nd</sup> fraction collector coupled

- Fraction collector
- General power cable
- Serial communication cable (Sub-09 male/female)
- Tubing with 2 fittings and 2 ferrules (l = 3m - 1/8" - ID = 1.6mm)

### 18. Ballasting kit for 3/16" and 1/8" tubing

### 19. Racks set

### 20. Stand Alone column holder for flash

## 4 > Operation

### 4.1 Control the system

When all the installations and connections are made, switch on the *PURIFLASH® 800*. After few seconds and the initialisations of different instruments, the control software “Interchim® Software” automatically starts. Note the good starting (no message “Communication Error”), is necessary for the correct working of the system.



Once the user interface opened, several graphic elements allow to see in real time the state of the instruments and the evolution of the working parameters. The application has been developed with the aim of being convivial and intuitive. Refer to the “Instruction Manual” supplied for more explanations on the software use.

Use the stylus on the tactile screen to create methods in the software and control the system.


## 4.2 Make a liquid injection



The liquid injection is done directly through the chromatography column by using a syringe.



In the panel "Settings" of the method edited, select "Equilibration" in the section "Injection Mode". According to the type of column chosen, a time and a flow rate are written by default, they depend on its parameters but you can modify them if needed.



After having launched the method (= "RUN"), the equilibration phase starts. That allows to condition the column with the solvents used. At the end of the time (or at a click on "START"), the pump stops and the following dialog box appears on the screen.

It is the moment to inject the sample :



- 1 -

Unlock the black ring of the fitting "Adapter PEEK Luer lock male to 1/4 - 28 Male" of the column holder.



- 2 -

Place the syringe filled with sample on the column and inject manually the desired quantity.



- 3 -

Remove the syringe.



- 4 -

Screw the fitting "Adapter PEEK Luer lock male to 1/4 - 28 Male" adapter. Answer "yes" to the question of the software. The method starts.




### 4.3 Make a solid injection

The solid injection is done via a dry-load in which the solid sample has been previously deposited. The dry-load is emplaced directly on the pre-column holder.

An optional separated holder is available to use greater pre-columns and have a support comparable to the one of the main column (P/N: PF4530).



In the panel "Settings" of the method edited, select "Equilibration" in the section "Injection Mode". According to the type of column chosen, a time and a flow rate are written by default, they depend on its parameters but you can modified them if need be.



After having launched the method (= "RUN"), the equilibration phase starts. That allows to condition the column with the solvents used. At the end of the time (or at a click on "START"), the pump stops and the following dialog box appears on the screen.

It is the moment to put the dry-load:



- 1 -

Unlock the black ring of the fitting "Adapter PEEK Luer lock male to 1/4 - 28 Male" of the pre-column holder.



- 2 -

Move up the large black POM holder on the vertical guidance rail and place the bottom fitting of the dry-load or pre-column on the "stainless steel luer lock female to 1/4"-28 male" on the stainless steel bulkhead union 90°.



- 3 -

Move down the large black POM holder and lock again the black ring of the fitting "adapter PEEK luer lock male to 1/4"-28 male" to the top fitting of the dry-load.

- 4 -

Answer "Yes" to the question of the software. The method starts.

## 5 > Safety regulations

This equipment is built according to the E.U. security standards. However, risks and dangers could remain if the system is used in a different way from which it is intended or if it is used by untrained personnel.

The **PURIFLASH® 800** should only be used by laboratory trained or experienced persons. The user should inform the distributor in case of problems and safety-related matters that occur during the use of the instrument.



*Alerts you to potential hazardous situations that could result in serious injury or damage of the system.*

This system is designed for chromatographic purposes. It must be operated using appropriate solvents and within specified ranges for pressure, flows and temperatures as described in this manual. If it is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

When handling potentially hazardous solvents and samples, it is important to follow standard laboratory safety procedures. Common chemicals used in the laboratory can be toxic, corrosive, flammable or a combination of these properties.

**A list of all safety aspects is outside the scope of this manual but the following points should be considered:**

- ✓ Ensure adequate ventilation,
- ✓ Wear eye protection,
- ✓ Wear gloves when handling chemicals,
- ✓ Ensure all liquid containers are secure,
- ✓ No smoking or naked flames,
- ✓ Clean-up all spillages...



*Eye damage could occur from particles or chemicals. Use proper eye protection.*

If a line ruptures or a valve opens accidentally under pressure, a potentially hazardous spray of liquid might be generated by the pump. Never open a solvent line or valve under pressure. Stop the pump first and let the pressure drop to zero. Wear protection glasses.



*Hazardous voltages are present inside the instrument. Disconnect from main power before removing screw-attached panels.*

Removal of protective panels may only be performed by an authorized person. When it is necessary to use a non-original power cord plug, make sure the replacement cord adheres to the color-coding and polarity described in the manual and local building safety codes. Replace blown fuses with fuses of the size and rating stipulated in the manual. Replace faulty or frayed power cords immediately with the same type and rating. Make sure that voltage sources and line voltage match the value for which the instrument is wired.

**To ensure your own safety, thus of your co-workers, and safe operation of the equipment, observe the following instructions:**


- ✓ Use a 3-wire mains socket with ground connection, grounding is necessary to ensure operator safety and proper operation.
- ✓ Take care of inflammable solvents.
- ✓ Store the system after having passed a clean solvent in the hydraulic circuit when not in use. Do not leave buffer in the system, it may cause blockages and damage the seals on start-up.
- ✓ Before moving the instrument, the external connections have to be disconnected, the hydraulic circuit rinsed and the inlets/outlets corked. The instrument must be carefully carried in its original box.

## 6 > Maintenance

### 6.1 Introduction

This chapter describes actions that should be performed on a routine basis to ensure long term safe and trouble free operation of the system. The frequency of the maintenance activities is dependent on the nature of the application (the solvents used, the quantity of the mobile phase delivered by the pump, the level of cleanliness of the facility, etc...). The pump has been designed for reliability and needs very little routine maintenance when operated correctly.

The following figure represents the liquid connections of the pump:



Reduced accuracy of the pump may be due to seal wear or to a fault in the check valves. This section provides information and procedures about how to replace consumable parts and how to maintain the instrument. The removal of protective panels may only be performed by an authorized person.

### 6.2 Removing a pump head

When disassembling or reassembling the pump, make sure that each component is clean and take care that the system is assembled in a clean environment.

**Remove it if you need to replace its seal (details hereafter).**

1. Purge the pump head with a suitable solvent or alcohol (1-propanol or 2-propanol).
2. Switch off the module and disconnect it from the electrical supply.
3. Remove the plastic cover to gain access to the pump heads.
4. Disconnect the Suction and Delivery tubing and fittings.
5. Remove the two screws (Allen wrench 3 mm). Carefully unscrew, alternating from one to the other.



6 - Carefully remove the pump head. The piston and the seal are now visible.

## 6.3 Replacing a seal

### Piston seal




There are in all four seals, one of every head and one of every cleaning disc. You can exchange them by removing the pump heads and cleaning disc as explained previously.


1. Remove the seal using a small screwdriver. Be sure to avoid scoring any part of the pump head and the cleaning disc as a scratch will prevent proper sealing.

*If the seat is scored, it is necessary to replace the head or the cleaning disc and these are not under warranty.*

### Head (back views)



### Cleaning disc (back views)



2. Thoroughly clean the seal seat with alcohol to remove possible seal particles and plunge the new seal in propanol-2 for lubrication.
3. Place the new seal into the pump head or the cleaning disc.
4. Note that the seal is well-placed in the pump head or in the cleaning disc.
5. Finally replace the cleaning disc and the head on the piston and put its two fixing screws with care, alternating from one to the other.
6. Connect the Suction and Delivery tubing and fittings.

### Note :

*The life time of the seals is dependent on the flow rate, pressure, type of liquids used and temperature but mostly on the cleanliness of the mobile phase and sample. The presence of micro-particles will cause accelerated wear and tear of the seals. Similarly any dried buffer particles on the piston will damage the seal.*

## 6.4 Cleaning / replacing a check valve

### Check valve



Abnormal pressure fluctuation or the observation of pump noise on the data system is usually due to check valve problems. The noise is regular in nature and is due to pressure and flow variation. Pump accuracy is also compromised. Details on check valve repairs are given below.

There are four check valves, two per head, one for each inlet (model 20mm, P/N: LO2620) and one for each outlet (model 10mm, P/N: PFS930). If they become dirty, they will not open and close correctly for a long time.

1. Purge the pump head with a suitable solvent or alcohol (1-propanol or 2-propanol).
2. Switch off the module and disconnect it from the electrical supply.
3. Remove the plastic cover to gain access to the pump heads.
4. Disconnect the Suction and Delivery tubing and fittings.
5. Unscrew the check valve housing using a spanner.
6. Clean the check valves or provide yourself with news.

The cleaning allows to remove the particles present in the valves by placing them in an ultrasonic bath with methanol or acetone during approximately five minutes.


**Ensure that the arrow is directed downwards during the cleaning. Else damage to the surface quality or the ball and the seat will occur.**

7. Replace the check valves. Whether it is for the inlet (20mm) or the outlet (10mm), **make sure the arrow is always directed upwards.**



### Inlet (Suction)

### Outlet (Delivery)



8. Hard-tighten the valve housing in the head until hearing a little crunch which means that the new check valve is enough tight and the imperviousness properly done. This noise corresponds to the collapsing of the top and bottom membranes.

9. Connect the Suction and Delivery tubing and fittings.




## 6.5 Replacing the detector deuterium lamp

### Deuterium lamp




1. Purge the system with a suitable solvent, with then alcohol, and with air.
2. Switch off the system and disconnect it from the electrical supply. It will be reconnected only when the panel will be installed back in its place.
3. Let the lamp cool down fifteen minutes after switching it off.
4. Remove tubing and connections from the detector module and remove it from the whole system (you must first remove the pump module).
5. Position the detector module on its side in order to view the bottom and unscrew the three M3 cylindrical-head screws indicated below using an Allen wrench (2.5mm). Take care not to lose M4 flat washers.



6. Put the detector in the upright position and take off the front plastic cover.
7. Then unscrew the two M4 domed-head screws located on the back of the module using an Allen wrench (3mm) and take off the top sheet cover.
8. Loosen the screw on the socket of light tube behind the cell compartment using T-20 Torx adapter. Turn the lamp socket to the horizontal position and tighten the same screw to fix it again.



9. Take the lamp power supply cables and disconnect the 3-pin white connector. Loosen and remove the two screws which hold the lamp using T-10 Torx adapter.





**10- Remove the lamp from the lamp socket and provide yourself with a new one. Never touch quartz bulb of the lamp with bare fingers to avoid damages causing faulty working of the detector.**



11. Carefully insert the new lamp assuring that it is correctly positioned.

**Should you touch the glass of the lamp accidentally, clean it thoroughly with a lint free cloth and alcohol.**

12. Screw it in securely and connect the plug in the 3-pin white socket.

13. Reassemble the unit by making the same steps but in the reverse order.

*After the lamp replacement, recalibration and zeroing of the counter of operating hours must be carried out (cf. Service Info#39).*


## 6.6 Cleaning the detector flow cell

### Flow cell



A "contaminated" cell in lowered light transmission may increase noise level and make difficulties at zeroing. The simplest cleaning method is flushing dismantled cell with suitable solvents. The cell must be taken out from the instrument before washing.

1. Follow the instructions 1 to 7 described in the previous chapter 6.5 "Replacing detector deuterium lamp" in order to properly disassemble the system.
2. Disconnect the outlet tubing from the cell (PEEK nuts).
3. Remove easily the cell from the system by tilting the locking lever, placing one hand under the cell to prevent it from falling to the bottom and tilting the locking lever with the other hand.
4. Select the solvent type according to character of contamination by a series of mutually miscible solvents. It is possible to use both organic and inorganic solvents and diluted solutions of acids.



5. Flush the cell with pure solvent and rinse extensively.

6. Reassemble the module by carrying out the same steps but in the reverse order.

*After the operation, check cell cleanness purging the system and verify there is no baseline drift.*

## 7 > Preventive measures

### Upkeep schedule

By considering the use of clean solvents and all precautions already mentioned in the chapter "6. Maintenance", here is a table which regroup the recommended actions in order to assuring a good working of the system:

Operation	Frequency
<b>Pump</b>	
Check the suitable tightening of fittings (strongly hand tight) and the proper state of tubing (not bended).	Weekly (or when inaccuracy flow rates, solvents leaks or air bubbles are observed)
Replace the tubing and associated fittings.	Yearly (or when damaged parts are observed)
Clean and rinse the gradient valve with suitable solution.	As often as possible (specially after having pumped aqueous solvents or strong buffers)
Replace seals and backup rings.	Yearly
Clean the check valves.	Quarterly
Replace the check valves.	Yearly (or when abnormal pressure / flow rate fluctuation or noisy baseline are observed)
Clean the piston rods.	Quarterly (or when pump head is disassembled)
<b>Detector</b>	
Check the efficiency of the deuterium lamp.	Quarterly (or when higher noise levels or decreased sensitivity are observed)
Replace the deuterium lamp.	When its efficiency is unreliable
Clean the flow cell.	Quarterly
<b>General</b>	
Verify the cleanliness of all liquid containers.	Daily
Check the suitable tightening of fittings (strongly hand tight) and the proper state of tubing (not bended).	Weekly (or when inaccuracy flow rates, solvents leaks or air bubbles are observed)
Clean the injection valves and all parts in contact with solvents or samples with suitable solutions. Rinse abundantly using a mixing of water and alcohol.	As often as possible (specially after having pumped aqueous solvents or strong buffers)
Check the tangential fan working to maintain a good air recycling inside the different modules.	Monthly

## 8 > Warranty

### 8.1 Generalities

Interchim guarantees the *PURIFLASH*<sup>®</sup>800 for one year for parts and labour at the discretion of Interchim in normal conditions of use and installation from the date of signing of the installation report by the service provider approved by Interchim to install the machine and by the customer.

The *PURIFLASH*<sup>®</sup>800 device is guaranteed against material and manufacturing defaults in normal conditions of use by approved professionals and within the technical characteristics compatible with the functions defined in the user manual.

#### I - Terms and conditions of the guarantee:

- ✓ The device is only guaranteed if installed by a service provider approved by Interchim.
- ✓ The guarantee covers the supply of parts found defective by Interchim free-of-charge, as a minimum, within the limits of the parts listed in part 6.2.
- ✓ The device must have been used in normal operating conditions and in accordance with the instructions.
- ✓ The device must have been used with the consumables recommended by Interchim.
- ✓ Exhaustive list of original parts covered by the guarantee. (part 6.2)

#### II - Guarantee exclusions:

- ✓ The guarantee will not cover devices installed by the customer or by a service provider not approved by Interchim.
- ✓ The guarantee will not cover equipment used in a manner which is non-compliant with the provisions in the instructions.
- ✓ The guarantee will not cover equipment subject to interventions, repairs or modifications by personnel without Interchim approval.
- ✓ The guarantee will not cover:
  - parts modified or changed by the customer or the service provider without Interchim's approval.
  - parts damaged by this modification,
  - parts not recommended by Interchim,
  - parts subject to wear and consumables.
- ✓ The guarantee will not cover any parts not listed in part 6.2.
- ✓ The guarantee will not cover electrical and/or electronic and/or IT incidents caused by external factors.
- ✓ The guarantee will not cover damage caused to the software or hardware due to contamination by an IT virus.
- ✓ The guarantee will not cover damage or failures caused by impact and/or bad weather.
- ✓ The guarantee will be cancelled in case of damage caused by abnormal mechanical forces applied to the device and exceeding the limits defined in the user guide.
- ✓ The guarantee will be cancelled in case of corrosion to the device due to solvent leakage or samples.
- ✓ The guarantee will be cancelled in case of corrosion to electronic components caused by highly corrosive gas.
- ✓ The guarantee will not cover damage or failures caused by assembly, dismantling, modification or transport after initial installation.
- ✓ The guarantee will be cancelled in case of failure or damage due to non-compliance with the closing and shutdown procedure for the device.
- ✓ The guarantee will be cancelled in case of failure or damage due to wrong installation
- ✓ The guarantee will be cancelled in case of failure or damage due to wrong AC power supply
- ✓ The guarantee will be cancelled in case of failure or damage due to mechanical force to the unit
- ✓ The guarantee will be cancelled in case of failure or damage due to acts of nature
- ✓ The guarantee will be cancelled in case of failure or damage due to non-compliance with safety procedures.
- ✓ The customer is liable for transport risks.

In case of damage during transport, the beneficiary must issue all reserves to the transport firm before accepting the delivery of the device.

- ✓ The guarantee will not cover damage caused in accidents, external events, contingencies or force majeure, due to negligence, or a lack of surveillance by the customer or due to non-compliance with safety rules.
- ✓ The guarantee will not cover damage if the maintenance procedures recommended by the manufacturer are not complied with.
- ✓ Interchim will not guarantee the results of the use of the *PURIFLASH*<sup>®</sup>800.

## 8.2 Limited warranty of consumable parts

Item	Part Number	Warranty period
<b>Fittings, unions &amp; tubing</b>		
TEFZEL ferrule with Stainless steel lock ring for 1/8" tubing	GV1690	<b>3 months</b>
PEEK 1/4" - 28 short nut for 1/8" tubing for GV1690	GV1700	
PEEK 1/4" - 28 long nut for 1/8" tubing for GV1690	GV1710	
PEEK 5/16" nut + ETFE ferrule for 3/16" tubing	JO8520	
Adapter Luer Lock Male to 1/4" - 28 Male PEEK	JO2630	
Adapter ETFE Luer Lock Female to 1/4" - 28 Male	PFS780	
Adapter Stainless Steel Luer Lock Female to 1/4" - 28 Male	GV1680	
PEEK tee union for 3/16" tubing	LV9960	
Stainless steel bulkhead union 90°	PFS370	
Suction tube (FEP - 1.50m - 3/16" - ID = 3.2mm) + 1 PEEK nut + 1 ETFE ferrule (x 2)	LV9830	
Suction tube (FEP - 1.50m - 3/16" - ID = 3.2mm) + 1 PEEK nut + 1 ETFE ferrule (x 4)	LV9840	
Tube Gradient valve to Tee pump inlet (FEP - 28cm - 3/16" ID=3.2mm) + 2 PEEK nuts + 2 ETFE Ferrules	LV9850	
Tube pump outlet to Pre-column holder inlet (ETFE - 60cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9860	
Tube Pre-column holder outlet to Column holder inlet (ETFE - 90cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9870	
Tube Column holder outlet to detector inlet (ETFE - 70cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9880	
Tube Detector outlet to collector inlet (ETFE - 30cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9890	
Waste outlet tube (ETFE - 1.50m - 1/8" - ID = 2.4mm) + PEEK nut + ETFE ferrule	LV9910	
Tube Collection valve to collection needle (ETFE - 1m - 1/8" - ID=2.4mm) + 1 PEEK nut + 1 TEFZEL Ferrule + spring + pipe	LV9900	
<p><b>&gt; Terms of use :</b></p> <ul style="list-style-type: none"> <li>- Normal tightening of the fittings</li> <li>- Use of clean solvents and clean containers</li> <li>- Follow up of the preventive maintenance actions</li> </ul> <p><b>&gt; Damages not covered :</b></p> <ul style="list-style-type: none"> <li>- Damaged thread or crushed internal part due to unsuitable tightening</li> <li>- Blocked passage due to uncleanliness of installation or unsuitable mobile phase</li> </ul>		





Item	Part Number	Warranty period
<b>Valves</b>		
Binary Gradient valve (CPU2010)	LV9550	<b>6 months</b>
Quaternary Gradient valve (CPU2010)	LV9570	
Collection 3-ways valve	HO7260	
<p><b>&gt; Terms of use :</b></p> <ul style="list-style-type: none"> <li>- Normal tightening of the fittings</li> <li>- Use of clean solvents and clean containers</li> <li>- Follow up of the preventive maintenance actions</li> </ul> <p><b>&gt; Damages not covered :</b></p> <ul style="list-style-type: none"> <li>- Damaged thread or crushed internal part due to unsuitable tightening</li> <li>- Blocked passage due to uncleanliness of installation or unsuitable mobile phase</li> </ul>		
<b>Pump</b>		
Check valve 10mm	PFS930	<b>6 months</b>
Check valve 20mm	LO2620	
Seal 14mm	LO2640	
Piston 14 x 90 + kit	LV9790	
<b>Detector</b>		
Flow cell (0.3mm)	HO7270	<b>1 year</b>
Deuterium lamp	PFS970	<b>6 months or 1000 hours</b>
<p><b>&gt; Terms of use :</b></p> <ul style="list-style-type: none"> <li>- Use of clean solvents and clean containers</li> <li>- Respect precisely the maintenance procedures</li> </ul> <p><b>&gt; Damages not covered :</b></p> <ul style="list-style-type: none"> <li>- Damaged piece due to wrong use or replacement or inappropriate solvent pumping</li> <li>- Blocked passage due to uncleanliness of installation or unsuitable mobile phase</li> </ul>		




## 9 > List of spare parts + Blow-up views

Item	Part Number	N°
<b>Fittings, unions &amp; tubing</b>		
TEFZEL ferrule with Stainless steel lock ring for 1/8" tubing	GV1690	
PEEK 1/4" - 28 short nut for 1/8" tubing for GV1690	GV1700	
PEEK 1/4" - 28 long nut for 1/8" tubing for GV1690	GV1710	
PEEK 5/16" nut + ETFE ferrule for 3/16" tubing	JO8520	
Adapter Luer Lock Male to 1/4" - 28 Male PEEK	JO2630	
Adapter ETFE Luer Lock Female to 1/4" - 28 Male	PFS780	
Adapter Stainless Steel Luer Lock Female to 1/4" - 28 Male	GV1680	
PEEK tee union for 3/16" tubing	LV9960	
Suction tube (FEP - 1.50m - 3/16" - ID = 3.2mm) + 1 PEEK nut + 1 ETFE ferrule (x 2)	LV9830	
Suction tube (FEP - 1.50m - 3/16" - ID = 3.2mm) + 1 PEEK nut + 1 ETFE ferrule (x 4)	LV9840	
Tube Gradient valve to Tee pump inlet (FEP - 28cm - 3/16" ID=3.2mm) + 2 PEEK nuts + 2 ETFE Ferrules	LV9850	
Tube pump outlet to Pre-column holder inlet (ETFE - 60cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9860	
Tube Pre-column holder outlet to Column holder inlet (ETFE - 90cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9870	
Tube Column holder outlet to detector inlet (ETFE - 70cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9880	
Tube Detector outlet to collector inlet (ETFE - 30cm - 1/8" - ID=2.4mm) + 2 PEEK nuts + 2 TEFZEL Ferrules	LV9890	
Waste outlet tube (ETFE - 1.50m - 1/8" - ID = 2.4mm) + PEEK nut + ETFE ferrule	LV9910	
Tube Collection valve to collection needle (ETFE - 1m - 1/8" - ID=2.4mm) + 1 PEEK nut + 1 TEFZEL Ferrule + spring + pipe	LV9900	
Spring for tubes 5cm	PFS800	
Spring for suction tubes 20cm	PFT030	
<b>Electronic boards</b>		
Electronic board CPU2010	IO2950	
Electronic board Collector	PFS500	
System "Flex" for collector connections	PFS490	

Item	Part Number	N°
<b>Pump</b>		
		
Check valve 10mm	PFS930	<b>4</b>
Check valve housing 10mm	PFS010	<b>6</b>
Check valve 20mm	LO2620	<b>3</b>
Check valve housing 20mm	LV9580	<b>2</b>
Seal 14mm	LO2640	<b>5</b>
Head	LV9590	<b>1</b>
		
Piston 14 x 90 + kit	LV9790	
Cleaning disc	LV9800	
Camshaft + 2 Cams	LV9820	
Compression spring	LV9810	
Pump opto-sensor (+ cable)	PFS180	
Opto-sensor disc	PFS190	

Item	Part Number	N°
------	-------------	----

**Pump (suite)**




Belt	LV9970	
Motor (Black)	IO1940	
Motor driver (Motor black)	IO1950	
Axial fan + kit	PFS320	
AC Power supply 230/24V 60W	PFS290	
Pump/Collector fuses 2.5 A-110V (x 2)	PFS990	
Pump/Collector fuses 1.6 A-220V (x 2)	PFS980	



Item	Part Number	N°
<b>Pump (suite)</b>		
Mixing chamber	LV9610	
Seal for mixing chamber	LV9620	
Frit 20µm for mixing chamber (x 2)	PFS260	
Pressure transducer seal	LV9630	<b>7</b>
Pressure transducer + cable	LV9640	<b>6</b>
Leakage sensor	PFS300	
<b>Detector</b>		
Flow cell (0.3mm)	HO7270	
Deuterium lamp	PFS970	
Axial Fan	PFT050	
Power Board	GV3320	
Detector fuse 500mA - 250V	LV7070	

Item	Part Number	N°
------	-------------	----

**Fraction collector**




Positioning guide for racks	PFS400	
Fraction collector motor X-axis + cable	PFS410	
Fraction collector motor Y-axis + cable	PFS420	<b>5</b>
X-axis belt	PFS430	<b>4</b>
Y-axis belt	PFS440	
Pulleys (1 motor + 1 receptor)	PFS450	
Collector opto-sensor X-axis	PFS460	
Collector opto-sensor Y-axis	PFS470	
ON/OFF switch with EMI filter	PFS330	
AC Power supply 230/24V 50W	PFS480	
Pump/Collector fuses 2.5A - 110V (x 2)	PFS990	
Pump/Collector fuses 1.6A - 220V (x 2)	PFS980	
Black POM for column holder fixed (x 2) + screws	IO2960	

Item	Part Number	N°
<b>Column holder</b>		
		
TEFZEL ferrule with Stainless steel lock ring for 1/8" tubing	GV1690	
PEEK 1/4 - 28 long nut for 1/8" tubing for GV1690	GV1710	
Adapter Luer Lock Male to 1/4 - 28 Male PEEK	JO2630	<b>4</b>
Guidance part for large black POM holder	FO4400	<b>3 + 5</b>
Stainless steel union for flash configuration	FV6560	<b>8</b>
Large black POM injection holder + grub screw	IO4960	<b>7 + 10</b>
		
Stainless Steel bulkhead union 90°	PFS370	<b>9</b>
Adapter ETFE Luer Lock Female to 1/4-28 Male or	PFS780	<b>2</b>
Adapter Stainless Steel Luer Lock Female to 1/4 - 28 Male	GV1680	<b>2</b>
PEEK Ring (Outlet column holder)	IO2690	<b>3</b>
Binding column	IO2700	<b>6</b>




Item	Part Number	N°
------	-------------	----

## Pre-column holder



TEFZEL ferrule with Stainless steel lock ring for 1/8" tubing	GV1690	
PEEK 1/4" - 28 long nut for 1/8" tubing for GV1690	GV1710	
Adapter Luer Lock Male to 1/4" - 28 Male PEEK	JO2630	<b>4</b>
Guidance part for large black POM holder	FO4400	<b>3 + 5</b>
Stainless steel union for flash configuration	FV6560	<b>7</b>
Large black POM holder + grub screw	IO4960	<b>6 + 9</b>



Stainless Steel bulkhead union 90°	PFS370	<b>5</b>
Adapter ETFE Luer Lock Female to 1/4" - 28 Male      or	PFS780	<b>2</b>
Adapter Stainless Steel Luer Lock Female to 1/4" - 28 Male	GV1680	<b>2</b>
PEEK Ring (Outlet column holder)	IO2690	<b>3</b>
Screw pre-column holder	IO9140	<b>6</b>

Item	Part Number	N°
<b>General</b>		
Central Processing Unit + Touch screen	PFS520	
USB hub	PFS530	
USB extension cord (l = 0.5cm)	PFS540	
Stylus	FK2400	
AC power supply 230/24V 150W	PFS550	
UC axial fan (V2)	LV9380	
Multi-stage performance EMI filter	PFS570	
General ON/OFF switch	PFS590	
General fuses (x 2) 3.15 A - 220V	PFT010	
General fuses (x 2) 6.3 A - 110V	PFT100	
Serial communication cable male/female	PFS600	
Serial communication cable female/female	PFS610	
Distribution power cable male/female (for main voltage)	PFS620	
Distribution power cable male/female (for 24V)	PFS630	
Power cord FR	PFS640	
Power cord UK	PFS650	
Power cord US	PFS660	

## 10 > Specifications

- ✓ **Model name:** *PURIFLASH*® 800 system
- ✓ **Serial number:** XXXX-X-XX
- ✓ **Flow rate range:** up to 800mL/min
- ✓ **Flow rate performance:** 1% accuracy (with H<sub>2</sub>O degassed at 20°C)  
Repeatability better than 0.5%
- ✓ **Gradient former:** Linear response from 2 to 98%
- ✓ **Injection:** Liquid (manual / loop) or Solid
- ✓ **System protection:** Software, Pressures Min/Max adjustable, Leakage detection
- ✓ **Maximum pressure:** 10 bar (145 psi)
- ✓ **Control:** via PC software "INTERCHIM SOFTWARE VX.X"
- ✓ **Storage:** up to 32 methods in software memory (plus 160Go of hard disk memory)
- ✓ **AC mains supply voltage:** 220-240VAC ; 50-60Hz or 110-120VAC ; 50-60Hz
- ✓ **Electrical protection:** Delayed action fuse 3.15A or 6.3A
- ✓ **Power distribution:** TT or TN power system only, transient over-voltages according to installation category II, Class 1 equipment, requires mains socket plug with protective earth terminal.
- ✓ **IP Code:** IP20
- ✓ **Dimensions:** 680 x 600 x 560mm (Width x Height x Depth without connection)
- ✓ **Weight:** 60kg
- ✓ **Temperature:** 5 to 40°C
- ✓ **Others:** Indoor use, Altitude up to 2000m, Pollution degree 2, Noise level < 35dB.

## ► Appendix 1 : Rack sets available

As mentioned in the chapter **3.5 Prepare the fractions collection**, several racks are usable with the fraction collector. Here is a recapitulative list indicating their designation, the features of tubes used and the total capacity of the rack set :

SHORT RACKS (a rack set includes 4 racks)					
Designation	Tubes size	Tubes volume	Tubes by rack	Tubes by rack set	Total volume
PF4390	18 x 150mm	25mL	44	176	4.4L
	18 x 180mm	32mL	44	176	5.6L
PF4370	13 x 100mm	9mL	60	240	2.1L
PF4380	16 x 150mm	22mL	44	176	3.8L
PF4410	21 x 150mm	35mL	27	108	3.7L
PF4420	25 x 150mm	60mL	16	64	3.8L
PF4430	28 x 150mm	75mL	14	56	4.2L
PF4440	29.5 x 200mm	110mL	14	56	6.1L



**PURIFLASH® 800**

800 ml/min @ up to 10 bar

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:



► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:

► **Service report** ◀

<input type="checkbox"/> Installation	<input type="checkbox"/> End of warranty visit
<input type="checkbox"/> Repair	<input type="checkbox"/> Maintenance

<input type="checkbox"/> Warranty
<input type="checkbox"/> Out of warranty

System description: .....

S/N: ..... Module: .....

Service action: .....  
.....  
.....

P/N and description of defective spare part: .....  
.....

Any further action required: .....  
.....

Date:..... Name:..... Signature:



211 bis, avenue J. F. Kennedy  
BP 1140 - 03100 Montluçon

Tél. +33 4 70 03 88 55  
Fax +33 4 70 03 82 60



Hotline : +33 4 70 03 73 01

E-mail : [interfine@interchim.fr](mailto:interfine@interchim.fr)