


Stimuli™



**Integrated Refraction System
Direct Access Remote Controller
Programmable Tests**

Based on International Standards

Version 4.0 Advanced

Accommodata Corporation
U.S. Patent Number 6,108,634
U.K. Patent Number 2329992
Japan Patent 3302697
Germany Patent 19781691
Canada Patent 2,251,080
Other Patents Pending

Setting Up Your System



Installing the Remote Unit

Simply connect the remote receiver to the supplied cable and plug the cable into any available USB port on the computer or the keyboard. Extension cables for the receiver are available for long examination rooms.

Place the receiver where it is most likely to pick up the signals from the hand controller. Batteries last for approximately 6 months and are available at most drugstores or electronics supply store. Batteries are also available directly from your dealer.

In the event you notice the unit is not always responding to button presses, or that you must point the controller directly at the remote receiver, install a fresh pair of batteries in the hand controller.




Installing the Security Key

To install the Security Key, simply plug the Key into any available USB port on the computer.

Once finished setting up and calibrating the unit, the Keyboard and mouse may be disconnected and stored. If you had installed your Security Key in a USB port on the keyboard, it now must be moved to the spare USB port on the computer.

Starting Up


Your Stimuli system from Accommodata Corporation can automatically turn on and off at preset times. This is adjusted through the “Schedule” option under “Energy Saver” of the System Preferences control panel. Please refer to Apple documentation. You may also start up manually at any time by pressing the Power button on the unit. The system will be available for use in about a minute.

Stimuli is set to launch when the computer is started. If it is not configured this way, or you have exited Stimuli, you may launch Stimuli by pressing the “Star”  button on the hand controller. (If a keyboard and mouse are connected, you may also launch Stimuli by double-clicking its icon on the desktop.)

When you see the Stimuli startup screen press the Next Test  button on the remote control to advance to the tests.

See “Setting Up the System” (above) for additional installation notes.

The Stimuli Remote






Designed for ease of use, the lightweight Stimuli Remote allows wireless control of all functions without the need for a bulky wired controller, keyboard or mouse. The buttons are luminescent for use in low light, while the geometric layout allows the practitioner to make selections without looking at the remote.

Logical groupings of functions make control of each test intuitive. The column of Direct Access Size buttons are also useful for directly accessing all major test suites as described below. Its range is over 7.6M (30 feet) with a wide dispersion angle.



A "Help" Screen and all Stimuli setup screens are available through the "Menu" button on the bottom center of the wireless remote. The Help option will show which buttons affect each particular test.

Controlling the Stimuli System


Control of the Stimuli System is reached through the Menu  button on the wireless remote. To return from any setup screen, press the Menu  button again.



Help	Displays functions of the Wireless Remote Control .
Calibration	Allows access to Calibration controls of the unit.
Configuration	Permits adjustment of Stimuli options & features.
Optotypes	Supports selection of optotypes presented.
Media	Adjusts all film clip display and volumes.
Quit	Allows you to leave the Stimuli program.
Shutdown	Quits and shuts down the computer.


Navigation of the control screen is affected by the Up  and Down  navigation buttons. Simply scroll to the desired feature and press the Select button to enter the chosen control screen.

The Help Screen








The help screen provides explanations of the different tests and the buttons used to control test features. The graphical help for each test is adjusted to the features which are included with the system. For example, here is a help screen for a DuoChrome test.

Calibrating the System



Your system should be delivered calibrated for your examination room from your local dealer. If for some reason you wish to change the calibration or if you have moved the unit and need to adjust the calibration, follow these instructions:


Press the Menu  button on the handheld remote control to reach the main Menu screen. Using the Up  and Down  navigation buttons on the remote, scroll down until Calibration is highlighted and press "Select". You can navigate the calibration screen with your Right  and Left  buttons and then use the Up/Down buttons to adjust that aspect.


The on-screen calibration calculator tells the user the measurable on-screen size. The maximum resolution setting for the computer must be used. (For older CRT displays that have adjustable rasters, the following procedure requires confirming the size of the test "E" by measuring it on the screen with a ruler.)


Calibration of the Stimuli System requires three steps.

Step #1. Adjust display to its maximum resolution.


The maximum resolution setting for the computer must be used.



The computer should arrive already set to its maximum resolution. To confirm this or to adjust the resolution of the screen, launch the System Preferences on top of Stimuli by pressing the Star  button.



Choose the “Displays” function in under “Hardware” to set to the maximum resolution. Please refer to the Apple documentation that came with your computer. If the mouse is not connected to the system simply plug it into a USB port and it will begin working.



ACTION:


Adjust the Resolution of the monitor to its maximum by choosing the setting. Close the Preferences window and click on the Stimuli screen.

OUTCOME:

Your system is now set for the highest resolution.

Step #2. Selecting the Platform.

ACTION:

Ensure that the correct Platform is selected by pressing the Right  Circle button to move the white Highlight Box to the “Platform” pane and select the correct “Platform” or computer you are using. For example, when using a typical iMac-based system choose “iMac_21.5”. The “Other” setting is for any other computer display size such as a larger iMac, or an older platform.






If using “Other” the horizontal and vertical resolutions must be entered in the fields provided. Then screen’s resolution may be entered as: dots-per-inch (DPI); pixels-per-millimeter; or simply by entering the screen’s actual width in inches or millimeters. The other terms will be automatically calculated.

OUTCOME:

The above adjustments have selected the type of computer you are using and provided the screen resolution.

Step #3. Setting the Refraction Distance

Measure the refraction distance from screen to patient's eyes.

Using the Left  and Right  Navigation buttons on the hand controller, move the Highlight Box to the "Foot" adjustment and select the refraction distance to the nearest foot using the Up  and Down  Circle buttons. Using the Right  Circle button, highlight the "Inches" adjustment and adjust the inches for the refraction distance as measured in the room. (If you prefer to work in meters then navigate to the Configuration screen and select Units: "Metric").




For an LCD flat-panel monitor:

Confirm the measurement of the calibration target "E". Measurement is best effected by closing one eye and looking directly at the edge of the letter while aligning a ruler laid on the screen. The indicated measurement should be within a fraction of a millimeter of the measured size.

For a CRT monitor:

To calibrate Cathode Ray Tube monitors which are separate from the computer, using the Height and Width size adjustments on the monitor itself (refer to monitor's User Manual), adjust the size of the calibration target "E" until it exactly equals the calculated size.





If the "Displays" Control Panel is in front of the letter and blocking your view, simply use the mouse to move the window by clicking and holding the title bar at the top of the window and dragging the window to a new position. Note: If you happen to click anywhere outside the Control Panel, then the Stimuli Calibration Setup screen will come to the front, occluding the Control Panel. Simply select the "System Preferences" function again. Close the "Displays" Control Panel. Click once on the Calibration Setup screen to ensure it is in the front to receive commands.


Now press the Menu  button once on the hand controller to get to the Menu screen and using the Down  Navigation button select "Quit" and press the Select  button. The system will save all settings then quit. Simply relaunch as described in "Starting Up" (above) and the Stimuli system is ready to use.

OUTCOME:

You have calibrated your system to a fine accuracy. You can easily check the calibration of your system at any time by going to the Calibration window and measuring the actual refraction distance and the size of the calibration target, and checking it against the denoted size and refraction distance.

Configuring The System

Your Stimuli Vision Station can be adjusted to suit a great number of individual preferences. Highlight a particular feature with the Right  or Left  Navigation buttons then use Up  and Down  buttons to select the option you desire, the system allows changing numerous features.



Calibration units

Optotype Units

Progression

Red/Green Extinction

View

Feedback

Legends

Max. Characters/Line

Max. Lines/Chart

Progressive

Master Volume

English or Metric units for refraction computer.

Snellen, logMAR or Decimal units for optotype legends.

American modified Snellen or true logarithmic chart progression.

Matches colors to Red/ Green Glasses.

Direct viewing or Reflected, (mirrored) examination lane.

Clicks when tests change with volume control, or Silent operation.


Allows selection of Far, Near , Desktop Readout or Audible.

Adjust number of characters on a line.

Allows selection number of lines on multi-line charts.

Equal or Filled charts.



Overall volume control.

In general, the Menu  button steps back up a level each time it is pressed all the way back to the current test chart.

Setup Optotypes


Many optotypes are available. Some are options, others are included and may be turned on or off at the clinician's preference.








Highlight an optotype using the Up  and Down  buttons.

Activate or deactivate optotypes with the Select  button.


The system is configured to allow customizing of all optotype sets.

The Right  Navigation button moves to a customization screen for each hi-lighted optotype.

Choices include: “All optotypes”, “Standard Selection”, or three “Custom” sets.

To make a custom set, use the black Up  and Down  buttons to select “Custom 1, 2 or 3”
Use the Left  and Right  Navigation buttons to make the cursor square traverse the full list of Optotypes, and the  button to select the optotypes desired.


Setup Media



Stimuli's fixation/educational videos set the standard.

Our exclusive Animal Parade is a pediatric favorite with interesting motorized animals that don't wear out or require replacement. The automatic sequencing feature moves the system on to a new Animal every few seconds without intervention by the practitioner. Both child and examiner can stay focused on the task at hand with Stimuli™ Animal Parade.


A set of cartoons is also included.

The Video Play  button is used to Play and Pause the selected video clip in the review window. Clips can be directly selected with the "Size" buttons. A number of Video Clip folders are available:

- The Pediatric Clips folder is included with the system.
- Sets of Multiple Clips may be added.

The Sets of Multiple Clips allow different varieties of video clips to be kept in different sets. This is detailed in "Video Clip Sets" on page 19.

Acuity Charts



The four main varieties of Acuity Chart are: Progressive Lines, Single Line, Single Character, and Multiple Lines of the same size. These are chosen with the Left (←) and Right (→) Navigation buttons.



Sizes are varied with both the Up (↑) and Down (↓) buttons, and the column of Direct Access buttons. Note: additional sizes between those represented by the Direct Access buttons are available, both larger and smaller, with variations in progression (see Setup).

In Chart mode, the black Up (↑) and Down (↓) buttons select the different optotypes available for testing (some of which are optional).



Option Buttons

The two feature buttons (⏏) and (⏏) directly above the navigation buttons provide features which act upon the selected chart. These features include a single column when in chart mode, isolation, contour bars, pointer, and an optional inverse (light on dark).


Pointer

A pointer is available when in Line mode. This helps with many younger, older and foreign speaking patients. It is easily controlled by the  and  buttons along with other test features. The pointer is used to indicate the desired test object in a given line of characters. The pointer automatically sizes to be easily viewable to the patient.



Isolation / Contour Bars


When operating in the **Line** or **Single Character** mode, **isolation** is available. The two feature buttons directly above the navigation buttons  and  selectively provide an isolation box along with the other features. When in “**single character**” mode, contour bars are also available.

Random Charts




For common letter charts, the first set in a size is a common selection of characters. Each time the Select  button is pressed, a new chart of optotypes is presented which eliminates memorization by the patient. The second and third sets of common optotypes are different, but consistent so that the practitioner can eventually memorize them. After the third chart, all charts are randomly generated. This randomly generated selection eliminates the problematic "Memory Effect" of conventional automatic and manual projector systems.

Legends

Legends are controlled by the star  button. Visible legends can be selected in the “Setup” menu. If a remote display panel is used, it will indicate the Chart, Mode, and test details. A second press of the star  button provides a second page of legends for multi-line optotype charts.

Audible Legends may also be selected which announce the size of the first optotype line, and then announce the identity of each optotype upon a second press of the star  button. Young children like this “Talking” test feature which can simplify some testing.

Optional Foot Pedals

Stimuli was first with foot pedal control of charts. This hands-free innovation is a boon to those working with younger patients or patients with special needs. This option is especially useful if the practitioner wishes to use a number of hand held devices such as an occluder. The optional foot pedal unit allows handsfree control of the chart size Up  and Down  and Select  buttons.

Direct Access to Tests

Stimuli has always been the easiest system to use. In addition to providing for direct access to different optotype sizes, and for direct access to film clips when playing, Stimuli supports direct access to tests, This feature is easier than 1-2-3.

1 Press the Pause button.



2 Press any direct access button.


Button numbers correspond to the different available tests:

Programmed Favorite Chart	200
Acuity Charts	160
Red / Green	125
Fixation	100
Astigmatic Dial	80
4-Dot Test	50
Contrast Sensitivity	30
<reserved>	20

Programmability


The initial test screen presented to the patient is progressive from 20/20 up to and including 20/80 (depending on screen size). However, Stimuli also supports a programmable chart starting point. To program your favorite chart, size, optotype, and presentation option into memory:

1. Navigate to your favorite test chart, such as the Pediatric progressive column of optotypes.
2. Press the Pause  button on the remote. The test screen will go into “Sleep” mode (black).
3. Press the bottom left button, Star  button twice, you will hear a Ding/Dong, your test is now programmed.

To return to the programmed favorite chart, from anywhere go into Pause  and then press the 200 direct access button. Stimuli will return to your programmed choice. This feature greatly speeds vision testing and helps keep the patient’s attention.

Note that the buttons both for programming and for navigating to the programmed chart are at the corners of the remote, making reprogramming and navigation a snap.


Sizable 4-Dot Test




Our unique sizable Four Dot chart sets the standard. The test is correctly calibrated for the refraction distance.

The intensity of of the circles may be adjusted for different Red / Green eyewear to assure maximum extinction (via the Four Dot Setup).

This test can be sized with the Up  and Down  Circle buttons.

The external and internal angles in degrees can be displayed or spoken via the  button.

Fixation Dot







Upon entering this test, you will be presented a white dot on a black background. Stimuli includes many features to make this the most flexible and useful fixation device.



The different features supported include:

1. Disappearing Dot Game.
2. White Dot on Black background.
3. Black Dot on White background.
4. Animated dot.
5. Flashing Dot.
6. Combinations.


Pressing the center button makes the dot disappear and reappear.




The  and  feature buttons Blink and Animate the dot respectively. For convenience, the Left  and Right  buttons duplicate the feature buttons.

The animated dot moves in a varying trajectory in the patient's central field.

The Up  and Down  buttons turn the background white or black respectively.

Astigmatic Dial






The dial can be presented and a confirmation dot can be displayed with the Left  and Right  Circle buttons. The indicators can be cleared with the Select  button.

The optional foot pedals duplicate the Left  and Right  circle buttons and the Select  button.


Playing Video




The Video Play  button controls playing and stopping of video clips. When a video clip is playing the Up  and Down  Circle buttons and the Direct Access buttons control selection of additional clips.

The black Up  and Down  buttons are used to adjust volume when a film clip is playing. The Mute  button mutes the audio while the clip continues to play.


Sleep Mode



If video or audio is playing, the Pause  button will pause the video or audio, otherwise the system will be put to sleep.

Video Clip Sets




In addition to the provided set of pediatric cartoons, other sets of clips may be added such as: a set of Patient Education clips, and a set of advertisements for the dispensary may be added. Any number of sets of video clips or still images may be customized for the practice. This is done simply by creating a folder of videos or stills, and placing it in the Videos folder in the Stimuli folder. Stimuli will automatically recognise the new folder and add it to the clip sets.




Sets of multiple video clips can be selected using the  and  feature buttons when a clip is playing.

When a Clip Set is selected and playing, the Direct Access buttons allow direct selection of the first eight clips in each set. The Up  and Down  buttons allow selection of every clip in the set.

Contrast Sensitivity



An Contrast chart, with similar features to the normal Acuity charts, has the ability to lower contrast of the optotypes.

Sizes are selected with the Up  and Down  Circle buttons.
A new line of optotypes is generated each time the Select  button is pressed.

The  and  buttons will adjust contrast level from 100% down to 0.6%.
The contrast level and size can be displayed via the Legend  button.

Software License Agreement

This is a legal agreement between you and Accommodata Corporation covering your use of Accommodata Corporation software. Be sure to read the following agreement before using the software. **IF YOU DO NOT AGREE TO THE TERMS OF THIS AGREEMENT, DO NOT USE THE SOFTWARE AND DESTROY ALL COPIES OF IT.**

The software (including any documentation, electronic keys, images, photographs, animations, video, audio, music, text, and fonts) accompanying this License whether on disk, in read only memory, on any other media or in any other form are licensed to you by Accommodata Corporation. Accommodata Corporation does not transfer title to the software to you.

This License allows you to install and use the software on a single computer at a time. You cannot use the software without obtaining an electronic key from the copyright owner. You may not decompile, reverse engineer, disassemble, modify, rent, lease, loan, sublicense, distribute or create derivative works based upon the software in whole or part or transmit the software over a network or from one computer to another. If you violate any part of this agreement, your right to use this software terminates automatically and you then must return all copies of this software, documentation, and electronic keys to Accommodata Corporation.

The software and related documentation are provided "AS IS" and without warranty of any kind and Accommodata Corporation expressly disclaims all other warranties, express or implied, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose and noninfringement of third party rights. Under no circumstances shall Accommodata Corporation be liable for any incidental, special or consequential damages that result from the use or inability to use the software or related documentation, even if Accommodata Corporation has been advised of the possibility of such damages. In no event shall Accommodata Corporation's liability exceed the license fee paid, if any.

This Agreement shall be governed by the laws of the State of Ohio. If for any reason a court of competent jurisdiction finds any provision of this Agreement, or portion thereof, to be unenforceable, that provision of the Agreement shall be enforced to the maximum extent permissible so as to affect the intent of the parties, and the remainder of this Agreement shall continue in full force and effect.

The software and documentation is provided with RESTRICTED RIGHTS. Use, duplication, or disclosure by the Government is subject to restrictions as set forth in subdivision (b)(3)(ii) of The Rights in Technical Data and Computer Software clause at 252.227-7013. Manufacturer is Accommodata Corporation, 20950 Edgecliff Road, Euclid OH 44123.

Accommodata Corporation
20950 Edgecliff Road
Euclid OH 44123

Under License from Accommodata Corporation
U.S. Patent Number 6,108,634
U.K. Patent Number 2329992
Japan Patent 3302697
Germany Patent 19781691
Canada Patent 2,251,080
Other Patents Pending