


PART V

TUTORIAL



TUTORIAL 1: LOOKING AT GELS

Introduction

This section describes a simple session with the **Melanie II** software using the gels provided on the distribution tape, in the directory `Gels`. You may read this section and at the same time perform the described steps. The `>` sign indicates the command prompt on the Sun workstation.

This tutorial explains how to open gels, then how to manipulate them in order to improve their visualization. *Tutorial 2: analyzing gels* on page 17-1 presents more advanced features of the Melanie II software.

Manipulating two gels.

PC / Macintosh

In Melanie II for **PC or Macintosh**:

- Go to the directory¹ `Gels`, in which the sample gels are located.
- Double-click on the **MelView** icon.

Unix

In Melanie II for **Unix**:

- Go to directory `Gels`, in which the sample gels are located. (for example: `/usr/local/MelanieII/Gels`):

1. On the Macintosh, a directory is called a folder.

```
> cd /usr/local/MelanieII/Gels
```

- Type the following command to run the **MelView** program:

```
> MelView &
```

- The & sign indicates that the MelView program has to run in the background. Figure 16-1 shows the screen that appears when running MelView.




Figure 16-1. The MelView main screen

- Click on the Melanie II logo to get rid of it.

Opening two 2-DE images

PC | In Melanie II for **PC**:

- Click on the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 16-2). Hold down the Control key and click on gels 94-0002 and




Figure 16-2. The Open window allows the selection of the gels that will be opened and loaded into MelView, PC release

94-0005 in the *Files* column. This selects the two gels. Click on the *OK* button. The gels will be opened.

Unix

In Melanie II for **Unix**:

Click on the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 16-3) Hold down the Control key and click on gels 94-0002 and




Figure 16-3. The Open window allows the selection of the gels that will be opened and loaded into MelView, Unix release

94-0005 in the *Files* column. This selects the two gels. Click on the *OK* button. The gels will be opened.

Macintosh

In Melanie II for the **Macintosh**:

- Click on the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 16-4). Click on gels 94-0002 and 94-0005, and each time click




Figure 16-4. The *Open* window allows the selection of the gels that will be opened and loaded into MelView, Macintosh release

on the *Add* button. This selects the two gels and adds them to the *Open gel* box. Click on the *Done* button. The gels will be opened.

- In the *Select* menu, choose the *Gels* → *Visible Gels* option. This will select all gels and draw a green line around each of them showing that they have been selected (Fig. 16-5). In MelView the green color always indicates selected objects.




Figure 16-5. The two opened gels.

Adjusting colors

- Select the Hand tool (the second tool from the top on the left hand side). Put the cursor in one of the gels. Hold down the right mouse button¹ and move the gels until you see a region with many spots.

1. On the Macintosh, hold down the Command key(**⌘**) and the mouse button.

- Now select the Region tool (the third tool from the top in the tool bar) and draw a small rectangle in one of the gels, using the mouse button. In the *Select* menu, choose *Select Mode* → *Adjust Colors*. Select a logarithmic color function and adjust its parameter to the value 4. The change will directly be reflected in the selected region (Fig. 16-6).




Figure 16-6. The *Adjust Colors* option. The modifications are directly reflected in the selected region.

- Click on the *OK* button. The colors will be adjusted in the two selected gels.

Zooming

- With the Select tool active (first tool in the tool bar), select only one gel. Then, in the *Show* menu, choose *Gels* → *Selected Gels*. Only the selected gel stays in the MelView window. The other gel is hidden. In menu *Select*, choose *Select Zoom* → *1/2*. The gel is zoomed out by a factor of 2 and a larger part of the gel can be seen (Fig. 16-7).




Figure 16-7. A gel after it has been zoomed out.

Stacking gels

- Choose *Show Gels* → *All Gels*, and *Select Gels* → *Visible Gels*, then *Stack* → *Stack Selected Gels*. The two gels are now displayed on top of each other. To switch from one gel to the other, choose *Stack* → *Front To Back*, or just Control¹-F. By pressing Control¹-F quickly several times, the differences between the two gels can be observed visually.

1. On the Macintosh, use ⌘ instead.

Aligning gels

We are now going to align the two gels in order to correct their relative distortion. The first step consists in setting a few landmarks in the two gels. Each landmark must correspond to the same feature in both gels. The landmarks must be well distributed over the gel. Follow the next steps:

- Put gel 94-0002 in front. Select the Landmark tool (fourth tool in the tool bar). Press Control¹-F to toggle.
- To add a landmark, double-click, in the center of the gel, on a spot that is also in the other gel. The Create Landmark window pops up. Enter a landmark name (for example L1). Click *OK*. The landmark now appears on the gel image (Fig. 16-8).




Figure 16-8. The L3 Landmark

- Press Control¹-F to toggle to gel 94-0005. To move the top gel so that a spot is exactly above the corresponding spot in the other gel, follow these steps: with the Landmark tool selected, put the cursor in the center of a feature (spot). Hold down the middle mouse button² and at the same time press Control¹-F once and move the cursor to the center of the corresponding spot in the other gel. Release the mouse button.
- Double-click on the same feature (spot) in gel 94-0005. The same landmark name should already be shown in the Create Landmark window. Click *OK*. The L1 landmark will be created on the second gel. Repeat this step four times, in order to add four landmarks around each of the four corners of the gels.
- Press Control¹-R to unstack the gels. Press Control¹-S to save the gels.
- Now that both gels contain related landmarks we are going to perform the alignment. Choose *Process* → *Align Gels* → *Align Gels*. In the *Align Gels* window, select gel 94-0002 as the reference gel. Click *OK*. Gel 94-0005 will then be aligned relative to 94-0002.
- Choose *Stack* → *Stack Selected Gels* to stack the gels again. To visually compare both gels, switch from one gel to the other by pressing Control¹-F quickly several times.

1. On the Macintosh, use ⌘ instead.

2. On the PC mouse with 2 buttons, hold down the Control key and the right mouse button. On the Macintosh, hold down the Option key (⌥) and the mouse button.

Magnifying glass

- Press Control¹-R to unstack the gels. Select the Magnify tool (last tool in the tool bar). Click the mouse button and move the cursor in one of the gels. A small area under the cursor will be magnified. Now hold down the right mouse button² and move the cursor. The same region in both gels will be magnified at the same time (Fig. 16-9).





Figure 16-9. The Magnify tool applied to aligned gels.

Exiting MelView

To exit MelView, do *File* → *Quit*. To try out more advanced features, such as detecting spots and matching gels, you may proceed to *Tutorial 2: analyzing gels* on page 17-1.

1. On the Macintosh, use ⌘ instead.

2. On the Macintosh, hold down the Command key(⌘) and the mouse button.



TUTORIAL 2:

ANALYZING GELS

Introduction

This section describes a typical session with the **Melanie II** software using the gels provided on the distribution tape, in directory¹ `Gels`. You may read this section and at the same time perform the described steps. The directory `Gels-analyzed` contains the same gels after all the analysis described in this tutorial has been performed. You may therefore check what the result of this analysis should be. To do the complete analysis, perform the following steps. The `>` sign indicates the command prompt on the Sun workstation.

Analyzing a set of three gels.

PC / Macintosh

In Melanie II for **PC or Macintosh**:

- Go to the directory `Gels`, in which the sample gels are located.
- Double-click on the **MelView** icon.

Unix

In Melanie II for **Unix**:

- Go to directory `Gels`, in which the sample gels are located. (for example: `/usr/local/MelanieII/Gels`):

1. On the Macintosh, a directory is called a folder.

```
> cd /usr/local/MelanieII/Gels
```

- Type the following command to run the **MelView** program:

```
> MelView &
```

The & sign indicates that the MelView program has to run in the background.

Figure 17-1 shows the screen that appears when running MelView.




Figure 17-1. The MelView main screen

- Click on the Melanie II logo to get rid of it.

Setting preferences

- This tutorial is also available on-line. To consult it on-line:

Unix

On **Unix**, choose *Preferences* in the *File* menu. Set the *Melanie Home* directory to the directory in which the *man* directory resides, usually the same as the Melanie II binaries.

PC / Macintosh

On the **PC** and **Macintosh**, make sure that the Help directory is located in the same directory as the MelView application.

Then, in the *Help* menu, choose *Tutorial*. On-line user manual and index are available too. To access the on-line tutorial and user manual, you must have the Netscape Navigator® World-Wide Web browser installed.

Opening 2-DE images

PC

In Melanie II for **PC**:

- Click in the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 17-2). Hold down the Control key and click on each of the three




Figure 17-2. The Open window allows the selection of the gels that will be opened and loaded into MelView, PC release

file names in the *Files* column. This selects the three gels. Click on the *OK* button. The gels will be opened.

Unix

In Melanie II for **Unix**:

- Click in the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 17-3). Hold down the Control key and click on each of the three




Figure 17-3. The *Open* window allows the selection of the gels that will be opened and loaded into MelView, Unix release

file names in the *Files* column. This selects the three gels. Click on the *OK* button. The gels will be opened.

Macintosh | In Melanie II for the **Macintosh**:

- Click in the *File* menu and select the *Open* option. The *Open* window will be displayed (Fig. 17-4). Click on each of the three gels and each time click on the




Figure 17-4. The *Open* window allows the selection of the gels that will be opened and loaded into MelView, Macintosh release

Add button. This selects the three gels and adds them to the *Open gel* box. Click on the *Done* button. The gels will be opened.

- In the *Select* menu, choose the *Gels → Visible Gels* option. This will select all gels and draw a green line around each of them showing that they have been selected (Fig. 17-5). In MelView the green color always indicates selected objects.




Figure 17-5. The three opened gels.

Adjusting colors

- The 94-0002 and 94-0005 gels are raw, not yet analyzed images. ECOLI is an already analyzed master 2-DE map of *Escherichia Coli* from the **SWISS-2DPAGE** database. In the *Show* menu, choose *Hide All* to hide the already detected and annotated features in the ECOLI image.
- Select the Hand tool (the second tool from the top in the left hand side tool bar). Put the cursor in one of the gels. Hold down the right mouse button¹ and move the gels until you see a region with many spots.
- Now select the Region tool (the third tool from the top in the tool bar) and draw a small rectangle in one of the gels, using the mouse and the button. In the *Select* menu, choose *Select Mode* → *Adjust Colors*. Select a logarithmic color function and adjust its parameter to the value 4. The change will directly be reflected in the selected region (Fig. 17-6).




Figure 17-6. The *Adjust Colors* option. The modifications are directly reflected in the selected region.

- Click on the *OK* button. The colors will be adjusted in all three selected gels.

1. On the Macintosh, hold down the Command key(**⌘**) and the mouse button.

Reporting on gels

- Choose *Show* → *Gels* → *Report On Gels*. A window pops up giving basic information about the selected gels, such as image size in pixels, minimum and maximum grey levels, number of detected features (spots), etc. At this point, only the ECOLI gel has detected features (Fig. 17-7).




Figure 17-7. The *Report On Gels* window.

Printing gels

- You first have to set the printer.

Unix

In Melanie II for **Unix**, select *File* → *Preferences*. Make sure to define an existing printer (example: lp).

PC / Macintosh

In Melanie II for the **PC** or **Macintosh**, select *File* → *Page Setup* and verify the settings.

Then select *File* → *Print* → *Print MelView Window*. The image of the whole MelView window is sent to the printer, it means to the Postscript® printer when using the Unix release.

- Using the Select tool, now select two gels. To achieve this, click on one gel, then hold down the shift key and click on another gel. Choose *File* → *Print* → *Print Selected Gels*.
- In the pop-up window select *Color Adjustment* and click on the *OK* button. The selected gels will be printed in full size.

Setting pseudo colors

- Choose *Select* → *Select Mode* → *Pseudo Colors*. Choose one of the color lookup table, for example *FiveRamps*. Click *OK*. The look-up table will be changed accordingly (Fig. 17-8). Choose *Select* → *Select Mode* → *Grey Levels* to go back to normal grey level display.




Figure 17-8. The gels displayed using the *FiveRamps* color lookup table.

Detecting and quantifying features

- Select the gels that have not yet been analyzed (94-0002 and 94-0005) using the *Select* tool and holding down the shift key. Choose *Process* → *Detect Features*. Click *OK*. Features (spots) will be detected using the default parameters.

Now choose *Process* → *Quantify Features* to compute the optical density, area and volume of the detected features (Fig. 17-9). Select *File* → *Save* → *Save Changes* to save the features and their associated quantification values.





Figure 17-9. Features have been detected. The *Quantify Features* option has been selected.

Reporting on features

- With the Feature tool (fifth on the tool bar) select a few spots: click on one spot. It will be displayed in green, showing that it has been selected. Then hold down the shift key and click on a few other spots. Choose *Show* → *Features* → *Report On Features*. This brings up a report window showing specific information on the selected features (Fig. 17-10).



Report On Features (EC)

| GelName | FeatureID | POL | NPOL | AREA | MSD | OD | pI | pH | M:SW | M:pl | M:W | M:AC |
|---------------------|-----------|-------|-------|-------|-------|-------|------|--------|-----------|------|--------|--------|
| Gels-analyzed-Ecoli | 344 | 0.009 | 0.003 | 0.376 | 0.131 | 0.140 | 4.94 | 104993 | 20-000002 | 4.94 | 104993 | |
| Gels-analyzed-Ecoli | 345 | 0.009 | 0.003 | 0.245 | 0.137 | 0.150 | 5.15 | 103399 | 20-000002 | 5.15 | 103399 | |
| Gels-analyzed-Ecoli | 346 | 0.045 | 0.016 | 0.613 | 0.363 | 0.400 | 5.37 | 103399 | 20-000002 | 5.17 | 103399 | |
| Gels-analyzed-Ecoli | 347 | 0.027 | 0.112 | 1.531 | 1.099 | 1.272 | 5.20 | 101760 | 20-000002 | 5.20 | 101760 | P07395 |
| Gels-analyzed-Ecoli | 348 | 0.250 | 0.075 | 1.011 | 1.039 | 1.179 | 5.69 | 103399 | 20-000002 | 5.69 | 103399 | |
| Gels-analyzed-Ecoli | 349 | 0.142 | 0.051 | 0.898 | 0.737 | 0.852 | 5.71 | 103399 | 20-000007 | 5.71 | 103399 | |
| Gels-analyzed-Ecoli | 350 | 2.899 | 1.047 | 6.860 | 1.832 | 2.119 | 5.75 | 98626 | 20-000009 | 5.75 | 98626 | P17240 |
| Gels-analyzed-Ecoli | 351 | 0.400 | 0.145 | 1.286 | 1.420 | 1.643 | 5.72 | 100495 | 20-000002 | 5.72 | 100495 | |
| Gels-analyzed-Ecoli | 352 | 0.025 | 0.009 | 0.490 | 0.325 | 0.268 | 5.41 | 101125 | 20-000002 | 5.41 | 101125 | |
| Gels-analyzed-Ecoli | 353 | 0.240 | 0.090 | 0.827 | 1.387 | 1.409 | 5.73 | 100495 | 20-000002 | 5.73 | 100495 | |
| Gels-analyzed-Ecoli | 354 | 0.267 | 0.133 | 1.206 | 1.390 | 1.503 | 5.79 | 100495 | 20-000003 | 5.79 | 100495 | P00479 |
| Gels-analyzed-Ecoli | 361 | 0.040 | 0.014 | 0.429 | 0.423 | 0.469 | 5.54 | 98626 | 20-000007 | 5.54 | 98626 | |

Buttons: [Dismiss] [Save] [Print]

Figure 17-10. Report on selected features.

Stacking gels

- Choose *Select* → *Gels* → *Visible Gels*, then *Show* → *Hide All* to hide the detected features. This also has the side effect of displaying the images faster.
- Select the ECOLI and 94-0002 gels. Choose *Show Gels* → *Selected Gels*, then *Stack* → *Stack Selected Gels*. The two gels are now displayed on top of each other. To switch from one gel to the other, choose *Stack* → *Front To Back*, or just Control¹-F. By pressing Control¹-F quickly several times, the differences between the two gels can be observed visually.

Aligning gels

This step is optional, but it is very useful when analyzing gels that are not completely superimposable. Gel alignment corrects gels pixelwise according to a reference gel, in order to make all gels superimposable.

1. On the Macintosh, use ⌘ instead.

The first step consists on setting a few landmarks in each of the three gels. Each landmark must correspond to the same feature in all gels. Start with the two currently stacked gels. The landmarks must be well distributed over the gel. Follow the next steps:

- Put the ECOLI gel in front. Select the Landmark tool (fourth tool in the tool bar). Press Control¹-F to toggle.
- To add a landmark, double-click, in the center of the ECOLI gel, on a spot that is also in the other two gels. The Create Landmark window pops up. Enter a landmark name (for example L1). Click *OK*. The landmark now appears on the gel image (Fig. 17-11).




Figure 17-11. The L3 Landmark

- Press Control¹-F to toggle to gel 94-0002. To move the top gel so that a spot is exactly above the corresponding spot in the other gel, follow these steps: with the Landmark tool selected, put the cursor in the center of a feature (spot). Hold down the middle mouse button² and at the same time press Control¹-F once and move the cursor to the center of the corresponding spot in the other gel. Release the mouse button.
- Double-click on the same feature (spot) in gel 94-0002. The same landmark name should already be shown in the Create Landmark window. Click *OK*. The L1 landmark will be created on the second gel. Repeat this step four times, in order to add four landmarks near each of the four corners of the gels.
- Press Control¹-R to unstack the gels. Choose *Show → Gels → All Gels*. Then, with the Select tool, select the ECOLI and 94-0005 gels. Press Control¹-Z and Control¹¹-T to display the two selected gels in stack mode. Add the same five landmarks to gel 94-0005.
- Choose *Show → Gels → All Gels*, then *Select → Gels → Visible Gels*. Press Control¹-S to save the gels.
- Now that all three gels contain related landmarks, we are going to perform the alignment. Choose *Process → Align Gels → Align Gels*. In the *Align Gels* window, select the ECOLI gel as the reference gel. Click *OK*. Gels will then be aligned relative to the ECOLI gel.
- Once gels have been aligned, many tools and actions that work on one gel also work on all gels at the same time. Instead of pressing the mouse button to use a tool on one gel, press the right mouse button³. The tool will be applied to all aligned gels.

1. On the Macintosh, use ⌘ instead.

2. On the PC mouse with 2 buttons, hold down the Control key and the right mouse button. On the Macintosh, hold down the Option key (⌥) and the mouse button.

3. On the Macintosh, hold down the Command key(⌘) and the mouse button.

Defining regions

- Select the Region tool (third tool in the tool bar). In the ECOLI gel, hold down the mouse button, and drag the mouse for about one inch. Release the button. A small rectangle will be drawn, defining a region on the ECOLI gel. Now hold down the right mouse button¹, drag the mouse. Release the mouse button. The same region will be defined in all gels.

Magnifying glass

- Select the Magnify tool (last tool in the tool bar). Click the mouse button and move the cursor in any of the gels. A small area under the cursor will be magnified. Now hold down the right mouse button¹ and move the cursor. The same region will be magnified in all aligned gels at the same time (Fig. 17-12).




Figure 17-12. The Magnify tool applied to aligned gels.

1. On the PC mouse with 2 buttons, hold down the Control key and the right mouse button. On the Macintosh, hold down the Option key (⌥) and the mouse button.

Matching gels

We will now compare the features (spots) in gels 94-0002 and 94-0005 with those in the ECOLI gel. This process is called **gel matching**. It is going to pair each feature in each of the two gels with the corresponding features in ECOLI. The matching algorithm in Melanie II is able to automatically match features without any manual intervention. Nevertheless, when gels are quite dissimilar, or just to speed up the matching process, you may manually define starting feature pairs. As we have already defined landmarks in all three gels, we will use these as starting pairs.

- Select all three gels, using the Select tool (or Control¹-B).
- Choose *Process* → *Match* → *Match Landmarks*. This will perform pairwise matches on the features containing landmarks (in our case landmarks L1 to L5).
- Choose *Process* → *Match* → *Match Gels*. The Match Gels window pops up, asking for the **reference** gel, i.e. the gel the other selected gels will be matched with. Choose ECOLI and click *OK*. This will now compare all features in gels 94-0002 and 94-0005 to the features in gel ECOLI. Depending on the type of computer you are using, this might take a few minutes.
- After matching has been successfully completed, press Control¹-S to save the pairs.
- There are various ways of looking at the result. First choose *Show* → *Gels* → *Report On Matches*. A window shows how many features have been paired across the gels.

1. On the Macintosh, use ⌘ instead.

- A second way of viewing results of the matching operation is to select a few features in the reference gel (ECOLI) using the Feature tool, then to choose *Select* → *Pairs* → *For Features*. This will select (and display) the paired features in the other two gels (Fig. 17-13).




Figure 17-13. Selected paired features across the three gels.

- Thirdly, you may select two gels, then select a few features on them, then *Show* → *Pairs* → *Report On Pairs* to get the feature IDs of the pairs.
- Finally, select gel 94-0005 and ECOLI, then Control¹-Z and Control¹-T. The two selected gels are now stacked. Choose *Show* → *Pairs* → *Show Pairs*. All the paired features will be shown by small blue vectors (Fig. 17-14).




Figure 17-14. Matched gels: blue vectors indicate paired features.

1. On the Macintosh, use ⌘ instead.

Creating a synthetic gel

Let us now create a synthetic gel by merging the three gels. Prior to the creation of a synthetic gel, pairwise matching has to be performed on all gels involved. Therefore, we still have to match 94-0002 with 94-0005.

- Choose *Show* → *Gels* → *All Gels*, then select gels 94-0002 and 94-0005. Choose *Process* → *Match* → *Match Landmarks*, then *Process* → *Match* → *Match Gels*. Wait for the matching to finish.
- Press Control¹-S to save the pairs.
- Choose *Image* → *Create Gels* → *Create Synthetic Gels*. Then choose ECOLI as the reference gel. This means that the new synthetic gel will be based on ECOLI, adding spots from the other gels.
- In the next pop-up window, give a name for the synthetic gel. We will just call it *Synthetic*. Once created, the synthetic gel will be displayed. You might display it in the full screen to have a better look at it (select the gel, then Control¹-Z). See Figure 17-15.




Figure 17-15. A part of the synthetic gel, displayed with GECOLOR.

Matching gels to the synthetic gel

We are now going to use the synthetic gel as a reference gel, because it contains most of the important features from the three other gels.

- Add the L1 to L5 landmarks to the synthetic gel (See *Aligning gels* on page 17-10).

1. On the Macintosh, use ⌘ instead.

- Match the three gels with Synthetic: Control¹-B, then *Process* → *Match* → *Match Landmarks*. Choose Synthetic as the reference gel. Then *Process* → *Match* → *Match Gels*, and choose Synthetic as reference gel again. The three gels should now be matched with Synthetic (Fig. 17-16).




Figure 17-16. Three selected pairs.

- Once the gels have been matched, you may unalign the gels, in order to work with the original images. All pairs will be preserved: *Process* → *Align Gels* → *Unalign Gels*.

Grouping features

When several gels have been matched with one reference gel, the paired features form so-called **groups**. A group is a set of features, all paired together. Figure 17-16 shows three such groups. Groups may be visualized in reports and in histograms. They may also be exported to files, in order to be used in external programs such as Excel.

- To get a report on groups, choose *Show* → *Groups* → *Report On Groups*. Select a gel as reference gel, for example ECOLI. If a master gel has already been set, the latter will be taken as reference gel (See *Using a master gel* on page 17-17). Select the value type, for example the optical density (OD) or the feature area.

1. On the Macintosh, use ⌘ instead.

- To get histograms on groups, choose *Show* → *Groups* → *Histograms On Groups*. Figure 17-17 shows histograms for two of the selected features. The




Figure 17-17. Histograms for two feature groups.

IDs of the selected groups are shown on the reference gel.

- To export the groups, choose *File* → *Export*. Select the export format. You may export to *Outgroups* format (this is the Melanie 1 table format), Excel tables and Excel lists.

Using a master gel

The ECOLI gel is a SWISS-2DPAGE master gel. This means that several proteins have been identified on the 2-DE map. Corresponding entries exist in the SWISS-2DPAGE database on the ExPASy World-Wide Web server. Each identified protein has a label associated with the corresponding feature.

- Choose *Database* → *Set* → *Master*. Select ECOLI. The ECOLI gel is now defined as the master gel. Its name is printed with a red background.
- Select ECOLI with the Select tool, then Control¹-Z, to have it alone on the screen.
- Choose *Show* → *Labels* → *Show Labels*. All existing labels will appear, showing the SWISS-2DPAGE accession numbers (ACs) of all identified proteins (Fig. 17-18).




Figure 17-18. Labels showing the proteins' accession numbers.

1. On the Macintosh, use ⌘ instead.

- Choose *Show* → *Labels* → *Show Name*. The SWISS-2DPAGE short protein names will be shown instead of the accession numbers (Fig. 17-19).




Figure 17-19. Labels showing the proteins' short names (ID).

- Choose *Show* → *Gels* → *All Gels* to display all gels again, then Control¹-B, then *Select* → *Features* → *All Features*. Finally, choose *Database* → *Labels* → *Add From Master*. Labels will be propagated to paired features. If necessary, choose *Show* → *Labels* → *Show Labels* (Fig. 17-20).



Figure 17-20. Labels propagated to other gels.

Using the Internet

If your computer is connected to the Internet, you may directly access the SWISS-2DPAGE or SWISS-PROT databases on the Expasy server, as well as other databases on various servers.

1. On the Macintosh, use ⌘ instead.

- Choose *Database* → *Set Master* and choose ECOLI. In the **Unix** release, choose *Database* → *Set Browser* and verify that the right name and path are given for Netscape. Choose *Database* → *Set Server*. Enter **www.expasy.ch** and finally, choose *Database* → *Set Database* and enter **SWISS-2DPAGE**.
- Using the Feature tool (fifth tool on the tool bar), select one feature corresponding to an identified protein. Then choose *Database* → *Query Server*. MelView will remotely control Netscape Navigator® and request the corresponding entry from the SWISS-2DPAGE database on Expasy.

Analyzing data

When gels have been matched with one reference gel, groups have been created, and statistical data analysis may be performed on sets of features.

Close all gels that you have been using so far by choosing *File* → *Close* and restart with another set of already analyzed gels of Plasma. Perform the following steps:

PC / Unix ■ In Melanie II for **PC** or **Unix**:

- Click in the *File* menu and select the *Open* option. The *Open* window will be displayed. Double-click in the *Gels Analyzed* directory to open it. Hold down the Control key and click on each of the six gel files called 94-0613, 94-0614, 94-0616, 94-0622, 94-0623 and 94-0624, in the *Files* column. This selects the six gels. Click on the *OK* button. The gels will be opened.

Macintosh ■ In Melanie II for the **Macintosh**:

- Click in the *File* menu and select the *Open* option. The *Open* window will be displayed. Select the *Gels Analyzed* directory by double-clicking in it. Click on each of the six gel files called 94-0613, 94-0614, 94-0616, 94-0622, 94-0623 and 94-0624 and each time click on the *Add* button. This selects the six gels and adds them to the *Open gel* box. Click on the *Done* button. The gels will be opened.
- All six gels have already been analyzed. Gel 94-0623 was used as the reference map.
- In the *Select* menu, choose the *Gels* → *Visible Gels* option. This will select all gels and draw a green line around each of them showing that they have been selected. In MelView the green color always indicates selected objects. In the *Show* menu, choose *Hide All* to hide the already detected and annotated features in the gels.
- Select the Hand tool (the second tool in the tool bar). Place the cursor in one of the gels. Hold down the right mouse button¹ and move the gels until you see a region with many spots.

1. On the Macintosh, hold down the Command key(**⌘**) and the mouse button.

Differential analysis with classes set

To use this analysis option, you need to create two classes of gels.

- To create the first class, pick the **Select** tool in the tool bar. Put the cursor in gel 94-0613 and click on the mouse button. Then hold the shift key on the keyboard and select gels 94-0614 and 94-0616 by clicking on them. This will select the three gels and draw a green line around each of them showing that they have been selected.
- Choose *Analyze → Set Class*. A pop-up window asks you to give a name to this class. Accept the example name, class A, by clicking *OK*.
- Now select gels 94-0622, 94-0623 and 94-0624, as described above.
- Choose *Analyze → Set Class* and accept the example name, class B, by clicking *OK*.
- Now select the **Region** tool (the third tool in the tool bar) and draw a small rectangle in one of the gels, using the mouse and the button.
- To select groups in this region, choose *Select → Groups → In Region*. Select gel 94-0623 as reference gel by choosing *Database → Set → Master* and picking up 94-0623 as reference gel.
- Choose *Analyze → With Classes Set → Differential Analysis*.
- Select the value type, for example the feature area.
- Accept the default ratio value (200) to be considered between the mean value of one class for a given group and the mean value for the other class on the same group (See *Differential analysis with classes set* on page 13-2).
- If there aren't any features selected, increase the rectangle size in one of the gels.
- To get a report on differential analysis results, choose *Show → Analysis → Report On Analysis*. You may print or save the report.

Student T Test with classes set

To visualize the Student T Test results, keep the classes set, classes A and B, and select all groups in the six visible gels.

- Choose *Select → Groups → All Groups*.
- Choose *Analyze → With Classes Set → Student T Test*.
- Select the value type, for example the feature area.
- Accept the default value of the T Test (0.95) to be considered between both classes (See *Student T Test with classes set* on page 13-3).
- To get a report on Student T Test analysis results, choose *Show → Analysis → Report On Analysis*. You may print or save the report.

Differential analysis without classes set

For this analysis method, you have to unset the existing classes.

- Choose *Analyze → Unset All Classes*.

- Now select the Region tool (the third tool in the tool bar) and draw a small rectangle in one of the gels, using the mouse button.
- To select groups in this region, choose *Select → Groups → In Region*.
- Choose *Analyze → Without Classes Set → Differential Analysis*.
- Select the value type, for example the feature area.
- Accept the default value of the ratio (200) to be considered between the highest and the lowest values in groups (See *Differential analysis without classes set* on page 13-4).
- To get a report on differential analysis results, choose *Show → Analysis → Report On Analysis*. You may print or save the report.

Correspondence analysis without classes set

- Keep the same region as in the preceding paragraph and select groups in this region, choose *Select → Groups → In Region*.
- Choose *Analyze → Without Classes Set → Correspondence Analysis*.
- Select the value type, for example the feature area.
- Accept the default value of the number of significant groups that will be displayed in the results (See *Correspondence analysis without classes set* on page 13-4).
- To get a report on correspondence analysis results, choose *Show → Analysis → Report On Analysis*. You may print or save the report.
- You can also visualize the results by displaying a projection of gels and features in the 2-D factorial space. Choose *Show → Analysis → Graphics on Analysis*. Choose to display gels and significant features.

Exiting MelView

To exit MelView, do *File → Quit*.

Other Melanie II features

Many other options and features are available in the MelView program, including:

- Importing/exporting images from and to TIFF, GIF, PPM, etc.;
- Filters for smoothing, contrast enhancement, background subtraction, etc.;
- Managing gel images (flipping, scaling, duplicating, etc.);
- Landmark, feature and label selection by various criteria;
- Landmark, feature, label, pair and group editing.

