

MACOTRONIC V4

USER MANUAL



Rue Lorthiois – 59420 Mouvaux – France
Tel: +33 (0)3 20 11 84 30 - Direct: + 33 (0)3 20 11 86 50 - Fax : +33 (0)3 20 11 84 04

CONTENTS

1	DESCRIPTION.....	1
1.1	GENERAL DESCRIPTION	1
1.2	TECHNICAL SPECIFICATIONS.....	1
1.2.1	Machine Specifications	1
1.2.2	PC Specification:	2
2	INSTALLATION.....	3
2.1	INSPECTION.....	3
2.2	UNPACKING	3
2.3	PACKING LIST	3
2.4	MACHINE INSTALLATION INSTRUCTIONS	4
2.5	WARNING / SECURITY	5
2.6	CONNECTING THE PC AND PERIPHERALS	7
2.6.1	Connecting the PC to the MACOTRONIC:	7
2.6.2	Connection of PC peripherals:	7
3	STARTUP	8
3.1	STARTING UP THE MACOTRONIC	8
3.2	STARTING UP THE PC.....	8
3.3	PROGRAMME LAUNCH	8
3.4	PLACEMENT OF PLASMA PACKS	9
4	SOFTWARE PARAMETERS	11
4.1	ACCESS TO "PARAMETERS" MENU	11
4.2	ILLUMINATION AND TEMPERATURE SETTINGS	12
4.3	SOFTWARE CONFIGURATION	13
4.3.1	User Identification	13
4.3.2	Illumination History Management.....	13
4.3.3	Creation of a text file (ASCII) for illuminations	14
4.3.4	Product Code.....	15
4.3.5	Printing an Illumination report	16
4.3.6	Printing Barcode labels	16
4.3.7	Non-authorized Barcode entry on keyboards	16
4.3.8	Definition / modification of password access to the "Parameters" menu.....	16
5	USE OF THE SOFTWARE	18
5.1	DEFINITION OF DATA RELATING TO AN ILLUMINATION	18
5.2	STARTING AN ILLUMINATION.....	19
5.3	PRINTING BARCODE LABELS	21
5.4	VERIFICATION OF BARCODE	ERREUR ! SIGNET NON DEFINI.
5.5	SAVING ILLUMINATION FILES	23
5.6	HOW TO EDIT THE COMMENTS BOX FOR AN ILLUMINATION	25
5.7	HOW TO START A NEW ILLUMINATION.....	26
5.8	HOW TO CONSULT PREVIOUS ILLUMINATION FILES.....	27

5.9	HOW TO VIEW ILLUMINATION TUBE STATUS	27
6	MAINTENANCE.....	28
7	ERROR MESSAGES.....ERREUR ! SIGNET NON DEFINI.	
8	GARANTEE	35

1 DESCRIPTION

1.1 GENERAL DESCRIPTION

The MACOTRONIC is used for the illumination of plasma following the THERAFLEX-MB Plasma procedure.

The MACOTRONIC is an illumination device with the following characteristics:

- Completely Computer-controlled illumination process,
- Capacity of 4 packs per cycle,
- Illumination System: Low-pressure Sodium Vapour lamps with U-shaped tube, (91 W.)
- Individual control of energy supplied by each lamp,
- Plasma cooled under continual agitation,
- Data Traceability (time, intensity, temperature, lot and donation numbers, lamp control),
- Lid Opening/Closing System

1.2 TECHNICAL SPECIFICATIONS

1.2.1 Machine Specifications

MACOTRONIC®

Illumination device for
methylene blue pathogen
reduction process with
THERAFLEX-MB Plasma



CE 0459

Ref. HMC4P000

Specifications

Dimensions :	Macotronic: 94cmx80cmx110cm (width x depth x height) IT Cart with material: 60cmx80cmx141cm (width x depth x height)
Weight :	280 kg
Capacity :	4 bags/cycle
Power supply :	230 volts ~ (5%)/50Hz/2800 watts
Fuses :	15 A ATM dim 10.3x38
Dissipated energy :	1480 watts
Electrical opening/closing system for the lid	



- **Illumination characteristics**

Light source :	8 sodium low pressure lamps (91 watts/lamp)
Delivered energy :	180J/cm2 (as per THERAFLEX procedure) with individual control of energy delivered by each lamp.
Wave length :	590nm

- **Temperature control**

Integrated cooling system to maintain the temperature inside the illumination shafts
Regulated temperature: +18°C
Programmable temperature alarm (22°C by default)
Bag mixing during illumination

- **Conditions of use**

Air conditionned room: at 21°C

- **Process control & Traceability**

Displayed data : Light intensity, energy delivered, temperature graphs & indicator of advancement
Plasma ID: Donor number, product code, lot number
Illumination report and label printing
Data back-up and export (under ASCII files format)

As an option, process control and traceability ensured by RFID (Radio Frequency Identification device).

- **Accessories (incl. in the ref.HMC4P000)**

Central unit, keyboard, screen, mouse, IT cart, barcode reader, label printer, report printer



- **Regulations :**

Conformity : MDD 93/42 EEC



MacoPharma
F - 59420 MOUVAUX

HMC4P000_A_B 28/12/2006

1.1.2 PC Specification:

Minimal Configuration:

- PC XT/AT compatible
- 1 parallel port
- 1 serial port (RS 232)
- 1 USB port
- WINDOWS 2000 Operating System

2 INSTALLATION

2.1 INSPECTION

MACOTRONIC devices are tested prior to consignment and delivery. It is recommended that the state of the delivered material is verified in order to point out any apparent damage from transporter. Setting up and use of the equipment should follow the instructions in this document stringently.

2.2 UNPACKING

The MACOTRONIC is packed on a special palette and covered in packing material. Remove the packing material as shown below and lower the front part, placing the ramps at the front of the machine. Remove the harness and unlock the brakes on the wheels in front of the machine, then carefully lower the machine from the palette.



2.3 PACKING LIST

The delivered MACOTRONIC system comprises 3 packages:

- MACOTRONIC (Dim: 108x108x130 cm; mass: 330 kg)
- Trolley (Dim: 65x65x135 cm; mass: 50 kg)
- PC unit – printer (Dim: 100x50x50 cm; mass : 30 kg)

The MACOTRONIC carton contains the following pieces:

- MACOTRONIC
- PC – MACOTRONIC connection cable
- Frontal wheel-cover to fit to the machine

The PC/Printer carton contains the following pieces:

- PC
- Screen
- Report printer

- Barcode label printer
- Barcode reader
- Documentation: Certificate of conformity for the machine & software, calibration certificates for light sensors.

2.4 MACHINE INSTALLATION INSTRUCTIONS

WARNING!



Work space

The MACOTRONIC must be used in a temperature-regulated environment, controlled to between 15°C and 21°C.




Electrical connections

Verify that the voltage of your equipment corresponds to that of the available supply (230v~50 Hz).

Connect to the electrical supply using an EARTHED socket.

Installation must be carried out by a qualified technician:

- Place the MACOTRONIC in the designated workspace and lock the wheel brakes.
- Install the frontal wheel-protector, as shown below:



CAUTION: Leave at least 30cm of free space at the rear of the machine in order to allow sufficient air circulation for correct operation of the cooling unit.

- Install the PC and Printer trolley and all associated parts.

2.5 WARNING / SECURITY

- The interrupter switch  puts the system into stand-by but is NOT a general power cut-off switch. In order to disconnect the unit from the mains supply, unplug at the wall socket. Never unplug by pulling on the wire.
- This apparatus must be connected to an earthed socket. Earthing the system provides essential protection for users and the system.
- In order to avoid risk of electrocution, do not expose the unit to rain or excess moisture.
- If a long period of non-use is planned, unplug the equipment at the wall socket. Unplug the cord from the machine using the adapter unit, never pull on the cord itself to unplug.
- Do not obstruct the air intake grill on this apparatus.
- Any use of the equipment contrary to that laid out in the instructions for use documentation may affect the quality of treated products and the safety of users.
- Any and all use of accessories not supplied by the manufacturer may affect the quality of treated products.
- Do not use the MACOTRONIC under dangerous atmospheric conditions or with dangerous materials which have not been designated for use with the MACOTRONIC.

Guide to symbols found on the machine:

ATTENTION



This symbol indicates the equipotential limit to which all parts of the machine are linked




Alternating current



This symbol indicates that the system is on standby



Refer to maintenance and user notes



Risk of electrocution



Moving/rotating part



Risk of crushing

OBSERVATIONS



Unplug the apparatus at the wall socket if the equipment is not to be used for a prolonged period. Unplug only at the adapter point, never by pulling on the cord.



Connect the electrical cord to the supply using an earthed wall socket.



Do not dismantle the machine or attempt repairs. Maintenance is strictly reserved for qualified personnel.

2.6 CONNECTING THE PC AND PERIPHERALS

2.6.1 Connecting the PC to the MACOTRONIC:



Using a serial cable, attach the SubD9b reference "RS 232" on the MACOTRONIC to the PC (serial port)



Important :

All equipment connected to the system must be certified to IEC 950 standard or IEC/ISO applicable to equipment.

2.6.2 Connection of PC peripherals:



The PC is delivered with pre-installed software as standard (installation and updates are carried out by a qualified technician).

3 STARTUP

3.1 STARTING UP THE MACOTRONIC

Switch on the machine by pressing on the interrupter switch, (rear right of the casing).

3.2 STARTING UP THE PC

Ensure that the PC has power from the wall socket.


Verify that all peripherals are powered (printers, screen).

Power up the PC by pressing on the interrupter button on the front panel.

Wait for automatic launch of the MACOTRONIC software. This is the default configuration.

3.3 PROGRAMME LAUNCH

The screen below will appear:



3.4 PLACEMENT OF PLASMA PACKS

Opening and closing the lid

On the front of the Macotronic are 4 buttons. (2 green for opening and 2 red for closure of the lid).

- In order to open the lid, push simultaneously on both green buttons, (during opening, the orange light on the lid will flash.)
- To close the lid, push simultaneously on both red buttons until the lid has completely closed.


CAUTION!



- When the lid is opening, or closing, verify that there is no-one close to the machine.
- Do not place hands / fingers underneath the lid.
- Check that there are no obstructing objects under the lid.


Placement of plasma packs

- Open the two glass panes, (see below):



- Place each of the 4 packs on its own printed area (cf. diagram below).

VIEW FROM ABOVE



Front face of MACOTRONIC


CAUTION

- Once the packs are correctly placed, close the two panes.
- In case of breakdown, the lid may be opened. This type of opening must only be carried out in an emergency, by a qualified technician (Cf. point 6).


4 SOFTWARE PARAMETERS

4.1 ACCESS TO “PARAMETERS” MENU

To access and modify MACOTRONIC software options select the “Software Parameters” menu from the “System” Tab. The following dialogue box (below) will appear:



4.2 ILLUMINATION AND TEMPERATURE SETTINGS



Definition of Settings:

The setting for illumination for the THERAFLEX MB Plasma process is fixed at **180J/cm²**

Verify that the value entered is 180J/cm².

The defined energy delivery is the sole parameter used by the system to validate or invalidate an illumination cycle.

Definition of maximum temperature:

The maximum temperature setting must be fixed at **22°C**.

This value corresponds to a warning level. Temperatures in excess of this value are noted by a cross-hatched area on the illumination graphs.

The operator is warned that the defined temperature has been exceeded during the illumination cycle. In all cases the cycle remains valid if the required energy is delivered.

The state of an illumination cycle is programmed around this warning temperature and will not begin until the temperature is less than 21.5°C (22°C – a variability of 0.5°C, defined by default). This variable is defined as an additional security margin against excess temperature during illumination.

4.3 SOFTWARE CONFIGURATION


4.3.1 User Identification

When this option is selected, user identification is obligatory in order to allow each illumination cycle.

This function provides a way to define multiple users who may obtain secure access to the MACOTRONIC software by means of a password system. At least one user must be defined in order for this option to become available.

To select this option from set-up, click as before on "User Management" and then go on to define a new user (a user name and a password.) The "User Identification" option will now become available.

The user name and the password are now required in order to access the MACOTRONIC software.




4.3.2 Illumination History Management

When the "History management" option is selected in "SOFTWARE PARAMETERS" DIALOGUE BOX, a history file is automatically created for each illumination.


All history data is stored in a file: "ILLUMINATIONS.DAT", created in the following location: C:\Program Files\MacoPharma\Macotronic.

At the end of an illumination cycle, click on "SYSTEM" then on "ILLUMINATION HISTORY" to view the history file.



NOTE: The text editor used to view the history files can read files of 32kb or less. If the file exceeds 32kb, only the latest illuminations are recorded.

4.3.3 Creation of a text file (ASCII) for illuminations



When the option "Export ASCII automatically" is selected in the "SOFTWARE PARAMETERS" dialogue box, an ASCII format data file is automatically created for each illumination.

If this option is not selected, an ASCII file may be created at the end of an illumination by clicking on "FILE", then on "EXPORT ASCII FILE".

Click on "SAVE" in order to save the ASCII file to a user-specified location.

The file is created with a ".TXT" file extension.


This file will be organized with “n” lines (from 1 to 4). For each bag, you will follow the next fields:

Lot number / Donor code / Product code / Starting illumination status / Ending illumination status / Macotronic Serial number / User's name / Illumination date / Illumination starting time / sequential number of cycle / programmed energy / Bag location / illumination duration / maximum temperature reached

The selection of the parameter « Full illumination ASCII export » allows the automatic creation of the ASCII file along with the points of temperature, energy and intensity curves of the illumination process (detailed format on request).

4.3.4 Product Code

When this option is selected, the line: “Product Code” is activated in “Bag Definition” window. This identification field may then be filled in.



4.3.5 Printing an Illumination report

If the “Print Illumination Report” is selected in the software parameters, the report will automatically be printed, as defined by default.

There are two printing modes: Standard (with graphs) and Simplified (without graphs.)

If this option is not selected, it is still possible to print a single report. In order to do so:

Select the Option “PRINT PREVIEW” in the “FILE” menu in order to view the file to be printed.

Click on the  icon to start printing the document, or select “PRINT” on the “FILE” menu.

4.3.6 Printing Barcode labels

If the “Label printing” option is selected in the software parameters, the following window will appear once illumination is finished.

The number of labels per bag has to be specified but (except specific case), this one must be “1”.

On this label, there are : the illumination status (“Invalid treatment” or “MB treated plasma”), the date and time of illumination, Macotronic serial number, lot number, donor code, bag location and file name .MTC.

When the option “automatic label printing” is selected, it will print directly the label without operator intervention.

With this mode, the labels are printed sequentially.

4.3.7 Barcode comparison

This option allows checking that the user has stuck the Macotronic Label on the correct bag.

If this option has been validated, the donor codes checking is requested for each registered bag at the previous step (bags identification phase).

4.3.8 Non-authorised Barcode entry on keyboards

When activated, this option prevents the entry of barcode information by use of the keyboard. The keyboard remains useable for all other applications: i.e. text editing for comments box, etc...

4.3.9 Definition / modification of password access to the “Parameters” menu

By default, the software is delivered without a password. Access to password definition is open t all at this point. Access can be controlled by password, which is recommended:

In the "SYSTEM" menu, select the "PARAMETERS" option

- Click on "Modify Password"
- Enter the new password in the appropriate field
- Click on "OK" to validate



- Confirm the password by repeating the same series of steps

From this point, the new password must be used each time that you wish to access the software parameters section.

In order to delete the system password, use a blank password by entering nothing into the dialogue box.

5 USE OF THE SOFTWARE


5.1 DEFINITION OF DATA RELATING TO AN ILLUMINATION

User Identification


Click on the “USER NAME” input field and enter the text desired using the keyboard.
If the option “User Identification” is selected, this field is automatically filled in (Cf. Section 4.3.1)

Identification of a plasma unit

In the “PACKS” zone, click on the square representing the pack to identify, the field opposite will be the field filled out:



Fill out the fields using the barcode reader and validate the entries by clicking on “OK” in the dialogue box:



It is possible to erase the data and to redefine the identification of a pack by once again on the tick-box representing the number of the pack to be redefined.

Note: The “Product Code” function is optional and may be deactivated in the Software Parameters section (See Section 4.3.4).

5.2 STARTING AN ILLUMINATION

An illumination cycle will not begin until the interior temperature of the system is correct.

Illumination is started by clicking on the “Illumination” button on the software screen (see below).



During the illumination phase, the data shown on screen are kept up to date automatically.

The four bar charts indicate, in real time, the total energy delivered compared to the requested energy (180J/cm²).

Each bar represents the energy supplied by a unit of two lamps.



Localisation of sources:

“B-t”: Unit at rear of machine (Back) – cover section (Top)


“B-b”: Unit at rear of machine (Back) – central section (Bottom)

“F-t”: Unit at front of machine (Front) – cover section (Top)


“F-b”: Unit at front of machine (Front) – central section (Bottom)

The "TEMPERATURE" field independently shows the current temperature inside the MACOTRONIC.

Illumination stops for each unit automatically, once the energy delivered by the tube unit equals the defined quantity ($180\text{J}/\text{cm}^2$). The cycle is finished once the last tube switches off. The following screen appears :



It is possible to visualise the graphs of each unit of tubes by clicking on the tab corresponding to the tube unit, (Source B-t, B-b, etc) :



Illumination data is automatically saved with a predetermined file name (Serial Number of machine_date_illumination number.MTC).
Cf. Section 5.4.

Abort Illumination:


Click on the "STOP" button to abandon an ongoing illumination cycle. An interrupted cycle may not be restarted. The treatment of an interrupted cycle is not valid. Any partially-treated plasma is unusable.

NOTE :

Temperature: If the temperature exceeds the maximum defined temperature during an illumination, the corresponding section of the graph is cross-hatched in red. Illumination is NOT stopped. A warning message appears on the screen and in the report. The labels will still indicate "Plasma treated with MB"
The user should analyse the amount by which the warning temperature was exceeded and determine the validity of the procedure by means of internal protocols.


Viewing instantaneous illumination data:

Place the mouse cursor on the graph: a black line appears in the position of the cursor and the appropriate values are shown on a legend on screen (see below).



5.3 PRINTING BARCODE LABELS & VERIFICATION OF BARCODE

Following to the illumination and if the function “Automatic label printing” is activated, the procedure is the following one :



If there is a comparison error (scan failure or bad printing), the user has the possibility to re-print a label by the set-up.


First comparison

The dialog box is titled 'Barcode matching / Bag n° 1'. It contains two input fields: 'Code for bag n° 1' and 'Code for macotronic label n° 1'.

If comparison failure




2nd comparison



Possibility to re-print the label by a password (Set-up)

3 comparison failures can be authorized. Up to 3, the next bag label will be printed. The following message appears :



Each failure during the barcode comparison will be mentioned on the illumination report :

References :


Macotronic serial number : 0610383
 WARNING : This illumination did not finish correctly.
 An error of communication occurred.
 File: C:\Program Files\Macopharm\Macotronic\Data\0610383_11_27_2007_001.MTC
 User :

Illumination date : 11/27/2007
15:34

	Bag N°1	Bag N°2	Bag N°3	Bag N°4
Donation number	AA12345			
Product code	AA12345			
Bag lot number	AA12345			

Bag 1 : The barcodes are different - Invalid comparison

NB: In the menu "Option", a fonction "re-print" allows to print the label followed by a new comparison window.





NB : All printing action will be mentioned and dated on the illumination report.
In the case of we use the option “label re-print” from the Option menu, the illumination report will be automatically updated and printed.

References :

Macotronic serial number : 0610383

Illumination date : 11/30/2007

WARNING : This illumination did not finish correctly.

11:48

An error of communication occurred.

File : C:\Program Files\MacoPharma\Macotronic\Data\0610383_11_30_2007_001.MTC

User :

	Bag N°1	Bag N°2	Bag N°3	Bag N°4
Donation number	AA12345			
Product code	AA12345			
Bag lot number	AA12345			

Bag 1 : The barcodes are different - Invalid comparison

Label n°1 reprinted the 11/30/07 at 11:48 - Invalid comparison

5.4 SAVING ILLUMINATION FILES

Once an illumination cycle is complete, a file is automatically created in the following way.

MACHINE SERIAL NUMBER _DATE_ ILLUMINATION NUMMBER FOR THE DAY.MTC

For example:

A000001_05_03_2000_007.MTC represents the 7th illumination, carried out on 5th March 2000, on the machine with serial N° A000001.

All saved files are placed within a folder in the following location by default:
C:\Program Files\MacoPharma\Macotronic\Data.

To update an existing file

It is not possible to update existing files.

Files written by the software are in “Read-Only” format.

An illumination file may be saved again by changing the filename, for instance in order to add comments.

5.5 HOW TO EDIT THE COMMENTS BOX FOR AN ILLUMINATION

Comments may be defined for each illumination.

Once comments are added through the comments editor, they may be seen by means of the comments window which is attached to all illumination data

Click on the  icon, or select the "COMMENTS EDITOR" in the "OPTIONS" menu, or via the "Comments" button at the base of the screen.




This editor is used in exactly the same way as other applications such as WordPad. In order to quit the window, select "Close" or "File" and "Exit".

5.6 HOW TO START A NEW ILLUMINATION

NOTE: It is only possible to manage one illumination at the same time.

In order to create a new file for a new illumination run, click on the  icon, or select the “NEW” option from the FILE” menu.

The following window will appear:




In order to start the illumination cycle, follow the steps as described in sections 5.1 to 5.3.

5.7 HOW TO CONSULT PREVIOUS ILLUMINATION FILES

Click on the  icon, or “OPEN” from the “FILE” menu

The various saved illumination files are shown:



Select the desired file and click on the “OPEN” button in the dialogue box


The data are automatically loaded and shown in the MACOTRONIC main window

5.8 HOW TO VIEW ILLUMINATION TUBE STATUS

It is possible to view the total illumination time for each of the tubes, and therefore their status.

In the “SYSTEM” menu, select the "SYSTEM INFORMATION" icon.

The following dialogue box opens:



6 ERROR & WARNING MESSAGES

6.1 FAILURE DURING ILLUMINATION CYCLE

Rappel :


The Macotronic system prints automatically validation labels and report that show if the illumination cycle is well done or not.




Only plasma bags with “MB treated plasma” mention can be considered as well treated according to the Theraflex process.

The plasma bags with « invalidated treatment / do not use » have to be thrown away. This mention means that the cycle has been interrupted / a failure during the cycle (cf. list below).


A treated bag without label has to be considered as an invalidated treatment (case or power outage).

6.1.1 Critical error : Invalid treatment




Message	Comment
<u>Agitating failure</u>	
 A screenshot of a Macotronic software window titled 'Macotronic'. It contains a yellow warning triangle icon and the text: 'WARNING! This illumination has not finished correctly. The agitator stopped during the illumination.' There is an 'OK' button at the bottom.	<u>Cause :</u> Failure of the agitating motion or sensor. <u>Consequence :</u> Illumination stop Data storage with error code (invalid treatment) <u>Corrective intervention :</u> Ask to a qualified technician
<u>Lid failure</u>	
 A screenshot of a Macotronic software window titled 'Macotronic'. It contains a yellow warning triangle icon and the text: 'The cover is not properly closed. The illumination is not valid and is going to be cancelled.' There is an 'OK' button at the bottom.	<u>Cause :</u> Lid failure. Failure of the lid sensor <u>Consequence :</u> Illumination stop Data storage with error code (invalid treatment) <u>Corrective intervention :</u> Ask to a qualified technician

Message	Comment
<u>Light failure :</u> Note: The lamps are controlled after a warning period regarding the 8 sensors (1 photodiode for each lamp).	
	<u>Cause :</u> Light blinking detection (ratio up to 20 between minimum and maximum illumination) <u>Consequence :</u> Illumination shutdown. The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> Ask to a qualified technician for lamp(s) replacement
	<u>Cause :</u> Dispersion between tubes for one shaft (maximum variance in term of illumination duration between 4 sources). At the shutdown of the first lamp, the maximum time is calculated according to the following formula. $\text{TimeMaxIllum} = (\text{TimeIllumTubeoff} - \text{TimeIWarmingTubeoff}) * 1.5$ <u>Consequence :</u> Illumination stops as soon as the calculated time is reached. The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> Ask to a qualified technician for cartridge(s) replacement
	<u>Cause :</u> No homogeneity for a light source (based on a maximum ratio between 2 lamps of a source). The test is active when delivered energy reaches half of the programmed energy. For each source, if the ration is up to 50%, the light source appears in red <u>Consequence :</u> Illumination stop. The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> Ask to a qualified technician for cartridge(s) replacement

	<u>Cause :</u> Maximum illumination time (based on a maximum duration for a cycle). The maximum duration is defined at 40 min to the max intensity (Without warming) <u>Consequence :</u> Illumination stop. The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> Ask to a qualified technician for cartridge(s) replacement
---	--

Message	Comment
<u>Illumination stop by the operator</u>	
	<u>Cause :</u> Click on the « stop » button <u>Consequence :</u> Illumination Stop The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> NA
<u>Communication loss</u>	
	<u>Cause :</u> Communication between Macotronic and PC <u>Consequence :</u> Illumination Stop The failed source appears in red (screen). Data back-up with error code (invalid treatment for every bags) <u>Corrective intervention :</u> Check the connections Check the serial port Check the status of the cable PC-Macotronic

6.1.2 Warnings

Message	Comment
<u>Lamp aging alarm</u>	
	<u>Cause :</u> The cycle duration goes beyond 30min <u>Consequence :</u> None. The cycle is validated <u>Corrective intervention :</u> Sources replacement (the warming period has to be take into account)
<u>Temperature alarm</u>	
	<u>Cause :</u> The temperature inside shaft(s) is up to 22°C <u>Consequence :</u> The overtaking appears on the curves of the interface. It is mentioned in the data file and on the report. Nevertheless, illuminated bags are considered as well treated plasma. Labels have the (MB treated plasma) mention <u>Corrective intervention :</u> Check that the machine is well located according to the user's manual indications
	<u>Cause :</u> The temperature inside the shaft is up to the defined limit (22 - 0,5 = 21.5°C) <u>Consequence :</u> The software waits the decrease of the temperature and the illumination will be started automatically as soon as the temperature is under the limit <u>Corrective intervention :</u> Check that the machine is well located according to the user's manual indications

6.2 OTHER ERROR MESSAGES

Error Message	Corrective Action
 Macotronic Invalid password. Try again or click on "Cancel". OK	Enter the correct password as defined for the user. (cf. section 4.3.1 and 4.3.9)
 Macotronic Barcode formats not specified- Please define them in the barcode definition module. OK	Define the appropriate barcode format in the "Definition of barcodes" section. (cf. User manual for barcode definition)
 Macotronic Communication error !. Check: - The connections between the computer and Macotronic unit. - The Macotronic unit is switched on. Click on "Retry" to restart the illumination. Réessayer Annuler	Communication problem between MACOTRONIC and PC : Verify connections, Verify serial port, Verify state of the PC – MACOTRONIC cable (cf. section 2.6)
 Macotronic Barcode formats not specified- Please define them in the barcode definition module. OK	Barcode format does not correspond to the predefined format.
 Macotronic The codes are different ! OK	The scanned barcodes are not equivalent. Verify that the label has been placed upon the correct pack.
 Macotronic The bags have not been removed from the shafts ! You cannot illuminate them twice ! Replace them by new ones and restart the illumination. OK	Message indicates opening of the lid of the MACOTRONIC. The lid was not correctly shut between the end of the last illumination and the start of the new cycle.

7 MAINTENANCE



When the equipment requires disassembly for repair, (to be carried out only by a qualified technician), unplug the equipment at the wall socket first. Unplug the power cable by pulling on the plug unit and never by pulling on the cable itself.

In order to ensure optimal system function, several recommendations should be followed:

Regular Maintenance:

Cleaning of the interior and glass parts

Use a non-greasy, non-abrasive cleaning or warm water. Dry with a lint-free, clean cloth.

Suggested Frequency: Once per week, minimum.

In order to gain access to the glass parts at the rear, (bottom of the interior); pull out the 2 plates by undoing the two nuts.

CAUTION: Handle the glass plates with care.

After cleaning the two lower glass panes, reposition the plates and retighten the two nuts.

Air Intake Grilles

Never obstruct the air intake grilles. Clean these frequently to avoid blockage.

Suggested Frequency: Once per month.

Laboratory temperature

The MACOTRONIC has a refrigeration unit allowing regulation of its interior temperature during the illumination cycle. However, the MACOTRONIC system must always be used within a temperature-controlled environment (from 15°C to 21°C.)

Maintenance / Preventative:

Replacement Lamps

The MACOTRONIC equipment includes a monitoring system for the function of each amp unit. This allows an illumination to be stopped in case of a large disparity between the energy delivered between lamps. In case of a problem, all 4 units must be replaced.

CAUTION :

Request the intervention of a qualified technician in order to replace the lamps.

Sensor Calibration

The MACOTRONIC is equipped with 8 sensors which detect and quantify the energy delivered by each lamp.

These sensors are integrated with each lamp unit. When the unit is replaced, the sensors are provided pre-calibrated.

The units (lamps + sensors) must be replaced every 24 months by a qualified technician.

Opening the lid (In case of Emergency)

In case of breakdown, it is possible to open the cover manually.

CAUTION :

Manual opening of the lid must not be carried out except in case of emergency except by a qualified technician.

Before any intervention, unplug the machine at the wall socket.

In order to open the lid manually, follow these instructions:

Remove the panel at the rear of the machine.

Loosen the nuts and remove the cover in order to lift the handle of the pneumatic pump.

Tip the handle towards the rear to allow access to the lifting arm.

It is now possible to open the cover of the MACOTRONIC.

CAUTION:

Once the opening procedure is finished, carefully and slowly lower the cover into a closed position once more.

Put back the lever and replace those parts previously removed, effectively carrying out the procedure in reverse.

For all questions regarding maintenance, please contact the after sales service by email at the following address: sav.equipment@macopharma.com or by telephone on: +33 (0)3 20 11 86 50.

8 GARANTEE

The information contained in this manual is validated for the use of the **MACOTRONIC**. MacoPharma will not be held responsible for any errors or omissions which might exist in this manual, or for any consequences arriving from them.

MacoPharma can under no circumstances be held responsible for and damages or problems arising from the use of the apparatus in question outside the specifications laid down in this manual.

MacoPharma will not be held responsible for damages nor problems arising from the use of parts or accessories other than those approved by **MacoPharma**.

Our equipment (not including lamp units) is guaranteed against all defects and errors of construction for one year from the date of delivery, except for incorrect use, negligence or incorrect maintenance and handling.

The return of materiel under guarantee will not be accepted except following written authorisation from the supplier.

CONFORMITY

The medical device for inactivation of viruses in plasma by photo-treatment with methylene (THERAFLEX MB Plasma procedure) conforms to the European Directive 93/42/CEE.



The MACOTRONIC illumination equipment conforms to the following directives and norms :

Directives: 89/336/CEE (CEM),
73/23/CEE (Low tension)
98/37/CE (Machines)

Norms: EN 61326-1 part 1997/A1 part 1997/A2 part 2000 (Measuring, Control and Laboratory Electrical systems)
EN 61010-1: 2001 - 2nd edition (Security rules for electrical apparatus for measuring, regulation and the laboratory.)