Content

| | | |
|---|--|----|
| Introduction | | 3 |
| I. User Login | | 4 |
| II. Operation | | 6 |
| 1. Demography | | 6 |
| 1.1. Basic Information | | 7 |
| 1.2. Life Style Questionnaire | | 8 |
| 2. Anthropometry | | 9 |
| 2.1. Disease and Injury Record | | 10 |
| 2.2. Health Assessment | | 11 |
| 2.2.1. Video Conferencing | | 12 |
| 2.2.2. Blood Pressure | | 13 |
| 2.2.3. Heart Sound | | 14 |
| 2.2.4. Breathing Sound | | 15 |
| 2.2.5. Blood Oxygen Saturation | | 17 |
| 2.2.6. Weight and Body Fat Percentage | | 18 |
| 2.2.7. Skin Image | | 19 |
| 2.3. Injury Survey | | 21 |
| 3. Nutrition | | 22 |
| 3.1. Diet Record | | 23 |
| 3.2. Add Food | | 24 |
| 4. Conditioning | | 25 |
| 4.1. Body Measurement | | 26 |
| 4.2. NSCA measure | | 27 |
| 4.3. Physiotherapy examination | | 28 |
| 4.4 NSCA basic | | 29 |
| 4.5. Flexibility | | 30 |
| 5. Report | | 31 |
| 5.1. Diet Report | | 32 |
| 5.2. Physical Training Report | | 33 |
| 5.2.1. Body Measurement Report | | 34 |
| 5.2.2. NSCA measure Report | | 36 |
| 5.2.3. Physiotherapy examination Report | | 37 |
| 5.2.4. NSCA basic Report | | 38 |
| 5.2.5. Flexibility Report | | 39 |
| 5.3. Injury Survey Report | | 40 |

| | | |
|------|--------------------|----|
| III. | Account Management | 41 |
| 1. | Create User | 42 |
| 2. | Change Login Name | 43 |
| 3. | Change Password | 44 |


Introduction

Welcome for using this Tele-Sports-Health system. Use this guide to find out about the health check empowerment and health record retrieval on this system and to learn more about how they work.

This system is a convenient tool enabling users to conduct their own health assessment and vital signs monitoring. Since it is a light and highly portable box, it is easy for storage and users can carry it anywhere for measuring vital signs such as heart rate, pulse, blood pressure and body mass index (BMI) in comfort, without embarrassment and at their own pace. The data can be stored, thus helping users to better monitor and manage their own health.

This system includes three main components:

- 1) Tele-Sports-Health box case
- 2) System Database
- 3) Scale

I. User Login

1. Click “Start” on the screen.

2. Login the system with **User ID** and **Password**.

3. Main Menu is accessed.

II. Operation

1. Demography

Click “Demography” in the Main Menu.

Two items, **Basic Information** and **Questionnaire**, are displayed.

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2. Anthropometry

Click “Anthropometry” in the Main Menu.

Three items, **Disease and Injury Record**, **Health Assessment** and **Injury Survey**, are displayed.

2.1.Disease and Injury Record

Click “Disease and Injury Record” for inputting information on illness history and injury surveillance.

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2.2. Health Assessment

Click “Health Assessment”.

Seven items, **Video Conferencing, Blood Pressure, Health Sound, Breathing Sound, Blood Oxygen Saturation, Weight and Body Fat Percentage and Skin Image**, are displayed.

2.2.1. Video Conferencing

The purpose of Video Conferencing is to provide the platform for interactive communication between the practitioner and patients.

a) Click “Video Conference” on the screen

- b) The system will connect to the designated practitioner or expert station in School of Nursing, The Hong Kong Polytechnic University, so you can have discussion or consultation at different locations.
- c) You can see the video displaying the image from the side of practitioner. At the same time, your image will also be transmitted to the practitioner's side through the camera attached on the Tele-Sports-Health box case.
- d) You can talk to the practitioner by using the microphone provided and you can hear the practitioner's voice from the speaker.

2.2.2. Blood Pressure

The purpose of blood pressure assessment is to determine whether the blood pressure is within an acceptable range

a) Click “Blood Pressure” on the screen.

b) Subject should sit down and rest his/her arm on a table such that the brachial artery is level with the heart.

c) Wrap the sphygmomanometer cuff around the subject's upper arm, just above the elbow and place a stethoscope on the hollow of the elbow, over the artery.

d) When you are ready for the test, click “Start” on the screen.

e) Wait until the record end.

f) After the completion of the test, your systolic blood pressure, diastolic blood pressure and pulse rate will be shown on the screen.

2.2.3. Heart Sound

The purpose of recording heart sound is to indicate blood flow turbulence.

a) Click “Heart Sound” on the screen.

b) Press the button with “Heart” icon on the electronic stethoscope and a red light can be seen.

c) Put the electronic stethoscope on the left side of the chest, between the 4th and 6th ribs, almost directly under the breast.

d) You can hear a series of magnified “lub-dub” sounds representing each contraction of the heart.

e) Click “Record” on the screen to start recording.

2.2.4. Breathing Sound

The purpose of recording breathing sound is to indicate pneumonia or pulmonary edema.

a) Click “Breathing Sound” on the screen.

- b) Press the button with “Lung” icon on the electronic stethoscope and a blue light will be seen.

- c) Put the electronic stethoscope on the chest.
d) Subject is asked to breath in and out deeply and slowly.
e) Normal magnified “hollow” sound will be heard.
f) Click “Record” on the screen to start recording.

2.2.5. Blood Oxygen Saturation

The purpose of pulse oximetry is to provide estimates of hemoglobin saturated by oxygen and you can have early sign of lung disease or heart disease.

a) Click “Blood Oxygen Saturation” on the main page.

b) Put a subject's finger into the non-invasive sensor clip.

c) When you are ready for the test, click “Start” on the screen.

d) Wait until the record end.

*Note: Do not remove the finger from clip during measurement.

e) After the completion of the test, your assessment results will be shown on the screen.

2.2.6. Weight and Body Fat Percentage

The purpose of measuring the weight is to calculate the body mass index (BMI) which indicating whether your weight is underweight, normal, overweight or even obese.

- Users are required to connect the “USB to Pin Socket Cable” to the computer and the electronic scale. Once connected, the control will be handled automatically by the system and the application.
- Click “Weight” on the screen.

- Subject is asked to step on the scale machine without wearing shoes and socks.
- Readings will be directly transferred to the system for monitoring.
- Height of the subject is loaded if the “Basic Information” in the Demography has been filled in previously.
- After the completion of the test, your body weight, BMI and body fat percentage will be shown on the screen.

2.2.7. Skin Image

The purpose of capturing skin images is to detect skin lesions or to record wounds.


a) Click “Skin Image” on the screen.

b) Pick up the skin image camera and turn it on.

- c) When the light comes on, locate the skin image camera on the site of mole or any wounds on the skin.
- d) Frame the subject to be photographed in the screen.
- e) Click “Capture” to capture the image.
- f) The captured image is photographed on the right frame.
- g) The other capturing image is framed on the left frame.

2.3. Injury Survey

Click “Injury Survey” on the screen for inputting information on injury.

Click on the region of body that is injured, and type in a value representing the level of pain.

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3. Nutrition

Click “Nutrition” on the Main Menu.

Two items, **Diet Record** and **Add Food**, are displayed.

3.1. Diet Record

Click “Diet Record” for inputting information on diet.

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3.2.Add Food

Click “Add Food” to add new food information in the system database.

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Browse a picture of food that you want to add into the system.

Fill in the nutrient profiles regarding to one portion of your added food in the form above.

Click “Send” and then this food item will be shown in the “Diet record” automatically.

If you do not have nutrient profiles of the food, you could send the food photo to Expert Station in School of Nursing, The Hong Kong Polytechnic University for help.

4. Conditioning

Click “Conditioning” on the Main Menu.

Five items, **Body measure**, **NSCA measure**, **Physiotherapy examination**, **NSCA basic** and **Flexibility**, are displayed.

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5. Report

Click “Report” on the Main Menu.

Three items, **Diet Report**, **Physical Training Report** and **Injury Survey Report**, are displayed.

5.1. Diet Report

Click “Diet Report”.

Daily diet record and nutrient profiles can be shown.

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5.2. Physical Training Report

Click “Physical Training Report”.

Five Conditioning reports can be chosen.

5.2.1. Body Measurement Report

Click “Body Measurement”.

Click on the hyperlinks for reports you would like to read.

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When a table is shown, select a parameter first.

Click the calendar icons and select a “Starting Date” and a “Finishing Date” to indicate the period you would like the information generated in your report.

A graph showing the selected parameter in different episodes will be generated.

5.2.3. Physiotherapy examination Report

Click “Physiotherapy examination”

Click on the hyperlinks for reports you would like to read.

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When a table is shown, select a parameter first.

Click the calendar icons and select a “Starting Date” and a “Finishing Date” to indicate the period you would like the information generated in your report. A graph showing the selected parameter in different episodes will be generated.

5.2.4. NSCA basic Report

Click “NSCA basic”

Click on the hyperlinks for reports you would like to read.

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
When a table is shown, select a parameter first.

Click the calendar icons and select a “Starting Date” and a “Finishing Date” to indicate the period you would like the information generated in your report. A graph showing the selected parameter in different episodes will be generated.

5.3. Injury Survey Report

Click “Injury Survey Report”.

Detail information on injury can be shown.

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III. Account Management

Click “Account Management” for personal account setting.

Three items, **Create User**, **Change Login Name** and **Change Password**, can be displayed.

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3. Change Password

Click “Change Password”

Fill in the old password, the new password, and reconfirm the new password.

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