
MA-900

Hall effect magnetic lens power supply

User's Manual



Microanalytical Research Centre
School of Physics
University of Melbourne
VICTORIA 3010
AUSTRALIA

Fax: + 61 (0)3 9347 4783

Ph: + 61 (0)3 8344 5376

Email: marc@physics.unimelb.edu.au

Web: <http://www.ph.unimelb.edu.au/marco>

February 2002

Preface

This manual contains a user's guide to the MA-900 magnetic lens power supply.

For further technical assistance please contact MARC via the following email address:

marc-support@physics.unimelb.edu.au

Limitation of Liability

Micro Analytical Research Centre does not assume any liability arising out of the use of the information contained within this manual. This document may contain or reference information and products protected by copyrights or patents and does not convey any license under the patent rights of Micro Analytical Research Centre, nor the rights of others. Micro Analytical Research Centre will not be liable for any defect in hardware or software or loss or inadequacy of data of any kind, or for any direct, indirect, incidental, or consequential damages in connections with or arising out of the performance or use of any of its products. The foregoing limitation of liability shall be equally applicable to any service provided by Micro Analytical Research Centre.

Note

No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying or otherwise, without the prior written permission of MARC.

Manual Version: 1.0

Manual Date: February 2002

Table of Contents

1. INTRODUCTION	3
2. PRINCIPLE OF OPERATION	5
3. INSTALLATION AND OPERATING INSTRUCTIONS	7
4. SPECIFICATIONS	8
4.1 GENERAL.....	8
4.2 KEPCO JQE UNITS.....	8
APPENDIX 1: JQE GENERAL SPECIFICATIONS	9

1. Introduction

Magnetic lenses in a nuclear microprobe system need a convenient to use, high stability lens power supply. The MA-900 Hall effect unit is designed for magnetic lenses that operate at currents from 10 to 90 Amperes.

The MA-900 magnetic lens power supply was developed to control the magnetic lenses of the Melbourne nuclear microprobe system that employs a Russian antisymmetric quadruplet of MARCO type 4 or type 5 magnetic quadrupole lenses.

MARCO type 2 magnetic lenses typically require currents of up to 90 A at voltages up to 6 V and long-term stability.

The MA-900 magnetic lens power supply consists of the following components:

- MARCO Hall effect lens stabilisation unit and hand controller
- 2 x Kepco JQE 0-6V 0-90A C173 power supply units
- Interconnection cables



Figure 1: The major components of the MA-900 system.

The layout of the major components in relation to the magnetic lens system of the nuclear microprobe is shown below in Fig. 2.




Figure 2: Layout of major system components. Fine line represents control cable, heavy line represents power cable.

The hand controller (Fig. 3) incorporates coarse and fine current controls for both channels. When not in use, following fine focusing of the lens system, it should be stowed in a safe place to prevent mechanical disturbance which may cause the lenses to drift out of focus.



Figure 3: The MA-900 hand controller

The MA-900 unit itself is shown in Fig. 4. From left to right, the items are: current in channel 1 (Amperes), current in channel 2 (Amperes), socket for hand controller and mains on/off switch with pilot light.



Figure 4: The MA-900 Hall effect stabilisation unit.

The rear of the MA-900 unit is shown in Fig 5.



Figure 5: The rear current connections for the MA-900 unit.

From left to right the connections from the rear are: 240 VAC input, channel 2 control signal, Channel 1 control signal, Channel 2 high current connections, Channel 1 high current connections.

The MA-900 contains a precision two channel Hall effect transducer clamped on the high current supply cable. The Hall effect transducer is temperature compensated and provides a precision feedback signal that does not require the use of a current sensing series resistor.

The unit is to be connected up as in Fig. 6.

From tests in the MARC laboratory, no noticeable drift in the beam spot size, above 1 micrometre, occurs over time-spans of 24 hours.




Figure 6: Wiring diagram for channel one (channel 2 is similar). Note the 7.5 k resistor and the links on the control block of the Kepco JQE unit.

2. Principle of operation

The current leads from the Kepco units pass through the core of an internal Hall effect current transducer that measures the current and provides a voltage feedback signal. This is compared to an internal reference, controlled by the handpiece, to generate a programming voltage to the Kepco units.

The Hall effect current transducer is a LOHET2 unit which is temperature compensated with an offset stability of $\pm 0.02\%/^{\circ}\text{C}$ and a gain stability of $\pm 0.02\%/^{\circ}\text{C}$. The LOHET2 is installed in a ferrite core transformer where it operates in the range of -400 to $+400$ Gauss to produce a voltage output of 5 mV/Gauss.

The zero current setting of the MA-900 may be adjusted by trim-pot R5 (see Fig. 7 for the location of this trim-pot on the feed back circuit boards). The procedure for this is as follows:

1. Connect the MA-900 unit to the lens system as described in the previous section.
2. With the mains power off and the cable disconnected, remove the top cover of the MA-900 unit.
3. Set the controls on the hand controller to zero.
4. **CAUTION:** Potentially hazardous voltages are present inside the MA-900 unit when powered. Observe all necessary precautions. Qualified service personnel only should perform the next steps.

5. Power on all units.
6. With an appropriate screwdriver, adjust R5 until the output current reads zero on the front panel digital displays.
7. Power off all units and reassemble.




Figure 7: The circuit board of one channel of the MA-900 Hall effect stabilisation unit showing the location of the R5 trimpot.

3. Installation and operating instructions

Please observe all normal precautions for electrical equipment and lifting heavy apparatus, in addition to the instructions here, when installing the MA-900 unit.

These instructions should be read in conjunction with the separate manuals for the Kepco JQE units themselves.

CAUTION: Ensure the correct mains voltage has been selected on each unit before connection to the mains.

To install and power-up the system the following procedure must be followed:

1. Ensure all power switches on all units are off. Mount the Kepco units in suitable 19" rack, together with the MA-900 stabilisation unit. Connect all units into the mains.
2. Connect the Kepco units to the MA-900 stabilisation unit and connect the loads as shown in figure 3. Be sure to observe the polarity of the lenses.
3. Plug the handpiece into the MA-900 Hall effect unit.
4. Ensure the potentiometers on the hand-piece are set to zero.
5. Turn the current limit control on the Kepco units to maximum.

CAUTION: It is essential that steps 6 and 7 are followed in the correct order to ensure safe operation!

6. Power on the MA-900 Hall effect unit.
7. Power on the Kepco units.
8. Adjust the controls on the hand-piece to set the desired current in the lenses.
9. Allow the system to warm-up for 30 minutes and re-adjust as necessary.

To power off the system:

1. Make a note of the settings of the controls of the hand-piece and the currents in each channel for future reference.
2. Turn the controls on the hand-piece to zero.

CAUTION: It is essential that steps 3 and 4 be followed in the correct order to ensure safe operation

3. Switch off the Kepco units.
4. Switch off the MA-900 unit.

4. Specifications

4.1 General

The MA-900 magnetic lens power supply is designed to operate with a MARCO nuclear microprobe system with type 2 lenses. The probe forming lens system consists of an MP2 Russian antisymmetric magnetic lens system or a high excitation quintuplet system of type 4 lenses.

The MA-900 power supply parameters are:

- Two independent channels of power supply output.
- Output: 0 – 90A, continuously adjustable by the coarse and fine control controls on the hand-piece, at 0-6V
- Stability: After the warm-up time, will maintain 1 μ m focus on any lens system with excitation aberration coefficients less than 260 μ m/(mrad %).

4.2 Kepco JQE Units

The power supply parameters are (independent of the MA-900 stabilisation unit):

- 10 turn voltage control for exceptional resolution.
- Analog output control by resistance: 1000W/Volt; or by a voltage delivering 0-1mA.
- Digital listen only control using SN-series digital interfaces.
- Current limited, front panel control (not programmable) 10%-105% I_o max.
- JQE can control current with an external current-sense resistor.

Appendix 1: JQE General Specifications


JQE GENERAL SPECIFICATIONS			
SPECIFICATION		RATING/DESCRIPTION	CONDITION
INPUT			
a-c Voltage		105-125, 210-250V a-c	User selectable
Current		See model table	Max load, 115V a-c
Frequency		47-65Hz	Range
OUTPUT			
d-c Output		Series pass	Transistor
Type of stabilizer		Voltage stabilizer	Current limited
Voltage		0 to 100% of rating	Adjustment range for temp 0-71°C
Current		0 to 100% of rating	
Error sense		0.5V per load wire	Static voltage allowance
Isolation voltage		500V d-c or peak	Output to ground
Leakage current		<5 microamperes	rms at 115V a-c
Output to ground		<50 microamperes	p-p at 115V a-c
Series connection		500V	Max voltage off ground
Parallel connection		Automatic	Use current mode limiting
		Current sharing	Use master-slave connection
		Redundancy type	External steering diodes
OVP		Not available	
CONTROL			
Type	Voltage	Fixed input, variable gain	
	Current	Differential comparison	
Voltage	Local	10-turn precision rheostat	
	Remote analog	1000 ohms per volt or 0 to 1mA control current	
	Remote digital	Use SN/SNR interface	12 bit Listen-only
Current	Local	Multiturn pot	
	Remote analog	Not provided	See Series ATE models
Dynamics	Normal (slow)	$dV/dt = I/C$	See tabulated value of C in the model table
	Fast mode	Not provided	See Series ATE models

JQE MODEL TABLE


MODEL	D-C OUTPUT RANGE		OUTPUT IMPEDANCE				MAX. INPUT AMPS @ 125V _{ac}
	VOLTS	AMPS	VOLTAGE MODE		CURRENT MODE		
			SERIES R	SERIES L ⁽¹⁾	SHUNT R	SHUNT C ⁽²⁾	
QUARTER-RACK							
JQE 6-10M	0-6	0-10	30μΩ	1μH	50kΩ	3kμF	2.0
JQE 15-6M	0-15	0-6	125μΩ	1μH	84kΩ	1kμF	2.1
JQE 25-4M	0-25	0-4	300μΩ	1μH	125kΩ	700μF	2.2
JQE 36-3M	0-36	0-3	600μΩ	1μH	165kΩ	400μF	2.2
JQE 55-2M	0-55	0-2	1.4mΩ	1μH	250kΩ	220μF	2.3
JQE 75-1.5M	0-75	0-1.5	2.5mΩ	1μH	330kΩ	160μF	2.3
JQE 100-1M	0-100	0-1	5mΩ	2μH	500kΩ	110μF	2.1
HALF-RACK							
JQE 6-22M	0-6	0-22	14μΩ	0.5μH	23kΩ	5.8kμF	4.2
JQE 6-45M	0-6	0-45	7μΩ	0.5μH	11kΩ	8kμF	9.0
JQE 15-12M	0-15	0-12	63μΩ	0.5μH	42kΩ	2.7kμF	4.0
JQE 15-25M	0-15	0-25	30μΩ	0.5μH	20kΩ	4.5kμF	8.4
JQE 25-10M	0-25	0-10	125μΩ	0.5μH	50kΩ	2.4kμF	5.3
JQE 25-20M	0-25	0-20	63μΩ	0.5μH	25kΩ	4.3kμF	10.5
JQE 36-8M	0-36	0-8	225μΩ	0.5μH	62.5kΩ	1.4kμF	6.0
JQE 36-15M	0-36	0-15	120μΩ	0.5μH	33kΩ	3.6kμF	9.5
JQE 55-5M	0-55	0-5	550μΩ	1μH	100kΩ	850μF	5.0
JQE 55-10M	0-55	0-10	275μΩ	1μH	50kΩ	2.1kμF	9.0
JQE 75-3M	0-75	0-3	1.25mΩ	1μH	165kΩ	850μF	4.0
JQE 75-8M	0-75	0-8	469μΩ	1μH	62.5kΩ	1.2kμF	10.0
JQE 100-2.5M	0-100	0-2.5	2mΩ	1μH	200kΩ	600μF	4.5
JQE 100-5M	0-100	0-5	1.25mΩ	1μH	100kΩ	600μF	8.4
JQE 150-1.5M	0-150	0-1.5	5mΩ	2μH	330kΩ	440μF	4.6
JQE 150-3.5M	0-150	0-3.5	2.2mΩ	2μH	140kΩ	440μF	8.7
FULL-RACK							
JQE 6-90M	0-6	0-90	3.5μΩ	0.5μH	3.5kΩ	17.6kμF	15.7
JQE 15-50M	0-15	0-50	15μΩ	0.5μH	10kΩ	12kμF	16.6
JQE 25-40M	0-25	0-40	31μΩ	0.5μH	12.5kΩ	14kμF	21.0
JQE 36-30M	0-36	0-30	60μΩ	0.5μH	16kΩ	11kμF	19.0
JQE 55-20M	0-55	0-20	138μΩ	1μH	25kΩ	7.3kμF	18.0
JQE 75-15M	0-75	0-15	250μΩ	1μH	33kΩ	4.2kμF	18.0
JQE 100-10M	0-100	0-10	0.62mΩ	1μH	50kΩ	2.2kμF	17.0
JQE 150-7M	0-150	0-7	1.1mΩ	2μH	72kΩ	1kμF	18.0

(1) For determining dynamic impedance in voltage mode.


(2) For determining dynamic impedance in current mode.




FRONT VIEW



BOTTOM VIEW



REAR VIEW



SIDE VIEW

JQE MECHANICAL SPECIFICATIONS		
SPECIFICATION	RATING/DESCRIPTION	CONDITION
Input connection	Detachable IEC type 3-wire (1)	¼ and ½ rack size
	Permanently wired (2)	Full rack size
Output connections	Front panel binding posts	Models under 15A
	Rear barrier strip	¼ and ½ size
	Rear compression studs	Full rack size
Meters	Two 1½" vertical 3%, analog	Front panel
Indicators	Neon	Pilot
Mounting (in std 19" racks)	Use RA 24 rack adapter	¼ and ½ size
	Mounting "ears" provided	Full rack
Cooling	Forced air	Exhaust to rear