

Procarta™ Transcription Factor Whole Cell Lysis Kit

User Manual

P/N 13939 Rev. A 061119



Panomics, Inc.

Procarta Transcription Factor Whole Cell Lysis Kit User Manual

Copyright

© Copyright 2006, Panomics, Inc. All rights reserved.

Trademarks

Procarta is a trademark of Panomics, Inc.

Citing Procarta in Publications

When describing a procedure for publication using this product, we would appreciate it if you would refer to it as the Procarta™ TF Whole Cell Lysis Kit.

If a paper cites a Procarta product and is published in a research journal, the lead author(s) may receive a travel stipend for use at a technology conference or tradeshow by sending a copy of the paper to our technical support group at techsupport@panomics.com or via fax at (510) 818-2610.

Disclaimer

Panomics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice.

Although this manual has been prepared with every precaution to ensure accuracy, Panomics, Inc. assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.

Contents

About the User Manual	5
Who Should Read this Manual.	5
What this Manual Covers	5
Safety Warnings and Precautions	5
For More Information	5
Procarta Transcription Factor Whole Cell Lysis Assay Kit	5
About the Transcription Factor Whole Cell Lysis Kit	5
Kit Contents and Storage	5
Required Materials and Equipment Not Provided	6
Materials and Equipment	6
Cell Preparation	7
Growing Cells	7
Whole Cell Lysis Procedure	7
Assay Guidelines	7
Preparing Working Reagent	7
Preparing Whole Cell Lysates From Adherent Cells	8
Preparing Whole Cell Lysates From Suspension Cells	9
Contacting Panomics	10
Technical Help	10
For Additional Services	10

About the User Manual

Who Should Read this Manual Anyone that has purchased a Procarta Whole Cell Lysis Kit from Panomics to prepare whole cell lysates for use in Panomics' Procarta Transcription Factor (TF) Assay Kits.

What this Manual Covers This manual provides the following:

- ◆ Kit contents
- ◆ Required materials and equipment
- ◆ Whole cell lysis procedure

Safety Warnings and Precautions **CAUTION** All chemicals should be considered potentially hazardous. We recommend that this product and its components be handled by those trained in laboratory techniques and be used according to the principles of good laboratory practice.

Note This product is intended for research use only.

For More Information For information about the Procarta products mentioned in this manual, visit our website at www.panomics.com.

Procarta Transcription Factor Whole Cell Lysis Assay Kit

About the Transcription Factor Whole Cell Lysis Kit The Procarta TF Whole Cell Lysis Kit contains reagents and procedures for the preparation of whole cell lysates for use in our Procarta TF Assay Kits. The Procarta TF Whole Cell Lysis Kit contains sufficient reagents for the preparation of 40 whole cell lysates from cultured cells grown in 6-well culture plates or 240 whole cell lysates from cultured cells grown in 96-well culture plates.

Kit Contents and Storage The Procarta Transcription Factor Whole Cell Lysis Assay Kit contains the following components.

Procarta TF Whole Cell Lysis Kit components:

Component	Quantity	Storage
Lysis Buffer I	12.0 mL	-20 °C
Lysis Buffer II	1.2 mL	-20 °C
DTT	120 µL	-20 °C
Protease inhibitor cocktail	120 µL	-20 °C

Required Materials and Equipment Not Provided

**Materials and
Equipment**

Item	Source
1X PBS	Invitrogen (P/N 14190-144)
Rocking platform	VWR, Rocking Platform, Model #100 or equivalent
Centrifuge	Eppendorf #5804R
Protein determination kit	Bio-Rad DC Protein Assay Kit (P/N 500-0112) or equivalent
Microcentrifuge tubes	Major laboratory supplier (MLS)
15 mL conical centrifuge tubes	MLS
Adjustable single and multi-channel precision pipettes	MLS
PCR plates	MLS

Cell Preparation

Growing Cells In all cases, cells are grown to about 90% confluence. The following table provides recommendations for the cell requirements for each culture vessel type. However, it is important to realize that cell types vary in size and actual numbers of cells/vessel may vary.

Use the table below as a guide.

Culture Vessel	Cell Number
100-mm culture dish	1 x 10 ⁷ cells/dish
6-well plate	1–5 x 10 ⁶ cells/well
24-well plate	2–5 x 10 ⁵ cells/well
96-well plate	2–5 x 10 ⁴ cells/well

Whole Cell Lysis Procedure

Assay Guidelines **IMPORTANT** All components and PBS must be kept on ice at all times. Lysis Buffer 1 Working Reagent must be kept on ice and should be used within 2 hours of preparation.

Preparing Working Reagent

To prepare working reagent:

Step	Action										
1	<p>For preparation of 1 mL Lysis Buffer 1 Working Reagent, combine the following reagents, then invert tube to mix:</p> <p>1 mL Lysis Buffer 1</p> <p>10 µL DTT</p> <p>10 µL Protease Inhibitor Cocktail</p> <p>Scale Lysis Buffer preparation according to experimental design. Use the table below as a guide.</p> <table border="1" data-bbox="591 1419 1390 1650"> <thead> <tr> <th>Vessel</th> <th>Quantity of Lysis Working Reagent</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>1 mL/dish</td> </tr> <tr> <td>6-well plate</td> <td>250 µL/well</td> </tr> <tr> <td>24-well plate</td> <td>100 µL/well</td> </tr> <tr> <td>96-well plate</td> <td>50 µL/well</td> </tr> </tbody> </table>	Vessel	Quantity of Lysis Working Reagent	100 mm culture dish	1 mL/dish	6-well plate	250 µL/well	24-well plate	100 µL/well	96-well plate	50 µL/well
Vessel	Quantity of Lysis Working Reagent										
100 mm culture dish	1 mL/dish										
6-well plate	250 µL/well										
24-well plate	100 µL/well										
96-well plate	50 µL/well										

Preparing Whole Cell Lysates From Adherent Cells

To prepare working reagent:

Step	Action										
1	Remove the culture media from all wells and wash cells twice with an appropriate volume of cold 1X PBS.										
	<table border="1"> <thead> <tr> <th>Vessel</th> <th>Quantity of PBS</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>10 mL/dish</td> </tr> <tr> <td>6-well plate</td> <td>1 mL/well</td> </tr> <tr> <td>24-well plate</td> <td>500 μL/well</td> </tr> <tr> <td>96-well plate</td> <td>200 μL/well</td> </tr> </tbody> </table>	Vessel	Quantity of PBS	100 mm culture dish	10 mL/dish	6-well plate	1 mL/well	24-well plate	500 μ L/well	96-well plate	200 μ L/well
	Vessel	Quantity of PBS									
	100 mm culture dish	10 mL/dish									
	6-well plate	1 mL/well									
24-well plate	500 μ L/well										
96-well plate	200 μ L/well										
2	Following the second wash, make sure the PBS is completely removed.										
3	Add the appropriate volume of Lysis Buffer I Working Reagent to the wells.										
	<table border="1"> <thead> <tr> <th>Vessel</th> <th>Quantity of Lysis Buffer I Working Reagent</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>1 mL/dish</td> </tr> <tr> <td>6-well plate</td> <td>250 μL/well</td> </tr> <tr> <td>24-well plate</td> <td>100 μL/well</td> </tr> <tr> <td>96-well plate</td> <td>50 μL/well</td> </tr> </tbody> </table>	Vessel	Quantity of Lysis Buffer I Working Reagent	100 mm culture dish	1 mL/dish	6-well plate	250 μ L/well	24-well plate	100 μ L/well	96-well plate	50 μ L/well
	Vessel	Quantity of Lysis Buffer I Working Reagent									
	100 mm culture dish	1 mL/dish									
	6-well plate	250 μ L/well									
24-well plate	100 μ L/well										
96-well plate	50 μ L/well										
4	Transfer culture vessel(s) to an ice bucket and transfer ice bucket to a rocking platform at 200 rpm for 10 minutes.										
5	Add the appropriate volume of Lysis Buffer II Working Reagent to the wells.										
	<table border="1"> <thead> <tr> <th>Vessel</th> <th>Quantity of Lysis Buffer II Working Reagent</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>100 μL/dish</td> </tr> <tr> <td>6-well plate</td> <td>25 μL/well</td> </tr> <tr> <td>24-well plate</td> <td>10 μL/well</td> </tr> <tr> <td>96-well plate</td> <td>5 μL/well</td> </tr> </tbody> </table>	Vessel	Quantity of Lysis Buffer II Working Reagent	100 mm culture dish	100 μ L/dish	6-well plate	25 μ L/well	24-well plate	10 μ L/well	96-well plate	5 μ L/well
	Vessel	Quantity of Lysis Buffer II Working Reagent									
	100 mm culture dish	100 μ L/dish									
	6-well plate	25 μ L/well									
24-well plate	10 μ L/well										
96-well plate	5 μ L/well										
6	Transfer culture vessel(s) to an ice bucket and transfer ice bucket to a rocking platform at 200 rpm for 1 hour.										
7	<p>Pipet up and down several times, then transfer each sample to a 1.5 mL microcentrifuge tube and centrifuge at 14,000 x g for 3 minutes at 4°C. Note the orientation of tube in centrifuge as pellets may not be visible.</p> <p>For 96-well plates, pipet up and down several times, then transfer well contents from the cell culture plate to a PCR plate and centrifuge the PCR plate at 2,250 x g for 5 minutes.</p>										
8	Transfer supernatant(s) to new microcentrifuge tubes or a new PCR plate, this is your whole cell lysate.										

To prepare working reagent: *(continued)*

Step	Action										
9	Measure the protein concentration of each sample using a protein quantitation assay (sold separately). Then prepare 5 μ L aliquots for each of the samples. Store samples at -80°C or use immediately in Procarta TF Assay Kit.										
	<table border="1"> <thead> <tr> <th>Vessel</th> <th>Typical Protein Yields</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>300–800 $\mu\text{g}/\text{dish}$</td> </tr> <tr> <td>6-well plate</td> <td>200–500 $\mu\text{g}/\text{well}$</td> </tr> <tr> <td>24-well plate</td> <td>40–80 $\mu\text{g}/\text{well}$</td> </tr> <tr> <td>96-well plate</td> <td>20–40 $\mu\text{g}/\text{well}$</td> </tr> </tbody> </table>	Vessel	Typical Protein Yields	100 mm culture dish	300–800 $\mu\text{g}/\text{dish}$	6-well plate	200–500 $\mu\text{g}/\text{well}$	24-well plate	40–80 $\mu\text{g}/\text{well}$	96-well plate	20–40 $\mu\text{g}/\text{well}$
	Vessel	Typical Protein Yields									
	100 mm culture dish	300–800 $\mu\text{g}/\text{dish}$									
	6-well plate	200–500 $\mu\text{g}/\text{well}$									
24-well plate	40–80 $\mu\text{g}/\text{well}$										
96-well plate	20–40 $\mu\text{g}/\text{well}$										

Preparing Whole Cell Lysates From Suspension Cells

To prepare whole cell lysates from suspension cells:

Step	Action										
1	Transfer cells to a 1.5 mL or 15 mL centrifuge tube as appropriate and centrifuge at 500 x g for 5 minutes. For 96-well plates, transfer cells to a PCR plate and centrifuge at 500 x g for 5 minutes										
2	Remove the culture media and wash cells by resuspending in 1 mL of cold 1X PBS followed by centrifugation at 500 x g for 5 minutes. Repeat wash step. For 96-well plates, remove the culture media and wash cells by resuspending in 200 μ L of cold 1X PBS followed by centrifugation at 500 x g for 5 minutes. Repeat wash step. Following the second wash step, ensure that the 1X PBS solution is completely removed from the cells. Note Before the second centrifugation, transfer contents from the 15 mL centrifuge tube to a 1.5 mL microcentrifuge tube.										
3	Immediately add the appropriate volume of Lysis Buffer 1 Working Reagent to cell pellets. Mix by pipetting up and down several times. <table border="1"> <thead> <tr> <th>Original Culture Vessel</th> <th>Lysis Buffer 1 Working Reagent</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>1 mL/dish</td> </tr> <tr> <td>6-well plate</td> <td>250 $\mu\text{L}/\text{well}$</td> </tr> <tr> <td>24-well plate</td> <td>100 $\mu\text{L}/\text{well}$</td> </tr> <tr> <td>96-well plate</td> <td>50 $\mu\text{L}/\text{well}$</td> </tr> </tbody> </table>	Original Culture Vessel	Lysis Buffer 1 Working Reagent	100 mm culture dish	1 mL/dish	6-well plate	250 $\mu\text{L}/\text{well}$	24-well plate	100 $\mu\text{L}/\text{well}$	96-well plate	50 $\mu\text{L}/\text{well}$
Original Culture Vessel	Lysis Buffer 1 Working Reagent										
100 mm culture dish	1 mL/dish										
6-well plate	250 $\mu\text{L}/\text{well}$										
24-well plate	100 $\mu\text{L}/\text{well}$										
96-well plate	50 $\mu\text{L}/\text{well}$										
4	Transfer tube(s) or culture plate to an ice bucket and place on a rocking platform at 200 rpm for 10 minutes.										

To prepare whole cell lysates from suspension cells: *(continued)*

Step	Action										
5	Add the appropriate volume of Lysis Buffer II to the cells.										
	<table border="1"> <thead> <tr> <th>Original Culture Vessel</th> <th>Quantity of Lysis Buffer II Working Reagent</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>100 µL/dish</td> </tr> <tr> <td>6-well plate</td> <td>25 µL/well</td> </tr> <tr> <td>24-well plate</td> <td>10 µL/well</td> </tr> <tr> <td>96-well plate</td> <td>5 µL/well</td> </tr> </tbody> </table>	Original Culture Vessel	Quantity of Lysis Buffer II Working Reagent	100 mm culture dish	100 µL/dish	6-well plate	25 µL/well	24-well plate	10 µL/well	96-well plate	5 µL/well
	Original Culture Vessel	Quantity of Lysis Buffer II Working Reagent									
	100 mm culture dish	100 µL/dish									
	6-well plate	25 µL/well									
24-well plate	10 µL/well										
96-well plate	5 µL/well										
6	Transfer tube(s) or culture plate to an ice bucket and place on a rocking platform at 200 rpm for 1 hour.										
7	For microcentrifuge tubes, pipet up and down several times and centrifuge at maximal speed (12,000 x g) for 3 minutes at 4°C. Note the orientation of tube in centrifuge as pellets may not be visible. For PCR plates, pipet up and down several times and centrifuge the PCR plate at 2,250 x g for 5 minutes.										
8	Transfer supernatant(s) to a new microcentrifuge tube or a new PCR plate, this is your whole cell lysate.										
9	Measure the protein concentration of each sample using a protein quantitation assay (sold separately). Then prepare 5 µL aliquots for each of the samples. Store samples at -80 °C or use immediately in Procarta TF Assay Kit.										
	<table border="1"> <thead> <tr> <th>Vessel</th> <th>Typical Protein Yield</th> </tr> </thead> <tbody> <tr> <td>100 mm culture dish</td> <td>300–800 µg/dish</td> </tr> <tr> <td>6-well plate</td> <td>200–500 µg/well</td> </tr> <tr> <td>24-well plate</td> <td>40–80 µg/well</td> </tr> <tr> <td>96-well plate</td> <td>20–40 µg/well</td> </tr> </tbody> </table>	Vessel	Typical Protein Yield	100 mm culture dish	300–800 µg/dish	6-well plate	200–500 µg/well	24-well plate	40–80 µg/well	96-well plate	20–40 µg/well
	Vessel	Typical Protein Yield									
	100 mm culture dish	300–800 µg/dish									
	6-well plate	200–500 µg/well									
24-well plate	40–80 µg/well										
96-well plate	20–40 µg/well										

Contacting Panomics

Technical Help For technical questions, contact our technical support group by telephone at 1-877-726-6642 option 3 or by email at techsupport@panomics.com (US and Canada) or techsupport_europe@panomics.com (Europe), or visit our website www.panomics.com for an updated list of FAQs and product support literature.

For Additional Services For information about Panomics products or for ordering information, contact your Regional Sales Manager, or visit our website at www.panomics.com.
