



# Intelect<sup>®</sup> Neo

Therapeutic excellence, at the touch of your finger

The Intelect Neo is the new standard in physical medicine modalities. Its intelligent design (the result of over 100,000 hours of intensive R&D) is clever in its features, usability and clinical technology. Delivering an exceptional patient and therapist experience. Every element of the Intelect Neo has been expertly crafted, empowering the clinician to provide patients with an unsurpassed level of rehabilitation.

Intelect Neo offers multimodality with 6 plug-and-play modules. Each unit can be assembled specifically for your customized clinical needs. And your unit can also adapt to your future needs by adding additional modules. The plug-and-play modules are easily installed in no time, without worrying about configuration and other settings.

Neo is also simple to use and a pleasure to navigate, thanks to the Clinical Protocol Setup™ (CPS), which leads you through the functions of the device and each therapy.

Additionally, Neo offers a stunning anatomic library that illustrates an array of pathologies, making it easier for you to communicate with patients about their condition and educate them on further treatment options.

































A completely modular setup allows you to choose among the components that best fit your practice. The unit has an integrated base with discreet, strong handles for carrying.

The optional high-quality cart is stable, height-adjustable, moves with ease, and includes three roomy storage drawers with sturdy pull tracks.

Clinical excellence and increased productivity, Intelect Neo gives you the freedom to have both.



# Electrotherapy Waveforms

Waveforms	Description	Output Intensity	Phase Duration	Frequency (Hz)
Interferential Current IFC (4-Pole)	 Interferential Current is distributed through two channels (four electrodes). The currents cross each other and interfere, resulting in a modulation of the intensity (the current intensity increases and decreases at a regular frequency).			2000-10000 (Carrier) 0 - 200 (Beat)
Premodulated IFC (2-Pole)	 Premodulated Current comes out of one channel (two electrodes). The current intensity is modulated: it increases and decreases at a regular frequency.			2000-10000 (Carrier) 1 - 200 (Beat)
VMS™	 A symmetrical biphasic waveform with a 100 µsec interphase interval. The short pulse has a low skin load, ideal for high intensities applications, such as muscle strengthening.	0-255 mA	20-1000 µsec	1-200
VMS™ Burst	 A burst version of VMS™	0-255 mA	20-1000 µsec	1-200
VMS™ FR	 A version of VMS where physiologically based channel interaction in which one channel stimulates the agonist & the other the antagonist of the muscle that is being exercised.	0-150 mA/ 0- 150 V	20-400 µsec	1-200
Asymmetric Biphasic TENS	 This waveform has a short pulse duration. It is capable of strong stimulation of the nerve fibers in the skin as well as of muscle tissue.	0-110 mA	20-1,000 µsec	1-250
Symmetric Biphasic TENS	 A waveform with a short pulse duration. It is capable of strong stimulation of the nerve fibers in the skin as well as of muscle tissue.	0-80 mA	20 - 1,000 µsec	0-250
Alternating Rectangular	 An interrupted biphasic current with a rectangular pulse shape. Commonly used as a pain management application.	0-100 mA	20-1,000 µsec	0-250
Microcurrent	 Microcurrent is a monophasic waveform of very low intensity. The literature reports beneficial effects of this waveform in the treatment of wounds.	0-1000 µA		
Han Stimulation	 This waveform provides optimal parameters with a precisely controlled sequence of Dense-and-Disperse (DD) modes of stimulation, simultaneously releasing all 3 kinds of opioid peptides leading to an enhanced analgesic effect.	0-80 mA/ 0-80 V	20-180 µsec	1-2
High Voltage Pulsed Current	 This monophasic waveform has a very brief pulse duration with two distinct peaks. The high voltage causes a decreased skin resistance making it comfortable and easy to tolerate.	0-500 V		10-120
Monophasic Rectangular TENS	 This waveform is an interrupted unidirectional current with a rectangular pulse shape. Commonly used with electro diagnostic testing & to stimulate Denervated muscle.	0-110 mA	20-1,000 µsec	1-200
Monophasic Rectangular: Pulsed	 An interrupted unidirectional current with a rectangular pulse shape.	0-80 mA	0.1-500 ms	
Monophasic Triangular: Pulsed	 An interrupted unidirectional current with a triangular pulse shape.	0-80 mA	0.1-500 ms	
Monophasic Rectangular: Surged	 A series of rectangular, monophasic pulses. The pulses surge to maximum power, hold and then decrease before the pause. This waveform well suited for muscle strengthening.	0-80 mA	0.2-5.0 ms	5-60
Monophasic Triangular: Surged	 As above but with a triangular shaped pulse.	0-80 mA	0.2-5.0 ms	5-60
Träbert (Ultrareiz)	 A monophasic waveform with a phase duration of 2 ms & a pause of 5 ms resulting in a frequency of ~ 143Hz.	0-80 mA	2 ms	~143
Galvanic: Continuous	 A direct current following in one direction only.	0-80 mA		
Galvanic: Interrupted	 A direct current following in one direction only.	0-80 mA		
Russian	 A sinusoidal waveform, delivered in bursts or series of pulses. Claimed to produce maximal muscle strengthening effects without significant discomfort to the patient.	0-100 mA		
Diadynamic Monophasic (MF)	 The Diadynamic waveforms are rectified alternating currents. Frequency of 50 Hz: phase duration of 10 ms followed by a pause of 10 ms.	0-80 mA	10 ms	50
Diadynamic Diphasic (DF)	 Frequency of 100 Hz: phase duration of 10 ms followed immediately by another identical phase of 10 ms.	0-80 mA	10 ms	100
Diadynamic CP	 1 second of monophasic followed by 1 second of diphasic.	0-80 mA		
Diadynamic CP-iso	 A combination of monophasic & diphasic waveforms.	0-80 mA		
Diadynamic CP-id	 A combination of monophasic & diphasic waveforms.	0-80 mA		
Diadynamic LP	 Rhythmical % uctuation between 2 monophasic currents.	0-80 mA		
Diadynamic MF + CP	 A period of MF followed by a period of CP.	0-80 mA		
Diadynamic MF + CP-id	 A period of MF followed by a period of CP-ID.	0-80 mA		
Diadynamic DF + LP	 A period of DF followed by a period of LP.	0-80 mA		
Diadynamic DF + CP	 A period of DF followed by a period of CP.	0-80 mA		
Iontophoresis	 The introduction of ionizable drugs through intact skin by the administration of continuous, direct electrical current.			
S/D Curve				

# Technical Specifications

## ULTRASOUND

Frequency:	1 MHz, ± 5%; 3.3 Mhz, ±5%
Duty Cycles:	10%, 20%, 50%, Continuous
Pulse Repetition Rate:	16, 48 or 100 Hz
Pulse Duration:	1 mSec, ±20%; 2 mSec, ±20%; 5 mSec, ±20%
Output Power:	Large Crystal: 0-15 W at 1 MHz, 0-10 W at 3.3 MHz • Medium Crystal: 0-6W @ 1 and 3.3 MHz • Small Crystal: 0-3 W @ 1 and 3.3 MHz
Amplitude:	0 to 2.5 w/cm <sup>2</sup> in continuous mode, 0-3 w/cm <sup>2</sup> in pulsed modes
Output accuracy:	± 20%, 10% of maximum
Temporal Peak to Average Ratio:	2:1, ± 20%, at 50% Duty Cycle • 5:1, ± 20%, at 20% Duty Cycle • 9:1, ± 20%, at 10% Duty Cycle
Beam Nonuniformity Ratio (BNR):	6:1 maximum Beam Type, Collimating
IPXX Rating for Unit:	IPX0
IPXX Rating for Applicator:	IPX7
Effective Radiating Areas (ERA):	Large Crystal: 5.0 cm <sup>2</sup> (minimum) • Medium Crystal: 2.0 cm <sup>2</sup> (minimum) • Small Crystal: 1.0 cm <sup>2</sup> (minimum)
Treatment Time:	1 to 30 min

### Head Warming Feature

The Head Warming feature of a Intellect® Neo Therapy System utilizes Ultrasound output, resulting in warming of the Applicator to increase patient comfort. With Head Warming enabled, ultrasound is emitted without pressing the Start button while an ultrasound treatment is being setup. The Applicator LED will not illuminate during the Head Warming period. US Channel will indicate "Head Warming".

Output:	0 - 50% Cycling of maximum power
Frequency:	3.3 MHz
Applicator Temperature:	29.4 °C - 43.3 °C (85 °F - 110 °F)

## LASER


















Output Type:	Infrared Lamp (Laser)
Laser Class:	3B

### Laser Technical Specifications

Pulse Frequencies:	8 Hz - 10000 Hz and continuous
Wavelengths:	670-950 nm (dependent on applicator)
Output:	100-1440 mW (dependent on applicator)
Output accuracy:	+/- 20% of nominal

### Description of Device Markings

The markings on the unit are assurance of its conformity to the highest applicable standards of medical equipment safety and electromagnetic compatibility. One or more of the following markings may appear on the device:

 Refer to Instructional Manual Booklet	 Equipment capable of delivering output values in excess of 10 mA r.m.s. or 10V r.m.s. averaged over any period of 5s						
 Testing Agency	 Dangerous Voltage	 Electrical Type B	 Electrical Type BF	 Laser	 Ultrasound	 Stim	
 Start	 Stop	 Pause	 Intensity	 Lock/Unlock	 ON/OFF	 Laser Stop Switch	
 This unit is considered to be a Class 3B laser product and thus emits visible and invisible laser radiation (lR). Avoid direct eye exposure to the Laser beam. The symbol to the right is located on the back of the applicator and indicates the active radiant surface (the area on the applicator that emits infrared laser energy and the direction of the beam of light)							

### Laser Applicator Technical Specifications

For all single diode and cluster laser and LED applicators, the expected increase in the measured quantities after manufacture combined with the values measured at the time of manufacture is ± 20%. The software incorporates a cooling function that forces the user to cool the laser clusters prior to the next treatment when treatment times exceed 3 minutes per application. The software will calculate the cooling time needed as follows:

- For a 3 minute treatment, it will force a 15 second cool down period
- For a 4 minute treatment, it will force a 2 minute cool down period
- The software extrapolates for times between 3 and 4 minutes

Note: A message will display for 5 seconds on the screen informing the user that the probe is cooling down and the time period required. If the user attempts to use the probe before the cool down period is completed, the message re-displays.

When cool down is complete, a message displays the unit is ready for use.

# Technical Specifications

## VACUUM

**Vacuum Range:** 0 to 600 mbar maximum +/- 5%  
**Vacuum Modes:** Continuous or Pulsed in 10 steps

**Continuous:**  
10 setting over vacuum range, 60mbar setting+10mbar to -10mbar per setting

**Pulsed Mode:**  
Maximum Vacuum settings 2 to 10, +10mbar to 10mbar per setting. Minimum Vacuum settings in 1 to 9, +10mbar to -10mbar per setting Hold Time in minimum & maximum vacuum settings, 0-20 seconds, in 1 second steps, +/-0.5 seconds

**Power:**  
Input: 20-25 Vdc  
4A Electrical Class: CLASS I  
Electrical: TYPE BF

## POWER (Combination and Electrotherapy Units)

**Mains:** 100 - 240 V AC- 1.0 A, 50/60 Hz  
**Electrical Class:** CLASS I  
**Mode of Operation:** Continuous

**Electrical Type (Degree of Protection):**  
Ultrasound: TYPE B  
Laser: TYPE B  
Electrotherapy: TYPE BF  
Electrotherapy & sEMG: TYPE BF  
Electrotherapy & Vacuum: TYPE BF  
Ultrasound & Electrotherapy: TYPE B

**Note:**  
All waveforms except High Voltage Pulsed Current (HVPC) have been designed with a 200 mA current limit. VMS™ Waveforms and all TENS waveform output intensities are measured, specified, and listed to peak, not peak to peak.

## GENERAL SYSTEM OPERATING & STORAGE TEMPERATURE

**Operating conditions:**

- Temperature: 10° C to 45° C
- Relative Humidity: 0% to 90%
- Atmospheric Pressure: 700hPa to 1060hPa

**Transport and storage conditions:**

- Temperature: Above 0° C freezing to +60°C
- Relative Humidity: max 95%
- Atmospheric Pressure: 700hPa to 1060hPa

## DIMENSIONS & WEIGHTS

	Dimensions (WxDXH)	Weight
<b>Vacuum Module</b>	26,39 cm x 13,92 cm x 28,12 cm	Vacuum weight is 2.22 Kg
<b>Modules, ultrasound, 2-channel stim, 2-channel Stim with EMG, Laser</b>	28,24 cm x 16,10 cm x 3,63 cm	2.45 Kg
<b>Head @ 45 degree with Base (Tabletop)</b>	40,36 cm x 40,36 cm x 56,01 cm	Tabletop weight is 9.38 Kg
<b>Cart Lowered (with casters)</b>	60.80 cm x 66.52 cm x 69.62 cm	13.33 Kg
<b>Cart Raised (with casters)</b>	60.80 cm x 66.52 cm x 76.58 cm	13.33 Kg
<b>Head and raised cart with screen @ 90deg</b>	60.80 cm x 66.52 cm x 134.23 cm	22.18 Kg

# Custom Made Modality

Build the precise combination for your needs with five slide-in, plug-and-play module options; Channel 1/2 stim, Channel 1/2 stim/ EMG, Channel 3/4 stim, Laser and Ultrasound. Plus the option of a vacuum electrode module. Each unit is assembled and shipped specifically according to your customized clinical needs.



2 channel tabletop stim with SEMG	2 channel tabletop combo with SEMG	2 channel combo with cart	2 channel combo with sEMG and cart	2 channel combo with sEMG, laser and cart	4 channel combo with cart	4 channel combo with sEMG and cart	4 channel combo with sEMG, laser and cart
<b>Required Part Numbers</b>							
6001 70004	6001 70002 70004	6001 70000 70001 70002	6001 70002 70001 70004	6001 70002 70001 70004 70005	6001 70000 70001 70002 70003	6001 70001 70002 70003 70004	6001 70001 70002 70003 70004 70005

# Ordering information

## Intellect Neo Base Unit + Optional Cart

Part Number	Description
<b>Intellect® Neo</b>	
6001	Intellect Neo Base Unit
<b>Standard Accessories</b>	
13-7651	Intellect Neo User Manual
13-7652	Intellect Neo User Manual on CD
<b>Optional Accessories</b>	
70001	Intellect Neo Therapy System Cart



## Intellect Neo Ultrasound Module

Part Number	Description
<b>Intellect® Neo Ultrasound Module</b>	
70002	Intellect Neo Ultrasound Module
<b>Standard Accessories</b>	
13-8911	Ultrasound User Manual
13-7718	Modules User Manual on CD
4248	Ultrasound Gel
<b>Optional Accessories</b>	
70001	Intellect Neo Therapy System Cart
70008	Intellect Neo Operator Remote
28900	Small Intellect Neo Applicator
28901	Medium Intellect Neo Applicator
28902	Large Intellect Neo Applicator



## Intellect Neo sEMG + Stim Module

Part Number	Description
<b>Intellect Neo sEMG and Channel 1/2 Stimulation Module</b>	
70004	Intellect Neo sEMG and Channel 1/2 Stimulation Module
<b>Standard Accessories</b>	
13-8905	Stim 1/2 Module User Manual
13-7718	Modules User Manual on CD
10648	Nylatex Wrap
79967	6 x 8 cm carbon electrodes (4x)
79970	6 x 8 cm sponges (4x)
42192	Dura-Stick® Plus 5 cm round electrodes (4x)
114.610EXP	Intravaginal Probe
27321	Channel 1 EMG Leadwire
27322	Channel 2 EMG Leadwire
<b>Optional Accessories</b>	
70001	Intellect Neo Therapy System Cart
70008	Intellect Neo Operator Remote
79977	HiVolt Probe
57007	Microcurrent Probe



## Intellect Neo Stim Module

Part Number	Description
<b>Intellect Neo Channel 1/2 Stimulation Module</b>	
70000	Intellect Neo Channel 1/2 Stimulation Module
<b>Standard Accessories</b>	
13-8905	Stim 1/2 Module User Manual
13-7718	Modules User Manual on CD
10648	Nylatex Wrap
79967	6 x 8 cm carbon electrodes (4x)
79970	6 x 8 cm sponges (4x)
42192	Dura-Stick® Plus 5 cm round electrodes (4x)
57007	Microcurrent Probe
<b>Optional Accessories</b>	
70001	Intellect Neo Therapy System Cart
70008	Intellect Neo Operator Remote
79977	HiVolt Probe
57007	Microcurrent Probe

**Note: No applicator included as standard**



## Intellect Neo Channel 3/4 Stim Module

Part Number	Description
<b>Intellect Neo Channel 3/4 Stimulation Module</b>	
70003	Intellect Neo Channel 3/4 Stimulation Module
<b>Standard Accessories</b>	
13-8893	Stim 3/4 Module User Manual
13-7718	Modules User Manual on CD
10648	Nylatex Wrap
79967	6 x 8 cm carbon electrodes (4x)
79970	6 x 8 cm sponges (4x)
42192	Dura-Stick® Plus 5 cm roundelectrodes (4x)
<b>Optional Accessories</b>	
70001	Intellect Neo Therapy System Cart
70008	Intellect Neo Operator Remote
79977	HiVolt Probe
57007	Microcurrent Probe





## Intelect Neo Vacuum Stimulation Module

Part Number	Description
<b>Intelect Neo Vacuum Module</b>	
70006	Intelect Neo Vacuum
<b>Standard Accessories</b>	
13-8809	Vacuum Module User Manual
13-7652	Intelect Neo User Manual on CD
56852	60 mm Vacuum Electrodes Stainless Steel (8x)
56806	60 mm Vacuum Sponges (8x)
56902	Ch 1/3 Vacuum Lead Hose (Black)
56901	Ch 2/4 Vacuum Lead Hose (Gray)
<b>Optional Accessories</b>	
70030	Leadhose Vacuum Ch 1/2 Kit
70031	Leadhose Vacuum Ch 3/4 Kit
70032	Leadhose Vacuum Ch 1/2 Kit XL
70033	Leadhose Vacuum Ch 3/4 Kit XL
70034	90 mm Vacuum Electrode Kit
70035	60 mm Vacuum Electrode Kit
70036	30 mm Vacuum Electrode Kit
70037	90 mm Vacuum Sponge Kit
70038	60 mm Vacuum Sponge Kit
70039	30 mm Vacuum Sponge Kit
70040	Vacuum Module Electrode/Leadhose Kit
70041	Vacuum Plug Kit



## Intelect Neo Laser Therapy Module

Part Number	Description
<b>Intelect Neo Laser Therapy Module</b>	
70005	Intelect Neo Laser Therapy Module
<b>Standard Accessories</b>	
13-8907	Laser Module User Manual
13-7718	Modules User Manual on CD
70008	Intelect Neo Operator Remote
27525	Laser Protection Glasses
<b>Optional Accessories</b>	
70001	Intelect Neo Therapy System Cart
<b>Laser Applicators</b>	
<b>Singles</b>	
27799	LED Diode 10mW
27803	Laser Diode 40mW
27840	Laser Diode 100mW
27804	Laser Diode 150mW
27841	Laser Diode 200mW
27805	Laser Diode 300mW
<b>9 Diode Cluster</b>	
27810	290mW Total: 5x50mW Lasers, 4x10mW LED
27811	540mW Total: 5x100mW Lasers, 4x10mW LED
27812	1040mW Total: 5x200mW Lasers, 4x10mW LED
<b>13 Diode Cluster</b>	
27813	Total: 3x50mW Lasers, 7x10mW LED, 3x15mW LED
27814	415mW Total: 3x100mW Lasers, 7x10mW LED, 3x15mW LED
27816	715mW Total: 3x200mW Lasers, 7x10mW LED, 3x15mW LED
<b>19 Diode Cluster</b>	
27815	325mW Total: 6x10mW LED, 7x25mW LED, 6x15mW LED
<b>33 Diode Cluster</b>	
27809	565mW Total: 12x10mW LED, 13x25mW LED, 8x15mW LED
27802	690mW Total: 5x50mW Lasers, 12x10mW LED, 8x25mW LED, 8x15mW LED
27807	940mW Total: 5x100mW Lasers, 12x10mW LED, 8x25mW LED, 8x15mW LED
27808	1440mW Total: 5x200mW Lasers, 12x10mW LED, 8x25mW LED, 8x15mW LED

**Note: No applicator included as standard**



## DJO GLOBAL

### AUSTRALIA:

**T:** +1300 66 77 30  
**F:** +1300 66 77 40  
**E:** service@djortho.com.au

### BENELUX:

**T:** Belgium 0800 18 246  
**T:** Netherlands 0800 0229442  
**T:** Luxemburg 8002 27 42  
**E:** benelux.orders@DJOglobal.com

### CANADA:

**T:** +1 1866 866 5031  
**F:** +1 1866 866 5032  
**E:** canada.orders@DJOglobal.com

### CHINA:

**T:** (8621) 6031 9989  
**F:** (8621) 6031 9709  
**E:** information\_china@DJOglobal.com

### DENMARK, FINLAND, NORWAY & SWEDEN:

**T:** Sweden 040 39 40 00  
**T:** Norway 8006 1052  
**T:** Finland 0800 114 582  
**T:** Denmark +46 40 39 40 00  
**E:** info.nordic@DJOglobal.com

### FRANCE:

**T:** +33 (0)5 59 52 86 90  
**F:** +33 (0)5 59 52 86 91  
**E:** sce.cial@DJOglobal.com

### GERMANY:

**T:** +49 761 4566 01  
**F:** +49 761 456655 01  
**E:** infoservice@DJOglobal.com

### ITALY:

**T:** +39 02 484 63386  
**F:** +39 02 484 09217  
**E:** it.info@DJOglobal.com

### INDIA:

**T:** +91 44 66915127  
**E:** customercare.india@DJOglobal.com

### SOUTH AFRICA:

**T:** +27 (0) 87 3102480  
**F:** +27 (0) 86 6098891  
**E:** info.southafrica@DJOglobal.com

### SPAIN:

**T:** +34 934 803 202  
**F:** +34 934 733 667  
**E:** ventas@DJOglobal.com

### SWITZERLAND:

**T:** +41 (0) 21 695 2360  
**F:** +41 (0) 21 695 2361  
**E:** info@compex.ch

### UK & IRELAND:

**T:** +44 (0)1483 459 659  
**F:** +44 (0)1483 459 470  
**E:** ukorders@DJOglobal.com

### UNITED STATES:

**T:** +1 800 336 6569  
**F:** +1 800 936 6569  
**E:** customercare@DJOglobal.com

## DJO GLOBAL, EXPORT CENTRES

### ASIA-PACIFIC:

DJO Asia-Pacific Limited  
Unit 1905, 19/F, Tower II  
Grand Central Plaza  
138 Shatin Rural Committee Road  
Shatin  
HONG KONG  
**T:** +852 3105 2237  
**F:** +852 3105 1444  
**E:** info.asia@DJOglobal.com

### EUROPE, MIDDLE EAST & AFRICA:

DJO Benelux  
Welvaartstraat 8  
2200 Herentals  
BELGIUM  
**T:** +32 (0) 14248350  
**F:** +32 (0) 14248358  
**E:** info.emea@DJOglobal.com

### LATIN AMERICA:

DJO Global, Inc  
1430 Decision Street  
Vista  
CA 92081-8553  
U.S.A.  
**T:** 1 800 336 6569  
**F:** 1 800 936 6569  
**E:** info.latam@DJOglobal.com



*Together in Motion™*