# Electrophoresis Power Supply EPS 3500



**User Manual** 



#### Important user information

## Reading this entire manual is necessary for full understanding and safe use of this product.



The exclamation mark within an equilateral triangle is intended to alert the user to the presence of important operating and maintenance instructions in the literature accompanying the instrument.



The lightning symbol within an equilateral triangle is intended to alert the user to the risk of exposure to high voltages.

Should you have any comments on this manual, we will be pleased to receive them at:

Pharmacia Biotech AB S-751 82 Uppsala Sweden

Pharmacia Biotech reserves the right to make changes in the specifications without prior notice.

#### Warranty and Liability

Pharmacia Biotech AB guarantees that the product delivered has been thoroughly tested to ensure that it meets its published specifications. The warranty included in the conditions of delivery is valid only if the product has been installed and used according to the instructions supplied by Pharmacia Biotech AB.

Pharmacia Biotech AB shall in no event be liable for incidental or consequential damages, including without limitation, lost profits, loss of income, loss of business opportunities, loss of use and other related exposures, however caused, arising from the faulty and incorrect use of the product.

### Copyright© 1994 Pharmacia Biotech AB

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form by any means, without permission in written form from the company.

### Contents

1.	Intro	ductionduction	3
2.	Safet	y information	4
	2.1	•	
	2.2	In-built safety features	4
3.	Unpa	acking and installation	
	Technical description		
	4.1	Front panel	6
		4.1.1 Display	6
		4.1.2 Keyboard	6
		4.1.3 Output sockets	8
	4.2	Rear panel	8
5.	Oper	ation	9
	5.1	Overview	9
	5.2	Programming	9
	5.3	Running a program	11
	5.4	Optional programming	13
	5.5	Choosing run parameters	14
	5.6	Short instructions	16
6.	Main	tenance	17
7.	Troul	ole shooting	18
8.	Techi	nical specifications	19
9.	Orde	ring information	20

### 1. Introduction

The Pharmacia Electrophoresis Power Supply EPS 3500 is a high quality, high precision and safe power supply for electrophoresis applications that require a maximum of 3500 volts.

EPS 3500 is primarily designed for the following techniques:

- DNA sequencing
- IEF (Isoelectric focusing)

EPS 3500 is also suitable for:

- SDS-PAGE (Polyacrylamide Gel Electrophoresis)
- Native PAGE
- Agarose electrophoresis
- Electroblotting
- DNA pulsed field electrophoresis

Electrophoresis separations can be controlled by voltage, current or power. The EPS 3500 automatically switches over the controlling parameter according to programmed limits and conductivity variations in the system.

Two electrophoresis units can be connected to the EPS 3500 and run with the same programmed method at one time. Three programs can be saved.

### 2. Safety information

### 2.1 Safety precautions





Extreme caution should be exercised in the operation of this instrument as it can develop sufficient voltage and current to produce a lethal shock.

To avoid any risk of injury, the instrument should only be operated by properly trained personnel and always in accordance with the instructions provided.

Read this entire manual before using this power supply.

- 1. This instrument is designed for indoor use only.
- 2. The instrument must always be used with the protective earth lead of the power cord correctly grounded to earth at the mains outlet.
- 3. To permit sufficient cooling, ensure that the vents in the rear and sides of the instrument are not covered.
- 4. Do not operate the instrument in extreme humidity (above 95%). Avoid condensation by letting the unit equilibrate to ambient temperature when taking the power supply from a colder to a warmer environment.
- 5. Keep the instrument as dry and clean as possible. Wipe regularly with a soft damp cloth. Let the power supply dry completely before use. If wetted, unplug the power supply until the instrument is dry.
- 6. Use only undamaged electrical wire and equipment specified for the voltages you will use. High voltage wires must meet the requirements of the IEC 1010-2-031: 1993 electrical standard.

Any electrophoresis equipment connected to the power supply should meet the requirements of the IEC 1010-1:1993 electrical standard.

7. Note that the output is connected to the chassis/reference earth.

The EPS 3500 is designed in accordance with the IEC 1010 electrical safety standard. The power supply also has several built-in safety functions. Error messages are also shown on the display:

1. Functional earth leakage.

Should the power supply be connected to an electrophoresis unit that has a leakage path to earth, the EPS 3500 will detect this fault and the high voltage is turned off.

2. Start current check.

To ensure that an electrophoresis unit is connected correctly, the power supply checks that the resistance is not higher than a specified limit at a low safety voltage (<40V). If this resistance is too high, the voltage is turned off. Too high a resistance can also be caused by using buffers with extremely low conductivity. The high voltage is also turned off in this case and an error message is shown. This function can be disabled to perform certain applications (see 5.4 Optional programming).

3. Sudden load change detection.

This function prevents accidents under running conditions due to a break in the electrical circuit such as a bad connection to the electrophoresis unit. The high voltage is turned off in such an event.

### 2.2 In-built safety features

# Unpacking and installation

#### Unpacking

Check the contents against the packing list supplied. Inspect for any damage that may have occurred during transit. Report any damage immediately to your local Pharmacia representative and to the transport company concerned.

#### **Mains connection**



Select the appropriate voltage range, 100-120 or 220-240 V, see Fig. 2.

Warning! If the power supply is connected to 220-240 V with the range set to

100-120 V, the instrument can be severely damaged.

Select the appropriate mains cable and connect one end to the mains socket on the EPS 3500 power supply, see Fig. 2, and the other end to an AC grounded outlet.

### Local regulation for Great Britain WARNING IMPORTANT

This appliance must be earthed.

The wires in the mains lead are coloured in accordance with the following code:

Green and yellow Blue Neutral Brown Live

If the plug provided is unsuitable for socket outlets, the plug must be cut off and suitable plug fitted. The cut-off plug should be disposed of and must not be inserted into any 13 amp socket as this can result in electric shock. The plug or adapter or the distribution panel should be provided with 13 amps fuse. As the colours of the wires in the mains lead of this appliance may not correspond with coloured markings identifying the terminals in your plug, proceed as follows:

The green and yellow wire must be connected to the terminal in the plug which is marked with the letter E or by the earth symbol, or coloured green, or green and yellow.

The blue wire must be connected to the terminal which is marked with the letter N or coloured black.

The brown wire must be connected to the terminal which is marked with the letter L or coloured red.

Note:

After replacing or changing a fuse, the fuse cover in the plug must be replaced with a fuse cover which corresponds to the colour of the insert in the base of the plug or the word that is embossed on the base of the plug, and the appliance must not be used without a fuse cover.

Only 13 Amps fuse approved to B.S. 1362 A.S.T.A. should be used.

Switch on the power. Each time the instrument is turned on a self diagnostic test is done. If an error is detected during the test a message will appear on the display and an alarm will sound.

#### Connection of the electrophoresis unit(s)

Connect the leads from the electrophoresis unit (red to red, and black or blue to blue), see Fig. 1. The red lead is the positive and black or blue is the negative.





Warning! Use only undamaged electrical wire and equipment approved for the voltage you will use.

Two electrophoresis units can be run simultaneously with the same program. Please remember to double the limiting current and power if two electrophoresis units are run at the same time. The voltage will be the same regardless of the number of units.

### 4. Technical description

#### 4.1 Front panel

The front panel consists of an alphanumeric display, a keyboard with 9 membrane keys, a light emitting diode (LED) that lights when voltage is applied (HV on) and connectors for two electrophoresis units.

#### 4.1.1 Display

A 32 digit alphanumeric display guides you through the programming, shows current parameter values during the electrophoresis and final parameter values afterwards. It also asks questions and shows error messages. The display has an upper and lower row.

Fig. 1 shows the display in the start position when power is switched on for the first time. The mode (in this case SET) is shown in the upper row on the left. The program number is in the center of the upper row. The number shown is that of the previously entered program. As no program has yet been entered, the first available program number, P1, is shown. The blinking figure, in this case "1", indicates that it can be changed by using the

keys. The upper right shows the programmed breakpoint. In this case it is off.

The lower left, middle and right positions show voltage, current and power respectively.




Fig. 1. The front panel of the EPS 3500.

#### 4.1.2 Keyboard

enter

#### Enter

Enters a value, locks the value if correct, and moves programming to the next field. Valid values are voltage 35-3500 V, current 1-150 mA, power 1-100 W, time 0:01-500 h, volthours 1-500 000 Vh.

set	Set
	Puts the instrument into SET, its programming and starting mode. In the
	RUN mode, pressing also shows the programmed parameters for
	the actual run. In addition, allows you to make changes in the
	program during a run after first pressing continue.
	Change up/Change down Changes the parameter, value or other choice in the field which is blinking. Numerical values are changed in an accelerating manner when a key is held down. Clicking a key changes the value in preset increments. Parameters or units (e.g. Vh) and YES/NO choices are changed with one key push. The
	keys can also be used to switch between time and volthours in RUN, PAUSE and END. The values scroll i.e. they automatically change from maximum to minimum value or vice versa.
[run	Run
run	Pressing run starts the run and puts the program into RUN mode. The
	program number, current values for voltage, current and power are shown on the display. The elapsed time or volthours are also displayed. Switch
	between time and volthours with .
	Pause
continue	Puts the instrument in PAUSE mode and switches off the voltage. The display shows the status of the run at the time the key was pressed.
	only operates in RUN mode. Time and integration values are retained.
	In the pause mode, set can be used to make changes in the program.
	Return to RUN mode by pressing pause or by pressing or by pressing .
STOP	Stop
	Stops the run and puts the instrument in END mode. The voltage is switched off and the end parameters are displayed. Switch between time and
	integrated voltage by pressing $ \bigcirc $ . A run cannot be continued after pressing $ \bigcirc $ .
	Continue by pressing run to run the same method again or pressing set
	to choose another program, program a new method or make changes in an
	existing method.

more

#### More

You can disable the start current check that otherwise detects if the resistance is too high and switches off the power. See Section 5.4.

exit

#### **Exit**

Stops the execution of an operation, such as the entry of a value. Only values/units that have already been confirmed by enter are retained.

Returns the instrument to the mode that was left or the start position in SET.  $\,$ 

#### 4.1.3 Output sockets

There are two sets of output sockets to allow two electrophoresis units to be connected and run at the same time, see Fig. 1. The voltage output is 0-3500 V. The negative output socket gives between 0 and -1750 V and the positive gives between 0 and +1750 V.

#### 4.2 Rear panel

The rear panel is shown in Fig. 2. On the rear panel there is:

- 1. A mains switch. Press in I to switch on the power to the power supply. Press **0** to switch off the power.
- 2. A socket for the mains cable.
- 3. A switch for voltage range. The left position corresponds to 100-120 V and the right to 220-240 V.
- 4. Fan vents.




Fig. 2. The rear panel of the EPS 3500

### 5. Operation

#### 5.1 Overview

The operation of the EPS 3500 is described in the following sequence.

- 1. Programming or editing a method.
- 2. Running a method.

When the method is running, it can be interrupted (put in PAUSE mode) to apply samples and/or change the program. During a run it is also possible to look at the program settings.

Programming and running are discussed in more detail in the following two sections.

Blinking characters are shown as **bold** characters.

#### 5.2 Programming

#### Start position

When the power supply is switched on, the display shows the start position in SET mode. The previous program set is shown and that program number is blinking, for example:

SET:	P <b>1</b>	OFF
3000 V	50mA	20 W
3000 V	AIIIUC	20 VV

(If the power supply is switched on for the first time, the breakpoint parameter and the alarm are off and the values for voltage, current, and power are 35 V, 150 mA and 100 W respectively).

Up to 3 programs can be entered. Press   to confirm the program
shown on the display or use keys to choose another.  Confirm with enter.  The programming procedure can be followed step by step in Fig. 3.
Setting voltage, current and power limits  The display will now flash for the set voltage limit. Using the
keys, select the maximum voltage desired for the run. Confirm
with $\begin{tabular}{ l l l l l l l l l l l l l l l l l l l$
power. Programmable values for voltage are 35-3500, current, 1-150 mA, power, 1-100 W, time, 0:01-500 h, volthours, 1-500 000 Vh.

#### Setting break-point

Choosing a program

Choose between manual break or automatic break. For automatic break, choose breakpoint units in either Volthours (Vh) or time (h). Select the correct unit or, for manual break, choose, "OFF" with




Fig. 3. Step-by-step summary of programming.

Confirm with enter.
If volthours or h (time) were chosen, the display will flash for the break
value. Set the value with and confirm with enter.
Choosing alarm When a breakpoint has been entered, the following question is shown:
Alarm at End? Y
Select YES or NO with and confirm with enter. If YES is
selected, a small "A" appears on the right of the program number, i.e. 2 <sup>A</sup> .
Back to start position
The program is now back to the start position in the SET mode with the program number flashing. It is possible to go back to this position at any
stage during programming in SET mode by pressing exit or Stop . Note
that the program is saved automatically when exit is pressed.
Disabling the current check See Section 5.4 if you want to use this feature.
See Section 3.1 if you want to use this feature.
Connect the leads from the electrophoresis unit (red to red, and black or blue to blue). Red is positive and black or blue negative. Up to two electrophoresis units can run at the same voltage at one time. Please remember to double the maximum current and power conditions if two units are to be run. Voltage will be the same regardless of how many units are run. The current should also be doubled if two gels are run on the same unit.
Choosing a program
Press set and select the program you wish to run by pressing
until the value after the "P" is correct. Confirm with enter.
(Omit this step if you have just programmed a method as described in section 5.2.)
Running a program is described schematically in Fig. 4.
Running
Press run to start the electrophoresis. The display will show current
values for voltage, current, power and elapsed time or volthours. You can
switch between showing the elapsed time or volthours by \[ \bigcup \].
The limiting parameter will attain the programmed limit. The unit
for that parameter is underlined. A light emitting diode shows when voltage is applied (HV on).

5.3 Running a program




Fig. 4. Step-by-step summary of run and pause function

If no current is displayed, please check the electrical connections to the electrophoresis equipment.

#### **Pausing**

You can interrupt the electrophoresis for sample loading and/or changing parameter values by pressing parameter values by pressed. The display shows the status of the run when pause pause pause pause pause pause pause pause again or pause pause again or pause pause again or pause pause again or p

entered you can make changes as described under programming. The P-SET mode is the same as the SET mode apart from restrictions in setting the breakpoint. Naturally, it is not possible to enter a time or volthour that is already passed.
Press exit or pause continue to go back from P-SET to PAUSE. Press or pause continue
to proceed with the electrophoresis.
View programmed values
It is also possible to view the programmed values during a run by pressing
set . Note that no values can be changed here. The display returns
automatically to show RUN values after 5s. Alternatively use exit or run.
Stopping the run and viewing end parameter values When the programmed time or volthour is attained, the program will enter
the END mode. It is also possible to break the run manually by pressing $^{\text{STOP}}$ .
In both cases, the voltage, current and power will go to zero as indicated by the HV on LED switching off. The end parameter values are displayed.
Switch between elapsed time and integrated voltage by $\left[\begin{array}{c c} \end{array}\right] \left[\begin{array}{c c} \end{array}\right]$ .
An alarm will sound if selected in the program. You can stop the alarm by pressing stop.
A run cannot be continued after pressing stop.
Disconnect the leads and proceed with post-electrophoretic techniques.
Since diffusion will begin as soon as the voltage is turned off, you should remove the gel and begin staining, blotting or autoradiography immediately.
You can disable the start current check that otherwise checks that the resistance is not higher than a specified limit. This disabling function is set in the MORE mode instead of the SET mode in order not to confuse users who do not need this feature, and for safety reasons.
The MORE mode is described in Fig. 5. After pressing , the program
asks you for the program number and if you want to keep the start current check.
Choose program number and change between "YES "and "NO" by
and confirm with center . Choosing "NO" for start current check
means that the function is disabled and high voltage can be applied despite a
very high resistance. Leave the MORE mode by exit.
Warning! By disabling the start current check, the power supply can
deliverhigh voltage even if it is not connected to
electrophoresis equipment.

5.4 Optional

programming

13




Fig. 5. Disabling the current check.

### 5.5 Choosing run parameters

EPS 3500 is an automatic cross-over power supply that allows the user to set limits for voltage, current and power. During electrophoresis, only one of the parameters is limiting at a time. The limiting parameter determines, together with the conductivity in the electrophoresis system, the values for the other two parameters. Voltage, current, power and conductivity are related by the following equations:

$$U = I / L \qquad (1)$$

$$P = U \times I \qquad (2)$$

Where U = Voltage, I = Current, P = Power and L = Conductivity

Equation (1) is more familiar if the conductivity is replaced by the reciprocal resistance (R):

$$U = R \times I$$
 (Ohm's law)

The electric field E, measured in V/cm, is the driving force behind electrophoresis.

$$E = U/d$$

where E = Electrical field strength, U = Voltage, d = distance

The electrical field strength is achieved by applying a voltage. The higher the voltage, the faster the electrophoresis. Fast electrophoresis is beneficial since it counteracts diffusion.

The disadvantage of increasing the voltage too much is that most of the generated electrical energy, the product of power and time, is transformed to heat. Therefore cooling of electrophoresis equipment is recommended.

Cooling will also reduce "smiling" effects which are caused by mobility differentials across an electrophoresis gel resulting from poor heat transfer. Since the cooling efficiency cannot be increased indefinitely, the power should be limited when programming the power supply.

The parameter that should be chosen as the constant and thus control the electrophoresis depends on the type of electrophoresis. In the case of homogeneous buffers throughout the system (same electrode and gel buffer), the conductivity is constant during the electrophoresis. If the conductivity is constant, the voltage will be proportional to the current and the power to the square of the current, according to (1) and (2). This means that the result of the electrophoresis will be the same, regardless which parameter is chosen for constant. For historical and practical reasons, voltage is most commonly used for regulation. Submarine gel electrophoresis of DNA/RNA and pulsed field electrophoresis are usually run at constant voltage. SDS-PAGE using continuous buffer systems is run at constant voltage or current.

For discontinuous buffer systems, the resistance will increase as the electrophoresis proceeds due to a decrease in conductivity. Running at constant voltage will result in decreasing current and power. Constant voltage will thus be "safe" in the respect that the power will not increase and produce more and more heat. On the other hand, the separation will slow down and impair resolution due to an increased time available for diffusion. Running at constant power would give a faster electrophoresis and controlled power, while running at constant current would, at the first sight, seem to be problematic because of increasing voltage and power.

During discontinuous electrophoresis, however, the voltage is not distributed evenly across the gel. These gels have a region with low ionic strength that causes a high electrical field strength. This region increases as the electrophoresis proceeds. This means that the main part of the voltage is spread over a longer and longer distance and a higher and higher power is tolerated. This is the reason why constant current is chosen for such applications. It is, however, recommended to also limit the power as a precaution against overheating the gel. The power supply will probably switch over to limiting power at the end of the run due to increased voltage.

The crossing-over between different parameters controlling the electrophoresis can be illustrated by IEF (isoelectric focusing) using carrier ampholytes. A graphical representation of the changes in power, voltage and current that may occur during a typical IEF run is given in Fig. 6. Since the pK values of the carrier ampholytes and the proteins are temperature dependent, IEF must be carried out at a constant temperature. Therefore cooling of electrophoresis equipment and controlling by power is recommended. The main part of the IEF is thus controlled by power (phase II). The conductivity is gradually decreasing because the carrier ampholytes and sample will loose their net charge during the build up of the pH gradient. Thus the current will decrease and the voltage increase at constant power. During the early stage of the formation of the pH gradient it is important to limit the current. Otherwise the gradient will be irregularly shaped (phase I). The conductivity is not constant along the gel so it is important to also limit the voltage. This means running at constant voltage for the last phase (phase III) to prevent local overheating.

For detailed information about parameter values, always follow the gel supplier's recommendations.





Fig. 6. Power fluctuations which occur during a typical IEF run.

### 5.6 Short instructions

This section summarizes the main points covered earlier in this chapter. Use it as a check list once you are familiar with the detailed programming and running instructions. Refer also to the separate schematic operating guide included with the power supply. We recommend you keep this schematic guide close by the instrument.



by a small A as in "1<sup>A</sup>" to the right of the program number in the

display.

- 7. When programming is completed, connect your electrophoresis unit to the outputs, and press run.
- 8. The actual voltage, current and power are shown on the display together with the elapsed time or volthours (switch with ).
- 9. pause can be used to interrupt and then continue the run.
- 10. The programmed parameters can be checked during the run by  $\ensuremath{\text{pressing}}$  pressing  $\ensuremath{\mbox{\tiny [set]}}$  .
- 11. The programmed run stops automatically or by pressing stop

### 6. Maintenance

Check regularly that the fan is working properly. Wipe the instrument regularly with a damp cloth. Let the instrument dry completely before use. Otherwise no user maintenance is necessary. All servicing should be entrusted to qualified personel only. Please contact your local Pharmacia representative for more service information.

### 7. Trouble shooting

If an error that can be corrected by the user occurs, either during a run or when switching on the power supply, the program enters the HALT mode and the output is switched off. Four different types of errors can cause HALT. The following list shows the error message on the display, the cause and the remedy.

#### **Error message**

#### Cause

#### Remedy

HALT: Low start current!

The current is less than the lower limit. This can be due to incorrect connection of the electrophoresis equipment or due to use of buffers with extremely low conductivity.

The current to ground leakage in the

1. Check connections and/or buffers.

Press , switch off the start current check in the MORE mode.

IMPORTANT! Read Section 5.4 first.



HALT: Ground leakage current

The mains voltage is too low s

Check the electrophoresis unit.

HALT: Mains Voltage too low!


The mains voltage is too low, see Specifications.

electrophoresis unit is too high.

Check voltage selector. Check mains voltage.

HALT: Mains fail Program stopped

Mains power failure for more than 7 s.



If a serious error occurs, the program enters the FAIL mode. The output is switched off and an error message is shown in the display.

FAIL:Code No: xxx Call service

Please read the error code number and contact your Pharmacia LKB Biotechnology representative.

### 8. Technical specifications

Regulation Maximum voltage, current and power with automatic cross-over

at preset limits

Output range Voltage: 0-3500 V DC,

Current: 0-150 mA Power: 0-100 W

Programming range Voltage: 35-3500 V DC

Current: 1-150 mA Power: 1-100 W Time: 00:01-500 h Volthour: 1-500 000 Vh

Output resolution Voltage: 1 V

Current: 100 µA, 0-14.9 mA

1 mA, 15-150 mA Power: 1 W

Programming resolution Voltage: 5 V

Current: 1 mA Power: 1 W

Time: 1 min, 00:01-99:59,

1 h, 100-500 h

Volthour: 1 Vh, 1-9999 Vh,

100 Vh, 10.0-99.9 kVh 1 kVh, 100-500 kVh

Accuracy Voltage: 1%, ± 5 V

Current: 1%, ± 1 mA Power: 2%, ± 1 W

Timer:  $0.1\% \pm 1 \text{ min}, 00:01-99:59 \text{ h}$ 

0.1%, ± 10 min, 100-500 h

Line regulation < 0.2%

Load regulation < 1% at load change 10-90% of maximum load

Ripple <1 % at 3500 V

Short term stability < 0.2%/10 h after warm up

Long term stability <1%/year

Start current check Resistance not greater than 1  $M\Omega$  at 40 V

(current less than 40 µA). Can be disabled

Sudden load change check

Ground leakage check Leakage not greater than 500 µA

Output protection Fully protected against any overload conditions

Recovery after power failure Duration <8 s, The program continues automatically

Duration  $\geq$  8s, The program continues after manual restart

Ambient operating temperature 4-40 °C Ambient operating humidity 0-95%

Ambient operating pressure 106-84 kPa (corresponds to 0-2000 m)

Mains requirements 100-120 V/220-240 V; 50/60 Hz

Power consumption Max 140 W

Dimensions (WxDxH) 250 x 385 x 95 mm

Weight 4.7 kg

## 9. Ordering information

Designation	Code No.
Electrophoresis Power Supply - EPS 3500	19-3500-00

