


| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 1 of 83 |

Introduction

ProSolv® CardioVascular is a web-based, software-only, vendor- and hardware neutral solution designed to provide multi-modality cardiovascular image and information management in a versatile and cost effective manner supporting all cardiovascular modalities. From a single workstation, with a single log-on at any location, a physician can review patient studies and create reports within all cardiovascular modalities. ProSolv® CardioVascular accepts images and information from any DICOM compliant acquisition device. In addition, our reporting tools have set the standard for customization and advanced analysis capabilities, with auto-interpretation, interactive diagrams; with auto-summary, and spell-check. Our ability to customize reports that mimic exacting requirements is unsurpassed in the industry. The ProSolv® CardioVascular software has been carefully designed to address each of the major issues of importance to the field of cardiovascular image and information management, including multimodality viewing, advanced analysis, customized workflow automation through HIS, ADT/Orders, DICOM Modality Worklists, Results and Billing integrations. ProSolv CardioVascular offers two different clinical user interfaces:

ProSolv CardioVascular Client is an installed application for Microsoft Windows. Much of the content of this User Manual is devoted to the functions of the Client. The Client gives full capability for image viewing and analysis, and on-line completion and distribution of clinical reports. When viewing images, the computer monitor screen is organized to maximize the available pixels for digital images rather than controls, buttons and icons. The default screen permits two side-by-side, high-resolution viewing windows. This layout highlights the ability of digital imaging to display two or more cardiac images simultaneously, which is one of the principle advantages of digital imaging over videotape or ciné film and has many clinical uses such as viewing serial studies or different views of the same study. To ensure that the highest resolution images are displayed, none of the lateral borders of the monitor are used for anything other than imaging. Almost all manipulation of images can be performed by clicking on toolbar buttons or pressing keyboard keys. The controls are designed to be intuitive, requiring relatively brief instruction.

ProSolv CardioVascular WebViewer is a lightweight, Web-based application that provides for efficient, distributed viewing of images and completed reports. The WebViewer is described more fully in a separate document, *Web Viewer User Manual* (document control: PMAN-0014755-A).

Intended Use

Use

ProSolv CardioVascular provides the functionality for viewing, measuring, analyzing, archiving, and reporting digital studies generated by digital-imaging instruments from a multitude of manufacturers.


Users

ProSolv CardioVascular is intended for use by doctors and non-physicians (e.g., sonographers, technologists, nurses) who are affiliated with medical labs that use images obtained from digital-imaging DICOM instruments (echo, vascular, nuclear cardiology, cardiac catheterization, MRI, etc.)

Administration of Server-based Applications

ProSolv® CardioVascular is a client-server type application with a significant component of the system's capability being provided by a set of server-based applications.

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 2 of 83 |

Software/Hardware

The ProSolv CardioVascular image management and reporting system is software that runs on standard PC equipment using Windows 2000 or later operating systems. The software can be used in a variety of network configurations:

- Stand-alone workstation
- Network of workstations connected through TCP/IP networking, utilizing the ProSolv CardioVascular Core Server components
- Mini-server Workstation providing both client and server applications

For more information, refer to the *Hardware Requirements* document (document control PSUP-0013745-B).

Disclaimer

No part of this manual may be reproduced without prior written permission.

The information contained in this manual may be subject to change without prior notice.

ProSolv CardioVascular, A FUJIFILM Company, shall not be liable for malfunctions and damages resulting from improper installation, relocation, remodeling, maintenance and repair performed by ProSolv CardioVascular, A FUJIFILM Company, and all other FUJIFILM Subsidiaries or Distributors.

ProSolv CardioVascular, A FUJIFILM Company, shall not be liable for malfunctions and damages of ProSolv CardioVascular, A FUJIFILM Company products due to products of other manufacturers not supplied by ProSolv CardioVascular, A FUJIFILM Company.

ProSolv CardioVascular, A FUJIFILM Company, shall not be liable for malfunctions and damages resulting from remodeling, maintenance and repair using repair parts other than those specified by ProSolv CardioVascular, A FUJIFILM Company.


ProSolv CardioVascular, A FUJIFILM Company, shall not be liable for malfunctions and damages resulting from negligence.

ProSolv CardioVascular, A FUJIFILM Company, shall not be liable for malfunctions and damages resulting from natural disasters.

ProSolv CardioVascular, A FUJIFILM Company, may have patents or pending patent applications, trademarks, copyrights, or other intellectual property rights covering various subject matters contained in this document. The furnishing of this document does not extended any rights to said patents, trademarks, copyrights, or other intellectual property rights, except as expressly set forth in any written license agreement from ProSolv CardioVascular, A FUJIFILM Company.

ProSolv CardioVascular is a computer program which is protected by copyright law and international treaties. Unauthorized reproduction or distribution of this program, or any portion of this program, is strictly prohibited.

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 3 of 83 |

Certain images/diagrams contained in this manual may not appear clearly on high resolution monitors.

ProSolv CardioVascular is a registered trademark of ProSolv CardioVascular, A FUJIFILM Company.

©Copyright 2008 ProSolv CardioVascular, A FUJIFILM Company. All rights reserved.

All brand, product and company names are trademarks or registered trademarks of their respective companies or organizations.

Safety Information

DANGER, WARNING and CAUTION Safety Statements

The DANGERS, WARNINGS and CAUTIONS that appear throughout this manual are not only preventative measures designed to uphold the safety of both the service person and operator, but also enhance equipment reliability.

The definitions for DANGER, WARNING and CAUTION comply with ANSI Z535.2, American National Standard for Environmental and Facility Safety Signs and ANSI Z535.4, Product Safety Signs and Labels, issued by the American National Standards Institute.

DANGER indicates an imminently hazardous situation, which if not avoided, will result in death or serious injury.

WARNING indicates an imminently hazardous situation, which if not avoided, could result in death or serious injury.

CAUTION indicates a potentially hazardous situation, which if not avoided, may result in minor to moderate injury. It may also be used to alert against unsafe practices.

Important indicates an emphasized note. It is something you should be particularly aware of; something not readily apparent. Important is typically used to prevent equipment damage.

Note indicates a note. It is something you should be particularly aware of; something not readily apparent. A note is typically used as a suggestion.

Electrical Power


Do not perform service with equipment energized. Always de-energize the equipment when performing service. Even when de-energized, equipment components can retain an electrical charge capable of causing electric shock. Always use caution and common sense when performing service.

When installing equipment electrical wiring must comply with the local electrical codes in your area governing the installation of electronic equipment.

Heavy Objects

Heavy objects are objects that weigh 50 pounds or more, or are objects that are too large or bulky to be

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 4 of 83 |

moved or handled by one service person. Installing, removing or servicing heavy objects may require two or more service personnel. When installing, removing or servicing heavy objects, always use suitable gear such as heavy gloves and eye protection. Always uses the appropriate support machinery when installing, removing or moving heavy objects. Always use caution and common sense when installing, removing or servicing heavy objects.

Safety Devices


Safety devices are protection components, which include, but not limited to, fuses, circuit breakers, interlock switches, equipment panels and covers. Do not disable, alter, modify or remove safety devices when performing service.

Confidential and Proprietary Information

Use, duplication, or distribution of this document without permission from ProSolv CardioVascular, A FUJIFILM Company, is strictly prohibited.

Table of Contents

| | | |
|---------|---|----|
| 1. | Service Manager | 7 |
| 1.1. | Integrations | 8 |
| 1.1.1. | WebLab Integration Configuration | 8 |
| 1.1.2. | Synapse Integration Configuration | 9 |
| 1.1.3. | OpenECG Integration Configuration | 11 |
| 1.1.4. | MacLab Integration Configuration | 22 |
| 1.1.5. | Siemens Sensis Integration Configuration | 25 |
| 1.2. | ProSolv® Server | 28 |
| 1.2.1. | Configuration | 28 |
| 1.3. | License Server | 33 |
| 1.3.1. | Configuration | 33 |
| 1.3.2. | Enterprise Licensing Methodology | 34 |
| 1.4. | Watchdog Service | 35 |
| 1.4.1. | Configuration | 35 |
| 1.5. | DICOM Server Service (DICOM Store SCP) | 37 |
| 1.5.1. | Configuration | 37 |
| 1.6. | DICOM Query/Retrieve SCP Service | 42 |
| 1.6.1. | Configuration | 42 |
| 1.7. | DICOM Store SCU Service | 43 |
| 1.7.1. | Configuration | 44 |
| 1.8. | Image Compression Service | 45 |
| 1.8.1. | Configuration | 45 |
| 1.9. | Archiver Service | 47 |
| 1.9.1. | Configuration | 47 |
| 1.9.2. | Archiver Job Wizard | 48 |
| 1.9.3. | Archive to DICOM PACS | 58 |
| 1.10. | Retriever Service | 60 |
| 1.10.1. | Configuration | 60 |
| 1.10.2. | Retrieve Queues | 62 |
| 1.11. | Measurement Importer Service | 63 |
| 1.11.1. | Configuration | 63 |
| 1.12. | HPDT Importer | 65 |
| 1.13. | HD11XE Importer | 66 |
| 1.14. | Burdick/Marquette/Quinton Treadmill Importers | 67 |
| 1.15. | Mennen Hemodynamic Import Service | 67 |
| 1.15.1. | Configuration | 67 |
| 1.16. | HL7 ADT Service | 68 |
| 1.16.1. | Configuration | 69 |
| 1.17. | HL7 Report Messenger Service | 71 |
| 1.17.1. | Configuration | 71 |
| 2. | DICOM Modality Worklist Application | 74 |
| 2.1. | Configuration | 74 |
| 2.2. | Starting the ProSolv® DICOM MWL Application | 76 |
| 2.3. | Exiting the ProSolv® DICOM MWL Application | 77 |
| 3. | DICOM Structured Reporting Application | 77 |
| 3.1. | Configuration | 77 |

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 6 of 83 |

| | | |
|------|--|----|
| 4. | Witt Hemodynamic Import Application..... | 80 |
| 4.1. | Installation..... | 80 |
| 4.2. | Configuration | 80 |


Preliminary

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

1. Service Manager

Windows service-mode applications that are integrated with ProSolv®'s Service Manager framework can run without a user interface and without having a local user being logged on. They share a common, tab-based user interface for managing their configuration and state.

The Service Manager provides a User Interface (UI) to configure and view the status of the ProSolv® Services. Each of the Services has a separate tab on the main window. If a service is not installed, a tab does not appear for that service. A shortcut to launch the Service Manager is installed on the user's Desktop. The first tab labeled, "All Services", lists the services that are installed along with the version numbers. Each Service tab has the same set of common controls described below.



Service Control Section
View service status
Start/stop/restart the service
Set the Service Log Level

Log Entries Section
View the last entries of the log file.
Clicking the **View Log** button opens the log file in a text editor.

1.1. Integrations

1.1.1. WebLab Integration Configuration

ProSolv.WebLab is a Web Service that allows ProSolv CardioVascular Client to communicate with ProSolv Application Server Service and access patient data via HTTP or HTTPS protocols. ProSolv.WebLab acts as a Client's proxy, using its own credentials (ProSolv Service User account) to access and route patient data from the server to the Client application.

To configure the WebLab integration:

1. Start the ProSolv® Service Manager.
2. Click the "WebLab" tab.

ProSolv Service Manager

DICOM QR Mennen DT Report Messenger Retriever Store SCU
All Services Archiver ProSolv Server License Server DICOM Server Watchdog **WebLab** Synapse ADT

Settings

| | |
|-------------------------|--------------------|
| PCode Server Port | 704 |
| PCode Server | 127.0.0.1 |
| WebLab AE Title | ProSolvWebLab |
| Image Directory | \\Image Path |
| Update Directory | \\Updates Path |
| HL7 Reports Directory | c:\\HL7\\Reports |
| GE Measurement Transfer | c:\\GE_Measurement |
| Fax Server | c:\\FaxServer |

GE Measurement Transfer
This is path of GE Measurement Transfer Directory.

Optional Save Directory Settings

| Optional Directory Name | Optional Directory Path |
|-------------------------|-------------------------|
| | |
| | |
| | |

Add Edit Delete Save

3. Modify settings by clicking in the Settings table and entering new information or browsing to required locations.

Note: The values modified in the GUI are read from and written to an xml file in \ProSolv CardioVascular\Services 4.0\ProSolv.WebLab\Web.Config. Values can be modified by the user in GUI itself. A snapshot of the xml file is provided here:

```

Web.Config
<?xml version="1.0"?>
<configuration>
  <configSections>
    <section name="FileLocations" type="System.Configuration.IgnoreSectionHandler"/>
  </configSections>
  <appSettings>
    <add key="pCodeServer" value="127.0.0.1"/>
    <add key="pCodeServerPort" value="704"/>
    <add key="WebLabTitle" value="ProSolvWebLab"/>
  </appSettings>
  <connectionStrings>
    <add name="ProsolvDB" connectionString="Provider=SQLEDB.1;Password=psci;Persist Security Info=Tr
      providerName="System.Data.OleDb" />
    <!--
    <add name="ProsolvDB"
      connectionString="Provider=Microsoft.Jet.OLEDB.4.0;Data source=c:\Testing\New Folder\db1-Passwor.mdl
      providerName="System.Data.OleDb" />
    --> </connectionStrings>
  <FileLocations>
    <FileLocation name="images" location="\\Image Path"/>
    <FileLocation name="updates2" location="\\Updates Path"/>
    <FileLocation name="HL7Reports" location="c:\HL7Reports" />
    <FileLocation name="gmeasurement_transfer" location="c:\GE_Measurement"/>
    <FileLocation name="Fax_server" location="c:\FaxServer" />
  </FileLocations>
  <system.web>
    <compilation defaultLanguage="c#" debug="true"/>
    <authentication mode="Windows" />
    <httpRuntime maxRequestLength="16384" />
    <customErrors mode="Off"/>
  </system.web>
</configuration>

```


4. Click the **Save** button to accept changes.
5. To update directory paths, click the **Browse** button to browse to a specific directory and select the required folder.
6. Click the **Save** button to accept any changes.

NOTE: The Optional Save Directory Settings should not be modified by the customer.

1.1.2. Synapse Integration Configuration

To configure the Synapse integration:

1. Start the ProSolv® Service Manager.
2. Click the "Synapse" tab. All available settings will be displayed in the Synapse Settings table.



3. Modify Settings by using the **Add**, **Delete** or **Settings** buttons to the right of the table.
4. To add a setting, click the **Add** button to the right of the Synapse Settings table.
5. At the Synapse Settings dialog box, enter the new information.

6. Click the **Save** button to accept the changes.
7. To edit a setting, click the **Settings** button to the right of the Synapse settings table and enter the updated information.
8. Click the **Save** button to accept the changes.
9. To delete a setting, click the **Delete** button to the right of the Synapse Settings table.
10. At the confirmation dialog box, confirm or cancel the operation. Click the **Yes** button to delete all related data from the database and refresh the Optional Save Directory Settings table. Click the **No** button to cancel the change.

1.1.3. OpenECG Integration Configuration

OpenECG is the integration with the Epiphany product. The modalities supported in the OpenECG integration are: Resting ECG and Holter. This integration allows physicians to view studies and modify and sign reports using the OpenECG software from an existing workstation configured with ProSolv Client. Studies captured via OpenECG will appear in a ProSolv study list, and when selected will launch OpenECG web client.


The ProSolv Server will import measurements and PDF report files from the OpenECG server, in fashion consistent with ProSolv report data and measurements, and will be completely searchable using SQL queries.

The ProSolv system will also act as a bi-directional “gateway” for ADT, Orders, Billing, and HL7 Result messages between a Hospital Information Systems and the OpenECG server. ProSolv will forward ADT and Orders messages from HIS to OpenECG and will also forward Billing and Result messages from OpenECG to HIS. This approach will consolidate and improve the system integrations.


All ProSolv users that are going to use the integration must have their ProSolv User Names added to the user list that is tied to the OpenECG server. The User Name is used for authentication between the two systems and hence must always match. Passwords on the other hand are not used for authentication and they don't have to match.

To configure the OpenECG integration:

1. Select Tools, then Preferences from the ProSolv menu.



2. Click on the **Machine Settings** button on the Preferences window.
3. On the Machine Settings Window, check the “Enable OpenECG features” checkbox and enter the URL of the OpenECG Web

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 13 of 83 |

Interface (Note: Typically this is: http://<ECG Management Server Machine Name>).

Enable OpenECG features

Server URL (HTML)

4. A short cut to the OpenECG URL needs to be placed on the desktop of all workstations. The short cut needs to point to the Server URL listed above in step '3'.
5. To configure client components, the Measurement Importer, ADT and Report Messenger Plug-ins must be installed.

ECG Management Module DT Service Configuration

The ECG Management Server will export an XML and a PDF report whenever a study is confirmed within the Web Interface. These files are stored on file shares on the ECG Management Server. The ProSolv Measurement Import Service needs to be configured to import these files into the system. This is done by configuring the ECG Management Module Plug-in that is part of the Measurement Import Service. Please note that in order for the Plug-in to be able to access files on the remote ECG Management Server server, the Measurement Import Service must be running under a Domain Account that has read access to those file share.

Importer Settings

Plugin Type: Epiphany DT

Database Settings

| | |
|------------------------------------|--|
| ProSolv Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security |
| Images Directory | I:\images |

File Importer Settings

| | |
|------------------------|----------------------------------|
| Search Pattern | *.xml |
| Input Directory | \\<OpenECG Server IP>\cs\hl7\xml |
| Monitor Subdirectories | False |

Importer Settings

| | |
|----------------|------------------|
| Machine Name | OpenECG Importer |
| Sleep Interval | 5 |

OpenECG Settings

| | |
|---------------------|-----------------------------|
| OpenECG Server Root | \\<OpenECG Server IP>\files |
| File Age Limit | 30 |
| Output Folder | I:\Epiphany\FromOpenECG |

Buttons: Save, Cancel

| Name | Description | Recommended Values |
|------------------------------------|---|--|
| Database Settings | | |
| ProSolv Database Connection String | Connection string to the ProSolv® database. Click the “...” button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=psci;Persist Security Info=True;UserID=prosolv;Initial Catalog=Prosolv;DataSource=localhost |
| Image Directory | Path to the Images Folder. | \\<ImageStore>\psdata\$ |
| File Importer Settings | | |
| Search Pattern | File pattern that the importer is going to be watching for. | *.xml |
| Input Directory | UNC path to the CardioServer export folder. | \\<CardioServer Server IP>\cs\hl7\xml |
| Monitor Subdirectories | An option to enable monitoring of subdirectories. The default settings of the CardioServer Server does not allow for exporting files to subdirectories. | False |
| Importer Settings | | |
| Machine Name | A unique name that describes the importer task. | OpenECG Importer |
| Sleep Interval | Number of seconds between polls of the Input Directory. | 5 |


THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| Name | Description | Recommended Values |
|-------------------------|---|---|
| OpenECG Settings | | |
| OpenECG Server Root | UNC path to the root folder of the CardioServer server. | \\<CardioServer Server IP>\files |
| File Age Limit | Number of days to keep incoming files before deleting them. | 30 |
| Output Folder | Path to the folder where the input files are stored after they have been processed. | \\<AppServer>\ProSolv Database\Epiphany\From CardioServer |

Note: Use recommended values unless Site requests otherwise.

HL7 Message Routing Configuration

All HL7 messages are routed through the ProSolv Server to/from the CardioServer. The ProSolv Report Messenger and ProSolv ADT Services are used in tandem to perform the message forwarding.



HL7 Message Routing Configuration

The ADT service needs to be configured to receive HL7 messages from the HIS on some predefined port number. This task will be responsible


for receiving all incoming HL7 from the HIS. These messages will be processed and entered into the ProSolv Database. The messages will also be saved to a folder so that they can be forwarded to the OpenECG server.

Create a directory: I:\HL7ADT

| ADT Delimiter Settings | |
|--|---|
| HL7 Delimiters | ^~\& |
| ADT Settings | |
| Create Study | False |
| Images Directory | I:\Images |
| Listening Port | 3333 |
| Max Stream File Size | 32 |
| Message File Age Limit | 30 |
| Message Store Directory | I:\HL7ADT |
| Patient ID as Main ID | False |
| Process Message | True |
| Send ACK | True |
| Store HL7 Messages | False |
| Study UID Prefix | 1.2.826.0.1.3680043.2.1331.3175947431. |
| Task Name | HL7 Input Listener |
| Database Settings | |
| ADT Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info |
| DICOM Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info |
| ProSolv Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info |
| Task Name Unique name to describe this ADT task. | |

| Name | Description | Recommended Values |
|-------------------------------|--|-------------------------|
| ADT Delimiter Settings | | |
| HL7 Delimiters | HL7 Delimiters can be added as needed | ^~\& |
| ADT Settings | | |
| Create Study | Changing this option to TRUE will result in studies getting created in the system for each order | False |
| Images Directory | Path to the images folder | \\<ImageStore>\psdata\$ |
| Listening Port | Port that the service is listening on | 3333 |
| Max Stream File Size | Number of MB to set for the maximum stream file size. | 32 |
| Message File Age Limit | Number of days to keep incoming files before deleting them. | 30 |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 17 of 83 |

| Name | Description | Recommended Values |
|------------------------------------|--|---|
| ADT Delimiter Settings | | |
| Message Store Directory | Path to the location of the incoming HL7 messages. | \\<AppServer>\ProSolv database\HL7ADT |
| Patient ID as Main ID | Boolean for setting Patient ID as the main identifier. | False |
| Process Message | Boolean for setting whether the message needs to be processed. Setting this to False results in just storing messages. | True |
| Send ACK | Boolean for turning on ACK messages, this might be required by the HIS system. | True |
| Store HL7 Messages | Boolean for storing HL7 messages. Setting this to False will not create individual message files. | True |
| Study UID Prefix | Study UID Prefix | 1.2.826.0.1.3680043.2.1331.3175947431. |
| Task Name | Desired task name | HL7 Input Listener |
| Database Settings | | |
| ADT Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |
| DICOM Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |
| ProSolv Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |

Note: Use recommended values unless Site requests otherwise.

Report Messenger Service Configuration

The Report Messenger service needs to be configured to forward HL7 messages over to the Open ECG Server. The service will monitor the output folder of the ADT service for HL7 messages. As it finds HL7 messages it will forward to the Open ECG Server that is preconfigured to list on a certain port number.

IP address or Machine Name of the OpenECG Server

Port that the OpenECG server is listening on.

Configure according to incoming HL7 messages (Sending facility / application)

| Name | Description | Recommended Values |
|-------------------|---|---------------------------------------|
| Machine Name | Desired machine name | CardioServer Incoming Forwarder |
| Seconds to Sleep | Number of seconds between polls of the Input Directory. | 5 |
| Watch Directory | Directory to watch for input files | \\<AppServer>\ProSolv Database\HL7ADT |
| Search Pattern | Pattern to search on | *.hl7 |
| Recurse Directory | Option to check subdirectories | Unchecked |
| IP Add.\Host Name | Host IP | <IP OpenECG Server> |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| Name | Description | Recommended Values |
|-----------------|--|--------------------------------|
| Watch Port | Host Port | <Port OpenECG is Listening on> |
| Wait Time | Number of milliseconds between polls of the Input Directory. | 10000 |
| Num Wait Cycles | Number of tries | 3 |
| Persistent | Checking for persistent connections | Unchecked |
| Remove LF | Removes Line Feed from the message | Unchecked |
| Seg Delim | Segment delimiter | 0A |
| Delimiters | Standard HL7 delimiters | MSH ^~\& |
| Send By | Sending facility / application | PROSOLV ;PSV |
| Return ACK | Option to enable return of ACK message | Checked |

Note: Use recommended values unless Site requests otherwise.

ADT Listening Task Configuration

An ADT Listening Task needs to be configured that will receive a message stream from the OpenECG Server on some predefined port number. These messages will not be processed and stored in the ProSolv Database. To accomplish this, the "Process Message" flag should be set to FALSE for this task. The messages will be saved to a folder so that they can be forwarded to the HIS. The ADT Task should be configured to store messages to a folder on the ProSolv Server.


Create a directory: I:\Epiphany\FromCardioServer.

ADT Task Settings

- ADT Delimiter Settings**
 - HL7 Delimiters | ^~\&|
- ADT Settings**
 - Create Study | False
 - Images Directory | I:\Images
 - Listening Port | 3333
 - Max Stream File Size | 32
 - Message File Age Limit | 30
 - Message Store Directory | I:\Epiphany\FromCardioServer
 - Patient ID as Main ID | False
 - Process Message | True
 - Send ACK | True
 - Store HL7 Messages | False
 - Study UID Prefix | 1.2.826.0.1.3680043.2.1331.3175947431.
 - Task Name | HL7 Output Listener
- Database Settings**
 - ADT Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info
 - DICOM Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info
 - ProSolv Database Connection String | Provider=SQLOLEDB.1;Password=psci;Persist Security Info

Task Name
Unique name to describe this ADT task.

Save Cancel

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 20 of 83 |

| Name | Description | Recommended Values |
|------------------------------------|--|---|
| ADT Delimiter Settings | | |
| HL7 Delimiters | HL7 Delimiters can be added as needed. | ^~\& |
| ADT Settings | | |
| Create Study | Changing this option to TRUE will result in studies getting created in the system for each order | False |
| Images Directory | Path to the images folder | \\<ImageStore>\psdata\$ |
| Listening Port | Port that the service is listening on | 3333 |
| Max Stream File Size | Number of MB to set for the maximum stream file size. | 32 |
| Message File Age Limit | Number of days to keep incoming files before deleting them. | 30 |
| Message Store Directory | Path to the location of the incoming HL7 messages. | \\<AppServer>\ProSolv Database\Epiphany\From CardioServer |
| Patient ID as Main ID | Boolean for setting Patient ID as the main identifier. | False |
| Process Message | Boolean for setting whether the message needs to be processed. Setting this to False results in just storing messages. | True |
| Send ACK | Boolean for turning on ACK messages, this might be required by the HIS system. | True |
| Store HL7 Messages | Boolean for storing HL7 messages. Setting this to False will not create individual message files. | True |
| Study UID Prefix | Study UID Prefix | 1.2.826.0.1.3680043.2.1331.3175947431. |
| Task Name | Desired task name | HL7 Input Listener |
| Database Settings | | |
| ADT Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |
| DICOM Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |
| ProSolv Database Connection String | Connection string to the ProSolv® database. Click the "... " button to bring up the Connection editor. | Provider=SQLOLEDB.1;Password=p sci;Persist Security Info=True;User ID=prosolv;Initial Catalog=Prosolv;Data Source=<sqldb> |

Note: Use recommended values unless Site requests otherwise.

Report Messenger Forward Task Configuration

A Report Messenger Forward Task needs to be configured to monitor the folder that the ADT Task stores messages. The Report Messenger

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

task will pull in the messages and send them to the HIS on a predefined port number.


The screenshot shows the 'Report Messenger Settings' dialog box with the following fields and callouts:

- Machine Name:** CardioServer Outgoing Forwarder
- Seconds To Sleep:** 5
- Directory Information:**
 - Watch Directory:** I:\Epiphany\From CardioServer
 - Search Pattern:** *.HL7
 - Recurse Directory:**
- Socket Info:**
 - IP Add.\Host Name:** HIS IP (Callout: IP address or Machine Name of the OpenECG Server)
 - Watch Port:** HIS Port (Callout: Port that the OpenECG server is listening on.)
 - Wait Time (ms):** 10000
 - Num Wait Cycles:** 3
 - Persistent:**
- File Processing:**
 - Remove LF:**
 - Seg Delim:** 0A
 - Delimiters:** MSH|^~\&|
 - Send By:** * XYZIABC (Callout: Configure according to incoming HL7 messages (Sending facility / application))
 - Return ACK:**

At the bottom, there is a note: "Separate with ','" and buttons for "Save" and "Cancel".

| Name | Description | Recommended Values |
|------------------|---|-------------------------------|
| Machine Name | Desired machine name | CardioServer Forwarder to HIS |
| Seconds to Sleep | Number of seconds between polls of the Input Directory. | 5 |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 22 of 83 |

| Name | Description | Recommended Values |
|-------------------|--|---|
| Watch Directory | Directory to watch for input files | \\<AppServer>\ProSolv Database\Epiphany\From CardioServer |
| Search Pattern | Pattern to search on | *.hl7 |
| Recurse Directory | Option to check subdirectories | Unchecked |
| IP Add.\Host Name | IP Address of the OpenECG server | <IP of HIS system> |
| Watch Port | OpenECG HL7 Port Number (Default is 2575) | <Port that HIS is listening on> |
| Wait Time | Number of milliseconds between polls of the Input Directory. | 10000 |
| Num Wait Cycles | Number of tries | 3 |
| Persistent | Checking for persistent connections | Unchecked |
| Remove LF | Removes Line Feed from the message | Unchecked |
| Seg Delim | Segment delimiter | 0A |
| Delimiters | Standard HL7 delimiters | MSH ^~\& |
| Send By | Sending facility / application | <Sending Facility> |
| Return ACK | Option to enable return of ACK message | Checked |

Note: Use recommended values unless Site requests otherwise.

1.1.4. MacLab Integration Configuration

The MacLab equipment exports data by depositing files in a directory that the ProSolv Measurement Importer monitors. These files end with the extension "h7c" for cardiac and "h7e" for electrophysiology. The cardiac file (not the EP) are retrieved and put into the database. Once a match is made on the client the values are added to variables and are available to be put on the report.


MacLab is also capable of exporting tiff image files that contain wave forms. MacLab's output may or may not contain these image files. The export with tiff image files may populate the ProSolv Importer input directory in two formats:

- 1) Deposited directly in the input folder with the names matching the "h7?" and ending with snapxxx.tif :
 - a. (0a9ac3e9-0981-41cc-aa0b-939d5d8c8acdsnap0001.tif)
 - b. 0a9ac3e9-0981-41cc-aa0b-939d5d8c8acd.h7c
- 2) A new subdirectory holding the data and tiff files put into the input directory with names that may not match.

These image tiff files are put into the studies image directory for viewing when a match takes place at the client

To configure the MacLab integration:


1. Start the ProSolv® Service Manager.
2. Select the "Measurement Importer" tab.




3. Select <Add> to configure a new MacLab Importer.

4. Make sure that MacLab Hemo Importer is selected in the window before configuring.

Select "MacLab Hemo Importer"



| Name | Description | Recommended Values |
|------------------------------------|---|---|
| Database Settings | | |
| Images Directory | Path to the Images Folder. | I:\Images |
| ProSolv Database Connection String | Connection string to the ProSolv® database. Click the "..." button to bring up the Connection editor. | Provider = SQLOLEDB.1 Password = psci Persist Security Info = True User ID = prosolv Initial Catalog = Prosolv Data Source = localhost |
| File Importer Settings | | |
| Input Directory | Path to the folder where the hemodynamic output files will be deposited by MacLab. | I:\GEMacLabIn\Input |
| Search Pattern | Filter pattern of files to import. MacLab = "*.h7?" This allows | *.h7? |

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 25 of 83 |

| Name | Description | Recommended Values |
|--------------------------|---|--------------------|
| | for h7c and h7e files | |
| Monitor Subdirectories | False = will not process files in subdirectories. | True |
| Importer Settings | | |
| Machine Name | A unique name that describes the importer task. | MacLab Importer |
| Sleep Interval | Number of seconds between polls of the Input Directory. | 30 |

Note: Use recommended values unless Site requests otherwise.

5. After completing the entries select <Save>, select the entry and start the service.


Note: Configure the above per your specific system. The MacLab system will also need to be configured so that it deposits files in the Input Directory that is specified in this configuration step.

1.1.5. Siemens Sensis Integration Configuration


The new generation of our recording solution is adapting to the user's workflow featuring customizable workflow support programs and new advanced measurement functions. As a truly integrated solution AXIOM Sensis XP combines both electrophysiology and hemodynamics in one system.

To configure the Siemens Sensis integration:

1. Start the ProSolv® Service Manager.
2. Select the "Measurement Importer" tab.



3. Select <Add> to configure a new Siemens Importer.
4. Make sure that Siemens Hemo Importer is selected in the window before configuring.




| Configuration Setting | Description | Recommended Value |
|------------------------------------|--|-------------------|
| Plugin Type | Drop down list of supported importers | |
| Image Directory | Directory location where the ProSolv Study images are located. | |
| ProSolv Database Connection String | Connection string pointing to the ProSolv database. | |
| Search Pattern | Filter pattern of files to import. Siemens = "*.his" | |
| Input Directory | Directory location where the hemodynamic output files will be deposited by the OEM and then processed. | |
| Monitor Subdirectories | False = will not process files in subdirectories. True = not supported | |
| Machine Name | Unique name to describe this importer task. | |
| Sleep Interval | Amount of time to sleep in between checks for new input files. | |

5. After completing the entries hit <Save>, select the entry and start the service.

Note: Configure the above per your specific system. The Siemens Sensis system will also need to be configured so that it deposits files in the Input Directory that is specified in this configuration step.

Note: The screen below displays the recommended configuration values.




1.2. ProSolv® Server

The ProSolv® Server application enables users to share centrally located databases, as well as images and reports. The ProSolv® Server is actually a background process that gets spawned by the ProSolv® Launcher Windows Service. This allows a site to have multiple ProSolv® Servers running on the same machine. The ProSolv® WebLab Web Service provides HTTP communications between the ProSolv® Client and Server.


1.2.1. Configuration

To configure the ProSolv® Server:


1. Start the ProSolv® Service Manager.
2. Click the “ProSolv® Server” tab.






3. To add a new ProSolv® Server process, click the **Add** button. This will bring up the Add ProSolv® Server screen. Alternatively, select existing ProSolv® Server and click **Settings**.
4. There needs to be a minimum of two entries under ProSolv® Server -> Databases list. **NOTE:** USERLOG and server are case-sensitive. USERLOG should be all uppercase and server must be all lowercase.
 - a. USERLOG – for user activity auditing



b. server – the actual ProSolv® server configuration object:




c. Complete the ProSolv® Server tab appears as follows:



| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Process Name | Unique name to describe this ProSolv® Server process | |
| Executable | Path to the ProSolv.exe to spawn | |
| Arguments | Optional arguments to pass to the ProSolv®.exe process, default is 'apps\server\server.pxe' | |
| Port Number | TCP Port Number | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.


| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 32 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Saved Files Directory | Folder that stores the saved P-Code files | |
| Database | List of databases that this server will use. The following names must be used: server - ProSolv® Database USERLOG – User log database | |

5. Click **OK** to save the settings.
6. To start the newly configured process, click the Green **Start** button to the right of the Process List.
7. To stop the process, click the Red **Stop** button.
8. To restart the process, click the **Restart** button.
9. To remove a process, select it in the list and click the **Delete** button.
10. To modify a process' settings, select it in the list and click the **Settings** button.
11. To disable a process from being spawned and pinged by Watchdog, click the **Disable** button.

To configure the Server Service – Non Service Manager:

1. **Server Files:** Navigate to the Services 4.0\apps\server directory and modify deletepaths.txt to reference the locations of the deletedfiles.mdb database and the deletetemp directory. This is typically under the \prosolv database directory.



2. **Deleted Files:** Navigate to the \ProSolv Database directory referenced above and verify the users group has full permissions to the deletedfiles.mdb database and the deletetemp folder.

1.3. License Server


The License Server is in charge of distributing concurrent licenses to client users. The License Server is a Web Service that handles license requests from the ProSolv® Server. The License Server determines the amount of licenses for each of the application types from the license file that is installed on the server.

The License Server Tab in the Service Manager displays a list of all users that are currently logged into the ProSolv® Server. It also displayed the number of licenses that are in use and how many are remaining for each application type.


1.3.1. Configuration

To configure the License Server:

1. Start the ProSolv® Service Manager.
2. Click the "License Server" tab.



- If this is a remote/secondary server that is requesting licenses from a primary license server put the address of the main license server on the URL line. Ex: <http://192.168.0.1/LicenseServer/Service.aspx> where the IP address is the IP address of the main license server.
- To modify the configuration, click the **Settings** button. This will bring up the License Server Security Settings screen.




| Configuration Setting | Description | Recommended Value |
|---------------------------|--|-------------------|
| Enforce License Limits | Boolean flag indicating if the License Server should reject a license request if all of the licenses are currently in use. | |
| Remote Server Access List | List of Remote ProSolv® Servers that can request a license from this License Server. | |

- If you are configuring a primary/main license server that will be sharing licenses with other license servers, you must include the IP address and Port of the requesting ProSolv® server on this screen.
- To save the settings, click the **OK** button.

1.3.2. Enterprise Licensing Methodology

In many cases, an Enterprise purchases a multi-user concurrent license for the entire enterprise comprised of a main license server and multiple local (or remote) license servers. Configure the license servers as follows:

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 35 of 83 |

1. On the primary/main License Server, create and add a license for the total number of concurrent licenses for the entire enterprise.
2. On the local/remote License Server(s) create and add a single user concurrent license. The URL will point to the main license server.
3. By default the "Enforce License Limits" flag on both should be unchecked.

1.4. Watchdog Service

The Watchdog Service monitors the ProSolv® Services and is capable of sending alert emails to administrative staff. The Watchdog Service will ping each service on a designated time interval. The default setting performs the ping every 5 minutes. If the Watchdog cannot ping the service, then it attempts to restart the service. It will continue to try to restart the service for the next 30 minutes. If the Watchdog cannot restart the service after 30 minutes, it flags the service to indicate that the service cannot be restarted. Watchdog can also monitor drive space and send an alert email if free space falls below a certain percentage.


The Watchdog Tab in the Service Manager shows the list of services that are currently being monitored. It displays an icon next to each task. The following icons can be displayed:

- Green Check Icon – indicates that the Service is running
- Yellow X Icon – indicates that the Service has failed a ping test and that Watchdog is trying to restart the service
- Red X Icon – indicates that Watchdog was not able to restart the service after 30 minutes


1.4.1. Configuration

To configure the Watchdog Service:

1. Start the ProSolv® Service Manager.
2. Click the "Watchdog" tab.




3. To modify the configuration, click the **Settings** button. This will bring up the Watchdog Settings screen.



| Configuration Setting | Description | Recommended Value |
|------------------------------|---|-------------------|
| Perform Tasks every | Time interval to perform the ping tests. | |
| Send warning and error email | Boolean flag indicating if Watchdog should sent alert emails. | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 37 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| SMTP Server | The name of the SMTP server to use for sending email. | |
| SMTP Port | The port number for SMTP | |
| SMTP Username | User account that has access to the SMTP server | |
| SMTP Password | SMTP user account password | |
| Enable SSL | Boolean flag indicating if SSL should be used when sending email | |
| Warnings | List of email addresses that will receive warning emails | |
| Errors | List of email addresses that will receive error emails | |
| Disk Space Monitoring | List of drives and percentage limits that Watchdog will monitor | |
| Service Monitoring | List of Services that Watchdog should monitor. Services that are not checked will not be pinged by Watchdog | |


1.5. DICOM Server Service (DICOM Store SCP)

The ProSolv® DICOM Server is intended to receive and convert DICOM images sent over a network from DICOM devices. Studies received by the DICOM Server are stored in the ProSolv® database. The DICOM Server is actually a background process that gets spawned by the ProSolv® Launcher Windows Service. This allows a site to have multiple DICOM Servers running on the same machine.


1.5.1. Configuration

To configure the DICOM Server Service:

1. Start the ProSolv® Service Manager.
2. Click the "DICOM Server" tab.



3. To add a new DICOM Server process, click the **Add** button. This opens the Add DICOM Server screen.



| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Process Name | Unique name to describe this DICOM Server process | |
| Executable | Path to the psdcm serv.exe to spawn | |
| Arguments | Optional arguments to pass to the psdcm serv.exe process | |
| AE Title | The AE Title of this Store SCU process | |
| Port Number | The port number for this Store SCU process | |
| Max PDU Length | The maximum length that a PDU can be | |
| ARTIM Timeout | ARTIM timeout value in milliseconds | |
| Read Timeout | Network read timeout in milliseconds | |
| Write Timeout | Network write timeout in milliseconds | |
| Image Path | Directory location for Image Store | |
| Store Commitment | Boolean flag indicating if Store Commitment should be included in association | |
| Auto Shuffle | Boolean flag indicating if the DICOM Server should attempt to autos shuffle the study images | |
| Show Debug Code | Boolean flag indicating if the DICOM | |


THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| | Server should log debug statements | |
| Database Settings | Database settings for the DICOM and ProSolv® databases | |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| ADT Enabled | Boolean flag indicating if the DICOM Server should perform ADT operations | |
| Database Type | Database settings for the ADT Database | |
| Accession Field Name | The Accession Field Name | |
| Accession Padding | The Accession padding | |
| MRN Padding | MRN Padding | |
| Operator Field Name | The name to use for the Operator DICOM tag | |
| Store SCU Path | The folder to store images for Store SCU | |
| Store SCU Path #2 | The second folder to store images for Store SCU | |
| SR Path | The folder where SR data is stored | |


THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Accept RLE | Boolean Flag indicating if Server should accept RLE images | |
| View Label | Label to use for the view field | |
| Default Exam Type | Value to use for the default exam field | |
| Webviewer IP Address | The URL of the ProSolv® webviewer | |
| Status Update Folder | The folder to store status update logs | |
| Ignore Accession | Boolean flag indicating if accession should be ignored | |



| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Misc Settings | Select which of the 3 misc items to configure | |
| Add | Add a new entry | |
| Delete | Remove an entry | |
| Settings | Edit an entry | |
| SOP | A DICOM SOP Class | |
| SOP Type | Exam type to use for that SOP Class | |
| Tag Name | DICOM Tag to match. This is a hex number without any leader (i.e. 00080018). | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 42 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Tag DB | The name of the column in the EXAM table to store the tag data | |
| Tag Length | The maximum length of data to store | |
| AE Title | The AE Title of an instrument used to determine an exam location. Leave this field blank to specify a default Exam Location if nothing else matches. The default location should be the last one in the list. | |
| Location | The Exam Location to use for that AE title | |

4. To start the newly configured process, click the Green **Start** button to the right of the Process List.
5. To stop the process, click the Red **Stop** button.
6. To restart the process, click the **Restart** button.
7. To remove a process, select it in the list and click the **Delete** button.
8. To modify a process' settings, select it in the list and click the **Settings** button.
9. To disable a process from being spawned and pinged by Watchdog, click the **Disable** button.


1.6. DICOM Query/Retrieve SCP Service

The QR SCP Service provides DICOM QR services to requesting QR SCU devices. The QR SCP service is multithreaded and is capable of running multiple "tasks" that can act as separate QR SCP devices.

1.6.1. Configuration

To configure the DICOM Query/Retrieve SCP Service:

1. Start the ProSolv® Service Manager.
2. Click the "DICOM QR" tab.
3. To configure the service, it must be started. Click the Green **Start** button in the Service Control group box to start the service.
4. To add a new task, click the **Add** button. This opens the DICOM QR Task Settings screen.



| Configuration Setting | Description | Recommended Value |
|----------------------------------|--|-------------------|
| DICOM Database Connection String | Connection string to the ProSolv® DICOM database. Click the “...” button to bring up the Connection editor | |
| Acceptable AE Titles | List of QR SCU AE Titles that the QR SCP Service will accept association from | |
| AE Title | The AE Title of this QR SCP task | |
| Destination Store SCP's | List of Store SCP's that can receive images from the QR SCP Task | |
| Task Name | Unique name to describe this QR SCP Task | |

5. Click the **Save** button to add the task after it has been configured.
6. To start the newly configured task, click the Green **Start** button to the right of the Task List.
7. To stop the task, click the Red **Stop** button.
8. To restart the task, click the **Restart** button.
9. To remove a task, select it in the list and click the **Delete** button.
10. To modify a task's settings, select it in the list and click the **Settings** button.

1.7. DICOM Store SCU Service

The Store SCU Service monitors an input folder for DICOM images and sends them to a designated Store SCP device. The service is multithreaded and is capable of running

multiple “tasks” that can monitor several input folders and forward images to several SCP devices.


1.7.1. Configuration

To configure the DICOM Store SCU Service:

1. Start the ProSolv® Service Manager.
2. Click the “Store SCU” tab.
3. To configure the service, it must be running. started. Click the Green **Start** button in the Service Control group box to start the service.
4. To add a new Task, click the **Add** button. This opens the Store Task Settings screen.

| Store Task Settings | |
|---|------------|
| Input Folder Settings | |
| Input Folder | |
| Search Pattern | *.* |
| Store SCP Settings | |
| SCP AE Title | PROSOLV |
| SCP Port Number | 104 |
| SCP Server | |
| Store SCU Settings | |
| Include SR | True |
| Lookup Interval | 200 |
| Read Timeout | 40 |
| SCU AE Title | PROSOLVSCU |
| Storage Commitment | False |
| Storage Commitment Port | 104 |
| Task Name | |
| Write Timeout | 200 |
| Include SR Include SR related SOP's on the same association. | |
| <input type="button" value="Save"/> <input type="button" value="Cancel"/> | |

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Input Folder | Directory that this task should monitor for input DICOM images | |
| Search Pattern | Search filter to reduce the file types that the task should monitor. *.* will include all files. | |
| SCP AE Title | The AE Title of the Store SCP device that is to receive the images | |
| SCP Port Number | The port number of the Store SCP device that is to receive the images | |
| SCP Server | The machine name or IP address of the Store SCP device that is to receive the images | |

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 45 of 83 |

| Configuration Setting | Description | Recommended Value |
|-------------------------|--|-------------------|
| Include SR | Boolean flag that indicates if SR related SOP's are included on the association | |
| Lookup Interval | Amount of time in seconds that the task will sleep in between checks for new image files | |
| Read Timeout | Network read timeout in milliseconds | |
| SCU AE Title | The AE Title of this Store SCU Task | |
| Storage Commitment | Boolean flag indicating whether or not to include Storage Commitment on the same association | |
| Storage Commitment Port | Port number to use for Storage Commitment | |
| Task Name | Unique name used to describe this Store SCU Task | |
| Write Timeout | Network writes timeout in milliseconds | |

5. Click the **Save** button to add the task after it has been configured.
6. To start the newly configured task, click the Green **Start** button to the right of the Task List.
7. To stop the task, click the Red **Stop** button.
8. To restart the task, click the **Restart** button.
9. To remove a task, select it in the list and click the **Delete** button.
10. To modify a task's settings, select it in the list and click the **Settings** button.


1.8. Image Compression Service

The Image Compression Service creates compressed versions of images in ProSolv® Studies. It monitors the file system for the creation of new ProSolv® Studies. When it detects a new study, it processes the contents of the study image folder and creates compressed versions of the images files. Note, Image Compression Service currently only compresses loops and only from modalities XA (Catheterization) or US (Ultrasound).

1.8.1. Configuration


To configure the Image Compression Service:

1. Start the ProSolv® Service Manager.
2. Click the "Compression" tab.



| Configuration Setting | Description | Recommended Value |
|-------------------------------------|---|-------------------|
| ProSolv® Database Connection String | Connection string pointing to the ProSolv® database | |
| Original | Check this box to create JPEG2K versions of the original images | |
| Clinical | Check this box to create Clinical versions of the images at the designated compression factor | |

3. Navigate to '\ProSolv CardioVascular\Services 4.0\Image Compression Service' folder and open ImageCompression.exe.config file using Notepad for editing.
4. Update configuration values as follows:
 - a. ImageDir: location of image files
 - b. ImageDir: MachineName: host name of the server running the Compression service, if cluster, use virtual host name defined in ProSolv® Apps cluster group



```
<?xml version="1.0" encoding="utf-8"?>
<configuration>
  <configSections>
    <section name="compressionConfigs" type="System.Configuration.IgnoreSectionHandler" />
  </configSections>
  <appSettings>
    <add key="RemotePort" value="9998" />
    <add key="imageDir" value="D:\Images" />
  </appSettings>
  <connectionStrings>
    <add name="ProsolvDB" connectionString="Provider=sqloledb;uid=prosolv;pwd=psci;data source=localhost;initial
catalog=Prosolv"
    providerName="System.Data.OleDb" />
  </connectionStrings>
  <compressionConfigs>
    <compressionConfig machineName="machine1" className="Prosolv.ImageCompression.Compressor.Compressor"
assemblyName="Prosolv.ImageCompression.Compressor" secondsToSleep="100" loCompression="40" medCompression="5"
compressLo="True" compressMed="False" compressOriginal="False" />
  </compressionConfigs>
</configuration>
```


1.9. Archiver Service

The Archiver Service archives ProSolv® Study images to free disk space on the image store. The Archiver can archive images to a DICOM PACS, FTP Site, or a Network Share. The Archiver can run multiple “jobs” that can be scheduled to run at different times.

1.9.1. Configuration

To configure the Archiver Service:


1. Start the ProSolv® Service Manager.
2. Click the “Archiver” tab.




3. To create a new Archive Job, click the **Add Job** button. This opens the Archive Job Wizard.
4. To delete a job, select the job in the list and click **Delete Job**.
5. To modify a job's settings, select the job in the list and click **Settings**.
6. To manually start an Archive Job, select the job in the list and click **Archive Now**.

1.9.2. Archiver Job Wizard

The Archiver Job Wizard provides an interface to allow users to more easily configure the settings for an Archive Job. It guides the user through the proper steps based on user input.




| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Archiver Job Name | Unique name to describe this Archive Job | |





| Configuration Setting | Description | Recommended Value |
|------------------------------|---|-------------------|
| Image Path Location | Directory where the ProSolv® Study images are located | |
| Delete study after archival? | Checking yes here will delete the local ProSolv® Study images after they have been archived. Checking No here will not delete the local copy. | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.




| Configuration Setting | Description | Recommended Value |
|------------------------|--|-------------------|
| FTP Server | The study will be archived to a designated FTP server. | |
| Local or Network Drive | The study will be archived to a designated local or remote file share. | |
| ProSolv® Server | The study will be archived to another ProSolv® Server. | |
| DICOM PACS | The study will be archived to a designated DICOM PACS. | |



| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 51 of 83 |


| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| IP Address | FTP Server IP Address or host name | |
| User | FTP account Username | |
| Password | FTP account Password | |
| Root Path | The root path of where to upload the images on the FTP server | |




| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Archive Label | Unique name used to label this file share | |
| Root Path | Directory path to the file share resource. NOTE: Remote files shares must be designated as an UNC path. | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.


| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| IP Address | IP Address or host name of the ProSolv® Server | |
| Port | Port number that the ProSolv® Server is running | |
| User | ProSolv® Server user account to use when archiving | |
| Password | ProSolv® Server user account password to use when archiving | |
| Root Path | Root file path of where to store the images on the ProSolv® Server | |




| Configuration Setting | Description | Recommended Value |
|-------------------------|---|-------------------|
| Store SCP | The Store SCP device that is to receive the archived images | |
| Store SCU | The Store SCU device that is to perform the archiving | |
| Move AE Title | The AE Title of the Device that is to perform the DICOM move | |
| Query SCP | The Query SCP device that should be used when querying studies | |
| Manage DICOM Components | Click this button to launch the DICOM Device Configuration Screen to add or remove DICOM components | |



| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Add | Click this button to define a new device | |
| Delete | Select the component to delete, then click this button to delete the selected device | |
| Modify | Select the component to modify in the list, then click this button to modify the selected device | |
| Save | Click this button to save the newly defined device | |
| Name | Unique name to describe this DICOM device | |
| IP Address | IP Address or hostname where the device resides | |
| Port | Port number that the device is running | |
| Service Type | Available options here are: QR SCP, Store SCU, or Store SCP | |


| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 55 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Maximum Connections | The Maximum number of connections allowed on this device (not required) | |
| Timeout | Network read/write timeout values (not required) | |




| Configuration Setting | Description | Recommended Value |
|----------------------------------|--|-------------------|
| All Studies | Check this option to archive all of the ProSolv® studies | |
| Only studies that | Check this option to define a set of criteria that will determine which studies should be archived | |
| Have preliminary reports | Archive studies that have preliminary reports | |
| Have finalized reports | Archive studies that have finalized reports | |
| Are of a specific type | Archive studies that are of a specific type | |
| Were read at a specific location | Archive studies that were read at a specific location | |
| Had a specific reader | Archive studies that had a specific reader | |
| Were performed after | Archive studies that were performed after a designated date | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.


| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 56 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Are older than | Archive studies that are older than a designated number of days | |
| Are smaller than | Archive studies that are smaller than a designated size | |
| No more than | Set a maximum number of studies that can be archived | |




| Configuration Setting | Description | Recommended Value |
|---------------------------|--|-------------------|
| Left list box | Lists the current study types that will be archived | |
| Right text box | Designate a new study type to be archived. Click New to add the new type. | |
| Load types from this file | Select a predefined type's file to load the study types from. | |


THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.



| Configuration Setting | Description | Recommended Value |
|---------------------------|---|-------------------|
| Left list box | Lists the current study locations that will be archived | |
| Right text box | Designate new study locations to be archived. Click New to add the new type. | |
| Load types from this file | Select a predefined locations file to load the study locations from. | |




| Configuration Setting | Description | Recommended Value |
|---------------------------|--|-------------------|
| Left list box | Lists the current study reading doctors that will be archived | |
| Right text box | Designate a new reading doctor to be archived. Click New to add the new type. | |
| Load types from this file | Select a predefined reading doctor's file to load the study reading doctors from. | |



| Configuration Setting | Description | Recommended Value |
|---------------------------------------|--|-------------------|
| Performed more than | Designate a number of days that a study must be in age to be removed after archiving | |
| Reloaded more than | Designate a number of days that a study must have been reloaded by to be removed after archiving. | |
| Verified by a third-party application | Check this box to remove a study after archiving if it has been verified by a third-party application. | |

1.9.3. Archive to DICOM PACS

Archiving to DICOM PACS involves six processes. First the Archiver Store SCU (1) sends images to the PACS Store SCP (2). When retrieving the images the Archiver Query/Retrieve(QR) SCU(3) tells the PACS QR SCP(4) to

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 59 of 83 |

send the images with the PACS Store SCU(5) back to our Dicom Server Store SCP(6).

From the PACS we need information about their Store SCP and QR SCP. The information needed is the AE Title, IP Address, and port number of each SCP. In some cases, this information is the same for both the Store and QR SCP.

Likewise the PACS needs to know information about the Archiver Store SCU, Archiver QR SCU, and the Dicom Server Store SCP. For the Store SCU and QR SCU, they need the AE Title and in some cases they also want the IP Address. For the Dicom Server Store SCP, they need the AE Title, IP Address, and port number. Note: A different ProSolv® Dicom Server instance may be configured for the retrieve but it is not necessary.

From the "Manage Dicom Components" you need to setup 3 different Dicom components: Archiver Store SCU, PACS Store SCP, and PACS QR SCP.


Store SCU Component

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Name | StoreSCU | |
| IP Address | 127.0.0.1 (this is not used but still need to enter data) | |
| AETitle | PROSOLV® (this is what is normally used and what you need to give to PACS) | |
| Port | 0 (this is also not used but still need to enter data) | |
| Service Type | Choose STORESCU | |
| Maximum Connections | 2 (consult with PACS if they have a limit of connections) | |
| Timeout | 3 (minutes) | |

Store SCP Component

| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Name | PACS Store SCP | |
| IP Address | xx.xxx.xxx.xxx (PACS Store SCP IP is supplied by PACS) | |
| AETitle | PACSAE (PACS Store SCP AE Title is supplied by PACS) | |
| Port | 104 (PACS Store SCP Port Number is supplied by PACS) | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 60 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Service Type | Choose STORESCP | |
| Maximum Connections | 2 (consult with PACS if they have a limit of connections) | |
| Timeout | 3 (minutes) | |

QRore SCP Component

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Name | PACS Store SCP | |
| IP Address | xx.xxx.xxx.xxx (PACS Store SCP IP is supplied by PACS) | |
| AETitle | PACSAE (PACS Store SCP AE Title is supplied by PACS) | |
| Port | 104 (PACS Store SCP Port Number is supplied by PACS) | |
| Service Type | Choose QRSCP | |
| Maximum Connections | 2 (consult with PACS if they have a limit of connections) | |
| Timeout | 3 (minutes) | |

After setting up these components, they must be selected from the drop down list boxes on the "Configure Dicom Access" dialog. If you have only setup these three they should already be selected.

For the Move AE Title, enter the AE Title of the Dicom Server SCP that you gave to PACS (the SCP that PACS sends the images to during its retrieval process).


1.10. Retriever Service

The Retriever Service retrieves ProSolv® Studies that have been archived from an external PACS.


1.10.1. Configuration

To configure the Retriever Service:

1. ADT needs to be configured to send MRN based requests to retrieve studies from external PACS. Please refer to the ADT configuration section for details on how to configure ADT to send requests to the Retriever.
2. Go to Archiver tab. Click **Add Job, Select Retriever**, and click **Next**.




3. Please use Manage DICOM Components as described in a previous section.
 - a. Store SCP
 - b. Query/Retrieve SCP



| Configuration Setting | Description | Recommended Value |
|-----------------------|--|-------------------|
| Retriever AETitle | Retriever's AE Title. By default, this is set to ProSolv@Retriever. | |
| Store SCP | AETitle for Move SCP. The server where the images need to be sent to store be stored | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.


| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 62 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Query SCP | AETitle for Query/Retrieve SCP | |
| Modality Filter | Modality filter used to query external PACS. In the above example, only "US" exams will be retrieved from QR SCP "PROSOLV_QRSCP" | |
| Date Range | Date range filter to query external PACS. If no date range is selected then all the studies in the external PACS will be queried. | |
| Previous Studies | By default, it is set to ALL. A value in this field indicates previous number of studies to be retrieved for a particular patient. For example, if value 10 is specified in this field then only the last 10 studies for a patient will be retrieved. | |

4. Click **Apply** to save the configuration.
5. Click **Finish** to close the configuration dialog.

1.10.2. Retrieve Queues

The Retriever Tab in the Service Manager displays status on its current tasks. The external Requests queue lists all the requests that need to be retrieved from external PACS based on Patients MRNs. The Internal Requests queue lists all the retrieve requests based on Study UIDs.




1.11. Measurement Importer Service

The Measurement Importer Service imports measurements from several different medical devices. The service consists of a set of “plugins”. Each plugin can handle importing a measurement from a different type of device (e. g. HPDT, Burdick, etc.). The General Measurement Import Service is multithreaded and can run each plugin in a separate thread. All US measurement interfaces not supported by the Measurement Importer Service are supported by stand-alone applications, such as Toshiba, ATL, and Cypress.

1.11.1. Configuration

To configure the Measurement Importer Service:

1. Start the ProSolv® Service Manager.
2. Click the “Measurement Importer” tab.



3. To add a new importer task, click the **Add** button. This opens the Importer settings screen and allows the user to configure the setting for the importer.
4. To start the importer, click the Green **Start** button.
5. To stop the importer, click the Red **Stop** button.
6. To restart the importer, click the **Restart** button.
7. To remove an importer, select it in the list and click the **Delete** button.
8. To modify an importer's settings, select it in the list and click the **Settings** button.


1.12. HPDT Importer

| Configuration Setting | Description | Recommended Value |
|-------------------------------------|--|-------------------|
| ProSolv® Database Connection String | Connection string pointing to the ProSolv® database | |
| Image Directory | Directory location where the ProSolv® Study images are located | |
| Sleep Interval | Amount of time to sleep in between checks for new input files | |
| Machine Name | Unique name to describe this importer task | |
| Port Number | TCP Port number to listen on for input data | |

1.13. HD11XE Importer

| Configuration Setting | Description | Recommended Value |
|-------------------------------------|---|-------------------|
| ProSolv® Database Connection String | Connection string pointing to the ProSolv® database | |
| Image Directory | Directory location where the ProSolv® Study images are located | |
| Search Pattern | Filter pattern of files to import. *.* will include all files | |
| Input Directory | Directory location where the Epiphany XML output files will be located | |
| Monitor Subdirectories | Boolean indicating if the service should monitor subdirectories for input files | |
| Machine Name | Unique name to describe this importer | |
| Sleep Interval | Amount of time to sleep in between checks for new data | |

1.14. Burdick/Marquette/Quinton Treadmill Importers



| Configuration Setting | Description | Recommended Value |
|-------------------------------------|---|-------------------|
| ProSolv® Database Connection String | Connection string pointing to the ProSolv® database | |
| Image Directory | Directory location where the ProSolv® Study images are located | |
| Search Pattern | Filter pattern of files to import. *.* will include all files | |
| Input Directory | Directory location where the Epiphany XML output files will be located. | |
| Monitor Subdirectories | Boolean indicating if the service should monitor subdirectories for input files | |
| Machine Name | Unique name to describe this importer task | |
| Sleep Interval | Amount of time to sleep in between checks for new input files | |
| PDF Off | Boolean flag indicating if PDF Conversion should be turned on/off | |

1.15. Mennen Hemodynamic Import Service


The Mennen Service imports measurements from the Mennen Hemodynamic device.

1.15.1. Configuration

To configure the Mennen Hemodynamic Import Service:

1. Start the ProSolv® Service Manager.

2. Click the "Mennen DT" tab.



| Configuration Setting | Description | Recommended Value |
|----------------------------|---|-------------------|
| Input Folder | Directory location to monitor for input data files | |
| Archive Folder | Directory location to store input data files after processing | |
| Date Time Format String | Date/Time format to use when add data to the database | |
| ProSolv® Database Settings | Connection string pointing to the ProSolv® database | |


1.16. HL7 ADT Service

The Admission Discharge Transfer (ADT) Service accepts and subsequently parses HL7 messages. The ADT Service retrieves information necessary to enable an efficient patient data and image management solution. The ADT Service is multithreaded and is capable of running multiple "tasks" that can serve multiple TCP input ports.


1.16.1. Configuration

To configure the HL7 ADT Service:

1. Start the ProSolv® Service Manager.
2. Click the "ADT" tab.




3. To configure the service, it must be running. Click the Green **Start** button in the Service Control group box to start the service.
4. To add a new Task, click the **Add** button. This opens the ADT Task Settings screen.



| Configuration Setting | Description | Recommended Value |
|-------------------------------------|---|-------------------|
| HL7 Delimiters | List of delimiters to use for the HL7 messages | |
| Create Study | Boolean flag indicating whether the service should create ProSolv® studies in the database | |
| Images Directory | Location of the ProSolv® study images | |
| Patient ID as Main ID | Boolean flag indicating if the Patient ID should be used as the main identifier | |
| Store HL7 Messages | Boolean flag indicating if HL7 messages should be written to log files | |
| Study UID Prefix | The prefix to use for Study UIDs | |
| Task Name | Unique name to describe this ADT task | |
| ADT Database Connection String | Connection string to the ADT database (by default, included into the ProSolv® database). Make sure the password is included in the connection string. | |
| DICOM Database Connection String | Connection string to the DICOM database. Make sure the password is included in the connection string. | |
| ProSolv® Database Connection String | Connection string to the ProSolv® database. Make sure the password is included in the connection string. | |

5. Click the **Save** button to add the task after it has been configured.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 71 of 83 |

6. To start the newly configured task, click the Green **Start** button to the right of the Task List.
7. To stop the task, click the Red **Stop** button.
8. To restart the task, click the **Restart** button.
9. To remove a task, select it in the list and click the **Delete** button.
10. To modify a task's settings, select it in the list and click the **Settings** button.


1.17. HL7 Report Messenger Service

The Report Messenger Service sends HL7 formatted messages from a designated folder to the designated receiving system. The Report Messenger Service is multithreaded and is capable of running multiple "tasks" that can monitor several input folders and forward to several receiving systems.

1.17.1. Configuration

To configure the HL7 Report Messenger Service:


1. Start the ProSolv® Service Manager.
2. Click the "Report Messenger" tab.



3. To configure the service, it must be started. Click the Green **Start** button in the Service Control group box to start the service.
4. To add a new task, click the **Add** button. This opens the Report Messenger Settings screen.

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Machine Name | Unique name to describe this HL7 task | |
| Seconds to Sleep | Amount of time to sleep between checks for new HL7 files | |
| Watch Directory | Input directory to monitor for new HL7 messages | |
| Search Pattern | Search filter to reduce the types of files that should be processed | |
| Recurse Directory | Boolean flag indicating if subdirectories should be monitored | |
| IP Add.\Host Name | IP Address or Host Name of the receiving application | |
| Watch Port | Port number of the receiving application | |
| Wait Time | Time in milliseconds to wait | |
| Num Wait Cycles | Numbers of cycles to perform when waiting | |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 74 of 83 |

| Configuration Setting | Description | Recommended Value |
|-----------------------|---|-------------------|
| Persistent | Boolean flag indicating TBD | |
| Remove LF | Boolean flag indicating if line feeds should be removed during processing | |
| Seg Delim | Designates the Segment Delimiters for the message | |
| Delimiters | Designates the message delimiters | |
| Send By | Sending device name | |
| Return ACK | Boolean flag indicating if a Return ACK should be sent | |

5. Click the **Save** button to add the task after it has been configured.
6. To start the newly configured task, click the Green **Start** button to the right of the Task List.
7. To stop the task, click the Red **Stop** button.
8. To restart the task, click the **Restart** button.
9. To remove a task, select it in the list and click the **Delete** button.
10. To modify a task's settings, select it in the list and click the **Settings** button.

2. DICOM Modality Worklist Application

The ProSolv® DMWL interface responds to the Worklist requests sent to the application by the different modalities.

2.1. Configuration

To configure the DICOM Modality Worklist Application:


1. Place the ProSolv® DMWL executable file (**ProSolv® WorkList Server.exe**) and the **psmwls.ini** file in the same application folder (\ProSolv CardioVascular\SR).
2. Open (double-click) the **psmwls.ini** file and specify the appropriate configuration information as shown below. The format of the file is dependent on the type of database engine used (SQL Server or Jet Access; only SQL Server is supported as of 4.0.)
3. The shortcut to ProSolv® DMWL should be added to the Startup folder so ProSolv® DMWL will launch on login.

psmwls.ini for a SQL Server database engine:

```


psmwls - Notepad
File Edit Format View Help
PROSOLVDMWL
112
45
1
<BLANK>
SQLServer
localhost
ADT_TMH
prosolv
psci
    
```

| Line Numbers | Explanation | Possible Values | Example |
|----------------------------------|---|--|--|
| 1 | This line indicates the AE title of the ProSolv DMWL application (the application will use this title to associate with the different modalities). | Printable <u>ASCII</u> characters that are supported by the Operating System | PROSOLVDMWL |
| 2 | This line indicates the port number used by the DMWL application to accept requests from the different modalities. | 1 – 65535 | 112 113 |
| 3 | This line indicates whether the ProSolv DMWL application should generate and send the study UID to the modality (in response to the study/order). | yes or no | no |
| 4 | This line indicates the number of AE titles following this line (if there are 2 acceptable AE titles that this application should use, enter 2 in this line). | 1 - infinity | 2 |
| 5 thru (4 + number at line 4) | This line indicates the acceptable AE titles for this application. The application will refuse associations to modalities with AE titles other than those specified in these lines. | Printable <u>ASCII</u> characters that are supported by the Operating System | <BLANK> (using this AE title would allow associations from all possible clients) |
| 6 (assuming line 4 contains "1") | This line indicates the type of the ADT database used by this application. | SQLServer or Access | SQLServer |

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 76 of 83 |

| Line Numbers | Explanation | Possible Values | Example |
|--------------|--|--|--|
| 7 | This line indicates the network address of the SQL Server (the IP address of the machine that has SQL Server running on it). | localhost w.x.y.z (where w,x,y,z are in range 0-255) | localhost 192.168.229.81 |
| 8 | This line indicates the name of the ADT database on the SQL Server. | Printable <u>ASCII</u> characters that are supported by the Operating System | ADT ADTTEST |
| 9 | This line indicates the username associated with the <u>ADT</u> database that would allow the ADT application to connect to the database directly. | Printable <u>ASCII</u> characters that are supported by the Operating System | ProSolv® adt sa |
| 10 | This line indicates the password associated with the ADT database for the user specified in line 9. | Printable <u>ASCII</u> characters that are supported by the Operating System | psci hl7231adt |

The following screen capture displays an example of how the **psmwls.ini** file should look when configuration is completed for a SQL Server database engine. It reflects a setup specifying one acceptable AE title (a **1** is placed in line 4 and there is only one line following it):




```

psmwls.ini - Notepad
File Edit Format View Help
PROSOLVMWL
112
yes
1
<BLANK>
SQLServer
localhost
ADT
prosolv
psci

```


2.2. Starting the ProSolv® DICOM MWL Application

To start the application, double-click the **ProSolv® Worklist** icon. The application will begin listening to requests automatically.

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 77 of 83 |

2.3. Exiting the ProSolv® DICOM MWL Application

To exit the application, click the **Exit** button (refer to the screen capture below):



3. DICOM Structured Reporting Application

The ProSolv® DICOM SR application processes DICOM Structured Reporting files from the specified folder. The application also stores the measurement information retrieved from the files into the ProSolv® database (to be subsequently displayed by the ProSolv® Client application).


3.1. Configuration

To configure the DICOM Modality SR application:

1. Place the ProSolv® DICOMSR executable file (ProSolv® SR EchoVascular.exe) and the SR.ini file into the application folder (\ProSolv CardioVascular\MWL).
2. Open the SR.ini file (double-click the icon) and specify the appropriate configuration information as shown below. The format of the file is dependent on the type of the database engine used (SQL Server or Jet Access).
3. The shortcut to ProSolv® DICOM SR should be added to the Startup folder so ProSolv® SR will launch on login. The first table is for SQL, the second is for Access database.


SQL

| Line number | Explanation | Possible Values | Example |
|-------------|-------------|-----------------|---------|
|-------------|-------------|-----------------|---------|

| | | |
|--|---|----------------------------|
| ProSolv CardioVascular  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 78 of 83 |

| Line number | Explanation | Possible Values | Example |
|-------------|--|--|---|
| 1 | The full path to the folder that contains the DICOM SR files to be processed by this application. | Windows formatted file path ending with a '**' symbol to indicate all files in that folder (directory path should be valid and existing) | C:\SRFolder* |
| 2 | How often DICOM SR application will check the directory for new DICOM SR objects. (milliseconds) | 5000 | 5000 |
| 3 | The type of ProSolv® database where retrieved file data is stored. | SQL Server or Access | SQL Server |
| 4 | The network address of the SQL Server (IP address of the machine that the SQL Server is running on). | localhost w.x.y.z (where w,x,y,z are in range 0-255) | localhost 192.168.229.81 |
| 5 | The name of the ProSolv® database on the SQL Server. | Printable <u>ASCII</u> characters that are supported by the Operating System | ProSolvDB ProSolvDBTEST |
| 6 | The username associated with the ProSolv® database that allows the DICOM SR application to connect to the database directly. | Printable <u>ASCII</u> characters that are supported by the Operating System | prosolv sa |
| 7 | The password associated with the ProSolv® database and the user specified in line 5. | Printable <u>ASCII</u> characters that are supported by the Operating System | psci |
| 8 | The type of DICOM Server database. | SQL Server or Access | SQL Server |
| 9 | The network address of the SQL Server (IP address of the machine that the SQL Server is running on). | localhost w.x.y.z (where w,x,y,z are in range 0-255) | localhost 192.168.229.81 |
| 10 | The name of the DICOM Server database on the SQL Server. | Printable <u>ASCII</u> characters that are supported by the Operating System | DCMServ DCMServTEST |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| | | |
|--|---|----------------------------|
| ProSolv CardioVascular  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 79 of 83 |

| Line number | Explanation | Possible Values | Example |
|-------------|---|--|-----------------------|
| 11 | The username associated with the DICOM Server database that allows the ADT application to connect to the database directly. | Printable <u>ASCII</u> characters that are supported by the Operating System | <i>prosolv sa</i> |
| 12 | The password associated with the DICOM Server database and the user specified in line 9. | Printable <u>ASCII</u> characters that are supported by the Operating System | <i>psci</i> |

Access

| Line number | Explanation | Possible Values | Example |
|-------------|--|---|--------------------------------|
| 1 | The full path to the folder that contains the DICOM SR files that are processed by this application. | Windows formatted file path ending with a '*' symbol to indicate all files in that folder (directory path should be valid and existing) | <i>C:\SRFolder*</i> |
| 2 | How often DICOM SR application will check the directory for new DICOM SR objects. (milliseconds) | <i>5000</i> | <i>5000</i> |
| 3 | The type of ProSolv® database where retrieved file data is stored. | <i>SQL Server</i> or <i>Access</i> | <i>Access</i> |
| 4 | The full path to the <i>.mdb</i> file on the ProSolv® database. | Windows formatted file path | <i>C:\DCMServ\prosolv.mdb</i> |
| 5 | The password associated with the ProSolv® database. | Printable <u>ASCII</u> characters that are supported by the Operating System | <i>stick</i> |
| 6 | The type of DICOM Server database. | <i>SQL Server</i> or <i>Access</i> | <i>Access</i> |
| 7 | The full path to the <i>.mdb</i> file on the DICOM Server database. | Windows formatted file path | <i>C:\DCMServ\dcmserve.mdb</i> |

THIS DOCUMENT CONTAINS CONFIDENTIAL AND PROPRIETARY INFORMATION OF PROBLEM SOLVING CONCEPTS, INC. D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY, AND SHALL NOT BE DUPLICATED OR DISTRIBUTED OR OTHERWISE TRANSFERRED TO A THIRD PARTY WITHOUT THE EXPRESS WRITTEN CONSENT OF PROBLEM SOLVING CONCEPTS, INC D/B/A PROSOLV CARDIOVASCULAR, A FUJIFILM COMPANY.

| Line number | Explanation | Possible Values | Example |
|-------------|---|--|---|
| 8 | The password associated with the DICOM Server database. | Printable <u>ASCII</u> characters that are supported by the Operating System | This line can be left blank if the database is not password protected |

4. Witt Hemodynamic Import Application

The Witt Service imports measurements from the Witt Hemodynamic device. It consists of two Windows applications that have to be run from the desktop.

4.1. Installation


To install the Witt Hemodynamic Import Application:

1. Copy the Witt files to the desired installation folder on the target machine (\\ProSolv CardioVascular\\Witt).

4.2. Configuration

To configure the Witt Hemodynamic Import Application:

1. Copy Witt.mdb file from the installation folder to main ProSolv Database folder.
2. Open the text file, pswitt.ini, from the installation folder in a text editor.
3. Edit the file to contain the information listed in the tables below. The first table is for SQL, the second is for Access database.




SQL

| Line number | Explanation | Possible Values | Example |
|-------------|-------------|-----------------|---------|
|-------------|-------------|-----------------|---------|

| Line number | Explanation | Possible Values | Example |
|-------------|---|---|---|
| 1 | The full path to the folder that contains the Witt files that are processed by this application. | Windows formatted file path ending with a '*' symbol to indicate all files in that folder (directory path should be valid and existing) | \\<AppServer>\From_Witt* |
| 2 | The type of Witt database where retrieved file data is stored. | SQLServer or Access | SQLServer |
| 3 | The network address of the SQL Server (IP address of the machine that the SQL Server is running on). | localhost w.x.y.z (where w,x,y,z are in range 0-255) | localhost 192.168.229.81 |
| 4 | The name of the Witt database on the SQL Server. | Printable <u>ASCII</u> characters that are supported by the Operating System | ProSolvWittDB |
| 5 | The username associated with the Witt database that allows the application to connect to the database directly. | Printable <u>ASCII</u> characters that are supported by the Operating System | prosolv sa |
| 6 | The password associated with the Witt database and the user specified in line 5. | Printable <u>ASCII</u> characters that are supported by the Operating System | psci |


Access

| Line number | Explanation | Possible Values | Example |
|-------------|--|---|--|
| 1 | The full path to the folder that contains the Witt files that are processed by this application. | Windows formatted file path ending with a '*' symbol to indicate all files in that folder (directory path should be valid and existing) | \\<AppServer>\From_Witt* |
| 2 | The type of Witt database where retrieved file data is stored. | SQLServer or Access | Access |
| 3 | The full path to the Witt .mdb file. | Windows formatted file path | \\<AppServer>\ProSolv_Database\Witt.mdb |

| | | |
|---|--|----------------------------|
|  | ProSolv® CardioVascular 4.0.2 Server Administrator's Manual | PMAN-0014756-A |
| | | Effective Date: 02/06/2009 |
| | | Page 82 of 83 |

| Line number | Explanation | Possible Values | Example |
|-------------|--|--|--|
| 4 | The full path to the images folder of the ProSolv® database. | Printable <u>ASCII</u> characters that are supported by the Operating System | \\<AppServer>\ProSolv Database\images |

- Open the text file, rt2f.ini, from the installation folder in a text editor.



- Edit the file to contain the information listed in the tables below.

| Line number | Explanation | Possible Values | Example |
|-------------|---|------------------------------------|--|
| 1 | Port number for HL7RealTime2File.exe application. | 1-65535 | 3336 |
| 2 | The full path to the folder where Witt files should be stored for processing. | SQL Server or Access | \\<AppServer>\From_Witt* |

- Launch HL7RealTime2File.exe. This application reads the HL7 data stream and writes out data files.
- Launch Witt.exe. This application processes the data files. **Note** that both applications must be running for the Witt interface.

