

**" Annexe  $P_E$  max -  $P_I$  max "**  
**SPIROMETER**  
**D A T O S P I R - 120**  
**511-8D0-MU2**  
**\* USER'S MANUAL \***

**2012-05 • Rev.: 2.01**

**SPIROMETER DATOSPIR-120**  
**USE MANUAL "Annexe P<sub>E</sub>max - P<sub>I</sub>max"**  
**E. &\$%**

**1. INSTRUCTIONS FOR USE AND INSTALLATION**

**1.1. INTRODUCTION**

**1.2. PROGRAM OF MAXIMAL PRESSURES**

**1.3. SETUP OF MAXIMAL PRESSURES**

**1.4. PROCEDURE FOR MAXIMAL PRESSURES TEST**

**1.4.1. INPUT OF TEST DATA**

**1.4.2. PERFORMANCE OF MAXIMAL PRESSURE TESTS**

**1.4.3. DISPLAY OF THE RESULTS**

**1.4.4. PRINTING AND/OR SAVING IN THE DATA BASE**

**1.4.5. TRANSFER OF TESTS**

**2. TECHNICAL SPECIFICATIONS**

**2.1. PARAMETERS**

**2.2. RANGES AND MEASUREMENTS**

**2.3. CONTROL**

**3. INTRODUCTION TO FUNCTIONING**

**3.1. SIGNALS ACQUISITION**

**3.2. CONTROL PROGRAM**


**4. TECHNIQUE OF MEASUREMENT AND REFERENCE VALUES**

**5. PREVENTIVE AND CORRECTIVE MAINTENANCE**

**6. MODIFICATIONS**

**Revised**  
**Date:** 2012-05

**Approved**  
**Date:** 2012-05



**José María Plana**  
**Technical Manager**



**Carlos Recio**  
**Sales Manager**



IN ACCORDANCE WITH  
93/42/EEC Medical Devices Directive  
II a Class

**SIBEL S.A.**

Rosselló, 500  
08026 Barcelona  
Spain

Tel. : +34 93 436 00 08 Fax : +34 93 436 16 11  
e-mail : [sibel@sibelmed.com](mailto:sibel@sibelmed.com) WEB : [www.sibelmed.com](http://www.sibelmed.com)

**AFTERSALES SERVICE**

Tel. : +34 93 433 54 50 FAX : +34 93 436 16 11  
e-mail : [sat@sibelmed.com](mailto:sat@sibelmed.com)

# 1. INSTRUCTIONS FOR USE AND INSTALLATION

## 1.1. INTRODUCTION

This manual is an annexe to the Spirometer **DATOSPIR-120** User's Manual, exclusively dedicated to the functioning of the maximal pressures. For consults and observations referred to the device as a whole, see chapter 1 of the general manual.

The maximal pressures module has been designed in collaboration with the "Pulmonary Function Laboratory" of the Hospital de la Santa Creu i Sant Pau of Barcelona. It is based on the criteria expressed by J.L. Clausen in the Thoracic Society of California.

It allows a measurement range of 300 cmH<sub>2</sub>O in both inspiratory and expiratory tests, and has several reference available values to be setup by the user.

## 1.2. PROGRAM OF MAXIMAL PRESSURES

Next the different options of the maximal pressures program are shown:

### Data of the test

- Patient
- Code
- Age, height, weight and sex
- Name and surname

### Change of patient

### Start of the manoeuvre

### Change from $P_E$ max to $P_I$ max or vice versa

### Display of the test data

- Display of all the manoeuvres
- Selection of the curve to be printed or memorized
- Save in the Data Base
- Performance of the Report

### Performance of the Report


### Display and selection of the manoeuvres

- Display of the manoeuvre in progress
- Display of the 3 best manoeuvres.
- Display of one manoeuvre


## 1.3. SETUP OF THE MAXIMAL PRESSURES

The **Maximal Pressures** module is an option included in the **DATOSPIR-120**. It is advisable for the user to set it up according to his/her needs.

To enter the Maximal Pressures setup menu, press key F2 in the main menu and next, press F8.




- F1 ESC, escapes this screen and moves back to the previous one
- F2 Setup of references
- F3 Selection of graphics and time lag



- F1 ESC, escapes this screen and moves back to the previous one
- F2 Selects references for adults
- F3 Selects references for children
- F5 & F6 Moves the cursor back or forwards.
- F7 Validates the input data and goes to the next screen.

Not all the references have values for adults and for children; the program will not allow to select non existent references.



- F1 ESC, escapes this screen and moves back to the previous one.
- F5 y F6 Moves the cursor back or forward .
- F7 Activates or deactivates the option.
- F8 Validates the input data and go to the next screen.

The time lag for the measurement has a **default value of 1.0 s** ,although the user can configurate it between 0.1s and 4.9 s.

This time lag has an effect on the calculation of the measurement. Therefore no value in this first second ( or time setup by the user) is taken in account.

The options for graphics allow to activate or deactivate the options for saving the graphics in the data base and print the curves of  $P_{E\max}$  or  $P_{I\max}$ .

## 1.4. PROCEDURE FOR MAXIMAL PRESSURE TEST

To have access to the Maximal Pressure test, press key F7 from the main menu.

### 1.4.1. INPUT OF THE TEST DATA

When entering the Maximal Pressures Program, first of all enter the test data.

The screenshot shows a terminal window titled "MAXIMAL RESPIRATORY PRESSURES". Inside, there is a sub-window titled "TEST DATA" with the following fields and values:

Code:	1234567890		
A(y):	43	Sex:	♂
H(cm):	176	W(Kg):	92
Name:	FRANCISCO A.		
Surn.:	SOPENA IBAÑEZ		

At the bottom of the screen, the date and time are displayed as "11.04.2000 12:58".

Navigation keys are indicated by icons and labels:

- F1: ESC key
- F2: Left arrow key
- F3: Erase key
- F4: Right arrow key
- F5: Up arrow key
- F6: Down arrow key
- F7: Left arrow key
- F8: Right arrow key


See section 1.8.1 of the general manual, for the input of data.

If a performed test is still in memory, the program will go directly to the next screen (section 1.4.2).

### 1.4.2. PERFORMANCE OF MAXIMAL PRESSURE TESTS

After validating the data with F7, enter the tests performance screen.

In this screen the axis appear, and then the **DATOSPIR 120** is ready to perform maximal pressures tests.



- F1 ESC, escapes this screen and goes back to the previous one
- F4 Enables to modify the Patient data
- F5 Starts the manoeuvre
- F6 Changes the test: from P<sub>E</sub>max to P<sub>I</sub>max or vice versa
- (F7) Has access to the Data screen (If there are performed manoeuvres)
- (F8) Performs the Report (If there are performed manoeuvres)


**It is convenient for the technician to know the usual procedure required for the patient to perform the test correctly. If not, it is advisable to see some bibliography on this ( See also chapter 4).**

In the performance of the test, bear in mind the following steps:


1 Connect the module in the connector no. 9. The device detects automatically that the module is connected. If the module is not connected, the device will show it on screen and it will not allow to enter the screen of tests performance. Make sure that the shutter Valve is connected to the module.

2 Train the patient about the test performance, as his/her collaboration is essential for the correct execution. Put him/her the nose clip.

3 Press key F5 and wait until an intermittent arrow appears on the screen.




The shutter valve must be in open position, in order to allow the patient to breath normally.



Once the patient has reached the position of Total Pulmonary Capacity, move the cursor to the closed position and perform the manoeuvre.

4 Once the manoeuvre is finished, repeat step 3 in order to perform more manoeuvres.

It is advisable a minimum of 3 manoeuvres. The best 3 manoeuvres should not differ more than 5 % between them, and the last manoeuvre should not be the best.



5 The above screen presents (or can present ) the following information:

- Date, patient code and type of test.
- Warnings:
  - >5%: indicates that the variation between the best 3 manoeuvres is over 5%.
- Pressure/time Graphic for the best 3 manoeuvres.
- Pressure/time Graphic for the selected manoeuvre (pressing one of the keys M1 to M5).
- Horizontal line of dots showing the reference value
- Vertical line of dots showing the maximal value for each manoeuvre.
- Values of maximal pressure of the performed manoeuvres, average and standard deviation of the best 3 ones.

The manoeuvres of maximal pressure are ordered chronologically according to their performance. This is made because it is important to follow the evolution of the patient manoeuvres.

Anyway, although the manoeuvres are saved in the memory in temporal order, the best 3 manoeuvres are displayed on screen ordered from the best to the worst (the manoeuvre of highest pressure value is considered as the best ).

It is also important to point out that more than 5 manoeuvres may be performed, despite the fact that the device has space only for 5 manoeuvres.


When performing the sixth manoeuvre, the first performed will be deleted (M1), if this is not the best. Should it be the best manoeuvre, the second one would be deleted (M2). Once the corresponding manoeuvre has been deleted, the system will reorder them,

being the performed manoeuvre the M5, the M5 will be the M4, the M4 will be M3, and so on.

**NOTA: Remember that with the backspace key F1 "ESC" is possible to move back in the menu without losing the information available up to that moment. Only if the patient is changed by entering a new code, or performing other operation, the data might be lost. Anyway, the device will always warn with a message on screen before deleting the information.**

### 1.4.3 DISPLAY OF THE RESULTS

Press key F7.



- When entering this option, the maximal pressure values are displayed, for the expiratory and inspiratory tests, as well as for the average and standard deviation of the best three manoeuvres.
- The reference value is also displayed, as well as the percentage with regard to the reference of:
  - \* the selected manoeuvre
  - \* the average of the best 3 manoeuvres
- The selected inspiratory and expiratory manoeuvres are displayed
- The selected references in the setup are displayed
- The Icons give access to:
  - F2: Deletes the selected manoeuvre of the selected test
  - F3: Changes from  $P_{E\max}$  to  $P_{I\max}$  and vice versa
  - F7: Prints the report of Maximal Pressures
  - F8: Saves the test in the Data Base

**WARNING:**

The system displays, saves or prints the values for all the manoeuvres, but only one curve for each manoeuvre is saved or printed in the data base.

Default, the system selects the curve of the best manoeuvre ( unlike the spirometry it may not be number 1 ). If the specialist prefers, he/she can select another curve, as it is explained next.

To select a manoeuvre different to the one selected default by the system, press the corresponding memory key (M1 to M5).

The upper right part of the screen indicates which test is selected, the expiratory (TEST PEmax) or the inspiratory (TEST PImax).

To change from PEmax to PImax and vice versa, press key F3.

**1.4.4. PRINTING AND/OR SAVING IN THE DATA BASE**

After performing the wanted manoeuvres, of both expiratory and inspiratory pressure, the following operations can be performed:

**A.- Printing the results**

A report displaying the same information of the data screen will be made, along with the patients data and the curves of the selected manoeuvres.

DATOSPIR 120 SIBEL S.A.  
DOS DE MAIG 290 08025 BARCELONA

**DATOSPIR 120 - SIBELMED**


Code: 1234567890 Date: 27.04.2000  
Name: FRANCISCO A.  
SOPENA IBANEZ  
Sex: Male Age(y): 43  
Hght.(cm): 176 Wght(Kg): 92  
Motive:  
Precedence:  
Techn. :  
References: P MORALES

REPORT P<sub>Imax</sub> - P<sub>E<sub>max</sub></sub>

	P <sub>Imax</sub>	P <sub>E<sub>max</sub></sub>
M1	144	171
M2	145	164
M3	143	153
M4	0	165
M5	0	157
Mean	143	167
Desv. Std.	0	0
Reference	143	206
(%)Manoeuvre	101	83
(%)Mean	100	81

Comments: .....

.....



## B.- Saving in the internal data base

The test is saved in the device internal data base, in order to display it afterwards, print it and/or transfer it to a computer.

### 1.4.5. TRANSFER OF TESTS

Like the spirometric tests, the maximal pressure tests saved in the data base can be transferred to a computer (See section 1.16.1 of the general manual).

## 2. TECHNICAL SPECIFICATIONS

### 2.1 PARAMETERS

In both expiratory and inspiratory test, the next values are measured:

- \* Maximal pressure of the 5 manoeuvres.
- \* Average of the best 3.
- \* Standard deviation of the best 3

### 2.2 RANGES AND MEASUREMENTS

Range of Measurement (cmH <sub>2</sub> O)	0 - 300
Resolution (cmH <sub>2</sub> O)	1
Accuracy	3 %
Sample Rate (Hz)	100

### 2.3 CONTROL

- *Number of manoeuvres.*

Up to 5 manoeuvres for each type can be performed (P<sub>E</sub>max and P<sub>I</sub>max)

- *Duration of the manoeuvres.*

The maximum duration is 8 seconds

- *Start of Manoeuvre.*

When the threshold of 3 cmH<sub>2</sub>O is surpassed

- *End of Manoeuvre.*

Variation lower than 1 cmH<sub>2</sub>O in the last 2 seconds.

- *Delay in the calculation of the maximum pressure value.*

Selectable in the setup between 0.1 and 4.9 seconds. 1 second as a default.

- *Internal Data Base*

The Spirometer shares the internal data base for all test types.

If only Maximal Pressure tests are saved, counting a duration of 4 seconds for the curves, it is possible to save:

	With Graphics		Without Graphics
	1 curve	2 curves	
<b>Data Base "L"</b>	150	150	150
<b>Data Base "H"</b>	400	200	1100

In this case, **Data Base "H"** for spirometry is reduced taking into account the next table (compare with table of 1.14):

<b>Data Base "H"</b>	With Graphics		Without Graphics
	Fleisch/Des	Turbina	
	300	540	900

### 3. INTRODUCTION TO FUNCTIONING

The **DATOSPIR-120** optionally has an external Module of Maximal Respiratory Pressures. This module is connected to the **DATOSPIR** through a phone connector.

#### 3.1 SIGNAL ACQUISITION

The module of maximal pressures includes basically a pressure sensor, a differential amplifier and a filter.

The pressure sensor covers the range from 0 to 300 cmH<sub>2</sub>O. The output of this sensor is amplified and filtered so that the analogical signal has the intended range and bandwidth.

Then, this signal is transformed into digital with a A/D converter. The converter used is the same as the one used in the spirometry. (See section 3.7 of the general manual).

#### 3.2 CONTROL PROGRAM

The control program performs the processing of the acquired signal. Among the important calculations, we must point out :

-Start and end of the manoeuvre:

The start of a manoeuvre is considered when the level of 3 cmH<sub>2</sub>O is surpassed. The end of the manoeuvre is considered when there is a pressure variation lower than 1 cmH<sub>2</sub>O for 2 s.

-Delay for the calculation:

For the calculation of the maximal pressure, to avoid artefacts, a initial period of time is discarded. This period is set to 1 s as a default but it can be set up between 0.1 and 4.9 s.

-Order of the manoeuvres:

The manoeuvres are saved in the memory in the same order as they have been performed.

-Automatic deleting of a remaining manoeuvre:

When the 5 memories are full, and a new manoeuvre is performed, the first manoeuvre is always deleted, if this is not the best. If the first is the best, then the second one is deleted.

#### 4. TECHNIQUE OF MEASUREMENT AND REFERENCE VALUES

The maximal inspiratory pressure ( $P_{I\max}$ ) depends directly on the strength developed by the inspiratory muscles. Its measurement is very useful in the diagnosis and follow up of patients with neuro-muscular illnesses, specific alterations of the respiratory muscles, different processes of the thoracic area, air trapping or modifications produced by breath depressive medicines, among other clinic situations. The maximal expiratory pressure ( $P_{E\max}$ ) also informs about the cough efficacy and the drainage of the bronchial secretions. In both cases, they are measures easy to perform and comfortable for the patient. They can be obtained near the patient and allow to follow his/her evolution.

The procedure to obtain the  $P_{I\max}$  and the  $P_{E\max}$  must complete the following aspects:

- Correct collaboration of the patient to get a maximum effort in the manoeuvres. The patient must keep his/her hands against the cheeks and leave them flaccid to attenuate the pressure generated by the facial muscles.
- To measure the  $P_{I\max}$ , the patient must perform a deep inspiration from the position of Residual Volume and maintain the effort for at least 3 seconds.
- To measure the  $P_{E\max}$ , the patient has to inspire previously up to the position of Total Pulmonary Capacity and from there blow with the maximum effort.
- It is advisable a minimum of three satisfactory manoeuvres for each pressure.

If measuring errors are suspected, we must make sure that there are no leaks in the system and the correct collaboration of the patient in the maximum effort manoeuvres.

The Module of Maximal Pressures is useful for Pneumology, Anaesthesiology, ICU and Neurology units, for isolated determinations or for the follow up of interned patients.

The device includes 4 different references, to be selected by the user. Next, the bibliographic reference is detailed.

##### A-) P. MORALES

*Presiones respiratorias estáticas máximas en adultos. Valores de referencia de una población caucásica mediterránea.*

P. Morales, J. Sanchís, P.J. Cordero y J.L. Díez.

ARCHIVOS DE BRONCONEUMOLOGÍA. Vol 33, num 5, 1997.

##### B-) SH WILSON

*Predicted normal values for maximal respiratory pressures in caucasian adults and children*

SH Wilson, NT Cooke, RHT Edwards, SG Spiro.  
THORAX 1984; 39:535-538.

**C-) BLACK & HYATT**

*Maximal Respiratory Pressures: Normal values and relationship to age and sex*  
Black LF, Hyatt RE  
AM REV RESPIR DIS. 99:696-702, 1969

**D-) MR CHARFI**

*Les pressions ventilatoires maximales à la bouche chez l'adulte: valeurs normales et variables explicatives*  
MR Charfi, R Matran, J Regnard, MO Richard, J Champeau, J Dall'ava, A Lockhart  
REV MAL RESP, 1991, 8, 367-374.

## **5. MAINTENANCE, PREVENTIVE AND CORRECTIVE MAINTENANCE**

### **5.1 MAINTENANCE**

The maintenance is the action directed to keep the device in correct functioning. The person in charge does not require any technical qualification, except the knowledge of the functions and handling of the device. Usually the user of the device must do it. The operations to be performed are the following:

#### **5.1.1 MODULE**

The manometer is cleaned softly with a dry or slightly wet with water cloth, drying the rests of humidity. Take care that no liquid goes into the inner space or in the connectors or connections.

Do not use abrasive or solvent substances.

#### **5.1.2 CLEANING OF THE SHUTTER PROBE**

The shutter probe, and specially the reusable mouthpiece, must be washed with soap or some disinfectant liquid which does not chemically attack the methacrylate or Teflon of the composition.

After the rinse and drying of the probe, lubricate the friction space between the probe and the cursor with Vaseline. Finally check that the components are not obstructed and that it works correctly.

### **5.2 PREVENTIVE MAINTENANCE**

The preventive maintenance consists of all the actions directed to keep the device in a good status of use.

There are two types of preventive maintenance:

1 The first type, that can be performed by the user, consists of a periodical supervision of the aspect of the different interconnections and the other external elements of the device. It will be verified that all the interconnections are perfectly connected, that the cables and/or connectors and the rest of elements do not have breakage or external damage.

In case of detecting any anomaly which the user cannot solve, inform the aftersales service of SIBEL S.A. or your Distributor, to proceed with the checking or repair.

2 The second type consist of a general technical verification of the safety systems,

adjustments, functions, etc, which configurate the device.

**THESE CHECK UPS WILL BE PERFORMED WITH ANNUAL PERIODICITY** according to the Verification and Adjustment Procedure of the MAXIMAL PRESSURES MODULE, available by the manufacturer. This type of operation must be performed by the qualified technical staff of the maintenance department or the technical service in the distributor or manufacturer.

In any case, SIBEL S.A. as manufacturer, must give a written authorization to the corresponding technical service , at least during the warranty period , so that they are able to perform such maintenance. Sibel is not responsible for any damage, malfunction etc. which can come as a consequence of a faulty maintenance made by persons not belonging to SIBEL S.A.

### 5.3 CORRECTIVE MAINTENANCE

The corrective maintenance consists of leaving the device in a good state of use, which has failed in its correct functioning by bad functioning or misuse and has to be repaired.

In the case of detecting a failure in the device, which prevents its normal utilization, disconnect the device from the mains and contact with the Aftersales Service of SIBEL S.A., specifying with detail the type of anomaly produced.

**6. MODIFICATIONS**