

Cryostat Microtome RD-2298

Operation Manual

UnRegistered

Helpful Suggestions:

- The User manual should be considered part of the instrument. The user should read this user manual carefully before first installation and before operating the machine, and place the manual close to the machine for reference.
- The cryostat microtome must be operated by a trained operator. Before operation, please carefully read the operating instructions, and familiarize yourself with the steps and processes.
- Unit must be shipped in the upright position. If the unit was shipped lying down or stored lying down then it MUST stand upright for 36 hours before use for all lubricating oils to enter the compressor
- Do not put instrument under extreme temperatures or a high humidity environment. Failure to follow this procedure will cause severe damage
- Do not place instrument directly under sunshine or near air conditioning vents
- No other obstruction should be around the instrument within a 30cm distance.
- Check that your power outlet is well grounded. Ensure that power supply is stable, constant and adequate. Compressor requires a start-up current between 45 and 50A.. Therefore, a professional electrical engineer should inspect the electrical circuit before installation in order to ensure the smooth operation of the instrument. Requirements:
 - The instrument should be switched on 24 hours every day.
 - Blade is sharp, exercise extreme caution when changing it, remove when cleaning the inside of the chamber
 - Do not place blade anywhere with the cutting edge facing upwards.
 - Before changing specimens, always lock the hand wheel and cover the blade edge with the blade guard.
 - In case of malfunction, contact a qualified service technician. Don't try to solve problem by yourself.

1. Introduction

The Semi-automatic Cryostat Microtome, is a national regulated 1st class medical device, is used for rapid freezing of pathological sections of human and animal body tissue. It can be widely used for pathological diagnosis, analysis and research in hospitals, medical colleges, and by legal medical experts and propagation institutes.

The cryostat microtome is the eighth generation of our cryostat. It utilizes advanced horizontal feeding, The specimen chamber area is much larger than previous models.

The instrument consists of 4 major parts:

- 1.1 The computer control panel is on the upper part. This area displays the temperature and the general working conditions of the instrument.
- 1.2 The center part is the constant low temperature chamber used to rapidly freeze the tissue and to conduct the cutting.
- 1.3 Behind the chamber is the mechanical transmission and motor drive.
- 1.4 The lower part is the compressor unit , cooling condenser and evaporator.

2. Scope of Application

Pathological sectioning in order to study the histology of plant, animal or human tissue.

3. Technical Parameters

Section thickness range: 1-90 μ m Adjustable , Increment: 1 μ m, Increment: \pm 20%。

Trimming thickness range: 10 ~ 400 μ m Adjustable , Increment: 10 μ m , Increment: \pm 20%。

Specimen retraction: 20 μ m

Horizontal stroke: 25mm

Vertical stroke: 59mm

Specimen orientation: 12° (X、Y、Z three directions)

Voltage: 220V \pm 10%

Frequency: 50Hz

Power: 650W

Largest start-up current (5 second) 50A

Chamber temperature: -10 $^{\circ}$ C ~ -35 $^{\circ}$ C Adjustable

Freezing shelf temperature: \leq -50 $^{\circ}$ C

Number of freezing station: 18

Peltier number: 2

Refrigerant: R404a,300g \pm 10g

Compressor oil: 0.6L EMKARATE RL-22S, ICI

Dimensions:

Width (No hand wheel) 650mm

Width (With hand wheel) 790mm

Depth: 700mm

Height: 1150mm

Weight: 140kg

4. Working Condition and Pre-adjustment

4.1 The instrument is a movable floor type with 4 caster wheels. Two front wheels can be up and down and are adjusted by turning knobs on the foot of the wheel. These front wheels only function when the instrument needs to move. After removing the packaging material, put the front wheels on the ground and move instrument to the working position and secure it, then disable the front wheels.

4.2 The instrument should only be used after 4 hours of standing still to ensure proper lubricant distribution.

4.3 Confirm the power supply is grounded (triple socket), then connect the instrument with the main power socket, turn on the power switch at the back of the instrument, the machine is now in operation mode.

4.4 Areas on both sides of the instrument should be clear to allow proper air flow to the cooling unit. There should be a minimum of 300mm space to facilitate ventilation and heat dispersion.

4.5 The instrument uses a single phase AC power supply, the voltage is 110-115V±11V.

4.6 The instrument working ambient temperature is +5°C~+28°C. If the temperature is above 26°C, air conditioning must be installed; otherwise the life time of the cooling compressor will be shortened.

4.7 Environment relative humidity should be no more than 80%.



High temperature or humidity will influence freezing effect and shorten compressor life.

5. Installation

5.1 Freezing chamber assembly

- Place left and right shelf in the freezing chamber outside of microtome;
- Place the waste tray in the freezing chamber;
- Install blade holder on the microtome base and lock it.
- Insert blade to precool;
- Place all tools needed for specimen preparation in the chamber;
- Close glass window to enclose chamber.


5.2 Connect Power



6. Operation

6.1 General overview

Check your power socket is well grounded. Ensure that power supply is stable, constant and adequate. Compressor requires a start-up current between 45 and 50A., Therefore, a professional electrical engineer must inspect the electric circuit before installation to ensure the smooth operation of the instrument.






Fig 2 Main control panel

Fig 3 Left side control panel




Fig 4 Right side control panel

6.2 Turning the Instrument On



Fig 5

Confirm the power supply is grounded (triple socket), then use attached power line to connect the instrument with the main power socket, turn on the power switch on the right of the instrument, the machine is in operation mode. Main power is shown as Fig 5 on the right of the instrument. The switch must be in the top position for "on" and in the bottom position for switching "off".

Turn on the power switch on the right side of the instrument, the machine is now in operation mode, the sectioning arm automatic reset. LED on the left side control panel

should now display “06” μm , that indicates default section thickness is “06” μm and present work mode is “sectioning” mode. LED on the right side control panel displays “----” to waiting for system reset. After reset buzzer, LED on the right control panel displays “0” to show total sectioning number.

6.3 Function of main control panel (Please refer to Fig 2)

Main control panel is separated into 4 parts. From left to right in sequential order are Freezing shelf area, Freezing chamber area, Time setting area and Special function area.

6.4 Peltier element





Fig 6 Peltier element


Under a normal state, the peltier element is off. LED shows freezing shelf temperature. Press  button , indicator light comes on and peltier element begins to work to cool down the freezing stations. LED displays “L15”, which means peltier element will work 15 minutes. The countdown of the remaining cooling time is now permanently displayed.


The Peltier element turns off automatically after 10 minutes.

Once the remaining cooling time displayed is 4 minutes, the figure 4 is followed by a point (“ 4. ”). At this stage the Peltier element may be deactivated by pressing peltier button again. Once deactivated, the display indication returns to ‘PE’. pressing the button during the peltier element working process will stop peltier freezing immediately.

6.5 Temperature setting

Compressor will start to work after turning on power switch allowing the freezing chamber to cool down. (depending on the interval between turning off and on, if switched within 3 minutes, the instrument requires 3 minutes for compressor to start again). Freezing chamber temperature is shown and set by temperature set area as Fig 6. Normally LED displays temperature of sensor in the freezing chamber.



Freezing chamber temperature display and setting: System default temperature is -22°C . When the instrument is under idle state, press “” to increase



temperature; press “ ”to decrease temperature.
 At present, the unit’s LED is flashing to display set value. When set to desired temperature, wait about 15 seconds, the instrument go back to normal state.

Fig 6 Freezing chamber temperature User can adjust temperature within the adjustable range.

6.6 Time setting




Fig 7 Time setting area

Displays buttons on the panel stand for:

- SUN: Sunday MON: Monday
- TUE: Tuesday WED: Wednesday
- THU: Thursday FRI: Friday
- SAT: Saturday
- AM: Morning PM: Afternoon
- W-DAY: Working day H: Hour
- M: Minute



6.6.1 Present time setting

After turning on system default time, if it is not correct, then reset it. When system is under idle state, press “D” button to adjust day, light will flash correspondingly. Press “H” to adjust hour, the left two digit place of LED will flash; Press “M” to adjust minute, the right two digit place of LED will flash. About 10 seconds after setting, the instrument will return normal state



Note: if setting time is over 10 seconds between two operations, system will automatically return to normal state. Below setting is the same.




6.6.2. Programming working time

When under idle state, Press  button once,  symbol lights and is ready to enter programming mode. Left side LED displays “F1.1”; Front F1 stands for first mode; the other 1 stands for switch on mode. The lights from Mon to Fri all light up concurrently. That means the instrument is on the working temperature From Monday to Friday; instrument goes into stand by mode on the other time. Under this condition, if “D” button is pressed, it does not function. Hour button “H” and minute button “M” are used to adjust hour and minute respectively. When in programming



mode, LED right digit will flash; Then press “ ” again to go into “switch off” mode. Left side LED displays “F1.0”; F1 this is the same function as above F1; 0 stands for switch off mode. Switch off programming method is same as above. 10 seconds after finishing programming, the instrument returns working state.


6.6.3 Programming the working time for each day

Press “  “ button three times continuously,  indicator lights up and unit goes into time programming. Left side LED displays “F2-1”; first 2 stands for second time programming mode, 1 at the end stands for switch on mode. User may program different time for every day within one week. If sigle light is on among Monday to Sunday, user may program that day working time. Press “ D” to change among days. Hour button “H” and minute button “M” are to adjust hour and minute respectively. When under programming, LED flash with right digit; user then needs to press “  “ again to go into switch off mode. Left side LED displays “F2-0”; front 2 same function as above; 0 stands for switch off mode. Switch off programming method is same as above. 10 seconds after finish programming, the instrument returns working state.




Note: If Sunday is out of work, program switch on and off time as same time.

6.6.4 Programming the defrost cycle

Press “  ” button that is in time control area. Defrost indicator light is on and enter defrost programming. Indicator lights of week on the top line go out currently. It represents that program same defrost time every day. ‘W’ button is not available. Hour button “H” and minute button “M” are to adjust hour and minute respectively. Programming time is defrost beginning time. About 10 seconds after programming, system return idle state.

6.7 Function buttons



There are fast defrost, stand by, UV sterilization lamp and illumination buttons in function area as shown in Fig 8.

Fig 8


6.7.1 Locking key




When user press locking key, time setting area is locked and could not be changed unless press it again.

6.7.2



Fast defrost button. Under normal freezing state, it is in unlit. Press “” to activate fast defrost function and indicator light is on, system goes into fast defrost mode. Press button again to cancel fast defrost, indicator light goes out simultaneously. Unit should then defrost withing a 6 minute time period.



If system automatically enter fast defrost mode,  button lights up too. Press it again to stop fast defrost and indicator light goes out.

6.7.3 UV sterilize lamp



is UV sterilization button. Press it to make UV lamp active; press once more to stop it. Indicator light is on or off correspondingly. It is better to start to sterilize with UV lamp under stand by situation.



Note: UV radiation is harmful for human body. When doing UV sterilize, do not forget to close glass window and keep hands away from chamber.




Note: Before doing UV sterilization, please remove all specimens and clean chamber first.

6.7.4 Stand by



is stand by button. Under normal state, it is dark. Press it one time, system goes into stand by mode. Press it again to return normal working state. Indicator light is on or off correspondingly. Under stand by state, freezing chamber temperature is maintained between $-4^{\circ}\text{C} \sim -9^{\circ}\text{C}$, to conserve electricity.



If system automatically enter stand by mode,  button lights up too. Press it again to stop stand by function and indicator light goes out.

6.7.5 Illumination



Illumination button. It is dual function button. First press lamp is on; second press lamp is off.

6.8 Use of the liquid waste container

Use the wasted liquid barrel to collect defrosted liquid, tissue thawing solution, flushing liquid and waste objects. In order to prevent human and environmental hazard from the waste, put 200ml 10% formalin or other sterile solution into the barrel in advance. The disposal of the waste shall follow all the related regulations of the hospital.





Note: Use the right and effective method to treat the waste.

6.9 Left side control panel (Refer to Fig 3)


6.9.1 Specimen forward and backward feed

Specimen forward and backward feed is for user fast move specimen position. If

user press  button, specimen goes to home position in normal speed; press  button, specimen

goes forward in normal speed. Press  button, specimen moves to home position in fast

speed; press  button moves to opposite direction in fast speed. Press

 button to cancel moving command.

6.9.2 Section/Trim thickness selection



Fig 9 is shows section/trim control. Above LED displays the thickness value of section/trim.



If "SECT" lights up, LED displays section thickness, otherwise, is trim thickness. Press "T/S" button per to change from section to trim one time, or conversely.

Fig 9 Section/Trim selection

and display

Section thickness range: 1-90 μ m (default value is 6 μ m)

- 1~20 μm , in 1 μm steps;
- 20~40 μm , in 2 μm steps;
- 20~60 μm , in 5 μm steps;

Trim thickness range: 10~400 μm (default value is 30 μm)

- 10~50 μm , in 5 μm steps;
- 50~100 μm , in 10 μm steps;
- 100~400 μm , in 50 μm steps;

6.9.3 Thickness selection

Under current working mode, adjust “” or “” to select section or trim thickness.

6.9.4 Section

Turn the hand wheel one rotation, the system will cut the one piece of section one time, and the thickness is the preset thickness.

Note: When adjusting specimen quick fro and back movements, the hand wheel must lock at the right position. Refer to Fig. 18.

two control sections.

6.10 Right side control panel (Refer to Fig 3)

6.10.1 Counter




Σn
 $\Sigma \mu\text{m}$



As left Fig 10, above LED displays number of total pieces and thickness. When Σn lights up, total pieces is shown in LED; on the contrast, $\Sigma \mu\text{m}$ lights up, total thickness is shown. Press “M” button to change from total pieces to total thickness one time, or conversely. Press “C” to clear all the numbers to 0.

Fig 10



6.10.2 Use of vacuum extractor

Vacuum extractor is designed for cleaning tissue

scrap. Press button on the right control panel to enable extractor motor. Indicator light on the upper corner of symbol is on. Move speed controller on the left side to adjust motor speed to choose suitable air volume. Press again to disable it; indicator light is off simultaneously.

7. Mechanical introduction

7.1 Freezing shelf




Fig 12 Freezing shelf description


Freezing shelf with 20 stations is located to the left of freezing chamber. The lowest temperature of freezing shelf may arrive to -40°C . Two project stations in the front of freezing shelf are peltier stations that the lowest temperature may arrive to -60°C .

Procedure of preparing specimen is as follows:

- Cut specimen into proper size;
- Put specimen onto freezing disc between layers of tissue freezing medium;
- Put freezing disc into station on the freezing shelf;
- Use hammer plat the specimen surface;
- The specimen is frozen by low temperature quickly;
- Move away hammer after it is frozen ;
- Insert freezing disc with specimen into hole of clamp on the microtome, then tighten locking handle;
- Adjust angle of knife holder and blade;

- Rotate handwheel to cut specimen into sections.

7.2 Microtome component introduction



1. Freezing disc locking handle
2. Orientation adjustment handle
3. Orientation lever
4. Extractor hole
5. Knife holder locking block
6. Specimen hole (Inserted freezing disc).

Fig 13 Microtome parts introduction

7.3 Inserting freezing disc



Note: Before inserting freezing disc, please lock handwheel well

- Release freezing disc locking handle;
- Insert the shaft of the freezing disc in the location hole of the specimen clamp. Make sure that the shaft of the freezing disc is fully inserted. The entire rear surface of the prism should contact well with the specimen clamp, there should be no gap between them;
- Lock freezing disc locking handle and fix it,
- Observe cutting angle, if necessary, adjust orientation;
- Procedure for adjusting orientation:
 - Loosen orientation adjustable handle;
 - Adjust orientation lever to choose best cutting angle;
 - Relock orientation adjustable handle.

7.4 Blade holder position adjustment




Fig 14 Blade holder introduction

Release (anti clockwise) blade holder lock handle on the left (④ Fig .13), manually move the holder forward and back and lock it at the desired position. When the holder moves to the terminal point the whole blade holder can be taken out and put it back by reversing procedure. Release (counter clockwise) blade holder lock handle on the right, manually move the holder left and right and lock it at the desired position.

Note: Remember to lock the blade holder to right position when doing forward or backward specimen adjustment.

7.5 Change the blade

Release clamping lever ① on the upper right to open blade clearance; insert blade from left to right; then lock the clamp lever (Fig 13).



Note: Blade is sharp, beware when you change it.

7.6 Adjust the angle of the blade

Loosen the lowest adjustment handle③ on the right side of blade holder. Refer to reference scale to turn rotator in order to choose best cutting angle of blade. Retighten the adjustable handle.

7.7 Adjust the anti rolling plate




Fig 15 Making parallel between edge of blade and slide

The anti roll plate is made of Plexiglas.

Adjustment includes:

- 1、 Make the edge of the anti rolling plate parallel to the edge of the blade. (Fig 15)
 - Insert blade and lock the clamping lever(①)
 - Lock two knobs⑥.
 - Use hexagonal wrench to make blade and anti rolling plate parallel
 - Fasten hexagonal screw⑦
 - Turn height turning knob⑤ to adjust height of anti rolling plate



Note: Please use ⑧Handle to replace scrap collector , do not touch other knobs unless you want to adjust parallel.

2. Adjust the anti rolling plate angle with the blade. (Fig 16)




Fig 16



The angle shall be no more than 10°

3. Make the edge of plate and edge of blade flush (Fig 17)




Fig17 Making flush between edge of blade and anti rolling plate




1. The upper part of the anti-roll plate should not contact the tissue specimen.
2. When the quality of the tissue section is not good, check the edge of the blade, check any segment of tissue on the plate,

3. Avoid touching the front of the plate, because high temperature will adhere the tissue section to plate.

7.8 Section reference

1. When using the cryostat microtome, controlling the blade speed and correctly adjusting the anti-roll plate are important factors in cutting high quality sections. The proper speed of the blade can only be acquired by practice and experience, along with skillful control of the hand wheel. The proper adjustment of the anti-roll plate requires even more practice. Sometimes they interfere with each other; and require careful adjustment.
2. When freezing living tissue, the water within the tissue will condense into ice, the tissue will become hard. The hardness changes with the temperature, the lower temperature the harder the tissue becomes. Different tissue reacts differently under different temperatures, achieving high quality sections, can only be obtained through practice. cutting temperatures of the specimen tissue without fat and formalin, should be between $-13^{\circ}\text{C}\sim-23^{\circ}\text{C}$.
3. For high quality sectioning, pay attention to the following:
 - (1) Proper selection of working temperature in the chamber
 - (2) Correct slicing operation and hand wheel speed
 - (3) Fine adjustment of the anti-rolling plate
 - (4) Sharp blade edge, correct blade angle
4. Frozen sections are just the opposite of paraffin slicing, the frozen sections should not be cut in wide sections, better to cut less wide sections so that the line of contact is shorter.

7.9 Hand wheel lock



To achieve operator convenience and safety, we designed a handwheel locking function for the cryostat. Lock handwheel at the top position not in use to avoid unnecessary harm. If the locking position “3-8 o'clock” are at the top while gearing handle is facing to operator as showing in the Fig 18.

Fig 18 Hand wheel locker

7.10 Installing and Changing the Section Extraction



Please take care of this precision instrument carefully. It needs daily cleaning and regular maintenance to achieve a long life for the instrument.





Fig 19 Extractor assembly

Connect assemblies of extractor according to Fig 19. Assembly ④ in Fig 19 is extractor install hole on the instrument. Screw assembly ③ into assembly ④ to



connect with main instrument. Press symbol  to activate extractor.◦



Filter should be cleaned and changed regularly and be kept clear.
Other assemblies such as hose should be cleaned regularly as well.



If extractor is not being used, please use rubber cap to cover extractor hole to keep the freezing chamber enclosed to avoid moisture entering chamber.

8 Troubleshooting

Problems	Causes/Remedies
Tissue cracking	Select thinner tissue, The freezing time is not quick enough
Tissue specimen off	Apply more embedding material The cutting tissue is soft The temperature of the specimen disc is low
Tissue moves across blade, but without slicing	Ensure that the blade is fastened and seated well. Fix the specimen clamp securely Check that the the specimen disc is securely fixed on the clamp- Increase the blade angle
Tissue section rolling	Temperature is not low enough, chooser lower temperature. Blade is not sharp, change to new one. Check the edge of anti-rolling plate if broken, change to new one . Anti-rollplate is not clean. Clean with alchohol. Check the gap between the anti-roll plate and the blade, adjust the angle of blade Increase the height of the anti-rolling plate Check the angle of the anti-roll plate and the blade
Tissue becomes soft when cutting	Prolong the refrigeration time of the anti-roll plate and the

	blade
Tissue section stick to the anti-rolling plate	Use short hair brush to clean the anti-roll plate Prolong the refrigeration time
Tissue section overlapping	Clean the blade with camel hair brush The clearance angle between the blade and the anti-roll plate is too narrow
Vertical crack of tissue section	Edge of the blade has a defect, Foreign matter on the blade
Tissue section narrows	Adjust the clearance angle of anti-roll plate and the blade evenly
Tissue section broken	Too low temperature of the tissue, too long time freezing, Blade is not sharp, change to new one. Blade sticking with some impurity, clean and brush
Vibration when slicing	Unfixed blade The angle of the blade is too big or too small Specimen did not arrive at slicing temperature, please continue to freeze it. Specimen clamp is unfixed The specimen disc is not fixed well Cutting speed is too fast
Uneven thickness of the tissue section	Temperature is not proper Microtome did not dry completely Dull blade Angle is too small The specimen disc is not fixed
Rotate handwheel, but no feeding	Microtome inside is not dry completely- allow to dry Microtome itself damaged. Please contact with technical service. Defective Photoelectric sensor. Defective stepper motor Defective stepper motor driver

No display after power on	Check the power line, Reconfirm power to line. Fuse blown, Change fuse.
Can't reset after turning on instrument	Defective Stepper motor. Defective Stepper motor driver. Lead screw or nut is broken, Upper limitation Photoelectric sensor is defective
No time and temperature display	Connection cable between main board and display board has problem Display board is defective
Compressor doesn't operate after turning on instrument more than 4 minutes	Check the solid state relay and change if necessary. Compressor is damaged
No defrost	Defective Solenoid valve. Defective Solenoid valve power supply,
Lamp does not come on	Lamp is broken, Lamp wire has problem,
Chamber temperature stays same.	Defective Temperature sensor.
Freezing shelf temperature stays same.	Defective Temperature sensor.

chamber temperature can't reach to the programming temperature (under running condition)	<p>Lack of refrigerant, please call to service to add some T404A. and check for leaks</p> <p>Room temperature is too high. keep more ventilation, or install air conditioning to regulate chamber temperature at 25°C.</p> <p>Bad ventilation on both sides of the instrument. Move things away from both sides, keep good ventilation.</p>
Compressor is working without freezing	<p>Capillary tube (pressure relief valve) blocked, call service.</p> <p>Leakage of the refrigerant call service.</p>
Fan working but compressor not	Defective Compressor.
Freezing shelf temperature is unstable	Water in the cooling distribution system, ice blockage. If the outlet of the compressor is sometimes hot and some times cold, indicates water causing ice blockage. Call to manufacturer to vacuume and recharge system

9. Safety, Maintenance and importance

9.1 Make sure that hand wheel is in the lock position when the instrument is not in use.

9.2 Change the broken lamp accordingly (110V 12W), pay attention to safety when doing the exchange.






9.3 Change the broken lamp accordingly (110V 8W), pay attention to safety when doing the exchange.

9.4 The elements of the electric system, such as control panel, electric control box, require expert maintenance.

9.5 The cooling system mainly consists of compressor, condenser, filter, and evaporator. The condenser (with fan) may accumulate dust after long usage; this can affect the cooling function. It is necessary to open the side ventilation board and use brush to clean the dust on the projecting vanes of condenser, or use high pressure air to blow it clean.

9.6 Frequently clean the constant temperature chamber, to keep it clean. Remove blade before cleaning

9.7 There are 5 fuse seats at the rear of the instrument, from left to right models and functions are as following:

F1	F2	F3	F4	F5
				
5A	1A	1A	3A	3A
power	UV lamp	lamp	extractor	peltier

If broken, please change them according to above model.

10. Clean and maintenance



These machines are precision instruments. Pay attention to cleanliness, maintenance and repair work. Regularly clean and disinfect.

10.1 Daily cleaning work for freezing chamber

Conduct the following steps before cleaning each time:

- Take out all specimens from freezing chamber;
- Turn up the specimen grip to the top and lock the hand wheel.
- Remove the blade from the blade holder and put it back to the blade box;
- Activate extractor to clean away scraps;
- Empty waste tray;
- Use alcohol or common disinfectants based on alcohol to clean baffle and brush shelf;
- Collect waste cleaning liquid into wasted liquid barrel, then dispose it.



Do not use organic solvents or any other aggressive substances to clean and sterilize the chamber! Only use the cleaning agents specified in this manual such as alcohol or common disinfectants based on alcohol.



Please do manual defrost before cleaning if chamber full of frost.



If necessary, the external painted surface can be cleaned with light-duty commercial housework cleaner. Use wet cloth rub it until dry. Do not apply xylene or alcoholic liquid (eg: glass cleaner)

10.3 Spray disinfection

The cryostat has to be disinfected after each daily use and cleaning. Steps are as following:

- Cool down freezing chamber temperature to -20°C ;
- Spray disinfectant after temperature decrease enough;
- Wait for 15 minutes;
- Wipe chamber off with tissue paper or dry cloth .
- Reprogram temperature according to slicing requirement after disinfection.

10.4 Blade holder cleaning

If blade holder sticks contains impurities and is sticking, dismount and clean it as below:

- Turn the adjustable handles in the lateral of rotator and draw it out from sideward.
- Push the knife clamp back which have knife clip and shift it out from the rotator.
- Turn the eccentric rod handle in the lateral of the knife clap and draw it out from sideward.
- Clean all the assemblies of blade holder;
- Dry off all the assemblies and apply a drop of low temperature lubricating oil on them;
- Reassembly blade holder and put into freezing chamber.



Before putting the blade holder back make sure all the assemblies are dry.

10.5 Lubrication



Do low temperature oil lubrication for the following parts weekly. (1~2 drop is well enough).

10.6 Removal of the microtome

After long periods of not being used, inside of microtome needs lubrication, please remove microtome out of freezing chamber. Steps are as following:

- Remove glass window;
- Take all the accessories out chamber such as brush shelf, left storage shelf, right storage shelf, freezing disc, waste tray, blade holder.
- Use Allen key to loosen screw;
- Take off plug of step motor;
- Remove temperature sensor;
- Slightly lift the microtome and pull it to the left to disengage the plastic coupling (5) connecting the axes.
- Take the microtome out of the freezing chamber;
- Loosen the screws on both sides of cover
- Take away cover and dry off microtome;
- You may place microtome into oven at 40°C-50°C temperature about 4 hour to dry it. But you must apply low temperature lubrication oil to cross roller bearings.

Please pay attention to below notice when user put microtome back to freezing chamber carrying on opposite procedure:

- Ensure that the clamp is in the lowest position, and then place the microtome slightly into the left side of freezing chamber;
- Apply a drop of low temperature lubrication oil on the plastic coupling so that it may connect well with shaft;
- Use right hand to put handle of handwheel at the lowest position, and universal clamp at the lowest position too;
- Use left hand push microtome to right side, turn the handwheel up and down to ensure proper alignment of plastic coupling and shaft.



Do not directly use heater to dry off freezing chamber. This can cause damage to the cooling system!