

# Pierce<sup>®</sup> NHS-Activated Magnetic Beads

**88826 88827**

2320.0

Number	Description
88826	<b>Pierce NHS-Activated Magnetic Beads</b> , 1mL, supplied at 10mg/mL in <i>N,N</i> -dimethylacetamide (DMAC)
88827	<b>Pierce NHS-Activated Magnetic Beads</b> , 5mL, supplied at 10mg/mL in <i>N,N</i> -dimethylacetamide (DMAC)

**Note:** Before using, refer to the product label for the expiration date.

**Storage:** Upon receipt store at 4°C. Product shipped with an ice pack.

## Table of Contents

Introduction .....	1
Important Product Information .....	1
Procedure for Protein Immobilization .....	2
Procedure for Immunoprecipitation.....	5
A. Manual Antigen Immunoprecipitation.....	5
B. Automated Antigen Immunoprecipitation .....	5
Troubleshooting.....	6
Additional Information .....	7
Related Thermo Scientific Products .....	7

## Introduction

The Thermo Scientific Pierce NHS-Activated Magnetic Beads covalently immobilize proteins for the affinity purification of antibodies, antigens and other biomolecules. The activated magnetic beads contain *N*-hydroxysuccinimide (NHS) functional groups, which react with primary amines on proteins or other molecules to form stable amide linkages. The coupling reaction is performed in an amine-free buffer at pH 7-9. The beads are manually removed from the solution using a magnetic stand or by automation using an instrument such as the Thermo Scientific KingFisher Flex System. Automated instruments are especially useful for large screening of multiple samples.

**Table 1. Characteristics of Thermo Scientific Pierce NHS-Activated Magnetic Beads.**

<b>Composition:</b> <i>N</i> -hydroxysuccinimide (NHS) functional groups on a blocked magnetic bead surface
<b>Magnetization:</b> Superparamagnetic (no magnetic memory)
<b>Mean Diameter:</b> 1µm (nominal)
<b>Density:</b> 2.0g/cm <sup>3</sup>
<b>Bead Concentration:</b> 10mg/mL in DMAC
<b>Binding Capacity:</b> ≥ 26µg of rabbit IgG/mg of beads

## Important Product Information

- Magnetic beads are moisture-sensitive. To protect the beads, cap the bottle immediately after removing the slurry and wrap lab film around the cap before storing at 4°C.
- Do not centrifuge, dry or freeze the magnetic beads. Bead aggregation and loss of binding activity can result from using these methods.

- Estimate the amount of protein coupled to the magnetic beads with a protein assay (e.g., Thermo Scientific Pierce 660nm Protein Assay, Product No. 22660 and 22662) and subtract the amount of flow-through protein from the loaded protein. To measure the amount of protein on the bead directly, use the Thermo Scientific Pierce Micro BCA Protein Assay (Product No. 23235, see Tech Tip #75 from our website).
- For coupling antibodies to magnetic beads, ensure the antibody storage solution does not contain a protein stabilizer (e.g., BSA, gelatin), which inhibits coupling of the antibody to the beads. Protein stabilizers can be removed using the Thermo Scientific Pierce Antibody Clean-up Kit (Product No. 44600, see Tech Tip #55 from our website). For best results, buffer exchange the antibody into 50mM borate, pH 8.5 (e.g., Thermo Scientific BupH Borate Buffer Packs, Product No. 28384).
- Primary amine-containing buffers (e.g., Tris and glycine) inhibit coupling of protein to the magnetic beads. Remove primary amine-containing buffer using dialysis (e.g., Thermo Scientific Slide-A-Lyzer G2 Dialysis Cassettes, 10K MWCO, 3mL, Product No. 87730) or desalting (e.g., Thermo Scientific Zeba Spin Desalting Columns, 7K MWCO, 0.5mL, Product No. 89882).
- The coupling efficiency of protein to the magnetic beads varies depending on the specific protein. Typically, 0.1-2.0mg/mL of protein produces optimal protein coupling; however, optimize the concentration for each specific protein. As a reference, binding capacities of different proteins are listed in Table 2.

**Table 2. Binding capacity of proteins with different molecular weights on the Thermo Scientific Pierce NHS-Activated Magnetic Beads.**

Protein	Molecular weight (kDa)	Binding capacity of NHS beads (µg/mg of bead)
IgG	150	50
Streptavidin	53	24
Protein A/G	50	21
Protein L	36	39

**Note:** Results will vary depending on the number of accessible primary and secondary amines.

## Procedure for Protein Immobilization

**Note:** The following coupling procedure is for 300µL of magnetic beads in a 1.5mL microcentrifuge tube; scale the procedure as needed.

### A. Additional Materials Required

- Wash Buffer A: ice-cold 1mM hydrochloric acid
- Coupling Buffer: 50mM borate, pH 8.5 (Product No. 28384) or other amine-free buffer, pH 7-9
- Protein Solution: 0.1-2.0mg/mL in Coupling Buffer. For proteins already in solution, completely remove primary amine-containing buffer (e.g., Tris or glycine) using desalting or dialysis.
- Quenching Buffer: 3M ethanolamine, pH 9.0
- Storage Buffer: Coupling Buffer with 0.05% sodium azide
- 1.5mL microcentrifuge tubes
- Ultrapure water
- Wash Buffer B: 0.1M glycine, pH 2.0

### For Manual Coupling:

- Magnetic stand (e.g., Thermo Scientific Magnabind Magnet for 6 × 1.5mL microcentrifuge tubes, Product No. 21359)

---

**For Automated Coupling:**

- KingFisher® Flex System with 96 deep well head (Product No. 5400630) or KingFisher 96 System (Product No. 5400500)
- Thermo Scientific Microtiter Deep Well 96 Plate, V-bottom, polypropylene (100-1000µL; Product No. 95040450)
- KingFisher Flex 96 Tip Comb for Deep Well Magnets (Product No. 97002534)

**B. Manual Procedure for Coupling and Blocking**

1. Equilibrate protein solution and magnetic beads to room temperature.

**Note:** To ensure homogeneity, mix the magnetic beads thoroughly before use by repeated inversion, gentle vortexing or using a rotating platform.

2. Place 300µL of magnetic beads into a 1.5mL microcentrifuge tube.
3. Place the tube into a magnetic stand, collect the beads and discard the supernatant.
4. Add 1mL of ice-cold Wash Buffer A into the tube and gently vortex for 15 seconds to mix.
5. Place the tube into a magnetic stand, collect the beads and discard the supernatant.

**Note:** Immediately proceed with adding the protein solution.

6. Add 300µL of protein solution into the tube and vortex for 30 seconds.
7. Incubate the tube for 1-2 hours at room temperature on a rotator. During the first 30 minutes of the incubation, vortex the tube for 15 seconds every 5 minutes. For the remaining time, vortex the tube for 15 seconds every 15 minutes until incubation is complete.

**Note:** If required, incubate overnight at 4°C.

8. Collect the beads with a magnetic stand and save the flow-through.
9. Add 1mL of Wash Buffer B to the beads and vortex the tube for 15 seconds.
10. Place the tube into a magnetic stand, collect the beads and discard the supernatant.
11. Repeat Steps 9 and 10 one time.
12. Add 1mL of ultrapure water to the beads and vortex for 15 seconds.
13. Place the tube into a magnetic stand, collect the beads and discard the supernatant.
14. Add 1mL of Quenching Buffer to the beads and vortex the tube for 30 seconds.
15. Incubate the tube for 2 hours at room temperature on a rotator.
16. Place the tube into a magnetic stand, collect the beads and discard the supernatant.
17. Add 1mL of purified water to the tube, mix well, collect the beads with a magnetic stand and discard the supernatant.
18. Add 1mL of Storage Buffer to the tube, mix well, collect the beads with a magnetic stand and discard the supernatant. Repeat this wash two additional times.
19. Add 300µL of Storage Buffer to the beads, mix well and store at 4°C until ready for use.

**Note:** The final concentration of the protein-coupled magnetic beads is 10mg/mL.

**C. Automated Procedure for Coupling and Blocking**

**Note:** The following protocol is designed for use with the KingFisher Flex or KingFisher 96 Instrument. The protocol can be modified according to your needs using the Thermo Scientific BindIt Software provided with the instrument.

1. Enter the “Protein Coupling” protocol from Table 3 into the BindIt® Software on an external computer.
2. Transfer the protocol to the KingFisher Flex or KingFisher 96 Instrument from an external computer. See the BindIt Software User Manual for detailed instructions on importing protocols.

3. Set up plates according to Table 3.

<b>Table 3. Pipetting instructions for the Coupling and Blocking protocol using the Microtiter Deep Well 96 Plates.</b>				
<b>Plate #</b>	<b>Plate Name</b>	<b>Content</b>	<b>Volume</b>	<b>Time/Speed</b>
1	Beads	NHS-Activated Beads	300µL	Collect Beads
2	Wash 1	Wash Buffer A	1mL	10 seconds/Slow
3	Coupling	Protein Sample in Coupling Buffer	300µL	1-2 hours/Slow
4	Wash 2	Wash Buffer B	1mL	10 seconds/Slow
5	Wash 3	Wash Buffer B	1mL	10 seconds/Slow
6	Wash 4	Purified Water	1mL	10 seconds/Slow
7	Quench	Quenching Buffer	1mL	2 hours/Slow
8	Wash 5	Purified Water	1mL	30 seconds/Slow
9	Wash 6	Storage Buffer	1mL	30 seconds/Slow
10	Storage	Storage Buffer	300µL	Release Beads
11	Tip Plate	KingFisher Flex 96 Tip Comb for Deep Well Magnets	-	-

4. Select the protocol using the arrow keys on the instrument keypad and press Start. See the KingFisher Flex or KingFisher 96 Instrument User Manual for detailed information.
5. Slide open the door of the instrument's protective cover.
6. Load the Tip Plate (plate #11) and press Start. The instrument places the Tip Comb onto the magnet head.
7. Remove the Tip Plate.
8. Load plates #1-8 into the instrument according to the protocol requests; place each plate in the same orientation. Confirm each action by pressing Start.
9. After sample processing through plate #8, the instrument will pause and instruct to remove each individual processed plate while simultaneously loading the remaining three plates, plates #9-11. For example, remove plate #1 and load plate #9 into that position. Confirm each action by pressing Start.
10. After sample processing, remove the plates as instructed by the instrument display. Press Start after each plate. Press Stop after removing all of the plates.

**Notes:**

- Load ice-cold 1mM Wash Buffer A in plate #2 immediately before instrument loading to ensure the buffer remains cold.
- If using fewer than 96 wells, fill the same wells in each plate. For example, if using wells A1 through A12, use these same wells in all plates.
- Combine the Tip Comb with a Deep Well 96 Plate. See the instrument user manual for detailed instructions.

---

## Procedure for Immunoprecipitation

### A. Additional Materials Required

- 1.5mL microcentrifuge tubes
- Binding/Wash Buffer: Tris-buffered saline (TBS, Product No. 28360) containing 0.05% Tween-20<sup>®</sup> Detergent and 0.5M NaCl
- Low-pH Elution Buffer: 0.1M glycine, pH 2-2.5
- Antibody for immunoprecipitation
- Antigen sample for immunoprecipitation (e.g., cell lysate)
- Neutralization Buffer: High-ionic strength alkaline buffer such as a 1M phosphate or 1M Tris; pH 7.5-9
- Optional: Protease inhibitor cocktail (e.g., Thermo Scientific Halt Protease Inhibitor Single-Use Cocktail EDTA-free, Product No. 78425 )

### For Manual IP:

- Magnetic stand (e.g., MagnaBind Magnet for 6 × 1.5mL Microcentrifuge Tubes, Product No. 21359)

### For Automated IP:

- KingFisher<sup>®</sup> Flex System with 96 deep well head (Product No. 5400630) or KingFisher 96 System (Product No. 5400500)
- Thermo Scientific Microtiter Deep Well 96 Plate, V-bottom, polypropylene (100-1000µL, Product No. 95040450)
- KingFisher Flex 96 Tip Comb for Deep Well Magnets (Product No. 97002534)

### B. Manual Antigen Immunoprecipitation

1. Add 25µL (0.25mg) of antibody-coupled magnetic beads to a 1.5mL centrifuge tube.
2. Place tube in a magnetic stand, collect the beads and discard the supernatant.
3. Dilute antigen sample for immunoprecipitation to 1-2mg/mL using Binding/Wash Buffer.
4. Add 500µL of diluted antigen sample to the tube containing antibody-coupled magnetic beads and incubate for 1-2 hours at room temperature on a rotator or mixer. Gently vortex the beads every 10-15 minutes during incubation to ensure the beads remain in suspension.
5. Collect the beads with a magnetic stand, remove the unbound sample and save for analysis.
6. Add 500µL of Binding/Wash Buffer to the tube and gently mix. Collect the beads on a magnetic stand and discard the supernatant. Repeat this step one time.
7. Add 500µL of ultrapure water to the tube and gently mix. Collect the beads on a magnetic stand and discard the supernatant.
8. Add 100µL of Elution Buffer to the tube. Incubate for 5 minutes at room temperature on a rotator or mixer. Magnetically separate the beads and save the supernatant containing the target antigen.  
**Note:** One elution may be sufficient; however, optimization is required for each system.
9. Repeat Step 8 and combine the two eluates. To neutralize the low pH of the solution, add 10µL of Neutralization Buffer for each 100µL of eluate.

### C. Automated Antigen Immunoprecipitation

**Note:** The following protocol is designed for use with the KingFisher Flex or KingFisher 96 Instruments. The protocol can be modified according to your needs using the BindIt Software provided with the instrument.

1. Enter the “Direct IP” protocol from Table 3 into the BindIt Software on an external computer.

- Transfer the protocol to the KingFisher Flex or KingFisher 96 Instrument from an external computer. See the BindIt Software User Manual for detailed instructions on importing protocols.
- Set up plates according to Table 4.

**Table 4. Pipetting instructions for the IP protocol using Microtiter Deep Well 96 Plates.**

Plate #	Plate Name	Content	Volume	Time/Speed
1	Bind	Beads	25 $\mu$ L	2 hour/Slow
		Antigen Sample for IP	500 $\mu$ L	
2	Wash 1	Binding/Wash Buffer	500 $\mu$ L	30 seconds/Slow
3	Wash 2	Binding/Wash Buffer	500 $\mu$ L	30 seconds/Slow
4	Wash 3	Ultrapure Water	500 $\mu$ L	30 seconds/Slow
5	Elution 1	Elution Buffer	100 $\mu$ L	5 minutes/Slow
6	Elution 2	Elution Buffer	100 $\mu$ L	5 minutes/Slow
7	Tip Plate	KingFisher Flex 96 Tip Comb for Deep Well Magnets	-	10 seconds/Fast

- Select the protocol using the arrow keys on the instrument keypad and press Start. See the KingFisher Flex or KingFisher 96 Instrument User Manual for detailed information.
- Slide open the door of the instrument's protective cover.
- Load plates into the instrument according to the protocol requests, placing each plate in the same orientation. Confirm each action by pressing Start.
- After sample processing, remove the plates as instructed by the instrument's display. Press Start after each plate. Press Stop after removing all of the plates.

**Notes:**

- The low-pH elutions must be neutralized by adding 10 $\mu$ L of Neutralization Buffer for each 100 $\mu$ L of eluate directly to each well immediately after incubation.
- If fewer than 96 wells are used, fill the same wells in each plate. For example, if using wells A1 through A12, use these same wells in all plates.
- Combine the Tip Comb with a Deep Well 96 Plate. See the instrument user manual for detailed instructions.

**Troubleshooting**

Problem	Possible Cause	Solution
Low coupling efficiency	Primary amine-containing buffer was not completely removed before coupling	Dialyze or desalt sample to completely remove Tris and glycine
	Protein addition was delayed	Immediately mix protein with beads after washing
Protein is not soluble in coupling buffer	Molecule was hydrophobic	Dissolve protein in coupling buffer containing up to 20% DMSO
Beads aggregate during the coupling process	Proteins on bead surface adhered to tube plastic	After blocking, add 0.05% detergent (e.g., Tween-20 Detergent) to the water wash and the Storage Buffer (Part C) <b>Note:</b> Do not use detergent in the coupling step

---

## Additional Information

Visit [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce) for additional information relating to this product, including the following:

- Frequently Asked Questions
- Tech Tip #43: Protein stability and storage
- Tech Tip #55: Remove BSA and gelatin from antibody solutions using Melon Gel
- Tech Tip #75: Measure protein bound to Pierce NHS-Activated Magnetic Beads
- Visit [www.thermoscientific.com/kingfisher](http://www.thermoscientific.com/kingfisher) for information on KingFisher Products.
- In the U.S.A., purchase KingFisher Supplies from VWR. Outside the U.S.A., contact your local Thermo Fisher Scientific office to purchase KingFisher Supplies.

## Related Thermo Scientific Products

88828	Pierce Direct Magnetic IP/Co-IP Kit
88802	Pierce Protein A/G Magnetic Beads, 1mL
88803	Pierce Protein A/G Magnetic Beads, 5mL
88804	Pierce Protein A/G Magnetic IP/Co-IP Kit
88805	Pierce Crosslink Magnetic IP/Co-IP Kit
88816	Pierce Streptavidin Magnetic Beads, 1mL
88817	Pierce Streptavidin Magnetic Beads, 5mL
88821	Pierce Glutathione Magnetic Beads, 4mL
88822	Pierce Glutathione Magnetic Beads, 20mL
23235	Micro BCA Protein Assay Kit
21030	Gentle Ag/Ab Binding and Elution Buffer Kit
44600	Pierce Antibody Clean-up Kit

Tween is a trademark of Croda International Plc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.