



# Agilent 1260 Infinity Manual Injector

User Manual



**Agilent Technologies**

# Notices

© Agilent Technologies, Inc. 2006, 2008, 2010

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

## Manual Part Number

G1328-90012

## Edition

06/2010

Printed in Germany

Agilent Technologies  
Hewlett-Packard-Strasse 8  
76337 Waldbronn

**This product may be used as a component of an in vitro diagnostic system if the system is registered with the appropriate authorities and complies with the relevant regulations. Otherwise, it is intended only for general laboratory use.**

## Warranty

**The material contained in this document is provided “as is,” and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.**

## Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

## Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as “Commercial computer software” as defined in DFAR 252.227-7014 (June 1995), or as a “commercial item” as defined in FAR 2.101(a) or as “Restricted computer software” as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies’ standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will

receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

## Safety Notices

### CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

### WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

# In This Guide

This manual covers the Agilent 1260 Infinity Manual Injector (G1328C)

## **1 Introduction**

Operation and mechanical hardware

## **2 Installing the Manual Injector**

Installation of the manual injector

## **3 Using the Manual Injector**

How to use the manual injector

## **4 Maintenance**

Instructions on simple, routine repair procedures

## **5 Parts and Materials for Maintenance**


Detailed illustrations and lists for identification of parts and materials

## **6 Appendix**

Additional information

# Contents

<b>1</b>	<b>Introduction</b>	<b>5</b>
	Introduction to the Manual Injector	6
<b>2</b>	<b>Installing the Manual Injector</b>	<b>7</b>
	Unpacking the Manual Injector	8
	Installing the Manual Injector	9
	Flow Connections	13
	Leak Drainage	15
<b>3</b>	<b>Using the Manual Injector</b>	<b>17</b>
	Solvent Information	18
	Choice of Injection Seal	19
	Needles	20
	Injecting Sample	21
<b>4</b>	<b>Maintenance</b>	<b>23</b>
	Overview of Maintenance	24
	Flushing the Manual Injector	25
	Cleaning the Manual Injector	26
	Stator Face	27
	Injection-Valve Seal	29
	Position-Sensing Switch	32
<b>5</b>	<b>Parts and Materials for Maintenance</b>	<b>35</b>
	Manual Injector	36
	Injection-Valve Assembly	38
<b>6</b>	<b>Appendix</b>	<b>41</b>
	Agilent Technologies on Internet	42



# 1 Introduction

Introduction to the Manual Injector 6

Operation and mechanical hardware




## Introduction to the Manual Injector

**NOTE**


See also the Rheodyne 7725i operating instructions supplied with the injection valve.

The Agilent 1260 Infinity Manual Injector uses a Rheodyne 7725i 7-port sample injection valve. Sample is loaded into the external 20- $\mu$ l sample loop through the injection port at the front of the valve. The valve has a ceramic stator and Vespel™ injection seal (for pH above 10, a Tefzel™ seal is available). A make-before-break passage in the stator ensures flow is not interrupted when the valve is switched between the INJECT and LOAD positions, and back again (see also “Needles” on page 20 and “Flow Connections” on page 13).

The valve is mounted on a steel mounting pole, and can be installed at the left- or right-hand side of the LC system.



**Figure 1** Rheodyne 7725i Injection Valve



## 2 Installing the Manual Injector

Unpacking the Manual Injector	8
Damaged Packaging	8
Delivery Checklist	8
Installing the Manual Injector	9
Flow Connections	13
Leak Drainage	15

Installation of the manual injector



## Unpacking the Manual Injector

### Damaged Packaging

Upon receipt of your manual injector, inspect the shipping containers for any signs of damage. If the containers or cushioning material are damaged, save them until the contents have been checked for completeness and the manual injector has been mechanically checked. If the shipping container or cushioning material is damaged, notify the carrier and save the shipping material for the carriers inspection.

### Delivery Checklist

Ensure all parts and materials have been delivered with the manual injector. The delivery checklist is shown in [Table 1](#) on page 8. To aid in parts identification, please see [“Parts and Materials for Maintenance”](#) on page 35. Please report missing or damaged parts to your local Agilent Technologies sales and service office.

**Table 1** Manual Injector Checklist

Description	Quantity
Manual injection valve (p/n 5063-6502) with start cable, <i>including</i> : operating instructions, needle port cleaner, vent tubes (×2) and fittings, 5/64 and 9/64-inch hex keys	1
Mounting pole (p/n 5001-3738)	1
Connection capillary, 0.17 mm id, 500 mm (p/n G1328-87600)	1
Base plate (p/n G1328-44111)	1
Organizer plate (p/n 5042-8553)	1
Catch tube cap (p/n 5042-8576)	1
Valve syringe, fixed needle, 50 µL (p/n 5182-9619)	1
User Manual (p/n G1328-90011)	1



## Installing the Manual Injector

### CAUTION

"Defective on arrival" problems


If there are signs of damage, please do not attempt to install the module. Inspection by Agilent is required to evaluate if the instrument is in good condition or damaged.

- Notify your Agilent sales and service office about the damage.
- An Agilent service representative will inspect the instrument at your site and initiate appropriate actions.

### NOTE

The manual injector can be installed at the left- or right-hand side of the instrument stack.

- 1 Place the baseplate on the bench.
- 2 Connect the two organizer plates to the base plate.




**Figure 2** Connecting the Organizer Plates

- 3 Screw the mounting pole into one of the two holes in the organizer plate.

## 2 Installing the Manual Injector


### Installing the Manual Injector

- Slide the manual injector onto the mounting pole (see [Figure 3](#) on page 10). Tighten the mounting screw.



**Figure 3** Installing the Mounting Pole and Manual Injector


- 5 Install other system modules on top of the manual injector baseplate (see [Figure 4](#) on page 11).



**Figure 4** Installing the System

## 2 Installing the Manual Injector

### Installing the Manual Injector



**Figure 5** Installing the Start Cable

---

1. See [Figure 6](#) on page 14

---

**6** Connect the capillaries to the manual injector (see “[Flow Connections](#)” on page 13).

## Flow Connections

### WARNING

**Toxic and hazardous solvents and flammable liquids**

**The handling of solvents and reagents can hold health risks.**

- When working with solvents observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents and flammable liquids are used.
- 

### CAUTION

*Prevent siphoning*


- The outlets of the two vent capillaries (ports 5 and 6) and the needle port must be at the same level to prevent siphoning (see [Figure 7](#) on page 14).
- 

- 1 Connect the pump outlet capillary to port 2.
- 2 Connect the column-compartment inlet capillary to port 3.
- 3 Connect the sample loop between ports 1 and 4.


## 2 Installing the Manual Injector

### Flow Connections

- 4 Connect one vent capillary (supplied with valve) to port 5 and one to port 6.



**Figure 6** Flow Connections



**Figure 7** Vent Capillaries


## Leak Drainage

**WARNING**

*Leaking injector fittings*

In the event of a leak, solvent will drop into the leak channel in the baseplate, from where it is channelled to the front and back of the baseplate.

→ Check the manual injector fittings periodically for signs of leakage.



**Figure 8** Leak Drainage

## **2** **Installing the Manual Injector** Leak Drainage





### 3 Using the Manual Injector

Solvent Information 18

Choice of Injection Seal 19

Needles 20

Injecting Sample 21

How to use the manual injector



## Solvent Information

Observe the following recommendations on the use of solvents.


### Flow Cell

Long term operation at pH > 11 should be avoided. Never leave strongly alkaline solutions in the flow cell without flow.

### Solvents

Always filter solvents through 0.4 µm filters, small particles can permanently block filters, frits and capillaries. Avoid the use of the following steel-corrosive solvents:

- Solutions of alkali halides and their respective acids (for example, lithium iodide, potassium chloride, and so on).
- High concentrations of inorganic acids like sulfuric acid, especially at higher temperatures (replace, if your chromatography method allows, by phosphoric acid or phosphate buffer which are less corrosive against stainless steel).
- Halogenated solvents or mixtures which form radicals and/or acids, for example:



This reaction, in which stainless steel probably acts as a catalyst, occurs quickly with dried chloroform if the drying process removes the stabilizing alcohol.

- Chromatographic grade ethers, which can contain peroxides (for example, THF, dioxane, di-isopropylether). Such ethers should be filtered through dry aluminium oxide which adsorbs the peroxides.
- Solutions of organic acids (acetic acid, formic acid, and so on) in organic solvents. For example, a 1-% solution of acetic acid in methanol may attack steel.
- Mixtures of carbon tetrachloride with 2-propanol or THF. dissolve stainless steel.

## Choice of Injection Seal

The manual injector is supplied with a Vespel™ injection seal as standard. Vespel™ is sensitive to alkaline attack, so avoid mobile phases with a pH of 10 or more. Use the Tefzel™ injection seal (see [“Injection-Valve Assembly”](#) on page 38) for mobile phases with a pH above 10.

## Needles

**CAUTION**

Needle can damage valve

→ Always use the correct size needle.

---

The manual injector is not supplied with syringes or needles.

Use needles with 0.028-inch outer diameter (22 gauge) × 2-inch long needle, without electro-taper, and with 90° point style (square tip).

## Injecting Sample

### WARNING


#### Ejection of mobile phase

When using sample loops larger than 100  $\mu\text{l}$ , mobile phase may be ejected from the needle port as the mobile phase in the sample loop decompresses.

- Please observe appropriate safety procedures (for example, goggles, safety gloves and protective clothing) as described in the material handling and safety data sheet supplied by the solvent vendor, especially when toxic or hazardous solvents are used.

### LOAD Position

In the LOAD position (see [Figure 9](#) on page 21), the pump is connected directly to the column (ports 2 and 3 connected), and the needle port is connected to the sample loop. At least 2 to 3 sample-loop volumes (more if better precision is required) of sample should be injected through the needle port to provide good precision. The sample fills the loop, and excess sample is expelled through the vent tube connected to port 6.




**Figure 9** LOAD Position


### 3 Using the Manual Injector Injecting Sample

#### INJECT Position

In the INJECT position (see [Figure 10](#) on page 22), the pump is connected to the sample loop (ports 1 and 2 connected). All of the sample is washed out of the loop onto the column. The needle port is connected to the vent tube (port 5).



**Figure 10** INJECT Position



## 4 Maintenance

Overview of Maintenance	24
Flushing the Manual Injector	25
Cleaning the Manual Injector	26
Stator Face	27
Injection-Valve Seal	29
Position-Sensing Switch	32

Instructions on simple, routine repair procedures



## Overview of Maintenance

**Table 2** Overview of Repair Procedures

<b>Procedure</b>	<b>Typical Frequency</b>	<b>Time Required</b>	<b>Notes</b>
Flushing the injector	After using aqueous buffers or salt solutions	5 minutes	See <a href="#">“Flushing the Manual Injector”</a> on page 25
Exchanging the stator face	When visibly scratched, or when the valve performance shows indication of leakage or wear	10 minutes	See <a href="#">“Stator Face”</a> on page 27
Exchanging the injection-valve seal	After approximately 10000 to 20000 injections, or when the valve performance shows indication of leakage or wear	10 minutes	See <a href="#">“Injection-Valve Seal”</a> on page 29
Exchanging the position-sensing switch	When cable damaged or when no start signal is sent when switching to the inject position	10 minutes	See <a href="#">“Position-Sensing Switch”</a> on page 32



## Flushing the Manual Injector


### CAUTION

#### *Damage through crystal formation*

The use of aqueous buffers or salt solutions can lead to crystal formation which may cause scratches on the injection seal.

→ Always rinse the valve with water after using aqueous buffers or salt solutions.

- 1 Switch the valve to the INJECT position.
- 2 Use the pump to flush the sample loop and seal grooves.
- 3 Use the needle-port cleaner (supplied with the valve) and syringe to flush the needle port and vent capillary.



**Figure 11** Needle-port Cleaner

## **Cleaning the Manual Injector**

The manual injector base should be kept clean. Cleaning should be done with a soft cloth slightly dampened with water or a solution of water and a mild detergent.


## Stator Face

- When**
- Poor injection-volume reproducibility
  - Leaking injection valve


**Tools required** Hex key, 9/64 inch (supplied with valve)

<b>Parts required</b>	<b>#</b>	<b>p/n</b>	<b>Description</b>
	1	0100-1859	Stator face

**1** Loosen the three stator screws.




**2** Remove the stator head and stator face.




## 4 Maintenance

### Stator Face


**3** Insert the new stator face onto the stator head.



**4** Install the stator head and stator face onto the valve. Ensure the pin in the stator ring is aligned with the hole in the stator head.



**5** Secure the stator face and stator head in place with the stator screws. Tighten each screw alternately  $\frac{1}{4}$ -turn until the stator head is secure.




# Injection-Valve Seal

- When**
- Poor injection-volume reproducibility
  - Leaking injection valve


**Tools required** Hex key, 9/64 inch (supplied with valve)

Parts required	#	p/n	Description
	1	0100-0623 (VespeI™), 0100-0620(Tefz eI™)	Rotor seal

**1** Loosen the three stator screws.




**2** Remove the stator head and stator face (see “Stator Face” on page 27).




## 4 Maintenance

### Injection-Valve Seal


**3** Remove the stator ring.




**4** Remove the seal.




**5** Install the new seal. Ensure the seal is positioned as shown.




**6** Install the stator ring. Ensure the pin in the stator ring is aligned with the hole in the valve body.




7 Insert the stator face onto the stator head.



8 Install the stator head and stator face onto the valve. Ensure the pin in the stator ring is aligned with the hole in the stator head.



9 Secure the stator face and stator head in place with the stator screws. Tighten each screw alternately ¼-turn until the stator head is secure.



## 4 Maintenance

### Position-Sensing Switch


# Position-Sensing Switch

**When** • No start signal when switching to the inject position


**Tools required** Hex key, 9/64 inch (supplied with valve)

Parts required	#	p/n	Description
	1	0490-1849	Position-sensing switch

**1** Loosen the three stator screws.




**2** Remove the stator head and stator face (see "Stator Face" on page 27).







**3** Remove the stator ring.




**4** Disconnect the sensor cable from the start cable. Pull the sensing switch out of the stator ring.



**5** Insert the new sensing switch into the stator ring.




**6** Install the stator ring. Ensure the pin in the stator ring is aligned with the hole in the valve body.




## 4 Maintenance

### Position-Sensing Switch


**7** Insert the stator face onto the stator head.




**8** Install the stator head and stator face onto the valve. Ensure the pin in the stator ring is aligned with the hole in the stator head.



**9** Secure the stator face and stator head in place with the stator screws. Tighten each screw alternately  $\frac{1}{4}$ -turn until the stator head is secure.



**10** Reconnect the sensor cable to the start cable.





## 5 Parts and Materials for Maintenance

Manual Injector 36

Injection-Valve Assembly 38


Detailed illustrations and lists for identification of parts and materials



## Manual Injector

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	5063-6502	Manual injection valve
2	5001-3738	Mounting pole
3	G1328-44111	Base plate
4	5042-8553	Organizer plate
	5042-8576	Catch tube cap
5	5042-8901	Name plate
	5182-9619	Valve syringe, fixed needle, 50 µL
	G1328-87600	Connection capillary, 0.17 mm id, 500 mm
	0100-1677	Start cable

Manual injector valve see [“Injection-Valve Assembly”](#) on page 38.




**Figure 12** Manual Injector

## Injection-Valve Assembly

Manual injection valve (p/n 5063-6502) with starts cable (complete assembly), including operating instructions, needle port cleaner, vent tubes (×2) and fittings, 5/64 and 9/64-inch hex keys. Includes items 1 – 8 .

<b>Item</b>	<b>p/n</b>	<b>Description</b>
1	0101-1254	Bearing ring — order rebuild kit instead
2	1535-4046	Isolation seal
3	0101-0623	Rotor seal (Vespe <sup>TM</sup> )
	0101-0620	Rotor seal (Tefzel <sup>TM</sup> )
	0101-1255	Rotor seal (PEEK <sup>TM</sup> )
4	0100-1859	Stator face
5	0100-1860	Stator head
6	1535-4857	Stator screws
7	8710-0060	Hex-key wrench, 9/64 inch
8	0490-1849	Position-sensing switch



**Figure 13** Injection-Valve Assembly

## 5 Parts and Materials for Maintenance

### Injection-Valve Assembly

#### Sample loops stainless steel

<b>p/n</b>	<b>Description</b>
0101-1248	Sample loop 5 µL
0100-1923	Sample loop 10 µL
0100-1922	Sample loop 20 µL
0100-1924	Sample loop 50 µL
0100-1921	Sample loop 100 µL
0101-1247	Sample loop 200 µL
0101-1246	Sample loop 500 µL
0101-1245	Sample loop 1 mL
0101-1244	Sample loop 2 mL
0101-1243	Sample loop 5 mL

#### Sample loops PEEK

<b>p/n</b>	<b>Description</b>
0101-1241	Sample loop 5 µL
0101-1240	Sample loop 10 µL
0101-1239	Sample loop 20 µL
0101-1238	Sample loop 50 µL
0101-1242	Sample loop 100 µL
0101-1227	Sample loop 200 µL
0101-1236	Sample loop 500 µL
0101-1235	Sample loop 1 mL
0101-1234	Sample loop 2 mL
0101-1230	Sample loop 5 mL





## 6 Appendix

Agilent Technologies on Internet 42

Additional information



## **Agilent Technologies on Internet**

For the latest information on products and services visit our worldwide web site on the internet at:

<http://www.chem.agilent.com>

# Index

## A

Agilent on internet 42  
APG-remote 12

## B

buffers 25

## C

capillary connections 13  
cleaning the manual injector 26

## D

damaged packaging 8  
delivery checklist 8  
description 6

## E

exchanging  
    injection seal 29  
    position-sensing switch 32  
    stator face 27

## F

flow cell  
    solvent information 18  
flow connections 13  
flushing the manual injector 25

## I

information  
    on internet 42  
INJECT 6, 22

injecting sample 6  
injection seal 29  
    tefel 19  
    vespel 19  
installation 9  
installing the manual injector 9  
internet 42, 42

## L

leak channel 15  
leak drainage 15  
leaks 13  
LOAD 6, 21

## M

make-before-break 6

## N

needle type 20  
needle-port cleaner 25  
needles 20

## P

position-sensing switch 32  
position-sensor cable 12  
precision 21

## R

repair procedures  
    injection seal 29  
    position-sensing switch 32  
    stator face 27  
Rheodyne 7725i 6

## S

salt solutions 25  
sample loops 6  
sample volume 21  
solvents  
    acids 18  
    buffers 25  
    ethers 18  
    rinsing 25  
    salt solutions 25  
    steel-corrosive 18  
start cable 12  
stator face 27

## U

unpacking  
    damaged packaging 8  
    delivery checklist 8

## W

worldwide web 42

[www.agilent.com](http://www.agilent.com)

## In This Book

This manual contains user information about the Agilent 1260 Infinity Manual Injector. The manual describes the following:

- introduction to the manual injector,
- installing the manual injector,
- using the manual injector
- maintenance of the manual injector,
- parts and materials, and
- additional information.

© Agilent Technologies 2006, 2008, 2010

Printed in Germany  
06/2010



G1328-90012



**Agilent Technologies**