Your **Power** for Health







CheckExtractor™

www.gbo.com/diagnostics



CheckExtractor™

Instructions for Use

For in-vitro diagnostic use by professional laboratory personnel only

Revision: BQ-395-01 / June 2015



Greiner Bio-One GmbH Maybachstr. 2 72636 Frickenhausen Germany ermany REF 865 070



1. GENERAL

These Instructions for Use are written for operators of the CheckExtractor[™] using software version 1.1 or higher. Due to software evolution slight changes in the dialog boxes displayed in this manual may occur. When a modification is made that may have adverse effects on the system using this manual, an update or amendment will become available.

Every effort has been made to ensure that the information contained in this manual is updated and accurate as of the date of issue or revision and it is consistent with the product described. However, no guarantee is given or implied that the document is error-free or accurate with regard to any specifications. The instrument manufacturer reserves the right to make changes to this manual without notification.

1.1 Intended Use

Automated System for the extraction of nucleic acids.

The system has been validated for the following applications:

- oCheck[®] DNA Extraction Kit CheckExtractor™ (REF 517070)
- PapilloCheck® for 288 analyses (REF 465088)
- PelvoCheck[®] for 288 analyses (REF 504288)

1.2 Copyright

© 2015, Greiner Bio-One. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from Greiner Bio-One.

The information in this guide is subject to change without notice. Greiner Bio-One reserves the right to change its products and services at any time to incorporate the latest technological developments. Although this guide has been prepared with every precaution to ensure accuracy, Greiner Bio-One assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information. Greiner Bio-One welcomes customer input on corrections and suggestions for improvement.

Microsoft[®], Windows[®] and Excel[®] are registered trademarks of the Microsoft[®] Corporation. Microlab[™] is a registered trademark of Hamilton Bonaduz AG. The Greiner Bio-One Logo, the CheckExtractor[™] Logo and oCheck® are property of Greiner Bio-One. All other company and product names might be trademarks of the respective companies with which they are associated.

Patent no: EP 1171240 / US 7033543 (Hamilton Bonaduz AG)

Approvals

This instrument is in compliance with the following directives of the Council of the European Community and the according standards: 98/97/EC, 2004/108/EC, ISO 13485:2007, DIN EN ISO 9001:2000, DIN EN 18113-3, ISO 14971.

1.3 Trademark

The used brand names and labels are the property of the respective owner.

1.4 Terms used in this Manual

GLOSSARY OF SYMBOLS

	Caution
4	Electronic equipment can be the source of electrical shocks. DO NOT attempt to gain access to parts of the instrument marked with this label. DO NOT attempt to repair any electrical parts. DO NOT open the instrument housing. Service and repair should only be carried out by authorized and qualified personnel of the manufacturer.
	Hot Surface DO NOT touch.
	Laser beam DO NOT stare into the beam.
	Biohazard
	Harmful
	UV-light
	Waste Electrical and Electronic Equipment: Collect separately ac- cording to national laws.
!	Note
Â	Warning



1.4.1 Abbreviations

Name	Meaning	
DNA	Deoxyribonucleic Acid	
GLP	Good Laboratory Practices	
GBO	Greiner Bio-One	
GUI	Graphical User Interface	
HHS	Heater Shaker	
ID	Identifier	
IFU	Instructions For Use	
IVD	In-Vitro-Diagnostic	
LED	Light Emitting Diode	
LIMS	Laboratory Information Management System	
.NET	Software Framework for running Windows-based applications	
PC	Personal Computer	
PCR	Polymerase Chain Reaction	
PDF	Portable Document Format	
RF	Radio frequency	
PK	Proteinase K	
RNA	Ribonucleic Acid	
rpm	Rotations per Minute	
SW	Software	
UPS	Uninterruptible Power Supply	
USB	Universal serial bus	
WEEE	Waste electrical and electronic equipment	

1.4.2 Name and Description of Components (Overview)

Name	Meaning
Controlling PC	The personal computer (PC) which runs the software that controls the CheckExtractor ^{TM} .
Eluate	Extracted DNA in elution buffer
Elution plate	Barcoded PCR plate into which the extracted DNA is eluted.
MASTERBLOCK®	Deep well plate in which the extraction steps are carried out.
Final MasterMix	The solution containing the oCheck [®] PCR MasterMix mixed with the enzymes DNA polymerase and UNG.
Final PCR reaction / PCR reaction mix	The mix of the final MasterMix and the eluate, ready for the PCR.
Heater shaker, HHS	Device on the deck that heats up the MASTERBLOCK [®] and also shakes it for mixing purposes.
Input file	Excel [®] file with samples to be processed.
Liquid waste	All liquid material that requires disposal.
Liquid waste bottle	Container which collects the liquid waste.
Loading deck	Platform at the front of the CheckExtractor [™] on which the carriers are placed to be loaded.
Magnetic separator	Plate with 24 magnets required for the separation steps. The magnetic separator must be placed in the CheckExtractor [™] on the position shown in Figure 5.
Method	Script or program which defines the steps the CheckExtractor [™] has to carry out and the amounts of liquid to pipet.
(oCheck [®]) PCR MasterMix and final PCR MasterMix	oCheck [®] PCR MasterMix is delivered in the oCheck [®] assay without the enzymes DNA polymerase and UNG. Therefore, the final PCR MasterMix contains oCheck [®] PCR MasterMix, HotStarTaq [®] Polymerase and Uracil-N-Glycosylase
PCR plate	Barcoded PCR plate which will contain the final PCR reaction.
PCR setup	Preparation of the PCR, including preparation of the final MasterMix and the PCR setup method.
PCR setup method	Method in which the final MasterMix is added to the eluate.
Pipetting tracking file	File in which all pipetting steps that were carried out in an extraction run are listed. Volumes and possible errors are stored here. See also Report file.
Plate carrier	Carrier for the MASTERBLOCK [®] , the bottle adapter, the elution plate and the PCR plate.
Report file	File in which all pipetting steps that were carried out during the run (extraction run or PCR setup run) are listed. Volumes and possible errors are stored here. The scheme of the plate is indicated.
Run Control	Part of the CheckExtractor [™] 's software that does the communication between the controlling PC and the CheckExtractor [™] itself. It therefore controls the sequence of the runs.
Sample carrier	Carrier for sample tubes.
Sample ID	Unique identifier for each sample.
Tip carrier	Carrier for the 300 μL tips and the 1000 μL tips.
Track	Position or rail on the loading desk. All carriers are loaded in at least one of these tracks. The tracks are numbered.
Trough carrier	Carrier for five reagent troughs.

Tube carrier	Carrier for the tubes for the final PCR MasterMix, the proteinase K and the magnetic beads (M0).
User ID	Unique identifier for each user. Following IVD and FDA formalities it is mandatory that every user, who works with the CheckExtractor [™] has to log-in with its own unique user ID.
UV cover	Protective cover with handles which is used during UV-light decontamination.
Worklist	List with samples to be processed

TABLE OF CONTENTS

1.	GE	NERAL		3
	1.1	Intend	ed Use	3
	1.2	Copyr	ight	3
	1.3	Trader	nark	3
	1.4	Terms	used in this Manual	4
		1.4.1	Abbreviations	6
		1.4.2	Name and Description of Components (Overview)	7
2.	SAI	FETY I	NSTRUCTIONS AND PRECAUTIONS	12
	2.1	Gener	al Safety Information	12
	2.2	Warniı	ngs and Precautions	12
		2.2.1	Contamination prevention	. 12
		2.2.2	General precautions	. 12
		2.2.3	Working safely	. 12
	2.3	Instru	nent	13
		2.3.1	Operating the Instrument	. 14
		2.3.2	Biological Safety	. 14
		2.3.3	Results and Data	. 14
3.	SHI	IPMEN	T AND STORAGE	15
4.	WA	STE D	ISPOSAL	16
	4.1	Reage	nts and samples	16
5.	SYS	STEM		
			DESCRIPTION	17
	5.1	Syster	DESCRIPTION	17 17
	5.1	Syster 5.1.1	DESCRIPTION n Overview Extraction	17 17 .18
	5.1	Syster 5.1.1 5.1.2	DESCRIPTION n Overview Extraction PCR Setup	17 17 .18 .18
	5.1 5.2	Syster 5.1.1 5.1.2 Hardw	DESCRIPTION n Overview Extraction PCR Setup	17 17 .18 .18 18
	5.1 5.2	Syster 5.1.1 5.1.2 Hardw 5.2.1	DESCRIPTION n Overview Extraction PCR Setup vare Deck	17 17 .18 .18 18 .18 .19
	5.1 5.2	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2	DESCRIPTION n Overview Extraction PCR Setup vare Deck Pipetting	17 17 .18 .18 18 .19 .19 .19
	5.1	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3	DESCRIPTION m Overview Extraction PCR Setup vare Deck Pipetting Gripper	17 17 18 18 18 18 19 . 19 . 21
	5.1	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4	DESCRIPTION m Overview Extraction PCR Setup vare Deck Pipetting Gripper Heater Shaker	17 17 .18 .18 .19 .19 .21
	5.1	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5	DESCRIPTION	17 17 .18 .18 .18 .19 .21 .21 .22
	5.1	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6	DESCRIPTION m Overview Extraction PCR Setup vare Deck Pipetting Gripper Heater Shaker Barcode Scanner UV-light	17 17 .18 .18 .19 .21 .21 .22 .22
	5.1 5.2 5.3	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa	DESCRIPTION n Overview Extraction PCR Setup rare Deck Pipetting Gripper Heater Shaker Barcode Scanner UV-light are (Graphical User Interface)	17 17 .18 .18 .19 .21 .21 .22 .22 22
	5.1 5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa Specia	DESCRIPTION	 17 17 .18 .18 .18 .19 .21 .21 .22 .24
	5.1 5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa 5.4.1	DESCRIPTION	 17 17 18 18 18 19 21 21 22 22 22 24 24
	5.1 5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa 5.4.1 5.4.1 5.4.2	DESCRIPTION	17 17 .18 .18 .19 .21 .22 .22 22 24 .24 .24
	5.1 5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa 5.4.1 5.4.2 5.4.3	DESCRIPTION	 17 18 18 19 21 21 22 22 22 22 22 24 24 24 28
	5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa 5.4.1 5.4.2 5.4.3 5.4.3 5.4.4	DESCRIPTION	17 .18 .18 .19 .21 .22 .22 22 .22 .22 .24 .24 .24 .28 .29
	5.1 5.2 5.3 5.4	Syster 5.1.1 5.1.2 Hardw 5.2.1 5.2.2 5.2.3 5.2.4 5.2.5 5.2.6 Softwa 5.4.1 5.4.2 5.4.3 5.4.3 5.4.4 5.4.5	DESCRIPTION	 17 18 18 19 21 21 22 23

6.	PR	OCED	URE	32
	6.1	Order	Information	32
	6.2	Requi	red Consumables and Equipment	32
	6.3	Prepa	ration Prior to Installation	33
	6.4	Bringi	ng into Operation	33
	6.5	Install	ation and configuration	34
	6.6	Prepa	ration Prior to Operation	34
		6.6.1	Samples	.34
		6.6.2	Filling the sample drawers	34
		6.6.3	Preparing Elution Plates	.35
		6.6.5	Reagent Barcodes	35
		6.6.6	Liquid Waste	36
7.	INS	TRUC	IONS FOR THE OCHECK [®] WORKFLOW	37
	7.1	Gener	al instructions	37
	7.2	Room	separation	37
	7.3	Sealin	g and storage of plates	38
		7.3.1	Sealing and storage of the elution plate	38
		7.3.2	Sealing of the PCR plate	38
		1.3.3	Reclosing of the PCR plate	39
8.	SH	UTDO	WN PROCEDURE	40
	8.1	Regul	ar Shutdown	40
	8.2	Decor	ntamination	40
	8.3	Using	UV-Light	40
9.	MA	INTEN	IANCE	43
	Q 1	Introd	uction	13
	9.1	Maint	anance Intervals	43
	0.2	Mater	ials required	13
	9.0	Maint	enance Procedures	43
	0.1	9.4.1	Daily Maintenance	.44
		9.4.2	Weekly Maintenance	.45
	9.5	Printin	ng a Maintenance process status Report	48
	9.6	If Mai	ntenance Fails	49
	9.7	Furthe	er Maintenance and Service	49
	9.8	Finish	ing the Maintenance and Cleaning	49
10.	TE	CHNIC	AL SUPPORT	50
11.	TR	OUBLE	ESHOOTING	51
	11.1	Error I	Vessages	51
	11.2	2 Error (Codes	52
	11.3	BCause	es of Trouble	55
	11.4	Assist	ance	55

12.	TECHNIC	AL SPECIFICATIONS	. 56
	12.1 Physic	cal Properties	. 56
	12.1.1	CheckExtractor™	56
	12.1.2	Distances	56
	12.1.3	Enviromental Conditions	56
	12.1.4	Operating Data	56
	12.1.5	Pipetting Specifications	57
	12.1.6	Connectivity	57
	12.1.7	Instrument Life Span	57
	12.2 Barco	des	. 58
	12.2.1	Barcode Specifications	58
	12.2.2	Sample Barcodes	58

2. SAFETY INSTRUCTIONS AND PRECAUTIONS

2.1 General Safety Information

Please make sure that you have read this user manual carefully and pay particular attention to the safety information before using the CheckExtractor[™]. The instructions and safety information in this user manual must be followed to ensure safe operation of the instrument and to maintain the instrument in a safe condition.

Although the CheckExtractor[™] is developed with all care and tested throughout, errors may occur due to various reasons. Greiner Bio-One does not guarantee the obtained results. Any diagnostic result generated, using the CheckExtractor[™] should be interpreted in conjunction with other clinical or laboratory findings.

2.2 Warnings and Precautions

2.2.1 Contamination prevention

- Lab coats must be worn throughout the procedure and different sets of lab coats are required for each laboratory room.
- Lab cleanness: The working place must be decontaminated with DNA-AWAY® or any other appropriate cleaning solution before and after work.
- Gloves must be worn during each step of the analysis and must be changed frequently, especially during DNA extraction.
- **Sample tubes:** Never touch the inside of a reaction tube cap. To avoid cross-contamination, open only one tube at a time.
- **Pipetting:** Appropriate micropipette filter tips with aerosol barriers must be used (sterile, free of DNase, RNase and human DNA). Pipette tips should always be changed between liquid transfers.

2.2.2 General precautions

- Upon arrival, check the kit components for damage. If one of the components is damaged (e.g. buffer bottles), contact your local Greiner Bio-One distributor. Do not use damaged kit components, since their use may lead to poor kit performance.
- Do not use kits after the expiry date.
- Do not use expired reagents.
- Do not mix reagents from different batches.
- Use only reagents/equipment provided with the kit and those recommended by the manufacturer.
- Regular calibration/maintenance should be performed for the CheckExtractor[™] (see chapter MAINTENANCE) and other equipment, e.g. pipettes.
- Pipetting of small amounts of liquid in the microliter range is a challenge. Therefore take care to pipette as accurately as possible.
- To avoid microbial contamination of the reagents, take care when removing aliquots from reagent tubes.
- Unused reagents and waste material must be disposed of in accordance with federal and state guidelines.

2.2.3 Working safely

- Always wear a suitable lab coat, disposable gloves and protective goggles!
- The kits are for in vitro diagnostic use only and should be exclusively used by personnel trained

in in vitro diagnostic laboratory practice.

- Take care whilst handling biological samples containing potential human infectious material. To minimise the risk of infection from potentially infectious material, it is recommended to work under laminar airflow conditions until sample lysis is completed. Handle and dispose all biological samples as if they were capable of transmitting infectious agents.
- Avoid direct contact with the biological samples as well as splashing or spraying. Always wear lab coat, gloves and goggles while working with human samples.
- Never pipette solutions by mouth.
- Do not eat, drink, smoke or apply cosmetic products in the work areas.
- Wash hands carefully after handling of samples and reagents.

2.3 Instrument

The CheckExtractor[™] conforms to European standards with regards to interference immunity. However, if the CheckExtractor[™] is subjected to electromagnetic RF fields or if static electricity is discharged directly onto the CheckExtractor[™], its Liquid Level Detection ability may be adversely affected. It is therefore recommended that the CheckExtractor[™] is kept away from any other equipment that emits electromagnetic RF fields in the laboratory, and that static electricity is minimised in its immediate environment.

During operation, the CheckExtractor[™] must be shielded from direct sunlight and intense artificial light. The instrument should be positioned in the laboratory such that personnel have access to the front and sides of the instrument, for operation, maintenance, opening and removal of protective covers. To calculate how much room is needed, consider the dimensions of the instrument (see chapter TECHNICAL SPECIFICATIONS) and sufficient room for a person to move and work comfortably around it.



Do not expose the CheckExtractor™ to direct sunlight.

Periodic maintenance (daily maintenance, weekly maintenance, six-monthly maintenance) is a mandatory part of the work routine. Discard used tips and do not reuse them. Do not empty the tip waste during a run. Do not try to open the front cover during a run as it may cause the run to abort. This will cause the loss of sample material and reagents, since a new run has to be started.

Never lift a fully installed instrument to transport it from one place to another. It must be re-installed in the new work location by an authorized service technician only. The instrument weights more than 150 kg. Necessary precautions should be taken when transporting the instrument.

Only certified technicians are authorized to perform mechanical maintenance on the CheckExtractor™.

For repair or shipment, all mechanical parts must be put in their rest positions. A CheckExtractor[™] sent away for repair must also be decontaminated (see chapter MAINTENANCE) if it was in a laboratory environment with infectious or hazardous materials. The CheckExtractor[™] must be repacked in the original shipping crate only by an authorised service technician (contact your local Greiner Bio-One representative). There must be no containers or tips on the CheckExtractor[™] during transportation.

The service technician and the laboratory share the responsibility for the installation qualification (IQ) and the operation qualification (OQ), i.e. verification and training. The process qualification (PQ) is the sole responsibility of the laboratory.

Only original Greiner Bio-One CheckExtractor[™] specific parts and tools may be used with the CheckExtractor[™], e.g. carriers, racks, CO-RE tips, and waste containers.

A breakdown of the power supply during a run may cause the loss of data. If data loss is unacceptable, use an independent power supply or an Uninterruptible Power Supply (UPS).

2.3.1 Operating the Instrument

When using the CheckExtractor[™], Good Laboratory Practices (GLP) must be observed. Suitable protective clothing, safety glasses and protective gloves must be worn, particularly when dealing with a malfunction of the instrument where the risk of contamination from spilled liquids exists. Any persons operating the CheckExtractor[™] and the controlling PC, which runs the CheckExtractor[™] software, must have attended a certified training course. Any deviation from the procedures given here could lead to erroneous results or CheckExtractor[™] malfunction.

2.3.2 Biological Safety

The CheckExtractor[™] is designed to prepare medical samples. This implies that the samples may contain infectious material especially viruses. Therefore it is mandatory and the responsibility of the laboratory's manager that only trained personnel work with the CheckExtractor[™] and the associated reagents and samples. All work has to be carried out in accordance to national regulations. Also the disposal of sample material, reagents and used laboratory equipment, like pipette tips and micro plates, needs to follow national regulations.

It is in the responsibility of the laboratory's manager to make sure that the safety and disposal standards are fulfilled.



The sample to be extracted and prepared for the analysis in the CheckExtractor™ are potentially infectious! Handle these samples with extreme care, to keep you and your environment save.

2.3.3 Results and Data

Data loss can be a common problem when dealing with computers. This may be due to technical problems or user errors. Although the CheckExtractor[™] is designed and developed in a way that reduces the risk of data loss to a minimum, data loss can never be fully excluded. It is therefore highly recommended that you back-up your data regularly, in particular if the CheckExtractor[™] is not connected to a LIMS (see chapter BACK-UP)). The back-up of data is in the sole responsibility of the laboratory's manager.



Back-up your data regularly.

3. SHIPMENT AND STORAGE

The CheckExtractor[™], when delivered to a customer, is set up by a technician from Greiner Bio-One or one of its representatives. Everything that is required for a successful installation of the CheckExtractor[™] is carried out by this technician.



Do not unpack and install the CheckExtractor[™] yourself. Special training is needed! Unpacking and installation of the CheckExtractor[™] by any unauthorised personnel will result in loss of warranty!

Nevertheless please check the outside of the packaging for any damage. If you find any damage, please contact your Greiner Bio-One representative immediately.

The CheckExtractor[™] is delivered on two pallets. Store them closely together, to have both readily available for the installation.



Figure 1: CheckExtractor™ when delivered

Parameter	Range
Temperature range	-25 °C – +70 °C
Relative humidity	10 % – 90 % (no condensation, indoors)
Indoor storage only	



Be aware that the CheckExtractor™ is delivered on two pallets with a weight of approx. 260 kg.

The pallets have approximately the dimensions given in the following table:

Dimension	Pallet I	Pallet II
Width	1280 mm	840 mm
Height	1300 mm	760 mm
Depth	900 mm	610 mm
Weight	220 kg	40 kg

4. WASTE DISPOSAL

The CheckExtractor[™] is a piece of electrical equipment and therefore has to be treated as such and ultimately disposed of in compliance with all local laws. Electronic waste has to be collected and separated following the regulations of the electronic waste act.

For users in the European Union the directive 2012/19/EU contributes to sustainable production and consumption by the re-use, recycling and other forms of recovery waste electrical and electronic equipment (WEEE). This reduces the unnecessary disposal of waste and contributes to the efficient use of resources and the retrieval of valuable secondary raw materials. To achieve this goal the members of the EU ensure the correct treatment of all collected WEEE and therefore WEEE is disposed of separately from other items.



Figure 2: WEEE



It is necessary to dispose of WEEE separately.

4.1 Reagents and samples

Guanidine hydrochloride may be in some buffers and can form highly reactive compounds when combined with bleach. Bleach or acidic solutions should not be added directly to the sample preparation waste.

Contamination of the liquid waste with residual infectious material is highly unlikely, but cannot be completely excluded. Thus, liquid waste must be considered infectious and be handled and discarded according to local safety regulations. Please carefully follow federal, state and local guidelines for waste disposal.



Liquid waste has to be considered as potentially infectious material! Dispose this material according local safety regulations!

Further information can be found in the instructions for use of the Greiner Bio-One oCheck[®] DNA Extraction kit - CheckExtractor[™] and the appropriate oCheck[®] kit.

5. SYSTEM DESCRIPTION

5.1 System Overview



Figure 3: The CheckExtractor™

The CheckExtractor[™] is a fully automated system for DNA extraction from patient samples and the subsequent preparation steps necessary for PCR-based diagnostic tests. Especially designed for the use with the oCheck[®] products from Greiner Bio-One, e.g. PapilloCheck[®] and PelvoCheck[®] CT/NG, it offers a maximum in reliability, a high throughput combined with a marginal hands-on time.

The CheckExtractor[™] is based on superior air displacement pipetting technology. This increases accuracy and repeatability while providing chain of custody with pipette condition monitoring and recording. The CheckExtractor[™] is configured with one arm consisting of an 8-pipetting-channel and one gripping device. The auto load option provides barcode tracking of samples, lab ware, racks and carriers.

The purpose of the instrument is to perform the following tasks:

- Sample transfer from primary tubes in a collection block (optional, can be purchased separately on request)
- Deoxyribonucleic acids (DNA) extraction from clinical samples with GBO reagents
- PCR setup of a microplate in association with GBO reagents

The three tasks are executed independently (sample transfer or extraction or PCR setup).

The CheckExtractor[™] is designed to process up to 96 samples per run and optimised for 48 and 96 samples, respectively to process in parallel.

Please refer to the IFU of the corresponding Greiner Bio-One Extraction kit or Greiner Bio-One oCheck[®] kit for further details relating to:

- recommended number of samples
- sample volume
- run time for the extraction and the PCR setup

5.1.1 Extraction

The CheckExtractor[™] is designed for the extraction of viral, bacterial and human genomic DNA from samples of human origin by professionals. By using the oCheck[®] DNA Extraction Kit – CheckExtractor[™] (REF 517 070) a chemical lysis supported by heat treatment is performed. By using the magnetic bead technology the DNA is bound and fixated. After well-defined washing steps to remove cell debris and proteins the purified DNA is eluted from the beads in a suitable solution for the following PCR setup.

5.1.2 PCR Setup

After a successful extraction, the CheckExtractor[™] prepares the PCR reaction mix by adding the DNA eluate to the final PCR MasterMix. Afterwards the prepared and sealed PCR plate can be placed in the PCR cycler to perform the DNA amplification.

5.2 Hardware

The CheckExtractor[™] comprises of two parts: The CheckExtractor[™] itself and the controlling PC. The PC with the graphical user interface (GUI) guides the user through all the steps that require manual intervention when setting up a run. The two devices communicate via a USB connection all necessary data and commands. The CheckExtractor[™] carries out the run and either confirms the successful completion to the PC or reports any errors.



Figure 4: CheckExtractor[™] and the controlling PC

The CheckExtractor[™] performs both pipetting operations and transport of plates placed on its work surface. Pipetting refers to the transfer of small quantities of liquid from one container to another. A pipetting operation is achieved by aspirating (drawing) liquid from a source container, then transferring and dispensing (dropping) it into a target container.

5.2.1 Deck

The CheckExtractor[™] internal work surface is known as the "deck". The deck is where all the carriers are housed. These carriers hold reagent containers, such as tubes, microplates, or other kinds of labware. The deck of the CheckExtractor[™] is divided into equal tracks (T) for loading carriers in predetermined positions. This eliminates the need for precise measurement of positions. The tracks are numbered from left to right, starting with -1.

Figure 5 shows the layout of the deck. It comprises from left to right

- two carriers of tip racks,
- one plate carrier for the PCR-plate, the elution plate, the MASTERBLOCK®, and 10 mL bottles,
- one reagent carrier for the reagent troughs,
- one tube carrier for the reagent tubes
- multiple sample carriers for the sample tubes,
- the position for the heater shaker and the magnetic separator and
- the waste position



Figure 5: Deck layout

The tip racks are divided into three tip regions A, B, and C.

- Region A: 300 µL tips
- Region B: 1000 µL tips, the CheckExtractor[™] will only pick up a complete column of 8 tips at a time from this region.
- Region C: 1000 µL tips, the CheckExtractor[™] pick up single tips here.

The positions on all carriers are numbered from the end that is far inside the CheckExtractor[™] to the front. In other words: The through or plate positions are numbered in the order in which the barcode scanner sees them.

5.2.2 Pipetting

The CheckExtractor[™] is equipped with eight 1000 µL pipetting channels. These pipetting channels work in parallel for simultaneous transfer of liquids and support pipetting with disposable tips. The Dynamic Positioning System (DPS) of the CheckExtractor[™] moves each pipetting channel independently on the y-Axis, as well as the z-Axis. Each channel uses its own high-precision motors

and electronics to reach any position on the deck.



Figure 6: Pipetting tips

The CheckExtractor[™] is based on the air displacement pipetting principle. Air displacement means that the liquid is aspirated into and dispensed from a disposable tip or needle by the movement of a plunger. Between the plunger and the liquid surface is air. No system liquid of any kind is involved in the pipetting procedure.



Figure 7: Air displacement technology

The first task for the CheckExtractor[™] Pipettor is to pick up a disposable tip. Due to the unique compression-induced O-ring expansion (CO-RE) technology, precise tip attachment and positioning is achieved. The system requires no vertical force for tip attachment or tip ejection, thus eliminating mechanical stress and improving overall system reliability along with pipetting speed and capability.



Figure 8: Figure 7: Co-Re technology

The principle of the CO-RE technology has the following advantages:

- Allowing different sizes of tips on the same pipetting head in the same run
- Picking up a gripper and other tools
- Eliminating aerosol production upon tip ejection



Do not use any tip types other than the Hamilton CO-RE tips. Other tips may not be mounted correctly and may lead to erroneous results.

5.2.3 Gripper

The CO-RE Gripper is the plate handling tool that is picked up by two pipetting channels during a run. The CO-RE Gripper transports the different plates (MASTERBLOCK[®], Collection Block and Elution plate) to and from positions on the deck of the CheckExtractor[™].



Figure 9: CO-RE Gripper

5.2.4 Heater Shaker

The Heater Shaker (HHS) is designed to heat and shake microplates that conform to the SBS format. Loading and unloading as well as the independent heating and shaking function of the HHS is under control of the CheckExtractor[™] software.

Before shaking, the microplates are locked into position in the center of the HHS. When the shaking process has finished, the microplates are unlocked and can be easily removed from the heater shaker using the CO-RE Gripper.



Figure 10: Heater Shaker (HHS)

5.2.5 Barcode Scanner

The red instrument on the right side in the housing is the barcode scanner. Every time the autoload loads a drawer into the CheckExtractor[™], the barcode scanner is moved into position, adjacent to the drawer. Since all disposables and samples in the drawer are labelled with a barcode, the barcode scanner registers all these components.

This allows the CheckExtractor™:

- Confirmation of the presence and correct positioning of all reagents.
- Confirmation of the number and position of all samples throughout the extraction process.
- Complete traceability of all the samples through the extraction process.
- Monitoring of the PCR and elution plates
- Mapping of the sample position from plate to plate.

5.2.6 UV-light

The CheckExtractor[™] is equipped with a UV-light. This UV-light helps to decontaminate the desk. The UV-light bulb is situated inside the housing. The housing itself is UV-resistant to protect the user and the environment against UV radiation. Please use the UV-light during shutdown (see chapter SHUTDOWN PROCEDURE) or maintenance.

5.3 Software (Graphical User Interface)

The Graphical User Interface (GUI) has been designed to guide the user through the entire process, is easy to understand and self-explanatory.

To start the GUI just a double click on the CheckExtractor[™] icon on the desktop is necessary.



Figure 11: CheckExtractor™ icon on desktop

After the double click on the desktop icon the GUI start screen appears where you can choose which program (Sample Transfer (optional), Extraction or PCR setup) the CheckExtractor[™] should perform. The GUI guides you through the information needed to start the run.



Figure 12: GUI start screen

In order to avoid false operating several error messages are part of the GUI and appear in case of an erroneous handling. For explanations and recommended steps please refer to chapter TROUBLESHOOTING.



Figure 13: Software overview

5.4 Special Functions

5.4.1 User Management

The onboard software of the CheckExtractor[™] allows access to features at different levels. For routine use of the CheckExtractor[™], the user level is defined as "Lab Operator". All other user levels are for service purposes only and only accessible by Greiner Bio-One representatives.

Group	User	Authorisation
Lab Operator	Routine User	Operators may run any method

The adjustment of the Windows[®] user settings on the controlling PC will be within the responsibility of the system administrator. Please ask your system administrator to add all user in the lab to the group "Lab Operator".

5.4.1.1 Protection of PCR config files

Users who have been assigned as "Lab Operator" will only have 'read' access to these files and will not be able to edit or delete them.

5.4.1.2 Protection of method and library files

Users who have been assigned as "Lab Operator" will only have 'read' access to these files and will not be able to edit or delete them.

5.4.1.3 Protection of output files

The output files (see chapter DATA MANAGEMENT) will be protected by Windows[®] access rights. Users who have been assigned as "Lab Operator" will only have 'read' access to these files and will not be able to edit or delete them.

5.4.2 Input files

For every run a worklist is required. This worklist contains all necessary information about the samples to be run on the CheckExtractor[™] and for the tracking of the samples e.g. with a LIMS.

5.4.2.1 Input file for Extraction Method

There are three common ways how the CheckExtractor[™] creates the worklist:

- 1. An input file is provided by the LIMS. The user just loads it from the LIMS directory using the "Import" button. The CheckExtractor™ converts the import file into a worklist. Please contact your LIMS provider to create a input file in the form, depicted below.
- 2. An input file is manually written into an Excel[®] file. This Excel[®] file is then loaded using the "Import" button. The CheckExtractor[™] converts the import file into a worklist.
- 3. It is also possible, although not recommended, to create a worklist at the beginning of the run. Instead of using the "Import" button, click in the first sample field and type in the sample ID and the barcode in the respective columns. The barcode can although be scanned with the handheld barcode scanner.

If the input file is created manually, please observe the following instructions. (Your LIMS has to follow the instructions as well.) A template of the Excel[®] file is already installed on the controlling PC. Nevertheless, in the case this template is missing, the following instructions guide to a functional input file.

All data is treated as text in the worklist spreadsheet in Excel[®] so ensure all cells are formatted as text. This is especially important, when sample ID or barcodes start with zeroes (e.g. "0012345").



Make sure, all columns in the Excel® file are formatted as "Text".

To enable this text formatting:

1. Select all cells by either clicking in the topmost left corner of the spreadsheet, i.e. the area above the row numbers and left to the column headings, or press on the keyboard simultaneously the "Ctrl" key and the "A" key. All cells are highlighted or coloured (see Figure 14).

🗶 🛃 🎝 🔹 (* 🗉 📼					
File Home In	nsert Pag	e Layout Fo	ormulas Dat	a Review	View
Cut	Calibri	* 11	· A A	= = =	≫~ 📑 Wr
Paste V Format Painte	B I	<u>u</u> • <u>m</u> •	<u>ð</u> - <u>A</u> -	₽ ₹ ₹	岸 📳 🔤 Me
Clipboard	r ₂	Font	G		Alignment
A1	• (*	f_{x}			
A 4	В	С	D	E	F
1					
2					
3					
4					
5					
6					
7					
8					
0					

Figure 14: Select all cells.

2. Click the right mouse button on the heading of the first column, i.e. the grey area marked "A". A menu appears. Select the item "Format Cells" (see Figure 15).



Figure 15: Select the item to format the cells as text.

3. When the item "Format cells" is selected a window appears. Activate the first tab named "Number" and select in the list "Category" the item "Text" (see Figure 16).

ormat Cell	s						8	X
Number	Alignment	Font	Border	Fill	Protection			
Category: General Number Currency Accountin Date Time Percenta Special Custom	ge	Text fc The cel	ie rmat cells a lis displaye	re treated d exactly a	as text even w s entered.	ihen a numbe	r is in the	cell.
					ſ	ОК	Can	cel

Figure 16: Format cells as "Text".

4. Finish the steps by clicking the "OK" button.

In the next step rename the first (and active) spreadsheet. It is mandatory to name this spreadsheet "Samples". To rename the spreadsheet double-click its name in bottom-left corner. The name, it is normally "Sheet1", is highlighted. Type in "Samples" and confirm your input by clicking in an arbitrary cell in the spreadsheet (Figure 17).

19	
20	
21	
22	
23	
24	
25	
26	
27	
I4 4	🕨 🕅 Samples / Sheet2 / Sheet3 / 💱 /

Figure 17: Rename the Spreadsheet into "Samples".

The spreadsheet is now correctly formatted for your input file. The first row of the input file has to contain the four entries (see Figure 18):

SampleID SampleBarcode SampleAtrix

Under these titles the information for each sample is given. While the columns "Description" and "SampleMatrix" are optional, the first two columns "SampleID" and "SampleBarcode" are compulsory, they have to be filled in.

X	9.6	(* #	2.91 W 2.91						-			Extraction	Vorklist48Sam	plesEXT - Mi	crosoft Excel										
- Fé	Home	Inse	ert Page Layout	Formula:	D	ata Revie	w Vie																	۵ 🕜 ۱	- 8 2
	X Cut		Calibri *	11 - /	۲.	= = =	\$>	🐨 Wrap Text	Text				Normal	Bad		Good	Neutral	Calcula	tion	🜁 ች		E AutoSum * A	A		
Past	Copy ·	Dainter	B / U · 🖽	- 3-	Α.		课课	Merge & Center		% · % 4	Condition	l Format	Check Ce	II Exp	lanatory	Input	Linked Cell	Note		Insert Delete	Format	Sort	& Find &		STATES .
	Clipboard	14	Font		- 14		Aligner	ent	15 N	lumber /	Formatting	· as rapie		The second	Styl	es				Cells		Editing	a. resea.		1000
	836		• (* fe 41	21300010	060																				*
	A		В			С		D	E	F	(6	Н	1	J	K	L	М	N	0	Р	Q	R	S	-
1	Samplel	D	SampleBarcode	2	Desi	ription		SampleMatri	x																
2	sample0	01	412130001002	22																					
3	sample0	02	412130001002	23																					
4	sample0	03	412130001002	24																					
5	sample0	04	412130001002	26																					
6	sample0	05	412130001002	27																					
7	sample0	06	412130001002	28																					
8	sample0	07	412130001002	29																					
9	sample0	08	412130001003	30																					
10	sample0	09	412130001003	31																					
11	sample0	10	412130001003	32																					
12	sample0	11	412130001003	33																					
13	sample0	12	412130001003	35																					
14	sample0	13	412130001003	36																					
15	sample0	14	412130001003	37																					
16	sample0	15	412130001003	38																					
17	sample0	16	412130001003	39																					
18	sample0	17	412130001004	11																					
19	sample0	18	412130001004	12																					
20	sample0	19	412130001004	13																					
21	sample0	20	412130001004	14																					
22	sample0	21	412130001004	15																					
23	sample0	22	412130001004	16																					
24	sample0	23	412130001004	17																					
25	sample0	24	412130001004	18																					
26	sample0	25	412130001004	19																					
27	sample0	26	412130001005	50																					
28	sample0	27	412130001005	52																					
29	sample0	28	412130001005	53																					
30	sample0	29	412130001005	54																					
31	sample0	30	412130001005	55																					
22	> H Shee	ti sh	eet2 / Sheet3	1												1							1		•
Read	v	State State						and the second second second													and the second	16	130%		

Figure 18: Example for a input file.

When finished save the input file on the controlling PC. It is necessary to save the input file in the Excel® 97-2003 format, i.e. with the extension ".xls". Do not use the Excel® 2010 format (.xlsx).



The input file needs to be saved in Excel® 97-2003 format (.xls)! Do NOT save input files in the Excel® 2010 format (.xlsx)!

Therefore always use the "Save As" option in the "File" menu. In the dialog click on the drop down menu named "Save as type". Choose the item "Excel 97-2003 Workbook".

Save As			
N	ew folder	· · · · · · · · · · · · · · · · · · ·	Search New folder
Organize 👻 Ne	sw folder		
K Microsoft Excel		No items match your search.	
🗙 Favorites			
💐 Libraries			
Computer			
Vetwork			
File name:	ExtractionWorklist48SamplesEXT		
Save as type:	Excel 97-2003 Workbook		
Authors	Excel Workbook		
Authors:	Excel Macro-Enabled Workbook		
Tala.	Excel Binary Workbook		
Title	Excel 97-2005 Workbook		
Subject	Single File Web Page		
Manager:	Web Page		
Company:	Excel Template		
Categories:	Excel Macro-Enabled Template		
Comments:	Excel 97-2003 Template		
	Text (Tab delimited)		
	VMI Spreadchest 2003		
	Microsoft Excel 5.0/95 Workbook		
	CSV (Comma delimited)		
Hide Folders	Formatted Text (Space delimited)		
	Text (Macintosh)		

Figure 19: Save the input file in the Excel® 97-2003 format.

It might be helpful to save data in the Excel[®] 97-2003 format by default. This can be achieved by: "File -> Options -> Save -> Save file in this format". For detailed information please refer to your Microsoft[®] Excel[®] manual.

In the combo box named "File name" type in the name of the input file. This name is up to you. The filename is only limited by the parameters set by the Microsoft[®] Windows[®] operating system¹ and must end with ".EXT" (starting with a dot "." and all letters in upper case).

 $^{^{1}}$ For example: The characters <, >, :, ,, /, \, |, ?, * must not be used.



The file name of the input file for the extraction method must end with ".EXT".

1

An example for the file name of a input file would be: ExtractionWorklist48Samples.EXT.xls

5.4.2.2 Input file for PCR setup method

A worklist is also needed for the PCR setup method. This worklist is created by the CheckExtractor[™] at the start of the PCR setup method. The worklist again is based on an input file, that is again a simple Excel[®] file. In contrast to the extraction method, the input file of the PCR setup method is created by the CheckExtractor[™] itself.



The file name of the input file for the PCR setup method must end with ".SMP".

I

An example for the file name of a input file would be: 2013-12-12_18-07-30-330_E13070BF000561_GBO_PCRSampleList.SMP.xls

At the end of every extraction run the input file for the corresponding PCR setup is stored in the output directory (see chapter DATA MANAGEMENT). Likewise the PCR setup method creates a input file for a subsequent PCR setup run, i.e. if you choose to re-run a elution plate the input file is already stored.

5.4.3 Tip Editor

During every tip loading step the tip editor window comes up (see Figure 20). For each tip region (A, B, and C) (see Figure 5 Deck Layout) the tip editor appears separately. The tips which are expected to be present on the deck are shown in brown. Normally this situation is correct, since the CheckExtractor[™] tracks which tips have been used. Nevertheless check the tip situation before every run to avoid an abort of the run due to insufficient tips.



Figure 20: Tip Editor

If the tip positions displayed on the screen do not correspond with the real situation on the deck, you can edit the CheckExtractor[™] software using one of the following methods:

- To add or remove tips click on the tip position in the tip editor. The colour of the tip toggles between brown and white, i.e. "tip in place" to "no tip in place".
- To toggle more tips at once position the mouse pointer besides the tips, press the left mouse button and drag the mouse pointer. Release the mouse button to toggle the marked tips.
- To remove all tips in the tip editor press the "Remove All" button.
- You can zoom in or out by turning the mouse wheel.

When the marked tips in the tip editor corresponds with the real situation, click the "OK" button.



Always load full tip racks only.

5.4.4 Manual Barcodes Entry

In the unlikely event that a barcode is destroyed and not readable anymore it is possible to manually type in the barcode. To do so position the cursor by clicking in the box where the barcode needs to be. A window comes up with two entry boxes (see Figure 21). Type in the barcode in the upper box and press the return key on your keyboard. For confirmation type in the barcode again in the lower box. While you type the barcode in the lower box, the barcode in the upper box will be hidden. This double blind entry makes sure that no typing errors will corrupt the database of the LIMS.

Position 1 Bar	Error code Error		Assigned recove	ery
sition description: Assigned barcode: "	Insert Barcode [Po Enter barcode Confirm barcode	sition 1]	Help	
Repeat Continue			Barcode Exclude	Unload carrier

Figure 21: Popup window for manual barcode entry

5.4.5 Data Management

For each run, the CheckExtractor[™] generates different output files, which contain all the information that relates to that particular run. All the files from one run are stored in one individual subdirectory. And all these subdirectories are stored collectively in the directory: C:\MethodOutput\

5.4.5.1 Extraction method

The subdirectories in this folder are named with the date and time of the run. Additionally the barcode of the elution plate is added.

For example 2013-11-21 10-26-25-067 E1012ELF000124

is a subdirectory of a run started on the 21st of November, 2013 at 10:26:25,067. The eluate can be found in the elution plate E1012ELF000124.

The directory contains files whose names begin with the same information as the directory they are contained within i.e. date, time and elution plate barcode. Each file will then have a different suffix based on the type of file and the information it relates to. The file names end with one of the following endings:

Name	Meaning
GBO_EXT_General.xls	General information about the extraction run (Excel® file)
GBO_EXT_PipettingTracking.xls	Information about all pipetting steps, i.e. volumes, source and target position etc. (Excel $^{\otimes}$ file)
GBO_EXT_Samples.xls	Description of the pipetted samples (Excel® file)
GBO_PCRSampleList.SMP.xls	The extraction run already creates the input file for the PCR setup method, it is stored in this file.
GBO_EXT_ConsolidatedReport.pdf	Information about sample position and result of extraction

5.4.5.2 PCR Setup Method

For each run of the PCR setup method, the CheckExtractor[™] creates output files. These files are stored in the same subdirectory as the files which have been created during the extraction method. So, for the above example this would again be the subdirectory:

C:\MethodOutput\2013-11-21_10-26-25-067_E1012ELF000124

Similar to the extraction method the file names created by the PCR setup method all start with the date and time of the run of the PCR setup method, followed by the barcode of the PCR plate. The file names end with one of the following endings:

Name	Meaning
GBO_PCR_General.xls	General information about the PCR run (Excel® file)
GBO_PCR_Report.xls	Information about all pipetting steps, i.e. volumes, source and target position etc. (Excel $^{\circ}$ file)
GBO_PCRSampleList.SMP.xls	The PCR setup method creates the input file for the next run of PCR setup method, it is stored in this file.
GBO_PCR_Worklist.APDWL.xls	Internal worklist
GBO_PCR_ConsolidatedReport.pdf	Information about sample position and result of PCR setup

5.4.5.3 Trace files

Additionally you will find trace files in the subdirectory. For every run, independent whether extraction method or PCR setup method, a new trace file is created. These trace files contain low level information about every step the CheckExtractor[™] carried out during the run. These files are stored for service purposes only. They do not contain any information necessary for daily work.

5.4.6 Back-up

For the safety of your data it is highly recommended to back-up your data regularly. The CheckExtractor[™] does not carry out data back-ups. It is the sole responsibility of the user or lab personnel to initialise the back-up process.



It is the sole responsibility of the user to back-up all data on a regular basis.

Please use the functionality of the Microsoft[®] Windows operating system or your LIMS for the backup. Save all the files discussed in chapter DATA MANAGEMENT for a complete back-up (including the trace files). I.e. back-up the complete directory: C:\MethodOutput\

6. PROCEDURE

6.1 Order Information

Equipment	Greiner Bio-One CatNo.	Quantity
CheckExtractor™	865 070	1
Line cord		2
USB cable		1
PC		1
Keyboard		1
Mouse		1
Hand held barcode scanner		1

6.2 Required Consumables and Equipment

Consumables	Greiner Bio-One CatNo.	Quantity/box
oCheck [®] DNA Extraction Kit for CheckExtractor ^{TM}	517 070	for 288 analyses
Papillo Check ®	465088	for 288 analyses
PelvoCheck [®] CT/NG	504288	for 288 analyses
1000 μL pipette tips	865 806	3840
300 μL pipette tips	865 807	5760
Sample tubes	481 478	1000
Caps for the sample tubes	481 479	1000
Reagent troughs 60 mL	865 808	28
MASTERBLOCK®	780 270	50
PCR and elution plate	652 290-CEX	40
Pierce Seal	865 804	100
MICROLAB™ Detergent & Disinfectant	865 809	1
MICROLAB™ Disinfectant Spray	865 810	1
Waste bag	865 815	25

Equipment	Greiner Bio-One CatNo.	Quantity/box
Magnetic separator	865 801	1
Heat sealer	865 802	1
Adapter for heat sealer	865 803	1
Plate centrifuge	865 805	1

6.3 Preparation Prior to Installation

Before the CheckExtractor[™] is delivered, make sure you consider the physical and operational properties of the instrument (see chapter TECHNICAL SPECIFICATIONS) and the requirements listed here. Choose a place for the CheckExtractor[™] in the lab according to these specifications.



Since the CheckExtractor[™] will be connected to a power line, make sure that an earthed electrical socket is within 1.5 m of the installation site.



Be sure to have a magnetic separator on hand when the CheckExtractor™ is set-up (see chapter ORDER INFORMATION).

The CheckExtractor[™] requires a minimum width of 800 mm to pass through any opening. So, any doors that the unit will pass through on delivery and installation will require this dimension.

6.4 Bringing into Operation

Liquid waste will arise during a run. Therefore the CheckExtractor[™] is equipped with a liquid waste bottle. This bottle is connected to the liquid waste via a hose. As no pump is fitted to the hose the liquid waste bottle has to be placed underneath the CheckExtractor[™] (see Figure 4).

Prior to each run make sure that:

- Only trained personnel are operating the CheckExtractor[™].
- For special training contact your local Greiner Bio-One representative.
- All maintenance steps have been carried out according to chapter MAINTENANCE.
- Only Greiner Bio-One Extraction kits, Greiner Bio-One oCheck[®] kits and approved accessories are loaded into the Check Extractor[™] (see chapter ORDER INFORMATION).
- Samples are put into the CheckExtractor[™] in the defined sample tubes (see chapter ORDER INFORMATION).

Follow the instructions of the relevant Greiner Bio-One Extraction kit and Greiner Bio-One oCheck® kit.



The CheckExtractor[™] does not include any controls. Positive and negative controls must be inserted in the runs like ordinary samples.

6.5 Installation and configuration

The CheckExtractor[™], when delivered to a customer, is set up by a technician from Greiner Bio-One or one of its representatives. Everything that is required for a successful installation of the CheckExtractor[™] is carried out by this technician.



Do not unpack and install the CheckExtractor[™] yourself. Special training is needed! Unpacking and installation of the CheckExtractor[™] by any unauthorised personnel will result in loss of warranty!

6.6 Preparation Prior to Operation

6.6.1 Samples

6.6.1.1 Sample Barcode

Each sample needs to be equipped with a unique identification code. It is recommended to print this unique code in form of a linear barcode. For the requirements of the barcode please refer to chapter BARCODES. Attach this barcode to the corresponding sample tube. The barcode must be aligned in axial direction, i.e. along the length of the tube.



The barcode is constricted to 20 alpha-numeric characters.

6.6.2 Filling the sample drawers

6.6.2.1 Preparing the samples

During the extraction method the CheckExtractor[™] pipettes the sample from the sample tube into the MASTERBLOCK[®]. It is therefore necessary that the samples are easy to pipette. Please follow the instructions in the oCheck[®] DNA Extraction Kit - CheckExtractor[™] IFU for preparing the samples prior use.



Avoid the transfer of cell-clots, since these may plug the pipette tips.

6.6.2.2 Setting-up the samples in the drawer

The system pipettes the samples starting with the leftmost sample drawer (track 20). It takes the samples from the rear end of the drawer, i.e. the sample deepest inside the CheckExtractor[™] is taken first.

Since the scanner is placed on the right side of the sample drawer, all sample barcodes must be turned to the right opening in the drawer positions.



Be aware that the CheckExtractor[™] is able to process any number of samples between 1 and 96. Nevertheless the oCheck[®] DNA Extraction kit - CheckExtractor[™] allows only 48 or 96 extractions per run. So numbers other than 48 or 96 will cause loss of reagent material.

6.6.3 Preparing Elution Plates

For the PCR setup method the CheckExtractor[™] takes an elution plate with DNA material extracted previously. If the extraction run was carried out right before the PCR setup, the elution plate may still be on the plate carrier. In this case no further steps are required.

Alternatively the DNA may have been extracted and stored at +4 °C or in a freezer. In this situation, make sure that the contents of the elution plate is fully thawed and spun down. Open the plate, no sealing film or cover must be left on the plate.



Be sure that the elution plate is fully thawed, spun down and open, before loading it into the CheckExtractor™.

6.6.4 Reagent Barcodes

All reagents have to be barcoded. For details please refer to the IFU of the corresponding extraction kit and oCheck[®] kit, respectively.

Attach the barcode labels to the reagent troughs before starting the run. To do so, place all needed troughs into the trough carrier. Hold the trough carrier in a way that the rounded edge of the troughs faces to the left side. Then place the barcode label vertically in the middle of the troughs. Be sure to attach the barcodes no closer than 20 mm and 10 mm from the edges of the trough as drawn in Figure 22.

Each trough corresponds to a defined reagent. Attach the labels in the right order. The right order is indicated on the screen, when loading the reagent troughs. You may also refer to the IFU of the extraction kit.







Make sure that the barcode is orientated vertically on the trough and is placed inside the borders indicated (see Figure 22).

6.6.5 Magnetic Separator

The magnetic separator must be in place before starting the run. Place the magnetic separator on the holder right in front of the heater shaker (HHS). Please refer to Figure 5 for the exact position of the magnetic separator.

6.6.6 Liquid Waste

Since liquid waste is generated during the run, it is necessary to connect the CheckExtractor[™] to the liquid waste bottle.

Empty the liquid waste bottle prior to a run, before you connect it with the CheckExtractor[™]. Follow the local regulations for disposing liquid waste (see also chapter Waste Disposal). The volume of the liquid waste bottle should last for a few runs. The CheckExtractor[™] reminds you to empty the bottle if not enough volume is left for the next run.



Liquid waste has to be considered as potentially infectious material! Dispose this material according local safety regulations!

Place the lid on the liquid waste bottle. Be sure that the hose and the sensor cable are properly attached to the lid and the CheckExtractor[™]. If there are doubts about the connections of the hose or the sensor cable call your Greiner Bio-One representative!

7. INSTRUCIONS FOR THE oCHECK® WORKFLOW

For the performance this instruction for use is designed to individually choose the applications you need or you prefer. The general part of the instruction for use ends here and in the following supplements the respective application procedures are described in detail.

7.1 General instructions

When implementing currently used state-of-the-art techniques in molecular biology into a laboratory, the following instructions must be considered to ensure both maximum safety for laboratory staff and high quality results.

Execution of molecular biology techniques such as DNA extraction, amplification and detection of the amplification products require appropriately qualified personnel. In addition, a clean and well structured workflow is required to prevent erroneous results, such as those occurring due to DNA degradation or contamination by amplification products. To ensure this, it is necessary to separate the areas of extraction, amplification and detection as described in Chapter 8.1.2.

Each area should be equipped with separate equipment, consumables, lab coats and gloves. Never transfer lab coats, gloves or equipment from one distinct area to another.

7.2 Room separation

It is absolutely necessary to separate the areas of DNA extraction, amplification, and hybridisation. Sample material must not be introduced in the area for setting up the final PCR MasterMix.



Figure 23 shows an example of how a laboratory may be separated into 3 distinct sections.

Figure 23: Room separation

Room 1 is designated for DNA extraction with the CheckExtractor[™] and the processing of the PCR plate. The setup of the PCR reaction mix should be carried out in room 2, the preparation is optimally performed in a protected surrounding, e.g. a PCR hood, to avoid reaction contamination. The prepared MasterMix is transferred back to the CheckExtractor[™] in room 1 for performing the PCR setup. In room 3 the hybridisation, washing, drying and the analysis of the processed oCheck[®] chip using the CheckScanner[™] and CheckReport[™]Software is performed.

Each room is used exclusively for the application or technique indicated to prevent sample contamination. The use of colour coding could be advantageous in avoiding the accidental exchange of equipment and consumables between areas.

7.3 Sealing and storage of plates

In this chapter the sealing of the plates is described in order to store them at \leq -20 °C.



Nevertheless, it is recommended to perform the PCR immediately after the extraction.

7.3.1 Sealing and storage of the elution plate

For storage at \leq -20 °C, the elution plate has to be closed after unloading from the CheckExtractorTM.

 Close the elution plate after removal from the device by attaching a Greiner Bio-One adhesive foil (Silver seal, 80.0 / 140 MM, REF 676 090). Press it firmly onto the surface of the plate. Store the plate at ≤ -20 °C.

It is recommended to use the Greiner Bio-One adhesive foil (Silver seal, REF 676 090) for storage of the elution plate as it can easily be removed from the plate if the plate has to be reused for a PCR setup using the CheckExtractor[™].

But it is also possible to use the sealing procedure described for the PCR plate to close the elution plate.

7.3.2 Sealing of the PCR plate

After the unloading procedure of the CheckExtractor[™], the PCR plate must be sealed before transferring it to the PCR cycler. For this purpose, a sealing foil (Pierce Seal, REF 865 804) is welded onto the PCR plate using a heat sealer (Model: 4S3, REF 865 802). This procedure is described in Figure 24.

- Turn on the heat sealer using a switch at the rear part of the housing. Turn it on about 5 minutes prior to use to ensure that the required sealing temperature of 170 °C can be reached.
- Check the programmed sealing temperature and duration of the sealing procedure, by pressing the "SET" key (1x: the sealing temperature is displayed; 2x: the duration of the sealing procedure is displayed). For the correct sealing of the plate, a sealing temperature of 170 °C and a sealing duration of t2.0 seconds must be set. If one of these values or both are not correct, please adjust the values using the "▲" and "▼" keys.



In order to adjust the sealing temperature, press the "SET" key in the top left corner of the front of the device. The currently set sealing temperature appears on the display and flashes. The required temperature can then be set using the " \blacktriangle " and " \blacktriangledown " keys. By pressing again, the last visible sealing temperature is confirmed and the currently set duration of the sealing procedure is displayed (also flashing). The duration can now also be adjusted using the " \bigstar " and " \blacktriangledown " keys. By pressing again the "SET" key, the duration time is confirmed and stored on the device.

- Wait until the device has reached a temperature of 170 °C (the actual temperature is constantly displayed). In addition to the display, this is also identifiable by a beep sound that is heard when the device is ready for use.
- Open the device by pressing the "OPERATE" key. The tray with insert for the PCR plate will open up automatically.
- Remove the metal frame and put the PCR plate onto the insert of the tray. Then place the sealing foil (Pierce Seal, REF 865 804) onto the PCR plate. Make sure that the blue line is facing side up

and that all wells are well covered. Place the metal frame on the top of the plate.

- Close the device by pressing again the "OPERATE" key. The tray will then close automatically and the sealing procedure starts.
- After the sealing procedure the tray opens again automatically. First, remove the metal frame and then the sealed PCR plate.

After removal of the PCR plate, the device can be closed by pressing the "CLOSE" key. To avoid bending the metal frame it is recommended to keep the frame inside the device. When closed, the device can be turned off anytime, using the power switch at the rear part of the housing. Turn on the heat sealer 5 minutes prior to use. Use sealing temperature of 170 °C and duration of t2.0 seconds. Press "OPERATE" to open device. Remove the metal frame and put the PCR plate onto the insert of the tray. Place sealing foil (Pierce Seal, REF 865 804) onto the PCR plate. Place metal frame on top of the plate. Press "OPERATE". Remove sealed plate, press "CLOSE". Turn off sealer.

Figure 24: Sealing procedure using a heat sealer

7.3.3 Reclosing of the PCR plate

After performing the hybridisation, the PCR plate can be resealed and stored again at \leq -20 °C.

 Close the PCR plate by attaching a Greiner Bio-One adhesive foil (Silver seal, REF 676 090). Press it firmly onto the surface of the plate. Store the plate at ≤ -20 °C.



The adhesive foil is stuck directly onto the sealing film. When removing the foil, it is possible that the underlying sealing foil may partially detach from the plate. However, this is not a problem for taking out PCR products from the PCR plate.

8. SHUTDOWN PROCEDURE

8.1 Regular Shutdown

At the end of each day it is recommended to shut down the CheckExtractor™.



Make sure all reagents, samples and plates are removed from the deck!

8.2 Decontamination

Following is the recommended procedure for decontaminating of the CheckExtractor™:

- Spray the front and side cover with MICROLAB[™] Disinfectant Spray (see chapter ORDER INFORMATION).
- Open the front cover and spray the deck with MICROLAB[™] Disinfectant Spray (see chapter ORDER INFORMATION) or wipe it clean with a soaked cloth.
- Remove the tip eject plate of the tip waste station and clean it.
- Spray MICROLAB™ Disinfectant Spray directly onto the surface of the tip waste station.
- Remove the frame that holds the plastic bag in place, and discard the plastic bag in the laboratory's contaminated waste. Put the tip eject plate back in place.
- Spray or wipe out the tipwaste drawer.
- Clean the tip eject sleeve (outer part of the pipetting channels) with a lint-free cloth soaked in MICROLAB[™] Disinfectant Spray (see chapter ORDER INFORMATION).
- Clean all carriers with MICROLAB[™] Disinfectant Spray (see chapter ORDER INFORMATION) and let them dry. If they are heavily soiled, soak them afterwards in a solution of MICROLAB[™] Detergent and Disinfectant (see the product data sheet for further information).
- Use the UV-light.

8.3 Using UV-Light

To assist the decontamination processes the CheckExtractor[™] is equipped with an UV-light, since UV-light destroys DNA. It therefore makes sense to decontaminate the deck by using UV-light at the end of each day. Make sure the CheckExtractor[™] is unloaded after the last run. The tips may stay on the deck.



UV-light may be harmful to the reagents. Therefore please remove all reagents before using the UV-light.

To start the decontamination of the deck using UV-light, close the front cover with the UV cover. Place it in the opening of the front cover by holding the UV cover horizontally by the handles and with the two catches towards the CheckExtractor[™]. Hook the catches under the front cover. Move the lower edge of the UV cover downwards and towards the CheckExtractor[™] until closes and seals the front of the unit.



UV-light is harmful to your eyes and skin! Never start the UV-light before the front cover is tightly closed with the UV cover.

1	Double-click the icon for the UV-light which lays on the desktop.						
2	A warning reminds you to close the UV cover and to wear suitable clothing. Follow the in- structions and press the "OK" button. To test whether the UV-light bulb works, the UV-light will turn on for approx. one second.						
	Warning!						
	During the check of the UV light: Avoid exposure to UV light. Exposure may cause severe eye or skin injury.						
	Wear suitable face shield, gloves and protective clothing. Make sure that the UV cover is placed in the correct position to completely close the instrument before switching on the UV lights.						
	Press OK to test the UV light!						
	OK Cancel						
3	Please confirm with "Yes" if the UV-light is working. If no press "Cancel" and call your Greiner Bio-One representative for further instructions.						
	Check UV light result confirmation.						
	Have you seen the UV light flashing?						
	Yes: UV light is working. No: Retest the UV light.						
	Cancel: Aborts the run.						
	Yes No Cancel						
4	You will then be asked for the duration of the decontamination cycle and if any delay should be applied. Insert these two times in minutes and press the "OK" button. A decontamination time of 30 min is recommended.						
	UV light timer: Please set the decontamination duration and delay timer for the UV light (in minutes).						
	Prompt Value Minimum Maximum Decontamination duration (minutes) 30 1 60 Time to wait before decontamination starts 0 0 10						
	ОК Нер						

5 ^B m T	efore s nove fre hen sta	tarting the UV-light make sure that the pipetting arm of the CheckExtractor™ can eely. Doublecheck the tightness of the UV cover and your personal UV protection. art the UV-light by clicking the "OK" button.
	Warning!	
		During the UV decontamination process the pipetting arm will move. Please ensure that no collision can take place. Avoid exposure to UV light. Exposure may cause severe eye or skin injury. Wear suitable face shield, gloves and protective clothing. Make sure that the UV cover is placed in the correct position to completely close the instrument before switching on the UV lights. OK Cancel

9. MAINTENANCE

9.1 Introduction

The CheckExtractor[™] is designed to require a minimum amount of maintenance by the user. However, to preserve the accuracy and reliability of the instrument, the following preventive maintenance procedures must be performed on a regular basis according to the maintenance schedules given below. Periodic maintenance routines need to be run in order to ensure safe and reliable operation of the CheckExtractor[™] and its accessories.

Since the CheckExtractor[™] is an IVD accessory the maintenance routines are mandatory.

9.2 Maintenance Intervals

In order to maintain the CheckExtractor[™] in good working condition, the following maintenance intervals should be observed:

Daily: Recommended before CheckExtractor™ start-up or shut-down

Weekly: Recommended at the end of the week before CheckExtractor™ shut-down Six-monthly: Preventive service maintenance carried out by a service technician



If the operator decides not to run either daily or weekly maintenance before shutdown, these routines must be executed at the next start-up.

If any parts of the instrument, carriers, or racks have become contaminated, the weekly maintenance procedure must be performed immediately.

9.3 Materials required

- Disposable latex gloves
- Protective glasses
- Lab coat
- Paper towels
- Lint-free cloths or Q-tips
- Ethanol (70 %), (Do not use to clean the cover!)
- De-ionized water
- MICROLAB[™] Detergent & Disinfectant (see chapter ORDER INFORMATION)
- MICROLAB[™] Disinfectant Spray (see chapter ORDER INFORMATION)



Use cleaning, disinfecting and decontaminating fluid in accordance with manufacturer's instructions. Do not use disinfecting materials which contain hypoclorite (Javel water, Chlorox) or bleaching fluids. Prepare disinfectant fluids according to manufacturers guidelines.

9.4 Maintenance Procedures

The operator will be guided by the CheckExtractor[™] software through the maintenance procedures. To execute a maintenance procedure, double-click the following icon on the desk-top:

In the Hamilton "Maintenance and Verification Run" window, the process status information view lists all maintenance and verification processes for the connected/selected instrument. It is from here that all the maintenance procedures can be initiated.

🖬 Hamilto	n Maintenance and Verification	Run				
File View	Process Tools Help					
😻 🖻						
	Instrument: Microlab® STAI	r IVD		Global mandatory flag:	mandatory	
	Serial Number: 2509					
Execute	Process	Status	Expired Date / Time	Processed Date / Time	Processed Status	Mandatory
	Daily maintenance	required	26.02.2010 17:10:00	25.02.2010 14:15:00	successful	mandatory
	Weekly maintenance	completed	05.03.2010 02:15:00	25.02.2010 14:15:00	successful	mandatory
	Verification	completed	20.07.2010 02:59:00	25.02.2010 14:15:00	successful	mandatory
·						
1						

Figure 25: "Maintenance and Verification Run" Window

Select the desired maintenance routine by clicking the specific check box and by pressing the "Run Process" button. The CheckExtractor™ software will issue on-screen instructions detailing all procedures required to perform the selected maintenance routine.



Always wear disposable gloves during maintenance.

Do not clean the instrument in the vicinity of open flames or devices which could create sparks. Do not use hot air blowers to dry the instrument. The liquids used for cleaning are flammable. When cleaning CO-RE gripper tool, be very careful of sharp edges on the CO-RE gripper's pincers. This manual provides indications regarding general disposal of waste. In addition, any regulations specific to the country of operation must be taken into account and observed.

Routine Completion

A maintenance routine is completed once the procedure has been fully executed and the results are within the specifications.

Aborting Maintenance Procedures

Aborting a maintenance procedure will lead to a 'failed' status, and maintenance will need to be started again.

9.4.1 Daily Maintenance

After clicking "**Yes**" the daily maintenance procedure will be started. Pressing "**No**" will abort the procedure.

The front cover (the hinged Plexiglas window that shields the instrument in front) can be opened for user intervention.

- 1. Once the maintenance procedure has been started, the pipetting arm moves to the left side. The operator now has access to the deck to check if cleaning is needed or not.
- If the deck is clean, continue with the daily maintenance.
- If the deck needs to be cleaned, the daily maintenance can be interrupted. Instead of the daily maintenance carry out the weekly maintenance.
- 2. Continuing the daily maintenance procedure will lead the user to the next maintenance task. The tip waste needs to be emptied. Discard tip waste with the rest of the laboratory's contaminated waste.



The tip waste has to be considered as contaminated.

3. For the next steps, the maintenance needles are required.



Figure 26: Maintenance needles

- 4. The procedure continues with the tightness check of the pipetting channels. The pipetting arm will travel to the right side to pick up the maintenance needles. Two checks are done with the pipetting channels, the over-pressure and the under-pressure check.
- 5. For the capacitive liquid level (cLLD) check, the needles are picked up again. One channel after the other is checked for the proper functioning of the cLLD.
- 6. The end of the daily maintenance is displayed.



Figure 27: "Daily Maintenance - successful" Window

The daily maintenance process status is saved on the instrument and a report file is created (see the section 12.5).



If any parts of the instrument, carriers or racks have become contaminated, the weekly maintenance procedure must be performed.

9.4.2 Weekly Maintenance

- The following tasks are carried out by the weekly maintenance:
- Cleaning the deck and carriers
- Checking the condition of the carriers
- Empting and cleaning of the tip waste
- Checking the tightness of the pipetting channel
- Verifying the functioning of the cLLD
- Cleaning of the pipetting head: stop disk, O-ring, tip eject sleeve
- Cleaning of the covers, Autoload protecting ribbon

After initialisation of the instrument, the operator will be asked to execute the weekly maintenance:



Figure 28: "Weekly Maintenance" Window

When the instrument is initialised, the weekly maintenance program advises the user to unload the deck manually. This step is carried out automatically by the Autoload.

1. Clean all carriers and with the MICROLAB[™] Disinfectant spray or a cloth, that is soaked with disinfectant, and leave them to dry. If they are heavily soiled, soak these carriers afterwards in a solution of MICROLAB[™] Detergent and Disinfectant (see the product data sheet for further information).

Examine each carrier for scratches on the barcode and any signs of damage. If damage is apparent, replace with new carriers.

2. Once the weekly maintenance program is resumed, the Autoload will be repositioned to the right hand side of the instrument.

Open the front cover and wipe the deck with a cloth saturated with MICROLAB[™] Disinfectant Spray. The slide blocks must be checked and cleaned if required. Close the front cover.



Do not directly spray at the Autoload unit or any electrical boards or connectors. Use a soaked cloth in the surroundings of electrical parts.

3. The next step of the maintenance procedure will reposition the Autoload to the left hand side of the instrument. At this point, the tip waste needs to be emptied and cleaned. Dispose of tip waste with the rest of the laboratory's contaminated waste.



The tip waste, the tip eject plate and the plastic bag should always be treated as if they are contaminated.

4. Remove the tip eject plate from the tip waste station and clean it: Spray MICROLAB[™] Disinfectant Spray directly onto the surface, and wipe. Remove the frame that holds the plastic bag in place, and discard the plastic in the laboratory's contaminated waste. Pull a new plastic bag over the frame and re-attach it. Put the clean tip eject plate back in place.



Do not spray the maintenance needles.

5. To prevent unreliable barcode reading, check the laser scanner window of the barcode reader and clean it with a lint-free cloth or Q-tips lightly soaked in Ethanol (70 %).



The laser scanner window must be completely dry and free from dust and fibres before the instrument can be reused.

6. For the next steps the maintenance needles are required.



Figure 29: Maintenance needles

- 7. The procedure continues with the tightness check of the pipetting channels. The pipetting arm will travel to the right hand side to pick up the maintenance needles. Two checks are carried out with the pipetting channels, the over-pressure and the under-pressure check.
- 8. For the capacitive liquid level (cLLD) check the needles are picked-up again.
- 9. The end of the weekly maintenance program is displayed:

Weekly maintenance - successful	
Weekly maintenance execute successfully completed	
ОК	

Figure 30: "Weekly Maintenance - successful" Window

The weekly maintenance process status is saved on the instrument and a report file is created (see the section 12.5).

Following the weekly maintenance, clean the rest of the instrument as follows:

10. Clean the tip eject sleeve (outer part of the pipetting channels) with a lint-free cloth soaked in MICROLAB™ Disinfectant Spray.



Figure 31: Cleaning of eject sleeve



Take care that no liquid gets inside the tip channel. Whenever it is necessary to move Channels on the x-Arm, move them gently by pushing close to their y-slide. Never force them as this may lead to damage. If possible, switch on the instrument as this will result in a smoother motion when Channels have to be moved on the x-Arm.

11. Clean the stop disk and the O-rings of the pipetting head (outer part of the pipetting channels) with a lint-free cloth soaked in MICROLAB[™] Disinfectant Spray.



Figure 32: Cleaning of O-rings



Avoid any liquid ingress into the tip channel.

12. Spray the front and side cover with MICROLAB[™] Disinfectant Spray and wipe dry.



Do not use highly concentrated alcoholic solutions, because they turn the cover panes dull.

13. Clean the Autoload protecting ribbon with a cloth soaked in MICROLAB[™] Disinfectant Spray, and wipe without exerting pressure.



Do not directly spray the Autoload unit or any electrical boards or connectors. Use a soaked cloth in the surroundings of electrical parts.

14. Clean the x-guide shaft behind the upper front cover with a dry cloth at least once a month.

9.5 Printing a Maintenance process status Report

The maintenance process report is printed as follows: :

- 1. Go to "File -> Open Report". All maintenance and verification processes which are found in the default "Report Path" are listed.
- 2. If necessary, change the report path using the browse button "..."
- 3. Select a report and press the "Open" button. The Report Viewer displays the selected report file
- 4. Go to "File -> Print" to print the report file

9.6 If Maintenance Fails

If an error is encountered during a maintenance procedure, try to resolve the problem and restart the maintenance procedure. If you cannot resolve the error yourself, call your local service representative.

9.7 Further Maintenance and Service

Because the CheckExtractor[™] is an IVD accessory, it is compulsory to keep the instrument always in optimal condition. This guarantees the highest reliability of the results of the CheckExtractor[™]. It is therefore mandatory to regularly service the instrument. A built in timer locks the instrument if service has not been performed on the instrument in the previous 200 days.



Make sure to service your CheckExtractor[™] on time. The CheckExtractor[™] will be unusable if it has not been serviced within the previous 200 days.

Please be sure to contact your local service representative on time for the next service.

The service includes:

- cleaning and lubricating the guides
- · checking parts that are subject to wear and tear
- checking of the pipetting and liquid system
- · checking of the functioning of the instrument
- · checking of the positioning of all moving parts
- adjustment, if required
- Installation of software and firmware updates within a version (if necessary)

9.8 Finishing the Maintenance and Cleaning



10. TECHNICAL SUPPORT

If you have any questions or experience any difficulties concerning one of the products, please do not hesitate to contact your local Greiner Bio-One distributor or directly the technical support department: **support.dx@gbo.com**.

Greiner Bio-One employs a technical service department staffed with experienced scientists with extensive practical and theoretical expertise in molecular biology and on the different oCheck[®] products.

11. TROUBLESHOOTING

11.1 Error Messages

Error message	Cause	Solution
Cannot start method xxx Daily/Weekly maintenance is required.	The maintenance process has not been carried out recently.	Start the maintenances method and follow its instructions (see section 10).
Please select at least one sample.	No input file was loaded or no sample is selected.	Make sure that a input file is loaded and at least one sample is checked.
You selected more than 96 samples. Please deselect some samples in order to start the run.	More than 96 samples are selected.	Make sure the number of selected samples is less than or equal to 96. Click on the green tick to deselect the sample.
Please specify a user id!	The white input box besides "User ID" is empty.	Please fill in your user ID.
Please specify an extraction kit barcode!	The white input box besides "Bar- code/Ref" is empty.	Position the cursor in this box and scan the kits barcode or enter it manually.
Expiry day of extraction kit must be in the future!	The kit is expired.	Use a new kit.
Please specify an extraction kit expiry date!	The white input box besides "Expiry Date" is empty.	Enter the expiry date of the kit in this box or move the cursor in the "Bar- code/Ref" box an scan the barcode.
Please specify an extraction kit lot number.	The white box besides the tack "Lot Number" is empty.	Enter the lot number of the kit in this box or move the cursor in the "Bar- code/Ref" box an scan the barcode.
Cover not closed.	The front cover is still open.	Close the front cover.
The liquid waste cap switch indicates that it is not fastened on the bottle or the cable is not plugged in. Please check the bottle and the cable.	The sensor at the liquid waste bottle does not give a signal.	Check that the bottle top is closed and the cable is connected on both ends.
The liquid waste bottle is full and has to be emptied before the run can start.	The liquid waste does not have enough volume left to take the liquid waste for the next run.	Please empty the liquid waste bottle. Be aware that the liquid waste is po- tentially infectious material (see 5.2)!
Barcode Error	A barcode is expected at the given position of the carrier, but cannot be read.	Make sure the barcode is undamaged and clean and faces to the right- hand side. If necessary re-enter the barcode manually.
Barcode Error	A barcode is expected at the given position of the carrier, but cannot be read.	Be sure the barcode fulfils the requi- rements given in 14.2, especially the 3 mm of the quiet zone. If necessary re-enter the barcode manually.
Barcode Error	The CheckExtractor™ finds a reagent barcode at a wrong position.	Make sure the reagents are at the indicated positions in the correct carrier. If necessary re-enter the barcode manually.
Tubes with the following barcodes are missing but have been defined in the worklist.	Some samples defined in the input file are not present.	Either load all the samples which are defined in the input file or change the input file to reflect the samples loaded.
Tubes with the following barcodes have been loaded but have not been defined in the worklist.	Samples have been loaded into the instrument, which are not in the worklist.	Please remove these samples or change the input file.

Error message	Cause	Solution
Not enough 1000 ul tips for less than 8 channels were loaded:	There are not enough tips in region C.	Load a sufficient number of tip racks into the instrument.
Carrier 'xxx' not found on the loading tray.	The carrier which is to be loaded is either not on the loading desk or is not inserted correctly.	Please slide the carrier carefully on the loading desk until your feel a mechanical stop.
Double blind entry failed.	When entering a barcode manually the two entries are different.	Re-enter the barcode in both fields.
The extraction plate is missing. A reload of the carrier is required.	While loading the plate carrier no extraction plate was detected.	Unload the carrier, insert the extrac- tion plate and reload the carrier again.
Wrong carrier type detected.	The loaded carrier is not the carrier expected.	Make sure to load each carrier at its specific track number.
Not enough buffer xx was loaded in trough xx.	Trough xx does not contain enough buffer for the run.	Refill the trough.
Too much Buffer M1 was loaded in trough 1. A reload of the carrier is required.	Since M1 is mixed with the protei- nase K it is compulsory to provide a precise amount of M1 to the system.	Reload the trough 1 with the exact amount of M1.
Error while creating PCR input file.	PCR input file is corrupted.	Call your Greiner Bio-One represen- tative.
Unable to update general file.	General file is corrupted.	Call your Greiner Bio-One represen- tative.
The elution plate on position 2 does not have the expected barcode.	A PCR setup method is set-up. But an elution plate with a barcode that does not match the input file was loaded.	Either place the correct elution plate into the instrument or change the worklist to the one that corresponds to the current elution plate.
An error occurred while running vector. The error description is: File 'xxx' must not be used in routine work! (File not validated) (xxx)	The method loaded is not validated and has been corrupted in some way.	Call your Greiner Bio-One representa- tive to reinstall a validated method.

11.2 Error Codes

All pipetting steps are monitored by the CheckExtractor[™]. You will find the results of the pipetting step of the Extraction method in the "MethodOutput" directory of the corresponding run. The file name ends with "_GBO_EXT_PipettingTracking.xls".

D	Natei Start Einfüg	en Seitenlayout Formeln Dat	en Überprüfen Ansi								
Ein	Ausschneiden	MS Sans Serif 10 A A		Zeilenumbruch	Standard -	Bedingte Als Tabelle	Standard Gut	Neutral Sesuchter Hy	Schlecht	Einfügen Löschen Format	Σ AutoSumme
	🗸 📝 Format übertrag	jen 🔭 🏧 🔍 👾 📥 🖊		verbinden und zentheren -	-3 . 10 000 ,00 -,0	Formatierung * formatieren	- Haspase	contracting according to the later	7		📿 Löschen 🐐
	Zwischenablage	G Schriftart	a Ausri	chtung i	S Zahl 5	i	Format	vorlagen		Zellen	
	A1 -	fx StepDescription									
1	StepDescription	SourceLabwareID	SourceLabwarePosition	MainErrorCodeAsp Reco	overyAsp VolumeCheck	TargetLabwareID	TargetLabwarePosition	MainErrorCodeDisp Recovery	Disp VolumeDisp	ChannelNumber DateAnd	Time
2	Add PK	MaN_VB_10x45_skirted_2200_0001	1	0	0 Not checked	ReagentTrough1	1	0	0 759	1 2015-05-2	0 11:49:13
3	Add PK	MaN_VB_10x45_skirted_2200_0001	0	0	0 Not checked	ReagentTrough1	5	0	0 759	1 2015-05-2	0 11:49:14
4	Add Buffer M1	ReagentTrough1	0	0	0 Not checked	ProcessingPlateHhs	A1	0	0 145	1 2015-05-2	0 11:50:36
5	Add Buffer M1	ReagentTrough1	2	0	0 Not checked	ProcessingPlateHhs	B1	0	0 145	2 2015-05-2	0 11:50:36
6	Add Buffer M1	ReagentTrough1	3	0	0 Not checked	ProcessingPlateHhs	C1	0	0 145	3 2015-05-2	0 11:50:36
7	Add Buffer M1	ReagentTrough1	4	0	0 Not checked	ProcessingPlateHhs	D1	0	0 145	4 2015-05-2	0 11:50:36
8	Add Buffer M1	ReagentTrough1	5	0	0 Not checked	ProcessingPlateHhs	E1	0	0 145	5 2015-05-2	0 11:50:36
9	Add Buffer M1	ReagentTrough1	6	0	0 Not checked	ProcessingPlateHhs	F1	0	0 145	6 2015-05-2	0 11:50:36
10	Add Buffer M1	ReagentTrough1	?	0	0 Not checked	ProcessingPlateHhs	G1	0	0 145	7 2015-05-2	0 11:50:36
11	Add Buffer M1	ReagentTrough1	8	0	0 Not checked	ProcessingPlateHhs	H1	0	0 145	8 2015-05-2	0 11:50:36
12	Add Buffer M1	ReagentTrough1	1	0	0 Not checked	ProcessingPlateHhs	A2	0	0 145	1 2015-05-2	0 11:50:37
13	Add Buffer M1	ReagentTrough1	2	0	0 Not checked	ProcessingPlateHhs	B2	0	0 145	2 2015-05-2	0 11:50:37
14	Add Buffer M1	ReagentTrough1	3	0	0 Not checked	ProcessingPlateHhs	C2	0	0 145	3 2015-05-2	0 11:50:37
15	Add Buffer M1	ReagentTrough1	4	0	0 Not checked	ProcessingPlateHhs	D2	0	0 145	4 2015-05-2	0 11:50:37
16	Add Buffer M1	ReagentTrough1	5	0	0 Not checked	ProcessingPlateHhs	E2	0	0 145	5 2015-05-2	0 11:50:37
17	Add Buffer M1	ReagentTrough1	6	0	0 Not checked	ProcessingPlateHhs	F2	0	0 145	6 2015-05-2	0 11:50:37
18	Add Buffer M1	ReagentTrough1	7	0	0 Not checked	ProcessingPlateHhs	G2	0	0 145	7 2015-05-2	0 11:50:37

Figure 33: Pipetting Tracking file

If an error occurs, open this file. In the 1st column the pipetting step is shown. In the 4th and 9th column error codes are given for aspiration and dispensing respectively. The following table lists all error codes that might appear.

For the PCR setup method similar data is saved in the file that ends with "_GBO_PCR_Report.xls".

Date	Start Ein	fügen Seitenla	iyout Formein	Daten Über	orüfen A	Ansicht												
Ê	🔏 Ausschneide	MS Sans S	Serif * 10 *	A A = =	_ >>-	Zeilenumbruc	h S	tandard	×	1 ≦5		Standard	Gut	Neutral	Schlecht		P 🏢	Σ AutoSumme
Einfüg	en 🍠 Format über	tragen F K	<u>u</u> • 🖽 • 🌺	• <u>A</u> • 🔳 🗉	= # #	Verbinden un	d zentrieren 👻 🔮	🤧 - % 000	•,0 ,00 ,00 •,0	Bedingte Formatierung *	Als Tabelle formatieren	Ausgabe	Berechnung	Besuchter Hy	Eingabe	Ţ Einfügen L	öschen Forma	t 🖉 Löschen -
	Zwischenablage	G	Schriftart	G	A	usrichtung	G.	Zahl	G				Formatvorlagen			2	lellen	
	P5	• () f _x																
																		1 0
																	Reco	wery Recovery
1 S	ourceSeqName	SourceLiquidID) SourceBarcode	SourceLabware	ID Source	eLabwarePosition	TargetBarcode	 TargetLat 	bwareID	TargetLabwar	ePosition \	//orklistVolumes	PipettedVolumes	PipettingStatus	MainErrorCodeAsp	MainErrorCode	eDisp Asp	Disp
2 s	eqElutionPlate	Sample001	E1406013001207	7 Gre_96_VB_PC	R_2A1		E140601300118	33 Gre_96_V	B_PCR_	A1		5	5	0k	(0	0 0
3 s	eqElutionPlate	Sample002	E1406013001207	7 Gre_96_VB_PC	R_2B1		E140601300118	33 Gre_96_V	B_PCR_	B1		5	5	0k	(0	0 0
4 s	eqElutionPlate	Sample003	E1406013001207	7 Gre_96_VB_PC	R_2C1		E140601300118	33 Gre_96_V	B_PCR_	C1		5	5	Ok	(0	0 0
5 s	eqElutionPlate	Sample004	E1406013001207	7 Gre_96_VB_PC	R_2D1		E140601300118	33 Gre_96_V	B_PCR_	D1		5	5	Ok	(0	0 0
6 s	eqElutionPlate	Sample005	E1406013001207	7 Gre_96_VB_PC	R_2E1		E140601300118	33 Gre_96_V	B_PCR_	E1		5	5	Ok	(0	0 0
7 s	eqElutionPlate	Sample006	E1406013001207	7 Gre_96_VB_PC	R_2F1		E140601300118	33 Gre_96_V	B_PCR_	F1		5	5	Ok	(0	0 0
8 s	eqElutionPlate	Sample007	E1406013001207	7 Gre_96_VB_PC	R_2G1		E140601300118	33 Gre_96_V	B_PCR_	G1		5	5	Ok	(0	0 0
9 s	eqElutionPlate	Sample008	E1406013001207	7 Gre_96_VB_PC	R_2H1		E140601300118	33 Gre_96_V	B_PCR_	H1		5	5	Ok	(0	0 0
10 s	eqElutionPlate	Sample009	E1406013001207	7 Gre_96_VB_PC	R_2A2		E140601300118	33 Gre_96_V	B_PCR_	A2		5	5	Ok	(0	0 0
and the second s		10 1 010																

Figure 34: PCR Report file

Here the error codes are given in the 12th and 13th column for aspiration and dispensing respectively.

Error Code	Туре	Description
0	No error	-
1	Syntax Error	There is a wrong set of parameters or parameter ranges.
2	Hardware Error	Steps lost on one or more hardware components, or component not initialized or not functioning.
3	Not Executed Error	There was an error in previous part command.
4	Clot Error	Clot detected.
5	Barcode Error	Barcode could not be read or is missing.
6	Insufficient Liquid Error	Not enough liquid available.
7	Tip Present Error	A tip has already been picked up.
8	No Tip Error	Tip is missing or not picked up.
9	No Carrier Error	No carrier present for loading.
10	Execution Error	A step or a part of a step could not be processed.
11	Pressure LLD Error	A dispense with pressure liquid level detection is not allowed.
12	Calibrate Error	No capacitive signal detected during carrier calibration proce- dure.
13	Unload Error	Not possible to unload the carrier due to occupied loading tray position.
14	Pressure LLD Error	Pressure liquid level detection in a consecutive aspiration is not allowed.
15	Parameter Error	Dispense in jet mode with pressure liquid level detection is not allowed.
16	Cover Open Error	Cover not closed or cannot be locked.
17	Improper Aspiration Error	The pressure-based aspiration/dispensation control reported an error (not enough liquid).
18	Wash Liquid Error	Waste full or no more wash liquid available.
19	Temperature Error	Incubator temperature out of range.
20	TADM overshot	Overshot of limits during aspirate or dispense. Note: On aspirate this error is returned as main error 17. On dispense this error is returned as main error 4.
21	Labware Error	Labware not available.
22	Labware Gripped Error	Labware already gripped.

Error Code	Туре	Description
23	Labware Lost Error	Labware lost during transport.
24	Illegal target plate position	Cannot place plate, plate was gripped in a wrong direction.
25	Illegal Intervention Error	Cover was opened or a carrier was removed manually.
26	TADM undershot	Undershot of limits during aspirate or dispense. Note: On dispense this error is returned as main error 17. On aspirate this error is returned as main error 4.
27	Position Error	The position is out of range.
28	Unexpected cLLD Error	The cLLD detected a liquid level above start height of liquid level search.
29	Area already occupied	Instrument region already reserved.
30	Imposible to occupy Area	An region on the instrument cannot be reserved
31	Anti drop control error	Anti-drop controlling out of tolerance
99	Slave Error	Slave error
100	Wrong Carrier Error	Wrong carrier barcode detected.
101	No Carrier Barcode Error	Carrier barcode could not be read or is missing.
102	Liquid Level Error	Liquid surface not detected.
103	Not Detected Error	Carrier not detected at deck end position.
104	Not Aspirated Error	Dispense volume exceeds the aspirated volume.
105	Improper Dispensation Error	The dispensed volume is out of tolerance.
106	No Labware Error	The labware to be loaded was not detected by Autoload module.
107	Unexpected Labware Error	The labware contains unexpected barcode (may only occur on a Reload Carrier step).
108	Wrong Labware Error	The labware to be reloaded contains wrong barcode.
109	Barcode Mask Error	The barcode read doesn't match the barcode mask defined.
110	Barcode Not Unique Error	The barcode read is not unique. Previously loaded labwa- re with same barcode was loaded without unique barcode check.
111	Barcode Already Used Error	The barcode read is already loaded as unique barcode (it's not possible to load the same barcode twice).
112	Kit Lot Expired Error	Kit Lot expired.
113	Delimiter Error	Barcode contains character which is used as delimiter in result string.

11.3 Causes of Trouble

Trouble	Cause	Solution
Barcodes are not read.	Direct sunlight interferes with the cor- rect action of the barcode-reader	Place the CheckExtractor™ in an area without direct sunlight. If necessary re-enter the barcode manually.
Barcodes are not read.	The barcodes do not fulfil the specifi- cations. Especially the quiet zone and the length of the barcode have to be observed.	Make sure the barcodes fulfil the specifications (see 12.2). If necessary re-enter the barcode manually.
Barcodes are not read.	Barcode are misaligned or misplaced.	Place the barcodes vertically on troughs and tube as shown in the corresponding figures. If necessary re- enter the barcode manually.
Carrier is not loa- ded.	The Autoload cannot catch the carrier.	Make sure the carrier is pushed in all the way to the first mechanical block.
CheckExtractor™ refuses to work at all.	The period of 200 days between two services has been exceeded.	A service is necessary.
Input file is not found.	The CheckExtractor [™] allows for the extraction method only input file which file names end with "EXT"	Rename you input file accordingly.
PCR setup errone- ous	Elution plate was not correctly pre- pared.	Be sure that the elution plate is fully thawed, spun down and open, before loading it into the CheckExtractor [™] .
Run stopped due to missing tips	The rip racks were not filled as shown in the tip editor.	Start a run only with enough tips. The number of required tips is denoted by the software.
Run stopped un- expectedly and all samples are marked with an error.	The front cover was opened during the run.	Do not open the front cover before you are asked to clean up the desk.
Sample ID not cor- rect in the worklist.	While registering the sample IDs in the Excel® spreadsheet leading zeroes are vanish.	Make sure, all columns in the Excel® file are formatted as "Text".
Sample was not extracted.	Not enough sample volume.	900 μl of sample are necessary to process the sample.
Sample was not extracted.	The liquid level detection found a wrong volume or the tip touched the bottom of the sample tube.	Be sure that the sample tube is inserted firmly and all the way down into the drawer
Uncontrolled or unpredictable be- haviour	Strong electromagnetic RF fields in- terfere with the CheckExtractor™. RF sources may be mobile phones.	Do not place the CheckExtractor™ close to equipment that emits electromagnetic RF fields.

11.4 Assistance

For further information please contact your local Greiner Bio-One representative.

12. TECHNICAL SPECIFICATIONS

12.1 Physical Properties

12.1.1 CheckExtractor™

Dimension	Value
Width	1124 mm
Height	903 mm
Depth	795 mm (without loading desk) 1010 mm (with loading desk)
Weight	approx. 140 kg

12.1.2 Distances

A minimum amount of space should be left on all sides of the CheckExtractor[™] to allow for servicing, operations and ventilation of the system. This space should be no less than:

Side	Distance
Left side	300 mm
Right side	300 mm
Above	500 mm
Behind	100 mm
In front	1000 mm
Below	600 mm

12.1.3 Enviromental Conditions

Parameter	Range
Tomporatura	Operating: 20–25 °C
Temperature	Storage: -25–70 °C
Luxidity.	Operating: 20-80 %, non-condensing
numidity	Storage: 10–90 %, non-condensing, indoor
Sunlight	No direct sunlight or other intensive light sources Indoor storage only
Altitude	Up to 2000 m above mean sea level
Autude	Storage: as required for air travel
Dust	Preferred dust-free environment

12.1.4 Operating Data

Parameter	Range
Maximum power consumption	≤ 600 VA
Voltage	115 V~/230 V~
Frequency	50/60 Hz ± 5 %
Delayed action fuse	115 V~: 6.3 A (T6.3AL250) 230 V~: 3.15 A (T3.15A250)
Installation	Ш
Pollution Degree	2
Temperature range	15–30 °C
Relative humidity	15–85 % (no condensation)
Noise level	< 65 dBA (regarding EN27779)
Altitude	Up to 2000 m above sea level
Barcode Reader	Laser class 2
UV Lamp	Wave length: 254 nm Power: 3 x 30 W
Indoor use only	

The CheckExtractor[™] fulfils the requirements of directive 2004/108/EC (EMC).

12.1.5 Pipetting Specifications

Volume [µl]	Tip Size 300µl	Tip Size 1000µl	Accuracy [%]	Precision [%]
10	Х		± 10	5.0
50	Х	Х	± 5	2.5
100	Х	Х	± 5	2.0
200	Х	Х	± 5	1.5
1000		Х	± 5	1.5

12.1.6 Connectivity

The CheckExtractor[™] itself and the controlling PC are equipped with standard USB 2.0 sockets for communication between the two devices. The PC has additional USB sockets, which can be used to connect external data storage devices (like USB sticks) for the back-up of data or the importing of input files.

The CheckExtractor[™] can be integrated into a LIMS via the ethernet (RJ45) connector on the PC.

12.1.7 Instrument Life Span

The instruments have been designed for a lifespan of 6 years under the assumption that the user maintenance and preventative maintenance is performed at the appropriate intervals.



The maximum lifespan of the CheckExtractor[™] is limited to 6 years! Working with an older instrument may lead to false results. Return the instrument to your supplier after 6 years!

12.2 Barcodes

12.2.1 Barcode Specifications

Length of string	Maximum 20 characters excluding start, stop and check characters, depending on the code length (see label dimensions).			
Code Density, Tolerance				
Check character	ISBT standard	One character		
	Code 128	One character		
	Code 39	None		
	Codabar	None		
	Code 2 of 5 Interleaved	None		
	UPC A	One character		
Quiet Zone	10 times the x dimension, but at least 3 mm.			
Print quality	The barcode print must be of a high quality. A printed barcode with an ANSI/CEN/ISO grade A or B is required. Offset, typographic, intaglio and flexographic printing are suitable. Mechanical dot matrix and thermo matrix printing are not suitable. The surface may be treated, sealed or plastic-coated.			

12.2.2 Sample Barcodes



Figure 35: Barcode dimensions

Dimension	Min.	Max.
A Label length	-	45 mm
B Code length	-	39 mm
C Quiet zone	3 mm	
D Label width	12 mm	-
E Code width	12 mm	-
F Distance from code to label edge	-	1 mm



Figure 36: Tube barcode

The label must be attached within a range of between 20 mm to 65 mm from the bottom of the tube. The label must fit tightly at an angle of approximately 90° to the tube's axis. The label must fit tightly over its whole length.