


OrisTM *Cell Invasion & Detection Assay Kit*
Product No.: CIA101DE & CIA200DE

96-well, 3-D Assay for Investigating
Cell Invasion of Adherent Cell Lines

Protocol & Instructions

Patent Pending



PLATYPUS TECHNOLOGIES, LLC
5520 NOBEL DRIVE, SUITE 100
MADISON, WI 53711
TOLL FREE: (866) 296-4455
PHONE: (608) 237-1270
FAX: (608) 237-1271

WWW.PLATYPUSTECH.COM

PLATYPUS
TECHNOLOGIES

*Bringing Science to the Surface*TM

RM0033.01


Oris™ CELL INVASION & DETECTION ASSAY

I. INTRODUCTION

The Oris™ Cell Invasion & Detection Assay is a reproducible, sensitive, and flexible assay that can be used to monitor cell invasion. Formatted for a 96-well plate, the assay utilizes Oris™ Cell Seeding Stoppers (made from a medical-grade silicone) to restrict cell seeding to the outer annular regions of the wells. Removal of the stoppers reveals a 2 mm diameter unseeded region in the center of each well, i.e., the invasion zone, into which the seeded cells may then invade. The Oris™ Detection Mask is applied to the plate bottom and restricts visualization to the invasion zone, thus allowing only invading cells to be detected (see Figure 1). The Oris™ Cell Invasion & Detection Assay is designed to be used with any commercially available stain or labeling technique, but comes complete with a Calcein AM Staining Reagent. Calcein AM is a fluorescent dye that passes through the membrane of live cells and is useful for short-term labeling of cells at the end of the invasion period. Calcein AM is processed by esterases inside the cell allowing the dye to bind calcium and be trapped inside the cell, and has an excitation/emission wavelength of 495/515 nm. Assay readout can be performed by microscopic examination or by using a plate reader. The Oris™ Cell Invasion & Detection Assay kit has been uniquely designed to detect cellular invasion *in vitro* within a 3-dimensional extracellular matrix comprised of a basement membrane extract (BME) purified from Murine Engelbreth-Holm-Swarm (EHS) tumor. The Oris™ Cell Invasion & Detection Assay system has been designed for use with adherent cell cultures. Performance of the Oris™ Cell Invasion & Detection Assay was optimized using the invasive HT-1080 fibrosarcoma and the non-invasive 3T3-Swiss albino fibroblast cell lines.

Using the Oris™ Cell Invasion & Detection Assay offers the following benefits:

- **Membrane-free Cell Invasion**- no cumbersome cell culture inserts or Transwell® membrane devices to limit cellular movement; there is no membrane that restricts the ability to image cells.
- **Preserve Cell Morphology**- the Oris™ Cell Invasion & Detection Assay provides a more native 3-D environment since cells are ensconced in an extracellular matrix.
- **Real-time Monitoring**- Invasion-induced changes in cell structure can be monitored in real-time with a microscope or digital imaging system.
- **Versatile**- analyze cells using multiple fluorescent probes, labels or colorimetric stains in a single well.
- **Flexible**- design kinetic or endpoint assays using a fluorescence plate reader.



II. ORIS™ PLATE DIMENSIONS (per well)

Diameter of Well	6.5 mm
Diameter of Stopper Space (Detection Zone)	2 mm
Suggested Media Volume per Well (populated with Stoppers)	100 µl
Effective Area of Outer Annular Region (seeding region) per Well	30.03 mm ²
Effective Area of Central Detection Zone per Well	3.14 mm ²

Important: Read Instructions Before Performing any Oris™ Assay.

III. MATERIALS PROVIDED

Product No.: CIA101DE

Component	Quantity	Storage
Oris™ 96-well Plate	1	Room Temperature
Oris™ Cell Seeding Stoppers	96	Room Temperature
Oris™ Detection Mask	1	Room Temperature
Oris™ Stopper Tool	1	Room Temperature
Oris™ Basement Membrane Extract (BME)	5 mL	≤ -20°C*
Calcein AM Reagent	20 µL	-20°C with dessicant**

Product No.: CIA200DE

Component	Quantity	Storage
Oris™ 96-well Plates	2	Room Temperature
Oris™ Cell Seeding Stoppers	2 x 96	Room Temperature
Oris™ Detection Mask	2	Room Temperature
Oris™ Stopper Tool	2	Room Temperature
Oris™ Basement Membrane Extract (BME)	2 x 5 mL	≤ -20°C*
Calcein AM Reagent	2 x 20 µL	-20°C with dessicant**

* Oris™ BME reagent can be stored at -20°C in a manual defrost freezer if kit will be used within 3 months of receipt. For long-term storage, store Oris™ BME reagent at -80°C.


** Calcein AM can be stored at -20°C in a manual defrost freezer with a dessicant for use within 6 months of receipt.

IV. MATERIALS REQUIRED

- Biological Cells
- Cell Culture Growth Medium with Fetal Bovine Serum
- Sterile PBS
- Hanks Balanced Salt Solution (HBSS)
- Serum-Free Cell Culture Medium
- Sterile Pipette Tips and Pipette or Multi-Channel Pipette
- Trypsin or Non-Enzymatic Cell Removal Reagent or Scraper
- Inverted Microscope (optional)
- Fluorescence Microplate Reader (optional)
- Cell Labeling Fluorescent Agent (e.g., CellTracker™ Green) - *required if performing staining in addition to or in place of Calcein AM.*

V. PRECAUTIONS AND RECOMMENDATIONS

- For Research Use Only. Not for use in diagnostic procedures.
- Oris™ Cell Seeding Stoppers are for single use only; Platypus can not guarantee the integrity or performance of the stopper material after a second sterilization procedure.
- Recommendations for Oris™ BME:
 - Thaw on ice (2-8°C) overnight.
 - The inherent variability of invasiveness between different cell lines may make it necessary to dilute the Oris™ BME prior to application. Less invasive cell lines may require a more permissive barrier. Do not dilute Oris™ BME lower than 9 mg/mL as this will inhibit gel formation.
 - Aliquot and freeze any remaining Oris™ BME. Avoid repeated freeze-thaw cycles.
- It should be appreciated that there will be lot-to-lot differences in the performance of biologically derived materials, such as BME and FBS that may impact the extent of invasion for a given cell line, as well as in the behavior of different cell lines with respect to invasive capabilities. To minimize the effects of biological variation in cell invasion and improve reproducibility, we recommend, 1) conducting invasion assays with cells lines within a narrow passage number [10-12 passages], and 2) using an invasion inhibitor (i.e., Latrunculin) with an established dose response curve. Controlling these two parameters will help establish a reference point when comparing data from different experiments.



V. ORIS™ CELL INVASION & DETECTION ASSAY PROTOCOL

The following steps should be performed in a biological hood using aseptic technique to prevent contamination.

1. If desired, cells can be starved by incubating for 18 - 24 hours in serum-free medium prior to assay (0.5% fetal bovine serum may be used if needed).
2. Mix 500 μ L of the thawed Oris™ BME with 1.5 mL of cold Hanks Balanced Salt Solution (HBSS) or serum-free media to create a BME Coating Solution.



NOTE: Oris™ BME will gel in 5-10 min above 15°C; therefore, keep the BME Coating Solution and the remaining Oris™ BME on ice until ready to use (Oris™ BME will be used again in Step 17). The use of chilled pipette tips/reservoirs may be beneficial.

3. Pipette 100 μ L of the BME Coating Solution into the wells of Column 1. Gel loading pipette tips may be helpful. Avoid bubble formation in the BME Coating Solution by not fully expelling all contents of the pipette. Immediately remove the BME Coating Solution from the wells of Column 1 and return contents back into the BME Coating Solution reservoir. Repeat procedure for the remainder of plate.
4. After the BME Coating Solution has been applied to all wells of the plate, check each well for excess solution. If excess solution is apparent, remove via pipette.
5. Incubate the plate in a humidified chamber (37°C, 5% CO₂) for 15 - 30 minutes.
6. Under sterile conditions, populate the 96-well plate with Oris™ Cell Seeding Stoppers:
 - Vertically position the tip ends of two, 4-stopper strips into one full column of 8 wells at a time (Figure 2A).
 - Gently press down on the strip backbone to partially insert the stoppers halfway into the well (Figure 2B).
 - When both stopper strips have been partially inserted in 1 column, ensure that the position of the stoppers is vertical with respect to the well wall, making any necessary adjustments (Figure 2C).
 - Using the Oris™ Stopper Tool, firmly press down on the strip backbone to fully insert the stoppers into each well (Figure 2D and 2E). Repeat for all remaining columns.



NOTE: It is extremely important to ensure that the stoppers are inserted perpendicular to the well bottom and fully engaged with the bottom of the well. Failure to do so will increase the CV of your data set.

7. Visually inspect the underside of the populated 96-well plate to ensure that the bottom of the Oris™ Cell Seeding Stoppers are firmly sealed against the well bottom. To inspect the stoppers, turn the plate over and examine the stoppers for sealing (see Figure 3). If incomplete sealing is observed, return the plate to the upright position and use a sterile instrument to gently push the stopper back into the well until sealing is observed.
8. Apply the Oris™ Detection Mask to the bottom of the 96-well plate.

First Time Users: In order to prevent splashing of well contents, familiarize yourself with the attachment and removal of the Detection Mask before any liquids are placed in the wells.

- Orient the chamfered corners of the mask with those of the 96-well plate, ensuring that the A1 corner of the mask is aligned with the A1 well of the plate.
- Align the holes in the attachment lugs with the bosses on the bottom of the 96-well plate and gently press the mask onto the bottom of the plate.




Figure 2. Stopper Insertion Process. A) Placement of Stoppers into Wells, B) Close-up of Stoppers Partially Inserted into Wells, C) Proper Placement of Stoppers, D) Pressing of Stoppers into Wells, and E) Fully Inserted Stoppers




Figure 3. Stoppers that are A) Partially Sealed, B) Unsealed, & C) Completely Sealed



NOTE: It may be necessary to wash the mask with ethanol to remove dust and debris since the mask is **not** sterile. The mask may be applied at any point during the assay. For kinetic assays, it is often most convenient to apply the mask at the beginning of the assay before any liquids are placed in the well. For endpoint assays, using fixed and stained cells, it is often most convenient to apply the mask just before reading assay results.

- If performing a kinetic analysis of Cell Invasion, pre-stain the cells with a fluorescent stain now.
- Collect cells and prepare a suspension that is 10-fold greater in density than the optimal seeding concentration using cell culture growth medium containing 10% FBS.

First Time Users: The optimum seeding density of cells must be determined as an integral part of the design of the cell invasion assay. Please see Appendix I for a discussion of this process.

- Pipette 100 μ l of suspended cells into each test well through either of the side ports of the Cell Seeding Stopper.



NOTE: For best results, add or extract media by placing the pipette tip along the wall of the well. Care should be taken not to disturb the Cell Seeding Stopper or the BME layer when introducing the pipette tip into the well. A gel loading tip may be useful.

- IMPORTANT:** Lightly tap the plate on your work surface to evenly distribute well contents (extreme tapping may result in splashing of well contents and lead to contamination).

- Incubate the seeded plate containing the Oris™ Cell Seeding Stoppers in a humidified chamber (37°C, 5% CO₂) for at least 4 hours (cell line dependent) to permit cell attachment.

- Remove plate from incubator.

- Use the Oris™ Stopper Tool to remove all stoppers. See Figure 4.
 - Secure the 96-well plate by holding it firmly against the deck of your work space. Slide the tines of the stopper tool under the backbone of the stopper strip, keeping the underside of the tool flush with the top surface of the plate.
 - Lift the stopper tool **vertically** to gently remove the stopper.



NOTE: DO NOT use the stopper tool as a lever to pry the stoppers from the well (see Figure 4E), as doing so may cause displacement of seeded cells.

- Remove media with a pipette and gently wash wells with 100 μ l PBS (or media) to remove any unattached cells. Do not aspirate using an in-house vacuum.

- Add 40 μ l of the thawed Oris™ BME to each well (supplements, such as FBS or growth factors, may be mixed with BME).



NOTE: Oris™ BME gels in 5-10 min above 15°C, therefore, you must keep the thawed Oris™ BME on ice until ready to use. In addition, the use of chilled pipette tips/reservoirs might be beneficial. Since different cell lines and different treatments can result in a wide range of invasive behavior, the BME overlay may be optimized to fit each experiment by diluting this material to permit more invasion. However, do NOT dilute the Oris™ BME below a concentration of 9 mg/mL.

- Incubate plate in a humidified chamber (37°C, 5% CO₂) for 30 - 60 minutes to permit polymerization of Oris™ BME.
- For reference purposes, designate several 'reference' wells (e.g., $n = 4$) to which a staining agent will be added to serve as $t=0$ pre-invasion controls. After staining of the reference wells, capture well data via microscope, plate reader or digital imaging system.
- Add 100 μ l of serum-free cell culture medium on top of the BME overlay. **Optional:** Invasion inhibitors or stimulants may be added to the media.
- Incubate plate in a humidified chamber (37°C, 5% CO₂) to permit cell invasion (length of incubation is cell line dependent). The Oris™ BME will remain gelled for up to 14 days. Refresh media or supplements, every 48 - 72 hours, as needed, for the duration of the invasion experiment.
- For an endpoint analysis of cell invasion, please see Appendix I for the Calcein AM staining procedure.




Figure 4. Removal of Stoppers. Panels A, B, and C) Position the Tines of the Stopper Tool between the Stopper Tips, D) Lift Vertically, and E) Do NOT Pry Stoppers

VI. DATA ACQUISITION

The readout of the Oris™ Cell Invasion & Detection Assay can be conducted at any time, allowing the user to perform a kinetic assay or an endpoint assay. The Oris™ Cell Invasion & Detection Assay is designed to be used with any commercially available stain or labeling technique, but is supplied complete with the Calcein AM Reagent. The readout can be performed by microscopic examination or by plate reader.

Microscopic Analysis


- Cell counting or image capture / analysis (using software, such as Image J freeware, available from NIH)

Plate Reader Analysis

- Setup on individual plate readers varies according to make and model. Consult your user manual for proper operation.
- The plate reader MUST be set to use the bottom probe read.

Sample Data Obtained via Microscopic Examination and Plate Reader are shown in Figure 5.

- Wells were seeded with 50,000 HT-1080 cells (i.e., 100 μ l of 5×10^5 cells/mL) that had been serum starved for 18 hours and the plate was then incubated for 4 hours. The stoppers were removed from the wells and the Oris™ BME (with 10% FBS or without FBS) was overlaid on the cells. The plate was incubated in a humidified chamber for 48 hours to permit cell invasion. Cells were labeled with Calcein AM and images were captured using a Zeiss Axiocvert microscope (5X magnification). Fluorescence in the analytic zone was quantified by using a plate reader. Each column represents the mean \pm SD of at least 4 wells. A non-invasive cell line, 3T3-Swiss albino, served as the negative control. The images below, captured without a detection mask in place, illustrate representative data from pre-invasion (t=0 hrs) and post-invasion (t=48 hrs) wells. The graph depicts the average RFU's in the invasion zones for each condition, confirming the invasion augmenting effect of FBS on serum starved HT-1080 cells and the lack of invasion by non-invasive 3T3-Swiss albino cells.
- As observed in Figure 5, HT-1080 cells formed invadopodial structures projecting into the central analytic zone (panel 1b) as compared to a time zero control (panel 1a).
- Invasion by serum-starved HT-1080 cells was augmented by including 10% v/v fetal bovine serum (FBS) in the BME overlay (panel 1c).
- In contrast, as observed in images 2b and 2c, the non-invasive 3T3-Swiss albino cells did not form any invadopodial structures and showed a nominal amount of FBS-independent migration into the central detection zone as compared to a time zero control (image 2a).



Immunostaining Analysis

- As determined by immunostaining (shown in Figure 6), the structures projecting into the BME formed by the invasive HT-1080 cells exhibited classic hallmarks of invadopodia, namely colocalization of F-actin with cortactin (panel a), and the expression of proteases such as MMP-9 and cathepsin B (panels b and c). These invadopodia could also be observed in 3-dimensions at different focal points along the z-axis. This is in contrast to any migratory activity observed by non-invasive cells that occurs only within the 2-dimensional space of the x and y axes (data not shown).
- Cells from the Oris™ Cell Invasion & Detection Assay were fixed, permeabilized, and pre-treated by sequential incubations in 3.7% formaldehyde, PBS, 0.5% Triton X-100, PBS, Image-iT™ FX signal enhancer, and PBS. Immunostaining was performed by incubating with primary antibodies (2 hours, 37°C) and Alexa Fluor® 488 conjugated secondary antibodies at 1:200 dilution (1 hour, 37°C) followed by an Alexa Fluor® 555 phalloidin counterstain (1 hour). Images were collected using a Nikon TE300 inverted microscope equipped with a Photometrics Coolsnap fx CCD camera, deconvolved using Slidebook™ v4.2 (Intelligent Imaging Innovations) and processed using Adobe® Photoshop® CS2 (Adobe Systems).




Figure 6. Confirmation of Cell Invasion using Immunostains


References:

1. Weaver, AM; Clin Exp Metastasis 2006. "Invadopodia: specialized cell structures for cancer invasion." 23:97-105.
2. Furmaniak-Kazmierczak et al; Circulation Research 2007. "Formation of extracellular matrix-digesting invadopodia by primary aortic smooth muscle cells." 100:1328-1336.
3. Hulkower, KI et al; Eur. J. Biochem 2000. "Fluorescent microplate assay for cancer cell-associated cathepsin B." 267:4165-4170.

Oris™ is a trademark of Platypus Technologies, LLC.

CellTracker™ Green is a trademark of Invitrogen Corporation.

Transwell® is a registered trademark of Corning, Inc.



VII. ORDERING INFORMATION

Product No.	Product Description	Package Size
CIA101DE	Oris™ Cell Invasion & Detection Assay, 1-pack: Oris™ 96-well plate, 1 Oris™ Basement Membrane Extract, 5 mL Oris™ Cell Seeding Stoppers, 96 Oris™ Detection Mask, 1 Oris™ Stopper Tool, 1 Calcein AM Reagent, 20 µL	1-pack
CIA200DE	Oris™ Cell Invasion & Detection Assay, 2-pack: Oris™ 96-well plates, 2 Oris™ Basement Membrane Extract, 2 x 5 mL Oris™ Cell Seeding Stoppers, 2 x 96 Oris™ Detection Mask, 2 Oris™ Stopper Tool, 2 Calcein AM Reagent, 2 x 20 µL	2-pack

To place an order, visit the Platypus Technologies website at: www.platypustech.com/order_main.html.
For technical assistance, contact Technical Support at (866) 296-4455 or techsupport@platypustech.com.

VIII. TERMS & CONDITIONS

Certain uses of these products may be covered by U.S. Pat. No. 6,284,197, No. 7018838, No. 10/597,118, No. 11/342,413, and No. 60/836,109, licensed to PLATYPUS. Certain applications of PLATYPUS products may require licenses from other parties. Determining the existence and scope of such third party intellectual property is the responsibility of the PURCHASER. Purchase of the product provides the PURCHASER with a limited non-transferable license under any PLATYPUS patents or patent applications to use the product for internal research unless there is a written limitation to this license in the product literature. PURCHASER is responsible for carefully reviewing the product literature and respecting any limitations to this license, e.g. limitations for commercial use or research by for-profit institutions. These products may not be resold, modified for resale, used to manufacture commercial products, or used to develop commercial products without the express written approval of PLATYPUS.


These products are intended for research or laboratory use only and are not to be used for any other purposes, including, but not limited to, unauthorized commercial purposes, *in vitro* diagnostic purposes, *ex vivo* or *in vivo* therapeutic purposes, investigational use, in foods, drugs, devices or cosmetics of any kind, or for consumption by or use in connection with or administration or application to humans or animals.

PLATYPUS warrants that its products shall conform substantially to the description of such goods as provided in product catalogues and literature accompanying the goods until their respective expiration dates or, if no expiration date is provided, for 6 months from the date of receipt of such goods. PLATYPUS will replace, free of charge, any product that does not conform to the specifications. This warranty limits PLATYPUS's liability only to the replacement of the nonconforming product.

THIS WARRANTY IS EXCLUSIVE AND PLATYPUS MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The stated express warranties, and the remedy provided for breach thereof, are in lieu of all other liability or obligations of PLATYPUS for any damages whatsoever arising out of or in connection with the delivery, use, misuse, performance, or the inability to use any of its products. IN NO EVENT SHALL PLATYPUS BE LIABLE UNDER ANY LEGAL THEORY (INCLUDING BUT NOT LIMITED TO CONTRACT, NEGLIGENCE, STRICT LIABILITY IN TORT, OR WARRANTY OF ANY KIND) FOR ANY INDIRECT, SPECIAL, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES (INCLUDING BUT NOT LIMITED TO LOST PROFITS) EVEN IF PLATYPUS HAD NOTICE OF THE POSSIBILITY OF SUCH DAMAGES. Without limiting the effect of the preceding sentence, PLATYPUS's maximum liability, if any, shall not exceed the purchase price paid by PURCHASER for the product.

This warranty shall not be effective if PLATYPUS determines, in its sole discretion that PURCHASER has altered or misused the goods or has failed to use or store them in accordance with instructions furnished by PLATYPUS. PLATYPUS's sole and exclusive liability and PURCHASER's exclusive remedy with respect to goods proved to PLATYPUS's satisfaction (applying analytical methods reasonably selected by PLATYPUS) to be defective or nonconforming shall be the replacement of such goods free of charge, upon the return of such goods in accordance with our instructions, although at its discretion, PLATYPUS may provide a credit or refund. If PLATYPUS manufactures custom goods for PURCHASER based on instructions, specifications, or other directions provided by PURCHASER, PLATYPUS shall not be liable for the lack of sufficiency, fitness for purpose or quality of the goods to the extent attributable to such instructions, specifications, or other directions. PLATYPUS shall not be liable for any loss, damage or penalty as a result of any delay in or failure to manufacture, deliver or otherwise perform hereunder due to any cause beyond PLATYPUS's reasonable control.

PLATYPUS shall not be liable for injury or damages resulting from the use or misuse of any of its products.



APPENDIX I: Determining Optimal Cell Seeding Concentration & Calcein AM Reagent Staining Protocol

This appendix is intended to assist in determining the cell seeding density needed to achieve confluency of your cell line when using the Oris™ Cell Invasion & Detection Assay. To that end, several dilutions of cell suspensions will be investigated.

NOTE: The Oris™ Detection Mask **MUST** be removed from the 96-well plate prior to the start of the following steps:

1. Collect cells by trypsinization or by non-enzymatic means (such as mechanical scraping) and calculate total number of cells.
2. Pellet cells by centrifugation and resuspend to a final concentration of 1,000,000 cells/mL in culture media.
3. Seed a 100 µl portion of cells, at 2-fold serial dilutions in the 96-well plate starting at 100,000 cells/well (a suggested starting amount), as shown below. Keep in mind that the cell seeding area of the well with the stopper in place is ~ 0.3 cm² and based on the typical seeding density of your cells, you can infer the appropriate cell number for your first serial dilution.

Column	1	2	3
Cells / well	100,000	50,000	25,000
Number of wells	6	6	6

4. Incubate the plate in a humidified chamber (37°C, 5% CO₂) for 16 hours (or ample time to allow for one doubling) with cell seeding stoppers in place.
5. Following cell attachment, remove the Oris™ Cell Seeding Stoppers from each well (see Figure 6) and **gently** wash the wells with PBS to remove non-adhered cells.
 - Secure the 96-well plate by holding it firmly against the deck of your work space. Slide the tines of the stopper tool under the backbone of the stopper strip, keeping the underside of the tool flush with the top surface of the plate.
 - Lift the stopper tool **vertically** to gently remove the stopper. Do not use the tool as a lever to pry the stoppers from the well as doing so may cause displacement of the seeded cells.
6. Visually inspect the cells by microscope to determine the cell seeding concentration that yields a confluent layer.



NOTE: If you plan to obtain the results of the Oris™ Cell Invasion & Detection Assay via colorimetric or microscopic analysis, you have successfully determined the optimal cell seeding concentration for your cell line. Proceed to Step 2 of the Cell Invasion Assay Protocol. If you plan to obtain the results of the Oris™ Cell Invasion & Detection Assay via a fluorescence plate reader, proceed with the following steps to optimize your plate reader settings.

7. The Oris™ Cell Invasion & Detection Assay has been designed to work with all types of fluorescent stains and staining techniques, but is supplied complete with the Calcein AM Reagent. The precise method for staining cells with fluorescent stains varies according to the nature of the individual stain. Please consult the manufacturer of your fluorescent stain for specific considerations.

Calcein AM Reagent Staining Protocol:

- a. To stain one fully-seeded 96-well plate, combine 5 µl of Calcein AM Reagent with 10 mL of phenol red-free and serum-free media or 1x PBS. Protect solution from light until ready to use.
 - b. Carefully remove culture medium from wells, being very careful not to disturb the BME overlay.
 - c. Wash wells with 100 µl of PBS (containing both Calcium and Magnesium).
 - d. Add 100 µl of diluted Calcein AM Reagent to each well.
 - e. Incubate plate at 37°C for 45 - 60 minutes.
 - f. Remove plate from incubator and remove staining solution. Do not aspirate using an in-house vacuum.
 - g. Wash wells with PBS, removing PBS at last step.
 - h. Attach mask and read in plate reader at 485/528 nm, sensitivity 55 nm.
8. Apply the Oris™ Detection Mask to the plate.
 9. Using the bottom probe of a fluorescence plate reader, obtain the total output from each well (adjust the gain settings to achieve optimal dynamic range). To determine optimal dynamic range, consider the following factors:
 - a) The gain setting that permits detection of the lowest concentration of cells.
 - b) The gain setting that permits discrimination between cell numbers at higher densities.



NOTE: When using a plate reader to analyze the Oris™ Cell Invasion & Detection Assay, it is important to stain cells using a fluorescent reagent that uniformly stains cells.

Once you have determined the optimal cell seeding concentration for your cell line, proceed to Step 2 of the Cell Invasion & Detection Assay Protocol.