WATER DISTILLER Stillo

USER MANUAL

UPDATING

The following table lists the updating history of the Operating Manual. The field "Description" shortly explains the modifications.

Edit.	Rev.	Date	Description
1	0	03-02	First issue
	1	03-03	Modified the e-mail address

TABLE OF CONTENTS

FOREWORDS	1
FOR YOUR SAFETY	2
PACKING CONTENTS DIMENSIONS AND WEIGHT LIST OF PACKING CONTENTS	3
PRODUCT DESCRIPTION	4
INSTRUCTIONS	5 5 5
APPENDIX A - TECHNICAL DATA	6
APPENDIX C - MAINTENANCE CLEANING THE OUTSIDE CLEANING THE BODY COOLING DUCT	7 7
APPENDIX D - TROUBLE SHOOTING	B
APPENDIX G - DECLARATION OF CONFORMITY	9
APPENDIX Z - CUSTOMER SERVICE	D

FOREWORDS thank you for the faith that you have given us by acquiring our water distiller **Stillo**, and we hope that the performances of this product result of your complete satisfaction.

In this manual you will find all the procedures for the correct use and the indications for the complete exploitation of the equipment performances.

We remain at your disposal for any more explanation, as well as suggestions turned to improve the product and the service.

WARNING SYMBOLS

PAY PARTICULAR ATTENTION TO THE PARAGRAPHS MARKED BY THE FINGER SYMBOL.



ATTENTION! THIS SYMBOL POINTS OUT A POTENTIAL DANGER FOR PEOPLE. PLEASE OPERATE ACCORDING TO THE SUITABLE PROCEDURES OF THE MANUAL IN ORDER TO PREVENT POSSIBLE DAMAGES TO THE USER AND/OR TO THIRD PARTIES.



CAUTION! This symbol points out a potential danger for property. Act according to the suitable procedures of the manual in order to prevent possible damages to material, equipment and/or properties.

USE OF THE EQUIPMENT

This manual is an integral part of the thermal sealer; keep it near to the machine in order to ensure easy and quick reference.

This equipment is designed for professional use only and must be used for the **production of distilled water**.

FOR YOUR SAFETY

The product should be used in compliance with the procedures described in the manual and never for purposes differing from the ones foreseen.



THE USER IS RESPONSIBLE FOR THE LEGAL FULFILLMENTS CONCERNING THE INSTALLATION AND USE OF THE PRODUCT.

IF THE PRODUCT IS NOT CORRECTLY INSTALLED AND USED, OR A SUITABLE MAINTENANCE IS NOT OPERATED, THE MANUFACTURER CANNOT BE CONSIDERED RESPONSIBLE FOR POSSIBLE BREAKS, MALFUNCTIONS, DAMAGES, LESIONS TO PROPERTY AND/OR PEOPLE.

In order to avoid any danger situations, with possible consequent damages or lesions to property and/or people, the following precautions should be followed:

- Place the distiller on a dry, flat and heat-tolerant surface and keep it away from any heat source.
- Because the cooling fan uses air to cool it, we suggest you place it in a well-ventilated space.
- Do not use a plastic bottle to store hot water.
- Unplug the power line if you wish to terminate the distillation process prematurely.





USE A POWER POINT WITH PROPER VOLTAGE AND DO NOT SHARE THE SAME OUTLET WITH OTHER ELECTRICAL APPLIANCES TO PREVENT OVERLOADING.



UNPLUG THE MAIN POWER CORD WHEN FILLING OR CLEANING THE DISTILLER OR WHEN NOT IN USE.

- Never immerse any part of the distiller in the water for the purpose of cleaning.
- Make sure that the electric plant is provided with earth connection in accordance with the current laws;
- Do not remove any label or plate from the equipment; in case call for new ones;
- **Use exclusively** original spare parts.



THE NON-OBSERVANCE OF WHAT ABOVE DESCRIBED MAKES ANY RESPONSIBILITY OF THE MANUFACTURER TO DECAY.



LIST OF PACKING The packing contains the following items: **CONTENTS** Water distiller; * ۵m Protective cover (Ref. 1) * * Water collection tank (Ref. 2); * Kit with 6 activated charcoal (Ref. 3) * Filter container (Ref. 4); * Residue cleaner bottle (Ref. 5); * Power cord (Ref. 6); User manual * STORE THE ORIGINAL PACKAGE FOR POSSIBLE FUTURE TRANSPORT OF THE EQUIPMENT. THE USE OF DIFFERENT PACKAGE COULD CAUSE DAMAGES TO THE PRODUCT DURING THE SHIPPING.

HOW IT WORKS



PRODUCTWater IDESCRIPTIONand effect

Water Distiller makes the purest and cleanest water in the most natural and effective way.

When the distiller is turned on, the water temperature rises gradually.
Chemicals, heavy metals and other contaminants, will be left as a residue at the bottom of the stainless steel tank. Impurities and unpleasant odours will be removed as steam collects on the surface of the cooling duct, which forms distilled water droplets and passes through the activated charcoal sachet. You then have sweet, pure water!



PLEASE DON'T DRINK THE FIRST 2000CC. OF WATER AS THIS MAY CONTAIN ANY RESIDUES AFTER MANUFACTURE.





Stillo

locom

APPENDIX A - TECHNICAL DATA

Equipment:	Water Distiller
Model:	Stillo
Manufacturer:	Megahome Corporation No. 24, Lane 44, Ping-Feng Road TAINAN 708 TAIWAN, R.O.C.
	for M.O.COM. S.r.I. Via delle Azalee, 1 20090 Buccinasco (MI) - Italy
Power supply:	220/230 V (other voltages on request)
Mains frequency:	50 Hz (other frequencies on request)
Fuse melting temperature:	142°C
External Dimensions (D x H):	200 mm x 360 mm
Nominal Power:	580 W
Insulation class:	Class I
Installation category:	Cat. II
Utilization:	Internal use
Environmental operating conditions:	Temperature:+15°C ÷ +35°CRelative humidity80% max (without condensation)
Net weight:	3,5 Kg
Capacity of the stainless steel reservoir:	4 liters

ENGLISH

APPENDIX C - MAINTENANCE

CLEANING THE OUTSIDE	Use household neutral detergent to clean the surface of the body and keep it dry (the rubber cushion within the upper cover is set with a specific direction; please do not remove it on your own).
CLEANING THE BODY	The interior of the stainless steel tank is smooth. Never use abrasive cleaners on your distiller. After each cycle rinse the inside of the boiling chamber with clean water. To remove heavy deposits, use the distiller cleaner as follows:
	1. Pour hot water into the stainless steel chamber until it covers the deposits.
	2. Pour one to four spoonfuls of cleaner (based on the amount of deposits) into the tank.
	3. Leave the upper cover off, plug in the main power plug, and turn on the power reset to heat the stainless steel tank until boiling.
	4. Discard the water fifteen to twenty minutes later, then wash the interior with clean water.
	5. Repeat the process until the chamber is clean.
	Small amounts of residue are harmless, and you can continue to use it safely.
COOLING DUCT	If the condensing coil is clogged with dust, the efficiency of distillation may be affected. To clean the condensing coil, remove all three screws from the underside of the cover. Use a tooth brush or a vacuum cleaner to remove the accumulated dust.
	TAKE CARE" DO NOT IMMERSE THE MOTOR IN THE WATER WHEN YOU



ARE CLEANING THE COVER. WE RECOMMEND YOU REGULARLY ADD A LITTLE OIL ON THE MOTOR SHAFT TO HELP MAINTAIN THE MOTOR IN GOOD WORDING ORDER.

APPENDIX D - TROUBLE SHOOTING

PROBLEMS

Not making any distilled water

- Main power cord not plugged into power point or lid cord not plugged into base of the distiller.
- No power to electrical outlet.
- Reset button not pressed.
- -

Not making enough distilled water during each cycle

- Water not filled to maximum fill line.
- Lid gasket not sealing properly or up the wrong way.
- Collection bottle not placed directly under the spout.

Water or steam leaks from spout

- Carbon sachets not properly placed in the nozzle.
- Collection bottle not correctly lined up with distiller.



CONTACT YOUR DISTRIBUTOR FOR FURTHER ADVICE.



APPENDIX G - DECLARATION OF CONFORMITY

Sistemi innovativi di sterilizzazione	DECLARATION OF CONFORMITY Application of the EEC Directive 89/336				
Name of the Manufactur	rer: Megahome Corporation				
	for M.O.COM. S.r.I. Via delle Azalee, 1 20090 Buccinasco (MI) – Italy				
Address of the Manufac	turer: No. 24, Lane 44, Ping-Feng Road TAINAN 708 TAIWAN, R.O.C.				
Product description:	Water Distiller				
Model:	STILLO				
Made in:	TAIWAN, R.O.C.				
The undersigned hereby declares that the above mentioned goods entirely					
	conform				
to the	EEC Directive 89/336 (and following updates).				
Reference standards:	EN 61000-3-3: 1995 EN 55014-2: 1997 EN 55014-1: 1993+A1+A2 EN 61000-3-2: 1995+A1+A2+A14/01.00				
01/01/2002	t. Qua				
Date	Signature				
	Alfio VILLA Name and Surname				
	Legal Agent Function				

APPENDIX Z - CUSTOMER SERVICE

The Customer Service by M.O.COM. Ltd. Co. is at complete disposal of the Customers for any request or technical information concerning the product, as well as for any suggestions on the steam sterilization procedures.

Please contact at the following address:

M.O.COM. Srl Customer Service

Via delle Azalee, 1 20090 Buccinasco (MI) ITALY

Tel.	(+39) 02-45701505
Fax	(+39) 02-48840514
e-mail	at@mocom.it

For improving the quality of the product and service, please send Your comments and/or suggestions to the followings **e-mail** addresses:

uc@mocom.it (marketing and sales)

Besides, You can send comments and/or suggestions by **letter or fax** to the above indicated address.

Thanks in advance for the valuable assistance you would supply to us.