

User manual Pico4, Pico8 og Pico16 – short version

Revisionlist

Date	Name	Description
22.12.14	Bent Sollid	First version.

Introduction

This is a short version of the user manual for Picomeds IR-transmitters Pico4, Pico8 and Pico16. For more information, please use the complete user manual.

Overview

The IR-transmitters have 4, 8 and 18 buttons and 2 pages. Functions:

- Learnable IR-signals.
- Built in IR-signals.
- Scanning, 1 switch.
- Integrated fingerguide.
- Possible to choose continious IR-signals or not when «press & hold» on keys.
- Programmable: sound on keypress, spasm-filter, parameters for scanning, backup.

The figure shows an overview of a Pico4. Pico8 and Pico16 have the same functions.



If **page 2 is be enabled** the Pico4 and Pico8 use the lower left key as a page-select key. In this case it will not transmit any IR-signals. Pico4 will have a maximum of 6 channels and Pico8 will have a maximum of 14 channels if page 2 is enabled. Pico16 has a dedicated page select key.


Lock-code key transmits IR-signals in the same way as the other keys. It is also able to transmit «lock-code» and «common code». It may be programmed as a page-independent key.

The battery is rechargeable and it will normally not be necessary to exchange it.

- Low battery: red light indicator flashes for 2 seconds when it transmits IR-signals.
- Charging: green light flashes slowly. It has a steady green light when it is full charged.
- Typical time of charging for an empty battery is 4-5 hours.



Programming – on screen



Android products uses an app downloaded from Google-play. Search for «Picomed» and install the app «Picomed Starter».

- Connect the «OTG» cable to Picomeds IR-transmitter and the USB cable to your unit.
 When a USB icon is shown on screen:
 Open the Picomed Starter app, choose «Enable Connection» and enable «USB-internet sharing», go one step backward in the app on your unit.
- Choose «Configure Product». The app will download a web page from Pico4/8/16.

Apple products will work without any installation on iPhone 4 or newer, and for iPad 3 or newer which has a SIM-card holder.

- You have to turn on mobile data on the unit to let it use the USB cable for downloading an internet page from the Picomed transmitters. Tap Settings Cellular Cellular data = on. If this alternative is not available it may be that your unit or your mobile supplier not has this option.
- Connect the «OTG» cable to Picomeds IR-transmitter and the USB cable to your unit. Then go to this adress in browser: picomed.local

Windows PC has to install a driver the first time you are going to use it. Download the driver from Picomeds web-site. It may be different drivers for different Windows versions.

- Use a web-browser and go to www.picomed.no.
- \bullet Choose: Products IR-transmitters Pico4/8/16 and follow instructions on screen. Then:
- Connect USB cable between computer and IR-remote, open a web-browser and go to: picomed.local



Programming – on IR-transmitters

Without a web-browser it is possible to copy IR-signals into it and to enable/disable commoncode.

- Enter programming mode: press and hold keys 1, 2 and 3 for 3 seconds \rightarrow long beep.
- **Leave programming mode**: press and hold key «1» for 3 seconds, or wait for timeout after 60 sec without any keystrokes.
- If low battery: red flashing light for 2 seconds, then it leaves programming mode.
- To enable/disable **commoncode**: press and hold lock-code key for 3 seconds. Commoncode off: a short signal from beeper, light in lock code key and the green light. Commoncode on: a long signal from beeper, light in lock code key and the green light.
- Norwegian version transmits IR = Picomed, usercode 1, channel 1 (also named: NA1-1).
- Swedish version transmits IR = Gewa, channel 16.
- Copying IR signals into Picomed IR-transmitters:
- 1. Position Pico4/8/16 and the other IR-transmitter against each other.





- 2. Press shortly on the actual key on Pico4/8/16: steady green light in actual key.
- 3. Press shortly on actual key on the other IR-remote, it transmits IR to Pico4/8/16.
- 4. When IR is received by Pico4/8/16, it will respond with single green flashing.
- 5. Press once again shortly on actual key on the other IR-remote.

 A successful IR-learning will be indicated by 3 green flashing, then light is turned off.
- 6. Go back to step 2 to copy an IR-signal to an other key.
- 7. Leave programming mode when finished.
- IR-learning will be interrupted after 20 seconds if no IR is received. It will then go to ordinary programming mode.
- To switch to page 2 in programming mode; press and hold page select key for 3 seconds.









Programming map

