

Biometra®

a Whatman company


TGGE System

230 V

Code-No. 024-000

115 V

Code-No. 024-090



Manual

March 1999



!! Warning !!

**Please read these instructions carefully
before using this apparatus!**



Biometra biomedizinische Analytik GmbH

Rudolf-Wissell-Straße 30, D-37079 Göttingen

P.O. Box 1544, D-37005 Göttingen

Tel: ++49 – (0)5 51 / 50 68 6-0 * Fax: ++49 – (0)5 51 / 50 68 6-66

e-mail: info@biometra.co.uk

internet: <http://www.biometra.de>

The **TGGE method is covered by patents** issued to Diagen (now QIAGEN GmbH).
The polymerase chain reaction (PCR) process is covered by patents issued to Hoffmann-La Roche.

Acryl-Glide is a trademark of Amresco Inc.
Biometra is a trademark of Biometra GmbH.
Whatman is a trademark of Whatman International Ltd.

The POLAND software service established by Gerhard Steger, Department of Biophysics, University of Duesseldorf, is available by internet <http://www.biophys.uni-duesseldorf.de/service/polandform.html>.

Table of Contents

1	The TGGE System	4
1.1	Introduction	4
1.2	Principle of method	4
1.3	Special features of the TGGE System	5
1.4	How to start with TGGE analysis	6
2	General Recommendations	8
2.1	Safety warnings	8
2.2	Notes for use.....	8
3	Components of the TGGE System	9
3.1	Electrophoresis unit	10
3.2	Controller unit.....	17
3.2.1	Instrument keys and ports	17
3.2.2	Programming of the TGGE Controller.....	18
4	Sample preparation	24
4.1	Purity of samples	24
4.2	Quantity and volumes of samples	25
5	Setting up polyacrylamide gels	25
5.1	Selecting Concentration of PAA gels	25
5.2	Setting up the gel solution.....	25
5.3	Some remarks corresponding to standard TGGE conditions	28
5.4	Assembling the gel sandwich.....	30
5.5	Disassembling the gel sandwich	32
6	Electrophoresis with the TGGE System	33
6.1	Electrophoresis conditions	33
6.2	Preparing the electrophoresis unit	33
6.3	Perpendicular TGGE.....	35
6.4	Parallel TGGE.....	39
6.5	Silver staining.....	44
6.6	Ethidium bromide-staining	48
6.7	Blotting	48
6.8	Autoradiography.....	48
6.9	Elution of DNA from the TGGE gel	48
7	TGGE in analysis of point mutations in dsDNA	49
7.1	 Theoretical background of a detection rate approximating 100% for point mutations - calculations with the POLAND program	49
7.2	The “old” Poland program (Old Server).....	51
7.2.1	About the Poland service.....	51
7.2.2	Program-specific information	51
7.2.3	How to use the “old” Poland program (Standard)	53
7.3	The “new” Poland program (New Server)	55
7.3.1	About the Poland service.....	55
7.3.2	Program-specific information	55
7.3.3	References for Poland Service	56
7.3.4	HELP for Poland Service:	56

7.3.5	How to use the “new” Poland program	65
7.4	The optimized DNA fragment.....	66
7.4.1	Asymmetric GC-clamps for PCR primers used for TGGE analysis.....	68
7.4.2	Chemical clamp with Psoralen (Furo[3,2-g]coumarin, C ₁₁ H ₆ O ₃).....	68
7.4.3	POLAND analysis of samples.....	69
8	Optimizing parallel TGGE by perpendicular TGGE.....	70
8.1	Check short DNA fragments for their melting behavior	70
8.2	From perpendicular to parallel TGGE	71
9	TGGE / SSCP.....	72
9.1	Running an SSCP on the TGGE.....	72
9.2	DNA sample preparation.....	73
9.3	Gel casting	73
9.4	Electrophoresis	75
9.5	Routine analysis.....	75
10	TGGE in RNA analysis	76
10.1	Completely double-stranded RNA.....	76
10.2	Partly double-stranded RNA, e.g. viroid RNA	76
10.3	Single-stranded RNA with single hairpin structures, m-RNA secondary structures	76
10.4	Staining.....	77
11	TGGE in protein analysis	78
11.1	Buffers	78
12	Trouble-shooting	80
13	TGGE Testkit.....	93
13.1	Introduction	93
13.2	Protocol.....	93
13.2.1	Gel composition:.....	93
13.2.2	Running buffer:	94
13.2.3	Electrophoresis parameters:.....	95
13.3	Gel images.....	95
14	Appendix	96
14.1	Technical Data	96
14.1.1	System.....	96
14.1.2	Electrophoresis Chamber	96
14.1.3	TGGE System Controller with integrated power pack.....	96
14.2	Buffers	97
14.3	Silver staining solutions:	99
15	References	101
16	Order information and spare parts.....	105
17	Instructions for return shipment	106

18 Equipment Decontamination Certificate.....	107
19 Warranty	108

Declaration of Conformity

1 The TGGE System

1.1 Introduction

Temperature Gradient Gel Electrophoresis is a new and powerful electrophoresis method for separation of nucleic acids like DNA or RNA or for analysis of proteins. The TGGE method, which is covered by patents, uses the temperature dependent change of conformation for separating molecules (for review see Reference 1).

Since the introduction of the first commercial available TGGE apparatus in 1989 temperature gradient gel electrophoresis has gained high interest in scientific and clinical laboratories due to the unprecedented resolution capability and easiness of analysis. The range of scientific publications using the TGGE method is broad and covers all disciplines which use molecular biology methods: e.g. Oncology²⁻⁴, Virology^{5,6}, Immunology^{7,8}, RNA Viroid Research⁹⁻¹², Prion Research¹³, Population Analysis¹⁴⁻¹⁵. The TGGE method has also been used for quantitative analysis in industry¹⁶⁻¹⁷ and for conformational analysis of proteins¹⁸⁻¹⁹.

1.2 Principle of method

Conventional protein or nucleic acid electrophoresis separates molecules mainly according to size or charge. **TGGE** adds a new parameter for separation, namely the **shape of the molecule**.

The shape is mostly determined by the secondary and tertiary structure of the molecule and can be changed by external influences like temperature, salt concentration, pH etc.

The conformation both of proteins and nucleic acids depend on weak binding forces like hydrogen bonds or van der Waals bonds. Increasing the temperature above a certain limit breaks down these bonds. The molecules will adopt a so called denatured conformation in contrast to the native one.

E.g. with DNA it is possible to determine the temperature which is necessary to break down hydrogen bonds along double stranded DNA. This temperature is called midpoint of transition (T_M) or melting temperature and characteristic for a certain stretch of DNA (see figure 1). TGGE uses the melting temperature to identify DNA which differs in sequence among a mixture of molecules of the same size.

TGGE therefore not only separates molecules but gives additional information about the sequence (DNA or RNA) or the stability of proteins.




Figure 1: Schematic drawing of different conformations of DNA during temperature gradient gel electrophoresis.

1.3 Special features of the TGGE System

The microprocessor controlled gradient block of the TGGE System allows strictly defined linear gradients with high resolution. A run distance of 2 mm which can easily be detected by eye corresponds to a maximum temperature difference of about 0.6°C. Therefore even slightest differences of molecules can be detected by the new TGGE System.

Because of the small amount of material used for separation DNA or RNA fragments appear as fine bands which can clearly be distinguished from each other. Even complex band patterns can be analyzed due to the high resolution capability of the gradient block.

Comparing the TGGE method with another screening method like SSCP shows superior performance of the TGGE method²⁰⁻²².

The controlled temperature conditions make repetition of experiments easy and lead to reproducible gel results. The small format of the gradient block has been optimized in order to reduce sample volume and especially to save experimental time.

Perpendicular and parallel TGGE are two different modes applicable with the Biometra TGGE System without need for specialized parts or equipment.

Whereas perpendicular TGGE is mostly used for defining optimal separation conditions, parallel TGGE allows the analysis of multiple samples (e.g. screening).

perpendicular TGGE: temperature gradient is perpendicular to the electrophoretic run direction

→ one sample is spread over a broad temperature range

parallel TGGE: temperature gradient is parallel to run direction

→ multiple samples are spread over a narrow temp. range




Figure 2: Schematic drawing of typical results after perpendicular TGGE (left panel) and parallel TGGE (right panel).

1.4 How to start with TGGE analysis

Starting with the gene of interest one develops the right combination of PCR primers for amplification of the desired gene fragment. The polymorphic site must be located inside the amplicon and not at the far end. The POLAND software helps one to identify suitable primers and the adequate fragment length of the amplicon. This program gives a rough estimation too, what temperature parameters will fit best to the desired separation²³⁻²⁵.

A revised version of the POLAND program can be found on the world wide web ([http:// www.biophys.uni-duesseldorf.de/service/polandform.html](http://www.biophys.uni-duesseldorf.de/service/polandform.html)).

For example the polymorphic site is represented by allele A and allele a. The two alleles can either exist as homoduplices (AA or aa) or heteroduplices (Aa or aA) (see figure 3).




Figure 3: Schematic drawing of double stranded DNA with polymorphism „A“ or „a“ (left panel) and corresponding DNA sequence (right panel). Each line represents double stranded DNA.

The TGGE System used as perpendicular TGGE (see above, figure 2) gives the possibility to identify the different alleles by their individual melting behavior. Samples with homoduplex „AA“ or „aa“ have a distinct melting temperature (e.g. T_{m1} and T_{m2} , see figure 4), at which double stranded DNA separates into branched DNA. At even higher temperature the branched DNA separates into individual strands.

After perpendicular TGGE, a heteroallele sample like „Aa“ normally shows four different T_m values. The PCR amplification of a heteroallele sample results in four different double stranded DNA types: The two homoduplices "AA" and "aa" as well as two heteroduplices "Aa" and "aA", which have a non-pairing base at the polymorphic site. This non-pairing base will lead to a shift of the T_m to lower values (T_{m3} and T_{m4}).


Perpendicular TGGE shows at which temperature the different DNA strands will separate. For future analysis by perpendicular or parallel TGGE a narrower temperature range which includes the T_M values of homo- and heteroduplices can be used.

Please remember always, that wild type and mutant DNA have to be mixed before PCR to get the 4 bands with different melting behavior. Sometimes the melting difference between heteroduplices cannot be resolved but remember that 3 bands on the TGGE gel are enough to detect a mutation.



*Figure 4: Schematic drawing of melting behavior of double stranded, heterozygotic DNA with allele type „Aa“ after **perpendicular TGGE**.*

Screening of multiple samples is performed by using parallel TGGE (see above, figure 2). Parallel TGGE looks like conventional SSCP-analysis, but has a higher probability to identify possible mutants.



*Figure 5: Schematic drawing of a screening for double stranded, heterozygotic DNA with allele type „Aa“ after **parallel TGGE**. (Sometimes the melting difference between heteroduplexes can not be resolved.)*

Different plates with pre-fixed slots for 8, 12 or 18 samples are available for screening purposes (see chapters 6.4 and 15).

2 General Recommendations

2.1 Safety warnings



Check the voltage of the power supply and of the control unit before use.

In any case of malfunction of the power supply, controller or electrophoresis chamber, do not open the case but contact Biometra or your local distributor.



First switch off the power switch of the power supply before opening the lid of the electrophoresis chamber if you want to interrupt or stop the electrophoresis run.

During electrophoresis don't touch the electrode wires or the buffer inside the electrophoresis chamber. High voltage. Danger of life!



Whenever polyacrylamide gels are handled pay attention to standard laboratory safety regulations, e.g. wear lab coat, protective gloves and eye shield. Polyacrylamide is neurotoxic.

2.2 Notes for use

- Do not scratch **the protective foil** of the gradient block. In case of a damaged foil contact Biometra or your local distributor for a replacement foil.
- Do not use strong acids or basic solutions or organic solvents for **cleaning** glass plates, the electrophoresis chamber or the gradient block.
- Do not incubate glass plates over night in **cleaning solution**.
- Wear **protective gloves** during all steps of the silver staining protocols to avoid staining artifacts due to the high sensitivity of the staining protocol.

3 Components of the TGGE System

The TGGE System contains all components which are necessary to get started. All kinds of TGGE applications (parallel or perpendicular TGGE, Constant Temperature GE, Time resolved TGGE) can be run with the System. For certain applications which need different numbers or sizes of sample slots the adequate parts are available and can be ordered from **Biometra** or your local distributor (see 6.3. Order Information).

The **TGGE System** (Order number: 024-000) consists of:

• <u>TGGE Controller</u> with integrated power supply, 100 program stores and control function of temperature and electrophoresis conditions		<u>024-001</u>
• <u>TGGE-Electrophoresis unit</u> with 2 removable buffer chambers, Peltier-element powered gradient block and control cable		<u>024-002</u>
• <u>TGGE Starter Kit</u> contains		<u>024-003</u>
3 plane „Bonding“ glass plates	024-021	
1 glass plate with 8 slots for parallel TGGE	024-022	
1 glass plate with 1 rectangular slot for perpendicular TGGE	024-023	
1 glass pate with 12 slots for parallel TGGE	024-025	
Electrophoresis wicks (100/pkg)	024-015	
Polybond film (25/pkg)	024-030	
1 Cover glass plate and 10 cover films	024-031	
1 Acryl-Glide™ (100 ml)	211-319	
• Manual		

When unpacking your System please check whether all mentioned parts are included. If individual parts are missing call **Biometra** or your local distributor.

3.1 Electrophoresis unit

The electrophoresis unit consists of 4 parts:

- ❶ 1 safety lid with 2 electric plugs (anode and cathode)
- ❷ 2 removable electrophoresis chambers each with platinum wires and electric connectors (volume: max. 250 ml)
- ❸ housing with Peltier-element powered gradient block
- ❹ 37 pin connecting cable to control unit




Figure 6: Parts of the electrophoresis unit





Figure 6b: Gradient block with temperature lines and marks for positioning the gel slots (top). Safety lid with arrow indicating the running direction of nucleic acids (middle). Gradient block covered with safety lid; setup for a perpendicular TGGE run (bottom, left); setup for a parallel TGGE run (bottom, right).

The gradient block is centrally positioned in the middle of the electrophoresis unit and protected by a layer of white foil. This foil is necessary to protect the electronic parts beneath it from liquid, buffer or other harmful chemicals.

If the protection foil has been scratched during use
stop working and
exchange the protection foil with a new one.

The two opposite sides of the gradient block are marked with lettering T_1 and T_2 . Beneath these symbols the Peltier-elements which build up the temperature gradient during electrophoresis can be found. Both sides of the gradient block can reach any preset temperature from 15°C - 80°C. The orientation of the temperature gradient, i.e. which side of the gradient shall be cold or hot, can be freely determined.

Between symbols T_1 and T_2 six thick lines (L1 to L6) and five thin lines (not coded) are marked on the block, which represent the entire linear range of the gradient block (see figure 7). The temperature difference between two lines is identical from line to line. E.g. if T_1 is 30°C and T_2 is 75°C, the temperature difference is 7.4°C between two thick lines or 3.7°C between a thick and a thin line.



When performing parallel TGGE the beginning and the end of the linear temperature gradient are represented by the first and the last line. Whereas when performing perpendicular TGGE the ends of each line represent the linear temperature range. When choosing a temperature gradient e.g. between 25°C and 65°C these two temperatures can actually be found at the first marked line (figure 7: 10mm distance to the block edge) and at the last marked line (figure 7: 50 mm distance) on the block. The block areas to the left and right of these lines are slightly hotter respectively cooler (see figure 7 and table 1 + 2).

Linearity of block





Figure 7: upper panel: *Temperature profile of gradient block measured by a micro sensor on top of the gradient plate every 1mm beginning from one edge of the block. 10 mm distance and 50 mm distance correspond to first (L1) and last marked thick line (L2) on the gradient block respectively.*


lower panel: *Schematic drawing of the block.*

The maximum temperature difference between the two sides of the gradient block (T_1 and T_2) during electrophoresis is 45 Kelvin. That means it is possible to build up a gradient between 35°C and 80°C or between 25°C and 70°C, just to give two examples.

Programming T_1 and T_2 , the actual temperature of L_1 to L_6 can be calculated by the following formula:

$$L_n = L_1 + (n - 1) \Delta\vartheta \quad (L_n = \text{Temperature of line } n; n = 1 \dots 6)$$

$$\Delta\vartheta = \frac{L_6 - L_1}{5} \quad (\Delta\vartheta = \text{Temperature difference between two thick lines})$$



When leaving out the gradient function the block can be cooled down to 4°C or heated up to 80°C. Although Peltier elements reach lower respectively higher temperature values, the surrounding plastic material does not permit the temperature range of the TGGE System to be extended.

The electrophoresis buffer chambers can freely be positioned around the gradient block. This makes it easy to switch between parallel or perpendicular TGGE. When the electrophoresis buffer chambers stay parallel to the lines of the gradient block (see figure 8 left panel) parallel TGGE applications can be run. By simply switching the chambers perpendicular to the lines of the gradient block (see figure 8 right panel) it is possible to run perpendicular TGGE.




Figure 8: Orientation of the two electrophoresis buffer chambers (- represents cathode, + represents anode) relative to the centrally positioned gradient block.

Programming L1 and L6 (°C):								
T1	L1	L2	L3	L4	L5	L6	T2	DELTA
15,94	20	27,4	34,8	42,2	49,6	57	61,06	7,4
16,94	21	28,4	35,8	43,2	50,6	58	62,06	7,4
17,94	22	29,4	36,8	44,2	51,6	59	63,06	7,4
18,94	23	30,4	37,8	45,2	52,6	60	64,06	7,4
19,94	24	31,4	38,8	46,2	53,6	61	65,06	7,4
20,94	25	32,4	39,8	47,2	54,6	62	66,06	7,4
21,94	26	33,4	40,8	48,2	55,6	63	67,06	7,4
22,94	27	34,4	41,8	49,2	56,6	64	68,06	7,4
23,94	28	35,4	42,8	50,2	57,6	65	69,06	7,4
24,94	29	36,4	43,8	51,2	58,6	66	70,06	7,4
25,94	30	37,4	44,8	52,2	59,6	67	71,06	7,4
26,94	31	38,4	45,8	53,2	60,6	68	72,06	7,4
27,94	32	39,4	46,8	54,2	61,6	69	73,06	7,4
28,94	33	40,4	47,8	55,2	62,6	70	74,06	7,4
29,94	34	41,4	48,8	56,2	63,6	71	75,06	7,4
30,94	35	42,4	49,8	57,2	64,6	72	76,06	7,4
31,94	36	43,4	50,8	58,2	65,6	73	77,06	7,4
32,94	37	44,4	51,8	59,2	66,6	74	78,06	7,4
33,94	38	45,4	52,8	60,2	67,6	75	79,06	7,4
34,94	39	46,4	53,8	61,2	68,6	76	80,06	7,4
35,94	40	47,4	54,8	62,2	69,6	77	81,06	7,4
36,94	41	48,4	55,8	63,2	70,6	78	82,06	7,4
37,94	42	49,4	56,8	64,2	71,6	79	83,06	7,4
38,94	43	50,4	57,8	65,2	72,6	80	84,06	7,4

Table 1: Examples for actual temperatures on the gradient block programming L1 and L6


Programming T1 and T2 (°C):								
T1	L1	L2	L3	L4	L5	L6	T2	DELTA
20	24,05	31,43	38,81	46,19	53,57	60,95	65	7,38
21	25,05	32,43	39,81	47,19	54,57	61,95	66	7,38
22	26,05	33,43	40,81	48,19	55,57	62,95	67	7,38
23	27,05	34,43	41,81	49,19	56,57	63,95	68	7,38
24	28,05	35,43	42,81	50,19	57,57	64,95	69	7,38
25	29,05	36,43	43,81	51,19	58,57	65,95	70	7,38
26	30,05	37,43	44,81	52,19	59,57	66,95	71	7,38
27	31,05	38,43	45,81	53,19	60,57	67,95	72	7,38
28	32,05	39,43	46,81	54,19	61,57	68,95	73	7,38
29	33,05	40,43	47,81	55,19	62,57	69,95	74	7,38
30	34,05	41,43	48,81	56,19	63,57	70,95	75	7,38
31	35,05	42,43	49,81	57,19	64,57	71,95	76	7,38
32	36,05	43,43	50,81	58,19	65,57	72,95	77	7,38
33	37,05	44,43	51,81	59,19	66,57	73,95	78	7,38
34	38,05	45,43	52,81	60,19	67,57	74,95	79	7,38
35	39,05	46,43	53,81	61,19	68,57	75,95	80	7,38
36	39,96	47,18	54,39	61,61	68,82	76,04	80	7,216
37	40,87	47,92	54,97	62,03	69,08	76,13	80	7,052
38	41,78	48,67	55,56	62,44	69,33	76,22	80	6,888
39	42,69	49,41	56,14	62,86	69,59	76,31	80	6,724
40	43,6	50,16	56,72	63,28	69,84	76,4	80	6,56

Table 1: Examples for actual temperatures on the gradient block programming T_1 and T_2

3.2 Controller unit


The controller is a highly integrated, micro processor driven unit for controlling the temperature, ramping time and ramping rate of the gradient block as well as supplying the power for the electrophoresis unit. For entering and storing run parameters the front panel of the controller offers 4 function keys and a full numerical key pad. During the run the display of the controller continuously shows the current parameters.

3.2.1 Instrument keys and ports



- ❶ Power on/off switch
- ❷ Display
- ❸ 4 function keys A, B, C, D
- ❹ Alphanumerical key pad

Figure 9: *Front panel of the TGGE System Controller.*



- ❶ Computer port (RS 232)
- ❷ Connectors for electrode cable:
red = anode, black = cathode
- ❸ Printer port
- ❹ Interface to electrophoresis unit
- ❺ Mains and fuses

Figure 10: *TGGE System Controller from the rear*

3.2.2 Programming of the TGGE Controller

After switching on the controller the display shows the instruments name and the software version. Immediately afterwards the main menu appears.

Main menu

T1: 22.0°C	T2: 22.0°C
block off	
A? B Elpho C programs D	
+	

At the bottom line of the display 4 possible options which can be retrieved

by the 4 functions keys ↵A, B, C, D are shown. These 4 options change

during programming relative to the chosen menu.

(Temperatures T1 and T2 shown in the display are dependent on the room

temperature.)

↵A “?”: comments or tips about the current program step

↵B “Elpho”: commands to load and start the program

↵C “programs”: commands to edit new programs, change or delete existing programs

↵D “+”: different options like printing of program stores or of running protocols, choice of language, choice of signal

In general select by scrolling, activate by pressing enter; except when selecting program numbers or temperature values!

Function key D: Options

1 print programs
2 signal
3 language
A ↑ B ↓ C quit D enter

1: Printing of program stores. A dot matrix printer can be connected to the controller by using the port at the rear (see figure 10).

2: Choice whether a beep signal can be heard at the end of a program or when the program has reached an infinite time step

3: Choice between ↵A “German” or ↵B “English”

4 standard mode
5 test mode
6 void
A ↑ B ↓ C quit D enter

4: not occupied

5: not occupied

6: not occupied

↵A ↑ and ↵B ↓ allow scrolling of display.

By ↵C “quit” you will return to the main menu.

Function key C: Editing of programs

program no: __

A list **B** del **C** quit **D** enter

After entering a number of a non occupied program store, the display shows:

Name: > _____ <
 ABCDEFGHIJKLMNOPQRST
 UVWXYZ – () α σ /, <> & +. %
A → **B** ABC **C** quit **D** enter

In the main menu pressing $\rightarrow$ C will offer the possibility to edit a new program. First choose a program store number.

$\rightarrow$ A “list” displays all program stores from 0 to 99 with names or the information <empty>. $\rightarrow$ B “del” deletes the last entry.

To give the program an individual name strike $\rightarrow$ B “ABC” to jump with the cursor into the letter field. Moving inside this field is possible by the keys $\rightarrow$ A “→” and $\rightarrow$ B “←”. Pressing $\rightarrow$ D “enter” the current high lighted letter is stored in the name field. This step can be repeated 8x times. Pressing $\rightarrow$ D “enter” two times leads to the next step.

1:L1: __ L6:
 alternative:
 T1: __ T2:
A ? **B** T1 **C** quit **D** →

The temperature gradient between T1 and T2 must not exceed 45 K.

It is now possible to enter for the first program step temperature values for both sides of the gradient block:
 Alternatively the temperatures for L1 and L6 (first and last thick lane on the gradient block) or for T1 and T2 (left and right edge of the gradient block) can be programmed. The change between programming L and T can be done by pressing $\rightarrow$ B “T1”. If a number has been entered at the field L1, you have to confirm by pressing $\rightarrow$ D “enter”. The cursor jumps to L6. After programming L6 you will be asked “ok ?”: Pressing $\rightarrow$ B “no→L1” allows new programming of L1 and L2. Pressing $\rightarrow$ D “yes” confirms the temperatures.

1:L1: 25.0°C L6:
 Alternative:
 T1: __ T2:
A ? **B** delete **C** quit **D** →

After entering a temperature for T1, T2, L1 or L6 the temperature can be Deleted by pressing $\rightarrow$ B “delete”.

1:L1: __ L6:
 Alternative:
 T1: __ T2:
A ? **B** L1 **C** quit **D** →

If programming of T1 and T2 (left and right edge of the gradient block) has been done but L1 and L6 are preferred, changing to L1 and L6 can be done by pressing $\rightarrow$ B “L1”.

1: L1: 25.0°C L6: 60.0°C
 T1: 21.1°C T2: 63.8°C
 Ok ?:
A ? **B** no→L1 **C** quit **D** →
Or

After entering L6 or T2 all four temperatures (L1, L6, T1, T2) are displayed.

“ok ?”: Pressing $\rightarrow$ D “yes” leads to the programming of the Electrophoresis parameters.

Pressing $\rightarrow$ B “no→L1” leads back to programming L1. Pressing $\rightarrow$ B

“T1” allows programming of T1 and T2.

Pressing $\rightarrow$ B “no→T1” leads back to programming T1. Pressing $\rightarrow$ B “L1” allows programming of L1 and L6.

1: L1: 25.0°C L6: 60.0°C
 T1: 21.1°C T2: 63.8°C
 Ok ?:
A ? **B** no→T1 **C** quit **D** →

Depending on the decision of programming L1 and L6 or T1 and T2 the programmed temperatures will be shown in the display. After programming L1 and L6 this temperatures will be shown in the display during the next programming steps.

1:L1: 25.0°C L6: 60.0°C
time: __0m 0s
El: 0V 500mA 30W
A ? B V*h C quit D →

Time: You can choose electrophoresis time. Standard values for a TGGE run are 30 – 45 min. If you enter 30 and confirm by pressing ↵D “enter” you get 30 s. If you enter 30 • you will get 30 min. If you enter 30 ** you will get 30 h.

1:L1: 25.0°C L6: 60.0°C
time: 30m 0s
El: 0V 500mA 30W
A ? B V*h C quit D →

V*h: Pressing ↵B “V*h” replaces times by the more precise Volt x hour integration.

1:L1: __ L6:
V*h: __0.00 Vh
El: 0V 500mA 30W
A ? B Time C quit D →

Pressing ↵B “Time” replaces V*H by time.

1:L1: 25.0°C L6: 60.0°C
time: 30m 0s
El: __ 0V 500mA 30W
A ? B special C quit D →

El: Three different electrophoresis parameters (voltage, current or wattage) can be set. Current and wattage are pre-set at max. values of 500 mA and 30W respectively. In the beginning we recommend to set only the Voltage. Depending on the resistance of the gel electrophoresis the controller will regulate the other two parameters automatically.

1:L1: 25.0°C L6: 60.0°C
time: 30m 0s
El: 250V 500mA 30W
A ? B special C quit D →

After confirming the voltage by pressing ↵D “enter” and pressing two times ↵D “→”, the following two parameters are not changed.

1: special functions
ramptime: __0m 0s
A ? B standard C quit D →

Ramptime = Ramping time
Pressing ↵B “special” gives you the choice to choose how fast the Gradient block is going to the established gradient. Normally you choose 1s which means maxi ramping speed. In this case enter “1” and confirm with ↵D “enter”
Pressing ↵B “standard” results in the standard display.

1:L1: 25.0°C L6: 60.0°C
Time: 30m 0s
El: 250V 500mA 30W
A ? B special C quit D →

2:L1: __ L2:
alternative:
T1: T2:
A ? B T1 C quit D →

Pressing ↵D “→” starts the programming of step 2 in this program.

program no:
pgm end:step(s)
runtime:h.....m.....s

By pressing ↵C “quit” any change will be saved and the following messages appear:

L1: 22.0°C L6: 22.0°C
block off
A? B Elpho C programs D +

After a few seconds the main menu appears.

Function key B: Start/Stop

Start block function and electrophoresis

Start program: __

A list **B** del **C** quit **D** enter

In the main menu pressing key $\rightarrow$ B “block” offers the possibility to start a program.

L1: 25.0 °C L6: 60.0°C

ramp: 1 time: 0m 1s
El: 250 V 20mA 20W
A list **B** block **C** program **D** +

After entering a program number or choosing a program from the list ($\rightarrow$ A “list”) the block starts to establish the gradient. The timer for the ramping starts immediately. The electrophoresis is started as soon as the gradient block reaches the programmed temperature (gradient. The limiting factor (const. V, mA or W) is indicated by an blinking arrow.)

L1: 25.0 °C L6: 60.0°C

hold: 1 0m 1s 0.00 Vh
El: 250 V 8mA 20W
A list **B** block **C** program **D** +

After establishing the gradient, line two of the display changes. The timer now starts again and counts the electrophoresis running time. Additionally the volt/hour integrator starts to count.

Pressing $\rightarrow$ C “program” during block run

Program no: 4

Pgm is active!

A copy **B** del **C** quit **D** display

It is possible to review a program during run. After pressing $\rightarrow$ C “program” and the corresponding store number a warning message appears.

program no: 4

name:

A copy **B**
C quit **D** display

Then it is possible to display the active program ($\rightarrow$ D “display”) or to copy it into a new store ($\rightarrow$ A “copy”). It’s not possible to change the currently running program
By pressing $\rightarrow$ C “quit” the main menu appears again.

Stop block function and electrophoresis

program: 0 TEST
pause?
stop?

A ? **B** pause **C** quit **D** stop

To stop a current running program you press $\rightarrow$ B “Elpho” and again $\rightarrow$ D “stop”.

By pressing $\rightarrow$ C “quit” you return to previous display without any changes.

When pressing $\rightarrow$ D “stop” the run will be aborted and you leave the actual running program.

L1: 25.0 °C L6: 60.0°C

hold: 1 pause 0.00 Vh
El: 250 V 8mA 20W
A ? **B** Elpho **C** programs **D** +

Pressing $\rightarrow$ B “pause” holds the actual situation of the gradient (stopping electrophoresis, holding the temperature gradient) “pause” is blinking and shown in the display alternatively with the time.

program: 0 TEST
continue?
stop?

A ? **B** contin **C** quit **D**
stop

Pressing $\rightarrow$ B “Elpho” (in pause status):

Pressing $\rightarrow$ B “contin” will continue the program.

Function key A "?"

T1: 22.0°C T2: 22.0°C
block off

A? **B** Elpho **C** programs **D** +

B: start/stop/pause

C: edit/delete/copy

D: special functions

A? **B** Elpho **C** programs **D**
+

Pressing $\rightarrow$ A "?" in the main menu results in the following display:

L1: 25.0 °C L6: 60.0°C
hold: 1 0m 1s 1.16Vh
El: 250 V 8mA 20W
A list **B** block **C** program **D** +

Pressing $\rightarrow$ A "?" in a running program results in the following display:

T1: 21.1 °C T2: 63.8°C
L1→3: 21.1 32.0 39.0
L4→6: 46.0 53.0 60.0
← = const rtime: h m

The actual temperatures of T1 and T2 as well as the actual temperatures of L1 to L6 are shown in the display.

The limiting factor (const. V, mA or W) is indicated by a blinking arrow (←).

The actual remaining electrophoresis time is shown on the bottom (right side of the display).

Error messages

TGGE - System check connection to thermoblock	TGGE connector cable is not connected to gradient block and / or system Controller. Check connections!!!!
warning: gradient too large! max. grad. T1→T2: 45°C A ? B no→L1 C quit D enter	Programmed temperature gradient too large.
or warning: gradient too large! max. grad. T1→T2: 45°C A ? B no→T1 C quit D enter	
program no: TEST pgm is active! A copy B del C quit D display	It is possible to review a program during run. After pressing  C “program” and the corresponding store number a warning message appears.
program no: __ name: not programmed! A ↑ B ↓ C quit D enter	This program number has not been programmed.
1:L1: __ L6: entry required T1: __ T2: A ? B L1 C quit D →	No temperature or time has been programmed.
or 1:L1: __ L6: entry required T1: __ T2: A ? B T1 C quit D →	

4 Sample preparation

4.1 Purity of samples

Due to the high sensitivity of the staining procedure after TGGE it is recommended to use purified DNA, RNA or protein samples. Any impurities might be misinterpreted after TGGE, thereby making the analysis of gels difficult. Nevertheless it is possible to use even crude mixtures for TGGE analysis.

PCR-amplified DNA fragments can usually be analyzed without purification. But note, that the presence of high amounts of nonspecific, secondary PCR products may result in difficulties with interpretation of band pattern, melting profile, etc. For example, in parallel TGGE, nonspecific bands with a higher molecular weight than the specific PCR product may be misinterpreted as heteroduplexes, or analogs with lower thermal stabilities. Therefore, before running a TGGE gel, please check the PCR product and, if necessary, purify the specific PCR product of interest, e.g., by agarose gel electrophoresis and subsequent gel extraction.

Sample preparation for direct DNA analysis

1 volume of DNA/RNA samples are dissolved with 1 volume of TBE or Na-TAE loading buffer or with 0.1 volume of the total loading volume ME loading buffer (see Appendix 2). The resulting mixture is loaded directly on to the polyacrylamide gels. Secure that the slots are filled up to maximum (if necessary add loading buffer to fill up the slots to maximum).

In case of low-concentration samples we recommend to prepare 5x conc. loading buffer. 0.2 volume of this concentrated loading buffer is mixed with the sample and loaded onto the polyacrylamide gel.

Denaturation/Renaturation for heteroduplex analysis of DNA

For heteroduplex analysis the samples are denatured and renatured prior to TGGE. Quantitative denaturation is accomplished by heating in 4 M urea. The following protocol is recommended for all DNA fragments with GC-contents of 50 - 70%. Depending on the buffer to be used for electrophoresis add one sample volume of corresponding DR buffer (denaturation/renaturation buffer) to the sample and mix. Heat at 95°C for 5 minutes (denaturation). Incubate at 50°C for 15 minutes (renaturation). The sample is then loaded directly to the gel. In order to achieve the recommended loading volumes for diagonal or perpendicular TGGE (refer to chapter 4.2), the samples should be filled up with running (or loading) buffer.

- **Renaturing at 50°C:**
 - higher temperatures (higher stringency) can be chosen for high GC contents to avoid artificial hybrids
 - lower temperatures (e.g. 37°C) are applied for expected hybrids with multiple mismatches or for sequences with very low GC content. The renaturation temperature should be approx. 10°C below the T_m of the desired hybrid.

4.2 Quantity and volumes of samples

Depending on the slot size of the gel Biometra recommends the following amounts of material:

glass plate with 1 rectangular slot + 2 marker slots	50 µl volume, approx. 50 ng DNA/RNA of interest 5 µl volume, 3 - 5 ng DNA/RNA of interest
glass plate with 8 slots:	5 µl volume, 3 - 5 ng DNA/RNA of interest
glass plate with 12 slots:	3 µl volume, 1 - 3 ng DNA/RNA of interest
glass plate with 18 slots:	2 µl volume, approx. 1 ng DNA/RNA of interest

If less volumes fill up the slots with running buffer or loading buffer to the volumes listed on top! This will create better results.

5 Setting up polyacrylamide gels

5.1 Selecting Concentration of PAA gels

The TGGE System represents a highly optimized system for performing flat bed polyacrylamide gel electrophoresis under defined temperature conditions. In addition to typical TGGE applications the system is ideally suited to run standard fragment separations without temperature gradient very quickly.

Depending on the molecular weight of the sample we recommend the following acrylamide/bisacrylamide concentrations:

Conc.	DNA fragment length
3%	> 1000 bp
5%	500 - 1000 bp
8%	< 500 pb

5.2 Setting up the gel solution

Each gel sandwich contains approx. 2.5 ml polyacrylamide solution.

We therefore recommend to prepare **10 ml solution to pour 3-4 gels** at the same time. Polymerized gels which are not immediately used must be stored at room temperature. To inhibit any gel drying we recommend to wrap the polymerized gels into saran foil or wet plastic bags. Wet towels can be used only for short time storage.

Keep in mind that polymerized polyacrylamide gels which include urea should not be used after 2 – 4 days of storage (depending on storage conditions)!

Recipe for 10 ml gel solution (3 – 4 gels) for TBE running buffer:

	3% Gel	5% Gel	8% Gel
Urea ($C_{\text{End}} = 7 \text{ M}$)	4.2 g	4.2g	4.2 g
Acrylamide/bis Acrylamide stock solution (30 : 0,8), 40% (w/v)	0.75 ml	1.25 ml	2.0 ml
10x conc.TBE ($C_{\text{End}} = 0.1 \times \text{conc.}$)	0.1 ml	0.1 ml	0.1 ml
50% Glycerol ($C_{\text{End}} = 2\%$)	0.5 ml	0.5 ml	0.5 ml
Water, distilled	3.5 ml	3 ml	2.5 ml
Make sure that the urea has been completely resolved. It is possible to heat up the urea containing solution slightly (40°C – 50°C) for a <u>short time</u> in order to improve the solubilization of urea. De-gas the solution under gentle vacuum for 3 - 5 min.			
Water, distilled	fill up to 10 ml		
TEMED	22.5 µl	22.5 µl	22.5 µl
APS (4%)	42 µl	42 µl	42 µl
Mix gently. Avoid air bubbles!			
Pour the gel solution into the glass plate sandwich immediately thereafter (see chapter 4.1.2) without air bubbles.			

Recipe for 10 ml gel solution (3 – 4 gels) for Na-TAE running buffer:

	3% Gel	5% Gel	8% Gel
Urea ($C_{\text{End}} = 8 \text{ M}$)	4.8 g	4.8g	4.8 g
Acrylamide/bis Acrylamide stock solution (30 : 0,8), 40% (w/v)	0.75 ml	1.25 ml	2.0 ml
10x conc. Na-TAE, pH 8.4 ($C_{\text{End}} = 0.2 \times \text{conc.}$)	0.2 ml	0.2 ml	0.2 ml
40% Glycerol ($C_{\text{End}} = 2\%$)	0.5 ml	0.5 ml	0.5 ml
Water, distilled	3.5 ml	3.0 ml	2.5 ml
Make sure that the urea has been completely resolved. It is possible to heat up the urea containing solution slightly (40°C – 50°C) for a <u>short time</u> in order to improve the solubilization of urea. De-gas the solution under gentle vacuum for 3 - 5 min.			
Water, distilled	fill up to 10 ml		
TEMED	14 µl	14 µl	14 µl
APS (4%)	45 µl	45 µl	45 µl
Mix gently. Avoid air bubbles!			
Pour the gel solution into the glass plate sandwich immediately thereafter (see chapter 4.1.2) without air bubbles.			

Recipe for 10 ml gel solution (3 – 4 gels) for ME (MOPS/EDTA) running buffer:

	3% Gel	5% Gel	8% Gel
Urea ($C_{\text{End}} = 8 \text{ M}$)	4.8 g	4.8g	4.8 g
Acrylamide/bis Acrylamide stock solution (30 : 0,8), 40% (w/v)	0.75 ml	1.25 ml	2.0 ml
50x conc. ME-buffer ($C_{\text{End}} = 1 \times \text{conc.}$)	0.2 ml	0.2 ml	0.2 ml
40% Glycerol ($C_{\text{End}} = 2\%$)	0.5 ml	0.5 ml	0.5 ml
Water, distilled	3.5 ml	3.0 ml	2.5 ml
Make sure that the urea has been completely resolved. It is possible to heat up the urea containing solution slightly (40°C – 50°C) for a short time in order to improve the solubilization of urea. De-gas the solution under gentle vacuum for 3 - 5 min.			
Water, distilled	fill up to 10 ml		
TEMED	17 μl	17 μl	17 μl
APS (4%)	76 μl	76 μl	76 μl
Mix gently. Avoid air bubbles!			
Pour the gel solution into the glass plate sandwich immediately thereafter (see chapter 4.1.2) without air bubbles.			

5.3 Some remarks corresponding to standard TGGE conditions

Electrophoresis buffer (running buffer):

- Always membrane filtrate (e.g. 0.45 μm pore size) the buffers before use!
- Running buffer: always use the concentration identical with the gel condition
- TBE is the most common used buffer system but the electrophoresis is not as fast as with Na-TAE buffer. It is possible to add up to 5 mM NaCl if a higher ionic strength is desired, for reversible melting processes which are required for parallel TGGE in multiple sample analysis. A higher NaCl concentration should not be used because it causes an unacceptable high electrical current.
- Na-TAE is the buffer for fastest electrophoresis.
- ME buffer meets all the requirements of a variety of TBE-buffers with different ionic strengths but is only stable for a very short time. (Stable for about 3 days. Do not use as the buffer becomes yellow.)
- ME buffer allows Na^+ concentrations up to 20 mM which greatly favors "reversible melting" and still allows short run times for TGGE electrophoresis. Mobile Cl^- ions which slow down the migration velocity of nucleic acids are avoided by using the sodium salt form of MOPS. Due to their reduced mobility the large MOPS anions keep the current low.

Gel conditions:

- 4 M urea can be used for low GC and high degree of mismatches. This concentration increases T_m instead of the standard concentration (8 M urea) by approx. 16 - 20°C using TBE buffer or 8 - 12°C using Na-TAE buffer.
- 10 M urea can be used for high GC and lowers T_m instead of the standard concentration (8 M urea) by approx. 8 - 10°C using TBE buffer or 4 - 6°C using Na-TAE buffer.
- Glycerol reduces the steepness of very cooperative transition curves, broadening the profile and expanding the temperature range for detecting small T_m differences of closely related nucleic acids:
 - 0% Glycerol increases cooperativity (>200 bp, narrower transitions)
 - >2% Glycerol lowers cooperativity (>200 bp, broader transitions)

TGGE conditions:

- Voltage can be raised to 400 V if the current is below 30 mA
- Current should not exceed 30 mA
- Running buffer:

	pre-run	run
0.2 x conc. Na-TAE	250 V, 10-12 mA, 2-5 min.	250V, 15-20 mA, 30-60 min.
0.1 x conc. TBE	250 V, 4- 5 mA, 2-5 min.	250 V, 9-12 mA, 30-60 min.
- $T_1 = 20^\circ\text{C}$: can be raised to obtain optimized resolution
- $T_1 = 20^\circ\text{C}$: lowering should be avoided by using 4M urea
- $T_2 = 60^\circ\text{C}$: higher temperatures should be avoided by using 10 M urea or/and 10 mM MOPS.
- $T_2 = 60^\circ\text{C}$: can be lowered to obtain optimized resolution
- $T_2 = 80^\circ\text{C}$: maximum temperature, can be used if the gel is carefully protected against evaporation.

Sample preparation:

- Denaturation / renaturation cycle:
 - Renaturing at 50°C : higher temperatures (higher stringency) can be chosen for high GC contents to avoid artificial hybrids
 - Renaturing at 50°C : lower temperatures (e.g. 37°C) are applied for expected hybrids with multiple mismatches or for sequences with very low GC content. The renaturation temperature should be approx. 10°C below the T_m of the desired hybrid.
- Samples should be dissolved in buffers with ionic strength identical to the ionic strength of the running buffer.
- If samples are dissolved in buffer different to the running buffer, the samples have to be equilibrated against the running buffer (e.g. using dialysis). Nucleic acids can be precipitated with ethanol and dissolved in denaturation/renaturation buffer or running buffer.

5.4 Assembling the gel sandwich

The thinness of the gel makes it necessary to cast polyacrylamide gels on a gel support film (Polybond, Order Number: 024-030). Each sandwich consists of four elements:

- Bonding glass plate without spacer
- Polybond film
- Polyacrylamide gel
- Glass plate with fixed spacer and fixed slot former (different types of slots available).

Glass plates

- Glass plates must be **dry and free of any dirt or dust**. Biometra recommends to wear **powder-free gloves** even during cleaning of glass plates in order to prevent any skin debris which might interfere with silver staining.
- **Do not** use strong acidic or basic solutions or organic solvents for cleaning the glass plates.
- **Do not** incubate glass plates over night in cleaning solutions.

Pretreatment of glass plate with spacer and slot former

- Glass plate with spacer and slot former must be carefully treated with **Acryl-Glide solution** (Order Number 211-319) or a similar hydrophobic solution. Drop about 0.5 ml of solution onto the plate and **especially between the slot former**. (This protection layer helps to withdraw the polyacrylamide gel from the sandwich after polymerization.) **Wait 2 – 3 min. and then polish the plate with soft paper** to remove any haze!
- This procedure should be repeated after each run.
- **Do not** drop Acryl-Glide onto the spacer of the glass plate! This possibly leads to leakage during polymerization.
- **Clean spacers with ethanol** before assembling the glass plate sandwich.
- If necessary treat the spacers with a small amounts of silicone grease (to protect leakage).

Polybond film

- Use only original pre-cut Polybond film which perfectly fits onto the gradient block.
- Biometra recommends **to use Polybond film only once**. Repeated usage and especially staining might weaken the strength of the Polybond film.
- The Polybond film has **two different sides**: one hydrophobic side which repels water and one hydrophilic side on which a water drop will adhere. You can test the different sides with a drop of water. (The protection paper is attached to the hydrophobic site.)
- Remove the protecting paper sheet before assembling the sandwich. Handle the Polybond film only with **powder-free gloves**.
- The **hydrophobic side** of the Polybond film must be **orientated to the Bonding glass plate**. On the hydrophilic side the polyacrylamide gel will polymerize and stick.
- Press the hydrophobic side of the Polybond film firmly to the Bonding glass plate by using your thumb, a rubber or a gel casting clip. This will prevent that any polyacrylamide solution running between the Bonding glass plate and the Polybond film.

Glass plate sandwich

- Assemble the Acryl-Glide treated glass plate with spacer, the Polybond film and the Bonding glass plate as indicated in figure 11.




Figure 11: Setting up a gel sandwich for PAGE.

Bonding glass plate (①), Polybond film (②) and glass plate with spacer and slot former (③) are assembled (left panel). The Polybond film is only visible at the inclined edge of the glass plate with spacer (middle panel). Fasten the clamps (④) above the spacer to increase the pressure and to ensure a leakage free sandwich.




Figure 12: Pouring the polyacrylamide gel.

Initially 1 ml of polyacrylamide gel solution must be poured. The gel sandwich must be held at an **angle of 45° when pouring** (left panel). The solution must run along one side of the plate sandwich to **avoid air bubbles**. The remaining 1.5 ml of solution is poured into the plate sandwich during which time the plate sandwich is slowly brought back into a vertical position (right panel).

- Fill up the plate sandwich as high as possible. The gel solution is overlayed with 200µl of **Isopropyl Alcohol** (2-Propanol, Isopropanol) or **Isobutyl Alcohol** (2-Methyl-1-propanol, Isobutanol) or **distilled water** to produce a horizontal surface of the gel.

- Polymerization of the gel must be at least for 0.5 h, better for 1 - 1.5 h or optional over night at room temperature. The sandwich should stand up vertically and must not be moved during polymerization.
- Gels may be **stored up to 4 days at room temperature** (wrapped in wet paper towels in a plastic bag). **Do not store at 4°C!**

5.5 Disassembling the gel sandwich

- Remove the clamps from the plate sandwich.
- **Remove the Bonding glass plate from the sandwich by sliding it smoothly!** The gel polymerized to the Polybond film will adhere to the other glass plate.
- Withdraw the Polybond film with the adhering polyacrylamide gel carefully from the other glass plate. In **the area of the slot former** remove the Polybond film **very carefully** to avoid any damage to the slots.
- If slots show distortion or wrinkles don't fill in samples because after electrophoresis bands in this lane will show distortion as well.

6 Electrophoresis with the TGGE System

The electrophoresis unit of the TGGE System has been designed to accommodate all TGGE and related applications like CTGE, TTGE and SSCP, without cumbersome changes. It's easy to switch between perpendicular, parallel or diagonal TGGE (for adequate accessories see chapter 6.4).

6.1 Electrophoresis conditions

The electrophoresis conditions depend on the
 kind of material to be separated, e.g. fragment size differences,
 kind of application, e.g. parallel or perpendicular TGGE,
 sample preparation, e.g. high salt or low salt preparation,
 buffer system.

Any recommendations can only be used as guidelines to start with. Further improvements to the analysis is easily possible by adjusting the run conditions to the individual needs.

Voltage:	100 V - 400 V;	usually 250 V
Current:	5 mA - 25 mA;	usually 10-20 mA
Run Time:	10 min - 2 h;	usually 30 min

6.2 Preparing the electrophoresis unit

- Use the leveling eye on the electrophoresis unit and the 4 leveling feet to adjust the unit.
- Remove the safety lid and fill in max. **250 ml of the desired running buffer per buffer chamber** (e.g. 0.1 x conc. TBE, see Appendix 13.2.). The same running buffer should only be used once.
- Soak the pre-cut electrode wicks (order number: 024-020) in the running buffer before use.
- **Drop 0.3 - 0.5 ml of thermal coupling solution like 0.1% Triton or 0.1% Tween 20 on the surface** of the gradient block (see figure 13 *left panel*). The thermal coupling solution will increase the adhesion of the Polybond film with the attached polyacrylamide gel and therefore supports temperature equilibration between gradient block and polyacrylamide gel. The whole block must be covered by the thermal coupling solution layer. No air bubbles must form.




Figure 13: Positioning of polyacrylamide gel on gradient block. A small volume of thermal coupling solution (2) is applied to the gradient block (1) (left panel). The Polybond film with the polyacrylamide gel on top (3) is put on the gradient block. Slightly bend the Polybond film (right panel) in order to spread the thermal coupling solution evenly.

- To position the polyacrylamide gel on the gradient block the Polybond film should be held between thumb and middle finger and slightly bended. This leads to an even distribution of thermal coupling solution beneath the Polybond film.
- If **air bubbles** are visible beneath the polyacrylamide gel, try to **squeeze them out** by moving the gel slightly back and forth. If this will not succeed completely remove the gel and repeat the aforementioned steps. Don't touch the polyacrylamide gel directly with your fingers or your gloves.
- **Excess thermal coupling solution** must be removed from the gradient block by using paper towels.

6.3 Perpendicular TGGE

During perpendicular TGGE a mixture of molecules is separated over a **wide temperature range**. The temperature gradient is perpendicular to the electrophoresis run direction.

Steps before TGGE

- sample preparation (refer to chapter 4.1)
- programming of temperature gradient and electrophoresis parameters (refer to chapter 3.2.2)

Prepare in advance:

- polyacrylamide gel attached to Polybond film
- running buffer (250 ml for each chamber)
- pre-cut and pre-soaked electrode wicks
- pre-cut cover film
- cover glass plate (treated on both sides with Acryl-Glide)

- Because of the fixed orientation of the temperature gradient the removable electrophoresis chambers must be positioned as indicated in figure 14.




Figure 14: Orientation of the two electrophoresis buffer chambers for perpendicular TGGE.

- Fill in the running buffer (e.g. 0.1 x conc. TBE) into each electrophoresis chamber.



Before you place the gel onto the gradient block be sure that the sample is ready for loading and cover film is available.

- To make full use of the linear range of the gradient block the polyacrylamide gel attached to the Polybond film should be positioned as indicated in figure 15. The rectangular and the marker slots of the gel (⊙) are positioned at the beginning of the gradient block (the marked lines represent the beginning of the linear range of the gradient block).




Figure 15: Orientation of the polyacrylamide gel attached to the Polybond film (1) on the gradient block (2) during perpendicular TGGE. See the position of the gel slots (3) relative to the marked lines of the gradient block.

- Two pre-cut and pre-soaked **electrode wicks** must be **positioned at the start and the end** of the polyacrylamide gel. Wick and gel have to overlap (see figure 16).




Figure 16: Side view of the polyacrylamide gel (2) on top of the gradient block during pre-run. Pay attention to the position of electrode wicks (4) on top of the polyacrylamide gel. (Polybond film (1), slot of polyacrylamide gel (3)).

- Avoid any contact between sample slots and electrode wicks.** Otherwise the samples will diffuse into the electrode wicks.
- Load the samples quickly** at room temperature without air bubbles. **Do not start the temperature gradient** (the temperature gradient is established after the samples have fully entered the polyacrylamide gel).

glass plate with 1 rectangular slot

+ 2 marker slots

50 μ l volume,
approx. 50 ng DNA/RNA of interest
5 μ l volume,
3 – 5 ng DNA/RNA of interest




If less volumes fill up the slots with running buffer or loading buffer to the volumes listed on top! This will create better results.



The time between mounting the gel onto the gradient block and loading the sample must not exceed 5 minutes.

- Close the safety lid of the electrophoresis chamber and start electrophoresis at 20°C or 25°C and 250 V for 2 - 5 min. Standard electrophoresis conditions are given in chapters 5.3 and 6.1.
- Make sure that the orientation of the gel and the safety lid is exact as indicated in the following:



- Wait at the electrophoresis chamber until the samples have fully entered the polyacrylamide gel (unlike with the former TGGE System of QIAGEN this process will only take 1-3 minutes) and have moved about 3 - 5 mm in the gel.
 - Stop the electrophoresis run, open the safety lid.
 - **Rinse** the - now empty - slots with 0.5 - 1 ml running buffer.
 - **Cover the polyacrylamide gel** including the slots with the 7 x 6 cm pre-cut cover film (see figure 17). A small buffer layer must remain between cover foil and gel.
- Avoid air bubbles!**



The cover film must be positioned with the long side parallel to the buffer chambers (= perpendicular to the arrow on the safety lid).

- **Soak** any excessive buffer from the side of the gel. The gel must not swim in buffer solution.




Figure 17: The polyacrylamide gel (2) has to be covered by a pre-cut hydrophobic cover film (4). A small buffer layer remains between gel and cover film. (Polybond film (1), wicks (4), slot of polyacrylamide gel (5))

- Bring the **electrode wicks** to an **overlap with the cover film**. The overlap between wick and cover film should be almost 2 cm (see figure 17). Avoid air bubbles.
- Be sure that the 2 **silicone barriers** are fixed to the cover glass plate before use:



- **Cover the sandwich with the Acryl-Glide treated cover glass plate** (see figure 18).



The silicone barriers have to be positioned perpendicular to the wicks and never on top of the wicks!




Figure 18: Side view of the polyacrylamide gel (2) on top of the gradient block. Sandwich of polyacrylamide gel (2), pre-cut cover foil (3), electrode wicks (4) and cover glass plate with silicone barriers (6) during perpendicular TGGE run (Polybond film (1), slot of polyacrylamide gel (5)).

- **Start the temperature gradient** and wait until the gradient has been established (usually 0.5 - 1 minute).
- **Start the electrophoresis run.**
- The Bromophenol blue dye only gives you an indication how far the samples have migrated until you have optimized the best run time.
- After the electrophoresis run, switch off the controller, open the safety lid of the electrophoresis unit, remove the polyacrylamide gel and proceed further for staining the gel (chapter 6.5). It is recommended to fix the gel immediately in order to improve the analysis.

6.4 Parallel TGGE

During parallel TGGE a mixture of molecules is separated over a **narrow temperature range** determined by perpendicular TGGE. The temperature gradient is parallel to the electrophoresis run direction.

Steps before TGGE

- sample preparation (refer to chapter 4.1)
- programming of temperature gradient and electrophoresis parameters (refer to chapter 3.2.2)

To prepare in advance:

- polyacrylamide gel attached to Polybond film
- running buffer (250 ml for each chamber)
- pre-cut and pre-soaked electrode wicks
- pre-cut cover foil
- glass plate (treated on both sides with Acryl-Glide)

- For parallel TGGE the removable electrophoresis buffer chambers must be positioned as indicated in figure 19.




Figure 19: Orientation of the electrophoresis chambers for parallel TGGE.

- Fill in the running buffer (e.g. 0.1 x conc. TBE) into each electrophoresis chamber.



Before you place the gel onto the gradient block be sure that the sample is ready for loading and cover film is available.

- The polyacrylamide gel attached to the Polybond film must be positioned as indicated in figure 20. The slots of the gel (❶) should be positioned at the beginning of the gradient block. The first marked line (L1) represents the beginning of the linear range of the gradient block.




Figure 20: Orientation of the polyacrylamide gel attached to the Polybond film (1) on the gradient block (2) during parallel TGGE. See the position of the gel slots (3) relative to the marked lines of the gradient block.

- Two pre-cut and pre-soaked **electrode wicks** must be **positioned at the start and the end** of the polyacrylamide gel. Wick and gel have to overlap (see figure 16).




Figure 16: Side view of the polyacrylamide gel (2) on top of the gradient block during pre-run. Pay attention to the position of electrode wicks (4) on top of the polyacrylamide gel. (Polybond film (1), pre-cut Polybond film (3), slot of polyacrylamide gel (5)).

- **Avoid any contact between sample slots and electrode wicks.** Otherwise the samples will diffuse into the electrode wicks.
- **Load the samples quickly** at room temperature without air bubbles. **Do not start the temperature gradient** (the temperature gradient is established after the samples have fully entered the polyacrylamide gel).

Depending on the slot size of the used gel Biometra recommends the following amounts of material:

glass plate with 8 slots:	5 μ l volume, 3 - 5 ng DNA/RNA of interest
glass plate with 12 slots:	3 μ l volume, 1 - 3 ng DNA/RNA of interest
glass plate with 18 slots:	2 μ l volume, approx. 1 ng DNA/RNA of interest




If less volumes fill up the slots with running buffer or loading buffer to the volumes listed on top! This will create better results.



The time between mounting the gel onto the gradient block and loading the sample must not exceed 5 minutes.

- Close the safety lid of the electrophoresis chamber and start electrophoresis at 20°C or 25°C and 250 V for 2 - 5 min. Standard electrophoresis conditions are given in chapters 5.3 and 6.1.
- Make sure that the orientation of the gel and the safety lid is exact as indicated in the following:



- Wait at the electrophoresis chamber until the samples have fully entered the polyacrylamide gel (unlike with the former TGGE System of QIAGEN this process will only take 1-3 minutes) and have moved about 3 - 5 mm in the gel.
- Stop the electrophoresis run, open the safety lid.
- **Rinse** the - now empty - **slots** with 0.5 - 1 ml running buffer.
- **Cover the polyacrylamide gel** including the slots with the 7 x 6 cm pre-cut cover film (see figure 17). A small buffer layer must remain between cover foil and gel.
Avoid air bubbles!



The cover film must be positioned with the long side parallel to the buffer chambers (= perpendicular to the arrow on the safety lid).

- **Soak** any excessive buffer from the side of the gel. The gel must not swim in buffer solution.




Figure 17: The polyacrylamide gel (2) has to be covered by a pre-cut hydrophobic cover film (3). A small buffer layer remains between gel and cover film. (Polybond film (1), wicks (4), slot of polyacrylamide gel (5))

- Bring the **electrode wicks to an overlap with the cover film**. The overlap between wick and cover film should be almost 2 cm (see figure 17). Avoid air bubbles.
- Be sure that the 2 **silicone barriers** are fixed to the cover glass plate before use:



- **Cover the sandwich with the Acryl-Glide treated cover glass plate** (see figure 18).



The silicone barriers have to be positioned perpendicular to the wicks and never on top of the wicks!




Figure 18: Side view of the polyacrylamide gel (2) on top of the gradient block. Sandwich of polyacrylamide gel (2), pre-cut cover foil (3), electrode wicks (4) and cover glass plate with silicone barriers (6) during perpendicular TGGE run (Polybond film (1), slot of polyacrylamide gel (5)).

- **Start the temperature gradient** and wait until the gradient has been established (usually 0.5 - 1 minute).
- **Start the electrophoresis run.**
- The Bromophenol blue dye only gives you an indication how far the samples have migrated until you have optimized the best run time.
- After the electrophoresis run, switch off the controller, open the safety lid of the electrophoresis unit, remove the polyacrylamide gel and proceed further for staining the gel (chapter 6.5). It is recommended to fix the gel immediately in order to improve the analysis.

6.5 Silver staining

Aside from autoradiography silver staining is the most sensitive method for detecting small amounts of DNA, RNA or proteins in polyacrylamide gels. Due to the low thickness of the gels (0.5 mm) the staining procedure takes no more than 35 minutes.

Other staining protocols may be used, but generally exhibit less sensitivity. This must be considered in relation to the amount of DNA loaded on the gel.

All incubation steps are done in small plastic containers which are agitated on a rocking platform (e.g. order number 042-400 or 042-500).

For handling several polyacrylamide gels simultaneously, Biometra offers a semi-automated instrument called Blot Processor (order number 015-000 or 015-090). Please contact Biometra or your local distributor to receive further information about the Blot Processor.



Wear non-powdered protective gloves during all steps of the silver staining protocol to avoid staining artifacts due to the high sensitivity of the staining protocol.

- Remove the protective plastic sheets from the gel.
- Carefully remove any residual thermal coupling solution from the back of the gel (Polybond film) prior to staining
- Put the polyacrylamide gel with the gel side upwards into the staining tray. Avoid air bubbles during all staining steps.
- It's recommended to prepare at least 100 ml solution for each incubation step.
- If NaCl has been added to the gel running buffer, incubate the TGGE gel for 15 min in Fixation solution to remove the NaCl.

Standard method:

Step	Time	Solutions*	Notes
Fixation	5 min	Fixation solution	
Silver Binding	10 min	AgNO ₃ -Solution	prepare freshly
Washing	3 x 1 min	Fresh ddH ₂ O	demineralised water may be ok
Developing	10 min	Developer	prepare freshly
Stopping	5 min	Stopping Solution	
Washing	10 min	Rinse under fresh ddH ₂ O	Demineralised water may be ok
Preparing for storage	1-5 h	50% glycerol	Not absolutely necessary!

Preparing the gel for storage (1 - 5 h at room temperature in 50% glycerol) is not necessary.

Staining solutions:

Fixation	10% EtOH 0.5% Acetic Acid 100 ml ethanol and 5 ml acetic acid are adjusted with distilled water to 1 liter.
Silver Binding	0.19% AgNO ₃ 1.9g AgNO ₃ is dissolved in 1 liter of distilled water. (Can be reused 5 times) Store dark!
Developing Solution	1.5% NaOH 0.08% NaBH ₄ 0.1% Formaldehyde Dissolve 15 g NaOH in 1 liter distilled water. Add 0.8g NaBH ₄ and 2.7 ml formaldehyde stock solution (37% in water). This buffer must be freshly prepared immediately before use!
Stopping Solution	0.75% Na ₂ CO ₃ Dissolve 7.5 g sodium carbonate in ddH ₂ O. Total volume: 1 liter

Quick method (for PCR products) (Sanguinetti et al):

Step	Time	Solutions*	Notes
Fixation	3 min	Fixation solution	prepare freshly
Silver Binding	5 min	AgNO ₃ -Solution	prepare freshly
Washing	3 x 1 min	Fresh ddH ₂ O	Demineralised water may be ok
Developing	5 min	Developing solution	prepare freshly
Stopping	5 min	Ethanol and acetic acid solution	
Washing	10 min	Rinse under fresh ddH ₂ O	Demineralised water may be ok
Preparing for storage	1-5 h	50% glycerol	Not absolutely necessary
Drying		Room temperature	

Preparing the gel for storage (1 - 5 h at room temperature in 50% glycerol) is not necessary.

Staining solutions:

Fixation:	10% EtOH 0.5% Glacial Acid 100 ml ethanol and 5 ml acetic acid are adjusted with double distilled water to 1 liter. Prepare freshly !
Silver Binding	0.2% AgNO ₃ 2.0 g AgNO ₃ is dissolved in 1 liter of distilled water. (Can be reused 5 times) Store dark!
Developing Solution	3.0 % NaOH 0.5% Formaldehyde Dissolve 3 g NaOH and 1.35 ml formaldehyde stock solution (37% in water) in 100 ml double distilled water. This buffer must be freshly prepared immediately before use!
Stopping Solution:	identical with Fixation solution (10% EtOH, 0.5% Glacial Acid)

**Quick method using the AMRESCO SilverPAGE™ staining kit
(Code No. 211-761)**

Step	Time	Solutions	Notes
Fixation	15 min	2 x 100 ml Fixation solution	
Sensibilisation	10 min	2 x 100 ml 30% ethanol	prepare freshly!
Washing	10 min	3 x 200 ml fresh ddH ₂ O	Demineralised water may be ok
Silver Binding	15 min	reconstituted Silver Binding Agent + Formaldehyde	prepare freshly!
Washing	0.5-1 min	<i>Rinse under fresh ddH₂O</i>	<i>!!!!</i>
Developing	1 - 2 min	reconstituted Developing solution + Formaldehyde	prepare freshly! Develop to desired level!
Stopping	5 min	7.5% Acetic acid	
Preparing for storage	1-5 h	50% glycerol	Not absolutely necessary
Drying		Room temperature	

Preparing the gel for storage (1 - 5 h at room temperature in 50% glycerol) is not necessary.

Staining solutions:

- Fixation: 30% EtOH
10% Acetic Acid
300 ml ethanol and 100 ml acetic acid are adjusted with double distilled water to 1 liter.
- Sensibilisation: 30% EtOH
Prepare freshly 60 ml ethanol in 140 ml double distilled water.
- Silver Binding: Prepare Silver Binding Agent by reconstituting contents of one pouch in 1 l of ddH₂O. **(This solution must be prepared fresh every time!)** Immediately before staining, add 0.7 ml of 37% Formaldehyde to 200 ml of reconstituted Silver Binding Agent.
- Developing Solution: Just prior to use, prepare developing solution by reconstituting contents of one pouch of Developer I and 15 mg of Developer II in 200 ml of ddH₂O. **(This solution must be prepared fresh every time!)** Immediately before developing, add 0.7 ml of 37% Formaldehyde to 200 ml of reconstituted developing solution.
- Stopping Solution: 7.5% Acetic Acid
75 ml acetic acid are adjusted with double distilled water to 1 l.

6.6 Ethidium bromide-staining

Incubate the gel in staining solution (0.5 µg/ml ethidium bromide in 1 x conc. TBE) for 30 - 45 min. Analyze under UV radiation (27).

6.7 Blotting

DNA from TGGE gels can be **blotted** onto a solid-state support either by electroblotting (Fastblot) or vacuum-blotting (Vacu-Blot System). If DNA is to be blotted after TGGE analysis, the TGGE gel must be poured onto the hydrophobic side of the gel support film (Polybond film). Otherwise, the gel cannot be detached from the gel support film!

6.8 Autoradiography

TGGE gels can also be directly exposed to x-ray films if radiolabeled samples are analyzed.

Direct exposure:

Incubate the TGGE gel for 15 min. in Fixation solution (see 6.5 Silver staining). Optional: Silver stain the gel.

Remove residual buffer from the gel. Expose to an x-ray film at room temperature.

Exposure of dried TGGE gels:

Incubate the TGGE gel for 15 min. in Fixation solution (see 6.5 Silver staining). Optional: Silver stain the gel.

Incubate the gel in 2 - 5% glycerol for 10 minutes to prevent the gel from cracking. Incubate an appropriate sheet of cellophane (no Saran wrap!!!!) in 2 - 5% glycerol. Layer the cellophane on the gel. Air dry at room temperature for one day or use a geldryer at 50°C for at least 3h. Exposure to an X-ray film.

6.9 Elution of DNA from the TGGE gel

DNA fragments which have been separated on TGGE, for example, different alleles of one gene, can be eluted from silver-stained TGGE gel and reamplified by PCR.

Using a Pasteur pipet, puncture the gel and extract a µl piece containing the particular DNA duplex. Incubate in 20 µl TE buffer overnight. Use a 1 µl aliquot for reamplification.

7 TGGE in analysis of point mutations in dsDNA

For analysis of point mutations in dsDNA, an extremely high detection rate of greater than 95% is routinely achieved when the experiment is carefully planned. The next two chapters provide information for optimizing detection of base substitutions.

7.1 Theoretical background of a detection rate approximating 100% for point mutations - calculations with the POLAND program

DNA does not melt by deannealing base pair by base pair from one end to the other, but by cooperative denaturation of long stretches, called melting domains. The length of a melting domain is 25 to several hundred base pairs. The midpoint melting temperature T_m and the length of a melting domain are mainly determined by the nucleotide sequence of the DNA. The T_m of DNA fragments differing by even small changes, such as point mutations, can differ by as much as 1.5°C. When heteroduplexes, hybrids of two species of DNA fragments differing in their base composition, have been formed, the mismatches lower the T_m value significantly. Thus the heteroduplex analysis is the preferable because of the additional resolution provides (1, 28).

The principle by which TGGE uses differences in T_m is that the DNA fragments are electrophoresed through a linear temperature gradient in the polyacrylamide gel. When the fragments reach the temperature at which the lowest melting domain starts to melt, they take on a branched, Y-shaped configuration, which slows down mobility in the TGGE gel matrix. The electrophoretic migration of fragments differing by single base changes is retarded by branching at different temperatures, thus they are resolved from one another during temperature gradient electrophoresis.

The denaturing behavior of any DNA fragment can be predicted, if its sequence is known. For this purpose the POLAND calculation can be used. The POLAND software is available in the internet:

<http://www.biophys.uni-duesseldorf.de/service/polandform.html>.

POLAND software predicts location and T_m values of melting domains for dsDNA and dsRNA as well as their perpendicular TGGE pattern. The ability to predict the melting behavior of particular DNA fragments enables one to construct DNA fragments with optimized melting behavior, resulting in a nearly 100% detection rate for point mutations inside of this fragment.

Since the end of February 1999 the POLAND program is available in two versions. Using the above internet address allows the user to select between the old POLAND request form (Standard), the old POLAND expert request form or the new POLAND request form. The following information is visible on the screen:

Poland Server: ANNOUNCEMENT

The WWW server for prediction of nucleic acid's thermal stability (called Poland server according to the author of the basic mathematics) will move during the near future to another computer. This is not a mere relocation of the program that you have used up to now, but the input/output procedure is completely rewritten for the new server. The new server produces better/nicer (?) plots and has a better/more elongated help file. But be aware of new bugs, which might be introduced during the rewriting and relocation.

The old server, both the standard and the expert form, are unchanged. Both forms will be available for the near future. However, that server is running on our DEC Alpha under OpenVMS, and we run into more and more trouble to support that machine.

Now, make a decision:

NEW Server	OLD Server
Poland request form	Poland request form
	Poland expert request form

Institut für Physikalische Biologie
(Department of Biophysics)
Heinrich Heine-Universität Düsseldorf, Germany
Feb. 26, 1999

[G. Steger](#) / M. Labensky / A. Jäger

7.2 The “old” Poland program (Old Server)

7.2.1 About the Poland service

The Poland program is an **experimental** service of the University of Düsseldorf Biophysics Department, and thus the whole set-up, access and service are subject to change.

The Poland server will calculate thermal denaturation profiles and temperature-dependent uv absorbance or gel mobility of double stranded RNA or DNA, based on sequence input and parameter settings in the [request](#) form. - Details [below](#) !

The program used in these calculations was developed by Gerhard [Steger](#), for comparing theoretical predictions to experimental data, mainly optical denaturation profiles, taken at 260 and 280 nm, and [TGGE](#) (temperature gradient gel electrophoresis) experiments. The original version was written in VAX Fortran (VMS), using the Graphics Kernel System GKS for data presentation.

7.2.2 Program-specific information

Calculation is based on D. [Poland](#)'s algorithm in the implementation described by Gerhard [Steger](#).

The Poland algorithm calculates the denaturation profile for double-stranded nucleic acid using nearest-neighbor stacking interactions and loop entropy functions described in the literature. An extension of the algorithm, the 'virtual stack' model, allows for the incorporation of specific mismatched sequence positions in the stability calculations, as described by Heinz [Werntges](#).

The input data required for calculation are:

- **the sequence ($\leq 1,000$ bp).** Use GCG-format or plain format without spacing. Plain format accepts only 180 characters per line!
- **mismatched positions (optional),**
- the strand concentration, affecting the dissociation temperature (use the programmed standard),
- parameter set selection (DNA/DNA low salt/RNA; oligo/long ds),
- **output format options (choose GIF format).**

Data sets predicted by the program comprise the following:

- A perspective view on the temperature-dependent denaturation profile, that is denaturation probability vs. sequence position vs. temperature.
- The temperature-dependent relative uv hypochromicities as they would be measured in optical melting, at wavelengths of 260 and 280 nm (282 nm in case of DNA), respectively (full hypochromaticity corresponds to approx. 30% of the OD at low temperature).
- The derivative form of above hypochromicities, showing the melting temperature(s) and corresponding half width(s) of the transitions(s), giving hints about the transition cooperativity.

- Predicted relative gel mobility as calculated according to Lerman *et al.*, as a graph vs. temperature for different values of the 'retardation length' parameter. This plot can be used for direct comparison with TGGE experiments; superpositions of plots generated with or without mismatched positions given are useful as a hint whether specific mismatched duplexes could be detected among homoduplexed DNA in a mixture of sample and reference double strands having undergone a denaturation-renaturation cycle, using either perpendicular or parallel TGGE.
- A 'half-denaturation temperature' plot showing the half-denaturation temperature for each base. This plot can also be used to estimate the destabilizing effect of mismatches on the surrounding part of the sequence: a temperature-shift of the TGGE transition can be expected if the lowest melting part of the sequence is directly affected by the mismatch!

Calculations can be done for oligonucleotides (>15 bases) or long double strands (>50 bases), respectively. In the case of oligonucleotide mode, a length-dependent correction for the strand dissociation process is applied; the temperature range is adapted as well. We do not have sufficient experimental results to stringently check for this mode to give valid results, but for the length range of about 20 nucleotides there is at least experimental evidence. Using 'oligo' mode with far longer sequences gives misleading results!

Graphics output is possible in Postscript, HPGL, GIF and PBM format, numeric results are available as well. All graphics results are directly sent to the WWW client, for GIF inline images and pbm images links are provided to retrieve a copy for the external viewer (or for download to disk).

7.2.3 How to use the “old” Poland program (Standard)

Poland service request form

This form is an **experimental** service of the Biophysics Department, further informations are available [here](#)!

The Poland server will calculate the thermal denaturation profile of double stranded RNA or DNA based on sequence input and parameter settings in this form.

Calculation is based on D. [Poland](#)'s algorithm in the implementation described by [G. Steger](#).

Calculations can be done for oligonucleotides (>15 bases) or long double strands (>50 bases), respectively. A form allowing for '[expert](#)' parameter settings is available, too. Graphics output is possible in Postscript, HPGL, GIF and PBM format, numeric results are available as well. Graphics results are directly sent to your WWW client.

For a (more or less) detailed description of the various parameters, you may read a [help page](#).

Sequence title line:

Sequence:

(plain format without spacing; max. 180 chars per line)

[Mismatched positions:](#)

(comma-separated)

[Strand concentration:](#)

(default 1.0e-6 M; don't give the unit)

[Thermodynamic parameters:](#)

[Sequence length:](#)

Output options: __

Click here to , or

click here to the form to defaults.

If you have comments or suggestions on this service, you can send us [mail](#) here !

BiophysWWW / [G. Steger](#) Oct. 1996

Working with the web-based **POLAND** program only need 4 steps:

1. Enter DNA sequence (≤ 1000 bases)
2. Enter mismatch position (**optional**)
3. Choose **GIF** format
4. Press submit

7.3 The “new” Poland program (New Server)

7.3.1 About the Poland service

The Poland server will calculate thermal denaturation profiles and temperature-dependent UV absorbance or gel mobility of double stranded RNA, DNA, or RNA/DNA-hybrids based on sequence input and parameter settings in the Poland request form. -- Details of the Poland program are given below.

The program used in these calculations was developed by Gerhard Steger for comparing theoretical predictions to experimental data, mainly optical denaturation profiles, taken at 260 and 280 nm, and TGGE (temperature gradient gel electrophoresis) experiments.

The original version was written in VAX Fortran (VMS), using the Graphics Kernel System GKS for data presentation.

7.3.2 Program-specific information

Calculation is based on D. Poland's algorithm including the modification by Fixman & Freire in the implementation described by Gerhard Steger. The Poland algorithm calculates the denaturation profile for double-stranded nucleic acid using nearest-neighbor stacking interactions and loop entropy functions described in the literature. The input data required for calculation are:

- the sequence, of course (and no default here!),
- optional mismatched positions,
- the strand concentration, affecting the dissociation temperature,
- the method to calculate the final dissociation into single strands,
- the thermodynamic parameter set (DNA/DNA low salt/RNA), and
- the temperature range in which the calculation is performed.

In case you need access to the full range of input options, more options are available to the experts. Data sets predicted and figures drawn by the program are described below; see also for OUTPUT.

- A perspective view on the temperature-dependent denaturation profile (denaturation probability vs. sequence position vs. temperature. This plot does not include the dissociation of dsNA into single-strands; thus it shows most clearly the relative stability of the different parts of the NA.
- The temperature-dependent relative UV hypochromicities as measured in optical melting, at wavelengths of 260 and 280 nm (282 nm in case of DNA), respectively (full hypochromicity corresponds to approx. 30% of the OD at low temperature).
- The derivative form of above hypochromicities, showing the melting temperature(s) and corresponding half width(s) of the transitions(s), giving hints about the transition cooperativity.

- Predicted relative gel mobility, calculated according to Lerman et al., vs. temperature for different values of the 'retardation length' parameter L_r . This plot can be used for direct comparison with TGGE experiments; superpositions of plots generated with or without mismatched positions given are useful as a hint whether specific mismatched duplexes could be detected among homoduplexed DNA in a mixture of sample and reference double strands having undergone a denaturation-renaturation cycle, using either perpendicular or parallel TGGE.
- A 'half-denaturation temperature' plot showing the temperature at which each base pair has a probability of 50% to be in the open state. Similar to the three-dimensional plot, this plot can be used to estimate the destabilizing effect of mismatches on the surrounding part of the sequence: a temperature-shift of the TGGE transition can be expected if the lowest melting part of the sequence is directly affected by the mismatch.

Calculations can be done for oligonucleotides (>15 bases) or long double strands (>50 bases), respectively. In the case of oligonucleotide mode, a length-dependent correction for the strand dissociation process is applied. We do not have sufficient experimental results to stringently check for this mode to give valid results, but for the length range of about 20 nucleotides there is at least experimental evidence. Using 'oligo' mode with far longer sequences gives misleading results!

7.3.3 References for Poland Service

Description of implemented programs

Steger, G. (1994). *Nucleic Acids Res.* **22**, 2760-2768.

Thermal denaturation of double-stranded nucleic acids: prediction of temperatures critical for gradient gel electrophoresis and polymerase chain reaction.

Original version of algorithm:

Poland, D. (1974). *Biopolymers* **13**, 1859-1871.

Recursion relation generation of probability profiles for specific-sequence macromolecules with long-range correlations.

Fixman & Freire (1977). *Biopolymers* **16**, 2693-2704.

7.3.4 HELP for Poland Service:

7.3.4.1 POLAND

The program POLAND simulates transition curves of double-stranded nucleic acids (DNA and RNA as well as DNA/RNA hybrids).

Additional information available:

OUTPUT, RELATED PROGRAMS, RESTRICTIONS, ALGORITHM,
SUGGESTIONS, PARAMETERS

7.3.4.2 OUTPUT

The program writes its output in numeric format, which is converted to graphics by Tk/Tcl.

Additional information available:

General description, resolution, 3Dplot, GelPlot, MeltPlot, Temp50%Plot

Resolution of graphics output

The primary graphics output is produced as PostScript® (vector format). That format is converted to GIF® (raster format); this is a format directly displayed by your WWW browser. The resolution of the GIF images is selectable: 72 dots per inch (72 dpi) is the standard screen resolution; 150 dpi or 300 dpi are nice for printing. But be aware on NanoWeak® systems: the higher resolutions need a lot of memory and tend to crash your system.

3DPlot

Probability of an open base-pair is plotted as a function of position in sequence and temperature.

GelPlot

Relative mobility is plotted as a function of temperature for the three different stiffness parameters.

MeltPlot

Relative hypochromicity and its derivative is plotted as a function of temperature at 260 nm and 282 nm (RNA 280 nm).

References:

for RNA:

Coutts, S.M. (1971). *Biochim. Biophys. Acta* **232**, 94-106. Thermodynamics and kinetics of GC base pairing in the isolated extra arm of serine-specific tRNA from yeast

for DNA:

Blake, R.D. & Haydock, P.V. (1979). *Biopolymers* **18**, 3089-3109. Effect of sodium ion on the high-resolution melting of lambda DNA

Temp50%Plot

Temperature is plotted at which the corresponding base stack has a probability of 50% to be in the open state. The two horizontal lines in the plot mark the temperature range of calculation; i.e., a curve coinciding with such a line is not valid.

7.3.4.3 RELATED PROGRAMS

The Poland program calculates the denaturation behavior of double-stranded NA. LinAll, RNAfold, and mFold calculate the secondary structure of single-stranded (R)NA; in addition LinAll and RNAfold allow the prediction of denaturation behavior of ssRNA.

7.3.4.4 Restrictions

The sequence has to be shorter than 1001 nucleotides but longer than 5 nucleotides.
Valid nucleotides are A, G, C, U, and T.

Calculation of asymmetric or bulge loops is not possible.

7.3.4.5 Algorithm

Calculation is based on Poland's algorithm including the modifications proposed by Fixman & Freire.

With the original algorithm of Poland computing time is proportional to the square of the sequence length.

With the modification according to Fixman & Freire computing time is proportional to 10 times the sequence length, but it works only with loop parameters according to Poland.

7.3.4.6 SUGGESTIONS

Hints for combination of parameters and their values.

Additional information available:

RNA

DNA

RNA/DNA

Ionic_strength_dependence

RNA

Thermodynamic values according to Turner et al. are ideally suited for calculation in 1 M ionic strength after correction of all DeltaS values by 1.021 and all DeltaSGC values by 0.961. These corrections are equivalent to a shift in T_m values of A:U stacks by -7 K or -2%, respectively and of G:C/G:C stacks by +7 K or +2%, respectively.

Optimal (?) parameter combination for Turner *et al.*:

-d 1.021 1.000 0.961	(DeltaS, DeltaS(A:U), and (DeltaS(G:C) factor)
-n 1.e-3	(Dissociation constant β)
-c 1.e-6	($\beta \cdot c_0 = 1E-8$ to $1E-10$)
-s 1.000	(loop parameter Sigma)
-l g	(internal loops according to Gralla & Crothers)
-t 90. 120. 0.5	(Temperature range and steps)

Optimal (?) parameter combination for Pörschke et al.:

-d 1.000 1.040 0.970 (DeltaS, DeltaS(A:U), and (DeltaS(G:C) factor)
 -n 1.e-3 (Dissociation constant β)
 -c 1.e-6 ($\beta \cdot c_0 = 1E-9$ to $1E-11$)
 -s 1.e-6 (loop parameter Sigma)
 -l p (internal loops according to Poland)
 -a f (algorithm according to Fixman & Freire)
 -t 90. 120. 0.5 (Temperature range and steps)

DNA

Thermodynamic values according to Gotoh et al. and Klump, both, are ideally suited for calculations. The parameter set of Breslauer et al. does not fit our experiments (?). The parameter set of Allawi & SantaLucia is based on a reevaluation of all known parameter sets for DNA; i.e., this set may be the optimal one. For references to the original thermodynamic parameter sets see here.

Additional information available:

Gotoh

Breslauer *et al.*

Klump

SantaLucia *et al.*

Allawi & SantaLucia

- **Gotoh**

Optimal (?) parameter combination for Gotoh:

-d 1.000 (DeltaS factor)
 -n 1.e-3 (Dissociation constant β)
 -c 1.e-6 ($\beta \cdot c_0 = 1E-9$ to $1E-11$)
 -s 1.e-3 (loop parameter Sigma)
 -l p (internal loops according to Poland)
 -a f (algorithm according to Fixman & Freire)
 -t 60. 80. 0.5 (Temperature range and steps)

- **Breslauer *et al.***

No optimal parameter combination for Breslauer et al.

- **Klump**

Optimal (?) parameter combination for Klump:

-d 1.000 (DeltaS factor)
 -n 1.e-3 (Dissociation constant β)
 -c 1.e-6 ($\beta \cdot c_0 = 1E-9$ to $1E-11$)
 -s 1.e-3 (loop parameter Sigma)
 -l p (internal loops according to Poland)
 -a f (algorithm according to Fixman & Freire)
 -t 70. 90. 0.5 (Temperature range and steps)

- **SantaLucia et al.**
- **Allawi & SantaLucia**

This is the "unified parameter set"!

RNA/DNA

Thermodynamic values according to Sugimoto et al. (1995) in 1 M NaCl. The top strand is RNA, the bottom strand is DNA (5'-r-3'/3'-d-5'); the input sequence is the RNA strand.

Ionic strength dependence

Following values may be used for correction of calculated T_m-values:

$$\frac{T_{m,2} - T_{m,1}}{\log(c_2/c_1)} = f(G:C) \cdot I(G:C) + (1-f(G:C)) \cdot I(A:U)$$

with T_m = transition (midpoint, melting) temperature
 c = ionic strength (=concentration of Na ions)
 f(G:C) = G:C content
 I(X:Y) = dependence of ionic strength of
 base pair type X:Y

DNA
 I(A:T) = 18.3 °C (Owen, Hill, & Lapage (1969).
 Biopolymers **7**, 503-516.)

 I(G:C) = 11.3 °C (Frank-Kamenetskii (1971).
 Biopolymers **10**, 2623-2624.)

RNA
 I(A:U) = 20.0 °C (Steger, Müller & Riesner (1980).
 I(G:C) = 8.4 °C *Biochim. Biophys. Acta* **606**, 274-284.)

7.3.4.7 PARAMETERS

Additional information available:

BASE_STACKING_(thermodynamic_parameters)

ENTROPY_CORRECTION_of_base_stacking

LOOP_PARAMETERS_(thermodynamic_parameters)

TEMPERATURE_RANGE_OF_CALCULATION

MISMATCHED_POSITIONS_in_original_sequence

CONCENTRATION_and DISSOCIATION_CONSTANT

STIFFNESS_of_nucleic_acid

THERMODYNAMIC PARAMETER SETS for BASE STACKING

You can select between five different thermodynamic parameter sets of base stacking (for loop parameters see below):

Additional information available:

RNA

DNA

RNA/DNA

- **RNA**

- for RNA in 1 M NaCl

Freier, S.M., Kierzek, R., Jaeger, J.A., Sugimoto, N., Caruthers, M.H., Neilson, T. & Turner, D.H. (1986). *Proc. Natl. Acad. Sci. USA* **83**, 9373-9377.

Improved free-energy parameters for predictions of RNA duplex stability.

- for RNA in 1 M NaCl

Pörschke, D., Uhlenbeck, O.C. & Martin, F.H. (1973). *Biopolymers* **12**, 1313-1335.

Thermodynamics and kinetics of the helix-coil transition of oligomers containing GC base pairs.

- **DNA**

- for DNA in 0.019 M NaCl

Gotoh, O. (1983). *Adv. Biophys.* **16**, 1-52.

Prediction of melting profiles and local helix stability for sequenced DNA. •for DNA

- in 1 M NaCl

Breslauer, K.J., Frank, R., Bloecker, H. & Marky, L.A. (1986). *Proc. Natl. Acad. Sci. USA* **83**, 3746-3750.

Predicting DNA duplex stability from the base sequence. •

- for DNA in 0.1 M NaCl

Klump, H.H. (1987). *Canad. J. Chem.* **66**, 804-809.

Energetics of order/order transitions in nucleic acids.

Klump, H. (1990). in Landolt-Börnstein, New Series, Group VII Biophysics, Vol. 1 Nucleic Acids, Subvol. c Spectroscopic and Kinetic Data, Physical Data I, (W. Saenger, ed.), Springer-Verlag Berlin, p. 244-245.

Calorimetric studies on DNAs and RNAs.

- for DNA in 1 M NaCl

SantaLucia, J. Jr., Allawi, H.T. & Seneviratne, P.A. (1996). *Biochemistry* **35**, 3555-3562.

Improved nearest-neighbor parameters for predicting DNA duplex stability.
for DNA in 1 M NaCl

Allawi, H.T. & SantaLucia, J. Jr. (1997). *Biochemistry* **36**, 10581-10594.

Thermodynamics and NMR of Internal G-T Mismatches in DNA.

- **RNA/DNA**

- for RNA/DNA hybrids in 1 M NaCl

Sugimoto, N., Nakano, S., Katoh, M., Matsumura, A., Nakamuta, H., Ohmichi, T., Yoneyama, M. & Sasaki, M. (1995). *Biochemistry* **34**, 11211-11216.

Thermodynamic parameters to predict stability of RNA/DNA hybrid duplexes.

- **ENTROPY CORRECTION of base stacking**

(Option not available by WWW)

DeltaS values of base stacking may be corrected by factors in order to simulate deviating ionic strengths.

The Delta S values of the thermodynamic parameters are multiplied with these factors. The first is used for correction of all Delta S values, the second only for A:U stacks, the third only for G:C stacks.

Different values are used as defaults in dependence on the chosen thermodynamic parameter set.

- **LOOP PARAMETERS (thermodynamic parameters)**

(Option not available by WWW; i.e., Sigma is fixed to 1.e-3, and loop entropy is calculated according to Poland.)

Loops which appear during denaturation by internal base stack opening may be calculated by three different methods:

- -l p ==> $\Delta S(\text{loop}) = \text{SIGMA} * (\text{loop} + 1)^{-1.75}$
(according to Poland or Fixman & Freire)
Use only with stacking parameters according to
Pörschke *et al.* (-p p),
Gotoh (-p g), or
Klump (-p k).
- -l g ==> $\Delta S(\text{loop}) = \text{SIGMA} * \Delta S(\text{loop})$
(according to Gralla & Crothers)
Use only with stacking parameters according to
Turner *et al.* (-p t).
- -l t ==> $\Delta S(\text{loop}) = \text{SIGMA} * \Delta S(\text{loop})$
(according to Turner *et al.*)
Use only with stacking parameters according to
Turner *et al.* (-p t).

Therefore, Sigma influences the cooperativity and the half width of each transition. With '-a f' you can change the default algorithm (only in case '-l p'). With the original algorithm of Poland (default), computing time is proportional to the square of the sequence length. With '-a f' the modified algorithm of Fixman & Freire is used which results in computing time proportional to 10 times the sequence length but works only with loop parameters according to Poland (-l p) up to a sequence length of 1000 base pairs.

References:

Poland, D. (1974) *Biopolymers* **13**, 1859-1871. Recursion Relation Generation of Probability Profiles for Specific-Sequence Macromolecules with Long-Range Correlations.

Fixman and Freire (1977) *Biopolymers* **16**, 2693-2704. Theory of DNA melting curves.

Gralla, J. & Crothers, D.M. (1973) *J. Mol. Biol.* **78**, 301-319. Free energy of imperfect nucleic acid helices. III. Small internal loops resulting from mismatches.

Freier, S.M., Kierzek, R., Jaeger, J.A., Sugimoto, N., Caruthers, M.H., Neilson, T. & Turner, D.H. (1986) *Proc. Natl. Acad. Sci. USA* **83**, 9373-9377. Improved free-energy parameters for predictions of RNA duplex stability.

• TEMPERATURE RANGE OF CALCULATION

The temperature range for calculations has to be adapted to the other parameters; thus see the topic Suggestions.

In principal not more than 110 temperature points are allowed for a single calculation.

• MISMATCHED POSITIONS in original sequence

Mismatches are given as a comma-separated list of sequence positions; f.e.

-m 2,3,111

specifies mismatched 'base pairs' at positions 2, 3, and 111. If the mismatch is longer than a 'base pair', the position of each base has to be given separately. The sequence position of the mismatched base pair(s) may be given in any order. Calculation of asymmetric or bulge loops is not possible; these have to be modeled by larger mismatches (internal loops).

• CONCENTRATION and DISSOCIATION CONSTANT

-c 1.e-6 Concentration of single strands C0

-n 1.e-3 Dissociation constant β

$\beta \cdot c_0$ influences temperature T_m and half width of the second order transition, i.e. the strand separation.

The dissociation constant β has to be in the range $1. \Rightarrow \beta \Rightarrow 1.E-5.$, and the strand concentration has to be in the range $1. \Rightarrow C_0 \Rightarrow 1.E-13.$

If case of short oligonucleotides, β might be calculated according to

Benight, A.S. & Wartell, R.M. (1983). *Biopolymers* **22**, 1409-1425. Influence of base-pair changes and cooperativity parameters on the melting curves of short DNAs.

and

Benight, A.S., Wartell, R.M. & Howell, D.K. (1981). *Nature* **289**, 203-205. Theory agrees with experimental thermal denaturation of short DNA restriction fragments.

- **STIFFNESS of nucleic_acid (Lr)**

Stiffness of nucleic acid or gel pore size is given f.e. with
-I 40 90 200.

References:

Lerman, L.S., Fischer, S.G., Hurley, I., Silverstein, K. & Lumelsky, N. (1984). *Ann. Rev. Biophys. Bioeng.* **13**, 399-423.

Fischer, S.G. & Lerman, L.S. (1982). *Proc. Natl. Acad. Sci. USA* **80**, 1579-1583.

Riesner, D., Henco, K. & Steger, G. (1991). In: *Advances in Electrophoresis*, Vol. 4 (Chrambach, A., Dunn, M.J. & Radola, B.J., eds.) VCH Verlagsgesellschaft, Weinheim, pp. 169-250.

Temperatur-gradient gel electrophoresis: A method for the analysis of conformational transitions and mutations in nucleic acids and proteins

7.3.5 How to use the “new” Poland program

Poland service request form

Sequence title line:

Sequence:

(plain format;
no numbers;
max. 1000 nts;
min. 5 nts)

Mismatched positions:

(comma-separated numbers)

Thermodynamic parameters:Dissociation constant β :Strand concentration:

(default: 1.0E-6 M)

Temperature range:Which graphics do you want:Graphics size:

(GIF format)

DNA (100 mM NaCl, Klump). ☐Oligonucleotide
(β is function of seq.length)Long double strand
(default: $\beta=1.0E-3/M$)Low temperature limit:
(default: 40.0°C)High temperature limit:
(default: 110.0°C)Temperature step
size: (default: 2.0°C)

$T_m(p=50\%)$ plot	3d plot	Mobility plot	Melting curve	Diff. melting curve
☒	☒	☒	☒	☒

72x72 dpi ☐Click here to , orClick here to the form to defaults.

7.4 The optimized DNA fragment

The optimized fragment for detection of point mutations in dsDNA is derived by PCR amplification and has a length of 200 / 300 to max. 800 / 900 bp. It consists of 1 - 2 melting domains derived from the native sequence plus a synthetic stabilizing region (**GC-clamp**). The GC clamp is highly stable because of 3 hydrogen bridges between G and C whereas there are only 2 hydrogen bridges between A and T. This clamp may either consist of a 40 bp artificial stretch of GC base pairs (29) or a covalent chemical clamp (**Psoralen** = Furo[3,2-g]coumarin, $C_{11}H_6O_3$). Psoralen is the better choice for temperature gradients at high temperatures because Psoralen intercalates with the double helix and after UV-treatment it links both strands covalently (irreversible binding). Both clamps are introduced into the DNA fragment by a 5' overhang of one of the PCR amplification primers. (For easier reading of the following text both kinds of clamps will be referred to as "GC-clamps" further on.)

The melting properties of a DNA fragment are best described by the two-dimensional "TempPlot diagram" (fig. 21). In this diagram the T_m value is given on the y-axis and the base pair number on the x-axis. The optimal fragment for TGGE analysis shows a "stair type" profile with 2 or 3 "steps", respectively (fig. 21, lower diagram). The highest "step" (the melting domain with the highest melting Temperature T_m) is the artificial GC-clamp. Length and midpoint melting temperature of the 1 - 2 lower "steps" (melting domains with lower T_m s) are determined by the original sequence of the DNA under study.




Figure 21: Constructing the optimized DNA fragment

When constructing an optimized fragment, start with the "TempPlot" of an approximately 1000 bp sequence. For PCR amplification, select a fragment consisting of 1 - 2 melting domains. Put the GC-clamp at the more stable end or at any end, if the fragment contains one melting domain. Figure 22 schematically demonstrates how fragments with the utmost longest part of the sequence could be selected according to the "TempPlot" diagram.

The primers used for PCR amplification have to meet the following rules:

- Use non complementary primer sequences. Do not allow base-pairings of the last 3 bases at the 3'-end, neither with any other bases in the primer itself, nor with the counterpart primer.
- Select primers 20 - 25 bp in length
- Be sure that there are no additional primer annealing sites in the DNA sequence.




Figure 22: Schematic "TempPlot" diagram of a DNA sequence

Fragments which have not been "clamped" may also be analyzed on TGGE, but they must contain at least two melting domains. In this case, the most stable melting domain may act as a "natural GC-clamp", provided that electrophoresis is terminated before this second domain reaches its respective T_m . Under these experimental conditions nucleotide changes within this highest melting domain will not lead to a shift of bands on the TGGE gel, and hence, will not be detectable. In conclusion, the absence of GC-clamps will still allow nearly 100% detection rate for mutations in the low melting temperature domain(s), but virtually no ability to detect mutations in the domain with the highest T_m .

7.4.1 Asymmetric GC-clamps for PCR primers used for TGGE analysis

The 5' end or the 3' end of the primer for the end of the segment at which a clamp is optimal must carry a GC-clamp. The length of the clamp depends on the sequence of the sample. The denaturing behavior of the modified sample can be tested using the POLAND software.

short GC-clamp (23 bp):

cccgc cgcgcc cccgc cgccc gcc

long GC-clamp (40 bp)⁴⁴:

cgccc gccgc gcccc gcgcc cggcc cgccg ccccc gcccg

long GC-clamp (39 bp)⁴⁵:

ccccg ccccc gccgc ccccc ccgcg cccgg cgccc ccgc

7.4.2 Chemical clamp with Psoralen (Furo[3,2-g]coumarin, $C_{11}H_6O_3$)


The 5' end of the primer for the end of the segment at which a clamp is optimal (POLAND program) must carry a appropriately linked psoralen moiety at the end adjacent to T or A, preceding the genomic sequence. The optimal primer sequence may be 5'(Pso)pTaPpnpnp.....3', given the preference of psoralen for binding between TpA and ApT pairs^{13,46,47}. Crosslinking of the PCR product is done e.g. in a flat-bottom microtiter plate using a 365 nm UV source. Working with small volumes it may be necessary to minimize evaporation by cross-linking at 4 – 10°C. The yield is not affected by temperature. The distance of the sample from the UV source affects the yield. 15 min at 0.5 cm distance of the sample from an 8 W UV lamp is sufficient.

7.4.3 POLAND analysis of samples

Unoptimized DNA fragment


Domains not distinctly different

Poland analysis of Withgc_18_30_43.Seq on September 23, 1998 18:30:47.54
Sequence vs. Probability vs. Temperature



Second order line is flat


Poland analysis of Withgc_18_30_43.Seq on September 23, 1998 18:30:57.16
Stack-index vs. Temp. of 50% probability



Optimized DNA fragment (GC clamp attached)


Domains distinctly different

Poland analysis of Withgc_18_28_34.Seq on September 23, 1998 18:28:37.77
Sequence vs. Probability vs. Temperature



Two distinct second order line plateaus


Poland analysis of Withgc_18_28_34.Seq on September 23, 1998 18:28:48.33
Stack-index vs. Temp. of 50% probability



Optimized DNA fragment with mismatch


Mismatch position

Poland analysis of Mismatch_18_29_34.Seq on September 23, 1998 18:29:38.55
Sequence vs. Probability vs. Temperature



Mismatch changes melting behavior

Poland analysis of Mismatch_18_29_34.Seq on September 23, 1998 18:29:49.73
Stack-index vs. Temp. of 50% probability



8 Optimizing parallel TGGE by perpendicular TGGE

8.1 Check short DNA fragments for their melting behavior

All short DNA fragments (100 - 150 bp) should be checked first by a perpendicular TGGE gel. This is not only a good place to start for practical optimization of parallel TGGE, but also verifies the reversible melting behavior of the DNA fragment (fig. 23). "Reversible melting" can only occur if the DNA fragment consists of at least two separate melting domains (fig 23 a, c, d). Reversible melting behavior must be verified since it is required for successful parallel TGGE analysis. Thus, be sure to check all fragments on a perpendicular TGGE gel which do not contain 40 bp GC-clamp or which have not been evaluated with the aid of the POLAND program.




Figure 23: Short DNA fragments without GC-clamp

8.2 From perpendicular to parallel TGGE

Perpendicular TGGE routinely uses a standard temperature gradient from 20 - 60°C in combination with buffers that contain minimum 8 M urea. Using this same standard gradient in parallel TGGE is possible, but time-consuming, since a longer running time is required to move the sample from the slot (top of the gel) to the effective range of separation (middle or lower part of the gel where the domains begin to melt). In other words, much of the time consumed by electrophoresis is nonproductive since melting will not begin to occur until the DNA fragment has migrated a considerable distance.

Parallel TGGE can be easily optimized by the information acquired from a preliminary perpendicular TGGE gel.




Figure 24: From perpendicular to parallel TGGE

Steps:

1. Determine the temperature range of effective separation. This temperature interval is defined by two temperatures, T_{low} and T_{high} . If possible, use two different DNA fragments, one wild-type and one mutated, and perform heteroduplex analysis. In this case, T_{low} is determined by the highest temperature where all duplices remain double-stranded (no retardation in the gel), T_{high} is determined by the melting temperature (T_m) of the most stable homoduplex (see fig. 24a). In the event that only a wild-type sequence is available, determine T_{high} by the melting temperature T_m of this sequence, and define $T_{low} = (T_{high} - 10^\circ\text{C})$.
Note: The temperature range of effective separation will have to be determined for each new DNA fragment analyzed.
2. Using the information from step 1, select a temperature gradient for parallel TGGE which overlaps the range of effective separation. Program the temperature of T_{low} for L1 (first thick lane on the gradient block, close to T_1) and program the temperature of $T_{high} + 5^\circ\text{C}$ for L6 (last thick lane on the gradient block, close to T_2) for the parallel TGGE run (fig. 24b).

9 TGGE / SSCP

9.1 Running an SSCP on the TGGE

TGGE can be used in combination with SSCP ("single strand conformation polymorphism") to improve (often dramatically) the frequency of detecting SSCP markers. TGGE/SSCP is non-radioactive, because it utilizes silver staining detection. SSCP relies upon the separation of single-stranded DNA or RNA which have formed hairpin secondary structures. Different conformations exhibit different electrophoretic mobilities. The conformation which a particular single-stranded molecule adopts is sequence-dependent, and mutations are detected by their influence upon the secondary structure, and hence, the altered electrophoretic mobility.




Figure 25: *Effect of conformational differences*

As indicated in figure 25, conformational differences between two different base-substituted fragments can be achieved only within a limited temperature range (T_{low} - T_{high}): At temperatures below T_{low} , both fragments adopt the characteristic hairpin structure. On the other hand, at temperatures higher than T_{high} , neither fragment will form the hairpin, and both will have the same mobility. Only at temperatures within the range of T_{low} - T_{high} will the base-substituted DNA be distinguishable from the wild-type DNA.

By virtue of the temperature gradient which TGGE imposes upon the gel, one particular area of the gel will provide the appropriate temperature range (i.e., T_{low} - T_{high}) to allow formation of the hairpin, and hence, visualization of differences in mobility.

9.2 DNA sample preparation

Add 3 μ l of 95% formamide/10 mM EDTA to a 3 μ l aliquot of the PCR-amplified sample. Heat to 90°C for 5 min. **Immediately chill on ice!!!!!!!**

Note: For denaturation, please refer to the above mentioned protocol. **Do not use NaOH!!!!!!**

9.3 Gel casting

Recipe for 10 ml gel solution (3 – 4 gels) for SSCP-ME running buffer:

	3% Gel	5% Gel	8% Gel
Acrylamide/bis Acrylamide Stock solution (30 : 0.8), 40% (w/v)	0.75 ml	1.25 ml	2.0 ml
50x conc. SSCP-ME-buffer ($C_{End} = 1 \times \text{conc.}$)	0.2 ml	0.2 ml	0.2 ml
40% Glycerol ($C_{End} = 2\%$)	0.5 ml	0.5 ml	0.5 ml
Water, distilled	2 ml	1.5 ml	1 ml
Make sure that the urea has been completely resolved. It is possible to heat up the urea containing solution slightly (40°C – 50°C) for a short time in order to improve the solubilization of urea. De-gas the solution under gentle vacuum for 3 - 5 min.			
Water, distilled	fill up to 10 ml		
TEMED	17 μ l	17 μ l	17 μ l
APS	76 μ l	76 μ l	76 μ l
Mix gently. Avoid air bubbles!			
Pour the gel solution into the glass plate sandwich immediately thereafter (see chapter 4.1.2) without air bubbles.			

Cast the TGGE gels according to the instructions given under "Setting up polyacrylamide gels".

9.4 Electrophoresis

Perpendicular TGGE with a 12 or 18 slot gel (as normally used for parallel TGGE) is the fastest approach to determine the optimal conditions for SSCP.

Proceed as described for parallel TGGE. Establish the temperature gradient from 5°C (cathode, black) to 30°C (anode, red) and use SSCP-ME buffer for the run. Run at 200 V for 2 - 5 hours. Covering the gel with a protective cellophane sheet is not absolutely necessary because the gel will not dessicate at the temperatures used in this protocol.

Silver stain the gel.

9.5 Routine analysis

When the temperature range of SSCP separation has been determined, routine analysis can be performed in a gel with constant temperature.

10 TGGE in RNA analysis

TGGE is a perfect tool to analyze RNA for secondary structures. Applications have been published on the differentiation of plant pathogen variants (1, 33, 34, 35, 36, 37), analysis of intermediates of plant pathogens (1, 11, 30, 31, 32, 33), and also on the analysis of hairpin structures in m-RNA (1, 38). The technique of choice is perpendicular TGGE. Depending on the species of RNA (m-RNA, r-RNA etc.) which is to be analyzed, standard protocols and conditions described in this manual are generally applicable, but may require minor modifications.

10.1 Completely double-stranded RNA

The T_m values of a dsRNA sequence will be about 20°C higher in comparison to the corresponding ds-DNA sequence. For dsRNA with low GC-content, start with the standard protocols using ME-buffer. For very GC-rich sequences, raise the temperature in the TGGE gel (40°C - 80°C) or lower the ionic strength in the electrophoresis buffer and gel by using another buffering system (electrophoresis buffer: 0.1 x conc. TBE: 8.9 mM Tris, 8.9 mM boric acid, 0.24 mM EDTA; gel: 0.1 x conc. TBE buffer, 5% polyacrylamide, 8 M urea; temperature gradient 35 - 60°C; published for analysis of dsCARN5 (1, 33, 34, 36, 37) and reovirus RNA (34)).

Note: The reduction in ionic strength will lower the T_m !!!

10.2 Partly double-stranded RNA, e.g. viroid RNA

Use the suggested protocols provided in this manual for ME-buffer or refer to the buffer and gel systems published in various papers (electrophoresis buffer: 0.2 x conc. TBE: 17.8 mM Tris, 17.8 mM boric acid, 0.4 mM EDTA; gel: 0.2 x conc. TBE, 5% polyacrylamide, no urea (1, 32, 33, 34, 37)).

10.3 Single-stranded RNA with single hairpin structures, m-RNA secondary structures

Use standard protocols given in this manual for ME-buffer or refer to other buffer and gel systems described in the literature (38) (electrophoresis buffer: 10 mM sodium phosphate, pH = 6.0, with or without 1 mM $MgCl_2$; gel: 8% polyacrylamide, 10 mM sodium phosphate, pH = 6.0, with or without 1 mM $MgCl_2$).

10.4 Staining

For detection of the RNA, silver-staining is recommended. For identification of double-stranded virus RNA from crude plant extracts, a protocol based on immunoblotting has been published (35).

11 TGGE in protein analysis

TGGE can be successfully applied to investigation of protein/Structural transitions, and also thermostability of protein-nucleic acid interactions. In comparison to conventional methods such as spectroscopy, hydrodynamics or calorimetry, TGGE analysis offers several advantages:

- Only minimal amounts of sample material are required.
- TGGE may be carried out by using crude protein extracts.
- The effect of additives that influence the protein stability may easily be investigated.

11.1 Buffers

In contrast to nucleic acids, which can generally all be analyzed with standard buffer conditions, each protein requires its own special buffer system. This buffer has to fulfill the same requirements as those for native gel electrophoresis of proteins:

- The protein has to be in its native conformation at low temperature (e.g. room temperature).
- The protein has to carry a net charge.
- The protein has to be soluble in the used buffer.
- The protein has to migrate as a homogeneous band under native conditions.

This single one prerequisite has to be tested on a non denaturing gel before testing the denaturation behavior on TGGE.

Additionally:

- The pH value of buffer should be virtually independent of the temperature.
- In order to avoid excessive current in TGGE gel above 100 mA, the ionic strength of the electrophoresis buffer should not exceed 30 mM.

Temperature dependence of the pH-value of different electrophoresis buffers:

Buffer	pH (20°C)*	$\Delta\text{pH} / \Delta T$ 50°C**
15 mM glycine / NaOH	11.9	- 1.75
30 mM H ₃ BO ₃ / NaOH	10.0	- 0.52
40 mM Borax + 20 mM H ₃ BO ₃	9.1	- 0.44
25 mM glycine / NaOH	9.4	- 1.66
30 mM Borax + 75 mM H ₃ BO ₃	8.6	- 0.33
89 mM Tris / H ₃ BO ₃	8.3	- 0.68
25 mM Tris / glycine	8.3	- 1.05
375 mM Tris / HCl	8.2	- 1.29
61 mM NaH ₂ PO ₄ / 10 mM Na ₂ HPO ₄	7.5	+ 0.08
25 mM Na ₂ HPO ₄ / 25 mM KH ₂ PO ₄	7.3	- 0.07
30 mM Na ₂ HPO ₄ / 8.7 mM KH ₂ PO ₄	6.8	- 0.04
125 mM Tris / HCl	6.8	- 1.39
23 mM Na ₂ HPO ₄ / 132 mM NaH ₂ PO ₄	6.0	+ 0.10
690 mM glycine + 240 mM H ₃ BO ₃	5.6	- 0.59
48 mM Tris / H ₃ PO ₄	5.5	- 0.18
48 mM KOH / acetic acid	4.8	+ 0.07
20 mM sodium acetate / acetic acid	4.5	+ 0.03
48 mM KOH / acetic acid	3.6	+ 0.09
690 mM glycine / acetic acid	3.5	- 0.17

* The pH-value of 20°C is that with the optimum buffer capacity.
 ** ΔpH is given for an increase in temperature of 50°C ($\Delta T = 50^\circ\text{C}$).

In the literature, different buffer systems have been reported for the following proteins:

Dehydrogenases (41), β -lactamase (1), tet-repressor from *E. coli* (33, 42), alpha-amylases (1, 34, 39), and serine proteases (40).

12 Trouble-shooting

The following trouble-shooting guide may be helpful in solving any problem that you may encounter. If you need further assistance, please do not hesitate to contact your local Biometra distributor or Biometra directly.

In any case where you recognize a failure which is marked in the list by an exclamation mark, please stop working with the instrument and call the local representative for replacing faulty parts.

Problem	Cause	Solutions
Preparing the gel solution		
Urea can not be dissolved in gel solution	1. Dissolving urea is an endothermic process and requires energy in form of heat.	1. Heat up the acrylamide/urea solution – but not more than 40°C – 50°C. Mix the solution.
Gel does not polymerize	1. Old chemicals. 2. Gel solution prepared incorrectly. 3. Too much oxygen in the gel solution.	1. Prepare acrylamide/bis-acrylamide solution freshly. Prepare 4% APS freshly and freeze in small aliquots. 2. Check all reagents that have been included in the gel solution and mix thoroughly. 3. Degas solution before adding TEMED and APS.
Gel polymerizes too fast	1. Too much TEMED and APS has been added to the gel solution. 2. The gel solution has been heated in order to dissolve the urea.	1. Check the concentrations of TEMED and APS. Use the amounts given in the standard protocol. 2. Allow the gel solution to cool down to room temperature before adding TEMED and APS. (Note: the gel solution should not be warmed up to more than 50°C.)

Preparing the gel sandwich		
Gel solution leaks out of sandwich	1. Gel sandwich has been set up incorrectly 2. Clips not correct positioned. 3. Scratches on the spacers or old clips.	1. Clean spacers with methanol. 2. Fasten the clips above the spacer to increase the pressure. 3. Use silicone grease along the glass spacer, but never on the sample slots!!!!
Polybond film with gel can not easily removed from the sandwich	1. Gel sticks to the glass plate with spacer	1. Glass plate with spacer must be treated with Acryl-Glide before use.
Gel does not stick to the Polybond film	1. Gel has been poured onto the hydrophobic side of the Polybond film.	1. Pour the gel onto the hydrophilic side if you want to link the gel covalently to the Polybond film. Hydrophobic face of the Polybond film must face the bonding plate. (Check the Polybond film with a drop of water for the hydrophilic and hydrophobic side.)
Front (top) of gel shows a zig zag line	1. Gel has not been overlaid with solution.	1. Overlay the gel solution in the sandwich with 200µl Isopropyl- or Isobutyl-Alcohol. (Alternatively bi-distilled water can be used.)
Air bubbles in the gel	1. The glass plate has not been cleaned carefully. 2. Sandwich was hold vertical during pouring. 3. Solution poured to quick or in the middle of the sandwich. 4. The glass plate needs treatment with Acryl-Glide	1. Clean the glass plate and slot formers with ddH₂O and EtOH before use. Avoid the intensive use of organic solvents. They will dissolve the glue of the spacer and slot formers and thus remove them from the glass plate. 2. Hold the sandwich at an angle of 45° during pouring. 3. Pour the gel solution slowly along one side of the glass plate. 4. Treat glass plate with spacer with Acryl-Glide
Air bubbles between slots	1. The glass plate needs treatment with Acryl-Glide	1. Treat glass plate with spacer with Acryl-Glide before each use. (Clean spacers with EtOH!)
Slots are distorted	1. Gel sticks to the slot former and/or glass plate. 2. Polybond film with gel have been removed beginning from the bottom or to quick.	1. Treat glass plate with spacer with Acryl-Glide before each use. 2. Remove slowly the Polybond film with gel from the glass plate with spacer beginning from the top.

Problem	Cause	Solutions
Electrophoresis unit and Controller		
Scratches in the white cover film of the gradient block	1. Cover film of the gradient block damaged.	1. Remove the cover film and replace by a new one.
Gel running		
No or minimal current (< 5 mA), Marker dyes stop in the gel.	1. Safety lid is not seated properly. 2. Assembly of the TGGE system is incorrect. 3. Gel is drying out. White opaque areas can be seen inside the gel. 4. Electrodes are dirty or damaged. 5. Wicks to dry and/or placed not correct. 6. Programmed voltage to low.	1. Position the safety lid correctly. 2. Check the assembly of the TGGE system and check plug connections. 3. Carefully protect the gel against evaporation. Take special care of the slots. After the sample has migrated into the gel, cover the gel with cover film and additionally with the special Cover glass plate (with 2 silicone barriers). 4. Check/clean the electrodes inside the buffer chambers. 5. Immerse the electrophoresis wicks in the buffer and place them properly on the gel. 6. Increase voltage.
Extremely high current (> 50mA)	1. High ionic strength in electrophoresis buffer, wrong buffer concentration.	1. Check the composition of electrophoresis buffer and gel.
Silver staining artifacts		
Cloudy yellow or brown staining		1. Wash extensively after silver binding step. 2. Use rocking table for staining protocols. 3. Place gel up side in the staining tray – use sufficient solution.
Brown spots in the gel	1. Un-dissolved crystals of urea remain in the gel 2. Contaminated chemicals 3. Gel casting glass plate or Cover film have not been cleaned properly.	1. Dissolve urea completely before gel pouring. 2. Use fresh stock solutions. Filter prior to use. Wear only non-powdered gloves during handling the gel. 3. Clean gel casting plate and the Cover film carefully before they come in direct contact with the gel. Wear gloves when handling the gel.

Problem	Cause	Solutions
Gel is stained completely black or looks like a "silver mirror"	1. Chloride ions are contaminating the staining solutions, the electrophoresis buffer or the gel solution. During the staining protocol silver chloride (AgCl) is precipitated in the gel (gel looks "milky" and is then reduced to elemental silver 2. High amounts of chloride ions contaminate the sample.	1. a) Check if tap water has been used for one of the buffers (tap water always contains chloride ions). Prepare and use only fresh solutions with ddH₂O. b) Check, if tap water has been used for washing the gel twice after incubation with 0.1% AgNO₃ in dest. H₂O. Use ddH₂O instead. c) If deionized water is used, its integrity should be checked (no chlorid ions): Take approximately 1 ml of the deionized water or the buffer you want to check and add some drops of 0.1% AgNO₃. If you see a milky precipitation (silver chloride, AgCl), the solution is contaminated with chloride ions. Use distilled water for preparing the buffers and washing gel. 2. Desalt the sample prior to TGGE.
Gel has a heavy background	1. Heavy background, caused by smearing of the sample: High amounts of proteins or polysaccharides (also stained by the silver) may contaminate the sample. 2. Old chemicals, especially acrylamide/bis solution 3. Silver-staining protocol has been carried out incorrectly.	1. Check the purity of the sample prior to TGGE. 2. Use fresh stock solutions. 3. Follow the silver staining protocol as exactly as described. Use an excess of freshly distilled water when washing the gel twice. Do not prolong the incubation in Developing solution.
No DNA bands are visible in the gel	1. No DNA, or amount of DNA sample is below the level of detection. 2. Too much DNA, inverse silver-staining.	1. Check the amounts of DNA. 2. If bands contain high amounts of DNA, the silver-staining may result in an inverse staining: the background is darker than the band itself. Reduce amounts of DNA.
Gel fades out	1. Stopping solution (0.75% Na₂CO₃) was not sufficient.	1. Incubate the gel in the staining buffer for 10 min.

Problem	Cause	Solutions
Band pattern in general		
Bands are diffuse	1. Sample volume is too large. 2. Sample is contaminated. 3. The DNA has undergone diffusion inside the gel, because the DNA has not been fixed after running the gel.	1. Reduce amount of DNA/RNA 2. Purify the sample in order to reduce contaminating proteins or polysaccharides 3. Proceed with silver staining protocol immediately after electrophoresis.
Interpretation of the TGGE band pattern		
DNA bands are diffuse	1. Sample volume is too large. 2. Too much DNA in the sample, gel is overloaded. 3. Temperature gradient has not been stable during electrophoresis. 4. The gel has been shifted during electrophoresis, thus the temperature gradient inside the gel has not been stable. 5. The DNA has undergone diffusion inside the gel, because of extremely prolonged electrophoresis time at low voltage. 6. The DNA has undergone diffusion inside the gel, because the DNA has not been fixed after running the gel. Only in parallel TGGE: 7. The band has migrated to a temperature which causes an irreversible transition of the DNA into single strands.	1. Load the correct sample volumes. 2. Check the amount of DNA. 3. Check the thermal coupling solution used. Use 0.1% Triton or Tween 20. 4. Check the volume of thermal coupling solution used. The gel should not change the position during the run. 5. Run the electrophoresis at 200 - 300 V (depending on the buffer used). 6. Incubate the gel directly after the run in Fixation solution (of the silver staining protocol, 10% EtOH/0.5% acetic acid). Silver stain. 7. Check the melting behavior of your DNA fragment in a perpendicular TGGE gel. If perpendicular TGGE shows a sigmoidal (S-shaped) curve, determine the effective range of separation. Set up new conditions for parallel TGGE. If parallel TGGE does not show a sigmoidal (S-shaped) curve, check all items listed under "No S-shaped curve in perpendicular TGGE".

Problem	Cause	Solutions
Band pattern is disturbed or distorted	1. Air bubbles in the gel. 2. Gel has been punctured (pipette tip) during loading of the sample. 3. The edges of the gel have dried out during electrophoresis. The front with the bands "smiles". 4. Wrong composition of gel and / or electrophoresis buffer. A "salt front" is in the gel. During electrophoresis run this salt front is indicated by an abnormal mobility of the marker dyes. The dye bands look extremely sharp, sometimes the bromophenolblue band moves at the same position as the xylene cyanol blue band. 5. High amounts of salt ions in the sample. Symptoms like described under 4.	1. Run a new gel without air bubbles. 2. Load the sample carefully. Don't touch the gel by the pipette tip. 3. Carefully protect the gel against evaporation. Cover the gel additionally with the special Cover glass plate (with 2 silicone barriers). 4. Check the composition of the gel and electrophoresis buffer. 5. Desalt the sample before loading.
No S-shaped curve in perpendicular TGGE (fig. 25)	1. Wrong buffer concentration, ionic strength in the gel is too high. DNA has been denatured. 2. Amount of urea in the gel is not sufficient. The DNA has not been denatured. 3. Ineffective renaturation of the DNA, only ssDNA has been loaded onto the gel. 4. DNA fragment is not stabilized by a GC-rich part of the sequence. Irreversible melting in a one-step transition into completely single-stranded DNA. 5. Unstable temperature gradient or no temperature gradient at all.	1. Check the composition of the gel and the electrophoresis buffer. 2. Check the amount of urea added to the gel. The standard protocol requires 8 M urea when dsDNA fragment (GC-content 55 - 75%) is analyzed on TGGE. 3. Check the denaturation / renaturation protocol used to form heteroduplexes. 4. Consider adding a stabilizing clamp to the fragment. Evaluate the optimized DNA fragment for TGGE analysis by calculating the melting pattern of the sequence with the POLAND program. 5. Check the gradient block. Purge any air bubbles from the gradient block.




Figure 25: No S-shaped curve in perpendicular TGGE

Problem	Cause	Solutions
PCR product derived from a putative heterozygous locus exhibits no heteroduplex bands in parallel TGGE (fig. 26 a and b).	1. Electrophoresis time is too short. 2. The temperature gradient range used in the experiment is not sufficient. The temperature at the cathode (-) is too high. The melting domain containing the point mutation has already been denatured, when the DNA enters the gel. 3. Masking of homoduplex bands. Only heteroduplex bands are visible in the gel. The DNA fragments have already passed the "effective range of separation". Homo- and heteroduplexes have been separated, but the fastest running homoduplex bands have also passed "T_{diss}", the temperature of total denaturation. Due to the irreversible transition into completely ssDNA, the homoduplex bands become diffuse, sometimes nearly invisible. 4. The sample loaded onto the gel only contains two different kinds of only homoduplexes with nearly identical T_m. 5. The point mutation is located in one of the most stable melting domains (parts) of the fragment. 6. Only one species of DNA has been loaded onto the gel. 7. Wrong composition of buffer and/or gel, unstable or nonexistent temperature gradient. No sigmoidal (S-shaped) curve at all in the	1. Run a perpendicular TGGE at first. Estimate the effective temperature range of separation and the running time required for the DNA fragment to reach this range in a parallel TGGE run. 2. See 1 3. See 1 4. Force heteroduplex formation by heating and reannealing the PCR sample prior to electrophoresis. 5. Calculate the melting map of your DNA fragment by the POLAND program. Construct a new fragment, which contains the site of mutation in one of the melting domains with the lowest T_m values. (See "Theoretical backbone of a detection rate approximately 100%".) 6. Check for the possibility that only one species of DNA was in your sample: ineffective or nonexistent PCR amplification of a particular allele. Loss of heterozygosity in a cell line, etc.... 7. Check all items listed under "No S-shaped curve in perpendicular TGGE".

	corresponding perpendicular TGGE.	
--	--------------------------------------	--




Figure 26 a: No heteroduplex bands in parallel TGGE




Figure 26 b: No heteroduplex bands in parallel TGGE.

Problem	Cause	Solutions
Gel visualizes more bands than expected (fig 27).	1. PCR artifacts: Nonspecific high molecular weight products are sometimes obtained by PCR amplification of genomic DNA samples. Due to the extremely sensitive silver-staining they become visible on a TGGE gel, although the PCR product seemed to be "clean" on an agarose gel stained by ethidium bromide. 2. ssDNA due to asymmetric PCR or ineffective renaturation during the formation of heteroduplexes is visualized on the gel by a band of orange to brown-red color. 3. DNA simply contains more species than expected. 4. The slot formers on the glass plates are damaged. "Ghost bands" are caused by these imperfect slot formers.	1. Recheck the sample by running a non-denaturing polyacrylamide gel or a perpendicular TGGE gel. Silver stain this gel. Hot start PCR. Adjust cycle parameters to finesse primer annealing. 2. Check on a perpendicular gel for ssDNA which does not show the sigmoidal inflection. Change the PCR conditions in order to achieve even amplification of both single strands and force heteroduplex formation by heating and reannealing the PCR sample prior to electrophoresis. 3. Check the number of sigmoidal (S-shaped) curves on a perpendicular TGGE. 4. Always use glass plates slot formers that are not damaged.




Figure 27: *More bands than expected.*

13 TGGE Testkit

13.1 Introduction

The TGGE Testkit contains samples for 4 parallel and 4 perpendicular TGGE runs with Na-TAE buffer. The samples contain a wild type mixture with one DNA double strand, a mutant mixture with one DNA double strand and a heteroduplex mixture with four different DNA double strands: two homoduplices and two heteroduplices.

Wild type: min. 20 µl
 Mutant: min. 20 µl
 Heteroduplex: min. 120 µl
 Loading buffer: min. 180 µl (with blue marker for Na-TAE buffer)

13.2 Protocol

13.2.1 Gel composition:

8% PAA
 8 M Urea
 0.2 x Na-TAE
 2% Glycerol

Recipe for 10 ml gel solution (3 – 4 gels) for Na-TAE running buffer:

	8% Gel
Urea ($C_{End} = 8\text{ M}$)	4.8 g
Acrylamide/bis Acrylamide stock solution (30 : 0,8), 40% (w/v)	2.0 ml
10x conc. Na-TAE, pH 8.4 ($C_{End} = 0.2 \times \text{conc.}$)	0.2 ml
40% Glycerol ($C_{End} = 2\%$)	0.5 ml
Water, distilled	2.5 ml
Make sure that the urea has been completely resolved. It is possible to heat up the urea containing solution slightly (40°C – 50°C) for a <u>short time</u> in order to improve the solubilization of urea. De-gas the solution under gentle vacuum for 3 - 5 min.	
Water, distilled	fill up to 10 ml
TEMED	14 µl
APS (4%)	45 µl
Mix gently. Avoid air bubbles!	
Pour the gel solution into the glass plate sandwich immediately thereafter without air bubbles.	

13.2.2 Running buffer:

0.2 x Na-TAE

Buffer composition:

Na-TAE, pH 8.4, 1 M Sodium acetate
Stock solution (10 x conc.) 10 mM EDTA
 400 mM TRIS
 pH = 8.4 (titrate with Acetic Acid, never use HCl!)

Na-TAE Running Buffer 0.2x conc. Na-TAE, pH 8.4

TBE buffer: (not recommended for the TGGE Testkit)

Using TBE as running buffer results in less sharp bands, longer running time and lower melting temperatures!

Perpendicular gel:

pre-run: 4 min + run: 40 min.

Using the 30°C - 70°C temperature gradient allows no optimization of parameters for parallel TGGE as the melting curve starts just before the DNA is completely molten into single strands.

Parallel gel:

pre-run: 4 min + run: 10 min.

Using the 40°C - 60°C temperature gradient gives no resolution, as all double strands are molten and only single strands exist. The temperature gradient has to be reduced to 20°C - 50°C. 20 min running time necessary.

13.2.3 Electrophoresis parameters:

perpendicular gel (figure 28):

Sample volume: 20µl Heteroduplex (Hd)
+ 5µl load. buffer + 25µl run. buffer
Temperature range (T_1 - T_2): 30°C - 70°C
Pre-run: 4 min., 250 V, 20°C
Run: 40 min., 250 V, 30°C - 70°C

parallel gel (figure 29):

Sample volume: 2 µl sample
+ 0.5 µl load. buffer + 2.5 µ run. buffer
Temperature range (T_1 - T_2): 40°C - 60°C
Pre-run: 4 min., 300 V, 25°C
Run: 10 min., 300 V, 40°C - 60°C

13.3 Gel images

This is what you should see on your gels if you proceed according to the protocol:




Figure 28: perpendicular gel




Figure 29: parallel gel

M = Marker
Wt = Wild type
Mt = Mutant
Hd = Heteroduplex

14 Appendix

14.1 Technical Data

14.1.1 System

Working temperature	4°C – 35°C
Humidity	10 % - 80%

14.1.2 Electrophoresis Chamber

Temperature range of gradient block	5°C – 80°C
Maximum linear temperature range (above 20°C, room temperature)	45 K, for example 30°C – 75°C or 25°C – 70°C
Start of linear gradient (parallel)	first marked line
Start of linear gradient (perpendicular)	end point of each line
Temperature Accuracy	± 0.3°C
Temperature Uniformity	± 0.3°C/ 2 mm
Gradient block size	6 cm x 6 cm
Gel size (W x L)	7.4 cm x 8.0 cm
Run distance	6.2 cm (parallel), 6.2 cm (perpendicular)
Size (LxWxH)	22.5 cm x 22.5 cm x 23 cm
Weight	4.2 kg

14.1.3 TGGE System Controller with integrated power pack

Program stores	100
Display	LCD
Languages	German, English
Mains voltage	115V / 230 V
Frequency	50-60 Hz
Wattage	max. 30 VA
Fuses	2x 3.15 AT (115 V) / 2x 1.6 AT (230 V)
Interfaces	1 parallel port (Centronics for printer) 1 serial port (RS232)
Size (LxWxH)	31 cm x 22 cm x 11.5 cm
Weight	3.8 kg
Integrated power pack	
Voltage	max. 400 V
Current	max. 500 mA
Power	max. 30 W

14.2 Buffers

Loading buffers:

Loading buffer TBE:	0.1x conc. TBE (up to 1x conc. TBE is possible) 0.1% Triton-X 100 0.01% Bromophenol Blue dye 0.01% Xylene Cyanol dye
Loading buffer Na-TAE:	0.2x conc. Na-TAE, pH 8.4 0.1% Triton-X 100 0.01% Bromophenol Blue dye 0.01% Xylene Cyanol dye
Loading buffer ME:	10 x conc. MOPS 10 mM EDTA 0.05% Bromophenol Blue dye 0.05% Xylene Cyanol dye pH = 8.0

Denaturation / Renaturation (DR) Loading buffers:

DR Loading buffer TBE:	0.1x conc. TBE (up to 1x conc. TBE is possible) 0.1% Triton-X 100 0.01% Bromophenol Blue dye 0.01% Xylene Cyanol dye 7 M Urea
DR Loading buffer Na-TAE:	0.2x conc. Na-TAE, pH 8.4 0.1% Triton-X 100 0.01% Bromophenol Blue dye 0.01% Xylene Cyanol dye 8 M Urea
DR Loading buffer ME:	20 x conc. MOPS 20 mM EDTA 0.01% Bromophenol Blue dye 0.01% Xylene Cyanol dye 8 M Urea pH = 8.0

Running buffers:

TBE	890 mM Boric Acid
Stock solution (10 x conc.)	20 mM EDTA 890 mM TRIS Do not titrate to adjust pH!
TBE Running buffer	0.1 x conc. TBE (up to 1x conc. TBE is possible)
Na-TAE, pH 8.4,	1 M Sodium acetate
Stock solution (10 x conc.)	10 mM EDTA 400 mM TRIS pH = 8.4 (titrate with Acetic Acid, never use HCl!)
Na-TAE Running Buffer	0.2x conc. Na-TAE, pH 8.4
ME (MOPS/EDTA)	1 M MOPS
Stock solution (50 x conc.)	50 mM EDTA pH = 8.0
ME-Running Buffer	1 x conc. ME, pH = 8.0

SSCP buffer:

SSCP-ME buffer	1 M MOPS
Stock solution (50 x conc)	250 mM EDTA (Free Acid) pH = 8.0
SSCP-ME Running buffer	1 x conc. SSCP-ME, pH = 8.0

Others:

TE buffer	10 mM Tris/HCl 0.1 mM EDTA pH = 8.0
TEMED	Solution of N,N,N',N'tetramethylethylenediamine
APS	4% Ammonium persulfate
Glycerol 40%	40% glycerol in water
Glycerol 50%	50% glycerol in water

14.3 Silver staining solutions:

Standard method:

Fixation	10% EtOH 0.5% Acetic Acid 100 ml ethanol and 5 ml acetic acid are adjusted with distilled water to 1 liter. Prepare freshly !
Silver Binding	0.19% AgNO ₃ 1.9g AgNO ₃ is dissolved in 1 liter of distilled water. (Can be reused for 5 gels.) Store dark!
Developing Solution	1.5% NaOH 0.08% NaBH ₄ 0.1% Formaldehyde Dissolve 15 g NaOH in 1 liter distilled water and add 0.8g NaBH ₄ . Immediately before developing add 2.7 ml formaldehyde stock solution (37% in water). This solution must be prepared fresh every time!
Stopping Solution	0.75% Na ₂ CO ₃ Dissolve 7.5 g sodium carbonate in ddH ₂ O. Total volume: 1 liter

Quick method (for PCR products):

Fixation:	10% EtOH 0.5% Glacial Acid 100 ml ethanol and 5 ml acetic acid are adjusted with double distilled water to 1 liter.
Silver Binding	0.2% AgNO ₃ 2.0 g AgNO ₃ is dissolved in 1 liter of distilled water. (Can be reused for 5 gels.) Store dark!
Developing Solution	3.0 % NaOH 0.5% Formaldehyde Dissolve 3 g NaOH and 1.35 ml formaldehyde stock solution (37% in water) in 100 ml double distilled water. This solution must be prepared fresh every time!
Stopping Solution:	identical with Fixation solution (10% EtOH, 0.5% Glacial Acid)

**Quick method using the AMRESCO SilverPAGE™ staining kit
(Code No. 211-761)**

Fixation:	30% EtOH 10% Acetic Acid 300 ml ethanol and 100 ml acetic acid are adjusted with double distilled water to 1 liter.
Sensibilisation:	30% EtOH Prepare freshly 60 ml ethanol in 140 ml double distilled water.
Silver Binding:	Prepare Silver Binding Agent by reconstituting contents of one pouch in 1 l of ddH ₂ O. <i>(This solution must be prepared fresh every time!)</i> Immediately before staining , add 0.7 ml of 37% Formaldehyde to 200 ml of reconstituted Silver Binding Agent.
Developing Solution:	Just prior to use, prepare developing solution by reconstituting contents of one pouch of Developer I and 15 mg of Developer II in 200 ml of ddH ₂ O. <i>(This solution must be prepared fresh every time!)</i> Immediately before developing , add 0.7 ml of 37% Formaldehyde to 200 ml of reconstituted developing solution.
Stopping Solution:	7.5% Acetic Acid 75 ml acetic acid are adjusted with double distilled water to 1 l.

15 References

1. Riesner, D., Henco, K. and Steger, G. (1990): Temperature-Gradient Gel Electrophoresis: A method for the analysis of conformational transitions and mutations in nucleic acids and protein. Page 169-250 In Chrambach, A., Dunn, M.J., Radola, B.J.: *Advances in Electrophoresis*, Vol. 4, VCH Verlagsgesellschaft Weinheim
2. Kappes, S. et al.(1995): p53 mutations in ovarian tumors, detected by temperature-gradient gel electrophoresis, direct sequencing and immunohistochemistry. *Int. J. Cancer* 64: 52-59
3. Milde-Langosch, K. et al. (1995): Presence and persistence of HPV and p53 mutation in cancer of the cervix uteri and the vulva. *Int. J. Cancer* 63: 639-645
4. Horn, D. et al.(1996): Three novel mutations of the NF1 gene detected by temperature gradient gel electrophoresis of exons 5 and 8. *Electrophoresis* 17: 1559-1563
5. Wieland, U. et al.(1996): Quantification of HIV-1 proviral DNA and analysis of genomic diversity by polymerase chain reaction and temperature gradient gel electrophoresis. *J. Virology Methods* 57: 127-139
6. Kuhn, J.E. et al. (1995): Quantitation of human cytomegalovirus genomes in the brain of AIDS patients. *Journal of Medical Virology* 47: 70-82
7. Linke, B. et. al. (1995): Identification and structural analysis of rearranged immunoglobulin heavy chain genes in lymphomas and leukemia. *Leukemia* 9: 840-847.
8. Menke M.A. et al. (1995): Temperature gradient gel electrophoresis for analysis of a polymerase chain reaction-based diagnostic clonality assay in the early stages of cutaneous T-cell lymphomas.
9. Hecker, R. et al. (1988): Analysis of RNA structure by temperature-gradient gel electrophoresis: viroid replication and processing. *Gene* 72: 59-74
10. Baumstark, T. and Riesner, D. (1995): Only one of four possible secondary structures of the central conserved region of potato spindle tuber viroid is a substrate for processing in a potato nuclear extract. *Nucleic Acids Research* 23: 4246-4254
11. Loss, P., Schmitz, M., Steger, G. and Riesner, D. (1991): Formation of a thermodynamically metastable structure containing hairpin II is critical for the potato spindle tuber viroid. *EMBO Journal* 10: 719-728
12. Riesner, D. (1998): Nucleic acid structures. In: *Antisense Technology. Practical Approach Series*. Oxford University Press. p1-24 (in press)

13. Wiese, U. et al. (1995): Scanning for mutations in the human prion protein open reading frame by temporal temperature gradient gel electrophoresis. *Electrophoresis* 16: 1851-1860
14. Nubel, U. et al. (1996): Sequence heterogenities of genes encoding 16S rRNAs in *Paenibacillus polymyxa* detected by temperature gradient gel electrophoresis.
15. Lessa, E.P. and Applebaum, G. (1993): Screening techniques for detecting allelic variation in DNA sequences. *Molecular Ecology* 2: 119-129
16. Richter, A., Plobner, L., Schumacher, J. 1997: Quantitatives PCR-Verfahren zur Bestimmung der Plasmidkopienzahl in rekombinanten Expressionssystemen. *BIOforum* 20: 545-547
17. Henco, K. and Heibey, M. (1990): Quantitative PCR – the determination of template copy numbers by temperature gradient gel electrophoresis. *Nucleic Acids Research* 18: 6733-6734
18. Birmes, A. et al. (1990): Analysis of the conformational transition of proteins by temperature-gradient gel electrophoresis. *Electrophoresis* 11: 795-801
19. Arakawa, T. et al. (1993): Analysis of the heat-induced denaturation of proteins using temperature gradient gel electrophoresis. *Analytical Biochemistry* 208: 255-259
20. Chen, X. et al. (1995): High resolution SSCP by optimization of the temperature by transverse TGGE. *Nucleic Acids Research* 23: 4524-4525
21. Scholz, R.B. et al. (1993): Rapid screening for Tp53 mutations by temperature gradient gel electrophoresis: a comparison with SSCP analysis. *Human Molecular Genetics* 2: 2155-2158
22. Elphinstone and Baverstock, P.R. (1997): Detecting mitochondrial genotypes by temperature gradient gel electrophoresis and heteroduplex analysis. *BioTechniques* 23: 982-986
23. Poland, D. (1974): Recursion relation generation of probability profiles for sequence-specific macromolecules with long-range correlations. *Biopolymers* 13:1859-1871
24. Lerman, L.S. and Silverstein, K. (1987): Computational simulation of DNA-melting and its application to denaturing gradient-gel electrophoresis. *Meth. Enzymol.* 155: 482-501
25. Steger, G. (1994): Thermal denaturation of double-stranded nucleic acids: prediction of temperatures critical for gradient electrophoresis and polymerase chain reaction.

-
26. Schumacher, J, Randels, J.W. and Riesner,D. (1983): A two dimensional electrophoretic technique for detection of circular viroids and virusoids. *Anal. Biochem.* 135, 288 - 295
 27. Sambrook, J., Fritsch, E.F. and Maniatis,T. (1989): Molecular cloning, Cold Spring Habor Laboratory press
 28. Steger, G. and Riesner, D. (1992): Temperaturgradienten-Gelelektrophorese: eine Methode zur Analyse von Konformationsübergängen und Mutationen in Nukleinsäuren und Proteinen. In Radola, B.J. (ed) *Handbuch der elektrophorese*, VCH Verlagsgesellschaft, Weinheim
 29. Sheffield, V.C., Cox, D.R. and Lerman, R.M. (1989): Attachment of a 40-base-pair G+C-rich sequence (GC-clamp) to genomic DNA fragments by the polymerase chain reactiob results in improved detection of single-base changes. *Proc. Natl. Acad. Sci. USA* 86, 232 - 236
 30. Hecker R., Wang Z., Steger G. and Riesner D. (1988): Analysis of RNA structure by temperature-gradient gel electrophoresis: viroid replication and processing. *Gene* 72, 59-74
 31. Jiang L., Chen W., Tain L.P. and Liu Y. (1991): Temperature-gradient gel electrophoresis of apple scar skin viroid. *Acta Microbiol. Sin.* 30, 278-283
 32. Riesner D., Hecker R. and Steger G. (1988): Structure of viroid replication intermediates as studied by thermodynamics and temperature-gradient gel electrophoresis. In Sarma R.H. and Sarma M.H. (eds.) *Structure & Expression*, Vol. I: From Proteins to Ribosomes, Adenine press, 261-285
 33. Riesner D., Steger G., Zimmat R., Owens R.A., Wagenhöfer M., Hillen W., Vollbach S. and Henco K. (1989): Temperature-gradient gel electrophoresis of nuleic acids: Analysis of confor-mational transitions, sequence variations, and protein-nucleic acid interactions. *Electrophoresis* 10, 377-389
 34. Rosenbaum V. and Riesner D. (1987): Temperature-gradient gel electrophoresis: thermodynamic analysis of nucleic acids and proteins in purified form and in cellular extract. *Biohys. Chem.* 26, 235-246
 35. Schönborn J., Oberstraß J., Breyel E., Tittgen J., Schumacher J., Lukacs N. (1991): Monoclonall antibodies to double-stranded RNA as probes of RNA structure in crude nucleic acid extracts. *Nucleic Acids Res.* 19, 2993-3000
 36. Po Tien, Steger G., Rosenbaum V., Kaper J. and Riesner D. (1987): Double-stranded cucumovirus associated RNA5: experimental analysis of necrogenic and non-necrogenic variants by temperature-gradient gel electrophoresis. *Nucleic Acids Res.* 15, 5069-5083

37. Zimmat R., Gruner R., Hecker R., Steger G. and Riesner D. (1991): Analysis of mutations in viroid RNA by non-denaturing and temperature-gradient gel electrophoresis. In R.H. Sarma and M.H. Sarma (eds.) Structure & Methods, Vol. 3:, DNA & RNA, Adenine Press, 339-357
38. Rosenbaum V., Klahn T., Lundberg Holmgren E., von Gabain A. and Riesner D. (1992): Co-existing structures of an mRNA stability determinat: The 5' region of the Escherichia coli and Serratia marcescens ompA mRNA, J.Mol.Biol., in press
39. Birmes A., Sättler A., Maurer S.O. and Riesner D. (1990): „Analysis of the conformational transitions of proteins by temperature-gradient gel electrophoresis“. Electrophoresis 11, 795-801
40. Sättler A., Kanka S., Schrörs W. and Riesner D. (1992): „Random mutagenesis of the weak calcium binding side in Subtilisin Carlsberg and screening for thermal stability by temperature-gradient gel electrophoresis“. Accepted for: 1st International Symposium of Subtilisin Enzymes, EMBL, Hamburg
41. Thatcher D. and Hodson B. (1981): „Denaturation of proteins and nucleic acids by thermal-gradient electrophoresis“. Biochem. J. 197, 105-109
42. Wagenhöfer M., Hansen D. and Hillen W. (1988): „Thermal denaturation of engineered tet repressor proteins and their complexes with tet operator and tetracycline studie by temperature-gradient gel electrophoresis“. Analytical Biochem. 175, 422-432
43. Sanguinetti C.J., Neto E.D. and Simpson A.J.G. (1994): BioTechniques 17, 915
44. Kappes, S., Milde-Langosch, K., Kressin, P., Passlack, B., Dockhorn-Dwornczak, B., Röhlke, P. and Löning, T. (1995): "p53 Mutations in ovarian tumors, detected by temperature-gradient gel electrophoresis, direct sequencing and immunohistochemistry". Int. J. Cancer 64, 52 - 59
45. Kluwe, L., MacCollin, M., Tatagiba, M., Thomas, S., Hazim, W., Haase, W. and Mautner, V.-F. (1998): "Phenotypic variability associated with 14 splice-site mutations in the NF2 gene". American Journal of Medical Genetics 77, 228 - 233
46. Lerman, L. S. and Beldjord, C. (1998). "Comprehensive mutation detection with denaturing gradient gel electrophoresis". In R.G.H. Cotton, E. Edkins and S. Forrest (eds.) Mutation Detection, A Practical Approach, Oxford University Press, 35 - 62
47. Gamper, H., Piette, J. and Hearst, J.E. (1984): Photochem. Photobiol. 40, 29 ff

16 Order information and spare parts

Biometra offers a wide range of accessories parts and consumables related to Temperature Gradient Gel Electrophoresis.

Item	Order No.
TGGE System ; 230 V, electrophoresis unit with Peltier element-powered gradient block and 2 removable electrophoresis buffer chambers, system controller with integrated power supply, Starter Kit and manual	024-000
TGGE System , 115 V, dito	024-090
TGGE System Controller with integrated power supply, 230 V	024-001
TGGE System Controller with integrated power supply, 115 V	024-091
TGGE electrophoresis unit with Peltier element powered gradient block, 2 removable electrophoresis buffer chambers and connector cable	024-002
TGGE connector cable (controller to electrophoresis unit)	024-033
TGGE Starter Kit with 3 Bonding glass plates, 3 types of glass plates with slots, pre-cut electrode wicks, pre-cut Polybond film and cover film	024-003
TGGE Testkit , samples for 4 parallel and 4 perpendicular TGGE runs (20 µl wild type, 20 µl mutant, 120 µl Heteroduplex, 180 µl loading buffer with blue-marker for Na-TAE buffer)	024-050
Accessories	
TGGE removable electrophoresis buffer chambers, 2 pcs	024-010
TGGE electrode wicks, pre-cut 8 x 7 cm, 100 pcs	024-015
TGGE bonding plate, 9 x 9 cm, w/o spacer	024-021
TGGE glas plate 9 x 9 cm, 8 slots 4x3x0.4 mm, approx. 5 µl	024-022
TGGE glas plate 9 x 9 cm, 1 slot (rectangular) 40x3x0.4 mm, approx. 50 µl	024-023
TGGE glas plate 9 x 9 cm, 1 slot (diagonal) 62x3x0.4 mm, approx. 75 µl	024-024
TGGE glas plate 9 x 9 cm, 12 slots 3x2x0.4 mm, approx. 3 µl	024-025
TGGE glas plate 9 x 9 cm, 18 slots 2x2x0.4 mm, approx. 2 µl	024-026
TGGE glas plate 9 x 9 cm, 0.5 mm spacer, no slots	024-027
TGGE Polybond film, pre-cut 8.8 x 8.8 cm, 25/pkg	024-030
TGGE Cover glass plate with silicone barriers + 10 cover films, pre-cut 7 x 6 cm	024-031
TGGE cover film, pre-cut 7 x 6 cm 25/pkg	024-032
TGGE Polybond film, pre-cut 8.8 x 8.8 cm, 100/pkg	024-034
TGGE cover film, pre-cut 7 x 6 cm 100/pkg	024-035
TGGE Slotformer („Slot forming units“), 10 x „multi“, 9 x „long“	024-121
Gel casting clips, 3/pkg	010-007
Consumables	
Acrylamide/bis Acrylamide, 40% (30:0.8), 500 ml	210-254
Ammonium Persulfate, 4x 25 g	210-486
EDTA (Tetrasodium Salt, Dihydrate), 500 g	210-245
Glycerol, 1 l	210-854
Sodium Acetate, anhydrous, 500 g	210-602
TEMED, 25 ml	210-761
Acryl-Glide, 100 ml	211-319
Silver PAGE™, Silver staining kit for 20 stains	211-761
TRIS, 1 kg	220-826

Consumables from Biometra are not available outside Germany. Please source from another supplier.

17 Instructions for return shipment



If you would like to send the unit back to us, please read the following return instructions.

Should you have any problems with the TGGE System, please contact your local **Biometra** dealer or our service department:

Biometra biomedizinische Analytik GmbH
Service Department
Rudolf-Wissell-Straße 30
D-37079 Göttingen
Phone: ++49 – (0)5 51 / 50 68 6-0
Fax: ++49 – (0)5 51 / 50 68 6-66

- Return only defective devices. For technical problems which are not definitively recognisable as device faults please contact the Technical Service Department at **Biometra**.
 - Use the original box or a similarly sturdy one.
 - Label the outside of the box with “CAUTION! SENSITIVE INSTRUMENT!”
 - Please enclose a **precise description of the fault**, which also reveals during which procedures the fault occurred, if possible.
 - **Important:** Clean all parts of the instrument from residues, and of biologically dangerous, chemical and radioactive contaminants. Please **include a written confirmation** (use the “**Equipment Decontamination Declaration**” following on the next page) that the device is free of biologically dangerous and chemical or radioactive contaminants in each shipment. If the device is contaminated, it is possible that **Biometra** will be forced to refuse to accept the device.
- 
- The **sender of the repair order will be held liable** for possible losses resulting from insufficient decontamination of the device.
 - Please enclose a note which contains the following:
 - a) Sender's name and address,
 - b) Name of a contact person for further inquiries with telephone number.

18 Equipment Decontamination Certificate

To enable us to comply with german law (i.e. §28 StrlSchV, §17 GefStoffV and §19 ChemG) and to avoid exposure to hazardous materials during handling or repair, will you please complete this form, prior to the equipment leaving your laboratory

COMPANY / INSTITUTE _____

ADDRESS _____

TEL NO _____ FAX NO _____

E-MAIL _____

EQUIPMENT

Model

Serial No

If on loan / evaluation

Start Date: _____ Finish Date _____

Hazardous materials used with this equipment

Has the equipment been cleaned and decontaminated? YES / NO (delete)

Method of cleaning / decontamination

NAME _____

POSITION _____

(HEAD OF DIV./ DEP./ INSTITUTE / COMPANY)

SIGNED _____

DATE _____

PLEASE RETURN THIS FORM TO BIOMETRA GMBH OR YOUR LOCAL BIOMETRA DISTRIBUTOR TOGETHER WITH THE EQUIPMENT.

PLEASE ATTACH THIS CERTIFICATE OUTSIDE THE PACKAGING. INSTRUMENTS WITHOUT THIS CERTIFICATE ATTACHED WILL BE RETURNED TO SENDER.

19 Warranty

This Biometra instrument has been carefully built, inspected and quality controlled before dispatch. Hereby Biometra warrants that this instrument conforms to the specifications given in this manual. This warranty covers defects in materials or workmanship for 12 month as described under the following conditions:

This warranty is valid for **12 month** from date of shipment to the customer from Biometra or an authorized distributor. This warranty will not be extended to a third party without a written agreement of Biometra.

This warranty covers only the instrument and all original accessories delivered with the instrument. This warranty is valid only if the instrument is operated as described in the manual.

Biometra will repair or replace each part which is returned and found to be defective.

This warranty does not apply to wear from normal use, failure to follow operating instructions, negligence or to parts altered or abused.