

# 40min Fast Plasmid Maxiprep Kit - Syringe Format

## Description

The Biotool endofree system uses a specially formulated buffer that extracts the endotoxin from bacterial lysate. The **final endotoxin level is < 0.3 EU/μg plasmid DNA**.

Purified endofree DNA is ready for high performance downstream applications, such as transfection of robust cells, endotoxin-sensitive cell lines, primary cultured cells or embryo microinjection.

## Components

Contents	Cat#:B22213	Cat#:B22216
Preps	25	100
BtBind DNA Columns	25	100
Filter Syringe (60 mL)	25	100
2.0 mL Microfuge tube	50	200
Plastic wrench	2	4
Buffer A1	270 mL	270 mL×4
RNase A (20 mg/mL)	27 mg (1.35 mL)	27 mg (1.35 mL) ×4
Buffer B1	270 mL	270 mL×4
Buffer N3	80 mL	80 mL×4
Buffer RET	250 mL×2	250 mL×8
DNA Wash Buffer	2×54 mL	10×54 mL
Endofree Elution Buffer	60 mL	60 mL×4
User Manual	1	1×4

### Notes for handling buffers

#### Buffer A1 & RNase A:

Spin down and add RNase A into Buffer A1, mix well and **store at 4°C**.

#### Buffer B1:

Buffer B1 precipitates below 10°C. **Warm up Buffer B1 at 50°C if precipitates before use. Keep the cap of Buffer B1 tightly closed after use.**

#### Buffer N3:

Buffer N3 may precipitate below 10°C. **Warm up Buffer N3 at 37°C if precipitates before use.**

#### DNA Wash Buffer:

**Add 216 mL 100% ethanol according to bottle label before use and mix well.**

#### Endofree Elution Buffer:

**Pre-warm the Elution Buffer at 60°C to increase plasmid DNA yield.**

## Storage

- **After added RNase A, Buffer A1 should be stored at 4°C.**
- All other materials can be stored at room temperature (20-26°C).
- The guaranteed **shelf life is 12 months** from the date of purchase.

## Notice

Below is **Steps where mistakes are easily made**.  
Please read the following carefully!

### Step 5:

When transferred into **Filter Syringe**, **let stand the lysate at room temperature for 1 to 2 min** until the white precipitates float to the top, then press the plunger.

- Optional: **Spin down the lysate at 3500 rpm for 1 to 2 min** to remove most of the precipitates, and then transfer the supernatant into the **Filter Syringe**.

### Step 6:

- **Buffer RET** is recommended to be pre-warmed at **50°C** to dissolve salt precipitates if any exists.
- The solution mixture **MUST be immediately filtered** to avoid salt precipitation after adding **Buffer RET** and **100% ethanol**.

### Step 7:

- Use a **vacuum pump** to suck the solution mixture through the column to absorb the plasmid DNA onto the membrane.
- Optional: Transfer the solution mixture into the **BtBind DNA** column and press the plunger to expel the solution to bind DNA onto the membrane.

### Step 11:

- **Endofree Elution Buffer or ddH<sub>2</sub>O must** be pre-warmed at **60°C** to effectively elute the plasmid DNA.

## Protocol

The kit is optimized for high copy number plasmid purification. Scale up both the volume of bacterial culture and buffers for low copy number plasmids.

The yield of plasmid DNA depends on the origin of the replication and the size of the plasmid.

Ensure **all the buffers are handled as described above (critical for efficient plasmid purification)** and provide materials listed below:

- 100% ethanol.
- High speed centrifuge (Equipped for 1.5 mL or 2 mL EP tubes).




#### Order & Inquiry

Tel: (713)732-2181 Fax: +1-866-747-4781  
E-mail: [order@biotool.com](mailto:order@biotool.com)



#### Order & Inquiry

Tel: +49-89-46148500 Fax: +49-89-461485022  
E-mail: [eu.order@biotool.com](mailto:eu.order@biotool.com)



**figure 1. Flow chart of protocol**

1. Harvest **150-250 mL** fresh bacterial **overnight culture** (at 37°C for 14-16 h by vigorous shaking with initial inoculation ratio of 1:500 to 1: 1000) by centrifugation at 5,000 g for 10 min at room temperature. Pour off the supernatant and blot the inverted tube on paper towels to completely remove residual medium.

**Note:**

- **Complete removal of residual medium is critical for bacteria lysis in the next steps.**
- **Do not use more than 250 mL culture.** If processing > 250 mL of bacterial culture, the buffers **MUST** be scaled up.

2. Add **10 mL Buffer A1** and **completely resuspend** bacterial pellet by vortexing or pipetting. (Ensure that **RNase A** has been added into Buffer A1 before use).

**Note:** Complete resuspension is critical for optimal yield.

3. Add **9 mL Buffer B1**, **mix thoroughly** by inverting 10 times with **gentle shaking**. **Incubate for 5 to 10 min at room temperature** to obtain a clear lysate. The mixture of completely lysed bacteria will appear transparent.

**Note:**

- Buffer B1 precipitates below 10°C, **warm up at 50°C** to dissolve **if precipitates**.
- **Shake gently** and **don't incubate over 10 min** to avoid genomic DNA fragmentation which greatly hampers plasmid DNA purity.

4. Add **3 mL Buffer N3**, **mix completely** by inverting the tube 10 times then **shaking vigorously** for 5 times. It is critical to mix the solution well. More shaking is required to completely neutralize the solution until the mixture appears congealed, brownish or viscous.

5. Set the **Filter Syringe** upon a new 50 mL conical tube and **pour the lysate into the barrel of the filter syringe**. **Let stand at room temperature for 1 to 2 min** until the white precipitates float to the top. Insert the plunger and expel the clear lysate (DNA remains in the clear lysate) into the tube. Stop when the syringe encounters strong resistance. Flip the syringe upside down and the precipitate will float to the top. Expel the volume again to get maximal clear lysate. Small amounts of lysate may remain in the precipitate and can be ignored.

Optional (Low speed centrifuge): After adding buffer N3, **spin down the lysate mixture at 3500 rpm for 1 to 2 min** in a 50 mL conical tube to remove most of the precipitate. The DNA will remain in the supernatant, which is transferred to the filter syringe. The plunger can then be pushed more easily to filter the lysate.

6. Add **0.7-1.0 volume Buffer RET** (e.g. **14 to 20 mL of Buffer RET to 20 mL** of clear lysate), and **10 mL 100% ethanol** into the lysate. Mix well by **vigorous shaking** for 5 times. The solution mixture **MUST be immediately filtered** through the **BtBind DNA Columns**.

**Note:** Immediate filtration is critical to avoid salt precipitation.

7. Use a vacuum pump to suck the lysate/ethanol mixture through the column to absorb the plasmid DNA onto the membrane.  
Optional: Set the assembled DNA column in an upright position upon a 50 mL conical tube (using the Plastic Wrench as the column holder). Add the lysate/ethanol mixture into a DNA column and use the plunger to expel the mixture through the column. Take the plunger out of the DNA column and add the remaining mixture. Expel the plunger until all of the mixture goes through the DNA binding membrane.

**Note:**

There is no need for disassembling the membrane column before taking the plunger out of the DNA column.



Order & Inquiry

Tel: (713)732-2181 Fax: +1-866-747-4781  
E-mail: order@biotool.com



Order & Inquiry

Tel: +49-89-46148500 Fax: +49-89-461485022  
E-mail: eu.order@biotool.com

8. Add **10 mL DNA Wash Buffer** into DNA column. Hold the plunger and invert the syringe upside down several times to wash the cavity. Expel wash buffer using the plunger or a vacuum pump.

9. Repeat Step 8.

10. Use the plastic wrench to detach the membrane from the DNA column and insert the membrane into a **2.0 mL Microfuge tube**. Spin the column at **~12,000 rpm for 1 min**, discard the flow and spin again at **~12,000 rpm for 2 min**. Air-dry (air-blast in fume hood preferred) the membrane at room temperature for 3-5 min.

11. Insert the membrane into a new sterile 1.5 mL or 2.0 mL **Microfuge tube**. Add **500 µL Endofree Elution Buffer or ddH<sub>2</sub>O** (Pre-warmed at 60°C) to the center of the membrane and incubate for 1 min at room temperature. Centrifuge at **~12,000 rpm for 1 min** to elute the DNA. Add the DNA eluate back onto the membrane for a second elution.

**Note:**

The first elution normally yields about 60-70% of the DNA and the second elution yields another 30% of the DNA.

**Optional:**

To obtain the maximum recovery, add another **300 µL Endofree Elution Buffer or ddH<sub>2</sub>O** to the center of the membrane and incubate for 1 min at room temperature. Elute the DNA by centrifugation at **12,000 rpm for 1 min**. The final DNA concentration will decrease with the increased volume.

## Trouble shooting

Problem	Possible Reason	Suggested Improvement
No DNA	Plasmid lost in Host E. coli	Prepare fresh culture.
Low Yield	Bacterial culture overgrown or not fresh	Grow bacterial 12-16 hours. Spin down culture and store at -20°C until use.
Low Yield	Low copy-number plasmid	Scale up culture volume and the volume of Buffers A1, B1, N3 and 100% ethanol proportionally.
Low Yield	Poor Cell lysis	<ul style="list-style-type: none"> <li>Resuspend pellet thoroughly by vortexing and pipetting prior adding Buffer B1.</li> <li>Make fresh Buffer B1 if the cap had not been closed tightly.</li> </ul>
Genomic DNA contamination	Improper manipulation after adding Buffer B1	<ul style="list-style-type: none"> <li>Don't vortex or mix aggressively.</li> <li>Don't incubate over 10 minutes.</li> </ul>
OD <sub>260/280</sub> >2.0	RNA contamination	Add RNase A into Buffer A1 and Store at 4°C.
Low ratio of OD <sub>260/230</sub>	Presence of contaminants (EDTA, phenol, salt etc.) which absorbs at 230 nm	<ul style="list-style-type: none"> <li>Pre-warm Buffer RET at 50°C for 15 min to avoid salt precipitation.</li> <li>Increase washing times.</li> <li>Air-dry the membrane longer or in the ventilation system.</li> </ul>
Plasmid DNA floats out of wells while running in agarose gel, DNA doesn't freeze or smell of ethanol	Ethanol traces were not completely removed from column	Make sure that no ethanol remains in the silicon membrane before eluting the plasmid DNA. Re-centrifuge or vacuum again if necessary.



**Order & Inquiry**

Tel: (713)732-2181 Fax: +1-866-747-4781  
 E-mail: order@biotool.com



**Order & Inquiry**

Tel: +49-89-46148500 Fax: +49-89-461485022  
 E-mail: eu.order@biotool.com