


For best results please review these guidelines.

Printing Requirements

- Have a volumetric DICOM data set.
- THE PRIMARY/AXIAL/ORIGINAL SERIES/SLICES IS REQUIRED TO PRINT.
- Archive the entire study in uncompressed DICOM format with one slice per file.

Technical Guidelines

- The accuracy of the patient-specific anatomical model is dependent upon the quality of the CT scan or ConeBeam scan.
- Motion artifact will affect the accuracy of the anatomical model. The file quality will be poor if it contains motion artifacts.
- If scanning, use a scanning algorithm that produces high-resolution images suitable for 3D reconstruction or stereotactic planning (consult with your scanner manufacturer to use the correct settings).
- High spatial resolution scans are critical. The series should be acquired in continuous slices at a thickness of 0.5mm – 1.0mm. Slice spacing and thickness should be the same.
- The scanning FOV must contain the entire region of interest with a 5mm margin.
- Please do not use gantry tilt. Position the patient with the occlusal plane parallel to the image plane.
- For CBCT, use a long scan time and a Voxel Size of .3 to .4mm with the largest FOV.



Getting Started uploading the DICOM data via Bespoke Modeling

Go to www.3dsystems.com/bespokemodeling

Select your subscription level from the “Subscribe” menu and then register as a new user

<http://www.3dsystems.com/bespoke/modeling/sign-in>

Ensure that [Microsoft Silverlight](#) is loaded

Download the [“upload”](#) application.

Download the [“user manual”](#)

View the video tutorials on the [“tour”](#) page

Upload the DICOM data set and select “view” and you will have a 3D medical model in a few seconds.

If you have any questions please don’t hesitate to contact us at support.bespoke@3DSystems.com