




ICP-RIE

(Oxford Plasmalab 100)
Basic User Manual



1) ICP-RIE System Overview

Gas Sensors

Gas Leak Alarm




ICP-RIE
Loadlock

Temperature
Control Box

Logon RIE Oxford Plasmalab on FOM System

Sign Log Book



WUSTL - FOM - Schedule

Notes from instrument manager
Contact Nathan Reed 24 hours in advance for RIE Start Up and availability.

Instrument Schedule: RIE Oxford Plasmalab

- RIE Oxford Plasmalab is now Available
- Your user level on this instrument is: Instrument Manager.

02/11 02/18 02/25 03/04 03/11 03/18 03/25 04/01

Today Apr 8, 2013

Mon 04/08	Tue 04/09	Wed 04/10	Thu 04/11	Fri 04/12
09:00 - 09:30	09:00 - 09:30	09:00 - 09:30	09:00 - 09:30	09:00 - 09:30
09:30 - 10:00	09:30 - 10:00	09:30 - 10:00	09:30 - 10:00	09:30 - 10:00
10:00 - 10:30	10:00 - 10:30	10:00 - 10:30	10:00 - 10:30	10:00 - 10:30
10:30 - 11:00	10:30 - 11:00	10:30 - 11:00	10:30 - 11:00	10:30 - 11:00
11:00 - 11:30	11:00 - 11:30	11:00 - 11:30	11:00 - 11:30	11:00 - 11:30
11:30 - 12:00	11:30 - 12:00	11:30 - 12:00	11:30 - 12:00	11:30 - 12:00
12:00 - 12:30	12:00 - 12:30	12:00 - 12:30	12:00 - 12:30	12:00 - 12:30
12:30 - 13:00	12:30 - 13:00	12:30 - 13:00	12:30 - 13:00	12:30 - 13:00
13:00 - 13:30	13:00 - 13:30	13:00 - 13:30	13:00 - 13:30	13:00 - 13:30
13:30 - 14:00	13:30 - 14:00	13:30 - 14:00	13:30 - 14:00	13:30 - 14:00
14:00 - 14:30	14:00 - 14:30	14:00 - 14:30	14:00 - 14:30	14:00 - 14:30
14:30 - 15:00	14:30 - 15:00	14:30 - 15:00	14:30 - 15:00	14:30 - 15:00
15:00 - 15:30	15:00 - 15:30	15:00 - 15:30	15:00 - 15:30	15:00 - 15:30
15:30 - 16:00	15:30 - 16:00	15:30 - 16:00	15:30 - 16:00	15:30 - 16:00
16:00 - 16:30	16:00 - 16:30	16:00 - 16:30	16:00 - 16:30	16:00 - 16:30
16:30 - 17:00	16:30 - 17:00	16:30 - 17:00	16:30 - 17:00	16:30 - 17:00
17:00 - 17:30	17:00 - 17:30	17:00 - 17:30	17:00 - 17:30	17:00 - 17:30
17:30 - 18:00	17:30 - 18:00	17:30 - 18:00	17:30 - 18:00	17:30 - 18:00
18:00 - 18:30	18:00 - 18:30	18:00 - 18:30	18:00 - 18:30	18:00 - 18:30
18:30 - 19:00	18:30 - 19:00	18:30 - 19:00	18:30 - 19:00	18:30 - 19:00
19:00 - 19:30	19:00 - 19:30	19:00 - 19:30	19:00 - 19:30	19:00 - 19:30
19:30 - 20:00	19:30 - 20:00	19:30 - 20:00	19:30 - 20:00	19:30 - 20:00
20:00 - 20:30	20:00 - 20:30	20:00 - 20:30	20:00 - 20:30	20:00 - 20:30
20:30 - 21:00	20:30 - 21:00	20:30 - 21:00	20:30 - 21:00	20:30 - 21:00
21:00 - 21:30	21:00 - 21:30	21:00 - 21:30	21:00 - 21:30	21:00 - 21:30
21:30 - 22:00	21:30 - 22:00	21:30 - 22:00	21:30 - 22:00	21:30 - 22:00
22:00 - 22:30	22:00 - 22:30	22:00 - 22:30	22:00 - 22:30	22:00 - 22:30
22:30 - 23:00	22:30 - 23:00	22:30 - 23:00	22:30 - 23:00	22:30 - 23:00
23:00 - 23:30	23:00 - 23:30	23:00 - 23:30	23:00 - 23:30	23:00 - 23:30
23:30 - 00:00	23:30 - 00:00	23:30 - 00:00	23:30 - 00:00	23:30 - 00:00


Instrument Logon

Previous user comment:
Service Work in PO 031513-01
Consumables:

Confirm instrument logon:

2) Vent Loadlock

- Under "Cycling Loadlock Pumping" click **Stop** then **Vent**
- Wait until "Vent Time Left" is 0 secs (150 sec total vent time)
- Open the Loadlock chamber door

The screenshot shows the 'PUMP CONTROL' interface with the following elements:



- Process Utilities** menu at the top left.
- User** and **PUMP CONTROL** tabs at the top right.
- Status Indicators:** ARM HOME (green), ARM EXTENDED (green), ARM FAULT (green), WAFFER LIFT (DOWN), WAFFER CLAMP (DOWN).
- Load complete, wafer is in Loadlock** message.
- N2 Flow Rate:** 63.1 sccm.
- at speed** indicator.
- Set Base Pressure** button.
- Base pressure reached** message.
- Cycling Loadlock Pumping** section (circled in red):

Lid	CLOSED
Pirani	3.28e-12 Torr
Vent Time Left	0 secs
- Base pressure reached** section:

Lid	CLOSED
Process interlock	OK
Penning	2.01e-06 Torr
Cm Gauge	0 mTorr
Vent Time Left	0 secs
- Buttons:** EVACUATE, STOP, and VENT are visible at the bottom of both sections.


3) Insert Wafer

- Add wafer to Loadlock with the flat edge between the pins and close lid
- Click **Stop** → Evacuate in the Loadlock panel
- Enter a Wafer Name in the “Load Wafer or pump loadlock” pop-up window → click **OK**



4) Recipe



- Input your recipe temperature in the temperature control box
- Push the T button → use the arrow keys to select a “Temperature Setpoint” → input a new temperature using the number keypad if needed → push OK
- In the PC2000 software go to “Process → Recipes”



Recipe cont.


- Click **Load** and choose a recipe from the list
- When the temperature control box reaches the desired temperature start the etch process by clicking **Run**

(recipes should include a pump step before (PTP) and after the etching step (PTB) to automatically pump down the Loadlock, move the wafer into the process chamber, pump down the process chamber after etching, and move the wafer back to the Loadlock)




Recipe cont.

- You may edit a recipe that has been loaded by clicking on a step in the list on the left and choosing "Edit Step"
- Make a new recipe by clicking **New** → in pop-up window "Clear current recipe?" click **OK** → then click in the recipe step box on the left and choose "Insert Step" or click and drag a preprogrammed step from the "Step Library" list on the right into the recipe box and choose "Edit Step"



Recipe cont.

- In the "Process Step Editor" input Step Name, Step Time, Process Gas flows, Chamber Pressure, Helium Backing Pressure, RF Generator Forward Power, and ICP Generator Forward Power
- Click **OK** when done or if button is grayed out press **Enter** on the keyboard (press Esc to cancel)
- In the "Recipe" screen give the recipe a new name and click **Save**




5) Etch Process

- During the etch process the chamber screen will display the recipe progress and all input gas, power, and pressure settings
- When the etch process is complete a "Yellow Alert" window will pop up to notify you of the process completion → click **Continue**

(It is normal for this alert to remain active in the right hand corner of the screen)

Etch Process cont.

- In the pop-up window "Wafer has finished processing. Ready to be removed." click **OK**
- Click **Stop** → **Vent** in the Loadlock panel and wait for the vent cycle to finish
- Remove your wafer and insert another → repeat steps starting on page 11



6) Shut down

- When done using the instrument remove your sample, close the Loadlock door, and click **Stop** → **Evacuate** in the Loadlock panel

