

# MULTI-CAPILLARY DS-31 (DYE SET D w/VIC) MATRIX STANDARD KIT



Product P/N 4345829  
Insert P/N 4362883 REV A  
Printed in USA

For Research Use Only.  
Not for use in diagnostic procedures.

The Matrix Standard Set DS-31 is used to generate the "multicomponent matrix" required when analyzing 6FAM -, VIC-, NED -, and ROX - labeled DNA fragments on the Applied Biosystems 3130 and 3100 Series Systems, and 3700 DNA Analyzer. The Data Collection Software for these instruments uses the multicomponent matrix to automatically analyze the four different colored fluorescent dye-labeled samples in a single capillary.

Matrix standards do not need to be run with every set of sample injections. The standard only needs to be run once in order to generate a matrix file which is then applied to samples run under similar conditions. For more information on the use of matrix standards, refer to the instrument User's Manual or Getting Started Guide.

The kit consists of one tube of matrix standard, which is sufficient for a minimum of eight array runs. The Matrix Standard Set DS-31 contains specific sizes of DNA fragments labeled with a unique fluorescent dye label. This standard is diluted in 1X TE buffer and is stable for one year when stored at 2°C to 8°C. (Do not freeze.)

## Preparing the Matrix Standard Set DS-31 for the Applied Biosystems 3130 and 3100 Series Systems:

**WARNING! CHEMICAL HAZARD. Hi-Di™ Formamide.** Exposure causes eye, skin, and respiratory tract irritation. It is a possible developmental and birth defect hazard. Read the MSDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

1. Thoroughly mix the contents of the tube and spin briefly in a microcentrifuge.
2. Prepare the matrix standard

For a 36 cm array, combine 10 µL of standard with 190 µL of Hi-Di™ Formamide (P/N 4311320) in a 1.5 mL microcentrifuge tube.

For a 50 cm array, combine 5 µL of standard with 195 µL of Hi-Di™ Formamide (P/N 4311320) in a 1.5 mL microcentrifuge tube.

3. Mix thoroughly and spin briefly in a microcentrifuge.
4. **If a GeneAmp® 9600/9700 thermal cycler is available for denaturation, follow steps A and B below.**
  - A) Dispense 10 µL of the matrix standard/Hi-Di™ formamide mixture into a 96-well microtiter plate. For 16 capillaries, dispense into two columns (e.g. wells A1-H1, A2-H2). For 4 capillaries, dispense into 4 wells (e.g. wells A1-D1).
  - B) Cover the plate and denature at 95°C for 5 minutes. Immediately place on ice.

**If a GeneAmp® 9600/9700 thermal cycler is not available, follow steps C and D below.**

- C) Heat the mixture at 95°C for 5 minutes to denature, and immediately place on ice.
  - D) Dispense 10 µL of the matrix standard/Hi-Di™ formamide mixture into a 96-well microtiter plate. For 16 capillaries, dispense into two columns (e.g. wells A1-H1, A2-H2). For 4 capillaries, dispense into 4 wells (e.g. wells A1-D1).
5. Place the 96-well microtiter plate on the plate deck of the instrument.
  6. For specifics on setting up a run, refer to your User's Manual or Getting Started Guide.

#### **Preparing the Matrix Standard Set DS-31 for the ABI PRISM 3700 DNA Analyzer:**

1. Thaw and thoroughly mix the contents of the tube and spin briefly in a microcentrifuge.
2. Prepare the matrix standard by combining 10  $\mu\text{L}$  of standard with 390  $\mu\text{L}$  of Hi-Di™ Formamide (P/N 4311320) in a 1.5 mL microcentrifuge tube. Mix thoroughly.
3. Heat the mixture at 95°C for 5 minutes to denature and immediately place on ice.
4. Divide the mixture into two 200  $\mu\text{L}$  MicroAmp® Reaction Tubes and place tubes at position 9 and 10 in the right 8-bar of the ABI PRISM 3700 DNA Analyzer plate deck.
5. For specifics on setting up a run, refer to the ABI PRISM 3700 DNA Analyzer User's Reference Manual.

**NOTE:** If the signal height of any one of the matrix standard fragments is saturated on the above instruments, reprepare the matrix standard mix with half the volume of matrix required in step 2. For example, in step 2 for the 3130 and 3100 Series Systems with a 36 cm array, reprepare the matrix standard mix with 5  $\mu\text{L}$  of matrix standard and 195  $\mu\text{L}$  of Hi-Di™ Formamide (P/N 4311320). Rerun the spectral calibration.

#### **Notice to Purchaser: Disclaimer of License**

This product is optimized for use in the DNA sequencing or fragment analysis methods covered by patents owned or licensable by Applied Biosystems. No license under these patents to use the DNA sequencing or fragment analysis methods is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the DNA sequencing or fragment analysis methods for certain research and development activities accompanies the purchase of certain Applied Biosystems reagents when used in conjunction with an authorized DNA sequencing machine, or is available from Applied Biosystems. Further information on purchasing licenses to practice the DNA sequencing or fragment analysis methods may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.



VENDOR INFORMATION,  
PRODUCT INSERT  
FULL SHEET

Insert P/N 4362883 REV A  
Printed in USA



GUIDELINE, INSERT FOLD INSTRUCTIONS

- 1 Specification is intended for Full Sheet, two sheets max.
- 2 Material: 50# Offset (White)
- 3 Dimensions: Full Sheet (8 1/2" x 11") folded size = 1.833" x 2.125" +/- .125"
- 4 Print on one sheet, front and back.
- 5 After final folding, upper left-hand corner of page 1 is to face out, with title and part number displayed.
- 6 Folded inserts to be packaged as follows:
  - Bundled in groups of 10ea.
  - Shrink-wrap 10ea groups for 100 inserts total.
  - Label shrink-wrapped bundles [Part# & Rev (font large enough to be legible), Qty, Manufacturer, and Barcode Part #128C].
  - Seal shrink-wrapped bundles in an appropriately sized corrugated shipper (not to exceed 25#, or have outside dimensions greater than 12" x 13 1/2" x 24").
  - Appropriately label shipper [Part# & Rev (font large enough to be legible at 8 feet), Description, P.O.#, Qty, Mfg Date, Manufacturer, Country of Origin, and Barcode Part # in 128C] on at least one end panel (12" x 13 1/2" side).

**NOTICE TO PRINTER**

Print page(s) prior  
to this page **ONLY**.

This page is  
folding instructions  
which is NOT and  
MUST NOT BE PRINTED  
as part of this insert.

8 1/2"			
First Page Header Title of the Insert (Up to 3 Lines of Information) Product P/N: 00000000 Insert P/N: 00000000 Printed in (Country)		Company Logo For Research Use Only. Not for use in diagnostic procedures.	
First Page Footer ©Copyright 2002. Applied Biosystems. All rights reserved. ... are registered trademarks and ... are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.			Product Insert DD Mon.YYY
11"			

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION  
OF APPLIED BIOSYSTEMS AND IS TO BE RETAINED  
IN CONFIDENCE, NOT TO BE REPRODUCED, COPIED,  
OR INCORPORATED IN ANY MANNER.

© APPLIED BIOSYSTEMS 2004