



CODONICS[®]

We bring the future into focus

Container Labeling System[™]

Manufactured for Codonics, Inc. by BIQ Consultores de Engenharia e Gestão Industrial, Lda.

User's Manual

August 29, 2013

Version 2.4.0

Codonics, Inc.
17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1 440.243.1198 Phone
+1 440.243.1334 Fax
Email info@codonics.com
www.codonics.com

Copyright © 2012–2013 by Codonics, Inc. All rights reserved, worldwide.

Part number 905-070-401

No part of this document may be copied or reproduced in any form by any means without prior written consent of Codonics, Inc., 17991 Englewood Dr., Middleburg Heights, Ohio 44130 USA.

Although every effort has been made to ensure the accuracy of this document, Codonics, Inc. assumes no responsibility for any errors that may appear. Codonics, Inc. makes no commitment to update nor to keep current the information contained in this document.

Codonics, the Codonics logo, and “We bring the future into focus” are registered trademarks, and Container Labeling System is a trademark of Codonics, Inc.

Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries. Intel Core is a trademark of Intel Corporation in the U.S. and other countries. Athlon is a trademark of AMD Corporation in the U.S. and other countries. Motorola is a trademark of Motorola Corporation in the U.S. and other countries. Zebra is a trademark of Zebra Corporation in the U.S. and other countries.

European Authorized Representative:



CEpartner4U
Esdoornlaan 13,
3951DB Maarn
The Netherlands
Tel.: +31 343 442 524

Regulatory & Agency:

- **IT class equipment meeting EN55022(B), EN61000-3-2, EN61000-3-0, EN55024 and EN 60950 standards**
- **CLS is intended for use in medication production and inventory management control environments.**
- **CLS is not intended for use in patient care areas and is not a medical device**
- **Components safety marked by original manufacturer (CE, CCC, etc.)**

**USER'S MANUAL
TABLE OF CONTENTS**

1. PURPOSE AND SCOPE	3
2. MANUFACTURER INFORMATION.....	3
3. PRODUCT INFORMATION.....	4
4. INDICATIONS FOR USE.....	4
5. CONTAINER LABELING SYSTEM COMPONENTS	4
6. SETTING UP YOUR SYSTEM.....	5
7. GLOSSARY OF TERMS	6
8. OVERVIEW OF CONTAINER ID FORMATTING FOR DRUG CONTAINERS	8
9. OVERVIEW OF DATABASE CREATION FOR MASTER CLS DATA.....	9
10. HOW TO UPGRADE CLS.....	11
10.1 Backup the Database	11
10.2 Uninstall CLS	12
10.3 Remove CLS folder.....	12
10.4 Install new CLS version	12
10.5 Update the license	12
10.5 Import Database	12
11. HOW TO INSTALL CLS.....	17
Installation.....	17
12. HOW TO RUN THE APPLICATION FOR THE FIRST TIME	20
License	20
13. HOW TO CONFIGURE SYSTEM SETTINGS.....	24
General Settings	24
File Path	24
Database Path.....	24
CLS Settings.....	25
Parameters of the System	26
Container ID Formatting	27
WARNING: Changing the Container ID Format and Database.....	27
WARNING: Overview of Editing the Container ID Format.....	27
System Settings	30
14. HOW THE CLS APPLICATION IS STRUCTURED.....	31
15. HOW TO USE THE APPLICATION FOR THE FIRST TIME	32
Adding Users.....	33
Edit Users.....	35
16. HOW TO ADD A NEW USER TO CLS.....	37
17. HOW TO EDIT USERS	39
18. HOW TO ADD CONCENTRATION VALUES.....	41
19. HOW TO EDIT CONCENTRATION VALUES.....	42
20. HOW TO DELETE CONCENTRATION VALUES	43
21. HOW TO ADD CONCENTRATION UNITS.....	44
22. HOW TO EDIT CONCENTRATION UNITS.....	45

23.	HOW TO DELETE CONCENTRATION UNITS	46
24.	HOW TO ADD ROUTES.....	47
25.	HOW TO EDIT ROUTES	48
26.	HOW TO DELETE ROUTES.....	49
27.	HOW TO ADD RULES	50
28.	HOW TO EDIT RULES.....	54
29.	HOW TO INACTIVATE (OBSOLETE) RULES	56
30.	HOW TO TEST BARCODE PARSING	58
31.	HOW TO ADD A PRINTER.....	59
	Choosing the label layout	60
	Flag	60
	Rectangle.....	60
	Round.....	60
32.	HOW TO EDIT A PRINTER.....	61
33.	HOW TO DELETE A PRINTER	62
34.	HOW TO CALIBRATE A PRINTER.....	63
35.	HOW TO EDIT YOUR USER ACCOUNT.....	64
36.	HOW TO CHANGE USER PASSWORD.....	65
37.	HOW TO ADD MANUFACTURERS.....	66
38.	HOW TO EDIT MANUFACTURERS.....	67
39.	HOW TO ADD A NEW DRUG TO CLS	71
40.	HOW TO EDIT A DRUG.....	72
41.	HOW TO INACTIVATE (OBSOLETE) A DRUG.....	74
42.	HOW TO ASSOCIATE DRUGS TO MANUFACTURERS	76
43.	HOW TO EDIT MANUFACTURER DRUGS ASSOCIATIONS	79
44.	HOW TO DELETE MANUFACTURER DRUGS ASSOCIATIONS.....	81
45.	HOW TO PRINT A CONTAINER LABEL	83
	Using Manufacturer Barcodes.....	84
	Exceptions.....	85
	When there isn't a Manufacturer Barcode.....	88
	Exceptions.....	91
46.	HOW TO VERIFY A CONTAINER CODE	92
47.	HOW TO VERIFY A MANUFACTURER CODE.....	93
48.	HOW TO ANALYZE THE LOGS	95
49.	HOW TO EXPORT DATA TO A CSV FILE	97
50.	HOW TO USE CLS ON A NETWORK.....	98
51.	HOW TO SECURE YOUR CLS	99
	Backup	99
	Protect against Malware and Intruders.....	99

1. PURPOSE AND SCOPE

Refer to this User's Manual for procedures on how to perform the Container Labeling System's user operations, including:

- Setting up the hardware and software
- Overview of terms and system configuration
- Entering configuration information
- Entering master data
- Printing labels
- Exporting data
- Maintaining the system
- Monitoring events and log files

2. MANUFACTURER INFORMATION

The following image is the Help | About screen that displays the manufacturer information of BIQ Consulting Inc. for Codonics as the Container Labeling System (CLS).



3. PRODUCT INFORMATION

For technical assistance with the Container Labeling System, call Codonics Technical Support at the following number:

Phone: +1.440.243.1198
Toll Free: 800.444.1198 (USA only)

Codonics Customer Service is available weekdays from 8:30 AM to 5:30 PM EST (U.S. holidays excluded). Technical Support is also available online via email and the Codonics web site:

Email: support@codonics.com
Web Site: www.codonics.com

General product information can also be requested by sending email to:

Email: info@codonics.com

Please include your postal mailing address and telephone number in the email message. Basic product information is returned via email. Complete product literature is sent by postal mail.

4. INDICATIONS FOR USE

The Codonics Container Labeling System (CLS) integrates software and hardware to provide site management over the labeling of medication vials for identification purposes under pharmacy control. CLS is a pharmacy level system/tool that aids in the "production/inventory/site identification" of parental vial medications. Site specified UDI are in addition to manufacturer vial labeling and do not obscure original vial labeling in any way. CLS labels added to vials assist in site tracking, inventory management, and identification of medications purchased and approved for distribution within the institution. CLS barcodes are configured to meet industry symbology standards for easy integration with devices that read machine readable codes. The product is not a medical device and the PC (computer), barcode scanner and printer components are marked by the original manufacturer for their intended use, installation environments and WEEE disposal instructions. CLS is not to be placed in a patient care area.

5. CONTAINER LABELING SYSTEM COMPONENTS

The Codonics Container Labeling System (CLS) is comprised of software, hardware and label media.

The software includes the following features, which are detailed throughout this manual.

- Integration with desktop printer
- Integration with barcode scanner.
- Container ID code formatting
- User group roles
- Manufacturer configuration
- Drug configuration
- Manufacturer drugs configuration
- Route configuration

-
- Concentration configuration
 - User configuration
 - Printer configuration
 - Verification workflow
 - Workflow approval
 - Lot and expiration date tracking
 - Container code verification
 - Manufacturer code verification
 - Label printing
 - Export of CLS data in CSV format

The hardware provided with your system includes the following items:

- An easy-to-use Zebra® desktop direct thermal printer
- An easy-to-use Motorola® barcode scanner

The label media (see NOTE below) provided with your system includes three unique label designs for use in pharmacy, which are detailed throughout this manual.

- Flag labels intended to be wrapped around the neck of drug containers.
- Round labels intended to be applied to the bottom of drug containers.
- Rectangular labels intended to be applied to the side of drug containers.

NOTE: Use only Codonics labels to ensure proper operation of the device and proper labeling of drug containers. Using unapproved labels could lead to unacceptable results, including poor print quality and poor label adhesion to the drug containers. Please contact Codonics Customer Support to order label media for your CLS system.

NOTE: Storage of Codonics labels is recommended to be in a dark place at a relative humidity between 45 percent and 65 percent and a temperature below 77°F (25°C) will assure satisfactory performance.

6. SETTING UP YOUR SYSTEM

It is recommended that the Container Labeling System be set up in the following order:















- 1.) Install the printer according to the Zebra installation guide and Zebra disc included.
 - a. Unpack.
 - b. Install the power supply.
 - c. Install connectivity (USB or Ethernet)
 - d. Run the Zebra set up software – ZebraSetupUtilities & ZebraDesigner
 - e. Install label media.
- 2.) Install your scanner according to Motorola quick start guide.
 - a. Unpack.
 - b. Install connectivity

-
- 3.) Configure your scanner for use with the Container Labeling System software. Please refer to the Codonics supplied Release Notes document – part number 901-448-001 – with your scanner named “Container Labeling System Barcode Scanner Configuration” to program your scanner.
 - 4.) Install the Container Labeling System software per the instructions provided in this manual.
 - 5.) Configure your printer in the Container Labeling System software per the instructions provided in this manual.
 - 6.) Perform a printer calibration in Container Labeling System software per the instructions provided in this manual.

7. GLOSSARY OF TERMS

User Groups

There are four user group types (administrator, manager, technician, operator). Each user account created is granted permissions in the system by being assigned to one of the four group types.

	System Configuration	Barcode Search	Manual Search	Workflow Approval	Print Container Labels
Administrator					
Manager					
Technician					
Operator					

Workflow Approval

An administrative or managerial workflow process of the CLS allows for the quality assurance review of confirming a technician's or operator's label selection prior to the labels being printed.

This optionally enabled verification as assigned at the Drug record will require authorized Workflow Approval users to review the containers to be labeled and the label intended to be printed. Authorized users with Workflow Approval will enter their user name and password to approve the label printing or cancel the pending system request if a selection error occurred by the technician or operator.

Verification Workflow

An administrative workflow process of the CLS which allows for quality assurance review for confirming the configuration of specific container codes to manufacturer codes.

This optionally enabled verification will provide a CLS to have a “double-check” confirmation process prior to releasing each specific association of a container code to a manufacturer code into an “in use (production)” state.

Manufacturers

The name of the companies or laboratories that will be managed in the CLS as providers of drugs used in the hospital's pharmacy.

Drugs

Drugs that will be tracked in the CLS. Associated data includes a container code, drug name, Tallman spelling, concentration value, concentration units and route.

Manufacturer Drugs

A CLS association of a specific drug as supplied by a specific manufacturer.

Manufacturer IDs

The value assigned to a manufacturer as defined by your hospital.

Manufacturer Codes (Barcodes)

The barcode value printed on a manufacturer's box of drug containers.

Container Code

The ID value entered or system assigned to uniquely identify a drug type by container.

Concentrations

The values and units that are defined to a container's drug.

Routes

The intended delivery method type of the container's drug to a patient.

Labels

Labels are printed labels by the CLS printer for application to drug containers. Labels are available in three die-cut shapes:

Flag (Neck) – intended to be wrapped around and adhered to the neck of a container.

Rectangular (Seam) – intended to be applied and adhered to the side of a container.

Round (Bottom) – intended to be applied and adhered to the bottom of a container.

Events

Details of actions and activities as performed in the CLS for an administrator to review.

8. OVERVIEW OF CONTAINER ID FORMATTING FOR DRUG CONTAINERS

The Container ID Formatting configuration should be carefully considered prior to the start of the utilization of the Container Labeling System. Evaluation of the identification numbering system for the drugs in your pharmacy is critical to ensure that each drug's container can be uniquely labeled. Considerations include identification methodology that are already in use in your hospital, national standards in place by your country's government or emerging national standards that are being implemented for use within your country, or with drugs imported from neighboring and/or partner countries.

If your hospital and/or country does not have any existing standards for an identification numbering system, a schema of parameters can be defined for use in the Container Labeling System to assign unique identifiers for each drug. The defined schema can then be used via the Container Code Auto-Assign functionality when each drug is entered into your CLS.

The Container ID (CID) parameters will be required to be set for the format of choice for your Container Labeling System upon system start up.

Application Identifier – this is site configurable, with a default value of “COD”. This value is a prefix intended to help uniquely identify the drug containers used by your hospital, i.e. your hospital network. This value can be configured from 0-4 characters with alpha-numeric values (A-Z, a-z and 0-9). Optionally, if no value is configured, there will be no Application Identifier included in the barcode for the CID.

It is highly recommended to utilize a value for the Application Identifier as this will allow a distinguishing value to be used in your printed barcode labels which can mitigate the possibility of conflicting with any existing or emerging standard. An example with the default configuration of the Application Identifier value of “COD” and the Variable Identifier (described below) of 13 digits would create a barcode value of “COD00000000000123”.

Fixed Identifier – this is site configurable, with no default value. This value is a secondary prefix intended to further help in uniquely identifying the drug containers used by your hospital network, i.e. a regionally based hospitals where you may want to distinguish drug containers in one region differently than those at another region. The value can be configured from 0-4 characters with numeric values only (0-9). If included in your configuration, this prefix value is included in the drug's overall CID barcode value. Optionally, if no value is configured, there will be no Fixed ID included in the barcode for the CID.

It is recommended to utilize a value for the Fixed Identifier as this will allow for further distinction of your barcode values for mitigation of conflicts with any existing or emerging standard. An example with the Application Identifier of “COD”, a Fixed ID of “7777” and a Variable Identifier of 9 digits (described below) would create a barcode value of “COD7777000000123”.

Variable Identifier – this is site configurable, with a default value of 13 digits. This value is required for uniquely identifying the drug containers used by your hospital. The value can be configured numerically from a minimum of a 6-digit length up to 15 digits in length. This value will be included in the barcode for the CID.

Please refer to the “Parameters of the System” section of this manual for the details on the configuration of the Container ID Formatting parameters for your Container Labeling System.

Please refer to the “Parameters of the System” section of this manual for the details on the configuration editing of the Container ID Formatting parameters for your Container Labeling System.

Please also see the WARNING section in Chapter 12 entitled “WARNING: Changing the Container ID Format and Database” to ensure that you fully understand and accept the implications for changing your Container ID parameters after your CLS has been configured and placed into use.

9. OVERVIEW OF DATABASE CREATION FOR MASTER CLS DATA

The database creation process for the Container Labeling System requires the entry of the following types of data based upon the drugs inventoried and utilized within your hospital’s pharmacy:

Manufacturers

The name of the company or laboratory that manufactured the drugs inventoried and utilized in your hospital’s pharmacy.

Drugs

The drugs inventoried and utilized in your hospital’s pharmacy.

Manufacturer Drugs

An association that connects a specific manufacturer to a specific drug for utilization in your hospital’s pharmacy.

It is recommended that the database entry process occur in the order of the terminology outlined above, which is Manufacturers followed by Drugs and completed with Manufacturer Drugs.

Database entry is only available to users that are granted with rights of the Administrative group. Please see the User Group permissions table.

Entering information into your CLS database for Manufacturers is the process for managing specific companies or laboratories that your hospital procures drugs from for your pharmacy. This can be accomplished by entering the data based upon lists of Manufacturers that you have available from your hospital’s information system or inventory management systems as well as the physical boxes of drugs inventoried in your pharmacy. Please refer to the “How to Add Manufacturers” section of this manual for the details of the data entry process.

Entering information in your CLS database for Drugs is the process for managing specific drugs that your hospital procures for your pharmacy. This can be accomplished by entering the data based upon lists of Drugs that you have available for your hospital’s information system or inventory management systems as well as the physical boxes of drugs inventoried in your pharmacy. Please refer to the “How to Add a Drug” section of this manual for the details of the data entry process.

Entering information in your CLS database for Manufacturer Drugs is the process for associating a specific Drug to a specific Manufacturer. Having the physical drug container boxes available on hand is required as you will scan the barcode provided by the drug’s manufacturer or laboratory. The physical boxes that have an available barcode to scan may either be on an individual box of drug containers or on a larger carton containing multiple smaller boxes of drug containers. Please refer to the “How to Add a Manufacturer Drug” section of this manual for the details of the data entry process.

To accomplish the entry of data into your system, an administrative user will perform the following steps to create CLS master data for utilization. An example data entry use case for CLS master data would be entered in the following order:

1. Manufacturers
2. Drugs
3. Manufacturer Drugs

The entry of the data would be accomplished in the following manner with this example data of 3 manufacturers, 10 drugs and 10 manufacturer drugs associations:

As an administrator of the system, enter your 3 manufacturers into the CLS database.

Manufacturer 1
Manufacturer 2
Manufacturer 3

Now that you have your 3 manufacturers entered, you will then enter your 10 drugs into the CLS database.

Drug 1
Drug 2
Drug 3
Drug 4
Drug 5
Drug 6
Drug 7
Drug 8
Drug 9
Drug 10

At this point, your CLS database includes 3 manufacturers and 10 drugs. The last step in setting up your data is to associate specific manufacturers to specific drugs.

Manufacturer 1
 Drug 1
 Drug 5
 Drug 8

Manufacturer 2
 Drug 2
 Drug 3
 Drug 9
 Drug 10

Manufacturer 3
 Drug 4
 Drug 6
 Drug 7

Your CLS database now includes 3 manufacturers, 10 drugs and 10 manufacturer drugs for use in printing barcode labels for your drug containers.

10. HOW TO UPGRADE CLS


This option is used to upgrade the CLS software

Note: This option is only available before the application is upgraded and while there is no data in the new CLS database.


10.1 Backup the Database

Please backup your CLS database. This database backup file will be used during the upgrade process to import the existing data from your current CLS software version to the upgrade version.

To find your Database file, click *Settings* and then *System*. The *System* configuration window is displayed.




Under the *General* tab you can find the *Database Path* field which defines the path to the database that the CLS is using. The default path is C:\CLS\BD\DB.mdb. Backup your database file in a secure location – out of CLS application folder.



Check the root CLS folder defined in the *File Path* field. The default root folder for the application is C:\CLS. You will need to remove the CLS folder later in this section of the manual.

10.2 Uninstall CLS

Go to Windows *Control Panel*. Run the *Add or Remove Programs* option.



The *Currently installed programs* list is displayed. Select the CLS application from the list and uninstall it.

10.3 Remove CLS folder

In Windows Explorer, remove the CLS root folder with all subfolders and files. The default root folder is C:\CLS.

10.4 Install new CLS version

For further information on how to install a new CLS version please see the *HOW TO INSTALL CLS* section of this manual.


10.5 Update the license

When starting a new version of CLS for the first time, a new license will be required to activate the CLS.


For further information on how to activate the product please see the *HOW TO RUN THE APPLICATION FOR THE FIRST TIME* section of this manual.

10.5 Import Database

When running the new version of CLS application for the first time after the license process, the *System* configuration window is displayed. Otherwise, click *Settings* and then *System*.



Click on the tab *Import*.




Click on the folder icon to select the database backup file. Select the backup file you created previously in this section.

Note: You need a valid backup file of the CLS database in order to proceed.


Press the *Verify* button to verify the selected file.

If the selected backup file is correct, a *Correct Database version* message is prompted.

Press the *OK* button after database version validation.




Press the *OK* button after database structure validation.




If validation process completes successfully, the *Update* button is available.


Press the *Update* button to update the application.




Press the *OK* button after update process is finished.



Press the *Save Configuration* button after update process is finished.



Note: If the wrong backup file is selected, “*Incorrect Database version*” message is prompted.



11. HOW TO INSTALL CLS

This chapter explains how to install the CLS on your computers.

The following minimum system standards are required for deployment on intended computers:

- Windows XP Professional SP3 or Windows 7 Professional SP1 (32-bit or 64-bit)
- .Net Framework 3.5
- 500MB available disk space
- 2GB RAM
- Intel Core 2 Processor or AMD Athlon64 X2


Installation

Copy the files from the CLS software installation set to your local operating system. Please note where you copy the files, as they will be required later in the configuration process.

To install CLS, you need to double click on the setup.exe icon.




Welcome to the CLS Setup Wizard, click on the *Next* button.




Choose the destination folder (the user should have writing permissions) and click *Next*.



On the next window, confirm the installation and click the *Next* button.



When the installation is complete, click on the *Close* button.



To access the application, the icon will be available on your computer's desktop and on the Start menu.

12. HOW TO RUN THE APPLICATION FOR THE FIRST TIME



To run the application, double click on the  icon on the desktop or select from the Start menu.

License

When starting the application for the first time, a license will be required to activate the CLS.


Note: the license that is provided to activate your CLS software is workstation-specific. Please ensure that the application is installed on the specific workstation or workstations that will utilize the Container Labeling System software. A license cannot be transferred to another workstation.

Enter the requested information in the fields displayed with the given information.

A screenshot of the CLS application's license request dialog box. The window has a title bar with the CLS logo and standard Windows window controls. Inside, there are several input fields: 'Distributor Code:', 'Site Name:', 'Language File:', and 'Language in Use:' (which is a dropdown menu). Below these is a large blue circular button labeled 'Get Key'. Underneath the button are two large text areas labeled 'Key:' and 'License:'. At the bottom, there is a block of text providing instructions on how to request a license via email, and two buttons labeled 'OK' and 'Cancel'.

In the Language File option, find the '*Language.xml*' file. It will be located with the CLS setup.exe where you placed the files on your local operating system from the CLS software installation set at the start of your installation process.

Note: The '*Language.xml*' file is not located within the CLS directory created during the installation process.




Upon selecting the *Language.xml* file, the default language in use will be English. You can select another supported language as provided in the drop-down list.

To generate a key for the application, click on *Get Key* button.

A screenshot of the CLS application window. It contains several input fields: 'Distributor Code:' with 'Codonics' entered, 'Site Name:' with 'Codonics' entered, 'Language File:' with 'C:\CLS_v2.3.0\Language_2012-02-08_173400.xml' entered, and 'Language in Use:' with a dropdown menu showing 'English'. A prominent 'Get Key' button is located below these fields. Below the button are two large text areas labeled 'Key:' and 'License:'. At the bottom, there is explanatory text about email support and license activation, followed by 'OK' and 'Cancel' buttons.


The application key will be generated and displayed in the *Key* field.



The screenshot shows the CLS application window. The 'Distributor Code' and 'Site Name' fields are both set to 'Codonics'. The 'Language File' is 'C:\CLS_v2.3.0\Language_2012-02-08_173400.xml' and 'Language in Use' is 'English'. The 'Key' field is populated with a long alphanumeric string: RvVwHJM9pVTSUV8QPFpKfNDmNBdWAomLhwJbM+<P6sIXVIgzwHpeHnpXXPUdcofOuaaZO2KScAMXJH2Q2Rn1.KPIAsUePIS3VivZc9dxbagB0HnL9r7VVqaCCEDvVr55rPeVQ9QJUA+5R.IAxR7bzmmawuccu.C2BvWCveH0tkireLHnH pVrgFucpq8GmLrkH9llie5UeTcnbDB8O:nHgOsOwllmCHIAKJpm:h5rIDtHpG6qMbCw/p+DtGNLn+8U:URAAQTQ==. The 'License' field is empty. Below the fields, there is a message: 'Please copy and paste your key into an email with the subject line of "CLS License Request" and include your site name and distributor name and send it to clssupport@codonics.com. A license will be emailed to you within the next 24 hours to be applied into the license field to activate your CLS.' At the bottom are 'OK' and 'Cancel' buttons.

To request the *License* for activating your system, please copy the key and paste it into an email with the subject line ‘*CLS License Request*’ and send it to clssupport@codonics.com. Please include your distributor name and site name in the email. A response will be sent within 24 hours.

When you receive the email with the License, copy and paste it into the *License* field.



The screenshot shows the CLS application window. The 'Distributor Code' and 'Site Name' fields are both set to 'Codonics'. The 'Language File' is 'C:\CLS_v2.3.0\Language_2012-02-08_173400.xml' and 'Language in Use' is 'English'. There is a 'Get Key' button above the 'Key' field. The 'Key' field is populated with the same alphanumeric string as in the previous screenshot. The 'License' field is now populated with a long alphanumeric string: DPFxIEf.S0du1MRcQJnCbupmmw35SzUyitz4ZHonp9AH0edYEDWnr7GnJN+ wVem85CX8LIX7HkA5GIMgBvBH97FwLVQ93.B+ ON6hhRjHyUj.o52xoCl94ltCnrfmInSj.NqlBICND5rTupnYGUrwUihE8aA/G8GyKtY5+y1CD3Q4I82JhG/cHeb5gZ4pH2e5E1xmW0EgzeYZZLIEHmBtLbRiVFFOI5HS GQ689u3erCOMI1dbd6Q8penh37d.CX6hmBAImhqXpFujTQJp0ik.O5tIAFcaQITQesropFEONoR29ehA16VPAGgrA3G2L6xh77dW0xH2op2u+ygLLHczbPdEzwJ/2a4EIJU8C/dIkOT0y. Below the fields, there is a message: 'Please copy and paste your key into an email with the subject line of "CLS License Request" and include your site name and distributor name and send it to clssupport@codonics.com. A license will be emailed to you within the next 24 hours to be applied into the license field to activate your CLS.' At the bottom are 'OK' and 'Cancel' buttons.

To activate your CLS, click on the *OK* button.

If any of the fields were selected incorrectly upon generating your system's key, the license activation will notify you with an error message.



If all fields were entered correctly upon generating your system's key, the license activation will be accepted and a message will be displayed welcoming you to the CLS.




Now you will be prompted to configure the system settings, explained in *HOW TO CONFIGURE SYSTEM SETTINGS*.


Upon completion of the system settings, the CLS is ready to be utilized.

13. HOW TO CONFIGURE SYSTEM SETTINGS

When running the CLS application for the first time after the license process, the *System* configuration window is displayed. Otherwise, click *Settings* and then *System*.



General Settings




File Path

The *File Path* field defines the destination folder for exported CSV files. The default destination folder is C:\CLS\Export. If a different folder is desired, select the browse folder icon to choose a different folder, which will then be displayed in the *File Path* section.

Database Path

The *Database Path* field defines the path to the database that the CLS will use. The default path is C:\CLS\BD\DB.mdb. If a different path is desired, select the browse folder icon and select the path to where the DB.mdb file was placed, which will then be displayed in the *Database Path* section.

After defining the database path, test the database connection by pressing the *Test Connection* button to check its availability. If *Test Connection* fails, the CLS will give an error message.




If the connection to the database is successful, the CLS will give a success message.




CLS Settings

CLS Description identifies the computer in use with the CLS system as more than one computer can be configured for use. This description will identify each computer being utilized by the CLS system.




Note: All the fields are required, except the *CLS Description* field.

If you leave *File Path* field empty, the CLS will provide an error message.




If CLS can't find a database on the given path , the application will provide an error message.



Parameters of the System

To change the application parameters click on the tab *Parameters*



Container ID Formatting

The Container ID Formatting for the Container Labeling System may be edited or deleted. Please ensure that you fully read and understand the information below prior to editing or deleting your container ID formatting.

WARNING: Changing the Container ID Format and Database

The CLS is designed to ensure that each new drug added to the system is associated with a unique Container ID. CLS will not allow duplicate Container IDs to exist in the database. This is an important feature for the safe operation of the CLS. Many devices and applications that use drugs labeled by the CLS rely on each drug having a unique Container ID.

Because having a unique Container ID associated with each drug is so important, this section of the User's Manual describes some HIGH RISK functions of CLS that can cause duplicate Container IDs to be assigned to drugs. The functions described in this section should only be used by experienced administrators who fully understand the operation of the CLS and the implications of changing the Container ID Format or deleting the database.

The CLS provides two special functions to Edit Container ID Format and Delete the Container ID database. These functions should only be used when no other options are available. Please study the information in the section carefully before using these functions.

Warning: Use of the functions described in this section can result in duplicate Container IDs.

WARNING: Overview of Editing the Container ID Format

The CLS can be set up to produce unique Container IDs that conform to a specific format. The Container ID Format includes an Application Identifier, a Fixed ID and a total Container ID Length. Once the Container ID Length is set, the CLS will not allow the length to be changed without deleting the entire Container ID database.


The Application Identifier and Fixed ID number can be changed, but should be done so using EXTREME CAUTION.

The Fixed ID is often used in Hospitals systems with multiple locations that are labeling drugs. Assigning a unique Fixed ID number to each location using the CLS allows Container IDs to be created that will not conflict with other sites. This can be useful if multiple locations are exchanging drugs that were labeled on different CLS systems. Changing the Fixed ID format at a location may cause a conflict with another location if the same Fixed ID is already in use.

EXAMPLE: If Location A is using Fixed ID 3333 and the current Container ID is 333300000001, and Location B is using Fixed ID 4444 and the current Container ID is 444400000001, then changing Location A from Fixed ID 3333 to 4444 will result in the same Container ID being issued by two different CLS systems. Different drugs labeled at each location will have the same Container ID.

It is HIGHLY RECOMMENDED that once the CLS is configured with a Container ID Format and drugs are being labeled by CLS, that Container ID Format not be changed.

The CLS will provide a warning and require confirmation when a user attempts to Edit the Container ID Format.



WARNING: Overview of Deleting the Container ID Database

The CLS stores all Container IDs in a database. Each time a new drug is added, a corresponding Container ID is created. CLS always verifies that each new Container ID is not a duplicate of another in the database. This guarantees that each new drug added to the CLS will have a unique Container ID assigned.


Deleting the Container ID database will erase all information about previous Container IDs from the system and prevent the CLS from verifying that new Container IDs are unique. When deleting the Container ID database, the CLS will continue to automatically generate new Container IDs using the last number that was created before the database was deleted.

Deleting the Container ID database on CLS can result in different drugs being labeled with the same Container ID. This applies mostly to Container IDs that manually entered into the CLS but can also apply to Container IDs that are generated automatically by the CLS.

EXAMPLE: If Container ID 0000000000123 was manually entered into CLS and originally assigned to the drug Propofol, deleting the CLS Container ID database will allow Container ID 0000000000123 to be manually re-entered and assigned to a different drug such as Penicillin. Scanning a Propofol drug vial with Container ID 0000000000123 would report Penicillin.


It is **HIGHLY RECOMMENDED** that once the CLS is configured with a Container ID Format and drugs are being labeled by CLS, that the CLS Container ID database not be deleted.

The CLS will provide a warning and require confirmation when a user attempts to Delete the Container ID database.



Container ID Formatting

The CLS is designed to ensure that each new drug added to the system is associated with a unique *Container ID*. The length of the *Container ID* can be set between 6 and 15 characters, sharing it with the *Fixed ID*. We can allocate up to 4 characters of the *Container ID* length to be a *Fixed ID*.




Keep in mind that the *Container ID Length* is shared with the *Fixed ID* and there's a minimum of 6 characters following the *Fixed ID*. For example, having the *Container ID Length* of 10 characters and the *Fixed ID* of 4 characters (i.e. 0001) leaves the minimum of 6 characters available as a *Container ID*.

System Settings

If the *Verification Workflow* field is checked then CLS will set *Manufacturer Code Status Pending* for new associations of Manufacturers to Container Codes. If not checked, then all associated Manufacturer Codes will be set to Approved/In Use.


If the *Box Quantity Editable on Printing* field is checked, the user can change the number of labels to print with each print job.

If the *Lot and Expiration Date Mandatory* field is checked, the user must specify the lot and expiration date prior to printing labels.




The screenshot shows the 'System' dialog box with the 'Parameters' tab selected. The 'Container ID Formatting' section contains three input fields: 'Container ID Length [6..15]:' with the value '11', 'Fixed ID [max length = 4]:' with the value '2222', and 'Application Identifier:' with the value 'COD'. To the right of the first two fields is a 'Delete' button, and to the right of the third is an 'Edit' button. Below this section is the 'System Settings' section, which contains three checked checkboxes: 'Verification Workflow', 'Box Quantity Editable on Printing', and 'Lot and Expiration Date Mandatory'. At the bottom of the dialog are two buttons: 'Save Configuration' and 'Cancel'. A blue circular highlight is placed over the 'Save Configuration' button.

To save the changes made, press the *Save Configuration* button and the following dialog will be presented confirming your changes.



14. HOW THE CLS APPLICATION IS STRUCTURED




Menu – Access to different features of the application. The menu is different for each type of user.

Window Controls – Allow to minimize, maximize and close the CLS application.

Username Bar – Show the currently logged in user.

15. HOW TO USE THE APPLICATION FOR THE FIRST TIME

The image shows the login interface for the CLS application. At the top left is the Codonics logo, which consists of a blue hexagon with a stylized 'C' inside, followed by the word 'CODONICS' in a bold, blue, sans-serif font. To the right of the logo, the text 'CLS' is displayed in a smaller, black, sans-serif font. Below the logo and text, there are two input fields. The first field is labeled 'Username' and the second is labeled 'Password'. Both fields are empty and have a light gray border. Below the password field, there are two buttons: a blue 'Login' button and a yellow 'Cancel' button.

The first time you logon to CLS, the database is empty, therefore no users are available yet. To logon to the application for the first time, you need to use the default Codonics username and password as shown below:

Username: admin
Password: DDMM*

*Where DD represents the day and MM represents the month of date that you are accessing the application.

Example: You are accessing the application on the September 26th, so the password will be 2609. This password is always 4 numeric characters long.


If you mistyped the username or password the following message will be displayed:



If the username and password are entered correctly, the application will be available.

Adding Users

In the main menu choose *Manage* then *Users*.



Enter the Name, Username, Password, Status and Group of the User and press the *Add User* button.


Note: See the “HOW TO ADD A NEW USER TO CLS” section for definition about application roles available for each CLS Group type.

If the *Perform Verification Workflow* field is checked then the user will have permission to Approve or Reject manufacturer’s codes that were associated to container codes. In case of approval, it will be available in system to be used, if rejected, it won’t be available for use to print labels.


This procedure is later explained in the section named “HOW TO VERIFY A MANUFACTURER CODE”.

All of these fields are required except the checkbox for *Perform Verification Workflow*, which is optional per user.


NOTE: No user will be able to log in until an Administrator User is enabled. The “Admin” Username is a reserved username that will no longer be available for log on (it is a one-time use only log on to initialize the system).

A screenshot of the 'Manage Users' dialog box. The 'Add User' tab is selected. The form contains the following fields: 'Name:', 'Username:', 'Password:', 'Status:' (with 'Active' selected), and 'Groups:' (a list box with 'Administrator', 'Manager', 'Technician', and 'Operator' options). Below these is a checkbox for 'Perform Verification Workflow'. At the bottom are 'Add' and 'Close' buttons. The 'Add' button is circled in blue.


If any of the required fields are not configured, the corresponding messages will be displayed:



It is not possible to add users with the same username. The following message will appear if a username is duplicated:




If no error(s) occurs, the user was entered successfully and the following message appears:




Note: After you add an Administrator user, the Codonics username and password is no longer available (it is a one-time use only log on). Please be sure to create at least one Administrative user account for ongoing maintenance and configuration of the CLS system.

Edit Users

The *Edit User* tab allows editing of the existing users in the database. To edit a user, double click on the desired user or select it and press the *Edit User* button.




Edit the user data and press the *Save* button.



If any of the fields are not configured, the corresponding messages will appear:




If no error occurs, the following message will appear:



Now you're ready to logout and start using CLS with the created users.















16. HOW TO ADD A NEW USER TO CLS

In the main menu choose *Manage* then *Users*.



Fill in the Name, Username, Password, Status and Group of the User and press the *Add User* button.

There are four user group types (administrator, manager, technician, operator). Each user account created is granted permissions in the system by being assigned to one of the four group types.

	System Configuration	Barcode Search	Manual Search	Workflow Approval	Print Container Labels
Administrator					
Manager					
Technician					
Operator					

If the *Perform Verification Workflow* field is checked then the user will have permission to Approve or Reject manufacturer's codes that were associated to container codes.

All of these fields are required except *Perform Verification Workflow*.

NOTE: No user will be able to log in until an Administrator User is enabled. The Admin Username is a reserved username that will no longer be available for log on.

Manage Users

Add User Edit User

Name:

Username:

Password:

Status:

☒ Active

☐ Inactive


Groups:

	Group
<input checked="" type="checkbox"/>	Administrator
<input type="checkbox"/>	Manager
<input type="checkbox"/>	Technician
<input type="checkbox"/>	Operator


☐ Perform Verification Workflow

Add Close


If any of the fields are not configured, the following messages will appear:



It is not possible to add users with the same username. If a username is duplicated, the following message will appear:




If no error occurs, the user was entered successfully and the following message appears:




17. HOW TO EDIT USERS

In the main menu choose *Manage* then *Users*.



Click on the *Edit User* tab, double click on the desired user or select it and press the *Edit User* button.

Note: A user may be inactivated from the CLS, but not deleted. This allows the Events logs to maintain the history performed in the system by all users (active or inactive). Additionally, an inactivated user may also be reactivated as needed by the administrator.



Edit the user data and press the *Save* button.


The 'Users' dialog box contains the following fields and options:

- Name:** Dr. John Smith
- Username:** john
- Password:** (masked with asterisks)
- Status:** ☒ Active, ☐ Inactive
- Group:** A list box with 'Administrator' selected (checked), and 'Manager', 'Technician', and 'Operator' as options.
- Perform Verification Workflow:** ☒
- Buttons:** 'Save' (circled) and 'Close'.

If any of the fields are not configured, the corresponding messages will appear:



If no error occurs, the following message will appear:




18. HOW TO ADD CONCENTRATION VALUES


In the main menu choose *Manage*, then *Concentration Units*.

Note: The CLS is pre-populated with a default set of concentration values.

Note: The default values are separated with periods. Values can be modified as needed including using commas rather than periods.




To add a new concentration value, fill in the *Concentration Value* field and press the *Add* button.




When you've finished adding values, press the *Close* button.

19. HOW TO EDIT CONCENTRATION VALUES

In the main menu choose *Manage*, then *Concentration Units*.




To edit an existing concentration value, double click the desired value and modify it as needed.




When you've finished editing values, press the *Close* button.

20. HOW TO DELETE CONCENTRATION VALUES

In the main menu choose *Manage*, then *Concentration Units*.



To delete a concentration value, click the desired value and press the *Delete* button.




When you've finished deleting values, press the *Close* button.


21. HOW TO ADD CONCENTRATION UNITS

In the main menu choose *Manage*, then *Concentration Units*.

Note: The CLS is pre-populated with a default set of concentration units.




To add a new concentration unit, select the *Concentration Units* tab, fill in the *Concentration Unit* field and press the *Add* button.




When you've finished adding units, press the *Close* button.

22. HOW TO EDIT CONCENTRATION UNITS

In the main menu choose *Manage*, then *Concentration Units*.




To edit an existing concentration unit, double click the desired unit and modify it as needed.




When you've finished editing units, press the *Close* button.

23. HOW TO DELETE CONCENTRATION UNITS

In the main menu choose *Manage*, then *Concentration Units*.




To delete a concentration unit, click the desired unit and press the *Delete* button.




When you've finished deleting units, press the *Close* button.

24. HOW TO ADD ROUTES

In the main menu choose *Manage*, then *Routes*.




To add a new route, fill in the *Route* field and press the *Add* button.




When you've finished adding *Routes*, press the *Close* button.

25. HOW TO EDIT ROUTES

In the main menu choose *Manage*, then *Routes*.




To edit an existing route, double click the desired route and modify it as needed.




When you've finished editing, press the *Close* button.

26. HOW TO DELETE ROUTES

In the main menu choose *Manage*, then *Routes*.




To delete, click the desired route and press the *Delete* button.



When you've finished deleting *Routes*, press the *Close* button.

27. HOW TO ADD BARCODE RULES

In the main menu choose *Manage*, then *Barcode Rules*.



The *Barcode Rules* option is used to define Barcode parsing for the barcodes scanned on manufacturer's drug boxes and CLS printed barcodes. In this configuration, the administrative user can define barcode parsing rules according to Barcode Symbology, Application Identifier and Barcode length.

Additionally, the configuration enables the ability to set up Lot and Expiration Date if those values are available on the barcodes provided by the manufacturers.

The administrative user can create multiple barcode rules.


Note: Only the *Manufacturer Code* field is mandatory for a barcode parsing rule.

1 – In the *Manufacturer Code* scan a Manufacturer Barcode on a selected drug box.


Note: your CLS hand scanner must be configured with the programmed barcode configuration provided in CLS Technical Briefs provided with your hand scanner. Please see the CLS software version 2.4.0 or higher Technical Brief for the DS9808 hand scanner or the CLS software version 2.4.0 or higher Technical Brief for the DS4208-HC hand scanner.

Note: the barcode symbology is automatically detected by the hand scanner and identified on the *Add Barcode Rules* screen as shown below. Additionally, the barcode symbology identifier is added as a prefix to the barcode value.

2- Using the mouse, highlight the characters that define the *Application Identifier* and press the F2 key on the keyboard.



3 - With the mouse, highlight the characters that define the *Manufacturer Code* and press the F3 key on the keyboard.



4 - With the mouse, highlight the characters that define the *Expiration Date* and press the F4 key on the keyboard.

The screenshot shows the 'Barcode Rules' window with the 'Add Barcode Rule' tab selected. The 'Barcode Mask' section contains the text 'AD135790L45670061215'. The 'Expiration Date' field, which contains '061215', is highlighted with a blue selection box. To the right of the text input fields is a table with three columns: 'Field', 'Start Position', and 'Length'. The rows are: (F2) Application Identifier (AD, Start Position 1, Length 2), (F3) Manufacturer Code (135790, Start Position 3, Length 6), (F4) Expiration Date (061215, Start Position 15, Length 6), and (F5) Lot (empty, Start Position empty, Length empty). A 'Format' dropdown menu is set to 'YYYYMM'. An inset image of a keyboard shows the F4 key highlighted. At the bottom of the window are two buttons: 'AddBarcode Rule' and 'Clear Mask'.

Field	Start Position	Length
(F2) Application Identifier	1	2
(F3) Manufacturer Code	3	6
(F4) Expiration Date	15	6
(F5) Lot		

5 – From the drop-down menu, select the proper format for the *Expiration Date*.

This screenshot is similar to the previous one, but the 'Format' dropdown menu for the 'Expiration Date' field is open. The dropdown list shows four options: 'YYYYMM' (highlighted in orange), 'YYMMDD', 'DDMMYY', and 'MMYYYY'. A blue circle is drawn around the 'DDMMYY' option. The rest of the interface, including the text input fields and the table, remains the same as in the previous screenshot.

6 - With the mouse, highlight the characters that define the *Lot* and press the F5 key on the keyboard.

The image shows a software window titled "Barcode Rules". It has two tabs: "Add Barcode Rule" (selected) and "Edit Barcode Rule".

Barcode Mask

Manufacturer Code (Scan Barcode on Box)
AD135790L45670061215 Code 39, Code 32

Field	Value	Start Position	Length
(F2) Application Identifier	AD	1	2
(F3) Manufacturer Code	135790	3	6
(F4) Expiration Date	061215	15	6
(F5) Lot	L45670	9	6


Format: YYYYMM

At the bottom, there are two buttons: "AddBarcode Rule" (circled in blue) and "Clear Mask". A keyboard overlay is visible on the right side of the dialog box.


Press the Add Barcode Rule button to save the rule.

28. HOW TO EDIT BARCODE RULES

In the main menu choose *Manage*, then *Barcode Rules*.



To edit an existing rule, go to the *Edit Barcode Rule* tab and double click the desired rule or select the rule from the grid and press the *Edit Barcode Rule* button.



The *Edit Barcode Rule* window will be opened. Modify the rule as needed and apply the changes by pressing the *Save* button.


Field	Start Position	Length
(F2) Application Identifier		
(F3) Manufacturer Code	1	10
(F4) Expiration Date	11	6
(F5) Lot	18	6

To clear the existing mask, press the *Clear Mask* button.


Note: For further information on how to edit the *Barcode Rule* please read the *HOW TO ADD RULES* section of the manual.

29. HOW TO INACTIVATE (OBSOLETE) BARCODE RULES

In the main menu choose *Manage*, then *Barcode Rules*.



To edit an existing rule, go to the *Edit Barcode Rule* tab and double click the desired rule or select the rule from the grid and press the *Edit Barcode Rule* button.



The *Edit Barcode Rule* window will be opened. Change the rule's state to *Inactive* and apply the changes by pressing the *Save* button.

Edit Barcode Rule

Barcode Mask

Manufacturer Code (Scan Barcode on Box)

1234567890L123456011020

(F2) Application Identifier

(F3) Manufacturer Code

1234567890

(F4) Expiration Date

L12345

(F5) Lot

011020

Start Position

Length

Format

YYYYMM

Barcode Rule Status

☐ Active

☒ Inactive


Save

Clear Mask

Cancel


30. HOW TO TEST BARCODE PARSING

In the main menu choose *Manage*, then *Barcode Rules*.




To test a barcode parsing rule, go to *Edit Barcode Rule* tab and scan a barcode in *Manufacturer Code* text box.

The barcode will be parsed and the respective text boxes will display the data per the configured barcode rule. The rule used for the parsing will be highlighted in the grid.



31. HOW TO ADD A PRINTER

In the main menu choose *Manage*, then *Printers*.



Windows Printer Name

Enter the same name of the printer configured in Microsoft Windows.

Printer Description

Enter a description clear enough to distinguish which printer and label is associated.

NOTE: You can add the multiple printers with different layouts. The printer description must be unique in the example as shown below:

Two physical printers & 3 virtual printers

- Windows Printer Names – Zebra1 and Zebra2
- CLS Printer Descriptions – Zebra1 Round, Zebra1 Rectangle and Zebra2 Flag
- Default Label Layouts – Round (Bottom), Rectangle (Seam) and Flag (Neck) respectively to CLS Printer Descriptions

Choosing the label layout

There are three kinds of label layouts available, enabling you to selectively choose the right one for each type of drug container. The human-readable portion of the label is for identification of the containers to be labeled and the die-cut is for the machine-readable barcode to be applied to the containers.

NOTE: Use only Codonics labels to ensure proper operation of the device and proper labeling of drug containers. Using unapproved labels could lead to unacceptable results, including poor print quality and poor label adhesion to the drug containers. Please contact Codonics Customer Support to order label media for your CLS system.

NOTE: Storage of Codonics labels is recommended to be in a dark place at a relative humidity between 45 percent and 65 percent and a temperature below 77°F (25°C) will assure satisfactory performance.


Flag

Flag labels have this structure. These are best suited for container necks that are ½ inch in diameter or less.




Rectangle

Rectangle labels have this structure.




Round

Round labels have this structure.



Fill in all the fields and press the *Add Printer* button.




The 'Printers' window contains the following fields and controls:


- Add Printer** and **Edit Printer** tabs at the top.
- Windows Printer Name:** A text input field.
- Printer Description:** A text input field.
- Default Label Layout:** A dropdown menu currently showing 'Flag (Neck)'.
- Add** button at the bottom right, circled in blue.

32. HOW TO EDIT A PRINTER

In the main menu choose *Manage*, then *Printers*.



Select the *Edit Printer* tab and double click on the desired printer in the list or select it and press the *Edit Printer* button.




The 'Printers' window shows the 'Edit Printer' tab selected. A table lists the printers, with the first one selected:

Windows Printer Name	Printer Description
ZDesign 5K420d (Copy 1)	Zebra (Rectangle)

At the bottom, there are **Edit** and **Remove** buttons. The 'Edit Printer' tab, the selected printer row, and the 'Edit' button are circled in blue.


The Edit Printer window will be displayed to edit the selected printer.




Make the necessary changes and click the *Save* button to save or the *Cancel* button to cancel the changes.

33. HOW TO DELETE A PRINTER

In the main menu choose *Manage*, then *Printers*.




To remove an existing printer, select the printer you want to delete and press the *Remove Printer* button.




34. HOW TO CALIBRATE A PRINTER


In the main menu choose *Manage*, then *Printers*.



Select the *Edit Printer* tab and double click the desired printer in the list or select it and press the *Edit Printer* button.




Press the *Calibrate* button, which will send a print calibration job to your configured printer.




Close the window by either selecting the *Save* button or the *Cancel* button.

35. HOW TO EDIT YOUR USER ACCOUNT

In the main menu choose *Manage*, then *User Password*.




You are able to change the name of the user in the *Name* field.




After the intended changes are made, press the *OK* button to save or press the *Close* button to cancel.

36. HOW TO CHANGE USER PASSWORD

In the main menu choose *Manage*, then *User Password*.




Enter your current password, then type your new password and confirm the new password in the respective fields. The changes only apply when the current password is correct, and the confirmation of the new password corresponds to the new password.

A screenshot of the 'Edit User' dialog box. It contains several input fields: 'Name' (with 'admin1' entered), 'Username' (with 'admin1' entered), 'Current Password', 'New Password', and 'Confirm Password'. At the bottom, there are two buttons: 'OK' and 'Close'. The 'OK' button is circled with a blue line.


After the intended changes are made, press the *Ok* button to save or press the *Close* button to cancel.

37. HOW TO ADD MANUFACTURERS


In the main menu choose *Manage*, then *Manufacturers*.



To add new manufacturers, enter the Manufacturer ID and Manufacturer Name and press the *Add Manufacturer* button.




Manufacturer IDs must be unique or the CLS will prompt the user with the following message:




Note: Leading or trailing spaces in the *Manufacturer Name* will be truncated.

38. HOW TO EDIT MANUFACTURERS

In the main menu choose *Manage*, then *Manufacturers*.



To edit a Manufacturer, select the manufacturer from the list and press the *Edit Manufacturer* button, or double click on the desired manufacturer.



When editing a manufacturer a new window is displayed.

Edit Manufacturer

Manufacturer Name

Manufacturer ID: 74291829

Manufacturer Name: Abbott Laboratories

Manufacturer Status

☒ Active
☐ Inactive

Barcode Mask for Lot and Expiration Date

Barcode:

	Start Position	Length
(F2) Lot		
(F3) Day		
(F4) Month		
(F5) Year		

Save Clear Mask Cancel

In this window you can change the Manufacturer Name, its Status and define a *Barcode Mask*. The *Barcode Mask* is optional, which is only utilized if the manufacturer includes lot and expiration date data in the box's barcode. The Barcode mask defines this manufacturer barcode structure, and is configured in the following manner:

- 1 – In the Barcode field, scan a Manufacturer Barcode.

Edit Manufacturer

Manufacturer Name

Manufacturer ID: 74291829

Manufacturer Name: Abbott Laboratories

Manufacturer Status

☒ Active
☐ Inactive


Barcode Mask for Lot and Expiration Date

Barcode: 300911GT&fd51


	Start Position	Length
(F2) Lot		
(F3) Day		
(F4) Month		
(F5) Year		

Save Clear Mask Cancel


2- With the mouse, click and drag to select the characters that define the Lot and press the F2 key on the Keyboard.




3- With the mouse, click and drag to select the characters that define the Day and press the F3 key on the keyboard.



4- With the mouse, click and drag to select the characters that define the Month and press the F4 key on the keyboard.



5- With the mouse, click and drag to select the characters that define the Year and press the F5 key on the keyboard.




To save the changes, click the *Save* button, and you will return to the *Edit Manufacturer* window.

Note: The manufacturer barcode parsing does not interact in any way with the barcode parsing based on the *Barcode Parsing Rules*.

The *Manufacturer Barcode Mask* is optional and is only used when the manufacturer provides a specific structure to include lot and expiration date data in the box's barcode.

39. HOW TO ADD A NEW DRUG TO CLS

In the main menu choose *Manage*, then *Drugs*.



To add a new drug, fill in the Container Code or click on the *Auto Assign* button icon to give a system-generated unique ID based upon the *System* settings for *Container ID Formatting*. The *Auto Assign* gets the last Container Code, auto assigned or manual assigned, and increments by one.

Note: If you manually enter a Container Code with the maximum permitted number of digits based upon your system setting configuration (See “How to Configure System Settings”), the *Auto Assign* will not have any additional numbers to automatically assign.


If desired, fill the *Tallman* field for the drug. Enter a description in the *Drug 1* field, select its *Concentration Value* and *Concentration Units*, *Route* and a *Default Printer*.

The Workflow Approval configuration requires administrative approval for each print job for the drug labels.

Only the *Container Code*, *Drug 1*, *Concentration Value*, and *Concentration Unit* fields are required.


Note: *Drug 1* field has a maximum length of 30 characters when *Drug 2* field is empty. *Drug 1* field has a maximum length of 15 characters when *Drug 2* field is not empty.

To save, click the *Add* button.


A screenshot of the 'Drugs' form in the CLS application. The form has a title bar 'Drugs' and two buttons: 'Add Drug' and 'Edit Drug'. The form fields include: 'Container Code' with an 'Auto Assign' button next to it; 'Tallman'; 'Drug 1'; 'Concentration of Drug 1' with a dropdown and a text input; 'Drug 2 (optional)'; 'Concentration of Drug 2 (optional)' with a dropdown and a text input; 'Route' with a dropdown; 'Default Printer' with a dropdown; and a checkbox for 'Workflow Approval'. At the bottom of the form is a large 'Add' button. The 'Auto Assign' button and the 'Add' button are circled in blue.

40. HOW TO EDIT A DRUG

In the main menu choose *Manage*, then *Drugs*.



To edit a drug, first select the drug from the list and click the *Edit Drug* button, or double click on the desired drug.




Make the necessary changes and press the *Save* button or the *Cancel* button to cancel.


The screenshot shows a software window titled "Edit Drug". It contains several input fields and dropdown menus. The "Container Code" field has the value "30". The "Tallman" field is empty. The "Drug 1" field has the value "Aspex". The "Concentration of Drug 1" dropdown is set to "2" and the unit dropdown is set to "%". The "Drug 2 (optional)" field has the value "Aspirin". The "Concentration of Drug 2" dropdown is set to "0,1" and the unit dropdown is set to "mg/L". The "Route" dropdown is set to "INJECTABLE". The "Default printer" dropdown is set to "Flag Printer". There is a checkbox for "Workflow Approval" which is currently unchecked. The "Drug Status" section has two radio buttons: "Active" (which is selected) and "Inactive (Obsolete)". At the bottom of the window, there are two buttons: "Save" and "Cancel". The "Save" button is circled in blue.

41. HOW TO INACTIVATE (OBSOLETE) A DRUG

In the main menu choose *Manage*, then *Drugs*.




To inactivate a drug, select the drug from the list and click the *Edit Drug* button, or double click on the desired drug.



On *Drug Status*, change from *Active* to *Inactive (Obsolete)* and press the *Save* button or the *Cancel* button to cancel.

Note: Drugs marked as *Inactive (Obsolete)* are not available to be selected in the system for a new Manufacturer Drug record. Drugs that are marked as *Inactive (Obsolete)*, remain unchanged as previously assigned to Manufacturer Drug records. If a drug needs to be inactivated for use, please delete the drug from the Manufacturer Drug records (please see the “HOW TO DELETE MANUFACTURER DRUG ASSOCIATION” section of this document).




The screenshot shows the 'Edit Drug' window with the following fields and values:


- Container Code:** 30
- Tallman:** (empty)
- Drug 1:** Aspx
- Concentration of Drug 1:** 2 %
- Drug 2 (optional):** Aspirin
- Concentration of Drug 2:** 0,1 mg/L
- Route:** INJECTABLE
- Default printer:** Flag Printer
- Workflow Approval:** ☐
- Drug Status:**
 - ☒ Active
 - ☐ Inactive (Obsolete)


The 'Inactive (Obsolete)' radio button and the 'Save' button are circled in blue.

42. HOW TO ASSOCIATE DRUGS TO MANUFACTURERS


In the main menu choose *Manage*, then *Manufacturer Drugs*.



To associate a drug to a manufacturer, first select the desired manufacturer with the  button.


A screenshot of the 'Manufacturer Drugs' window. The window has a title bar 'Manufacturer Drugs' and standard window controls. It contains two main sections. The top section, titled 'Manufacturer', has input fields for 'Manufacturer ID:' and 'Manufacturer Name:'. To the right of these fields is a circular button with a magnifying glass icon and the text 'Select Manufacturer'. The bottom section, titled 'Manufacturer Drugs', contains a table with the following columns: 'Manufacturer Code', 'Container Code', 'Tallman', 'Drug', 'Concentration', and 'Verification Status'. Each column has a dropdown arrow. Below the table are three buttons: 'Add', 'Edit', and 'Delete'.

Double click on the desired manufacturer or select it and press the *Select* button.



Now click on the *Add* button to bring the *Add Manufacturer Drugs* window to select the drug to be associated to the select manufacturer.

Note: The drugs already associated to this manufacturer appear in the list below.



Choose the drug from the list and enter a Manufacturer Code by barcode scanning the manufacturer box's barcode. Enter the appropriate number of drug containers that are supplied in the manufacturer's box in the *Box Quantity* field.

The 'Add Manufacturer Drugs' window contains the following data:


Container Code	Tallman	Drug	Concentration
10011		Succinylcholine	20 mg/ml
10012		Imipenem	250 mg
10013		Ketamine	10 mg/ml
10015	ePHEDrine	Ephedrine	5 mg/ml
10016		Atropine	0.4 mg/ml
10017		Rocuronium	10 mg/ml
10018		Esmolol	10 mg/ml
10019		Diazepam	0.08 g/ml

The 'Manufacturer' section includes the following fields:

- Manufacturer ID: 74290811
- Manufacturer Name: B Braun
- Manufacturer Code (Scan Barcode on Box): 124851876
- Box Quantity: 10


The 'Add' button is circled in blue.

To save, click on the *Add* button.



43. HOW TO EDIT MANUFACTURER DRUGS ASSOCIATIONS


In the main menu choose *Manage*, then *Manufacturer Drugs*.




To edit associations of drugs to a manufacturer, first select the desired manufacturer with the *Select Manufacturer* button.

A screenshot of the 'Manufacturer Drugs' window. The window has a title bar 'Manufacturer Drugs' and a search icon labeled 'Select Manufacturer'. Below the search bar are input fields for 'Manufacturer ID' and 'Manufacturer Name'. A table titled 'Manufacturer Drugs' is displayed with columns: 'Manufacturer Code', 'Container Code', 'Tallman', 'Drug', 'Concentration', and 'Verification Status'. Each column has a dropdown arrow. At the bottom of the window are three buttons: 'Add', 'Edit', and 'Delete'.


Double click the desired manufacturer.



Now, double click the desired drug or select it and press the *Edit* button.





Make the necessary changes and press the *Save* button or the *Close* button to cancel.




44. HOW TO DELETE MANUFACTURER DRUGS ASSOCIATIONS


In the main menu choose *Manage*, then *Manufacturer Drugs*.




To delete associations of drugs associated to a manufacturer, first select the desired manufacturer with the *Select Manufacturer*  button.



Double click the desired manufacturer.




Now, select it and press the *Delete* button.




45. HOW TO PRINT A CONTAINER LABEL

In the main menu choose *Label*, then *Print Container Labels*.




You can print labels by scanning the Manufacturer Barcode on the box or by a manual search of the drug description.

Note: Manual Searching is only available to users assigned to the Administrator, Manager or Technician groups. Operators are not able to perform a Manual Search.



Using Manufacturer Barcodes

With the Barcode Search field selected, you can scan the barcode of the manufacturer box to fill the Scan Manufacturer Barcode on Box field to obtain information of the scanned drug.




Print Container Label

Manufacturer Code (Barcode on Box)
2323124344

OK

Barcode Search
Manual Search

Confirm the information displayed and press the *Confirm and Print Container Labels* button.



Print Container Label

Manufacturer Code (Barcode on Box)
4030539043144

OK

Barcode Search
Manual Search

Manufacturer ID:
74290911

Manufacturer Name:
B Braun

Container Code: 10020

Tallman:

Drug 1: **Heparin**

Concentration of Drug 1: **5000 Units**

Drug 2:

Concentration of Drug 2:

Route: IV

☒ Lot and Expiration Date

Lot:
384221


Expiration Date:
July, 2013

Sun	Mon	Tue	Wed	Thu	Fri	Sat
30	1	2	3	4	5	6
7	8	9	10	11	12	13
14	15	16	17	18	19	20
21	22	23	24	25	26	27
28	29	30	31	1	2	3
4	5	6	7	8	9	10

Today: 7/18/2013

Confirm and Print Container Labels


Finally, choose the printer and the quantity of labels that you want to print.




Note: Each print job submitted will include one blank label to be included at the end to allow for ease of separating the printed labels from the label stock. It is recommended that the printed labels be separated by using the perforation on the label media and leaving the blank label connected to the label stock.

Exceptions

If the Manufacturer Code scanned is pending for verification a message box will be displayed indicating that it is not yet approved for use. If the user has user rights for *Perform Verification Workflow*, they will be able to choose to verify the Manufacturer Code.




If the user decides to verify it, a new window will appear asking for a username and password with permissions to *Perform Verification Manufacturers Codes*.



A dialog box titled "Confirm" with a blue header bar. It contains two text input fields: "Username" and "Password". Below the fields are two buttons: "OK" and "Cancel".

If a valid username is entered with Perform Verification Workflow user rights, then it will be possible to Approve or Reject the introduced Manufacturer Code.

In the event the user chooses to reject the Manufacturer Code, it will not be available to print labels. If the user Approves the Manufacturer Code then the window will close and the process will continue.



A dialog box titled "Verify" with a blue header bar. It displays a list of fields with their corresponding values:

Manufacturer Code:	1234567890
Manufacturer ID:	1
Manufacturer Name:	Abbot
Container ID:	34
Tallman:	Wrox
Internal Item:	Wrox
Concentration:	0.05 g
Internal Item 2:	
Concentration 2:	g
Route:	Intramuscular
Created By:	alex
Created On:	31-10-2011 10:09:39
Verification Status:	Select One
Approved/In Use	
Rejected	

Below the list, there is a "Confirm Verification Status" button. A blue circle highlights the "Select One" dropdown menu.

If the manufacturer of this container has Lot and Expiration Date configured in a Barcode Mask, the barcode can be scanned in this screen to verify the information.

Print Container Label

Manufacturer Code (Barcode on Box): 458712 OK Barcode Search Manual Search

Manufacturer ID: 74291029

Manufacturer Name: Abbott Laboratories

Container Code: 10013

Tallman:

Drug 1:

Concentration of Drug 1:

Drug 2:

Concentration of Drug 2:

Route:

☒ Lot and Expiration Date

Lot:

Expiration Date: July, 2013

Sun Mon Tue Wed Thu Fri Sat

1 2 3 4 5 6

7 8 9 10 11 12 13

14 15 16 17 18 19 20

21 22 23 24 25 26 27

28 29 30 31 1 2 3


4 5 6 7 8 9 10

Today: 7/18/2013

Close

Confirm and Print Container Labels

This will check the expiration date of the container. If the container date is expired, the following message will be displayed.



If the code for Lot and Expiration Date is correct, it will display the information for this container.

If this container has workflow approval configured a window for approval will appear, asking a user with permissions to enter their username and password to verify the labels to be printed for this container.

Confirm

Username:

Password:

Approve Reject

With the granted permission, the Print Preview window is displayed.

When there is no Manufacturer Barcode

With the Manual Search field selected, you can search by drug description.

Note: Manual Searching is only available to users assigned to the Administrator, Manager or Technician groups. Operators are not able to perform Manual Searching.

Print Container Label

Enter Drug: 232244522

Barcode Search

Manual Search

OK

Enter part of or the full name of the drug, press the *OK* button and a list of drugs will be displayed in a new window.


Drugs

Add Drug Edit Drug

Container Code	Tallman	Drug	Concentration	Drug 2
10016		Atropine	0.4 mg/mL	
10019		Cilastatin	0.08 g/mL	Diazoxide Meglumine
10015	ePHEDrine	Ephedrine	5 mg/mL	
10007	EPINEPHrine	Epinephrine	0.005 mg/mL	Lidocaine
10018		Esmolol	10 mg/mL	
10006		Fentanyl	50 mg/mL	
10020		Heparin	5000 Units	
10012		Imipenem	250 mg	Cilastatin
10013		Ketamine	10 mg/mL	
10009		Lidocaine	1 %	
10010		Morphine	1 mg/mL	
10008		Propofol	1 %	
10017		Rocuronium	10 mg/mL	
10011		Succinylcholine	20 mg/mL	
10014		Vecuronium	1 mg/mL	

Edit


Double click the desired drug to select it.



The 'Print Container Label' dialog box contains the following fields and controls:

- Enter Drug:** A text box with 'Morphine' and an 'OK' button.
- Search Options:** Radio buttons for 'Barcode Search' and 'Manual Search' (selected).
- Manufacturer ID:** An empty text box.
- Manufacturer Name:** An empty text box.
- Container Code:** A text box with '10010'.
- Tallman:** A text box.
- Drug 1:** 'Morphine'.
- Concentration of Drug 1:** '1 mg/mL'.
- Drug 2:** An empty text box.
- Concentration of Drug 2:** An empty text box.
- Route:** 'INJECTABLE'.
- Lot and Expiration Date:** A checked checkbox.
- Lot:** An empty text box.
- Expiration Date:** A calendar widget showing July 2013, with the 18th selected. Below the calendar, it says 'Today: 7/18/2013'.
- Find Manufacturer:** A magnifying glass icon with the text 'Find Manufacturer'.
- Confirm and Print Container Labels:** A button at the bottom.

You can associate, at this point, the selected drug to a manufacturer, by clicking on the *Find Manufacturer* button where the list of manufacturers appears.




The 'Manufacturers' dialog box displays a table of manufacturers with the following data:

Manufacturer ID	Manufacturer Name	Active
2	Bayer	<input checked="" type="checkbox"/>
10	Dispers	<input checked="" type="checkbox"/>
40	JHC	<input checked="" type="checkbox"/>
3	Pfizer Ltd	<input checked="" type="checkbox"/>
7	Solvay Healthcare	<input checked="" type="checkbox"/>
9	Space Medicals	<input checked="" type="checkbox"/>
8	Univar	<input checked="" type="checkbox"/>

At the bottom of the dialog are two buttons: 'Select' and 'Cancel'. A blue circle highlights the 'Bayer' row in the table.

Select the desired manufacturer and press the Select button or double click it.

Confirm the information displayed and press the *Confirm and Print Container Labels* button.



The 'Print Container Label' dialog box contains the following fields and controls:

- Enter Drug:** Morphine (with an OK button and radio buttons for Barcode Search and Manual Search).
- Manufacturer ID:** 74293311
- Manufacturer Name:** AstraZeneca (with a Select Manufacturer button).
- Container Code:** 10010
- Tallman:**
- Drug 1:** Morphine
- Concentration of Drug 1:** 1 mg/mL
- Drug 2:**
- Concentration of Drug 2:**
- Route:** INJECTABLE
- Lot and Expiration Date:** (checked)
 - Lot:** (empty field)
 - Expiration Date:** July, 2013 (with a calendar view showing the date 7/18/2013 highlighted). Below the calendar, it says 'Today: 7/18/2013'.
- Confirm and Print Container Labels** (button, circled in blue).

Finally, choose the printer and the quantity of labels that you want to print.




The 'Printing' dialog box contains the following fields and controls:

- Container Code:** 10010
- Tallman:**
- Label Preview:** Morphine, 1 mg/mL, Lot: 347123, Exp: 2013/07/18. Below the text are two small images of the label and a circular dashed line.
- Printer:** Zebra (Rectangle) (dropdown menu)
- Quantity:** 1 (spinner control)
- Print Labels** (button)
- Calibrate** (button)
- Cancel** (button)


Note: Each print job submitted will include one blank label to be included at the end to allow for ease of separating the printed labels from the label stock. It is recommended that the printed labels be separated by using the perforation on the label media and leaving the blank label connected to the label stock.

Exceptions

If the container date is expired the next message will show up.




If this container has workflow approval configured a window for approval will appear, asking a user with permissions to enter their username and password to verify the labels to be printed for this container.




46. HOW TO VERIFY A CONTAINER CODE


To check all information on a printed label, go to the main menu choose *Label* then *Verify Container Code*.




Using your barcode scanner, scan the barcode on the CLS printed label.



Check the data displayed on screen.




47. HOW TO VERIFY A MANUFACTURER CODE



If the Verification Workflow is selected in the System Settings, this option allows the user to scan a Manufacturer Code that is pending for Approval, and for user convenience there is a list of Manufacturer Codes that are pending Approval.


The list of Containers is for reference only as you are not able to manually select a container to verify. To start the verification process, please scan a manufacturer's barcode from one of the listed pending codes.

Note: You should have the manufacturer box and container available for review in order to properly verify the manufacturer code and container.



After scanning a Manufacturer Barcode, the following window appears and the user will have the option to Approve or Reject.

To confirm the option made, click on the button Confirm Verification Status.



The Verify window displays the following information:


- Manufacturer Code: 1234567890
- Manufacturer ID: 1
- Manufacturer Name: Abbot
- Container ID: 34
- Tallman: **Wrox**
- Internal Item: **Wrox**
- Concentration: **0.05 g**
- Internal Item 2:
- Concentration 2: **g**
- Route: Intramuscular
- Created By: alex
- Created On: 31-10-2011 10:09:39
- Verification Status: Select One

The dropdown menu for Verification Status is open, showing the following options:


- Select One
- Approved/In Use
- Rejected

The 'Confirm Verification Status' button is highlighted.

CLS informs that the action was successful.



Approving or Rejecting the Manufacturer Code will update the pending list on the Verify Manufacturer Barcode Window.



The Verify Manufacturer Barcode Window displays the following information:


- Manufacturer Code (Barcode on Box):
- OK
- Pending:

Manufacturer ID	Manufacturer	Manufacturer Code	Container Code	Tallman	Drug	Concentration	Created By	Edited By
-----------------	--------------	-------------------	----------------	---------	------	---------------	------------	-----------

If approved, you can now print the labels.

48. HOW TO ANALYZE THE LOGS

In the main menu choose *Logs*, then *Events*.



This allows you to view the events of adding manufacturers, printing labels, resetting the system parameters, adding and editing barcode rules and approving drugs performed by all users.

Add Barcode Rule

Presents a list of barcode rules that were added to CLS, followed by the responsible user, symbology, barcode rule and date/time of the procedure.

Add Manufacturers

Presents a list of manufacturers that were added to CLS, followed by the responsible user and date/time of the procedure.

Edit Barcode Rule

Presents a list of barcode rules that were edited, followed by the responsible user, symbology, barcode rule and date/time of the procedure.

Print Labels

Presents a list of labels that were printed by CLS, containing all information about the label, the number of copies, the responsible user and date/time of the procedure.

Reset

Presents a list of the users and date/time that reset the system parameters.


Drugs

Presents a list of containers that were associated with a lot and validation for printing, followed by the responsible user and date/time of the procedure.

Verification

Presents a list of manufacturer codes that were Approved or Rejected, the manufacturer information, the responsible user and date/time of the procedure.

To see an Event, first select the Procedure, the date's interval and press the *Apply Filter* button.




Pressing the *Export to Excel* button exports the list to an Excel file.

49. HOW TO EXPORT DATA TO A CSV FILE


When a user logs out, if drugs were added or edited, CLS will export this data into a CSV file to the folder specified in *System Settings*.

If you want to export in a middle of a session, go to the main menu and choose *Export* then *Export CLS Data to CSV File*.



All added and edited drugs in this session will show in this window. To export this data, select the intended folder in *File Path* field (the *File Path* from *System Settings* will default) and press the *Export* button.


If all of the information in your CLS database is required to be exported, please contact Codonics Technical Support for use of the *Export All* button.




50. HOW TO USE CLS ON A NETWORK

In Windows Explorer, create in a shared folder on the network, a folder for the exports and a folder to accommodate the database. Copy the database you want to use to this folder.

In CLS, in the main menu choose *Settings*, then *System*.



Here you should change the *File Path* to the export folder you just created and the *Database Path* to the database you just copied.



Repeat the same process in all the CLS clients that are going to share this database.

51. HOW TO SECURE YOUR CLS

Backup

Your CLS data should be backed up as often as possible, ideally after every major session on your PC. Put your CLS folder (defined in the installation) in your backup list. In the event of a serious system crash, you can usually restore the operating system and programs that you use, but your CLS data may be irreplaceable.

There are dozens of programs available to help you to make backups. They vary in price but most have the same basic features - you create a list of the files and folders you want to back up and decide whether you want to save the registry entries. You can either backup your files and folders straight away or set the program to perform the backup automatically later. A basic but useable backup program comes with all Windows operating systems. You can find the Microsoft Backup Utility in the Start menu under Accessories > System Tools.

Never store your backups on your main hard drive. If you do and your hard drive crashes, your backup files will be lost along with everything else. Ideally, backups of your most important data should not be stored in the same room or even on the same site as the original.

Protect against Malware and Intruders

While there is no method that guarantees your computer will not be infected with malware, there are a number of steps that you can take to lessen the probability of this happening.

- Install antivirus and antispyware programs from a trusted source
- Update software regularly
- Use strong passwords and keep them secret
- Never turn off your firewall
- Use flash drives cautiously

To protect your computer, you can install and run an up-to-date antivirus product such as Microsoft Security Essentials, a free solution from Microsoft available for small businesses with up to 10 PCs. If your business has more than 10 PCs, you can protect them with Microsoft Forefront Endpoint Protection. This kind of software provides real-time protection from viruses, spyware, and potentially unwanted software.

It is best practice to run a scan with your antivirus/antispyware product on a regular basis.