

> Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: info@bpsbioscience.com

## **Data Sheet**

# PRMT1 Chemiluminescent Assay Kit

Catalog # 52004L Size: 96 reactions

**DESCRIPTION:** The *PRMT1 Chemiluminescent Assay kit* is designed to measure PRMT1 activity for screening and profiling applications. The *PRMT1 Chemiluminescent Assay Kit* comes in a convenient format, with wells precoated with histone H4 peptide substrate, the antibody against methylated arginine residue of histone H4, the secondary HRP-labeled antibody, S-adenosylmethionine, methyltransferase assay buffer, and purified PRMT1 enzyme for 96 enzyme reactions. The key to the *PRMT1 Chemiluminescent Assay Kit* is a highly specific antibody that recognizes methylated R3 residue of Histone H4. With this kit, only three simple steps are required for methyltransferase detection. First, S-adenosylmethionine is incubated with a sample containing assay buffer and methyltransferase enzyme. Next, primary antibody is added. Finally, the strip plates are treated with an HRP-labeled secondary antibody followed by addition of the HRP substrate to produce chemiluminescence that can then be measured using a chemiluminescence reader.

#### **COMPONENTS:**

Catalog #	Component	Amount	Sto	rage
51040	PRMT1 human recombinant enzyme	10 μg	-80℃	
52120	20 μM S-adenosylmethionine	250 μΙ	-80℃	
52150	Primary antibody 4	100 μΙ	-80℃	
52131H	Secondary HRP-labeled antibody 2	10 μΙ	-80℃	
52170	4x HMT assay buffer 2	3 ml	-20℃	(Avoid
52100	Blocking buffer	50 ml	+4℃	freeze/
	HRP chemiluminescent substrate A	6 ml	+4℃	thaw
	(transparent bottle)			cycles!)
	HRP chemiluminescent substrate B	6 ml	+4℃	
	(brown bottle)			
	96-well plate precoated with histone	1 plate	+4℃	
	substrate			

## MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

TBST buffer (1 x TBS, pH 8.0, containing 0.05% Tween20) Luminometer or fluorescent microplate reader capable of reading chemiluminescence Adjustable micropipettor and sterile tips Rotating or rocker platform

**APPLICATIONS:** Great for studying enzyme kinetics and HTS applications.



Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: info@bpsbioscience.com

**CONTRAINDICATIONS:** DMSO >1%, strong acids or bases, ionic detergents, high salt

**STABILITY:** One year from date of receipt when stored as directed.

REFERENCE: Dillon SC, Zhang X, Trievel RC, Cheng X. Genome Biology 2005; 6:227.

#### ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

#### Step 1:

- 1) Rehydrate the microwells by adding 150 µl of TBST buffer (1 x TBS, pH 8.0, containing 0.05% Tween-20) to every well. Incubate 15 minutes at room temperature. Tap the strip plate onto clean paper towels to remove liquid.
- 2) Thaw **S-adenosylmethionine** on ice. Upon first thaw, briefly spin tube containing **S-adenosylmethionine** to recover full contents of the tube. Aliquot **S-adenosylmethionine** into single use aliquots and store at -80°C. Note: **S-adenosylmethionine** is very sensitive to freeze/thaw cycles. Avoid multiple freeze/thaw cycles.
- 3) Prepare the master mixture: N wells × (7.5 μl **4X HMT assay buffer 2** + 2.5 μl **20** μ**M S-adenosylmethionine** + 15 μl water). Add 25 μl of master mixture to all wells labeled "Positive Control", "Test Sample" and "Blank". For wells labeled "Substrate control", add 7.5 μl **4X HMT assay buffer 2** + 17.5 μl water.

	Blank	Substrate Control	Positive Control	Test Sample
4X HMT assay buffer 2	7.5 µl	7.5 μl	7.5 μl	7.5 µl
20 μM S-adenosylmethionine	2.5 μΙ	-	2.5 μΙ	2.5 μl
H <sub>2</sub> O	15 μl	17.5 μl	15 μl	15 μl
Test Inhibitor	ı	-	1	5 μl
Inhibitor buffer (no inhibitor)	5 μΙ	5 μΙ	5 μΙ	-
1X HMT assay buffer 2	20 μΙ	ı	ı	1
Diluted PRMT1 (0.1-0.5 ng/μl)	-	20 μΙ	20 μΙ	20 μΙ
Total	50 μl	50 μl	50 μl	50 μl

- 4) Add 5 µl of inhibitor solution of each well designated "Test Inhibitor."
- 5) For the "Positive Control," "Substrate Control," and "Blank," add 5 μl of the same solution without inhibitors.



Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: <u>info@bpsbioscience.com</u>

- 6) Thaw **PRMT1** enzyme on ice. Upon first thaw, briefly spin tube containing enzyme to recover full contents of the tube. Aliquot PRMT1 enzyme into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C. Note: PRMT1 enzyme is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.
- 7) Dilute PRMT1 enzyme in **1X HMT assay buffer 2** to 0.1-0.5 ng/µl (2-10 ng/20 µl). Keep diluted enzyme on ice until use. Discard any unused diluted enzyme after use. *Note: optimal enzyme concentration may vary with the specific activity of the enzyme.*
- 8) Add 20 µl of **1X HMT assay buffer 2** to the wells designated "Blank".
- 9) Initiate reaction by adding 20 µl of diluted **PRMT1** enzyme to the wells designated "Positive Control", "Substrate Control", and "Test Sample ". Incubate at room temperature for 20 minutes.
- 10) Remove the supernatant from the wells and wash the strip three times with 200 µl TBST buffer. Blot dry onto clean paper towels.
- 11) Add 100 µl of **Blocking buffer** to every well. Shake on a rotating platform for 10 minutes. Remove supernatant as described above.

#### Step 2:

- 1) Dilute "Primary antibody 4" 100-fold with Blocking buffer.
- 2) Add 100 µl per well. Incubate 1 hour at room temperature with slow shaking.
- 3) Remove the supernatant from the wells and wash the strip three times with 200  $\mu$ l TBST buffer and incubate in **Blocking buffer** as described in steps 1-10 and 1-11.

#### Step 3:

- 1) Dilute "Secondary HRP-labeled antibody 2" 1,000-fold with Blocking buffer.
- 2) Add 100 µl per well. Incubate for 30 minutes at room temperature with slow shaking.
- 3) Remove the supernatant from the wells and wash the strip three times with 200  $\mu$ l TBST buffer and incubate in **Blocking buffer** as described in steps 1-10 and 1-11.
- 4) Just before use, mix on ice 50 μl **HRP chemiluminescent substrate A** and 50 μl **HRP chemiluminescent substrate B** and add 100 μl per well. Discard any unused chemiluminescent reagent after use.



Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: info@bpsbioscience.com


5) Immediately read sample in a luminometer or microtiter-plate reader capable of reading chemiluminescence. "Blank" value is subtracted from all other values.

## **Reading Chemiluminescence:**

Chemiluminescence is the emission of light (luminescence) which results from a chemical reaction. The detection of chemiluminescence requires no wavenlength selection because the method used is emission photometry and is not emission spectrophotometry.

To properly read chemiluminescence, make sure the plate reader is set for LUMINESCENCE mode. Typical integration time is 1 second, delay after plate movement is 100 msec. Do not use a filter when measuring light emission. Typical settings for the Synergy 2 BioTek plate reader are: use the "hole" position on the filter wheel; Optics position: Top; Read type: endpoint. Sensitivity may be adjusted based on the luminescence of a control assay without enzyme (typically we set this value as 100).

## **Example of Assay Results:**



PRMT1 enzyme activity, measured using the PRMT1 Chemiluminescent Assay Kit, BPS Bioscience #52004L. Luminescence was measured using a Bio-Tek fluorescent microplate reader. Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.



Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: info@bpsbioscience.com

### **RELATED PRODUCTS**

<b>Product Name</b>	Catalog #	<u>Size</u>
PRMT1 (expressed in E. coli)	51040	50 μg
PRMT1 (expressed in Sf9 cells)	51041	20 μg
PRMT3 (expressed in E. coli)	51043	50 μg
PRMT4 (expressed in HEK293)	51047	20 μg
PRMT4 (expressed in Sf9 cells)	51044	20 μg
PRMT5 (expressed in HEK293)	51045	20 μg
PRMT5 (expressed in Sf9 cells)	51048	20 μg
PRMT6 (expressed in HEK293)	51046	20 μg
PRMT8 (expressed in Sf9 cells)	51052	20 μg
PRMT3 Chemiluminescent Assay Kit	52005	96 reactions
PRMT4 Chemiluminescent Assay Kit	52041L	96 reactions
PRMT5 Chemiluminescent Assay Kit	52002	96 reactions
PRMT6 Chemiluminescent Assay Kit	52046	96 reactions
Histone H4(R3) Universal Assay Kit	52074	96 reactions
PRMT1 Homogeneous Assay Kit	52052	384 reactions
PRMT3 Homogeneous Assay Kit	52055	384 reactions
PRMT5 Homogeneous Assay Kit	52054	384 reactions
PRMT6 Homogeneous Assay Kit	52056	384 reactions
PRMT8 Homogeneous Assay Kit	52058	384 reactions



Tel: 1.858.829.3082 Fax: 1.858.481.8694

Email: info@bpsbioscience.com

#### TROUBLESHOOTING GUIDE

Problem	Possible Cause	Solution		
Luminescence signal of positive control reaction is weak	PRMT1 enzyme has lost activity	Enzyme loses activity upon repeated freeze/thaw cycles. Use fresh enzyme (PRMT1, BPS Bioscience #51040). Store enzyme in single-use aliquots. Increase time of enzyme incubation. Increase enzyme concentration.		
	Antibody reaction is insufficient	Increase time for primary antibody incubation. Avoid freeze/thaw cycles of antibodies.		
	Incorrect settings on instruments	Refer to instrument instructions for settings to increase sensitivity of light detection. See section on "Reading Chemiluminescence" above.		
	Chemiluminescent reagents mixed too soon	Chemiluminescent solution should be used within 15 minutes of mixing. Ensure both reagents are properly mixed.		
Luminescent signal is erratic or varies widely among wells	Inaccurate pipetting/technique	Run duplicates of all reactions. Use a multichannel pipettor. Use master mixes to minimize errors.		
	Bubbles in wells	Pipette slowly to avoid bubble formation. Tap strip lightly to disperse bubbles; be careful not to splash between wells.		
Background (signal to noise ratio) is high	Insufficient washes	Be sure to include blocking steps after wash steps. Increase number of washes. Increase wash volume. Increase Tween-20 concentration to 0.1% in TBST.		
	Sample solvent is inhibiting the enzyme	Run negative control assay including solvent. Maintain DMSO level at <1% Increase time of enzyme incubation.		
	Results are outside the linear range of the assay	Use different concentrations of enzyme (PRMT1, BPS Bioscience #51040) to create a standard curve.		