# RayBio®Phospho Stat1 (Ser727) and Pan Stat1 ELISA Kit

For Measuring Phospho-Stat1 (Ser727) and Pan Stat1 in Human and Mouse Cell Lysates

User Manual (Revised Mar 1, 2012)

RayBio<sup>®</sup> Phospho-Stat1 (Ser727) and Pan Stat1 ELISA Kit Protocol

(Cat#: PEL-Stat1-S727-T)



We Provide You With Excellent Protein Array System And Service

Tel:(Toll Free)1-888-494-8555 or 770-729-2992; Fax:770-206-2393;

Web: www.raybiotech.com Email: info@raybiotech.com



# RayBio® Phospho-Stat1 (Ser727) and Pan Stat1 ELISA Kit Protocol

# **TABLE OF CONTENTS**

I.	Introduction	2
II.	Material Provided	3
III.	Storage	3
IV.	Additional Materials Required	4
V.	Sample Preparation	4
VI.	Reagent Preparation	5
VII.	Assay Procedure.	7
VIII	. Assay Procedure Summary	9
IX.	Typical Data	10
	i. Positive Control.	10
i	i. Sensitivity	11
ii	i. Recombinant Human EGF Stimulation of	
	A431 Cell Lines	12
X.	References.	13
XI.	Troubleshooting Guide	14

#### I. INTRODUCTION

RayBio® Phospho-Stat1 (Ser727) and Pan Stat1 ELISA (Enzyme-Linked Immunosorbent Assay) kit is a very rapid, convenient and sensitive assay kit that can monitor the activation or function of important biological pathways in cell lysates. By determining phosphorylated Stat1 protein in your experimental model system, you can verify pathway activation in your cell lysates. You can simultaneously measure numerous different cell lysates without spending excess time and effort in performing a Western Blot analysis.

This Sandwich ELISA kit is an in vitro enzyme-linked immunosorbent assay for the measurement of human phospho-Stat1 (Ser727) and pan Stat1 (help normalize the results of phospho-Stat1 from different cell lysate being compared). An anti-Stat1 (Ser727) (half plate, red marker on left side) and anti-pan Stat1 antibody (half plate, black marker on right side) has been coated onto a 96-well plate. Samples are pipetted into the wells and phosphorylated (left side) and pan (right side) Stat1 present in a sample is bound to the wells by the immobilized antibody. The wells are washed and biotinylated Stat1 is used to detect phosphorylated or pan Stat1. washing away unbound antibody, HRP-conjugated After Streptavidin is pipetted to the wells. The wells are again washed, a TMB substrate solution is added to the wells and color develops in proportion to the amount of Stat1 (Ser727) or pan Stat1 bound. The Stop Solution changes the color from blue to yellow, and the intensity of the color is measured at 450 nm.

#### II. MATERIAL PROVIDED

- 1. Stat1 Microplate (Item A): 96 wells (12 strips x 8 wells) coated with anti-phospho-Stat1 (Ser727) (half plate, red marker on left side) and anti-Stat1 antibody (half plate, black marker on right side).
- 2. Wash Buffer Concentrate (20x) (Item B): 25 ml of 20x concentrated solution.
- 3. Assay Diluent (Item E2): 15 ml of 5x concentrated buffer. For diluting cell lysate sample, detection antibody (Item C) and HRP-Streptavidin Concentrate (Item G).
- 4. Detection Antibody Stat1 (Item C): 2 vial of biotinylated anti-Stat1 (each vial is enough to assay half microplate).
- 5. HRP-Streptavidin Concentrate (Item G): 2 vials, 200 µl/vial, HRP-conjugated streptavidin concentrate.
- 6. TMB One-Step Substrate Reagent (Item H): 12 ml of 3,3',5,5'-tetramethylbenzidine (TMB) in buffered solution.
- 7. Stop Solution (Item I): 8 ml of 0.2 M sulfuric acid.
- 8. Cell Lysate Buffer (Item J): 5 ml 2x cell lysis buffer (not including protease and phosphatase inhibitors).
- 9. Positive Control A431S002-1 (Item K): 1 vial of lyophilized powder from A431 cell lysate.

### III. STORAGE

Upon receipt, the kit should be stored at -20 °C. Please use within 6 months from the date of shipment. After initial use, Wash Buffer Concentrate (Item B), Assay Diluent (Item E2), TMB One-Step Substrate Reagent (Item H), HRP-Streptavidin (Item G), Stop

Solution (Item I) and Cell Lysate Buffer (Item J) should be stored at 4 °C to avoid repeated freeze-thaw cycles. Return unused wells to the pouch containing desiccant pack, reseal along entire edge and store at -20 °C. Reconstituted Positive Control (Item K) should be stored at -70 °C.

# IV. ADDITIONAL MATERIALS REQUIRED

- 1 Microplate reader capable of measuring absorbance at 450 nm.
- 2 Protease and Phosphatase inhibitors.
- 3 Shaker.
- 4 Precision pipettes to deliver 2 µl to 1 ml volumes.
- 5 Adjustable 1-25 ml pipettes for reagent preparation.
- 6 100 ml and 1 liter graduated cylinders.
- 7 Distilled or deionized water.
- 8 Tubes to prepare sample dilutions.

# V. SAMPLE PREPARATION

Cell lysates - Rinse cells with PBS, making sure to remove any remaining PBS before adding the Cell Lysate Buffer. Solubilize cells at 4 x 10<sup>7</sup> cells/ml in 1x Cell Lysate Buffer (we recommend adding protease and phosphatase inhibitors to Cell Lysate Buffer prior to sample preparation). Pipette up and down to resuspend and incubate the lysates with shaking at 2 - 8° C for 30 minutes. Microcentrifuge at 13,000 rpm for 10 minutes at 2 - 8° C, and transfer the supernates into a clean test tube. Lysates should be used immediately or aliquoted and stored at -70 °C. Avoid repeated

freeze-thaw cycles. Thawed lysates should be kept on ice prior to use.

For the initial experiment, we recommend to do a serial dilution testing such as 5-fold and 100-fold dilution for your cell lysates with Assay Diluent (Item E2) before use.


Note: The fold dilution of sample used depends on the abundance of phosphorylated proteins and should be determined empiricallys. More of the sample can be used if signals are too weak. If signals are too strong, the sample can be diluted further.

Cell Lysate Buffer should be diluted 2-fold with deionized or distilled water before use (recommend to add protease and phosphatase inhibitors).

# VI. REAGENT PREPARATION

- 1. Bring all reagents and samples to room temperature (18 25°C) before use.
- 2. Item E2, Assay Diluent should be diluted 5-fold with deionized or distilled water before use.
- 3. Preparation of Positive Control: Briefly spin the Positive Control vial of Item K. Add 700 µl 1x Assay Diluent (Item E2, Assay Diluent should be diluted 5-fold with deionized or distilled water before use) into Item K vial to prepare a Positive Control (P-1).

Dissolve the powder thoroughly by a gentle mix (it can be removed by centrifuge if any precipitate in the solution is found). Pipette 260 µl 1x Assay Diluent into each tube. Use the P-1 to produce a dilution series (shown below). Mix each tube thoroughly before the next transfer. 1x Assay Diluent serves as the background.



- 4. If the Wash Concentrate (20x) (Item B) contains visible crystals, warm to room temperature and mix gently until dissolved. Dilute 20 ml of Wash Buffer Concentrate into deionized or distilled water to yield 400 ml of 1x Wash Buffer.
- 5. Briefly spin the biotinylated antibody (Item C) before use. Add 100 µl of 1x Assay Diluent into the vial to prepare a biotinylated antibody concentrate. Pipette up and down to mix gently (the concentrate can be stored at 4°C for 5 days or at 80°C for one month). The biotinylated Stat1 antibody should be diluted 55-fold with 1x Assay Diuent and used in step 4 of Part VII Assay Procedure.
- 6. Briefly spin the HRP-Streptavidin concentrate vial (Item G) and pipette up and down to mix gently before use. HRP-Streptavidin

concentrate should be diluted 40 fold or 200 fold (see "step 6" in page 8 for detail) with 1x Assay Diluent.

For example: Briefly spin the vial (Item G) and pipette up and down to mix gently. Add 150µl of HRP-Streptavidin concentrate into a tube with 6 ml 1x Assay Diluent to prepare a 40-fold diluted HRP-Streptavidin solution (don't store the diluted solution for next day use). Mix well.


7. Cell Lysate Buffer should be diluted 2-fold with deionized or distilled water before use (recommend to add protease and phosphatase inhibitors).

#### VII. ASSAY PROCEDURE:

- 1. Bring all reagents to room temperature (18 25°C) before use. It is recommended that all samples or Positive Control should be run at least in duplicate.
- 2. Add 100 µl of each sample or positive control into appropriate wells (see the following 96 well microplate formate). Cover well with plate holder and incubate for 2.5 hours at room temperature or over night at 4°C with shaking.

96 well microplate coated with phosphorylated and pan antibodies:

# Anti-Stat1 (Ser 727) Anti-pan Stat1



- 3. Discard the solution and wash 4 times with 1x Wash Solution. Wash by filling each well with Wash Buffer (300 µl) using a multi-channel pipette or autowasher. Complete removal of liquid at each step is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating or decanting. Invert the plate and blot it against clean paper towels.
- 4. Add 100 µl of prepared 1x biotinylated Stat1 antibody (Reagent Preparation step 5) to each well. Incubate for 1 hour at room temperature with shaking.
- 5. Discard the solution. Repeat the wash as in step 3.
- 6. Add 100 µl of 40 fold diluted HRP-Streptavidin solution (see Reagent Preparation step 7) to each well coated with anti-Stat1 (Ser 727) on the left side (red marker, see Assay Procedure step1). Add 100 µl of 200 fold diluted HRP-

Streptavidin solution to each well coated with anti-pan Stat1 on the right side (black marker). Incubate for 1 hour at room temperature with shaking.

- 7. Discard the solution. Repeat the wash as in step 3.
- 8. Add 100 µl of TMB One-Step Substrate Reagent (Item H) to each well. Incubate for 30 minutes at room temperature in the dark with shaking.
- 9. Add 50 µl of Stop Solution (Item I) to each well. Read at 450 nm immediately.

#### VIII. ASSAY PROCEDURE SUMMARY

1. Prepare all reagents, samples and standards as instructed.

 $\int$ 

2. Add 100 μl sample or positive control to each well. Incubate 2.5 hours at room temperature or over night at 4°C.

 $\int$ 

3. Add 100 µl prepared biotinylated primary antibody to each well. Incubate 1 hours at room temperature.

 $\int$ 

4. Add 100 μl prepared 1X HRP-Streptavidin solution. Incubate 1 hour at room temperature.

IJ

5. Add 100 µl TMB One-Step Substrate Reagent to each well. Incubate 30 minutes at room temperature.

 $\int$ 


6. Add 50 µl Stop Solution to each well. Read at 450 nm immediately.

#### IX. TYPICAL DATA

ELISA data analysis: Average the duplicate readings for each sample or positive control.

#### i. Positive Control

A431 cells were treated with recombinant human EGF at 37°C for 20 min. Solubilize cells at 4 x 10<sup>7</sup> cells/ml in lysis buffer. Serial dilutions of lysates were analyzed in this ELISA. Please see step 3 of Part VI. Reagent Preparation for detail.




Positive control dilution series


# ii. Recombinant Human EGF Stimulation of A431 Cell Lines

A431 cells were treated or untreated with 100 ng/ml recombinant human EGF for 10 min. Cell lysates were analyzed using this phosphoELISA and Western Blot.

# A). ELISA




# B). Western-Blot Analysis




#### iii. SENSITIVITY

The A431 cells were treated with 100 ng/mL recombinant human EGF for 20 minutes to induce phosphorylation of Stat1. Serial dilutions of lysates were analyzed in this ELISA and by Western blot. Immunoblots were incubated with anti-phospho-Stat1 (Ser 727).

# A). ELISA



# B). Western-Blot Analysis



# X. REFERENCES:

- 1. Wen, Z. et al. (1995) Cell. 82: 241-250.
- 2. Zhang, X., et al. (1995) Science. 267: 1990-1994.
- 3. Hackel, P.O. et al. (1999) Curr. Opin. Cell Biol. 11, 184-189.
- 4. Cooper, J.A. and Howell, B. (1993) Cell 73, 1051-1054.

# XI. TROUBLESHOOTING GUIDE

Problem	Cause	Solution
1. Sample signals:	Cause	Solution
a. Too low	Sample concentration is too low	a. Increasing sample concentration
b. Too high	<ul><li>b. Sample concentration is too high</li></ul>	b. Reducing sample concentration
2. Large CV	a. Inaccurate pipetting	a. Check pipettes
3. High background	a. Plate is insufficiently washed	a. Review the manual for proper washing. If using an automated plate washer, check that all ports are unobstructed.
	<ul><li>b. Contaminated wash buffer</li></ul>	<ul><li>b. Make fresh wash buffer</li></ul>
4. Positive Control: Low signal	a. Improper storage of the ELISA kit	<ul> <li>a. Upon receipt, the kit should be stored at -20 °C. Store the positive control at -70°C after reconstitution.</li> </ul>
	b. Stop solution	<ul> <li>b. Stop solution should be added to each well before measurement and read OD immediately.</li> </ul>
	<ul><li>c. Improper primary or secondary antibody dilution</li></ul>	c. Ensure correct dilution

# RayBio® ELISA kits:

Over 200 ELISA kits, custom ELISA kit choose from over 300 list visit <a href="https://www.raybiotech.com">www.raybiotech.com</a> for details.

RayBiotech, Inc., the protein array pioneer company, strives to research and develop new products to meet demands of the biomedical community. RayBio's patent-pending technology allows detection of over 180 cytokines, chemokines and other proteins in a single experiment. Our format is simple, sensitive, reliable and cost effective. Products include: Cytokine Arrays, Chemokine Arrays, ELISA kits, Phosphotyrosine kits, Recombinant Proteins, Antibodies, and custom services.

#### **Antibody Array**

Cytokine Antibody Array: Simultaneous detection up to 200 proteins (cytokine, chemokine, growth factor, adipokine, angiogenic factor, protease) in one experiment

#### **Phosphorylation Antibody Array**

- RTK antibody array
- EGFR phosphorylation antibody arrays

Label based antibody array: Simultaneous detection more than 500 proteins in one experiment

Quantibody Array: Quantitative measurement of multiple protein levels

### **Protein Array**

**ELISA** 

**Cell-Based Phosphorylation ELISA** 

**Tissue MicroArray** 

Protein: Cytokine, Chemokine, Adiplokine, Angiogenic factor, Virus, bacteria and infectious disease protein, hormone, Enzyme, other

#### **Peptide**

Antibody: Cytokine, Adipokine, Angiogenic factor, Signal transduction,
Transcription factor, Receptor, Adhesion molecule, Virus, bacteria and other
infectious agents, Secondary antibody, Tag antibody, Immunoglobulin,
Hormone, Cell surface, Protease, other

Antibody array, Protein array, Peptide array, ELISA, Phosphorylation assay Tissue array

Assay service: just simply send your samples and get data in 1 to 2 weeks.

Antibody array, Protein array, ELISA, Quantibody array

Antibody production: highest quality with very competitive price
Monoclonal antibody, Recombinant antibody, Polyclonal antibody, Phase
display, Antibody angineering, Antibody conjugation

**Recombinant protein production** 

**Assay development** 

**Array printing** 

Contact and non-contact arrayers. All kinds of substrates of your choice including glass slides, membranes and plates.

This product is for research use only.

