

USER MANUAL



Supported products

Software	PistonXP version 1.40
PDD-301/sh	Spirometer
PDD-301/r	Rhinomanometer
PDD-301/sco	Breath CO monitor and spirometer
PDD-301/rco	Breath CO monitor and rhinomanometer
PDT-111/p	Body Plethysmograph
PDT-111/d	CO-Diffusion
PDT-111/pd	Body Plethysmograph and CO-Diffusion

Piston Ltd.
1033 Budapest, Szőlőkert u. 4/b

VOLUMES INCLUDED

This User Manual includes the following volumes

PC Software

DB

Installation
Settings
Maintenance
Patient's database
Common tasks

Spirometry

SR

Installation
Use
Maintenance
Troubleshooting

TGV and CO-Diffusion test

PD

Installation
Use
Maintenance
Troubleshooting

Appendix

AX

Certificates
Reference tables

DB

PC SOFTWARE



Supported product

PistonXP version 1.40

Piston Ltd.
1033 Budapest, Szőlőkert u. 4/b

CE 1979

Rev.: DB-EN-04.00

Date modified: 04/07/2012

TABLE OF CONTENT

Installation.....	3
Electric shock protection	3
Minimum PC configuration.....	4
Compatibility	5
Software installation	5
User interface.....	10
Icons.....	10
Settings	13
Institute data	13
Doctor's data	13
Language selection.....	14
Patient identification format.....	14
Graph settings.....	14
Curve magnification.....	15
Reference values	15
Displayed parameters.....	16
Automatic backup	16
Patient database.....	17
User interface overview	17
Patient's personal data	19
Finding a patient in the database.....	19
Viewing previous measurements	20
Comment field for patients.....	21
Measurements.....	22
Patient selection	22
Preparations	22
Calibration	22
Measurement	22
Enter comment	23
Store	23
Printing	23
Export report into PDF and other graphical formats	24
Interface to information systems.....	26
PRE/POST.....	27
The PRE/POST measurement.....	27
Retrieve measurement.....	27
Printing	27
Troubleshooting.....	28
Possible problems	28


Electric shock protection



The electric shock protection instructions in this section must be followed!

Only Piston Ltd., as manufacturer, or its authorized distributor's personnel, or the distributor's representatives may install the medical device. The above mentioned companies only accept responsibility for systems installed by them.

Before installing the medical devices the personnel must make sure the computer, the monitor and the printer installed as medical electronic devices comply with the standards, for the given country or the user declares concerning this with responsibility.



Information exchange with the computer goes through USB connection. For low leakage current relating to medical devices standards this connection is optically isolated inside the device.

Before shipping we check the device's leakage current. The operator has the opportunity to have the leakage current checked periodically, if he / she finds it necessary.

The system must be installed so the examined person is at least 1.5m away from those devices that are electrically connected to the computer equipment.

Parts of the system (computer, monitor, printer) can only be replaced in case of failure, modification or for any other reason, if the part to be installed has the same electric shock protection conditions as the original one.

The personnel installing the device will train the operator concerning operation electric shock protection. This training includes the contents of this section. The operator verifies the training in official written form.

Minimum PC configuration

The operation of the lung diagnostics system requires a personal computer with the following minimum configuration:

Item	Minimum*	Recommended**
Operating system	Windows XP SP3	
Processor (for PDD-301, PDD-401 device family)	800 MHz clock frequency	Intel Pentium 3 / Celeron / Core 2 család AMD Athlon család
Processor (for PDT-111 device family)	1 GHz clock frequency	Intel Celeron / Pentium 4 / Dual Core / Core 2 család AMD Athlon család
Processor (for PRE-101 device family)	1.5 GHz clock frequency	Intel Core 2 család, Core i3 / i5 / i7 AMD Sempron / Phenom / Athlon család
Memory	256 MB	1 GB for Windows XP 2 GB for Windows 7
Memory (for PRE-101 device family)	512 MB	
Hard drive	1.5 GB free space for software system and database management system	4 GB more space for user data
Screen resolution	1024×768	1280×1024
Screen resolution (for PRE-101 device family)	1280×1024, two monitors required	
Printer	Windows compatible	Color for PRE-101 device family
Internet connection		for Software update

* Minimal configuration is the theoretically required minimum for running the operating system and associated services

** Recommended configuration is the required minimum for fluent daily work. Before purchasing a workstation is recommended to consider increasing performance need of later expansion of operating system and other system services.

Compatibility

The medical diagnostic software suite is a PC based Microsoft Windows compatible system that was tested under the following versions of Windows:

32 bit

- Windows XP Home, Professional SP2, SP3
- Windows Vista Home Basic
- Windows Vista Home Premium
- Windows Vista Business
- Windows Vista Enterprise
- Windows Vista Ultimate
- Windows 7 Home Basic
- Windows 7 Home Premium
- Windows 7 Professional SP1
- Windows 7 Enterprise
- Windows 7 Ultimate SP1

64 bit


- Windows XP Professional x64 Edition SP2
- Windows Vista Home Basic
- Windows Vista Home Premium
- Windows Vista Business
- Windows Vista Enterprise
- Windows Vista Ultimate
- Windows 7 Home Basic
- Windows 7 Home Premium
- Windows 7 Professional SP1
- Windows 7 Enterprise
- Windows 7 Ultimate SP1
- Windows Server 2008 R2 SP1

Software installation


Perform the installation from the included CD.

The most up-to-date version is available from our website:

<http://www.pistonmedical.com> In the Downloads / Software section.



Click on the **Start** menu and select **Run**




Click on the **Browse** button and select the install program.

When installing from the CD select the CD drive.

Find the pxp_setup.exe file in the Programs folder.


Click **OK**




The install program starts

Select the preferred language for the setup and the installed software

Click **OK**




A welcome screen appears, just click **Next**




Carefully read the License Agreement, click **I accept the agreement** and click **Next**

If you do not accept the agreement, exit the installation




You can specify the install destination.

Click **Next**




You can select which part of the program to install (experienced users)

Click **Next**




You can enter the name the program appears under in the Start menu (experienced users)

Click **Next**




You can select whether a PistonXP icon should be created on the desktop (experienced users)


Click **Next**




An install summary window appears, and if all settings are acceptable
Click **Install**



The install process begins
Please wait until it finishes
installing the software
After installing the software,
external components will be
installed also



The installation of the Oracle XE
Database Server / Client runs in
background and the process can
take several minutes
Please wait until it finishes
installing the software




Finally USB Drivers are being
installed
Please wait until it finishes
installing the software

**In case of WindowsXP**

If the windows generates a warning click on the [**Continue Anyway**] button

**In case of Windows Vista and 7**

If the windows generates a warning click on the [**Install this driver software anyway**] button



A window indicates the end of the installation

Click **Finish** to close the install wizard

This concludes installation

Start the program

The program automatically detects the connected devices

Icons

Main window



Open patient database



Open expertise editor



Submenu for Hospital information systems



Open report editor and printing



Exit software

Patient database



Clear patient quick search fields



Enter new patient



Modify patient data



Store entered / modified data



Cancel changes



Load all measurements from the selected meeting(s)



Load selected measurements

Settings



Set institute data



Doctor records



Devices' settings connected to the PC



Program operation related settings



Display graphs and other program parts



Maintenance, safety backup related settings



Reference value calculating algorithms



List of parameters to be displayed



Service panels



Enter new doctor



Modify doctor's data



Store entered / modified data



Cancel changes



Accept changes and close Options panel

Report editor



Print preview for lung function tests



Print preview for Compliance test



Print preview for Rhinomanometry



Print preview for Provocation test



Print preview for Ergospirometry



Print preview for Audiometry



Print preview for Calibration



Print selected measurement results



Store a report as PDF document or image




Close Report editor

Settings

The **Setup/Options** menu item allows customization of the system.

Settings that can be changed during measurement are also available in the **Setup** tab of the measurement windows.

Program settings appear grouped on the left side.



Institute data

You can enter the following information at the **Setup/Options/Institute setup** menu item:

Institute name, Site address, Mailing address, Phone number, Fax number, Web page, E-mail address.

This data appears in the header of the printed report.

Doctor's data

The doctor's data can be entered at the **Setup/Options/Doctors** menu item.

New doctor

Press the **[New Doctor]** button to enter data for a new doctor.

Complete the fields.

Make sure that two doctors cannot have the same identifier.

Press the **[Save]** button to store the entered data.

Modify data

Select the doctor from the **[Doctor's name]** drop down list whose data you would like to modify.

Click the **[Modify]** button.

Change the desired fields.

When done, press the **[Save]** button.

You will see feedback about the success of the data storage.

If you do not wish to store the entered data, press the **[Discard]** button.

About deleting ...

To preserve consistency and for future searches, it is not possible to delete from the database.

All diagnosis has traces in the database.

Language selection

You can select the program's language in the **Setup/Options/Operation** menu item.

All supported languages are displayed in English and in the specific language as well.

Select the language you would like to use.

Patient identification format

You can enter the patient identification format in the **Setup/Options/Operation** menu item.

Format descriptions may be found in the Appendix I. section.

Graph settings

Graph displays may be set in the **Setup/Options/Display** menu item.

Graph scheme

You can select the graph colour settings:

- Dark background, bright lines
- Bright background, dark lines
- Same as printed (white background)

Raster

The grid may be enabled or disabled on the graph

Show curves

It can be selected for several same type measurements:

- The diagrams appear in one coordinate system.
- All the diagrams appear in different coordinate systems.

Visible part of the curve

For easier overview curve sections unrelated to the evaluation can be hidden.

OnFly Analysis

When this function is enabled, the program monitors the patient's breathing during measurement, separates normal breathing from deep exhalations and inhalations.

Active curves after measure

In the **Setup/Options/Operation** menu item those curves can be selected which will be automatically indicated as active ones after each measurement:

- Just the best measure
- First three
- All measurements

Curve magnification

Click on any graph with the right mouse button.

Select the required size from the menu that appears.

The following magnifications are available:

Resistance measurement: 0.5×-, 1×-, 2×-

Other measurements: 0.5×, 1×, 2×, 3×, 5×

Reference values

The desired algorithm may be selected in the **Setup/Options/Prediction** menu item:

- ECCS – European Community for Coal and Steel
- Knudson
- Cotton & Dust
- Crapo-HSU
- Österreichisch

To turn it off, select the **No reference** values option.

Displayed parameters


In the **Setup/Options/Parameters** menu you can enter which parameter to display on the screen and which one to print.

Automatic backup

In the **Setup/Options/Maintenance** menu you can select the frequency and the location of automatic backup of the patient database.

User interface overview

Main window



Quick search

Helps find a patient.

Patient list

A list of patients meeting the search criteria.

Details

Displays the selected patient's most important parameters for the selected measurement.

Control panel

Basic database operations: enter new patient, modify patient data, store.

Visits

Dates of previous visits.

Measurements

A list of measurements for the selected date or measurement type.

Load

Control buttons to display the selected measurements.

Measurement selection

Measurements may be listed according to measurement type as well.

Data input form

Content of the Data input form can be set in the menu **Setup / Options / Display / Contents of Patient's Datasheet***

The screenshot shows the 'Patient database' data input form. On the left, there is a vertical toolbar with buttons: 'Diagnose', 'Spirometry', 'Audiometry', 'Print...', 'LINK', and 'Exit'. The main form area is divided into several sections. 'Identification data' points to the top section containing fields for Title, Last name, First name, Sex (Female/Male), Born (date), ID, and Mother's name. 'Contact information' points to the section below containing Postal code, City, Street, Phone, Mobile, and Mail. 'Body mass index' points to the section containing Height (cm) and Weight (kg), with a calculated 'Body Mass Index' and 'Age' (58 years) displayed. 'List of incomplete fields' points to a list box showing 'First name', 'Last name', 'Mother's name', 'Birth date', and 'Sex'. 'Control panel' points to the bottom section containing 'New Patient', 'Modify', 'Save', and 'Cancel' buttons. At the very bottom, there is a status bar showing 'IVC 4.36' and 'FVC 3.99'.

Identification data

Group of data essentially identifying the patient: Name, date of birth, social security number, sex, etc.

Contact information

Patient's accessibility: Address, phone numbers, e-mail address.

Body mass index (calculated value)

The patient's current body weight index: square of the height of the patient in meter divided by body weight

List of incomplete fields

A list of fields those both have to be completed and are still empty, or that have been filled out incorrectly.

Control panel

Basic database operations: new patient, modification, store.

*You have to close and re-open Patient's Database to apply changes

Patient's personal data

The program can store an arbitrary number of patients.

Pink fields indicate fields that have to be completed.

Anthropometrics data

You have to enter the patient's body mass and height

These data are required to calculate reference value

The database stores the body mass and height of the patient for each visit, so changes may be followed in time.

New patient

To enter a new patient, press the **[New Patient]** button. Complete the fields and make sure that two patients cannot have the same identifier.

To store the patient, press the **[Save]** button.

You will receive feedback about the success of data storage.

If you do not wish to save the data, press the **[Cancel]** button.

Modify data

Select the patient to modify

Click on the **[Modify]** button

After modification press the **[Save]** button

You will receive feedback about the success of data storage

If you do not wish to save the modified data, press the **[Cancel]** button

About deleting ...

To preserve consistency and for future searches, it is not possible to delete from the database. All diagnosis has traces in the database.

Finding a patient in the database

The top section of the patient database window is the search block.

You can search based on several criteria. When those criteria change, the program automatically lists the patients meeting the updated criteria

Normal search

Search only based on the patient's family and surname.

Enter the patient's name or part of it.

Detailed search

Click the **[Detailed search]** button.

You can refine the search criteria in the window:

- patient's sex
- date of birth with interval
- address or part of it
- doctor
- identifier (social security number)

Viewing previous measurements

All previous measurements can be reloaded, so reports can be printed at anytime.

Viewing previous measurements

To reload previous measurements:

- Select the patient
- Select the visit by date
- If you only wish to view the results of certain measurement mode select the one from the list
- Select required measurements

If you wish to see all measurement results of a selected visit, click the **[Open]** button

If only certain measurements are important; click them while holding the CTRL button down

After selection click on the **[Open]** button!

If you wish to include further measurements to the report click the **[Patient database]** button to reopen the Patient database.

Select further measurements and click the **[Add measure]** button to include them to the report.

**WARNING:**

You can only simultaneously load eight measurements of the same mode.

For this reason, the [Load all] button is not always available.

PRE/POST evaluation

To load the data for all previous visits, check the **[All measurements]** checkbox.

This displays a patient's all previous measurements sorted according to the following:

- Date
- Measurement mode
- Measurement results quality

Select the results of at least two identical measurement mode, for example two FVC measurements.

Load the data as mentioned earlier.

PRE/POST measurements are detailed in the PRE/POST section.

Comment field for patients

Comments may be entered about the patients even for every visit. All comments are stored separately in the database and may be retrieved individually.

To enter a comment:

- Open the Patient database
- Selected the desired patient
- Click the [Diagnose] button to open the text editor window
- Select the [Patient] operating mode from the list
- Enter the comment
- Press the [Store] button to store the comment

Previous diagnosis

All previous comments about the patient may be retrieved from the **[History]** list.

The currently entered text is not lost when viewing a previous diagnosis.

To display the currently entered text again, select the **[Patient]** option from the list of operating modes again.

Patient selection

Before starting the measurement it is necessary to enter patient data using one of the following methods:

- Enter new patient
- Search for patient already in the database

Preparations

Device

Connection

Make sure that the device you wish to use is connected to the computer.

If not, connect the device as detailed in the Installation section.

Selection

Select the device you wish to use from the **[Device selection]** list, because basic lung function tests (FVC, IVC, MVV) can be performed with any of the devices.

Calibration

Regular calibration ensures maximum accuracy.

Automatic warnings

Warning time interval can be entered for all device types. When this expires the device warns the user to perform calibration again. In this case it is recommended to perform the calibration.

Calibration is detailed in the Calibration section of the Spirometer specific volume.

Measurement

The individual measurement operating modes are detailed in the Measurement Modes section in each device specific volumes.

Enter comment

A separate comment may be entered for all measurement modes.

All comments are stored separately in the database.

- Click on the **[Diagnose]** icon to open Diagnose composer
- Select the measurement mode or the Patient mode from the list to which you would like to add a comment
- Enter the comment
- Press the **[Store]** button to store the comment in the database attached to the measurement

Previous comments

Previously created comments for the given measurement mode can be viewed anytime in the **[History]** list. The currently entered comment is not lost when viewing a previous or another measurement modes' comments. To display the comment select the measurement mode you would like to edit from the list.

Store

Pulmonary function test and Audiometry

To store curves marked **Visible** and **Questionable**, press the **[Store]** button.

Successful data storage returns a feedback.

Ergospirometry and Resting metabolic test

To store the actual measurement push the **[Store]** button.

Successful data storage returns a feedback.

Printing

Printable data is divided into several groups:

- Complex report on lung function tests: FCV, IVC, MVV, TGV, DLCO, PMAX
- Compliance
- Rhinomanometer
- Ergospirometry and Resting metabolic test
- Provocation test
- Audiometry

PRE/POST

The system can print two types of reports:

- Normal report: Three measurements' results simultaneously.
- PRE/POST report: Two measurements' results simultaneously and their difference in absolute and percentage format.

Customized reports

The printed report has the following parts:

- Header
- Parameter table
- Graphs
- Comment

The header is the only fixed part of the header, the other three may be turned on and off arbitrarily, only the desired parts make it into the report.

Highlight rows

To highlight odd rows at colour or greyscale printing select the **[Highlight odd rows]** checkbox

Simplified report

To print only the best measurements of all modes select the **[Just the best measure]** checkbox

Printing

Before printing measurement results have to be stored so the printed reports can be followed up.

- Click the **[Print...]** icon in the main menu
- Select the graphs, tables and manual diagnosis you would like to print
- Select the report language
- Select the report type: normal or PRE/POST
- Click on a button in the **[Print preview]** section to view the print preview
- After making the necessary settings, click the **[Print]** button

During printing graph display is similar to on-screen display:

- Complete curve or only the representative curve section
- One or more graphs

Export report into PDF and other graphical formats

This feature provides export of the printed report into the commonly used graphical formats. Exported reports can be stored and for example sent as an attachment to an e-mail.

Supported formats

- PDF, Adobe Acrobat document
- GIF picture
- JPEG picture
- BMP Windows Bitmap picture
- EMF and WMF vector graphics

Settings

The Export function is in the Report edition window

- Click on the [Print...] button at the main menu

Contents and the format of the exported report are fully identical to the printed version. More information can be found in the chapter Pulmonary function test and Audiometry

To store curves marked **Visible** and **Questionable**, press the [Store] button.

Successful data storage returns a feedback.

Ergospirometry and Resting metabolic test

To store the actual measurement push the [Store] button.

Successful data storage returns a feedback.

Printing.

Export procedure

Prior to printing and exporting results of measurements have to be stored in order to provide reliable traceability

- Click on the [Print...] button at the main menu
- Select graphs, tables and text fields to be exported
- Select the language of the report
- Select the type of the report PRE/POST
- For previewing the report click on any button at the [Print preview] section
- After setting click on the [Store] button

Exported graphs are fully identical to the graphs shown on the screen:

- Full curve or only the important part
- One or more graphs

Interface to information systems

Interface to frame systems

Our system provides communication according to the more commonly used protocols:

- Health Level Seven (HL7)
- Geräte Daten Träger (GDT)

These protocols provide exchange of the patient data and measured results between the lung diagnostics equipment and the frame systems. These protocols are predefined by the System administrator consequently can not be modified by the user.

Receiving the request for tests

Click on the [LINK...] button in the main menu and open the Import/Export window

According to your frame system type click one of the [HL7] or [GDT] buttons in the Import section in order to receive a Request for tests

If a Request for test is available the system automatically acquires it and lists all the requested tests

Exporting

Click on the [LINK...] button in the main menu and open the Import/Export window

According to your frame system type click one of the [HL7] or [GDT] buttons in the Export section. The system automatically exports the results of the tests.

Filling special forms

The system provides filling customer defined forms. Templates of the forms can be compiled in any ASCII format (HTML, XML, CSV etc.). Compilation of the form is the competence of the System administrator.

Filling a form

User may select a form from the preinstalled templates.

Click on the [LINK...] button in the main menu to open the Import/Export window

Click on the [HTML] button to open the Custom Report window

Select the desired template from the a [Templates] list

Click on the [Select] button and the form is automatically filled out

The PRE/POST measurement

The system supports measurement comparison – previous measurements may be compared against measurements made later:

- Select the patient
- Select and load the PRE (or previous) measurements
- Measure the current, POST values with the patient
- Select the two measurement to be compared
- Print the PRE/POST report

Retrieve measurement

Perform the steps described in the Patient Database section:

- Open the database
- Select the patient
- Select one or more measurements
- Load the measurement results

Max. 8 measurements may be displayed simultaneously, so if you loaded 6 measurements, you can perform 2 more measurements.

Notice

The program also makes it possible to print the PRE/POST report from the actually performed measurements.

Printing

Printing is similar to normal report printing:

- Select the PRE/POST option in the report edit window
- Select the parts of the report you would like to print

During printing the graphs are displayed similar to the screen:


- Complete curve or only the representative curve section
- One or more graphs

Possible problems

Software		
Problem	Diagnosis	Solution
Cannot find a patient.	Too many search criteria.	Enter less search criteria.
The patient's data cannot be loaded with the [Selected] button.	Not a single measurement has been selected.	If there is only one measurement in the list, use the [Load all] button.
When making a PRE/POST report, the program only prints previous data loaded from the database.	The new measurement has not been stored.	The measurement must be stored with the [Store] button before printing.
Starting up the program the connection to the database server is unsuccessful	After starting up the Windows the database server is not started yet	Wait a few minutes while all the processes of the operational system are fully running
It is impossible to access to the database connected via network	Network connection is interrupted The problem might be caused by the database server	Check the connection with the database server Consult with the system administrator or with the installer
It is impossible to access the local database even after 5 minutes as the Windows started up	The problem might be caused by the database server	Make a cold restart of the PC
It is permanently impossible to access the local database permanently even after re-starting the PC	The problem might be caused by the database server or by the malfunction of the operational system	Consult with the system administrator or with the installer

SR

SPIROMETRY



Supported devices

PDD-301/sh Spirometer
PDD-301/r Rhinomanometer
PDD-301/sco Breath CO monitor and spirometer
PDD-301/rco Breath CO monitor and rhinomanometer

Piston Ltd.
1033 Budapest, Szőlőkert u. 4/b

CE 1979

Rev.: SR-EN-04.00

Date modified: 04/07/2012

TABLE OF CONTENT

Introduction	4
Devices	4
Symbol annotation	5
Technical overview	5
Installation.....	6
PDD 301/s Spirometer installation	6
PDD 301/r Rhinomanometer installation	8
PDD 301/sco and PDD 301/rco Breath CO monitor installation.....	10
User interface	11
Settings	11
Icons.....	11
User interface general design	14
Calibration	16
Flow meter calibration	16
Checking calibration results	17
Viewing previous calibration data	17
Measurements.....	18
Daily routine – overview	18
System overview table	19
Available examinations	20
Entering environmental data	23
Zero setting.....	23
Preparations	24
Connection between the patient and the device	25
Measurement evaluation – Pulmonary function test.....	25
Compilation of PRE/POST reports	27
Measurement modes	29
Forced Vital Capacity	29
Inspiratory Vital Capacity	31
Maximal voluntary ventilation.....	33
Rhinomanometry.....	35
Breath carbon monoxide monitoring	37
Maintenance	39
Device maintenance	39
Flow meter maintenance.....	39
Breath CO monitor maintenance.....	40
Single-use parts	40
Reusable parts and accessories.....	40
Troubleshooting.....	41

Possible problems	41
Technical data.....	42
Warranty	42
Limited liability	42
Safety instructions	42
Informing values	43
Electrical data	43
Mechanical data	44
Guaranteed values	44
List of accessories	45
EMC Guidance and manufacturer's declaration	46

Devices

Piston Ltd.'s respiratory diagnostics product family contains the following members:

PDD-301/s Spirometer

Measurement operating modes

- Forced inspiration and expiration
- Static vital capacity
- Maximum voluntary ventilation

Design

- Flow meter: PPF-18 PinkFlow, symmetric Pitot tube
- USB computer connection
- Portable design

PDD-301/r Rhinomanometer and Spirometer

Measurement operating modes

- Forced inhalation and exhalation
- Static vital capacity
- Maximum voluntary ventilation
- Nasal respiratory resistance measurement with active anterior and posterior methods

Design

- Flow meter: PPF-18 PinkFlow, symmetric Pitot tube
- USB computer connection
- Portable design

PDD-301/sco and PDD-301/rco Breath carbon monoxide monitor

Measurement operating modes

- Forced inhalation and exhalation
- Static vital capacity
- Maximum voluntary ventilation
- Breath carbon monoxide concentration
- Version PDD-301/rco additionally provides nasal respiratory resistance measurement with active anterior and posterior methods

Design

- Flow meter: PPF-18 PinkFlow, symmetric Pitot tube
- Side stream gas analysis for higher accuracy and shorter response time

- USB computer connection
- Portable design

Symbol annotation

The following symbols indicate which descriptions apply to which device.



Spirometer



Rhinomanometer



Breath carbon monoxide monitor

Technical overview

Lung diagnostic device family main parts description:

Flow meter PDD-301/s, PDD-301/r PDD-301/sco, PDD-301/rco




PPF-18 PinkFlow, symmetric Pitot tube flow meter, which provides pressure difference in proportion with the flow speed.

A differentiate pressure sensor converts the pressure difference to electric signal.

PDD 301/s Spirometer installation

S



Connect the USB cable to the PC



For the connection push the PinkFlow flow meter into the quick connector
For removal push the metal button



To remove the PinkFlow flowmeter push the metal button at the back side of the quick connector and pull off the flowmeter

The PinkFlow flow meter can be used without bacterial filter as well. In this case a clean PinkFlow meter should be installed prior to each patient.



Connect one MPA-30 mouthpiece to bigger diameter end of the PinkFlow flow meter

If there is no possibility to provide a clean PinkFlow flow meter prior to each patient you have to use a bacterial and viral filter to avoid cross contamination.


PBF-100MG bacterial and viral filter can be used.



Connect one PBF-100MG bacterial and viral filter to bigger diameter end of the PinkFlow flow meter

PDD 301/r Rhinomanometer installation

R




Connect the USB cable to the device and to the PC



Connect the blue connector of the twin tubing to the blue coded socket indicated with the Flow meter label and the white connector to the white coded socket




Connect the green connector of the nasal pressure tubing to the green coded socket
For the proper connection the lock should be turned 180 degrees clock wise



The disc filter of the pressure port prevents the device from the contamination
When the disc filter gets dirty it has to be re-placed



For the connection push the PinkFlow flow meter into the quick connector
For removal push the metal button



Select a proper size Nasal probe and lace its tubing thru the PinkFlow flow meter. The plug of the Nasal probe has to face the patient side of the flow meter namely it has to be on the opposite side to the release button of the pneumatic quick connector.



Connect the tubing of the Nasal probe to the barbed fitting of the disc bacterial filter



Warning: When a tubing of the Nasal probe is laced thru the PinkFlow flow meter the sensitivity of the flow meter is modified. This modification is automatically corrected in the Rhinomanometer mode.

Do not use the Nasal probe during any other measurement!

Only the Piston made Nasal probes can be used with the device.



The PinkFlow flow meter can be used only without bacterial and viral filter in Rhinomanometer mode consequently a clean PinkFlow flow meter has to be connected prior to each measurement.



Lace the plug of the Nasal probe thru the adapter of the facial mask and connect the PinkFlow flow meter to the adapter.

PDD 301/sco and PDD 301/rco Breath CO monitor installation

C

Connect the USB cable to the device and to the PC




Connect the blue connector of the twin tubing to the blue coded socket indicated with the Flow meter label and the white connector to the white coded socket



Connect the yellow connector indicated with the Gas sample label to the yellow coded socket
For the proper connection the lock should be turned 180 degrees clock wise



For the connection push the PinkFlow flow meter into the quick connector
For removal push the metal button



The disc filter of the gas sampling tubing prevents the device from the contamination
When the disc filter gets dirty it has to be replaced

Settings

Part of the parameters are for system data that can seriously effect measurement accuracy.

You can view these parameters through the user interface, but they cannot be modified.

Only professionals can modify these data in the PistonXP.ini file.

Calibration Syringe

You can set the calibrating pump volume at the **Setup/Options/Devices/Calibration Syringe** menu

Number of calibrating cycles

The number of calibrating cycles with the calibration pump may be set in the **Setup/Options/Devices/Calibration Syringe** menu:

- Minimum: 2
- Maximum: 20
- Recommended: 10

Spirometer



Select the **Spirometer** group in the **Setup/Options/Devices** menu.

The system senses the spirometer connection to the USB port in 2 seconds.

Calibration time interval

You can set how often the device should remind you of the need for calibration

Rhinomanometer



Select the **Rhinomanometer** group in the **Setup/Options/Devices** menu.

The system senses the rhinomanometer connection to the USB port in 2 seconds.

Calibration time interval

You can set how often the device should remind you of the need for calibration.

Icons

Main window – Pulmonary function test



Open calibration measurement window. The checkmark indicates that all connected devices are calibrated.



Open calibration window. The exclamation point indicates that one or more connected device needs to be calibrated.



Forced vital capacity (FVC)



Inspiratory vital capacity (IVC)



Maximal voluntary ventilation (MVV)



Rhinomanometry



Breath carbon monoxide monitoring



Manual

Calibration



Start calibration



Skip specific phase during plethysmograph calibration, continue from the next step



Abort calibration



Store measured results



Print calibration report

Measurement windows – Pulmonary function test



Start measurement in at FVC, IVC, MVV, TGV, TLC measurement




Measure left nostril resistance








Measure right nostril resistance














Finish measurement (in case of successful measurement)

	Abort measurement (partial results are lost)
---	--

	The measurement is technically correct
	The measurement is most likely technically incorrect


	The icon in the summary table indicates the active curves, the specific measurement's curve is also displayed
	Indicates questionable curves The curves appear dashed
	The program does not store curves marked like this and they do not appear on the graph either

	Store active and questionable curves, measurement and their parameters
	Mark all curves as active
	Mark the three best curves as active, hide all other curves
	Effectively delete the selected curve
	List of Lung function Parameters
	Instructions
	PRE/POST
	Warnings
	Setup
	Animation
	Miller Quadrant

User interface general design



The following image shows the general design of the measurement screens. The individual measurement windows may differ from each other but the main controls are identical.



Device selector

Select the device to be used from the drop down list.

This is necessary if, for example, you own a Plethysmograph and a Spirometer, and would like to perform IVC measurement.

Zero setting

Runs manual Zero setting of the selected device

Without manual Zero setting the system automatically sets zero before all measurements

Menu

The program's general main menu, which contains the grouped basic functions.

Navigator

Controls that group the basic phases of daily routine.

Patient data

Contains the most important measurement data for the patient selected from the database.

BTPS data

These are the environment data measured by members of the PDT-111 device family. If you only own the PDD-301 device, this is where you can set the individual values manually.

Complex curves

The more complex curves of the individual measurement operating modes. For example, in case of FVC measurement, the flow-volume loops, in case of Plethysmograph measurement Resistance and TGV loops.

Graph settings

This is where you can set graph display modes. These settings are also available on the Options panel, details may be found in the **Hiba! A hivatkozási forrás nem található.** (page **Hiba! A könyvjelző nem létezik.**) section.

Control

This filed contains the basic control functions during the measurement. The appropriate Function buttons are shown in square brackets:

- Start measurement [F3]
- Start special measurement section [F4]
- Finish measurement after a successful measurement [F5]
- Stop measurement, abort measurement (for example, in case of malfunction) [ESC]
- Store, print

Spirogram

Volume – time graph, which monitors the patient's breathing during the measurement.

Information panel

This section contains information, settings and functions:

- Current measurement parameter list
- PRE/POST settings and parameters
- Measurement related warnings, error messages
- Measurement instructions

After turning the device on and entering the BTPS data, calibration is recommended for maximum measurement accuracy.

Calibration is recommended when starting a new shift, after flow sensor disinfections or replacement.



IMPORTANT

If work environment conditions (temperature, air pressure, humidity) change significantly, re-calibration is recommended.

Flow meter calibration



The flow meter volume calibration ensures maximum accuracy and is an efficient way to check the proper operation.

It is possible to perform measurements without calibration but at least 5% additional error must be taken into account.

Connecting the flow meter

Connect the patient side, the bigger diameter side of the PinkFlow flow meter of the Spirometer or the Rhinomanometer directly to the calibrating pump.

Calibration process

Spirometry/Calibration

In case of several connected devices, select the one to be calibrated from the **[Device selection]** list.

The calibration should be performed in two steps. At first the peak flow should be at about 1,0 l/s and afterwards at about 5,0 l/s

Press the **[Start]** button to start calibration.

Push the fully drawn out calibrating pump with uniform speed all the way in, then pull it out all the way.

Horizontal lines on the loop curve indicate optimal flow limits. During calibration make sure the calibration curve peaks are within these lines.

The number of calibration cycles may be set as described in the **Hiba! A hivatkozási forrás nem található.** section (page **Hiba! A könyvjelző nem létezik.**). The number of recommended cycles is 10.

The first part of calibration should be done with the peak flow at about 1,0 l/s (red curves)

The second part of calibration should be done with the peak flow at about 5,0 l/s (green curves)

After the calibration process the system automatically calculates calibration factors for the different flow values.

The following values appear in the calibration result table:

- Param – name of the measured parameter
- Pred – reference value
- Meas – the measured value during calibration
- % – difference of measured value from the reference

Possible error messages

Calibration must contain at least 10 exhalations and inhalations.
--

Calibration was not performed properly:

- There were less calibration cycles than prescribed
- The flow meter slipped out of the calibrating pump during calibration

Asymmetry error

In this case either calibration was performed incorrectly or an error occurred in the system:

- You did not pull out or push in the calibrating pump all the way
- Check pneumatic connections
- Check flow meter assembly
- Check that the twin tube is not broken or punctured
- Check that there is no liquid in the flow sensor or the twin tube

Flowmeter error out of allowed range

If during calibration the device measures the calibration volume with greater than 20% error, there is a chance for hardware problems.

Checking calibration results

It is recommended to store calibration results as the tendency over time can help draw conclusions concerning device stability and possible aging.

Click the **[Store]** button to store calibration results.


Click the **[Print]** button to print calibration results.

Viewing previous calibration data

Select the **[Result]** tab on the calibration window Information panel.

You can search for previous calibration results from the **[Reload calibration data]** time-sorted list.

Daily routine – overview



System overview table

	PDD-301/spf	PDD-301/rpf	PDD-301/sco	PDD-301/o	PDT-111/p	PDT-111/d	PRE-101
	Spirometer	Rhinomanometer	Breath CO meter	Oscillometer	Body Plethysmograph	CO-Diffusion	Ergospirometer
PinkFlow* flowmeter	+	+	+	+	+	+	+
Forced Vital Capacity	+	+	+	+	+	+	+
Static Vital Capacity	+	+	+	+	+	+	+
Maximal occlusion pressure	+	+	+	+	+	+	+
Rhinomanometry		+	optional	+	optional		
Breath CO			+				
Impulse oscillometry				+			
Thoracic gas volume					+		
Airways resistance					+		
Légzési munka					+		
Compliance					optional		
CO-Diffusion					optional	+	
Ergospirometry							+
EKG with 12 leads							+
Weight	220 g	220 g	420 g	2,5 kg	200 kg	5,5 kg	4,5 kg
Dimensions	150 * 82 * 45 mm	150 * 82 * 45 mm	190 * 138 * 68 mm	260 * 155 * 160 mm	1680 * 925 * 790 mm	320 * 200 * 240 mm	320 * 200 * 140 mm
Power source	USB port	USB port	USB port	90 – 260 VAC 50/60 Hz	90 – 260 VAC 50/60 Hz	90 – 260 VAC 50/60 Hz	90 – 260 VAC 50/60 Hz

Available examinations

Parallel measurements

The program makes it possible to perform eight different measurements in all measurement modes. All eight measurements' data can be stored and reloaded later.

Forced exhalation and inhalation



The most widely applied method for dynamic pulmonary function test. Detailed description may be found in the Measurement modes (page 29) section.

In this operating mode the device measures the following parameters:

FVC **Forced Vital Capacity**

Expired volume after full inspiration at the highest possible flow

FEV*0,5 **Forced Expiratory Volume 0,5 sec**

The amount of air exhaled in the first 0.5s during forced exhalation

FEV*1,0 **Forced Expiratory Volume 1,0 sec**

The amount of air exhaled in the first 1.0s during forced exhalation

FEV*0,5/IVC

The ratio of FEV*0.5 and the static vital capacity

FEV*0,5/FVC

The ratio of FEV*0.5 and the forced vital capacity

FEV*1,0/IVC

The ratio of FEV*1.0 and the static vital capacity

FEV*1,0/FVC

The ratio of FEV*1.0 and the forced vital capacity

PEF **Peak Expiratory Flow rate**

Highest flow during forced exhalation

FEF*25-75% **Forced mid-Expiratory Flow rate**

The average volume-flow speed calculated for the middle half of forced exhalation

MEF*75% **Forced Expiratory Flow at 75% lung volume**

Flow when 75% of the forced vital capacity is still in the lung

MEF*50% **Forced Expiratory Flow at 50% lung volume**

Flow when 50% of the forced vital capacity is still in the lung

MEF*25% **Forced Expiratory Flow at 25% lung volume**

Flow when 25% of the forced vital capacity is still in the lung

FET **Forced Expiratory Time**

The duration of forced exhalation.

MTT **Mean Transit Time**

The average leaving time from the lung of gas molecules during forced expiration

FIVC **Forced Inspiratory Vital Capacity**

Inspired volume after full expiration at the highest possible flow

FIV*0,5 **Forced Inspiratory Volume 0,5 sec**

The amount of air inhaled during the first 0.5 seconds of forced inhalation.

FIV*1,0 **Forced Inspiratory Volume 1,0 sec**

The amount of air inhaled during the first 1.0 seconds of forced inhalation.

PIF **Peak Inspiratory Flow rate**

Highest inhalation flow speed during forced inhalation

FIF*25-75% **Forced mid-Inspiratory Flow rate**

The average flow calculated for the middle half of the forced inhalation.

The following graphs are displayed during measurement:

- Volume/time curve
- Flow/volume loop

Static vital capacity

The most widely used method for the static lung function test.
Detailed description may be found in the Measurement modes (page 29) section

In this operating mode the device measures the following parameters:

IVC **Inspiratory Vital Capacity**

Total inspired volume after a full expiration

IRV **Inspiratory Reserve Volume**

The inspiration reserve volume is volume, what the patient can inhale from the average inhalation endpoints of quiet breathings

ERV **Expiratory Reserve Volume**

The expiration reserve volume is volume, what the patient can exhale from the average exhalation endpoints of quiet breathings

TV **Tidal Volume**

The average volume moved during quiet breathing

SVC **Slow Vital Capacity**

Total expired volume after full inspiration

The following graphs are displayed during measurement:

- Volume/time curve
- Flow/volume loop

Maximal voluntary ventilation



A rarely used dynamic lung function test.

Details may be found in the Measurement modes (page 29) section.

In this operating mode the device measures the following parameters:

MVV

Maximal Voluntary Ventilation

The maximum respiratory volume measured during voluntary respiration, calculated for one minute

MVV*f

Maximal Voluntary Ventilation Frequency

The maximum respiratory frequency measured during voluntary respiration, projected for one minute

The following graphs are displayed during measurement:

- Volume/time curve
- Flow/volume loop

Breath CO measurement



The device insures the measurement of breath carbon monoxide concentration. It is inevitable in the smoking cessation program.

Details may be found in the Measurement modes (page 29) section.

In this operating mode the device measures the following parameters:

CO ppm

Breath CO concentration

%COHb

Carboxyhemoglobin %

SVC

Slow Vital Capacity

The Slow Vital Capacity is the volume which was exhaled slowly by the patient after a total inspiration.

The following graphs are displayed during measurement:

- Volume/time curve

Rhinomanometry



Measuring the nasal airway resistance.

Detailed information may be found in the Measurement modes (page 29) section.

In this operating mode the device measures the following parameters:

Flow (50 Pa)

The flow speed at 50 Pa drive pressure.

Flow (75 Pa)

The flow speed at 75 Pa drive pressure.

Flow (100 Pa)

The flow speed at 100 Pa drive pressure.

Flow (150 Pa)

The flow speed at 150 Pa drive pressure.

Flow (300 Pa)

The flow speed at 300 Pa drive pressure.

The following graphs are displayed during measurement:

- Volume/time curve
- Flow/P_{Ch} Resistance loop

Entering environmental data



Entering exact environmental data is necessary for proper BTPS correction.

If the temperature, humidity or air pressure changes, the data must be re-entered.

BTPS

The top right parts of the individual measurement windows contain the BTPS data panel where you can enter the environmental data.

Automatic BTPS parameter measurement

The following Piston devices provide automatic measurement of ambient temperature, humidity and pressure:

- PDD-301/shm Integrated BTPS module
- PAM-201 Individual BTPS module with USB connection

These devices automatically measure the environmental data and display them in the BTPS panel.

Warning!

Entering incorrect environmental data may cause even 15% inaccuracy

Zero setting



For exact volume measurement the zero setting of flow meter channel must be performed immediately before the measurement.

Preparation

There cannot be any flow through the flow meter during zero setting, so the patient cannot take the connected mouthpiece into the mouth.

Notice

In case of Plethysmograph and the Diffusion capacity test pneumatic valves detach the flow meter from the pressure transducer, so zero setting occurs automatically in the background.

Patient may continue breathing thru the flow meter.

Zero setting process

The program automatically starts the zero setting process immediately before each measurement.

The system evaluates the data measured during the zero setting process, and displays an error message and repeats the zero setting process if a zero error is encountered.

Manual zero setting

You can reset the currently selected device anytime with the **[Zero]** button next to the **[Device selection]** list in the program header.

Notice

Zero setting is automatically performed before calibration.

Preparations ---

Device

Patient circuit

To avoid cross contamination connect a clean PinkFlow flow meter or a new disposable bacterial and viral filter before each patient measurement.

Patient

This chapter describes those general considerations which are necessary to inform the patient and get the patient prepared for the tests

Recommended body position

- Sitting on a chair
- Straight back
- Level head
- Tight clothing or jewels must not prevent free breathing

Directions

Respiratory examination requires patient cooperation so patient preparation and instructions are important for the measurement:

- Let the patient know the measurement process and goal
- Show the patient how to take in the mouthpiece, especially in case of the bite-grip mouthpiece used with Plethysmograph and Diffusion capacity test
- Prepare the patient for any unusual and uncomfortable events, such as shutter closing or breath holding.
- In case of the Diffusion capacity meter it is possible to practice without inhaling the gas mixture.
- The basic pulmonary function test including forced expiration and vital capacity tests may take about 15 minutes
- A complex pulmonary function test including may take about 30 minutes

Connection between the patient and the device

The proper usage of the patient circuit is necessary for the accurate measurement

Nose clip

In order to avoid any nasal breathing during tests apply nose clip. Even at a good cooperation level there can be leakage thru nostrils without nose clip.

Of course at the measurements with facial mask usage of the nose clip is senseless.

Mouth piece

The leakage free connection between the patient and the device is a key issue for the accurate measurements.

We provide the following options:

- MPA-30 Anatomically shaped mouthpiece, for the basic pulmonary function tests
- PMP-30 Mouthpiece with bite-on grip guarantees the perfect sealing for the most demanding tests like whole body Plethysmograph and diffusion capacity test
- PBF-100M-G Bacterial and viral filter with elliptic patient side

Measurement evaluation – Pulmonary function test

You can simultaneously perform and display max. 8 measurements.

Measurements deemed not appropriate can be deleted and a new one can be performed.

The system selects the best measurements based on different aspects for each measurement operating mode.

Sort order aspects

Forced Vital Capacity

Decreasing order based on **FVC+FEV*1.0** sum

Larger values are better

Inspiratory Vital Capacity

Decreasing order based on **IVC** value

Larger values are better

Hyperventilation

Decreasing order based on **MVV** value

Larger values are better

Rhinomanometry

Increasing order based on **RES [75]** value

Smaller values are better

Pairing

R

In case of the Rhinomanometer the two nostrils' resistance measurement may differ in time, so before storing them in the database the two sides' measurements must be paired.

The Rhinomanometer's measurement screen has a separate summary table for right and left nostril measurements.

The program automatically pairs the right and left nostril measurement pairs based on quality or measurement time.

Normal mode

In case of a simple measurement it is recommended to sort the measurements based on airway resistance.

PRE/POST

In case of comparison measurement it is recommended to sort the measurements based on measurement time, so the first right side measurement is paired with the first left side measurement.

Measurement selection

Measurement management requires the selection of individual measurements:

- Click on a single point on a curve with the mouse button
- or
- Click on the line in the summary table belonging to the curve.

The selected curve appears on the graph with a dotted line, the summary table's appropriate line is light blue.

Measurement selection for storage

3	✓	📈
2	✓	📈
6	✓	📈
4	✓	?
8	✓	✗
7	✓	✗
1	⚠	✗
5	✓	✗

All measurements are displayed in the measurement summary table

Colour of serial number is identical to the colour of the curve

There are two icons next to their number:

The first icon indicates the measurement's technical quality:



The measurement is technically correct.



The measurement is technically incorrect.

The second icon indicates the given measurement's status:

Visible curve



The measurement appears on the graph with a solid line.
The system can store and print the measurement.

Questionable curve



The measurement appears on the graph with a dashed line.
The system can store and print the measurement.
However, they are easily distinguishable on the graph for the user.

Turned off curve



The system does not store or print this measurement.
Unlike when deleted, the curves can be displayed again at any-time.

Changing status:

- Select the measurement and the required curve will be dashed.
- Clicking on the selected line again will rotate the curve status.

Selecting the best measurement

Press the **[Best]** button to have the program automatically display the three best measurements and turn off the rest.

Press the **[All]** button to display all performed measurements.

Click the **[Report]** button to select just one from a couple of measurements. Only the selected curve will be visible all the others will be switched off.

Delete measurement

It is possible to delete measurements not already stored:

- Select the measurement and the required curve will be dashed.
- Press the **[Delete]** button.

Compilation of PRE/POST reports

Measurement selection

Quick keys

Parameter list

<div> <div>Warnings</div> <div>Setup</div> <div>Animation</div> <div>Miller Quadrat</div> <div>Parameters</div> <div>Instructions</div> <div>PRE/POST</div> </div>				
PRE	4. 2007. 03. 05. 13:27:48			
POST	2. 2007. 03. 05. 13:33:35			
PARAM	REF	PRE	POST	DIF
FVC	5.24	5.68	6.17	-0.49
IVC	5.48	5.62	5.85	-0.23
FEV*0.5		3.04	3.21	-0.16
FEV*1.0	4.41	4.23	4.49	-0.26

PRE/POST report compilation:

- Select the PRE/POST tab on the measurement window information panel
- Select the two measurements to compare from the measurement selection list
- You can also use the [Quick keys] to select the measurement, use the mouse to select the desired measurement – either from the summary table or directly on the graph
- Selecting the two curves automatically refreshes the parameter table.

Forced Vital Capacity



Measurement goal

The goal of the measurement is to get parameters of the volume forced expiration and inspiration



Measurement process


Instruct the patient to perform the following manoeuvres:

- Put on the nasal clip so he / she can only breath through the nose.
- Take at least three quiet breathings
- Take as deep as possible inspiration
- Take as fast and deep as possible expiration
- Take as fast and deep as possible inspiration

The patient has 60 second to perform the FVC manoeuvre

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct FVC manoeuvre

Phases: quiet breathing, deep inspiration, forced exhalation, forced inhalation, return to normal breathing.

Miller Quadrant


The Miller Quadrant an effective graphical tool which helps making the lung function diagnoses.

The vertical axle shows the ratio of FVC measured value and the reference value

The horizontal axle shows the ratio of FEV*1,0 measured value and the reference value

The diagram is divided into four quadrants:

- Normal
- Restrictive
- Obstructive
- Combined



Animation

Animation with blowing away dandelions helps with motivating children to reach their maximum effort during FVC manoeuvre

Inspiratory Vital Capacity



Measurement goal

The goal of the measurement is to get the parameters of the maximal inspiration.



Measurement process


Instruct the patient to perform the following manoeuvres:

- Put on the nasal clip so he / she can only breath through the nose
- At least three quiet breaths
- As deep expiration as possible
- As deep inspiration as possible
- Optionally as deep expiration as possible to get the SVC (Slow Vital Capacity) parameter

The patient has 60 second to perform the IVC manoeuvre

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct IVC manoeuvre


Phases: Quiet breathing, complete deep expiration, complete deep inspiration, return to normal breathing.

Maximal voluntary ventilation



Measurement goal

The goal of the measurement is to get the amount of volume the patient can move in a given time.



Measurement process


Instruct the patient to perform the following manoeuvres:

- Put on the nose clip so he / she can only breath through the nose.
- Move as much air as possible.
- The patient can choose the breath frequency as convenient
- To prevent hypocapnia it is not recommended to continue this measurement for more than 15 seconds.

The patient has 60 seconds to perform MVV manoeuvre.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct MVV measurement


The patient breaths quickly, evenly during the measurement.

Rhinomanometry



Measurement goal

The goal of the measurement is to get the patient's nasal airway resistance.



Measurement process


The following series of manoeuvres must be performed:

- The patient must clean the nasal canals
- Place the appropriate size nasal plug into the side opposite the measured one
 So if you would like to measure the resistance of the right nasal canal, place the nasal plug into the left nostril, and vice versa
- Have the patient hold the appropriate size facial mask against his / her face to prevent leaking
- The patient can only breathe through the free nostril, not through the mouth

The patient has 60 seconds to perform the manoeuvres, but usually a few even respiratory cycles are enough.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct Rhinomanometry measurement

The patient breathes relaxed during the measurement.

Selecting the loop curve

The system records several respiratory cycles during the measurement and by default displays their average.

However, you have the option to view the curves one-by-one:

- Select the measurement you would like to modify from the measurement summary table
- Click with the mouse on the spirogram – the marker jumps to the next cycle
- The system automatically recalculates the parameters

Breath carbon monoxide monitoring




Measurement goal

The goal of the measurement to follow up the smoking habits of the patient and assist the smoking cessation program.



Attention

Avoid measuring patients whose exhalation may contain alcohol because the CO sensor of the device may get wrong



Measurement process


The following series of manoeuvres must be performed:

- Place the Nasal clip to avoid the breathing thru the nostrils
- Perform at least 3 quiet breathings
- Breath hold et least for 10 seconds
- Slow and even full expiration

The patient has 60 second to perform the manoeuvres.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct Breath CO measurement

Phases: Quiet breathing, complete deep inspiration, breath holding, complete slow expiration

Evaluation of the results

The following evaluations are given by the system after the measurement:

- Non Smoker
- High value for non-smoker
- Smoker
- Frequent smoker
- Addicted smoker
- Heavily addicted smoker
- Dangerously addicted smoker

Device maintenance

Our lung diagnostics devices do not require special maintenance. For continuous reliable operation take care of the following:

- To prevent device contamination and patient cross-contamination, use a new disposable bacterial and viral filter for all patient measurements
- The flow sensor must be contamination free
- The filter elements must be replaced according to instructions
- The tubes must always be dry and cannot be broken

Flow meter maintenance

The flow meter condition and cleanliness affects measurement accuracy.

Cleaning measurement head main parts

The individual patient circuit type installations are described in section Installation.

The plastic parts may be disinfected with cold water and appropriate chemicals (for example, Sekusept), and may be used after rinsing and drying.

Cleaning the flow meter


- Disconnect the PinkFlow flow meter from the docker
- Clean the flow meter in a cold disinfecting solvent
- After it is completely dried, reassemble the flow meter

Cleaning the pneumatic twin-tubes


- Disconnect the twin-tube from the device and the flow meter
- Rinse the tube
- After it is completely dried, reconnect the tube

Facial mask maintenance

The facial mask's pneumatic cushion may deflate with time. For appropriate fitting the escaped air must be replaced



Fill a LUER cone-shaped syringe with air



Fit the syringe into the facial mask valve opening, push it in all the way to open the valve. Push the air in



Repeat the previous two steps until the facial mask is properly inflated. Never over inflate as it will not properly fit the face

Breath CO monitor maintenance

The expected life time of the CO sensor is 2 years.

For replacement of CO sensor please contact the local responsible of the Manufacturer!

Single-use parts



It is strictly prohibited to clean and/or reuse the single-use parts

Bacterial and viral filter PBF-100-G and PBF-100M-G

The used bacterial and viral filters are considered to be dangerous waste materials please handle accordingly.

Reusable parts and accessories

The following plastic parts can be cleaned in a cold disinfecting water-based solution (for example: Glutaraldehyde, Sekusept, Cidex e.t.c.)

Type	Description	Material
PPF-18	PinkFlow flow meter	Polystyrol
MPA-30	Mouthpiece	Polypropylene
PMP-30	Mouthpiece with bite-on grip	Thermoplastic elastomer

Some discoloration may occur to the plastic parts after frequent disinfection.

Possible problems

Spirometry		
Problem	Diagnosis	Solution
During quiet breathing the Volume(time) curve drifts up or down	<p>After several quiet breaths have the patient remove the mouthpiece</p> <p>The program continues to display the curve.</p> <p>Watch the spirogram for at least 10 seconds</p>	<p>Set Zero again and repeat the measurement.</p> <p>Check that liquid did not get into the flow sensor or the twin-tube leading to it.</p>
The measured values deviate from expected to a considerable extent	<p>The device has to be recalibrated</p> <p>Environmental data must be checked</p>	<p>If the situation does not get better even after recalibration, clean the pneumatic twin tubes and check the flow sensor according to Flow meter maintenance (page 39.) section</p>
Rhinomanometer		
The resistance curves are too steep	The pressure meter's or the nasal plug's pneumatic tube is not connected appropriately, it maybe punctured	<p>Check the pressure tubes</p> <p>The nasal plug, the filter, or the pressure release tube is clogged</p>
The resistance curves are too flat	The device measures the drive pressure to be too high	<p>Check the pressure tubes</p> <p>The nasal plug, the filter, or the pressure release tube is clogged</p>

Warranty

The device complies with the effective Technical Specifications.

The manufacturer guarantees the product according to the terms of the Installation/Delivery protocol.

The warranty does not cover post-delivery careless shipping, unprofessional storage, violent damaging, abnormal operation, unprofessional operation, inefficient protection against external effects, natural disasters, or not following the contents of the User Manual.



**Check package condition after delivery.
If packaging is damaged, notify the carrier and Piston Ltd., or its representative.**

Use of any broken or otherwise damaged products (devices, accessories etc) is dangerous and forbidden!

Limited liability

Piston Ltd. and its carriers, according to the valid laws, do not accept any responsibility for any individual, unforeseeable, direct or indirect damages (including loss of business profit, interruption of business activity, loss of business data, or any other damages due to financial loss), resulting from the use or non-usefulness of the product.

Safety instructions

To avoid possible damages and accidents, please pay attention to the following safety instructions:

- Make sure the mains voltage is the same as that on the product label
- Make sure the connection cable is not damaged
- Take care of your device according to the maintenance section
- Only use the device according to the manual
- Do not use any accessories not recommended for the device
- Store the device in a dry place
- Keep the cable away from heat source, sharp objects, rough surfaces and check the cable's good condition
- Do not expose the device to direct sunlight or strong light (more than 1500 lux)
- Do not use the device in a highly dusty environment
- Do not use the device in a highly vibrating environment
- Take care to ensure the current environmental conditions

The equipment complies with the applicable requirements of laws and standards.

Shipping conditions

Air temperature:-30 °C ÷ +60 °C
 Relative humidity: 10% ÷ 100%
 Atmospheric pressure: 500 ÷ 1060 mbar

Storage conditions

Air temperature: 0 °C ÷ +50 °C
 Relative humidity: 10% ÷ 85%
 Atmospheric pressure: 500 ÷ 1060 mbar

Operating conditions

Air temperature:+10 °C ÷ +40 °C
 Relative humidity: 30% ÷ 75%
 Atmospheric pressure: 700 ÷ 1060 mbar

Informing values

Expected lifetime

Devices 8 years
 Measurement head lifetime 2 years

Forced inhalation and exhalation

Measurement duration 60 s
 Volume measurement limit 15 l

Vital capacity measurement

Measurement duration 60 s
 Volume measurement limit 15 l

Maximal voluntary ventilation

Measurement duration 60 s
 Volume measurement limit 250 l/min

Sampling frequency

PDD-301 device family 250 Hz

Other data

Analog-digital converter resolution 16 bit

Electrical data

The connected computer's and printer's electrical data is found in the respective manufacturer provided specifications.

The following values apply only to the Piston Ltd. manufactured devices:

PDD-301/sh – Spirometer

PC connectionUSB 1.1
Power..... Does not require external power

PDD-301/r –Rhinomanometer

PC connectionUSB 1.1
Power..... Does not require external power

PDD-301/sco and PDD-301/rco – Breath CO monitor

PC connectionUSB 1.1
Power..... Does not require external power

Mechanical data

PDD-301/sh – Spirometer

Flow meter PPF-18 PinkFlow
Dimensions..... H 150 * W 82 * D 45 mm
Weight255 g

PDD-301/r – Rhinomanometer

Flow meter PPF-18 PinkFlow
Dimensions..... H 150 * W 82 * D 45 mm
Weight210 g

PDD-301/sco and PDD-301/rco – Breath CO monitor

Flow meter PPF-18 PinkFlow
Dimensions.....H 185 * W 140 * M 60 mm
Weight650 g

Guaranteed values

PPF-18 – PinkFlow flow meter

Type PPF-18 PinkFlow
Principle of operation Symmetric Pitot tube
Flow range ± 18 l/s
Dead space 36 ml
Resistance..... 60 Pa/l/s @ 15 l/s
Weight34 g

PDD-301 – Spirometer and Rhinomanometer

Flow meter PPF-18 PinkFlow
Flow measurement range ± 18 l/s
Flow measurement accuracy..... $\pm 2\%$ or ± 10 ml/s

Flow measurement resolution 10 ml/s
 Linearity $\pm 2\%$
 Volume measurement range ± 20 l
 Volume measurement accuracy $\pm 2\%$ or ± 10 ml
 Volume measurement resolution 10 ml

PDD-301/r – Rhinomanometer

Flow meter PPF-18 PinkFlow
 Pressure measurement range ± 2 kPa
 Pressure measurement accuracy $\pm 3\%$ or ± 15 Pa
 Resistance measurement accuracy $\pm 3\%$ or ± 30 Pa/l/s

PDD-301/sco – Breath CO monitor and spirometer

Accuracy of CO measurement $\pm 2\%$ or ± 2 ppm
 All other measurement parameters are the same as of the PDD-301/sh spirometer

PDD-301/rco – Breath CO monitor and Rhinomanometer

Accuracy of CO measurement $\pm 2\%$ or ± 2 ppm
 All other measurement parameters are the same as of the PDD-301/r rhinomanometer

List of accessories

Included accessories

The current Shipping contract contains the list of accessories included in the purchase price.

Optionally purchased accessories

The following information must be provided when ordering accessories and disposables:

- Description
- Type
- Part number
- Device type and serial number for which the accessories are used

EMC GUIDANCE AND MANUFACTURER'S DECLARATION

Guidance and manufacturer's declaration – electromagnetic emissions


The PDD 301/r Rhinomanometer (particular implementation PDD 301/s Spirometer) is intended for use in the electromagnetic environment specified below. The customer or the user of the PDD 301/r Rhinomanometer should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PDD 301/r Rhinomanometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PDD 301/r Rhinomanometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

<p>The PDD 301/r Rhinomanometer (particular implementation PDD 301/s Spirometer) is intended for use in the electromagnetic environment specified below. The customer or the user of PDD 301/r Rhinomanometer should assure that it is used in such an environment.</p>			
IMMUNITY test	IEC 60601 test level	Compliance Level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<p><5 % U_T (>95 % dip in U_T) for 0,5 cycle</p> <p>40 % U_T (60 % dip in U_T) for 5 cycles</p> <p>70 % U_T (30 % dip in U_T) for 25 cycles</p> <p><5 % U_T (>95 % dip in U_T) for 5 s</p>	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the PDD 301/r Rhinomanometer requires continued operation during power mains interruptions, it is recommended that the PDD 301/r Rhinomanometer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The **PDD 301/r** Rhinomanometer (particular implementation **PDD 301/s** Spirometer) is intended for use in the electromagnetic environment specified below. The customer or the user of the **PDD 301/r** Rhinomanometer should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 V_{rms} 0,15-80 MHz</p> <p>3 V/m 80 MHz – 2,5GHz</p>	<p>3 V_{rms} 0,15-80 MHz</p> <p>3 V/m 80MHz – 2,5GHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the PDD 301/r Rhinomanometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1,17\sqrt{P}$ $d = 1,17\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,33\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PDD 301/r Rhinomanometer is used exceeds the applicable RF compliance level above, the PDD 301/r Rhinomanometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the PDD 301/r Rhinomanometer.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Recommended separation distances between portable and mobile RF communications equipment and the PDD 301/r Rhinomanometer

The **PDD 301/r** Rhinomanometer (particular implementation **PDD 301/s** Spirometer) is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of **PDD 301/r** Rhinomanometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **PDD 301/r** Rhinomanometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz – 80 MHz $d = 1,17\sqrt{P}$	80 MHz – 800 MHz $d = 1,17\sqrt{P}$	800 MHz – 2,5 GHz $d = 2,33\sqrt{P}$
0,01	0,12	0,12	0,24
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,7	3,7	7,38
100	11,7	11,7	23,33


For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PD

Body Plethysmograph CO-Diffusion



Supported devices

PDT-111/p	Body Plethysmograph
PDT-111/d	CO Diffusion
PDT-111/pd	Body Plethysmograph with CO-Diffusion

Piston Ltd.
1033 Budapest, Szőlőkert u. 4/b

CE 1979

Version: PD-EN-04.10

Revision date: 2013.07.25.

TABLE OF CONTENT

Introduction	4
Devices	4
Symbol annotation	5
Technical overview	5
Installation.....	7
PDT-111/pd Body-Plethysmograph and CO-Diffusion installation.....	7
PDT-111/d és PDT-111/pd Patient circuit assemble	8
PDT-111/p and PDT-111/pd Cabin placement	10
Connecting the gas cylinder to the PDT-111/pd	11
User Interface.....	12
Icons.....	12
Settings	13
Calibration	18
Flow meter calibration	18
Plethysmograph calibration	18
Diffusion capacity test calibration	19
Measurements.....	20
Available examinations	20
Zero setting.....	25
Preparations	25
Connection between the patient and the device	25
Measurement evaluation – Pulmonary function test.....	26
Measurement modes	27
Thoracic Gas Volume	27
Diffusion capacity test	30
Single-Breath Diffusion capacity test	31
Intra-Breath Diffusion capacity test.....	36
Compliance.....	40
Maximum inspiratory and expiratory pressure	47
Maintenance	49
Device maintenance	49
Flow meter maintenance.....	49
Diffusion capacity test	49
Single-use parts	50
Reusable parts and accessories.....	50
Refilling of the test gas	51
Troubleshooting.....	52
Escape from the Plethysmograph cabin.....	52
Possible problems	52

Technical data.....	54
Warranty	54
Limited liability	54
Safety instructions	54
Informing values	55
Electrical data	55
Mechanical data	56
Guaranteed values	57
List of accessories	57
EMC Guidance and manufacturer's declaration	58

Devices

Piston Ltd.'s respiratory diagnostics product family PDT-111 contains the following members:

PDT-111/p Whole-body plethysmograph

Measurement operating modes

- Forced inhalation and exhalation
- Static vital capacity
- Maximum voluntary ventilation
- Thoracic gas volume measurement
- Residual volume measurement
- Respiratory resistance measurement
- Respiration work measurement
- Nasal respiratory resistance measurement with active anterior and posterior methods

Design

- Heated flow meter with stainless steel screen
- Hermetically sealed cabin
- USB computer connection

PDT-111/d Diffusion capacity measurement

Measurement operating modes

- Forced inhalation and exhalation
- Static vital capacity
- Maximum voluntary ventilation
- Residual volume measurement
- Transfer factor measurement

Design

- Heated flow meter with stainless steel screen
- Measurement gases: CO and CH₄
- USB computer connection

PDT-111/pd Whole-body Plethysmograph and diffusion capacity test

Measurement operating modes

- The combined device's measurement operating modes are identical to the PDT-111/p and PDT-111/d devices' measurement operating modes

Symbol annotation

The following symbols indicate which descriptions apply to which device.



Plethysmograph



Diffusion capacity test

Technical overview

Lung diagnostic device family main parts description:

Flow meter PDT-111/p, PDT-111/d és PDT-111/pd



PPF-18 PinkFlow, symmetric Pitot tube flow meter, which provides pressure difference in proportion with the flow speed.

A differentiate pressure sensor converts the pressure difference to electric signal.

Patient circuit PDT-111/p



The patient circuit ensures device-to-patient connection and contains the following parts:

- Lilly-type heated screen flow sensor
- Shutter magnetic valve
- Metronome

Patient circuit PDT-111/d and PDT-111/pd



The patient circuit ensures device-to-patient connection and contains the following parts:

- Lilly-type heated screen flow sensor
- Shutter magnetic valve
- Metronome
- Demand valve
- Gas sampling line

Plethysmograph cabin PDT-111/p and PDT-111/pd



The closed cabin makes it possible to measure alveolar pressure by non-invasive method through two transfers.

- The cabin has two types of leakage time constants, user selectable depending on measurement method
- The cabin door can only be locked from the out, from the inside it can only be opened
- In case of locking mechanism malfunction, the cabin can be opened from the outside by undoing a few screws

Gas supplying system PDT-111/d and PDT-111/pd



To determine diffusion capacity the patient must inhale an air mixture containing CO 0.3% and CH₄ 0.3% gases.

Parts of the gas supplying system:

- High-pressure gas cylinder
- Pressure reductor
- Main valve
- Pressure limiting safety valve
- Demand valve, which ensures gas amount required by the patient's inhalation

Gas analyzer PDT-111/d and PDT-111/pd



The patient's exhaled has mixture must be sampled and analyzed to determine diffusion capacity.

Parts of the gas analyzer:

- Gas sampling capillary
- PermaPure capillary moisture exchanger to normalize the humidity of the gas samples
- Sampling pump, vibrating membrane design
- NDIR (Non Dispersive Infra Red) multi-channel quick gas analyzer

Environment status measurement module

PDT-111/p, PDT-111/d and PDT-111/pd



BTPS correction requires the measurement of the following environmental data:

- Environment temperature
- Environment relative humidity
- Atmospheric pressure

Power supply PDT-111/p, PDT-111/d and PDT-111/pd



Medical design switching power supply, which enables operating the device from any mains voltage:

- Mains voltage: 90~264 VAC
- Mains frequency: 50~60 Hz

PDT-111/pd Body-Plethysmograph and CO-Diffusion installation




Figure: PDT-111/d CO-Diffusion




Figure: PDT-111/pd Plethysmograph equipped with CO-Diffusion

Part number	Description	Comment
PPC30/P	Patient circuit holder	Only for standalone Plethysmograph cabin
PPC30/D	Patient circuit holder with demand valve	Only for standalone diffusion capacity test
PPC30/PD	Patient circuit holder with demand valve including demand valve blocker	For plethysmograph cabin with integrated diffusion capacity test

PCV-30	Check valve	
PDY-30	Y-tube	For diffusion capacity test
PPF-18	PinkFlow flow sensor	
VM-MP-2000	Gumi szájsutóra	
SV-36	Shutter valve plate	
PPF-18/C	PinkFlow pneumatic docker	
PPF-18/TT	PinkFlow twin tube	Blue and white connectors
BPF220	Gas sampling capillary fast release connector	
VM-3435	Gas sampling capillary disk filter	
PPBF-900	Gas sampling moisture exchanger capillary	
PPT-300	Extension tube of the gas sapling capillary	Yellow connector


Connection

Connect the power supply cables and pneumatic connectors of the patient circuit according to the types and colour codes of the connectors.

Rhinomanometer and Compliance measurement

Prior to the rhinomanometer and Compliance measurement remove the PinkFlow flow sensor from the Y-tube and make the measurement this way!

PDT-111/d és PDT-111/pd Patient circuit assemble



Item	Part number	Description	Comment
1	PDY-30	Y-tube	
2	PCV-30	Check valve	

3	SV-36	Valve plate	
4	PPF-18	PinkFlow flow sensor	
5	PPF-18A	PinkFlow flow sensor docker	
6/A		Patient circuit holder with demand valve	For standalone diffusion capacity test
6/B		Patient circuit holder with demand valve including demand valve blocker	For plethysmograph cabin with integrated diffusion capacity test



In the case of standalone diffusion capacity test fasten the vice of the patient circuit holder to the table


The photo shows the patient circuit holder with demand valve



Push the valve plate (SV-36) to the rod of the pneumatic valve



Insert the check valve to the Y-tube in a way that the black rubber valve faces inside



Connect the Y-tube to the demand valve push it firmly till it finds its final position



Connect a clean PinkFlow flow sensor to the Y-tube



Connect the pneumatic dock to the PinkFlow flow sensor



Connect the gas sampling capillary to the pneumatic dock



Connect a clean bite-on type mouth piece to the PinkFlow flowsensor

PDT-111/p and PDT-111/pd Cabin placement



Place the plethysmograph warily. Since extremely small pressure differences have to be measured in the plethysmograph cabin, external interference adversely may affect measurement accuracy.

The following have to be taken into the consideration during installation and operation:

- Do not install the device in a draughty place
- There must not be a radiating heat source (radiator, sunny window) within 1 meter
- Air-conditioning unit's orifice within 1 meter
- Open window or door during measurement

The temperature in the laboratory has to be equalized refrain from any rapid change of the temperature!

Connecting the gas cylinder to the PDT-111/pd

D

Fix the gas cylinder into position



If necessary, replace the gas cylinder sealing ring. Connect the pressure meter reductor to the gas cylinder



Connect the high-pressure gas tube to the device's pneumatic quick connector



After opening the main valve of the gas cylinder set the secondary pressure to 6 bar

**WARNING!**

Always turn off the gas cylinder at the end of every shift and in all situations when diffusion capacity measurements are not made for longer periods

Icons

Main window – Pulmonary function test



Thoracic gas volume and Resistance measurement



Compliance measurement



Diffusion capacity test



Maximum inspiratory and expiratory pressure

Measurement windows – Pulmonary function test



Start gas mixture inhalation



Close shutter for respiratoric pressure measurement



Start Resistance loop recording (if the measurement already started and the patient breaths quietly)



Start Dynamic Compliance loop recording (if the measurement already started and the patient breaths quietly)



Prepare Static Compliance loop recording (if the measurement already started and the patient breaths quietly)



Start Static Compliance recording

Settings

Parts of the parameters are for system data that can seriously effect measurement accuracy. You can view these parameters through the user interface, but they cannot be modified. Only professionals can modify these data in the PistonXP.ini file.

Measuring the environmental status

If you have any member of the PDT-111 device family, the device automatically measures the environmental data necessary for BTPS correction:

- air pressure
- temperature
- humidity

You can select in the **Setup/Options/Operation/Source of environment** info menu which device's environmental data the system should use:

- USB diffusion
- USB Plethysmograph
- Automatic

In case of automatic option

- The device calculates with the environment data measured in the cabin during Resistance and TGV measurement.
- During all other lung function test the system uses the environment data measured by the diffusion capacity meter outside the cabin.

Plethysmograph

P

Select the **Plethysmograph** group in the **Setup/Options/Devices** menu.

The following device operation related settings are available:

Enable Metronome

You can turn the light and sound signal of metronome on or off.

Low breathing rate

The plethysmograph cabin has two selectable leakage time constants.

Low breath frequency means less stress for the patient, but thermal compensation takes longer.

High breathing rate

High breath frequency means more stress for the patient, but thermal compensation happens sooner.

AutoBTPS

If enabled, resistance loops are automatically closed

Efficiency of BTPS

If AutoBTPS option is off, this is where you can set the theoretical lung model BTPS correction effectiveness.

Default value: 45%

Balance time

This is where you can set the waiting time after closing the cabin's door to the beginning of the very first TGV measurement in order to reach the adequate thermal balance in the cabin

Number of Resistance loops

This is where you can set the number of Resistance loops the device should record in one measurement.

Number of TGV loops


This is where you can set the number of breaths the shutter should close.

Barometer calibration

It is possible to exactly set the barometer measuring environmental pressure. Enter the exact environmental pressure value and click the **[OK]** button.


Resistance Calculation

This is where you can set the algorithm used to calculate Resistance loop slope.




Method by Matthys

Resistance loop is intersected with ± 0.5 l/s flow value and place the steepness indicating line on the geometric bisecting points of the horizontal intersections



Method of Peak pressure

We place the steepness indicating line on the peak value points of the pressure measured in the cabin

**Method of Maximum flow**

We place the steepness indicating line on the maximum flow points

Calibration time interval

You can set how often the device should remind you of the need for calibration.

Diffusion capacity meter



Select the **Diffusion** group in the **Setup/Options/Devices** menu.

The following device operation related settings are available:

CH₄ concentration

Original methane concentration of the test gas mixture

CO concentration

Original carbon monoxide concentration of the test gas mixture

Wash out volume (WOV)


The volume the device releases from the start of expiration until the beginning of gas analysis.

Breath hold time (BHt)

The period of breath holding, or closing the shutter.

BHt calculation method

Since diffusion in the lung starts from the beginning of inhalation and lasts till the end of expiration, the system offers several algorithms to calculate the effective diffusion time:



Ogilvie method

Period start: Start inhalation (1)

Period end: Start gas sampling (4)

Jones and Meade method

Period start: Start inhalation 1/3 (2)

Period end: Middle of gas sampling period (5)

Epidemiologic Standardisation Project method

Period start: Middle of inhalation period (3)

Period end: Start gas sampling (4)

Barometer calibration

It is possible to exactly set the environment pressure measuring barometer. Enter the exact environmental pressure value and press the **[OK]** button.

Calibration before all measurement

You can select to have automatic Zero setting and automatic calibration of gas analyzer before each test

Recommended option.

Using the shutter

You can select to have the shutter closed during breath holding.

If the shutter is closed, the patient has to make sure to hold the breath and not to press the shutter, because this can increase alveolar pressure and change the level of diffusion.

If the shutter is open, the patient has to make sure not to inhale nor exhale when holding the breath, because this changes the gas concentration in the lung and effects the measurement accuracy.

Limit values during breath holding

If you select the **Shutter close** option, you can enter the limit value for mouth pressure fluctuation, beyond which the device emits a warning.

If you select the **Shutter open** option, the inhaled volume limit value during breath holding cannot exceed ± 200 ml.

Calibration interval

You can set how often the device should remind you of the necessity to calibrate.

After turning the device on and entering the BTPS data, calibration is recommended for maximum measurement accuracy.

Calibration is recommended when starting a new shift, after flow sensor disinfections or replacement.



IMPORTANT

If work environment conditions (temperature, air pressure, humidity) change significantly, re-calibration is recommended.

Flow meter calibration

Connect the patient side of the flow meter of Plethysmograph or Diffusion capacity test directly or with the included corrugated tube to the calibrating pump.

Proceed the calibration process according to the Spirometry volume

Plethysmograph calibration

P

Diffusion capacity meter calibration is performed in the

Spirometry/Calibration menu

In case of several connected devices, select **Plethysmograph** from the **[Device selection]** list.

Plethysmograph calibration has to be started with the flow meter volume calibration.

The rest of calibration has to be performed with close cabin door:

Please close the cabin door!

After closing the cabin door, you have to wait for complete thermal compensation, which is 60s. The automatic calibration process requires the following steps:

Check cabin leakage test

During leakage test both cabin time constants are checked. In this case the built-in pump inflates the cabin to 15Pa pressure, and discharge like curves indicating cabin leaking appear on the screen.

If the time constants deviate from the prescribed values, the following error message appears on the screen:

Cabin leakage out of range

In this case there is something wrong with the cabin tightness:

- Check that the cabin door is closed properly.
- Check that there is no foreign object between the cabin door and the sealing.

- Check that the cabin door sealing is intact.

Cabin pressure meter calibration

In this case the built-in sine pump creates a 15Pa pressure difference in the cabin. The change in cabin pressure is visible on the screen.

If the cabin pressure meter calibration constant falls outside the 0.5 and 2.0 range, the following error message appears:

Cabin pressure out of range

In this case there is something wrong with the cabin tightness:

- Check that the cabin door is closed properly.
- Check that there is no foreign object between the cabin door and the sealing.
- Check that the cabin door sealing is intact.

Diffusion capacity test calibration



Diffusion capacity meter calibration is performed in the **Spirometry/Calibration** menu

In case of several connected devices, select **Diffusion** from the **[Device selection]** list.

Before calibrating the diffusion capacity meter, you have to ensure measurement gas supply:

- Open the gas cylinder main valve
- Check that the secondary pressure is set to 6 bar
- Disconnect the gas sampling capillary from the pneumatic docker of the flow sensor and reconnect to the calibration outlet of the device

Flow sensor volume calibration must be performed as detailed in the chapter Flow sensor Calibration with the following additions:

The calibration of the gas analyzer is performed automatically during the zeroing phase

Available examinations

The PDT-111 device family offers the following basic pulmonary function tests:

- Forced in- and expiration
- Static vital capacity
- Hyperventilation
- Rhinomanometry (optional)

Guide to these measurement modes can be found in the Spirometry volume of this user manual.

Parallel measurements

The program makes it possible to perform eight different measurements in all measurement modes. All eight measurements' data can be stored and reloaded later.

Plethysmograph measurement

P

Complex lung function test to measure the mechanical parameters of the respiratory system.

Details may be found in the Measurement modes (page 27) section.

In this operating mode the device measures the following parameters:

TLC **Total Lung Capacity**

Total lung capacity

TGV **Thoracic Gas Volume**

Thoracic gas volume

RV **Residual Volume**

Residual volume

RV/TLC **Residual Volume/Total Lung Capacity**

The ratio of residual volume and total lung capacity

Raw **Resistance of Airways**

Airway resistance

Rin **Resistance of Airways at Inspiration**

Airway resistance during inhalation

Rex **Resistance of Airways at Expiration**

Airway resistance during exhalation

Req **Equivalent Resistance**

Equivalent airway resistance

sRaw **Specific Resistance of Airways**

Specific airway resistance

sRin **Specific Resistance of Airways at Inspiration**
Specific airway resistance during inhalation

sRex **Specific Resistance of Airways at Expiration**
Specific airway resistance during exhalation

Gaw **Conductance of Airways**
Airway conductance

Gin **Conductance of Airways at Inspiration**
Airway conductance during inhalation

Gex **Conductance of Airways at Expiration**
Airway conductance during exhalation

sGaw **Specific Conductance of Airways**
Specific airway conductance

sGin **Specific Conductance of Airways at Inspiration**
Specific airway conductance during inhalation

sGex **Specific Conductance of Airways at Expiration**
Specific airway conductance during exhalation.

BF **Frequency of Breathing at Resistance Measurement**
Breathing frequency during resistance measurement.

W **Work of Breathing at Resistance Measurement**
Breathing work during resistance measurement.

The following graphs are displayed during measurement:

- Volume/time curve
- Flow/ P_{alv} Resistance loop
- P_{mouth}/P_{cab} TGV loop

Diffusion capacity test



Complex lung function test for the measurement of oxygen bounding capacity of the lung.

Detailed information may be found in the Measurement modes (page 27) section.

In this operating mode the device measures the following parameters:

TLC **Total Lung Capacity**
Total lung capacity

IVC **Inspiratory Vital Capacity**
Total inspired volume after a full expiration

IRV **Inspiratory Reserve Volume**
The inspiration reserve volume is volume, what the patient can inhale from the average inhalation endpoints of quiet breathings

ERV	Expiratory Reserve Volume
The expiration reserve volume is volume, what the patient can exhale from the average exhalation endpoints of quiet breathings	
TV	Tidal Volume
The average volume moved during quiet breathing	
RV	Residual Volume
Residual volume	
RV/TLC	Residual Volume/Total Lung Capacity
The ratio of residual volume and total lung capacity	
FRC	Functional Residual Capacity
Function residual capacity	
FRC/TLC	
The ration of Function residual capacity and total lung capacity	
Single-Breath:	
Tlco	Transfer factor of the lung for CO
The lung's transfer factor for carbon monoxide	
The SI unit of measurement: mmol/min/Pa	
Dlco	Transfer factor of the lung for CO
The lung's transfer factor for carbon monoxide.	
The imperial unit of measurement: ml/min/mmHg	
Klco	Transfer coefficient of the lung for CO
Transfer coefficient.	
BHt	Breath hold time
Effective time breath hold time	
Intra-Breath:	
Tlco IB	Transfer factor of the lung for CO
The lung's transfer factor for carbon monoxide	
The SI unit of measurement: mmol/min/Pa	
Dlco IB	Transfer factor of the lung for CO
The lung's transfer factor for carbon monoxide	
The imperial unit of measurement: ml/min/mmHg	
Klco IB	Transfer coefficient of the lung for CO
Transfer coefficient	
Auxiliary parameters supporting evaluation of measurements	
WOV	Wash-Out Volume
Wash-out volume	
FACH4	Expiratory concentration CH4
Exhaled CH4 concentration	

FACO	Expiratory concentration CO
Exhaled CO concentration	
CCCH4	Calibration constant of CH4 channel
CH4 channel calibration constant	
CCCO	Calibration constant of CO channel
CO channel calibration constant	
GSL	Gas Sample Lag
Gas sampling time displacement	
CH4 IB L	CH4 concentration by IB on Left side
Exhaled CH4 concentration at the left side of gas sampling window during Intra-Breath measurement	
CH4 IB R	CH4 concentration by IB on Right side
Exhaled CH4 concentration at the right side of gas sampling window during Intra-Breath measurement	
CO IB L	CO concentration by IB on Left side
Exhaled CO concentration at the left side of gas sampling window during Intra-Breath measurement	
CO IB R	CO concentration by IB on Right side
Exhaled CO concentration at the right side of gas sampling window during Intra-Breath measurement	

The following graphs are displayed during measurement:

- Volume / time curve
with CH₄, CO, CO₂ gas concentration functions
- Flow/volume loop

Compliance



The measurement of compliance of the lungs

Detailed information may be found in the Measurement modes (27. page) section.

In this operating mode the device measures the following parameters:

Cdyn	Dynamic Compliance
Dynamic Compliance of lungs	
Edyn	Dynamic Elastance
Reciprocal ratio of the Dynamic Compliance	
Cstat	Static Compliance
Quasi static compliance of lungs	
Estat	Static Elastance
Reciprocal ratio of quasi Static Compliance	
Wcomp	Work of Breathing at Cdyn measurement
Viscous work of breathing during Dynamic Compliance measurement	

ReqComp**Equivalent Resistance**

Equivalent resistance

Parameters listed below are calculated only when the value of the TLC was measured before of the Compliance measurement at TGV or Diffusion capacity test mode:

Cdyn/TLC**Dynamic Compliance/Total Lung Capacity**

Ratio of Dynamic Compliance and TLC

Cstat/TLC**Static Compliance /Total Lung Capacity**

Ratio of Static Compliance and TLC

Cdyn/FRC**Dynamic Compliance/ Functional Residual Capacity**

Ratio of Dynamic Compliance and FRC

Cstat/FRC**Static Compliance/ Functional Residual Capacity**

Ratio of Static Compliance and FRC

Cdyn/TGV**Dynamic Compliance/Thoracic Gas Volume**

Ratio of Dynamic Compliance and TGV

Cstat/TGV**Static Compliance/Thoracic Gas Volume**

Ratio of Static Compliance and TGV

Maximal respiratory pressure

The measurement of the respiration muscles' maximum strength.

Detailed information may be found in the Measurement modes (page 27) section.

In this operating mode the device measures the following parameters:

PE_{max} Maximal expiratory pressure

Maximum expiratory pressure.

PI_{max} Maximal inspiratory pressure

Maximum inspiratory pressure.

The following graphs are displayed during measurement:

- Volume / time curve

Automatic BTPS parameter measurement

Piston PDT-111 family devices are measuring BTPS parameters automatically.

**Attention!**

The Plethysmograph measures the environmental conditions inside the cabin. Data measured in a long closed, overheated cabin is not suitable for measurements with a separate device.

Zero setting

For exact volume measurement the zero setting of flow meter channel must be performed immediately before the measurement.

In case of Plethysmograph and the Diffusion capacity test pneumatic valves detach the flow meter from the pressure transducer, so zero setting occurs automatically in the background.

Patient may continue breathing thru the flow meter.

Zero setting process

The program automatically starts the zero setting process immediately before each measurement.

The system evaluates the data measured during the zero setting process, and displays an error message and repeats the zero setting process if a zero error is encountered.

Manual zero setting

You can reset the currently selected device anytime with the **[Zero]** button next to the **[Device selection]** list in the program header.

Notice

Zero setting is automatically performed before calibration.

Preparations

Device

Patient circuit

To prevent cross contamination a new disposable bacterial and viral filter must be connected before each patient measurement.

The height of patient circuit of the Plethysmograph and Diffusion capacity test is freely adjustable according to patient comfort.

Patient

This chapter describes those subjects which are inevitable in the information and preparation of the patient.

Directions

Inform the patient about the goal and process of the pulmonary function tests. The detailed description can be found in the Spirometry volume.

Connection between the patient and the device

The proper usage of the patient circuit is necessary for the accurate measurement

Nose clip

In order to avoid any nasal breathing during tests apply nose clip. Even at a good cooperation level there can be leakage thru nostrils without nose clip.

Of course at the measurements with facial mask usage of the nose clip is senseless.

Mouth piece

The leakage free connection between the patient and the device is a key issue for the accurate measurements.

We provide the following options:

- MPA-30 Anatomically shaped mouthpiece, for the basic pulmonary function tests
- PMP-30 Mouthpiece with bite-on grip guarantees the perfect sealing for the most demanding tests like whole body Plethysmograph and diffusion capacity test
- PBF-100M-G Bacterial and viral filter with elliptic patient side

Esophageal Balloon

The esophageal balloon is used during the compliance test. Insertion of the esophageal balloon can be monitored by the system. It is fully described in the Compliance test chapter

Measurement evaluation – Pulmonary function test

The system selects the best measurements based on different aspects for each measurement operating mode.

Sort order aspects

TGV and Resistance

Decreasing order based on **TLC** value: Larger values are better

or

Increasing order based on **Req** value: Smaller values are better

Use the switches above the summary table to change the sort order

Diffusion capacity test

Decreasing order based on **TLC** value

Larger values are better

Thoracic Gas Volume

P

Measurement goal

The goal of the measurement is to get the mechanical parameters of the patient's respiratory system:

- Thoracic gas volume
- Total lung capacity
- Airway resistance and its components



Premises

To get the IVC parameter value required to determine Thoracic Gas Volume measurement parameters:

- Perform an IVC measurement before the TGV measurement. The advantage to this is that the TGV manoeuvre becomes simpler.
- Retrieve an IVCh measurement from the database that looks recent.
- The IVC measurement can be performed even during the TGV measurement, after shutter opening. The advantage to this is that

you get all the parameters in one measurement; the disadvantage is that the measurement is more complicated

The system always selects the most recent IVC measurement result.

Preparation

TGV is a complex measurement, requiring considerable cooperation from the patient:

- Seat the patient in the cabin
- Set the chair height
- Set the height of the patient circuit
- Place the nasal clip so the patient can only breath through the mouth
- Explain the exact measurement process to the patient
- Prepare the patient for unexpected and unusual events, such as shutter closing
- Close the cabin door
- Set the speaker volume


Measurement process

The following series of maneuvers must be performed:

- The metronome turns on after the measurement starts.
- The countdown to thermal compensation begins.
- The patient can practice breathing to the metronome during the waiting period, the current respiratory frequency is continuously displayed in the top part of the measurement window.
- After the waiting period tell the patient to evenly breath according to the metronome. The current respiratory frequency is continuously displayed in the top part of the measurement window, in the **[Breathing frequency]** field. In case of respiratory frequency that is extremely different from the settings, the field changes to red.
- When the patient is breathing evenly according to the metronome, press the **[Start]** button to start Resistance, then TGV loop recording.
- During shutter closing have the patient try to perform even expiration and inspiration manoeuvres. The patient does not have to exert too much force, but has to suck on it and has to push against the shutter.
- After shutter opening let the patient rest, then optionally instruct him / her to perform a complete IVC manoeuvre – a complete expiration followed by a complete inspiration.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct TGV measurement

Phases: quiet breathing, quiet breathing during Resistance loop recording, shutter closing during TGV loop recording, complete IVC manoeuvre, return to normal breathing.

Selecting resistance loop

Several (5 by default) Resistance loops are recorded during the measurement. From these the system automatically selects the loop closest to the average.

If you would like to select the Resistance loop into the report manually, follow these steps:

- Select the measurement to be modified in the measurement summary table.
- Click with the mouse on the spirogram – the marker jumps to the next cycle.
- The system automatically recalculates the parameters.

Diffusion capacity test

The system provides two alternative methods for measuring diffusion capacity:

- Single-Breath – Breath holding method

The Breath holding method is the standard way for diffusion capacity measurement. The patient has to make a deep inspiration from the gas mixture and has to hold the breath for a certain period of time and afterwards has to exhale slowly and evenly

- Intra-Breath – Without breath holding method

An alternative method for diffusion capacity measurement for the poorly cooperating patients who are unable for the breath holding manoeuvre

It is enough for the patient to make a deep inspiration from the gas mixture and afterwards the patient may start the slow and even expiration immediately



After changing the gas cylinder enter the actual gas concentrations according to Certificate issued by the Filling station as described in the Installation (page 7) section

Single-Breath Diffusion capacity test



The Single-Breath method is the standard way for diffusion capacity test

Measurement goal


The primary goal of the measurement is to get the oxygen binding capacity of the patient's lung:

- Transfer factor calculated for carbon monoxide
- Functional residual capacity
- Total lung capacity
- Residual lung capacity



IMPORTANT!

There must be at least 5 minutes between measurements so the measurement gas mixture can completely clear out from the patient's lung.



Preparation

Diffusion capacity measurement is a complex process, requiring considerable cooperation from the patient:

- Seat the patient.
- Set chair height.
- Set measurement head height.
- Place the nasal clip so the patient can only breath through the mouth.
- Tell the patient the exact measurement process.
- Prepare the patient for unexpected and unusual events, such as gas mixture inhalation through the Demand valve, and holding the breath.

Measurement process

The following series of manoeuvres must be performed:

- The patient must breath evenly after the measurement started.
- The [Gas] button is enabled after the third relaxed breath.
- Instruct the patient to exhale relaxed, deep, completely.
- Press the [Gas] button when the patient started the complete deep exhalation.
- During the next inhalation the patient inhales the measurement gas mixture and must be instructed to inhale completely, deeply.
- Breath holding begins after inhalation, the remaining time appears on the screen.
- Warn the patient to hold the breath, do not press it onto the shutter and do not try to inhale.
- The metronome indicates the end of breath holding with a visible and audible signal.
- When the shutter opens, instruct the patient to exhale relaxed; the patient cannot inhale until gas sampling is performed.
- The metronome indicates the end of gas sampling with a visible and audible signal.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.

Practice operating mode

Turning on the **[Practicing mode]** option lets the patient practice the complete measurement process without inhaling the measurement gas mixture.

Use the practice operating mode for a few measurements with a badly cooperating patient.

Using the shutter

You can select whether the shutter is closed or opened during breath holding:

If the shutter is closed, the patient must make sure that he / she holds the breath and does not press it onto the shutter, because this increases the alveolar pressure and change the level of diffusion.

If the shutter is open, the patient must make sure not to inhale or exhale during breath holding, because this changes the gas concentration in the patient's lung and effects measurement accuracy.

Measurement instructions

During measurement the system continuously indicates the upcoming manoeuvre.

It is important to note that the system does not indicate the immediately performed task, but the next in line.

For example "Complete, deep inhalation" appears towards the end of the complete deep exhalation, however it is obvious that inhalation comes after the completely finished exhalation and not immediately.

Possible error messages during closed-shutter breath holding

Mouth pressure out of allowed range during breath hold time. Repeat the measurement if possible!

Due to the mouth pressure variation the measured Tlco, Dlco and Klco values can be used only for orientation.

Valsalva- or Miller effect.

The patient must take care to hold his / her breath and not press it onto the shutter and not try to inhale.

This can increase or decrease the alveolar pressure, and change the diffusion amount.

Possible error messages during open-shutter breath holding

Expiration during inspiration of gas mixture. Repeat the measurement if possible!

If the shutter is open, the patient must take care not to inhale or exhale during breath holding.

This can change the gas concentration in the lung and effect measurement accuracy.

The patient did not breath in, the measured RV value is correct, but Tlco, Dlco and Klco values are invalid.

Since methane quickly mixes into the lung, the breath holding time has less effect on the RV and TLC measurement result. However, carbon monoxide diffusion greatly depends on breath holding time, so transfer factor calculation is not possible.

Patient made breathing during the breath hold time, the measured RV, TLC, Tlco, Dlco and Klco values are invalid.

As a result of inspiration the gas mixture in the lung is diluted to an unknown degree, so measurement evaluation is not possible, and the measurement must be repeated!

Possible error messages during and after gas sampling**Inspiration during expiring gas mixture**

The patient breathed back during gas sampling, so the results are inaccurate, the measurement must be repeated.

Too short expiration time


The patient did not exhale long enough for reliable gas sampling.

If the gas concentration curves and the measured values appear correct, you can keep the measurement.

Measurement refinement

The system automatically optimizes the gas sampling position, but manual refining may be needed:


- Select the measurement to be refined in the measurement summary table
- The gas sampling marker appears on the spirogram
- Use the arrows on the gas concentration panel to move the marker left and right.
- Use the methane (CH₄) curve (blue graph) to find the optimal marker position, where the curve is almost horizontal, but is close enough to the curve inflexion point
- When manually setting the marker you have to be careful, because it significantly affects the measurement results.




Correct Diffusion capacity measurement

Phases: quiet breathing, complete deep expiration, even gas mixture inspiration, breath holding, even expiration during gas sampling

Incorrect marker positions



The marker is **too much to the left**, at the beginning of the expiration curve



The marker is **too much to the right**, at the end of the expiration curve

Intra-Breath Diffusion capacity test



An alternative method for diffusion capacity measurement for the poorly cooperating patients who are unable for the breath holding manoeuvre

Measurement goal


The primary goal of the measurement is to get the oxygen binding capacity of the patient's lung:

- Transfer factor calculated for carbon monoxide
- Functional residual capacity
- Total lung capacity
- Residual lung capacity



IMPORTANT!

There must be at least 5 minutes between measurements so the measurement gas mixture can completely clear out from the patient's lung.



Preparation

Diffusion capacity measurement is a complex process, requiring considerable cooperation from the patient:

- Seat the patient.
- Set chair height.
- Set measurement head height.
- Place the nasal clip so the patient can only breathe through the mouth.
- Tell the patient the exact measurement process.
- Prepare the patient for unexpected and unusual events, such as gas mixture inhalation through the Demand valve, and holding the breath.

Measurement process

The following series of manoeuvres must be performed:

- The patient must breathe evenly after the measurement started.
- The [Gas] button is enabled after the third relaxed breath.
- Instruct the patient to exhale relaxed, deep, completely.
- Press the [Gas] button when the patient started the complete deep exhalation.
- During the next inhalation the patient inhales the measurement gas mixture and must be instructed to inhale completely, deeply.
- When the shutter opens instruct the patient to exhale slowly and evenly at approximately 0,5 litre/sec flow rate
- The metronome indicates the end of gas sampling with a visible and audible signal.

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.

Practice operating mode

Turning on the **[Practicing mode]** option lets the patient practice the complete measurement process without inhaling the measurement gas mixture.

Use the practice operating mode for a few measurements with a badly cooperating patient.

Using the shutter

If the patient circuit is equipped with the optional exhalation orifice and you are going to use it select the Shutter option.

If the patient circuit has no exhalation orifice or you are not going to use it unselect the Shutter option.

Measurement instructions

During measurement the system continuously indicates the upcoming manoeuvre.

It is important to note that the system does not indicate the immediately performed task, but the next in line.

For example “Complete, deep inhalation” appears towards the end of the complete deep exhalation, however it is obvious that inhalation comes after the completely finished exhalation and not immediately.

Possible error messages during and after gas sampling

Too short expiration time


The patient did not exhale long enough for reliable gas sampling.

If the gas concentration curves and the measured values appear correct, you can keep the measurement.

Measurement refinement

The system automatically optimizes the gas sampling position, but manual refining may be needed:


- Select the measurement to be refined in the measurement summary table
- The gas sampling marker appears on the spirogram
- Use the arrows on the gas concentration panel to move the marker left and right.
- Use the carbon monoxide (CO) curve (purple graph) to find the optimal marker position where part of the curve inside of the markers almost straight but it is close enough to the curve inflexion point
- When manually setting the marker you have to be careful, because it significantly affects the measurement results.




Correct Diffusion capacity measurement

Phases: quiet breathing, complete deep expiration, even gas mixture inspiration, slow and even expiration during gas sampling

Incorrect marker positions



The marker is **too much to the left**, at the beginning of the expiration curve



The marker is **too much to the right**, and the CO curve is not straight yet

Compliance



Measurement goal

The primary goal of the measurement is to get the compliance of the lungs or with other words to determine the change of the lungs volume as a function of the pressure against the wall of the lungs.

There is a possibility to measure the Dynamic and Static Compliance sequentially during one manoeuvre.

To get the value of the Transpulmonary pressure we have to measure the pressure of the intrapleural cavity with the help of oesophageal balloon:

$P_{TP} = P_{AV} - P_{PL}$ where P_{TP} - Transpulmonary pressure

P_{AV} - Alveolar pressure

P_{PL} - Pleural pressure (oesophageal pressure)



Oesophageal balloon


We propose the following oesophageal balloon for the Compliance measurement:

Manufacturer: Cooper Surgical
Trumbull, CT06611, USA
Tel: +1-203-601-5200

Description: Esophageal Balloon Catheter Set

Type: 47-9005

Of course any other type of oesophageal balloon can be used but the current description is relevant only for the specified type.



Parts of the Esophageal Balloon Catheter Set:

1. Esophageal Balloon Catheter
2. Stylet Y connector
3. Extension tube
4. Syringe, Glass 5 cc
(The set does not contain it)
5. Three way Stopcock
6. Pressure transducer tubing

Preparation

1. Have available for use a nemesias basin, tissues, a protective drape, topical anaesthetic, water –soluble lubricant and a 1 to 5 cc glass syringe, a glass of water and a straw
2. Select a naris with the best airflow for catheter insertion
3. If necessary apply a suitable topical anaesthetic (e.g. 2-4% Lidocaine Spray) to the patient's nasal passage and throat
4. Remove the sterile radiopaque catheter with the stylet from its protective sleeve from the catheter and discard
5. Remove the yellow protective sleeve from the catheter (1) and discard
6. Apply water-soluble lubricant to the distal tip of the catheter (1)
7. With the patient's head in a neutral position or flexed slightly forward, slowly insert the catheter (1) through the naris and hypo pharynx using a gentle advancing motion. The insertion may be easier if the patient is sipping water thru the straw.
8. Avoid placement of the catheter (1) in the trachea. Tracheal placement can be identified by patient choking or airway obstruction causing an increase in airway resistance and pressure.


9. To estimate depth in which to place catheter calculate the product of the patient's height x 0.288.
For example the patient is 175 cm height the principal depth is 50 cm.
10. Advance the catheter (1) to the calculated depth mark. (If the catheter meets obstruction, DO NOT FORCE THE CATHETER. Remove it and insert it through the other naris). At this depth, the balloon will be entering the thoracic cavity.
11. Attach the extension tubing (3) to the „Y” connector of the stylet (2), a syringe (4) and an isolated physiologic transducer to the 3-way stopcock (5)

Catheter Placement

12. Turn the 3-way stopcock (5) open to the syringe (4) and extension tube (3). Evacuate all the air from the balloon (1) by pulling back on the syringe plunger (4) and then allowing the plunger to return to a nonvacuum position. Use of a glass syringe avoids creating a vacuum in the balloon catheter.
13. Turn the 3-way stopcock (5) off to the extension tube (3), remove the syringe (4) and fill the syringe (4) with 1 cc of air.
14. Attach the syringe (4) to the 3-way stopcock (5) , open the 3-way stopcock (5) from the syringe (4) to the extension tube (3). Introduce 1 cc of air into the balloon (1). The balloon (1) will now be semi inflated. An incorrect amount of air in the balloon will adversely affect pressure wave performance (see Trouble Shooting Guide below).
15. After 1 ml of air is injected, turn the stopcock (5) off to the syringe (4), and open from the extension tube (3) to the transducer to read pressure from the catheter.
If no or a damped pressure signal is seen, the catheter (1) may need to be advanced further into the thoracic cavity or may be kinked on itself and needs to be withdrawn.
In the absence of diaphragmatic paralysis, the pressure recorded should be negative on inspiration. A positive inspiratory pressure recording may indicate gastric placement and the catheter (1) should be pulled back.
Pressures taken through the „Y” connector (2) are for balloon placement only. The stylet, luer cap and „Y” connector (2) are bonded into a single unit and must be removed from the catheter before taking pressure measurements for clinical purposes.
16. Once the catheter (1) has been positioned properly, disconnect the extension tube (3) from the stylet „Y” connector (2) and remove the stylet assembly from the catheter. Excessive curvature of the catheter may cause the stylet to bind in the catheter making removal difficult. If the stylet binds in the catheter during removal, instruct the patient to raise their head to straighten the catheter.

Pressure Data Acquisition

17. After removing the stylet assembly (2), reattach the extension tube (3) to the luer of the catheter (1) and repeat steps 12 through 15.
18. Radiographic placement of the balloon (1) is recommended to verify proper balloon placement.
19. When the catheter (1) is properly positioned it can be secured with tape to prevent extubation or movement.
20. Take pressure measurement.
21. Upon completion of the pressure measurements, deflate the catheter (1) prior to removal.



Premises

For determination some parameters of the Compliance measurement there is a need to have the TLC and/or TGV and/or FRC measured previously. The system provides a couple opportunities for that:

- Prior to Compliance measurement make a TGV measurement
- Prior to Compliance measurement make a diffusion capacity test
- Retrieve from the database a TGV or diffusion capacity test result as an Actual measurement

If there are more TLC values the system always selects the largest one.

Preparation

Compliance measurement is a complex process, requiring considerable cooperation from the patient:

- Inform patient on the process of the measurement
- Prepare patient for the unexpected and unusual events like closing of the shutter during Static Compliance measurement
- Pushing the [Monitor] button it is possible to monitor continuously the breathing of the patient. The spirogram and the pressure of oesophageal balloon are displayed simultaneously.
- As it is described before insert the oesophageal balloon into the proper depth.

The pressure of the oesophageal balloon displayed over the spirogram with a light blue line.

The pressure of the oesophageal balloon is displayed during the preparation and the measurement furthermore during browsing when a given measurement is selected.


Measurement process

The following series of manoeuvres must be performed:

- Starting the [Monitor] mode the metronome starts as well
- Patient may practice breathing according to the metronome in [Monitor] mode. The actual breathing frequency is displayed in the middle of the upper part of the screen.
- When the oesophageal balloon is at the proper position the oesophageal pressure is in counter phase with the spirogram
- When the patient is breathing quietly and evenly push the [Cdyn] button to start the measurement
- After three quiet breathing cycles the Dynamic Compliance loops are recorded
- Pushing the [Cstat] button the Dynamic Compliance mode is finished and the a [Shutter] button becomes active
- Instruct the patient for the full deep inspiration and afterwards for a quiet, even and very slow exhalation (200-500 ml/sec)
- Push the [Shutter] button as the patient starts the full deep inspiration
- The Static Compliance measurement starts automatically as the patient starts expiration.
During the expiration the Shutter closes automatically after certain exhaled volumes and interrupts expiration. When the Shutter is closed the Mouth pressure is recorded.
- The Static Compliance manoeuvre can be finished reaching the FRC level or it is finished automatically when the expiration flow less than 10 ml/s
- Bad dynamic or static maneuvers can be repeated within the 1 minute period but parameters will always be calculated from the last dynamic or static maneuver
- The system provides possibility for recording 8 consecutive Dynamic and Static Compliance measurement pairs

Pushing the [Done] button the measurement can be finished.

Pushing the [Discard] button the measurement can be deleted.



Correct Compliance measurement

Phases: Quiet breathing, Quiet breathing during recording of Dynamic Compliance loops, Preparation for Static Compliance with deep inspiration, Static Compliance measurement during slow and even expiration, Returning to the relaxed breathing

Selection of Dynamic Compliance loop

During the Dynamic Compliance measurement more loops are recorded.

The system selects automatically the loop which is closest to the average of the all loops.

However there is a possibility to select the desired loop manually as well:

- Select the measurement to be refined in the measurement summary table
- Click on the spirogram and marker jumps to the next cycle
- All parameters will be recalculated automatically

Manual correction

The system calculates automatically the steepness of the Dynamic and Static Compliance loops.

However there is a possibility to modify the steepness manually as well:

- Select the measurement to be corrected
- The steepness of Dynamic and Static Compliance loops can be corrected separately in the [Setup] panel at the a [Correction] filed
- All parameters will be recalculated automatically

Selection of TLC or TLC% display mode

There are two possibilities for interpretation of Compliance loops which can be selected at [Setup] panel in [Scale] section:

- Litre – The gradation of vertical axle is in absolute volume
- TLC % – The gradation of vertical axle is in percentage of TLC value

Warning! If there is no available TLC or RV value it is impossible to display the Compliance loops in TLC% mode
--


Maximum inspiratory and expiratory pressure



Measurement goal

The goal of the measurement is to measure the respiratory muscle strength:

- Maximum inspiration pressure
- Maximum expiration pressure




Measurement process

The following series of manoeuvres must be performed:

- Place the nasal clip so the patient can only breathe through the mouth
- Have the patient breathe evenly, relaxed
- Instruct the patient to take a complete, deep breath
- When the patient began the deep inhalation, press the [Shutter] button
- Have the patient try to exert as great force as possible exhaling against the shutter; the current mouth pressure value appears on the screen
- Pressing the [Shutter] button again opens the shutter immediately
- Have the patient breathe evenly, relaxed again
- Instruct the patient to exhale completely, deep
- When the patient began the deep exhalation, press the [Shutter] button
- Have the patient try to exert as great force as possible inhaling against the shutter; the current mouth pressure value appears on the screen
- Pressing the [Shutter] button again opens the shutter immediately

Push the **[Done]** button to stop the measurement.

Push the **[Discard]** button to delete the measurement.



Correct PImax / PEmax measurement

Phases: quiet breathing, complete deep inspiration, PEmax measurement
quiet breathing, complete deep expiration, PImax measurement, quiet breathing

Device maintenance

Our lung diagnostics devices do not require special maintenance.

For continuous reliable operation take care of the following:

- To prevent device contamination and patient cross-contamination, use a new disposable bacterial and viral filter for all patient measurements
- The flow sensor must be contamination free
- The filter elements must be replaced according to instructions
- The PermaPure moisture exchange capillary must be replaced according to instructions
- The tubes must always be dry and cannot be broken

Flow meter maintenance

The flow meter condition and cleanliness affects measurement accuracy.

Cleaning measurement head main parts

The individual patient circuit type installations are described in section Installation.

The plastic parts may be disinfected with cold water and appropriate chemicals (for example, Sekusept), and may be used after rinsing and drying.

Cleaning the pneumatic twin-tubes

- Disconnect the twin-tube from the device and the flow meter
- Rinse the tube
- After it is completely dried, reconnect the tube

Diffusion capacity test

It is advisable to validate and check the linearity of the non-dispersive infrared gas analyzer of the diffusion capacity test at least once a year.

☞ Request your local dealer for this service!

Single-use parts



It is strictly prohibited to clean and/or reuse the single-use parts

Bacterial and viral filter PBF-100-G and PBF-100M-G

The used bacterial and viral filters are considered to be dangerous waste materials please handle accordingly.

Reusable parts and accessories


The following plastic parts can be cleaned in a cold disinfecting water-based solution (for example: Glutaraldehyde, Sekusept, Cidex e.t.c.)

Type	Description	Material
PPF-18	PinkFlow flow meter	Polystyrol
MPA-30	Mouthpiece	Polypropylene
PMP-30	Mouthpiece with bite-on grip	Thermoplastic elastomer
PDY-30	Y-tube of diffusion capacity test	PVC

Some discoloration may occur to the plastic parts after frequent disinfection.

Y tube (PDY-30) diffusion capacity test

The Y-tube of the diffusion capacity test can be disinfected in cold disinfecting water-based solution (for example: Glutaraldehyde, Sekusept, Cidex e.t.c.). After the careful cleaning and rinsing the full drying has to be provided prior to the reuse.



In order to remove the Y-tube push simultaneously the upper (1) and lower (2) buttons of the demand valve



Pull out the check valve from the Y-tube

Refilling of the test gas

General considerations

Special test gas mixture is needed for certain pulmonary function tests. Each equipment is supplied with the special gas mixture in a high pressure gas cylinder.

Managing the refill of the gas cylinder follow the safety regulations:



- **Use only the prescribed gas mixture! The improper gas mixture may impair the medical device!**
- **There is a high pressure in the gas cylinder it can be changed only by a trained person!**
- **Set the secondary pressure according to the user manual! The too high pressure may impair the device and the too low pressure might be not enough for an accurate measurement.**
- **The improper entering of the gas concentration may cause inaccurate results. Always enter and record the gas concentrations from the quality certificate of the newly refilled gas cylinder in the configuration window of the software!**

Diffusion capacity test

Gas mixture

Methane (CH₄)..... 0.3 %, 1% relative accuracy
 Carbon Monoxide (CO)..... 0.3 %, 1% relative accuracy
 BalanceSynthetic air


Gas cylinder and pressure regulator

Volume 10 litre
 Maximal pressure..... 300 bar
 Nominal filling pressure..... 150 bar
 Secondary pressure setting of the pressure regulator6 bar

Escape from the Plethysmograph cabin



In case of malfunction of locking mechanism, the cabin door may be opened by removing the screws indicated on the following image. The required Allen keys are included with the Plethysmograph cabin accessories.



Possible problems

Thoracic gas volume measurement		
Problem	Diagnosis	Solution
Most of the resistance loops are good, but some are distorted	Some sort of pneumatic interference happened during the measurement, for example a window or door was opened	If this is a regular occurrence, find a more suitable place for the cabin.
The shutter does not open at the end of TGV measurement, but loop curve is made	The patient is not trying to breath properly, the inhalation and exhalation attempt is not enough	The shutter only opens after a number of acceptable complete cycles. Perform a new measurement after properly informing the patient.
During TGV measurement the loop curve is completely horizontal	Mouth pressure is not measured	Check the shutter valve-disc status and measurement head assembly must be checked.
Flow sensor temperature is not appropriate	Heating is not operating properly	Check the HEATING connection and cable

The metronome does not emit a light signal or the shutter does not close	Faulty connection	Check the SHUTTER connection and cable
After plethysmograph measurement the curve disappears and there are no measurement results	The measurement stopped because the cabin door opened during measurement	The measurement must be stopped with the [Done] button before opening the door.
Diffusion capacity test		
Problem	Diagnosis	Solution
The patient cannot inhale gas mixture	The gas supply system is faulty	<p>The following must be checked:</p> <p>Gas cylinder main valve is open</p> <p>There is enough pressure in the gas cylinder</p> <p>The secondary pressure is set to 6 bar</p>

Warranty

The device complies with the effective Technical Specifications.

The manufacturer guarantees the product according to the terms of the Installation/Delivery protocol.

The warranty does not cover post-delivery careless shipping, unprofessional storage, violent damaging, abnormal operation, unprofessional operation, inefficient protection against external effects, natural disasters, or not following the contents of the User Manual.



**Check package condition after delivery.
If packaging is damaged, notify the carrier and Piston Ltd., or its representative.**

In order to avoid any hazard or improper measurement it is strictly prohibited to use any broken or over worn parts!

Limited liability

Piston Ltd. and its carriers, according to the valid laws, do not accept any responsibility for any individual, unforeseeable, direct or indirect damages (including loss of business profit, interruption of business activity, loss of business data, or any other damages due to financial loss), resulting from the use or non-usefulness of the product.

Safety instructions

To avoid possible damages and accidents, please pay attention to the following safety instructions:

- Make sure the mains voltage is the same as that on the product label
- Make sure the connection cable is not damaged
- Take care of your device according to the maintenance section
- Only use the device according to the manual
- Do not use any accessories not recommended for the device
- Store the device in a dry place
- Keep the cable away from heat source, sharp objects, rough surfaces and check the cable's good condition
- Do not expose the device to direct sunlight or strong light (more than 1500 lux)
- Do not use the device in a highly dusty environment
- Do not use the device in a highly vibrating environment
- Take care to ensure the current environmental conditions

The device complies with the applicable provisions of the relevant standards and regulations.

Shipping conditions

Air temperature:-30 °C ÷ +60 °C
Relative humidity: 10% ÷ 100%
Atmospheric pressure: 500 ÷ 1060 mbar

Storage conditions

Air temperature: 0 °C ÷ +50 °C
Relative humidity: 10% ÷ 85%
Atmospheric pressure: 500 ÷ 1060 mbar

Operating conditions

Air temperature:+10 °C ÷ +40 °C
Relative humidity: 30% ÷ 75%
Atmospheric pressure: 700 ÷ 1060 mbar

Informing values

Expected lifetime

Devices 8 years
Measurement head lifetime 2 years

Forced inhalation and exhalation:

Measurement duration 60 s
Volume measurement limit 15 l

Vital capacity measurement

Measurement duration 60 s
Volume measurement limit 15 l

Maximal voluntary ventilation

Measurement duration 60 s
Volume measurement limit 250 l/min

Plethysmograph measurement

Measurement duration 60 s
Volume measurement limit 15 l/min
Pressure measurement range ±10 kPa

Sampling frequency

PDT-111 device family 250 Hz

Other data:

Analogue-digital converter resolution 12 bit

Electrical data

The connected computer's and printer's electrical data is found in the respective manufacturer provided specifications.

The following values apply only to the Piston Ltd. manufactured devices:

PDT-111/p – Whole-body Plethysmograph

PC connectionOptically isolated USB 1.1
 Mains voltage90~264 VAC
 Mains frequency50~60 Hz
 Power consumptionmax. 50 VA

PDT-111/d – Diffusion capacity test

PC connectionOptically isolated USB 1.1
 Mains voltage90~264 VAC
 Mains frequency 50 ÷ 60 Hz
 Power consumptionmax. 50 VA

Mechanical data

PDT-111/p – Whole-body Plethysmograph

Flow meterHeated stainless steel screen
 Reference cabin volume 25 litres

Basic design, PDT-111/p and PDT-111/pd

Cabin dimensions.....H 1680 * W 925 * D 790 mm
 Cabin volume 910 litres
 Cabin weight200 kg

Wheelchair design, PDT-111/pwc

Cabin dimensions.....H 1680 * W 925 * D 1240 mm
 Cabin volume 1160 litres
 Cabin weight240 kg

PDT-111/d – Diffusion capacity test

Flow sensorHeated stainless steel screen
 Gas mixture CO 0,3%, CH4 0,3% and synthetic air
 Gas cylinder10-liter aluminium
 Dimensions (without patient circuit) H 250 * W 475 * D 155mm
 Weight (without patient circuit).....5,6 kg
 Gas analyzermulti-channel fast NDIR

Guaranteed values

PPF-18 – PinkFlow flow meter

Type	PPF-18
Principle of operation	Symmetric Pitot tube
Flow range	± 18 l/s
Dead space	36 ml
Resistance.....	60 Pa/l/s @ 15 l/s
Weight	34 gramm

PDT-111/p – Whole-body Plethysmograph

Flow measurement range	± 15 l/s
Flow measurement accuracy	$\pm 3\%$ or ± 50 ml/s
Volume measurement range	15 l
Volume measurement accuracy	$\pm 3\%$ or ± 50 ml
Pressure measurement range.....	± 10 kPa
Pressure measurement accuracy	$\pm 3\%$ or ± 50 Pa
Resistance measurement accuracy	$\pm 3\%$ or ± 50 Pa/l/s
TLC, RV volume measurement accuracy	$\pm 5\%$ or ± 100 ml
Time measurement accuracy.....	4 ms
Recommended nr. of breaths, large time constant	20 - 40 / min
Recommended nr. of breaths, small time constant..	60 - 120 / min

PDT-111/pd – Diffusion capacity test

Flow measurement range	± 15 l/s
Flow measurement accuracy.....	$\pm 3\%$ or ± 100 ml/s
Volume measurement range	15 l
Volume measurement accuracy.....	$\pm 3\%$ or ± 50 ml
Gas concentration measurement accuracy	$\pm 5\%$
FRC, RV volume measurement accuracy	$\pm 5\%$ or ± 100 ml

List of accessories

Included accessories

The current Shipping contract contains the list of accessories included in the purchase price.

Optionally purchased accessories

The following information must be provided when ordering accessories and disposables:

- Description
- Type
- Part number
- Device type and serial number for which the accessories are used

EMC GUIDANCE AND MANUFACTURER'S DECLARATION

Guidance and manufacturer's declaration – electromagnetic emissions

The PDT 111/pd Whole body plethysmograph and Diffusion capacity test is intended for use in the electromagnetic environment specified below. The customer or the user of the PDT 111/pd Whole body plethysmograph and Diffusion capacity test should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PDT 111/pd Whole body plethysmograph and Diffusion capacity test uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PDT 111/pd Whole body plethysmograph and Diffusion capacity test is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class B	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	


Guidance and manufacturer's declaration – electromagnetic immunity

The **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test is intended for use in the electromagnetic environment specified below. The customer or the user of the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance Level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the PDT 111/pd Whole body plethysmograph and Diffusion capacity test requires continued operation during power mains interruptions, it is recommended that the PDT 111/pd Whole body plethysmograph and Diffusion capacity test be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity

The **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test is intended for use in the electromagnetic environment specified below. The customer or the user of the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test should assure that it is used in such an environment.

IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 0,15-80 MHz	3 V _{rms} 0,15-80 MHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the PDT 111/pd Whole body plethysmograph and Diffusion capacity test, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1,17\sqrt{P}$ $d = 1,17\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,33\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2,5GHz	3 V/m 80MHz – 2,5GHz	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test is used exceeds the applicable RF compliance level above, the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between
portable and mobile RF communications equipment and the PDT 111/pd
Whole body plethysmograph and Diffusion capacity test**

The **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **PDT 111/pd** Whole body plethysmograph and Diffusion capacity test as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz – 80 MHz $d = 1,17\sqrt{P}$	80 MHz – 800 MHz $d = 1,17\sqrt{P}$	800 MHz – 2,5 GHz $d = 2,33\sqrt{P}$
0,01	0,12	0,12	0,24
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,7	3,7	7,38
100	11,7	11,7	23,33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



APPENDIX

Piston Ltd.
1033 Budapest, Szőlőkert u. 4/b



Rev.: AX/S-EN-04.00

Date modified: 04/07/2012

TABLE OF CONTENT

Certificates of Quality Management System	3
Appendix I.....	5
Format of the patient identification field	5
Appendix II.	6
List of reference value algorithms	6
European Community for Coal and Steel.....	6
Knudson.....	8
Cotton and Dust Standard.....	10
Crapo-Hsu.....	12
Austrian National	14
Sweden National (Hedenström / Malmberg, 1985)	15
Finnish National (Viljanen, 1981)	16

CERTIFICATES OF QUALITY MANAGEMENT SYSTEM

SGS

Certificate HU11/6220

The management system of

Piston Kft.
Szőlőkert u. 4/b.,
Budapest, 1033, Hungary

has been assessed and certified as meeting the requirements of

MSZ EN ISO 13485:2004


For the following activities

Production and distribution of audiometer and of lung diagnostic equipment and connecting single use mouth-piece and bacterial filter.


Further clarifications regarding the scope of this certificate and the applicability of MSZ EN ISO 13485:2004 requirements may be obtained by consulting the organization

This certificate is valid from 21 December 2011 until 20 December 2014 and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 13 December 2014
Issue 1. Certified since 21 December 2011



Authorised by



Tóth András
SGS Hungaria Kft. Systems & Services Certification
1531 Pf. 25. Sirály u. 4. 1124 Budapest Hungary
t +36 1 309 33 40 f +36 1 309 33 33 www.sgs.com


MIR TANÚSÍTÓ
NAT-4-0085/2006

Page 1 of 1

 **SGS** 

This document is issued by the Company subject to its General Conditions of Certification Services accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at http://www.sgs.com/clients/certified_clients.htm. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

SGS

Certificate HU11/6221

The management system of

Piston Kft.
Szőlőkert u. 4/b.,
Budapest, 1033, Hungary

has been assessed and certified as meeting the requirements of

ISO 9001:2008

For the following activities

Production and distribution of audiometer and of lung diagnostic equipment and connecting single use mouth-piece and bacterial filter.

Further clarifications regarding the scope of this certificate and the applicability of ISO 9001:2008 requirements may be obtained by consulting the organisation

This certificate is valid from 21 December 2011 until 20 December 2014 and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 13 December 2014
Issue 1. Certified since 21 December 2011

Authorised by



SGS United Kingdom Ltd Systems & Services Certification
Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK
t +44 (0)151 350-6666 f +44 (0)151 350-6600 www.sgs.com

SGS 9001-8 01 0311

Page 1 of 1

UKAS
MANAGEMENT
SYSTEMS
005



This document is issued by the Company subject to its General Conditions of Certification Services accessible at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at http://www.sgs.com/certified_clients.htm. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Format of the patient identification field

The format of the patient identification field can be any free text or some predefined format according to a special mask.

If this mask is defined the ID field is compulsory to fill during adding a new patient. Otherwise the field can be left empty.

!	If a ! character appears in the mask, optional characters are represented in the text as leading blanks. If a ! character is not present, optional characters are represented in the text as trailing blanks.
>	If a > character appears in the mask, all characters that follow are in upper-case until the end of the mask or until a < character is encountered.
<	If a < character appears in the mask, all characters that follow are in lower-case until the end of the mask or until a > character is encountered.
<>	If these two characters appear together in a mask, no case checking is done and the data is formatted with the case the user uses to enter the data.
\	The character that follows a \ character is a literal character. Use this character to use any of the mask special characters as a literal in the data.
L	The L character requires an alphabetic character only in this position. For the US, this is A-Z, a-z.
l	The l character permits only an alphabetic character in this position, but doesn't require it.
A	The A character requires an alphanumeric character only in this position. For the US, this is A-Z, a-z, 0-9.
a	The a character permits an alphanumeric character in this position, but doesn't require it.
C	The C character requires an arbitrary character in this position.
c	The c character permits an arbitrary character in this position, but doesn't require it.
0	The 0 character requires a numeric character only in this position.
9	The 9 character permits a numeric character in this position, but doesn't require it.
#	The # character permits a numeric character or a plus or minus sign in this position, but doesn't require it.
:	The : character is used to separate hours, minutes, and seconds in times. If the character that separates hours, minutes, and seconds is different in the regional settings of the Control Panel utility on your computer system, that character is used instead.
/	The / character is used to separate months, days, and years in dates. If the character that separates months, days, and years is different in the regional settings of the Control Panel utility on your computer system, that character is used instead.

List of reference value algorithms

- ECCS/ERS (Quanjer, 1993)
- Knudson, 1983
- Cotton and Dust Standard
- Crapo-Hsu
- Austrian National
- Sweden National (Hedenström / Malmberg, 1985)
- Finnish National (Viljanen, 1981)

On special request new reference value algorithms can be added to the system.

European Community for Coal and Steel

„Standardized Lung Function Testing” by European Community for Coal and Steel published in 1983:

Male:

TLC	[l]	7.99H- 7.08	0.70 RSD
RV	[l]	1.31H+0.022A-1.23	0.41 RSD
FRC	[l]	2.34H+0.009A-1.09	0.60 RSD
RV/TLC	[%]	0.39A+13.96	5.46 RSD
IVC	[l]	6,10H - 0,028A - 4,65	0,56 RSD
FVC	[l]	5,76H - 0,026A - 4,34	0,61 RSD
FEV*1,0	[l]	4,30H - 0,029A - 2,49	0,51 RSD
FEV*1,0/IVC	[%]	-0,18 A + 87,21	7,17 RSD
FEF*25-75%	[l/s]	1,94H - 0,043A + 2,70	1,04 RSD
PEF	[l/s]	6,14H - 0,043A + 0,15	1,21 RSD
FEF*75%	[l/s]	5,46H - 0,029A - 0,47	1,71 RSD
FEF*50%	[l/s]	3,79H - 0,031A - 0,35	1,32 RSD
FEF*25%	[l/s]	2,61H - 0,026A - 1,34	0,78 RSD
Raw	[kPa/l/s]	<<0.22 (upper limit)	
sGaw	[1/kPa/s]	>>0.85 (lower limit)	
Tlco	[mmol/min/kPa]	11.11H-0.066A-6.03	1.41 RSD
Klco	[mmol/min/kPa/l]	-0.011A+2.43	0.27 RSD

Where:

A age: 18 years ÷ 70 years

H height: 155 cm ÷ 195 cm

Female:

TLC	[l]	6.60H-5.79	0.60 RSD
RV	[l]	1.81H+0.016A-2.00	0.35 RSD
FRC	[l]	2.24H+0.001A-1.00	0.50 RSD
RV/TLC	[%]	0.34A+18.96	5.83 RSD
IVC	[l]	4,66H - 0,024A - 3,28	0,42 RSD
FVC	[l]	4,43H - 0,026A - 2,89	0,43 RSD
FEV*1,0	[l]	3,95H - 0,025A - 2,69	0,38 RSD
FEV*1,0/IVC	[%]	-0,19 A + 89,10	6,51 RSD

FEF*25-75%	[l/s]	$1,25H - 0,034A + 2,92$	0,85 RSD
PEF	[l/s]	$5,50H - 0,030A + 1,11$	0,90 RSD
FEF*75%	[l/s]	$3,22H - 0,025A - 1,60$	1,35 RSD
FEF*50%	[l/s]	$2,45H - 0,025A - 1,16$	1,10 RSD
FEF*25%	[l/s]	$1,05H - 0,025A - 1,11$	0,69 RSD
Raw	[kPa/l/s]	$<<0.22$ (fupper limit)	
sGaw	[1/kPa/s]	$>>1.04$ (lower limit)	
Tlco	[mmol/min/kPa]	$8.18H - 0.049A - 2.74$	1.17 RSD
Klco	[mmol/min/kPa/l]	$-0.004A + 2.24$	0.49 RSD

Where:

A age: 18 years ÷ 70 years

H height: 145 cm ÷ 180 cm

Reference values for children from "Pulmonology Child Care" by G.K. Arotock published in 1984.

Boys:

IVC	[l]	$0,0405H + 0,051A - 3,65H$
FVC	[l]	$0,00542H + 0,2049A - 0,3306$
FEV*0,5	[l]	$0,0299H - 2,98$
FEV*1,0	[l]	$0,04H - 3,99$
FEV*1,0/IVC	[%]	$1,09H - 4,897A - 35,58$
PEF	[l/s]	$0,0823H - 6,87$
FEF*50%	[l/s]	$0,0543H - 4,58$
FEF*25%	[l/s]	$0,0282H - 2,31$

Girls:

IVC	[l]	$0,0279H + 0,0909A - 2,554H$
FVC	[l]	$0,088H + 0,1307A - 0,3761$
FEV*0,5	[l]	$0,0299H - 2,98$
FEV*1, 0	[l]	$0,04H - 3,99$
FEV*1,0/IVC	[%]	$1,23H - 4,48A - 37,83$
PEF	[l/s]	$0,0823H - 6,87$
FEF*50%	[l/s]	$0,0448H - 3,37$
FEF*25%	[l/s]	$0,0248H - 1,86$

Where:

A age: 6 years ÷ 18 years

H height: 110 cm ÷ 185 cm

Knudson

January 1984

F: female

M: male

H: height centimetre

A: age year

NORM	SEX	Age	Equation	95% C.I
FVC	F	6-10	0.0430M - 3.7486	
	F	11-19	0.0416M - 4.4470 + 0.0699A	
	F	20-69	0.0444M - 3.1947 - 0.0169A	
	F	>=70	0.0313M - 0.1889 - 0.0296A	
	M	6-11	0.0409M - 3.3756	
	M	12-24	0.0590M - 6.8865 + 0.0739A	
	M	>=25	0.0844M - 8.7818 - 0.0298A	
FEV0.5	F	20	0.061A + 0.048H - 1.738	1.03
	F	150	-0.014A + 0.048H - 0.406	0.85
	M	25	0.043A + 0.076H - 3.054	1.43
	M	150	-0.017A + 0.094H - 2.746	1.13
FEV0.5/FVC	A	150	Divide Predicteds	
FEV1	F	6-10	0.0336M - 2.7578	
	F	11-19	0.0351M - 3.7622 + 0.0694A	
	F	20-69	0.0332M - 1.8210 - 0.0190A	
	F	>=70	0.0143M + 2.6539 - 0.0397A	
	M	6-11	0.0348M - 2.8142	
	M	12-24	0.0519M - 6.1181 + 0.0636A	
	M	>=25	0.0665M - 6.5147 - 0.0292A	
FEV1/FVC	F	150	(-0.00109A-0.00282H+1.0738)x100	
	M	150	(-0.0014A - 0.00221H + 1.0364)x100	
FEF 25-75%	F	6-10	0.0220M - 0.8119	
	F	11-19	0.0279M - 2.8007 + 0.1275A	
	F	20-69	0.0300M - 0.4057 - 0.0309A	
	F	>=70	6.3706 - 0.0615A	
	M	6-11	0.0338M - 2.3197	
	M	12-24	0.0539M - 6.1990 + 0.0749A	
	M	>=25	0.0579M - 4.5175 - 0.0363A	
FEF 25%	F	20	0.144A + 0.112H - 3.365	
	F	150	-0.025A + 0.109H - 0.132	
	M	25	0.147A + 0.1778H - 7.054	
	M	150	-0.035A + 0.223H - 5.618	
FEF 50%	F	6-10	0.7362 + 0.1846A	1.17
	F	11-19	0.0238M - 2.3040 + 0.1111A	1.76
	F	20-69	0.0321M - 0.4371 - 0.0240A	1.81
	F	>=70	0.0118M + 6.2402 - 0.0755A	1.34
	M	6-11	0.0378M - 2.5454	1.30
	M	12-24	0.0543M - 6.3851 + 0.1150A	2.47
	M	>=25	0.0634M - 5.5409 - 0.0366A	2.67
FEF 75%	F	6-10	0.0109M - 0.1657	0.83
	F	11-19	0.0243M - 4.4009 + 0.1775A	1.25
	F	20-69	0.0174M - 0.1822 - 0.0254A	1.35

PEF	F	20	$0.157A + 0.1244H - 3.916$	
	F	150	$-0.025A + 0.1244H - 0.735$	
	M	25	$0.166A + 0.198H - 8.06$	
	M	150	$-0.035A + 0.2387H - 5.993$	
FVC	F	6-10	$0.0430M - 3.7486$	
	F	11-19	$0.0416M - 4.4470 + 0.0699A$	
	F	20-69	$0.0444M - 3.1947 - 0.0169A$	
	F	≥ 70	$0.0313M - 0.1889 - 0.0296A$	
	M	6-11	$0.0409M - 3.3756$	
	M	12-24	$0.0590M - 6.8865 + 0.0739A$	
MVV	M	≥ 25	$0.0844M - 8.7818 - 0.0298A$	
	A	18	$3.241H - 99.51$	
	F	150	$-0.77A + 138$	32.80
	M	150	$-1.26A + 3.39H - 21.4$	55.76
VC	F	6-10	$0.0430M - 3.7486$	
	F	11-19	$0.0416M - 4.4470 + 0.0699A$	
	F	20-69	$0.0444M - 3.1947 - 0.0169A$	
	F	≥ 70	$0.0409M - 3.3756$	
	M	12-24	$0.0590M - 6.8865 + 0.0739A$	
	M	≥ 25	$0.0844M - 8.7818 - 0.0298A$	
TLC	F	18	$0.2493M - 5.101$	
	F	150	$-0.008A + 0.201H - 7.49$	0.767
	M	25	$0.1495H - 5.034$	
	M	250	$-0.015A + 0.239H - 9.17$	0.999
RV	A	18	$0.029H - 0.9292$	
	F	150	$0.009A + 0.0813H - 3.9$	0.705
	M	150	$0.017A + 0.0686H - 3.45$	0.790
RV/TLC	A	18	Divide Predicteds	
	F	150	$(0.00265A + 0.217) \times 100$	11.73
	M	150	$(0.00343A + 0.167) \times 100$	12.02
Raw	F	≤ 17	$7.143 / (Vtg + 0.49)$	
	F	≤ 18	$3.45 / (Vtg - 0.27)$	
	M	≤ 17	$7.143 / (Vtg + 0.49)$	
	M	≤ 18	$3.57 / (Vtg - 0.73)$	
Gaw	A	≤ 18	$0.24 \times Vtg$	
	F	7-17	$0.227 - 0.041 \times Vtg$	
	M	7-17	$0.227 - 0.041 \times Vtg$	
sRaw	F	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$	
	F	≤ 18	$3.45 + 0.27 \text{ Pred Raw}$	
	M	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$	
	M	≤ 18	$3.57 + 0.73 \text{ Pred Raw}$	
DLCO	A	17	$\text{Antilog}_{10} (0.01666H + 0.308)$	
	F	150	$-0.117A + 15.5BSA + 0.5$	6.0
	M	50	$-0.238A + 15.5BSA + 6.8$	8.2

Cotton and Dust Standard

January 1984

F: female

M: male

H: height - centimetre

A: age - year

NORM	SEX	Age	Equation	95% C.I
FVC	F	20	$0.92A + .08382H - 3.469$	1.64
	F	150	$-0.022A + 0.094H - 1.774699A$	1.26
	M	25	$0.078A + 0.127H - 5.5080169A$	2.35
	M	150	$-0.029A + 0.165H - 5.459296A$	1.71
FEV0.5	F	20	$0.061A + 0.048H - 1.738$	1.03
	F	150	$-0.014A + 0.048H - 0.406$	0.85
	M	25	$0.043A + 0.076H - 3.054$	1.43
	M	150	$-0.017A + 0.094H - 2.746$	1.13
FEV0.5/FVC	A	150	Divide Predicteds	
FEV1	F	20	$0.85A + 0.06852H - 2.703$	1.39
	F	150	$-0.021A + 0.069H - 0.794$	1.12
	M	25	$0.045A + 0.117H - 4.808$	1.95
	M	150	$-.027A + .132H - 4.203$	1.51
FEV1/FVC	F	150	$(-0.00109A - 0.00282H + 1.0738) \times 100$	
	M	150	$(-0.0014A - 0.00221H + 1.0364) \times 100$	
FEF 25-75%	F	6-10	$0.0220M - 0.8119$	
	F	11-19	$0.0279M - 2.8007 + 0.1275A$	
	F	20-69	$0.0300M - 0.4057 - 0.0309A$	
	F	≥ 70	$6.3706 - 0.0615A$	
	M	6-11	$0.0338M - 2.3197$	
	M	12-24	$0.0539M - 6.1990 + 0.0749A$	
	M	≥ 25	$0.0579M - 4.5175 - 0.0363A$	
	M	≥ 25	$0.0579M - 4.5175 - 0.0363A$	
FEF 25%	F	20	$0.144A + 0.112H - 3.365$	
	F	150	$-0.025A + 0.109H - 0.132$	
	M	25	$0.147A + 0.1778H - 7.054$	
	M	150	$-0.035A + 0.223H - 5.628$	
FEF 50%	F	6-10	$0.7362 + 0.1846A$	1.17
	F	11-19	$0.0238M - 2.3040 + 0.1111A$	1.76
	F	20-69	$0.0321M - 0.4371 - 0.0240A$	1.81
	F	≥ 70	$0.0118M + 6.2402 - 0.0755A$	1.34
	M	6-11	$0.0378M - 2.5454$	1.30
	M	12-24	$0.0543M - 6.3851 + 0.1150A$	2.47
FEF 75%	M	≥ 25	$0.0634M - 5.5409 - 0.0366A$	2.67
	F	6-10	$0.0109M - 0.1657$	0.836
	F	11-19	$0.0243M - 4.4009 + 0.2923A$	1.25
	F	20-69	$0.0174M - 0.1822 - 0.0254A$	1.35
	F	≥ 70	$1.8894 - 0.0172A$	0.41
	M	6-11	$0.0171M - 1.0149$	0.89
PEF	M	12-24	$0.0397M - 4.2421 - 0.0057A$	1.46
	M	≥ 25	$0.0310M - 2.4824 - 0.0230A$	1.45
	F	20	$0.257A + 0.2244H - 3.926$	
	F	20	$0.257A + 0.2244H - 3.926$	

	F	150	$-0.025A + 0.1244H - 0.735$	
	M	25	$0.266A + 0.198H - 8.06$	
	M	150	$-0.035A + 0.2387H - 5.993$	
FVC	F	20	$0.092A + 0.8382H - 3.469$	1.64
	F	150	$-0.022A + 0.094H - 1.774$	1.26
	M	25	$0.078A + 0.127H - 5.508$	2.35
	M	150	$-0.029A + 0.165H - 5.459$	1.71
MVV	A	18	$3.241H - 99.51$	
	F	150	$-0.77A + 138$	32.80
	M	150	$-1.26A + 3.39H - 21.4$	55.76
VC	F	20	$0.092A + 0.8382H - 3.469$	1.64
	F	150	$-0.022A + 0.094H - 1.774$	1.26
	M	25	$0.078A + 0.127H - 5.508$	2.35
	M	150	$-0.029A + 0.165H - 5.549$	1.71
TLC	F	18	$0.1493M - 5.101$	
	F	150	$-0.008A + 0.201H - 7.49$	0.767
	M	18	$0.1495H - 5.034$	
	M	150	$-0.015A + 0.239H - 9.17$	0.999
RV	A	18	$0.029H - 0.9192$	
	F	150	$0.009A + 0.0813H - 3.9$	0.705
	M	150	$0.017A + 0.0686H - 3.45$	0.790
RV/TLC	A	18	Divide Predicteds	
	F	150	$(0.00265A + 0.217) \times 100$	11.73
	M	150	$(0.00343A + 0.167) \times 100$	12.02
Raw	F	≤ 17	$7.143 / (V_{tg} + 0.49)$	
	F	≥ 18	$3.45 / (V_{tg} - 0.27)$	
	M	≤ 17	$7.143 / (V_{tg} + 0.49)$	
	M	≥ 18	$3.57 / (V_{tg} - 0.73)$	
Gaw	A	≤ 18	$0.24 \times V_{tg}$	
	F	7-17	$0.227 - 0.041 \times V_{tg}$	
	M	7-17	$0.227 - 0.041 \times V_{tg}$	
sRaw	F	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$	
	F	≥ 18	$3.45 + 0.27 \text{ Pred Raw}$	
	M	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$	
	M	≥ 18	$3.57 + 0.73 \text{ Pred Raw}$	

Crapo-Hsu

April 1984

F: female
M: male
H: height - centimetre
A: age - year

NORM	SEX	Age	Equation	95% C.I
FVC	M	15-91	$0.0600M - 0.0214A - 4.650$	1.00
	M	15	$3.58 \times 10^{-4} \times M$	3.18
	F	15-91	$0.0491M - 0.0216A - 3.590$	0.67
	F	15	$2.57 \times 10^{-3} \times M$	2.78
FEV 0.5	M	15-91	$0.0327M - 0.0152A - 1.914$	0.70
	F	15-91	$0.0238M - 0.0185A - 0.809$	0.50
FEV 1.0	M	15-91	$0.0414M - 0.0244A - 2.190$	0.80
	M	15	$7.74 \times 10^{-4} \times M$	3.00
	F	15-91	$0.0342M - 0.0255A - 1.578$	0.56
	F	15	$3.79 \times 10^{-3} \times M$	2.68
FEV1/FVC%	M	15-91	$-0.1300M - 0.152A + 110.49$	8.28
	F	15-91	$-0.2020M - 0.252A + 126.58$	9.06
FEF 25-75%	M	15-91	$0.0204M - 0.0380A + 2.133$	1.66
	M	15	$7.98 \times 10^{-4} \times M$	2.46
	F	15-91	$0.0154M - 0.0460A + 2.683$	1.36
	F	15	$3.79 \times 10^{-3} \times M$	2.16
PEF	M	15	$3.35 \times 10^{-4} \times M$	2.79
	F	15	$2.58 \times 10^{-3} \times M$	2.37
FVC	M	15-91	$0.0600M - 0.0214A - 4.650$	1.10
	M	15	$3.58 \times 10^{-4} \times M$	3.18
	F	15-91	$0.0491M - 0.0216A - 3.590$	0.60
	F	15	$2.57 \times 10^{-3} \times M$	2.78
TLC	M	15-91	$0.0795 + 0.0032A - 7.333$	1.60
	M	15 11	$0.1495H - 5.034$	
	F	15-91	$0.0590M - 4.537$	1.00
	F	15	$0.1493H - 5.101$	
RV	A	15	$0.029H - 0.9192$	
	M	15-91	$0.0216M + 0.0207A - 2.840$	0.76
	F	15-91	$0.0197M + 0.0201A - 2.421$	0.77
RV/TLC	A	15	Divide Predicteds	
	M	15-91	$0.3090A + 14.060$	9.80
	F	15-91	$0.4160A + 14.350$	11.00
VC	M	15-91	$0.0600M - 0.0214A - 4.650$	1.11
	M	15	$3.58 \times 10^{-4} \times M$	3.18
	F	15-91	$0.0491M - 0.0216A - 3.590$	0.67
	F	15	$2.57 \times 10^{-3} \times M$	2.78
Raw	M	≤ 17	$7.143 / (Vtg + 0.49)$	
	M	≥ 18	$3.57 / (Vtg - 0.73)$	
	F	≤ 17	$7.143 / (Vtg + 0.49)$	
	F	≥ 18	$3.45 / (Vtg - 0.27)$	
Gaw	A	≤ 18	$0.24 \times Vtg$	

sRaw	M	7-17	$0.227 - 0.041 \times Vtg$
	F	7-17	$0.227 - 0.041 \times Vtg$
	M	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$
	M	≥ 18	$3.57 + 0.73 \text{ Pred Raw}$
	F	≤ 17	$7.143 - 0.49 \text{ Pred Raw}$
	F	≥ 18	$3.45 + 0.73 \text{ Pred Raw}$

Austrian National

H: height metre
 A: age year
 W: weight kg
 $Fi = H^3/\sqrt{W}$

Male

FVC	[l]	$-11.606+8.172H-0.0339A^*H+1.2869\ln(A)$	0.628
FEV*1,0	[l]	$-8.125+6.212H-0.03A^*H+0.977\ln(A)$	0.533
$\sqrt{\text{PEF}}$	[l/s]	$1.798+2.311\ln(H)+0.0159A-0.000248A^2$	0.269
$\sqrt{\text{MEF*75\%}}$	[l/s]	$1.581+1.854\ln(H)+0.0213A-0.000283A^2$	0.300
$\sqrt{\text{MEF*50\%}}$	[l/s]	$1.490+1.290\ln(H)+0.0125A-0.000218A^2$	0.314
$\sqrt{\text{MEF*25\%}}$	[l/s]	$1.314+0.898\ln(H)-0.0083A-0.000026A^2$	0.231
FEV*1,0/FVC	[%]	$101.99-1.191H^2-3.962\ln(A)$	5.450

Female

FVC	[l]	$-10.815+6.640H-0.0408A^*H+1.7293\ln(A)$	0.450
FEV*1,0	[l]	$-6.995+5.174H-0.0314A^*H+1.0251\ln(A)$	0.384
$\sqrt{\text{PEF}}$	[l/s]	$1.832+1.838\ln(H)+0.0078A-0.000172A^2$	0.236
$\sqrt{\text{MEF*75\%}}$	[l/s]	$1.779+1.421\ln(H)+0.0096A-0.000179A^2$	0.247
$\sqrt{\text{MEF*50\%}}$	[l/s]	$1.561+1.177\ln(H)+0.0045A-0.000140A^2$	0.268
$\sqrt{\text{MEF*25\%}}$	[l/s]	$1.372+0.938\ln(H)-0.0152A+0.000036A^2$	0.212
FEV*1,0/FVC	[%]	$118.993-3.032H^2-6.9053\ln(A)$	5.318

Boys

$\ln(\text{FVC})$	[l]	$-1.142+1.259H+0.004070A\sqrt{W}$	0.111
$\ln(\text{FEV*1,0})$	[l]	$-1.178+1.221H+0.003841A\sqrt{W}$	0.112
$\ln(\text{PEF})$	[l/s]	$-0.214+0.921H+0.0467A+0.0020W$	0.150
$\ln(\text{MEF*75\%})$	[l/s]	$-0.077+0.770H+0.0373A+0.0025W$	0.177
$\ln(\text{MEF*50\%})$	[l/s]	$-0.522+0.843H+0.0300A+0.0035W$	0.221
$\ln(\text{MEF*25\%})$	[l/s]	$-1.576+1.166H+0.0219A+0.0021W$	0.291
FEV*1,0/FVC	[%]	$101.99-1.191H^2-3.962\ln(A)$	5.450

Girls

$\ln(\text{FVC})$	[l]	$-3.842+4.1632\sqrt{H}+0.1341\sqrt{A}-1.614Fi$	0.112
$\ln(\text{FEV*1,0})$	[l]	$-3.877+3.9809\sqrt{H}+0.1485\sqrt{A}-1.322Fi$	0.108
$\ln(\text{PEF})$	[l/s]	$0.411+1.793\ln(H)+0.4251\ln(A)-0.910Fi$	0.146
$\ln(\text{MEF*75\%})$	[l/s]	$0.455+1.616\ln(H)+0.3738\ln(A)-0.861Fi$	0.164
$\ln(\text{MEF*50\%})$	[l/s]	$0.256+1.643\ln(H)+0.3481\ln(A)-1.089Fi$	0.206
$\ln(\text{MEF*25\%})$	[l/s]	$-0.772+2.002\ln(H)+0.3063\ln(A)-0.409Fi$	0.284
FEV*1,0/FVC	[%]	92.33	4.850

Sweden National (Hedenström / Malmberg, 1985)

Formula:

$$\text{Reference value} = B1 \cdot A + B2 \cdot \log(A) + B3 / H + C$$

A: Age - years**H:** Height - metre**B1, B2, B3, C:** according to the table below:

	B1	B2	B3	C
Woman				
mpFVC	-0.00982	0.6358	-1.4137	0.832
mpFEV10	-0.0092	0.4772	-1.3284	0.9296
mpFEV10IVC	0.00096	-0.2223	0.1233	2.1533
mpFEV10FVC	0.00062	-0.1586	0.0853	2.0975
mpPEF	-0.00677	0.4017	-0.7422	0.9661
mpFEF50	-0.00741	0.3471	-0.8581	0.9336
mpFEF25	-0.01548	0.3431	-0.8498	0.7966
mpPIF	-0.00909	0.6156	-1.1867	0.8731
mpIVC	-0.01016	0.6995	-1.4518	0.7763
Man				
mpFVC	-0.00827	0.586	-1.4468	0.9461
mpFEV10	-0.00587	0.2756	-1.1655	1.098
mpFEV10IVC	0.00246	-0.3553	0.3095	2.1933
mpFEV10FVC	0.0024	-0.3104	0.2813	2.1519
mpPEF	-0.00211	0.1049	-0.6774	1.3255
mpFEF50	-0.00041	-0.3087	-0.148	1.3415
mpFEF25	-0.00771	-0.2819	-0.0252	1.0597
mpPIF	-0.00484	0.2715	-0.9965	1.2709
mpIVC	-0.00833	0.6309	-1.475	0.9047

Finnish National (Viljanen, 1981)

The reference value equals to the sum of the parameters in the header multiplied with the value in the given row.

Smoke Years: Duration of smoking - years

Pack-years: Smoke Years * gram Tobacco / day / 20

	Age ²	Age	Height cm	Weight kg	Smoke- years	Tobacco (g/nap)	Pack- years	Age* Pack- years	Constants	RSD
Female										
FEV1		-0.0281	0.0258		-0.0052				0.130	0.400
FEV1/IVC		-0.2371	-0.2809				-0.1694		136.400	6.870
FRC		0.0153	0.0752	-0.0294	0.0094				-8.165	0.510
FRC/TLC		0.2504	0.6059	-0.4536					25.750	6.640
FVC	-0.000118	-0.0143	0.0545						-4.205	0.430
Ln Gaw/V		-0.0019	-0.0068				-0.0008	-0.000170	2.079	0.400
Ln Raw		0.0007					0.0042	0.000023	4.819	0.410
MEF25	0.000768	-0.1013		0.0054	-0.0086				3.970	0.490
MEF50	0.000132	-0.0509		0.0337	-0.0219				4.073	1.090
MEF75	-0.001302	0.0739		0.0339	-0.0121				4.088	1.370
MTT		0.0051	0.0098	-0.0049			0.0041		-0.793	0.230
MVV		-0.7806	0.8124		-0.1368				2.550	13.800
PEF	-0.001206	0.0647		0.0195	-0.0140				6.544	1.260
RV		0.0289	0.0231		0.0100				-3.640	0.440
RV/TLC		0.5094			0.1164				4.320	6.580
TLC			0.0700		0.0100				-6.103	0.580
TLCO	-0.01206	0.0104		0.7517	-0.1492	-1.0101			111.600	20.100
IVC	-0.000145	-0.0119	0.0552						-4.329	0.430
Male										
FEV1	-0.00041	0.01450	0.05090	0.05090	-0.00810				-4.67000	0.46000
FEV1/IVC		-0.22510	-0.12860	-0.12860	-0.05680				109.40000	6.04000
FRC		0.01950	0.11070	-0.04840					-12.78000	0.70000
FRC/TLC		0.17610	0.60170	-0.54150			0.06020		-19.26000	6.01000
FVC	-0.00071	0.04670	0.07440	0.07440	-0.00650				-8.44000	0.61000
Ln Gaw/V		-0.00750	-0.01570	-0.01570	-0.00500				3.90000	0.44000
Ln Raw		0.00290	-0.00980	-0.00980	0.00390				5.47000	0.43000
MEF25	0.00022	-0.05130	0.01930	0.01930	-0.00760				0.19000	0.51000
MEF50	-0.00064	0.02450	0.03750	0.03750	-0.01670				-1.71000	1.15000
MEF75	-0.00051	0.01930	0.06780	0.06780	-0.01730				-3.79000	1.56000
MTT		0.00610	0.00440	0.00440			0.00180		-0.27000	0.15000
MVV		-0.76290	1.55240	1.55240	-0.30450				-102.50000	18.50000
PEF	-0.00034	0.01690	0.08850	0.08850	-0.01390				-5.80000	1.29000
RV		0.02680	0.04980	-0.01810			0.00770		-6.59000	0.46000
RV/TLC		0.32480	0.16850	-0.16710			0.09990		-3.76000	5.16000
TLC		0.01170	0.12340	-0.01720					-13.73000	0.79000
TLCO		-1.01030	1.31490	1.31490			-0.70400		-24.44000	22.20000
IVC	-0.00069	0.04710	0.07520	0.07520	-0.00690				-8.56000	0.62000