MEDICAL LCD MONITOR

Before connecting, operating or adjusting this product, please read this instruction booklet carefully and completely.

MEDICAL LCD MONITOR FS-L190*D (19") FS-L190*DT (19") FS-L240*D (24") FS-L240*DT (24") FS-L260*D (26") FS-L320*D (32") FS-L420*D (42")

Model definition



- 1. FS : Monitor manufacturer.
- 2. L : Panel manufacturer.
- 3. XXX : Display size.
- 4. * : Signal input option.
- 5. D : Medical grade.
- 6. T : Touch screen installed.

Table of Contents

Symbol Definitions	2
Safety Instructions	4
Caution	7
FCC Information	10
Parts	15
FS-L190*D/FS-L190*DT	
FS-L240*D/FS-L240*DT	16
FS-L260*D	17
FS-L320*D	18
FS-L420*D	19
Connector	
FS-L190*D/FS-L240*D/FS-L260*D/FS-L320*D	20
FS-L190*DT / FS-L240*DT	21
FS-L420*D	22
Mechanical Product Drawing	23
19" FS-L190*D / FS-L190*DT	
24" FS-L240*D / FS-L240*DT	24
26" FS-L260*D	25
32" FS-L320*D	26
42" FS-L420*D	27
Control	
Power management	32
OSD	33
Standard Signal table	48
Signal connector Pin Assignments	49
Specification	
FS-L190*D/FS-L190*DT	52
FS-L240*D/FS-L240*DT	53
FS-L260*D	
FS-L320*D	55
FS-L420*D	56

Symbol Definitions

The following symbols appear on the product, its labels, or the product package. Each symbol carries a special definition, as defined below.



Dangerous : High Voltage.



Consult accompanying documents.



Direct Current.



Indicates protective earth ground.



DC Power control switch.



Serial Number.



Top-Bottom.



Fragile.



Do not get wet.



Maximum Stacking.(19"/24"/26"/32")



3

Maximum Stacking.(42")



Indicates proof of conformity to applicable European Economic Community Council directives and to harmonized standards published in the official journal of the European Communities.



Medical LCD monitor is in accordance with UL 60601-1 and CAN/CSA C22.2 No.601.1 in regards to electric shock, fire hazards, and mechanical hazard.



Tested to comply with FCC Class B standard.



This symbol indicates that the waste of medical LCD monitor must not be disposed as unsorted municipal waste and must be collected separately. Please contact the manufacturer or other authorized disposal company to decommission your medical LCD monitor.

Language : English

Note : A manual with a proper legal language will be provided for a member of EU. Foreseeson Custom Display Inc.language. If you find an English version instead of your language version, call your local distributor or Foreseeson Custom Display Inc. will include the manual written in the local . for the proper version. This declaration is applied to the member of EU where this medical LCD monitor is sold legally.

Safety Instructions

On Safety

- 1. Before connecting the AC power cord to the DC adapter outlet make sure the voltage designation of the DC adapter corresponds to the local electrical supply.
- 2. Never insert anything metallic into the cabinet openings of the medical LCD monitor. Doing so may create the danger of electric shock.
- To reduce the risk of electric shock, do not remove cover. No user-serviceable parts inside. Only a qualified technician should open the case of the medical LCD monitor.
- 4. Never use your medical LCD monitor if the power cord has been damaged. Do not allow anything to rest on the power cord, and keep the cord away from areas where people can trip over it.
- 5. Be sure to hold the plug, not the cord, when disconnecting the medical LCD monitor power cord from an electric socket.
- 6. Unplug your medical LCD monitor power cord when it is going to be left unused for an extended period of time.
- 7. Unplug your medical LCD monitor power cord from the AC outlet before any service.
- 8. If your medical LCD monitor does not operate normally, in particular, if there are any unusual sounds or smells coming from it, unplug it immediately and contact an authorized dealer or service center.
- 9. Please contact the manufacturer if the set should be installed in an inaccessible area.

Warning

Do not touch input or output connectors and the patient simultaneously.

Warning

This medical LCD monitor is intended for connection to input/output signals and other connectors that comply with relevant IEC standard (e.g.,IEC60950 for IT equipment and IEC60601 series for medical electrical equipment). In addition, all such combination-system shall comply with the standard IEC 60601-1-1, safety requirements for medical electrical systems. Any person who has formed a combination-system is responsible for the system to comply with the requirements of IEC 60601-1-1.

If in doubt, contact qualified technician or your local representative.

On installation

- Openings in the medical LCD monitor cabinet are provided for ventilation. To prevent overheating, these openings should not be blocked or covered. If you put the medical LCD monitor in a bookcase or some other enclosed space, be sure to provide adequate ventilation.
- 2. Put your medical LCD monitor in a location with low humidity and a minimum of dust.
- Do not expose the medical LCD monitor to rain or use it near water (in kitchens, near swimming pools, etc.). If the medical LCD monitor accidentally gets wet, unplug it and contact an authorized dealer immediately. You can clean the medical LCD monitor with a damp cloth if necessary, but be sure to unplug the medical LCD monitor first.
- 4. Place your medical LCD monitor near an easily accessible AC outlet.
- High temperature can cause problems. Don't use your medical LCD monitor in direct sunlight and keep it away from heaters, stoves, fireplaces, and sources of heat.
- 6. Don't place your medical LCD Monitor on an unstable stand, Medical LCD monitor may malfunction or fall.
- 7. This medical LCD monitor should not topple over when tilted at a 5° angle, in any position, during NORMAL USE, excluding transport.
- 8. In the position specified for transport, medical LCD monitor shall not overbalance when tilted at a 10 degree angle.

Environmental Conditions for operation and Storage

- Temperature range within 0°C to 40°C(operation), -20°C to 60°C(storage)
- Relative humidity range 10% to 85% Atmospheric pressure range within 500 to 1060hPa.

Intended Use

- This Medical LCD Monitor is an accessory intended for use with Medical Equipment to display alphabetical, numerical and graphical data.

CAUTION





This symbol alerts the user that important literature concerning the operation of this unit has been included. Therefore, it should be read carefully in order to avoid potential problems.



This symbol warns user that un-insulated voltage within the unit may have sufficient magnitude to cause electrical shock. Therefore, it is dangerous to make contact with any part inside the unit. To reduce the risk of electric shock, DO NOT remove cover (or back). There are no user serviceable parts inside. Refer servicing to qualified service personnel.

To prevent fire or shock hazards, do not expose this unit to rain or moisture. Also, do not use this unit's polarized plug with an extension cord receptacle or other outlets unless the prongs can be fully inserted. The display is designed to meet the medical safety requirements for a patient vicinity device.

This device may not be used in connection with life support equipment.



Underwriters Laboratories (UL) Classification:

UL safety Compliance:

This medical LCD monitor is U.L. Classified WITH RESPECT TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH UL 60601-1/CAN/CSA C22.2 NO. 601.1



EEC Safety Compliance

This medical LCD monitor unit meets the requirements of EN-60601-1 so as to conform to the Medical Device Directive 93/42/EEC (general safety information). Use 120V rating 5-15P type plug only in the U.S

This medical LCD monitor complies to the above standards only when used with the supplied medical grade power supply.

19"(FS-L190*D / FS-L190*DT) JMW190KA1200F02(BRIDGE POWER CORP.) 24"/26"(FS-L240*D / FS-L240*DT / FS-L260*D) JMW1150KA2400F04(BRIDGE POWER CORP.) 32"(FS-L320*D) JMW1180KA2400F01(BRIDGE POWER CORP.)

Caution: Make sure the power cord is the correct type that is required in your area. This medical LCD monitor has a universal power supply that allows operation in either 100-120V AC or 200-240V AC voltage areas (no user adjustment is required).

Use the proper power cord with correct attachment plug type. If the power source is 120 V AC, use a power cord which is a Hospital Grade Power Cord with NEMA 5-15 style plug, labeled for 125 volts AC with UL and C-UL approvals. If the power source is a 240 V AC supply, use the tandem (T blade) type attachment plug with ground conductor power cord that meets the respective European country's safety regulations.

The hospital-grade plug for medical products intended for use in Denmark has DEMKO approval and is rated 13 amps at 250Vac. Plug is recommended for use in medical applications and specifications are being added to the standard SB 107-2-D1. Plug mates with maker's Danish hospital-grade socket. Hospital sockets have slightly different shaped openings allowing only the hospital plug, not the standard Danish plug, to be inserted, to protect the ac circuit in specific medical settings.

Recycling



Follow local governing ordinances and recycling plans regarding the recycling or disposal of this equipment.

Cleaning Instructions



Follow your hospital protocol for the handling of blood and body fluids. Clean the display with a diluted mixture of mild detergent and water. Use a soft towel or swab.

Use of certain detergents may cause degradation to the labels and plastic components of the product.

Consult cleanser manufacturer to see if agent is compatible. Do not allow liquid to enter the display.

Servicing

Do not attempt to service the medical LCD monitor yourself, as opening or removing covers may expose you to dangerous voltages or other hazards, and will void the warranty. Refer all servicing to qualified service personnel.

Unplug the medical LCD monitor from its power source and refer servicing to qualified personnel under the following conditions:

- If the power cord or plug is damaged or frayed.
- If liquid has been spilled into the medical LCD monitor.
- If objects have fallen into the medical LCD monitor.
- If the medical LCD monitor has been exposed to rain or moisture.
- If the medical LCD monitor has been subjected to excessive shock by being dropped.
- If the cabinet has been damaged.
- If the medical LCD monitor seems to be overheated.
- If the medical LCD monitor emits smoke or abnormal odor.
- If the medical LCD monitor fails to operate in accordance with the operating instructions.

Accessories

Use only accessories specified by the manufacturer, or sold with the medical LCD monitor.

Classification

- Protection against electric shock : Class I including AC/DC Adapter
- Applied Parts : No Applied Parts
- Degree of safety in the presence of flammable anesthetics mixture with air or with oxygen or with nitrous oxide.

Not suitable for use in the presence of a flammable anesthetics mixture with oxygen or with nitrous oxide.

- Mode of operation : Continuous.

FCC Information

This medical LCD monitor unit has been tested and found to comply with the limits of a Class B digital device, pursuant to Part 15 of the FCC rules.

These limits are designed to provide reasonable protection against interference. This monitor can radiate radio frequency energy and, if not installed and used in accordance with the instructions, it may interfere with other radio communications equipment. There is no guarantee that interference will not occur in a particular installation.

If this equipment is found to cause harmful interference to radio or television reception, the user is encouraged to try to correct the interference by carrying out one or more of the following measures:

- 1. Reorient or relocate the receiving antenna.
- 2. Increase the distance between the medical LCD monitor and the subject of interference.
- 3. Plug the monitor into an outlet on a different electrical circuit than that to which the subject of interference is connected.
- 4. Consult the dealer or an experienced radio/TV technician for help.

NOTICES TO USER

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC WARNING

This medical LCD monitor generates or uses radio frequency energy. Changes or modifications to this medical LCD monitor may cause harmful interference unless the modifications are expressly approved in the instruction manual. The user could lose authority to operate this equipment if an unauthorized change or modification is made.

1. Guidance and manufacturer's declaration - electromagnetic emissions

The medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the medical LCD monitor should assure that it is used in such an environment.				
Emission test Compliance Electromagnetic environment -guidance				
RF Emissions CISPR 11 Group 1 function. Theref		The medical LCD monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment		
RF Emissions CISPR 11	Class B	The medical LCD monitor is suitable for use in all establishments.		
Harmonic emissions IEC 61000-3-2	D	including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes		
Voltage fluctuations IEC 61000-3-3	Complies			

2. Guidance and manufacturer's declaration - electromagnetic immunity

This medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the medical LCD monitor should assure that it is used in such an environment.					
Immunity test	Immunity test IEC 60601 Test level Compliance level Electromagnetic environment -guidance				
Electrostatic discharge(ESD) IEC61000-4-2	6kV Contact 8kV air	tact 6kV Contact Floors should be wood,concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%			
Electrical fast transient/burst IEC 61000-4-4	2kV for power supply lines 1kV for input/output lines	2kV for power supply lines 1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.		

Surge	1kV differential mode	1kV differential mode	Mains power quality should be that of a
IEC 61000-4-5	2kV common mode	2kV common mode	typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the medical LCD monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance : d $d = [\frac{3.5}{V_{f}}]_{\sqrt{P}}$
			where <i>P</i> is the maximum output power rating of the transmitter in watts (W)

3. Guidance and manufacturer's declaration - electromagnetic immunity

This medical LCD monitor is intended for use in the electromagnetic environment specified below. The customer or the user of monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3.0A/m	3.0A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycle 70 % UT (30 % dip in UT) for 25 cycle <5 % UT (<95 % dip in UT) for 5 sec.	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycle 70 % UT (30 % dip in UT) for 25 cycle <5 % UT (<95 % dip in UT) for 5 sec.	Main power quality should be that of a typical commercial or hospital environment. If the user of monitor requires continued operation during power mains interruptions, it is recommended that monitor be powered from an uninterruptible power supply or a battery. NOTE : UT the A.C. mains voltage prior to application of the test level.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment -guidance
Radiated RF IEC 61000-4-3	3 V/m 80.0 MHz to 2.5 GHz	3 V/m 80.0 MHz to 2.5 GHz	Recommended separation distance 80MHz to 800MHz $d = \left[\frac{3.5}{E_{1}}\right]_{\sqrt{P}}$ 80MHz to 2.5GHz $d = \left[\frac{7}{E_{1}}\right]_{\sqrt{P}}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, It should be less than the compliance level in each frequency range.

4. Recommended separation distances between portable and mobile RF communications equipment and this medical LCD monitor.

- The medical LCD monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the medical LCD monitor as recommended below, according to the maximum output power of the communications equipment.

	Separation dista	nce according to frequenc	y of transmitter[m]
Rated maximum output power of transmitter [W]	$d = \left[\frac{3.5}{V_1}\right] \sqrt{P}$	80MHz to 800MHz $d = \left[\frac{3,5}{E_1}\right] \sqrt{P}$	800MHz to 2.5GHz $d = \left[\frac{7}{E_1}\right] \sqrt{P}$
	V1=3Vrms	E1=3V/m	E1=3V/m

0.01	0.116	0.116	0.2333
0.1	0.368	0.3687	0.7378
1	1.166	1.1660	0.2333
10	3.687	3.6872	0.7375
100	11.660	11.6600	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m)can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1) At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



19" FS-L190*D / FS-L190*DT Monitor



> Accessories



D-SUB Cable

User Manual



DVI Cable



BNC Cable



AC-DC Adaptor (JMW190KA1200F02)



AC Power cord (Hospital Grade)



S-Video (Y/C) Cable (Option)



SCREW FH M3X6

SCREW BH M4X10



24" FS-L240*D / FS-L240*DT Monitor







D-SUB Cable



AC-DC Adaptor (JMW1150KA2400F04)



(Hospital Grade)

AC Power cord



S-Video (Y/C) Cable (Option)

圄

DC Cable Terminal Male / female (Option)

User Manual





DVI Cable



SCREW FH M3X6

SCREW BH M4X10





Accessories



User Manual



D-SUB Cable

BNC Cable



AC-DC Adaptor (JMW1150KA2400F04)



AC Power cord

(Hospital Grade)

S-Video (Y/C) Cable (Option)

DC Cable Terminal Male / female (Option)



DVI Cable



SCREW FH M3X6 SCREW BH M4X10



32" FS-L320*D Monitor



User Manual

DVI Cable

員



D-SUB Cable

BNC Cable





AC-DC Adaptor (JMW1180KA2400F01)



AC Power cord (Hospital Grade)



DC Cable Terminal Male / female (Option)





Ì SCREW FH M3X6 SCREW BH M4X10





> 42" FS-L420*D Monitor



Accessories









User Manual

D-SUB Cable

S-Video (Y/C) Cable (Option)

DVI Cable



BNC Cable



AC Power cord (Hospital Grade)

8



Remote controller E

Battery(AAA)



Connector

SFS-L190*D / FS-L240*D / FS-L260*D / FS-L320*D Input connector





FS-L190*DT/FS-L240*DT input connector



Dual DVI input







Mechanical Product Drawing

3 19" FS-L190*D / FS-L190*DT Dimension





Rear view (Without I/O Cover)

Rear view (Installed I/O Cover)



Unit : mm





24" FS-L240*D / FS-L240*DT Dimension

Rear view (Without I/O Cover)

Rear view (Installed I/O Cover)







Unit : mm





Unit : mm

English - 25

26" FS-L260*D Dimension



Front view



Rear view (Without I/O Cover)



Rear view (Installed I/O Cover)



Unit : mm





Rear view





Unit : mm

Control

SOSD Button



An 8 button keypad, located in bottom right corner on the front of the display, allows the user to make adjustments to various display parameters using the On Screen Display (OSD) system.

Power Indicating LED

Normal mode (ON): Green			
Standby mode : Blinking Green			
Off mode	: Monitor Off		

- Note 1 : LED normal and off mode sign can be changed according to the requirement of the customer whether normal mode green is on or off
- Note 2 : The main AC power switch, on the back panel, should be in the ON position. The DC power button is used to turn on the monitor

On-Screen Display (OSD) Function Button

- 1. POWER : Soft power turns the monitor ON or OFF
- 2. MENU : Used to activate to the OSD menu and exit from main menu or sub menu.
- 3. PIP : Enables PIP(picture in Picture) function. Selects PIP,PBP1,PBP2.
- 4. UP () : With the OSD deactivated, it is a Hot Key for increasing brightness. With the OSD activated, moves the cursor upward.
- 5. DOWN (): With theOSD deactivated, it is a Hot Key for decreasing brightness. With the OSD activated, moves the cursor downward.

6. PLUS (+)	: With the OSD deactivated, it is a Hot key for increasing contrast.
	With the OSD activated, enter sub menu and increases the
	adjustment of the selected function.
7. MINUS (-)	: With the OSD deactivated, it is a Hot key for increasing contrast.
	With theOSD activated, it decreases the adjustment of the selected
	function.
8. INPUT	: With the OSD deactivated, if pressed down for over 1 sec., it is a
	Hot Key for auto-adjustment control under DSUB ANALOG / RGBs
	signals.
	With theOSD activated, it changes the displayed signal source.

Remote button function

<Notice> Remote control is available in FS-L320*D,FS-L420*D models only.





- 1. SOURCE : Changes the display signal source.
- 2. POWER : Soft power turns the monitor ON or OFF.
- 3. UP (): With the OSD deactivated, it is a Hot key for increasing brightness. With the OSD activated, moves the cursor upward.
- 4. DOWN (): With the OSD deactivated, it is a Hot key for decreasing brightness. With the OSD activated, moves the cursor downward.
- 5. MINUS (-) : With the OSD deactivated, it is a Hot key for decreasing contrast. With the OSD activated, it decreases the adjustment of the selected function.
- 6. PLUS (+) : With the OSD deactivated, it is a Hot key for increasing contrast. With the OSD activated, it decreases the adjustment of the adjustment of the Selected function.
- 7. MENU : With the OSD deactivated, activates the OSD menu. With the OSD activated, exits from main menu or sub menu.
- 8. AUTO : Fits to the most appropriate screen on the D-SUB Analog signal.
- 9. PIP : Enables PIP(picture in Picture) function. Selects PIP,PBP1,PBP2
- 10. MUTE : Sound muted.
- 11. SWAP : Swaps the position of the Primary and Secondary images.

GPIO

There are four pins on the RJ69 GPIO connector. Each pin has a preprogrammed function assigned to it. The function is initiated when the pin is grounded.

Pin 1. Primary and Secondary Swap. Grounding this pin will swap the primary and secondary image. Pin 3. Record Indicator The record indicator is displayed in the top left corner when the pin is ground to pin 4. The indicator will vanish when the contact is opened.



Pin 2. PIP,PBP1,PBP2 Single Continuously grounding this pin causes the position and size choices to cycle.



- REC
- Pin 4. Connector Ground This is the common ground location.



Power management

This monitor does not adhere to the VESA DPMS standard when no signal is present on the video input.

Model	Status	LED sign	Power Consumption
FS-L190*D	Normal mode	Green on	<60W
	Standby mode	Green blinking	<20W
FS-L190*DT	Normal mode	Green on	<60W
	Standby mode	Green blinking	<20W
FS-L240*D	Normal mode	Green on	<100W
	Standby mode	Green blinking	<20W
FS-L240*DT	Normal mode	Green on	<100W
	Standby mode	Green blinking	<20W
FS-L260*D	Normal mode	Green on	<130W
	Standby mode	Green blinking	<20W
FS-L320*D	Normal mode	Green on	<150W
	Standby mode	Green blinking	<20W
FS-L420*D	Normal mode	Green on	<260W
	Standby mode	Green blinking	<20W

OSD

DSUB ANALOG / RGBS input source

ADJUST

	ADJUST			
	BRIGHTNESS			50
t 🗌	CONTRAST			50
	CLOCK			50
	PHASE			65
2	BACK LIGHT			100
	AUTO ADJUST			
		ELECT	🔶 MOVE	
	•		•	
DSUB	ANALOG	1	280X1024/80	(x75Hz

- 1. BRIGHTNESS
- Increases or decreases the brightness. (Range : 0~100) 2. CONTRAST
- Increases or decreases the contrast. (Range : 0~100) 3. CLOCK
- Increases or decreases the sampling frequency. (Range : 0~100) 4. PHASE
- Increases or decreases the Phase level. (Range : 0~100) 5. BACKLIGHT
- Adjusts backlight dimming level. (Range : 0~100)
- 6. AUTO ADJUST

Fits to the most appropriate screen on the D-SUB Analog / RGBs signal.

COLOR TEMP

COLO	COLOR TEMP				
MODI	C C	1			
tr 🕝					
5					
M EXIT	+ SELECT	♣ MOVE			
DSUB ANALO	G	1280X1024/80K x75Hz			

- 1. MODE
- Changes the color mode (C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)
- 2. RED
- Red balance. (Only works with USER Mode) (Range : 0~100) 3. GREEN
- Green balance. (Only works with USER Mode) (Range : 0~100) 4. BLUE
- Blue balance.(Only works with USER Mode) (Range : 0~100)

IMAGE

	IMAGE										
	IMAGE SIZE		FI	LL							
	H POSITION										50
	V POSITION							_	_		50
CT 🔜	GRMMA		2.	0							
5	FILTER		N	DRI	MAL	L					
	OVER SCAN		0	1	2	3	4	5	6	7	8
	IMAGE SETT	ING				١	101	DE	1		
M EXIT		+ SELEC	T				¢	M	0 V I		
DSUB	ANALOG				128	30 X	102	4/8	10 K	x7	5Hz

- 1. IMAGE SIZE
- Changes the image size. (Full, Fill aspect, 1:1, Normal)
- 2. H POSITION
- Adjusts the horizontal position of the displayed source image. (Range : 0~100) 3. V POSITION
- Adjusts the vertical position of the displayed source image. (Range : 0~100) 4. GAMMA
- Adjusts GAMMA value (VIDEO,BYPASS,1.8,2.0,2.2,2.4,2.6,PACS) 5. FILTER
- Sets the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest) 6. OVER SCAN
 - Adjusts the displayed size. (0~8)
- 7. IMAGE SETTING
- Changes the image setting.(Preset 1,2 / User 1,2,3)
- 8. ZOOM / PAN
- Enlarges the image, moves images left and right.
- 9. FREEZE FRAME

Keeps the image still.

SETUP

5	SETUP					
	LANGUAGE	ENGLISH				
	OSD COLOR	50				
	OSD POSITION	0				
	DURATION	120 SEC				
C D	RESET SETTINGS					
	AUTO SOURCE SELECT ON					
	USER NAME ENTRY					
M EXIT	.+ SELE	CT 🔷 MOVE				
DSUB	ANALOG	1280X1024/80K x75Hz				

- 1. LANGUAGE
- Changes the OSD language (8 language)
- 2. OSD COLOR
- Adjusts the OSD background from white opaque to translucent. 3 .OSD POSITION
 - Changes the osd position. (9 Positions)
- 4. DURATION
- Adjusts the length of time the OSD menu is present on the screen. (5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
- 5. RESET SETTING
- Changes all the OSD values make to factory default. 6. AUTO SOURCE SELECT
 - Disables or enable auto source select.
 - ON: Searches through all possible input source untill an active video
 - source is found.
 - OFF: Video input is manually selected.
- 7. USER NAME ENTRY
- Creates or changes the user name when powered on.

PIP

	PIP	S	INGLE
M EXIT		+ SELECT	🗢 MOVE
DSUB	ANALOG		1280X1024/80K x75Hz

- 1. LAYOUT
- Changes the OSD layout. (Single, PIP, PBP1, PBP2)
- 2. SOURCE
- Changes the secondary source.
- 3. SIZE
 - Changes the PIP size(Small, Large).
- 4. POSITION
- Changes the secondary source.
- 5. SWAP

Swaps the position and size of the Primary and Secondary image.
DVI OPTICAL / DVI DIGITAL input source

ADJUST

	ADJUST	
	BRIGHTNESS	50
17	CONTRAST	50
\bigcirc	BACK LIGHT	100
5		
		-
M EXIT	+ SELEC	T 🌩 MOVE
DVI DI	GITAL	1280X1024/80K x75Hz

1. BRIGHTNESS

- Increases or decrease the brightness. (Range : 0~100) 2. CONTRAST
- Increases or decreases the Contrast. (Range : 0~100) 3. BACKLIGHT
- Adjusts backlight dimming level. (Range : 0~100)

COLOR TEMP

COLOR	TEMP	
MODE	C1	
50		
S		
M EXIT	+ SELECT	♦ MOVE
DVI DIGITAL		1280X1024/80K x75Hz

1. MODE

- Changes the color temperature mode. (C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)
- 2. RED
- Red balance. (Only works with USER Mode) (Range : $0 \sim 100$) 3. GREEN
- Green balance. (Only works with USER Mode) (Range : 0~100) 4. BLUE
- Blue balance. (Only works with USER Mode) (Range : 0~100)

IMAGE



- 1. IMAGE SIZE
- Changes the image size. (Full, Fill aspect**, 1:1**, Norma) **Only in DVI Optical
- 2. GAMMA
- Adjusts GAMMA value. (VIDEO,BYPASS,1.8,2.0,2.2,2.4,2.6,PACS) 3. FILTER
- Sets the sharpness of image. (Softest, Soft, Normal, Sharp, Sharpest) 4. OVER SCAN
- Adjusts the displayed size. (0~8)
- 5. IMAGE SETTING
- Changes the image setting. (Preset 1,2 / User 1,2,3) 6. ZOOM/PAN
- Enlarges the image, moves images left and right.
- 7. FREEZE FRAME
- Keeps the image still.

SETUP



- 1. LANGUAGE
- Changes the OSD language (8 language)
- 2. OSD COLOR
 - Adjusts the OSD background from white opaque to translucent.
- 3 .OSD POSITION Changes the OSD position. (9 Positions)
- 4. DURATION
- Adjusts the length of time the OSD menu is present on the screen. (5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
- 5. RESET SETTING
- Changes all the OSD values to factory outgoing state. 6. AUTO SOURCE SELECT
- Disables or enables auto source select.
- ON: Searches through all possible input source until an active video source is found.
- OFF: Video input is manually selected.
- 7. USER NAME ENTRY
 - Creates or changes the user name when powered on.

PIP

	PIP			2
	LAYOUT		SINGLE	
\bigcirc				
4				
S				
Ú 🗆				
M EXIT		+ SELECT		MOVE
DSUB /	ANALOG		1280X1	024/80K x75Hz

- 1. LAYOUT
- Changes the OSD layout. (Single, PIP, PBP1, PBP2)
- 2. SOURCE
- Changes the secondary source.
- 3. SIZE Changes the PIP size(Small, Large).
- 4. POSITION
- Changes the secondary source.
- 5. SWAP

Swaps the position and size of the Primary and Secondary image.

YPbPr input source

ADJUST

	ADJUST		
	BRIGHTNESS		50
¢	CONTRAST		50
	SHARPNESS		50
	SATURATION		65
S	COLOR		_ 0
	BACKLIGHT		100
	PHASE		65
M EXIT	+ SELECT	🔷 MOVE	
YPbPr		1920X1080/80	IK x60Hz

1. BRIGHTNESS

- Increases or decreases the brightness. (Range : 0~100)
- 2. CONTRAST
- Increases or decreases the Contrast. (Range : 0~100)
- 3. SHARPNESS
- Adjusts the sharpness of video image. (Range : 0~100)
- 4. SATURATION
- Changes the tone of color. (Range : 0~100)
- 5. COLOR
- Changes the richness of color. (Range : Greenish 0~50, Redish 0~50) 6. BACKLIGHT
- Adjusts backlight dimming level. (Range : 0~100)
- 7. CLOCK

USER(7200K)

1. MODE

2. RED

3. GREEN

- Increases or decreases the sampling. (Range : 0~100) 8. PHASE
- Increases or decreases the Phase level. (Range : 0~100)

COLOR TEMP



Green balance. (Only works with USER Mode) (Range : 0~100) 4. BLUE

Red balance. (Only works with USER Mode) (Range : 0~100)

Changes the color temperature mode. (C1(Reddish,6500K), C2(Bluish,9300K),

Blue balance.(Only works with USER Mode) (Range : 0~100)

IMAGE

	IMAGE									
	IMAGE SIZE	FI	LL							
	V POSITON						_	_		50
3	H POSITON									50
C A	GRMMA	2.	0							
50	FILTER	N	DRI	MA	L					
	OVER SCAN	0	1	2	3	4	5	6	7	8
	IMAGE SETTING				ł	PRE	SE	T1		
M EXIT	+ SELECT					\$	M	0 V I		
YPbPr				192	20 X	108	10/8	10 K	X 6	OHz

- 1. IMAGE SIZE
- Changes the image size.(Full, Fill aspect, 1:1, Normal, Anamorphic) 2. H POSITION
- Adjusts the horizontal position of the displayed source image.(Range : 0~100) 3. V POSITION
- Adjusts the vertical position of the displayed source image. (Range : 0~100) 4. GAMMA
- Adjusts GAMMA value.(VIDEO, BYPASS, 1.8, 2.0, 2.2, 2.4, 2.6, PACS))
- 5. FILTER
- Sets the sharpness of image.(Softest, Soft, Normal, Sharp, Sharpest) 6. OVER SCAN
- Adjusts the displayed size. (0~8)
- 7. IMAGE SETTING
- Changes the image setting.(Preset 1,2 / User 1,2,3)
- 8. ZOOM / PAN
- Enlarges the image, moves images left and right.
- 9. FREEZE FRAME
- Keeps the image still.

SETUP

5	SETUP	
	LANGUAGE	ENGLISH
	OSD COLOR	50
	OSD POSITION	0
	DURATION	120 SEC
17 E	RESET SETTINGS	
	AUTO SOURCE SELECT	r on
	USER NAME ENTRY	
M EXIT	+ SELECT	♠ MOVE
YPbPr		1920X1080/80K x60Hz

- 1. LANGUAGE
- Changes the OSD language (8 language)
- 2. OSD COLOR
 - Adjusts the OSD background from white opaque to translucent.
- 3 .OSD POSITION
 - Changes the osd position. (9 Positions)
- 4. DURATION
- Adjusts the length of time the OSD menu is present on the screen. (5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
- 5. RESET SETTING
- Changes all the OSD values to factory outgoing state. 6. AUTO SOURCE SELECT
- Disables or enables auto source select.
- ON: Searches through all possible input source until an active video source is found.
- OFF: Video input is manually selected.
- 7. USER NAME ENTRY
 - Creates or changes the user name when powered on.

PIP

	PIP	
	LAYOUT	SINGLE
5		
ċ I		
M EXIT	+ SELECT	♠ MOVE
YPbPr		1920X1080/80K x60Hz

- 1. LAYOUT
- Changes the OSD layout. (Single, PIP, PBP1, PBP2)
- 2. SOURCE
 - Changes the secondary source.
- 3. SIZE
 - Changes the PIP size(Small, Large).
- 4. POSITION Changes the secondary source.
- 5. SWAP
 - Swaps the position and size of the Primary and Secondary image.

SVIDEO / CVIDEO input source

ADJUST

	ADJUST			
	BRIGHTNESS			50
<u>i</u>	CONTRAST			50
\bigcirc	SHARPNESS			50
	SATURATION			65
5	COLOR			0
	BACKLIGHT			100
M EXIT	+ s	ELECT	🔶 MOVE	
CVIDE	0		NTSC 60Hz / 4.	43MHz

1. BRIGHTNESS

- Increases or decreases the brightness. (Range : 0~100)
- 2. CONTRAST
- Increases or decreases the Contrast. (Range : 0~100)
- 3. SHARPNESS
- Adjusts the sharpness of video image. (Range : 0~100) 4. SATURATION
- Changes the tone of color. (Range : 0~100)
- 5. COLOR
- Changes the richness of color. (Range : Greenish 0~50, Redish 0~50) 6. BACKLIGHT
- Adjusts backlight dimming level. (Range : 0~100)

COLOR TEMP

	COLOR TEMP		
	MODE	C1	
17 📀			
S			
M EXIT	+	SELECT	♠ MOVE
CVIDE	0		NTSC 60Hz / 4.43MHz

1. MODE

Changes the color temperature mode. (C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)

- 2. RED
- Red balance.(Only works with USER Mode) (Range : 0~100) 3. GREEN
- Green balance.(Only works with USER Mode) (Range : 0~100) 4 BLUE
- Blue balance.(Only works with USER Mode) (Range : 0~100)

IMAGE

	IMAGE	
	IMAGE SIZE	FILL
	V POSITON	50
\bigcirc	H POSITON	50
	GRMMA	2.0
5	FILTER	NORMAL
	OVER SCAN	0 1 2 3 4 5 6 7 8
	IMAGE SETTING	PRESET1
M EXIT	+ SELECT	T 🔶 MOVE
CVIDE	0	NTSC 60Hz / 4.43MHz

- 1.IMAGE SIZE
- Changes the image size.(Full, Fill aspect, 1:1, Normal, Anamorphic)
- 2. H POSITION
- Adjusts the horizontal position of the displayed source image. (Range : 0~100) 3. V POSITION
- Adjusts the vertical position of the displayed source image. (Range : 0~100) 4. GAMMA
- Adjusts GAMMA value (VIDEO,BYPASS,1.8,2.0,2.2,2.4,2.6,PACS) 5. FILTER
- Sets the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest) 6. OVER SCAN
 - Adjusts the displayed size. (0~8)
- 7. IMAGE SETTING
- Changes the image setting.(Preset 1,2 / User 1,2,3)
- 8. ZOOM / PAN
- Enlarges the image, moves images left and right.
- 9. FREEZE FRAME
- Keeps the image still.

SETUP

<u> </u>	SETUP	
	LANGUAGE	ENGLISH
	OSD COLOR	50
\bigcirc	OSD POSITION	0
	DURATION	120 SEC
1 D	RESET SETTINGS	
	AUTO SOURCE SELEC	T ON
	USER NAME ENTRY	
M EXIT	+ SELEC	T 🔶 MOVE
CVIDE	0	NTSC 60Hz / 4.43MHz

- 1. LANGUAGE
- Changes the OSD language (8 language)
- 2. OSD COLOR
- Adjusts the OSD background from white opaque to translucent.
- 3 .OSD POSITION
- Changes the osd position. (9 Positions) 4. DURATION
- Adjusts the length of time the OSD menu is present on the screen. (5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
- 5. RESET SETTING
 - Changes all the OSD values to factory outgoing state.
- 6. AUTO SOURCE SELECT
 - Disables or enables auto source select.
 - ON: Searches through all possible input source until an active video source is found.
 - OFF: Video input is manually selected.
- 7. USER NAME ENTRY
- Creates or changes the user name when powered on.

PIP



- 1. LAYOUT
- Changes the OSD layout. (Single, PIP, PBP1, PBP2) 2. SOURCE
- Changes the secondary source.
- 3. SIZE
- Changes the PIP size(Small, Large).
- 4. POSITION
- Changes the secondary source.
- 5. SWAP
 - Swaps the position and size of the Primary and Secondary image.

40 - English

SDI input source

ADJUST

	ADJUST			
	BRIGHTNESS			50
17 🛄	CONTRAST			50
3	SHARPNESS			50
	SATURATION			65
5	COLOR			0
	BACKLIGHT			100
M EXIT	+ SELEC	T	🔶 MOVE	
SDI			1920X1080i/	30Hz

1. BRIGHTNESS

- Increases or decreases the brightness. (Range : 0~100)
- 2. CONTRAST
- Increases or decreases the Contrast. (Range : 0~100) 3. SHARPNESS
- Adjusts the sharpness of video image. (Range : 0~100) 4. SATURATION
- Changes the tone of color. (Range : 0~100)
- 5. COLOR
- Changes the richness of color. (Range : Greenish 0~50,Redish 0~50) 6. BACKLIGHT
- Adjusts backlight dimming level. (Range : 0~100)

COLOR TEMP



1. MODE

Changes the color temperature mode. (C1(Reddish,6500K), C2(Bluish,9300K), USER(7200K)

- 2. RED
- Red balance. (Only works with USER Mode) (Range : 0~100) 3. GREEN
- Green balance. (Only works with USER Mode) (Range : 0~100) 4. BLUE
- Blue balance.(Only works with USER Mode) (Range : 0~100)

IMAGE

	IMAGE	
	IMAGE SIZE	FILL
	GRMMA	2.0
8	FILTER	NORMAL
to 🔍	OVER SCREN	0 1 2 3 4 5 6 7 8
2	IMAGE SETTING	MODE1
	ZOOM/PAN	
	FREEZE FRAME	0 N
M EXIT	+ SELEC	T 🔷 MOVE
SDI		1920X1080i / 30Hz

- 1. IMAGE SIZE
- Changes the image size.(Full, Fill aspect, 1:1, Normal, Anamorphic) 2. GAMMA
- Adjusts GAMMA value (VIDEO,BYPASS,1.8,2.0,2.2,2.4,2.6,PACS) 3. FILTER
- Sets the sharpness of image (Softest, Soft, Normal, Sharp, Sharpest) 4. OVER SCAN
 - Adjusts the displayed size. (0~8)
- 5. IMAGE SETTING
- Changes the image setting.(Preset 1,2 / User 1,2,3)
- 6. ZOOM / PAN
- Enlarges the image, moves images left and right.
- 7. FREEZE FRAME
- Keeps the image still.

SETUP

19m	SETUP	
	LANGUAGE	ENGLISH
	OSD COLOR	50
	OSD POSITION	0
	DURATION	120 SEC
1 D	RESET SETTINGS	
	AUTO SOURCE SELECT	ON
	USER NAME ENTRY	
M EXIT	+ SELECT	🗢 MOVE
SDI		1920X1080i / 30Hz

- 1. LANGUAGE
- Changes the OSD language (8 language)
- 2. OSD COLOR
- Adjusts the OSD background from white opaque to translucent.
- 3 .OSD POSITION
- Changes the osd position. (9 Positions)
- 4. DURATION
- Adjusts the length of time the OSD menu is present on the screen. (5, 10, 20, 30, 60, 90, 120, 180, 240 seconds)
- 5. RESET SETTING
 - Changes all the OSD values to factory outgoing state.
- 6. AUTO SOURCE SELECT
- Disables or enables auto source select.
- ON: Searches through all possible input source until an active video source is found.
- OFF: Video input is manually selected.
- 7. USER NAME ENTRY
- Creates or changes the user name when powered on.

PIP



- 1. LAYOUT
- Changes the OSD layout. (Single, PIP, PBP1, PBP2) 2. SOURCE
- Changes the secondary source.
- 3. SIZE
- Changes the PIP size(Small, Large).
- 4. POSITION
- Changes the secondary source.
- 5. SWAP
 - Swaps the position and size of the Primary and Secondary image.

OSD System overview

Menus	Function Descriptions		
BRIGHTNESS	Press the BRIGHTNESS botton to display the ADJUST menu or UP,DOWN hot key. Setting the brightness too high or too low will decrease the amount of visible grayscales.		
CONTRAST	Press the CONTRAST button to display the ADJUST menu or +,- hot key.		
	Setting the Contrast too high or too low will cause loss of some grascales.		
CLOCK	Do not adjust. It will adjust automatically after auto adjustment. When frequency value is wrong, the horizontal image will have a wrong size or noise.		
PHASE	Do not adjust . It will adjust automatically after auto adjustment. When frequency value is wrong, the image will have a noise.		
BACKLIGHT	Adjusts backlight dimming level. Setting the backlight too low will cause dark image and too high will decrease the backlight life time.		
AUTO ADJUST	Fits to the most appropriate screen on the D-SUB Analog/RGBS signal		
SHARPNESS	Adjusts the sharpness of video image		
SATURATION	Changes the tone of color		
COLOR	Changes the richness of color(Range Greenish 0~50,Redish 0~50)		
COLOR TEMP C1	Default 6500K color setting		
COLOR TEMP C2	Default 9300K color setting		
COLOR TEMP USER	Default 7200K color setting,but it's changeable value by user		
IMAGE SIZE	DSUB/ DVI OPTICAL / DVI DIGITAL input source FULL FILL ASPECT		
	CAUTION : FILL ASPECT, NORMAL		
	Size depends on input size ratio		
	YPbPr / RGBS / SDI / CVIDEO / SVIDEO input source FULL FILL ASPECT 1:1		
	NORMAL ANAMORPHIC		

English - 43

Menus	Function Descriptions
H POSITION	Adjusts the Horizontal position of the image. It will return to the default state when executing AUTO ADJUST or RESET SETTINGS.
V POSITION	Adjusts the Verticall position of the image. It will return to the default state when executing AUTO ADJUST or RESET SETTINGS.
GAMMA	Adjusts the gamma curve of video image.
	1.8 2.0 2.2 2.4 2.6 PACS VIDEO Image: Constraint of the panel specification.
FILTER	Adjusts the gamma curve of video image. ABCDEFGHIJKLMN0P0RSTUVW
OVER SCAN	Enables a 10% over scan of orignal input image.
	O: 1: 2: ABCDEFGHIJKLIMNOPORSTU ABCDEFGHIJK
	3: 5: ICDEFGHIJKLMNOPQRS: CDEFGHIJKLMNOPQR CDEFGHIJKLMNOPQR ICDEFGHIJKLMNOPQRS: CDEFGHIJKLMNOPQR
	6: 7: 2: 2: 2: 2: 2: 2: 2: 2: 2: 2





Menus			Fund	tion Des	criptions				
OSD DURATION	Adjusts the time the OSD menu remains on the screen.								
RESET SETTING	Resets the mor	Resets the monitor back to factory settings							
AUTO SOURCE SELECT	When AUTO So input Source.	OURCE SE	LECT is on	, the monit	tor will auto	omatically	search for	an active	MAIN
USER NAME ENTRY	Enables the us	ser to enter	informatio	n desired	to be disp	played du	ring powe	r on.	
PIP LAY OUT	Selets on of fo Note : The as	•				,			
PIP SOURCE		Selects a PIP source input. You can change other sub windows through OSD pip menu. Below chart is PIP							
					Sub wir				
	Input source	DVI OPTICAL DVI DIGITAL 2	DVI DIGITAL	DSUB ANALOG	SDI	YPbPr	RGBS	CVIDEO	SVIDEO
	DVI OPTICAL DVI DIGITAL 2	Х	0	0	0	0	0	0	0
	DVI DIGITAL	0	Х	Х	0	Х	Х	0	0
	DSUB ANALOG	0	X	Х	O ¹	Х	Х	0	0
	SDI	0	0	O ¹	Х	O ¹	O ¹	Х	Х
	YPbPr	0	Х	Х	O ¹	Х	Х	0	0
	RGBS	0	X	Х	O ¹	Х	Х	0	0
	CVIDEO	0	0	0	Х	0	0	Х	Х
	SVIDEO	0	0	0	Х	0	0	Х	Х
	O ¹ : Support u	p to UXGA	,60Hz (162	2MHz) : O	ption to h	ave DVI 2			
PIP SIZE	Changes PIP m	node sub wi	ndow size						
	SMALL		LARG		1				
		Sub		Sub					
	Main		Main						
	25% / Panel	size	33% / Pa	nel size					
PIP POSITION	Changes PIP m	node sub wi	ndow Posit	on					
	$\begin{array}{c c} & & & \\ \hline \\ \hline$								
	Î	Main (L	Main Sub		Main) (> _ M	ain		



Standard Signal table

PC Supported Mode

Resolution	Horizontal Frequency (KHz)	Vertical Frequency (Hz)	Clock Frequency (MHz)				
640 X 350 @70Hz	31.50	70.00	25.17		▲	. ♠	
720 X 400 @70Hz	31.50	70.00	28.32				
640 X 480 @60Hz	31.50	60.00	25.18				
640 X 480 @75Hz	37.50	75.00	31.50				
800 X 600 @60Hz	37.90	60.00	40.00	F			
800 X 600 @75Hz	46.90	75.00	49.50	FS-L190*D(T	0*D	Ģ	
1024 X 768 @60Hz	48.40	60.00	65.00	-L19	/ FS-L420*D	260*	
1024 X 768 @75Hz	60.00	75.00	78.75	Ч. Ч.	/ FS	I-S-L	
1152 X 864 @60Hz	54.35	60.00	80.000		0*D		
1152 X 864 @75Hz	67.50	75.00	108.00		S-L320*D	0*D(
1280 X 720 @60Hz	44.70	60.00	74.40		LS LS	FS-L240*D(T) / FS-L260*D	
1280 X 1024@60Hz	64.00	60.00	108.00			S. L	
1280 X 1024@75Hz	80.00	75.00	135.00	↓			
1360 X 768@75Hz	47.70	60.00	108.75				
1600 X 1200@60Hz	75.00	60.00	162.00				
1920 X 1080@60Hz	67.50	60.00	148.50		V		
1920 X 1200@60Hz	74.00	60.00	154.12			↓	

SDI Video format

Output Signal	Description
SMPTE-274M	1080i (60 / 59.94 / 50) 1080p (30 /29.97 / 25 / 24 / 24sF / 23.98 / 23.98sF)
SMPTE-296M	720p (60 / 59.94 / 50)
SMPTE-260M	1035i (60 / 59.94)
SMPTE-125M	480i (59.94)
ITU-R BT.656	576i (50)

Signal connector Pin Assignments

VGA (15Pin D-SUB)

Pin No.	Assignment	Pin No.	Assignment
1	Red	9	No Connection
2	Green	10	Ground-Sync
3	Blue	11	Ground
4	Ground	12	DDC Data
5	DDC 5V Standby	13	H.Sync
	Cable Connection check	14	V.Sync
6	Ground-Red	15	DDC Clock
7	Ground-Green		
8	Ground-Blue		



DVI In,Out (24Pin DVI-D)

Pin No.	Assignment	Pin No.	Assignment
1	T.M.D.S. Data2-	13	No Connection
2	T.M.D.S. Data2+	14	+5V Power
3	T.M.D.S. Data2 Shield	15	Ground
4	No Connection	16	Hot Plug Detect
5	No Connection	17	T.M.D.S. Data0-
6	DDC Clock	18	T.M.D.S. Data0+
7	DDC Data	19	T.M.D.S. Data0 Shield
8	No Connection	20	No Connection
9	T.M.D.S. Data1-	21	No Connection
10	T.M.D.S. Data1+	22	T.M.D.S. Clock Shield
11	T.M.D.S. Data1 Shield	23	T.M.D.S. Clock+
12	No Connection	24	T.M.D.S. Clock-



C-Video (BNC)

Pin No.	Assignment
1	Composite
2	Ground



S-Video (BNC)

Pin No.	Assignment	
1	S-VIDEO/Y (Luma)	
2	S-VIDEO/C (Chroma)	
3	Ground	



RS232C (D-SUB 9Pin)

Pin No.	Assignment	
1	No Connection	
2	TXD	
3	RXD	
4	No Connection	
5	Ground	
6	No Connection	
7	No Connection	
8	No Connection	
9	No Connection	



SDI (BNC)

Pin No.	Assignment
1	SDI IN
2	SDI OUT
3	Ground



RGBHV/RGBS/YPbPr (BNC)

Pin No.	Assignment		
PIN NO.	RGBS	Y Pb Pr	
1	Red	Pr	
2	Green	Y	
3	Blue	Pb	
4	H-Sync / C-Sync	No Connection	
5	V-Sync	No Connection	
6	Ground		



OPTICAL

Pin No.	Assignment
1	OPICAL Clock
2	OPICAL Blue
3	OPICAL Green
4	OPICAL Red

1	2	3	4
0	 	G	R

GPIO

Pin No.	Assignment		
1	P,S Swap		
2	PIP,PBP1,PBP2 Select		
3	Record Indicator		
4	Ground		



Specification

SFS-L190*D / FS-L190*DT

Model		FS-L190*D / FS-L190*DT	
	Туре	19" TFT-LCD	
	Screen Size	19 inch	
	Maximum Resolution	1280 X 1024 @ 60Hz	
Optical	Pixel Pitch	0.294(H) mm X 0.294(V) mm	
Characteristics	Display Colors	16.7M	
	Contrast Ratio(Typical)	500:1	
	Viewing Angle	85° / 85° / 85° / 85°	
	Response Time	14 msec(Gray to Gray)	
	Luminance(Typical)	210cd/m ²	
Touch Specification	Touch Panel	ELO 5wire resistive touch screen	
(F-L190*DT only)	Interface	USB & SERIAL	
Front Filter (FS-L190*D only)	Acrylic	Double side Anti-Reflection coating	
Resolution	Signal input	VGA~SXGA	
Power	Maximum	Max 60W	
Fower	Standby Mode	Max 20W	
Control Key	Front Side	INPUT, -, +, , , PIP, MENU, POWER	
Input Signal	Video	1xDVI, 1x Optical DVI 1(DVI 2 option), 1x D-SUB, 1xBNC (CVBS)Video, 2 x BNC (SVHS Y/C), 1 x BNC(SDI), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)	
Output Signal	Video	1xDVI, 1xBNC (SDI)	
Input Power	DC 12V, 7A Max		
Dimension	Size and Weight	423(W)X351.5(H)X76.5(D) (mm) 7.3Kg -Without stand 16.653(W)X13.838(H)X3.011(D) (inch) 16.09bs-Without stand	

FS-L240*D / FS-L240*DT

Model		FS-L240*D / FS-L240*DT
-	Туре	24" TFT-LCD
	Screen Size	24 inch
	Maximum Resolution	1920 X 1200 @ 60Hz
Optical	Pixel Pitch	0.270(H) mm X 0.270(V) mm
Characteristics	Display Colors	16.7M
	Contrast Ratio(Typical)	700:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	320cd/m
Touch Specification	Touch Panel	ELO 5wire resistive touch screen
(F-L240*DT only)	Interface	USB & SERIAL
Front Filter (FS-L240*D only)	Acrylic	Double side Anti-Reflection coating
Resolution	Signal input	VGA~WUXGA
Davia	Maximum	Max 100W
Power	Standby Mode	Max 20W
Control Key	Front Side	INPUT, -,+, , , PIP, MENU, POWER
Input Signal	Video	1xDVI, 1x Optical DVI 1(DVI 2 option), 1x D-SUB, 1xBNC (CVBS)Video, 2 x BNC (SVHS Y/C), 1 x BNC(SDI), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1xDVI, 1xBNC (SDI)
Input Power	DC 24V, 6.25A Max	
Dimension	Size and Weight	580(W)X386(H)X95(D) (mm) 7.5Kg -Without stand 22.834(W)X15.196(H)X3.740(D) (inch) 16.53 lbs-Without stand

SFS-L260*D

Model		FS-L260*D
	Туре	26" TFT-LCD
	Screen Size	26 inch
	Maximum Resolution	1920X 1200 @ 60Hz
Optical	Pixel Pitch	0.2685(H) mm X 0.2685(V) mm
Characteristics	Display Colors	16.7M
	Contrast Ratio(Typical)	700:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	320cd/m ²
Resolution	Video input	VGA~WUXGA
Front Filter	Acrylic	Double side Anti-Reflection coating
Powr	Maximum	Max 130W
1 Own	Standby Mode	Max 20W
Control Key	Front Side	INPUT, -, +, , , PIP, MENU, POWER
Input Signal	Video	1xDVI, 1x Optical DVI 1(DVI 2 option), 1x D-SUB, 1xBNC (CVBS)Video, 2 x BNC (SVHS Y/C), 1 x BNC(SDI), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1xDVI, 1xBNC (SDI)
Input Power	DC 24V, 6.25A Max	
Dimension	Size and Weight	618(W)X412(H)X99.5(D) (mm) 8.8Kg -Without stand 24.330(W)X16.220(H)X3.917(D) (inch) 194.bs-Without stand

FS-L320*D

Model		FS-L320*D	
Ortical	Туре	32" TFT-LCD	
	Screen Size	32 inch	
	Maximum Resolution	1920X 1080 @ 60Hz	
	Pixel Pitch	0.363(H) mm X 0.363(V) mm	
Optical Characteristics	Display Colors	1.06 billion	
	Contrast Ratio(Typical)	850:1	
	Viewing Angle	89° / 89° / 89° / 89°	
	Response Time	6msec(Gray to Gray)	
	Luminance(Typical)	380cd/m ²	
Resolution	Video input	VGA~1080P	
Front Filter	Glass	Double side Anti-Reflection coating	
Powr	Maximum	Max 150W	
TOWI	Standby Mode	Max 20W	
Control Key	Front Side	INPUT, -, +, , , PIP, MENU, POWER	
Input Signal	Video	1xDVI, 1xOptical DVI 1(DVI 2 option), 1x D-SUB, 1xBNC (CVBS)Video, 2 x BNC (SVHS Y/C), 1 x BNC(SDI), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)	
Output Signal	Video	1xDVI, 1xBNC (SDI)	
Input Power	DC 24V, 7.5A Max	C 24V, 7.5A Max	
Dimension	Size and Weight	795(W)X485(H)X100(D) (mm) 15.8Kg -Without stand 31.299(W)X19.094(H)X3.937(D) (inch) 34.8bs-Without stand	

SFS-L420*D

Model		FS-L420*D
Ortical	Туре	42" TFT-LCD
	Screen Size	42 inch
	Maximum Resolution	1920X 1080 @ 60Hz
	Pixel Pitch	0.4845(H) mm X 0.4845(V) mm
Optical Characteristics	Display Colors	1.06 billion
	Contrast Ratio(Typical)	650:1
	Viewing Angle	89° / 89° / 89° / 89°
	Response Time	5 msec(Gray to Gray)
	Luminance(Typical)	400cd/m
Resolution	Video input	VGA~1080P
Front Filter	Glass	Double side Anti-Reflection coating
Powr	Maximum	Max 260W
FOWI	Standby Mode	Max 20W
Control Key	Front Side	INPUT, -, +, , , PIP, MENU, POWER
Input Signal	Video	1xDVI, 1xOptical DVI 1(DVI 2 option), 1xD-SUB, 1xBNC (CVBS)Video, 2x BNC (SVHS Y/C), 1 x BNC(SDI), 5 x BNC (Component Y/G, Pb/B, Pr/R, H/CS, VS)
Output Signal	Video	1xDVI, 1xBNC (SDI)
Input Power	AC 100~230V, 50~60Hz,3A Max	
Dimension	Size and Weight	1024.6(W)X617.4(H)X111.1(D) (mm) 28Kg -Without stand 40.338(W)X24.307(H)X4.374(D) (inch) 61.7bs-Without stand

Thank you for choosing our product.

Service

Please contact our customer service if you need any information or help with our products.

Warranty

One year, parts and labor

Customer service contact

D&T Inc.

Daedeok Valley, 59-9, Jang Dong, Yuseong Gu, Daejeon, Korea, 305-343 Tel.+82(42)360-8000 Fax. +82(42)360-8005 http://www.dntinc.co.kr

USA

FORESEESON CUSTOM DISPLAY INC

2210 E. Winston Road, Anaheim, CA 92806 USA http://www.foreseesonusa.com Tel.1-714-300-0540 Fax. 1-714-300-0546

EC Representative Foreseeson GMBH

Benzstr.9 61352. Bad, Homburg, Germany Tel.+49-(0)6172-18531-030 Fax.+49-(0)6172-18531-011

> Part No. 949494060008-02DT