Nano**EnTek**



FREND™ Vitamin D Total 25-Hydroxyvitamin D

Intended use

The FREND™ Vitamin D test is a fluorescent nanoparticle immunoassay designed for *in vitro* quantitative measurement of total 25-hydroxy(OH) vitamin D and related hydroxylated metabolites in human serum and plasma(EDTA, Lithium-heparin and Citrate). The FREND™ Vitamin D microfluidic flow cartridge is designed for use in the FREND™ System to aid in the assessment of vitamin D sufficiency.

Summary and explanation of test

Vitamin D is a fat-soluble prohormone known for its role in regulating calcium and phosphorus levels in bone mineralization¹². Sunlight exposure produces vitamin D via photochemical conversion of 7-dehydrocholesterol in the epidermis and is the primary source of vitamin D¹². Seasonal changes, amount of exposure, sunscreen use, and skin pigmentation can cause variation in the amount of vitamin D produced in the body. A minor source of vitamin D can be absorbed from food and vitamin supplements, with an estimated 10-20% absorbed by the body in this manner¹².

In circulation, 25-OH vitamin D is bound to vitamin D binding protein (VDBP) or albumin at 1000 times higher concentrations than the active form 1,25-(OH)2-vitamin D". Additionally, the 25-OH form has a half-life of 2-3 weeks, as compared to the less stable 1,25-(OH)2 form, which has a half-life of a few hours". 25-OH vitamin D exists as D2 (ergocalciferol) and D3 (cholecalciferol) isomers, with supplements available for both. Often, total 25-OH vitamin D is measured to assess the sufficiency in a patient and make appropriate clinical decision".

Principle of the assay

The FREND™ Vitamin D test cartridge is a one-time-use rapid "competitive" immunoassay utilizing fluorescent nanoparticle in microfluidic flow to capture and quantify total 25-OH vitamin D in serum and plasma specimens. The drop of 70 µL patient sample is placed in the FREND™ Vitamin D pretreatment tube, where the sample interacts with a proprietary mix of pretreatment solution. Initially, patient sample is mixed with vitamin D antibody-labeled particles, forming immune complexes with total 25-OH vitamin D in the patient sample, then, incubated for 10 minutes at 49°C. The drop of 35 µL mixture is added to the FREND™ Vitamin D cartridge, interacting with vitamin D conjugated fluorescent nanoparticles. The mixture moves via Capillary action to the detection region, where fluorescent nanoparticle complexes are grabbed. The fluorescence intensities from the complexes are measured and total 25-OH vitamin D concentration is calculated by the FREND™ System.

Material provided

One reagent kit (box) contains: Catalogue number

- FREND™ Vitamin D cartridge(s)
 FREND™ Vitamin D pretreatment tube(s)
- 20 FREND™ Vitamin D dilution tube(s)
- 30 Disposable pipette tip(s)
- 01 FREND™ Vitamin D Code chip
- 01 FREND™ Vitamin D Package Insert

FRVD 020

One cartridge contains:

Monoclonal mouse anti-vitamin D	8±0.8 ng
Monoclonal anti-testosterone antibody	76±7.6 ng
Fluorescent particles	4.8±0.48 ng

One pretreatment tube contains:

Gold nano particle conjugation antibody 5±0.5 μg

One dilution tube contains:

Perfluorohexanoic Acid 1.4±0.14 mg

Materials required but not provided

The FREND[™] System Micro-pipette capable of delivering 35 and 70 μ L Heating block for tube and cartridge incubation at 37 $^{\circ}$ C Timer for incubation step Personal protective equipment and biohazard waste

Warning and Precautions

- The FREND™ Vitamin D cartridges are intended for in vitro diagnostic use only.
- Vitamin D cartridges are only to be used on the FREND™ System.
- Vitamin D cartridges are disposable, single use devices. Do not reuse them under any circumstances.
- Allow sealed cartridges to come to room temperature for 15~30 minutes prior to use.
- · Reagent kit should not be frozen.
- Assure the humidity in the laboratory is in the 10~80% range when tests are run.

- Avoid cross-contamination between samples by using a new pipette tip for each new specimen.
- · Avoid high humidity, direct sunlight or heat in the area used for cartridge storage.
- · Inaccurate results are possible if the sample used is contaminated in any way.
- Using specimens containing clotted fibrin could result in erroneous results.
- · Over or under loading the cartridge with sample may result in inaccurate results.
- Human specimens are not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the lab may be derived from human materials. Please use Universal Precautions when handling all specimens and controls.
- Do not use the reagent components beyond the expiration date on the pouch.
- Do not use the reagent components if the pouch is damaged or the seal is broken.
- Perform testing as specified in the Package Insert and User Manual.
- · Keep the reagent components sealed in the pouch until just ready for use.
- Use the reagent component immediately after opening the pouch.
- Wear disposable gloves when handling the reagent components and the samples.
- · Wash hands thoroughly and often after handling reagent cartridges or samples.
- Do not ingest the silica gel package found in the cartridge pouch.
- Vitamin D has been designed so that the high dose "hook effect" does not affect the vast majority of samples.

Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature. Reagent stability has been demonstrated for twelve months from the date of manufacture.

The expiration date is clearly indicated on the product box and the cartridges.

Materials Catalogue number

FRVD 020

None

Refrigerator temperature storage (2~8 ℃)

FREND™ Vitamin D cartridges
FREND™ Vitamin D pretreatment tubes

FREND™ Vitamin D dilution tubes None

Specimen collection and handling

Human serum and plasma (Lithium-heparin, EDTA and Citrate) samples are suitable for use with FREND™ Vitamin D cartridges.

Follow instructions detailed in this package insert as well as the specimen collection tube manufacturer's instructions for specimen collection and preparation (including manufacturer's instructions for centrifugation time and speed.)

For serum, a blood sample is collected aseptically without additives by venous puncture. After allowing the sample to clot for 30 minutes at room temperature, the collection tube should be centrifuged for 10 minutes at 3,000 rpm.

For plasma (Lithium-heparin, EDTA and Citrate), a venous blood sample is collected aseptically with the designated additive. After allowing the specimen to sufficiently mix with anticoagulant at room temperature, the sample tube can be centrifuged for 10 minutes at 3,000 rpm.

Samples may be stored at 2^{-8} °C for up to 6 hours prior to analysis. If the analysis is scheduled to be done at some later time, the sample should be stored frozen at -20 °C or below for future use.

Repeated freeze-thaw cycles should be avoided. Prior to assay, slowly bring frozen samples to room temperature (18~25 $^{\circ}$ C) and mix gently but thoroughly before test.

For optimal results, avoid grossly hemolytic, lipidic, or turbid specimens. Specimens should be free of aggregated fibrin, red blood cells, or other particulate matter.

When pipetting into the FREND™ Vitamin D cartridge sample inlet, ensure that bubbles in the sample are avoided. Bubbles may restrict flow and result in an incomplete or erroneous test result.

Procedure

Calibration

There is no need for calibration to be performed by the end user as is generally required on other automated laboratory equipment. All calibration statistics and information have been electronically stored on the FREND™ Vitamin D Code chip included in each box of FREND™ Vitamin D cartridges. The FREND™ Vitamin D Code chip is specific for each manufactured lot of FREND™ Vitamin D cartridges.

Calibration information should always be checked by running external quality control samples to verify that the results obtained for Vitamin D on the FREND System using the FREND Vitamin D cartridges of a specific lot meet the laboratory criterion for acceptability.

Code chip installation

Please refer to the FREND™ System user manual for more detailed instructions relative to the Code chip installation. Abbreviated instructions follow here:

- Insert the FREND™ System electrical cord into an appropriate outlet.
- Insert the Code chip into the Code chip slot at the rear of the FREND™ System following the arrows.
- 3) Press the 'Setup' button on the 'Main' screen.
- 4) Press the 'Code chip' button on the 'Setup' screen.
- 5) The information embedded on the FREND™ Vitamin D Code chip is automatically saved on the FREND™ System.
- 6) When the Code chip installation is completed, press the 'OK' button to go to the 'Setup' screen.
- 7) Press the 'Item' button on the 'Setup' screen.
- Check the FREND™ Vitamin D cartridge lot number and the installation date of the Code chip.
- Press the 'Home' button to go to the 'Main' screen to begin running external quality control and patient samples.

Quality control

FREND™ System OC cartridges

FREND™ QC Cartridge contains multiple controls to check optic part of the system. By testing QC Cartridge, part of analytical components of the system of (1) laser power, (2) alignment, and (3) mechanical integrity are confirmed.

For each day of patient testing perform QC Cartridge testing. Refer to the quality control procedures section in the User Manual of FREND™ System. In brief, perform QC Cartridge testing for the following conditions:

- (1) Upon initial setup of the system,
- (2) Each day of patient testing,
- (3) When the system has been transported or moved,
- (4) Whenever there is uncertainty about the performance of the system,
- (5) Whenever required by your laboratory's quality control requirements.

Internal procedural controls

FREND™ Vitamin D test cartridge contains built-in control feature. Fluorescence signal in the reference zone of each cartridge shows: (1) that enough volume is added, (2) that proper flow is obtained, and (3) that the antibody is reactive. If this reference zone signal is missing or lower than threshold, the FREND™ System consider it as an incorrect or failed test, not producing a test result but an error message. In addition, with each cartridge run, the system monitors, in part, for (1) flow of sample, (2) speed of sample flow, (3) shelf-life of cartridge components, (4) function of internal barcode scanner, and (5) function of scanner's mechanical components.

External quality control testing

Commercially available controls from a variety of manufacturers are available that contain 25-hydroxy vitamin D as a measured analyte. It is recommended that a minimum of two (2) levels of controls be run at least once per month or once for each new lot, whichever comes earlier. However, Controls should be run with a minimum frequency, depending on number of tests run in the laboratory. Each laboratory should establish its own criteria based on the following parameters:

- (1) Each new lot,
- (2) Each new shipment (even if from the same lot previously received),

- (4) Monthly, as a continued check on storage conditions.
- (5) Whenever problems (storage, operator, or other) are identified.
- (6) Or other times as required by your laboratory's standard QC procedures.

Individual laboratory policy will dictate exactly which control materials and lot numbers should be run, the frequency with which controls are to be tested, criteria for acceptance of the results and required corrective action to be taken if results do not meet laboratory criteria. If any external quality control sample values are out of the acceptable range, it will be necessary to investigate the problem before reporting patient results to assure there is not an instrument or software malfunction. Do not assay patient samples on the FREND™ System using FREND™ Vitamin D if quality control on how to determine acceptability of external control material results. Each laboratory operates under a different set of regulations. Every laboratory must follow the standardized procedures acceptable to the regulatory agencies to which the laboratory is responsible.

Specimen processing

Preparation

Remove sufficient cartridges and pretreatment tubes of FREND™ Vitamin D from the refrigerator to test the number of patient samples and required external quality materials. Allow the tubes and the sealed pouches containing the cartridges to come to room temperature for 15~30 minutes prior to the start of the testing sequence. Heating block provided with FREND™ System should be turned on 7~8 minutes before use.

If using refrigerated patient samples, remove those from the refrigerator and allow to them to come to room temperature prior to testing. If frozen samples will be utilized, be sure these are removed from the freezer, thawed naturally and then mixed gently but thoroughly prior to testing. Testing should not begin on these previously frozen samples until they have reached room temperature.

There are no other reagents or sample preparations necessary.

Assay procedure

 Prepare the FREND™ Vitamin D cartridge, pretreatment tube, dilution tube and specimen. Open the pouch and place the FREND™ Vitamin D cartridge on the heating block.

- 2) Transfer 35 uL of specimen to the dilution tube and mix well.
- 3) Transfer 70 μ L of diluted sample to the pretreatment tube and mix well.
- △Caution: Pellet in the tube bottom should be completely dissolved.
- 4) Put the tube in the hole of heating block and incubate for 30 minutes.
- 5) Pipette 35 µL of the incubated sample into the sample inlet on the cartridge using a suitable micro-pipette equipped with a fresh pipette tip.
- 6) Press the 'Test' button on the 'Main' screen of the FREND™ System.
- 7) The system moves to the Patient ID screen automatically.
- 8) Type the Patient ID and press the 'Enter' button to begin the test.
- 9) Insert the cartridge into the cartridge slot using the cartridge arrow as a guide.
- ⚠Caution: Please check the direction of the cartridge before insertion and assure the insertion is complete. It is recommended to insert the cartridge after the sample injection after 30 seconds elapsed in less than 5 minutes to obtain the optimal result of test.
- 10) When the reaction in the cartridge is completed, the FREND™ System will automatically begin the reading process.
- When the measurements are completed, the cartridge will automatically be expelled and the results displayed.
- $\hat{\mathbb{L}}$ Caution: Do not disconnect power cord or shut off power from the FRENDTM System while a cartridge is in the reading chamber. This may cause a system error.
- 12) If the FREND™ System is connected to the optional printer, press the 'Print' button and the results will be output on the printer paper.

Procedural notes

If a specimen Vitamin D concentration is found to be greater than the linearity limit of the assay of 110.00 ng/mL and a definitive result is required, the specimen should be diluted with low concentration sample that has been previously measured on the FREND™ Vitamin D and then re-assayed according to the 'Assay Procedure'. The recommended dilution for samples with an initial result of >110.00 ng/mL is 1:2. Dilutions must be made manually and the final result on the diluted sample calculated manually by multiplying the result obtained on the diluted sample by the dilution factor.

** Concentration (high conc. sample) = Concentration (1:2 dilution sample) X 2 - Concentration (low conc. sample)

Calculation of results

The FREND™ System performs all sample and reagent handling operations automatically within the cartridge once the sample has been manually loaded to the sample inlet in the cartridge and the cartridge placed into the FREND™ System. The rate of fluorescence produced by the reaction is read at various intervals during the analysis process, blank reading are subtracted after which the net rate is automatically converted to Vitamin D concentration in ng/mL based upon information stored on the FREND™ Vitamin D Code chip.

This result is then output on the screen and to the optional printer. It is also stored in memory on the FREND $^{\rm m}$ System.

Screen displayed for various concentration scenarios

Displayed result	Description
General Section 10 Process Conference 10 Pro	Vitamin D concentration Less than 10.00 ng/mL
Common Principal State	Vitamin D concentration Not less than 10.00 ng/mL and not higher than 110.00 ng/mL
	Vitamin D concentration Higher than 110.00 ng/mL

^{¾ 100ng/dL=1ng/mL}

Limitations of the procedure

- When used for diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g., symptoms, results of other tests, clinical impressions, medical history, therapy, etc.)
- 2) The FREND™ System paired with a FREND™ Vitamin D cartridge, is programmed to report 110.00 ng/mL as the highest concentration of Vitamin D measurable without dilution. The lowest measurable concentration is 10.00 ng/mL - the assay limit of detection.
- 3) Specimens from patients with heterophilic antibodies, such as anti-mouse (HAMA), anti-goat (HAGA), or anti-rabbit (HARA) antibodies, maybe show falsely elevated or depressed values or may result in an incomplete test[11, 12]. Patients routinely exposed to animals or animal serum products can be prone to these types of heterophilic interferences.
- 4) Certain medications may interfere with assay performance. All results should be interpreted with respect to the clinical picture of the patient.
- 5) Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- 6) Lipemia has an insignificant effect on the assay except in the case of gross lipemia where interference with the lateral flow of the sample in the cartridge may occur.
- 7) The concentration of Vitamin D in a given sample determined with assays from different manufacturers can vary due to differences in assay methods, calibration, and antibody specificity.
- 8) Please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this insert sheet.
- FREND™ Vitamin D has not been validated in point-of-care settings.
- 10) FREND™ Vitamin D is to be used in licensed clinical laboratories with trained technologies.

Performance characteristics

Performance characteristics were evaluated for the FREND™ Vitamin D as follows:


Precision

Clinical and Laboratory Standards Institute (CLSI) document EP5-A2 was utilized as a guidance for precision studies on FREND[™] Vitamin D. Five human serum based panels were assayed in two replicates, twice per day over a period of 20 days.

sample		Repea	Repeatability		Between-run		Between-day		Within-laboratory	
ID	Value (ng/mL)	SD	CV%	SD	CV%	SD	CV%	SD	CV%	
В	18.00	1.443	8.0	0.348	1.9	0.282	1.6	1.511	8.4	
C	30.00	1.643	5.5	0.655	2.2	0.166	0.6	1.777	5.9	
D	60.00	3.112	5.1	0.893	1.5	1.142	1.9	3.433	5.7	
E	90.00	4.853	5.4	2.201	2.5	0.787	0.9	5.386	6.0	
F	120.00	5.486	4.6	2.606	2.2	1.333	1.1	6.218	5.2	


Dilution linearity

The range of linearity in serum was established using CLSI document EP6-A. FREND™ Vitamin D was found to have linearity within the reportable range of 8.00 ng/mL ~ 130.00 ng/mL with a mean recovery of 100±10%. The dilution linearity study was performed by diluting a high concentration Vitamin D specimen with a low concentration Vitamin D specimen. Results are summarized in the graphs below, but may vary in individual laboratories.



Comparative analysis

FREND™ Vitamin D was compared to the predicate device using guidelines outlined in CLSI document EP09-A2-IR. Samples (n=120) were measured in duplicate on both systems. Linear regression analysis was utilized to demonstrate a correlation coefficient (R2) of ≥0.92.



Specificity

FREND™ Testosterone demonstrates cross-reactivity(%) with the following related substances at the concentrations indicated bellows. CLSI document EP7-A2 was used to analyze and calculate cross-reactivity using serum samples at three levels of Testosterone within the normal range.

		Vitamin D concentration level			
Substance	Concentration	Low	Median	High	
		Cro:	s reactivity	(%)	
Vitamin D2	500 ng/mL	0.0	0.8	0.8	
Vitamin D3	500 ng/mL	0.2	-0.8	-0.5	
1,25-(OH)2-Vitamin D2	100 ng/mL	-0.7	-1.0	2.2	
1,25-(OH)2-Vitamin D3	100 ng/mL	0.0	0.9	0.3	
3-epi-25(OH) Vitamin D3	400 ng/mL	-0.8	-1.2	-1.8	
25(OH) Vitamin D2	25 ng/mL	103.2	106.5	96.7	
25(OH) Vitamin D3	25 ng/mL	113.5	99.4	97.7	

Analytical sensitivity

The limit of blank (LoB) and limit of detection (LoD) was determined using guidelines found in CLSI document EP17-A. LoB was determined from 60 replicate measurements using a calibrator A(vitamin D depleted serum). LoD was determined using 12 replicates measurements of five low level patient samples.

LoB (ng/mL)	LoD (ng/mL)	LoQ (ng/mL)
3.84	7.07	7.07

Interference

FREND™ Vitamin D demonstrates ≤15% interference with the following substances at the concentrations indicated below: Bilirubin, Biotin, Triglycerides, Human albumin, Hemoglobin. CLSI document EP7-A2 was used to analyze and calculate percent interference using serum samples at three levels of Vitamin D within the normal range.

Substance	Concentration
Bilirubin	2.5 mg/dL
Biotin	30.0 ng/mL
Triglyceride	370.0 mg/dL
Human albumin	10.7 g/dL
Hemoglobin	200.0 g/dL

References

- Holick, MF. High prevalence of vitamin D inadequacy and implications for health. Mayo Clinic Proceedings. 2006, 81 (3): 353-373.
- Bolland, MJ et al. The effect of vitamin D supplementation on skeletal, vascular, and cancer outcomes: a trial sequential meta-analysis. The Lancet Diabetes & Endocrinology, 2014, 2 (4): 307-320.

- 3) Wolf, G. The discovery of vitamin D: the contribution of Adolf Windaus. Journal of Nutrition, 2004, 134 (6): 1299-1302.
- Holick, MF. Vitamin D: importance in the prevention of cancers, type 1 diabetes, heart disease, and osteoporosis. The American Journal of Clinical Nutrition. 2004. 79 (3): 362-371.
- Calvo, MS et al. Vitamin D intake: a global perspective of current status. Journal of Nutrition, 2005, 135 (2): 310-316.
- 6) Ross, AC et al. Dietary Reference Intakes for Calcium and Vitamin D. Washington DC: National Academies Press 2011, p. 435.
- 7) Svasti, J et al. Molecular basis for the three major forms of human serum vitamin D binding protein. Biochemistry 1979, 18 (8): 1611-1617.
- 8) Holick, MF et al. Vitamin D2 is as effective as vitamin D3 in maintaining circulating concentrations of 25-hydroxyvitamin D. The Journal of Clinical Endocrinology and Metabolism, 2008, 93 (3): 677-681.
- 9) Maxmen, A. Nutrition advice: the vitamin D-lemma. Nature 2011, 475 (7354): 23-25.
- Chao, E.L.; Henshaw, J.L. Occupational Safety and Health Administration: Model Plans and Programs for the OSHA Bloodborne Pathogens and Hazard Communications Standards. OSHA 3186-06R, 2003.
- Schroff, R.W.; Foon, K.A.; Beatty, S.M.; Oldham, R.K.; Morgan, A.C. Human Anti-Murine Immunoglobulin Responses in Patients Receiving Monoclonal Antibody Therapy. Cancer Research. 1985, 45: 879-885.
- Boscato, L.M.; Stuart, M.C. Heterophilic Antibodies: A Problem for All Immunoassays. Clinical Chemistry, 1988, 34 (1): 27-33.

Glossary of symbols

8	Do not reuse
Exp. Dute (YYYY-MM-DD)	Use by YYYY-MM-DD
LOT	Lot number
REF	Catalog number
\triangle	Warning or Caution
***	Manufactured by
EC REP	Authorized representative in the Europe Community
IVD	In vitro diagnostic medical device
215-815	Temperature limitation
)n'	Contains sufficient for <n> tests</n>

Nano**EnTek**

sales@nanoentek.com

Manufactured by :

NanoEnTek Inc. (HQ)

12F, 5, Digital-ro 26-gil, Guro-gu, Seoul, 08389, Korea

EC REP

MT Promedt Consulting GmbH

Altenhofstrasse 80, 66386 St. Ingbert, Germany

US Branch

NanoEnTek USA Inc.

5627 Stoneridge Drive Suite 304, Pleasanton, CA 94588, USA Tel:+1-925-225-0108. +1-888-988-0108 (Toll free) Fax:+1-925-225-0109

China Branch

SK Medical (Beijing) Co., Ltd.

26F, SK Tower No.6 Jia, Jianguomenwai Avenue, Chaoyang District, Beijing, 100022 P.R. China Tel:+86-10-5920-7844. Fax:+86-10-5920-5697