**CHATTANOOGA ionto™**

**Dual Channel Iontophoresis System**

**Instructions For Use**

---

---

# **Chattanooga Ionto™ Dual Channel Iontophoresis System**

## **Instruction for Use Manual**

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## INTRODUCTION

The Chattanooga Ionto™ Dual Channel Iontophoresis System operates on the principle that electrically-like charges repel each other. A positively charged ion is forced away from a positively charged electrode.

## SAFETY AND CLINICAL USE

### **Indications**

Iontophoresis is an alternative to hypodermic injection for the administration of ionic solutions.

### **Contraindications**

- Known adverse reactions to the application of electrical current.
- Cardiac pacemakers or other electrically sensitive implanted devices.
- Known sensitivity to the drugs to be administered.

### **Warnings**

- The safety of electrical stimulation during pregnancy or delivery has not been established.
- Do not apply electrodes across the temporal region as safety has not been established.
- Avoid operation in close proximity (< 1m) to transmitting cellular (wireless) telephones or two-way radios. This equipment may produce sudden, unexpected changes in the device output.
- Do not apply electrodes over broken or compromised skin (e.g. sunburn, cuts, acne) due to increased risk of skin reactions.
- Keep out of the reach of children.
- Lead wires and electrodes should be removed before using equipment that generates high frequency or high energy electromagnetic radiation such as welding equipment, surgical electrocautery units, short-wave and microwave therapy. Operation in close proximity (< 3m) to this equipment may startle the user by causing output instability or improper operation of the stimulator. Simultaneous connection to the patient may result in burns and possible damage to the stimulator.
- Do not plug lead wires into power outlets such as wall sockets and line cord receptacles. Doing so could result in severe shock or burns whether or not the lead wires are attached to the stimulator.
- Do not use any other electrodes with this system. The safety of using electrodes other than the Chattanooga Ionto Buffered Iontophoresis Electrodes with the Chattanooga Ionto Iontophoresis Device has not been established.

- Do not apply external pressure, heat or cold to the electrode sites during treatment.

### **Precautions**

- Consult the directions for the drug being administered for indications, contraindications and warnings concerning its use before giving an iontophoresis treatment.
- The current setting should never exceed a comfortable level for the patient. Advise the patient to report any sensation of pain or burning during the treatment.
- Avoid exposure to water. If accidentally immersed, remove the battery and drain excess water. Air dry the device for 48 hours before operating.
- The maximum recommended dosage is 80mA/min.
- Patients with known skin sensitivities or allergies should be treated with lower current settings and longer treatment times to reduce the risk of skin reactions.
- During iontophoresis, electrical discharges (sparks) are possible that could ignite flammable material or gases.

### **Adverse Reactions**

Iontophoresis can cause transient erythematic reaction, inflammation of the skin, skin irritation or skin burns. This mild histamine response, characterized by red, warm, raised skin with occasional small white bumps the size of a pinhead, will usually disappear within 12 hours of the treatment. Advise the patient of this possibility before starting treatment. If a skin reaction does occur, discontinue treatment and consult the prescribing physician. If a skin reaction occurs, the patient should be evaluated for possible drug sensitivity.

## PRODUCT DESCRIPTION

The dual channel Chattanooga lonto lontophoresis System delivers ionic drug solutions directly to the body site being treated and offers:


- Two channels for treatment of two different sites simultaneously.
- Automatic calculation of time based on set dose and amplitude.
- Self diagnosing capability that alerts the user when the circuit is not complete.

System components and accessories may include:

- Chattanooga lonto lontophoresis Device
- User's Manual
- Chattanooga lonto Buffered lontophoresis Electrodes
- Chattanooga lonto lontophoresis Lead Wires
- 9-Volt Batteries

## Controls And Features

1. Channel Output Jacks connect lead wires to the device. Channels 1 and 2 are independently controlled by Current Knobs 1 and 2.
2. ON/OFF Indicator Lights illuminate continuously when the corresponding channel is delivering current. They flash when there is an interruption of current flow on the corresponding channel.
3. Current Knobs control the current output of Channel 1 and 2 between 0.5 and 4mA.
4. Treatment Mode Switch has four settings.
  - CH1 Set Up Mode, (S1), allows the target dose for CH1 to be set.
  - CH2 Set Up Mode, (S2), allows the target dose for CH2 to be set.
  - Pause Mode, (P), pauses the output during a treatment.
  - Run Mode, (R), initiates or resumes treatment.
5. Power Switch turns the device ON or OFF.
6. Display Button selects which parameter is shown in the display window, either Time (T), Dose (D) or Current (C).
7. Dosage Buttons ( $\Delta$  or  $\nabla$ ) set the dosage (mA/min) for each channel. The Treatment Mode Switch must be set to (S1) or (S2) to set dosage for CH1 or CH2 respectively.

8. Battery Compartment is located under the Front Panel Cover.


9. Display window shows the treatment settings and has indicators for low battery or circuit problems. Components of the display include:

a) CH1 and CH2 indicators

- a pointing icon  indicating which channel is having problems
- a numeric display showing the parameter settings

b) TIME, DOSE or CURRENT indicating which parameter is being displayed

- the units of the parameter being displayed
- a battery icon indicating low battery status

## OPERATING INSTRUCTIONS

Read the Instruction Manual thoroughly for indications, contraindications, warnings and precautions, and methods of operation before initiating treatment.

**STEP 1:** Assess Skin Condition at electrode sites.

**STEP 2:** Prepare the Skin by cleaning and drying it thoroughly. Do not use rubbing alcohol unless skin is excessively oily. If necessary, trim excess body hair at application site before applying electrodes. DO NOT shave application sites 24 hours prior to treatment.

**STEP 3:** Connect the Lead Wires to the Electrodes. Pay careful attention to the polarity required for the treatment electrode. Always match the polarity of the drug ion with the polarity of the drug delivery electrode.

**STEP 4:** Prepare the Electrodes following the instructions accompanying the electrodes.

**STEP 5:** Place the Electrodes on skin.

**STEP 6:** Connect the Lead Wires to the Device. Caution: Be sure the device is OFF before connecting the lead wires.

**STEP 7:** Set the Power Switch to ON.

**STEP 8:** Set the Desired Treatment Dosage for CH1 using the (S1) mode. The display will start with a dosage of 40mA/min.

**STEP 9:** Adjust the dosage for CH1, by pressing the  $\Delta$  or  $\nabla$  Dosage Button until the correct dosage is displayed.

- 
- STEP 10: Adjust the dosage for CH2, using the (S2) mode, and repeat Step 9.
- STEP 11: Set the Current Level for each channel using the (P) mode. The current level will automatically appear in the display. Adjust to the desired level using the current knobs. If the desired level is unknown, initiate treatment at a lower level (< 2mA).
- STEP 12: Initiate Treatment by setting the Treatment Mode Switch to (R). Current will ramp up to set intensity level in 30 seconds.
- STEP 13: At the End of Treatment current will ramp down to zero, then the device will beep and shut down. To turn the device OFF, the Power Switch must be in the OFF position. Disconnect the lead wires and remove the electrodes.

### **Operational Considerations**

---

- During treatment, (P) or (R) mode, the following parameters can be viewed by using the Display Button:
  - Time remaining [min]
  - Dose delivered [mA/min]
  - Current level [mA]
- During setup, (S1) or (S2) mode, Dose represents the total dose to be delivered and is the only treatment parameter which is displayed.
- The current setting can be changed at any time during the treatment without going into the (P) mode. When current is changed, the time remaining is automatically recalculated, ensuring the desired dosage is delivered.
- Pausing the treatment (P) does not affect the actual treatment length or dosage. For safety and comfort reasons, current will automatically ramp down when treatment is paused (P) and ramp up when treatment is resumed (R).
- When using both channels, treatment will start simultaneously.
- Treatment dosage can be retained for use at the next treatment by turning the Treatment Mode Switch to (P), (S1) or (S2). If not used within 20 minutes, the device will beep until power is turned OFF. The CH1 default dosage at start up is 40mA/min.

## **TROUBLESHOOTING AND REPAIR**

---

### **Troubleshooting**

---

- If the device does not function, make sure the battery is fresh, properly installed (check polarity markings) and the battery contacts are clean.
- If the  appears on the display, replace the battery. If it begins blinking, the device will shutdown soon. Record how much of the treatment was completed because this information will be lost when the battery is replaced.
- If an electrode becomes disconnected or if a high impedance situation occurs, the treatment is automatically paused, the ON/OFF indicator light for the impacted channel blinks and a beep sounds. The  will also indicate which

channel requires correction. To correct the situation, you must select (P) and identify and fix the fault in the electrical circuit. Select (R) to resume the treatment.

- If the power is not turned off after a treatment, the device will begin to sound a reminder beep after 20 minutes. The beep will continue every 30 seconds until the power is turned off or the battery is drained.

## **Repair**

There are no user serviceable parts inside the stimulator. If the unit appears to be non-functional, contact your Authorized Chattanooga Group Distributor, or contact Chattanooga Group directly.

## **MAINTENANCE, CLEANING, STORAGE AND DISPOSAL**

### **Maintenance**

Under normal conditions, the stimulator does not require periodic maintenance or calibration; however, all the information necessary for a qualified technician to perform periodic safety checks is included in the Technical Information Section.

### **Changing the Battery**

#### **⚠ CAUTION:**

Turn the device off and disconnect the lead wires before inserting a fresh battery. Change the battery when the  appears on the display or if the device will not turn on. Be sure the polarity ("+" and "-") markings on the battery match with the markings on the device.

#### **⚠ CAUTION:**

Inserting the battery incorrectly may cause it to rupture or generate intense heat, causing irreversible damage to the battery. If there are signs of this type of damage, discard or recycle the battery and order a replacement.

### **Cleaning**

Use a cloth moistened with soap and water to clean the exterior of the device and the lead wires. Use of other cleaning solutions may damage these items. Never immerse the device in liquids.

### **Storage**

To store the stimulator for an extended period of time, i.e. 90 days or more, remove the battery and store the stimulator in a cool, dry location.

### **Disposal**

Contact your local distributor or Chattanooga Group for information on disposal of the stimulator.

**WARRANTY REPAIR/OUT OF WARRANTY REPAIR**

1. Written statement containing the following information:

- RA Number - Obtain from Factory
- Unit Model Number
- Unit Serial Number
- Contact person with Phone and Fax Numbers
- Billing Address (for Out of Warranty Repair)
- Shipping Address (Where to Ship Unit after Repair)
- Detailed Description of Problem or Symptoms

2. Copy of original invoice issued at purchase of the unit.

3. Ship unit to Factory in the original container with all accessories and information as required in item one (1) above to:

Chattanooga Group  
4717 Adams Road  
Hixson, TN 37343  
Phone USA: 800-592-7329  
Outside USA: +1 (423) 870-7200  
FAX: (423) 875-5497

Service to these units should be performed only by Service Technicians certified by Chattanooga Group.

---

---

## TECHNICAL INFORMATION

<b>Accessories</b>	
Battery Type	9V Alkaline battery.
Electrodes	Any Chattanooga Iontophoresis electrode with a minimum active area of 8cm <sup>2</sup> .
Lead Wires	 Any Chattanooga Iontophoresis lead wire with the standard device connection as pictured.
<b>Standard Measurement Conditions</b>	77° F (25° C); 8.5V d.c. supply voltage
Output Specifications	± 15% unless noted
Waveform	Constant direct current
Current Range (100Ω to 12kΩ load)	0.5 to 4mA ±10%
Root Mean Square/ Effective Value (1kΩload)	<10mA
Dose Range (1mA/min increments)	0 to 160 mA/min
Dose Accuracy	± 2 mA/min
Ramp Up/ Ramp Down	30s
<b>Battery Information</b>	
Supply Voltage Range	5.6V to 10V
Low Voltage Indicator Threshold	6.5V
Low Voltage Shutdown	5.3V
Expected Battery Life (one channel @ 4mA)	7h (1700 mA/min)
<b>IEC 601 Standard Classification</b>	Type BF Applied Part. Internally powered only. Ordinary protection against entry of liquids. Continuous operation. Not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
<b>Physical Dimensions</b>	
Size (H x W x D)	3.7in H x 2.4in W x 0.9in D (9.4cm H x 6.2cm W x 2.2cm D)
Approx. Weight (without battery)	2.45oz (77g)
<b>Environmental Conditions:</b>	
Transportation and Storage	-13° F to 176° F (-25° C to +80° C); 90% humidity @ 104° F (40° C); normal atmospheric pressure. Store in a cool, dry place.
  Operation	50° F to 104° F (10° C to +40° C)

## DESCRIPTION OF DEVICE MARKINGS

One or more of the following markings may appear on your device.

The Notified Body, Intertek Semko AB 0413, has granted Chattanooga Group an EC Certificate according to Annex II, Clause 3 of Council Directive 93/42/EEC. (Medical Device Directive.)

CSA C22.2 No. 125-M1984  
Electromedical equipment, Canadian Electrical Code Part II: Safety Standards for Electrical Equipment. CSA Risk Class 2.

EN60601-1-2: 1993 (EMC Directive)  
In accordance with Article 10(1) of Council Directive 89/336/EEC. Medical electrical equipment. Part 1: General requirements for safety. Part 2.  
Collateral standard: Electromagnetic compatibility: Requirements and tests.

Classified by Underwriters Laboratories Inc<sup>®</sup> with respect to electric shock, fire and mechanical hazards only in accordance with UL 60601-1, and CAN/CSA C22.2 No. 601.1-M90.

## LIMITED WARRANTY

Chattanooga Group ("Company") warrants that the Chattanooga Ionto™ Dual Channel Iontophoresis System ("Product") is free of defects in material and workmanship. This warranty shall remain in effect for one year (12 months) from the date of original consumer purchase. If this Product fails to function during the one year warranty period due to a defect in material or workmanship, Company or the selling dealer will repair or replace this Product without charge within a period of thirty (30) days from the date on which the Product is returned to the Company or the dealer.

All repairs to the Product must be performed by a service center authorized by the Company. Any modifications or repairs performed by unauthorized centers or groups will void this warranty.

To participate in warranty coverage, this Product's warranty registration card (included with Product) must be filled out and returned to the Company by the original owner within ten (10) business days of purchase.

### **This Warranty Does Not Cover:**

Accessories including, but not limited to, chargers, rechargeable batteries, electrodes, lead wires, are excluded from the Warranty and are sold "AS IS" because their structure is such that they may be easily damaged before or during use.

Replacement parts or labor furnished by anyone other than the Company, the selling dealer or a certified Company service agent.

Defects or damage caused by labor furnished by someone other than Company, the selling dealer or a certified Company service agent.

Any malfunction or failure in the Product caused by product misuse, including, but not limited to, the failure to provide reasonable and necessary maintenance or any use that is inconsistent with the Product User Manual.

### **COMPANY SHALL NOT BE LIABLE IN ANY EVENT FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES.**

Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

To obtain service from Company or the selling dealer under this warranty:

1. A written claim must be made within the warranty period to the Company or the selling dealer. Written claims made to the Company should be sent to:

4717 Adams Road  
P.O. Box 489  
Hixson, TN 37343 USA  
Telephone: 423-870-2281  
Outside USA: +1 (423) 870-7200  
Facsimile: (423) 875-5497

2. The Product must be returned to the Company or the selling dealer by the owner.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state.

The Company does not authorize any person or representative to create for it any other obligation or liability in connection with the sale of the Product. Any representation or agreement not contained in the warranty shall be void and of no effect.

**THE FOREGOING WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED,  
INCLUDING ANY WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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© 2005 Encore Medical, L.P.
Hixson, TN 37343 U.S.A.
1-800-592-7329 U.S.A.
+ 1-423-870-2046 OUTSIDE U.S.A. FAX