

# **CAM-14 Acquisition Module**

## **field service manual**

PN 421315-001


Revision A



# CAM-14 Acquisition Module

## field service manual

PN 421315-001    Revision A



**marquette**

*A GE Medical Systems Company*

## NOTE

Due to continuing product innovation, specifications in this manual are subject to change without notice.

MD1320-005

© Copyright GE Marquette Medical Systems, Inc. 1998. All rights reserved.

Trademarked names appear throughout this document. Rather than list the names and entities that own the trademarks or insert a trademark symbol with each mention of the trademarked name, the publisher states that it is using the names only for editorial purposes and to the benefit of the trademark owner with no intention of improperly using the trademark.

900 SC, ACCUSKETCH, AccuVision, APEX, AQUA-KNOT, ARCHIVIST, Autoseq, BABY MAC, C Qwik Connect, CardioServ, CardioSmart, CardioSys, CardioWindow, CASE, CD TELEMETRY, CENTRA, CHART GUARD, CINE 35, COROLAN, CORO, COROMETRICS, Corometrics Sensor Tip, CRG PLUS, Digistore, Digital DATAQ, E for M, EAGLE, Event-Link, FMS 101B, FMS 111, HELLIGE, IMAGE STORE, INTELLIMOTION, LASER SXP, MAC, MAC-LAB, MACTRODE, MARQUETTE, MARQUETTE MAC, MARQUETTE MEDICAL SYSTEMS, MARQUETTE UNITY NETWORK, MARS, MAX, MEDITEL, MEI, MEI in the circle logo, MEMOPORT, MEMOPORT C, MINISTORE, MINNOWS, Monarch 8000, MULTI-LINK, MULTISRIPTOR, MUSE, MUSE CV, Neo-Trak, NEUROSCRIPT, OnlineABG, OXYMONITOR, Pres-R-Cuff, PRESSURE-Scribe, QMI, QS, Quantitative Medicine, Quantitative Sentinel, RAMS, RSVP, SAM, SEER, SILVERTRACE, SOLAR, SOLARVIEW, Spectra 400, Spectra-Overview, Spectra-Tel, ST GUARD, TRAM, TRAM-NET, TRAM-RAC, TRAMSCOPE, TRIM KNOB, Trimline, UNITY logo, UNITY NETWORK, Vari-X, Vari-X Cardiomatic, VariCath, VARIDEX, VAS, Vision Care Filter, are trademarks of GE Marquette Medical Systems, Inc., registered in the United States Patent and Trademark Office.

2SL, 15SL, Access, AccuSpeak, ADVANTAGE, BAM, BODYTRODE, Cardiomatic, CardioSpeak, CD TELEMETRY® -LAN, CENTRALSCOPE, Corolation, DASH, EK-Pro, EDIC, Event-Link Cumulus, Event-Link Cirrus, Event-Link Nimbus, HI-RES, ICMMS, IMAGE VAULT, IMPACT.wf, INTER-LEAD, IQA, LIFEWATCH, Managed Use, MARQUETTE PRISM, MARQUETTE® RESPONDER, MENTOR, MicroSmart, MMS, MRT, MUSE CardioWindow, NST PRO, NAUTILUS, OCTANET, O2 SENSOR, OMRS, PHi-Res, Premium, Prism, QUIK CONNECT V. QUICK CONNECT, QT Guard, RAC, SMARTLOOK, SMART-PAC, Spiral Lok, Sweetheart, UNITY, Universal, Waterfall, Walkmom are trademarks of GE Marquette Medical Systems, Inc.

# Contents

## 1

|                                          |            |
|------------------------------------------|------------|
| <b>Introduction .....</b>                | <b>1-1</b> |
| <b>Manual Information .....</b>          | <b>1-3</b> |
| Revision History .....                   | 1-3        |
| Manual Purpose .....                     | 1-3        |
| Intended Audience .....                  | 1-3        |
| <b>Safety Information .....</b>          | <b>1-4</b> |
| Definitions .....                        | 1-4        |
| Messages .....                           | 1-4        |
| Responsibility of the Manufacturer ..... | 1-5        |
| Intended Use .....                       | 1-5        |
| General .....                            | 1-5        |
| Equipment Symbols .....                  | 1-6        |
| <b>Service Information .....</b>         | <b>1-7</b> |
| Service Requirements .....               | 1-7        |
| Equipment Identification .....           | 1-7        |

## 2

|                                        |            |
|----------------------------------------|------------|
| <b>Equipment Overview .....</b>        | <b>2-1</b> |
| <b>Technical Characteristics .....</b> | <b>2-3</b> |
| General Description .....              | 2-3        |
| Power Requirements .....               | 2-4        |
| <b>Operation .....</b>                 | <b>2-5</b> |
| Operating Controls .....               | 2-5        |
| Leadwire Attachments .....             | 2-5        |
| Leadwire Adapters .....                | 2-6        |
| Lead Configurations .....              | 2-7        |

# 3

## Maintenance ..... 3-1

### Recommended Maintenance ..... 3-3

|                                 |     |
|---------------------------------|-----|
| General .....                   | 3-3 |
| Cleaning .....                  | 3-3 |
| Exterior Surfaces .....         | 3-3 |
| Electrodes .....                | 3-3 |
| Visual Checking .....           | 3-4 |
| Built-In Diagnostic Tests ..... | 3-4 |

### Domestic Electrical Safety Tests ..... 3-5

|                                                                     |     |
|---------------------------------------------------------------------|-----|
| AC Line Voltage Test .....                                          | 3-5 |
| 120 VAC, 50/60 Hz .....                                             | 3-5 |
| 240 VAC, 50/60 Hz .....                                             | 3-5 |
| Leakage Tests .....                                                 | 3-6 |
| Leakage Test Diagrams .....                                         | 3-7 |
| Test #1 Ground-wire-leakage-to-ground .....                         | 3-7 |
| Test #2 Chassis-leakage-to-ground .....                             | 3-7 |
| Test #3 Patient-cable-leakage-to-ground .....                       | 3-8 |
| Test #4 Patient-cable-leakage-into-patient Leads-from 120 VAC ..... | 3-8 |
| Ground Continuity .....                                             | 3-9 |

### Disassembly ..... 3-10

# 4

## Parts Lists and Drawings ..... 4-1

### Ordering Parts ..... 4-3

|                    |     |
|--------------------|-----|
| Introduction ..... | 4-3 |
|--------------------|-----|

### 900995-001A Acquisition Module Assembly ..... 4-4

### 420101-001B 14 Leadwire Kit ..... 4-7

### 417483-9XXA Leadwire, Multi-Link, Universal ..... 4-8

### 900179-201 Leadwire Adapter Kit - Banana ..... 4-10

### 900179-202 Leadwire Adapter Kit - Mactrode ..... 4-11

### 900179-203 Leadwire Adapter Kit - Grabber ..... 4-12

# 5

|                                            |            |
|--------------------------------------------|------------|
| <b>PCB Assemblies .....</b>                | <b>5-1</b> |
| 8001280-001A Data Acquisition Module ..... | 5-3        |
| SD8001280-001A Schematic, Main Board ..... | 5-8        |
| <b>Appendix A: Abbreviations .....</b>     | <b>A-1</b> |
| Standard Abbreviations .....               | A-3        |



# 1

## Introduction

|                                          |          |
|------------------------------------------|----------|
| <b>Manual Information .....</b>          | <b>3</b> |
| Revision History .....                   | 3        |
| Manual Purpose .....                     | 3        |
| Intended Audience .....                  | 3        |
| <b>Safety Information .....</b>          | <b>4</b> |
| Definitions .....                        | 4        |
| Messages .....                           | 4        |
| Responsibility of the Manufacturer ..... | 5        |
| Intended Use .....                       | 5        |
| General .....                            | 5        |
| Equipment Symbols .....                  | 6        |
| <b>Service Information .....</b>         | <b>7</b> |
| Service Requirements .....               | 7        |
| Equipment Identification .....           | 7        |



# Manual Information

## Revision History

Each page of the document has the document part number followed by a revision letter located at the bottom of the page. This letter identifies the document’s update level. The latest letter of the alphabet corresponds to the most current revision of the document.

The revision history of this document is summarized in the table below.

| Table 2-1. Revision History |                  |                  |
|-----------------------------|------------------|------------------|
| Revision                    | Date             | Comment          |
| A                           | 18 December 1998 | Initial release. |

## Manual Purpose

This manual supplies technical information for service representative and technical personnel so they can maintain the equipment to the assembly level. Use it as a guide for maintenance and electrical repairs considered field repairable. Where necessary the manual identifies additional sources of relevant information and or technical assistance.

See the host operator manual for instructions necessary to operate the equipment safely in accordance with its function and intended use.

## Intended Audience

This manual is intended for the person who uses, maintains, or troubleshoots this equipment.

# Safety Information

## Definitions



Indicates an imminently hazardous situation which, if not avoided, **WILL** result in death or serious injury.







Indicates a potentially hazardous situation which, if not avoided, **COULD** result in death or serious injury.



Indicates a potentially hazardous situation which, if not avoided may result in minor or moderate injury.

## Messages

Additional safety messages may be found throughout this manual that provide appropriate safe operation information.



## Responsibility of the Manufacturer

GE Marquette Medical Systems is responsible for the effects of safety, reliability, and performance only if:

- Assembly operations, extensions, readjustments, modifications, or repairs are carried out by persons authorized by Marquette.
- The electrical installation of the relevant room complies with the requirements of the appropriate regulations.
- The equipment is used in accordance with the instructions for use.

## Intended Use

The CAM-14 acquisition module is intended to acquire analog ECG signal, digitize it and transmit the signal to a host unit. The circuitry is designed to protect the host unit against the effects of cardiac defibrillator discharge to ensure recovery.

This device is intended for use under the direct supervision of a licensed health care practitioner.

## General

This equipment is protected against the effects of cardiac defibrillator discharge to ensure recovery, as required by test standards.

This equipment will not cause abnormal operation of the patient's cardiac pacemaker or other electrical stimulator.

This device uses a computerized ECG analysis program which can be used as a tool in ECG tracing interpretation. This computerized interpretation is only significant when used in conjunction with clinical findings. All computer-generated tracings should be overread by a qualified physician.

To ensure accuracy, use only computer-generated tracings and not the display for physician interpretation.

To ensure patient safety, use only parts and accessories manufactured or recommended by GE GE Marquette Medical Systems.

Contact GE Marquette Medical Systems for information before connecting any devices to this equipment that are not recommended in this manual.

If the installation of this equipment, in the USA, will use 240 V rather than 120 V, the source must be a center-tapped, 240 V, single-phase circuit.

Parts and accessories used must meet the requirements of the applicable IEC 601 series safety standards, and/or the system configuration must meet the requirements of the IEC 601-1-1 medical electrical systems standard.

The use of ACCESSORY equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of

safety of the resulting system. Consideration relating to the choice shall include:

- use of the accessory in the PATIENT VICINITY; and
- evidence that the safety certification of the ACCESSORY has been performed in accordance to the appropriate IEC 601-1 and/or IEC 601-1-1 harmonized national standard.

## Equipment Symbols

The following symbols appear on the equipment.



M13932

Type BF equipment. Type BF equipment is suitable for intentional external and internal application to the patient, excluding a direct conductive connection to the patient's heart. Type BF equipment has an F-type applied part. The paddles indicate that the equipment is defibrillator proof.

# Service Information


## Service Requirements

Refer equipment servicing to GE Marquette Medical Systems' authorized service personnel only. Any unauthorized attempt to repair equipment under warranty voids that warranty.

It is the user's responsibility to report the need for service to GE Marquette Medical Systems or to one of their authorized agents.

## Equipment Identification

Every GE Marquette Medical Systems device has a unique serial number for identification. The serial number appears on the product label on the base of each unit.



MD1113-022B

**Table 2-2. Equipment Identifications**

| Item | Name                    | Description                                                                                                                                                                                                                       |
|------|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| A    | name of device          | AM-114 Acquisition Module                                                                                                                                                                                                         |
| B    | manufacturer            | GE Marquette Medical Systems, Inc.                                                                                                                                                                                                |
| C    | serial number           | Unique identifier                                                                                                                                                                                                                 |
| D    | device characteristics  | One or two letters that further describe the unit, for example: P = prototype not conforming to marketing specification; R = refurbished equipment; S = special product documented under Specials part numbers; U = upgraded unit |
| E    | division                | F = Cardiology      G = Monitoring      J = GW Labs                                                                                                                                                                               |
| F    | product sequence number | Manufacturing number (of total units manufactured)                                                                                                                                                                                |
| G    | product code            | Two-character product descriptor      MF = CAM-14 Acquisition Module                                                                                                                                                              |
| H    | year manufactured       | 7 = 1997, 8 = 1998, 9 = 1999, (and so on)                                                                                                                                                                                         |
| I    | month manufactured      | A = January, B = February, C = March, D = April, E = May, F = June, G = July, H = August, J = September, K = October, L = November, M = December                                                                                  |



# 2

## Equipment Overview

|                                        |          |
|----------------------------------------|----------|
| <b>Technical Characteristics .....</b> | <b>3</b> |
| General Description .....              | 3        |
| Power Requirements .....               | 4        |
| <b>Operation .....</b>                 | <b>5</b> |
| Operating Controls .....               | 5        |
| Leadwire Attachments .....             | 5        |
| Leadwire Adapters .....                | 6        |
| Lead Configurations .....              | 7        |



# Technical Characteristics

## General Description

The acquisition module performs high resolution ECG data acquisition for use with host equipment (resting ECG analysis systems and exercise systems). The acquisition module has the following features:

- ac leadfail bias,
- lead off detection,
- 4000 Hz sampling rate,
- patient isolation,
- software updates via floppy diskette, and
- Hi-Res analysis on 14-lead units
- function key control of host equipment functions

The acquisition module provides patient electrical isolation for the host equipment. The minute ECG signals from the patient's skin are received by the electrodes and sent to the acquisition module via leadwires. The acquisition module then amplifies, digitizes, and performs some processing on the signals.

This whole process is controlled by a microprocessor in the acquisition module. The acquisition module then sends the serial, digitized, ECG data to the host equipment in 16 bit "words". The host equipment communicates to the acquisition module using this same serial line.

**Table 2-1. Safety**

| Item                                                  | Description                                                                                                                                                                                                                                                                                                                   |
|-------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Certification                                         | CE Marking for Council Directive 93/42/EEC                                                                                                                                                                                                                                                                                    |
| Type of Protection Against Electrical Shock           | Not applicable                                                                                                                                                                                                                                                                                                                |
| Degree of Protection Against Ingress of Liquids       | Ordinary                                                                                                                                                                                                                                                                                                                      |
| Handling of Disposable Supplies and Other Consumables | <ul style="list-style-type: none"> <li>■ Use only as manufactured or recommended by Marquette.</li> <li>■ Follow manufacturer's instructions for use for disposable/consumable product.</li> <li>■ Follow local environmental guidelines concerning the disposal of hazardous materials (e.g. lead acid batteries)</li> </ul> |
| Patient Mode of Operation                             | Continuous                                                                                                                                                                                                                                                                                                                    |
| Patient Leakage Current                               | less than 10μA                                                                                                                                                                                                                                                                                                                |
| Degree of Protection Against Electrical Shock         | type BF applied parts                                                                                                                                                                                                                                                                                                         |

**Table 2-1. Safety (Continued)**

| Item                  | Description                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Maintenance Frequency | <ul style="list-style-type: none"> <li>■ Recommended user daily visual inspection and cleaning.</li> <li>■ Recommended six-month routine maintenance checks and test procedures performed by qualified technical personnel.</li> </ul>                                                                                                                                                                                         |
| Repair Guidelines     | <p>Calibration instructions, equipment descriptions, and all other service information to repair those parts of the equipment designated as field repairable by qualified technical personnel is available in the service manual.</p> <p>Upon request, GE Marquette Medical Systems will provide circuit diagrams and component parts lists for printed circuit boards deemed repairable by qualified technical personnel.</p> |

**Table 2-2. Environmental**

| Item                   | Description                                       |
|------------------------|---------------------------------------------------|
| Operating Instructions |                                                   |
| Temperature            | 0°C to 50°C (32°F to 154°F)                       |
| Relative Humidity      | 20% to 95% noncondensing                          |
| Atmosphere Pressure    | 70 to 106 KPa (PRELIMINARY pending final testing) |
| Storage Conditions     |                                                   |
| Temperature            | –20°C to 60°C (–4°F to 172°F)                     |
| Relative Humidity      | 5% to 95% noncondensing                           |
| Atmosphere Pressure    | 50 to 106 KPa (PRELIMINARY pending final testing) |

## Power Requirements

The acquisition module draws its power, 12 V dc, from the host unit. (See the host equipment's field service manual for the host equipment's power requirements.) The acquisition module provides the patient electrical isolation for the host equipment.

# Operation

## Operating Controls


There are three controls on the acquisition module that are programmed by the host to record an ECG, record a rhythm, or stop the writer. See the host equipment operator's manual for the for specific information on using the acquisition module.



MD1320-002

## Leadwire Attachments


The leadwires are shown with 2 ends. The figure below shows how to use both ends of the leadwires. There are various options for leadwire adapters, the 4 mm pin, stress grabber, or MACTRODE clip. See chapter 5, "Parts Lists and Drawings" for part numbers.



MD1320-003B


### Leadwire Adapters

The acquisition module has universal leadwires which can be made into any lead of 3 basic types by using the adapters on the ends as shown in the figure below. The acquisition module leadwire adapters have a variety of configurations as shown below. (See chapter 4, “Parts Lists and Drawings” for part numbers for the various individual leadwire adapters.)



Refer to the figure on the next page for various labels that are affixed to the acquisition module for these configurations.

**Lead Configurations**      The various lead configurations are shown in the figure below. The acquisition module may be configured in any one of these ways. See chapter 5, “Parts Lists and Drawings” for part numbers.



MD1320-004



# 3

## Maintenance

|                                                                       |           |
|-----------------------------------------------------------------------|-----------|
| <b>Recommended Maintenance .....</b>                                  | <b>3</b>  |
| General .....                                                         | 3         |
| Cleaning .....                                                        | 3         |
| Exterior Surfaces .....                                               | 3         |
| Electrodes .....                                                      | 3         |
| Visual Checking .....                                                 | 4         |
| Built-In Diagnostic Tests .....                                       | 4         |
| <b>Domestic Electrical Safety Tests .....</b>                         | <b>5</b>  |
| AC Line Voltage Test .....                                            | 5         |
| 120 VAC, 50/60 Hz .....                                               | 5         |
| 240 VAC, 50/60 Hz .....                                               | 5         |
| Leakage Tests .....                                                   | 6         |
| Leakage Test Diagrams .....                                           | 7         |
| Test #1   Ground-wire-leakage-to-ground .....                         | 7         |
| Test #2   Chassis-leakage-to-ground .....                             | 7         |
| Test #3   Patient-cable-leakage-to-ground .....                       | 8         |
| Test #4   Patient-cable-leakage-into-patient Leads-from 120 VAC ..... | 8         |
| Ground Continuity .....                                               | 9         |
| <b>Disassembly .....</b>                                              | <b>10</b> |



## Recommended Maintenance

---

### General

Other than daily cleaning and occasional visual checking, the acquisition module requires no maintenance.

Only qualified service personnel should attempt repairing components and assemblies internal to the acquisition module. For service and repair, return the acquisition module to the factory for 48-hour turn around service. Call Tech Support for assistance. (See the “How to Reach Us” page in the front of the manual.)

#### NOTE

Unless you have an Equipment Maintenance Contract, GE Marquette Medical Systems does not in any manner assume the responsibility for performing the recommended maintenance procedures. The sole responsibility rests with the individual or institution using the equipment. GE Marquette Medical Systems service personnel may, at their discretion, follow the procedures provided in this manual as a guide during visits to the equipment site.

### Cleaning

#### Exterior Surfaces

Disconnect the acquisition module interface cable from the host equipment.

- Use a soft cloth moistened with water and a mild detergent. Wipe the exterior of the unit, the leadwires, and the acquisition module interface cable with the damp cloth. Dry all surfaces with a clean, soft cloth or paper towel.
- DO NOT allow any excess water to get inside the acquisition module or onto the leadwires or interface cable.
- Do not immerse acquisition module in water.
- Do not use alcohols, organic solvents, or abrasive cleaning agents.

#### Electrodes

- After each use, wipe reusable electrodes with a tissue or damp cloth to clean them of electrode paste. At the end of each day, wash reusable electrodes thoroughly with soap and water and dried.
- For suction electrodes, use a toothbrush to clean out the cups.

## Visual Checking

Inspect the acquisition module each time you clean it or if you suspect a problem.

- Check the leadwires and leadwire adapters for wear and loose connections. Replace these parts at the first sign of wear.
- Check the pins that the leadwires plug into. They should not be bent or loose. Contact Tech Support for any repairs needed. (See the “How to Reach Us” page in the front of the manual.)
- Check the acquisition module plastic case for any damage. Contact Tech Support for any repairs needed. (See the “How to Reach Us” page in the front of the manual.)

## Built-In Diagnostic Tests

The host equipment generally contains built-in diagnostic tests to verify the operation of the electrocardiograph.

These built-in diagnostic tests verify the operation of the acquisition module, as well. For example, there is a test that records raw ECG data on the thermal paper. This test checks for noise and gain in the acquisition module. Another test is a serial link test that determines if the microprocessor in the host equipment is communicating with the acquisition module.

For details on using these tests, see the field service manual for the host equipment.

# Domestic Electrical Safety Tests


## AC Line Voltage Test

This test verifies that the domestic wall outlet supplying power to the equipment is properly wired. For international wiring tests, refer to the internal standards agencies of that particular country.

120 VAC, 50/60 Hz

Use a digital voltmeter to check the voltages of the 120-volt AC wall outlet (dedicated circuit recommended). If the measurements are significantly out of range, have a qualified electrician repair the outlet. The voltage measurements should be as follows:

1. 120 VAC ( $\pm 10$  VAC) between the line contact and neutral and between the line contact and ground.
2. Less than 3 VAC between neutral and ground.




MD1287-006

240 VAC, 50/60 Hz

Use a digital voltmeter, set to measure at least 300 VAC, to check the voltages of the NEMA 6-20R, AC wall outlet (dedicated circuit recommended). If the measurements are significantly out of range, have a qualified electrician repair the outlet. The voltage measurements should be as follows:

1. 120 VAC ( $\pm 10$  VAC) between either “hot” contact and ground.
2. 210 to 230 VAC between the two “hot” contacts.



MD1287-007

## Leakage Tests

The leakage tests are safety tests to ensure that the equipment poses no electrical health hazards. Use the table below to determine which tests apply to the unit under test and the maximum allowable leakage currents. For international leakage limits, refer to the internal standards agencies of that particular country.

If the unit under test fails the leakage tests, do not allow the customer to use the equipment. Call Tech Support for assistance. (See the “How to Reach Us” page in the front of the manual.)

GE Marquette Medical Systems recommends that you perform these tests:

- Before applying power for the first time
- Every 6 months as part of routine maintenance
- Whenever internal assemblies are serviced

You need a leakage tester to perform the leakage tests.

### NOTE

The accuracy of the leakage tests depends on a properly-wired wall outlet. Do not proceed until you verify the integrity of the power source.

### ⚠ WARNING



**Total system leakage current must not exceed 100 microamperes.**

M15287-9D

**Table 3-1. Leakage Tests and Maximum Allowable Leakage Currents**


| Test                                                     | Applies To                         | Maximum Current ( $\mu$ A) |
|----------------------------------------------------------|------------------------------------|----------------------------|
| 1 Ground-wire-leakage-to-ground                          | Host equipment                     | 100                        |
| 2 Chassis-leakage-to-ground                              | Host equipment                     | 100                        |
| 3 Patient-cable-leakage-to-ground                        | Acquisition module (patient cable) | 10                         |
| 4 Patient-cable-leakage-into-patient-leads-from-120 V ac | Acquisition module (patient cable) | 20                         |

# Leakage Test Diagrams

These diagrams show only a representation of how a typical leakage current tester functions. Follow the instructions provided with the leakage current tester that you use.

## Test #1

Ground-wire-leakage-to-ground




MD1287-008

## Test #2

Chassis-leakage-to-ground


Apply line voltage to the UUT chassis for this test.



MD1287-009

## Test #3

## Patient-cable-leakage-to-ground




MD1287-010

# Test #4

Patient-cable-leakage-into-patient  
Leads-from 120 VAC

During this test, line voltage is applied to the patient cable connectors. To prevent erroneous readings, do not allow the leadwires to contact conductive materials such as metal handles, and do not place the leadwires on the floor.



MD1287-011

## Ground Continuity

This test verifies that there is continuity (less than 100 m $\Omega$  resistance) between all the exposed metal surfaces on the host equipment, which have the potential to become energized, and the ground prong on the mains AC power cord. If the metal surfaces are anodized or painted, scrape off a small area in an inconspicuous area for the probe to make direct contact with the metal.

Use a digital multimeter to check all the metal surfaces of the host equipment. Make adjustments for any resistance in the test leads. (See the host equipment field service manual for detailed information.)


If the measurements are significantly out of range, check for breaks in the power cord or in the internal connections within the unit.

# Disassembly

Strict antistatic precautions should be followed during disassembly/assembly of the CAM-14.

## NOTE

Use the proper ESD grounding techniques when handling components. Wear an antistatic wrist strap and an ESD protected mat. Store ESD sensitive components in antistatic bags before placing them on any surface.



MD1320-006

Follow the steps below when disassembling the CAM-14.

1. Remove the four Torx screws holding the top and bottom of the CAM-14 unit together using a #6 TORX driver.
2. Lift the top only high enough to allow access to the switch flex cable from the switch assembly.

3. To avoid damage to the flex cable, gently pry the cable off the switch at the main pcb with a small screwdriver.

**NOTE**

Damage to the switch flex cable connector will result in failure of one or more of the switch functions.

- **DO NOT** attempt to disconnect the switch flex cable by pulling on it.
  - **DO NOT** let the main pcb hang by the cable.
4. When the main pcb has been separated from the switch flex cable, disconnect the host cable interface flex from the main pcb by gently prying the connectors apart.



# 4

## Parts Lists and Drawings

|                                                             |           |
|-------------------------------------------------------------|-----------|
| <b>Ordering Parts .....</b>                                 | <b>3</b>  |
| Introduction .....                                          | 3         |
| <b>900995-001A    Acquisition Module Assembly .....</b>     | <b>4</b>  |
| <b>420101-001B    14 Leadwire Kit .....</b>                 | <b>7</b>  |
| <b>417483-9XXA    Leadwire, Multi-Link, Universal .....</b> | <b>8</b>  |
| <b>900179-201    Leadwire Adapter Kit - Banana .....</b>    | <b>10</b> |
| <b>900179-202    Leadwire Adapter Kit - Mactrode .....</b>  | <b>11</b> |
| <b>900179-203    Leadwire Adapter Kit - Grabber .....</b>   | <b>12</b> |



# Ordering Parts

---

## Introduction

This chapter provides upper-level drawings for the standard configurations of the acquisition module and any kits which are available. It includes enough detail to provide part numbers for field-serviceable assemblies in the equipment.

Subassembly drawings follow the upper-level drawings and are arranged in ascending numerical order. The assembly drawings generally include both a parts lists and an exploded view.


To order parts, contact Service Parts at the address or telephone number on the “How to Reach Us...” page at the beginning of this manual.

## 900995-001A

## Acquisition Module Assembly


| Item | Description                          | Part Number | Qty. |
|------|--------------------------------------|-------------|------|
| 1    | PCB CAM 14 MAIN BRD                  | 801280-001  | 1    |
| 2    | PCB CAM14 CABLE INTFC                | 801502-001  | 1    |
| 3    | ASSY TOP W/SWITCH CAM 14             | 422032-001  | 1    |
| 4    | ASSY WELD BTM BELT CLIP CAM 14       | 420148-001  | 1    |
| 5    | LABEL NEW ACQ MOD LD CODES           | 420152-001  | 1    |
| 6    | KIT LDWR AM11X SET 14                | 420101-001  | 1    |
| 7    | LABEL CLEAR OVERLAMINATE             | 413608-001  | 1    |
| 8    | LABEL CAM 14 BOTTOM                  | 419979-001  | 1    |
| 9    | SCREW 1 X .375 TORX T-6 FHPH         | 417866-002  | 4    |
| 10   | CARTON, MAILER 11X5X3 CAM 14         | 421927-001  | 1    |
| 11   | INSERT, DIE CUT CAM 14               | 421927-002  | 1    |
| 12   | I-STATIC, 6.00W x 8.00 L (not shown) | 9976-005    | 1    |
| 13   | PTR 14 SET AHA BANANA                | 900179-201  | 1    |

1



Notes:

1. Model/Serial # Label to be Marked:  
CAM 14  
(Serial # TBD)




420101-001B

14 Leadwire Kit

---



| Item | Description                   | Part Number | Qty |
|------|-------------------------------|-------------|-----|
| 1    | BAG, ANTISTATIC 6 X 8         | 9976-005    | 1   |
| 2    | LEADWIRE UNIV/MLINK AM11X 26" | 417483-905  | 10  |
| 3    | LEADWIRE UNIV/MLINK AM11X 36" | 417483-906  | 2   |
| 4    | LEADWIRE UNIV/MLINK AM11X 40" | 417483-903  | 2   |
| 5    | INSERT/REORDER CARD           | 70357-001   | 1   |




## 417483-9XXA

## Leadwire, Multi-Link, Universal

| Item | Description                          | Part Number | Qty  |
|------|--------------------------------------|-------------|------|
| 1    | TUBING FIT-221 3/32 BLK              | 4882-103    | 0.04 |
| 2    | RESIN ABS NATURAL                    | 5421-001    | A/R  |
| 3    | RESIN PVC 60 DURO NATURAL            | 5423-001    | A/R  |
| 4    | RESIN PVC 85 DURO NATURAL            | 5423-002    | A/R  |
| 5    | MOLD UNIVERSAL PREMOLD AM4           | 58973-058   | 1    |
| 6    | MOLD STRAIN RELIEF AM4 WIRE          | 58973-059   | 1    |
| 7    | INSERT MOLDED AM4 LEADWIRE           | 58973-060   | 1    |
| 8    | LDWR PLUGFACE INDV MULTILINK DGRAY   | 411191-004  | 1    |
| 9    | PLUG PREMOLD INDV LDWR MULTILINK     | 411192-900  | 1    |
| 10   | PLUG OVERMOLD INDV LDWR MULTILINK    | 411193-001  | 1    |
| 11   | LEADWIRE CONTACT                     | 411197-001  | 3    |
| 12   | BAND SPLICE                          | 412201-001  | 1    |
| 13   | COLORANT MULTIPURPOSE DARK GRAY      | 412793-028  | A/R  |
| 14   | CABLE COAX 1 COND 26 AWG. 125 OD PVC | 700136-007  | A/R  |



| LENGTH CHART |        |                                   |
|--------------|--------|-----------------------------------|
| REV          | SUFFIX | ITEM 16<br>DIM "A"<br>INCHES/FEET |
| A            | -901   | 20"/1.70 FT                       |
| A            | -902   | 29"/2.45 FT                       |
| A            | -903   | 40"/3.35 FT                       |
| A            | -904   | 51"/4.25 FT                       |
| A            | -905   | 26"/2.17 FT                       |
| A            | -906   | 36"/3.0 FT                        |




Notes:

1. For Overmold Mix PVC Resins Items 3 & 4 together to obtain  $72 \pm 5$  Durometer.  
Color: Gray, Munsell N7.
2. Month Lot Code: One letter in alphabetical order for each month where A = January, B = February . . . Skip letter "I".  
Year Lot Code: One digit, the last digit of the year where 7 = 1997, 8 = 1998 . . .
3. Performance Specification:  
Hi-Pot: Wire must withstand 5000 VDC 1 mA for 1 second between center Conductor and Shield.  
Continuity Test: Leadwire to be tested end to end between universal End & Center Conductor.
4. Finished Assembly to be free of dirt. Discolorations or other signs of poor workmanship. Molded Components to have Flash and Gate Marks trimmed flush to within .010 of surface.

## 900179-201


## Leadwire Adapter Kit - Banana

| Item | Description           | Part Number | Qty |
|------|-----------------------|-------------|-----|
| 1    | ASSY BANNANA RL GREEN | 406554-007  | 1   |
| 2    | ASSY BANNANA RA WHITE | 406554-001  | 1   |
| 3    | ASSY BANANA LL RED    | 406554-005  | 1   |
| 4    | ASSY BANANA LA BLACK  | 406554-003  | 1   |
| 5    | ASSY BANANA V1 BROWN  | 406554-009  | 1   |
| 6    | ASSY BANANA V2 BROWN  | 406554-011  | 1   |
| 7    | ASSY BANANA V3 BROWN  | 406554-013  | 1   |
| 8    | ASSY BANANA V4 BROWN  | 406554-015  | 1   |
| 9    | ASY BANANA V5 BROWN   | 406554-017  | 1   |
| 10   | ASSY BANANA V6 BROWN  | 406554-019  | 1   |
| 11   | ASSY BANANA E ORANGE  | 406554-021  | 1   |
| 12   | ASSY BANANA H ORANGE  | 406554-023  | 1   |
| 13   | ASSY BANANA I ORANG   | 406554-025  | 1   |
| 14   | ASSY BANANA M ORANGE  | 406554-027  | 1   |



**900179-202****Leadwire Adapter Kit - Mactrode**


| Item | Description              | Part Number | Qty |
|------|--------------------------|-------------|-----|
| 1    | ASSY MACTRODE RL GREEN   | 406551-007  | 1   |
| 2    | ASSY MACTRODE RA WHITE   | 406551-001  | 1   |
| 3    | ASSY MACTRODE LL RED     | 406551-005  | 1   |
| 4    | ASSY MACTRODE LA BLACK   | 406551-003  | 1   |
| 5    | ASSY MACTRODE V1 BROWN   | 406551-009  | 1   |
| 6    | ASSY MACTRODE V2 BROWN   | 406551-011  | 1   |
| 7    | ASSY MACTRODE V3 BROWN   | 406551-013  | 1   |
| 8    | ASSY MACTRODE V4 BROWN   | 406551-015  | 1   |
| 9    | ASY MACTRODE V5 BROWN    | 406551-017  | 1   |
| 10   | ASSY MACTRODE V6 BROWN   | 406551-019  | 1   |
| 11   | ASSY MACTRODE E ORANGE   | 406551-021  | 1   |
| 12   | ASSY MACTRODE H ORANGE   | 406551-023  | 1   |
| 13   | ASSY MACTRODE I ORANG    | 406551-025  | 1   |
| 14   | ASSY v MACTRODE M ORANGE | 406551-027  | 1   |



## 900179-203

## Leadwire Adapter Kit - Grabber

| Item | Description             | Part Number | Qty |
|------|-------------------------|-------------|-----|
| 1    | ASSY GRABBER RL GREEN   | 406552-007  | 1   |
| 2    | ASSY GRABBER RA WHITE   | 406552-001  | 1   |
| 3    | ASSY GRABBER LL RED     | 406552-005  | 1   |
| 4    | ASSY GRABBER LA BLACK   | 406552-003  | 1   |
| 5    | ASSY GRABBER V1 BROWN   | 406552-009  | 1   |
| 6    | ASSY GRABBER V2 BROWN   | 406552-011  | 1   |
| 7    | ASSY GRABBER V3 BROWN   | 406552-013  | 1   |
| 8    | ASSY GRABBER V4 BROWN   | 406552-015  | 1   |
| 9    | ASSY GRABBER V5 BROWN   | 406552-017  | 1   |
| 10   | ASSY GRABBER V6 BROWN   | 406552-019  | 1   |
| 11   | ASSY GRABBER E ORANGE   | 406552-021  | 1   |
| 12   | ASSY GRABBER H ORANGE   | 406552-023  | 1   |
| 13   | ASSY GRABBER I ORANG    | 406552-025  | 1   |
| 14   | ASSY v GRABBER M ORANGE | 406552-027  | 1   |



# 5

## PCB Assemblies

|                |                               |   |
|----------------|-------------------------------|---|
| 8001280-001A   | Data Acquisition Module ..... | 3 |
| SD8001280-001A | Schematic, Main Board .....   | 8 |



**8001280-001A****Data Acquisition Module**


| Item | Description                    | Part Number | Qty. |
|------|--------------------------------|-------------|------|
| 1    | RES COMP 100M 5% 1/4W          | 1001-107    | 2    |
| 2    | RES SM CER 100 5% 1/8W         | 1081-101    | 1    |
| 3    | RES SM CHIP 5% 1/8W 51 OHMS    | 1081-472    | 1    |
| 4    | RES SM CER 26.7 1% 1/8W        | 1082-909    | 1    |
| 5    | CAP SM CER COG 150PF 5% 50V    | 1181-151    | 13   |
| 6    | CAP SM CER X7R .1UF 50V        | 1187-104    | 2    |
| 7    | DIODE SM SERIES PR D7000       | 2013-201    | 4    |
| 8    | DIODE SM SCHTKY BARRIER 2800   | 2412-001    | 1    |
| 9    | DIODE SM LL FDSO1503           | 2414-001    | 14   |
| 10   | IC SM HC 74HC161               | 3038-161    | 1    |
| 11   | IC SM SSOP VHC04               | 3057-004    | 1    |
| 12   | IC SM SSOP VHC573              | 3057-573    | 1    |
| 13   | DIODE SM ZENER 5234 6.2V       | 401502-001  | 2    |
| 14   | CAP SM TANT 1.0UF 10% 16V      | 406883-002  | 3    |
| 15   | CAP SM TANT 10UF 10% 35V       | 406884-013  | 3    |
| 16   | IC SM DUAL N-CHAN MFET SI9955  | 407547-001  | 1    |
| 17   | BEAD SM 1206 600@100 200MA     | 408746-001  | 3    |
| 18   | IC SM HCMOS TC7S02F            | 408795-001  | 2    |
| 19   | IC SM OPTOCOUPLER CNW2611      | 408825-001  | 2    |
| 20   | IC SM HCMOS AND TC7S08F        | 409035-001  | 2    |
| 21   | RES SM 0603 100 1% 1/16W       | 410334-003  | 1    |
| 22   | RES SM CER 0603 10K 1% 1/16W   | 410334-013  | 22   |
| 23   | RES SM 0603 11.5K 1% 1/16W     | 410334-014  | 1    |
| 24   | RES SM 0603 3.01K 1% 1/16W     | 410334-018  | 13   |
| 25   | RES SM 0603 100K 1% 1/16W      | 410334-019  | 7    |
| 26   | RES SM CER 0603 0 OHM JUMPER   | 410334-027  | 1    |
| 27   | RES SM 0603 511 OHM 1% 1/16W   | 410334-034  | 1    |
| 28   | RES SM 0603 47.5K 1% 1/16W     | 410334-036  | 1    |
| 29   | RES SM 0603 200 OHM 1% 1/16W   | 410334-044  | 2    |
| 30   | RES SM 0603 7.68K 1% 1/16W     | 410334-049  | 1    |
| 31   | RES SM 0603 20.0K 1% 1/16W     | 410334-053  | 6    |
| 32   | RES SM CER 0603 4.99K 1% 1/16W | 410334-066  | 1    |
| 33   | RES SM CER 0603 35.7K 1% 1/16W | 410334-091  | 2    |
| 34   | RES SM 0603 7.87K 1%           | 410334-098  | 3    |
| 35   | RES SM 0603 750 1%             | 410334-099  | 3    |
| 36   | RES SM 0603 1.24K 1%           | 410334-101  | 1    |
| 37   | RES SM 0603 30.1 1% 1/16W      | 410334-109  | 1    |
| 38   | RES SM 0603 2.49K 1% 1/16W     | 410334-110  | 7    |
| 39   | RES SM CER 0603 8.06K 1% 1/16W | 410334-113  | 2    |

| Item | Description                    | Part Number | Qty. |
|------|--------------------------------|-------------|------|
| 40   | RES SM 0603 1.0M 1% 1/16W      | 410334-120  | 3    |
| 41   | RES SM 0603 33.2K 1% 1/16W     | 410334-165  | 1    |
| 42   | RES SM 0603 249K 1% 1/16W      | 410334-186  | 2    |
| 43   | RES SM 0603 15.8K 1% 1/16      | 410334-225  | 12   |
| 44   | CONN BLK,CAM 14 14POS W/10K R  | 411305-003  | 1    |
| 45   | CAP SM X7R 0603 .1UF 10% 16V   | 411575-002  | 80   |
| 46   | CAP SM X7R 0603 1000PF 5% 50V  | 411575-003  | 2    |
| 47   | CAP SM X7R 0603 390PF 5% 50V   | 411575-009  | 1    |
| 48   | CAP SM X7R 0603 2700PF 5% 50V  | 411575-010  | 2    |
| 49   | CAP SM X7R 0603 0.01UF 5% 50V  | 411575-012  | 2    |
| 50   | CAP SM X7R 0603 .022UF 16V 5%  | 411575-014  | 4    |
| 51   | CAP SM X7R 0603 220PF 5% 50V   | 411575-022  | 4    |
| 52   | CAP SM X7R 0603 3300PF 5% 50V  | 411575-024  | 1    |
| 53   | CAP SM NPO 0603 100PF 5% 50V   | 411575-025  | 2    |
| 54   | CAP SM X7R 0603 2200PF-5%      | 411575-033  | 3    |
| 55   | CAP SM NPO 0603 33PF 5% 50V    | 411576-004  | 1    |
| 56   | CAP SM NPO 0603 5.6PF +/-25PF  | 411576-017  | 13   |
| 57   | CAP SM NPO 0603 470PF 5% 50V   | 411576-018  | 3    |
| 58   | CAP SM X7R 0805 .047UF 5% 25V  | 411587-004  | 1    |
| 59   | IC SM DUAL OP AMP TLC272CD     | 411710-001  | 1    |
| 60   | RES SM 0603 10M 5% 1/16W       | 411723-003  | 1    |
| 61   | SM RESONATOR 4MHZ TYPE KBR     | 412010-001  | 1    |
| 62   | IC SM OPTOCOUPLER CNW136       | 412017-001  | 1    |
| 63   | CAP SM TANT 220UF 20% 10V      | 412662-002  | 2    |
| 64   | CAP SM TANT 33UF 20% 25V       | 412662-003  | 1    |
| 65   | IC SM OP AMP AD822             | 412753-001  | 1    |
| 66   | IC SM OP AMP LOW POWER AD820   | 412799-001  | 1    |
| 67   | CAP SM TANT 3.3UF 10% 10V      | 413119-001  | 1    |
| 68   | DIODE SM SCHTKY RECT MBR0520   | 413970-001  | 2    |
| 69   | SM FERRITE BEAD 0603 BLM11A601 | 414061-001  | 5    |
| 70   | CAP SM TANT 2.2UF 10V 10%      | 414084-001  | 2    |
| 71   | IC SM PWM CONTROLLER UCC3806DW | 414247-001  | 1    |
| 72   | IC SM 4.85V REG MIC5201        | 416998-001  | 1    |
| 73   | IC RS485/422 XCVR MAX490 SO8   | 417913-001  | 1    |
| 74   | CONN SCKT SM B/B 2MM HGT 10POS | 417955-001  | 2    |
| 75   | TRANSFORMER TNI ISO T8175      | 418434-001  | 1    |
| 76   | IC SM DSP TMS320F206 PQFP      | 419051-001  | 1    |
| 77   | IC SM D FLIP TC7W74FU          | 419131-001  | 1    |
| 78   | IC SM AV9170-02                | 419133-001  | 1    |
| 79   | IC SM OP AMP TLV2221 SOT23-5P  | 419134-001  | 13   |
| 80   | IC SM TPS 5904                 | 419135-001  | 1    |
| 81   | IC SM MAX4051                  | 419283-001  | 2    |

| Item | Description                   | Part Number  | Qty. |
|------|-------------------------------|--------------|------|
| 82   | IC SM INSTR AMP INA126E SSOP  | 419315-002   | 12   |
| 83   | IC SM QUAD OP-AMP 34184 TSSOP | 419383-001   | 1    |
| 84   | RES SM CHIP 0.2 5% 0.2W       | 419604-001   | 1    |
| 85   | IC SM PWR SPLY MONMAX823      | 419699-002   | 1    |
| 86   | IC SM ADS7809 16 BIT A/D      | 419952-001   | 1    |
| 87   | POST SM 1MM MALE 2X7          | 420007-007   | 1    |
| 88   | FRMWR COMMAND INTERFACE V001A | 421920-001   | 1    |
| 89   | CKT BD CAM 14 MAIN BD         | 801281-001   | 1    |
| 90   | SCHEM CAM14 MAIN BRD          | SD801280-001 | 0    |

## Drawing notes:

- Unless otherwise specified, all resistors are 10K, 1%; all capacitors are 0.1 $\mu$ F, 10%; all diodes are FDS01503; all ICs are INA126E.
- The bar shown on all polarized capacitors denotes the positive terminal.
- Parts not installed (solder side): R29, R30.
- Trim leads of R224 and R225 in the indicated area so they do not extend beyond edge of solder pad.





Component (Top) View




**SD8001280-001A****Schematic, Main Board**


1 of 13







2 of 13






3 of 13








4 of 13






5 of 13







|      |          |      |      |    |   |
|------|----------|------|------|----|---|
| U47  | 74VHC573 | 20   | +5VD | 10 | D |
| C171 | 0.1      | +5VD | D    |    |   |


6 of 13






## 7 of 13






8 of 13










10 of 13



D ground plane tie to A ground plane here

|     |        |   |      |   |      |
|-----|--------|---|------|---|------|
| U45 | TLC272 | 8 | +5VA | 4 | -5VA |
| BC2 | 0.1    | 2 | +5VA | 1 | A    |
| BC1 | 0.1    | 2 | -5VA | 1 | A    |

**11 of 13**

|     |         |      |  |      |  |
|-----|---------|------|--|------|--|
| U1  | INA126E | 4    |  | -5VA |  |
|     |         | 7    |  | +5VA |  |
| BC3 | 0.1     | +5VA |  | A    |  |
| BC4 | 0.1     | -5VA |  | A    |  |

|     |         |      |  |      |  |
|-----|---------|------|--|------|--|
| U2  | INA126E | 4    |  | -5VA |  |
|     |         | 7    |  | +5VA |  |
| BC5 | 0.1     | +5VA |  | A    |  |
| BC6 | 0.1     | -5VA |  | A    |  |

|     |         |      |  |      |  |
|-----|---------|------|--|------|--|
| U3  | INA126E | 4    |  | -5VA |  |
|     |         | 7    |  | +5VA |  |
| BC7 | 0.1     | +5VA |  | A    |  |
| BC8 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U4   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC9  | 0.1     | +5VA |  | A    |  |
| BC10 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U5   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC11 | 0.1     | +5VA |  | A    |  |
| BC12 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U6   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC13 | 0.1     | +5VA |  | A    |  |
| BC14 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U7   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC15 | 0.1     | +5VA |  | A    |  |
| BC16 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U8   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC17 | 0.1     | +5VA |  | A    |  |
| BC18 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U9   | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC19 | 0.1     | +5VA |  | A    |  |
| BC20 | 0.1     | -5VA |  | A    |  |

|      |         |      |  |      |  |
|------|---------|------|--|------|--|
| U10  | INA126E | 4    |  | -5VA |  |
|      |         | 7    |  | +5VA |  |
| BC21 | 0.1     | +5VA |  | A    |  |
| BC22 | 0.1     | -5VA |  | A    |  |

# PCB Assemblies: SD8001280-001A Schematic, Main Board

|      |         |      |      |
|------|---------|------|------|
| U11  | INA126E | 4    | -5VA |
|      |         | 7    | +5VA |
| BC25 | 0.1     | +5VA | A    |
| BC26 | 0.1     | -5VA | A    |

|      |         |      |      |
|------|---------|------|------|
| U12  | INA126E | 4    | -5VA |
|      |         | 7    | +5VA |
| BC27 | 0.1     | +5VA | A    |
| BC28 | 0.1     | -5VA | A    |

|      |         |      |      |   |      |
|------|---------|------|------|---|------|
| U13  | TLV2221 | 5    | +5VA | 2 | -5VA |
| BC45 | 0.1     | -5VA | A    |   |      |
| BC46 | 0.1     | +5VA | A    |   |      |

|      |       |      |      |   |      |
|------|-------|------|------|---|------|
| U17  | AD820 | 7    | +5VA | 4 | -5VA |
| BC47 | 0.1   | +5VA | A    |   |      |
| BC48 | 0.1   | -5VA | A    |   |      |

|      |          |      |      |   |   |
|------|----------|------|------|---|---|
| U41  | TC7W74FU | 8    | +5VD | 4 | D |
| BC49 | 0.1      | +5VD | D    |   |   |

|      |           |      |      |   |   |
|------|-----------|------|------|---|---|
| U38  | AV9170-02 | 7    | +5VD | 3 | D |
| BC50 | 0.1       | +5VD | D    |   |   |

|     |            |    |    |    |    |      |    |    |    |    |   |
|-----|------------|----|----|----|----|------|----|----|----|----|---|
| U19 | TMS320F206 | 11 | 16 | 4  | 7  | +5VD | 21 | 25 | 30 | 14 | D |
|     |            | 35 | 50 | 63 | 75 | +5VD | 37 | 42 | 48 | 54 | D |
|     |            | 91 |    |    |    | +5VD | 59 | 65 | 70 | 83 | D |
|     |            |    |    |    |    |      | 88 | 94 |    |    | D |

|      |     |      |   |      |     |      |   |
|------|-----|------|---|------|-----|------|---|
| BC51 | 0.1 | +5VD | D | BC58 | 0.1 | +5VD | D |
| BC52 | 0.1 | +5VD | D | BC59 | 0.1 | +5VD | D |
| BC53 | 0.1 | +5VD | D |      |     |      |   |
| BC54 | 0.1 | +5VD | D |      |     |      |   |
| BC55 | 0.1 | +5VD | D |      |     |      |   |
| BC56 | 0.1 | +5VD | D |      |     |      |   |
| BC57 | 0.1 | +5VD | D |      |     |      |   |

12 of 13

|      |         |      |    |      |    |   |
|------|---------|------|----|------|----|---|
| U18  | ADS7809 |      | 19 | +5VA | 10 | D |
|      |         |      | 20 | +5VD |    |   |
| BC38 | 0.1     | +5VA | A  |      |    |   |
| BC39 | 0.1     | +5VD | D  |      |    |   |

|      |         |      |    |   |      |  |   |  |   |
|------|---------|------|----|---|------|--|---|--|---|
| U20  | 74VHC04 |      | 14 |   | +5VD |  | 7 |  | D |
| BC37 | 0.1     | +5VD |    | D |      |  |   |  |   |

|      |         |      |   |   |      |  |   |  |   |
|------|---------|------|---|---|------|--|---|--|---|
| U21  | TC7S08F |      | 5 |   | +5VD |  | 3 |  | D |
| BC24 | 0.1     | +5VD |   | D |      |  |   |  |   |

|      |         |      |   |   |      |  |   |  |   |
|------|---------|------|---|---|------|--|---|--|---|
| U43  | TC7S02F |      | 5 |   | +5VD |  | 3 |  | D |
| BC34 | 0.1     | +5VD |   | D |      |  |   |  |   |

|      |         |      |   |   |      |  |   |  |   |
|------|---------|------|---|---|------|--|---|--|---|
| U46  | TC7S02F |      | 5 |   | +5VD |  | 3 |  | D |
| BC23 | 0.1     | +5VD |   | D |      |  |   |  |   |

|      |         |      |      |   |      |
|------|---------|------|------|---|------|
| U14  | MC34184 | 11   | -5VA | 4 | +5VA |
| BC36 | 0.1     | -5VA | A    |   |      |
| BC35 | 0.1     | +5VA | A    |   |      |

|      |          |      |      |   |   |
|------|----------|------|------|---|---|
| U15  | MAX 4051 | 16   | +5VA | 8 | A |
|      |          | 7    | -5VA |   |   |
| BC29 | 0.1      | +5VA | A    |   |   |
| BC30 | 0.1      | -5VA | A    |   |   |

|      |          |      |      |   |   |
|------|----------|------|------|---|---|
| U16  | MAX 4051 | 16   | +5VA | 8 | A |
|      |          | 7    | -5VA |   |   |
| BC31 | 0.1      | +5VA | A    |   |   |
| BC32 | 0.1      | -5VA | A    |   |   |

13 of 13

|      |           |    |     |    |     |
|------|-----------|----|-----|----|-----|
| U36  | 68HC705P9 | 28 | +5  | 14 | GND |
| BC40 | 0.1       | +5 | GND |    |     |

|      |       |    |     |   |     |
|------|-------|----|-----|---|-----|
| U44  | AD822 | 8  | +5  | 4 | GND |
| BC41 | 0.1   | +5 | GND |   |     |

|      |         |    |     |   |     |
|------|---------|----|-----|---|-----|
| U42  | 74HC161 | 16 | +5  | 8 | GND |
| BC42 | 0.1     | +5 | GND |   |     |

|      |        |    |     |   |     |
|------|--------|----|-----|---|-----|
| U34  | MAX490 | 1  | +5  | 4 | GND |
| BC43 | 0.1    | +5 | GND |   |     |

|      |         |    |     |   |     |
|------|---------|----|-----|---|-----|
| U35  | TC7S08F | 5  | +5  | 3 | GND |
| BC44 | 0.1     | +5 | GND |   |     |

|      |        |    |     |   |     |
|------|--------|----|-----|---|-----|
| U37  | MAX823 | 5  | +5  | 2 | GND |
| BC33 | 0.1    | +5 | GND |   |     |

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U48 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U49 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U50 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U51 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U52 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U53 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U54 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U55 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U56 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U57 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U58 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

|     |         |   |       |   |       |
|-----|---------|---|-------|---|-------|
| U59 | TLV2221 | 5 | +5V A | 2 | -5V A |
|-----|---------|---|-------|---|-------|

# Appendix A: Abbreviations

|                              |   |
|------------------------------|---|
| Standard Abbreviations ..... | 3 |
|------------------------------|---|



\_\_\_\_\_

|            |                                         |     |
|------------|-----------------------------------------|-----|
| Revision A | CAM-14 Acquisition Module<br>421315-001 | A-3 |
|------------|-----------------------------------------|-----|

|               |                                                   |                |                                           |
|---------------|---------------------------------------------------|----------------|-------------------------------------------|
| CRC           | cyclic redundancy check                           | EPLD           | electrically programmable logic device    |
| CRD           | cord                                              | EPROM          | erasable, programmable, read-only memory  |
| crt, CRT      | cathode ray tube                                  |                |                                           |
| CSA           | Canadian Standards Association                    | ESD            | electrostatic discharge                   |
| CTRL          | control                                           | etc, etc.      | et cetera                                 |
|               | <b>D</b>                                          | EURO           | Europe, European                          |
|               |                                                   | EXP            | Expanded                                  |
| D/A           | digital to analog                                 |                | <b>F</b>                                  |
| DA            | damping relay                                     |                |                                           |
| dac, DAC      | digital-to-analog converter                       | F              | fuse, Farad, female                       |
| DAN           | Danish                                            | F1-F5          | function keys 1 through 5                 |
| Dat/Tim       | date/time                                         | Fax            | facsimile                                 |
| dBm           | decibel (referenced to 1 milliwatt into 600 ohms) | FCC            | Federal Communications Commission         |
| dc, DC        | direct current                                    | FE             | front end                                 |
| DD            | double density, day                               | FILH           | fillister head                            |
| DDD           | Digital Diagnostic Diskette                       | FLH            | flat head                                 |
| DEC           | Digital Equipment Corporation, December           | FLRAM          | flash RAM                                 |
| Del           | delete                                            | FR             | French                                    |
| DEMO          | demonstration                                     | FrntEnd        | front end                                 |
| DES           | designation                                       | FSK            | frequency shift keying                    |
| DevId         | device identification                             | ft             | foot, feet                                |
| Diag          | diagnostic                                        | Furosem        | Furosemide                                |
| Digital       | Digitalis                                         |                | <b>G</b>                                  |
| Digitox       | Digitoxin                                         | g              | gram, acceleration due to gravity         |
| Digox         | digoxin                                           | GB             | Great Britain                             |
| Digoxin       | Digoxin-Lanoxin                                   | GERM           | German, Germany                           |
| DIP           | dual in-line package                              | GND            | ground, digital ground (dc common)        |
| Dirctry       | directory                                         | GRN            | green                                     |
| Diurt         | diuretics                                         | GRY            | gray                                      |
| DOB           | date of birth                                     |                | <b>H</b>                                  |
| DOS           | disk operating system                             |                |                                           |
| DP            | diametral pitch                                   | H              | high, vector electrode site, vector lead  |
| DPST          | double-pole, single-throw                         | HDLC           | high-level data link control              |
| DRAM          | dynamic RAM                                       | Hex, HEX       | hexagon, hexadecimal                      |
| DR/DT         | digital recording/digital transmission            | HH             | hour                                      |
| DSKTP         | desktop                                           | HiRes          | high-resolution                           |
| Dysopyr       | Dysopyramide                                      | Hr             | hour                                      |
|               | <b>E</b>                                          | Hydral         | Hydralazine                               |
|               |                                                   | Hz             | Hertz (cycles per second)                 |
| E             | enable, vector electrode site, vector lead        |                | <b>I</b>                                  |
| ecg, Ecg, ECG | electrocardiogram                                 |                |                                           |
| ECO           | Engineering Change Order                          | I              | on, input, vector electrode site          |
| EDIC          | Electrocardiograph Digital Information Center     | I, II, III     | limb leads                                |
|               |                                                   | II             | vector lead                               |
| EEPROM        | electrically erasable programmable ROM            | IC             | integrated circuit                        |
| e.g.          | for example                                       | ID             | identification                            |
| EGA           | enhanced graphics adapter                         | i.e.           | that is                                   |
| EMF           | electromotive force                               | IEC            | International Electrotechnical Commission |
| EMI           | electromagnetic interference                      | in             | inch                                      |
| ENG           | English                                           | IN             | input                                     |
| EOF           | end of file                                       | inc, inc., INC | incorporated                              |
| EPIC          | Electronic Patient Information Chart              | Info           | information                               |

|               |                                                                                                                                                                   |         |                                                       |
|---------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------|-------------------------------------------------------|
| Ins           | insert                                                                                                                                                            | min     | minutes, minimum                                      |
| I/O           | input/output                                                                                                                                                      | Misc    | miscellaneous                                         |
| I/P           | input                                                                                                                                                             | mm      | millimeter                                            |
| ISA           | industry standard architecture                                                                                                                                    | MM      | minute                                                |
| Isosorb       | Isosorbide                                                                                                                                                        | MMM     | month                                                 |
| IT            | Italian, Italy                                                                                                                                                    | mm/mV   | millimeter per millivolt                              |
|               | <b>J</b>                                                                                                                                                          | mm/s    | millimeter per second                                 |
| JAN           | January                                                                                                                                                           | MMS     | Marquette Medical Systems                             |
| JIS           | Japan Industrial Standards                                                                                                                                        | Modem   | modulator/demodulator                                 |
|               | <b>K</b>                                                                                                                                                          | MOS     | metal oxide semiconductor                             |
| k, K          | kilo, 1000, 1024                                                                                                                                                  | MPE     | metallized polycarbonate expitaxial                   |
| Kb, KB        | kilobyte                                                                                                                                                          | ms      | milliseconds                                          |
| kg, Kg        | kilogram                                                                                                                                                          | MS-DOS  | Microsoft Disk Operating System                       |
| kHz, KHz      | kilohertz                                                                                                                                                         | MTBF    | mean time between failures                            |
| kV, KV        | kilovolt                                                                                                                                                          | mtg     | mounting                                              |
| Kyb           | keyboard                                                                                                                                                          | MTR     | MOTOR                                                 |
|               | <b>L</b>                                                                                                                                                          | MUSE    | Marquette Universal System for<br>Electrocardiography |
| L             | line                                                                                                                                                              | mux     | multiplexer                                           |
| L1            | level one                                                                                                                                                         | mV      | millivolt                                             |
| L2            | level two                                                                                                                                                         | mVR     | minus (inverted) aVR                                  |
| LA            | left arm                                                                                                                                                          |         | <b>N</b>                                              |
| lb            | pound                                                                                                                                                             | N       | neutral                                               |
| LCD           | liquid crystal display                                                                                                                                            | n/a     | not available                                         |
| Lcl Line      | local line                                                                                                                                                        | NA      | not applicable                                        |
| Ld Grps       | lead groups                                                                                                                                                       | NC      | no connection                                         |
| LED           | light-emitting diode                                                                                                                                              | Nitrate | nitrates                                              |
| LH            | left hand                                                                                                                                                         | NLQ     | near letter quality                                   |
| Lidoca        | Lidocaine                                                                                                                                                         | NMI     | non-maskable interrupt                                |
| LL            | left leg                                                                                                                                                          | NMOS    | N-channel metal-oxide semiconductor                   |
| Loc           | location                                                                                                                                                          | No      | number                                                |
| LocPc         | Local MAC PC                                                                                                                                                      | NO      | normally open                                         |
| LogRetry      | log retry                                                                                                                                                         | norm    | normal                                                |
| Ltd           | limited                                                                                                                                                           | nS      | nanoseconds                                           |
|               | <b>M</b>                                                                                                                                                          | NSR     | Normal Sinus Rhythm                                   |
| m             | meter                                                                                                                                                             |         | <b>O</b>                                              |
| M             | megabyte, metric, vector electrode site,<br>vector lead, male                                                                                                     | O       | off, original                                         |
| mA            | milliamperes                                                                                                                                                      | OE      | other errors                                          |
| MAC           | Microcomputer Augmented Cardiograph                                                                                                                               | OEM     | original equipment manufacturer                       |
| mains voltage | voltage of a supply mains between 2 line<br>conductors of a polyphase system or<br>voltage between the line conductor and<br>the neutral of a single-phase system | OH      | off-hook relay                                        |
| max           | maximum                                                                                                                                                           | OneSec  | one sector                                            |
| Measure       | measurements                                                                                                                                                      | ORG     | orange                                                |
| Med           | medications                                                                                                                                                       | Orig    | original                                              |
| MEM           | memory                                                                                                                                                            | OUT     | output                                                |
| MF            | metal film                                                                                                                                                        | oz      | ounce                                                 |
| MHz           | megahertz                                                                                                                                                         |         | <b>P</b>                                              |
|               |                                                                                                                                                                   | P       | P wave (section of the ECG waveform)                  |
|               |                                                                                                                                                                   | p-p     | peak-to-peak                                          |
|               |                                                                                                                                                                   | PA      | P wave amplitude                                      |
|               |                                                                                                                                                                   | Params  | parameters                                            |

|         |                                                                               |                  |                                                              |
|---------|-------------------------------------------------------------------------------|------------------|--------------------------------------------------------------|
| Passwds | passwords                                                                     | RC               | resistor capacitor                                           |
| PatData | patient data                                                                  | RD               | R wave duration                                              |
| PatInfo | patient information                                                           | Ref              | reference, refresh                                           |
| PATN    | patient                                                                       | REN              | Ringer Equivalence Number                                    |
| PC      | printed circuit, personal computer                                            | Reserp           | Reserpine                                                    |
| PCB     | printed circuit board                                                         | REV              | revision                                                     |
| pF      | picofarad                                                                     | RevdBy           | reviewed by                                                  |
| Pgm     | program                                                                       | RevXmit          | reverse transmission                                         |
| PgmId   | program identification                                                        | rf               | radio frequency                                              |
| Phenoth | Phenothiazide                                                                 | RFI              | radio frequency interference                                 |
| Phenytn | Phenytoin                                                                     | RGB              | red, green, blue                                             |
| PID     | patient identification digit                                                  | RI               | ring indicate                                                |
| PLCC    | plastic leadless chip carrier                                                 | RL               | right leg                                                    |
| PM      | power module                                                                  | RMR              | Rhythm and Morphology Report                                 |
| pm, PM  | post meridiem, preventive maintenance                                         | ROM              | read only memory                                             |
| PM-2    | Power Module-2                                                                | RPA              | R wave amplitude                                             |
| PM-3    | Power Module-3                                                                | RPD              | R wave duration                                              |
| pn, PN  | part number                                                                   | rpt, Rpt         | report                                                       |
| PNH     | pan head                                                                      | RTC              | real time clock                                              |
| PPA     | P wave amplitude                                                              | RTI              | relative to patient input                                    |
| PR      | ECG signal interval                                                           | RTN              | return                                                       |
| Pro-Off | progressive offset                                                            | RVS              | reverse                                                      |
| Procain | Procainamide                                                                  | R/W              | read/write                                                   |
| PROM    | programmable read-only memory                                                 |                  |                                                              |
| Propran | Propranolol                                                                   |                  | <b>S</b>                                                     |
| PSK     | phase shift keying                                                            | 12SL             | 12 simultaneous leads                                        |
| PSU     | power supply unit                                                             | s, S             | second, select, switch                                       |
| Psych   | psychotropic                                                                  | SA               | s wave amplitude                                             |
| PUP     | pull-up signal                                                                | SB               | slow-blow                                                    |
| PVC     | polyvinyl chloride                                                            | SCL              | safe current limits                                          |
| PWM     | pulse-width modulation                                                        | SD               | schematic diagram, S wave duration                           |
| PWR     | power                                                                         | SE               | serial input/output errors                                   |
| PWR CRD | power cord                                                                    | sec              | second                                                       |
|         | <b>Q</b>                                                                      | sec.s            | seconds                                                      |
| Q       | transistor                                                                    | SEER             | Solid-state Electronic ECG Recorder                          |
| QA      | quality assurance, Q wave amplitude                                           | SING             | Singapore                                                    |
| QAD     | Quality Assurance Deviation                                                   | SP               | Spanish                                                      |
| QAM     | quadrature amplitude modulation (phase and amplitude modulation)              | SPA              | S wave amplitude                                             |
| QC      | quality control                                                               | SPDT             | single-pole, double-throw                                    |
| QD      | Q wave duration                                                               | SRAM             | static RAM                                                   |
| QRS     | QRS complex (portion of ECG waveform), interval of ventricular depolarization | ST-T             | ST-T wave (section of the ECG waveform)                      |
| QT      | QRS interval                                                                  | standrd, Standrd | standard                                                     |
| QTC     | QRS interval                                                                  | STD              | standard                                                     |
| QTY     | quantity                                                                      | STE              | ST segment displacement at the end                           |
| Quinid  | Quinidine                                                                     | STJ              | ST segment displacement at the J point                       |
|         | <b>R</b>                                                                      | STM              | ST segment displacement at the mid-point between STJ and STE |
| R       | resistor, red, reset                                                          | stmts, Stmts     | statements                                                   |
| RA      | right angle, right arm or R wave amplitude                                    | SumRam           | some RAM                                                     |
| RAM     | random access memory                                                          | supply mains     | permanently installed power source                           |
|         |                                                                               | SVT              | power cord type; 300 V                                       |
|         |                                                                               | sw, SW           | switch, software                                             |
|         |                                                                               | SW               | Swedish, Sweden                                              |

| <b>T</b> |                                                               | <b>X</b> |                                                                                                                                                                                                                           |
|----------|---------------------------------------------------------------|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| T Tone   | touch tone                                                    | x        | by (as in "8-1/2 x 11")                                                                                                                                                                                                   |
| TA       | T wave amplitude                                              | XCV      | transceiver                                                                                                                                                                                                               |
| Tant     | tantalum                                                      | XYZ      | orthogonal leads                                                                                                                                                                                                          |
| TDML     | treadmill                                                     |          |                                                                                                                                                                                                                           |
| TE       | timeout errors                                                |          | <b>Y</b>                                                                                                                                                                                                                  |
| Tech     | technical                                                     | Y        | year, yellow                                                                                                                                                                                                              |
| Thiazid  | Thiazide                                                      | yr       | year                                                                                                                                                                                                                      |
| TM       | trademark                                                     | yrs      | years                                                                                                                                                                                                                     |
| Tot      | total number of errors                                        | YY       | year                                                                                                                                                                                                                      |
| TP       | test point                                                    |          |                                                                                                                                                                                                                           |
| TPA      | T' wave amplitude                                             |          | <b>Symbols</b>                                                                                                                                                                                                            |
| TRAM     | Transport Remote Acquisition Monitor                          | ↑        | SHIFTed or alternate function                                                                                                                                                                                             |
| Tricyli  | Tricylic antidepressant                                       | μ        | micro                                                                                                                                                                                                                     |
| TTL      | transistor-transistor logic, TTL levels                       | μF       | microfarad                                                                                                                                                                                                                |
| TVS      | transient voltage suppressor                                  | μs, μsec | microsecond                                                                                                                                                                                                               |
| <b>U</b> |                                                               | 68K      | 68000                                                                                                                                                                                                                     |
| UE       | undefined errors                                              | &        | and                                                                                                                                                                                                                       |
| uF       | microfarad                                                    | #        | number                                                                                                                                                                                                                    |
| UL       | Underwriters' Laboratory, Inc                                 | °C       | degrees Celsius                                                                                                                                                                                                           |
| Unconf   | unconfirmed                                                   | °F       | degrees Fahrenheit                                                                                                                                                                                                        |
| UUT      | unit-under-test                                               | Ω        | Ohm, ohm                                                                                                                                                                                                                  |
| <b>V</b> |                                                               | %        | percent                                                                                                                                                                                                                   |
| v, V     | volt, volts                                                   | ®        | registered                                                                                                                                                                                                                |
| V1-V6    | precordial leads                                              | >        | greater than                                                                                                                                                                                                              |
| V123     | V1, V2, V3                                                    | <        | less than                                                                                                                                                                                                                 |
| V3R      | precordial lead                                               | ±        | plus or minus                                                                                                                                                                                                             |
| V456     | V4, V5, V6                                                    | *        | An asterisk after a signal name indicates the signal is active at its relatively lower potential, or "active-low." Signals without the asterisk suffix are active at their relatively higher potential, or "active-high." |
| V4R      | precordial lead                                               |          |                                                                                                                                                                                                                           |
| V ac     | volts, alternating current                                    |          |                                                                                                                                                                                                                           |
| V dc     | voltage, direct current                                       |          |                                                                                                                                                                                                                           |
| VA       | volt-amperes                                                  | 12SL     | 12 simultaneous leads                                                                                                                                                                                                     |
| Var      | variable                                                      |          |                                                                                                                                                                                                                           |
| VDE      | Verband Deutscher Elektrotechniker (German regulatory agency) |          |                                                                                                                                                                                                                           |
| Vent.    | ventricular                                                   |          |                                                                                                                                                                                                                           |
| VF       | ventricular fibrillation                                      |          |                                                                                                                                                                                                                           |
| VGA      | video graphics array                                          |          |                                                                                                                                                                                                                           |
| VIA      | versatile interface adapter                                   |          |                                                                                                                                                                                                                           |
| VIO      | violet                                                        |          |                                                                                                                                                                                                                           |
| Volt     | voltage                                                       |          |                                                                                                                                                                                                                           |
| VRAM     | video RAM                                                     |          |                                                                                                                                                                                                                           |
| vs       | versus                                                        |          |                                                                                                                                                                                                                           |
| <b>W</b> |                                                               |          |                                                                                                                                                                                                                           |
| w/       | with                                                          |          |                                                                                                                                                                                                                           |
| W        | watt                                                          |          |                                                                                                                                                                                                                           |
| Warfar   | Warfarin                                                      |          |                                                                                                                                                                                                                           |
| WHT      | white                                                         |          |                                                                                                                                                                                                                           |
| WI       | Wisconsin                                                     |          |                                                                                                                                                                                                                           |

