
# Service Manual

LAPARO CO<sub>2</sub> PNEU

Automatic Insufflator

2233

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## **Important general instructions for use**

Ensure that this product is only used as intended and described in the instruction manual by adequately trained and qualified personnel, and that maintenance and repair is only carried out by authorized specialized technicians.

Operate this product only in the combinations and with the accessories and spare parts listed in the instruction manual. Use other combinations, accessories and spare parts only if they are expressly intended for this use and if the performance and safety requirements are met.

Reprocess the products before every application and before returning them for repair as required by the instruction manual in order to protect the patient, user or third parties.



**Subject to technical changes!**

Due to continuous development of our products, illustrations and technical data may deviate slightly from the data in this manual.

CAUTION - USA only:

Federal law restricts this unit to be used or sold, except under the supervision of a medical doctor.

### Safety instructions and levels of danger

Symbol	Level of danger
	<b>WARNING!</b> Failure to observe can result in death or severe injury.
	<b>CAUTION!</b> Failure to observe can result in slight injury or damage to the product.
	<b>IMPORTANT!</b> Failure to observe can result in damage to the product or surrounding.
	<b>NOTE!</b> Tips for optimum use and other useful information.

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# ENGLISH

# Contents

<b>1</b>	<b>General information</b> .....	<b>1</b>
1.1	Ordering spare parts .....	1
1.2	Disassembly of housing cover .....	2
<b>2</b>	<b>Maintenance</b> .....	<b>3</b>
2.1	Important notes .....	3
2.2	Maintenance of device and accessories .....	3
2.3	Visual check .....	4
2.4	Electrical safety test to EN / IEC 60601-1 .....	4
2.5	Measuring Instrumentation and other Equipment .....	5
2.5.1	Measuring set-up for overpressure simulation .....	5
2.6	Function check .....	6
2.6.1	Alarm / safety functions .....	6
2.6.2	Error codes .....	8
<b>3</b>	<b>Modules</b> .....	<b>9</b>
3.1	Switching power supply (item 0160) .....	9
3.1.1	Disassembly and assembly of modules .....	9
3.2	Valve block (item 0170) .....	9
3.2.1	Disassembly and assembly of module .....	9
3.2.2	Function test .....	9
3.3	Pressure reducing valve unit (item 0180) .....	9
3.4	Control PCB (item 0210) .....	9
3.4.1	Disassembly and assembly of module .....	9
3.4.2	Function tests .....	9
3.5	Control PCB (item 0220) .....	10
3.5.1	Measuring points and supply voltages .....	10
3.5.2	Disassembly and assembly of modules .....	10
3.5.3	Function tests .....	10
<b>4</b>	<b>Anhang /Annex</b> .....	<b>11</b>
4.1	Reparaturteile / Repair parts .....	11
4.2	Explosionszeichnung / Exploded view .....	12
4.3	Lageplan Schaltnetzteil / Component drawing of Switching Power Supply (Pos./Item 0160) .....	13
4.4	Lageplan Bedien E-Karte / Component drawing Operation PCB (Pos./Item 0210) .....	14
4.5	Lageplan Steuer-E-Karte / Component drawing of Control PCB (Pos./Item 0270) .....	15
4.6	Blockschaltbild / Block diagram (64 157.108) .....	16
4.7	Verdrahtungsplan / Wiring diagram (64 156.251) .....	17
4.8	Pneumatikplan / Pneumatic diagram (64 280.021) .....	18
<b>5</b>	<b>Protokolle / Reports</b> .....	<b>19</b>
5.1	Prüfprotokoll / Test Report .....	19
5.2	Wartungsprotokoll / Maintenance report .....	20

# 1 General information

## 1.1 Ordering spare parts

The part number required for ordering a spare part is listed in the repair parts list under the item number used in the service manual.

 **IMPORTANT!**

*Specify the following numbers when ordering spare parts:*


- Part number of the spare part.
- Model/type number of the device.
- Serial number of the device.

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
#### Exploded view

## 1.2 Disassembly of housing cover

► Required tools: Phillips screwdriver.

- ◇ Unscrew the 4 housing feet (1.2), **fig.1**.
- ◇ Slide trim plates (1.3) from the rails (1.4), **fig.2**.
- ◇ Remove the 4 bottom Phillips screws (1.5), **fig.2**.
- ◇ Pull the housing cover (1.0) as indicated by the arrow (a) until the tongues (1.8) are disengaged from the corresponding grooves in the front plate (1.1), **fig. 3**.
- ◇ Lift off the housing cover (1.0) as indicated by the arrow (b), **fig. 3**.
- ◇ Disconnect the earthing/ground cable (1.9), **fig. 3**.

## 2 Maintenance

### 2.1 Important notes

This service manual describes the external service measures defined for the product.

**IMPORTANT!**

*For carrying out the service measures the product instruction manual must be observed.*

### 2.2 Maintenance of device and accessories

To protect the technician carrying out the inspection, the maintenance jobs and inspections must be carried out in the order specified below:

- ◇ Visual check
- ◇ Electrical safety test
- ◇ Function check

The inspection and measuring procedures listed in the section entitled "Electrical safety test" relate to the test in accordance with EN/IEC 60601-1.

Alternatively the test can be carried out in accordance with DIN VDE 751. The responsibility for complying with the current limit values and the documentation of any values measured for the first time lies in the hands of the operator.

After repair all measurement values (specified output values) of the device must be checked as described in the service manual and readjusted if deviations are observed.

**CAUTION!**

***Do not use this product if the specified measurement values or functions are not fulfilled.***

**NOTE!**

*Any maintenance or inspection job on the device or its accessories must be documented.*

## 2.3 Visual check

Item to be checked	Check for
Device and accessories	<ul style="list-style-type: none"> <li>◆ Safety-relevant soiling/contamination and general cleanliness</li> <li>◆ Mechanical damage</li> <li>◆ Loose or missing parts</li> </ul>
Controls	◆ Mechanical function and easy operation
Lettering / symbols	◆ Completeness and legibility
Safety-relevant labeling (e.g. warnings)	◆ Completeness and legibility
Device fuses	◆ The values specified by the manufacturer on the identification plate (Nominal current and melting characteristics)
Cables	◆ For perfect condition (position, insulation and brittleness)
PCBs	◆ Corrosion or other damage
Lamp plate (light sources only)	◆ In place and easy to read
	◆ Replace damaged parts immediately!

## 2.4 Electrical safety test to EN / IEC 60601-1

Item to be checked	Maximum admissible measurement values	Limit values to EN / IEC 60601-1
<b>Protective ground (earth) connection</b> Without power cord	◆ $\leq 0.1 \text{ Ohm}$	◆ $\leq 0.1 \text{ Ohm}$
<b>Protective ground (earth) connection</b> With power cord	◆ $\leq 0.2 \text{ Ohm}$	◆ $\leq 0.2 \text{ Ohm}$  <b>Test conditions:</b> $I_{\text{meas}} 25 \text{ A} \pm 10\%$ , $V_0 \leq 6 \text{ V}$ , $t_{\text{test}} 5 \text{ s to } 10 \text{ s}$ , $50 \text{ Hz / } 60 \text{ Hz}$  <b>To be tested:</b> The resistance between protective ground (earth) contact or protective grounding (earth) pin in the power plug and other exposed (bare) metal parts connected to the protective ground (earth) conductor.
<b>Ground (earth) leakage current</b>	◆ $\leq 170 \mu\text{A}$	◆ $\leq 500 \mu\text{A}$
<b>Patient leakage current</b>	◆ $\leq 10 \mu\text{A}$  <b>Measuring point:</b> On the insufflation connection	◆ (N.C.) Normal condition for BF type: $\leq 100 \mu\text{A}$ ◆ (N.C.) Normal condition for CF type: $\leq 10 \mu\text{A}$  <b>Test condition:</b> Measuring device (MD) and measuring setup to EN / IEC 60601-1  <b>To be tested:</b> The leakage current which can flow through the protective ground (earth) conductor or from the applied part via the patient to ground (earth).
<b>Additional information</b>	Disconnect the insufflation tube.	


## 2.5 Measuring Instrumentation and other Equipment

In order to carry out the maintenance, repair and test routines described in this Service Manual, the measuring instrumentation and other equipment listed in the following are compulsory. Before you start a repair job, make sure that the measuring equipment is properly calibrated and in perfect working condition.

Please observe the following: The test configurations indicated have been defined for testing purposes exclusively and do not represent a recommendation for use.

Designation	Additional data
Multimeter	True RMS
Flow tube CO <sub>2</sub>	1-10 l/min, e.g. by Platon Flow Management Jays Close, Viables, Internet: www. platon.co.uk
Digital pressure meter	Required resolution: 1 mm Hg (in pressure range below 50 mmHg)  e.g.: "VDM 100" by ASF, Munich, or "GDH 14 AN" by Greisinger electronic, D-93128 Regenstauf
CO <sub>2</sub> cylinder	0.75 kg, 2046.961
CO <sub>2</sub> pressure hose	74 021.029
Pressure reducing valve CO <sub>2</sub>	74 007.054
Abdomen model	must be able to withstand 30 mm Hg, e.g. an empty 3-liter infusion bag can be used
Hygiene filter	4171.111
Insufflation tube	8170.101
Connector	64 202.147
Syringe	50 ml
Silicon tube	approx. 1 m, 4 x 2 mm
Y-fitting	for silicon tubes
Manual	GA-A 223

### 2.5.1 Measuring set-up for overpressure simulation

## 2.6 Function check

### 2.6.1 Alarm / safety functions

Designation	Required checks
Preparation	<p>Connect the gas cylinder with pressure reducing valve and open the stop valve.</p> <p>During the test, the cylinder pressure should be in the "normal" range. The green LED of the cylinder pressure indicator is lit. The test must be stopped and the cylinder changed, if the cylinder pressure indicator changes from continuous light (green) to blinking light during the test.</p>
Switch-on condition (self-test)	<p>While the self-test is in process:</p> <p>The cylinder pressure indicator is lit and the flow indicator counts up. The self-test takes may 1 minute. The end of the self-test is acknowledged by a brief sound signal.</p> <p>The digital actual pressure indicator displays the current software version (2 digits).</p> <p>Upon completion of the self-test:</p> <p>After successful completion of the self-test, the digital displays show the following readings: gas volume consumed "0.0", flow 0 to 0.5 l/min and intraabdominal pressure <math>0 \pm 1</math> mmHg. There must be no error message EO1 to E14 on the "gas volume consumed" indicator. Every 10 minutes, a cyclic test is carried out, releasing CO<sub>2</sub>-gas through the relief valve.</p>
Flow control	<p>Connect hygiene filter with insufflation tube to the patient output. Pre-select a nominal pressure of 25 mm Hg on the LAPARO CO<sub>2</sub> PNEU automatic insufflator.</p> <p>Actuate the "INSUFFLATION" button, the button must light up. Read the current flow on the digital display..</p> <p>◆ Current flow: 1.1 l/min to 1.5 l/min</p> <p>Connect the flow tube CO<sub>2</sub> 1-10 l/min to the insufflation tube. Actuate the "high flow" button, the button must light up. Pre-select a nominal flow of 5 liters per minute.</p> <p>◆ Current flow in accordance with digital display 4 l/min to 6 l/min ◆ Current flow in accordance with flow tube: 4 l/min to 7 l/min</p> <p>Pre-select a nominal flow of 8 liters per minute. Read the current flow on the digital display and on the flow tube.</p> <p>◆ Current flow in accordance with digital display: 7 l/min to 9 l/min ◆ Current flow in accordance with flow tube: 7 l/min to 10 l/min</p> <p>Switch off insufflation</p>
Consumption indicator	<p>Actuate the "RESET" button. The digital display for the gas volume consumed reads "0.0".</p>

Designation	Required checks
Pressure control	<p>Pre-select a nominal pressure of 15 mm Hg and a nominal flow of 10 liters per minute. Connect the pressure gauge and the insufflation tube via the Y-fitting to the abdomen model.</p> <p>Actuate the "INSUFFLATION" button and then the "HIGH FLOW" button. After a short time, the pre-selected pressure of 15 mm Hg prevails in the abdomen model. It is permissible that the pressure is briefly exceeded and then released.</p> <ul style="list-style-type: none"> <li>◆ Tolerance of digital display: <math>\pm 3</math> mmHg</li> <li>◆ Deviation from reading of digital pressure gauge: <math>\pm 3</math> mmHg</li> </ul> <p>Switch off insufflation.</p> <p>Pre-select a nominal pressure of 20 mmHg. Actuate the "INSUFFLATION" button. After a short time, the pre-selected pressure of 20 mmHg prevails in the abdomen model. It is permissible that the pressure is briefly exceeded and then released.</p> <ul style="list-style-type: none"> <li>◆ Tolerance of the digital display: <math>\pm 3</math> mmHg</li> <li>◆ Deviation from reading of digital pressure gauge: <math>\pm 3</math> mmHg</li> </ul>
Gas supply monitoring	<p>Actuate "RESET". Close the stop valve on the gas cylinder. Pre-select a nominal pressure of 20 mmHg and a nominal flow of 5 liters per minute. Switch insufflation to "High Flow".</p> <p>After a short time, the gas supply indicator switches from green continuous light to blinking light. Medium alarm is sounded.</p>
Relief function in the "HIGH FLOW" mode	<p>Connect a pressure source (e.g. syringe) to the insufflation tube. The LAPARO-CO<sub>2</sub>-PNEU automatic insufflator must release the excessive pressure without triggering the overpressure alarm (medium alarm).</p> <p>Switch off insufflation.</p>
Alarm, overpressure, Overpressure switch	<p>Connect pressure source (e.g. syringe) to insufflation tube. Pre-select a nominal pressure of 15 mmHg. Connect the release opening in the housing bottom panel. Actuate the "Insufflation" button.</p> <p>If the nominal pressure is exceeded by 4 mmHg, the overpressure alarm is sounded.</p> <ul style="list-style-type: none"> <li>- Medium alarm is sounded.</li> </ul> <p>If the pressure in the insufflation tube exceeds 30 mmHg, the overpressure alarm &gt; 30 mmHg is sounded with a delay of 5 seconds.</p> <ul style="list-style-type: none"> <li>- High alarm is sounded.</li> <li>- The alarm LED is permanently lit.</li> <li>- The digital display for the gas volume consumed shows error message E11 and the digital display for pressure blinks.</li> </ul> <p>The pressure in the insufflation tube increases to approx. 50 mmHg. If the pressure in the insufflation tube exceeds 45 mmHg, the overpressure switch is triggered with some delay.</p> <ul style="list-style-type: none"> <li>- The high alarm changes to a continuous alarm.</li> <li>- The alarm LED is permanently lit.</li> <li>- Valve 2 switches off (switching sound)</li> <li>- Remove the pressure source.</li> <li>- Switch insufflation on.</li> </ul>

Designation	Required checks
Brief interruption	Switch off the device for approx. 10 seconds. When switching on again, all the pre-selected functions and settings must still be active.
Basic factory setting	Settings on the device: Nominal flow: 5 l/min Nominal pressure: 15 mm Hg

### 2.6.2 Error codes

<b>Error codes (are displayed on the “gas volume consumed” display).</b>
Error 01: CO <sub>2</sub> cylinder closed or empty
Error 02: Button is pressed when switching on
Error 03: Power supply out of tolerance
Error 04: V1 defective; Flow measurement defective; V3 defective
Error 05: V2 defective, insufflation tube closed
Error 06: Pressure measurement defective
Error 07: Flow measurement defective
Error 08: V3 defective; V2 defective; pressure measurement defective; flow measurement defective
Error 09: Invalid nominal value
Error 10: Not applicable
Error 11: Current pressure > 30 mm Hg for more than 5 seconds
Error 12: General hardware error; pressure measurement defective; V1 defective; flow measurement defective
Error 13: “Input pressure” pressure switch defective
Error 14: Electronic pressure switch defective
Error 15: Wrong DIL switch coding (electric pressure monitoring fast / slow) or Pressure switch defective

## 3 Modules

**WARNING !**

*Danger of electric shock.*

*Danger of touching live parts when working on open devices which are connected to the power supply / mains.*

*Any work on the open device while the device is connected to power, may only be carried out provided the valid safety regulations and protective measures are observed.*

### 3.1 Switching power supply (item 0160)

#### 3.1.1 Disassembly and assembly of modules

- ◇ Install pre-set card (PCB) and check complete function (in accordance with Section 2.6 and GA-A 223).

### 3.2 Valve block (item 0170)

#### 3.2.1 Disassembly and assembly of module

- ◇ Open the union nut / locking collar and pull off the outlet tube. Unscrew the additional screw connection to the housing bottom plate from underneath the device. The valves cannot be replaced individually.

#### 3.2.2 Function test

- ◇ Check the full range of functions after replacing the modules (in accordance with Section 2.6 and GA-A 223).

### 3.3 Pressure reducing valve unit (item 0180)

- ◇ The unit is supplied in pre-set condition.

### 3.4 Control PCB (item 0210)

#### 3.4.1 Disassembly and assembly of module

- ◇ Remove the front frame assembly and replace the PCB.

#### 3.4.2 Function tests

- ◇ Install pre-set card (PCB) and check full range of functions (in accordance with Section 2.6 and GA-A 223).

### 3.5 Control PCB (item 0220)

#### 3.5.1 Measuring points and supply voltages

- ◇ Plug X1:
  - Pin 1 = PFD
  - Pin 2 =  $-12V \pm 0.4V$
  - Pin 3 = GND
  - Pin 4 = + 24V
  - Pin 5 = GND
  - Pin 6 =  $+ 12V \pm 0.4V$
  - Pin 7 = GND
  - Pin 8 = + 5V
  - Pin 9 = NC
  - Pin 10 = + 5V
  
- Stecker X3:
  - Pin 1 = DS 2 (low)
  - Pin 2 = DS 1 (high)
  - Pin 3 =  $2V \pm 0.4V$
  
- Stecker X5:
  - Pin 1 = GND
  - Pin 2 = + 5V
  - Pin 3 = GND
  - Pin 4 = + 5V
  - Pin 5 = GND
  - Pin 6 = Heating
  - Pin 7 = GND
  - Pin 8 = Flow > 2 l / min
  - Pin 9 = GND
  - Pin 10 = SDA
  - Pin 11 = GND
  - Pin 12 = SCL
  - Pin 13 = GND

#### 3.5.2 Disassembly and assembly of modules

- ◇ Remove the front frame assembly and replace the PCB.

#### 3.5.3 Function tests

- ◇ Install the pre-set card (PCB) and check the full range of functions (in accordance with Section 2.6 and GA-A 223).

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**RICHARD WOLF Belgium S.A.**
**RICHARD WOLF France S.A.R.L.**