

High-Performance Oven 220 - 240 Voltage



Installation and Operation Manual

SMO5HP-2

Previously designated as

HF4-2

Shel Lab Oven High Performance 220 - 240 Voltage

Installation and Operation Manual

Part Number (Manual): 4861706

Revision: September 22, 2015

This oven requires permanent connect wiring (also known as hardwiring). It does not connect to a wall power source using a conventional plug-in power cord.

Purge Valve

The SMO5HP-2 can be built with a purge valve located on the back of the unit. The valve can be connected to a cylinder of nitrogen (N₂) or other neutral gas to purge the oven chamber. This option must be ordered prior to the construction of the oven.

These units are TÜV listed as high performance ovens for professional, industrial, or educational use where the preparation or testing of materials is done at an ambient air pressure range of 22.14 – 31.3 inHg (75 – 106 kPa), with no flammable, volatile, or combustible materials being heated.

CAN/CSA C22.2 No. 61010-1:2012

CAN/CSA C22.2 No. 61010-2-010:2004 Reaffirmed: 2014-07

UL 61010-1:2012-05

UL 61010A-2-010:2002-03

EN 61010-1:2010

EN 61010-2-010:2014

Supplemented by: UL 61010-2-010:2015

TABLE OF CONTENTS

INTRODUCTION.....	4
General Safety Considerations	4
Engineering Improvements	5
Contacting Assistance.....	5
Temperature Reference Sensor Device	5
RECEIVING YOUR OVEN.....	6
Inspecting the Shipment.....	6
Recording Data Plate Information	6
Orientation.....	7
INSTALLATION	9
Required Ambient Conditions	9
Power Source Requirements	10
Power Feed Wiring.....	10
Lifting and Handling	11
Leveling.....	11
Install the Oven	11
Initial Cleaning.....	11
Shelving Installation	12
GRAPHIC SYMBOLS.....	13
CONTROL PANEL OVERVIEW.....	14
OPERATION	15
Operating Precautions	15
Theory of Operations	16
Preparing the Oven for Use	18
Set the Over Temperature Limit.....	19
Set the Temperature Controller Set Point.....	19
First Use Burn In	20
Temperature Accuracy Verification	21
Closed Vents.....	23
Over Temperature Limit Active	23
Programmed Operations.....	24
Positive Venting of Exhaust	24
USER MAINTENANCE.....	26
Cleaning and Disinfecting	26
Door Latch and Electronics	27
Door Gaskets	27
Removing the Chamber Liner	28
Calibrate the Temperature display.....	30
UNIT SPECIFICATIONS.....	34
Weight	34
Dimensions.....	34
Capacity	34
Shelf Capacity by Weight	34
Temperature.....	35
Airflow Performance	35
Power	35
REPLACEMENT PARTS.....	36

INTRODUCTION

Thank you for purchasing a Sheldon Manufacturing High Performance Oven. We know that in today's competitive marketplace, customers have many choices when it comes to constant temperature equipment. We appreciate you choosing ours. Our continued reputation as a leading laboratory product manufacturer rests with your satisfaction. Sheldon Manufacturing, Inc. stands behind our products, and we will be here if you need us.

These ovens are intended for professional, industrial, and educational applications requiring horizontal shelf space airflows. The ovens are not designed for use at hazardous or household locations.

Before you use your oven read this entire manual carefully to understand how to install, operate, and maintain the oven in a safe manner. Keep this manual available for use by all oven operators. Ensure that all operators are given appropriate training before the oven begins service.

GENERAL SAFETY CONSIDERATIONS

Note: Failure to follow the guidelines and instructions in this manual may create a protection impairment by disabling or interfering with the unit's safety features. This can result in injury or death.

Your oven and its recommended accessories have been designed and tested to meet strict safety requirements. For continued safe operation of your oven, always follow basic safety precautions including:

- Follow all local or regional ordinances in your area regarding the use of this unit. If you have any questions about local regulations, please contact the appropriate agency.
- Use only approved accessories. Do not modify system components. Any alterations or modifications to your oven may be dangerous and will void your warranty.
- Always hardwire the unit's power feed to an earth grounded electrical source that conforms to national and local electrical codes. If the unit is not grounded, parts such as, knobs and controls may conduct electricity and cause serious injury.
- Avoid damaging the power feed. Do not bend it excessively, step on it, place heavy objects on it. A damaged power feed can easily become a shock or fire hazard. Never use a power feed after it has been damaged.
- Do not position the equipment in such a manner as to make it difficult to disconnect or uncouple the power feed.
- Do not attempt to move the unit while in operation or before the unit has cooled.

INTRODUCTION (CONTINUED)

ENGINEERING IMPROVEMENTS

Sheldon Manufacturing continually improves all of its products. As a result, engineering changes and improvements are made from time to time. Therefore, some changes, modifications, and improvements may not be covered in this manual. If your unit's operating characteristics or appearance differs from those described in this manual, please contact your Shel Lab dealer or customer service representative for assistance.

CONTACTING ASSISTANCE

If you are unable to resolve a technical issue with the oven, please contact Sheldon Technical Support. Phone hours for Sheldon Technical Support are 6am – 4:30pm Pacific Coast Time (west coast of the United States, UTC -8).

Please have the following information ready when calling or emailing Technical Support: the **model number** and the **serial number**. These will be found on the unit's data plate, which is located on the back of the unit at the top right, next to the power supply as mandated by regulatory requirement.

EMAIL: tech@shellab.com PHONE: 1-800-322-4897 extension 4, or (503) 640-3000 FAX: (503) 640-1366

TEMPERATURE REFERENCE SENSOR DEVICE

The oven does not come with a temperature reference device. A reference sensor device must be purchased separately for performing accuracy verifications and calibrations of the oven temperature display.

The reference device must be accurate to at least 2°C, and should be regularly calibrated, preferably by a third party.

For best results, use a digital device with wired-connected temperature sensing probe. For example: a wire thermocouple probe that can be introduced into the oven chamber, leaving the device outside. Readings taken from outside the unit avoid chamber door openings during verifications and calibrations, and eliminate subsequent waits for the chamber temperature to re-stabilize before proceeding.

Select a probe suitable for the application temperature you will be calibrating or verifying the display accuracy at.

Alcohol thermometers are insufficient for conducting accurate verifications and calibrations. Do not use a mercury thermometer. **Never place a mercury thermometer in the oven chamber.**

RECEIVING YOUR OVEN


Before leaving our factory, all High Performance Ovens are packaged in high-quality shipping materials to provide protection from transportation-related damage. When an oven leaves our factory, safe delivery becomes the responsibility of the carrier. Damage sustained during transit is not covered by the oven's warranty. When you receive your High Performance Oven, inspect it for concealed loss or damage to its interior and exterior. If you find any damage to the oven follow the carrier's procedure for claiming damage or loss.

INSPECTING THE SHIPMENT

Carefully inspect the shipping carton for damage. Report any damage to the carrier service that delivered the oven. If the carton is not damaged, open the carton and remove the contents. The unit should come with an Installation and Operation Manual, warranty card, a Certificate of Compliance, and the Programming Guide – Watlow EZ-Zone Heating Profiles.

Verify that the correct number of shelves, shelf slides, and leveling feet have been included (see the following table for quantities).

Included accessories

Shelves (2)	Shelf Slides (4)	Leveling Feet (4)
		

Carefully check all packaging before discarding. Save the shipping carton until you are certain that the unit and its accessories function properly.

RECORDING DATA PLATE INFORMATION

Locate the data plate on the back of the oven next to the power inlet. The data plate contains the oven model number and serial number. Enter this information below for future reference.

Date Plate Information

Model Number	
Serial Number	

RECEIVING (CONTINUED)

ORIENTATION




Figure 1: SMO5HP-2 Oven

GRAPHIC SYMBOLS (CONTINUED)

This page left blank

INSTALLATION

Installation of an SMO5HP-2 requires permanent connect wiring (commonly known as hardwiring), and **must be performed by a qualified electrical technician**. All other installation steps can be performed by the end user.

REQUIRED AMBIENT CONDITIONS

This oven is intended for use indoors, at room temperatures between **15°C and 40°C (59°F and 104°F)**, at no greater than **80% Relative Humidity** (at 25°C / 77°F).

Allow a minimum of **12 inches (30cm)** between the oven and any walls or partitions, and **24 inches (60cm)** of vertical headspace clearance above the top of the oven for unobstructed airflow and cooling. If the oven exhaust will be vented from the workspace through a duct or other channeling, **12 inches (30cm)** of vertical clearance will suffice. Make sure the intake and exhaust vent remain unobstructed. Do not place objects on top of oven, aside from the power exhaust accessory described in this manual.

Operating the unit outside these conditions may adversely affect the oven's temperature range, uniformity, and stability. For conditions outside of those listed above, please contact your distributor or Shel Lab Sales to explore other ovens suited to your production of laboratory environment.

When selecting a location to install the oven, consider all conditions that can affect temperature performance. For example:

- Ovens, autoclaves, and any device that produces significant radiant heat
- Heating and cooling ducts, or other sources of fast moving air currents
- High-traffic areas
- Direct sunlight

INSTALLATION (CONTINUED)

POWER SOURCE REQUIREMENTS

When selecting a location for the oven, check that **each** of the following requirements are satisfied:

- The oven must be positioned so that all operators have access to the power feed disconnect in case of emergencies.
 - The Disconnect must be in close proximity to the equipment and within easy reach of the operator.
 - The Disconnect must be marked as the disconnecting device for the equipment.
- The power source must be **single (1) phase** and **protective earth grounded**.
- The power source must match the voltage and ampere requirements listed on the oven data plate. This oven is intended for **220 - 240 volt, 50/60 Hz** applications at **12 Amps**.
- The power source for the oven must conform to all national and local electrical codes.
- **Supplied voltage must not vary more than 10% from the data plate rating.** Damage to the oven may result if supplied voltage varies more than 10%. A separate circuit is recommended to prevent possible loss of product due to overloading or failure of other equipment on the same circuit.
- **A switch or circuit-breaker must be used in the building installation.** The required circuit-breaker is **20 Amps**

Fuses: The oven is provided with (2) **two 20 amp fuses**. These are located in holders adjacent to the power feed braid on the back power panel of the unit. The fuses protect against overcurrent conditions related to the operation of the oven. If one fuse blows, the oven will shut down. For safe operation, the cause of a blown fuse should be determined prior to replacing it.

The oven is also provide with a pair (2) of 2 amp fuses installed adjacent to the external power outlet. Please see the [Orientation photo](#) on page 7. The fuses protect against overcurrent conditions related to the operation of the outlet and an attached power exhaust blower. If one fuse blows, the outlet will depower. The cause of a blown fuse should be determined prior to replacing it.

POWER FEED WIRING

The oven comes provided with an integral 6 inch (15 cm) wire braid of two 14 gauge high-temperature (300°C) hot wires and a 14 gauge earth ground. See the [Orientation Picture](#) on page 7.

The wires for power source connection should be Green/Yellow – Earth; Black – Hot; Red – Hot.

The oven must be earth grounded using the protective conductor terminal (green with yellow stripe wire). Do not remove the protective conductor (earth connection). Removing the protective conductor will negate the oven's protections against potentially dangerous electric shocks and create a possible fire hazard.

INSTALLATION (CONTINUED)

LIFTING AND HANDLING

The oven is heavy. Use appropriate lifting devices that are sufficiently rated for these loads. Follow these guidelines when lifting the oven:

- Lift the oven only from its bottom surface.
- Doors, handles, and knobs are not adequate for lifting or stabilization.
- Restrain the oven completely while lifting or transporting so it cannot tip.
- Remove all moving parts, such as shelves and trays, and secure the door in the closed position during transfers to prevent shifting and damage.

LEVELING

The oven must be level and stable for safe operation. Each Clean Room Oven ships with four leveling feet. Insert one leveling foot into each of the four holes in the bottom corners of the oven. Stand the oven upright. Then, adjust the foot at each corner until the oven stands level and solid without rocking. To raise a foot, turn it in a counterclockwise direction; to lower a foot, turn it in a clockwise direction.

Note: To prevent damage when moving the oven, turn each of the four leveling feet completely clockwise.

INSTALL THE OVEN

Install the unit in a workspace location that meets the criteria discussed in the previous entries of the Installation section. After wiring the oven to its power source, only energize the unit for brief periods as necessary to test the wiring power feed and source.

INITIAL CLEANING

The unit was cleaned at the factory, but not sterilized. Additionally, it may have been exposed to contaminants en route during shipping. See the [Cleaning and Disinfecting](#) topic in the User Maintenance section (page 26) for more information on how to clean and disinfect the oven chamber.

Remove all wrappings and coverings from shelving prior to cleaning and installation.

INSTALLATION (CONTINUED)

SHELVING INSTALLATION

The horizontal airflow within the chamber moves from the air duct on the right-hand side of the chamber, to the large holes on the left side. Place the shelves so they do not obstruct the duct holes in order to maximize airflow across the shelf space.

1. Install shelf slide hangers on opposites of the oven by inserting the slider's tabs into the chamber's mounting slots, then pushing down gently to secure the slider.
2. Hang the shelf from the two installed sliders.





Figure 3: Installed Sliding Shelf

GRAPHIC SYMBOLS

The oven is provided with multiple graphic symbols on its external and internal surfaces. The symbols identify hazards and the functions of the adjustable components, as well as important notes found in the user manual.

Symbol	Definition
	Indicates that you should consult your user manual for further instructions. Indique que l'opérateur doit consulter le manuel d'utilisation pour y trouver les instructions complémentaires.
	Indicates Adjustable Temperature Indique température réglable
	Indicates AC Power Repère le courant alternative
	Indicates I/ON and O/OFF I repère de la position MARCHÉ de l'interrupteur d'alimentation O repère de la position ARRÊT de l'interrupteur d'alimentation
	Indicates protective earth ground Repère terre électrique
	Indicates UP and DOWN respectively Touches de déplacements respectifs vers le HAUT et le BAS
	Indicates a Potential Shock Hazard Signale danger électrique
	Indicates the unit should be recycled (Not disposed of in land-fill) Indique l'appareil doit être recyclé (Ne pas jeter dans une décharge)
	Indicates: Caution hot surface Indique: Avertissement symbole de surface chaude

CONTROL PANEL OVERVIEW



Figure 3: Control Panel

Power Switch



The green, light up power switch controls all power to the oven. It must be in the illuminated I (ON) position before any systems or components are operational.

Temperature, Profile, and OTL Controller



The temperature controller for the oven is labeled TEMPERATURE AND HIGH LIMIT CONTROL and consists of a Watlow EZ-Zone Controller. The controller is provided with a digital display and Up and Down arrow keys. The keys are used to change the oven temperature set point, perform calibrations, and program heating profiles. The arrow keys also control the operations of the Over Temperature Limit system (OTL).



Figure 4: Temperature Controller



The EZ1 button is used to quick launch Profile 1 or abort any running heating profile (see [the Programmed Operations section](#) on page 24). The EZ2 button does not have a function in the SMO5HP-2.



The green advance button is used to scroll forward through menus and parameters lists.



The reset button is used to scroll the display back to a previous page or menu while programing heating profiles. Pressing the reset button repeatedly returns the display to the home page.

Note: On some older Watlow Controllers the Reset button may be labeled with an infinite ∞ symbol rather than RESET.

A small flashing “2” toward the center bottom of the display indicates when the controller is powering the oven’s heating elements. A small flashing “1” indicates internal communications functions within the controller.

Vent Valves

The SMO5HP-2 comes provided with two vent valve controls on its front control panel. These open or close the intake and exhaust valves located on the top of the unit.

OPERATION

OPERATING PRECAUTIONS



Warning: This SMO oven is not an explosion-proof unit!

Avertissement: Ce sont des fours pas résistants aux explosions.

1. The oven is not designed to safely contain combustible gasses.
2. Do not place explosive, combustible, or flammable materials into the chamber.
3. The bottom surface of the chamber should **not** be used as a work surface.
4. Do not operate the oven in an environment with noxious fumes.
5. The SMO oven is provided with a dampened exhaust vent. For safe and efficient oven operation follow these precautions:
 - a. In normal heating the damper is closed.
 - b. Outgassed byproducts may be hazardous to or noxious for operating personal. If either is the case, the oven exhaust should be positively ventilated to a location outside workspace in accordance with national and local regulations.
6. Do not place sealed or filled containers in the oven.
7. Do not place mercury thermometers in the oven.
8. The SMO oven is not designed for use in Class I, II, or III locations as defined by the US National Electric Code.



Warning: The vent dampers may be hot to the touch. These areas are marked with Hot Surface labels. Proper PPE should be employed to minimize burn risk.

Avertissement: Les clapets d'aération peuvent être chauds au toucher. Ces zones sont marqués avec des étiquettes de Surface chaude. Les EPI approprié devraient être employée pour réduire au minimum le risque de brûler.

OPERATION (CONTINUED)

THEORY OF OPERATIONS

To achieve the oven's stated operation specs, the vents must be fully closed. Operating the oven with the vents open impacts the oven chamber temperature uniformity and stability. This may or may not speed drying rates, depending on the nature of your samples or product. During normal operations the vents are intended to be opened only while the product or samples are cooling. Running the oven with an active power exhaust blower attached to the exhaust vent will reduce the oven chamber temperature.

Heating

The oven runs in either manual mode or an automatic heating profile. The oven comes from the factory set to run in manual mode. When powered and running in manual mode, the SMO oven chamber heats to and then indefinitely maintains the temperature set point shown on the home page of the controller. The manual mode set point can be adjusted by the end-user using the Up and Down temperature buttons. The oven will not heat if the temperature set point is set below its operational temperature range.

When the oven is powered, the chamber air temperature cannot go below the ambient room temperature **plus** the internal waste heat of the oven. Waste heat is generated primarily by the operation of the air-circulating blower fan motor and the resulting air compression in the duct spaces. In practice the temperature floor is **ambient +15°C**.

The oven heater controller (a Watlow EZ digital controller) stores the user-selected temperature set point, along with any user-programmed heating profiles. The controller monitors the oven chamber air temperature using a solid-state probe located in the airstream on the right wall of the chamber. When the processor detects that the chamber temperature has dropped below the current temperature set point, it pulses power to a heating element in a recirculation air duct located above the oven chamber.

The controller employs proportional-integral-derivative analytical feedback-loop functions when measuring and controlling the chamber air temperature levels. PID-controlled heating pulse intensities and lengths are proportional to the difference between the measured chamber temperature and the current set point. The frequency of pulses are derived from the rate of change in the difference. The integral function slows the rate of pulses when the temperature nears the set point to avoid overshooting.

SMO ovens rely on natural heat radiation for cooling. The oven depends on the operation of its blower fan to circulate air and maintain temperature uniformity and stability in the chamber.

The heating rates given in unit specification section of this manual are for a 25°C environment. The ambient temperature of the workspace around the oven will affect heating and cooling.

OPERATION (CONTINUED)

Heating Profiles

The Watlow controller can store up to four heating profiles, each with ten steps. Please see the [Programmed Operations](#) section on page 24, and the *Watlow EZ-Zone Profile Programing Guide* for more details. A profile can be launched whenever the unit is running in manual mode. Heating profiles should allow adequate time for the oven chamber to heat to the first temperature in the profile, prior to continuing to the second step of the profile.

Airflow and Applications

During normal operations, the internal blower fan motor pulls air out of the oven chamber through perforations on the left chamber wall. The captured air is pulled up and blown across a heating element in the recirculation duct immediately above the chamber, then forced down through a duct on the right side of the chamber. The air exits the duct through small holes on the right chamber wall. Upon leaving these holes, the heated air blows horizontally across the shelf space, helping to create temperature uniformity and enhance drying. After crossing the shelves, the air enters a duct on the left side of the chamber, and is then pulled back into the heater space above the oven chamber.

If the intake and exhaust vents are opened, some air will exit through the open exhaust vent on the top right of the unit. Lost air will be replaced through the intake vent. The rate of atmosphere exchange depends on several variable factors. These include: ambient conditions, position of the vent valves, the volume of material in the oven chamber, and if a power exhaust unit is connected to the exhaust vent.

Over Temperature Limit

The Watlow controller contains a backup temperature control system with independent circuitry. This system connects to a redundant temperature sensor probe inside the oven chamber. The over temperature heating cutoff system is intended to be set 5°C above your application or process temperature set point. In the event the oven chamber temperature exceeds the over temperature limit setting, the backup system will depower the heating element. This provides redundant heating control the event of a failure of the primary circuitry or the main temperature sensor probe.

OPERATION (CONTINUED)



Home Page

27°C Chamber air
temperature, 0°C
set point

PREPARING THE OVEN FOR USE

Perform the following steps to prepare the oven for use after installing it in a new workspace environment.

Optional: If you will be performing a temperature display accuracy verification as part of your set up, Sheldon Manufacturing recommends introducing the temperature probe of your temperature reference device into the oven chamber now, prior to powering the oven. Please see the procedure (page 21) for probe placement instructions. Placing the probe now offers significant time savings during the procedure. The probe must be rated for the temperature you will be running the oven at during the burn-in (see page 20).

1. Open the intake and the exhaust vents on the oven to prepare for the initial burn-in procedure.
2. Place the Power Switch in the ON position. The Temperature Display will illuminate, then show the controller home page.
3. Set up the oven for use by carrying out the following procedures found in the Operation section:
 - **Set the Over Temperature Limit** to 325°C (page 19) to prevent the OTL heating cutoff from interfering with the following three procedures.
 - **First User Burn-In** (page 20)
 - **Set the Temperature Controller Set Point** to the highest temperature you will use the oven at to run the Burn-In procedure manually (page 19).
 - Or –
 - **Program a heating profile** to automatically run the Burn-In. See pages 20 and 24, as well as the *Programming Guide – Watlow EZ-Zone Controller Heating Profiles*, which came with the oven.
 - **Optional: Perform a Temperature Display Accuracy Verification** (page. 21)
 - **Set the Temperature Controller Set Point** on the home page to your application temperature, if you will be running the oven in its manual mode (page 19). Otherwise, **program a heating profile** (page 24).
 - **Set the Over Temperature Limit** (page. 19)
 - This time set the control to 5°C above your application's manual mode temperature or the highest temperature of your heating profile. This prepares the oven for normal use.
 - Load the oven with your product or samples.

End of procedure

OPERATION (CONTINUED)

SET THE OVER TEMPERATURE LIMIT

This procedure sets the Over Temperature Limit backup system. Normally the system is set 5°C above the highest set point temperature you will be using during your current application or process. Always set at least 5°C above the hottest temperature of a heating profile.

Note: If this is the first use of the oven, set the Over Temperature Limit to 325°C or higher to prepare to help ensure the OTL system will not interfere with the Burn-In procedure.

1. Press the green Advance button until "Lh.S1" (Limit High Set Point) shows in the green-colored mid-level display.
2. Use the Up and Down arrows to enter a temperature limit in the red top display.
3. After entering the temperature limit, press the Advance button. The red display will show "SAFE", indicating that the temperature limit has been set.
4. Press the Reset button to scroll the display back to the home page.



Steps 1 and 2:
Limit set at 325°C




Step 3

SET THE TEMPERATURE CONTROLLER SET POINT

This procedure sets a new temperature set point for the oven's manual mode.

Note: Set the oven to the highest temperature you will run the oven at, if manually running the First Use Filter Burn-In procedure (see page 20). This procedure may also be run by programing and launching it as a heating profile.

1. On the home page, use the Up and Down Arrow keys to adjust the set point in the green display to your baking application or process temperature set point.
 - a. Holding down an arrow key will cause the temperature to advance in increments of ten (10) rather one (1).
2. Release the Arrow Button once you have reached your required set point. There may be a brief pause before the oven starts heating.
 - a. A small flashing 2 in the Watlow display indicates the temperature controller is calling for heat.



Set Point at 0°C,
Chamber temp
27°C



Step 2: Set Point
adjusted to 280°C

End procedure

OPERATION (CONTINUED)

FIRST USE BURN IN

The steps below help eliminate smoking from protective coatings on the heating elements, when running above 150°C. Perform the burn-in at the highest temperature the oven will be used at.

1. Open the intake and exhaust vents. Oven exhaust should be positively vented outside of the oven's workspace during the Burn-In.
2. If you have not already done so, set the Over Temperature Limit to 325°C, and the temperature set point the highest operating temperature the oven will be used at.
3. Run the oven for a minimum of one (1) hour at your highest operational temperature, until smoke from the heating elements dissipates.

OPERATION (CONTINUED)

TEMPERATURE ACCURACY VERIFICATION

Note: Performing an accuracy verification of the temperature display requires a temperature reference device. Please see the [Reference Sensor Devices entry](#) on page 5 for device requirements.

Optional: A verification of the temperature display accuracy may be carried out when preparing the oven for use, if required by your laboratory or production protocol. The verification compares the displayed temperature of the oven chamber with the actual temperature of the chamber air as provided by a reference sensor device.



Figure 4: Damper in Closed Position




Figure 5: Probe End 2 inches (5cm) From Shelf Surface

If a difference between the actual and displayed temperatures is discovered, perform a temperature calibration. Please see the [Calibrate Temperature Display procedure](#) on page 30 in the User Maintenance section.

Close the Vents

Both the intake and exhaust vents must be closed in order to perform an accurate temperature verification. Chamber temperature uniformity and stability will be too unstable when operating with the vents open to verify the accuracy of the oven temperature display.

Probes

A reference device sensing probe may be introduced through the through the chamber door space. Use heat-resistant, non-marking tape to secure the wires to the body of the incubator, and seal any gaps in the door space. The door must close and latch fully.

Place the sensor probe of the temperature reference device as close as possible to the geometric center of the oven chamber. A thermocouple sensor probe sleeve may be taped to the shelving, as long as the exposed copper end is 2 inches (5cm) away from the shelf (see Figure 5). An exposed sensor probe in direct contact with the shelving may experience heat sinking, which can result in an inaccurate temperature reading.

Pre-Heating

Verify that the oven is set to your application or baking procedures set point. Let the oven heat to the set point temperature and run for at least 1 hour undisturbed in order to stabilize.


Temperature Stability

The oven chamber is thermally stabilized when no fluctuations of $\pm 2.0^{\circ}\text{C}$ or greater have been detected with the verification reference device **for a minimum of 30 minutes**. Failure to wait for stabilization will result in a failure to verify the accuracy of the oven temperature display reading.



OPERATION (CONTINUED)

Note: If not already set to maximum, set the OTL limit to its max setting (325°C) to prevent the OTL cutoff from interfering with the verification process. Please see the [Set The Over Temperature](#) procedure on page 19.

Verifying the Temperature Display Accuracy	
1. Once the chamber temperature has stabilized, compare the reference device and the oven display temperature readings. a. If the readings are the same, or the difference between the two (2) falls within the acceptable range of your protocol, the display is accurately showing the oven chamber air temperature. The Temperature Verification procedure is now complete. a. See step 2 if a difference falls outside the acceptable range of your protocol.	Reference Device  
2. Perform a temperature calibration to match the controller display to the actual chamber temperature if a difference falls outside your protocol range. a. Please see 30 in the User Maintenance section.	Reference Device  

End of procedure

OPERATION (CONTINUED)

CLOSED VENTS

The unit is intended to be run with both vents closed in order to achieve its stated operating temperature specifications. The unit can be run with one or both vents open, however, the oven may not achieve its stated chamber uniformity, stability, or heating or recovery times. Drying times are dependent on the nature of your product or samples, as well as the volume placed in the oven chamber.

OVER TEMPERATURE LIMIT ACTIVE

When the oven chamber air temperature exceeds the user-programed Over Temperature Limit, the OTL system will depower the heating elements. This is normally accompanied by a mechanical click. A small “4” will show in the bottom display, indicating that the controller is rerouting power away from the heating elements. Two alternating alert screens will flash, indicating a failure and that the Limit High One (LiH1) alert has been activated.

The Over Temperature Limit will activate if one of two events happens:

- You have set the Over Temperature Limit below the temperature set point you are using for your application or process.
- The main controller circuitry or sensor probe has failed, and must be replaced in order to maintain safe and accurate oven operations.

If you have set the Over Temperature Limit too low, perform the following steps to take the unit out of shutdown:



1. Press the advance key until Ignore “i9nr” shows in the top display, and Limit High “Lih1” in the green display.



2. Press the green advance key again. Limit High Set Point “LhS1” will appear in the green display, and the Over Temperature Limit setting in the red top display.



3. Adjust the Over Temperature Limit to above your set point using the Up Arrow Button.



4. Push the Reset Button to save the new OTL setting.

If the Over Temperature Limit has activated and the red oven temperature is higher than your green temperature set point, the main controller circuitry or temperature sensor may have failed. **Cease using the oven and contact Sheldon Technical Support for assistance.**



Alternating OTL Alert Screens



Step 1



Step 3

OPERATION (CONTINUED)

PROGRAMMED OPERATIONS

The Watlow controller can hold four (4) ten-step heating profiles. Additionally, profiles may be combined by programming profile steps sequentially to run as a single profile of between two (2) and forty (40) steps.

Please see the **Programing Guide – Watlow EZ-Zone Controller Heating Profiles** document that comes included with the oven for instructions on how to program heating profiles. The guide is intended as simplified explanation for all major heating profile functions and programming steps. The guide is also available on the Shel Lab website in PDF format.

Refer to Chapter 7 of the “EZ-Zone PM User’s Manual” for highly detailed instructions on how to program the EZ Watlow Controller.

POSITIVE VENTING OF EXHAUST

Exhaust ducting can be connected to the oven’s exhaust port to channel or positively vent exhaust out of the workspace around the oven. Additionally, Sheldon Manufacturing sells a 220 - 240 volt High Performance Power Exhaust that can be mounted directly on the exhaust vent, and powered by the oven. An exhaust duct should not extend straight up from the oven, but should include a steep bend sufficient to stop condensation within the ducting from sliding down into the oven.



Figure 5: High Performance Power Exhaust (PN 9990741)

Mounting the Power Exhaust

1. Remove the screws on the cover of the exhaust vent assembly on the top of the oven, but leave the assembly in place.
2. Mount the power exhaust blower on the exhaust vent cover assembly, aligning the blower and the assembly’s screw holes. The open side of the blower’s mounting body should go over the handle of the exhaust port’s sliding vent cover.
3. Reinstall the screws.
4. Plug in the power exhaust into the 220 - 240 volt outlet on the back of the oven.

OPERATION (CONTINUED)

Activating and Deactivating the Power Exhaust

To activate the power exhaust as part of a heating profile step, do the following:



1. While selecting the parameters for the step, choose Event 1 On ("Ent1" "on").
 - a. The Event 1 parameter is located after the Time parameter when scrolling through a step menu using the green Advance Key.
 - b. Use either of the Arrow Keys to toggle "oFF" to "on".
2. Press the Advance key to scroll to the next event and save Event 1 as On.
3. Next time the profile runs, the 220 - 240 volt outlet on the back of the oven will power up when the controller reaches this step.
4. The plugged in power exhaust will run until the next step in which Event 1 Off is selected ("Ent1" "oFF"), or until the oven shuts down.
5. By default the Event 1 parameter is set to Off. It must be switched to On **for each successive step** that you want the blower to run during the execution of the profile.



USER MAINTENANCE



Warning: Prior to maintenance or service on this unit, disconnect the power feed from the power supply.

Avertissement: Avant d'effectuer toute maintenance ou entretien de cet appareil, débrancher le cordon secteur de la source d'alimentation.

CLEANING AND DISINFECTING

If a **hazardous material or substance** has spilled in the oven, immediately initiate your site's Hazardous Material Spill Containment protocol. Contact your local Site Safety Officer and follow instructions per the site policy and procedures.

The unit should be cleaned prior to first use.

Do not use spray-on cleaners or disinfectants. These can leak through openings and coat electrical components. Do not use cleaners or disinfectants that contain solvents capable of harming paint coatings or stainless steel surfaces. **Do not use chlorine-based bleaches or abrasives; these will damage the chamber liner.** Consult with the manufacturer or their agent if you have any doubts about the compatibility of decontamination or cleaning agents with the parts of the equipment or with material contained in it.



Warning: Never clean the oven with alcohol or flammable cleaners.

Avertissement: Ne jamais nettoyer l'appareil à l'alcool ou avec des nettoyeurs inflammables.

Cleaning

Perform the steps below to clean the oven:

1. Remove all removable interior components such as, shelving and accessories.
2. Clean the oven with a mild soap and water solution, including all corners. **Do not** use an abrasive cleaner that will damage stainless steel surfaces. **Do not use deionized water to rinse or clean with.**
3. Rinse with distilled water and wipe dry with a soft cloth.
4. Take special care when cleaning around the temperature sensor probes at the back of the chamber to prevent damage.

Disinfecting

Carry out the following steps to disinfect the oven, if required by your laboratory or production protocol:

1. Turn the unit off. Carry out your disinfection protocol.
2. If allowed by your protocol, remove all shelving. Disinfect all corners. Take special care not to damage the temperature sensor probes when disinfecting.
3. Disinfect the unit using commercially available disinfectants that are non-corrosive, non-abrasive, and suitable for use on stainless steel surfaces. Contact your Site Safety Officer for detailed information on the disinfectants compatible with your study or production protocol.

USER MAINTENANCE (CONTINUED)

DOOR LATCH AND ELECTRONICS

Periodically, inspect the door latch, trim, catch, and gasket for signs of deterioration. Failure to maintain the integrity of the door system shortens the life span of the Oven.

Electrical components do not require maintenance. If the oven fails to operate as specified, please contact your Shel Lab distributor or Sheldon Manufacturing Technical Support for assistance.

DOOR GASKETS

The Sheldon High Performance Oven is provided with a high temperature silicon door gasket. See the [Parts Section](#) on page 29 for ordering replacement gaskets.





USER MAINTENANCE (CONTINUED)

REMOVING THE CHAMBER LINER

This procedure removes the ceiling liner and the side wall air ducts. The liner and ducts should be removed periodically and cleaned, and the surfaces beneath them cleaned.

Note: Door removed in illustrations for clarity. Do not remove the door from the oven.

1. Remove all shelves and shelf sliders from the oven.
2. Unscrew the nuts located at the top of the back wall of the oven chamber.
3. Remove the chamber gasket so that there is room to pull the ceiling liner out through the door space. The gasket may be fully or partially removed.
4. Remove the chamber ceiling liner by pulling it out through the oven door space. The liner may require some side-to-side motion to free up.
5. Remove the right wall air duct by pulling it out through the door space.
6. Remove the left wall air duct by pulling it out through the door.



USER MAINTENANCE (CONTINUED)

This page left blank

USER MAINTENANCE (CONTINUED)

CALIBRATE THE TEMPERATURE DISPLAY

Note: Performing a temperature display calibration requires a temperature reference device. Please see the [Reference Sensor Devices entry](#) on page 5 for device requirements.

Temperature calibrations are performed to match the temperature display to the actual air temperature inside the oven chamber. The actual air temperature is supplied by a reference sensor device. Calibrations compensate for drifts in the controller as well as those caused by the natural material evolution of the sensor probe under temperature in the chamber space. Calibrate as often as required by your laboratory or production protocol, or regulatory compliance schedule.



Figure 6: Vent Damper in Closed Position

Vents

The vent dampeners must be close in order to perform an accurate display calibration. Chamber temperature uniformity and stability are too disrupted when operating with the vents open to verify the accuracy of the oven temperature display.

Probes

A reference device sensing probe may be introduced through the through the chamber door space. Use heat-resistant, non-marking tape to secure the wires to the body exterior and seal any gaps. The door must close and latch fully.




Figure 7: Probe End 2 inches (5cm) From Shelf Surface

Secure the wiring and cover any gaps in the door space using heat-resistant, non-stick tape.

Place the sensor probe of the temperature reference device as close as possible to the geometric center of the oven chamber. A thermocouple sensor probe sleeve may be taped to the shelving, as long as the exposed copper end is 2 inches (5cm) away from the shelf (see Figure 12). An exposed sensor probe in direct contact with the shelving may experience heat sinking, which can result in an inaccurate temperature reading.

Pre-Heating

Verify that the oven is set to your application or baking procedures set point. Let the oven heat to the set point temperature and run for at least 1 hour undisturbed in order to stabilize.

Temperature Stability

The oven chamber is thermally stabilized when no fluctuations of $\pm 2.0^{\circ}\text{C}$ or greater have been detected with the verification reference device **for a minimum of 30 minutes**. Failure to wait for stabilization will result in an inaccurate calibration and oven temperature display reading.




Procedure continued on next page

USER MAINTENANCE (CONTINUED)

Temperature Calibration

1. Once the chamber has stabilized with no fluctuations of 0.2°C or greater, compare the reference temperature device and chamber temperature display readings.
 - a. If the readings are the same, or the difference between the two (2) falls within the acceptable range of your protocol, the display is accurately showing the chamber temperature. **The Temperature Verification procedure is now complete.**
 - b. See step 2 if a difference falls outside the acceptable range of your protocol.

Reference Device




2. The display requires calibration.
 - a. The difference (also known as an error) between the reference device and the display is an **offset**.
 - b. Examples of offset values:

Reference Sensor Reading	Oven Temp. Display	Offset Value
152.0°C	150°C	2
149.1°C	150°C	-0.9
148.0°C	150°C	-2

- c. Note the offset value.

Reference Device



3. Place the controller display in the Operations selection menu.



BOTH




- a. Press and hold both the **Up** and **Down** Arrow buttons simultaneously for approximately 3 seconds.
 - b. Release the buttons when "A1" appears in the Upper Display and "Oper" appears in the Lower Display.



Operations Selection Menu



Continued on next page

USER MAINTENANCE (CONTINUED)

Temperature Calibration (Continued)	
4. Advance through the Operations menu to the Temperature Calibration page.  a. Push the green Advance button repeatedly until “i.CA” appears in the green middle display and a number value in the red top display.	 Temperature Calibration Page
  5. Adjust the number value in the top display to match the offset value, using the arrow buttons.	
6. Save the calibration offset and return to the Home Page.  a. Press the Reset Button repeatedly until the display shows the home page. b. The oven will begin to heat or passively cool to reach the current set point with the offset display value.	 Home Page
7. Allow the oven to stabilize after achieving the temperature set point using the offset display value.	

Continued on next page

USER MAINTENANCE (CONTINUED)

Temperature Calibration (Continued)	
8. Once the chamber has stabilized (no fluctuations of 0.2°C or greater) for a half hour, compare the reference temperature device and oven temperature display readings. - If the readings are the same, or the difference between the two (2) falls within the acceptable range of your protocol, the display is now accurately showing the chamber temperature. The Temperature Verification procedure is now complete. - See step 9 if a difference falls outside the acceptable range of your protocol again.	Reference Device  
9. Repeat steps 2 – 8, up to two more times. Three calibration attempts may be required to successfully calibrate ovens that are more than $\pm 2^{\circ}\text{C}$ out of calibration.	Reference Device  
10. If the temperature readings difference still falls outside your protocol after three calibration attempts, contact your distributor or Sheldon Technical Support for assistance.	

End of procedure

UNIT SPECIFICATIONS

The SMO5HP-2 is a 220 - 240 voltage unit. Please refer to the oven's data plate for individual electrical specifications.

Technical data specified applies to ovens with standard equipment at an ambient temperature of 25°C and a voltage fluctuation of $\pm 10\%$. The temperatures specified are determined in accordance to factory standard following DIN 12880 respecting the recommended wall clearances of 10% of the height, width, and depth of the inner chamber. All indications are average values, typical for units produced in the series. We reserve the right to alter technical specifications at all times.

WEIGHT

Shipping	Net Weight
425 lbs	325 lbs

DIMENSIONS

By Inches

Exterior W × D × H	Interior W × D × H
35 x 30.5 x 37.5	20.5 x 20.7 x 20

By centimeters

Exterior W × D × H	Interior W × D × H
88.9 x 77.5 x 95.3	52 x 52.7 x 50.8

CAPACITY

Cubic Feet	Liters
4.9	139

SHELF CAPACITY BY WEIGHT

Pounds	Kilograms
50 lbs. per shelf	22.8 kg per shelf

UNIT SPECIFICATIONS (CONTINUED)

TEMPERATURE

All temperature performance data for 25°C ambient conditions

Range	Uniformity	Stability
Ambient + 15 to 300°C	$\pm 1.7^{\circ}\text{C}$ @ 110°	0.1°C

Heat Up Time to 110°C	Heat Up Time to 180°C*
10 minutes	35 minutes

Recovery to 110°C*
4 minutes after closing the door

AIRFLOW PERFORMANCE


Airflow performance data is for both the intake and exhaust vents fully open.

Cubic Feet Per Minute	Liters Per Minute	Air Exchanges Per Minute
6.45	182.64	1.6

POWER

AC Voltage	Amperage	Frequency
220 - 240	12	50/60 Hz


REPLACEMENT PARTS

Description	Parts Number		Description	Parts Number
Adjustable Foot	 2700512		Shelf Slide	 5121189
Chamber Gasket Silicone (unit of sales is per foot, requires 11 feet)	 3450546		Shelf SMO5HP-2	 5121195
Door Gasket SMO5HP-2 Silicone (unit of sales is per foot, requires 11 feet)	 3450587			
Fuse 20 Amp	 3300538		Optional: Power Exhaust Blower Unit 220 - 240 Volt	 9990741

Ordering Parts and Consumables

If you have the Part Number for an item, you may order it directly from Sheldon Manufacturing by calling 1-800-322-4897 extension 3. If you are uncertain that you have the correct Part Number, or if you need that specific item, please contact Sheldon Technical Support for help at 1-800-322-4897 extension 4 or (503) 640-3000. Please have the **model number** and **serial number** of the oven ready, as Tech Support will need this information to match your unit with its correct part.





Sheldon Manufacturing Inc.

P.O. Box 627

Cornelius, Oregon 97113

EMAIL: tech@Shellab.com

INTERNET: <http://www.Shellab.com>

PHONE: 1-800-322-4897 (503) 640-3000

FAX: (503) 640-1366

