
Service Manual

Laparo CO₂ Pneu

2232

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⚠ **Important general instructions for use** ⚠

Ensure that this product is only used as intended and described in the instruction manual by adequately trained and qualified personnel, and that maintenance and repair is only carried out by authorized specialized technicians.

Operate this product only in the combinations and with the accessories and spare parts listed in the instruction manual. Use other combinations, accessories and spare parts only if they are expressly intended for this use and if the performance and safety requirements are met.

Reprocess the products before every application and before returning them for repair as required by the instruction manual in order to protect the patient, user or third parties.

Subject to technical changes!

Due to continuous development of our products, illustrations and technical data may deviate slightly from the data in this manual.

CAUTION – USA only:

Federal law restricts this unit to be used or sold, except under the supervision of a medical doctor.

Safety instructions and levels of danger

Symbol	Level of danger
	WARNING! Failure to observe can result in death or severe injury.
	CAUTION! Failure to observe can result in slight injury or damage to the product.
	IMPORTANT! Failure to observe can result in damage to the product or surrounding.
	NOTE! Tips for optimum use and other useful information.

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ENGLISH

Contents

1	General information	1
1.1	Ordering spare parts	1
1.2	Disassembly of housing cover	2
2	Maintenance	3
2.1	Important notes	3
2.2	Maintenance of device and accessories	3
2.3	Visual check	4
2.4	Electrical safety test	4
2.5	Measuring devices and other equipment	5
2.6	Functional check	6
3	Modules	10
3.1	Video PCB 2232 (Pos. 140)	10
3.1.1	Disassembly and assembly of module	10
3.2	Power supply (Pos. 160)	10
3.3	Heating PCB 2232 (Pos. 170)	10
3.3.1	Measuring points and supply voltages	10
3.3.2	Disassembly and assembly of module	10
3.4	Valve unit (Pos. 180)	10
3.4.1	Disassembly and assembly of module	10
3.4.2	Functional checks	10
3.5	Pressure reducer unit (Pos. 030 / Pos. 360)	10
3.6	Operation PCB 2232 (Pos. 260)	11
3.6.1	Disassembly and assembly of module	11
3.6.2	Functional checks	11
3.7	Control PCB 2232 (Pos. 270)	11
3.7.1	Measuring points and supply voltages	11
3.7.2	Disassembly and assembly of module	11
3.7.3	Functional checks	11
3.8	Filter PCB 2232 (Pos. 280)	11
3.8.1	Measuring points and supply voltages	11
3.8.2	Disassembly and assembly of module	12
3.8.3	Functional checks	12
3.9	Pressure reducer DM2 (Pos. 190)	12
4	Anhang / Annex	14
4.1	Reparaturteile / Repair parts	14
4.2	Explosionszeichnung Type 2232.031 / Exploded View Model 2232.031	16
4.3	Lageplan der E-Karte Video 2232 (Pos. 140) / Component drawing of PCB Video 2232 (Pos. 140)	17
4.4	Lageplan des Schaltnetzteils (Pos. 160) / Component drawing of Power Supply (Pos. 160)	18
4.5	Lageplan der E-Karte Heizung 2232 (Pos. 170) / Component drawing of PCB Heating 2232 (Pos. 170)	19
4.6	Lageplan der Bedien-E-Karte 2232 (Pos. 260) / Component drawing of PCB Operation 2232 (Pos. 260)	20
4.7	Lageplan der Steuer-E-Karte 2232 (Pos. 270) / Component drawing of PCB Control 2232 (Pos. 270)	21
4.8	Lageplan der Filter-E-Karte 2232 (Pos. 280) / Component drawing of PCB Filter 2232 (Pos. 280)	22

4.9	Blockschaltbild (2232.031) / Block diagram (2232.031)	23
4.10	Verdrahtungsplan 2232.031 / Wiring diagram 2232.031	24
4.11	Pneumatikplan (2232.031) / Pneumatic diagram (2232.031)	25
5	Protokolle / Reports	26
5.1	Prüfprotokoll / Test Report	26
5.2	Wartungsprotokoll / Maintenance report	27

1 General information

1.1 Ordering spare parts

The part number required for ordering a spare part is listed in the repair parts list under the item number used in the service manual.

IMPORTANT!

Specify the following numbers when ordering spare parts:


- *Part number of the spare part.*
- *Model/type number of the device.*
- *Serial number of the device.*

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
Exploded view

1.2 Disassembly of housing cover

◆ Required tools: Phillips screwdriver.

- ◇ Unscrew the 4 housing feet (1.2), fig.1.
- ◇ Slide trim plates (1.3) from the rails (1.4), fig.2.
- ◇ Remove the 4 bottom Phillips screws (1.5), fig.2.
- ◇ Pull the housing cover (1.0) as indicated by the arrow (a) until the tongues (1.8) are disengaged from the corresponding grooves in the front plate (1.1), fig. 3.
- ◇ Lift off the housing cover (1.0) as indicated by the arrow (b), fig. 3.
- ◇ Disconnect the earthing/ground cable (1.9), fig. 3.

2 Maintenance

2.1 Important notes

This Service Manual describes the external service measures defined for the product.

IMPORTANT!

For carrying out the service measures the product instruction manual must be observed.

2.2 Maintenance of device and accessories

To protect the technician who carries out the test, the maintenance jobs and tests must be carried out in the order specified below.

- ◇ Visual check
- ◇ Electrical safety check
- ◇ Functional check

Alternatively the electrical safety test can be carried out in accordance with DIN VDE 751, part 1.

After repair all measurement values (specified output values) of the device must be tested as described in the Service Manual, and readjusted if deviations are observed.

NOTE!

Any maintenance or testing job on the device or accessories must be documented.

2.3 Visual check

Item to be checked	Check for
Device and accessories	◆ Safety-relevant soiling / contamination and general cleanliness. ◆ Mechanical damage ◆ Loose or missing parts
Controls	◆ Mechanical function and easy operation
Labeling / Symbols	◆ Completeness and legibility
Safety-relevant labeling (e.g warnings)	◆ Completeness and legibility
Device fuses	◆ the values specified by the manufacturer on the identification plate (Nominal current and melting characteristics)
Cables	◆ for perfect condition (position, insulation and brittleness)
PCBs/cards	◆ Corrosion or other damage
	◆ Replace damaged parts immediately

2.4 Electrical safety test

Item to be checked	Check for
Protective earth connection to EN 60601-1 / IEC 601-1	Measured values: ◆ without power/mains cable: $\leq 0.1 \text{ Ohm}$ ◆ with firmly attached power/mains cable: $\leq 0.2 \text{ Ohm}$ Test conditions: $I_{\text{meas}} 25 \text{ A} \pm 10\%$, $V_0 \leq 6 \text{ V}$, $t_{\text{test}} 5 \text{ s to } 10 \text{ s}$, $50 \text{ Hz} / 60 \text{ Hz}$ To be tested: The resistance between protective earth contact or protective earth pin in the power plug and any other exposed (bare) metal part connected to the protective earth conductor.
Leakage current to EN 60601-1 / IEC 601-1	Measured values: ◆ earth leakage current: $\leq 170 \text{ } \mu\text{A}$ - limit value to EN 60601-1 / IEC 601-1, normal condition: $\leq 500 \text{ } \mu\text{A}$ ◆ patient leakage current: $\leq 10 \text{ } \mu\text{A}$ - limit value to EN 60601-1 / IEC 601-1: normal condition for BF type: $\leq 100 \text{ } \mu\text{A}$ normal condition for CF type: $\leq 10 \text{ } \mu\text{A}$ Test conditions: Measurement arrangement (MD) and measurement setup to EN 60601-1 / IEC 601-1 To be tested: The leakage current which can flow from each pole of the power supply unit through or via the insulation through the protective earth conductor or from the applied part via the patient to earth.
Additional information	Disconnect the video connection cables between the monitor and the camera controller before making the measurement.

2.5 Measuring devices and other equipment

To carry out the maintenance, repair and test routines described in this service manual, the following equipment is imperative. Before you start a repair job, ensure that the measuring equipment is calibrated and in perfect working condition.

Designation	Additional information
DC source	24 V DC power must be available
RMS multimeter	
CCD camera	with Y/C and video input
Video monitor	with Y/C- and video input (BNC and S-video)
Hygiene filter	4171.151
High flow insufflation tube	8170.232 (contains connector 64202.147)
Flow tube for CO ₂	4-40 l/min, e.g. by Platon Flow Management Jays Close, Viabes, Basingstoke RG22 4 BS UK internet: www.platon.co.uk
Digital measuring device	required resolution: 1 mbar (in pressure range below 0.5 bar) e.g.: "VDM 100" by ASF, Munich, Germany, or "GDH 14 AN" by Greisinger electronic, D-93128 Regenstauf, Germany
CO ₂ Cylinder	0.75 kg, 2046.961 For the test the cylinder pressure must be in the green range. The green LED of the cylinder pressure indicator is lit. If during the test the cylinder pressure indicator changes from green to yellow the test can still be completed without cylinder change.
High pressure connection hose	8170.801
Tummy model	must be able to withstand 30 mm HG, e.g. an empty 3 l infusion bag can be used
Syringe	50 ml
Y-connector	for silicone tubing
Silicone tubing	1 m, 4 x 2 mm


2.6 Functional check

Designation	Required checks
Preparation	Connect the gas cylinder and open the cylinder valve.
Switch-on state (self test)	During the self test: The cylinder pressure indicator blinks alternately in yellow and green and, upon successful completion of the self test the indicator light switches to green light (continuous light if cylinder is full). All level indicators show increased levels while the self test is proceeding, i.e. an increasing number of segments are lit. The self test takes max. 2 minutes (only in case of a defect). Then all LEDs are lit. At the end of the self test a brief signal is sounded. The digital display for the actual flow reads the current software version (2 digits). Upon completion of self test: Upon successful completion of the self test, the digital displays for gas volume consumed, flow and intraabdominal pressure read "0.0". The display for the gas volume consumed must not indicate any error message. The level indicators for pressure preselection and flow preselection must indicate the value last selected.
Flow control	Connect the hygiene filter with insufflation tube to the patient outlet. Preselect a nominal pressure of 25 mm Hg at the Laparo CO₂ Pneu. Press the "INSUFFLATION" button. The button must light up. Read the current flow on the digital display. Measurements: ♦ 1.1 l / min < Current (actual) flow < 1.5 l / min Connect the CO₂ flow tube (description in section 2.5) to the insufflation tube. Press the "HIGH FLOW" button. The button must light up. Preselect a nominal flow of 15 l/min . Read the current (actual) flow on the digital display and on the flow tube. Measurements: ♦ 14 l / min ≤ Current flow on digital display ≤ 16 l / min ♦ 12 l / min ≤ Current flow on flow tube. ≤ 18 l / min Preselect a nominal flow of 30 l/min. Read the current flow on the digital display and on the flow tube. Measurements: ♦ 28 l / min ≤ Current flow on digital display ≤ 32 l / min ♦ 24 l / min ≤ Current flow on flow tube ≤ 33 l / min Switch off insufflation.
Consumption indicator	Press, the "RESET" button. The digital display for the gas volume consumed reads "00.0".

Designation	Required checks
Pressure control	Preselect a nominal pressure of 15 mm Hg and a nominal flow of 10 l/min. Connect the pressure gauge (see section 2.5) and the insufflation tube via the Y-connector to the abdomen model . Actuate the "INSUFFLATION" button and then the "HIGH FLOW" button. After a short time the preselected pressure of 15 mm Hg will be observed in the abdomen model. Brief overshooting of the pressure value with subsequent pressure relief is permissible. Tolerances: ◆ Level indicators for nominal and actual pressures: + 2 segments ◆ Digital display: ± 3 mm Hg switch off insufflation briefly. Preselect a nominal pressure of 20 mm Hg and push the "INSUFFLATION" button . After a short time the preselected pressure of 20 mm Hg can be observed in the abdomen model. Brief overshooting of the pressure value with subsequent pressure relief is permissible. Tolerances: ◆ Level indicators for nominal and actual pressures: + 2 segments ◆ Digital display: ± 3 mm Hg
On-screen display (Only device versions 2232.011 and 2232.031)	Connect the camera and the monitor to the Laparo CO2 Pneu automatic insufflator. Press the "VIDEO ON-SCREEN DISPLAY" button. The button must light up. On the monitor the "Intraabdominal pressure", "flow" and "gas volume consumed" data on-screen display appears. By pressing the "VIDEO ON-SCREEN DISPLAY" button on-screen display is switched on and off. On-screen display must be checked through both a BNC cable (FBAS signal) and with an S-video cable (Y/C signal).
Gas supply monitoring	Press "RESET". Close the cylinder valve. Select a nominal pressure of 20 mm Hg and a nominal flow of 5 l/min. Switch on "HIGH FLOW" insufflation. After a short time the gas supply indicator switches from continuous green light to continuous yellow light. The low alarm is sounded. Then the gas supply indicator switches from yellow continuous light to yellow blinking light. "Medium alarm" is sounded.

Designation	Required checks
Alarm, overpressure , overpressure switch	Plug the filter and the connector (64202.147) of the high-flow insufflation tube 8170.232 onto the pressure fitting (insufflation opening) on the front panel, as shown in the drawing below. Instead of the insufflation tube which you pull off, connect the Y-shaped connector described in section 2.5 via a silicone tube (4 x 2 mm). Connect the syringe (50 ml) and the digital pressure gauge to the other outputs of the Y-shaped connector. Select a nominal pressure of 15 mm Hg. Close the pressure release opening in the lower housing panel and press the "INSUFFLATION" button. The pressure in the measurements setup now increases to 15 mm Hg. If the nominal pressure is exceeded by 4 mm Hg when actuating the syringe, the overpressure alarm is sounded. "Medium alarm" is issued. If the pressure in the measurement setup exceeds 30 mm Hg, after 5 seconds the overpressure > 30 mm Hg alarm will be sounded. At the same time "High Alarm" is sounded. The alarm LED is permanently lit. On the digital display for the gas volume consumed, error message E11 is issued and the digital display for pressure blinks. If the pressure in the measurement setup exceeds 45 mm Hg, the overpressure switch is triggered after a delay period. The "high alarm" sound then switches to continuous sound, the alarm LED is permanently lit and valve 2 switches off (switching sound).
Short-term interruption	Switch off the device for approx. 10 seconds. After switching on again any function previously preselected as well as the settings must still be active.
Factory setting	After the functional checks please return to the factory settings (basic settings). Nominal flow: 10 l / min Nominal pressure : 15 mm Hg On-screen display: ON

Diagram on overpressure simulation:

Error codes (are displayed in the 'gas volume consumed' display)

Error 01: CO₂-cylinder closed or empty

Error 02: Button depressed when switching on

Error 03: Power supply not within tolerance

Error 04: V1 defective ; flow measurement defective; V3 defective

Error 05: V2 defective, Insufflation tube connected

Error 06: Pressure measurement defective

Error 07: Flow measurement defective

Error 08: V3 defective ; V2 defective ; pressure measurement defective; Flow measurement defective

Error 09: nominal value invalid

Error 10: Helium

Error 11: Current pressure > 30 mm Hg for more than 5 Seconds

Error 12: General hardware error; pressure measurement defective; V1 defective;
flow measurement defective

Error 13: Pressure switch (input pressure) defective

Error 14: Electronic pressure switch defective

Error 15: Wrong DIL switch coding (electric pressure monitoring fast / slow) or
Pressure switch defective

3 Modules

3.1 Video PCB 2232 (Pos. 140)

3.1.1 Disassembly and assembly of module

- ◇ Install the preset card and check its complete function (in accordance with GA-A 153 and section 2.6 of this service manual).

3.2 Power supply (Pos. 160)

3.3 Heating PCB 2232 (Pos. 170)

3.3.1 Measuring points and supply voltages

- ◇ Plug X1:
 - Pin 1 = 12.5 V
 - Pin 2 = 0 V
 - Pin 3 = 12.5 V
 - Pin 4 = 0 V
- Plug X3:
 - Pin 1 = 120 V
 - Pin 2 = 230 V AC
 - Pin 3 = 120 V
 - Pin 4 = NC
 - Pin 5 = NC
 - Pin 6 = 0 V
 - Pin 7 = 230 V AC
 - Pin 8 = 0 V

3.3.2 Disassembly and assembly of module

- ◇ Install the preset card and check its complete function (in accordance with GA-A 153 and section 2.6 of this service manual).

3.4 Valve unit (Pos. 180)

3.4.1 Disassembly and assembly of module

- ◇ Open locking collar and pull off outlet tube. Release additional screw connection to the bottom panel of the housing from underneath the housing. The valves can not be replaced separately.

3.4.2 Functional checks

- ◇ After changing the module check the entire function (in accordance with GA-A 153 and section 2.6 of this service-manual).

3.5 Pressure reducer unit (Pos. 030 / Pos. 360)

◇ is delivered in preset condition.

3.6 Operation PCB 2232 (Pos. 260)

3.6.1 Disassembly and assembly of module

◇ Remove entire front frame and change PCB .

3.6.2 Functional checks

◇ Install the preset PCB and check all functions (in accordance with GA-A 153 and section 2.6 of this service manual).

3.7 Control PCB 2232 (Pos. 270)

3.7.1 Measuring points and supply voltages

◇ Plug X1:

- Pin 1 = PFD
- Pin 2 = - 12 V $\pm$ 0.4 V
- Pin 3 = GND
- Pin 4 = + 24 V
- Pin 5 = GND
- Pin 6 = + 12 V $\pm$ 0.4 V
- Pin 7 = GND
- Pin 8 = + 5 V
- Pin 9 = NC
- Pin 10 = + 5 V

Plug X3:

- Pin 1 = DS 2 (low)
- Pin 2 = DS 1 (high)
- Pin 3 = 2 V $\pm$ 0.4 V

Plug X5:

- Pin 1 = GND
- Pin 2 = + 5 V
- Pin 3 = GND
- Pin 4 = + 5 V
- Pin 5 = GND
- Pin 6 = Heating
- Pin 7 = GND
- Pin 8 = Flow > 2 l / min
- Pin 9 = GND
- Pin 10 = SDA
- Pin 11 = GND
- Pin 12 = SCL
- Pin 13 = GND

3.7.2 Disassembly and assembly of module

◇ Remove complete front frame and change PCB.

3.7.3 Functional checks

◇ Install the preset PCB and check all functions (in accordance with GA-A 153 and section 2.6 of this service manual).

3.8 Filter PCB 2232 (Pos. 280)

3.8.1 Measuring points and supply voltages

- ◇ Plug X3:
 - Pin 1 = PT100 (thermal resistor)
 - Pin 2 = PT100
 - Pin 3 = Heating
 - Pin 4 = PT100
 - Pin 5 = PT100
 - Pin 6 = Heating

3.8.2 Disassembly and assembly of module


- ◇ Remove entire front frame and replace PCB.

3.8.3 Functional checks

- ◇ Install the preset PCB and check all functions (in accordance with GA-A 153 and section 2.6 of this service manual).

3.9 Pressure reducer DM2 (Pos. 190)

- ◇ Laparo CO₂-Pneu 2232 must be switched off. Connect the device to the CO₂ gas cylinder by means of the high pressure hose.
- ◇ Open the main valve of the gas cylinder to allow CO₂ under pressure to flow to the CO₂-Pneu automatic insufflator.
- ◇ The valve block contains three valves (seen from the front of the device to the rear, in the following order: V2, V3, V1. Disconnect the two upper connectors on "V1", the lower cable (protective earth) can remain plugged.
- ◇ Pull off both silicone tubes 2 x 1.5 mm at V3 for the pressure sensors from connection **A** and connect digital pressure gauge (see 2.5) to the free connections via the Y-shaped connector. The Y-shaped connector is shown below:

- ◇ Pull off the front silicone tube 4 x 2 mm from connection **B** and close safely. Caution : if you do not close the tube the flow sensor may become damaged during the following test.
- ◇ Apply 24 V DC from a DC power source to the free contacts of V1 (Polarity has no importance).
- ◇ Close the front connection **B** (previously freed from the silicone tube 4 x 2 mm) with your finger or some other means.
- ◇ Read pressure gauge. Turn plastic bolt on the top of DM 2 until a pressure of 93 ± 2 mbar is built up in a stable manner.
- ◇ Open and close the front connection **B** several times. Check that the pressure is readjusted correctly.
- ◇ Disconnect DC power from V1. Only then (!) should the closed tube to the flow sensor be opened again.
- ◇ Reconnect all tubes to their original connectors.

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RICHARD WOLF Belgium S.A.
RICHARD WOLF France S.A.R.L.