


ARMEN SKID LC SYSTEM



Custom solution for preparative and industrial HPLC

Armen Instrument : ZI de Kermelin - 16, rue Ampère - 56890 St-Avé - France
Phone : +33(0)2 97 61 84 00 – Fax : +33 (0)2 97 61 85 00
contact@armen-instrument.com; www.armen-instrument.com

Example of Armen Skid LC system

AP-HPG2-500-125-SKID




Configuration



- Binary High Pressure gradient pump, 500 ml/min, 125 bar max
- Manual injector with loop
- Injection through pump unit
- Fraction collector 20 ways and 1 waste
- Priming valve
- Back flush valve unit
- Recycling valve unit
- UV/VIS with prep cell
- RI detector
- Armen Glider Prep software
- Computer 17" flat screen.
- Power supply cabinet including components

Specifications/Performances


- Flow rate range : 1 to 500 mL/min
- Flow rate performance : 1% accuracy (with H₂O degassed at 20 °C),
- Repeatability better than 0.5%.
- Gradient former : Linear response from 0.5 to 99.5%
- System protection : Pressures Min and Max adjustable, Software
- Control : via PC & software
- Dimensions : (Width x Height x Depth -in millimetres- without connections)
- Trolley : 1270 x 1281 x 950
- Weight : 250 kg
- Temperature : 5 °C to 40 °C



Linear and step gradient at 200 mL/min



Linear and step gradient at 500 mL/min



Armen Instrument : ZI de Kermelin - 16, rue Ampère - 56890 St-Avé - France

Phone : +33(0)2 97 61 84 00 – Fax : +33 (0)2 97 61 85 00

contact@armen-instrument.com; www.armen-instrument.com

AP-HPG2-TRIX-5000-70-SKID




Configuration



- AP-HPG2 binary High Pressure gradient pump, 5000 ml/min, 70 bar max
- Priming valve
- Injection through pump unit
- Fraction collector 4 ways and 1 waste
- Back flush valve unit
- Recycling valve unit
- UV/VIS with prep cell
- Armen Glider Prep software
- Computer 17" flat screen. Soft installed and tested.
- Power supply cabinet including components
- Flow rate range :

Specifications/Performances


- Flow rate range : 10 to 5 000 mL/min
- Flow rate performance : 1% accuracy (with H₂O degassed at 20 °C),
- Repeatability better than 0.5%.
- Gradient former : Linear response from 0.5 to 99.5%
- System protection : Pressures Min and Max adjustable, Software
- Control : via PC & software
- AC mains supply voltage : 220/240V – 20A max - 50 Hz
- Dimensions : (Width x Height x Depth -in millimetres- without connections)
- Trolley : 1270 x 1281 x 950
- Weight : 250 kg



Linear and step gradient at 2000 mL/min



Linear and step gradient at 5000 mL/min



Armen Instrument : ZI de Kermelin - 16, rue Ampère - 56890 St-Avé - France

Phone : +33(0)2 97 61 84 00 – Fax : +33 (0)2 97 61 85 00

contact@armen-instrument.com; www.armen-instrument.com

Custom system on request :

Armen SKID LC systems could be set up according to your needs and application. Large range of pumping system are available from 500 mL/min to 15 L/min, isocratic, low pressure or high pressure gradient and all kind of detectors, fraction collectors and hydraulic valves could be installed. All those peripherals are installing in one box and controlled through our dedicated Armen Glider Software for GMP compliance and convenient purification.

Design your system by indicated which peripherals you would like for a quote.

Pumping system :

AP-TRIX range : Isocratic or low pressure gradient mixing

500 mL/min pump - 250 bar - Isocratic	<input type="checkbox"/>
1000 mL/min pump - 200 bar - Isocratic	<input type="checkbox"/>
2000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
3000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
5000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
10000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
15000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
Quaternary gradient for AP-TRIX	<input type="checkbox"/>

AP-HPG2-TRIX range : High pressure binary gradient mixing

2 x 500 mL/min pump - 125 bar	<input type="checkbox"/>
2 x 1000 mL/min pump - 125 bar	<input type="checkbox"/>
2 x 2000 mL/min pump - 100 bar	<input type="checkbox"/>
2 x 3000 mL/min pump - 70 bar	<input type="checkbox"/>
2 x 5000 mL/min pump - 100 bar	<input type="checkbox"/>
2 x 6000 mL/min pump - 100 bar	<input type="checkbox"/>

Hydraulic devices :

Manual solvent selectors	<input type="checkbox"/>
Manual purge valve	<input type="checkbox"/>
Recycling valve	<input type="checkbox"/>
Backflush valve	<input type="checkbox"/>
Column selection valve (for 2 columns connection)	<input type="checkbox"/>

Injection unit :

Injection through external pump

50 mL/min pump - 300 bar - Isocratic	<input type="checkbox"/>
100 mL/min pump - 250 bar - Isocratic	<input type="checkbox"/>
250 mL/min pump - 200 bar - Isocratic	<input type="checkbox"/>
500 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
1000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
2000 mL/min pump - 100 bar - Isocratic	<input type="checkbox"/>
3000 mL/min pump - 70 bar - Isocratic	<input type="checkbox"/>

Loop injection through manual valve

Injector for AP-HPG2-250/500/1000	<input type="checkbox"/>
1 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
2 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
5 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
10 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
20 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
30 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
40 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>
50 mL stainless steel Loop, 1/8" OD	<input type="checkbox"/>

Injection through pump unit with AP-HPG2 pumps

Injection unit for AP-HPG2-250	<input type="checkbox"/>
Injection unit for AP-HPG2-500	<input type="checkbox"/>

Detection :

UV 06S SINGLE, standard filter 254 nm, other filter on request	<input type="checkbox"/>
UV 06S DUAL, standard filters 254 and 280 nm, other filters on request	<input type="checkbox"/>
UV 06S DAD 400, 200-400 nm, four wavelengths simultaneously, scan of spectrum	<input type="checkbox"/>
UV 06S DAD 600, 200-600 nm, four wavelengths simultaneously, scan of spectrum	<input type="checkbox"/>
UV 06S DAD 800, 200-800 nm, four wavelengths simultaneously, scan of spectrum	<input type="checkbox"/>
RI Detector KNAUER S2300	<input type="checkbox"/>
Evaporative Light Scattering Detector "ELSD" SEDERE SEDEX Model 80LT	<input type="checkbox"/>

Fraction collector :

Number of ways needed	<input type="checkbox"/>
-----------------------	--------------------------

General features :

Hardware interfaces

PC to "Armen" instruments - By RS232
Between "Armen" instruments - Bus Can
PC to UV detector - By RS232

Software interfaces


Controlling "Armen Glider Prep" software (AGP)

A supervision screen allows following and controlling all operating parameters :

- Pump flow rate, pressure, high and low pressure safety
- Sample volume or injection time
- Gradient former and all switching valves
- Detectors and data acquisition
- Fractions collector drove by time, level and derivate peak detection
- Level of the solvent, fractions and waste reservoirs.
- Methods and data acquisitions are saved

Recording CFR 21 part 11 compliant "Armen Glider Prep" list :

- Access control : Software access is limited to authorized only persons through a login/password mechanism
- Access to each screen is restricted following the user profile scheme
- Each event below is recorded (times tamped in a protected database)
 - o User Login/Logout
 - o Method Ran/Stopped
 - o Device communication errors
 - o Times stamped records of every event: RUN, STOP and PAUSE...
 - o Data records with all parameters every second (AU, flow rate, gradient, pressure, collector position, wavelength...)
 - o Collection information (timestamp, seal number, peak number, volume, area)
- Three files are generated while running methods (ASCII files suitable for EXCEL)
- Software versioning
- Documentation (User Requirements definition. User manual...)
- Automated installation procedure
- Software engineering practices : Integrated testing and debugging procedures



Regulation

- CE
- Directive C. E. M. n° 89/336/CEE modified by directive n° 92/31/CEE and 93/68/CEE.
- Directive Low voltage n° 73/23/CEE modified by Directive n° 93/68/CEE and according to the following norms : Norm CEM : EN55011 amendment 1 & 2, EN61000-4-2 amendment 1, EN61000-4-3 amendment 1 & EN61000-4-4.
- Norm electrical safety : EN61010-1 amendment 2 – 07/95.
- FDA
- "GMP" - Good Manufacturing Practice (added note) :
 - o The conception of this instrument was done to avoid any retention volume in the liquid circuit.
 - o A user manual and a service manual are delivered with the instrument, giving all necessary information to use the instrument, to replace the seals of the pumps and check valves. It's explained in detail how to clean every parts of the instrument.
 - o A practice training will be done during the IQ / OQ.
 - o The materials in contact with the liquids are : Stainless steel 316L, PTFE bonded with Glass fibers, PTFE bonded with carbon fibers, FEP, Kel'F, PVDF, PEEK, Kalrez, Sapphire, Rubies, Hastelloy, Zirconium

Armen Instrument : ZI de Kermelin - 16, rue Ampère - 56890 St-Avé - France

Phone : +33(0)2 97 61 84 00 – Fax : +33 (0)2 97 61 85 00

contact@armen-instrument.com; www.armen-instrument.com