GE Healthcare

Optima XR200amx/XR220amx Wireless Handswitch Installation

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WARNING

X-RAY EQUIPMENT IS DANGEROUS TO BOTH PATIENT AND OPERATOR UNLESS MEASURES OF PROTECTION ARE STRICTLY OBSERVED.

Though this equipment is built to the highest standards of electrical and mechanical safety, the useful x-ray beam becomes a source of danger in the hands of the unauthorized or unqualified operator. Excessive exposure to x-radiation causes damage to human tissue.

Therefore, adequate precautions must be taken to prevent unauthorized or unqualified persons from operating this equipment or exposing themselves or others to its radiation.

Before operation, persons qualified and authorized to operate this equipment should be familiar with the Recommendations of the International Commission on Radiological Protection (ICRP), contained in Annals Number 26 of the ICRP, and with applicable national standards.

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Legal notes

TRADEMARKS

Optima XR200amx™ and Optima XR220amx™ are trademarks of GE Healthcare. WPA and WPA2 are trademarks of the Wi-Fi Alliance. All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All material copyright © 2014 by General Electric Company, Inc. All rights reserved. The material presented and contained herein may not be reproduced in any form or manner, without the written permission of General Electric Company, Inc.

Important precautions

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report, call 1-800-548-3366. Select the option for "Install Support Services for FOA and MIS." Contact your local service coordinator for more information on this process.

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Healthcare personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. GE Healthcare will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements (NCRP) as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that GE Healthcare, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment. Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENT

DANGER - Risk of Explosion

Danger of explosion if battery is incorrectly replaced.

Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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TECHNICAL MANUAL UPDATES

When operating or servicing GE Healthcare products, please contact your GE representative for the latest revision of product documentation. Product documentation may also be available on-line at the GE Healthcare support documentation library.

OMISSIONS AND ERRORS

Customers, please contact your GE Healthcare sales or service representatives.

GE personnel, please use the GE Healthcare complaint handling process to report all omissions, errors, and defects in this publication.

Revision history

Revision	Date	Reason for change
1	24MAR2014	Initial release of document.
2	25APR2014	Validation update: Swapped and edited steps 3 and 4 in Chapter 2.

Preface Publication Conventions

Standardized conventions for representing information is a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

It is important that you read and understand hazard statements, and not just ignore them.

Safety & Hazard Information

Proper product safety labeling allows a person to safely use or service a product. The format and style for safety communications reflected in this publication represents the harmonization of IEC/ISO 3864 and ANSI Z535 standards.

Within this publication, different paragraph and character styles are used to indicated potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that are applicable to each specific safety statement.

Hazard Messages

Any action that will, could or potentially cause personal injury will be preceded by the safety alert symbol and an appropriate signal word. The safety alert symbol is the triangle with an exclamation mark within it. It is always used next to the signal word to indicate the severity of the hazard. Together, they are used to indicate a hazard exists.

Signal words describe the severity of possible human injuries that may be encountered. The alert symbol and signal word are placed immediately before any paragraph they affect. Safety information includes:

- 1.) Signal Word - The seriousness level of the hazard.
- 2.) Symbol or Pictorial - The consequence of interaction with the hazard.
- 3.) Word Message:
 - a.) The nature of the hazard (i.e. the type of hazard)
 - b.) How to avoid the hazard.

The safety alert symbol is not used when an action can only cause equipment damage.

Text Format of Signal Words

DANGER - INDICATES AN IMMINENTLY HAZARDOUS SITUATION WHICH, IF NOT AVOIDED, WILL RESULT IN DEATH OR SERIOUS INJURY. THIS SIGNAL WORD IS LIMITED TO THE MOST EXTREME SITUATIONS.

WARNING - INDICATES A POTENTIALLY HAZARDOUS SITUATION WHICH, IF NOT AVOIDED, COULD RESULT IN DEATH OR SERIOUS INJURY.

Caution - Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices.

NOTICE - Indicates information or a company policy that relates directly or indirectly to the safety of personnel or protection of property. This signal word is associated directly with a hazard or hazardous situation and is used in place of 'DANGER,' 'WARNING,' or 'CAUTION.' It can include:

- Destruction of a disk drive
- Potential for internal mechanical damage, such as to a X-ray tube

Symbols and Pictorials Used

The following Symbols and Pictorials are be used in this publication. These graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

NOTICE	CAUTION	WARNING	DANGER	
 keep_up	 magnetic	 biohazard	 compressgas	 ppe-hearing
 fragile	 impact	 corrosive	 heavyobject	 ppe-2people
 static_elec	 heat	 general	 laser	 ppe-respiratory
 keep_dry	 pinch	 radiation	 poisongas	 ppe-loto
 general	 explosive	 electrical	 flammable	 ppe-eye
 torque	 crush/mechanical	 tipping	 Read Manual	 ppe-gloves
 ce	 instuction	 poisonmatl	 entanglement	 instuction

Publication Conventions

General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as Purpose, Example, Comment or Note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION:

Purpose: Introduces and provides meaning as to the information contained within the chapter, section or subsection (such as used at the beginning this chapter, for example).

Note: Conveys information that should be considered important to the reader.

Example: Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

Comment: Represents “additional” information that may or may not be relevant to your situation.

Page Layout

Publication Part Number & Revision Number

Publication Title

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document part number and revision number appears in every header on every page. Odd numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well as the current page number.

Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono-spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operator "greater-than" and "less-than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example: Fixed
Output

This paragraph's font represents computer generated screen "fixed" output. Its output is fixed from the sense that it does not vary from application to application. It is the most commonly used style used to indicate filenames, paths and text that do not change from system to system. The character style used is a fixed width such as courier.

Example:
Variable Output

This paragraph's font represents computer screen output that is "variable". It is used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater-than and less-than operators for clarification. For example: <variable_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.

Example: Fixed
Input

This paragraph's font represents fixed input. It is computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p

Example:
Variable Input

This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as:

```
ypcat hosts | grep 3.45.120.3
```

without the greater-than and less-than operators.

Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example: Hard
Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example: Soft
Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

Table of contents

Chapter 1	Safety	3
	Section 1.1 Energy sources	3
	Section 1.2 Lock Out/Tag Out (LOTO) procedure for electrical power	4
	1.2.1 Preparing for LOTO.....	4
	1.2.2 Performing LOTO.....	6
	1.2.3 Returning system to service.....	9
	Section 1.3 Electrostatic discharge (ESD)	10
	1.3.1 Important ESD considerations when working on a mobile system.....	10
	1.3.2 Generating static.....	10
	1.3.3 Personal grounding methods and equipment	11
	1.3.4 Grounding the work area.....	11
	1.3.5 Recommended materials and equipment	12
Chapter 2	Installation procedure	13

Chapter 1 Safety

Section 1.1 Energy sources

Table 1-1

Energy source	Yes/No	Location of energy isolating means	Magnitude of energy
Electrical	Yes	AC power plug	120 VAC/60Hz 220 VAC/50Hz
Electrical	Yes	Circuit breaker	156VDC
Electrical (batteries)	Yes	Battery board connector J1	156VDC
Pneumatic	No	n/a	n/a
Hydraulic	No	n/a	n/a
Gas/water/steam	No	n/a	n/a
Chemical	No	n/a	n/a
Mechanical motion	No	n/a	n/a
Gravity	No	n/a	n/a
Rear suspension spring	No	n/a	n/a
Column spring	Yes	Scroll lock at top of column	610 N-m (5400 in-lbs)
Thermal	No	n/a	n/a
Stored energy (generator capacitors)	Yes	Time discharge and DVM verification	15kW generator: 2200J 30kW generator: 8800J
Air under pressure	No	n/a	n/a
Oil under pressure	No	n/a	n/a
Water under pressure	No	n/a	n/a
Gas under pressure	No	n/a	n/a
Steam	No	n/a	n/a
Other	No	n/a	n/a

Type of equipment and/or method selected to dissipate or isolate stored energy:

- Allow three minutes for stored energy to dissipate. Stored energy is indicated by safety lamps on the generator Aux boxes.
- After waiting period, verify dissipation with volt meter.

Type of equipment and/or method used to ensure disconnections:

- One lock and tag for each electrical power source
- Multiple locking device

Section 1.2 Lock Out/Tag Out (LOTO) procedure for electrical power

Name of equipment: Optima XR200amx/Optima XR220amx

Number of locks: 3

Title of employees authorized to perform LOTO: Those trained in LOTO

Title of affected employees and how to notify: Hospital personnel, notified by verbal communication

1.2.1 Preparing for LOTO

1. Have this equipment on hand:
 - Safety glasses
 - Voltage meter
 - LOTO kit 5421892, contains (see [Figure 1-1](#)):
 - * 3 brass padlocks and identification labels
 - * 3 LOTO tag packages (25 tags per package)
 - * 1 electrical plug lockout device 110V - for AC power plug
 - * 2 electrical plug lockout devices 160VDC - for battery connectors
 - * 1 nylon carry bag

Figure 1-1 LOTO kit contents

2. Locate the AC power plug (E1), system circuit breaker (E2), and battery connectors (E3). See [Figure 1-2](#).


Figure 1-2 Electrical energy lockout locations

1.2.2 Performing LOTO

1. Notify all affected personnel working in the area that LOTO is being performed.
2. Exit all system software.
3. Turn off the Optima XR200amx/Optima XR220amx system by pressing the power button and holding it for several seconds until the power button light ring turns blue, indicating that the system is shutting down.
4. When the power light ring begins to blink blue, indicating that the system is in standby, switch the circuit breaker (E2) to the OFF position. See [Figure 1-3](#).

Figure 1-3 Circuit breaker (E2) in OFF position

5. Unplug the power cord from the wall outlet.
6. Apply an approved AC power plug locking device to the AC power plug (E1), then apply your personal red lock and tag. See [Figure 1-4](#).



Figure 1-4 AC power plug with locking device applied

7. Wait three minutes for discharge of stored energy.


8. Disconnect the two green battery power connectors J3 and J4 and the two sense connectors J13 and J14 from the Cricket board. Apply a locking device, then apply your personal red lock and tag. See [Figure 1-5](#).

Item	Description
1	Battery power connectors (J3, J4)
2	Battery sense connectors (J13, J14)

Figure 1-5 Battery connectors J3/J4/J13/J14 and locking device

9. Disconnect the two green battery power connectors J1 and J2 and the two sense connectors J11 and J12 from the Cricket board. Apply a locking device, then apply your personal red lock and tag. See Figure 1-6.


Item	Description
1	Battery power connectors (J1, J2)
2	Battery sense connectors (J11, J12)

Figure 1-6 Battery connectors J1/J2/J11/J12 and locking device

10. Press the power button and verify that the system does not turn on.
11. Verify that the voltage meter is working properly by testing it on a known live voltage source (such as a wall outlet).

12. See [Figure 1-7](#). Verify that electrical power has been removed from the system by applying the DC voltmeter (+) lead to connector J5 pin 1 (Item 1, left-side screw terminal) on the Cricket board and the meter (-) lead to chassis ground (Item 2).

Item	Description
1	Connect DC voltmeter (+) to J5 pin 1 (left-side screw terminal)
2	Connect DC voltmeter (-) to chassis ground

Figure 1-7 Cricket board test points

**ELECTRICAL SHOCK HAZARD
IF THE SYSTEM CANNOT BE LOCKED OUT, OR IF THE SYSTEM FAILS
VERIFICATION, CONTACT YOUR SUPERVISOR.**

13. Perform the necessary service or repair.

1.2.3 Returning system to service

1. Notify affected personnel that LOTO devices are being removed and equipment is being re-energized.
2. Verify that the area is safe to re-energize the equipment.
3. Verify that any safety guards or devices have been reinstalled, and replace all guards and covers as necessary.
4. Remove tools and all non-essential equipment.
5. Verify that all personnel are clear and in a safe position.
6. Verify that all controls and switches are in the off or neutral position.
7. Remove locks and tags as appropriate.

8. Re-energize equipment and follow safe startup procedures.
9. Notify affected persons that energy has been restored.

Section 1.3 Electrostatic discharge (ESD)

A sudden discharge of static electricity from your finger or other conductor can destroy static sensitive devices or microcircuitry. Often the spark is neither felt nor heard, but damage occurs. An electronic device exposed to electrostatic discharge (ESD) might not appear to be affected at all and can work perfectly throughout a normal cycle. The device can function normally for a while, but it has been degraded in the internal layers, reducing its life expectancy.

Networks built into many integrated circuits provide some protection, but in many cases, the discharge contains enough power to alter device parameters or melt silicon junctions.

1.3.1 Important ESD considerations when working on a mobile system

NOTICE

Static Discharge Risk to Components

If you are using a static mat or wrist strap, be sure to connect the mat and strap to machine ground, not earth ground. This is necessary to isolate yourself from earth ground and equalize your potential with the machine.

ESD damage prevention includes these steps:

1. Perform LOTO.
2. Verify that the AC plug is not connected to an energy source/outlet.
3. Connect ground mat and wrist strap assembly to frame connection point of the mobile unit.

1.3.2 Generating static

The following table shows that different activities generate different amounts of static electricity, and that static electricity increases as humidity decreases.

Table 1-2

Event	Relative humidity		
	55%	40%	10%
Walking across carpet	7,500 V	15,000 V	35,000 V
Walking across vinyl floor	3,000 V	5,000 V	12,000 V
Motions of bench worker	400 V	800 V	6,000 V
Removing bubble pack from PCB	7,000 V	20,000 V	26,500 V
Packing PCBs in foam-lined box	5,000 V	11,000 V	21,000 V

Many electronic components are sensitive to ESD. Circuitry design and structure determine the degree of sensitivity. The following packaging and grounding precautions are necessary to prevent damage to electric components and accessories.

- Transport products in static-safe containers to avoid hand contact.

- Protect all electrostatic parts and assemblies with conductive or approved containers or packaging.
- Keep electrostatic sensitive parts in their containers until they arrive at static-free stations.
- Place items on a grounded surface before removing them from their container.
- When handling or touching a sensitive component or assembly, ground yourself by touching the chassis.
- Avoid contact with pins, leads, or circuitry.
- Place reusable electrostatic-sensitive parts from assemblies in protective packaging or conductive foam.

1.3.3 Personal grounding methods and equipment

⚠ WARNING

ELECTRICAL SHOCK HAZARD

This section defines a workzone that is not connected to the system under service. The mobile system is a battery powered unit that can deliver lethal current. Ensure you are not providing an earth ground path when practicing ESD procedures.

Use the following equipment to prevent static electricity damage to equipment:

- Wrist straps are flexible straps with a maximum of one-megohm +/- 10% resistance in the ground cords. To provide a proper ground, wear the strap against bare skin. The ground cord must connect snugly into the banana plug connector on the grounding mat or workstation.
- Heel, toe, and boot straps can be used at standing workstations and are compatible with most types of shoes or boots. On conductive floors or dissipative floor mats, use them on both feet with a maximum of one-megohm +/- 10% resistance between the operator and ground.

1.3.3.1 Static shielding protection levels

Table 1-3

Method	Voltage
Antistatic plastic	1,500
Carbon-loaded plastic	7,500
Metalized laminate	15,000

1.3.4 Grounding the work area

⚠ WARNING

ELECTRICAL SHOCK HAZARD

This section defines a workzone that is not connected to the system under service. The mobile system is a battery powered unit that can deliver lethal current. Ensure you are not providing an earth ground path when practicing ESD procedures.

To prevent static damage of components and parts in the work area:

- Cover the work surface with approved static-dissipative material. Provide a wrist strap connected to the work surface and properly grounded tools and equipment.
- Use static-dissipative mats, foot straps, or air ionizers to give added protection.
- Handle electrostatic sensitive components, parts, and assemblies by case or PCB laminate.
- Handle components only at static-free work areas.
- Turn off power and input signals before inserting and removing connectors or test equipment.
- Use static-safe fixtures when fixtures must directly contact dissipative surfaces.
- Keep work area free of non-conductive materials, like plastic assembly aids and Styrofoam.

1.3.5 Recommended materials and equipment

ELECTRICAL SHOCK HAZARD

This section defines a workzone that is not connected to the system under service. The mobile system is a battery powered unit that can deliver lethal current. Ensure you are not providing an earth ground path when practicing ESD procedures.

Materials and equipment that are recommended for use in preventing static electricity include:

- Anti-static tape, smocks, aprons, or sleeve protectors
- Conductive bins, foam, tabletop workstations with ground cord of one-megohm +/-10% resistance, and other assembly or soldering aids
- Static-dissipative table or floor mats with hard tie to ground
- ESD field service kits
- Static awareness labels
- Wrist straps and footwear straps providing one-megohm +/-10% resistance
- Material handling packages:
 - Conductive plastic bags, plastic tubes and tote boxes
 - Opaque shielding bags
 - Transparent metalized shielding bags and shielding tubes

Chapter 2 Installation procedure

Contact a GE service representative to install the wireless handswitch.
Leave the wired handswitch on the system as a backup.

1. Remove both side covers and the top cover. See System manual 5336122-1EN, Chapter 8, Section 8.1.2 - Top cover removal.
2. Identify the components of the kit shown below ([Table 2-1](#), [Figure 2-1](#), [Figure 2-2](#)):

Table 2-1 Wireless handswitch kit 5503600

Part Number	Description
5503600:	Wireless Handswitch kit, contains the following items:
5503601	Wireless Handswitch Collector (see Figure 2-1): PCB Module (with cable) Spare transmitter Remote Finder/Transmitter Cradle Receiver (with cable)
5486733	Receiver Mount (see Figure 2-2)
5497064	Self-tapping Screw, KA35 X 10mm, Phillips Pan Head (not shown)
5497362	Drill Bit, 0.5 inch with 0.375 inch shank (not shown)
5497364	Drill Bit, #22 (not shown)
5498227-1EN	Wireless Handswitch Installation Manual (this document)
46-302200P7	Rating Plate (not shown)


Figure 2-1 Wireless Handswitch Collector - 5503601


Figure 2-2 Receiver mount

3. Remove the existing nut and washer located behind where the PCB Module is to be placed, then attach the ground wire from the PCB Module to the system chassis stud and re-install the washer and nut. See [Figure 2-3](#).


Figure 2-3 Attach ground wire

4. Mount the PCB Module in the bottom of the Thorax as shown below ([Figure 2-4](#)) by removing the adhesive liner and applying pressure to activate the adhesive.


Figure 2-4 Mount PCB module in bottom of Thorax

5. Route the cable attached to the PCB Module up along the main cable bundle between the Spyder and Locust boards. See [Figure 2-5](#).


Figure 2-5 Route cable

6. Disconnect the J5 plug from the Spyder board and plug into End C of the main cable from the PCB Module.
Plug End B of the main cable into Spyder J5.
See [Figure 2-6](#) to see how connections appear after re-wiring is complete.


Figure 2-6 Cable connection to Spyder board

7. Remove the top-left mounting fastener on the Spyder board and mount the ground terminal, End D, of the main cable underneath the head of the fastener. See [Figure 2-7](#).


Figure 2-7 Ground connection

8. Continue routing the main cable (twisted black and red wires) up through the Thorax ([Figure 2-8](#)) to the LVLE2 and plug End A of the main cable into J6 on the LVLE2 ([Figure 2-9](#)).


Figure 2-8 Route red/black wires


Figure 2-9 Connect cable to LVLE2 J6

9. Position the Receiver Mount on the top cover. The mount is curved and will fit the contour of the top cover when it is located in the proper position. See [Figure 2-10](#) for approximate location.

Figure 2-10 Receiver mount location

10. While holding the Receiver Mount securely in position, mark the locations for the 3 mounting holes and the cable access hole. See [Figure 2-11](#).


Figure 2-11 Mark holes

11. Drill the holes in the top cover using the drill bits provided in the kit. See [Figure 2-12](#).

Figure 2-12 Drill holes

12. Temporarily remove the two nylon nuts from the standoffs on the Receiver. Run the cable of the Receiver through the $\frac{1}{2}$ " diameter hole in the Receiver Mount, then place the Receiver onto the Receiver Mount and reinstall the nylon nuts to secure it to the mount. Do not over-tighten the plastic nuts. See [Figure 2-13](#).

Figure 2-13 Assemble receiver and receiver mount

13. Use a pair of cutters to remove the excessive length of the standoffs so that they are flush with the surface of the Receiver Mount. See [Figure 2-14](#).


Figure 2-14 Trim standoffs

14. Run the cable of the Receiver through the 1/2" diameter hole in the top cover.
15. Attach the Receiver/Mount assembly to the top cover using the 3 self-tapping screws provided in the kit. See [Figure 2-15](#). Do not over-tighten the screws.


Figure 2-15 Attach Receiver/Mount assembly to top cover

16. Route the cable along the power switch cable of the top cover. See [Figure 2-16](#). Zip tie the cables as necessary.


Figure 2-16 Route the cable

17. Continue routing the cable along the side of the PC and down through the hole in the sheet metal of the Thorax and between the Locust and Spyder boards. Plug the cable into the connector on the PCB Module. See [Figure 2-17](#). Zip tie the cables as necessary.


Figure 2-17 Route cable to PCB Module

18. Re-install the top cover and the side covers.
19. Mount the cradle (in the customer's desired location) by removing the adhesive liner and applying pressure to activate the adhesive.
20. Power up the system.
21. Press the collimator light button on the wireless handswitch.
If the collimator light on the system goes on, the wireless handswitch is ready for use.
22. Confirm that the system will take exposures using the wireless handswitch.
23. Apply the rating plate:
 - a. Apply the rating plate label to front cover, next to the power cord area. See [Figure 2-18](#).


Figure 2-18 Apply rating plate

- b. Fill out and submit the product locator card(s).

© 2014 General Electric Company.
GE Medical Systems, a General Electric Company, going to market as GE Healthcare.
3000 N. Grandview Boulevard
Waukesha, Wisconsin 53188
USA
www.gehealthcare.com