

# SERVICE MANUAL

for **Colenta**<sup>®</sup>

## MEDIPHOT 900E

# X-Ray-Filmprocessor

Warnung: Diese Anleitung ist nur für qualifizierte Service Techniker bestimmt.  
Warning: For the use of qualified service personnel only.  
Avertissement: Réservé au personnel de service qualifié.



**Colenta**<sup>®</sup>

MMP89 V2.4 and up



02/2006 FR

# Colenta<sup>®</sup> MEDIPHOT 900E

## X-Ray-Filmprocessor

### Technical Specifications


Processing applications:	Cut sheets of all commonly used X-Ray films incl. Mammography	
Film thickness:	min. 0,10 mm	
Material width:	min. 10 cm (4") - max. 36 cm (14")	
Material length:	min. 10 cm (4")	
Time in Developer:	min. <b>15</b> sec	max. <b>84</b> sec
Intake Speed:	76cm/min (31inch/min)	14cm/min (14inch/min)
Cycle Time:	82 sec	459 sec
	(see the table on one of the next page)	

		<b>MP 900 E</b>
<b>Tank capacity -</b>	Developer	2,6 L
	Fixer:	2,5 L
	Wash water:	3 L
Solution heating(Fix and DEV):	variable in a range of 18°C - 43°C (350W Inline Heater)	
Dryer:	warm air	
Replenishment:	fully automatic. replenishment is microprocessor controlled and calculated from information received from sensors measuring the width and length of material entering the processor. Replenishment cycles are variable.	
power supply:	1 / N / PE~ 230V (+6% / -10%), 50Hz, 13A, 1.7KW (alternatives are available on request)	
water supply:	2-way magnetic valve, with 3/4" hose connection by using a DVGW-approved system-separating device or pipe-separating device.	
Wash water flow rate:	1,5 ltr/min when film is in process	
Wash water supply pressure:	3 - 10 bar	
Wash water supply:	filtered at a temperature of 10°C - 15°C	
		<b>MP 900E</b>
<u>Weight:</u>	Empty	58 kg
	With solution	66 kg

Technical specification subject to change without notice.

**INPUT SPEED / DEV-TIME / CYCLE-TIME:**


**Colenta**<sup>®</sup>  
**MEDIPHOT**  
**900E**



Dev. time [ s ]	Intake speed [ cm/min ]	Cycle time [ s ]
14.6	78.75	80
15.5	74.12	85
16.5	70.00	90
17.4	66.32	95
18.3	63.00	100
19.2	60.00	105
20.1	57.27	110
21.0	54.78	115
21.9	52.50	120
23.8	48.46	130
25.6	45.00	140
27.4	42.00	150
29.3	39.38	160
31.1	37.06	170
32.9	35.00	180

# MEDIPHOT 900E

## Dimensions of MEDIPHOT 900E:



## FACTORY SETTINGS: MP900E

<b>PROCESSOR: COLENTA MP 900 E</b>	
<b>with PDB</b>	<b>MINI 04 2004 P89C668</b>
<b>Displayanzeige bei Initialisierung:</b>	P89c668 MM_V2.4

<b>Program</b>	<b>P1</b>	<b>P2</b>	<b>RANGE</b>
Tank1 time [sec]	023	018	015 - 084
T1 [°C]	33,0	33,0	18,0 - 43,0
R1 [ml/sqm]	0400	0400	0040 - 2000
	<b>P3</b>	<b>P4</b>	<b>P5</b>
Tank1 time [sec]	030	030	030
T1 [°C]	33,0	33,0	33,0
R1 [ml/sqm]	0400	0400	0400
<b>Setup [Werte für HA - Z12 MINI - MEDICAL]</b>		<b>(Setup Bridge necessary)</b>	
Tank1 [pls]	0042		0010 - 1000
Machine [pls]	0248		0030 - 9999
Gear [pls/m]	219		050 - 999
Pump [ml/s]	02,6		00,5 - 99,9
Sensor distance [mm]	060		000 - 300
Replenish after each [sq.m.]	0,125		0,100 - 0,999
Ready after film left [pls.]	0030		0000 - 1000
Add before start [pulses]	14		00 - 99
Additional run [cm]	001		000 - 200
<b>Standby</b>			
SB replenishment [ml]	100		000 - 999
each [h]	6		0 - 9
SB self-cleaning [cm]	20		05 - 40
each [min]	10		01 - 20
Switch ON Dryer for [min]	05		00 - 20
DR On each [min]	10		00 - 50
for [min.]	01		01 - 10

## MEDIPHOT 900E

<b>Service</b>		
<b>Inputs</b>		
Cov=	0 = Cover closed 1 = Cover open	Cover switch
Lev=	0 = closed 1 = open	Option level switch
SP=	1 = Setup disabled 0 = Setup enabled (wire in)	Setup Bridge
Th= Id=		Motor pulse output for speed regulation
Vd=		output for speed regulation
<b>Outputs</b>		
H1=	1 = off 0 = on	Heater tanks
Cir=	1 = off 0 = on	Circulation pump
Rep=	1 = off 0 = on	Replenishment pump
Col=	1 = off 0 = on	Cooling magnetic valve
Wsh=	1 = off 0 = on	Wash magnetic valve
Hd=	1 = off 0 = on	Dryer heater and fan
<b>LED</b>		
Testing LED.Look at little disp.		Display led check
<b>Mem</b>		
Load default val    Medical NDT Exit	Medical	Factory settings
<b>Motor</b>		
Speed:	0 + -	for Motor test Range 0 - 60
Ph_Del:	x x	for motor service
T=	xxx	measured Motor speed
<b>SB</b>		
Sensor: S-bar:	0 + - - - - - -	for sensorbar test

## INDEX

1. INTRODUCTION	1
2. GENERAL SAFETY INSTRUCTIONS	2-6
3. PREINSTALLATION	7
4. CHEMISTRY- and WATERCIRCUTS	8,9
5. REPLENISHMENT	10
6. CHEMISTRY DRAINS / WATER DRAINS	11,12
7. INSERT THE TRANSPORT RACKS	13-15
8. THE FIRST STEPS	16
9. WORKING WITH THE FILMPROCESSOR	17,18
10. THE DISPLAY	19 - 25
11. MAINTENANCE	26
12. RECOM. MAINTENANCE EVERY 2-4 MOUNTH	27
13. SERVICE OF THE PROCESSOR	28-44
14.SOFTWARE UPDATE	45-46

## 1. INTRODUCTION

Congratulations upon your decision to buy a

### **COLENTA X-RAY FILMPROCESSOR.**

Your purchase has been designed to meet the highest technical standards.

Some outstanding design features are:

- \* ) compact, space-saving design
- \* ) full automatic processing cycle
- \* ) smooth roller transport system
- \* ) low tank volumes
- \* ) electronically controlled temperature system
- \* ) automatic replenishment
- \* ) low energy consumption

This manual is an instruction for routine use of your:

### **COLENTA X-RAY FILMPROCESSOR.**


**Colenta**<sup>®</sup>  
**MEDIPHOT 900E**





## 2. GENERAL SAFETY INSTRUCTIONS

- \* ) Staff in charge of maintenance of the processor need to be thoroughly familiar with the equipment.
- \* ) Only the Top Cover of the Filmprocessor may be removed by the operator (see picture below)
- \* ) The film processor must be separated from mains prior to carrying out any maintenance. To do so, switch the mains switch of the machine to „0“ position. Always wear safety goggles and protective clothing when handling chemicals.



- \* ) The filmprocessor may not be in operation without supervision.
- \* ) Make sure that clothing or other objects cannot get entrapped in gear drives or similar of the film processor.
- \* ) The installation, service, repair as well as the initial operation of the machine may be carried out by qualified and trained service personnel only.
- \* ) Built-in safety devices may not be eluded or made inoperative. Use only original COLENTA spare parts when exchanging failed electrical components.
- \* ) Do not wear any loose necklaces or bracelets!

### 2.1 CHEMICAL HANDLING

- \* ) Observe all safety technical regulations of the chemical manufacturer.
- \* ) Wear safety goggles and protective clothing when handling chemicals.
- \* ) Ensure proper ventilation of the room in which the chemicals are prepared.
- \* ) Spillage or overflow of chemicals has to be removed instantly.
- \* ) In case of contact with the eyes flush with plenty of cold water for approximately fifteen minutes and consult a physician.
- \* ) Chemical disposal has to be in accordance with the local environmental codes. Contact your local water treatment and sewer district authorities for more information.

## 2.2 IMPORTANT WARNING AND SAFETY INSTRUCTIONS (Please read these instructions carefully.)

### Processor Operation

Make sure that no long hair, loose clothing or jewellery can get entrapped in the gear drives or in the media transport area.

The „**Service Manual for Colenta Film Processor**“ is for the use of qualified service personnel only.

**The racks must be cleaned with running water outside the Film Processor.**

Do not clean the processor with running water.

### Electrical and Mechanical Hazards

Observe all safety warnings to minimize the risk of electrical shock, burns or equipment damage. Photographic Film Processors are complex machines with many electrical and mechanical parts as well as with a considerable amount of chemicals.

### Fire Prevention

The area around the processor must be kept clean at all times. Keep dust, wood shavings, scrap paper or other inflammable materials out of the dryer compartment.

Fire extinguishers must be available in the room where the processor is installed.

### Chemical Handling and Accident Prevention

Misuse of almost any chemical may create a hazard of some type. Generally photo chemicals are not any more dangerous as most of the regular cleaning agents. However, there is always a residual risk. When handling chemicals observe the procedures below.

1). Never sniff a container or open bottle to determine its contents. A cautious sniff of the cap or lid is safer.

2). When handling chemicals wear protective clothing, safety goggles and rubber gloves .

3). Label storage containers properly. Avoid storing hazardous chemicals on high shelves or in unprotected glass containers. Keep chemicals away from children. Do not store chemicals in a refrigerator used for food because they may contaminate food or be mistaken for edibles.

4). Ensure proper ventilation in the area where chemicals are used or stored.

5). Observe the manufacturer's recommendations for using and mixing chemicals.

Overexposure to photographic chemistry may cause skin irritation to certain individuals.

### PHOTOGRAPHIC CHEMICAL EMERGENCIES AND FIRST AID PROCEDURES:

- **SKIN** - Rinse thoroughly with water.
- **EYES** - Rinse thoroughly with water and consult a physician.
- **INGESTION** - Consult a physician immediately.

## Service Manual for MEDIPHOT 900E

### Chemical Disposal

Waste from photographic processing normally contains diluted chemicals. These chemicals should be collected and disposed in accordance with local environmental codes. Dumping chemicals into a drain system could lead to a pollution problem. Contact your local water treatment and sewer district authorities before disposing chemicals.

All plumbing must comply with local and national codes. The DRAIN must be made of chemical resistant and non-corrosive material. Use PVC or equivalent

### Exhaust, Temperature and Humidity

It is necessary to ensure proper ventilation in order to receive good processing results.

Make sure that the exhaust hose of the built-in exhaust fan is properly connected to the exhaust air socket (Picture 1).

The built-in exhaust fan exhausts the fumes from the film processor.

These chemical fumes are corrosive. The top cover of the film processor should be removed over night.

Room temperatures between 18-26 °C (65-80 °F) with a relative humidity between 35% and 75% are ideal for photographic processing and working.

The film processor is a complex machine with moving parts such as the gear train and media transfer components. It uses photo processing chemicals which are irritating to eyes, lungs and skin. High voltage is used to power the film processor.

The dryer compartment produces heat.

- High voltage may cause electric shocks, burns or even death.

- Hands or fingers may be pinched or injured by moving parts or when handling heavy parts.

- Dryer compartment heat may ignite combustible materials and cause fires.

- Eyes, skin and lungs may be irritated by photo chemicals. Before using photo processing chemicals always read the Material Safety Data Sheets (MSDSs) for information about the hazards of the particular chemicals and how to use them safely.

- Do not operate the film processor after consuming alcohol or taking strong medication.

- Do not wear jewellery or loose clothing when operating the processor.

### Electronical and Electrical Hazards

HAZARDOUS VOLTAGE CAN CAUSE ELECTRIC SHOCK, BURNS OR EVEN DEATH.

Qualified service personnel must verify during installation that the processor is permanently and reliably grounded according to standards in the National Electrical Code.

## Service Manual for MEDIPHOT 900E

Carry out the following steps prior removal of the top cover:

**1. Train operators of the filmprocessor.**

**2. Switch off the main power switch (,,0“-position) and secure against restart by locking with a padlock (see picture below).**



**Main switch of the film processor**

### Fire Hazards

DRYER COMPARTMENT PRODUCES HEAT -PAPER OR OTHER COMBUSTIBLES CAN BE IGNITED

· Keep the area within 10 feet of the processor clean. Do not store combustible materials, including paper, close to the filmprocessor.

· Verify that a functional 10 lb. ABC fire extinguisher is located close to the processor.

### Burn Hazard

DRYER COMPARTMENT PRODUCES HEAT -DRYER PANELS AND GUARDS GET HOT

· Therefore do not touch dryer panels or guards when dryer in operation

### Corrosive Liquids

CHEMICALS MAY IRRITATE EYES, LUNGS, SKIN AND DIGESTIVE TRACT

# Wear safety goggles, protective gloves and chemical aprons as indicated on Material Safety Data Sheets (MSDSs) when handling chemistry.

# Drain tanks carefully, avoid splashing. Always drain the system thoroughly before working on any of the external hose systems.

# Read the MSDSs for more information regarding the proper safety procedures for working with photo processing chemicals.

# Do not allow untrained personnel to handle photo processing chemicals or to operate the filmprocessor.

# To avoid hazardous situations, keep floors and floor coverings around the processor and associated drains clean and dry at all times. Any accumulation of fluids outside the film processor, should be cleaned up immediately

## Service Manual for MEDIPHOT 900E

In the event of an accumulation of liquid due to backup, overflow or other malfunctions of the drain associated with the filmprocessor call a plumber or other contractor to correct the problem with the drain. Colenta assume any responsibility or liability whatsoever for the service ability of any drain connected to the filmprocessor. Such drains are the sole responsibility of the customer.

**DRAINS** must be made of chemically resistant and non-corrosive material.

### Chemical Hazards

Chemicals can be the source for errors, contaminate the waste water, irritate skin or eyes. Spills must be cleaned up immediately as follows:

1. Prevent the spilled chemicals from entering a waste water drain.
2. Clean up the spill with a moist mop or rag.

### CAUTION!

**When handling chemicals, especially fixer, wear protective clothing, safety goggles and rubber gloves.**

If filmprocessor chemicals make contact with the eyes, rinse them thoroughly with lots of water. If irritation persists, consult a physician.

3. Dispose cleaning materials and collected waste water from the clean up according to environmental regulations and local codes.
4. Avoid any inhalation of chemicals as this is dangerous to health.
5. Observe all environmental regulations for storage and disposal of waste chemicals.
6. Use this manual together with the instructions for chemicals.  
When handling chemicals wear protective clothing, safety goggles and rubber gloves,

### Corrosive Vapours

CHEMICAL VAPOURS MAY IRRITATE EYES, LUNGS AND SKIN IF ALLOWED TO ACCUMULATE IN WORK AREA

Assure an adequate supply of fresh outdoor air through natural or mechanical ventilation.

- Make sure that qualified service personnel is checking the external exhaust system at regular intervals.
- Read the Material Safety Data Sheets (MSDSs) for more information regarding the proper safety procedures for working with photo processing chemicals.

### 3. PRE- INSTALLATION

- \*) Site preparation, e.g., water supply, drainage electrical supply must be completed prior installation.

#### 3.1. LOCATION

- \*) Processor can be installed "through-the-wall" or completely in the darkroom. Required measurements can be taken from the processor specification sheet. For "through-the-wall", a purpose built panel is required (optional accessory).

#### 3.2 ELECTRICAL SUPPLY

- \*) All electrical connections must meet national safety requirements. Correct fuses and electrical requirement can be taken from the processor specification sheet.

**Main switch of  
the film processor**

**Power cord**



#### 3.3 WATER SUPPLY


- \*) The processor must be connected to the local water supply by using a DVGW-approved system separating device or pipe separating device.
- \*) The cold water supply pipe must have a stopcock fitted connection to the processor and should be done by using the 3/4" hose connector, supplied. Easy access to the stopcock should be provided as it has to be opened and closed daily.
- \*) A built in magnetic valve reduces water consumption to a maximum of 2,5 ltr./minute using pressure and quantity control.
- \*) It is recommended to run a second cold water supply with 2.5 meters of hosing to facilitate easy cleaning of the racks and tanks (water supply kit - optional accessory).

**Water connection und  
drains in front of the  
processor.**



## 4. Chemistry and water circuits

Following illustrations/drawings and pictures confirms the circulation and Washwater system of the processor:




**DEVELOPER  
TopView**

in detail....

Valve to adjust the flowrate of the intermediate rinse


Temp.-probe DEV



DEV Replenishment inlet

**FIXER  
TopView**


in detail....



FIX Replenishment inlet

**WATER  
TopView**

in detail....





Valve to adjust the flowrate  
of the intermediate rinse




## 5. REPLENISHMENT

When operating a processor which uses chemicals for the continuous production of plate/film it is very important that the chemicals within the machine are kept in good working order so as to provide consistent processing quality. To achieve this consistency we use replenishment solutions, which are formulated by the chemical manufacturer and injected into the processor precisely for the area of material being produced. Replenishment of the chemical tanks is done automatically using infra red sensors located at the entrance to the processor. These sensors accurately monitor the width of material entering the processor, this information is then used by the microprocessor (CPU) control software to calculate the surface area for each plate loaded into the processor. Each sensor, when covered, will generate a pulse, which is then recorded on a decoder and counted – the more sensors that are covered then the faster the count. When the count reaches the programmed value of pulse counts it triggers the start of a replenishment cycle. During each replenishment cycle the replenishment pumps inject fresh solution from small storage bottles/ tank and into the corresponding "working" tank solutions for a pre-set time.




External Developer replenisher Tank




USE FLOATING LID TO PROTECT DEVELOPER AGAINST OXYDATION

External Fixier replenisher Tank



### WARNING


- \*) Do not use brass or copper in the drainage system.
- \*) Chemistry disposal must be in accordance with local environmental regulations.
- \*) To avoid back pressure in the drain, the hoses should be free of bends and with a constant downward gradient.



## 6. CHEMISTRY DRAINS / WATER DRAINS

To drain the filmprocessor (Developer, Fixer and the Wash) just open the drain tubes according the illustrations below. Take care that all the mentioned drain taps are close during re-fill up.


**IMPORTANT:** Used Developer and used Fixer has to be collected in suitable containers separately.



WATER Supply Line

DRAINS for:  
WATER      DEV      FIX

Drain Tubes



Lift to drain the WASH  
(No-Code)

Lift to drain the FIX  
(Blue-Code)


Lift to drain the DEV  
(Red-Code)

Check regularly that the rubber O-Ring is in correct position.




**6.1 INFRARED REPLENISHMENT SENSORBAR**

The automatic replenishment system is using an infrared-sensor-bar to detect the incoming film area. With that information the CPU of the filmprocessor will calculate the replenishment rate which will be need.



**IMPORTANT:** Special care must be taken to ensure that the processor entrance rollers are always clean and dry – any spillage of chemical or water onto the feed tray / feed rollers or sensor bar must be avoided. Any spillage must be cleaned immediately.

**6.2 REPLENISHMENT CONTAINERS**



The filters, shown above, should be checked monthly and be cleaned or replaced if necessary.

**WARNING:** Separate the Film Processor from mains. To do so, switch the main power switch of the processor to „0“ position. Wear safety goggles, protection gloves and clothing.

## 7. INSERT THE TRANSPORT RACKS


After cleaning or changing the chemistry, the transport racks has to be re-installed to the processor. To do that, follow the instruction and illustration mentioned below:

**WARNING: Separate the Film Processor from mains. To do so, switch the main power switch of the Filmprocessor to „0“ position. Wear safety goggles, protection gloves and clothing.**


- # Wash the tanks with water, remove any deposits in the chemistry tanks.
- # Remove any possible algas from the washtank.
- # Close the the drain-tube (DEV-FIX-WASH) as mentioned on the page before.
- # Fill the chemistry tanks (first the FIXER, than the DEVELOPER!!) to the red maker. (Fill carefully and slowly - prevent splashes)
- # Fill the Washtank to the red marker.
- # Insert the racks according the reference number in follqwing the sequence below. \*
- # Fix the rack according the next page.
- # Close the Topcover and switch on the processor.

**\*Sequence to insert the racks:**

1. DevRack into the DevTank
  2. FixRack into the FixTank
  3. Wash&Dry Rack into the WashTank
  4. COS-A between DEV & FIX
  5. COS-B between FIX & WASH
- COS...CrossOverSegment





**INSERT THE TRANSPORT RACKS**





**IMPORTANT:**


check that the drive gears of the transport racks (DEV / FIX / WASH) are in correct position.



**INSERT THE TRANSPORT RACKS**



After all the racks and Cross-OverSegment are installed, fix them by closing the rotary joints.



## 8. THE FIRST STEPS

**WARNING:** Separate the Film Processor from mains. To do so, switch the main power switch of the Filmprocessor to „0“ position. Wear safety goggles, protection gloves and clothing.

**8.1 Never switch on the processor (main switch position "1") when the tanks (DEV/FIX/WASH) are empty. This would create serious technical problems!!**

### 8.2 USING THE CHEMISTRIES

- \*) Only use chemistry suitable for roller transport systems.
- \*) Follow instructions of chemistry manufacturers.

#### FIXER BATH:

- \*) Empty fixer tank by opening the fix drain tap.
- \*) Remove the Fixer-rack.
- \*) Check fixer tank is free of alien material.
- \*) Close fix drain tap.
- \*) Fill fixer tank with ready-to-use-fixer solution to the red marker on the tank wall. Insert the Fixer-rack very carefully and slowly, add hardener solution if advised by the chemistry manufacturer.

#### DEVELOPER BATH:

- \*) Empty developer tank by opening dev drain tap.
- \*) Remove the Developer-rack.
- \*) Check developer tank is free of alien material.
- \*) Close dev drain tap.
- \*) Fill developer tank with ready-to-use-developer solution to the red marker on the tank wall. Insert the developer-rack very carefully and slowly. Replenishment tanks may be used to mix the chemistry. Any remaining can be used for replenishment.

<p><b>CAUTION:</b> Even the smallest quantity of fixer could contaminate the developer solution. Therefore, always fill with fixer first. When removing the fixer rack, always cover the developer tank. For removing the fixer rack use rack carrier tray (optional accessory)</p>
---

## 9. WORKING WITH THE FILMPROCESSOR

### IN THE MORNING

- \*) Turn on water supply.
- \*) Check the levels of the replenishment containers (DEV&FIX)
- \*) Switch on the Filmprocessor with the Filmprocessor main switch (position „1“).
- \*) Wait for the "READY" of the processor

### STARTING WORK

- \*) Check level of the replenishment containers (DEV&FIX)
- \*) Check level of the waste containers (DEV&FIX)
- \*) Select programme
- \*) Feed through one or two of cleaning films (optional item).
- \*) During feeding films, always check the free-signal, given form the display.
- \*) Ensure first rollers pull material.
- \*) Feed large format films in straight.

### IN THE EVENING

- \*) Turn off water supply.
- \*) Switch off the main power switch of the Filmprocessor. (Main switch in position „0“)
- \*) Open water drain tap to prevent algae growths in water tank.

\*) **Lift the top cover to prevent condensation !!**






## 9.1 PROCESSOR FUNCTIONS

PROGRAMMING:	Automatic processing parameters, e.g., temperature, speed and replenishment rates, can be stored in 3 different programmes.
WARMING-UP:	Once programmed, temperature settings are accurately controlled. Heating commences with switching on at the mains. Constant solution temperatures are maintained in the processing tanks. Temperatures tolerances +/- 0,2 °C are achieved by the microprocessor control unit while the solutions are circulated by circulation pumps. When temperature has reached PRE-SET levels, the filmprocessor enters STANDBY mode and is ready for use.
STANDBY:	In case no material is processed - after a programmable periode of time, since the last media has exited the filmprocessor transport, dryer and water supply is switched off automatically. The filmprocessor goes in standby mode and is ready for work.
ANTICRYSTALLI-ZATION CYCLE	During STANDBY mode - within a programmable cycle periode - transport and intermediate rinse bath water supply is activated - this prevent cristallization build up on crossover rollers.
ANTI-OXIDATION CYCLE	During STANDBY mode - and no material is processed during set time - an preprogrammable ANTI OXIDATION cycle (replenishment cycle) is available. The additional replenishment compensates the impact of airoxidation of the chemistry during standby mode und tops up the chemistry levels in the tanks, compensating evaporation of the water in the solutions during standby.
AUTO REPLENISHMENT:	The filmprocessor comes equipped with a film area measuring facility. Infrared sensors scan the film area touchless and when the preprogrammed amount of film (area) entered the filmprocessor, a replenish-cycle is activated.
AUTOMATIC START-STOP:	Infrared sensors also automatically control the startcycle of the filmprocessor. The filmprocessor changes from STANDBY to RUN once a film has interrupted the light barrier. As the rollers turn, water is supplied to the wash tank and to the intermediate rinse bath system. Once the last film has passed through, the filmprocessor reverts to STANDBY. The film can be taken out of the receiving basket or top cover lid.

## 10.THE DISPLAY

Number of programs	5
Temperature range, developer and fixer	18.0 ÷ 43.0°C
Temperature control tolerances	±0.2°C
Temperature measurement resolution	0.03°C
Developing time tolerances at max. speed	±0,2%
· Motor speed is quartz-stabilized and controlled by a separate microprocessor	



- manual operation
- move cursor
- select menu item/  
change value

### 10.1 After switching on the processor

After the processor is installed and with all the components and service in place, switch on the processor. The controller display will illuminate and confirm the software revision momentarily.



u24

...in the display

Never switch on the processor (main switch position "1") when the tanks (DEV/FIX/WASH) are empty. This would create serious technical problems!!

After the initialization-cycle is complete (it takes 1-2 sec.) the processor will start in program1 (P1)



21.1

...in the display


Light2

...will be on


The "21.1" indicate that the current Developer temp. The flashing **Light2** informs about that an Error is existing. This is a normal status after switching on the processor, especially in the morning and depends on the room temp. where the unit is installed. The processor will start now to activate the heating element for the DEV, this of course to heat the DEV to value which is programmed in the microprocessor. As soon as this value is reached, the **Light2** will stop to illuminate - that confirms that no error is pending anymore. In that case the processor is ready to process films.

NOTE: the processor will deactivate the error code 1°C before the programmed value of the DEV is reached - but of course the heating element will still heat until the required temp. is reached - during that the **Light2** will flash, now it is possible to load film. 0,3 C° before the programmed DEV Temp., the processor will stop heating, the light2 is now OFF - the processor is now ready.

**10.2 Error codes**

As mentioned before, when **Light2** is active (on), an error is existing. To see which error, use the button  on the display you will see:

**Er1**

use  again to check wether there are more errors - until on the display will show

**Er-**

this confirms that no more errors are pending any more. Following error codes are possible:.

Error Code	Action
Er1 – Cover opened.	Close the top cover of the processor
Er2 – Developer too cold.	Normal message after start-up, when active during normal working
Er3 – Developer too warm.	Check wether the water-tap is open. when active during normal working -
Err – Motor overloaded.	Check the drive system
Er- – no more errors.	
E99 – Setup invalid;	reprogramm and check the Setup
E98 – Standby invalid;	reprogramm and check the Standby
E97 – Program invalid	reprogramm and check the Program (P1 to P3)
??? – Temp.-probe problem	Call service personal


Use  to jump back to main page (actual Dev-Temp.).

**10.3 PROGRAMPARAMETERS**

The software offers the possibility to use 3 different processing-programs. They are called P1, P2 until P5. In each of that mentioned programes, the following parameters can be preprogrammed

- # DeveloperTemp
- # Replenishment amount / m<sup>2</sup>
- # DeveloperTime

To see and to check the current parameters, push  1 time, you will get:



### 10.4 Changing the program


As mentioned before, the software offers 5 different programs which can be chosen by using the display: to do that, follow the instruction/illustrations below:

Push  1 time, you will get:



...in the display

**Hnd**




Use  1 time you will get:




...in the display

**P1**

....or P2 until P5

Use  to select the program (P1,P2 until P5) you want to use. When the program is shown, push  1 time you will get **USE** confirm  1 time. After the processor finished the re-initialization cycle, the processor is now working according the new programmparameters.

Step by step:



10.5 Reprograming the DEVELOPER temperature


In each of the five programs, the developer temp. can be reprogrammed, this according the specs given by the film/chemistry supplier.

To do that, use  1 time, you will get:



**Hnd**

...in the display

To do that, use  1 time, you will get:



**Px**

...in the display

x for P1,P2  
until P5




1 time



**USE**

...in the display

Use again  1 time, you will get:




**33.0**

actual Dev.  
Temp. for Px


x for P1,P2  
until P5

## Service Manual for MEDIPHOT 900E

Now, push again **+✓** 1 time, the first digit will start flashing:



Again with **+✓**, change the number of the first digit to the values you wish. Use **➡** to move to the second digit and change the number with **+✓** in the same way. Use again **➡** to move now to the last digit and again **+✓** to change the number you need.




Now, when the new value is shown, the last digit is still flashing, push **➡** until (approx. 3 sec) you can see: **PrG** This confirms that the new value is stored.

At least the processor will start one re-initialization cycle:

**u24**

To change the DEV-temp. for other programs (P1 - P2....P5) in the same way: change the program in the way as disreided on page 19 - 7.4.




## 10.6 Manual Mode

It might be necessary to run the transport by hand, this for example to process film regardless the current DEV-temp. To start the dive motor by hand use  1 time you will get:



...in the display

**Hnd**

Now push  1 time, you can hear the running. To stop the motor push again . To exit this mode, use .


## 10.7 Distance between films

To prevent film jams, some minimum distance between the films is needed. This needed distance is fixed in the processor setup.

# Load film: as soon as the sensorbar recognized that film is there, the Light2 will be on.

# During loading in film, the Light2 will be on to confirm this. At this status, never load any films additional!!

# When Light2 is off (there will be a beep as well) the processor is ready to load film again.





# 11. MAINTENANCE

The filmprocessor is designed to produce consistent high quality production with the minimum of maintenance.

Regular maintenance minimizes the chances for equipment failure and loss of processing quality. A well trained person has to be responsible for performing the maintenance of the filmprocessor and must be familiar with the operation and function of the processor as well.


## 1. Daily maintenance

- \* ) Check levels of the external replenishment tanks - If necessary mix fresh solution.
- \* ) Cleaning feed tray.
- \* ) Cleaning spray-bar-guide for the fixer
- \* ) Before starting production we advise to feed some cleaning films to remove any overnight residue.



## 2. Weekly maintenance

- \* ) Wipe external surfaces of film processor / enclosures / panels with a wet cloth to remove any chemical / dirt deposits.
- \* ) Inspect and clean the wash tank and intermediate water rins drains. If algae present then they should be removed, in such a case we suggest to use a proven algae control system
- \* ) Check the COS (CrossOverSegments) remove any deposits.
- \* ) Take out the guide of the COS and clean them as well.



## 12. RECOMMENDED MAINTENANCE EVERY 2-4 MONTHS.

**(Period is subject to filmprocessor usage.)**

Good processing quality and the reliable operation of a filmprocessor is dependent upon regular and careful cleaning. Every 3-6 months, the chemicals in the tanks should be drained. A chemical cleaning of the processing tanks and wash tank is recommended. Always follow safety warnings as described in section 1 when cleaning the filmprocessor.

Prior to **carrying out any maintenance**, switch off the power at the main power switch (position "0") ensuring it cannot be accidentally switched back on.

- \*) Switch off the main power switch of the film processor first (position „0“).
- \*) Remove the top cover of the filmprocessor.
- \*) Drain individual tanks by open the draitaps in front of the filmprocessor.
- \*) Remove rack assemblies (water / DEV / FIX, see item 2.1) and put them aside.
- \*) Close taps and fill all tanks with water or better with suitable cleaning solution until the red mark inside the tanks are reached.
- \*) Put the racks back into the tanks of the filmprocessor and close the top cover.
- \*) Switch on the filmprocessor and start some replenishment cycles. The hoses will be cleaned with water as well. Also start the transport of the filmprocessor, the racks has to be in. Let the filmprocessor run for 10 to 15 minutes.
- \*) Switch off („0“ position)the main power switch of the filmprocessor and drain the filmprocessor tanks again.
- \*) NOTE: Use cleaning solution according to the manufacturer´s instructions.
- \*) After tank cleaning, the developer- and wash-tank should be filled twice with fresh water (eventually use neutralizer recommended by manufacturer). Let the filmprocessor run for approximately 10 minutes again. Check all external (outside of the filmprocessor) hose connectors (outside of filmprocessor) and fittings for leaks.
- \*) Drain all tanks.
- \*) Remove the water / DEV / FIX Racks and check for:
  - worn gears
  - damaged or worn bearings
  - loose screws
  - scratched or bent film guides
  - plastic flat springs in developer bottom underturn.
  
- \*) All repairs must be carried out by qualified service personnel.
- \*) Check the inside of the tank for contamination and alien substances.
- \*) Clean the rollers well.
- \*) Close the drain taps of all 3 tanks.
- \*) Fill developer and fixer tanks with fresh chemicals to the required level (1st fixer, 2nd developer)
- \*) Fill wash tank.
- \*) Re-install the racks carefully. Take care of correct sequence of the racks is followed and make sure the gears are in the right position. Secure the racks.
- \*) Insert the respective suction pipe to the correct external replenisher tank.  
Re-install the top cover and switch on the filmprocessor.
- \*) Process test films.

## 13. SERVICE OF THE PROCESSOR


Warnung: Diese Anleitung ist nur für qualifizierte Service Techniker bestimmt.  
Warning: For the use of qualified service personnel only.  
Avertissement: Réservé au personnel de service qualifié.

The COLENTA MP900E X-Ray Filmprocessor is equipped with a small 7 segment OperatorDisplay - this display allows the operator, limited in a small range, to set some different processing parameters - described in the Instruction Manual. For service, a separate Display has to be installed - this only for reprogramming or service and trouble shooting.

Operator Display  
(Standard)




Display for Service  
(Optional)



Switch off the power at the main power switch (position "0") ensuring it cannot be accidentally switched back on.

100cm cable




Open the Cover by opening the 2 fixing screws.


Cover E-Box

13.1 INSTALLATION OF THE SERVICE DISPLAY


Install the Service display following the illustrations below:



Operator Display (Standard)



Display for Service (Optional)



Sensorbar

Operator Display (Standard)

Temp. Probe DEV

Display for Service (Optional)


When the service display is installed as mentioned before, switch on the processor.





NOTE: At this moment, the cover of the E-box is open so there is a risk of electrical shocks, ensure that nobody is getting contact with components and/or wirings inside the E-Box.

13.2 WORKING WITH THE NEW DISPLAY





13.2.1 AUTOMATIC MODE

The processor is designed to work without operator assistance. Under normal circumstances the operator will use the front panel only to check the process parameters and progress.

To scroll through the pages, press . - Press  to jump back to main page.



-  back light ON/OFF
-  check errors /alarm shutdown
-  setup mode
-  back to top menu

-  manual operation
-  move cursor
-  select menu item/  
change value
-  scroll page down

## Service Manual for MEDIPHOT 900E


1) **55% Done** is the progress indicator. It means that **55%** of the developing process are done. When **100%** are reached, the processor will go to standby.

**T1=32.4°C** confirms the actually measured temperature in the DevTank. The value shown in the brackets is the programmed temperature.

2)

**Rep1 400ml** confirms the amount of replenishment for the Developer .

Note: To supply fresh chemistry to the Fixer as well, a double head replenishment pump is installed. Read more on this issue on page....



**Speed:** Linear speed of the Film "going" through the processor in cm/min.

3)

**Tank1 time** is the time the media needs to pass through **tank 1** (Developer)

**Dry to dry** is the length of the complete processing cycle (leading edge to heading edge)

4)


**Add 14 pls** confirms the additional added pulses for the input section of the DEV


**Rdy 10 pls** confirms the safety distance between loading films



## Service Manual for MEDIPHOT 900E

### 13.2.2 PROGRAMMING PROCEDURES


Switch on the processor it will start in work mode. Make sure that no media is being processed, re-programming is enabled only during standby.

Press . The programming menu will appear:



Use , to move the cursor to Program and confirm with .

You will get:






## Service Manual for MEDIPHOT 900E

Press  $\uparrow$  to change the number of the program you wish to modify (P1 to P5).

Use  $\rightarrow$ , to move the cursor to Modify and use  $\uparrow$  to confirm.

The programme (P1,P2 or P3) consists of three pages:



### To change the programmed values:

# use  $\rightarrow$  to move the cursor to the digit you want to modify.

# use  $\uparrow$  to add "+1" to the value (..8-9-0-1-2..)

# use  $\downarrow$  to jump to the next page (..8-9-0-1-2..)

# save / exit


NOTE: It's not possible to go back to the previous page. In case you programmed something wrong, start again form page one.

Use  $\downarrow$  to scroll through the pages. Use  $\uparrow$  to confirm your changes or move the to Cancel and confirm with  $\uparrow$  as well to abort.




## Service Manual for MEDIPHOT 900E

During saving, you might will get the message "Out of range" which confirms that one or more inputs are not allowed - even ot of range.



After 2 seconds the message will disappear and you will be taken back to re-programme the values. A parameter that was too high will be automatically reset to the maximum possible value. A parameter that was too low will be reset to the minimum possible value.

This can be used if you want to program extreme values - for instance you want to use the shortest developing time possible, but you don't remember the value. In this case just programme 000. After the "Out of range" message, the developing time will be reset to the minimum. Just select Save once again.

To go back to work mode, press  .


During programing by using the service display, on the operator display will be shown:


**Pro**






### 13.2.3 CHANING THE PROGRAM


To use another programme, use  you will get:



With , move the cursor under P1.

Press  to change the programme number.

With , move the cursor to Use and confirm with .

Press  to jump back to main page.

**Note:** a cyclic redundancy check is used to verify the data being read from the non-volatile memory. If some damage occurred to the program data, or the programme was never set up properly, you will get an error message *Program invalid*. The solution is to go to programming mode and re-program the data. This error will occur also if the EEPROM chip has been replaced in which case it contains random data.


## Service Manual for MEDIPHOT 900E

### 13.2.4 MANUAL START / STOP

The manual start/stop is possible only when no media is being processed.

During the processing the corresponding menu items are not selectable - you can't move the cursor there.

To run the motor manually:

Press 


With , move the cursor under Start and select it with .


This will run the motor. The menu item Start changes to **Stop**.

You can stop the motor by selecting **Stop**.

When you start the motor manually, this will be indicated on the main page as **M1 instead of P1**.

### 13.2.5 DISPLAY ILLUMINATION ON/OFF


In a dark room, it might be necessary to switch off the display backlight to prevent exposure. The  button toggles the backlight on/off.

When the backlight is off, all the buttons except  are disabled.

This is done to prevent pressing buttons by accident in a dark room. Switching the display off is a good idea if the processor is left unattended. This will reduce the chances for unauthorized people to operate the machine.

### 13.2.6 AUTOMATIC START


The processor will start automatically when media is fed, except in case the developer is too cold - more than 1°C below the programmed. In this case, feeding the media will not start the processor. Instead you'll get the message,





which will disappear after 2 seconds.


If you need to feed a film regardless of the low developer temperature, run the machine in manual mode.

13.2.7 ERRORS


If an error occurs, the indication P1 (or M1) will alternate with Er. If this happens Press the  button. This will stop the beeper and bring you to the error menu, so you can check what the processor is telling you.

In case more than one error occurred, press  to scan the rest of them. Press  to jump back to main page.

When the processor switched on at the beginning of the working hours ,it is expected to have low temperatures in the tanks.For this reason, the **Er** indication will be present, but without alarm. If, however, the temperature drops during normal work, the alarm will be activated.



Error Code	Action
<b>Er1</b> – Motor overloaded.	Only appears when the re-pulses given from the motor are not the same as the CPU of the processor are expected them. The reason could be a blocked drive system. In that case for the first step, take out all CrossOverSegments and try again. 2nd take out all the racks/dryer and try again. When after that 2 tests the message is gone, check the racks, eventually replace them. In case the message is still active, check the drive gears below the tankbody. If they are OK, replace the motor and/or main board.
<b>Er2</b> – Cover opened.	Close the top cover of the processor. In case the message is still present, check the magnets (2!) connected to the main board. Eventually replace them.
<b>Er3</b> – Developer too cold.	Normal message after start-up, when active during normal working, check the circulation pump, the heating element for the DEV (over temp. fuse!! installed on the heatexchanger housing) Check wether all outputs are driven - from the main borad. (Service Mode)
<b>Er4</b> – Developer too warm.	Check wether the water-tap is open. when active during normal working - Check wether the output on the main board is driven (Service Mode). In case yes, replace the 2 way Water valve.
<b>Er-</b> – no more errors.	no more errors



Drive system

Drive Motor



Circulation pump



2 way water valve



Temp. fuse



Heatexchanger

**E99** – Setup invalid;

Check all the SETUP values, reprogram them and save again. Check the WEB for software update.

**E98** – Standby invalid;

Check all the STANDBY values, reprogram them and save again. Check the WEB for software update.

**E97** – Program invalid

Check all the STANDBY values, reprogram them and save again. Check the WEB for software update.

**Err** – Problem with Temp. probe DEV


Replace the temperature probe of the DEV

13.2.8 STANDBY OPTIONS

The Processor supports an anti-oxidation and a anti-crystallization cycle.

That means, when in standby, the processor will start the transport of the rollers and the wash on regular intervals in order to prevent any crystallization or build ups on the rollers (Anti-crystallization). The anti-oxidation cycle activat, in free programmable time intervals, replenishment cycles. This will prevent oxidation of the the chemistry.

During such a cleaning cycle, the display will look like this.



During such a cycle the processor will accept media. It's not necessary to wait to the end of the cycle.

Press , use  to move the courser to "Standby" use  to confirm; use  to scroll through the pages:

1) **SB replenishment  
100ml each 6h**

The processor will activat a 100ml replenishment cycle each 6 hours.(Anti-oxidation-cycle)



2) **SB self-cleaning  
20cm each 10min**

The processor will activat 3 roller turns (1 roller turn @ 8 cm) each 3 min.(Anti-crystallization-cycle)



3) **Switch ON Dryer  
for 5 min.**

After switching on the processor (Main switch) the dryer will start heating for 5 min. countinually this heat up the dryer and the dryer components.



4) **DR On each 10min  
for 01 min**

When the the processor is on, but not in use for longer time, the dryersystem would loose his heat. To prevent this, each 10 min the dryer (heater&fan) will start work.

**Save Cancel**

13.2.9 MONITOR MODE

The feature can be used to check current status of the processor. That mens it confirms wether different outputs are driven or not, the sensorbar is occupied or not..... see the table below for more detail:

Motor OFF  
Spd=34  
v=00  
P=000  
t=000



S-Bar: - - - - -  
Area=0,0000 sqm



P89c668 MM\_V24  
ID=000-0000  
r=1



H1=1/0  
Cir=1/0  
Rep=1/0  
Col=1/0  
Wsh=1/0  
Hd=1/0

Main drive Motor is working or not  
Speed of the Main drive motor  
needed for internal calculation  
needed for internal calculation  
needed for internal calculation

Confirms, which sensor is currently active  
Confirms the ingoing filmarea

Confirms the software code  
Option for Identification No.  
needed for internal calculation

Output Heater driven or not...1=on, O=off  
Output Circulation pump driven or not...1=on, O=off  
Output Replenishment pump driven or not...1=on, O=off  
Output Cooling Valve driven or not...1=on, O=off  
Output Wash Valve driven or not...1=on, O=off  
Output Fan for Dry driven or not...1=on, O=off

13.2.10 SETUP OF THE PROCESSOR

Press , use  to move the cursor to "SETUP" use  to confirm; use  to scroll through the pages:

1) **Tank1 0042 pls**  
**Machine 0248 pls**

Start Point of the processor for counting pulses.  
Total pulses trough the processor.



2) **Gear 219 pls**  
**Pump 02.6 ml/s**

Gear is the number of motor pulses corresponding to 1m advance of the material. It is needed to calculate the processed film area (for the replenishment) and the processing speed in cm/min.

Pump is the number of milliliters per second of the replenishment pump. Needed to allow programming the replenishment in milliliters.



3) **Sensor distance 060 mm**

Sensor distance is the distance between two sensors of the sensorbar - needed for the film area calculation.



4) **Replenish after each 0.125 sqm**

Replenish after. This is that area, after one replenishment cycle has to be activated.

5) **Ready 0030 pls after film left**

Defines the safety distance between loading film (free signal for the next fim).


6) **Add 14 pulses before start**

Due the fact that the optical sensorbar is installed before the film will enter to the DEV, a calculation factor has to be added.

7) **Additional Run 001cm**

For safety reasons, the moto will run 1cm additional after the exited the dryer rollers

**Save Cancel**





**NOTE:** The shaft of the main drive motor has a tachometer installed . When the drive motor starts the tacho will rotate to generate a speed counter which sends pulses to the Processor CPU. By using these pulses the CPU is able to accurately control the transport speed/timing sequences



## Service Manual for MEDIPHOT 900E




It is possible to see data and to change working data in the SETUP menu but to save any new data entered into the SETUP menu then it will be necessary to insert a jumper on the main board of the processor:



Don't forget to remove the jumper when the reprogrammong is done.

When selecting Save in the setup menu, the program will refuse to save the data unless this wire jumper is present.

13.2.10 SERVICE / Electronic check (Inputs-Outputs-Components)

Press , use  to move the cursor to "SERVICE" use  to confirm you will get:

Inputs	Outputs
Mem	LED Motor SB

**Inputs** 

Cov=0	1=On and 0=Off	Input for Coverswitch
Lev=0	1=On and 0=Off	Input for Level Switch (Option,Not in use)
SP=0	1=On and 0=Off	Input for Setup jumper
Th=0	Re-pulses given from the motor disk	
ID=0	Input for speed regulation	
VD=0	Output for speed regulation	

---

**Inputs** 

H1=0	1=On and 0=Off	Output for HeaterH1
Cir1=0	1=On and 0=Off	Output for Circulation Pump
Rep1=0	1=On and 0=Off	Output for Replenishment Pump
Col=0	1=On and 0=Off	Output for Magnetic Valve: Cooling
Wsh=0	1=On and 0=Off	Output for Magnetic Valve: Wash
Hd=0	1=On and 0=Off	Output for DryerFan

**Mem** 

Load default val	
Medical NDT Exit	..... <u>Medical</u> when the processor is used for <u>medical</u> X-Ray application
	..... <u>NDT</u> when the processor is used for <u>industrial</u> X-Ray application

**LED** 

Use to check the LEDs of the operator Display

**Motor** 

Speed:	0 +-	Used to check the motor needed for internal calculation needed for internal calculation
Ph_Del:	0	
T=0		

**SB** 

Speed:	0 +-	no. of sensor
S-bar:	- - - - -	Sensor check (optical sensorbar)

## 14. SOFTWARE UPDATE



Warnung: Diese Anleitung ist nur für qualifizierte Service Techniker bestimmt.  
Warning: For the use of qualified service personnel only.  
Avertissement: Réservé au personnel de service qualifié.

Whenever a Software update available, we are announcing that on our Homepage. We recommend to check this side regularly to see whether a new software is available. In case yes, you can download the update free of charge.

- # [http://www.colenta.de /X-ray/X-Ray Processors/MP900E](http://www.colenta.de/X-ray/X-Ray Processors/MP900E)
- # click to Software update and enter the Password: **MINI56xf**
- # Download the software and store the file (hex) on your hard disc
- # Check whether some additional informations are available


Follow now the illustrations below to update the processor.

Firstly you need to install the program "Flash Magic" on your system. To do that, open the folder MINI\_FlashMagic supplied on that CD or download from WEB as well. Follow up the installation routine. When the installation was done successfully, start the program, open the pull-down menu "Options", confirm "Advanced Options" open the register card "Timeouts" and set the values as shown below:




## Service Manual for MEDIPHOT 900E

- # Switch off the Processor
- # Place the jumper on ST15 as shown.



JUMPER ST15


- # Connect the download cable
  - on your PC / RS232 - 9pin
  - on the MainBoard of the processor on RS232



JUMPER ST15


- # Switch on the Processor (the service display will confirm that the CPU is now in program mode by showing black lines)
- # Start the program "FlashMagic"

Don't forget to remove the jumper when the reprogramming is done.



- add1) Compare the values shown on the image before with the values mentioned on your FlashMagic program.
- add2) Make sure that the box "Erase all Flash+Security" is chosen.
- add3) Browse the new Download file
- add4) make sure that nothing is chosen
- add5) Push the button Start

**NOTE:** In case an Error is indicated right after starting, restart the processor, restart the program and try again.

**IMPORTANT:** After successfully updating the software, remove the update jumper (mentioned on the page before), put the setup jumper in position and enter the setup by pushing  (service display!). Choose Service, move over to MEM you will get :

Load default val  
Medical NDT Exit

According the application: choose Medical or NDT

**NOTES:**

---


**NOTES:**

---

**NOTES:**

---






			Datum	28-11-2005	RDM8-2_01	1	Colenta LABORTECHNIK A-2700 Wr. Neustadt	STROMLAUFPLAN RDM8-2 WIRING DIAGRAM RDM8-2	Blatt 1 von 3 Bl.
			Bearb.	HANDLER					
			Gepr.	fh					
Zust.	Aenderung	Datum	Name	Norm	DIN EN 60950	Z.Nr.	Zust.	Ers. fuer	


L (Page 1/12) →


N (Page 1/12) →

PE (Page 1/12) →




DETAIL DRAWING  
DRYER HEATER





		Datum	28-11-2005			Colenta	STROMLAUFPLAN RDM8-2			
		Bearb.	HANDLER	RDM8-2_03	0				LABORTECHNIK	
		Gepr.	fh			A-2700 Wr. Neustadt		WIRING DIAGRAMM RDM8-2		Blatt 3
Zust.	Aenderung	Datum	Name	Norm	DIN EN 60950	Z.Nr.	Zust.	Ers. fuer	von 3 Bl.	



			Datum	28-11-2005	PDB_MINI_04_2004 _P89C668	0	Colenta LABORTECHNIK A-2700 Wr. Neustadt	SICHERUNGSBESTÜCKUNG RDM 8-2		Blatt 1
		Bearb.	HANDLER	Zust.						
Zust.	Aenderung	Datum	Name	Norm				DIN EN 60950	Z.Nr.	