

Deep Temperature Freezer

DW-40W100•DW-40W255•DW-40W380



Effective models

This service manual is effective for following models

| Model name | Product code | Voltage(V) | Frequency(Hz) | Plug-type |
|------------|--------------|------------|---------------|-----------|
| DW-40W100 | | 220 | 50 | All |
| DW-40W255 | | 220 | 50 | All |
| DW-40W380 | | 220 | 50 | All |
| | | 220 | 50 | All |


Content

| | |
|---|----|
| Product Features..... | 4 |
| Structure Drawing of Product Appearance | 5 |
| Technical Data..... | 6 |
| Circuit Diagram..... | 7 |
| Task Control Principle & Parameters..... | 8 |
| Temperature Controller Parameters | 11 |
| Compressor Parameters..... | 13 |
| Trouble shootings..... | 14 |
| FAQ | 17 |
| Question & Answer | 18 |
| Structure Schematic Diagram | 19 |
| Wiring Diagram..... | 20 |

Product Features

1. Temperature range inside the freezer is $-20\sim-40^{\circ}\text{C}$ (adjustable) (with the environmental temperature of $10\sim25^{\circ}\text{C}$);
2. Microprocessor control with digital display of the temperature;
3. High-density insulation layer & low power consumption;
4. High/low temperature alarm controller, with certain temperature alarm point to be set as needed;
5. With famous-brand Compressor, flouid-free & environmentally friendly refrigerating medium and two alarm modes (buzzer alarm & blinking indicator alarm);
6. Security door lock design to prevent arbitrary opening of the door;
7. Stepped anti-corrosion liner design for location of food baskets in different types
8. Wide voltage band, suitable to the $187\sim242\text{V}$ voltage.

Structure Drawing of Product Appearance




Technical Data

| Model | DW-40W100 | DW-40W255 | DW-40W380 |
|---|-------------|--------------|--------------|
| Climate Type | N | N | N |
| Protection type to prevent electric shock | I | I | I |
| Usable Capacity (L) | 100 | 255 | 380 |
| Rated Power Supply | AC220V/50Hz | AC220V/50Hz | AC220V/50Hz |
| Rated Power (W) | 260 | 410 | 550 |
| Power Consumption (kW·h/24h) | | | |
| Refrigerating medium R600a (g) | | | |
| Boxed Dimension (WxDxH) (mm) | 680*600*805 | 633*1243*838 | 633*1550*838 |
| Net Weight/Gross Weight (kg) | 43/46 | 69/81 | 82/89 |
| Compressor | SC12CL | SC18CL | SC21CL |


Circuit Diagram

1) Applied model: DW-40W100



M1- Compressor M2-Cooling blower P - Electronic temperature controller
 KA-Starter relay C₁-Starter capacitor XP-Power plug

2) Applied model: DW-40W255/380




M1- Compressor M2-Cooling blower P - Electronic temperature controller
 KA-Starter relay C₁ Starter capacitor XP-Power plug
 C₂-Running capacitor

Task Control Principle & Parameters

1. Control Principle & Functions:

This freezer adopts an electronically temperature control system controlled by the Carel temperature controller.

Temperature controller panel & wiring diagram



Temperature Controller Panel

Temperature Controller Wiring Diagram

Debug Methods of the temperature controller:

| Serial No. | Button Operation | Display |
|------------|---|---|
| 1 | Press the button "Set" for about 2 seconds | To display the initially set temperature of -30 |
| 2 | If "-30" is not displayed, press "Up▲" / "Down▼" for resetting. | To display: -30 |
| 3 | Press the button "Set" | To display the current internal temperature |
| 4 | After check of the set temperature, it is still necessary to check the setting of other parameters. | |
| 5 | | To display the current internal temperature |
| 6 | Press the button "Set" for at least 5 seconds | To display: a flashing "PS" |
| 7 | Press the button "Set" | To display: 0 |
| 8 | ▲ By pressing the button "Up▲" | To display: 22 |
| 9 | Press the button "Set" | To display: PS |
| 10 | By pressing the button "Up▲" | To display: /C1 |
| 11 | Press the button "Set" | To display: 0 |
| 12 | If "0" is not displayed, press "Up▲" / "Down▼" for resetting. | To display: 0 |
| 13 | Press the button "Set" | To display: /C1 |
| 14 | By pressing the button "Up▲" | To display: /4 |
| 15 | Press the button "Set" | To display: 1 |
| 16 | If "1" it is not displayed, press "Up▲" / "Down▼" for resetting. | To display: 1 |

| | | |
|----|------------------------|----------------|
| 17 | Press the button "Set" | To display: /4 |
|----|------------------------|----------------|

| | | |
|----|--|-----------------|
| 18 | By pressing the button “Up▲” | To display: /5 |
| 19 | Press the button “Set” | To display: 0 |
| 20 | If “0” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 0 |
| 21 | Press the button “Set” | To display: /5 |
| 22 | By pressing the button “Up▲” | To display: /6 |
| 23 | Press the button “Set” | To display: 1 |
| 24 | If “1” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 1 |
| 25 | Press the button “Set” | To display: /6 |
| 26 | By pressing the button “Up▲” | To display: rd |
| 27 | Press the button “Set” | To display: 2 |
| 28 | If “2” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 2 |
| 29 | Press the button “Set” | To display: rd |
| 30 | By pressing the button “Up▲” | To display: r1 |
| 31 | Set Press the button “Set” | To display: -40 |
| 32 | If “-30” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: -40 |
| 33 | Press the button “Set” | To display: r1 |
| 34 | By pressing the button “Up▲” | To display: r2 |
| 35 | Press the button “Set” | To display: -20 |
| 36 | If “-10” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: -20 |
| 37 | Press the button “Set” | To display: r2 |
| 38 | By pressing the button “Up▲” | To display: c4 |
| 39 | Press the button “Set” | To display: 100 |
| 40 | If “100” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 100 |
| 41 | Press the button “Set” | To display: c4 |
| 42 | By pressing the button “Up▲” | To display: cc |
| 43 | 按 Set Press the button “Set” | To display: 0 |
| 44 | If “0” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 0 |
| 45 | Press the button “Set” | To display: cc |
| 46 | By pressing the button “Up▲” | To display: D1 |
| 47 | Press the button “Set” | To display: 8 |
| 48 | By pressing the button “Down▼” | To display: 0 |
| 49 | Press the button “Set” | To display: D1 |
| 50 | By pressing the button “Up▲” | To display: AL |
| 51 | Press the button “Set” | To display: 5 |
| 52 | If “5” is not displayed, press “Up▲” / “Down▼”for resetting. | To display: 5 |
| 53 | Press the button “Set” | To display: AL |
| 54 | By pressing the button “Up▲” | To display: AH |

| | | |
|----|---|---|
| 55 | SET Press the button "Set" | To display: 5 |
| 56 | If "5" is not displayed, press "Up▲" / "Down▼"for resetting. | To display: 5 |
| 57 | Press the button "Set" | To display: AH |
| 58 | By pressing the button "Up▲" | To display: Ad |
| 59 | Press the button "Set" | To display: 0 |
| 60 | If "0" is not displayed, press "Up▲" / "Down▼"for resetting. | To display: 0 |
| 61 | Press the button "Set" | To display: Ad |
| 62 | By pressing the button "Up▲" | To display: H1 |
| 63 | Press the button "Set" | To display: 0 |
| 64 | If "0" is not displayed, press "Up▲" / "Down▼"for resetting. | To display: 0 |
| 65 | Press the button "Set" | To display: H1 |
| 66 | By pressing the button "Up▲" | To display: H2 |
| 67 | Press the button "Set" | To display: 1 |
| 68 | If "1" is not displayed, press "Up▲" / "Down▼"for resetting. | To display: 1 |
| 69 | Press the button "Set" | To display: H2 |
| 70 | Press the button "Set" for over 5 seconds to save the set value | To display the current internal temperature value |

Note: Error signals of this temperature controller are "E0" and "E1". When there comes a short circuit, breakage, fault or incorrect connection for the temperature sensor, the display screen shall show "E0". When there are wrong settings for the parameter /4 (default 1) or the parameter D0 (default 2), the display screen shall show "E1" with the buzzer alarming and the alarm indicator blinking. If it is the case, it is necessary to check the wiring and parameters' setup.

Temperature Controller Parameters

| Symbol | Symbol Name | Unit | Set Value | |
|-----------|---|-------|-------------------|-----------------------|
| | | | DW-40W100/255/380 | DW-25W198/300/388/518 |
| | Set Temperature | °C/°F | -30 | -20 |
| PS | Password | — | 22 | 22 |
| / | Sensor Parameters | | | |
| /C | Temperature Calibration for the Environmental Sensor | °C/°F | 0 | 0 |
| /2 | Stability Measurement | ---- | 4 | 4 |
| /4 | Sensor Probe Setup | — | 0 | 0 |
| /5 | Choose °C/°F (0=°C 1=°F) | — | 0 | 0 |
| /6 | Decimal Point Forbidden | — | 1 | 1 |
| r | Control Parameters | | | |
| r1 | Permitted Minimum Temperature Set Value | °C/°F | -40 | -30 |
| r2 | Permitted Maximum Temperature Set Value | °C/°F | -20 | -10 |
| rd | Control Thermal Difference | °C/°F | 2 | 2 |
| A | Alarm Parameters | | | |
| A0 | Alarm & Blower Control Difference | °C/°F | 0 | 0 |
| Ad | Lag Time of Temperature Alarm | Min. | 0 | 0 |
| AL | Offset of Low Temperature Alarm Limit Value to the Set Value | °C/°F | 5 | 5 |
| AH | Offset of High Temperature Alarm Limit Value to the Set Value | °C/°F | 5 | 5 |
| c | Compressor Parameters | | | |

| | | | | |
|-----------|--------------------------------|------|---|---|
| c0 | Compressor Lag Time at Startup | Min. | 0 | 0 |
|-----------|--------------------------------|------|---|---|

| | | | | |
|-----------|--|------|-----|-----|
| | While the Controller is Working | | | |
| C1 | The Minimum Interval between the Adjacent Two Operations of the Compressor | Min. | 0 | 0 |
| c2 | Shortest Down Time of the Compressor | Min. | 0 | 0 |
| c3 | Shortest Up Time of the Compressor | Min. | 0 | 0 |
| c4 | "Duty setting" Function Valid or Not (0=off 100=constantly on) | Min. | 100 | 100 |
| cc | Duration of Continuous Operation | Hour | 0 | 0 |
| c6 | Alarm Lag Time after Continuous Operation | Hour | 2 | 2 |
| H | Other Parameters | | | |
| H1 | Special Allocation (for Definition of Defrosting or Not) 0=no 1=yes | — | 0 | 0 |
| H2 | Button Lock Function (Locked or Not) 0=yes 1=no | — | 1 | 1 |
| H4 | Buzzer Valid or Not 0=yes 1=no | — | 0 | 0 |

Compressor Parameters

| | | | |
|-------------------------|-----------|----------------------|--------------|
| Product Model | DW-40W100 | | |
| Compressor Model | SC12CL | Rated Frequency (Hz) | 50 |
| Refrigerating Ouput (W) | 253 | Motor Type | CSIR |
| Input Power (W) | 260 | Refrigerating Mode | Wind Cooling |
| Rated Voltage (V) | 198~254 | Oil Filling (cc) | 550 |

| | | | |
|-------------------------|-----------|----------------------|--------------|
| Product Model | DW-40W255 | | |
| Compressor Model | SC18CL | Rated Frequency (Hz) | 50 |
| Refrigerating Ouput (W) | 313 | Motor Type | CSR |
| Input Power (W) | 410 | Refrigerating Mode | Wind Cooling |
| Rated Voltage (V) | 198~254 | Oil Filling (cc) | 550 |

| | | | |
|-------------------------|-----------|----------------------|--------------|
| Product Model | DW-40W380 | | |
| Compressor Model | SC21CL | Rated Frequency (Hz) | 50 |
| Refrigerating Ouput (W) | 378 | Motor Type | CSR |
| Input Power (W) | 550 | Refrigerating Mode | Wind Cooling |

| | | | |
|-------------------|---------|------------------|-----|
| Rated Voltage (V) | 198~254 | Oil Filling (cc) | 600 |
|-------------------|---------|------------------|-----|

Trouble shootings

| Product Model | Serial No. | Parts Details | Specification or Model | Private Part No. |
|---------------|------------|------------------------|---|------------------|
| DW-40W100 | 1 | Compressor | Zanuxi HVT86AA | 0074090643 |
| | 2 | Temperature Controller | Carel Electronic Temperature Controller | 0074091216 |
| | 3 | Power Cord | Power Cord | 0075180003 |
| | 4 | Bonding Wires | Bonding Wires for Electric Connection | 0074090759A |
| | 5 | Evaporator | Evaporator $\Phi 7.8 \times 0.6$ | 0070700258 |
| | 6 | Condenser | Internal-Mounted Condenser $\Phi 4.56 \times 0.6$ | 0070700256 |
| | 7 | Compartment Shield | Compartment Shield | 0070103022 |
| | 8 | Filter | Drier Filter of Single-Chamber Vacuum Pumping | 0060702792 |

| Product Model | Serial No. | Parts Details | Specification or Model | Private Part No. |
|---------------|------------|------------------------|---|------------------|
| DW-40W255 | 1 | Compressor | Jiaxipera ZBU1122CY | 0060702534 |
| | 2 | Temperature Controller | Carrel Electronic Temperature Controller | 0074091216 |
| | 3 | Blower | BD-219H (Nigerian) Blower 220V/50HZ | 0074091180 |
| | 4 | Condenser | BD-319H Domestic/European Internal-Mounted Condenser with 7 Loops | 0070701804 |
| | 5 | Return Pipe | BD-319H Return Pipe Group | 0070701809 |
| | 6 | Evaporator | HBC-268 Evaporator | 0270700031 |
| | 7 | Connecting Wires | Temperature Controller wire 1 | 0075180057 |
| | 8 | Power Cord | Power Cord | 0075180003 |
| | 9 | Connecting Wires | Compressor & Junction Box Connection | 0075180002 |
| | 10 | Connecting Wires | Temperature Controller wires 3RV1*0.75MM2 | 0075160026 |
| | 11 | Connecting Wires | Temperature Controller wires 3RV1*0.75MM2 | 0075160025 |
| | 12 | Connecting Wires | ROHS-BD-60GPB Compressor Earth Lead | 0070401132 |
| | 13 | Condenser | BD-379H (Nigerian) Fiber Tube Condenser | 0070803664B |
| | 14 | Filter | (20G Molecular Sieve) Common Double-Chamber Vacuum Pumping Filter | 0074180002A |
| | 15 | Junction Box | Junction Box | 0072040045 |

| Product Model | Serial No. | Parts Details | Specification or Model | Private Part No. |
|---------------|------------|------------------------|---|------------------|
| DW-40W380 | 1 | Compressor | Jiaxipera ZBU1122CY | 0060702534 |
| | 2 | Temperature Controller | Carrel Electronic Temperature Controller | 0074091216 |
| | 3 | Condensate Fan | 220V/50HZ | 0074091180 |
| | 4 | | | |
| | 6 | Power Cord | SC-180 Power Cord Assembly | 0070401234A |
| | 7 | Bonding Wires | SC-280 Temperature Controller Wires (A) | 0070401561C |
| | 8 | Lighting | SC-137SCP Lighting Tube (11W Osram U-Type Lighting Tube) | 0071400593 |
| | 9 | Lamp Shade | BC-110 Lamp Shade with transparent PS | 0070200919 |
| | 10 | Evaporator | Clad Plate Evaporator | 0074010035 |
| | 11 | Condenser | SC-280 Compartment Condenser | 0070701642 |
| | 12 | Compartment Shield | SC-280 Compartment Shield | 0070203048A |
| | 13 | filter | Common Double-Chamber Vacuum Pumping Filter (20G Molecular Sieve) | 0074180002A |
| | 14 | Shelf | SC-280 Food Shelf | 0070101179A |
| | 15 | Shelf (bottom) | SC-280 Food Shelf (bottom) | 0070101179B |
| | 18 | Door | SC-280 Glass Door Assembly | 0070810664 |
| 19 | Lock | BC-50G Lock Assembly | 0070804947 | |


FAQ

| Main Problem | Cause Analysis | Maintenance Measures |
|--|---|--|
| The compressor cannot be started. | 1. Fuse burn-out | Replace the fuse. |
| | 2. Connection Damage of the socket connector(s) for the compartment wires | Replace the socket connector(s). |
| | 3. Incorrect connection of wires | Check and make correct installation. |
| | 4. Fault in the controller | Replace the temperature controller. |
| | 5. Damaged starter or thermal protector | Replace the starter or the thermal protector. |
| | 6. Fault in the compressor | Replace the compressor. |
| The temperature controller displays "E0". | Fault in the sensor. Replace the temperature controller or the temperature sensor. | |
| Big noise from the refrigerator | 1. Uneven position. | Change the refrigerator's position. |
| | 2. Resonance caused among the tubes or among the refrigerator cabinets during the compressor's operation. | Sort out the tubes to avoid resonance. |
| | 3. Loosening of the compressor's connecting bolts | Tighten the connecting bolts |
| Refrigerating difference of the refrigerator | 1. Serious leakage of refrigerating medium | Check the leakage and refill the refrigerating medium. |
| | 2. Dirty block and oil block of capillary tubes or of the system | Clean the capillary tubes or replace the filter. |
| No alarm | 1. Wrong setup of the parameters | Check parameter setup of the temperature controller. |
| | 2. Fault in the temperature controller | Replace the temperature controller. |


Question & Answer

| Main Problem | Answer to the problem |
|--|---|
| The internal temperature is too high or too low. | It is a bit relevant to the surrounding environmental temperature. When the internal temperature is too high or too low, we can solve this problem by adjusting the controlling temperature. |
| There is a sweeper phenomenon within the refrigerator after a certain duration of the its operation. | A little sweeper can be drained off periodically and be wiped off by a towel. If there is a lot more sweeper, we can check whether the drainpipe or other drainage are blocked and may get it though with a piece of fine wire. |

Structure Schematic Diagram



Wiring Diagram



Haier

Inspired living

Haier Medical & Laboratory Products Co., Ltd.

Room 403D, Brand Building, Haier Industrial Park, No.1 Haier Road

Qingdao China

Website: www.haiermedical.com