320 BB

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Identification – Identifikation

Service manual 320, Precision balances, s/n 5000000 and higher
350-8100 e1
June 2012
Subject to technical modifications

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www.precisa.com

Copyright

No reproduction in any form of this document, in whole or in part, may be made without written authorization from Precisa Gravimetrics AG.

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1	Components – Bestandteile	6
1.1	About spare parts – Über Ersatzteile	6
1.2	Balance groups – Waagengruppen	7
1.3	Housing's bottom – Gehäuseunterteil	8
1.4	Housing's top – Gehäuseoberteil	9
1.5	Weighing cell – Wägezelle	11
1.5.1	Overall view – Übersicht	11
1.5.2	Flexure repair sets – Lager-Reparatursets	13
1.5.3	Particular spare flexures – Jeweilige Ersatz-Lager	13
1.6	Weighing pan D80 Group A2 – Waagschale D80 Gruppe A2	14
1.7	Weighing pan 135, Group M – Waagschale 135, Gruppe M	14
1.8	Weighing pan 170, Group C – Waagschale 170, Gruppe C	15
1.9	Weighing pan 200, Group D – Waagschale 200, Gruppe D	16
1.10	Draftshield, XT – Windschutz, XT	17
1.11	Draftshield 180, XB – Windschutz 180, XB	19
1.12	Draftshield 260, XB – Windschutz 260, XB	21
1.13	Draftshield easy – Windschutz einfach	23
1.14	Draft deflector M – Windabweiser M	24
1.15	Twin Reference Weight Device – Zwilling's-Referenzgewichts-Modul	25
1.16	Main board – Hauptprint	26
1.17	Connector board – Steckerprint	26
1.18	Display board VFD – Anzeigeprint VFD	27
1.19	Universal switching power adapter – Universal-Schaltnetzteil	28
1.20	Universal switching power adapter IP65 – Universal-Schaltnetzteil IP65	28
2	Repairs	29
2.1	Before you get started	29
2.1.1	Proceeding	29
2.1.2	About these instructions	29
2.1.3	Required tools	30
2.1.4	Preparing the tools	31
2.1.5	Spare parts	31
2.2	Opening the balance	32
2.2.1	Opening the balance	32
2.2.2	Opening the IP65 balance	32
2.3	Removing the weighing cell from the housing's bottom	33
2.4	Dismantling the weighing cell	34
2.4.1	Removing the coupling	34
2.4.2	Removing the 2 flexure holders	34
2.4.3	Removing the coil	34
2.4.4	Removing the support piece	34
2.4.5	Removing the balance arm	35
2.5	Exchanging the horizontal flexures	36
2.5.1	Removing the old horizontal flexures	36
2.5.2	Installing new horizontal flexures	36
2.6	Cleaning	37
2.6.1	Cleaning the magnet pot	37
2.6.2	Cleaning the coil	37
2.6.3	Checking the cleanliness again	37
2.7	Assembling the weighing cell	38
2.7.1	Assembling the balance arm	38
2.7.2	Assembling the coil, the magnet cover and the transport safety device	38
2.7.3	Assembling the floating sensor and the stopper	39
2.7.4	Assembling new vertical flexures	40

2.7.5	Assembling the support piece and the inner cone	41
2.7.6	Assembling the flexure holders	42
2.7.7	Assembling a new coupling	43
2.7.8	Dismantling the assembly jigs	43
2.8	Installing the weighing cell into the housing's bottom	44
3	Adjustments.....	45
3.1	Before you get started	45
3.1.1	Table of specifications	45
3.1.2	Placing loads	45
3.2	Setting the balance to the factory mode (the balance is open)	45
3.3	Adjusting the symmetry	46
3.4	Checking the pre-load	46
3.5	Adjusting the corner load	47
3.5.1	Determining the corner load values	47
3.5.2	Adjusting the corner load (Group M/C/D)	48
3.5.3	Adjusting the corner load (Group A)	49
3.6	Checking the lowering process of the internal reference weights (SCS only)	52
3.6.1	Checking the lowering process of the internal reference weights (10-Key version only)	52
3.6.2	Checking the lowering process of the internal reference weights (5-Key version only)	52
3.7	Closing the balance	53
3.8	Setting the balance to the service mode (the balance is closed)	54
3.8.1	Setting the balance to the service mode (10-Key version only)	54
3.8.2	Setting the balance to the service mode (5-Key version only)	54
3.9	Adjusting the linearity	55
3.9.1	Adjusting the linearity (10-Key version only)	55
3.9.2	Adjusting the linearity (5-Key version only)	55
3.9.3	Checking the linearity	56
3.10	Checking the adjustments of the closed balance	56
3.10.1	Checking the repeatability (with SCS)	56
3.10.2	Checking the repeatability (without SCS)	56
3.10.3	Checking the corner load	57
3.10.4	Checking the linearity	57
3.11	Adjusting the S-correction	58
3.11.1	Measuring the S-curve	58
3.11.2	Calculating the S-correction factors	58
3.11.3	Entering the S-correction factors (10-Key version only)	58
3.11.4	Entering the S-correction factors (5-Key version only)	59
3.12	Adjusting the internal reference weight (SCS only)	59
3.13	Error messages	60
3.13.1	Operating errors; Error disappears by correct use or application	60
3.13.2	Fatal errors; Balance stays still, error remediable in service mode	60
3.13.3	Hardware errors; Balance stays still	60
3.14	Discussion with the manufacturer	61
4	Reparatur	62
4.1	Bevor Sie beginnen	62
4.1.1	Vorgehen	62
4.1.2	Über diese Anleitung	62
4.1.3	Benötigtes Werkzeug	63
4.1.4	Vorbereiten des Werkzeugs	64
4.1.5	Ersatzteile	64
4.2	Die Waage öffnen	65
4.2.1	Die Waage öffnen	65
4.2.2	Die IP65-Waage öffnen	65
4.3	Die Wägezelle aus dem Gehäuseunterteil herausnehmen	66
4.4	Die Wägezelle zerlegen	67

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- Certain spare parts (e.g. flexures) have to be determined more precisely because of the various capacities of balances. Therefore refer to the following table and find out to which group a balance belongs.

Example: A 620M belongs to group M1. Thus, in group specific spare parts select the ones of group M1.

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



1.5.2 Flexure repair sets – Lager-Reparatursets

A repair set consists of the required spare flexures for one balance, i.e. 1 coupling, 2 vertical flexures and 8 horizontal flexures.

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












2.1 Before you get started

2.1.1 Proceeding


- **Important:** In case of a broken flexure, all flexures must be replaced.
- It's safest to put all parts in their pristine place when re-assembling, particularly the screws!
- Enough space and a sorted way of working are best conditions.

2.1.2 About these instructions

- Numbers in round brackets refer to components shown in the overall view of chapter "1.5.1 Overall view – Übersicht" on page 11. Have a printout of this overview always at disposal!
- Numbers in square brackets refer to tools shown in chapter "2.1.3.1. Tool kit 350-8537" on page 30. Have a printout of this chapter always at disposal, too!

2.1.3 Required tools

2.1.3.1 Tool kit 350-8537

Pos.	Pcs.	Article number	Article
1	1	W 41-1688	Mounting plate
2	2	W 41-1688-5	Support post long
3	2	W 41-1688-6	Support post short
4	1	350-6203	Screwdriver Torx M4
5	1	240-7124	Screwdriver small No.1
6	2	PW 29.43.300	Pin for flexure holders
7	1	W 41-1700	Assembly jig with slot
8	1	W 41-1701	Assembly jig with hole
9	4	W 41-1677	Assembly pin for assembly jig
10	2	W 41-1691	Centring spigot for flexures
11	1	W 41-1774	Centring spigot for coil, 360/321
12	1	W 41-1675	Centring spigot for coil, 320/330
13	1	W 41-1674	Coupling anchor
14	2	PN 1100-172	Screw cyl. Torx M4x12
15	1	W 41-1699	Alignment jig
16	1	W 41-1578-8	Connecting cable
17	1	W 41-1681	Centring bush

2.1.3.2 Instruments

- Soldering iron
- Voltmeter

2.1.3.3 Standard tools:

- Box spanner Nos.4, 5.5 and 7
- Flat-bladed screwdriver Nos.1, 2 and 4
- Phillips screwdriver No.2
- Hexagon socket wrenches SW 1.5mm, SW 2.5mm
- Tweezers
- Lens

2.1.4 Preparing the tools

- Screw the 2 long [2] and 2 short [3] support posts on the assembly plate as shown in the figure of chapter "2.1.3.1 Tool kit 350-8537" on page 30.
- Preheat the soldering iron to 310 +/- 5°C (590 +/- 10°F).

2.1.5 Spare parts


- Have all required spare parts at disposal, the fitting flexures (3), (6) and (8), according to the instructions in chapter "1.5.2 Flexure repair sets – Lager-Reparatursets" on page 13.

2.2 Opening the balance

2.2.1 Opening the balance


1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 screws on the balance's top using a No. 4 screwdriver (2 screws cyl. Phillips M4x25 above the display and 1 screw cyl. slot M4x12 nearby the bubble level).
4. Slightly raise the housing's top.
5. Disconnect the membrane keypad ribbon cable from the main board.

2.2.2 Opening the IP65 balance


1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the windshield, the protective ring, the weighing pan and the pan holder.
3. Unscrew the 3 screws IP56 with a No. 1 screwdriver and remove the clamping ring.
4. Unscrew the 3 screws on the balance's top (2 screws cyl. Phillips M4x25 above the display and 1 screw cyl. slot M4x12 nearby the bubble level) and the 2 screws cyl. countersunk M4x12 on the balance's bottom using a No. 4 screwdriver.
5. Slightly raise the housing's top. Take care of the sealing!
6. Disconnect the membrane keypad ribbon cable from the main board.

2.3 Removing the weighing cell from the housing's bottom

1. **SCS only:** Remove the 2 screws of the left counterholder and remove it with the reference weight.
2. Do so with the right part of the reference weight module, too.

3. Unplug the floating sensor (10) from the main board using a small screwdriver.
4. Turn over the housing's bottom and unscrew the 3 nuts (30) on the underside of the balance's bottom.
5. Hold the weighing cell only on the chassis and lift it out of the housing's bottom.
6. The chassis can now be placed on the assembly plate for easier handling.
7. To protect the weighing cell, fix both assembly jigs [7] and [8] with the 4 assembly pins [9].

2.4 Dismantling the weighing cell

2.4.1 Removing the coupling

1. Unscrew the 2 screws (28), remove the coupling (3), its round bases (34) and spring washers (32) (**Group A:** only the XR-bases (35)). Do not separate the screw components!
2. Screw the coupling anchor [13] tightly on with its 2 screws [14] in place of the coupling.

2.4.2 Removing the 2 flexure holders

1. Fix the upper flexure holder (5) into the chassis, using the 2 pins for flexure holders [6].
2. Unscrew the 4 screws (26) of the flexure holder and remove them together with its 4 spring washers (32), 2 on the chassis, 2 on the support piece (**Group A:** with its 4 XR-bases (35)).
3. Remove the flexure holder together with the flexures which are still fixed on it, then remove the pins from the holder.
4. Do so with the lower flexure holder, too.


2.4.3 Removing the coil

1. Unsolder the 2 connecting wires (13) carefully from the connection board (11) using the tweezers.
2. Unscrew the 3 screws (25) from the magnet cover (16). Remove the cover together with the transport safety device (17) and the spring washer (32).
3. Unscrew the 2 brass screws (23) from the coil (4) on the balance arm (2) using a No.2 screwdriver. Remove the coil carefully from the magnet pot, together with the washers (31).

2.4.4 Removing the support piece

1. Unscrew the 3 grub screws (29) of the inner cone (9) and remove it.
2. **SCS only:** Remove the left (19) and the right (20) reference weight bearing.
3. Unscrew the upper screw (balance arm side) of the coupling anchor [13] and pull off the support piece together with the coupling anchor over the assembly pins [9].
4. Untighten the lower screw of the coupling anchor for roughly 2 turns. The support piece and the coupling anchor are loosely held together now.

2.4.5 Removing the balance arm

1. Unscrew the brass screw from the stopper, turn it inwards ca. 90° and gently re-tighten the brass screw.
2. Unscrew the 4 screws (28) from the vertical flexures (6) and remove them together with its round bases (36) and spring washers (32).

Note: Since the lug of the floating sensor is firmly attached to the balance arm, the floating sensor must be unscrewed to dismantle the balance arm and gently tilted to the rear:


3. Unscrew the screw (21) from the connection board (11).
4. Unscrew the 3 screws (24) of the floating sensor and gently tilt the floating sensor to the rear.
5. Remove the two inner assembly pins [9], i.e. unscrew them from the balance arm.
6. Remove the two outer assembly pins [9], i.e. unscrew them from the chassis.
7. Remove the two assembly jigs [7] and [8].
8. Carefully tilt the balance arm over the pot and the floating sensor and take it out.

2.5 Exchanging the horizontal flexures

2.5.1 Removing the old horizontal flexures

1. Release all screws from the flexure holders and remove the old flexures together with their spring washers.

2.5.2 Installing new horizontal flexures


1. Place a flexure holder on the assembly plate [1]. The mark has to be upside-down, i.e. assembly plate sided!
2. Place 4 new horizontal flexures (8) carefully on the assembly plate, observe the alignment of their cut edges!
3. Place the 4 spring washers (32) (**Group A: 4 XR-bases (35)**) on the flexures and fix them tightly to the holder with their screws (24).
4. Take the complete flexure holder carefully off the assembly plate. Do not touch the flexures!
5. Prepare the other flexure holder the same way.

2.6 Cleaning


2.6.1 Cleaning the magnet pot

1. Wrap a double-sided adhesive tape around a flat stick.

2. Dab with the stick on the ground of the pot.

3. Prepare a thin and round stick the same way and spin it through the gap between the pot's wall and the magnet.
4. Dab over the top of the magnet pot with a piece of adhesive tape.

2.6.2 Cleaning the coil


1. Using a lint-free cloth, remove any contamination from the coil. **Caution:** Do not clean the coil with adhesive tape; risk of injury!

2.6.3 Checking the cleanliness again

1. Inspect the annular gap and the coil for residual contamination and clean if necessary.

2.7 Assembling the weighing cell

2.7.1 Assembling the balance arm

1. Place the assembly jigs [7] + [8] with 2 assembly pins [9] on the chassis, position them horizontally using the alignment jig [15] and tighten the two assembly pins [9] with a No. 4 box spanner.
2. Carefully introduce the balance arm (2) over the magnet pot and pass it to the floating sensor.
3. Assemble the balance arm with 2 assembly pins [9] loosely behind the assembly jigs. The balance arm must remain moveable.
4. Firmly tighten first the one, then the other inner assembly pin [9] of the assembly jigs [7] and [8].

2.7.2 Assembling the coil, the magnet cover and the transport safety device


1. Place the coil onto the balance arm, centre it with the centring spigot [11] and fix it firmly with 2 brass screws (23) and 2 washers (31). The centring spigot must be smoothly removable!
2. Slightly pre-fix the magnet cover (16) with its left and right screw (25).
3. Gently pre-assemble the transport safety device (17) with the third screw (25) and spring washer (32).


4. The centring bush [17] shall be introduced between the balance arm-pilot and the transport safety device to assure a constant annular gap. Tighten the screw (25), remove the centring bush and tighten the other two pre-fixed magnet cover screws.

2.7.3 Assembling the floating sensor and the stopper

1. Press the floating sensor carefully into the cutout of the chassis in the direction of the pot and fix it to the chassis with 3 screws (24).
2. Fix the connection board (11) onto the chassis with the screw (21).

3. Solder the two connection wires (1) parallelly and strain-freely to the connection board.


4. Turn back the stopper by 90° and tighten it with its brass screw.


5. The stopper has to be free of play and lie centrally in the symmetry screw (18). If necessary, adjust the symmetry screw.

2.7.4 Assembling new vertical flexures


1. To facilitate the assembly, the chassis can be placed on the assembly plate. Take care of the floating sensor!
2. Bring a new flexure (6) into position, its slot has to be on the balance arm side (bottom).
3. Centre the flexure with a centring spigot [10] first on the chassis side (top), then on the balance arm side (bottom).
4. Tighten the first (upper) centring spigot, then remove the second (lower) centring spigot. Do not touch the flexure and fix it there with its round base (34) and spring washer (32) (**Group A:** with its XR-base (35) only) and its screw (27).
5. Remove the remained (upper) centring spigot and fix the flexure the same way on the chassis side.
6. Do so with the other flexure, too.

2.7.5 Assembling the support piece and the inner cone


1. Introduce the support piece with the loosely held coupling anchor over the two inner assembly pins [9]. It must be possible to push the support piece smoothly over the two inner assembly pins. If not, release and tighten them again.
2. Screw the support piece and the balance arm together using the coupling anchor [13] and its two screws [14].
3. Assemble the inner cone (12) with the 3 grub screws (29) to the support piece.
4. **SCS only:** Re-assemble the two reference weight bearings (19) and (20) with the 2 screws (27) and spring washers (32/**Group A:** 35).

2.7.6 Assembling the flexure holders


1. Introduce the two pins for flexure holders [6] in the first pre-assembled holder the way that the shoulder of the pins is opposite the screw heads on the holder.
2. Bring the two pins with the holder into position in the chassis.
3. Lower the holder carefully over the two pins on the chassis, until the flexures just touch the chassis and the support piece.
4. Fix the holder with its 4 spring washers (32) (**Group A:** with its 4 XR-bases (35)) and screws (26).
5. Push the two pins through, up to the installed horizontal flexures holder.

6. Introduce and fix on the opposite side the other pre-assembled flexure holder the same way.
7. Gently pull the two pins out of the flexure holders. They must be smoothly removeable!

2.7.7 Assembling a new coupling


1. Remove the coupling anchor [13].
2. Bring the new coupling into position, its slot has to be on the support piece side.
3. Centre the coupling first on the balance arm side, then on the support piece side, using a centring spigot [10].
4. Tighten the first (upper) centring spigot, then carefully remove the second (lower) centring spigot. Do not touch the coupling and fasten it on the support piece side with its round base (34) and spring washer (32) (**Group A:** with its XR-base (35) only) and screw (28).
5. Remove the second (upper) centring spigot and fasten the coupling on the balance arm side the same way.

2.7.8 Dismantling the assembly jigs

1. Remove the 2 inner assembly pins [9] carefully.
2. Remove the 2 outer assembly pins [9] carefully.
3. Take out the two assembly jigs [7] and [8] by sliding them sideways.
4. Remove the two assembly pins for flexure holders [6] carefully.

2.8 Installing the weighing cell into the housing's bottom

1. Carry out a further movement check with the lens on the balance arm. The balance arm must not touch the transport safety device (17)!
2. Place the weighing cell carefully into the housing's bottom. Hold the cell only on its chassis!
3. Tilt the housing's bottom to the side and fix the weighing cell from below with the 3 nuts (30).
4. Connect the floating sensor (10) to the main board.

5. **SCS only:** Insert the internal reference weights and mount the counterholders. The counterholders have to rest fully on the reference weights, impeding them to move.

3 Adjustments

3.1 Before you get started

3.1.1 Table of specifications

In carrying out the following adjustments, it is necessary to know the specifications of the balance in service. You may download the specifications for servicing from our website:


<http://www.precisa.com>

The following instructions refer often to that specifications, have a printout of them always at disposal!

3.1.2 Placing loads

Place all loads always carefully and, unless otherwise asked, in the center of the weighing pan!


3.2 Setting the balance to the factory mode (the balance is open)


1. Remove the jumper J1 from the main board.
2. Level the balance with the levelling screws and connect it to the mains with the power adapter. The balance starts up in the factory mode.

3.5 Adjusting the corner load

3.5.1 Determining the corner load values


1. Ensure the balance is levelled and set to the factory mode, see chapter "3.2 Setting the balance to the factory mode (the balance is open)" on page 45.
2. Place the upper part of the housing on the right alongside the lower part of the housing and connect the membrane keypad to the pin header S4.
3. Put the pan holder and the weighing pan on the balance.
4. Have, according to the specifications, a suitable corner load reference weight at disposal.

5. Place the weight on Z and tare the balance.
6. Place the weight on A, B, C and D. Every single value has to stay, according to the specifications, within the tolerated range!

3.5.2 Adjusting the corner load (Group M/C/D)


Note: Skip this chapter if the determined corner load values stay within the tolerated range.

- The difference of the diagonal A-C is adjusted with the left corner load screw **1**.
- The difference of the diagonal B-D is adjusted with the right corner load screw **2**.
- Turning the corner load screw clockwise  slides the difference of the diagonal in negative direction. The displayed value relatively decreases.
- Turning the corner load screw counterclockwise  slides the difference of the diagonal in positive direction. The displayed value relatively increases.

Procedure

1. Ensure the balance is levelled and set to the factory mode, see chapter "3.2 Setting the balance to the factory mode (the balance is open)" on page 45.

2. Place the weight on C, tare the balance, shift the weight to A and calculate the difference: $\Delta AC = A - C = A - 0 = A$ [d]; (C = tared = 0); **Example:** $\Delta AC = A = +12d$.
3. Place the weight on D, tare the balance, shift the weight to B and calculate the difference: $\Delta BD = B - D = B - 0 = B$ [d]; (D = tared = 0); **Example:** $\Delta BD = B = -7d$.

Important: Start adjusting the diagonal with the larger difference!


4. The difference of the diagonal A-C is larger and positive (+12d). Turning the corner load screw 1 clockwise  slides the difference in negative direction (+12 .. +11 .. +10 ..).
5. The difference of the diagonal B-D is smaller and negative (-7d). Turning the corner load screw 2 counterclockwise slides the difference in positive direction (-7 .. -6 .. -5 ..).
6. Keep determining the diagonal differences as described in step 2 and 3 and reduce them, until all 4 corner load deviations stay within the tolerated range, see previous chapter "3.5.1 Determining the corner load values" from step 4 on.

3.5.3 Adjusting the corner load (Group A)

Note: Skip this chapter if the determined corner load values stay within the tolerated range.

- The difference of the diagonal A-C is adjusted with the left outer and inner corner load screw **1**.
- The difference of the diagonal B-D is adjusted with the right outer and inner corner load screw **2**.

Rough adjustment: the outer corner load screws have the same turning sense as in group M/C/D


- Turning the outer corner load screw clockwise  slides the difference of the diagonal roughly in negative direction. The displayed value relatively decreases.
- Turning the outer corner load screw counterclockwise  slides the difference of the diagonal roughly in positive direction. The displayed value relatively increases.

Fine adjustment: Caution, the inner and outer corner load screws have contrary turning senses!

- Turning the inner corner load screw counterclockwise  slides the difference of the diagonal finely in negative direction. The displayed value relatively decreases.
- Turning the inner corner load screw clockwise  slides the difference of the diagonal finely in positive direction. The displayed value relatively increases.

Procedure

1. Ensure the balance is levelled and set to the factory mode, see chapter "3.2 Setting the balance to the factory mode (the balance is open)" on page 45.

Rough adjustment having diagonal differences $\geq \pm 35d$:

2. Place the weight on C, tare the balance, shift the weight to A and calculate the difference:
 $\Delta AC = A - C = A - 0 = A [d]$; (C = tared = 0); **Example:** $\Delta AC = A = +68d$.
3. Place the weight on D, tare the balance, shift the weight to B and calculate the difference:
 $\Delta BD = B - D = B - 0 = B [d]$; (D = tared = 0); **Example:** $\Delta BD = B = -57d$.
4. Release the holding grub screws of the outer corner load screws by means of a hexagon socket wrench SW 1.5mm:

Important: Start adjusting the diagonal with the larger difference!

5. The difference of the diagonal A-C is larger and positive ($+68d$). Turning the outer corner load screw 1 clockwise  slides the difference in negative direction ($+68 \dots +67 \dots +66 \dots$). **Caution:** Co-turn the inner corner load screw, ensuring its thread keeps staying until its upper end into the outer corner load screw!

6. The difference of the diagonal B-D is smaller and negative ($-57d$). Turning the outer corner load screw 2 counterclockwise $\curvearrowright$ slides the difference in positive direction ($-57 \dots -56 \dots -55 \dots$).
Important: Co-turn the inner corner load screw, ensuring its thread keeps staying until its upper end into the outer corner load screw!
7. Keep determining and reducing the diagonal differences, until they fall below $\pm 35d$.
8. Fix the holding grub screws of the outer corner load screws by means of a hexagon socket wrench SW 1.5mm again.

Fine adjustment having diagonal differences $< \pm 35d$: Caution, the inner and outer corner load screws have contrary turning senses!


9. Place the weight on C, tare the balance, shift the weight to A and calculate the difference:
 $\Delta AC = A - C = A - 0 = A [d]$; (C = tared = 0); **Example:** $\Delta AC = A = +22d$.
10. Place the weight on D, tare the balance, shift the weight to B and calculate the difference:
 $\Delta BD = B - D = B - 0 = B [d]$; (D = tared = 0); **Example:** $\Delta BD = B = -19d$.

Important: Start adjusting the diagonal with the larger difference!

11. The difference of the diagonal A-C is larger and positive ($+22d$). Turning the corner load screw 1 counterclockwise $\curvearrowright$ slides the difference in negative direction ($+22 \dots +21 \dots +20 \dots$).
12. The difference of the diagonal B-D is smaller and negative ($-19d$). Turning the corner load screw 2 clockwise $\curvearrowleft$ slides the difference in positive direction ($-19 \dots -18 \dots -17 \dots$).
13. Keep determining the diagonal differences as described in step 9 and 10 and reduce them, until all 4 corner load deviations stay within the tolerated range, see chapter "3.5.1 Determining the corner load values" from step 4 on.

3.6 Checking the lowering process of the internal reference weights (SCS only)

1. Ensure the balance is levelled and set to the factory mode, see chapter "3.2 Setting the balance to the factory mode (the balance is open)" on page 45.

Pin header S4
for keypad

2. Place the upper part of the housing on the right alongside the lower part of the housing and connect the membrane keypad to the pin header S4.

3.6.1 Checking the lowering process of the internal reference weights (10-Key version only)


1. Hold  until "TESTPROG 1" appears.
2. Lower the weights with **{DOWN}** and lift them with **{UP}**.
3. If the weights don't run cleanly into the weighing cell's bearings, re-align the reference weight device.

3.6.2 Checking the lowering process of the internal reference weights (5-Key version only)


1. Hold **MODE** until "TESTPROG 1" appears.
2. Lower the weights with **{DOWN}** (Hold **REF** until "FKT-KEY 4" appears) and lift them with **{UP}** (hold **REF** until "FKT-KEY 3" appears).
3. If the weights don't run cleanly into the weighing cell's bearings, re-align the reference weight device.

3.7 Closing the balance

1. Disconnect the balance from the mains (unplug the power adapter).
2. Remove the weighing pan and the pan holder.


3. Ensure the Jumper J1 is set and the membrane keypad is connected to the pin header S4.
4. Ensure that both, the display and the glass of the housing's top are clean.
5. Place the housing's top on its bottom and fasten together with the 3 screws (2 screws cyl. Phillips M4x25 above the display and 1 screw cyl. slot M4x12 nearby the bubble level).
6. **IP65 only:** Fix the housing with the 2 screws cyl. countersunk M4x14. Fix the membrane carefully with the clamping ring and its 3 screws IP65 without causing any wrinkles on the membrane.

7. Put the pan holder, the weighing pan and, if present, the draftshield on the balance.

3.8 Setting the balance to the service mode (the balance is closed)

3.8.1 Setting the balance to the service mode (10-Key version only)

1. Level the balance with the levelling screws and connect it to the mains with the power adapter.
2. Press the service switch on the rear of the housing's bottom with the small screwdriver [5] ⇒ the display flashes once.
3. Press  ⇒ "CODE 0" appears.
4. Press  ⇒ "CODE 0000" appears.
5. Enter the service access code 5236 with ,  and  and confirm with . The balance is now in the service mode.

3.8.2 Setting the balance to the service mode (5-Key version only)

1. Level the balance with the levelling screws and connect it to the mains with the power adapter.
2. Press the service switch on the rear of the housing's bottom with the small screwdriver [5] ⇒ the display flashes once.
3. Press **MODE** twice ⇒ "CODE 0000" appears.
4. Pressing **MODE** changes the value of the figure, holding changes to the next figure.
5. Enter the service access code 5236 and confirm with **PRINT**. The balance is now in the service mode.

3.9 Adjusting the linearity

Caution: Adjusting the linearity erases all S-correction factors, too. This chapter has therefore to be processed entirely!

3.9.1 Adjusting the linearity (10-Key version only)

1. Ensure the balance has reached its operating temperature, is levelled and set to the service mode, see chapter "3.8 Setting the balance to the service mode (the balance is closed)" on page 54.
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Hold  until "LINEARITY" appears.
4. Have the weighing pan empty and press **{0}** ⇒ the balance catches the zero point and the display flashes.
5. Wait until the flashing stops, place the half load and press **{1/2}** ⇒ the balance catches the half load value and the display flashes.
6. Wait until the flashing stops, place the full load and press **{1/1}** ⇒ the balance catches the full load value and the display flashes.
7. Wait until the flashing stops and press **{SET}** ⇒ the balance sets the linearisation factors and displays the full weight value.
8. Empty the weighing pan, tare and check whether at quarter load and three quarter load the linearity stays within the tolerated range, too. If it doesn't, process chapter "3.11 Adjusting the S-correction" on page 58.
9. **With SCS:** Process chapter "3.12 Adjusting the internal reference weight (SCS only)" on page 59.
Without SCS: Switch the balance off and on (by pressing **ON/OFF** twice) ⇒ the balance quits the service mode.

3.9.2 Adjusting the linearity (5-Key version only)

1. Ensure the balance has reached its operating temperature, is levelled and set to the service mode, see chapter "3.8 Setting the balance to the service mode (the balance is closed)" on page 54.
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Hold **MODE** until "LINEARITY" appears.
4. Have the weighing pan empty and press **{0}** (hold **REF** until "FKT-KEY 1" appears) ⇒ the balance catches the zero point and the display flashes.
5. Wait until the flashing stops, place the half load and press **{1/2}** (hold **REF** until "FKT-KEY 2" appears) ⇒ the balance catches the half load value and the display flashes.
6. Wait until the flashing stops, place the full load and press **{1/1}** (hold **REF** until "FKT-KEY 3" appears) ⇒ the balance catches the full load value and the display flashes.
7. Wait until the flashing stops and press **{SET}** (hold **REF** until "FKT-KEY 4" appears) ⇒ the balance sets the linearisation factors and displays the full weight value.
8. Empty the weighing pan, tare and check whether at quarter load and three quarter load the linearity stays within the tolerated range. If it doesn't, process chapter 3.11 Adjusting the S-correction" on page 58.
9. **With SCS:** Process chapter "3.12 Adjusting the internal reference weight (SCS only)" on page 59.
Without SCS: Switch the balance off and on (by pressing **ON/OFF** twice) ⇒ the balance quits the service mode.

3.9.3 Checking the linearity

1. Ensure the balance has reached its operating temperature, is levelled and in weighing mode (switch it off and on by pressing **ON/OFF** twice).
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 2x half load, 1x quarter load.
3. Tare the balance and check whether the linearity stays at half, full, quarter and three quarter load within the tolerated range. If it doesn't at half or full load, adjust the linearity again. If it doesn't at quarter or three quarter load, process chapter "3.11 Adjusting the S-correction" on page 58.

3.10 Checking the adjustments of the closed balance

3.10.1 Checking the repeatability (with SCS)


1. Process chapter "3.12 Adjusting the internal reference weight (SCS only)" on page 59.
2. Ensure the balance has reached its operating temperature, is levelled and in weighing mode (switch it off and on by pressing **ON/OFF** twice).
3. Have, according to the specifications, a suitable full load reference weight at disposal.
4. Tare the balance and check whether the repeatability stays within the tolerated range. If it doesn't, process firstly chapter "3.12 Adjusting the internal reference weight (SCS only)" on page 59, then this chapter completely again.

3.10.2 Checking the repeatability (without SCS)

1. Ensure the balance has reached its operating temperature, is levelled and in weighing mode (switch it off and on by pressing **ON/OFF** twice).
2. Have, according to the specifications, a suitable full load reference weight at disposal.
3. Ensure by the configuration menu (see operating instructions) the calibration/adjustment mode is set to "EXTERNAL".
4. Adjust the balance manually with the ready full load weight by holding **TARE** until "CALIBRATION"/"ADJUSTMENT" appears and following the displayed instructions, see also the operating instructions.
5. Tare the balance and check whether the repeatability stays within the tolerated range. If it doesn't, re-adjust the balance (see previous step) and process this chapter completely again.

3.10.3 Checking the corner load

1. Ensure the balance has reached its operating temperature, is levelled and in weighing mode (switch it off and on by pressing **ON/OFF** twice).
2. Have, according to the specifications, a suitable corner load reference weight at disposal.

3. Place the weight on Z and tare the balance.
4. Place the weight on A, B, C and D. Every single value has to stay, according to the specifications, within the tolerated range! If they don't, re-adjust the corner load as described in chapter "3.5.2 Adjusting the corner load (Group M/C/D)" on page 48 or "3.5.3 Adjusting the corner load (Group A)" on page 49, close the balance and process this chapter completely again.


3.10.4 Checking the linearity

1. Process chapter "3.9.3 Checking the linearity" on page 56.

3.11 Adjusting the S-correction

Important: The S-correction only has to be adjusted in case the linearity at half and full load stays within and at quarter and three quarter load stays out of the tolerated range:

Once the linearity is adjusted the balance response behaviour may follow an S-curve, i.e. the linearity at quarter and three quarter load may stay out of the tolerated range and has to be adjusted:

3.11.1 Measuring the S-curve

1. Ensure the balance has reached its operating temperature and is in weighing mode (turn it off and on by pressing **ON/OFF** twice).
2. Have, according to the specifications, suitable linearisation reference weights at disposal: 1x half load, 1x quarter load.
3. Tare the balance, place the quarter load and note the displayed value.
4. Place the three quarter load and note the displayed value.
5. Calculate the differences:
 $\Delta a = (\text{displayed quarter load}) - (\text{placed quarter load})$
 $\Delta b = (\text{displayed three quarter load}) - (\text{placed three quarter load})$

Example:

$$\Delta a = 50.0004g - 50.0000g = 0.0004g$$

$$\Delta b = 149.9998g - 150.0000g = -0.0002g$$

3.11.2 Calculating the S-correction factors

1. Indicate the differences in digits and change their signs:
 $Sa = -\Delta a [d]$
 $Sb = -\Delta b [d]$

Example:

$$\Delta a = 0.0004g = +4 d \Rightarrow Sa = -\Delta a = -4 d$$

$$\Delta b = -0.0002g = -2 d \Rightarrow Sb = -\Delta b = +2 d$$

3.11.3 Entering the S-correction factors (10-Key version only)

1. Ensure the balance has reached its operating temperature, is levelled and set to the service mode, see chapter "3.8 Setting the balance to the service mode (the balance is closed)" on page 54.
2. Hold  until "SCORRECTION", followed by "SET SECTOR-A +0" appears.
3. Enter the S-correction factor Sa with $\uparrow$ and $\downarrow$ and set it with **{SET}** (press $\Leftarrow$).
4. Switch to sector b with **{SECTOR}** (press $\Rightarrow$). "SET SECTOR-B +0" appears.
5. Enter the S-correction factor Sb with $\uparrow$ and $\downarrow$ and set it with **{SET}** (press $\Leftarrow$).
6. Hold  until "BALANCING" appears.
7. Process chapter "3.10 Checking the adjustments of the closed balance" on page 56.

3.11.4 Entering the S-correction factors (5-Key version only)

1. Ensure the balance has reached its operating temperature, is levelled and set to the service mode, see chapter "3.8 Setting the balance to the service mode (the balance is closed)" on page 54.
2. Hold **MODE** until "SCORRECTION", followed by "SET SECTOR-A +0" appears.
3. Enter the S-correction factor Sa
with {UP} (hold **REF** until "FKT-KEY 3" appears)
and {DOWN} (hold **REF** until "FKT-KEY 4" appears)
and set it with {SET} (hold **REF** until "FKT-KEY 1" appears).
4. Switch to sector b with {SECTOR} (hold **REF** until "FKT-KEY 2" appears). "SET SECTOR-B +0" appears.
5. Enter the S-correction factor Sb
with {UP} (hold **REF** until "FKT-KEY 3" appears)
and {DOWN} (hold **REF** until "FKT-KEY 4" appears)
and set it with {SET} (hold **REF** until "FKT-KEY 1" appears).
6. Hold **MODE** until "BALANCING" appears.
7. Process chapter "3.10 Checking the adjustments of the closed balance" on page 56.

3.12 Adjusting the internal reference weight (SCS only)

Important: The internal reference weight has to be adjusted, after the weighing cell has been repaired, or the linearity or the S-correction has been adjusted!

It is possible to vary the full load reference weight within particular step sizes which may be, according to your actual balance, 10, 50, 100 or 500g. This allows an adjustment either with various reference weights at disposal or close to the usual operating range of the balance.

Important: Ensure the reference weight lays within the permitted range of $\pm 0.3d$ of its nominal value! *Example:* A 4200C of a customer, who usually weighs loads of 2.8 kg, asks for 2000g. Its reference weight step size is 500g. A reference weight of 3000g ($2000g + 2 \text{ steps of } 500g$) is the closest to the usual operating range. Place such a reference weight if its mass lays within the range of 2999.997g .. 3000.003g (maximum deviation = $\pm 0.3d = \pm 0.003g$ at $d=0.01g$, see specifications).

1. Ensure the balance has reached its operating temperature, is levelled and set to the service mode, see chapter "3.8 Setting the balance to the service mode (the balance is closed)" on page 54.
2. Have, according to the specifications, a suitable full load reference weight at disposal.
3. Hold **TARE** until "CALIBRATION"/"ADJUSTMENT" appears. The balance starts carrying out an internal adjustment without setting the obtained correction factors, starts thereupon with the external adjustment and asks, as soon as the zero point is measured, for the reference weight. Its value is displayed flashing.
4. Place the asked weight, or better, the one placed at disposal. The balance catches its value and the display flashes. As soon as the flashing stops, the balance asks for the zero point. This is displayed flashing: "--0000 g".
5. Remove the weight. The balance continues with the external adjustment, determines newly the correction factors for the internal reference weight, carries out immediately an internal adjustment and finally turns itself off and on, in order to quit the service mode.

3.13 Error messages

3.13.1 Operating errors; Error disappears by correct use or application

Code	Cause	Remedy
1	Invalid firmware	Load right firmware
3	Load too small	Use higher load
4	Wrong anti theft code	Enter right code
5	7 times consecutively wrong anti theft code	Enter unlock code, contact service
8	Zero point out of range	Repeat linearisation
9.0	Bad adjustment factor	Repeat adjustment
9.1	Zero point deviation exceeded	Repeat adjustment
9.2	Deviation of internal reference weight exceeded	Repeat adjustment
9.3	Bad adjustment factor	Repeat adjustment
30	No internal reference weight	Acknowledge with TARE Adjust internal reference weight
31	No external reference weight	Acknowledge with TARE Adjust external reference weight
47	Too much values	Record less values
48	Not the pristine unit	Set to the pristine unit
49	Too few values	Record more values
53	Values out of the +/- 50% range	Place loads within the range
60	Value being below the converter range	
61	Value exceeding the converter range	

3.13.2 Fatal errors; Balance stays still, error remediable in service mode

Code	Cause	Remedy
10	No linearity values	Repeat linearisation
11	No adjustment factors	Repeat adjustment
21	No SCS temperature compensation values	Contact service
22	No temperature compensation values	Contact service
23	No SCS temperature non-linearity compensation values	Contact service

3.13.3 Hardware errors; Balance stays still

Code	Cause	Remedy
14	Compatibility conflict	Contact service
16	Broken internal RAM	Exchange the μ -controller
17	Broken external RAM	Exchange the RAM
18.1	Wrong firmware checksum	Load the right firmware Exchange the flash memory
18.2	No approved firmware in an approved balance	Load the right firmware Exchange the flash memory
18.9	Wrong μ -controller firmware checksum	Exchange the μ -controller
19	Compatibility conflict	Contact service
20	Erroneous temperature measurement	Check floating sensor connections

3.14 Discussion with the manufacturer

Important: Have at disposal the following information in discussions with the manufacturer:

- Balance model: *e.g. XT 220A*
- Order number: *e.g. 320-9201 (label on the rear of the balance)*
- Serial number: *e.g. 5001057 (label on the rear of the balance)*
- Firmware ID: *e.g. E01-0000.P01 (appears on the display when plugging the balance to the mains)*
- A precise description of the failure
- Weighing results

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