

TROUBLE SHOOTING MANUAL

SINGLE PATIENT DIALYSIS MACHINE

SURDIAL

Manufacturer

NIPRO CORPORATION

3-9-3, Honjo-Nishi, Kita-Ku, OSAKA 531-8510, JAPAN

PHONE +81-6-6372-2331 / FAX +81-6-6371-7422

Version No. 0 R 3 1 0 8 – 0 2 0 8

1. GENERAL DESCRIPTION

<p>When TROUBLE is detected, check symptom and condition and perform</p> <p>TROUBLE SHOOTING</p> <p>to check them for finding trouble cause to determine action to be taken.</p>	<p>Condition, check method, cause and action to be taken for each trouble are described in this Manual.</p> <hr/> <p>⚠ WARNING</p> <hr/> <ul style="list-style-type: none">● Perform check after reading Cautions for Safety thoroughly contained in Service Manual.● When power ON of SURDIAL is required to perform check, take care of electric shock. <hr/>
--	---



<p>In accordance with procedure specified in</p> <p>SERVICE MANUAL,</p> <p>take action and check the result.</p>	<p>Construction drawings and adjusting procedures of each Assy and component, list of components and the like are contained in Service Manual.</p>
---	--

2. ALARM/ERROR INDICATIONS AND DESCRIPTION

NO.	ALARM/ERROR INDICATION	NO.	ALARM/ERROR INDICATION
E01	RECOVER FROM EMERGENCY STOP	E41	VCa CLOSED CIRCUIT LEAK ERROR
E02	T1 TEMPERATURE SENSOR ERROR	E42	VCb CLOSED CIRCUIT LEAK ERROR
E03	T2 TEMPERATURE SENSOR ERROR	E43	VCb CLOSED CIRCUIT LEAK ERROR
E04	T3 TEMPERATURE SENSOR ERROR	E44	AIR ELIMINATION TANK ERROR
E05	DIALYSATE TEMPERATURE MAX ALARM	E45	BLOWN-OUT FUSE DETECTION
E06	DIALYSATE TEMPERATURE MIN ALARM	E46	WATCH-DOG ERROR
E07	DIALYSATE PRESSURE MAXLIMIT	E47	RAM ERROR
E08	DIALYSATE PRESSURE MINLIMIT	E48	EEPROM ERROR
E09	DAIALYSATE PRESSUREMAX ALARM	E49	EPROM ERROR
E10	DIALYSATE PRESSUREMIN ALARM	E50	RAM INITIALIZE
E11	VENOUS PRESSURE MAX LIMIT	E51	EEPROM INITIALIZE
E12	VENOUS PRESSURE MIN LIMIT	E52	SWITCH ERROR
E13	VENOUS PRESSURE MAX ALARM	E53	I/O PORT ERROR
E14	VENOUS PRESSURE MIN ALARM	E54	CONNECTION ERROR SBY0151-SBY0156
E15	ARTERIAL PRESSURE MAX LIMIT	E55	BP BATTERY ERROR
E16	ARTERIAL PRESSURE MIN LIMIT	E56	COUPLER SW ON ALARM(DIALYSIS)
E17	ARTERIAL PRESSURE MAX ALARM	E57	COUPLER SW ON ALARM(PREP.)
E18	ARTERIAL PRESSURE MIN ALARM	E58	COUPLER SW OFF ALARM
E19	TMP ALARM	E59	UF MAX ALARM
E20	DIALYSATE CONDUCTIVITY MAX ALARM	E60	SYRINGE INFUSION COMPLETION
E21	DIALYSATE CONDUCTIVITY MIN ALARM	E61	COMPLETION OF SYRINGE.P INFUSION
E22	BLOOD LEAK SENSOR ERROR	E62	PERFORMED PRIMING
E23	BLOOD LEAK ALARM	E63	PREFORMED TARGET UF.
E24	AIR BUBBLES ARE DETECTED	E64	NO WATER RINSING
E25	HEATER OVER HEAT ALARM	E65	PREPARE DIALYSATE CONC.
E26	SHORT WATER SUPPLY	E66	CONC CONNECTOR TO RINSING PORT.
E27	SHORT DIALYSATE FLOW	E67	RESET SW
E28	DIALYSATE FLOW STOP	E68	UF DIALYSE SW
E29	SN SWITCH TIME ALARM	E69	BYPASS SW
E30	BLOOD PUMP ERROR		ALARM/ERROR DESCRIPTION
E31	BLOOD PUMP STOP ERROR		Disabled LCD indication
E32	CLV POSITION ERROR		Disabled SW function
E33	VISCOUS PUMP DRIVING ERROR(VPa)		Disabled syringe pump operation
E34	VISCOUS PUMP DRIVING ERROR(VPb)		Blood pump error
E35	VP DRAWING VOLUME IS AT MAX LIM.		Level change in venous drip chamber
E36	VCa CLOSED CIRCUIT LEAK ERROR		Air inclusion in dialysate
E37	VCa CLOSED CIRCUIT LEAK ERROR		Unstable dialysate flow rate
E38	VCb CLOSED CIRCUIT LEAK ERROR		Disabled discharge by pressing drain SW
E39	VCb CLOSED CIRCUIT LEAK ERROR		UFR error exceeding allowable upper/lower limits
E40	VCa CLOSED CIRCUIT LEAK ERROR		Date and time are initialized.

● PREFACE

Thank you for purchase of NIPRO Single Patient Dialysis Machine “SURDIAL”.

Information for serviceman or technician required to perform check and maintenance of SURDIAL is given by this Trouble Shooting Manual. All operations described in this Manual should be performed by only qualified serviceman or technician who should be familiar with replacement of components of microprocessor based electronic devices and should have basic knowledge concerning conventional dialysate line of dialysis machine. The serviceman or technician should be familiar with all operations contained in this Manual by reading it thoroughly to full understand them before starting trouble shooting.

Use this Manual together with Operation and Service Manuals of SURDIAL. Same descriptions contained in them are minimal.

Generally, information required to operate SURDIAL is contained in Operation Manual of SURDIAL though technical description is also contained partially. Information related to maintenance, component replacement and various settings is contained in Service Manual of SURDIAL. In this Manual, information related to method of trouble shooting is contained.

Description concerning the trouble shooting contained in this Manual is protected by copyright owned by the manufacturer. Partial or whole copy, duplication or preparation of version in any foreign language of this Manual should not be allowed without obtaining prior written approval of the manufacturer.

E01 RECOVER FROM EMERGENCY STOP					
Detection	Recovery to normal from emergency machine stop due to power failure caused during execution of sequence				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse •Preparation •Prep. end •Dialysis •Recovery	○ (3 times)	Execution of sequence before error detection	—	—	—
Reset	Automatic recovery to normal		Note	Becoming error free of upper/lower alarm limits of venous, arterial and dialysate pressures by error detection	

Condition or cause	Check and action
Recovery from power failure	Not required
Intentional OFF-ON of power breaker	Not required
Intentional OFF-ON of power SW	Not required
OFF of breaker due to leakage or short circuit of AC230V (AC110V) line	<p>Check : Identify position(s) of leakage or short-circuit of AC230V (AC110V) line.</p> <p>Disconnect components one by one from the line to identify damaged one(s) by testing or operate components in the line one by one manually to do so.</p> <p style="text-align: center;">⚠ CAUTION</p> <p><u>Check correct earth connection of SURDIAL and disconnect components after disconnecting power plug from receptacle.</u></p> <p>Action: Repair or replace damaged one(s).</p>
Incorrect power plug connection	Action: Connect the plug to receptacle securely.
Incorrect power line connection	<p>Check: Pull lightly each of harnesses connecting POWER PLUG— SWITCH GEAR— TERMINAL BLOCK2— POWER SWITCH— NOISE FILTER— TERMINAL BLOCK1— POWER SUPPLY— SUB-BOARD CN412— SUB-BOARD CN413— MAIN BOARD CN1 in series to check loosening, incorrect contact and disconnection.</p> <p>Action : If the loosening is detected, tighten securely. If the incorrect contact and/or disconnection are/is detected, replace relevant harness(es) with new one(s).</p>

Condition or cause	Check and action
<p>Repetition of error of recovery from emergency stop</p>	<p>Check: Identify whether error cause(s) is (are) main board, sub-board and/or DC+5V/±12V switching power supply. Check that actual output voltage is $\geq 4.8\text{V}$ (nominal: 5VDC) when measured at test points on main and sub-board.</p> <p>Action :</p> <p>When $\geq 4.8\text{V}$ is not outputted to the main board and sub-board:</p> <ul style="list-style-type: none"> • Replace sub-board with new one or the switching power supply with new one. <p>When it is outputted to sub-board, while it is not outputted to main board:</p> <ul style="list-style-type: none"> • Replace main board with new one. <p>When it is outputted to the main board and sub-board:</p> <ul style="list-style-type: none"> • Replace main board with new one.

E02 T1 TEMPERATURE SENSOR ERROR

Detection	T1 temperature exceeding measurable limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ •↑	•Continuance •↑ •↑ •↑ •↑ •Dialysate line stop •↑ •↑ •↑ •↑	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Indication error of temperature and other analog system	<p>Check: Check whether indications of all temperature sensors and other sensors such as dialysate pressure sensor are incorrect. Check stable output of voltage of $\geq 4.9V$ (nominal: 5Vref) at test point on main board.</p> <p>Action: Replace main board with new one.</p>
Incorrect [T 1] terminal connection	<p>Check: Check no incorrect wiring, loosening, corrosion and/or disconnection at TM1 T1 terminal of main board.</p> <p>Action: If loosened, tighten securely. If incorrectly wired, correct. If corroded or disconnected, replace [T1] with new one.</p>
Internal disconnection or short circuit of [T1]	<p>Check: Remove TM1 T1 terminal from main board and measure resistance between Red and White to find short circuit, disconnection and fluctuation.</p> <p>Action: Replace [T1] with new one.</p>
Introduction of high temperature water from water supply source	<p>Action: Overhaul dialysate line to correct potential damage.</p>
Incorrect indicated temperature calibration of main board	<p>Action: Calibrate temperatures indicated by [T1], [T2] and [T3] circuits of main board using temperature calibrator. If the calibration is impossible, replace main board with new one.</p>

E03 T2 TEMPERATURE SENSOR ERROR

Detection	T2 temperature exceeding measurable limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
• Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery	• ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑	• Continuance • ↑ • ↑ • ↑ • ↑ • Dialysate line stop • ↑ • ↑ • ↑ • ↑	—	—	○
Reset	Press [RESET] SW.		Note	After Ver 3, measurable extent of upper limit while rinse isn't detectable.	

Condition or cause	Check and action
Indication error of temperature and other analog system	Check: Check whether indications of all temperature sensors and other sensors such as dialysate pressure sensor are incorrect. Check stable output of voltage of $\geq 4.9V$ (nominal: 5Vref) at test point on main board. Action: Replace main board with new one.
Incorrect [T2] terminal connection	Check: Check no incorrect wiring, loosening, corrosion and/or disconnection at TM1 T2 terminal of main board. Action: If loosened, tighten securely. If incorrectly wired, correct. If corroded or disconnected, replace [T2] with new one
Internal disconnection or short circuit of [T2]	Check: Remove TM1 T2 terminal from main board and measure resistance between Red and White to find short circuit, disconnection and fluctuation. Action: Replace [T2] with new one.
Incorrect indicated temperature calibration of main board	Action: Calibrate temperatures indicated by [T1], [T2] and [T3] circuits of main board using temperature calibrator. If the calibration is impossible, replace main board with new one.

E04 T3 TEMPERATURE SENSOR ERROR

Detection	T3 temperature exceeding measurable limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ •↑	•Continuance •↑ •↑ •↑ •↑ •Dialysate line stop •↑ •↑ •↑ •↑	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Indication error of temperature and other analog system	<p>Check: Check whether indications of all temperature sensors and other sensors such as dialysate pressure sensor are incorrect. Check stable output of voltage of $\geq 4.9V$ (nominal: 5Vref) at test point on main board.</p> <p>Action: Replace main board with new one.</p>
Incorrect [T3] terminal connection	<p>Check: Check no incorrect wiring, loosening, corrosion and/or disconnection at TM1 T3 terminal of main board.</p> <p>Action: If loosened, tighten securely. If incorrectly wired, correct. If corroded or disconnected, replace [T2] with new one</p>
Internal disconnection or short circuit of [T3]	<p>Check: Remove TM1 T3 terminal from main board and measure resistance between Red and White to find short circuit, disconnection and fluctuation.</p> <p>Action: Replace [T3] with new one.</p>
Incorrect indicated temperature calibration of main board	<p>Action: Calibrate temperatures indicated by [T1], [T2] and [T3] circuits of main board using temperature calibrator. If the calibration is impossible, replace main board with new one.</p>

E05 DIALYSATE TEMPERATURE MAX ALARM

Detection	Dialysate temperature T3 of \geq upper alarm limit during heater ON				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Preparation •Prep. end •Dialysis •Recovery	○ (Continuous)	Auto bypass Heater OFF	—	—	○
Reset	Press [RESET] SW to reset. Auto reset is caused after release from alarm condition for 10 seconds.		Note		

Condition or cause	Check and action
Actually higher dialysate temperature ○ Failure of main board ○ Incorrect calibration of indicated temperature of main board ○ Higher supply water temperature ○ Lower priming flow rate	<p>Check: Check whether actual dialysate temperature is higher than upper alarm limit by touching dialysate supply line with hand when the line is connected to dialyzer or by measurement of temperature of dialysate discharged from coupler outlet using thermometer when patient is not connected to the machine.</p> <p>Action: Replace main board with new one.</p> <p>Action: Calibrate temperatures indicated by [T1], [T2] and [T3] circuits of main board using temperature calibrator. If the calibration is impossible, replace main board with new one.</p> <p>Action: Adjust supply water temperature as specified.</p> <ul style="list-style-type: none"> • When heat exchanger is used : 5 to 30°C • When heat exchanger is not used : 17 to 30°C <p>Action: Adjust priming flow rate to 700 to 750mL/min.</p>

<p>Actually acceptable dialysate temperature</p> <p>○ Failure of T3 temperature sensor</p> <p>○ Incorrect calibration of main board</p>	<p>Check: Check whether actual dialysate temperature is acceptable by touching dialysate supply line with hand when the line is connected to dialyzer or by measurement of temperature of dialysate discharged from coupler outlet using thermometer when patient is not connected to the machine.</p> <p>Check: Check whether indicated temperature is fluctuated when harness of T3 temperature sensor is shaken. Check also whether indicated T3 temperature is deviated when temperature calibrator is connected.</p> <p>Action: If check result is unacceptable, replace T3 temperature sensor with new one.</p> <p>Check: Check whether indicated T3 temperature is deviated when temperature calibrator is connected.</p> <p>Action: Calibrate [T3] circuit of main board using temperature calibrator. If the calibration is impossible, replace main board with new one.</p>
---	--

E06 DIALYSATE TEMPERATURE MIN ALARM					
Detection	Dialysate temperature T3 of \leq lower alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (3 times)	Auto bypass	—	—	○
Reset	Press [RESET] SW to reset. Auto reset is caused after release from alarm condition for 10 seconds.		Note		

Condition or cause	Check and action
--------------------	------------------

<p>Actually lower dialysate temperature</p> <ul style="list-style-type: none"> ○ Failure of main board ○ Incorrect calibration of indicated temperature of main board ○ Lower supply water temperature ○ Failure of heater ○ Higher priming flow rate 	<p>Check: Check whether actual dialysate temperature is lower than lower alarm limit by touching dialysate supply line with hand when the line is connected to dialyzer or by measurement of temperature of dialysate discharged from coupler outlet using thermometer when patient is not connected to the machine.</p> <p>Action: Replace main board with new one.</p> <p>Action: Calibrate temperatures indicated by [T1], [T2] and [T3] circuits of main board using temperature calibrator.</p> <p>If the calibration is impossible, replace main board with new one.</p> <p>Action: Adjust supply water temperature as specified.</p> <ul style="list-style-type: none"> • When heat exchanger is used : 5 to 30°C • When heat exchanger is not used : 17 to 30°C <p>Check: Measure resistance between both ends of heater at room temperature after power OFF.</p> <p>Action: Replace heater if the resistance is out of the below ranges:</p> <p>230V spec. : 38.3 to 44.1 Ω</p> <p>110V spec. : 8.7 to 10.0 Ω</p> <p>Action: Adjust priming flow rate to 700 to 750ml/min.</p>
<p>Actually acceptable dialysate temperature</p> <ul style="list-style-type: none"> ○ Failure of T3 temperature sensor ○ Incorrect calibration of main board 	<p>Check: Check whether actual dialysate temperature is acceptable by touching dialysate supply line with hand when the line is connected to dialyzer or by measurement of temperature of dialysate discharged from coupler outlet using thermometer when patient is not connected to the machine.</p> <p>Check: Check whether indicated temperature is fluctuated when harness of T3 temperature sensor is shaken.</p> <p>Check also whether indicated T3 temperature is deviated when temperature calibrator is connected.</p> <p>Action: If check result is unacceptable, replace T3 temperature sensor with new one.</p> <p>Check: Check whether indicated T3 temperature is deviated when temperature calibrator is connected.</p> <p>Action: Calibrate [T3] circuit of main board using temperature calibrator.</p> <p>If the calibration is impossible, replace main board with new one.</p>

E07 DIALYSATE PRESSURE MAX LIMIT

Detection	Dialysate pressure of $\geq +400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑ 	UF stop ON of V4 and V10 for 2 sec.	2 sec.	2 sec.	○
Reset	Auto reset of BP, SP and CLV Press [RESET] SW to reset UF stop.		Note		

Condition or cause	Check and action
Pressure rise due to failure of [PR1]	Action: Adjustment of [PR1] In adjustment is impossible, repair or replace [PR1].
Leak from [V10]	Action: Repair or replacement of [V10]
Error occurrence at gas purge ○ Occlusion of [R2] ○ Excess tightening of [PC2] ○ Impossible opening of [V4]	Action: Change [R2] so that dialysate pressure becomes 0 to 100mmHg during gas purge. Action: Adjust [PC2] so that dialysate pressure becomes 0 to 100mmHg during gas purge. Check: Check whether [V4] functions when it is turned on or off in maintenance mode 4 screen. If not functions, remove connector from [V4] and check whether 24V is outputted to the connector when [V4] is ON in maintenance mode 4 screen. Action: <ul style="list-style-type: none"> • If 24V is not outputted, replace main board with new one. • If it is outputted, replace [V4] with new one.
Pressure rise due to folding of blood line	Action: Correct folding.

<p>Failure of dialysate pressure sensor board or the sensor</p>	<p>Check: Disconnect hose from the sensor Assy to expose to atmosphere for checking whether indicated pressure is nearly 0mmHg.</p> <p>Action: If the indicated pressure is not nearly 0mmHg, calibrate the sensor board.</p> <p>If calibration is impossible, replace the board or sensor with new one.</p>
---	--

E07 DIALYSATE PRESSURE MAX LIMIT

Condition or cause	Check and action
<p>Disconnection of harness from the sensor Assy or insufficient contact of connector of it</p>	<p>Check: Pull lead wire (harness) lightly to check whether it is separated or loosened.</p> <p>Action: Repair</p>
<p>Drain tube is blocked.</p> <p>○ No folding/kinking of Drain tube</p> <p>○ Drain tube tip is placed too higher position.</p> <p>○ Drain tube is too long.</p>	<p>Check : Remove drain tube from drain port and check drain liquid is drained from port correctly or not.</p> <p>Action : Remove the trouble cause.</p> <p>Action : Drain tube tip must be placed above the machine floor by 0 to 60 cm.</p> <p>Action : Use tube of an inter diameter of 8 mm and a length of 2 m or shorter to connect to the port.</p>

E08 DIALYSATE PRESSURE MIN LIMIT

Detection	Dialysate pressure of $\leq -400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑ 	UF stop ON of V4 and V10 for 2 sec.	2 ses.	2 sec.	○
Reset	Auto reset of BP, SP and CLV Press [RESET] SW to reset UF stop.		Note		

Condition or cause	Check and action
Occlusion of dialysate line to dialyzer	Check: Check whether the line is folded or clamped by forceps visually. Action: Correct folding or clamp.
Excess occlusion in dialyzer	Action: Replace dialyzer with new one.
Leak from [V4]	Check: Check stop of leak by clamping [V4] line with forceps. Action: Clean [V4] by rinsing after disassembling or replace it with new one.
Opening of [V4] due to malfunction of [FSW]	Action: Clean air-elimination tank [AS2] by rinsing after disassembling or replace it with new one.
Occlusion of [RV2]	Action: Adjust pressure in closed line. If adjustment is impossible, Clean [RV2] by rinsing after disassembling or replace it with new one.
Error occurrence during gas purge ○ Failure of [V10]	Check: Check whether [V10] functions when it is turned on or off in maintenance mode 4 screen. If not functions, disconnect connector from [V10] and turn on it in maintenance mode 4 screen to check whether 24V is outputted to the connector. Action: If 24V is not outputted, replace main board with new one. If it is outputted, replace [V10] with new one.

Failure of dialysate pressure sensor board or the sensor	<p>Check: Disconnect hose from the sensor Assy to expose to atmosphere for checking whether indicated pressure is nearly 0mmHg.</p> <p>Action: If the indicated pressure is not nearly 0mmHg, calibrate the sensor board. If calibration is impossible, replace the board or sensor with new one.</p>
Disconnection of harness from the sensor Assy or insufficient contact of connector of it	<p>Check: Pull lead wire (harness) lightly to check whether it is separated or loosened.</p> <p>Action: Repair</p>
Low water supply flow	Action : Increase water supply flow rate more than 750 mL/min.
No water supply	Check/Action : See E26 SHORT WATER SUPPLY.
Low water supply pressure	Check/Action : See E26 SHORT WATER SUPPLY.
Alarm while gaspurge	Check/Action : See E26 SHORT WATER SUPPLY.
On the acid rinse without an acid rinse option.	Check/Action : See E26 SHORT WATER SUPPLY.

E09 DIALYSATE PRESSURE MAX ALARM

Detection	Dialysate pressure of \geq set upper alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop ON of V4 and V10 for 2 sec.	2 sec.	2 sec.	○
Reset	Auto reset of BP, SP and CLV Press [RESET] SW to reset UF stop.		Note	By the alarm actuation, upper/lower alarm limits of venous pressure become error free.	

Condition or cause	Check and action
Simultaneous actuation of this alarm with venous pressure max alarm	Check: See E13 VENOUS PRESSURE MAX ALARM. Action: See E13 VENOUS PRESSURE MAX ALARM.
Dialysate pressure rise due to UF stop	Action: Remove cause of UF stop.
Dialysate pressure rise due to failure of [PR1]	Action: Adjust [PR1]. If adjustment is impossible, repair [PR1] or replace it with new one.
Leak from [V10]	Action: Repair [V10] or replace it with new one.
Leak from [V5] and [V7]	Action: Repair [V5] and [V7] or replace them with new ones.

E10 DIALYSATE PRESSURE MIN ALARM

Detection	Dialysate pressure of \leq set lower alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop ON of V4 and V10 for 2 sec.	2 sec.	2 sec.	○
Reset	Auto reset of BP, SP and CLV Press [RESET] SW to reset UF stop.		Note	By the alarm actuation, upper/lower alarm limits of venous pressure become error free.	

Condition or cause	Check and action
Simultaneous actuation of this alarm with venous pressure min alarm	Check: See E14 Venous pressure min alarm. Action: See E14 Venous pressure min alarm.
Occlusion of [RV2]	Action: Adjust pressure in closed line. If adjustment is impossible, clean [RV2] by rinsing after disassembling or replace it with new one.
Leak from [V4]	Check: Check no dialysate pressure change by clamping [V4] line with forceps. Action: In the case of no change, clean [V4] by rinsing after disassembling or replace it with new one.
Opening of [V4] due to failure of [FSW] The opening may be simultaneously caused with E44 AIR ELIMINATION TANK ERROR	Check: See E44 AIR ELIMINATION TANK ERROR Action: See E44 AIR ELIMINATION TANK ERROR
Leak from [V6] and [V8]	Action: Repair [V6] and [V8] or replace them with new ones.
Excess occlusion in dialyzer	Action: Replace dialyzer with new one if required.

E11 VENOUS PRESSURE MAX LIMIT

Detection	Venous pressure of $\geq +400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ 	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note		

Condition or cause	Check and action
Folding or occlusion of blood line between Venous Drip Chamber and patient or occlusion of the line with Robert clamp	Action: Remove the trouble cause.
Thrombus formed at mesh of venous drip chamber or at shunt connector	Action: Remove thrombus or replace blood line with new one.
Error occurrence during single needle dialysis	Action: Set upper pressure change limit of SN correctly.
Failure of venous pressure sensor board	Check: Disconnect blood line from venous pressure inlet port to expose to atmosphere and check indicated venous pressure is nearly 0mmHg. Action: If nearly 0mmHg is not indicated, calibrate the board. If calibration is impossible, replace the board with new one.

E12 VENOUS PRESSURE MIN LIMIT

Detection	Venous pressure of $\leq -400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ 	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note		

Condition or cause	Check and action
Simultaneous occurrence of this error with dialysate related error	Action: Remove cause of dialysate related error.
Simultaneous occurrence of this error with arterial pressure related error	Action: Remove cause of arterial pressure related error
Press of DRAIN SW without removing coupler of dialyzer under occluding venous line with Robert clamp	Note: DRAIN SW operation-miss
Failure of venous pressure sensor board	Check: Disconnect blood line from venous pressure inlet port to expose to atmosphere and check indicated venous pressure is nearly 0mmHg. Action: If nearly 0mmHg is not indicated, calibrate the board. If calibration is impossible, replace the board with new one.

E13 VENOUS PRESSURE MAX ALARM

Detection	Venous pressure of \geq set upper alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note	By the alarm actuation, upper/lower alarm limits of dialysate pressure become error free.	

Condition or cause	Check and action
Occlusion of venous line.	Check: Check shunt occlusion or blood line folding. Action: Remove the trouble.
Occlusion of arterial line	Check: Check shunt occlusion or blood line folding. Action: Remove the trouble.
Venous pressure change due to patient movement	Action: If the change is significant, raise somewhat upper alarm limit of venous pressure at discretion of doctor.
Significant discharge pressure change of blood pump	Check: Check whether discharge flow rate, pressure applied to rolling tube and rotor rotation of blood pump are normal. Action: Adjust the occlusion or the pump rpm correctly. If adjustment is impossible, replace the pump or sub-board with new one.

E14 VENOUS PRESSURE MIN ALARM

Detection	Venous pressure of \leq set lower alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note	By the alarm actuation, upper/lower alarm limits of dialysate pressure become error free.	

Condition or cause	Check and action
Blood pressure drop of patient	Action: Stop UF and monitor patient condition. Correct if set lower alarm limit of venous pressure is incorrect.
Dislodgment of venous line from venous pressure inlet port	Check: Check no introduction of blood in venous pressure sensor line. Action: When blood is not introduced in the line, connect venous line to the port. If introduced, replace the sensor line with new one.
Occlusion of venous line	Check: Check shunt occlusion or venous line folding. Action: Remove the trouble.
Occlusion of arterial line.	Check: Check shunt occlusion or arterial line folding. Action: Remove the trouble.
Venous pressure change due to patient movement	Action: If the change is significant, lower somewhat lower alarm limit of venous pressure at discretion of doctor.
Significant discharge pressure change of blood pump	Check: Check whether discharge flow rate, pressure applied to rolling tube and rotor rotation of blood pump are normal. Action: Adjust the occlusion or the pump rpm correctly. If adjustment is impossible, replace the pump or sub-board with new one.

E15 ARTERIAL PRESSURE MAX LIMIT

Detection	Arterial pressure of $\geq +400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ 	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note		

Condition or cause	Check and action
Simultaneous actuation of this error with dialysate pressure related error	Action: Remove cause of dialysate pressure related error.
Simultaneous actuation of this error with venous pressure related error	Action: Remove cause of venous pressure related error.
Folding or occlusion of blood line between arterial Drip chamber and dialyzer or occlusion of the line with Robert clamp	Action: Remove the trouble cause.
Folding or occlusion of blood line between venous Drip chamber and patient or occlusion of the line with Robert clamp	Action: Remove the trouble cause.
Thrombus formed at mesh of arterial drip chamber or	Action: Remove thrombus or replace blood line with new one.
Failure of arterial pressure sensor board	Check: Disconnect blood line from arterial pressure inlet port to expose to atmosphere and check indicated arterial pressure is nearly 0mmHg. Action: If nearly 0mmHg is not indicated, calibrate the board. If calibration is impossible, replace the board with new one.

E16 ARTERIAL PRESSURE MIN LIMIT

Detection	Arterial pressure of $\leq -400\text{mmHg}$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ 	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note		

Condition or cause	Check and action
Simultaneous actuation of this error with dialysate pressure related error	Action: Remove cause of dialysate pressure related error.
Thrombus formed at arterial shunt connector	Action: Remove thrombus or replace blood line with new one.
Folding or occlusion of blood line between arterial drip chamber and patient or occlusion of the line with Robert clamp	Action: Remove the trouble cause.
Failure of arterial pressure sensor board	Check: Disconnect blood line from arterial pressure inlet port to expose to atmosphere and check indicated arterial pressure is nearly 0mmHg. Action: If nearly 0mmHg is not indicated, calibrate the board. If calibration is impossible, replace the board with new one.

E17 ARTERIAL PRESSURE MAX ALARM					
Detection	Arterial pressure of \geq set upper alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note	By the alarm actuation, set upper/lower alarm limits of dialysate pressure become error free.	

Condition or cause	Check and action
Occlusion of venous line.	Check: Check shunt occlusion and blood line folding. Action: Correct the trouble.
Occlusion of arterial line.	Check: Check shunt occlusion and blood line folding. Action: Correct the trouble.
Arterial pressure change due to patient movement	Action: If the change is significant, raise set upper alarm limit of arterial pressure somewhat at discretion of doctor.
Significant discharge pressure change of blood pump	Check: Check whether discharge flow rate, pressure applied to rolling tube and rotor rotation of blood pump are normal. Action: Adjust the occlusion or the pump rpm correctly. If adjustment is impossible, replace the pump or sub-board with new one.

E18 ARTERIAL PRESSURE MIN ALARM

Detection	Arterial pressure of \leq set lower alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	UF stop	○	—	○
Reset	Press [RESET] SW to recover.		Note	By the alarm actuation, set upper/lower alarm limits of dialysate pressure become error free.	

Condition or cause	Check and action
Blood pressure drop of patient	Action: Stop UF and monitor patient condition. Correct if set lower alarm limit of venous pressure is incorrect.
Dislodgment of blood line form arterial pressure inlet port	Check: Check no introduction of blood in arterial pressure sensor line. Action: When blood is not introduced in the line, connect arterial line to the port. If introduced, replace the arterial line with new one.
Occlusion of venous line.	Check: Check shunt occlusion and blood line folding. Action: Correct the trouble.
Occlusion of arterial line.	Check: Check shunt occlusion and blood line folding. Action: Correct the trouble.
Arterial pressure change due to patient movement	Action: If the change is significant, raise set upper alarm limit of arterial pressure somewhat at discretion of doctor.
Significant discharge pressure change of blood pump	Check: Check whether discharge flow rate, pressure applied to rolling tube and rotor rotation of blood pump are normal. Action: Adjust the occlusion or the pump rpm correctly. If adjustment is impossible, replace the pump or sub-board with new one

E19 TMP ALARM					
Detection	TMP (VP - DP) of \geq set alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Dialysis •Recovery	○ (Continuous)	UF stop	○	○	○
Reset	Press [RESET] SW to recover.		Note		

Condition or cause	Check and action
Excess occlusion in dialyzer	Action: Replace dialyzer with new one at discretion of doctor.
Lower setting of TMP alarm limit	Action: Set TMP to appropriate value.
E09 DIALYSATE PRESSURE MAX ALARM E10 DIALYSATE PRESSURE MIN ALARM E13 VENOUS PRESSURE MAX ALARM Simultaneous actuation of this alarm with E14 VENOUS PRESSURE MIN ALARM	Action: See E09 DIALYSATE PRESSURE MAX ALARM See E10 DIALYSATE PRESSURE MIN ALARM See E13 VENOUS PRESSURE MAX ALARM See E14 VENOUS PRESSURE MIN ALARM
E07 DIALYSATE PRESSURE MAX LIMIT E08 DIALYSATE PRESSURE MIN LIMIT E11 VENOUS PRESSURE MAX LIMIT Simultaneous actuation of this alarm with E12 VENOUS PRESSURE MIN LIMIT	Action: See E07 DIALYSATE PRESSURE MAX LIMIT. See E08 DIALYSATE PRESSURE MIN LIMIT. See E11 VENOUS PRESSURE MAX LIMIT. See E12 VENOUS PRESSURE MIN LIMIT.

E20 DIALYSATE CONDUCTIVITY MAX ALARM

Detection	Dialysate conductivity of \geq upper alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	<input type="radio"/> (Continuous)	Auto bypass	—	—	—
Reset	Press [RESET]SW. Auto reset after release from alarm condition for 10 sec.		Note		

Condition or cause	Check and action
Abnormal actual conductivity	Check: Take dialysate sample from [SL] to measure actual conductivity by specified method.
<input type="radio"/> No adjustment of conductivity	Action: Set conductivity correctly and calibrate indicated conductivity in maintenance mode 5 screen.
<input type="radio"/> Failure of main board	Check: See E33 VISCOUS PUMP DRIVING ERROR(VPa) and E34 VISCOUS PUMP DRIVING ERROR(VPb). Action: See E33 VISCOUS PUMP DRIVING ERROR(VPa) and E34 VISCOUS PUMP DRIVING ERROR(VPb).
<input type="radio"/> Failure of [V13] or [V14]	Action: Repair [V13] or [V14] (V15)) or replace it with new one.
<input type="radio"/> Failure of [V2] or [V3]	Action Repair [V2] or [V3] or replace it with new one.
<input type="radio"/> Higher or lower priming flow rate	Action: Adjust the flow rate to 700-750mL/min.
<input type="radio"/> Failure of [C2]	Action: Replace [C2] with new one.
<input type="radio"/> Folding, occlusion or collapsing of concentrate line tube	Action: Correct the trouble.

E20 DIALYSATE CONDUCTIVITY MAX ALARM

Condition or cause	Check and action
<p>Normal actual conductivity</p> <p>○[CD] electrode contaminated with dirt such as CaCO_3 or leak from conductivity sensor</p> <p>○Incorrect [CD] terminal connection</p> <p>○Internal disconnection or short circuit of [CD] harness</p> <p>○Failure of main board</p>	<p>Check: Take dialysate sample from [SL] to measure conductivity by specified method.</p> <p>Action: Rinse the electrode with acetic acid followed by water rinsing sufficiently. If the dirt can not be removed, use specified detergent such as CARON or SYUNMA to remove it after disassembling conductivity sensor.</p> <p>Check: Check no miss-wiring, loosening, corrosion and disconnection of TM1 CD terminal of main board.</p> <p>Action: If loosened, tighten the connection securely. If miss-wired, correct. If corroded or disconnected, replace harness with new one.</p> <p>Check: After disconnecting TM1 CD terminal from main board, measure resistance between [CD] and [GND] to check no short circuit, disconnection and fluctuation of measured resistance.</p> <p>Action: Replace [CD] harness with new one.</p> <p>Check: Change ambient temperature of main board to check whether conductivity and temperature indications are changed.</p> <p>Action: If incorrectly indicated, replace main board with new one.</p>
<p>When Na infusion unit (option) is used:</p> <p>○Higher or lower concentration of Na concentrate</p> <p>○ Use of other chemical added concentrate</p> <p>○Use of other concentrate</p>	<p>Action: Use 10% Na concentrate.</p> <p>Action: Set proper upper alarm limit of conductivity by considering electrolyte concentration.</p> <p>Action: Use 10% Na concentrate.</p>

E21 DIALYSATE CONDUCTIVITY MIN ALARM					
Detection	Dialysate conductivity of \leq lower alarm limit				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	Auto bypass	—	—	—
Reset	Press [RESET]SW. Auto reset after release from alarm condition for 10 sec.		Note		

Condition or cause	Check and action
See E20 DIALYSATE CONDUCTIVITY MAX ALARM.	See E20 DIALYSATE CONDUCTIVITY MAX ALARM.

E22 BLOOD LEAK SENSOR ERROR					TYPE I (RK-GR1)
Detection	<ul style="list-style-type: none"> •Absolute blood leak REAL value of ≥ 201 in PREP. or PREP. END process •Absolute blood leak REAL value of < 201 in DIALYSIS or RECOVERY process •Signal of RINSE alarm 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	Line function stop	—	—	○
Reset	Press [RESET]SW.		Note	BLD SENS VALUE :No7 in maintenance mode 1 Cancellation of detection alarm is available with ROM Ver3.30 at maintenance mode 2.	

E22 BLOOD LEAK SENSOR ERROR					TYPE II (RK-GR2)
Detection	<ul style="list-style-type: none"> •Absolute blood leak REAL(BLD AUTO ZERO) value of ≥ 201 in PREP. or PREP. END process •Signal of RINSE alarm 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Preparation •Prep. end •Dialysis •Recovery 	○ (Continuous)	Line function stop	—	—	○
Reset	Press [RESET]SW		Note	•BLD AUTO ZERO:No6 in maintenance mode1 •Cancellation of detection alarm is available with ROM Ver3.30 at maintenance mode 2.	

E22 BLOOD LEAK SENSOR ERROR					TYPE II (BLD-01)
Detection	<ul style="list-style-type: none"> •Absolute blood leak REAL(BLD AUTO ZERO) value of ≥ 201 in PREP. or PREP. END process 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Preparation •Prep. end 	○ (Continuous)	Line function stop	—	—	○
Reset	Press [RESET]SW		Note	•BLD AUTO ZERO:No6 in maintenance mode1 •Cancellation of detection alarm is available with ROM Ver3.30 at maintenance mode 2.	

Condition or cause	Check and action
Excess leak from dialyzer	<p>Action: Replace dialyzer with new one.</p> <p style="text-align: center;">⚠CAUTION</p> <p><u>After setting to error free condition, error is continued when leaked blood is retained in drain side of viscous chamber.</u> <u>In such case, release the alarm to operate the machine until the blood is discharged. If required, clean leak sensor window by wiping with swab or the like after opening the sensor cover.</u></p>
Excess dirt of blood leak sensor window (RK-GR1,RK-GR2,BLD-01)	<p>Action: Clean the window by wiping with swab or the like after opening the sensor cover.</p>
Inclusion of large amount of air in blood leak sensor	<p>Check: Take dialysate sample from drain port of the machine or dialysate outlet of dialyzer to check blood in the sample visually or by use of specified tester.</p> <p>Check: Check visually whether fine air bubbles impossible to be removed by [AS1] and [AS2] remain in them and whether negative pressure line (example: pump inlet) in the machine is leaked.</p> <p>Action: Repair damaged part or replace damaged component with new one.</p>
Incorrect connection of blood leak sensor connector	<p>Check: Pull harness lightly to find disconnection and loose connection and check defect of connector.</p> <p>Action: Replace harness with new one.</p>
Deterioration or failure of LED of blood leak sensor	<p>Action: Readjust the sensor. If adjustment is impossible, replace the sensor Assy with new one.</p>
Leak of silicone oil from chamber	<p>Action. Replace chamber Assy with new one.</p> <p>Clean chamber inside with detergent to remove dirt formed by contact with the oil from there.</p> <p>Replace silicone oil line and/or connectors contaminated by the oil leak with new one(s) when required.</p>

E23 BLOOD LEAK ALARM					
Detection	BLD REAL VALUE of \geq set alarm limit (BLD REAL VALUE: No.5 of maintenance mode 1, Alarm limit: No.12 of maintenance mode 2)				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	UF stop Auto bypass	○	—	○
Reset	Press [RESET]SW.		Note		

Condition or cause	Check and action
Blood leak	<p>Check: Take dialysate sample from drain port of the machine or dialysate outlet of dialyzer to check blood in the sample visually or by use of specified tester.</p> <p>Action: If leaked from dialyzer: After shifting to RECOVERY process, press BYPASS SW and replace dialyzer with new one followed by priming it.</p>

<p>No blood leak</p> <p>○ Air inclusion</p> <p>○ Dirt of blood leak sensor window (Deposition of protein and/or CaCO₃)</p> <p>○ Deterioration or failure of LED of blood leak sensor</p> <p>○ Failure of main board</p> <p>○ Leak of silicone oil</p>	<p>Check: Take dialysate sample from drain port of the machine or dialysate outlet of dialyzer to check blood in the sample visually or by use of specified tester.</p> <p>Check: Check visually whether fine air bubbles impossible to be removed by [AS1] and [AS2] remain in them and whether negative pressure line (example: pump inlet) in the machine is leaked.</p> <p>Action: Repair damaged part or replace damaged component with new one.</p> <p>Action: Clean the window.</p> <p>Action: Readjust the sensor. If adjustment is impossible, replace the sensor Assy with new one.</p> <p>Action: Replace the sensor with new one and check correct function of it.. If the replacement is ineffective, replace main board with new one.</p> <p>Check: Remove cover from the sensor and check whether dialysate is turbid or is mixed with the oil.</p> <p>Action. Replace chamber Assy with new one. Clean chamber inside with detergent to remove dirt formed by contact with the oil from there. Replace silicone oil line and/or connectors contaminated by the oil leak with new one(s) when required.</p>
--	--

E24 AIR BUBBLES ARE DETECTED

Detection	Air bubbles when BUBBLE [ON/OFF] SW is ON				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑ 	UF stop	<ul style="list-style-type: none"> ○ SW-OFF	<ul style="list-style-type: none"> ○ 	<ul style="list-style-type: none"> ○
Reset	Press [RESET] SW to reset UF stop and CLV. Reset BP/SP by SW-ON.		Note	Detection of air bubbles only when BUBBLE [ON/OFF] SW is ON.	

Condition or cause	Check and action
Air bubbles are detected: ○ Air introduction from blood line	Action: Correct air introduction site and collect air in drip chamber.
Air bubbles are not detected: ○ Accidental dislodgment of blood line from air bubble sensor ○ Malfunction due to external noise ○ Disconnection of the sensor harness ○ Failure of air bubble sensor board	Action: Return the line to the sensor. Action: Remove source of noise or place the source at distant position. Action: Replace the sensor with new one. Action: Adjust the board. If adjustment is impossible, replace the board or the sensor with new one.

E25 HEATER OVER HEAT ALARM

Detection	Functioning of [LS1] placed at outer wall of heater				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Preparation •Prep. End •Dialysis •Recovery	○ (Continuous)	Heater OFF	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Heater overheat ○ Failure of [T1], [T2] or [T3] ○ Failure of sub-board ○ Failure of SSR ○ Failure of RELAY	<p>Check: Measure voltage at heater side of TERMINAL BLOCK 3 using tester to find whether 230V (or 110V) is not applied always.</p> <p style="text-align: center;">⚠ CAUTION</p> <p><u>DO NOT touch heater with hand because heater temperature exceeds 60°C when heater SW is ON.</u></p> <p>Action: See E02 T1 TEMPERATURE ERROR; E03 T2 TEMPERATURE ERROR; and E04 T3 TEMPERATURE ERROR</p> <p>Action: Replace sub-board with new one.</p> <p>Action: Replace SSR with new one.</p> <p>Action: Replace RELAY with new one.</p>
No heater overheat ○ Failure of [LS1] ○ Disconnection of lead wires connected to [LS1] or incorrect connector connection ○ Failure of main board	<p>Check: Check closed contact of [LS1] at room temperature. If possible, check opening of the contact when heater surface temperature is about 60°C.</p> <p>Action: Replace [LS1] with new one if malfunction of [LS1] is detected.</p> <p>Check: Pull lightly lead wires of CN19 and those of connector of heater side to check disconnection or loosening of the wires.</p> <p>Action: Repair.</p> <p>Action, If [LS1] is normal, replace main board with new one.</p>

E26 SHORT WATER SUPPLY

Detection	Detect alarm depends on the time of Water priming complete. <ul style="list-style-type: none"> • Disabled detection of PS (Pressure switch) within 160 sec after starts the Preparation. • Disabled detection of PS within 60 sec after suction of B liquid on Fresh Dialysate mixing cycle. (After suction of A liquid in case of Acetate Dialysis) • Disabled detection of PS within 65 sec after complete the UF leak check procedure. 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •○ (Continuous) •↑ •↑ •↑	Line function stop	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
No water supply ○Occlusion of water supply filter ○Folding or occlusion of water supply line ○Leak from water supply system ○Water supply stop	Check: Remove water supply port to check normal water flow. Action: Clean the filter or replace it with new one. Action: Correct the trouble. Action: Repair leaked site. Action: Operate water supply equipment.

<p>Low water supply pressure</p> <p>○ Low water supply pressure</p> <p>○ Incorrect adjustment or failure of [PR]</p> <p>○ Incorrect adjustment or failure of [RV1]</p> <p>○ Incorrect adjustment or failure of [PS]</p> <p>○ Failure of main board</p> <p>○ Alarm while gaspurge</p> <p>○ On the acid rinse without an acid rinse option.</p>	<p>Check: Remove water supply port to check normal water flow.</p> <p>Action: Raise the pressure to $\geq 0.049\text{MPa}(0.5\text{kgf/cm}^2)$</p> <p>Action: Adjust [PR] correctly. If adjustment is impossible, repair [PR].</p> <p>Action: Adjust priming completion pressure correctly by operating [RV1]. If adjustment is impossible, repair [RV1].</p> <p>Action: Adjust [PS] correctly. If adjustment is impossible, repair [PS].</p> <p>Check: Remove connector from [PS] and check ON and OFF of [PS] by short circuiting and opening connector pins respectively in maintenance mode 4 screen.</p> <p>Action: If ON and OFF are not executed, replace main board with new one</p> <p>Check: Check the capacity of a RO supply against number of Dialysis machine.</p> <p>Action: Adjust the ratio of RO supply capacity to a number of Dialysis machine.</p> <p>Check: Check a mode of acid rinse option .</p> <p>Action: Set a mode of acid rinse option to off.</p>
--	--

E27 SHORT DIALYSATE FLOW

Detection	No functioning of [FS] after 85 seconds following chamber changeover				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Preparation •Prep. end •Dialysis •Recovery	○ (Continuous)	—	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
--------------------	------------------

<p>Normal dialysate flow (at set rate)</p> <p>○ Low dialysate flow rate</p> <p>○ Disconnection of lead wires connected to [FS] or incorrect connector contact</p> <p>○ Failure of [FS]</p> <p>○ Failure of main board</p>	<p>○ Adjust the flow rate to 300 to 600ml/min.</p> <p>Check: Pull lightly lead wires of CN15 placed on main board and those of connector of [FS] side to check disconnection or loosening of the wires.</p> <p>Action: Repair.</p> <p>Check: Measure 24V power supply and output voltages of [FS] and pin voltage of flow sensor.</p> <ul style="list-style-type: none"> • Power supply voltage between + and -: about 24V • Output voltage between OUT and -: About 24V when float is UP condition About 0V when float is DOWN condition <p>Action: Adjust positions of [FS] and [FM]. If adjustment is impossible, replace [FS] or flow meter Assy with new one.</p> <p>Check: Measure 24V power supply and output voltages of [FS] and pin voltage of flow sensor</p> <p>Action: If voltage between + and - is not about 24V, replace main board with new one.</p> <p>Check: If [FS] is normal, check indication of [FS] in maintenance mode 4 screen or FS function indicating lamp window.</p> <ul style="list-style-type: none"> • Float UP condition: OFF in maintenance mode 4 screen OFF of FS function indicating lamp window • Float DOWN condition: ON in maintenance mode 4 screen ON of FS function indicating lamp window <p>Action: Replace main board with new one.</p>
E27 SHORT DIALYSATE FLOW	

Condition or cause	Check and action
<p>No dialysate flow</p> <p>○ Occlusion in closed circuit</p> <p>○ Failure of [P2]</p> <p>○ Failure of one or more of [V5] to [V9]</p>	<p>Check: Check no dialysate flow by [FM] placed at front side of the machine.</p> <p>Action: Repair occluded component or replace it with new one.</p> <p>Action: Repair [P2] or replace it with new one.</p> <p>Action: Repair failed valve or replace it with new one.</p>

E28 DIALYSATE FLOW STOP					
Detection	<ul style="list-style-type: none"> •Functioning of [FS] within 20 seconds after chamber changeover •Functioning of [FS] for 170 seconds in DIALYSIS process excluding SQHD 				
Detection process	Audible Alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Rinse (water rinse) •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •○ (Continuous) •↑ •↑ •↑ 	Line function stop	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
--------------------	------------------

<p>Normal dialysate flow</p> <p>○ High dialysate flow rate</p> <p>○ Low dialysate flow rate</p> <p>○ This alarm occur just after gas purge function.</p> <p>○ Disconnection of lead wires connected to [FS] or incorrect connector contact</p> <p>○ Failure of [FS]</p> <p>○ Failure of main board</p>	<p>Action: Adjust the flow rate to 300 to 600mL/min.</p> <p>Check: Float of Flow meter cover the FS (Flow sensor) or not by checking maintenance mode 4.</p> <p>Action: Adjust the flow rate to 300 to 600mL/min.</p> <p>Check: Detect the FS by decreasing dialysate flow to less than 300 mL/min by gas purge function under low dialysate flow condition.</p> <p>Action: Adjust the flow rate to 300 to 600mL/min.</p> <p>Check: Pull lightly lead wires of CN15 placed on main board and those of connector of [FS] side to check disconnection or loosening of the wires.</p> <p>Action: Repair.</p> <p>Check: Measure 24V power supply and output voltages of [FS] and pin voltage of flow sensor.</p> <ul style="list-style-type: none"> • Power supply voltage between + and -: about 24V • Output voltage between OUT and -: <ul style="list-style-type: none"> About 24V when float is UP condition About 0V when float is DOWN condition <p>Action: Replace [FS] or flow meter Assy.</p> <p>Check: Measure 24V power supply and output voltages of [FS] and pin voltage of flow sensor</p> <p>Action: If voltage between + and - is not about 24V, replace main board with new one.</p> <p>Check: If [FS] is normal, check indication of [FS] in maintenance mode 4 screen or FS function indicating lamp window.</p> <ul style="list-style-type: none"> • Float UP condition: OFF in maintenance mode 4 screen OFF of FS function indicating lamp window • Float DOWN condition: ON in maintenance mode 4 screen ON of FS function indicating lamp window <p>Action: Replace main board with new one.</p>
---	---

E28 DIALYSATE FLOW STOP

Condition or cause	Check and action
Connector is connected to Dialyzer on preparation process then Dialysate flow is auto-bypassed.	Action: Connect the connector to the Dialysate on preparation end process.
After connector connected to Dialyzer ,Turn BP on and unbypass.	Action: BP should be turned off while priming a Dialysis.
Too much of flow at [V10]	Check: There is [R3] at [V10]out line tube or not. Action: Put [R3] and change the ROM later than Ver.3 Then input D.PRESS offsetting.
Wrong dialyzer in use	Action: Change to the correct dialyzer.
No dialysate flow	Check: Check no dialysate flow by [FM] placed at front side of the machine.
<input type="radio"/> Reduced chamber capacity	Check: Measure chamber capacity. Action: If the capacity is not 345 to 375mL, add silicone oil or replace the chamber Assy with new one.
<input type="radio"/> Inclusion of large amount of air in closed circuit	Check: Check air inclusion in the line visually. Action: Repair air included component or replace it with new one.
<input type="radio"/> Low priming flow rate	Action: Adjust the flow rate to 700 to 750mL/min.
<input type="radio"/> Incorrect adjustment of [PS]	Action: Adjust [PS]. If adjustment is impossible, replace [PS] with new one.
<input type="radio"/> Occlusion in closed circuit	Action: Repair occluded component or replace it with new one.
<input type="radio"/> Failure of [P2]	Action: Repair [P2] or replace it with new one.
<input type="radio"/> Failure of one or more of [V5] to [V9]	Action: Repair failed valve or replace it with new one.

E29 SN SWITCH TIME ALARM

Detection	No arrival at SN switching pressure after set alerting time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	UF stop	○	○	○
Reset	Press [RESET] SW.		Note	Enabled alerting only when SN SW is ON	

Condition or cause	Check and action
Slow venous pressure rise ○ Large gap between rollers and housing of blood pump ○ Folding or fatigue of rolling tube ○ Large clamp gap ○ Incorrect connection of [CN18] clamp connector	Check: Check whether the pressure rise is normal by operating blood pump. Action: Adjust the gap so that rolling tube is occluded properly. Action: Replace blood line including rolling tube Action: Adjust clamp gap. Check: Remove connector to check corrosion and disconnection. Action: Replace harness with new one.
No functioning of clamp	Check: See E32 CLV POSITION ERROR. Action: See E32 CLV POSITION ERROR.
Normal venous pressure rise/lowering ○ Incorrect blood flow rate, SN operation range or set SN alerting time	Check: Observe venous pressure change during SN operation. Check: Measure SN switching time. Action: Set properly.

E30 BLOOD PUMP ERROR					
Detection	<ul style="list-style-type: none"> •Overload •Abnormal rotating direction •Abnormal rpm 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> •Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery 	<ul style="list-style-type: none"> •○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑ 	—	○ SW-OFF	—	—
Reset	Press [RESET] SW. Turn on [BP/SP] SW.		Note		

Condition or cause	Check and action
Motor shaft lock ○ Catch of foreign matters by rotor of blood pump ○ Damage of motor gear head ○ Failure of sub-board	Action: Turn off power to remove foreign matters. Action: Replace blood pump Assy with new one. Action: Replace sub-board with new one.
No motor rotation ○ Disconnection of lead wires connected to [BP] or incorrect connector contact ○ failure of blood pump motor ○ Failure of sub-board	Check: Pull lightly lead wires of CN415 placed on sub-board and those of connector of [BP] side to check disconnection and loosening of the wires. Action: Connect or tighten securely. Action: Replace blood pump Assy with new one. Action: Replace sub-board with new one.

E31 BLOOD PUMP STOP ERROR

Detection	<ul style="list-style-type: none"> •Opening of blood pump cover •OFF of [ON/OFF] SW of blood pump •0 setting of flow rate setting volume of blood pump. 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	UF stop	SP stop	—	—
Reset	Press [RESET] SW. Auto reset by release from error		Note		

Condition or cause	Check and action
Opening of blood pump cover <ul style="list-style-type: none"> ○Opening of the cover ○ Incorrect setting of rolling tube to blood pump ○ Impossible closing of blood pump cover due to deformation of hinge shaft. ○No press of SW by pin of blood pump cover due to loose holding of hinge shaft 	<p>Action: Close the cover.</p> <p>Action: Set properly.</p> <p>Action: Correct the shaft deformation or replace the shaft with new one.</p> <p>Action: Tighten screws used to hold hinge shaft securely.</p>
Closing of blood pump cover <ul style="list-style-type: none"> ○Failure of [BPC] ○ Disconnection of lead wires connected to [BPC] or incorrect connector contact ○Failure of sub-board 	<p>Check: Disconnect CN403 from sub-board to measure resistance between harness pins.</p> <ul style="list-style-type: none"> • When the cover is closed: short circuit • When the cover is opened: open circuit <p>Check: Check whether resistance between harness pins is open condition when the cover is closed.</p> <p>Action: Adjust position of [BPC]. If adjustment is impossible, replace harness with new one.</p> <p>Check: Pull lightly lead wires of CN403 placed on sub-board to check disconnection and loosening of the wires.</p> <p>Action: Replace harness with new one.</p> <p>Action: Replace sub-board with new one.</p>

E32 CLV POSITION ERROR					
Detection	Abnormal clamp positioning				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	○ SW-OFF	—	—
Reset	Press [RESET] SW. Turn on [BP/SP] SW.		Note		

Condition or cause	Check and action
Clamped condition	Check: Check whether [CLV] is OFF in maintenance mode 4 screen and also check clamp head positioning visually.
○Failure of main board	Check: Measure output voltage between 1 and 2 of CN18 placed on main board. Action: If about 24V is outputted, replace main board with new one.
○Failure of solenoid	Action: If about 24V is not outputted, replace solenoid with new one.

<p>Non-clamped condition</p> <p>○ Incorrect clamp positioning</p> <p>○ Incorrect positioning of [LS5] or screw for ON of [LS5]</p> <p>○ Failure of [LS5]</p> <p>○ Failure of main board</p>	<p>Check: Check whether [CLV] is ON in maintenance mode 4 screen and also check clamp head positioning visually.</p> <p>Check: Check clamp position visually.</p> <p>Action: Correct the position.</p> <p>Check: Check whether [LS5] is turned on by screw for ON visually.</p> <p>Action: Adjust position of screw for ON.</p> <p>Check: Disconnect CN18 from main board to measure resistance between pins 5 and 6.</p> <p>Turn on and off [LS5] manually.</p> <ul style="list-style-type: none"> · When [LS5] is ON : short circuit · When [LS5] is OFF : open circuit <p>Action: Replace harness with new one.</p> <p>Action: Replace main board with new one.</p>
--	--

E33 VISCOUS PUMP DRIVING ERROR(VPa)

Detection	<ul style="list-style-type: none"> • Disable to change OFF to ON of PH2 (Rotation pulse signal) during a specified time (depend on UFR) at the UF circle. • Disable to change OFF to ON of PH1 (Original position signal) during 20 sec at the procedure of Original position. • Detect opposite rotation during VP_a working. 				
Detection process	Audible Alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ○ (Continuous) • ↑ • ↑ • ↑ 	Line function stop	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
--------------------	------------------

<p>Normal rotation of viscous pump with stop at error detection</p> <p>○ Failure of rotation detector of sensor circuit board</p> <p>○ Failure of original point detector of sensor circuit board</p> <p>○ Disconnection of [PHA] harness or incorrect connector contact</p> <p>○ Failure of main board</p> <p>○ Occur with E35 VP DRAWING VOLUME IS AT MAX LIM.</p>	<p>Check: Disconnect connector from viscous pump and check whether lower photo-sensor responds to manual slow rotation of position detection plate under monitoring PH2 of maintenance mode 4 screen.</p> <p>If the sensor does not respond, disconnect connector from the board and measure power supply voltage between pins 1 and 4 to find whether it is about 12V.</p> <p>Connect the connector and measure voltage between test points V and G of the board. When the voltage is about 12V, check whether Hi - Lo cycle is repeated during the pump rotation by the sensor output between test points OUT2 and G.</p> <p>Action: If voltage between test points V and G of the board is not about 12V or that between test points OUT2 and G of the board is abnormal, replace the board with new one.</p> <p>Check: Check whether upper photo-sensor responds to manual slow rotation of position detection plate under monitoring PH1 of maintenance mode 4 screen.</p> <p>If the sensor does not respond, disconnect connector from the board and measure power supply voltage between pins 1 and 4 to find whether it is about 12V.</p> <p>Connect the connector and measure voltage between test points V and G of the board. If the voltage is about 12V, check whether the sensor output between test points OUT1 and G is Lo at detection of original point and Hi at no detection of the point.</p> <p>Action: If voltage between test points V and G of the board is not about 12V or the sensor output between test points OUT1 and G of the board is abnormal, replace the board with new one.</p> <p>Check: Pull lightly lead wires of CN7 placed on main board to check disconnection and loosening of the wires.</p> <p>Action: Replace the harness with new one.</p> <p>Action: If the sensor circuit board is normal, replace main board with new one.</p> <p>Action: Refer to E35. then push the reset button two or three times.</p>
--	---

E33 VISCOUS PUMP DRIVING ERROR(VPa)

Condition or cause	Check and action
<p>Disabled motor rotation or uneven rotation</p> <p><input type="radio"/> Insufficient motor torque</p> <p><input type="radio"/> High water supply pressure</p> <p><input type="radio"/> High secondary pressure</p> <p><input type="radio"/> High priming pressure</p> <p><input type="radio"/> Disconnection of lead wired connected to [Vpa] or incorrect connector contact</p> <p><input type="radio"/> Failure of sub-board</p> <p><input type="radio"/> Failure of VPa</p> <p><input type="radio"/> Sensor Guide Plate crosses over PH1 by some causes.</p>	<p>Check: Check whether motor torque is sufficient. Press position detection plate lightly by hand during motor rotation or at motor start to find whether motor stop is caused.</p> <p>Action: Replace chamber Assy.</p> <p>Action: Set the pressure to 0.5 to 7.5MPa(kgf/cm²).</p> <p>Action: Adjust [PR]. If adjustment is impossible, repair [PR] or replace it with new one.</p> <p>Action: Adjust [RV1]. If adjustment is impossible, replace [RV1] with new one.</p> <p>Check: Pull lightly lead wires of CN409 placed on sub-board and those of connector of Vpa side to check disconnection and loosening of the wires.</p> <p>Action: Correct the disconnection or loosening.</p> <p>Action: Replace sub-board with new one.</p> <p>Action: Replace [Vpa] or chamber Assy with new one.</p> <p>Action: Disconnect connector from Viscous pump and rotate Position Detection Plate manually lower than detect position of PH1.</p>

E34 VISCOUS PUMP DRIVING ERROR(VPb)

Detection	<ul style="list-style-type: none"> • Disable to change OFF to ON of PH4 (Rotation pulse signal) during a specified time (depend on UFR) at the UF circle. • Disable to change OFF to ON of PH3 (Original position signal) during 20 sec at the procedure of Original position. • Detect opposite rotation during VPb working. 				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ○ (Continuous) • ↑ • ↑ • ↑ 	Line function stop	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
--------------------	------------------

<p>Normal rotation of viscous pump with stop at error detection</p> <ul style="list-style-type: none"> ○ Failure of rotation detector of sensor circuit board ○ Failure of original point detector of sensor circuit board ○ Disconnection of [PHS] harness or incorrect connector contact ○ Failure of main board ○ Occur with E35 VP DRAWIN VOLUME IS AT MAX LIM 	<p>Check: Disconnect connector from viscous pump and check whether lower photo-sensor responds to manual slow rotation of position detection plate under monitoring PH4 of maintenance mode 4 screen.</p> <p>If the sensor does not respond, disconnect connector from the board and measure power supply voltage between pins 1 and 4 to find whether it is about 12V.</p> <p>Connect the connector and measure voltage between test points V and G of the board. When the voltage is about 12V, check whether Hi - Lo cycle is repeated during the pump rotation by the sensor output between test points OUT2 and G.</p> <p>Action: If voltage between test points V and G of the board is not about 12V or that between test points OUT2 and G of the board is abnormal, replace the board with new one.</p> <p>Check: Check whether upper photo-sensor responds to manual slow rotation of position detection plate under monitoring PH3 of maintenance mode 4 screen.</p> <p>If the sensor does not respond, disconnect connector from the board and measure power supply voltage between pins 1 and 4 to find whether it is about 12V.</p> <p>Connect the connector and measure voltage between test points V and G of the board. If the voltage is about 12V, check whether the sensor output between test points OUT1 and G is Lo at detection of original point and Hi at no detection of the point.</p> <p>Action: If voltage between test points V and G of the board is not about 12V or the sensor output between test points OUT1 and G of the board is abnormal, replace the board with new one.</p> <p>Check: Pull lightly lead wires of CN8 placed on main board to check disconnection and loosening of the wires.</p> <p>Action: Replace the harness with new one.</p> <p>Action: If the sensor circuit board is normal, replace main board with new one.</p> <p>Action: Refer to E35. then push the reset button two or three times.</p>
---	---

E34 VISCOUS PUMP DRIVING ERROR(VPb)

Condition or cause	Check and action
<p>Disabled motor rotation or uneven rotation</p> <p>○ Insufficient motor torque</p> <p>○ High water supply pressure</p> <p>○ High secondary pressure</p> <p>○ High priming pressure</p> <p>○ Disconnection of lead wired connected to [Vpb] or incorrect connector contact</p> <p>○ Failure of sub-board</p> <p>○ Failure of VPb</p> <p>○ Sensor Guide Plate crosses over PH3 by some causes</p>	<p>Check: Check whether motor torque is sufficient. Press position detection plate lightly by hand during motor rotation or at motor start to find whether motor stop is caused.</p> <p>Action: Replace chamber Assy.</p> <p>Action: Set the pressure to 0.5 to 7.5MPa(kgf/cm²).</p> <p>Action: Adjust [PR]. If adjustment is impossible, repair [PR] or replace it with new one.</p> <p>Action: Adjust [RV1]. If adjustment is impossible, replace [RV1] with new one.</p> <p>Check: Pull lightly lead wires of CN409 placed on sub-board and those of connector of Vpb side to check disconnection and loosening of the wires.</p> <p>Action: Correct the disconnection or loosening.</p> <p>Action: Replace sub-board with new one.</p> <p>Action: Replace [Vpb] or chamber Assy with new one.</p> <p>Action: Disconnect connector from Viscous pump and rotate Position Detection Plate manually lower than detect position of PH3.</p>

E35 VP DRAWING VOLUME IS AT MAX LIM.

Detection	Out of operable range of viscous pump				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •○ (Continuous) •↑ •↑ •↑	Line function stop	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Execution of high UF at low dialysate flow rate	Check: Measure chamber switching time using stopwatch under monitoring flow meter indication. Action: ROM Ver 3.60 and older Adjust UFR or dialysate flow rate so that $[UFR (L/h)] \times [Chamber Switching Time]$ becomes < 364 . ROM Ver 3.70 and later Adjust UFR or dialysate flow rate so that $[UFR (L/h)] \times [Chamber Switching Time]$ becomes < 260 .
Significant deviation of input values of viscous pump flow rates	Check: Check whether values of [11.VPa FLOW] and [12.VPb FLOW] are correctly inputted in maintenance mode 6 screen.
Significant deviation of input values of chamber capacities	Check: Check whether values of [9.VCa VOLUME] and [10.VCb VOLUME] are correctly inputted in maintenance mode 6 screen.
No chamber switching for long time due to failure of [FS]	Check: See E27 SHORT DIALYSATE FLOW. Action: See E27 SHORT DIALYSATE FLOW.
Refer E33 VISCOUS PUMP DRIVING ERROR(VPa) E34 VISCOUS PUMP DRIVING ERROR(VPb)	Check: See E33 VISCOUS PUMP DRIVING ERROR(VPa) and E34 VISCOUS PUMP DRIVING ERROR(VPb). Action: See E33 VISCOUS PUMP DRIVING ERROR(VPa) and E34 VISCOUS PUMP DRIVING ERROR(VPb). After removing the above error causes, press [RESET] SW twice or thrice.

E36 VCa CLOSED CIRCUIT LEAK ERROR

Detection	Insufficient negative pressure in A chamber No arrival at set pressure within set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
Leak from closed circuit	Check: Check leak from closed circuit visually or by contact the circuit with finger. Action: Replace leaked component with new one.
Leak from NC side of [V5] or [V6]	Action: Repair leaked [V5] or [V6] or replace it with new one.
Leak from NO side of [V7] or [V8]	Action: Repair leaked [V5] or [V6] or replace it with new one.
Detection of leak error in both A and B chambers	Action: Identify leaked site(s). Operate the machine in RINSE or DIALYSIS PREP. process to apply negative pressure to C2 line. Shift to DIALYSIS WAIT process. Turn on and off each of valves [V5] to [V8] manually in maintenance mode 4 screen and operate viscous pump manually after removing connector from the pump. Clamp IN or OUT side of dialysate pressure sensor with forceps to find whether leak is caused from V5,7 or V6,8 side. Repeat clamp position change to find other leaked sites.
○Leak from [V3] or [V4]	Action: Repair leaked [V3] or [V4] or replace it with new one.
○Leak from coupler	Action: Replace O-ring of coupler or coupler including O-ring with new one.
○ Insufficient deaeration performance of water supply	Check: Measure dissolved oxygen % in water. Action: Improve deaeration performance of the means or change set closed circuit leak detection value properly.
○Gas purge pressure setting is wrong .	Action: After input D.PRESS offset appropriately, adjust [PC2]to make gas purge pressure 0 ~ 100 mmHg .

E37 VCa CLOSED CIRCUIT LEAK ERROR					
Detection	Negative pressure change in A chamber Pressure change out of set range during set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)		—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode2.	

Condition or cause	Check and action
See E36 VCa CLOSED CIRCUIT LEAK ERROR	See E36 Vca CLOSED CIRCUIT LEAK ERROR.

E38 VCb CLOSED CIRCUIT LEAK ERROR

Detection	Insufficient negative pressure in B chamber No arrival at set pressure within set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
Leak from closed circuit	Check: Check leak from closed circuit visually or by contact the circuit with finger. Action: Replace leaked component with new one.
Leak from NC side of [V7] or [V8]	Action: Repair leaked [V7] or [V8] or replace it with new one.
Leak from NO side of [V5] or [V6]	Action: Repair leaked [V5] or [V6] or replace it with new one.
Detection of leak error in both A and B chambers	Action: Identify leaked site(s). Operate the machine in RINSE or DIALYSIS PREP. process to apply negative pressure to C2 line. Shift to DIALYSIS WAIT process. Turn on and off each of valves [V5] to [V8] manually in maintenance mode 4 screen and operate viscous pump manually after removing connector from the pump. Clamp IN or OUT side of dialysate pressure sensor with forceps to find whether leak is caused from V5,7 or V6,8 side. Repeat clamp position change to find other leaked sites.
○Leak from [V3] or [V4]	Action: Repair leaked [V3] or [V4] or replace it with new one.
○Leak from coupler	Action: Replace O-ring of coupler or coupler including O-ring with new one.
○Insufficient deaeration performance of water supply	Check: Measure dissolved oxygen % in water. Action: Improve deaeration performance of the means or change set closed circuit leak detection value properly.

E39 VCb CLOSED CIRCUIT LEAK ERROR					
Detection	Negative pressure change in B chamber Pressure change out of set range during set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
See E38 VCb CLOSED CIRCUIT LEAK ERROR.	See E38 VCb CLOSED CIRCUIT LEAK ERROR.

E40 VCa CLOSED CIRCUIT LEAK ERROR					
Detection	Insufficient positive pressure in A chamber No arrival at set pressure within set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
Occlusion of water supply filter See E36 VCa CLOSED CIRCUIT LEAK ERROR.	Action: Wash up filter or change to a new filter See E36 VCa CLOSED CIRCUIT LEAK ERROR.

E41 VCa CLOSED CIRCUIT LEAK ERROR					
Detection	Positive pressure change in A chamber Pressure change out of set range during set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
See E36 VCa CLOSED CIRCUIT LEAK ERROR.	See E36 Vca CLOSED CIRCUIT LEAK ERROR.

E42 VCb CLOSED CIRCUIT LEAK ERROR					
Detection	Insufficient positive pressure in B chamber No arrival at set pressure within set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
Occlusion of water supply filter See E38 VCb CLOSED CIRCUIT LEAK ERROR.	Action: Wash up filter or change to a new filter See E38 VCb CLOSED CIRCUIT LEAK ERROR.

E43 VCb CLOSED CIRCUIT LEAK ERROR					
Detection		Positive pressure change in B chamber Pressure change out of set range during set time			
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse (Water rinse)	○ (7 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from this error in maintenance mode 2 .	

Condition or cause	Check and action
See E38 VCb CLOSED CIRCUIT LEAK ERROR.	See E38 VCb CLOSED CIRCUIT LEAK ERROR.

E44 AIR ELIMINATION TANK ERROR

Detection	Functioning of air elimination tank [FSW] for ≥ 15 seconds in DIALYSIS process or for ≥ 30 seconds in any process other than DIALYSIS process				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ○ (Continuous) • ↑ • ↑ • ↑ 	<ul style="list-style-type: none"> • — • ↑ • ↑ • Line function stop • — 	—	—	○
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Lowered float positioning	Check: Check function of [FSW] in maintenance mode 4 screen. When ○ is indicated, float is lowered position.
○ Increased float density due to water absorption	Check: Measure float weight. Float weight exceeding 8g is unacceptable. Action: Replace float with new one.
○ Lowered float positioning due to attachment of foreign matters to float and float shaft	Action: Clean float and float shaft to remove foreign matters.
○ No air discharge due to failure of [V4]	Action: Repair [V4] or replace it with new one.
○ Occlusion of [R2] by foreign matters	Check: Measure gas purge pressure which shall be 0 to 100mmHg. Action: Clean [R2] by rinsing or replace it with new one.
○ No air discharge due to failure of electric circuit to drive [V4]	Check: Turn on [V4] manually to check whether 24V is outputted to [V4] in maintenance mode 4 screen. Action: If 24V is not outputted, replace main board with new one.

<p>○ Wrong dialyzer in use</p> <p>Raised float positioning</p> <p>○ Failure of float SW due to short circuit</p> <p>○ Failure of input circuit of float SW due to short circuit</p>	<p>Change to the correct dialyzer.</p> <p>Check: Check function of [FSW] in maintenance mode 4 screen. Failure of electric system is indicated by continuous indication of ○.</p> <p>Check: Measure resistance of float side after removing connector from float.</p> <p>Action: If short circuited when float is raised position, replace float SW with new one.</p> <p>Action: If float SW is normal, replace main board with new one.</p>
---	--

E45 BLOWN-OUT FUSE DETECTION

Detection	Blown out of pump fuse				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
• Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery	• ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑	Line function stop	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
No functioning of [P1] or [P2] ○ Blown out of pump fuse ○ Incorrect contact of pump driving line ○ Failure of sub-board	<p>Check: Observe [FUSE 402] (P) and [FUSE 404] (P2) placed on sub-board visually or check using tester to find whether fuses are blown out.</p> <p>Action: If blown out, replace blown out fuse with new one.</p> <p>Check: Pull lightly each of harnesses connecting POWER PLUG – SWITCH GEAR – TERMINAL BLOCK2 – SUB-BOARD CN404 – SUB-BOARD CN401 – PUMP in series to check disconnection and loosening of the harnesses.</p> <p>Action: Tighten harness connection securely if loosening is detected. Replace harness with new one if incorrect contact or disconnection is detected</p> <p>Action: Replace sub-board with new one.</p>

<p>Functioning of [P1] and [P2]</p> <p><input type="radio"/> Failure of sub-board</p> <p><input type="radio"/> Failure of main board</p> <p><input type="radio"/> Effect of noise</p>	<p>Action: Replace sub-board with new one.</p> <p>Action: Replace main board with new one.</p> <p>Check: Identify noise source such as hyperthermia device, electro-cardiograph, laser knife or clean tower</p> <p>Action: Place the machine at a position distant from the source. Remove noise by any means such as use of insulating transformer.</p>
---	--

E46 WATCH-DOG ERROR

Detection	Reset of system due to system error detected by watchdog				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑ 	—	—	—	—
Reset	Press [RESET] SW.		Note	Initialization of RAM after error detection	

Condition or cause	Check and action
Reset of system by the machine when system error is detected by watchdog	
Incorrect function of software by effect of external noise.	<p>Check: Identify noise source such as hyperthermia device, electro-cardiograph, laser knife or clean tower</p> <p>Action: Place the machine at a position distant from the source. Remove noise by any means such as use of insulating transformer.</p>
Failure of main board	Action: Replace main board with new one.

E47 RAM ERROR					
Detection	RAM error				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
<ul style="list-style-type: none"> • Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery 	<ul style="list-style-type: none"> • ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑ 	—	—	—	—
Reset	Press [RESET] SW.		Note	Error detection at OFF→ON of POWER SW or return to normal after power failure Initialization of RAM after error detection	

Condition or cause	Check and action
Detection of RAM error by self diagnosis after operation start of the machine	Action: Replace main board with new one.

E48 EEPROM ERROR					
Detection	EEPROM error				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	—	—	—
Reset	Press [RESET] SW.		Note	Error detection at OFF→ON of POWER SW or return to normal after power failure	

Condition or cause	Check and action
Detection of EEPROM error by self diagnosis after operation start of the machine	Action: Replace main board with new one.

E49 EPROM ERROR					
Detection	EPROM error				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	—	—	—
Reset	Press [RESET] SW.		Note	Error detection at OFF→ON of POWER SW or return to normal after power failure Initialization of RAM after error detection	

Condition or cause	Check and action
Detection of EPROM error by self diagnosis after operation start of the machine	Action: Replace main board with new one.

E50 RAM INITIALIZE					
Detection	Initialization of RAM				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	—	—	—
Reset	Press [RESET] SW.		Note	Error detection at OFF→ON of POWER SW or return to normal after power failure Automatic shift to DIALYSIS WAIT process after error detection	

Condition or cause	Check and action
Detection of different version of ROM by check of RAM version	Normal
Detection of error as E46 WATCHDOG ERROR.	Normal
Detection of error as E47 RAM ERROR	Normal
Detection of error as E49 EPROM ERROR	Normal
Discharged memory back up battery on the Main board.	Check: Check if you have not turned on power breaker of the machine for long term (about 7 to 10 days). Turn on power breaker when you install the machine. Action: Set date, time and display of SETTING MODE 1 because it was initialized. Turn on power breaker of the machine to full charge battery for 22 hours
Failure of main board	Action: Replace main board with new one
Failure of sub-board	Action: Replace sub-board with new one.

E51 EEPROM INITIALIZE					
Detection	Initialization of EEPROM				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	—	—	—
Reset	Press [RESET] SW.		Note	Error detection at OFF→ON of POWER SW or return to normal after power failure Automatic shift to DIALYSIS WAIT process after error detection	

Condition or cause	Check and action
Detection of different version of EEPROM by check of EEPROM version	Normal.
Failure of main board	Action: Replace main board with new one.

E52 SWITCH ERROR					
Detection	ON of any SW on operation panel for 20 sec continuously at POWER SW ON				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	—	—	—
Reset	Press [RESET] SW.		Note	●ROM Ver~3.93 Error detection at OFF→ON of POWER SW or return to normal after power failure ●ROM Ver 4.01~ Error detection at anytime	

Condition or cause	Check and action
Press of any SW on operation panel	Check: Press SW's on operation panel one by one to find whether any SW remains to be pressed. Action: If any SW remains to be pressed, replace operation panel with new one.
No press of any SW on operation panel	
○Failure of operation panel	Action: Replace operation panel with new one.
○Failure of operation panel board	Action: Replace operation panel board with new one.
○Failure of main board	Action: Replace main board with new one.

E53 I/O PORT ERROR

Detection	Error of I/O port on board				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	Line function stop	○	—	—
Reset	Turn on [POWER] SW.		Note	Disabled all machine and operation panel SW functions after error detection	

Condition or cause	Check and action
Latch up of IC (71054) on main board	Check: Check whether the machine function becomes normal by turning on power breaker.
○Application of abnormal voltage to IC	Check: Check whether abnormal voltage is applied to IC accidentally during repair or maintenance of the machine. Action: Turn on power breaker.
○Effect of noise	Check: Identify noise source such as hyperthermia device, electro-cardiograph, laser knife or clean tower Action: Place the machine at a position distant from the source. Remove noise by any means such as use of insulating transformer.
Failure of main board	Check: Check whether error is alerted at each time when power breaker is turned on frequently. Action: Replace main board with new one.

E54 CONNECTION ERROR SBY0151-SBY0156

Detection	Disconnection of harness used to connect main board to operation panel board				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	• ○ (7 tones) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑	Line function stop	○	—	—
Reset	Turn on [POWER] SW.		Note	Disabled all machine and operation panel SW functions after error detection	

Condition or cause	Check and action
Disconnection of [CN2] from operation panel board and that of [CN10] from main board	Check: Check whether [CN2] and [CN10] are disconnected from respective boards. Action: If disconnected, connect correctly.
Connection of [CN2] to operation panel board and that of [CN10] to main board ○ Disconnection of [OPERATION 2] harness or incorrect connector contact .	Check: Pull lightly lead wires of [CN2] of operation panel board and those of [CN10] of main board to check disconnection and loosening of the wires. Action: Replace harness with new one.
○ Failure of operation panel board	Action: Replace operation panel board with new one.
○ Failure of main board	Action: Replace main board with new one.

E55 BP BATTERY ERROR

Detection	Voltage drop of blood pump backup battery to $\leq 8.5V$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
• Rinse wait • Rinse • Rinse end • Air purge • Dialysis wait • Preparation • Prep. end • Dialysis • Recovery	• ○ (7 times) • ↑ • ↑ • ↑ • ○ (Continuous) • ↑ • ↑ • ↑ • ↑	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from error in maintenance mode 7 screen	

Condition or cause	Check and action
Discharged battery ○ Long time battery operation of blood pump or insufficient battery charge time ○ Revere connector connection ○ Battery life ○ Failure of battery charge circuit	Check: Turn off power breaker of the machine and operate blood pump by battery to check whether the pump is operated normally. If the operation is abnormal, the cause is discharge or damage of battery. Action: Turn on power breaker of the machine to full charge battery for 48 hours. Action: Correct. Action: Replace battery with new one. Action: Replace sub-board with new one.
Charged battery	Check: Measure battery voltage using tester under loaded condition. Action: If voltage is $\geq 8.5V$, replace main board or sub-board.

E56 COUPLER SW ON ALARM(DIALYSIS)

Detection	ON of [LS2]				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	Line function stop	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Functioning of [LS2] ○ ON of [LS2] in dialysis process	Check: Operation miss Action: No
No functioning of [LS2] ○ Failure of [LS2]	Check: Turn on and off [LS2] manually and check its function in maintenance mode 4 screen. Check: Disconnect connector from [LS2] side and measure resistance between pins 1 and 2. • When [LS2] is ON: the pins shall be short circuit condition. • When [LS2] is OFF: the pins shall be open circuit condition. Action: If malfunction of [LS2] is detected, replace it with new one.
○ Disconnection of harness connected to [LS2] or incorrect connector contact	Check: Pull lightly lead wires of [CN15] placed on main board and those of connector of [LS2] side to check disconnection or loosening of the wires. Action: Replace harness with new one.
○ Failure of main board	Action: Replace main board with new one.

E57 COUPLER SW ON ALARM(PREP.)

Detection	Press of [UF DIALYSE] SW under [LS2] being ON				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Prep. end •Recovery	○ (3 times)	—	—	—	—
Reset	Press [RESET] SW. Auto reset by release from error		Note		

Condition or cause	Check and action
<p>┐ Functioning of [LS2]</p> <p>○ Press of [UF DIALYSE] SW under setting of coupler to its receiver</p>	<p>Check: Operation miss</p> <p>Action: Set coupler to dialyzer.</p>
<p>No functioning of [LS2]</p> <p>○ Failure of [LS2]</p>	<p>Check: Turn on and off [LS2] manually to check its function in maintenance mode 4 screen.</p> <p>Check: Disconnect connector from [LS2] side and measure resistance between pins 1 and 2.</p> <ul style="list-style-type: none"> • When [LS2] is ON: the pins shall be short circuit condition. • When [LS2] is OFF: the pins shall be open circuit condition. <p>Action: If malfunction of [LS2] is detected, replace it with new one.</p>
<p>○ Disconnection of harness connected to [LS2] or incorrect connector contact</p>	<p>Check: Pull lightly lead wires of [CN15] placed on main board and those of connector of [LS2] side to check disconnection or loosening of the wires.</p> <p>Action: Replace harness with new one.</p>
<p>○ Failure of main board</p>	<p>Action: Replace main board with new one.</p>

E58 COUPLER SW OFF ALARM

Detection	OFF of [LS2]				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse	○ (3 times)	Line function stop	—	—	—
Reset	Press [RESET] SW. Auto reset by release from error		Note		

Condition or cause	Check and action
No setting of coupler to its receiver ○ Shift to RINSE process without setting coupler to its receiver or dislodgment of coupler from its receiver during RINSE process	Check: Operation miss Action: Set coupler correctly.
Setting of coupler to its receiver ○ Damage of buttons of [LS2] ○ Deviated [LS2] position due to loosening of set screws ○ Deviated coupler receiver position due to loosening of set screws ○ Failure of [LS2] ○ Disconnection of lead wires connected to [LS2] or incorrect connector contact ○ Failure of main board	Check: Turn on and off [LS2] manually and check its function in maintenance mode 4 screen. Action: Replace [LS2] with new one. Action: Correct the position and set [LS2] securely. Action: Correct the position and set the receiver securely. Check: Disconnect connector from [LS2] side and measure resistance between pins 1 and 2. · When [LS2] is ON: the pins shall be short circuit condition. · When [LS2] is OFF: the pins shall be open circuit condition. Action: If malfunction of [LS2] is detected, replace it with new one. Check: Pull lightly lead wires of [CN15] placed on main board and those of connector of [LS2] side to check disconnection or loosening of the wires. Action: Replace harness with new one. Action: Replace main board with new one.

E59 UF MAX ALARM

Detection	Non-conformance with requirement of $[UFR] \leq [Blood\ Flow\ Rate] \times 60 \times [Set\ Upper\ UF\ Alarm\ Limit]$				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	UF stop	—	—	—
Reset	Press [RESET] SW. Auto reset by release from error		Note		

Condition or cause	Check and action
Non-conformance with requirement of $[UFR] \leq [Blood\ Flow\ Rate] \times 60 \times [Set\ Upper\ UF\ Alarm\ Limit]$	
○ Decrease of blood pump flow rate	Action: Increase blood pump flow rate.
○ Increase of UFR	Action: Decrease UFR.
[UF Alarm Limit] improperly set in maintenance mode 2 screen	Action: Set [UF Alarm Limit] properly.

E60 SYRINGE INFUSION COMPLETION					
Detection	Output of syringe pump overload signal [LS3] during ON of SP SW				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Rinse wait •Rinse •Rinse end •Air purge •Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	•○ (7 times) •↑ •↑ •↑ •○ (Continuous) •↑ •↑ •↑ •↑	—	SP stop	—	—
Reset	Press [RESET] SW. Auto reset by release from alerting		Note	Alerting during ON of SP SW	

Condition or cause	Check and action
Completion of fluid drug infusion	Check: Observe no fluid drug in syringe visually. Action: Turn off [ON/OFF] SW of SP.
Catch of syringe by any object	Action: Turn off power to release the catch.

<p>Continuous indication of ○ (overload) at [LS3] in maintenance mode 4 screen</p>	
<p>○ Incorrect positioning of [LS3] or screw for ON of it</p>	<p>Check: Observe whether screw for ON is positioned to allow ON of [LS3] visually.</p>
	<p>Action: Adjust the screw position correctly.</p>
<p>○ failure of [LS3]</p>	<p>Check: Disconnect [CN16] from main board and measure resistance between pins 1 and 3.</p>
	<p>Turn on and off [LS3] manually.</p>
	<ul style="list-style-type: none"> • When [LS3] is ON (normal): the pins shall be short circuit condition. • When [LS3] is OFF (overload): the pins shall be open circuit condition.
	<p>Action: Replace harness with new one.</p>
<p>○ Disconnection of [LS3] harness or incorrect connector contact</p>	<p>Check: Pull lightly lead wires of [CN16] placed on main board to check disconnection and loosening of the wires.</p>
	<p>Action: Correct the disconnection and/or loosening.</p>
<p>○ Failure of main board</p>	<p>Check: Monitor indication of [LS3] in maintenance mode 4 screen.</p>
	<ul style="list-style-type: none"> • Live condition : no indication • Overload condition: indication of ○
	<p>Action: Replace main board with new one.</p>

E61 COMPLETION OF SYRINGE.P INFUSION

Detection	Arrival at set syringe infusion time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (3 times)	—	SP stop	—	—
Reset	Press [RESET] SW. Auto reset by release from alerting		Note		

Condition or cause	Check and action
Arrival of set syringe infusion time	<p>Check: Not required (Normal)</p> <p>Action: To re-infuse, set infusion rate by adjusting syringe pump flow rate setting SW. If not re-infused, turn off [ON/OFF] SW of SP.</p>
Setting of syringe stop time to left dialysis time in maintenance mode 2 screen	<p>Check: Not required(Normal)</p> <p>Action: To re-infuse, set infusion rate by adjusting syringe pump flow rate setting SW. If not re-infused, turn off [ON/OFF] SW of SP.</p>
Change of UFR and UF goal and setting of dialysis left time to syringe stop time	<p>Check: Not required(Normal)</p> <p>Action: To re-infuse, set infusion rate by adjusting syringe pump flow rate setting SW. If not re-infused, turn off [ON/OFF] SW of SP.</p>

E62 PERFORMED PRIMING					
Detection	Completion of priming of target volume				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Priming	○ (3 times)	—	BP stop SW-OFF	○	—
Reset	Press [RESET] SW. Auto reset by release from alerting		Note		

Condition or cause	Check and action
Completion of priming	Check: Not required Action: To re-prime, set target priming volume by operating [F·1] SW. If not to re-prime, terminate priming function by operating [F·4] SW.
Reduction of target priming volume to currently primed volume	Check: Not required Action: To re-prime, set target priming volume by operating [F·1] SW. If not to re-prime, terminate priming function by operating [F·4] SW.

E63 PERFORMED TARGET UF.

Detection	Completion of UF of target volume				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (Continuous)	UF stop	—	—	—
Reset	Press [RESET] SW. Auto reset by release from alerting		Note	Becoming error free of set upper alarm limits of venous and dialysate pressures by alerting	

Condition or cause	Check and action
Completion of UF of target volume	Check: Not required. Action: To perform UF again, increase target UF volume. If not to perform UF, terminate dialysis.
Setting of target UF volume lower than total UF volume	Check: Not required. Action: To perform UF again, increase target UF volume. If not to perform UF, terminate dialysis.

E64 NO WATER RINSING

Detection	Shift to DIALYSIS WAIT process without performing water rinsing for a set time after acid rinsing, disinfection, hotwater rinsing, and hotdisinfection acid rinsing				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis wait	○ (3 times)	—	—	—	—
Reset	Press [RESET] SW. Reset by performing water rinsing for a set time		Note	Hotwater rinsing, and hot disinfection are not available with ROM Ver.2.40 and older	

Condition or cause	Check and action
No completion of water rinsing after acid rinsing, disinfection, hotwater rinsing, and hotdisinfection.	Action: Perform water rinsing for a set time.

E65 PREPARE DIALYSATE CONC.

Detection	No connection to concentrate tank due to no detection of pressure rise in A or B concentrate line by [PS]				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Preparation	○ (Continuous)	Line function stop	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
No dislodgment of concentrate connector from rinse port during PREPARATION process	Action: Check the both A and B lines if connect to correct port or not
No rise of priming end pressure	Dislodge concentrate connector from rinse port.
○ No water supply	Action: Supply water.
○ Incorrect adjustment or failure of [PR]	Action: Adjust [PR]. If adjustment is impossible, repair [PR] or replace it with new one.
○ Incorrect adjustment or failure of [RV1]	Action: Adjust priming end pressure by operating [RV1]. If adjustment is impossible, repair [RV1] or replace it with new one.
○ Incorrect adjustment or failure of [RV3]	Action: Repair [RV3] or replace it with new one.
○ Failure of [P1]	Action: Repair [P1] or replace it with new one.
○ Failure of [C2]	Action: Replace [C2] with new one.
○ Leak from rinse port	Check: Observe visually whether leak is caused from rinse port. Action: Replace rinse port with new one.

<p>Rise of priming end pressure</p> <p><input type="radio"/> Incorrect adjustment or failure of [PS]</p> <p><input type="radio"/> Failure of main board</p>	<p>Check: Measure priming end pressure.</p> <p>Action: Adjust [PS].</p> <p>If adjustment is impossible, repair [PS] or replace it with new one.</p> <p>Check: Disconnect connector from [PS] and check ON and OFF of [PS] in maintenance mode 4 screen by short circuiting and opening the connector pins respectively.</p> <p>Action: If ON/OFF function of [PS] is incorrect, replace main board with new one.</p>
--	--

E66 CONC CONNECTOR TO RINSING PORT.

Detection	No connection to rinse port due to detection of pressure rise in A or B concentrate line by [PS]				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Rinse	○ (3 times)	Line function stop	—	—	—
Reset	Press [RESET] SW.		Note		

Condition or cause	Check and action
Acetate dialysis	Check: Check if the undiluted solution is correct or not Action: Connect line A to the undiluted solution. Connect line B to the B rinse port
No connection of concentrate connector to rinse port while on preparation process.	Action: Connect concentrate connector to rinse port.

<p>Connection of concentrate connector to rinse port</p> <p>○ Incorrect adjustment or failure of [PR]</p> <p>○ Incorrect adjustment or failure of [RV1]</p> <p>○ Incorrect adjustment or failure of [RV3]</p> <p>○ Incorrect adjustment or failure of [PS]</p> <p>○ Failure of main board</p> <p>○ Failure of [V11], [V13], [V14], [V5], [V6], [V7] and/or [V8]. (Failure of [V15] when Na option is used.)</p> <p>○ Failure of main board</p>	<p>Action: Adjust [PR]. If adjustment is impossible, repair [PR] or replace it with new one.</p> <p>Action: Adjust priming end pressure by operating [RV1]. If adjustment is impossible, repair [RV1] or replace it with new one.</p> <p>Action: Repair [RV3] or replace it with new one.</p> <p>Action: Adjust [PS]. If adjustment is impossible, repair [PS] or replace it with new one.</p> <p>Check: Disconnect connector from [PS] and check ON and OFF of [PS] in maintenance mode 4 screen by short circuiting and opening the connector pins respectively.</p> <p>Action: If ON/OFF function of [PS] is incorrect, replace main board with new one.</p> <p>Check: Check ON and OFF of each valve in maintenance mode 4 screen. If ON/OFF function of each valve is incorrect, check whether output voltage is 24V.</p> <p>Action: If 24V is outputted, replace relevant valve with new one. If 24V is not outputted, replace main board with new one.</p>
---	---

E67 RESET SW					
Detection	No press of [RESET] SW for a set time after alerting				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
•Dialysis wait •Preparation •Prep. end •Dialysis •Recovery	○ (Continuous)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from alerting in maintenance mode 2 screen	

Condition or cause	Check and action
No press of [RESET] SW	Check: Not required. Action: Press [RESET] SW.

E68 UF DIALYSE SW					
Detection	No press of [UF DIALYSE] SW for a set time after alerting				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Prep. end	○ (3 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from alerting in maintenance mode 2 screen	

Condition or cause	Check and action
No press of [UF DIALYSE] SW	Check: Not required Action: Press [UF DIALYSE] SW to shift to DIALYSIS process.

E69 BYPASS SW					
Detection	No release from manual ON condition of [BYPASS] SW for a set time				
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Dialysis	○ (3 times)	—	—	—	○
Reset	Press [RESET] SW.		Note	Enabled release from alerting in maintenance mode 2 screen	

Condition or cause	Check and action
No press of [BYPASS] SW	Check: Not required. Action: Press [BYPASS] SW to release from bypass condition.

NO INDICATION IN LCD					
Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Simultaneous alerting of this error with E54 CONNECTION ERROR SBY0151-SBY0156	Check: See E54 CONNECTION ERROR SBY0151-SBY0156. Action: See E54 CONNECTION ERROR SBY0151-SBY0156.
Screen disappearance during RINSE WAIT process ○ Functioning of auto OFF mode	Check: Not required. Note: Back-light is OFF after no SW operation for ≥ 30 minutes during “indication of rinse execution screen and Pre-Rinse Wait process” or “indication of rinse end screen”. Action: Press any SW on operation panel to turn on back-light.

<p>OFF of LCD</p> <p>○ Incorrect contrast adjustment</p> <p>○ Disconnection of [OPERATION PART 1] harness or incorrect connector contact</p> <p>○ Disconnection of LCD harness or incorrect connector contact</p> <p>○ Failure of LCD</p> <p>○ Failure of main board</p>	<p>Action: Adjust contrast by turning volume placed on operation panel board.</p> <p>Check: Pull lightly lead wires of [CN1] of operation panel board and [CN10] of main board to check disconnection and loosening of the wires.</p> <p>Action: Replace harness with new one followed by turning on power SW.</p> <p>Check: Pull lightly lead wires of [CN4] and [CN5] of operation panel board to check disconnection and loosening of the wires.</p> <p>Action: Replace lead wires (harness) with new ones.</p> <p>Action: Replace LCD with new one.</p> <p>Action: Replace main board with new one.</p>
--	--

NO INDICATION IN LCD

Condition or cause	Check and action
<p>Deviated indication</p> <p>○ Indication deviated by external noise</p> <p>○ Disconnection of [OPERATION PART 1] harness or incorrect connector contact</p> <p>○ Disconnection of LCD harness or incorrect connector contact</p> <p>○ Failure of LCD</p> <p>○ Failure of main board</p>	<p>Check: Recover by turn on power SW again.</p> <p>Identify noise source such as hyperthermia device, electro-cardiograph, laser knife or clean tower.</p> <p>Action: Move the source to a position distant from the machine.</p> <p>Use noise reducing means such as insulating transformer.</p> <p>Check: Pull lightly lead wires of [CN1] of operation panel board and [CN10] of main board to check disconnection and loosening of the wires.</p> <p>Action: Replace harness with new one followed by turning on power SW.</p> <p>Check: Pull lightly lead wires of [CN4] and [CN5] of operation panel board to check disconnection and loosening of the wires.</p> <p>Action: Replace lead wires (harness) with new ones.</p> <p>Action: Replace LCD with new one.</p> <p>Action: Replace main board with new one.</p>

DISABLED SW FUNCTION					
Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Simultaneous alerting of this error with E54 CONNECTION ERROR SBY0151-SBY0156	Check: See E54 CONNECTION ERROR SBY0151-SBY0156. Action: See E54 CONNECTION ERROR SBY0151-SBY0156.
Disabled SW function; or Enabled SW function with error (function of any SW different from that of pressed SW)	
○ Disconnection of operation panel sheet cable or incorrect connector contact	Check: Pull lightly flexible cable of [CN3] placed on operation panel board to check disconnection and loosening of the cable. Action: Replace flexible cable with new one.
○ Disconnection of [OPERATION PART 1] harness or incorrect connector contact	Check: Pull lightly lead wires of [CN1] of operation panel board and [CN10] of main board to check disconnection and loosening of the wires. Action: Replace harness with new one.
○ Failure of operation panel sheet	Action: Replace operation panel sheet with new one.
○ Failure of operation panel board	Action: Replace operation panel board with new one.
○ Failure of main board	Action: Replace main board with new one.

DISABLED SYRINGE PUMP OPERATION

Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Disabled motor operation <ul style="list-style-type: none"> <input type="radio"/> No selection of SYRINGE <input type="radio"/> Disconnection of [SP] harness or incorrect connector contact <input type="radio"/> Failure of sub-board <input type="radio"/> Failure of motor 	<p>Action: Select ON of [18.SYRINGE OPTION] in maintenance mode 7 screen.</p> <p>Check: Pull lightly lead wires of [CN410] placed on sub-board and those of connector of [SP] side to check disconnection and loosening of the wires.</p> <p>Action: Replace.</p> <p>Action: Replace.</p> <p>Action: Replace.</p>
Enabled motor operation <ul style="list-style-type: none"> <input type="radio"/> Slip of half nut <input type="radio"/> Idle run of pulley <input type="radio"/> Breakage of belt <input type="radio"/> Loosening of belt 	<p>Action: Replace.</p> <p>Check: Check whether pulley is loosely fixed to shaft.</p> <p>Action: Tighten.</p> <p>Action: Replace.</p> <p>Action: Tighten.</p>

BLOOD PUMP ERROR					
Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Uneven rotation <input type="radio"/> Improper gap between rotors and housing	Action: Adjust the gap.
Restricted top rotation speed. <input type="radio"/> Failure of sub-board or motor	Check: See E30 BLOOD PUMP ERROR. Action: See E30 BLOOD PUMP ERROR.
Higher or lower rotation speed <input type="radio"/> Incorrect BP flow rate adjustment <input type="radio"/> Failure of sub-board or motor	Action: Adjust BP flow rate correctly. Check: See E30 BLOOD PUMP ERROR. Action: See E30 BLOOD PUMP ERROR.
Heavy manual rotation <input type="radio"/> Disabled OFF of BP cover SW <input type="radio"/> Failure of sub-board or motor	Check: See E31 BLOOD PUMP STOP ERROR. Action: See E31 BLOOD PUMP STOP ERROR. Check: See E30 BLOOD PUMP ERROR. Action: See E30 BLOOD PUMP ERROR.

Noise <input type="radio"/> Wearing of gear head <input type="radio"/> Wearing of rotor bearing	 Action: Replace blood pump Assy. Action: Replace the bearing.
--	--

CHANGE OF LEVEL IN VENOUS DRIP CHAMBER

Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Rise of level in chamber <ul style="list-style-type: none"> ○ Leak from a site between venous drip chamber and pressure inlet port, in particular, from filter junction ○ Leak from a site between pressure inlet port and pressure sensor ○ Failure of pressure sensor ○ Folding or occlusion of venous line downstream of venous drip chamber 	<p>Check: Clamp line between the chamber and port by forceps to find leaked site.</p> <p>Action: Remove the leak cause.</p> <p>Check: Clamp line between the port and sensor by forceps to find leaked site.</p> <p>Action: Remove the leak cause.</p> <p>Action: Replace venous pressure sensor board with new one.</p> <p>Action: Correct.</p>
Lowering of level in chamber <ul style="list-style-type: none"> ○ Air inclusion due to insufficient blood supply ○ Air inclusion through arterial line ○ Air introduction from dialyzer 	<p>Action: Correct.</p> <p>Action: Correct.</p> <p>Action: Remove the air introduction cause.</p>

ABNORMAL DIALYSATE DISCHARGE BY PRESS OF DRAIN SW					
Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
Slow discharge flow rate <ul style="list-style-type: none"> ○ Deformation of diaphragm of electro-magnetic valve [V4] or catch of plunger of it by any object ○ Occlusion of orifice [R2] ○ Performance reduction of dialysate pump [P2] ○ Folding of tube downstream of air elimination tank [AS2] ○ High gas purge pressure 	<p>Action: Repair [V4] or replace it with new one.</p> <p>Action: Repair [R2] or replace it with new one.</p> <p>Action: Repair [P2] or replace it with new one.</p> <p>Action: Correct folding.</p> <p>Action: Adjust the gas purge pressure to 0 ~ 100mmHg.</p>

<p>No discharge</p> <p>○ Failure of electromagnetic valve [V4]</p> <p>○ Disconnection of [VALVE 2] harness or incorrect connector contact</p> <p>○ Failure of main board</p> <p>○ Failure of DRAIN SW</p>	<p>Action: Repair [V4] or replace it with new one.</p> <p>Check: Pull lightly lead wires of [CN12] placed on main board and those of connector of [V4] side to check disconnection and loosening of the wires.</p> <p>Action: Replace harness (lead wires) with new ones.</p> <p>Check: Turn on [V4] manually to check whether 24V is outputted to [V4] in maintenance mode 4 screen.</p> <p>Action: If 24V is not outputted, replace main board with new one.</p> <p>Action: Replace operation panel sheet or operation panel board with new one.</p>
--	--

UFR ERROR					
Detection					
Detection process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset				Note	

Condition or cause	Check and action
Leak from closed circuit	<p>Check: See E36~E43 VC CLOSED CIRCUIT ERROR.</p> <p>Action: See E36~E43 VC CLOSED CIRCUIT ERROR.</p>
Incorrect viscous pump rotation	<p>Check: See E33~E34 VISCOUS PUMP DRIVING ERROR.</p> <p>Action: See E33~E34 VISCOUS PUMP DRIVING ERROR.</p>
Incorrectly set UFR calibration coefficients: a, c and/or d	<p>Action: Set [13.UF CORR.COEFFICIENT a], [14.UF CORR.COEFFICIENT c] and/or [15.UF CORR.COEFFICIENT d] indicated in maintenance mode 6 screen to a, c and/or d shown in set value seal.</p>
Incorrect (deviated) VP flow rate	<p>Check: Check values of [11.VPa FLOW] and [12.VPb FLOW] in maintenance mode 6 screen.</p> <p>Action: Re-input normal values if VP flow rate is higher than normal (lower UFR) or lower than it (higher UFR).</p>
<p>Air inclusion from dialysate line of priming circuit .</p> <p>○Air inclusion from [AS1]</p>	<p>Check: Observe visually whether [AS1] is positioned slantwise.</p> <p>Action: If slanted, correct to horizontal placement. If not so, replace [AS1] with new one.</p>

Incorrect pressure adjustment of dialysate line	
○ Incorrect priming end pressure	
Incorrect adjustment or failure of [PR]	Action: Adjust [PR] correctly. If adjustment is impossible, repair [PR] or replace it with new one.
Incorrect adjustment or failure of [RV1]	Action: Adjust priming end pressure by operating [RV1]. If adjustment is impossible, repair [RV1] or replace it with new one.
Incorrect adjustment or failure of [RV3]	
Failure of [P1]	Action: Repair [RV3] or replace it with new one.
Failure of [C2]	Action: Repair [P1] or replace it with new one.
Leak from rinse port	Action: Replace [C2] with new one. Check: Observe visually whether leak is caused from rinse port.
○ Incorrect adjustment of pressure in closed circuit	Action: If leak is detected, replace rinse port.
Incorrect adjustment or failure of [RV2]	Action: Adjust [RV2]. If adjustment is impossible, repair [RV2] or replace it with new one.

DATE AND TIME ARE INITIALIZED					
Detection					
Detection Process	Audible alarm	Dialysate line function	BP/SP stop	CLV function	Reset
Reset			Note		

Condition or cause	Check and action
<p>Date and time are indicated to 1997 years 6 months 15 dates 12:00 time.</p> <p>○ Discharged memory back up battery on the Main board.</p>	<p>Check: Check if you have not turned on power breaker of the machine for long term (about 7 to 10 days). Turn on power breaker when you install the machine.</p> <p>Action: Set date and time. Turn on power breaker of the machine to full charge battery for 22 hours.</p>