

# **APC** ErgoTouch Pro

*Airborne Particle Counter*

*ETP-1.0 en  
January 2009*

**Operation and Service Manual  
(English)**





# **APC** ErgoTouch Pro

## *Airborne Particle Counter*

---

*ETP-1.0 en  
January 2009*

**Operation and Service Manual  
(English)**

**SHIP TO/MAIL TO:**

Biotest AG  
Technical Center Microbiology  
Wernher von Braun Straße 9  
D-69214 Eppelheim  
German

**E-mail address:**

[Hycon@biotest.de](mailto:Hycon@biotest.de)

**Website:**

<http://www.biotest.de>

**Technical Support:**

+49.6221.726.5130

**Fax:**

+49.6221.726.5139

**Local Technical Support:**



## **Manual History**

The following is a manual history of the  **APC** ErgoTouch Pro Handheld Airborne Particle Counter:

| <b>Revision</b> | <b>Date</b>   |
|-----------------|---------------|
| A               | November 2008 |
| ETP-1.0 en      | January 2009  |
|                 |               |

## Warranty

|   |   |
|---|---|
| <b>Revision</b>                             | ETP-1.0 en / January 2009   |
| <b>Copyright</b>                            | ©Biotest AG / 2009 / All rights reserved.   |
| <b>Address</b>                              | Biotest AG, Technical Center Microbiology, Wernher von Braun Straße 9, 69214 Eppelheim, Germany   |
| <b>E-mail Address</b>                       | <a href="mailto:HYCON@biotest.de">HYCON@biotest.de</a>  |
| <b>Limitation of Warranty and Liability</b> | <p>Seller warrants the goods sold hereunder, under normal use and service as described in the operator's manual, shall be free from defects in workmanship and material for (24) months, or the length of time specified in the operator's manual, from the date of shipment to the customer. This warranty period is inclusive of any statutory warranty. This limited warranty is subject to the following exclusions:</p> <ol style="list-style-type: none"><li>Parts repaired or replaced as a result of repair services are warranted to be free from defects in workmanship and material, under normal use, for 90 days from the date of shipment.</li><li>Seller does not provide any warranty on finished goods manufactured by others or on any fuses, batteries or other consumable materials. Only the original manufacturer's warranty applies.</li><li>Unless specifically authorized in a separate writing by Seller, Seller makes no warranty with respect to, and shall have no liability in connection with, goods which are incorporated into other products or equipment, or which are modified by any person other than Seller.</li></ol> <p>The foregoing is IN LIEU OF all other warranties and is subject to the LIMITATIONS stated herein. <b>NO OTHER EXPRESS OR IMPLIED WARRANTY OF FITNESS FOR PARTICULAR PURPOSE OR MERCHANTABILITY IS MADE.</b></p> <p>TO THE EXTENT PERMITTED BY LAW, THE EXCLUSIVE REMEDY OF THE USER OR BUYER, AND THE LIMIT OF SELLER'S LIABILITY FOR ANY AND ALL LOSSES, INJURIES, OR DAMAGES CONCERNING THE GOODS (INCLUDING CLAIMS BASED ON CONTRACT, NEGLIGENCE, TORT, STRICT LIABILITY OR OTHERWISE) SHALL BE THE RETURN OF GOODS TO SELLER AND THE REFUND OF THE PURCHASE PRICE, OR, AT THE OPTION OF SELLER, THE REPAIR OR REPLACEMENT OF THE GOODS. IN NO EVENT SHALL SELLER BE LIABLE FOR ANY SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES. SELLER SHALL NOT BE RESPONSIBLE FOR INSTALLATION, DISMANTLING OR REINSTALLATION COSTS OR CHARGES.</p> <p>No Action, regardless of form, may be brought against Seller more than 12 months after a cause of action has accrued. The goods returned under warranty to Seller's factory shall be at Buyer's risk of loss, and will be returned, if at all, at Seller's risk of loss.</p> <p>Buyer and all users are deemed to have accepted this LIMITATION OF WARRANTY AND LIABILITY, which contains the complete and exclusive limited warranty of Seller. This LIMITATION OF WARRANTY AND LIABILITY may not be amended, modified or its terms waived, except by writing signed by an Officer of Seller.</p> |
| <b>Service Policy</b>                       | Knowing that inoperative or defective instruments are as detrimental to Biotest as they are to our customers, our service policy is designed to give prompt attention to any problems. If any malfunction is discovered, please contact your nearest sales office or representative, or call Biotest's Customer Service department at +49.6103.801.590  |
| <b>Trademarks</b>                           | APC, ERGOTOUCH, ERGOTOUCH PRO and ERGOTRAK are trademarks of Biotest AG. Biotest and the Biotest logo are registered trademarks of Biotest AG. Microsoft and Excel are registered trademarks of Microsoft Corporation.  |




## Contents

|  |            |
|--|------------|
| <b>Manual History .....</b>  | <b>ii</b>  |
| <b>Warranty .....</b>  | <b>iii</b> |
| <b>Safety Information .....</b>  | <b>vii</b> |
| Laser Safety.....  | vii        |
| Labels .....   | viii       |
| Description of Caution / Warning Symbols .....                           | viii       |
| Caution and Warning Symbols .....  | viii       |
| Getting Help.....  | ix         |
| <b>Chapter 1 Introduction and Unpacking.....</b>                         | <b>1</b>   |
| Unpacking the APC ERGOTOUCH PRO Handheld Airborne Particle Counter ..... | 1          |
| Optional Accessories .....   |            |
| <b>Chapter 2 Getting Started .....</b>                                   | <b>4</b>   |
| Using the Instrument Stand and Stylus.....                               | 4          |
| Providing Power.....   | 5          |
| To Install the Lithium-ion Battery .....                                 | 5          |
| To Use AC Power .....  | 6          |
| Performing a Daily Zero Check.....                                       | 6          |
| To Perform a Zero Check .....  | 6          |
| Installing an Isokinetic Inlet.....                                      | 7          |
| Installing a Temperature / Relative Humidity Probe.....                  | 7          |
| <b>Chapter 3 Operation .....</b>   | <b>8</b>   |
| Screen Layout and Functionality.....                                     | 8          |
| Main Tab.....  | 8          |
| Zoomed Data Screen.....  | 11         |
| Setup Tab .....  | 11         |
| Software Input Panel (Keyboard).....                                     | 12         |
| Data Setup Screen.....   | 12         |
| Count Units Screen.....  | 12         |
| Clear Samples Screen .....   | 13         |
| System Setup Screen .....  | 13         |
| Change Power On Password Screen .....                                    | 14         |
| Change Setup Password Screen.....  | 14         |
| System Configuration Screen .....  | 15         |
| Print Setting Screen .....   | 15         |
| Print Schedule Screen .....  | 16         |
| Device Setup Screen .....  | 16         |
| Date and Time Screen .....   | 17         |
| Display Screen.....  | 17         |

|   |           |
|---|-----------|
| Diagnostics Screen .....  | 17        |
| Sampling Screen.....  | 18        |
| Channels Screen .....   | 18        |
| Sample Timing Screen.....   | 19        |
| Count Mode Screen .....   | 20        |
| Environment Screen .....  | 20        |
| Locations Screen .....  | 21        |
| Alarms Screen .....   | 22        |
| Recipe Screen .....   | 22        |
| Data Tab .....  | 23        |
| Reports Tab .....   | 25        |
| Room Definition Screen .....  | 26        |
| Generate Screen.....  | 27        |
| <b>Chapter 4 Data Handling .....</b>  | <b>28</b> |
| USB Computer Communication.....   | 28        |
| Installing Software .....   | 28        |
| Procedure .....   | 28        |
| Download Data .....   | 30        |
| Delete Data.....  | 32        |
| <b>Chapter 5 Maintenance.....</b>   | <b>33</b> |
| Maintenance Schedule .....  | 33        |
| Daily Zero Check .....  | 33        |
| Cleaning the Instrument Enclosure.....  | 33        |
| <b>Chapter 6 Troubleshooting .....</b>  | <b>34</b> |
| <b>Chapter 7 Contacting Customer Service .....</b>                                  | <b>36</b> |
| Contacting Customer Service .....   | 36        |
| Technical Contacts .....  | 36        |
| Returning the APC ERGOTOUCH PRO Handheld Airborne Particle Counter for Service..... | 36        |
| <b>Appendix A Specifications.....</b>   | <b>37</b> |
| APC ERGOTOUCH PRO Handheld Airborne Particle Counter .....                          | 37        |
| Temperature / Relative Humidity Probe Specifications (optional accessory) .....     | 38        |



## Safety Information

This section gives instructions to promote safe and proper handling of the  **APC ErgoTouch Pro** Handheld Airborne Particle Counters.

### IMPORTANT

There are no user-serviceable parts inside the instrument. Refer all repair and maintenance to a qualified factory-authorized technician. All maintenance and repair information in this manual is included for use by a qualified factory-authorized technician.

Any attempt to open or service the instrument by non-authorized agencies or personnel will void the warranty.

## Laser Safety

The APC ERGOTOUCH PRO Airborne Particle Counter is a handheld Class I laser- based instrument.

- During normal operation, the user will **not** be exposed to laser radiation.
- Precaution should be taken to avoid exposure to hazardous radiation in the form of intense, focused, visible light.
- Exposure to this light may cause blindness.

### Take these precautions:

- **DO NOT** remove any parts from the particle counter unless you are specifically told to do so in this manual.
- **DO NOT** remove the housing or covers. There are no user-serviceable components inside the housing.








### WARNING

The use of controls, adjustments, or procedures other than those specified in this manual may result in exposure to hazardous optical radiation.

## Labels

Advisory labels and identification labels are attached to the outside of the particle counter housing and to the optics housing on the inside of the instrument:

|  |  |
|--|--|
| <b>1. Serial Number Label</b><br>(back panel)                                |  <b>9306-01BT</b> <b>S/N: 93060849001</b> <br><br><b>12 VDC</b><br><br><small>COMPLIES WITH 21 CFR 1040.10 AND 1040.11<br/>         U.S. Patent 6167107<br/>         Werner von Braun StraÙe 9, D-69214 Eppelheim, Germany      HYCONService@biotest.de</small> |
| <b>2. Laser Radiation Label</b><br>(internal)                                | <div style="border: 1px solid black; padding: 10px; text-align: center;"> <p><b>DANGER!</b></p> <p>VISIBLE LASER RADIATION WHEN<br/>         OPEN. AVOID DIRECT<br/>         EXPOSURE TO BEAM</p> <p>WARNING: NO USER SERVICABLE<br/>         PARTS INSIDE. REFER SERVICING<br/>         TO QUALIFIED PERSONNEL</p> </div>   |
| <b>3. European symbol for non-disposable item.</b><br>Item must be recycled. |   |

## Description of Caution/Warning Symbols

Appropriate caution/warning statements are used throughout the manual and on the instrument that require you to take cautionary measures when working with the instrument:

### Caution



| <b>C a u t i o n</b>   |
|--|
| <p>Failure to follow the procedures prescribed in this manual might result in irreparable equipment damage. Important information about the operation and maintenance of this instrument is included in this manual.</p> |





### Warning



| <b>W A R N I N G</b>  |
|---|
| <p>Warning means that unsafe use of the instrument could result in serious injury to you or cause damage to the instrument. Follow the procedures prescribed.</p> |

## Caution or Warning Symbols

The following symbols may accompany cautions and warnings to indicate the nature and consequences of hazards:

|   |  |
|---|--|
|  | <p>Warns that uninsulated voltage within the instrument may have sufficient magnitude to cause electric shock. Therefore, it is dangerous to make contact with any part inside the instrument.</p> |
|  | <p>Warns that the instrument contains a laser and that important information about its safe operation and maintenance is included in the manual.</p>   |
|  | <p>Warns that the instrument is susceptible to electro-static dissipation (ESD) and ESD protection procedures should be followed to avoid damage.</p>  |
|  | <p>Indicates the connector is connected to earth ground and cabinet ground.</p>  |

## Getting Help

To obtain assistance with this product or to submit suggestions, please contact Customer Service:

Biotest AG  
 Technical Center Microbiology  
 Wernher von Braun Straße 9  
 69214 Eppelheim  
 Germany  
 Fax: +49.6221.726.5139  
 Telephone: +49.6221.726.5130  
 E-mail Address: [HYCON@biotest.de](mailto:HYCON@biotest.de)  
 Web site: [www.biotest.com](http://www.biotest.com)

## Chapter 1

### **Introduction and Unpacking**






The APC ERGOTOUCH PRO Airborne Particle Counter is a lightweight, handheld particle counter with a touch-screen interface. It operates on the included lithium-ion battery or AC power.




The device has a 0.1 ft<sup>3</sup> / min (2.83 L / min) flow rate and counts bin sizes from 0.3 to 10 µm. Up to 10,000 data sets can be downloaded for analysis and reporting using the ERGOTRAK Lite Data Download Software included with the device.

Typical applications for this particle counter include clean room monitoring, research, exposure assessment, indoor air quality, filter testing, clearance testing, quality assurance, and contaminant migration studies. All APC ERGOTOUCH airborne particle counters meet the JIS B 9921: 1997 standard.

### **Unpacking the APC ErgoTouch Pro Handheld Airborne Particle Counter**

Carefully unpack the APC ERGOTOUCH PRO Airborne Particle Counter from the shipping container and verify that all the items shown in the photos below and listed in the following tables are present. Contact Biotest immediately if items are missing or broken:









| Qty. | Item Description                            | Part/Model             | Reference Picture   |
|------|---|------------------------|---|
| 1    | APC ERGOTOUCH PRO Airborne Particle Counter | 942 550                |  |
| 1    | Power Cord                                  | (US)<br>(UK)<br>(EURO) |  |
| 1    | AC power adapter                            |                        |  |
| 1    | Isokinetic inlet                            |                        |  |
| 1    | Battery pack                                |                        |  |

| Qty. | Item Description   | Part/Model | Reference Picture   |
|------|--|------------|---|
| 1    | Computer cable (2 m), USB A to B                         |            |  |
| 1    | Stylus   |            |  |
| 1    | HEPA zero filter assembly                                |            |  |
| 1    | ERGOTRAK LITE data download utility CD (includes manual) |            | (not shown)   |
| 1    | Operation and Service Manual                             |            | (installed on CD)   |
| 1    | Calibration certificate                                  |            | (not shown)   |
| 1    | Quick Start Guide  |            | (installed on CD)   |

## Optional Accessories

The following photos and table list optional accessories. If you ordered optional accessories, make certain they have been received and are in working order:

### APC ERGOTOUCH PRO Airborne Particle Counter Optional Accessories

| Item Description   | Part/Model | Reference Picture   |
|--|------------|---|
| External battery charger with AC adapter and power cord    |            |    |
| External Printer   |            |    |
| Heavy Duty Carry case                                      |            |    |
| Printer paper (5 rolls)                                    |            | (not shown)   |
| Temperature/humidity probe                                 |            |  |
| Stainless Steel Isokinetic inlet                           |            |  |
| Stainless Steel Isokinetic probe (used with tubing)        |            |  |
| 0.1 cfm Barb Inlet Fitting                                 |            |  |
| Tubing, Superthane 1/8-inch ID x 1/4-inch OD, Clear 100 ft |            |  |

## Chapter 2

### Getting Started

This chapter provides information to help you use the APC ERGO TOUCH PRO Airborne Particle Counter including:

- Using the Instrument Stand and Stylus
- Providing Power
- Performing a Daily Check
- Installing an Isokinetic Inlet
- Installing a Temperature/Relative Humidity Probe

#### **Using the Instrument Stand and Stylus**

The APC ERGO TOUCH PRO is equipped with an integral instrument support stand. To open the stand, grasp it by the large finger hole and pull it out until it locks into place. Be careful not to overextend the stand. To store the stand out of the way when not in use, simply push the stand back until it snaps into place.



The APC ERGO TOUCH PRO is also equipped with a plastic stylus for use with the touch screen interface. The stylus locks into place in the case near the top of the unit when not in use.



## Providing Power

The APC ERGOTOUCH PRO may be powered using a removable rechargeable lithium-ion battery, or through an AC power cord.

### Notes:

1. When using AC power, the battery (if installed) charges when the instrument is on, but not while actively sampling.
2. Removing / changing the lithium-ion battery or disconnecting the AC power does not cause loss of data.
3. A new battery is installed as part of routine service, if necessary.

## To Install the Lithium-ion Battery

1. Remove the battery cover from the back of the instrument by lightly depressing the textured tab on the cover located on the lower left.
2. Slide the lithium-ion battery into the slot, press down lightly and slide it forward (toward the top of the unit) until it locks into place.
3. Replace the battery cover and slide it in place until you hear a click.



### WARNING

The battery supplied by Biotest has built in protection against explosion and fire hazard. Do **not** use a substitute.




### WARNING

Do **not** use non-rechargeable batteries in this instrument. Fire, explosions, or other hazards may result.



## To Use AC Power

1. Connect the AC power adapter to the power cord.
2. Insert the AC power adapter into the bottom of the APC ERGOTOUCH PRO.
3. Connect the power cord to an outlet.
4. Press the on/off button  (located in the center of the front of the instrument).
5. After a splash screen displays the Biotest logo, a brief start-up sequence begins as the Windows® CE operating system boots up.



## Performing a Daily Zero Check

A zero check should be performed at least once a day. It should also be performed before conducting any important testing or certification.

### To Perform a Zero Check

Turn on the instrument and wait until the main menu appears.

1. Remove the Isokinetic inlet if attached. The daily zero check cannot be performed when the isokinetic inlet is attached to the instrument.
2. Attach the zero filter to the inlet nozzle located on the top of the instrument.
3. Press the **Start** button and allow the instrument to purge for 2 minutes.
4. After the 2-minute purge, continue to sample. In accordance with JIS standards, there should be no more than 1 particle counted at any size in 5 minutes.

**Note:** If the instrument does not go to zero (1 particle in 5 min is considered zero), refer to Chapter 6, [Troubleshooting](#), for additional information.



5. Remove the zero filter and put the isokinetic inlet back on; the instrument is now ready for operation.

### ***Installing an Isokinetic Inlet***

The Isokinetic inlet smoothly accelerates air into the inlet of the instrument. To install, simply thread the inlet directly onto the inlet nozzle until finger tight. The inlet seals over an o-ring so it doesn't have to be very tight to seal.



### ***Installing a Temperature/Relative Humidity Probe***

To install the optional temperature/relative humidity probe:

1. Align the small red dot at the base of the probe to the corresponding red dot on the socket.
2. Press the probe into the socket until it clicks.
3. Temperature and relative humidity are automatically displayed in the upper-left corner.
4. Remove the probe by pulling straight up.




## Chapter 3

### Operation

The APC ERGOTOUCH PRO Airborne Particle Counter is controlled using a touch screen display. Use the plastic stylus or your finger tip. **DO NOT** use sharp objects (such as a pen point) that may damage the screen overlay.

To turn on the instrument, press the **on/off** button (located in the center of the front of the instrument). After a splash screen displays the Biotest logo, a brief start-up sequence begins as the Windows<sup>®</sup> CE operating system boots up.

The instrument is ready for operation when the main tab (shown right) appears. If an optional temperature/humidity probe is attached, those values will be shown in the upper-left corner also.

|               |          | 11/20/08<br>11:43:03 AM   |      |
|---------------|----------|---|------|
|               |          | Automatic  |      |
| $\mu\text{m}$ | $\Sigma$ | Location  |      |
| 0.3           | 0        | ▼   |      |
| 0.5           | 0        | Delay: 00:00:05   |      |
| 1.0           | 0        | Count: 0/10   |      |
| 3.0           | 0        | Time: 00:01:00  |      |
| 5.0           | 0        | Hold: 00:00:10  |      |
| 10.0          | 0        | Recs: 123/10000   |      |
| Main          |          | Setup   | Data |
|               |          | Reports   |      |

### Screen Layout and Functionality

There are four main screens (tabs): Main, Setup, Data, and Reports. The operation of each of these screens, the information displayed on them, and the operations you can perform from each are described in the remainder of this chapter.

Some screens require or allow you to enter information. To enter information, tap on the screen and an on-screen keyboard appears.

### Main Tab


The Main Tab is the default screen. The left side of the screen summarizes the concentrations for the currently selected location. Tap on the size and count portion of the screen to enable Zoomed Data Screen (see [Setup Tab](#)).









The display shows:

- Temperature\*
- Relative humidity\*
- Bin sizes
- Particle count/concentration


\*Temperature and Humidity are displayed only if the optional T/H probe is installed.



The status bar at the top of the screen shows the current time and date settings (see the [Setup Tab](#)) and indicates:

|               |          | 11/20/08<br>11:43:03 AM   |              |
|---------------|----------|---|--------------|
|               |          | Automatic  |              |
| $\mu\text{m}$ | $\Sigma$ | Location  |              |
| 0.3           | 0        | ▼   |              |
| 0.5           | 0        | Delay: 00:00:05   |              |
| 1.0           | 0        | Count: 0/10   |              |
| 3.0           | 0        | Time: 00:01:00  |              |
| 5.0           | 0        | Hold: 00:00:10  |              |
| 10.0          | 0        | Recs: 123/10000   |              |
| Main          |          | Setup   | Data Reports |

| Icon  | Description  |
|---|--|
|    | Laser requires service                                   |
|   | Sufficient flow through the APC ERGOTOUCH PRO            |
|  | Insufficient flow through the APC ERGOTOUCH PRO          |
|  | Operating on AC power, no battery installed              |
|  | Operating on AC power, battery is installed and charging |
|  | Battery charged  |
|  | Low battery  |
|  | Battery must be charged                                  |

The right side of the Main Tab shows locations and other information (delay, cycles, and so on). These can be configured using the Setup Tab.



| Field   | Description   |
|---|---|
| Location  | Use this dropdown box to display information about any of the available locations.  |
| Delay   | The initial delay between the time the Start button is pressed and the instrument begins sampling. Valid only when Automatic mode is selected.  |
| Count   | The number of samples that have been taken/the total number of samples. Valid only when Automatic mode is selected.   |
| Time  | The time for each sample. Valid only when Automatic mode is selected.   |
| Hold  | The time between samples. Valid only when Automatic mode is selected.   |
| Recs  | The total number of records in the database/10000 (maximum number of records).  |
| Manual/Automatic/Beep   | Mode Indicator refers to the "Data Count Mode" (see section below).   |
|  | Press the Start/Stop button the begin sampling in the configured mode. Start/Stop may also be entered using the triangle-shaped button above the power button on the front of the instrument. |
|  | Prints the current sample to the optional printer.  |


## Zoomed Data Screen



The Zoomed Data screen is entered by touching in the size and count part of the main tab display. The bottom portion of the screen summarizes the concentrations for the currently selected location. Tap the size and count portion of the display to switch back to the Main Tab display.

The display shows:

- Temperature \*
- Relative humidity \*
- Bin sizes
- Particle count/concentration

\* Temperature and Humidity are displayed only if the optional T/H probe is installed.




| Field   | Description  |
|---|--|
|  | Prints the current sample to the optional printer.                     |
| Location  | Label that displays information about the currently selected location. |
|  | Press the Start/Stop button the begin sampling in the configured mode. |

## Setup Tab

The setup tab provides access to the following:


|                       |  |
|-----------------------|--|
| <b>Data Setup</b>     | View Count Units and Clear Samples.  |
| <b>System Setup</b>   | Change Power On Password, Setup Password, System Configuration, Print Settings, and Print Schedule.              |
| <b>Device Setup</b>   | Set Visual, Date and Time, and Diagnostics.  |
| <b>Sampling Setup</b> | Set up Particle Channels, Sample Timing, Sample Count Mode, Environment, Locations, and Particle Channel Alarms. |
| <b>Recipes</b>        | Save a group of settings (recipes) that you use over and over so you don't have to reset individual settings.    |






## Software Input Panel (Keyboard)

Throughout the setup screens, a keyboard will appear on the screen when text may be entered. Data may be entered using this keyboard. When the entry is complete, press either the ↵ (**Enter**) or **Esc** keys. The keyboard will then be hidden until another text entry box is selected.



## Data Setup Screen


This screen lets you access the Count Units screen and the Clear Samples screen:



## Count Units Screen


This screen lets you set the way in which particle concentration information is displayed:

| Field         | Description   |
|---------------|---|
| Differential  | Select to display particle concentration as a differential $\Delta$ (the total number of counts is the number of particles <i>between</i> bin sizes). |
| Cumulative    | Select to display particle concentration as cumulative $\Sigma$ (the total number of counts includes all particles larger than the bin size).         |
| Concentration | Display concentration in $\text{ft}^3$ or $\text{m}^3$ . If Beep mode is selected, display of concentration values is not allowed.                    |




## Clear Samples Screen

The Clear Samples screen lets you clear all samples from the internal database. Select **Yes** to clear all samples. Select **No** to return to the Data Setup screen:



## System Setup Screen

From the System Setup screen you can select (or change) the power on password, set up a password, select system configuration parameters, select print settings, and schedule printing:






## Change Power On Password Screen

If a Power On password has been previously set, you must enter that password before being allowed to change the Power On password. If a Power On password is set, then on instrument startup a password screen will ask for the password before the instrument can be used. A blank password is regarded as no password and if set as the new password, will not prompt you for a password on system startup.

Tap on the screen to display the on-screen keyboard and enter the required information.

**Note:** Keep the password in a safe place. It is very difficult to reset the password and requires contacting the factory.




| Field                | Description   |
|----------------------|---|
| Old Password         | Enter your existing password (if one has already been set).   |
| New Password         | Enter a new password. The password can be any length and use any characters.                              |
| Confirm New Password | Retype the new password then press <b>OK</b> . A confirmation message appears if the password is changed. |

## Change Setup Password Screen

If a Setup password has been previously set, you must enter that password before being allowed to change the Setup password. If a Setup password is set, clicking on the setup tab at the bottom of the main screen brings up a password screen. That password must be entered in order to change instrument settings.

Tap on the screen to display the on-screen keyboard and enter the required information.

**Note:** Entering a blank password will turn off password protection.




| Field                | Description   |
|----------------------|---|
| Old Password         | Enter your existing password. (if one has already been set).  |
| New Password         | Enter a new password The password can be any length and use any characters.                               |
| Confirm New Password | Retype the new password then press <b>OK</b> . A confirmation message appears if the password is changed. |

## System Configuration Screen


Use this screen to set system configuration parameters. Press **OK** when finished:

| Field                         | Description   |
|-------------------------------|---|
| $\Delta$ and $\Sigma$ on Zoom | Select to zoom in on both cumulative ( $\Sigma$ ) and differential ( $\Delta$ ) counts on the Main Tab. To zoom the Main Tab, click on the left side of the Main Tab. (It takes a moment for the screen to update.) Click on the screen again to return to normal view. |

**System**

 10/24/08  
7:40:33 AM

**$\Delta$  and  $\Sigma$  on Zoom**



Defaults

OK


Cancel

## Print Settings Screen


A hard copy of a sample set or statistics can be printed from the instrument using an optional thermal printer. Use this screen to set print parameters. Press **OK** when finished.

| Field         | Description  |
|---------------|--|
| Serial Number | Indicates that the serial number of the particle counter used to collect the data will be printed.               |
| Model Name    | Indicates that the model number of the particle counter used to collect the data will be printed.                |
| Separator     | Indicates a line separator will be printed after the Model Name and Serial Number in the header of all printouts |
| Differential  | Indicates that the differential value of the data will be printed.   |
| Cumulative    | Indicates that the cumulative value of the data will be printed.   |

**Print Settings**

 10/24/08  
7:41:42 AM

**Serial Number**  
 **Model Name**  
 **Separator**  
 **Differential**  
 **Cumulative**



Defaults

OK


Cancel

**Note:** Printer paper has a colored strip printed on the last meter of each roll to indicate time to change the paper roll.

## Print Schedule Screen


Use this screen to schedule automatic printing. You can choose to either print when an alarm occurs or print whenever a sample is complete:

| Field              | Description                                |
|--------------------|--|
| Automatic Printing | Enables automatic printing                 |
| On Alarm           | Print data when an alarm condition occurs. |
| On Sample          | Print data whenever a sample completes.    |



## Device Setup Screen


Use this screen to access screens that let you set or change the date and time, make audio visual selections, and run diagnostics:



## Date and Time Screen


This screen lets you set the current date and time and set the date format. Press **OK** when finished. You can select options using the arrows or tapping on the screen.


| Field       | Description  |
|-------------|--|
| Date        | Press the down arrow to display a calendar then select the date from the calendar.   |
| Time        | Select the time component you want to change (hours; minutes; seconds) and then use the left and right arrows to adjust to the current time. |
| Date Format | Highlight the date format you want to use from the list.   |
| 24 Hour     | Time display is in 24 hour format.   |



## Display Screen


This screen lets you set or change visual parameters

| Field   | Description   |
|---|---|
| Screen Alignment<br> | Press this item to reset the screen alignment, and follow the directions on the alignment screen. |




## Diagnostics Screen

This screen lets you view the system's model, serial number, copyright, manufacture date, calibration date, next calibration date, firmware version, and USB IP address. Press **Close** when finished.



## Sampling Screen


Use this screen to access screens that let you set up how sampling is displayed and handled. You can select which channels to use, the sample timing, the count mode, environment, sampling locations, and alarm thresholds.



## Channels Screen

This screen lets you choose the channels that are enabled and set their particle size. Press **OK** when finished.

| Field  | Description   |
|--------|---|
| Enable | Select the channels you want to view on the main display.   |
| Size   | If the instrument model allows variable bin sizes, this box allows for changing the default size for any channel. Highlight the size information and use the on-screen keyboard to change its value. Channels cannot be set below 0.3 or above 10.0 $\mu\text{m}$ and they may not overlap one another. |



## Sample Timing Screen

This screen lets you select parameters for sampling. Use the up and down arrows or the on-screen keyboard to change or enter information. These parameters are only valid when the Model 9306 is running in Automatic mode. Press **OK** when finished.

**Sample Timing**
10/20/08  
2:50:04 PM

**Count:**


**Delay:**

**Hold:**

**Time:**

**Volume:**


Cubic Feet
Cubic Meter



| Field  | Description  |
|--------|--|
| Count  | Count is the total number of samples you want collected. In Automatic mode, a Count value of 0 will cause the instrument to count continuously using the settings for Delay, Time, and Hold until the Start/Stop button is pressed again. Use the up and down arrows or the on-screen keyboard to set the count.   |
| Delay  | Delay indicates how long it will be before the first sample is taken. Remember, it takes approximately 6 seconds for the pump to reach the flow set point; taking a measurement before the pump is functioning properly may result in a data error.<br><br>Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value. |
| Hold   | Hold indicates how long the instrument pauses between samples. Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value.   |
| Time   | Time indicates the duration of each sample run (count particles). Highlight the time component you want to change (hours, minutes, seconds) and use the up and down arrows or the on-screen keyboard to change the value.  |
| Volume | Volume sets the volume of air that will pass through the instrument for each sample. If you select volume, you must select Cubic Feet, Cubic Meters or Liters for measurement using the arrows.  |

## Count Mode Screen


Use this screen to set the sample count mode. Press **OK** when finished.



| Field     | Description   |
|-----------|---|
| Automatic | If you select this mode, the Model 9306 starts counting in automatic mode when you press the start button according to the setting on the <a href="#">Sample Timing Screen</a> .  |
| Manual    | If you select this mode, the Model 9306 starts sampling immediately when you press the start button and stops at the end of the sample time, which is configured on the Sample Timing Screen.   |
| Beep      | If you select this mode, the Model 9306 starts sampling data immediately and beeps whenever the threshold for the smallest bin is reached, as specified in Alarms Screen. This can be very useful when searching for leaks, especially around filters. If this mode is selected, Display mode is set to Particle Counts while in Beep mode. |

## Environment Screen

Use this screen to set the units for temperature, which is displayed on the Main Tab, and the printouts when a humidity and temperature probe is hooked up to the instrument.




| Field | Description                                |
|-------|--|
| °F    | Display temperature in degrees Fahrenheit. |
| °C    | Display temperature in degrees Celsius.    |




## Locations Screen


Associating collected samples with labeled locations can help keep your data organized. The APC ERGOTOUCH PRO allows you to create up to 250 labeled locations (up to 10 characters in length). Use this screen to add, remove, or modify a location names to the list of locations.



To modify a location name, highlight the name in the list, then click the **Edit...** button. In the “Enter Location” screen click the edit box in the middle and use the on-screen keyboard to modify a location name. (You cannot edit the empty location). Click **OK** when finished.



To add a location, click on the **Add...** button. In the “Add Location” screen click in the edit box in the middle and use the on-screen keyboard to add a location name. Click **OK** when finished.



To remove a location, click on location to be removed and click the **Remove** button.


Back in the main Locations screen, after all editing has been completed, press **OK** when finished.




## Alarms Screen

Use this screen to set the alarm threshold for each channel. Press **OK** when finished.

| Field     | Description  |
|-----------|--|
| Enable    | Select the channels on which you want to enable alarms.  |
| Threshold | To change the threshold for any channel, click the up and down arrows for that channel or use the on-screen keyboard to change its value. The threshold value units use the current display Count Units (see <a href="#">Count Units Screen</a> ). |




When a channel value exceeds the threshold value you set, the channel data is highlighted in red on the Main tab, an audible alarm sounds, and the alarm icon appears on the Main tab.

To clear the alarm, click the alarm icon . In addition, the record is printed if you have selected that option on the [Print Schedule Screen](#).

## Recipe Screen


Use this screen to load and save recipes. Recipes let you save a group of settings (recipe) that you use over and over so you don't have to reset individual settings. There may be up to 100 recipes stored in the unit.



| Field   | Description   |
|---------|---|
| Save    | <p>When you select <b>Save</b>, a new window opens that lets you enter a name for the recipe you want to save. <b>The settings/parameters that are saved include:</b></p> <p>For each channel (1-6):</p> <ul style="list-style-type: none"> <li>• Alarm setting (on/off)</li> <li>• Alarm threshold (value)</li> <li>• Channel setting (enabled/disabled)</li> <li>• Channel threshold (value)</li> </ul> <p>Sample Timing settings</p> <ul style="list-style-type: none"> <li>• Count mode</li> <li>• Count total</li> <li>• Start delay (in secs)</li> <li>• Hold delay (in secs)</li> <li>• Sample time (in secs)</li> </ul> <p>Count Mode/Units Settings</p> <ul style="list-style-type: none"> <li>• Display normalized</li> <li>• Units (count, ft<sup>3</sup> or m<sup>3</sup>)</li> <li>• Cumulative/Differential</li> <li>• Volume units</li> </ul> <p>Printing settings</p> <ul style="list-style-type: none"> <li>• Auto print and mode</li> <li>• Print cumulative/differential</li> <li>• Print reverse setting (if supported)</li> <li>• Print model, separator, serial number</li> </ul> |
| Save As | When you select <b>Save As</b> , a new window opens that lets you enter a name for the recipe you want to save.   |
| Load    | Highlight the recipe you want to load and press <b>Load</b> . The settings/parameters are reset to the values of that recipe.   |
| Delete  | Highlight the recipe you want to delete and press <b>Delete</b> . The recipe is deleted.  |


## Data Tab

The Data tab lets you preview data that has been collected. Use the elevator (slide) on the right to scroll through the records. The record number is displayed at the bottom of the tab. As each record displays, its data and relevant parameters are displayed.


| 11/18/08<br>02:07:02 PM  |                       |                      |                         |
|--|-----------------------|----------------------|-------------------------|
| m <sup>3</sup>   | Size µm               | Δ m <sup>3</sup>     | Σ m <sup>3</sup>        |
|  | 0.3                   | 275695               | 286863                  |
|  | 0.5                   | 9733                 | 11168                   |
|  | 1.0                   | 1214                 | 1435                    |
|   | 3.0                   | 53                   | 221                     |
|  | 5.0                   | 145                  | 168                     |
| <b>Print</b>   | 10.0                  | 23                   | 23                      |
| Location: Sample 003 Laser: OK Alarm: YES<br>Sample: 00:01:00 Vol: 131.0 L Flow: ALRM<br>Date: 11/18/08 Temperature: 24.4 °C<br>Time: 11:59:56 AM Humidity: 53 % |                       |                      |                         |
| <b>Record: 20 Records: 123 / 10000</b>   |                       |                      |                         |
| <a href="#">Main</a>   | <a href="#">Setup</a> | <a href="#">Data</a> | <a href="#">Reports</a> |

| Field                               | Description  |
|-------------------------------------|--|
| #, ft <sup>3</sup> , m <sup>3</sup> | Button used to change between counts and concentration displays.                                   |
| Size µm                             | Channel size.  |
| Δ                                   | Differential concentration.  |
| Σ                                   | Cumulative concentration.  |
| Location                            | Location where the data was collected.   |
| Sample                              | Duration of the sampling period.   |
| Date                                | Date on which the data was collected.  |
| Time                                | Time at which data was collected.  |
| Temperature                         | Temperature at the end of the time the data was collected (if probe connected during sampling).    |
| Humidity                            | Humidity level at the end of the time the data was collected (if probe connected during sampling). |
| Flow                                | Status of the flow.  |
| Alarm                               | Alarm threshold was triggered (YES) or not (NONE).   |
| Laser                               | Status of the laser.   |

The print button will allow a range of sample data to be printed using the optional external printer.





The print data screen will show progress on the current selected range of sample data to be printed. Press the **Cancel Printing** button to cancel the rest of the print job.





## Reports Tab


Use this screen to select various standard reports for viewing and printing

 10/03/08  
9:40:55 AM

  
**Federal Standard FT3**


  
**Federal Standard M3**



  
**ISO 14644-1**

  
**EC GMP**


Main
Setup
Data
Reports



The standard reports are shown below:

**Fed Std 209E FT**
 11/20/08  
11:55:25 AM


|                 |                       |  |   |
|-----------------|-----------------------|--|---|
| Room Area:      | 96.88 ft <sup>2</sup> |  |  |
| Class Level:    | 100                   |  |   |
| Air Flow:       | Unidirectional        |  |   |
| Room Status:    | Operational           |  |   |
| Min Locations:  | 4                     |  |   |
| Min Samples:    | 5                     |  |   |
|                 |                       |  |   |
| Allowable Sizes | Min Volume            |  |  |
| 0.2 um          | 1.00E-001 f3          |  |   |
| 0.3 um          | 1.00E-001 f3          |  |   |
| 0.5 um          | 2.00E-001 f3          |  |   |



Close

**Fed Std 209E M**
 11/20/08  
11:56:35 AM


|                 |                     |  |   |
|-----------------|---------------------|--|---|
| Room Area:      | 9.00 m <sup>2</sup> |  |  |
| Class Level:    | M2                  |  |   |
| Air Flow:       | Unidirectional      |  |   |
| Room Status:    | Operational         |  |   |
| Min Locations:  | 4                   |  |   |
| Min Samples:    | 5                   |  |   |
|                 |                     |  |   |
| Allowable Sizes | Min Volume          |  |  |
| 0.1 um          | 5.71E+000 L         |  |   |
| 0.2 um          | 2.64E+001 L         |  |   |
| 0.3 um          | 6.47E+001 L         |  |   |
| 0.5 um          | 2.00E+002 L         |  |   |



Close

**EC-GMP**
 11/20/08  
09:01:59 AM

|                 |                     |  |   |
|-----------------|---------------------|--|---|
| Room Area:      | 9.00 m <sup>2</sup> |  |  |
| Class Level:    | C                   |  |   |
| Air Flow:       | Unidirectional      |  |   |
| Room Status:    | At-Rest             |  |   |
| Min Locations:  | 3                   |  |   |
| Min Samples:    | 3                   |  |   |
|                 |                     |  |   |
| Allowable Sizes | Min Volume          |  |  |
| 0.5 um          | 2.00E+000 L         |  |   |
| 5.0 um          | 6.90E+000 L         |  |   |

Close

**ISO 14644-1**
 11/20/08  
09:00:51 AM

|                 |                     |  |   |
|-----------------|---------------------|--|---|
| Room Area:      | 9.00 m <sup>2</sup> |  |  |
| Class Level:    | 3                   |  |   |
| Air Flow:       | Unidirectional      |  |   |
| Room Status:    | At-Rest             |  |   |
| Min Locations:  | 3                   |  |   |
| Min Samples:    | 3                   |  |   |
|                 |                     |  |   |
| Allowable Sizes | Min Volume          |  |  |
| 0.1 um          | 2.00E+001 L         |  |   |
| 0.2 um          | 8.44E+001 L         |  |   |
| 0.3 um          | 1.96E+002 L         |  |   |
| 0.5 um          | 5.71E+002 L         |  |   |
| 1.0 um          | 2.50E+003 L         |  |   |

Close

| Field           | Description  |
|-----------------|--|
| Room Area       | Displays the area of the room in ft <sup>2</sup> or m <sup>2</sup> .                             |
| Class Level     | Depends on the report definition, see below.   |
| Air Flow        | Displays the airflow characteristics of the room.  |
| Room Status     | Displays the status of the room. See Room Definition below.                                      |
| Min Locations   | Displays the minimum number of locations that must be sampled in the room.                       |
| Min Samples     | Displays the minimum number of samples that must be taken at each location.                      |
| Allowable Sizes | Allowable channel sizes for the selected Class Level for that Standard.                          |
| Minimum Volume  | Displays the minimum volume (in cubic feet or meters) that must be sampled on each channel.      |
| Room Definition | Press to set definitions for the room. (See <a href="#">Room Definition Screen</a> below.)       |
| Generate        | Select to print a single record or a range of records. (See <a href="#">Print Screen</a> below.) |

### Room Definition Screen

Use this screen to define specific values for the room. Press **OK** when finished:


**Room Definition**
10/03/08  
9:45:27 AM

**Room Status**

**Air Flow**

**Class**


**Area**  
  ft<sup>2</sup>  m<sup>2</sup>



| Field       | Description  |
|-------------|--|
| Room Status | Select the room status: As Built, At Rest, or Operational.   |
| Air Flow    | Select the air flow: Unidirectional or Non-unidirectional.   |
| Class       | Select the class of the room: The class is dependent on the standard:<br>FED FT3: 1, 20, 100,1000,10000, 100000<br>FED M3: M1.0, M1.5, M2.0, M3.0, M3.5, M4.0, M4.5, M5.0, M5.5, M6.0, M6.5, M7.0<br>ISO14644-1: 1, 2, 3, 4, 5, 6, 7, 8, 9<br>EC GMP: A, B, C, D |
| Area        | Use the on-screen keyboard to enter the area of the room in ft <sup>2</sup> or m <sup>3</sup> .  |

## Generate Screen

This screen lets you generate the report using either a single record or a range of records. Press the **Generate...** button to generate the selected report.



**Generate...** 10/03/08 9:47:47 AM


Select range used for reports.

All  
 Records

1 – 107

Generate...  
Close

The generated report will be displayed on the screen and may be viewed on the screen or printed (optional printer must be attached) by pressing the **Print** button.



**Fed Std 209E M** 10/03/08 9:52:23 AM

\*\*\*\*\*

Federal Standard 209E\_M  
\*\* Model 9306 \*\*  
Serial Number: xxxxxxxxx  
Targeted Class: M7.0  
Room Area: 100m<sup>2</sup>  
Room Status: Operational  
Air Flow: Unidirectional  
Min Locations: 3  
Min Samples/Room: 5

Print Close

## Chapter 4 Data Handling

### USB Computer Communication

The APC ERGOTOUCH PRO Airborne Particle Counter is equipped with a USB compatible cable for uploading and downloading information to a PC. The cable plugs into the right side of the instrument.



### Installing Software

The ERGOTRAK Lite Data Transfer utility comes on a CD that loads software and communications drivers for the particle counter. To install the software, insert the CD into your computer drive and follow the instructions. Installation consists of two parts:


- Installation of ERGOTRAK Lite software.  
Run “setup.exe” from the provided CD and follow on-screen instructions.
- Installation of USB NDIS driver. This installation is executed transparently during the setup process and does not require user input. Once installation is finished, drivers are ready for use. When the APC ERGOTOUCH PRO particle counter is connected for the first time, system will automatically detect the device and will start driver installation process.

## Procedure

1. When asked if Windows update should be used to download necessary software, select “No, not this time” and click **Next**.




2. Select “Install the software automatically (Recommended)” and click **Next**.



3. Hardware Wizard will search for the driver and locate it in \System32\drivers directory.

Once that is done, the following screen will appear:






- Depending on your system setup you may see a warning message:



- Click **Continue Anyway** and the installation will proceed.
- Once everything is completed, click **Finish**.

This procedure is required only on first connection, all subsequent devices will automatically locate the necessary drivers and install without requiring user input.




## Download Data

To transfer data from an instrument to a computer via the USB connection for further analysis and report generation:


Make sure that ERGOTOUCH PRO particle counter is attached to the computer and turned on.

Start the application. If the ERGOTOUCH PRO particle counter is not connected or discovered by the application, the following message appears:




- Make sure that the ERGOTOUCH PRO particle counter is connected, turned on, and functioning properly. Restart the application.

- If communication with ERGOTOUCH PRO particle counter has been established, the following screen appears:



- Once data is successfully downloaded, the main application screen will appear:




1. There are two options for downloading data:

- Download only user selected records:  
Hold down the **CTRL** key and use the mouse to click on the records you want to retrieve. When you have selected the records, press the **Get Selected** button to retrieve only the selected records from the device.
- Download all records:  
Press the **Download All** button to retrieve all the records from the device.

2. After you press either the **Get Selected** or **Download All** buttons, the following dialog appears to allow you to select the folder where data will be saved:


- To cancel the data transfer, select **Cancel**
- To accept the data transfer, enter the file name under which data will be saved and select **Save**.
- Data is stored in a .CSV file that can be opened by most spreadsheet programs such as Microsoft® Excel® spreadsheet software.



## Delete Data

In order to delete data from the device click **Delete Records** button. The following warning will appear:

If **Yes** is selected, ERGOTRAK Lite software will erase data from the device and also from application memory.



**W A R N I N G**

Deleting data is an irreversible operation. Download and save data *before* deleting in order to have a copy for future use.

## Chapter 5

### Maintenance

#### Note

There are no user-serviceable parts inside this instrument. Opening the instrument case may void the warranty. Biotest recommends that you return the APC ERGOTOUCH PRO Airborne Particle Counter to the factory for any required maintenance or service not described in this manual.

### Maintenance Schedule

Biotest recommends annual factory cleaning and calibration for the APC ERGOTOUCH PRO Airborne Particle Counter. See [Chapter 7, "Contacting Customer Service"](#) for service/calibration.

### Recommended Field Maintenance Schedule

| Item                              | Frequency                   |
|-----------------------------------|-----------------------------|
| Daily zero check                  | Daily (or before each use). |
| Factory cleaning and calibration  | Annually.                   |
| Cleaning the instrument enclosure | As needed                   |

### Daily Zero Check

The daily zero check ensures that the instrument is properly assembled and free from leaks, residual particles and electronic noise. Please see Chapter 2, "[Getting Started](#)" for detailed instructions on performing the zero check.




### Cleaning the Instrument Enclosure

To clean the enclosure, dampen a lint-free cloth and gently wipe the surface until surface contamination is removed.

## Chapter 6

### Troubleshooting

| Symptom  | Possible Cause  | Corrective Action  |
|--|---|--|
| Counts are too low   | <ul style="list-style-type: none"> <li>Instrument is being operated outside temperature or relative humidity specifications.</li> <li>Internal parts have been damaged because instrument was stored at a temperature greater than 122 °F (50 °C).</li> <li>Instrument has contamination on the optics due to condensation or excessive loading</li> <li>Laser or pump control is damaged</li> <li>Unit is due for calibration</li> </ul>   | <ul style="list-style-type: none"> <li>Operate instrument within specifications</li> <li>Return to factory for service</li> <li>Return to factory for service</li> <li>Return to factory for service</li> <li>Return to factory for service</li> </ul>   |
| Instrument does not turn on  | <ul style="list-style-type: none"> <li>Battery is dead</li> <li>AC cord is not plugged into unit</li> </ul>   | <ul style="list-style-type: none"> <li>Recharge battery or connect to AC power</li> <li>Connect AC cord</li> </ul>   |
| Instrument does not meet zero count specification (<1 particle/5 mins) | <ul style="list-style-type: none"> <li>HEPA filter is not connected properly and room air is leaking into the HEPA filter assembly</li> <li>Residual particles from previous samples are shedding off internal parts and into the optics</li> <li>An internal component has been damaged due to operation outside of temperature specifications or one or more excessive bumps or jolts, and electronic noise is inducing false counts</li> <li>A leak has developed in the aerosol flow path</li> <li>Internal optics have become dirty</li> </ul> | <ul style="list-style-type: none"> <li>Check that the HEPA filter has been tightly connected to the inlet. Check that rubber o-ring (black) on the inlet is in place</li> <li>Purge instrument by running the instrument for 10–15 minutes before attempting zero count test</li> <li>Return to factory for service</li> <li>Return to factory for service</li> <li>Return to factory for service</li> </ul> |
| Battery does not charge  | <ul style="list-style-type: none"> <li>The unit must be turned on but not in sampling mode for the battery to charge</li> </ul>   | <ul style="list-style-type: none"> <li>Turn on unit.</li> </ul>  |

| Symptom  | Possible Cause  | Corrective Action  |
|--|---|--|
| LOW BATTERY ERROR<br>     | <ul style="list-style-type: none"> <li>Low battery</li> </ul>   | <ul style="list-style-type: none"> <li>Recharge battery or connect AC cord</li> </ul>  |
| PHOTODETECTOR ERROR  | <ul style="list-style-type: none"> <li>Direct light is entering the aerosol inlet</li> <li>Laser has become misaligned due to excessive bumps or jolts</li> <li>Internal optics have become dirty</li> </ul>                              | <ul style="list-style-type: none"> <li>Remove instrument from direct light</li> <li>Return to factory for service</li> <li>Return to factory for service</li> </ul>  |
| SYSTEM ERROR   | <ul style="list-style-type: none"> <li>Information is not being read properly by microprocessor</li> </ul>  | <ul style="list-style-type: none"> <li>Restart instrument. If problem persists, contact Biotest technical support</li> </ul>   |
| TEMPERATURE HUMIDITY PROBE ERROR   | <ul style="list-style-type: none"> <li>Temperature/RH probe was not recognized</li> </ul>   | <ul style="list-style-type: none"> <li>Detach and reconnect probe. If problem persists, contact Biotest technical support</li> </ul>   |
| FLOW ERROR<br>          | <ul style="list-style-type: none"> <li>Instrument was unable to control flow rate (if any tubing is connected to particle counter)</li> <li>Pressure drop across inlet may be too large</li> <li>Inlet not at ambient pressure</li> </ul> | <ul style="list-style-type: none"> <li>Restart measurement</li> <li>Lessen pressure drop across inlet by using larger diameter tubing, less tubing, and/or adding a bleed valve</li> <li>Do <b>not</b> subject the unit to other than ambient pressure conditions</li> </ul> |
| LASER POWER WARNING<br> | <ul style="list-style-type: none"> <li>Laser power has fallen outside of specification</li> </ul>   | <ul style="list-style-type: none"> <li>Return to factory for service</li> </ul>  |

## Chapter 7

### **Contacting Customer Service**

This chapter gives directions for contacting people at Biotest AG for technical information and directions for returning the APC ERGOTOUCH PRO Airborne Particle Counter for service.

#### **Technical Contacts**

- If you have any difficulty setting up or operating the APC ERGOTOUCH PRO, or if you have technical or application questions about this system, contact an applications engineer at Biotest AG, +49.6221.726.5130 or e-mail [hycon@biotest.de](mailto:hycon@biotest.de)
- If the APC ERGOTOUCH PRO, does not operate properly, or if you are returning the instrument for service, visit our website at <http://www.biotest.com>, or contact Biotest Customer Service at +49.6103.801.496.
- Call your local Biotest Subsidiary or Representative.

#### **Returning the APC ERGOTOUCH PRO Airborne Particle Counter for Service**

Visit our website at <http://www.biotest.com> or call Biotest Technical Support at +49.6221.726.5130 or Biotest Customer Service at +49.6103.801.590 for specific return instructions. Customer Service will need this information when you call:

- The instrument model number
- The instrument serial number
- A purchase order number (unless under warranty)
- A billing address
- A shipping address

Use the original packing material to return the instrument to Biotest. If you no longer have the original packing material, seal off any ports to prevent debris from entering the instrument and ensure that the display and the connectors on the instrument front and back panels are protected.

## APPENDIX A

### Specifications

All specifications meet or exceed JIS B 9921. They are subject to change without notice.

|                        |  |
|------------------------|--|
| Size Range             | 0.3–10 µm  |
| Channel Sizes          | Standard: 0.3, 0.5, 1.0, 3.0, 5.0, 10.0 µm                                     |
| Counting Efficiency    | 50% @ 0.3 µm; 100% for particles > 0.45 µm (per JIS)                           |
| Concentration Limits   | 2,000,000 particles/ft <sup>3</sup> @ 5% coincidence loss                      |
| Light Source           | Laser diode  |
| Zero Count Level       | <1 count/5 minutes (<2 particles/ft <sup>3</sup> ) Meets JIS B9921             |
| Flow Rate              | 0.1 CFM (2.83 L/min) with ±5% accuracy   |
| Flow Control           | Automatic  |
| Calibration            | NIST traceable   |
| Sample Probe/Tubing    | Isokinetic sampling probe  |
| Sampling Modes         | Manual, automatic, beep, concentration and cumulative/differential             |
| Sampling Time          | 1 second to 99 min 59 sec  |
| Sampling Frequency     | 1 to 999 cycles or continuous  |
| Sample Output          | Internal HEPA filter   |
| Vacuum Source          | Internal pump  |
| Communication Mode     | RS-232, USB  |
| Data Storage           | 10,000 samples   |
| Data Security          | Password protected   |
| External Alarm         | Internal audible alarm. Indicators for counts, low battery, and sensor failure |
| Environmental Sensors  | Optional temperature/RH probe supported  |
| Display                | QVGA 3.7-inch touch screen with Windows CE                                     |
| Languages              | English  |
| Reports                | ISO-14644-1, FS-209E & EC GMP  |
| Printer                | Optional external printer supported  |
| External Surface       | High impact injection molded plastic   |
| Power                  | 110 to 240 VAC 50 to 60 Hz Universal in-line power supply                      |
| Battery                | Removable Li-Ion   |
| Battery Life           | Up to 3 hours of continuous use  |
| Recharge Time          | 3 hours  |
| Dimensions (L x W x H) | 25.4 x 11.4 x 7.6 cm (10 x 4.5 x 3 in.)  |
| Weight                 | 1.0 kg (2.2 lbs) with battery  |
| Warranty               | 2 years. Extended warranties available   |



|                      |  |
|----------------------|--|
| Operating Conditions | 5°C to 35°C; 20% to 95% non-condensing relative humidity   |
| Storage Conditions   | 0°C to 50°C; Up to 98% non-condensing relative humidity  |
| Included Accessories | Power supply, power cord, battery, isokinetic inlet, stylus, purge filter, ERGOTRAK Lite data download software, operational manual on CD, computer cable, calibration certificate, and Quick Start Guide. |
| Optional Accessories | Temp R/H probe, stainless steel isokinetic inlet and probe, tubing, barbed inlet fitting, printer, printer paper, carrying case and external battery charger   |

***Temperature / RH Probe Specifications (optional accessory)***

|   |                         |
|---|-------------------------|
| <b>Temperature</b><br>Range .....       | 32 to 115°F (0 to 45°C) |
| Accuracy .....                          | ±2°F (±1°C)             |
| <b>Relative Humidity</b><br>Range ..... | 10 to 90% RH            |
| Accuracy .....                          | ±5% RH                  |