SERVICE MANUAL FOR ANALOG CONTROLLER MODELS

Information is reserved for trained, qualified technical personnel. It includes schematics & wiring diagrams. Any repair or service to a unit is to be performed by a certified technician ONLY. Attempting repairs outlined in this guide may result in voiding your warranty





Models 611, 71, 911

TABLE OF CONTENTS

- Section 1. Troubleshooting Procedure
- Section 2. Field Calibration Check
- Section 3. Testing Temperature Probe
- Section 4. Replacing Control Board
- Section 5. Spare Parts List
- Section 6. Schematics & Wiring Diagram

SECTION 1. TROUBLESHOOTING PROCEDURE

Note: If the unit seems to function normally but fails to control or hold temperature steady, first perform instructions as in "Conditions and Steps before Troubleshooting the Unit"

CONDITIONS AND STEPS BEFORE TROUBLESHOOTING THE UNIT

- 1. Set the HI/LO pump speed switch on the rear panel to HI.
- 2. Set up unit.
 - a. For a closed loop application on the model 911 circulator, connect a 1/4 inch diameter or greater hose (a few feet long in length) to the Inlet and then Outlet.
 - b. For the model 71 immersion circulator, clamp to the side of a bath container.
- 3. Use the same fluid and the same setpoint temperature as your application. Always use enough fluid to cover pump housing and heater.
- 4. Unit should visibly pump the fluid inside the reservoir. If not, check that the pump shaft spins freely and impeller(s) are in place. If bearings are worn out, replace the entire pump motor.
- 5. Check the fluid or fluid mixture:
 - a. Watery consistancy, not thick like maple syrup. Viscosity <25 centistokes.
 - b. Not boiling on the surface or freezing at the bottom of the reservoir or to the cooling coils.
 - c. Not within 5°C of the freezing or boiling point of the fluid.
- 6. On the model 911, if used for an external application:
 - a. Check the heat load does not exceed the rated cooling power of the refrigeration, or the external heat loss does not exceed the heater's 1000 watts. (Fluids within 40°C of their boiling point can cause so much heat loss as to require more than the 1000 watts of heat provided by the unit, especially in containers with a large open surface area.)
 - b. Cover/uncover and add/remove insulation of both the internal reservoir and external application and tubing. Insulate the top of the fluid with hollow floating plastic or Teflon balls.

PROBLEM	Corrective Action
#1 Unit does not heat	a) Check power supply voltage shown on circuit board diagram. If low, check R8(and R9 on 240V units) for proper resistance. If still low or high, replace U1. If still low, check C1 for shorted or leaky condition. Also check R2 for proper resistance.
or intermittent heat.	b)Check all test voltages on the circuit board diagram.
	c) Check heater H1 resistance. $120V = 14.4$ ohms Ω , $240V = 58$ ohms Ω .
	d)Check continuity of and to heater H1 and triac Q1.
	e)Try replacing U1, Q1 or see Replacing Main Control Board.
#2	a)Measure for correct line voltage while the bath is heating.
Unit won't heat to set point.	b) Check for DC voltage on the heater. If more than 30 VDC, replace triac Q1. Also check H1 heater resistance.
#3 Unit continuously	a) Remove U1 IC and test again. If still heating, replace Q1. If stopped heating, Q1 is okay. Re-insert IC
heats.	b)Check all test voltages on the circuit board diagram.
	c)Check if Q1 is properly heat sinked with heat sink compound.

#4 Unit does not pump.	 a)Place rear speed switch in HI position. Unit should visibly pump the fluid inside of the reservoir. b)Check that the pump shaft spins freely and pump impeller(s) are in place. If bearings are worn out, replace entire pump motor.
#5 Unit very slowly rises above the setpoint and won't come back down.	 a)Without refrigeration, the unit's pumping and stirring will generate heat. Check heat rise specification above room ambient. Try low pump speed, uncovering bath(s), removing insulation, use low viscosity fluid, add tap water cooling or a refrigeration unit. b)On model 911 turn cooling on. Set to 20°C. To check that refrigeration is working properly, feel for warm air exiting out of the rear vent.
#6 Controls only near 40°C.	P1 or R6 are not making contact or are defective.
#7 Setpoint knob inaccurate.	See Field Calibration Check, section 2.

SECTION 2. FIELD CALIBRATION CHECK

There is no calibration adjustment or procedure. The setpoint dial accuracy of these units are rated to be within $\pm 10^{\circ}$ C of NIST actual temperature over the entire temperature range.

- Check calibration against a NIST traceable thermometer of known accuracy of ± 1°C or better.
- Make sure the setpoint knob is mechanically centered across the dial range.
- Check all test voltages on the circuit board diagram.
- Check temperature probe as below.

SECTION 3. TESTING TEMPERATURE PROBE

Remove the probe wires from the board. Check that the resistance of the probe is within the tolerance shown on the schematic & wiring diagram on page 7. If okay, then check electrical readings from the probe metal sheath to the two probe leads on the following two items:

- Check for any resistance reading. If there is any, the probe is defective.
- On a digital voltmeter with an impedance of one megohm or more, check for any DC voltage present. If there is some, the probe is defective.

SECTION 4. REPLACING MAIN CONTROL BOARD

When removing/replacing the control board, always remove/replace the bath temperature probe included with the board.

Temperature probes are checked at the factory to be within tolerance only when used with that board.

HOW TO INSTALL A PC BOARD:

To remove the existing PC board:

- Unplug the unit from the power line.
- Remove the two screws holding the cover in place. Once the cover is loose, disconnect the ground wire from the lug terminal.

- On a sheet of paper, note the colors and locations of the wires on the PC Board. Or use the wiring diagram on pg. 7..
- Remove all wires from their terminal connectors on the PC board.
- Loosen the probe clip from its mounting bracket.
- Remove the knob insert, knob and nut holding the front heater control (potentiometer).
- Remove the three PC board mounting screws.
- Slide the probe up and away from the probe clip.

To install the replacement PC board (Follow the steps above in reverse order):

- Slide the probe back into the probe clip.
- Install the three PC board mounting screws.
- Install the knob insert, knob and nut holding front heater control.
- Tighten probe clip into its mounting bracket.
- Attach all wires to terminal connectors according to your diagram.
- Attach the ground wire to the lug terminal.
- Attach the cover using two screws.
- Plug power line into unit.

SECTION 5. SPARE PARTS LIST

REF#	DESCRIPTION	PART#
Q1	TRIAC, ISOLATED TAB, TECCOR, Q4015L5, 15A, 400V	200-058
M1	MOTOR, PUMP, MODEL 911, 120V, SHAFT LENGTH 6-11/16 INCHES FROM MOUNTING SURFACE	215-235
M1	MOTOR,PUMP,MODEL 911,240V,SHAFT LENGTH 6-11/16 INCHES FROM MOUNTING SURFACE	215-236
M1	MOTOR,PUMP,MODEL 71,120V,SHAFT LENGTH 7-11/16 INCHES FROM MOUNTING SURFACE	215-237
M1	MOTOR, PUMP, MODEL 71, 240V, SHAFT LENGTH 7-11/16 INCHES FROM MOUNTING SURFACE	<u>215-238</u>
PCB1	71-911 PC BOARD ASM., WITH RTD PROBE, TESTED, 120V	500-107
PCB1	71-911 PC BOARD ASM., WITH RTD PROBE, TESTED, 240V	500-108
H1	HEATER,TUBULAR,3 1/2" O.D. COIL,1KW,120V	215-067
H1	HEATER,TUBULAR,3 1/2" O.D. COIL,1KW,240V	215-068
NE1,NE2	LAMP,NEON,GREEN,RIGHT ANGLE,120V	215-161
NE1,NE2	LAMP,NEON,GREEN,RIGHT ANGLE,240V	215-139
M2	FAN FOR REFRIGERATION, AXIAL, 4" INCH DIAMETER, 120V	215-296
M2	FAN FOR REFRIGERATION, AXIAL, 4" INCH DIAMETER, 240V	215-297

S1	TEMPERATURE PROBE,7" X 1/4" O.D.,RTD,1855 ohm $\Omega @ 0^\circ\text{C}$	200-078
CB1	CIRCUIT BREAKER, MAGNETIC, 1-POLE, 12A, 120V	215-041
CB2	CIRCUIT BREAKER, MAGNETIC, 2-POLE, 15A, 240V	215-157
SW1,SW2	SWITCH, GREEN INDICATOR, MAIN POWER OR COOLING POWER	235-026
OTP	OTP SAFETY HYDRAULIC THERMOSTAT, IMMERSION, 120V	215-122
OTP	OTP SAFETY HYDRAULIC THERMOSTAT, IMMERSION, 240V	215-167
SW3	SWITCH,SLIDE,FOR PUMP HIGH/LOW SPEED	235-008

GRAM

 \triangleleft

σ

111

P

שמ

TPOINT / MIN 3-8.9V MAX AT 40°C BATH TEMP 4.15V MIN 4.45V MAX BOLTPL BOLTPL M H 911 TEST VOLTAGES. WHEN MEASURING. ENTIRE CIRCUIT IS AT LINE VOLTAGE POTENTIAL ß 6 ভি TRIANGLE WAVE SWINGS 1V TO GV 6 SUPPLY VOLTAGE 8.3V MIN 8.9V MAX 6 $\langle 0 \rangle$ (0)6 ଚ 0 (c S1 0 0 ୕୲ 01 ک -108/240SETPUINT SETPUINT 3.9V MIN 4.7V MAX 700mV A C $\langle 0 \rangle$ \mathbb{R}^2 6 HTR O 500-107/120/ Di С E 2 N N L INE VDL TAGE 0 \odot METER REFERENCE POINT 0 0 PolyScience \odot NEI 6 0 NE1 0 NE2 0 R8-R9 MTR (و_ 0 HTR 0 DLK BRN SW2 0 0 0 0 OTP BLU 5 0 0 0 2 0 CB/SW1 0 0 COMP 0 0 MTR 0 ANC 0 0 0 Q1 CAU С SW2 MODEL TAKE (٥ \bigcirc



Schematic - Circuit Board



Schematic and Wiring Diagram - Electrical